



# ANALYTICAL REPORT

## PREPARED FOR

Attn: Rita Schmon-Stasik  
Leidos, Inc.  
Picatinny Arsenal  
356 Ninth Avenue  
Suite 106  
Dover NJ 07801  
Generated 5/31/2023 9:08 PM

## JOB DESCRIPTION

RVAAP FWGW

## JOB NUMBER

280-176674-1

# Eurofins Denver

## Job Notes

This report may not be reproduced except in full, and with written approval from the laboratory. The results relate only to the samples tested. For questions please contact the Project Manager at the e-mail address or telephone number listed on this page.

The test results in this report relate only to the samples as received by the laboratory and will meet all requirements of the methodology, with any exceptions noted. This report shall not be reproduced except in full, without the express written approval of the laboratory. All questions should be directed to the Eurofins TestAmerica Project Manager.

## Authorization



Generated  
5/31/2023 9:08 PM

Authorized for release by  
Patrick J McEntee, Client Service Manager  
[Patrick.McEntee@et.eurofinsus.com](mailto:Patrick.McEntee@et.eurofinsus.com)  
303 736-0107



# Table of Contents

|                                      |     |
|--------------------------------------|-----|
| Cover Title Page . . . . .           | 1   |
| Data Summaries . . . . .             | 5   |
| Definitions . . . . .                | 5   |
| Case Narrative . . . . .             | 6   |
| Detection Summary . . . . .          | 8   |
| Client Sample Results . . . . .      | 10  |
| Default Detection Limits . . . . .   | 14  |
| Surrogate Summary . . . . .          | 15  |
| QC Sample Results . . . . .          | 16  |
| QC Association . . . . .             | 24  |
| Chronicle . . . . .                  | 27  |
| Certification Summary . . . . .      | 29  |
| Method Summary . . . . .             | 30  |
| Sample Summary . . . . .             | 31  |
| Manual Integration Summary . . . . . | 32  |
| Reagent Traceability . . . . .       | 58  |
| COAs . . . . .                       | 67  |
| Organic Sample Data . . . . .        | 364 |
| HPLC/IC . . . . .                    | 364 |
| 8330B_DOD5 . . . . .                 | 364 |
| 8330B_DOD5 QC Summary . . . . .      | 365 |
| 8330B_DOD5 Sample Data . . . . .     | 377 |
| Standards Data . . . . .             | 431 |
| 8330B_DOD5 ICAL Data . . . . .       | 431 |
| 8330B_DOD5 CCAL Data . . . . .       | 655 |
| Raw QC Data . . . . .                | 782 |

# Table of Contents

|                                            |      |
|--------------------------------------------|------|
| 8330B_DOD5 Blank Data . . . . .            | 782  |
| 8330B_DOD5 LCS/LCSD Data . . . . .         | 789  |
| 8330B_DOD5 MS/MSD Data . . . . .           | 815  |
| 8330B_DOD5 Run Logs . . . . .              | 846  |
| 8330B_DOD5 Prep Data . . . . .             | 853  |
| Inorganic Sample Data . . . . .            | 862  |
| General Chemistry Data . . . . .           | 862  |
| Gen Chem Cover Page . . . . .              | 863  |
| Gen Chem Sample Data . . . . .             | 864  |
| Gen Chem QC Data . . . . .                 | 870  |
| Gen Chem ICV/CCV . . . . .                 | 870  |
| Gen Chem Blanks . . . . .                  | 876  |
| Gen Chem MS/MSD/PDS . . . . .              | 877  |
| Gen Chem LCS/LCSD . . . . .                | 879  |
| Gen Chem Calibration . . . . .             | 882  |
| Gen Chem MDL . . . . .                     | 888  |
| Gen Chem Preparation Log . . . . .         | 900  |
| Gen Chem Analysis Run Log . . . . .        | 901  |
| Gen Chem Prep Data . . . . .               | 922  |
| Gen Chem Raw Data . . . . .                | 939  |
| Shipping and Receiving Documents . . . . . | 1195 |
| Client Chain of Custody . . . . .          | 1196 |
| Sample Receipt Checklist . . . . .         | 1207 |

# Definitions/Glossary

Client: Leidos, Inc.  
Project/Site: RVAAP FWGW

Job ID: 280-176674-1

## Qualifiers

### HPLC/IC

| Qualifier | Qualifier Description                                                                                                           |
|-----------|---------------------------------------------------------------------------------------------------------------------------------|
| J         | Estimated: The analyte was positively identified; the quantitation is an estimation                                             |
| J1        | Estimated: The quantitation is an estimation due to discrepancies in meeting certain analyte-specific quality control criteria. |
| M         | Manual integrated compound.                                                                                                     |
| Q         | One or more quality control criteria failed.                                                                                    |
| U         | Undetected at the Limit of Detection.                                                                                           |

### General Chemistry

| Qualifier | Qualifier Description                                                                                         |
|-----------|---------------------------------------------------------------------------------------------------------------|
| D         | The reported value is from a dilution.                                                                        |
| H         | Sample was prepped or analyzed beyond the specified holding time. This does not meet regulatory requirements. |
| J         | Estimated: The analyte was positively identified; the quantitation is an estimation                           |
| M         | Manual integrated compound.                                                                                   |
| U         | Undetected at the Limit of Detection.                                                                         |

## Glossary

| Abbreviation   | These commonly used abbreviations may or may not be present in this report.                                 |
|----------------|-------------------------------------------------------------------------------------------------------------|
| ▫              | Listed under the "D" column to designate that the result is reported on a dry weight basis                  |
| %R             | Percent Recovery                                                                                            |
| CFL            | Contains Free Liquid                                                                                        |
| CFU            | Colony Forming Unit                                                                                         |
| CNF            | Contains No Free Liquid                                                                                     |
| DER            | Duplicate Error Ratio (normalized absolute difference)                                                      |
| Dil Fac        | Dilution Factor                                                                                             |
| DL             | Detection Limit (DoD/DOE)                                                                                   |
| DL, RA, RE, IN | Indicates a Dilution, Re-analysis, Re-extraction, or additional Initial metals/anion analysis of the sample |
| DLC            | Decision Level Concentration (Radiochemistry)                                                               |
| EDL            | Estimated Detection Limit (Dioxin)                                                                          |
| LOD            | Limit of Detection (DoD/DOE)                                                                                |
| LOQ            | Limit of Quantitation (DoD/DOE)                                                                             |
| MCL            | EPA recommended "Maximum Contaminant Level"                                                                 |
| MDA            | Minimum Detectable Activity (Radiochemistry)                                                                |
| MDC            | Minimum Detectable Concentration (Radiochemistry)                                                           |
| MDL            | Method Detection Limit                                                                                      |
| ML             | Minimum Level (Dioxin)                                                                                      |
| MPN            | Most Probable Number                                                                                        |
| MQL            | Method Quantitation Limit                                                                                   |
| NC             | Not Calculated                                                                                              |
| ND             | Not Detected at the reporting limit (or MDL or EDL if shown)                                                |
| NEG            | Negative / Absent                                                                                           |
| POS            | Positive / Present                                                                                          |
| PQL            | Practical Quantitation Limit                                                                                |
| PRES           | Presumptive                                                                                                 |
| QC             | Quality Control                                                                                             |
| RER            | Relative Error Ratio (Radiochemistry)                                                                       |
| RL             | Reporting Limit or Requested Limit (Radiochemistry)                                                         |
| RPD            | Relative Percent Difference, a measure of the relative difference between two points                        |
| TEF            | Toxicity Equivalent Factor (Dioxin)                                                                         |
| TEQ            | Toxicity Equivalent Quotient (Dioxin)                                                                       |
| TNTC           | Too Numerous To Count                                                                                       |

## CASE NARRATIVE

**Client: Leidos, Inc.**

**Project: RVAAP FWGW**

**Report Number: 280-176674-1**

With the exceptions noted as flags or footnotes, standard analytical protocols were followed in the analysis of the samples and no problems were encountered or anomalies observed. In addition all laboratory quality control samples were within established control limits, with any exceptions noted below. Each sample was analyzed to achieve the lowest possible reporting limit within the constraints of the method. In some cases, due to interference or analytes present at high concentrations, samples were diluted. For diluted samples, the reporting limits are adjusted relative to the dilution required.

Calculations are performed before rounding to avoid round-off errors in calculated results.

All holding times were met and proper preservation noted for the methods performed on these samples, unless otherwise detailed in the individual sections below.

### **RECEIPT**

The samples were received on 5/18/2023 9:45 AM. Unless otherwise noted below, the samples arrived in good condition, and where required, properly preserved and on ice. The temperatures of the 3 coolers at receipt time were 0.2° C, 0.4° C and 3.3° C.

The results for sample FBQmw-180-230301-GW (280-176674-1) requested on a rush 5 BD TAT are reported under SDG 280-176674-2.

### **EXPLOSIVES (HPLC)**

Samples FWGmw-015-230401-GW (280-176674-2), LL3mw-241-230401-GW (280-176674-4), FBQmw-173-230401-GW-R (280-176674-7), FBQmw-175-230401-GW-R (280-176674-10) and FBQmw-175-230402-GW-R (280-176674-11) were analyzed for Explosives (HPLC) in accordance with 8330B. The samples were prepared on 05/19/2023 and analyzed on 05/20/2023, 05/21/2023 and 05/24/2023.

The matrix spike / matrix spike duplicate (MS/MSD) recoveries for preparation batch 280-613095 and analytical batch 280-613168 were outside control limits for one or more analytes. See QC Sample Results for detail. Sample matrix interference and/or non-homogeneity are suspected because the associated laboratory control sample (LCS) recovery is within acceptance limits. The following samples are impacted: FWGmw-015-230401-GW (280-176674-2), LL3mw-241-230401-GW (280-176674-4), FBQmw-173-230401-GW-R (280-176674-7), FBQmw-175-230401-GW-R (280-176674-10), FBQmw-175-230401-GW-R (280-176674-10[MS]), FBQmw-175-230401-GW-R (280-176674-10[MSD]) and FBQmw-175-230402-GW-R (280-176674-11).

The %RPD between the primary and confirmation column exceeded 40% for 2,4-Dinitrotoluene and RDX for the following sample: LL3mw-241-230401-GW (280-176674-4) in preparation batch 280-613095 and analytical batch 280-613168. The results from both columns has been qualified and reported in accordance with the laboratory's QAS.

In preparation batch 280-613095, the following samples: LL3mw-241-230401-GW (280-176674-4) and FBQmw-173-230401-GW-R (280-176674-7) were decanted prior to preparation.

No additional analytical or quality issues were noted, other than those described above or in the Definitions/Glossary page.

### **ALKALINITY**

Samples FBQmw-173-230401-GW-R (280-176674-7) and FBQmw-173-230402-GW-R (280-176674-8) were analyzed for Alkalinity in accordance with SM20 2320B. The samples were analyzed on 05/19/2023.

Alkalinity was detected in method blank MB 280-613263/33 at a level that was below one half the LOQ.

No additional analytical or quality issues were noted, other than those described above or in the Definitions/Glossary page.

### **AMMONIA**

Samples LL12mw-187-230401-GW (280-176674-3), LL12mw-185-230401-GW (280-176674-5), LL12mw-185-230402-GW (280-176674-6) and LL12mw-245-230401-GW (280-176674-9) were analyzed for ammonia in accordance with EPA Method 350.1. The samples were analyzed on 05/23/2023.

Sample LL12mw-187-230401-GW (280-176674-3)[200X] required dilution prior to analysis. The reporting limits have been adjusted accordingly.

No additional analytical or quality issues were noted, other than those described above or in the Definitions/Glossary page.

### **SULFIDE**

Samples FBQmw-173-230401-GW-R (280-176674-7) and FBQmw-173-230402-GW-R (280-176674-8) were analyzed for sulfide in accordance with EPA SW-846 Method 9034. The samples were prepared and analyzed on 05/19/2023.

No analytical or quality issues were noted, other than those described above or in the Definitions/Glossary page.

**ANIONS, ION CHROMATOGRAPHY**

Samples LL12mw-187-230401-GW (280-176674-3), LL12mw-185-230401-GW (280-176674-5), LL12mw-185-230402-GW (280-176674-6), FBQmw-173-230401-GW-R (280-176674-7), FBQmw-173-230402-GW-R (280-176674-8) and LL12mw-245-230401-GW (280-176674-9) were analyzed for Anions, Ion Chromatography in accordance with 9056A (48 Hours). The samples were analyzed on 05/18/2023, 05/19/2023 and 05/22/2023.

Samples LL12mw-187-230401-GW (280-176674-3)[2000X], LL12mw-185-230401-GW (280-176674-5)[20X] and LL12mw-185-230402-GW (280-176674-6)[20X] required dilution prior to analysis. The reporting limits have been adjusted accordingly.

Reanalysis of the following samples was performed outside of the analytical holding time due to the sample requiring a higher dilution for nitrate: LL12mw-187-230401-GW (280-176674-3).

No additional analytical or quality issues were noted, other than those described above or in the Definitions/Glossary page.

**ANIONS (28 DAYS)**

Samples FBQmw-173-230401-GW-R (280-176674-7) and FBQmw-173-230402-GW-R (280-176674-8) were analyzed for anions (28 days) in accordance with 9056A (28 Days). The samples were analyzed on 05/18/2023.

No analytical or quality issues were noted, other than those described above or in the Definitions/Glossary page.

**TOTAL ORGANIC CARBON**

Sample FBQmw-173-230401-GW-R (280-176674-7) was analyzed for total organic carbon in accordance with EPA SW-846 Method 9060. The samples were analyzed on 05/19/2023.

No analytical or quality issues were noted, other than those described above or in the Definitions/Glossary page.

# Detection Summary

Client: Leidos, Inc.  
Project/Site: RVAAP FWGW

Job ID: 280-176674-1

## Client Sample ID: FWGmw-015-230401-GW

## Lab Sample ID: 280-176674-2

No Detections.

## Client Sample ID: LL12mw-187-230401-GW

## Lab Sample ID: 280-176674-3

| Analyte      | Result | Qualifier | LOQ  | LOD | DL  | Unit | Dil Fac | D | Method | Prep Type |
|--------------|--------|-----------|------|-----|-----|------|---------|---|--------|-----------|
| Ammonia as N | 730    |           | 20   | 10  | 5.8 | mg/L | 200     |   | 350.1  | Total/NA  |
| Nitrate as N | 1500   | H D       | 1000 | 400 | 180 | mg/L | 2000    |   | 9056   | Total/NA  |

## Client Sample ID: LL3mw-241-230401-GW

## Lab Sample ID: 280-176674-4

| Analyte                    | Result | Qualifier | LOQ  | LOD   | DL    | Unit | Dil Fac | D | Method | Prep Type |
|----------------------------|--------|-----------|------|-------|-------|------|---------|---|--------|-----------|
| 1,3,5-Trinitrobenzene      | 4.6    |           | 0.21 | 0.20  | 0.086 | ug/L | 1       |   | 8330B  | Total/NA  |
| 2,4,6-Trinitrotoluene      | 2.0    |           | 0.11 | 0.10  | 0.046 | ug/L | 1       |   | 8330B  | Total/NA  |
| 2,4-Dinitrotoluene         | 0.12   | J1        | 0.10 | 0.082 | 0.028 | ug/L | 1       |   | 8330B  | Total/NA  |
| 2,4-Dinitrotoluene         | 0.045  | J J1      | 0.10 | 0.082 | 0.028 | ug/L | 1       |   | 8330B  | Total/NA  |
| 3,5-Dinitroaniline         | 0.36   | J         | 0.41 | 0.31  | 0.13  | ug/L | 1       |   | 8330B  | Total/NA  |
| 4-Amino-2,6-dinitrotoluene | 1.7    |           | 0.15 | 0.12  | 0.059 | ug/L | 1       |   | 8330B  | Total/NA  |
| HMX                        | 0.14   | J M       | 0.21 | 0.20  | 0.090 | ug/L | 1       |   | 8330B  | Total/NA  |
| RDX                        | 1.0    | J1        | 0.21 | 0.20  | 0.053 | ug/L | 1       |   | 8330B  | Total/NA  |
| RDX                        | 0.49   | J1        | 0.21 | 0.20  | 0.053 | ug/L | 1       |   | 8330B  | Total/NA  |

## Client Sample ID: LL12mw-185-230401-GW

## Lab Sample ID: 280-176674-5

| Analyte      | Result | Qualifier | LOQ | LOD | DL  | Unit | Dil Fac | D | Method | Prep Type |
|--------------|--------|-----------|-----|-----|-----|------|---------|---|--------|-----------|
| Nitrate as N | 67     | D M       | 10  | 4.0 | 1.8 | mg/L | 20      |   | 9056   | Total/NA  |

## Client Sample ID: LL12mw-185-230402-GW

## Lab Sample ID: 280-176674-6

| Analyte      | Result | Qualifier | LOQ | LOD | DL  | Unit | Dil Fac | D | Method | Prep Type |
|--------------|--------|-----------|-----|-----|-----|------|---------|---|--------|-----------|
| Nitrate as N | 74     | D M       | 10  | 4.0 | 1.8 | mg/L | 20      |   | 9056   | Total/NA  |

## Client Sample ID: FBQmw-173-230401-GW-R

## Lab Sample ID: 280-176674-7

| Analyte                        | Result | Qualifier | LOQ  | LOD  | DL    | Unit | Dil Fac | D | Method   | Prep Type |
|--------------------------------|--------|-----------|------|------|-------|------|---------|---|----------|-----------|
| 2,4,6-Trinitrotoluene          | 0.11   |           | 0.11 | 0.10 | 0.047 | ug/L | 1       |   | 8330B    | Total/NA  |
| 2-Amino-4,6-dinitrotoluene     | 1.2    |           | 0.11 | 0.10 | 0.052 | ug/L | 1       |   | 8330B    | Total/NA  |
| 3,5-Dinitroaniline             | 0.15   | J         | 0.41 | 0.31 | 0.14  | ug/L | 1       |   | 8330B    | Total/NA  |
| 4-Amino-2,6-dinitrotoluene     | 1.5    |           | 0.16 | 0.12 | 0.060 | ug/L | 1       |   | 8330B    | Total/NA  |
| Nitrate as N                   | 0.20   | J M       | 0.50 | 0.20 | 0.090 | mg/L | 1       |   | 9056     | Total/NA  |
| Sulfate                        | 42     | M         | 5.0  | 2.5  | 1.0   | mg/L | 1       |   | 9056A    | Total/NA  |
| Total Organic Carbon - Average | 0.85   | J         | 1.0  | 0.80 | 0.35  | mg/L | 1       | ☼ | 9060     | Total/NA  |
| Alkalinity                     | 18     |           | 10   | 6.4  | 3.1   | mg/L | 1       |   | SM 2320B | Total/NA  |

## Client Sample ID: FBQmw-173-230402-GW-R

## Lab Sample ID: 280-176674-8

| Analyte      | Result | Qualifier | LOQ  | LOD  | DL    | Unit | Dil Fac | D | Method   | Prep Type |
|--------------|--------|-----------|------|------|-------|------|---------|---|----------|-----------|
| Nitrate as N | 0.21   | J M       | 0.50 | 0.20 | 0.090 | mg/L | 1       |   | 9056     | Total/NA  |
| Sulfate      | 45     |           | 5.0  | 2.5  | 1.0   | mg/L | 1       |   | 9056A    | Total/NA  |
| Alkalinity   | 18     |           | 10   | 6.4  | 3.1   | mg/L | 1       |   | SM 2320B | Total/NA  |

## Client Sample ID: LL12mw-245-230401-GW

## Lab Sample ID: 280-176674-9

No Detections.

This Detection Summary does not include radiochemical test results.

Eurofins Denver

## Detection Summary

Client: Leidos, Inc.  
Project/Site: RVAAP FWGW

Job ID: 280-176674-1

**Client Sample ID: FBQmw-175-230401-GW-R**

**Lab Sample ID: 280-176674-10**

☐ No Detections.

**Client Sample ID: FBQmw-175-230402-GW-R**

**Lab Sample ID: 280-176674-11**

☐ No Detections.

This Detection Summary does not include radiochemical test results.

# Client Sample Results

Client: Leidos, Inc.  
Project/Site: RVAAP FWGW

Job ID: 280-176674-1

## Method: EPA 8330B - Nitroaromatics and Nitramines (HPLC)

Client Sample ID: FWGmw-015-230401-GW

Date Collected: 05/17/23 14:06

Date Received: 05/18/23 09:45

Lab Sample ID: 280-176674-2

Matrix: Water

| Analyte                    | Result | Qualifier | LOQ  | LOD   | DL    | Unit | D | Analyzed       | Dil Fac |
|----------------------------|--------|-----------|------|-------|-------|------|---|----------------|---------|
| 1,3,5-Trinitrobenzene      | 0.21   | U         | 0.22 | 0.21  | 0.087 | ug/L |   | 05/20/23 02:53 | 1       |
| 1,3-Dinitrobenzene         | 0.10   | U         | 0.11 | 0.10  | 0.038 | ug/L |   | 05/20/23 02:53 | 1       |
| 2,4,6-Trinitrotoluene      | 0.10   | U         | 0.11 | 0.10  | 0.047 | ug/L |   | 05/20/23 02:53 | 1       |
| 2,4-Dinitrotoluene         | 0.083  | U         | 0.10 | 0.083 | 0.028 | ug/L |   | 05/20/23 02:53 | 1       |
| 2,6-Dinitrotoluene         | 0.083  | U         | 0.10 | 0.083 | 0.042 | ug/L |   | 05/20/23 02:53 | 1       |
| 2-Amino-4,6-dinitrotoluene | 0.10   | U         | 0.11 | 0.10  | 0.053 | ug/L |   | 05/20/23 02:53 | 1       |
| 2-Nitrotoluene             | 0.21   | U         | 0.22 | 0.21  | 0.089 | ug/L |   | 05/20/23 02:53 | 1       |
| 3-Nitrotoluene             | 0.36   | U         | 0.42 | 0.36  | 0.20  | ug/L |   | 05/20/23 02:53 | 1       |
| 4-Amino-2,6-dinitrotoluene | 0.12   | U         | 0.16 | 0.12  | 0.060 | ug/L |   | 05/20/23 02:53 | 1       |
| 4-Nitrotoluene             | 0.42   | U         | 0.43 | 0.42  | 0.10  | ug/L |   | 05/20/23 02:53 | 1       |
| HMX                        | 0.21   | U M       | 0.22 | 0.21  | 0.091 | ug/L |   | 05/20/23 02:53 | 1       |
| Nitrobenzene               | 0.21   | U         | 0.22 | 0.21  | 0.094 | ug/L |   | 05/20/23 02:53 | 1       |
| Nitroglycerin              | 2.1    | U         | 2.2  | 2.1   | 0.96  | ug/L |   | 05/20/23 02:53 | 1       |
| PETN                       | 1.0    | U         | 1.1  | 1.0   | 0.46  | ug/L |   | 05/20/23 02:53 | 1       |
| RDX                        | 0.21   | U         | 0.22 | 0.21  | 0.053 | ug/L |   | 05/20/23 02:53 | 1       |
| Tetryl                     | 0.10   | U         | 0.11 | 0.10  | 0.033 | ug/L |   | 05/20/23 02:53 | 1       |

| Surrogate          | %Recovery | Qualifier | Limits   | Prepared       | Analyzed       | Dil Fac |
|--------------------|-----------|-----------|----------|----------------|----------------|---------|
| 1,2-Dinitrobenzene | 85        | M         | 83 - 119 | 05/19/23 13:39 | 05/20/23 02:53 | 1       |

Client Sample ID: LL3mw-241-230401-GW

Date Collected: 05/17/23 15:00

Date Received: 05/18/23 09:45

Lab Sample ID: 280-176674-4

Matrix: Water

| Analyte                    | Result | Qualifier | LOQ  | LOD   | DL    | Unit | D | Analyzed       | Dil Fac |
|----------------------------|--------|-----------|------|-------|-------|------|---|----------------|---------|
| 1,3,5-Trinitrobenzene      | 4.6    |           | 0.21 | 0.20  | 0.086 | ug/L |   | 05/20/23 04:02 | 1       |
| 1,3-Dinitrobenzene         | 0.10   | U         | 0.11 | 0.10  | 0.038 | ug/L |   | 05/24/23 19:48 | 1       |
| 2,4,6-Trinitrotoluene      | 2.0    |           | 0.11 | 0.10  | 0.046 | ug/L |   | 05/20/23 04:02 | 1       |
| 2,4-diamino-6-nitrotoluene | 0.92   | U         | 1.0  | 0.92  | 0.44  | ug/L |   | 05/20/23 04:02 | 1       |
| 2,4-Dinitrotoluene         | 0.12   | J1        | 0.10 | 0.082 | 0.028 | ug/L |   | 05/20/23 04:02 | 1       |
| 2,4-Dinitrotoluene         | 0.045  | J J1      | 0.10 | 0.082 | 0.028 | ug/L |   | 05/24/23 19:48 | 1       |
| 2,6-diamino-4-nitrotoluene | 0.92   | U M       | 1.0  | 0.92  | 0.23  | ug/L |   | 05/20/23 04:02 | 1       |
| 2,6-Dinitrotoluene         | 0.082  | U         | 0.10 | 0.082 | 0.041 | ug/L |   | 05/20/23 04:02 | 1       |
| 2-Amino-4,6-dinitrotoluene | 0.10   | U         | 0.11 | 0.10  | 0.052 | ug/L |   | 05/24/23 19:48 | 1       |
| 2-Nitrotoluene             | 0.20   | U         | 0.21 | 0.20  | 0.088 | ug/L |   | 05/20/23 04:02 | 1       |
| 3,5-Dinitroaniline         | 0.36   | J         | 0.41 | 0.31  | 0.13  | ug/L |   | 05/20/23 04:02 | 1       |
| 3-Nitrotoluene             | 0.36   | U         | 0.41 | 0.36  | 0.20  | ug/L |   | 05/20/23 04:02 | 1       |
| 4-Amino-2,6-dinitrotoluene | 1.7    |           | 0.15 | 0.12  | 0.059 | ug/L |   | 05/20/23 04:02 | 1       |
| 4-Nitrotoluene             | 0.41   | U M       | 0.42 | 0.41  | 0.10  | ug/L |   | 05/20/23 04:02 | 1       |
| DNX                        | 0.26   | U M       | 0.51 | 0.26  | 0.10  | ug/L |   | 05/20/23 04:02 | 1       |
| HMX                        | 0.14   | J M       | 0.21 | 0.20  | 0.090 | ug/L |   | 05/20/23 04:02 | 1       |
| MNX                        | 0.30   | U M       | 0.51 | 0.30  | 0.095 | ug/L |   | 05/24/23 19:48 | 1       |
| Nitrobenzene               | 0.20   | U M       | 0.21 | 0.20  | 0.093 | ug/L |   | 05/20/23 04:02 | 1       |
| Nitroglycerin              | 2.0    | U         | 2.1  | 2.0   | 0.94  | ug/L |   | 05/20/23 04:02 | 1       |
| PETN                       | 1.0    | U         | 1.1  | 1.0   | 0.46  | ug/L |   | 05/20/23 04:02 | 1       |
| RDX                        | 1.0    | J1        | 0.21 | 0.20  | 0.053 | ug/L |   | 05/20/23 04:02 | 1       |
| RDX                        | 0.49   | J1        | 0.21 | 0.20  | 0.053 | ug/L |   | 05/24/23 19:48 | 1       |
| Tetryl                     | 0.10   | U M       | 0.11 | 0.10  | 0.033 | ug/L |   | 05/20/23 04:02 | 1       |
| TNX                        | 0.26   | U         | 0.51 | 0.26  | 0.081 | ug/L |   | 05/20/23 04:02 | 1       |

Eurofins Denver



# Client Sample Results

Client: Leidos, Inc.  
Project/Site: RVAAP FWGW

Job ID: 280-176674-1

## Method: EPA 8330B - Nitroaromatics and Nitramines (HPLC) (Continued)

| Surrogate          | %Recovery | Qualifier | Limits   | Prepared       | Analyzed       | Dil Fac |
|--------------------|-----------|-----------|----------|----------------|----------------|---------|
| 1,2-Dinitrobenzene | 93        |           | 83 - 119 | 05/19/23 13:39 | 05/20/23 04:02 | 1       |
| 1,2-Dinitrobenzene | 88        |           | 83 - 119 | 05/19/23 13:39 | 05/24/23 19:48 | 1       |

Client Sample ID: FBQmw-173-230401-GW-R

Date Collected: 05/17/23 11:47

Date Received: 05/18/23 09:45

Lab Sample ID: 280-176674-7

Matrix: Water

| Analyte                    | Result | Qualifier | LOQ  | LOD   | DL    | Unit | D | Analyzed       | Dil Fac |
|----------------------------|--------|-----------|------|-------|-------|------|---|----------------|---------|
| 1,3,5-Trinitrobenzene      | 0.21   | U         | 0.22 | 0.21  | 0.087 | ug/L |   | 05/20/23 04:25 | 1       |
| 1,3-Dinitrobenzene         | 0.10   | U         | 0.11 | 0.10  | 0.038 | ug/L |   | 05/21/23 00:31 | 1       |
| 2,4,6-Trinitrotoluene      | 0.11   |           | 0.11 | 0.10  | 0.047 | ug/L |   | 05/20/23 04:25 | 1       |
| 2,4-diamino-6-nitrotoluene | 0.93   | U         | 1.0  | 0.93  | 0.45  | ug/L |   | 05/20/23 04:25 | 1       |
| 2,4-Dinitrotoluene         | 0.083  | U         | 0.10 | 0.083 | 0.028 | ug/L |   | 05/20/23 04:25 | 1       |
| 2,6-diamino-4-nitrotoluene | 0.93   | U         | 1.0  | 0.93  | 0.23  | ug/L |   | 05/20/23 04:25 | 1       |
| 2,6-Dinitrotoluene         | 0.083  | U         | 0.10 | 0.083 | 0.041 | ug/L |   | 05/20/23 04:25 | 1       |
| 2-Amino-4,6-dinitrotoluene | 1.2    |           | 0.11 | 0.10  | 0.052 | ug/L |   | 05/20/23 04:25 | 1       |
| 2-Nitrotoluene             | 0.21   | U         | 0.22 | 0.21  | 0.088 | ug/L |   | 05/20/23 04:25 | 1       |
| 3,5-Dinitroaniline         | 0.15   | J         | 0.41 | 0.31  | 0.14  | ug/L |   | 05/20/23 04:25 | 1       |
| 3-Nitrotoluene             | 0.36   | U         | 0.41 | 0.36  | 0.20  | ug/L |   | 05/20/23 04:25 | 1       |
| 4-Amino-2,6-dinitrotoluene | 1.5    |           | 0.16 | 0.12  | 0.060 | ug/L |   | 05/20/23 04:25 | 1       |
| 4-Nitrotoluene             | 0.41   | U         | 0.42 | 0.41  | 0.10  | ug/L |   | 05/20/23 04:25 | 1       |
| DNX                        | 0.26   | U         | 0.52 | 0.26  | 0.10  | ug/L |   | 05/20/23 04:25 | 1       |
| HMX                        | 0.21   | U         | 0.22 | 0.21  | 0.091 | ug/L |   | 05/20/23 04:25 | 1       |
| MNX                        | 0.30   | U         | 0.52 | 0.30  | 0.096 | ug/L |   | 05/20/23 04:25 | 1       |
| Nitrobenzene               | 0.21   | U         | 0.22 | 0.21  | 0.094 | ug/L |   | 05/20/23 04:25 | 1       |
| Nitroglycerin              | 2.1    | U         | 2.2  | 2.1   | 0.95  | ug/L |   | 05/20/23 04:25 | 1       |
| PETN                       | 1.0    | U         | 1.1  | 1.0   | 0.46  | ug/L |   | 05/20/23 04:25 | 1       |
| RDX                        | 0.21   | U M       | 0.22 | 0.21  | 0.053 | ug/L |   | 05/20/23 04:25 | 1       |
| Tetryl                     | 0.10   | U         | 0.11 | 0.10  | 0.033 | ug/L |   | 05/20/23 04:25 | 1       |
| TNX                        | 0.26   | U M       | 0.52 | 0.26  | 0.082 | ug/L |   | 05/20/23 04:25 | 1       |

| Surrogate          | %Recovery | Qualifier | Limits   | Prepared       | Analyzed       | Dil Fac |
|--------------------|-----------|-----------|----------|----------------|----------------|---------|
| 1,2-Dinitrobenzene | 91        | M         | 83 - 119 | 05/19/23 13:39 | 05/20/23 04:25 | 1       |
| 1,2-Dinitrobenzene | 85        |           | 83 - 119 | 05/19/23 13:39 | 05/21/23 00:31 | 1       |

Client Sample ID: FBQmw-175-230401-GW-R

Date Collected: 05/17/23 08:58

Date Received: 05/18/23 09:45

Lab Sample ID: 280-176674-10

Matrix: Water

| Analyte                    | Result | Qualifier | LOQ  | LOD   | DL    | Unit | D | Analyzed       | Dil Fac |
|----------------------------|--------|-----------|------|-------|-------|------|---|----------------|---------|
| 1,3,5-Trinitrobenzene      | 0.20   | U J1      | 0.21 | 0.20  | 0.085 | ug/L |   | 05/20/23 04:48 | 1       |
| 1,3-Dinitrobenzene         | 0.10   | U J1      | 0.11 | 0.10  | 0.037 | ug/L |   | 05/20/23 04:48 | 1       |
| 2,4,6-Trinitrotoluene      | 0.10   | U J1      | 0.11 | 0.10  | 0.045 | ug/L |   | 05/20/23 04:48 | 1       |
| 2,4-diamino-6-nitrotoluene | 0.90   | U J1      | 1.0  | 0.90  | 0.44  | ug/L |   | 05/20/23 04:48 | 1       |
| 2,4-Dinitrotoluene         | 0.080  | U J1      | 0.10 | 0.080 | 0.028 | ug/L |   | 05/20/23 04:48 | 1       |
| 2,6-diamino-4-nitrotoluene | 0.90   | U J1      | 1.0  | 0.90  | 0.22  | ug/L |   | 05/20/23 04:48 | 1       |
| 2,6-Dinitrotoluene         | 0.080  | U J1      | 0.10 | 0.080 | 0.040 | ug/L |   | 05/20/23 04:48 | 1       |
| 2-Amino-4,6-dinitrotoluene | 0.10   | U J1      | 0.11 | 0.10  | 0.051 | ug/L |   | 05/20/23 04:48 | 1       |
| 2-Nitrotoluene             | 0.20   | U J1      | 0.21 | 0.20  | 0.086 | ug/L |   | 05/20/23 04:48 | 1       |
| 3,5-Dinitroaniline         | 0.30   | U J1      | 0.40 | 0.30  | 0.13  | ug/L |   | 05/20/23 04:48 | 1       |
| 3-Nitrotoluene             | 0.35   | U J1      | 0.40 | 0.35  | 0.20  | ug/L |   | 05/20/23 04:48 | 1       |
| 4-Amino-2,6-dinitrotoluene | 0.12   | U J1      | 0.15 | 0.12  | 0.058 | ug/L |   | 05/20/23 04:48 | 1       |
| 4-Nitrotoluene             | 0.40   | U J1      | 0.41 | 0.40  | 0.10  | ug/L |   | 05/20/23 04:48 | 1       |
| DNX                        | 0.25   | U J1      | 0.50 | 0.25  | 0.098 | ug/L |   | 05/20/23 04:48 | 1       |
| HMX                        | 0.20   | U         | 0.21 | 0.20  | 0.088 | ug/L |   | 05/20/23 04:48 | 1       |

Eurofins Denver

# Client Sample Results

Client: Leidos, Inc.  
Project/Site: RVAAP FWGW

Job ID: 280-176674-1

## Method: EPA 8330B - Nitroaromatics and Nitramines (HPLC) (Continued)

Client Sample ID: FBQmw-175-230401-GW-R

Date Collected: 05/17/23 08:58

Date Received: 05/18/23 09:45

Lab Sample ID: 280-176674-10

Matrix: Water

| Analyte       | Result | Qualifier | LOQ  | LOD  | DL    | Unit | D | Analyzed       | Dil Fac |
|---------------|--------|-----------|------|------|-------|------|---|----------------|---------|
| MNX           | 0.29   | U         | 0.50 | 0.29 | 0.093 | ug/L |   | 05/20/23 04:48 | 1       |
| Nitrobenzene  | 0.20   | U J1      | 0.21 | 0.20 | 0.091 | ug/L |   | 05/20/23 04:48 | 1       |
| Nitroglycerin | 2.0    | U         | 2.1  | 2.0  | 0.93  | ug/L |   | 05/20/23 04:48 | 1       |
| PETN          | 1.0    | U J1      | 1.1  | 1.0  | 0.45  | ug/L |   | 05/20/23 04:48 | 1       |
| RDX           | 0.20   | U M J1    | 0.21 | 0.20 | 0.052 | ug/L |   | 05/20/23 04:48 | 1       |
| Tetryl        | 0.10   | U J1      | 0.11 | 0.10 | 0.032 | ug/L |   | 05/20/23 04:48 | 1       |
| TNX           | 0.25   | U M       | 0.50 | 0.25 | 0.080 | ug/L |   | 05/20/23 04:48 | 1       |

| Surrogate          | %Recovery | Qualifier | Limits   | Prepared       | Analyzed       | Dil Fac |
|--------------------|-----------|-----------|----------|----------------|----------------|---------|
| 1,2-Dinitrobenzene | 95        | M         | 83 - 119 | 05/19/23 13:39 | 05/20/23 04:48 | 1       |

Client Sample ID: FBQmw-175-230402-GW-R

Date Collected: 05/17/23 08:58

Date Received: 05/18/23 09:45

Lab Sample ID: 280-176674-11

Matrix: Water

| Analyte                    | Result | Qualifier | LOQ  | LOD   | DL    | Unit | D | Analyzed       | Dil Fac |
|----------------------------|--------|-----------|------|-------|-------|------|---|----------------|---------|
| 1,3,5-Trinitrobenzene      | 0.20   | U         | 0.21 | 0.20  | 0.085 | ug/L |   | 05/20/23 06:43 | 1       |
| 1,3-Dinitrobenzene         | 0.10   | U         | 0.11 | 0.10  | 0.037 | ug/L |   | 05/20/23 06:43 | 1       |
| 2,4,6-Trinitrotoluene      | 0.10   | U         | 0.11 | 0.10  | 0.045 | ug/L |   | 05/20/23 06:43 | 1       |
| 2,4-diamino-6-nitrotoluene | 0.91   | U         | 1.0  | 0.91  | 0.44  | ug/L |   | 05/20/23 06:43 | 1       |
| 2,4-Dinitrotoluene         | 0.080  | U         | 0.10 | 0.080 | 0.028 | ug/L |   | 05/20/23 06:43 | 1       |
| 2,6-diamino-4-nitrotoluene | 0.91   | U         | 1.0  | 0.91  | 0.22  | ug/L |   | 05/20/23 06:43 | 1       |
| 2,6-Dinitrotoluene         | 0.080  | U         | 0.10 | 0.080 | 0.040 | ug/L |   | 05/20/23 06:43 | 1       |
| 2-Amino-4,6-dinitrotoluene | 0.10   | U         | 0.11 | 0.10  | 0.051 | ug/L |   | 05/20/23 06:43 | 1       |
| 2-Nitrotoluene             | 0.20   | U M       | 0.21 | 0.20  | 0.086 | ug/L |   | 05/20/23 06:43 | 1       |
| 3,5-Dinitroaniline         | 0.30   | U         | 0.40 | 0.30  | 0.13  | ug/L |   | 05/20/23 06:43 | 1       |
| 3-Nitrotoluene             | 0.35   | U         | 0.40 | 0.35  | 0.20  | ug/L |   | 05/20/23 06:43 | 1       |
| 4-Amino-2,6-dinitrotoluene | 0.12   | U         | 0.15 | 0.12  | 0.058 | ug/L |   | 05/20/23 06:43 | 1       |
| 4-Nitrotoluene             | 0.40   | U         | 0.41 | 0.40  | 0.10  | ug/L |   | 05/20/23 06:43 | 1       |
| DNX                        | 0.25   | U         | 0.50 | 0.25  | 0.098 | ug/L |   | 05/20/23 06:43 | 1       |
| HMX                        | 0.20   | U         | 0.21 | 0.20  | 0.088 | ug/L |   | 05/20/23 06:43 | 1       |
| MNX                        | 0.29   | U         | 0.50 | 0.29  | 0.093 | ug/L |   | 05/20/23 06:43 | 1       |
| Nitrobenzene               | 0.20   | U         | 0.21 | 0.20  | 0.092 | ug/L |   | 05/20/23 06:43 | 1       |
| Nitroglycerin              | 2.0    | U         | 2.1  | 2.0   | 0.93  | ug/L |   | 05/20/23 06:43 | 1       |
| PETN                       | 1.0    | U         | 1.1  | 1.0   | 0.45  | ug/L |   | 05/20/23 06:43 | 1       |
| RDX                        | 0.20   | U M       | 0.21 | 0.20  | 0.052 | ug/L |   | 05/20/23 06:43 | 1       |
| Tetryl                     | 0.10   | U         | 0.11 | 0.10  | 0.032 | ug/L |   | 05/20/23 06:43 | 1       |
| TNX                        | 0.25   | U M       | 0.50 | 0.25  | 0.080 | ug/L |   | 05/20/23 06:43 | 1       |

| Surrogate          | %Recovery | Qualifier | Limits   | Prepared       | Analyzed       | Dil Fac |
|--------------------|-----------|-----------|----------|----------------|----------------|---------|
| 1,2-Dinitrobenzene | 93        | M         | 83 - 119 | 05/19/23 13:39 | 05/20/23 06:43 | 1       |

## General Chemistry

Client Sample ID: LL12mw-187-230401-GW

Date Collected: 05/17/23 13:40

Date Received: 05/18/23 09:45

Lab Sample ID: 280-176674-3

Matrix: Water

| Analyte                   | Result | Qualifier | LOQ  | LOD | DL  | Unit | D | Analyzed       | Dil Fac |
|---------------------------|--------|-----------|------|-----|-----|------|---|----------------|---------|
| Ammonia as N (EPA 350.1)  | 730    |           | 20   | 10  | 5.8 | mg/L |   | 05/23/23 17:04 | 200     |
| Nitrate as N (SW846 9056) | 1500   | H D       | 1000 | 400 | 180 | mg/L |   | 05/22/23 18:59 | 2000    |

Eurofins Denver

# Client Sample Results

Client: Leidos, Inc.  
Project/Site: RVAAP FWGW

Job ID: 280-176674-1

## General Chemistry

Client Sample ID: LL12mw-185-230401-GW

Date Collected: 05/17/23 09:45

Date Received: 05/18/23 09:45

Lab Sample ID: 280-176674-5

Matrix: Water

| Analyte                   | Result | Qualifier | LOQ  | LOD   | DL    | Unit | D | Analyzed       | Dil Fac |
|---------------------------|--------|-----------|------|-------|-------|------|---|----------------|---------|
| Ammonia as N (EPA 350.1)  | 0.050  | U         | 0.10 | 0.050 | 0.029 | mg/L |   | 05/23/23 13:01 | 1       |
| Nitrate as N (SW846 9056) | 67     | D M       | 10   | 4.0   | 1.8   | mg/L |   | 05/18/23 16:26 | 20      |

Client Sample ID: LL12mw-185-230402-GW

Date Collected: 05/17/23 09:45

Date Received: 05/18/23 09:45

Lab Sample ID: 280-176674-6

Matrix: Water

| Analyte                   | Result | Qualifier | LOQ  | LOD   | DL    | Unit | D | Analyzed       | Dil Fac |
|---------------------------|--------|-----------|------|-------|-------|------|---|----------------|---------|
| Ammonia as N (EPA 350.1)  | 0.050  | U         | 0.10 | 0.050 | 0.029 | mg/L |   | 05/23/23 13:25 | 1       |
| Nitrate as N (SW846 9056) | 74     | D M       | 10   | 4.0   | 1.8   | mg/L |   | 05/18/23 16:41 | 20      |

Client Sample ID: FBQmw-173-230401-GW-R

Date Collected: 05/17/23 11:47

Date Received: 05/18/23 09:45

Lab Sample ID: 280-176674-7

Matrix: Water

| Analyte                                     | Result | Qualifier | LOQ  | LOD  | DL    | Unit | D | Analyzed       | Dil Fac |
|---------------------------------------------|--------|-----------|------|------|-------|------|---|----------------|---------|
| Sulfide (SW846 9034)                        | 2.0    | U         | 4.0  | 2.0  | 1.6   | mg/L |   | 05/19/23 17:27 | 1       |
| Nitrate as N (SW846 9056)                   | 0.20   | J M       | 0.50 | 0.20 | 0.090 | mg/L |   | 05/18/23 19:26 | 1       |
| Nitrite as N (SW846 9056)                   | 0.10   | U         | 0.50 | 0.10 | 0.049 | mg/L |   | 05/18/23 19:26 | 1       |
| Sulfate (SW846 9056A)                       | 42     | M         | 5.0  | 2.5  | 1.0   | mg/L |   | 05/18/23 19:26 | 1       |
| Total Organic Carbon - Average (SW846 9060) | 0.85   | J         | 1.0  | 0.80 | 0.35  | mg/L | ✱ | 05/19/23 19:17 | 1       |
| Alkalinity (SM 2320B)                       | 18     |           | 10   | 6.4  | 3.1   | mg/L |   | 05/19/23 17:16 | 1       |

Client Sample ID: FBQmw-173-230402-GW-R

Date Collected: 05/17/23 11:47

Date Received: 05/18/23 09:45

Lab Sample ID: 280-176674-8

Matrix: Water

| Analyte                   | Result | Qualifier | LOQ  | LOD  | DL    | Unit | D | Analyzed       | Dil Fac |
|---------------------------|--------|-----------|------|------|-------|------|---|----------------|---------|
| Sulfide (SW846 9034)      | 2.0    | U         | 4.0  | 2.0  | 1.6   | mg/L |   | 05/19/23 17:27 | 1       |
| Nitrate as N (SW846 9056) | 0.21   | J M       | 0.50 | 0.20 | 0.090 | mg/L |   | 05/18/23 19:41 | 1       |
| Nitrite as N (SW846 9056) | 0.10   | U         | 0.50 | 0.10 | 0.049 | mg/L |   | 05/18/23 19:41 | 1       |
| Sulfate (SW846 9056A)     | 45     |           | 5.0  | 2.5  | 1.0   | mg/L |   | 05/18/23 19:41 | 1       |
| Alkalinity (SM 2320B)     | 18     |           | 10   | 6.4  | 3.1   | mg/L |   | 05/19/23 17:22 | 1       |

Client Sample ID: LL12mw-245-230401-GW

Date Collected: 05/17/23 14:26

Date Received: 05/18/23 09:45

Lab Sample ID: 280-176674-9

Matrix: Water

| Analyte                   | Result | Qualifier | LOQ  | LOD   | DL    | Unit | D | Analyzed       | Dil Fac |
|---------------------------|--------|-----------|------|-------|-------|------|---|----------------|---------|
| Ammonia as N (EPA 350.1)  | 0.050  | U         | 0.10 | 0.050 | 0.029 | mg/L |   | 05/23/23 13:28 | 1       |
| Nitrate as N (SW846 9056) | 0.20   | U M       | 0.50 | 0.20  | 0.090 | mg/L |   | 05/19/23 06:42 | 1       |

## Default Detection Limits

Client: Leidos, Inc.  
Project/Site: RVAAP FWGW

Job ID: 280-176674-1

### Method: 8330B - Nitroaromatics and Nitramines (HPLC)

Prep: 3535

| Analyte                    | LOQ  | DL    | Units |
|----------------------------|------|-------|-------|
| 1,3,5-Trinitrobenzene      | 0.21 | 0.084 | ug/L  |
| 1,3-Dinitrobenzene         | 0.11 | 0.037 | ug/L  |
| 2,4,6-Trinitrotoluene      | 0.11 | 0.045 | ug/L  |
| 2,4-diamino-6-nitrotoluene | 1.0  | 0.43  | ug/L  |
| 2,4-Dinitrotoluene         | 0.10 | 0.027 | ug/L  |
| 2,6-diamino-4-nitrotoluene | 1.0  | 0.22  | ug/L  |
| 2,6-Dinitrotoluene         | 0.10 | 0.040 | ug/L  |
| 2-Amino-4,6-dinitrotoluene | 0.11 | 0.051 | ug/L  |
| 2-Nitrotoluene             | 0.21 | 0.086 | ug/L  |
| 3,5-Dinitroaniline         | 0.40 | 0.13  | ug/L  |
| 3-Nitrotoluene             | 0.40 | 0.20  | ug/L  |
| 4-Amino-2,6-dinitrotoluene | 0.15 | 0.058 | ug/L  |
| 4-Nitrotoluene             | 0.41 | 0.10  | ug/L  |
| DNX                        | 0.50 | 0.097 | ug/L  |
| HMX                        | 0.21 | 0.088 | ug/L  |
| MNX                        | 0.50 | 0.093 | ug/L  |
| Nitrobenzene               | 0.21 | 0.091 | ug/L  |
| Nitroglycerin              | 2.1  | 0.92  | ug/L  |
| PETN                       | 1.1  | 0.45  | ug/L  |
| RDX                        | 0.21 | 0.052 | ug/L  |
| Tetryl                     | 0.11 | 0.032 | ug/L  |
| TNX                        | 0.50 | 0.079 | ug/L  |

### General Chemistry

| Analyte                        | LOQ  | DL    | Units |
|--------------------------------|------|-------|-------|
| Ammonia as N                   | 0.10 | 0.029 | mg/L  |
| Nitrate as N                   | 0.50 | 0.090 | mg/L  |
| Nitrite as N                   | 0.50 | 0.049 | mg/L  |
| Sulfate                        | 5.0  | 1.0   | mg/L  |
| Total Organic Carbon - Average | 1.0  | 0.35  | mg/L  |
| Alkalinity                     | 10   | 3.1   | mg/L  |

### General Chemistry

Prep: 9030B

| Analyte | LOQ | DL  | Units |
|---------|-----|-----|-------|
| Sulfide | 4.0 | 1.6 | mg/L  |

# Surrogate Summary

Client: Leidos, Inc.  
Project/Site: RVAAP FWGW

Job ID: 280-176674-1

## Method: 8330B - Nitroaromatics and Nitramines (HPLC)

Matrix: Water

Prep Type: Total/NA

### Percent Surrogate Recovery (Acceptance Limits)

| Lab Sample ID       | Client Sample ID       | 12DNB1<br>(83-119) |
|---------------------|------------------------|--------------------|
| 280-176674-2        | FWGmw-015-230401-GW    | 85 M               |
| 280-176674-4        | LL3mw-241-230401-GW    | 93                 |
| 280-176674-7        | FBQmw-173-230401-GW-R  | 91 M               |
| 280-176674-10       | FBQmw-175-230401-GW-R  | 95 M               |
| 280-176674-10 MS    | FBQmw-175-230401-GW-R  | 82 Q               |
| 280-176674-10 MS    | FBQmw-175-230401-GW-R  | 95 M               |
| 280-176674-10 MSD   | FBQmw-175-230401-GW-R  | 98                 |
| 280-176674-10 MSD   | FBQmw-175-230401-GW-R  | 109                |
| 280-176674-11       | FBQmw-175-230402-GW-R  | 93 M               |
| LCS 280-613095/2-A  | Lab Control Sample     | 93                 |
| LCS 280-613095/4-A  | Lab Control Sample     | 98                 |
| LCSD 280-613095/3-A | Lab Control Sample Dup | 94                 |
| LCSD 280-613095/5-A | Lab Control Sample Dup | 98                 |
| MB 280-613095/1-A   | Method Blank           | 102                |

#### Surrogate Legend

12DNB = 1,2-Dinitrobenzene

## Method: 8330B - Nitroaromatics and Nitramines (HPLC)

Matrix: Water

Prep Type: Total/NA

### Percent Surrogate Recovery (Acceptance Limits)

| Lab Sample ID | Client Sample ID      | 12DNB2<br>(83-119) |
|---------------|-----------------------|--------------------|
| 280-176674-4  | LL3mw-241-230401-GW   | 88                 |
| 280-176674-7  | FBQmw-173-230401-GW-R | 85                 |

#### Surrogate Legend

12DNB = 1,2-Dinitrobenzene

# QC Sample Results

Client: Leidos, Inc.  
Project/Site: RVAAP FWGW

Job ID: 280-176674-1

## Method: 8330B - Nitroaromatics and Nitramines (HPLC)

Lab Sample ID: MB 280-613095/1-A  
Matrix: Water  
Analysis Batch: 613168

Client Sample ID: Method Blank  
Prep Type: Total/NA  
Prep Batch: 613095

| Analyte                    | MB Result | MB Qualifier | LOQ  | LOD   | DL    | Unit | D | Analyzed       | Dil Fac |
|----------------------------|-----------|--------------|------|-------|-------|------|---|----------------|---------|
| 1,3,5-Trinitrobenzene      | 0.20      | U M          | 0.21 | 0.20  | 0.084 | ug/L |   | 05/19/23 23:27 | 1       |
| 1,3-Dinitrobenzene         | 0.10      | U            | 0.11 | 0.10  | 0.037 | ug/L |   | 05/19/23 23:27 | 1       |
| 2,4,6-Trinitrotoluene      | 0.10      | U            | 0.11 | 0.10  | 0.045 | ug/L |   | 05/19/23 23:27 | 1       |
| 2,4-diamino-6-nitrotoluene | 0.90      | U            | 1.0  | 0.90  | 0.43  | ug/L |   | 05/19/23 23:27 | 1       |
| 2,4-Dinitrotoluene         | 0.080     | U            | 0.10 | 0.080 | 0.027 | ug/L |   | 05/19/23 23:27 | 1       |
| 2,6-diamino-4-nitrotoluene | 0.90      | U M          | 1.0  | 0.90  | 0.22  | ug/L |   | 05/19/23 23:27 | 1       |
| 2,6-Dinitrotoluene         | 0.080     | U            | 0.10 | 0.080 | 0.040 | ug/L |   | 05/19/23 23:27 | 1       |
| 2-Amino-4,6-dinitrotoluene | 0.10      | U            | 0.11 | 0.10  | 0.051 | ug/L |   | 05/19/23 23:27 | 1       |
| 2-Nitrotoluene             | 0.20      | U            | 0.21 | 0.20  | 0.086 | ug/L |   | 05/19/23 23:27 | 1       |
| 3,5-Dinitroaniline         | 0.30      | U            | 0.40 | 0.30  | 0.13  | ug/L |   | 05/19/23 23:27 | 1       |
| 3-Nitrotoluene             | 0.35      | U            | 0.40 | 0.35  | 0.20  | ug/L |   | 05/19/23 23:27 | 1       |
| 4-Amino-2,6-dinitrotoluene | 0.12      | U            | 0.15 | 0.12  | 0.058 | ug/L |   | 05/19/23 23:27 | 1       |
| 4-Nitrotoluene             | 0.40      | U            | 0.41 | 0.40  | 0.10  | ug/L |   | 05/19/23 23:27 | 1       |
| DNX                        | 0.25      | U            | 0.50 | 0.25  | 0.097 | ug/L |   | 05/19/23 23:27 | 1       |
| HMX                        | 0.20      | U            | 0.21 | 0.20  | 0.088 | ug/L |   | 05/19/23 23:27 | 1       |
| MNX                        | 0.29      | U            | 0.50 | 0.29  | 0.093 | ug/L |   | 05/19/23 23:27 | 1       |
| Nitrobenzene               | 0.20      | U            | 0.21 | 0.20  | 0.091 | ug/L |   | 05/19/23 23:27 | 1       |
| Nitroglycerin              | 2.0       | U            | 2.1  | 2.0   | 0.92  | ug/L |   | 05/19/23 23:27 | 1       |
| PETN                       | 1.0       | U            | 1.1  | 1.0   | 0.45  | ug/L |   | 05/19/23 23:27 | 1       |
| RDX                        | 0.20      | U            | 0.21 | 0.20  | 0.052 | ug/L |   | 05/19/23 23:27 | 1       |
| Tetryl                     | 0.10      | U            | 0.11 | 0.10  | 0.032 | ug/L |   | 05/19/23 23:27 | 1       |
| TNX                        | 0.25      | U            | 0.50 | 0.25  | 0.079 | ug/L |   | 05/19/23 23:27 | 1       |

| Surrogate          | MB %Recovery | MB Qualifier | Limits   | Prepared       | Analyzed       | Dil Fac |
|--------------------|--------------|--------------|----------|----------------|----------------|---------|
| 1,2-Dinitrobenzene | 102          |              | 83 - 119 | 05/19/23 13:39 | 05/19/23 23:27 | 1       |

Lab Sample ID: LCS 280-613095/2-A  
Matrix: Water  
Analysis Batch: 613168

Client Sample ID: Lab Control Sample  
Prep Type: Total/NA  
Prep Batch: 613095

| Analyte                    | Spike Added | LCS Result | LCS Qualifier | Unit | D | %Rec | %Rec Limits |
|----------------------------|-------------|------------|---------------|------|---|------|-------------|
| 1,3,5-Trinitrobenzene      | 2.00        | 1.98       |               | ug/L |   | 99   | 73 - 125    |
| 1,3-Dinitrobenzene         | 2.00        | 1.90       |               | ug/L |   | 95   | 78 - 120    |
| 2,4,6-Trinitrotoluene      | 2.00        | 1.83       |               | ug/L |   | 91   | 71 - 123    |
| 2,4-Dinitrotoluene         | 2.00        | 1.82       |               | ug/L |   | 91   | 78 - 120    |
| 2,6-Dinitrotoluene         | 2.00        | 1.84       |               | ug/L |   | 92   | 77 - 127    |
| 2-Amino-4,6-dinitrotoluene | 2.00        | 1.80       |               | ug/L |   | 90   | 79 - 120    |
| 2-Nitrotoluene             | 2.00        | 1.56       |               | ug/L |   | 78   | 70 - 127    |
| 3,5-Dinitroaniline         | 2.00        | 1.73       |               | ug/L |   | 87   | 71 - 117    |
| 3-Nitrotoluene             | 2.00        | 1.54       |               | ug/L |   | 77   | 73 - 125    |
| 4-Amino-2,6-dinitrotoluene | 2.00        | 1.78       |               | ug/L |   | 89   | 76 - 125    |
| 4-Nitrotoluene             | 2.00        | 1.56       |               | ug/L |   | 78   | 71 - 127    |
| DNX                        | 2.00        | 1.93       | M             | ug/L |   | 96   | 66 - 119    |
| HMX                        | 2.00        | 1.62       | M             | ug/L |   | 81   | 65 - 135    |
| MNX                        | 2.33        | 2.45       |               | ug/L |   | 105  | 57 - 132    |
| Nitrobenzene               | 2.00        | 1.74       |               | ug/L |   | 87   | 65 - 134    |
| Nitroglycerin              | 20.0        | 18.9       |               | ug/L |   | 95   | 74 - 127    |
| PETN                       | 20.0        | 20.5       |               | ug/L |   | 102  | 73 - 127    |

Eurofins Denver

# QC Sample Results

Client: Leidos, Inc.  
Project/Site: RVAAP FWGW

Job ID: 280-176674-1

## Method: 8330B - Nitroaromatics and Nitramines (HPLC) (Continued)

Lab Sample ID: LCS 280-613095/2-A

Matrix: Water

Analysis Batch: 613168

Client Sample ID: Lab Control Sample

Prep Type: Total/NA

Prep Batch: 613095

| Analyte | Spike Added | LCS Result | LCS Qualifier | Unit | D | %Rec | %Rec Limits |
|---------|-------------|------------|---------------|------|---|------|-------------|
| RDX     | 2.00        | 1.79       |               | ug/L |   | 89   | 68 - 130    |
| Tetryl  | 2.00        | 2.07       |               | ug/L |   | 103  | 64 - 128    |
| TNX     | 2.01        | 1.88       | M             | ug/L |   | 94   | 50 - 150    |

| Surrogate          | LCS %Recovery | LCS Qualifier | Limits   |
|--------------------|---------------|---------------|----------|
| 1,2-Dinitrobenzene | 93            |               | 83 - 119 |

Lab Sample ID: LCS 280-613095/4-A

Matrix: Water

Analysis Batch: 613168

Client Sample ID: Lab Control Sample

Prep Type: Total/NA

Prep Batch: 613095

| Analyte                    | Spike Added | LCS Result | LCS Qualifier | Unit | D | %Rec | %Rec Limits |
|----------------------------|-------------|------------|---------------|------|---|------|-------------|
| 2,4-diamino-6-nitrotoluene | 2.00        | 1.61       | M             | ug/L |   | 81   | 68 - 122    |
| 2,6-diamino-4-nitrotoluene | 2.00        | 1.68       | M             | ug/L |   | 84   | 72 - 122    |

| Surrogate          | LCS %Recovery | LCS Qualifier | Limits   |
|--------------------|---------------|---------------|----------|
| 1,2-Dinitrobenzene | 98            |               | 83 - 119 |

Lab Sample ID: LCSD 280-613095/3-A

Matrix: Water

Analysis Batch: 613168

Client Sample ID: Lab Control Sample Dup

Prep Type: Total/NA

Prep Batch: 613095

| Analyte                    | Spike Added | LCSD Result | LCSD Qualifier | Unit | D | %Rec | %Rec Limits | RPD | RPD Limit |
|----------------------------|-------------|-------------|----------------|------|---|------|-------------|-----|-----------|
| 1,3,5-Trinitrobenzene      | 2.00        | 2.04        |                | ug/L |   | 102  | 73 - 125    | 3   | 20        |
| 1,3-Dinitrobenzene         | 2.00        | 1.96        |                | ug/L |   | 98   | 78 - 120    | 3   | 20        |
| 2,4,6-Trinitrotoluene      | 2.00        | 1.88        |                | ug/L |   | 94   | 71 - 123    | 3   | 20        |
| 2,4-Dinitrotoluene         | 2.00        | 1.87        |                | ug/L |   | 93   | 78 - 120    | 3   | 20        |
| 2,6-Dinitrotoluene         | 2.00        | 1.95        |                | ug/L |   | 98   | 77 - 127    | 6   | 20        |
| 2-Amino-4,6-dinitrotoluene | 2.00        | 1.84        |                | ug/L |   | 92   | 79 - 120    | 2   | 20        |
| 2-Nitrotoluene             | 2.00        | 1.64        |                | ug/L |   | 82   | 70 - 127    | 5   | 20        |
| 3,5-Dinitroaniline         | 2.00        | 1.83        |                | ug/L |   | 91   | 71 - 117    | 5   | 20        |
| 3-Nitrotoluene             | 2.00        | 1.62        |                | ug/L |   | 81   | 73 - 125    | 5   | 20        |
| 4-Amino-2,6-dinitrotoluene | 2.00        | 1.82        |                | ug/L |   | 91   | 76 - 125    | 2   | 20        |
| 4-Nitrotoluene             | 2.00        | 1.61        |                | ug/L |   | 80   | 71 - 127    | 3   | 20        |
| DNX                        | 2.00        | 1.95        | M              | ug/L |   | 97   | 66 - 119    | 1   | 20        |
| HMX                        | 2.00        | 1.71        | M              | ug/L |   | 86   | 65 - 135    | 6   | 20        |
| MNX                        | 2.33        | 2.47        |                | ug/L |   | 106  | 57 - 132    | 1   | 20        |
| Nitrobenzene               | 2.00        | 1.81        |                | ug/L |   | 91   | 65 - 134    | 4   | 20        |
| Nitroglycerin              | 20.0        | 19.5        |                | ug/L |   | 97   | 74 - 127    | 3   | 20        |
| PETN                       | 20.0        | 21.1        |                | ug/L |   | 106  | 73 - 127    | 3   | 20        |
| RDX                        | 2.00        | 1.84        |                | ug/L |   | 92   | 68 - 130    | 3   | 20        |
| Tetryl                     | 2.00        | 2.06        |                | ug/L |   | 103  | 64 - 128    | 0   | 20        |
| TNX                        | 2.01        | 1.91        | M              | ug/L |   | 95   | 50 - 150    | 2   | 20        |

| Surrogate          | LCSD %Recovery | LCSD Qualifier | Limits   |
|--------------------|----------------|----------------|----------|
| 1,2-Dinitrobenzene | 94             |                | 83 - 119 |

Eurofins Denver

# QC Sample Results

Client: Leidos, Inc.  
Project/Site: RVAAP FWGW

Job ID: 280-176674-1

## Method: 8330B - Nitroaromatics and Nitramines (HPLC) (Continued)

Lab Sample ID: LCSD 280-613095/5-A

Matrix: Water

Analysis Batch: 613168

Client Sample ID: Lab Control Sample Dup

Prep Type: Total/NA

Prep Batch: 613095

| Analyte                    | Spike Added | LCSD Result | LCSD Qualifier | Unit | D | %Rec | %Rec Limits | RPD | RPD Limit |
|----------------------------|-------------|-------------|----------------|------|---|------|-------------|-----|-----------|
| 2,4-diamino-6-nitrotoluene | 2.00        | 1.47        | M              | ug/L |   | 74   | 68 - 122    | 9   | 20        |
| 2,6-diamino-4-nitrotoluene | 2.00        | 1.56        | M              | ug/L |   | 78   | 72 - 122    | 8   | 20        |

| Surrogate          | LCSD %Recovery | LCSD Qualifier | Limits   |
|--------------------|----------------|----------------|----------|
| 1,2-Dinitrobenzene | 98             |                | 83 - 119 |

Lab Sample ID: 280-176674-10 MS

Matrix: Water

Analysis Batch: 613168

Client Sample ID: FBQmw-175-230401-GW-R

Prep Type: Total/NA

Prep Batch: 613095

| Analyte                    | Sample Result | Sample Qualifier | Spike Added | MS Result | MS Qualifier | Unit | D | %Rec | %Rec Limits |
|----------------------------|---------------|------------------|-------------|-----------|--------------|------|---|------|-------------|
| 1,3,5-Trinitrobenzene      | 0.20          | U J1             | 2.00        | 1.75      | Q            | ug/L |   | 87   | 73 - 125    |
| 1,3-Dinitrobenzene         | 0.10          | U J1             | 2.00        | 1.66      | Q            | ug/L |   | 83   | 78 - 120    |
| 2,4,6-Trinitrotoluene      | 0.10          | U J1             | 2.00        | 1.56      | Q            | ug/L |   | 78   | 71 - 123    |
| 2,4-Dinitrotoluene         | 0.080         | U J1             | 2.00        | 1.50      | Q J1         | ug/L |   | 75   | 78 - 120    |
| 2,6-Dinitrotoluene         | 0.080         | U J1             | 2.00        | 1.58      | Q            | ug/L |   | 79   | 77 - 127    |
| 2-Amino-4,6-dinitrotoluene | 0.10          | U J1             | 2.00        | 1.46      | Q J1         | ug/L |   | 73   | 79 - 120    |
| 2-Nitrotoluene             | 0.20          | U J1             | 2.00        | 1.34      | Q J1         | ug/L |   | 67   | 70 - 127    |
| 3,5-Dinitroaniline         | 0.30          | U J1             | 2.00        | 1.46      | Q            | ug/L |   | 73   | 71 - 117    |
| 3-Nitrotoluene             | 0.35          | U J1             | 2.00        | 1.23      | Q J1         | ug/L |   | 61   | 73 - 125    |
| 4-Amino-2,6-dinitrotoluene | 0.12          | U J1             | 2.00        | 1.44      | Q J1         | ug/L |   | 72   | 76 - 125    |
| 4-Nitrotoluene             | 0.40          | U J1             | 2.00        | 1.30      | Q J1         | ug/L |   | 65   | 71 - 127    |
| DNX                        | 0.25          | U J1             | 2.00        | 1.56      | M Q          | ug/L |   | 78   | 66 - 119    |
| HMX                        | 0.20          | U                | 2.00        | 1.46      | M Q          | ug/L |   | 73   | 65 - 135    |
| MNX                        | 0.29          | U                | 2.34        | 1.98      | M Q          | ug/L |   | 85   | 57 - 132    |
| Nitrobenzene               | 0.20          | U J1             | 2.00        | 1.47      | Q            | ug/L |   | 73   | 65 - 134    |
| Nitroglycerin              | 2.0           | U                | 20.0        | 17.5      | Q            | ug/L |   | 87   | 74 - 127    |
| PETN                       | 1.0           | U J1             | 20.0        | 18.1      | Q            | ug/L |   | 90   | 73 - 127    |
| RDX                        | 0.20          | U M J1           | 2.00        | 1.49      | M Q          | ug/L |   | 74   | 68 - 130    |
| Tetryl                     | 0.10          | U J1             | 2.00        | 1.72      | Q            | ug/L |   | 86   | 64 - 128    |
| TNX                        | 0.25          | U M              | 2.01        | 1.70      | M Q          | ug/L |   | 85   | 50 - 150    |

| Surrogate          | MS %Recovery | MS Qualifier | Limits   |
|--------------------|--------------|--------------|----------|
| 1,2-Dinitrobenzene | 82           | Q            | 83 - 119 |

Lab Sample ID: 280-176674-10 MS

Matrix: Water

Analysis Batch: 613168

Client Sample ID: FBQmw-175-230401-GW-R

Prep Type: Total/NA

Prep Batch: 613095

| Analyte                    | Sample Result | Sample Qualifier | Spike Added | MS Result | MS Qualifier | Unit | D | %Rec | %Rec Limits |
|----------------------------|---------------|------------------|-------------|-----------|--------------|------|---|------|-------------|
| 2,4-diamino-6-nitrotoluene | 0.90          | U J1             | 2.01        | 1.25      | M J1         | ug/L |   | 62   | 68 - 122    |
| 2,6-diamino-4-nitrotoluene | 0.90          | U J1             | 2.01        | 1.60      | M            | ug/L |   | 80   | 72 - 122    |

| Surrogate          | MS %Recovery | MS Qualifier | Limits   |
|--------------------|--------------|--------------|----------|
| 1,2-Dinitrobenzene | 95           | M            | 83 - 119 |



# QC Sample Results

Client: Leidos, Inc.  
Project/Site: RVAAP FWGW

Job ID: 280-176674-1

## Method: 8330B - Nitroaromatics and Nitramines (HPLC) (Continued)

Lab Sample ID: 280-176674-10 MSD

Matrix: Water

Analysis Batch: 613168

Client Sample ID: FBQmw-175-230401-GW-R

Prep Type: Total/NA

Prep Batch: 613095

| Analyte                    | Sample Result | Sample Qualifier | Spike Added | MSD Result | MSD Qualifier | Unit | D | %Rec | %Rec Limits | RPD | RPD Limit |
|----------------------------|---------------|------------------|-------------|------------|---------------|------|---|------|-------------|-----|-----------|
| 1,3,5-Trinitrobenzene      | 0.20          | U J1             | 2.19        | 2.21       | J1            | ug/L |   | 101  | 73 - 125    | 23  | 20        |
| 1,3-Dinitrobenzene         | 0.10          | U J1             | 2.19        | 2.13       | J1            | ug/L |   | 97   | 78 - 120    | 25  | 20        |
| 2,4,6-Trinitrotoluene      | 0.10          | U J1             | 2.19        | 2.00       | J1            | ug/L |   | 92   | 71 - 123    | 25  | 20        |
| 2,4-Dinitrotoluene         | 0.080         | U J1             | 2.19        | 2.00       | J1            | ug/L |   | 92   | 78 - 120    | 29  | 20        |
| 2,6-Dinitrotoluene         | 0.080         | U J1             | 2.19        | 2.02       | J1            | ug/L |   | 93   | 77 - 127    | 25  | 20        |
| 2-Amino-4,6-dinitrotoluene | 0.10          | U J1             | 2.19        | 1.96       | J1            | ug/L |   | 90   | 79 - 120    | 29  | 20        |
| 2-Nitrotoluene             | 0.20          | U J1             | 2.19        | 1.78       | J1            | ug/L |   | 82   | 70 - 127    | 28  | 20        |
| 3,5-Dinitroaniline         | 0.30          | U J1             | 2.19        | 1.86       | J1            | ug/L |   | 85   | 71 - 117    | 24  | 20        |
| 3-Nitrotoluene             | 0.35          | U J1             | 2.19        | 1.67       | J1            | ug/L |   | 76   | 73 - 125    | 30  | 20        |
| 4-Amino-2,6-dinitrotoluene | 0.12          | U J1             | 2.19        | 1.87       | J1            | ug/L |   | 86   | 76 - 125    | 26  | 20        |
| 4-Nitrotoluene             | 0.40          | U J1             | 2.19        | 1.75       | J1            | ug/L |   | 80   | 71 - 127    | 29  | 20        |
| DNX                        | 0.25          | U J1             | 2.19        | 1.93       | M J1          | ug/L |   | 88   | 66 - 119    | 21  | 20        |
| HMX                        | 0.20          | U                | 2.19        | 1.77       | M             | ug/L |   | 81   | 65 - 135    | 19  | 20        |
| MNX                        | 0.29          | U                | 2.55        | 2.37       | M             | ug/L |   | 93   | 57 - 132    | 18  | 20        |
| Nitrobenzene               | 0.20          | U J1             | 2.19        | 1.94       | J1            | ug/L |   | 89   | 65 - 134    | 28  | 20        |
| Nitroglycerin              | 2.0           | U                | 21.9        | 21.5       |               | ug/L |   | 98   | 74 - 127    | 20  | 20        |
| PETN                       | 1.0           | U J1             | 21.9        | 22.4       | J1            | ug/L |   | 102  | 73 - 127    | 21  | 20        |
| RDX                        | 0.20          | U M J1           | 2.19        | 1.83       | M J1          | ug/L |   | 84   | 68 - 130    | 21  | 20        |
| Tetryl                     | 0.10          | U J1             | 2.19        | 2.23       | J1            | ug/L |   | 102  | 64 - 128    | 26  | 20        |
| TNX                        | 0.25          | U M              | 2.19        | 1.98       | M             | ug/L |   | 90   | 50 - 150    | 15  | 20        |

| Surrogate          | MSD %Recovery | MSD Qualifier | Limits   |
|--------------------|---------------|---------------|----------|
| 1,2-Dinitrobenzene | 98            |               | 83 - 119 |

Lab Sample ID: 280-176674-10 MSD

Matrix: Water

Analysis Batch: 613168

Client Sample ID: FBQmw-175-230401-GW-R

Prep Type: Total/NA

Prep Batch: 613095

| Analyte                    | Sample Result | Sample Qualifier | Spike Added | MSD Result | MSD Qualifier | Unit | D | %Rec | %Rec Limits | RPD | RPD Limit |
|----------------------------|---------------|------------------|-------------|------------|---------------|------|---|------|-------------|-----|-----------|
| 2,4-diamino-6-nitrotoluene | 0.90          | U J1             | 2.00        | 1.03       | M J1          | ug/L |   | 51   | 68 - 122    | 19  | 20        |
| 2,6-diamino-4-nitrotoluene | 0.90          | U J1             | 2.00        | 1.32       | M J1          | ug/L |   | 66   | 72 - 122    | 19  | 20        |

| Surrogate          | MSD %Recovery | MSD Qualifier | Limits   |
|--------------------|---------------|---------------|----------|
| 1,2-Dinitrobenzene | 109           |               | 83 - 119 |

## Method: 350.1 - Nitrogen, Ammonia

Lab Sample ID: MB 280-613561/90

Matrix: Water

Analysis Batch: 613561

Client Sample ID: Method Blank

Prep Type: Total/NA

| Analyte      | MB Result | MB Qualifier | LOQ  | LOD   | DL    | Unit | D | Analyzed       | Dil Fac |
|--------------|-----------|--------------|------|-------|-------|------|---|----------------|---------|
| Ammonia as N | 0.050     | U            | 0.10 | 0.050 | 0.029 | mg/L |   | 05/23/23 12:13 | 1       |

Eurofins Denver

# QC Sample Results

Client: Leidos, Inc.  
Project/Site: RVAAP FWGW

Job ID: 280-176674-1

## Method: 350.1 - Nitrogen, Ammonia (Continued)

Lab Sample ID: LCS 280-613561/89  
Matrix: Water  
Analysis Batch: 613561

Client Sample ID: Lab Control Sample  
Prep Type: Total/NA

| Analyte      | Spike Added | LCS Result | LCS Qualifier | Unit | D | %Rec | %Rec Limits |
|--------------|-------------|------------|---------------|------|---|------|-------------|
| Ammonia as N | 2.50        | 2.58       |               | mg/L |   | 103  | 90 - 110    |

Lab Sample ID: 280-176674-5 MS  
Matrix: Water  
Analysis Batch: 613561

Client Sample ID: LL12mw-185-230401-GW  
Prep Type: Total/NA

| Analyte      | Sample Result | Sample Qualifier | Spike Added | MS Result | MS Qualifier | Unit | D | %Rec | %Rec Limits |
|--------------|---------------|------------------|-------------|-----------|--------------|------|---|------|-------------|
| Ammonia as N | 0.050         | U                | 1.00        | 1.05      |              | mg/L |   | 105  | 90 - 110    |

Lab Sample ID: 280-176674-5 MSD  
Matrix: Water  
Analysis Batch: 613561

Client Sample ID: LL12mw-185-230401-GW  
Prep Type: Total/NA

| Analyte      | Sample Result | Sample Qualifier | Spike Added | MSD Result | MSD Qualifier | Unit | D | %Rec | %Rec Limits | RPD | RPD Limit |
|--------------|---------------|------------------|-------------|------------|---------------|------|---|------|-------------|-----|-----------|
| Ammonia as N | 0.050         | U                | 1.00        | 1.03       |               | mg/L |   | 103  | 90 - 110    | 1   | 10        |

## Method: 9034 - Sulfide, Acid Soluble and Insoluble (Titrimetric)

Lab Sample ID: MB 280-613197/2-A  
Matrix: Water  
Analysis Batch: 613208

Client Sample ID: Method Blank  
Prep Type: Total/NA  
Prep Batch: 613197

| Analyte | MB Result | MB Qualifier | LOQ | LOD | DL  | Unit | D | Analyzed       | Dil Fac |
|---------|-----------|--------------|-----|-----|-----|------|---|----------------|---------|
| Sulfide | 2.0       | U            | 4.0 | 2.0 | 1.6 | mg/L |   | 05/19/23 17:27 | 1       |

Lab Sample ID: LCS 280-613197/1-A  
Matrix: Water  
Analysis Batch: 613208

Client Sample ID: Lab Control Sample  
Prep Type: Total/NA  
Prep Batch: 613197

| Analyte | Spike Added | LCS Result | LCS Qualifier | Unit | D | %Rec | %Rec Limits |
|---------|-------------|------------|---------------|------|---|------|-------------|
| Sulfide | 19.0        | 16.0       |               | mg/L |   | 84   | 44 - 110    |

## Method: 9056 - Anions, Ion Chromatography

Lab Sample ID: MB 280-612961/13  
Matrix: Water  
Analysis Batch: 612961

Client Sample ID: Method Blank  
Prep Type: Total/NA

| Analyte      | MB Result | MB Qualifier | LOQ  | LOD  | DL    | Unit | D | Analyzed       | Dil Fac |
|--------------|-----------|--------------|------|------|-------|------|---|----------------|---------|
| Nitrate as N | 0.20      | U            | 0.50 | 0.20 | 0.090 | mg/L |   | 05/18/23 14:55 | 1       |
| Nitrite as N | 0.10      | U            | 0.50 | 0.10 | 0.049 | mg/L |   | 05/18/23 14:55 | 1       |

Lab Sample ID: MB 280-612961/55  
Matrix: Water  
Analysis Batch: 612961

Client Sample ID: Method Blank  
Prep Type: Total/NA

| Analyte      | MB Result | MB Qualifier | LOQ  | LOD  | DL    | Unit | D | Analyzed       | Dil Fac |
|--------------|-----------|--------------|------|------|-------|------|---|----------------|---------|
| Nitrate as N | 0.20      | U M          | 0.50 | 0.20 | 0.090 | mg/L |   | 05/19/23 01:26 | 1       |

# QC Sample Results

Client: Leidos, Inc.  
Project/Site: RVAAP FWGW

Job ID: 280-176674-1

## Method: 9056 - Anions, Ion Chromatography (Continued)

Lab Sample ID: LCS 280-612961/11  
Matrix: Water  
Analysis Batch: 612961

Client Sample ID: Lab Control Sample  
Prep Type: Total/NA

| Analyte      | Spike Added | LCS Result | LCS Qualifier | Unit | D | %Rec | %Rec Limits |
|--------------|-------------|------------|---------------|------|---|------|-------------|
| Nitrate as N | 2.50        | 2.46       |               | mg/L |   | 98   | 88 - 111    |
| Nitrite as N | 2.50        | 2.42       |               | mg/L |   | 97   | 87 - 111    |

Lab Sample ID: LCS 280-612961/53  
Matrix: Water  
Analysis Batch: 612961

Client Sample ID: Lab Control Sample  
Prep Type: Total/NA

| Analyte      | Spike Added | LCS Result | LCS Qualifier | Unit | D | %Rec | %Rec Limits |
|--------------|-------------|------------|---------------|------|---|------|-------------|
| Nitrate as N | 2.50        | 2.50       |               | mg/L |   | 100  | 88 - 111    |

Lab Sample ID: LCSD 280-612961/12  
Matrix: Water  
Analysis Batch: 612961

Client Sample ID: Lab Control Sample Dup  
Prep Type: Total/NA

| Analyte      | Spike Added | LCSD Result | LCSD Qualifier | Unit | D | %Rec | %Rec Limits | RPD | RPD Limit |
|--------------|-------------|-------------|----------------|------|---|------|-------------|-----|-----------|
| Nitrate as N | 2.50        | 2.45        |                | mg/L |   | 98   | 88 - 111    | 0   | 10        |
| Nitrite as N | 2.50        | 2.41        |                | mg/L |   | 96   | 87 - 111    | 1   | 10        |

Lab Sample ID: LCSD 280-612961/54  
Matrix: Water  
Analysis Batch: 612961

Client Sample ID: Lab Control Sample Dup  
Prep Type: Total/NA

| Analyte      | Spike Added | LCSD Result | LCSD Qualifier | Unit | D | %Rec | %Rec Limits | RPD | RPD Limit |
|--------------|-------------|-------------|----------------|------|---|------|-------------|-----|-----------|
| Nitrate as N | 2.50        | 2.45        |                | mg/L |   | 98   | 88 - 111    | 2   | 10        |

Lab Sample ID: MRL 280-612961/10  
Matrix: Water  
Analysis Batch: 612961

Client Sample ID: Lab Control Sample  
Prep Type: Total/NA

| Analyte      | Spike Added | MRL Result | MRL Qualifier | Unit | D | %Rec | %Rec Limits |
|--------------|-------------|------------|---------------|------|---|------|-------------|
| Nitrate as N | 0.250       | 0.248      | J             | mg/L |   | 99   | 50 - 150    |
| Nitrite as N | 0.250       | 0.221      | J             | mg/L |   | 89   | 50 - 150    |

Lab Sample ID: MB 280-613351/6  
Matrix: Water  
Analysis Batch: 613351

Client Sample ID: Method Blank  
Prep Type: Total/NA

| Analyte      | MB Result | MB Qualifier | LOQ  | LOD  | DL    | Unit | D | Analyzed       | Dil Fac |
|--------------|-----------|--------------|------|------|-------|------|---|----------------|---------|
| Nitrate as N | 0.20      | U            | 0.50 | 0.20 | 0.090 | mg/L |   | 05/22/23 17:14 | 1       |

Lab Sample ID: LCS 280-613351/4  
Matrix: Water  
Analysis Batch: 613351

Client Sample ID: Lab Control Sample  
Prep Type: Total/NA

| Analyte      | Spike Added | LCS Result | LCS Qualifier | Unit | D | %Rec | %Rec Limits |
|--------------|-------------|------------|---------------|------|---|------|-------------|
| Nitrate as N | 2.50        | 2.41       |               | mg/L |   | 96   | 88 - 111    |

# QC Sample Results

Client: Leidos, Inc.  
Project/Site: RVAAP FWGW

Job ID: 280-176674-1

## Method: 9056 - Anions, Ion Chromatography (Continued)

Lab Sample ID: LCSD 280-613351/5  
Matrix: Water  
Analysis Batch: 613351

Client Sample ID: Lab Control Sample Dup  
Prep Type: Total/NA

| Analyte      | Spike Added | LCSD Result | LCSD Qualifier | Unit | D | %Rec | %Rec Limits | RPD | RPD Limit |
|--------------|-------------|-------------|----------------|------|---|------|-------------|-----|-----------|
| Nitrate as N | 2.50        | 2.42        |                | mg/L |   | 97   | 88 - 111    | 0   | 10        |

Lab Sample ID: MRL 280-613351/3  
Matrix: Water  
Analysis Batch: 613351

Client Sample ID: Lab Control Sample  
Prep Type: Total/NA

| Analyte      | Spike Added | MRL Result | MRL Qualifier | Unit | D | %Rec | %Rec Limits |  |  |
|--------------|-------------|------------|---------------|------|---|------|-------------|--|--|
| Nitrate as N | 0.250       | 0.239      | J             | mg/L |   | 95   | 50 - 150    |  |  |

## Method: 9056A - Anions, Ion Chromatography

Lab Sample ID: MB 280-612962/13  
Matrix: Water  
Analysis Batch: 612962

Client Sample ID: Method Blank  
Prep Type: Total/NA

| Analyte | MB Result | MB Qualifier | LOQ | LOD | DL  | Unit | D | Analyzed       | Dil Fac |
|---------|-----------|--------------|-----|-----|-----|------|---|----------------|---------|
| Sulfate | 2.5       | U M          | 5.0 | 2.5 | 1.0 | mg/L |   | 05/18/23 14:55 | 1       |

Lab Sample ID: LCS 280-612962/11  
Matrix: Water  
Analysis Batch: 612962

Client Sample ID: Lab Control Sample  
Prep Type: Total/NA

| Analyte | Spike Added | LCS Result | LCS Qualifier | Unit | D | %Rec | %Rec Limits |  |  |
|---------|-------------|------------|---------------|------|---|------|-------------|--|--|
| Sulfate | 50.0        | 50.5       |               | mg/L |   | 101  | 87 - 112    |  |  |

Lab Sample ID: LCSD 280-612962/12  
Matrix: Water  
Analysis Batch: 612962

Client Sample ID: Lab Control Sample Dup  
Prep Type: Total/NA

| Analyte | Spike Added | LCSD Result | LCSD Qualifier | Unit | D | %Rec | %Rec Limits | RPD | RPD Limit |
|---------|-------------|-------------|----------------|------|---|------|-------------|-----|-----------|
| Sulfate | 50.0        | 50.0        |                | mg/L |   | 100  | 87 - 112    | 1   | 10        |

Lab Sample ID: MRL 280-612962/10  
Matrix: Water  
Analysis Batch: 612962

Client Sample ID: Lab Control Sample  
Prep Type: Total/NA

| Analyte | Spike Added | MRL Result | MRL Qualifier | Unit | D | %Rec | %Rec Limits |  |  |
|---------|-------------|------------|---------------|------|---|------|-------------|--|--|
| Sulfate | 2.50        | 2.98       | J             | mg/L |   | 119  | 50 - 150    |  |  |

## Method: 9060 - Organic Carbon, Total (TOC)

Lab Sample ID: MB 280-613268/4  
Matrix: Water  
Analysis Batch: 613268

Client Sample ID: Method Blank  
Prep Type: Total/NA

| Analyte                        | MB Result | MB Qualifier | LOQ | LOD  | DL   | Unit | D | Analyzed       | Dil Fac |
|--------------------------------|-----------|--------------|-----|------|------|------|---|----------------|---------|
| Total Organic Carbon - Average | 0.80      | U            | 1.0 | 0.80 | 0.35 | mg/L | ☼ | 05/19/23 18:46 | 1       |

# QC Sample Results

Client: Leidos, Inc.  
Project/Site: RVAAP FWGW

Job ID: 280-176674-1

## Method: 9060 - Organic Carbon, Total (TOC) (Continued)

Lab Sample ID: LCS 280-613268/3  
Matrix: Water  
Analysis Batch: 613268

Client Sample ID: Lab Control Sample  
Prep Type: Total/NA

| Analyte                        | Spike Added | LCS Result | LCS Qualifier | Unit | D | %Rec | %Rec Limits |
|--------------------------------|-------------|------------|---------------|------|---|------|-------------|
| Total Organic Carbon - Average | 25.0        | 24.8       |               | mg/L | ☼ | 99   | 88 - 112    |

Lab Sample ID: 280-176674-7 MS  
Matrix: Water  
Analysis Batch: 613268

Client Sample ID: FBQmw-173-230401-GW-R  
Prep Type: Total/NA

| Analyte                        | Sample Result | Sample Qualifier | Spike Added | MS Result | MS Qualifier | Unit | D | %Rec | %Rec Limits |
|--------------------------------|---------------|------------------|-------------|-----------|--------------|------|---|------|-------------|
| Total Organic Carbon - Average | 0.85          | J                | 25.0        | 25.0      |              | mg/L | ☼ | 97   | 88 - 112    |

Lab Sample ID: 280-176674-7 MSD  
Matrix: Water  
Analysis Batch: 613268

Client Sample ID: FBQmw-173-230401-GW-R  
Prep Type: Total/NA

| Analyte                        | Sample Result | Sample Qualifier | Spike Added | MSD Result | MSD Qualifier | Unit | D | %Rec | %Rec Limits | RPD | RPD Limit |
|--------------------------------|---------------|------------------|-------------|------------|---------------|------|---|------|-------------|-----|-----------|
| Total Organic Carbon - Average | 0.85          | J                | 25.0        | 25.1       |               | mg/L | ☼ | 97   | 88 - 112    | 0   | 15        |

## Method: SM 2320B - Alkalinity

Lab Sample ID: MB 280-613263/33  
Matrix: Water  
Analysis Batch: 613263

Client Sample ID: Method Blank  
Prep Type: Total/NA

| Analyte    | MB Result | MB Qualifier | LOQ | LOD | DL  | Unit | D | Analyzed       | Dil Fac |
|------------|-----------|--------------|-----|-----|-----|------|---|----------------|---------|
| Alkalinity | 3.58      | J            | 10  | 6.4 | 3.1 | mg/L |   | 05/19/23 16:11 | 1       |

Lab Sample ID: LCS 280-613263/31  
Matrix: Water  
Analysis Batch: 613263

Client Sample ID: Lab Control Sample  
Prep Type: Total/NA

| Analyte    | Spike Added | LCS Result | LCS Qualifier | Unit | D | %Rec | %Rec Limits |
|------------|-------------|------------|---------------|------|---|------|-------------|
| Alkalinity | 200         | 181        |               | mg/L |   | 91   | 89 - 110    |

Lab Sample ID: LCSD 280-613263/32  
Matrix: Water  
Analysis Batch: 613263

Client Sample ID: Lab Control Sample Dup  
Prep Type: Total/NA

| Analyte    | Spike Added | LCSD Result | LCSD Qualifier | Unit | D | %Rec | %Rec Limits | RPD | RPD Limit |
|------------|-------------|-------------|----------------|------|---|------|-------------|-----|-----------|
| Alkalinity | 200         | 183         |                | mg/L |   | 92   | 89 - 110    | 1   | 10        |

# QC Association Summary

Client: Leidos, Inc.  
Project/Site: RVAAP FWGW

Job ID: 280-176674-1

## HPLC/IC

### Prep Batch: 613095

| Lab Sample ID       | Client Sample ID       | Prep Type | Matrix | Method | Prep Batch |
|---------------------|------------------------|-----------|--------|--------|------------|
| 280-176674-2        | FWGmw-015-230401-GW    | Total/NA  | Water  | 3535   |            |
| 280-176674-4        | LL3mw-241-230401-GW    | Total/NA  | Water  | 3535   |            |
| 280-176674-7        | FBQmw-173-230401-GW-R  | Total/NA  | Water  | 3535   |            |
| 280-176674-10       | FBQmw-175-230401-GW-R  | Total/NA  | Water  | 3535   |            |
| 280-176674-11       | FBQmw-175-230402-GW-R  | Total/NA  | Water  | 3535   |            |
| MB 280-613095/1-A   | Method Blank           | Total/NA  | Water  | 3535   |            |
| LCS 280-613095/2-A  | Lab Control Sample     | Total/NA  | Water  | 3535   |            |
| LCS 280-613095/4-A  | Lab Control Sample     | Total/NA  | Water  | 3535   |            |
| LCSD 280-613095/3-A | Lab Control Sample Dup | Total/NA  | Water  | 3535   |            |
| LCSD 280-613095/5-A | Lab Control Sample Dup | Total/NA  | Water  | 3535   |            |
| 280-176674-10 MS    | FBQmw-175-230401-GW-R  | Total/NA  | Water  | 3535   |            |
| 280-176674-10 MS    | FBQmw-175-230401-GW-R  | Total/NA  | Water  | 3535   |            |
| 280-176674-10 MSD   | FBQmw-175-230401-GW-R  | Total/NA  | Water  | 3535   |            |
| 280-176674-10 MSD   | FBQmw-175-230401-GW-R  | Total/NA  | Water  | 3535   |            |

### Analysis Batch: 613168

| Lab Sample ID       | Client Sample ID       | Prep Type | Matrix | Method | Prep Batch |
|---------------------|------------------------|-----------|--------|--------|------------|
| 280-176674-2        | FWGmw-015-230401-GW    | Total/NA  | Water  | 8330B  | 613095     |
| 280-176674-4        | LL3mw-241-230401-GW    | Total/NA  | Water  | 8330B  | 613095     |
| 280-176674-7        | FBQmw-173-230401-GW-R  | Total/NA  | Water  | 8330B  | 613095     |
| 280-176674-10       | FBQmw-175-230401-GW-R  | Total/NA  | Water  | 8330B  | 613095     |
| 280-176674-11       | FBQmw-175-230402-GW-R  | Total/NA  | Water  | 8330B  | 613095     |
| MB 280-613095/1-A   | Method Blank           | Total/NA  | Water  | 8330B  | 613095     |
| LCS 280-613095/2-A  | Lab Control Sample     | Total/NA  | Water  | 8330B  | 613095     |
| LCS 280-613095/4-A  | Lab Control Sample     | Total/NA  | Water  | 8330B  | 613095     |
| LCSD 280-613095/3-A | Lab Control Sample Dup | Total/NA  | Water  | 8330B  | 613095     |
| LCSD 280-613095/5-A | Lab Control Sample Dup | Total/NA  | Water  | 8330B  | 613095     |
| 280-176674-10 MS    | FBQmw-175-230401-GW-R  | Total/NA  | Water  | 8330B  | 613095     |
| 280-176674-10 MS    | FBQmw-175-230401-GW-R  | Total/NA  | Water  | 8330B  | 613095     |
| 280-176674-10 MSD   | FBQmw-175-230401-GW-R  | Total/NA  | Water  | 8330B  | 613095     |
| 280-176674-10 MSD   | FBQmw-175-230401-GW-R  | Total/NA  | Water  | 8330B  | 613095     |

### Analysis Batch: 613238

| Lab Sample ID | Client Sample ID      | Prep Type | Matrix | Method | Prep Batch |
|---------------|-----------------------|-----------|--------|--------|------------|
| 280-176674-7  | FBQmw-173-230401-GW-R | Total/NA  | Water  | 8330B  | 613095     |

### Analysis Batch: 613683

| Lab Sample ID | Client Sample ID    | Prep Type | Matrix | Method | Prep Batch |
|---------------|---------------------|-----------|--------|--------|------------|
| 280-176674-4  | LL3mw-241-230401-GW | Total/NA  | Water  | 8330B  | 613095     |

## General Chemistry

### Analysis Batch: 612961

| Lab Sample ID    | Client Sample ID      | Prep Type | Matrix | Method | Prep Batch |
|------------------|-----------------------|-----------|--------|--------|------------|
| 280-176674-5     | LL12mw-185-230401-GW  | Total/NA  | Water  | 9056   |            |
| 280-176674-6     | LL12mw-185-230402-GW  | Total/NA  | Water  | 9056   |            |
| 280-176674-7     | FBQmw-173-230401-GW-R | Total/NA  | Water  | 9056   |            |
| 280-176674-8     | FBQmw-173-230402-GW-R | Total/NA  | Water  | 9056   |            |
| 280-176674-9     | LL12mw-245-230401-GW  | Total/NA  | Water  | 9056   |            |
| MB 280-612961/13 | Method Blank          | Total/NA  | Water  | 9056   |            |
| MB 280-612961/55 | Method Blank          | Total/NA  | Water  | 9056   |            |

Eurofins Denver

# QC Association Summary

Client: Leidos, Inc.  
Project/Site: RVAAP FWGW

Job ID: 280-176674-1

## General Chemistry (Continued)

### Analysis Batch: 612961 (Continued)

| Lab Sample ID      | Client Sample ID       | Prep Type | Matrix | Method | Prep Batch |
|--------------------|------------------------|-----------|--------|--------|------------|
| LCS 280-612961/11  | Lab Control Sample     | Total/NA  | Water  | 9056   |            |
| LCS 280-612961/53  | Lab Control Sample     | Total/NA  | Water  | 9056   |            |
| LCSD 280-612961/12 | Lab Control Sample Dup | Total/NA  | Water  | 9056   |            |
| LCSD 280-612961/54 | Lab Control Sample Dup | Total/NA  | Water  | 9056   |            |
| MRL 280-612961/10  | Lab Control Sample     | Total/NA  | Water  | 9056   |            |

### Analysis Batch: 612962

| Lab Sample ID      | Client Sample ID       | Prep Type | Matrix | Method | Prep Batch |
|--------------------|------------------------|-----------|--------|--------|------------|
| 280-176674-7       | FBQmw-173-230401-GW-R  | Total/NA  | Water  | 9056A  |            |
| 280-176674-8       | FBQmw-173-230402-GW-R  | Total/NA  | Water  | 9056A  |            |
| MB 280-612962/13   | Method Blank           | Total/NA  | Water  | 9056A  |            |
| LCS 280-612962/11  | Lab Control Sample     | Total/NA  | Water  | 9056A  |            |
| LCSD 280-612962/12 | Lab Control Sample Dup | Total/NA  | Water  | 9056A  |            |
| MRL 280-612962/10  | Lab Control Sample     | Total/NA  | Water  | 9056A  |            |

### Prep Batch: 613197

| Lab Sample ID      | Client Sample ID      | Prep Type | Matrix | Method | Prep Batch |
|--------------------|-----------------------|-----------|--------|--------|------------|
| 280-176674-7       | FBQmw-173-230401-GW-R | Total/NA  | Water  | 9030B  |            |
| 280-176674-8       | FBQmw-173-230402-GW-R | Total/NA  | Water  | 9030B  |            |
| MB 280-613197/2-A  | Method Blank          | Total/NA  | Water  | 9030B  |            |
| LCS 280-613197/1-A | Lab Control Sample    | Total/NA  | Water  | 9030B  |            |

### Analysis Batch: 613208

| Lab Sample ID      | Client Sample ID      | Prep Type | Matrix | Method | Prep Batch |
|--------------------|-----------------------|-----------|--------|--------|------------|
| 280-176674-7       | FBQmw-173-230401-GW-R | Total/NA  | Water  | 9034   | 613197     |
| 280-176674-8       | FBQmw-173-230402-GW-R | Total/NA  | Water  | 9034   | 613197     |
| MB 280-613197/2-A  | Method Blank          | Total/NA  | Water  | 9034   | 613197     |
| LCS 280-613197/1-A | Lab Control Sample    | Total/NA  | Water  | 9034   | 613197     |

### Analysis Batch: 613263

| Lab Sample ID      | Client Sample ID       | Prep Type | Matrix | Method   | Prep Batch |
|--------------------|------------------------|-----------|--------|----------|------------|
| 280-176674-7       | FBQmw-173-230401-GW-R  | Total/NA  | Water  | SM 2320B |            |
| 280-176674-8       | FBQmw-173-230402-GW-R  | Total/NA  | Water  | SM 2320B |            |
| MB 280-613263/33   | Method Blank           | Total/NA  | Water  | SM 2320B |            |
| LCS 280-613263/31  | Lab Control Sample     | Total/NA  | Water  | SM 2320B |            |
| LCSD 280-613263/32 | Lab Control Sample Dup | Total/NA  | Water  | SM 2320B |            |

### Analysis Batch: 613268

| Lab Sample ID    | Client Sample ID      | Prep Type | Matrix | Method | Prep Batch |
|------------------|-----------------------|-----------|--------|--------|------------|
| 280-176674-7     | FBQmw-173-230401-GW-R | Total/NA  | Water  | 9060   |            |
| MB 280-613268/4  | Method Blank          | Total/NA  | Water  | 9060   |            |
| LCS 280-613268/3 | Lab Control Sample    | Total/NA  | Water  | 9060   |            |
| 280-176674-7 MS  | FBQmw-173-230401-GW-R | Total/NA  | Water  | 9060   |            |
| 280-176674-7 MSD | FBQmw-173-230401-GW-R | Total/NA  | Water  | 9060   |            |

### Analysis Batch: 613351

| Lab Sample ID     | Client Sample ID       | Prep Type | Matrix | Method | Prep Batch |
|-------------------|------------------------|-----------|--------|--------|------------|
| 280-176674-3      | LL 12mw-187-230401-GW  | Total/NA  | Water  | 9056   |            |
| MB 280-613351/6   | Method Blank           | Total/NA  | Water  | 9056   |            |
| LCS 280-613351/4  | Lab Control Sample     | Total/NA  | Water  | 9056   |            |
| LCSD 280-613351/5 | Lab Control Sample Dup | Total/NA  | Water  | 9056   |            |

# QC Association Summary

Client: Leidos, Inc.  
Project/Site: RVAAP FWGW

Job ID: 280-176674-1

## General Chemistry (Continued)

### Analysis Batch: 613351 (Continued)

| Lab Sample ID    | Client Sample ID   | Prep Type | Matrix | Method | Prep Batch |
|------------------|--------------------|-----------|--------|--------|------------|
| MRL 280-613351/3 | Lab Control Sample | Total/NA  | Water  | 9056   |            |

### Analysis Batch: 613561

| Lab Sample ID     | Client Sample ID     | Prep Type | Matrix | Method | Prep Batch |
|-------------------|----------------------|-----------|--------|--------|------------|
| 280-176674-3      | LL12mw-187-230401-GW | Total/NA  | Water  | 350.1  |            |
| 280-176674-5      | LL12mw-185-230401-GW | Total/NA  | Water  | 350.1  |            |
| 280-176674-6      | LL12mw-185-230402-GW | Total/NA  | Water  | 350.1  |            |
| 280-176674-9      | LL12mw-245-230401-GW | Total/NA  | Water  | 350.1  |            |
| MB 280-613561/90  | Method Blank         | Total/NA  | Water  | 350.1  |            |
| LCS 280-613561/89 | Lab Control Sample   | Total/NA  | Water  | 350.1  |            |
| 280-176674-5 MS   | LL12mw-185-230401-GW | Total/NA  | Water  | 350.1  |            |
| 280-176674-5 MSD  | LL12mw-185-230401-GW | Total/NA  | Water  | 350.1  |            |



# Lab Chronicle

Client: Leidos, Inc.  
Project/Site: RVAAP FWGW

Job ID: 280-176674-1

**Client Sample ID: FWGmw-015-230401-GW**

**Lab Sample ID: 280-176674-2**

**Date Collected: 05/17/23 14:06**

**Matrix: Water**

**Date Received: 05/18/23 09:45**

| Prep Type | Batch Type | Batch Method | Run | Dil Factor | Initial Amount | Final Amount | Batch Number | Prepared or Analyzed | Analyst | Lab     |
|-----------|------------|--------------|-----|------------|----------------|--------------|--------------|----------------------|---------|---------|
| Total/NA  | Prep       | 3535         |     |            | 481.7 mL       | 5 mL         | 613095       | 05/19/23 13:39       | EH      | EET DEN |
| Total/NA  | Analysis   | 8330B        |     | 1          | 1 mL           | 1 mL         | 613168       | 05/20/23 02:53       | JZ      | EET DEN |

**Client Sample ID: LL12mw-187-230401-GW**

**Lab Sample ID: 280-176674-3**

**Date Collected: 05/17/23 13:40**

**Matrix: Water**

**Date Received: 05/18/23 09:45**

| Prep Type | Batch Type | Batch Method | Run | Dil Factor | Initial Amount | Final Amount | Batch Number | Prepared or Analyzed | Analyst | Lab     |
|-----------|------------|--------------|-----|------------|----------------|--------------|--------------|----------------------|---------|---------|
| Total/NA  | Analysis   | 350.1        |     | 200        | 10 mL          | 10 mL        | 613561       | 05/23/23 17:04       | MMP     | EET DEN |
| Total/NA  | Analysis   | 9056         |     | 2000       | 10 mL          | 10 mL        | 613351       | 05/22/23 18:59       | EJS     | EET DEN |

**Client Sample ID: LL3mw-241-230401-GW**

**Lab Sample ID: 280-176674-4**

**Date Collected: 05/17/23 15:00**

**Matrix: Water**

**Date Received: 05/18/23 09:45**

| Prep Type | Batch Type | Batch Method | Run | Dil Factor | Initial Amount | Final Amount | Batch Number | Prepared or Analyzed | Analyst | Lab     |
|-----------|------------|--------------|-----|------------|----------------|--------------|--------------|----------------------|---------|---------|
| Total/NA  | Prep       | 3535         |     |            | 488.4 mL       | 5 mL         | 613095       | 05/19/23 13:39       | EH      | EET DEN |
| Total/NA  | Analysis   | 8330B        |     | 1          | 1 mL           | 1 mL         | 613168       | 05/20/23 04:02       | JZ      | EET DEN |
| Total/NA  | Prep       | 3535         |     |            | 488.4 mL       | 5 mL         | 613095       | 05/19/23 13:39       | EH      | EET DEN |
| Total/NA  | Analysis   | 8330B        |     | 1          | 1 mL           | 1 mL         | 613683       | 05/24/23 19:48       | JZ      | EET DEN |

**Client Sample ID: LL12mw-185-230401-GW**

**Lab Sample ID: 280-176674-5**

**Date Collected: 05/17/23 09:45**

**Matrix: Water**

**Date Received: 05/18/23 09:45**

| Prep Type | Batch Type | Batch Method | Run | Dil Factor | Initial Amount | Final Amount | Batch Number | Prepared or Analyzed | Analyst | Lab     |
|-----------|------------|--------------|-----|------------|----------------|--------------|--------------|----------------------|---------|---------|
| Total/NA  | Analysis   | 350.1        |     | 1          | 10 mL          | 10 mL        | 613561       | 05/23/23 13:01       | MMP     | EET DEN |
| Total/NA  | Analysis   | 9056         |     | 20         | 10 mL          | 10 mL        | 612961       | 05/18/23 16:26       | MEC     | EET DEN |

**Client Sample ID: LL12mw-185-230402-GW**

**Lab Sample ID: 280-176674-6**

**Date Collected: 05/17/23 09:45**

**Matrix: Water**

**Date Received: 05/18/23 09:45**

| Prep Type | Batch Type | Batch Method | Run | Dil Factor | Initial Amount | Final Amount | Batch Number | Prepared or Analyzed | Analyst | Lab     |
|-----------|------------|--------------|-----|------------|----------------|--------------|--------------|----------------------|---------|---------|
| Total/NA  | Analysis   | 350.1        |     | 1          | 10 mL          | 10 mL        | 613561       | 05/23/23 13:25       | MMP     | EET DEN |
| Total/NA  | Analysis   | 9056         |     | 20         | 10 mL          | 10 mL        | 612961       | 05/18/23 16:41       | MEC     | EET DEN |

**Client Sample ID: FBQmw-173-230401-GW-R**

**Lab Sample ID: 280-176674-7**

**Date Collected: 05/17/23 11:47**

**Matrix: Water**

**Date Received: 05/18/23 09:45**

| Prep Type | Batch Type | Batch Method | Run | Dil Factor | Initial Amount | Final Amount | Batch Number | Prepared or Analyzed | Analyst | Lab     |
|-----------|------------|--------------|-----|------------|----------------|--------------|--------------|----------------------|---------|---------|
| Total/NA  | Prep       | 3535         |     |            | 483.2 mL       | 5 mL         | 613095       | 05/19/23 13:39       | EH      | EET DEN |
| Total/NA  | Analysis   | 8330B        |     | 1          | 1 mL           | 1 mL         | 613168       | 05/20/23 04:25       | JZ      | EET DEN |
| Total/NA  | Prep       | 3535         |     |            | 483.2 mL       | 5 mL         | 613095       | 05/19/23 13:39       | EH      | EET DEN |
| Total/NA  | Analysis   | 8330B        |     | 1          | 1 mL           | 1 mL         | 613238       | 05/21/23 00:31       | JZ      | EET DEN |

Eurofins Denver

# Lab Chronicle

Client: Leidos, Inc.  
Project/Site: RVAAP FWGW

Job ID: 280-176674-1

**Client Sample ID: FBQmw-173-230401-GW-R**

**Lab Sample ID: 280-176674-7**

**Date Collected: 05/17/23 11:47**

**Matrix: Water**

**Date Received: 05/18/23 09:45**

| Prep Type | Batch Type | Batch Method | Run | Dil Factor | Initial Amount | Final Amount | Batch Number | Prepared or Analyzed | Analyst | Lab     |
|-----------|------------|--------------|-----|------------|----------------|--------------|--------------|----------------------|---------|---------|
| Total/NA  | Prep       | 9030B        |     |            | 50 mL          | 50 mL        | 613197       | 05/19/23 15:58       | CAI     | EET DEN |
| Total/NA  | Analysis   | 9034         |     | 1          |                |              | 613208       | 05/19/23 17:27       | CAI     | EET DEN |
| Total/NA  | Analysis   | 9056         |     | 1          | 10 mL          | 10 mL        | 612961       | 05/18/23 19:26       | MEC     | EET DEN |
| Total/NA  | Analysis   | 9056A        |     | 1          | 10 mL          | 10 mL        | 612962       | 05/18/23 19:26       | MEC     | EET DEN |
| Total/NA  | Analysis   | 9060         |     | 1          | 20 mL          | 20 mL        | 613268       | 05/19/23 19:17       | ABW     | EET DEN |
| Total/NA  | Analysis   | SM 2320B     |     | 1          |                |              | 613263       | 05/19/23 17:16       | KEG     | EET DEN |

**Client Sample ID: FBQmw-173-230402-GW-R**

**Lab Sample ID: 280-176674-8**

**Date Collected: 05/17/23 11:47**

**Matrix: Water**

**Date Received: 05/18/23 09:45**

| Prep Type | Batch Type | Batch Method | Run | Dil Factor | Initial Amount | Final Amount | Batch Number | Prepared or Analyzed | Analyst | Lab     |
|-----------|------------|--------------|-----|------------|----------------|--------------|--------------|----------------------|---------|---------|
| Total/NA  | Prep       | 9030B        |     |            | 50 mL          | 50 mL        | 613197       | 05/19/23 15:58       | CAI     | EET DEN |
| Total/NA  | Analysis   | 9034         |     | 1          |                |              | 613208       | 05/19/23 17:27       | CAI     | EET DEN |
| Total/NA  | Analysis   | 9056         |     | 1          | 10 mL          | 10 mL        | 612961       | 05/18/23 19:41       | MEC     | EET DEN |
| Total/NA  | Analysis   | 9056A        |     | 1          | 10 mL          | 10 mL        | 612962       | 05/18/23 19:41       | MEC     | EET DEN |
| Total/NA  | Analysis   | SM 2320B     |     | 1          |                |              | 613263       | 05/19/23 17:22       | KEG     | EET DEN |

**Client Sample ID: LL12mw-245-230401-GW**

**Lab Sample ID: 280-176674-9**

**Date Collected: 05/17/23 14:26**

**Matrix: Water**

**Date Received: 05/18/23 09:45**

| Prep Type | Batch Type | Batch Method | Run | Dil Factor | Initial Amount | Final Amount | Batch Number | Prepared or Analyzed | Analyst | Lab     |
|-----------|------------|--------------|-----|------------|----------------|--------------|--------------|----------------------|---------|---------|
| Total/NA  | Analysis   | 350.1        |     | 1          | 10 mL          | 10 mL        | 613561       | 05/23/23 13:28       | MMP     | EET DEN |
| Total/NA  | Analysis   | 9056         |     | 1          | 10 mL          | 10 mL        | 612961       | 05/19/23 06:42       | MEC     | EET DEN |

**Client Sample ID: FBQmw-175-230401-GW-R**

**Lab Sample ID: 280-176674-10**

**Date Collected: 05/17/23 08:58**

**Matrix: Water**

**Date Received: 05/18/23 09:45**

| Prep Type | Batch Type | Batch Method | Run | Dil Factor | Initial Amount | Final Amount | Batch Number | Prepared or Analyzed | Analyst | Lab     |
|-----------|------------|--------------|-----|------------|----------------|--------------|--------------|----------------------|---------|---------|
| Total/NA  | Prep       | 3535         |     |            | 497.5 mL       | 5 mL         | 613095       | 05/19/23 13:39       | EH      | EET DEN |
| Total/NA  | Analysis   | 8330B        |     | 1          | 1 mL           | 1 mL         | 613168       | 05/20/23 04:48       | JZ      | EET DEN |

**Client Sample ID: FBQmw-175-230402-GW-R**

**Lab Sample ID: 280-176674-11**

**Date Collected: 05/17/23 08:58**

**Matrix: Water**

**Date Received: 05/18/23 09:45**

| Prep Type | Batch Type | Batch Method | Run | Dil Factor | Initial Amount | Final Amount | Batch Number | Prepared or Analyzed | Analyst | Lab     |
|-----------|------------|--------------|-----|------------|----------------|--------------|--------------|----------------------|---------|---------|
| Total/NA  | Prep       | 3535         |     |            | 497 mL         | 5 mL         | 613095       | 05/19/23 13:39       | EH      | EET DEN |
| Total/NA  | Analysis   | 8330B        |     | 1          | 1 mL           | 1 mL         | 613168       | 05/20/23 06:43       | JZ      | EET DEN |

## Laboratory References:

EET DEN = Eurofins Denver, 4955 Yarrow Street, Arvada, CO 80002, TEL (303)736-0100

Eurofins Denver

# Accreditation/Certification Summary

Client: Leidos, Inc.  
Project/Site: RVAAP FWGW

Job ID: 280-176674-1

## Laboratory: Eurofins Denver

The accreditations/certifications listed below are applicable to this report.

| Authority | Program               | Identification Number | Expiration Date |
|-----------|-----------------------|-----------------------|-----------------|
| A2LA      | Dept. of Defense ELAP | 2907.01               | 10-31-23        |

## Method Summary

Client: Leidos, Inc.  
Project/Site: RVAAP FWGW

Job ID: 280-176674-1

| Method   | Method Description                                 | Protocol | Laboratory |
|----------|----------------------------------------------------|----------|------------|
| 8330B    | Nitroaromatics and Nitramines (HPLC)               | EPA      | EET DEN    |
| 350.1    | Nitrogen, Ammonia                                  | EPA      | EET DEN    |
| 9034     | Sulfide, Acid Soluble and Insoluble (Titrimetric)  | SW846    | EET DEN    |
| 9056     | Anions, Ion Chromatography                         | SW846    | EET DEN    |
| 9056A    | Anions, Ion Chromatography                         | SW846    | EET DEN    |
| 9060     | Organic Carbon, Total (TOC)                        | SW846    | EET DEN    |
| SM 2320B | Alkalinity                                         | SM       | EET DEN    |
| 3535     | Solid-Phase Extraction (SPE)                       | SW846    | EET DEN    |
| 9030B    | Sulfide, Distillation (Acid Soluble and Insoluble) | SW846    | EET DEN    |

### Protocol References:

EPA = US Environmental Protection Agency

SM = "Standard Methods For The Examination Of Water And Wastewater"

SW846 = "Test Methods For Evaluating Solid Waste, Physical/Chemical Methods", Third Edition, November 1986 And Its Updates.

### Laboratory References:

EET DEN = Eurofins Denver, 4955 Yarrow Street, Arvada, CO 80002, TEL (303)736-0100

# Sample Summary

Client: Leidos, Inc.  
Project/Site: RVAAP FWGW

Job ID: 280-176674-1

| Lab Sample ID | Client Sample ID      | Matrix | Collected      | Received       |
|---------------|-----------------------|--------|----------------|----------------|
| 280-176674-2  | FWGmw-015-230401-GW   | Water  | 05/17/23 14:06 | 05/18/23 09:45 |
| 280-176674-3  | LL12mw-187-230401-GW  | Water  | 05/17/23 13:40 | 05/18/23 09:45 |
| 280-176674-4  | LL3mw-241-230401-GW   | Water  | 05/17/23 15:00 | 05/18/23 09:45 |
| 280-176674-5  | LL12mw-185-230401-GW  | Water  | 05/17/23 09:45 | 05/18/23 09:45 |
| 280-176674-6  | LL12mw-185-230402-GW  | Water  | 05/17/23 09:45 | 05/18/23 09:45 |
| 280-176674-7  | FBQmw-173-230401-GW-R | Water  | 05/17/23 11:47 | 05/18/23 09:45 |
| 280-176674-8  | FBQmw-173-230402-GW-R | Water  | 05/17/23 11:47 | 05/18/23 09:45 |
| 280-176674-9  | LL12mw-245-230401-GW  | Water  | 05/17/23 14:26 | 05/18/23 09:45 |
| 280-176674-10 | FBQmw-175-230401-GW-R | Water  | 05/17/23 08:58 | 05/18/23 09:45 |
| 280-176674-11 | FBQmw-175-230402-GW-R | Water  | 05/17/23 08:58 | 05/18/23 09:45 |

## HPLC/IC MANUAL INTEGRATION SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Instrument ID: CHHPLC\_X3 Analysis Batch Number: 599427Lab Sample ID: IC 280-599427/11 Client Sample ID: \_\_\_\_\_Date Analyzed: 01/17/23 00:33 Lab File ID: 01160039.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME              | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|----------------------------|-------------------|--------------------|---------|----------------|
|                            |                   | REASON             | ANALYST | DATE           |
| 2,6-diamino-4-nitrotoluene | 6.43              | Baseline           | LV5D    | 01/17/23 12:22 |
| 2,4-diamino-6-nitrotoluene | 6.62              | Baseline           | LV5D    | 01/17/23 12:22 |

Lab Sample ID: IC 280-599427/12 Client Sample ID: \_\_\_\_\_Date Analyzed: 01/17/23 00:56 Lab File ID: 01160040.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME              | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|----------------------------|-------------------|--------------------|---------|----------------|
|                            |                   | REASON             | ANALYST | DATE           |
| 2,6-diamino-4-nitrotoluene | 6.46              | Baseline           | LV5D    | 01/17/23 12:23 |
| 2,4-diamino-6-nitrotoluene | 6.64              | Baseline           | LV5D    | 01/17/23 12:23 |

Lab Sample ID: IC 280-599427/13 Client Sample ID: \_\_\_\_\_Date Analyzed: 01/17/23 01:19 Lab File ID: 01160041.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME              | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|----------------------------|-------------------|--------------------|---------|----------------|
|                            |                   | REASON             | ANALYST | DATE           |
| 2,6-diamino-4-nitrotoluene | 6.46              | Baseline           | LV5D    | 01/17/23 12:23 |
| 2,4-diamino-6-nitrotoluene | 6.64              | Baseline           | LV5D    | 01/17/23 12:23 |

Lab Sample ID: IC 280-599427/14 Client Sample ID: \_\_\_\_\_Date Analyzed: 01/17/23 01:43 Lab File ID: 01160042.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME              | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|----------------------------|-------------------|--------------------|---------|----------------|
|                            |                   | REASON             | ANALYST | DATE           |
| 2,6-diamino-4-nitrotoluene | 6.46              | Baseline           | LV5D    | 01/17/23 12:30 |
| 2,4-diamino-6-nitrotoluene | 6.64              | Baseline           | LV5D    | 01/17/23 12:30 |

## HPLC/IC MANUAL INTEGRATION SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Instrument ID: CHHPLC\_X3 Analysis Batch Number: 599427Lab Sample ID: IC 280-599427/15 Client Sample ID: \_\_\_\_\_Date Analyzed: 01/17/23 02:06 Lab File ID: 01160043.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME              | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|----------------------------|-------------------|--------------------|---------|----------------|
|                            |                   | REASON             | ANALYST | DATE           |
| 2,6-diamino-4-nitrotoluene | 6.46              | Baseline           | LV5D    | 01/17/23 12:31 |
| 2,4-diamino-6-nitrotoluene | 6.64              | Baseline           | LV5D    | 01/17/23 12:31 |

Lab Sample ID: IC 280-599427/16 Client Sample ID: \_\_\_\_\_Date Analyzed: 01/17/23 02:30 Lab File ID: 01160044.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME              | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|----------------------------|-------------------|--------------------|---------|----------------|
|                            |                   | REASON             | ANALYST | DATE           |
| 2,6-diamino-4-nitrotoluene | 6.46              | Baseline           | LV5D    | 01/17/23 12:31 |
| 2,4-diamino-6-nitrotoluene | 6.63              | Baseline           | LV5D    | 01/17/23 12:31 |

Lab Sample ID: IC 280-599427/17 Client Sample ID: \_\_\_\_\_Date Analyzed: 01/17/23 02:53 Lab File ID: 01160045.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME              | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|----------------------------|-------------------|--------------------|---------|----------------|
|                            |                   | REASON             | ANALYST | DATE           |
| 2,6-diamino-4-nitrotoluene | 6.46              | Baseline           | LV5D    | 01/17/23 12:31 |
| 2,4-diamino-6-nitrotoluene | 6.64              | Baseline           | LV5D    | 01/17/23 12:31 |

Lab Sample ID: IC 280-599427/18 Client Sample ID: \_\_\_\_\_Date Analyzed: 01/17/23 03:16 Lab File ID: 01160046.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME              | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|----------------------------|-------------------|--------------------|---------|----------------|
|                            |                   | REASON             | ANALYST | DATE           |
| 2,6-diamino-4-nitrotoluene | 6.46              | Baseline           | LV5D    | 01/17/23 12:31 |
| 2,4-diamino-6-nitrotoluene | 6.64              | Baseline           | LV5D    | 01/17/23 12:31 |

## HPLC/IC MANUAL INTEGRATION SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Instrument ID: CHHPLC\_X3 Analysis Batch Number: 599427Lab Sample ID: ICV 280-599427/19 Client Sample ID: \_\_\_\_\_Date Analyzed: 01/17/23 03:40 Lab File ID: 01160047.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME              | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|----------------------------|-------------------|--------------------|---------|----------------|
|                            |                   | REASON             | ANALYST | DATE           |
| 2,6-diamino-4-nitrotoluene | 6.45              | Baseline           | LV5D    | 01/17/23 12:32 |
| 2,4-diamino-6-nitrotoluene | 6.63              | Baseline           | LV5D    | 01/17/23 12:32 |



## HPLC/IC MANUAL INTEGRATION SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Instrument ID: CHHPLC\_X3 Analysis Batch Number: 601664Lab Sample ID: IC 280-601664/11 Client Sample ID: \_\_\_\_\_Date Analyzed: 02/08/23 15:38 Lab File ID: 02080011.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|---------------|-------------------|--------------------|---------|----------------|
|               |                   | REASON             | ANALYST | DATE           |
| HMX           | 6.55              | Baseline Smoothing | LV5D    | 02/08/23 16:28 |

Lab Sample ID: IC 280-601664/12 Client Sample ID: \_\_\_\_\_Date Analyzed: 02/08/23 16:01 Lab File ID: 02080012.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|---------------|-------------------|--------------------|---------|----------------|
|               |                   | REASON             | ANALYST | DATE           |
| HMX           | 6.55              | Baseline Smoothing | LV5D    | 02/08/23 16:28 |

Lab Sample ID: IC 280-601664/13 Client Sample ID: \_\_\_\_\_Date Analyzed: 02/08/23 16:24 Lab File ID: 02080013.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|---------------|-------------------|--------------------|---------|----------------|
|               |                   | REASON             | ANALYST | DATE           |
| HMX           | 6.54              | Baseline Smoothing | LV5D    | 02/08/23 16:55 |
| RDX           | 7.56              | Baseline Smoothing | LV5D    | 02/08/23 19:09 |

Lab Sample ID: IC 280-601664/14 Client Sample ID: \_\_\_\_\_Date Analyzed: 02/08/23 16:47 Lab File ID: 02080014.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|---------------|-------------------|--------------------|---------|----------------|
|               |                   | REASON             | ANALYST | DATE           |
| HMX           | 6.55              | Baseline Smoothing | LV5D    | 02/08/23 17:11 |
| RDX           | 7.56              | Baseline Smoothing | LV5D    | 02/08/23 19:09 |

Lab Sample ID: IC 280-601664/15 Client Sample ID: \_\_\_\_\_Date Analyzed: 02/08/23 17:10 Lab File ID: 02080015.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|---------------|-------------------|--------------------|---------|----------------|
|               |                   | REASON             | ANALYST | DATE           |
| RDX           | 7.57              | Baseline Smoothing | LV5D    | 02/08/23 19:09 |

8330B

## HPLC/IC MANUAL INTEGRATION SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Instrument ID: CHHPLC\_X3 Analysis Batch Number: 601664Lab Sample ID: IC 280-601664/16 Client Sample ID: \_\_\_\_\_Date Analyzed: 02/08/23 17:33 Lab File ID: 02080016.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|---------------|-------------------|--------------------|---------|----------------|
|               |                   | REASON             | ANALYST | DATE           |
| RDX           | 7.56              | Baseline Smoothing | LV5D    | 02/08/23 19:09 |
| PETN          | 14.55             | Baseline Smoothing | LV5D    | 02/08/23 19:08 |

Lab Sample ID: IC 280-601664/17 Client Sample ID: \_\_\_\_\_Date Analyzed: 02/08/23 17:56 Lab File ID: 02080017.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|---------------|-------------------|--------------------|---------|----------------|
|               |                   | REASON             | ANALYST | DATE           |
| HMX           | 6.55              | Baseline Smoothing | LV5D    | 02/08/23 19:07 |
| RDX           | 7.57              | Baseline Smoothing | LV5D    | 02/08/23 19:09 |
| PETN          | 14.58             | Baseline Smoothing | LV5D    | 02/08/23 19:08 |

Lab Sample ID: IC 280-601664/18 Client Sample ID: \_\_\_\_\_Date Analyzed: 02/08/23 18:19 Lab File ID: 02080018.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|---------------|-------------------|--------------------|---------|----------------|
|               |                   | REASON             | ANALYST | DATE           |
| HMX           | 6.55              | Baseline Smoothing | LV5D    | 02/08/23 19:09 |
| RDX           | 7.56              | Baseline Smoothing | LV5D    | 02/08/23 19:09 |
| PETN          | 14.59             | Baseline Smoothing | LV5D    | 02/08/23 19:08 |

Lab Sample ID: IC 280-601664/19 Client Sample ID: \_\_\_\_\_Date Analyzed: 02/08/23 18:42 Lab File ID: 02080019.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|---------------|-------------------|--------------------|---------|----------------|
|               |                   | REASON             | ANALYST | DATE           |
| HMX           | 6.55              | Baseline Smoothing | LV5D    | 02/08/23 19:09 |
| RDX           | 7.56              | Baseline Smoothing | LV5D    | 02/08/23 19:08 |
| PETN          | 14.58             | Unspecified        |         |                |

## HPLC/IC MANUAL INTEGRATION SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Instrument ID: CHHPLC\_X3 Analysis Batch Number: 601664Lab Sample ID: ICV 280-601664/20 Client Sample ID: \_\_\_\_\_Date Analyzed: 02/08/23 19:05 Lab File ID: 02080020.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|---------------|-------------------|--------------------|---------|----------------|
|               |                   | REASON             | ANALYST | DATE           |
| HMX           | 6.55              | Baseline Smoothing | LV5D    | 02/08/23 19:31 |
| RDX           | 7.56              | Baseline Smoothing | LV5D    | 02/08/23 19:32 |

## HPLC/IC MANUAL INTEGRATION SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Instrument ID: CHHPLC\_X3 Analysis Batch Number: 603284Lab Sample ID: IC 280-603284/11 Client Sample ID: \_\_\_\_\_Date Analyzed: 02/24/23 15:00 Lab File ID: 02240011.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|---------------|-------------------|--------------------|---------|----------------|
|               |                   | REASON             | ANALYST | DATE           |
| TNX           | 6.44              | Baseline           | LV5D    | 02/24/23 16:41 |
| DNX           | 6.76              | Baseline           | LV5D    | 02/24/23 16:41 |
| MNX           | 7.19              | Baseline           | LV5D    | 02/24/23 16:41 |

Lab Sample ID: IC 280-603284/12 Client Sample ID: \_\_\_\_\_Date Analyzed: 02/24/23 15:23 Lab File ID: 02240012.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|---------------|-------------------|--------------------|---------|----------------|
|               |                   | REASON             | ANALYST | DATE           |
| TNX           | 6.44              | Baseline           | LV5D    | 02/24/23 16:41 |
| DNX           | 6.76              | Baseline           | LV5D    | 02/24/23 16:41 |
| MNX           | 7.19              | Baseline           | LV5D    | 02/24/23 16:41 |

Lab Sample ID: IC 280-603284/13 Client Sample ID: \_\_\_\_\_Date Analyzed: 02/24/23 15:46 Lab File ID: 02240013.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|---------------|-------------------|--------------------|---------|----------------|
|               |                   | REASON             | ANALYST | DATE           |
| TNX           | 6.44              | Baseline           | LV5D    | 02/24/23 16:41 |
| DNX           | 6.77              | Baseline           | LV5D    | 02/24/23 16:41 |
| MNX           | 7.19              | Baseline           | LV5D    | 02/24/23 16:41 |

Lab Sample ID: IC 280-603284/14 Client Sample ID: \_\_\_\_\_Date Analyzed: 02/24/23 16:09 Lab File ID: 02240014.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|---------------|-------------------|--------------------|---------|----------------|
|               |                   | REASON             | ANALYST | DATE           |
| TNX           | 6.44              | Baseline           | LV5D    | 02/24/23 16:41 |
| DNX           | 6.77              | Baseline           | LV5D    | 02/24/23 16:41 |
| MNX           | 7.20              | Baseline           | LV5D    | 02/24/23 16:40 |

## HPLC/IC MANUAL INTEGRATION SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Instrument ID: CHHPLC\_X3 Analysis Batch Number: 603284Lab Sample ID: IC 280-603284/15 Client Sample ID: \_\_\_\_\_Date Analyzed: 02/24/23 16:32 Lab File ID: 02240015.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|---------------|-------------------|--------------------|---------|----------------|
|               |                   | REASON             | ANALYST | DATE           |
| TNX           | 6.44              | Baseline           | LV5D    | 02/24/23 17:15 |
| DNX           | 6.77              | Baseline           | LV5D    | 02/24/23 17:15 |

Lab Sample ID: IC 280-603284/16 Client Sample ID: \_\_\_\_\_Date Analyzed: 02/24/23 16:55 Lab File ID: 02240016.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|---------------|-------------------|--------------------|---------|----------------|
|               |                   | REASON             | ANALYST | DATE           |
| TNX           | 6.45              | Baseline           | LV5D    | 02/24/23 17:30 |
| DNX           | 6.77              | Baseline           | LV5D    | 02/24/23 17:30 |
| MNX           | 7.20              | Baseline           | LV5D    | 02/24/23 17:30 |

Lab Sample ID: IC 280-603284/17 Client Sample ID: \_\_\_\_\_Date Analyzed: 02/24/23 17:18 Lab File ID: 02240017.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|---------------|-------------------|--------------------|---------|----------------|
|               |                   | REASON             | ANALYST | DATE           |
| TNX           | 6.45              | Baseline           | LV5D    | 02/24/23 17:51 |
| DNX           | 6.77              | Baseline           | LV5D    | 02/24/23 17:51 |

Lab Sample ID: IC 280-603284/18 Client Sample ID: \_\_\_\_\_Date Analyzed: 02/24/23 17:41 Lab File ID: 02240018.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|---------------|-------------------|--------------------|---------|----------------|
|               |                   | REASON             | ANALYST | DATE           |
| TNX           | 6.44              | Baseline           | LV5D    | 02/24/23 18:10 |
| DNX           | 6.77              | Baseline           | LV5D    | 02/24/23 18:10 |

## HPLC/IC MANUAL INTEGRATION SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Instrument ID: CHHPLC\_X3 Analysis Batch Number: 603284Lab Sample ID: ICV 280-603284/19 Client Sample ID: \_\_\_\_\_Date Analyzed: 02/24/23 18:04 Lab File ID: 02240019.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|---------------|-------------------|--------------------|---------|----------------|
|               |                   | REASON             | ANALYST | DATE           |
| TNX           | 6.45              | Baseline           | LV5D    | 02/28/23 11:45 |
| DNX           | 6.77              | Baseline           | LV5D    | 02/28/23 11:45 |

## HPLC/IC MANUAL INTEGRATION SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Instrument ID: CHHPLC\_X3 Analysis Batch Number: 613168Lab Sample ID: CCV 280-613168/32 Client Sample ID: \_\_\_\_\_Date Analyzed: 05/19/23 22:41 Lab File ID: 05190032.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|---------------|-------------------|--------------------|---------|----------------|
|               |                   | REASON             | ANALYST | DATE           |
| TNX           | 6.43              | Baseline           | LV5D    | 05/20/23 10:26 |
| HMX           | 6.55              | Baseline           | LV5D    | 05/20/23 10:26 |
| DNX           | 6.76              | Baseline           | LV5D    | 05/20/23 10:26 |

Lab Sample ID: CCV 280-613168/33 Client Sample ID: \_\_\_\_\_Date Analyzed: 05/19/23 23:04 Lab File ID: 05190033.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME              | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|----------------------------|-------------------|--------------------|---------|----------------|
|                            |                   | REASON             | ANALYST | DATE           |
| 2,6-diamino-4-nitrotoluene | 6.43              | Baseline           | LV5D    | 05/20/23 10:26 |
| 2,4-diamino-6-nitrotoluene | 6.62              | Baseline           | LV5D    | 05/20/23 10:26 |

Lab Sample ID: MB 280-613095/1-A Client Sample ID: \_\_\_\_\_Date Analyzed: 05/19/23 23:27 Lab File ID: 05190034.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME              | RETENTION<br>TIME | MANUAL INTEGRATION  |         |                |
|----------------------------|-------------------|---------------------|---------|----------------|
|                            |                   | REASON              | ANALYST | DATE           |
| 1,3,5-Trinitrobenzene      |                   | Invalid Compound ID | LV5D    | 05/20/23 10:26 |
| 2,6-diamino-4-nitrotoluene |                   | Invalid Compound ID | LV5D    | 05/20/23 10:26 |

Lab Sample ID: LCS 280-613095/2-A Client Sample ID: \_\_\_\_\_Date Analyzed: 05/19/23 23:50 Lab File ID: 05190035.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|---------------|-------------------|--------------------|---------|----------------|
|               |                   | REASON             | ANALYST | DATE           |
| TNX           | 6.43              | Baseline           | LV5D    | 05/20/23 10:26 |
| HMX           | 6.55              | Baseline           | LV5D    | 05/20/23 10:26 |
| DNX           | 6.76              | Baseline           | LV5D    | 05/20/23 10:26 |

## HPLC/IC MANUAL INTEGRATION SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Instrument ID: CHHPLC\_X3 Analysis Batch Number: 613168Lab Sample ID: LCSD 280-613095/3-A Client Sample ID: \_\_\_\_\_Date Analyzed: 05/20/23 00:13 Lab File ID: 05190036.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|---------------|-------------------|--------------------|---------|----------------|
|               |                   | REASON             | ANALYST | DATE           |
| TNX           | 6.43              | Baseline           | LV5D    | 05/20/23 10:26 |
| HMX           | 6.55              | Baseline           | LV5D    | 05/20/23 10:26 |
| DNX           | 6.76              | Baseline           | LV5D    | 05/20/23 10:26 |

Lab Sample ID: LCS 280-613095/4-A Client Sample ID: \_\_\_\_\_Date Analyzed: 05/20/23 00:36 Lab File ID: 05190037.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME              | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|----------------------------|-------------------|--------------------|---------|----------------|
|                            |                   | REASON             | ANALYST | DATE           |
| 2,6-diamino-4-nitrotoluene | 6.42              | Baseline           | LV5D    | 05/20/23 10:27 |
| 2,4-diamino-6-nitrotoluene | 6.61              | Baseline           | LV5D    | 05/20/23 10:27 |

Lab Sample ID: LCSD 280-613095/5-A Client Sample ID: \_\_\_\_\_Date Analyzed: 05/20/23 00:59 Lab File ID: 05190038.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME              | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|----------------------------|-------------------|--------------------|---------|----------------|
|                            |                   | REASON             | ANALYST | DATE           |
| 2,6-diamino-4-nitrotoluene | 6.42              | Baseline           | LV5D    | 05/20/23 10:27 |
| 2,4-diamino-6-nitrotoluene | 6.61              | Baseline           | LV5D    | 05/20/23 10:27 |

Lab Sample ID: 280-176674-2 Client Sample ID: FWGmw-015-230401-GWDate Analyzed: 05/20/23 02:53 Lab File ID: 05190043.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME      | RETENTION<br>TIME | MANUAL INTEGRATION  |         |                |
|--------------------|-------------------|---------------------|---------|----------------|
|                    |                   | REASON              | ANALYST | DATE           |
| 1,2-Dinitrobenzene | 8.53              | Baseline            | LV5D    | 05/20/23 11:05 |
| HMX                |                   | Invalid Compound ID | LV5D    | 05/20/23 11:05 |



## HPLC/IC MANUAL INTEGRATION SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Instrument ID: CHHPLC\_X3 Analysis Batch Number: 613168Lab Sample ID: CCV 280-613168/44 Client Sample ID: \_\_\_\_\_Date Analyzed: 05/20/23 03:16 Lab File ID: 05190044.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|---------------|-------------------|--------------------|---------|----------------|
|               |                   | REASON             | ANALYST | DATE           |
| TNX           | 6.43              | Baseline           | LV5D    | 05/20/23 11:05 |
| HMX           | 6.55              | Baseline           | LV5D    | 05/20/23 11:05 |
| DNX           | 6.76              | Baseline           | LV5D    | 05/20/23 11:05 |

Lab Sample ID: CCV 280-613168/45 Client Sample ID: \_\_\_\_\_Date Analyzed: 05/20/23 03:39 Lab File ID: 05190045.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME              | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|----------------------------|-------------------|--------------------|---------|----------------|
|                            |                   | REASON             | ANALYST | DATE           |
| 2,6-diamino-4-nitrotoluene | 6.43              | Baseline           | LV5D    | 05/20/23 11:05 |
| 2,4-diamino-6-nitrotoluene | 6.62              | Baseline           | LV5D    | 05/20/23 11:05 |

Lab Sample ID: 280-176674-4 Client Sample ID: LL3mw-241-230401-GWDate Analyzed: 05/20/23 04:02 Lab File ID: 05190046.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME              | RETENTION<br>TIME | MANUAL INTEGRATION  |         |                |
|----------------------------|-------------------|---------------------|---------|----------------|
|                            |                   | REASON              | ANALYST | DATE           |
| HMX                        | 6.55              | Baseline            | LV5D    | 05/20/23 11:05 |
| 2,6-diamino-4-nitrotoluene |                   | Invalid Compound ID | LV5D    | 05/20/23 11:05 |
| 4-Nitrotoluene             |                   | Invalid Compound ID | LV5D    | 05/20/23 11:06 |
| DNX                        |                   | Invalid Compound ID | LV5D    | 05/20/23 11:05 |
| Nitrobenzene               |                   | Invalid Compound ID | LV5D    | 05/20/23 11:05 |
| Tetryl                     |                   | Invalid Compound ID | LV5D    | 05/20/23 11:05 |

## HPLC/IC MANUAL INTEGRATION SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Instrument ID: CHHPLC\_X3 Analysis Batch Number: 613168Lab Sample ID: 280-176674-7 Client Sample ID: FBQmw-173-230401-GW-RDate Analyzed: 05/20/23 04:25 Lab File ID: 05190047.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME      | RETENTION TIME | MANUAL INTEGRATION  |         |                |
|--------------------|----------------|---------------------|---------|----------------|
|                    |                | REASON              | ANALYST | DATE           |
| 1,2-Dinitrobenzene | 8.53           | Baseline            | LV5D    | 05/20/23 11:06 |
| 1,3-Dinitrobenzene | 9.27           | Baseline            | LV5D    | 05/20/23 11:06 |
| RDX                |                | Invalid Compound ID | LV5D    | 05/20/23 11:06 |
| TNX                |                | Invalid Compound ID | LV5D    | 05/20/23 11:06 |

Lab Sample ID: 280-176674-10 Client Sample ID: FBQmw-175-230401-GW-RDate Analyzed: 05/20/23 04:48 Lab File ID: 05190048.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME      | RETENTION TIME | MANUAL INTEGRATION  |         |                |
|--------------------|----------------|---------------------|---------|----------------|
|                    |                | REASON              | ANALYST | DATE           |
| 1,2-Dinitrobenzene | 8.53           | Baseline            | LV5D    | 05/20/23 11:14 |
| RDX                |                | Invalid Compound ID | LV5D    | 05/20/23 11:14 |
| TNX                |                | Invalid Compound ID | LV5D    | 05/20/23 11:14 |

Lab Sample ID: 280-176674-10 MS Client Sample ID: FBQmw-175-230401-GW-R MSDate Analyzed: 05/20/23 05:11 Lab File ID: 05190049.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME | RETENTION TIME | MANUAL INTEGRATION |         |                |
|---------------|----------------|--------------------|---------|----------------|
|               |                | REASON             | ANALYST | DATE           |
| TNX           | 6.43           | Baseline           | LV5D    | 05/20/23 11:14 |
| HMX           | 6.55           | Baseline           | LV5D    | 05/20/23 11:14 |
| DNX           | 6.76           | Baseline           | LV5D    | 05/20/23 11:14 |
| MNX           | 7.18           | Baseline           | LV5D    | 05/20/23 11:14 |
| RDX           | 7.57           | Baseline           | LV5D    | 05/20/23 11:15 |

## HPLC/IC MANUAL INTEGRATION SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Instrument ID: CHHPLC\_X3 Analysis Batch Number: 613168Lab Sample ID: 280-176674-10 MSD Client Sample ID: FBQmw-175-230401-GW-R MSDDate Analyzed: 05/20/23 05:34 Lab File ID: 05190050.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|---------------|-------------------|--------------------|---------|----------------|
|               |                   | REASON             | ANALYST | DATE           |
| TNX           | 6.43              | Baseline           | LV5D    | 05/20/23 11:15 |
| HMX           | 6.55              | Baseline           | LV5D    | 05/20/23 11:15 |
| DNX           | 6.75              | Baseline           | LV5D    | 05/20/23 11:15 |
| MNX           | 7.18              | Baseline           | LV5D    | 05/20/23 11:15 |
| RDX           | 7.57              | Baseline           | LV5D    | 05/20/23 11:15 |

Lab Sample ID: 280-176674-10 MS Client Sample ID: FBQmw-175-230401-GW-R MSDate Analyzed: 05/20/23 05:57 Lab File ID: 05190051.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME              | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|----------------------------|-------------------|--------------------|---------|----------------|
|                            |                   | REASON             | ANALYST | DATE           |
| 2,6-diamino-4-nitrotoluene | 6.42              | Baseline           | LV5D    | 05/20/23 11:15 |
| 2,4-diamino-6-nitrotoluene | 6.60              | Baseline           | LV5D    | 05/20/23 11:15 |
| 1,2-Dinitrobenzene         | 8.52              | Baseline           | LV5D    | 05/20/23 11:15 |

Lab Sample ID: 280-176674-10 MSD Client Sample ID: FBQmw-175-230401-GW-R MSDDate Analyzed: 05/20/23 06:20 Lab File ID: 05190052.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME              | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|----------------------------|-------------------|--------------------|---------|----------------|
|                            |                   | REASON             | ANALYST | DATE           |
| 2,6-diamino-4-nitrotoluene | 6.42              | Baseline           | LV5D    | 05/20/23 11:15 |
| 2,4-diamino-6-nitrotoluene | 6.60              | Baseline           | LV5D    | 05/20/23 11:15 |

## HPLC/IC MANUAL INTEGRATION SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Instrument ID: CHHPLC\_X3 Analysis Batch Number: 613168Lab Sample ID: 280-176674-11 Client Sample ID: FBQmw-175-230402-GW-RDate Analyzed: 05/20/23 06:43 Lab File ID: 05190053.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME      | RETENTION<br>TIME | MANUAL INTEGRATION  |         |                |
|--------------------|-------------------|---------------------|---------|----------------|
|                    |                   | REASON              | ANALYST | DATE           |
| 1,2-Dinitrobenzene | 8.53              | Baseline            | LV5D    | 05/20/23 11:16 |
| RDX                |                   | Invalid Compound ID | LV5D    | 05/20/23 11:16 |
| TNX                |                   | Invalid Compound ID | LV5D    | 05/20/23 11:16 |
| 2-Nitrotoluene     | 12.50             | Baseline            | LV5D    | 05/20/23 11:16 |

Lab Sample ID: CCV 280-613168/56 Client Sample ID: \_\_\_\_\_Date Analyzed: 05/20/23 07:52 Lab File ID: 05190056.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|---------------|-------------------|--------------------|---------|----------------|
|               |                   | REASON             | ANALYST | DATE           |
| TNX           | 6.43              | Baseline           | LV5D    | 05/20/23 11:16 |
| HMX           | 6.55              | Baseline           | LV5D    | 05/20/23 11:16 |
| DNX           | 6.76              | Baseline           | LV5D    | 05/20/23 11:16 |

Lab Sample ID: CCV 280-613168/57 Client Sample ID: \_\_\_\_\_Date Analyzed: 05/20/23 08:15 Lab File ID: 05190057.D GC Column: UltraCarb5uOD ID: 4.6 (mm)

| COMPOUND NAME              | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|----------------------------|-------------------|--------------------|---------|----------------|
|                            |                   | REASON             | ANALYST | DATE           |
| 2,6-diamino-4-nitrotoluene | 6.43              | Baseline           | LV5D    | 05/20/23 11:16 |
| 2,4-diamino-6-nitrotoluene | 6.62              | Baseline           | LV5D    | 05/20/23 11:16 |

## HPLC/IC MANUAL INTEGRATION SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Instrument ID: CHHPLC\_X5 Analysis Batch Number: 610603Lab Sample ID: IC 280-610603/10 Client Sample ID: \_\_\_\_\_Date Analyzed: 04/28/23 18:24 Lab File ID: 04280010.D GC Column: Luna-phenylhe ID: 4.6 (mm)

| COMPOUND NAME              | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|----------------------------|-------------------|--------------------|---------|----------------|
|                            |                   | REASON             | ANALYST | DATE           |
| 3,5-Dinitroaniline         | 14.37             | Baseline Smoothing | LV5D    | 04/28/23 19:01 |
| 1,3-Dinitrobenzene         | 14.53             | Baseline Smoothing | LV5D    | 04/28/23 19:01 |
| 2-Nitrotoluene             | 15.66             | Baseline Smoothing | LV5D    | 04/28/23 19:01 |
| 4-Nitrotoluene             | 15.83             | Baseline Smoothing | LV5D    | 04/28/23 19:01 |
| 4-Amino-2,6-dinitrotoluene | 16.44             | Baseline Smoothing | LV5D    | 04/28/23 19:01 |
| 3-Nitrotoluene             | 16.73             | Baseline Smoothing | LV5D    | 04/28/23 19:01 |
| 1,3,5-Trinitrobenzene      | 17.37             | Baseline Smoothing | LV5D    | 04/28/23 19:01 |
| 2-Amino-4,6-dinitrotoluene | 17.37             | Baseline Smoothing | LV5D    | 04/28/23 19:01 |
| 2,6-Dinitrotoluene         | 18.56             | Baseline Smoothing | LV5D    | 04/28/23 19:01 |
| 2,4-Dinitrotoluene         | 19.02             | Baseline Smoothing | LV5D    | 04/28/23 19:01 |

Lab Sample ID: IC 280-610603/11 Client Sample ID: \_\_\_\_\_Date Analyzed: 04/28/23 18:59 Lab File ID: 04280011.D GC Column: Luna-phenylhe ID: 4.6 (mm)

| COMPOUND NAME              | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|----------------------------|-------------------|--------------------|---------|----------------|
|                            |                   | REASON             | ANALYST | DATE           |
| 2-Amino-4,6-dinitrotoluene | 17.40             | Baseline Smoothing | LV5D    | 04/28/23 20:43 |
| 1,3,5-Trinitrobenzene      | 17.44             | Baseline Smoothing | LV5D    | 04/28/23 20:43 |

Lab Sample ID: IC 280-610603/12 Client Sample ID: \_\_\_\_\_Date Analyzed: 04/28/23 19:34 Lab File ID: 04280012.D GC Column: Luna-phenylhe ID: 4.6 (mm)

| COMPOUND NAME              | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|----------------------------|-------------------|--------------------|---------|----------------|
|                            |                   | REASON             | ANALYST | DATE           |
| 2-Amino-4,6-dinitrotoluene | 17.33             | Baseline Smoothing | LV5D    | 04/28/23 20:43 |
| 1,3,5-Trinitrobenzene      | 17.39             | Baseline Smoothing | LV5D    | 04/28/23 20:43 |

## HPLC/IC MANUAL INTEGRATION SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Instrument ID: CHHPLC\_X5 Analysis Batch Number: 610603Lab Sample ID: IC 280-610603/13 Client Sample ID: \_\_\_\_\_Date Analyzed: 04/28/23 20:08 Lab File ID: 04280013.D GC Column: Luna-phenylhe ID: 4.6 (mm)

| COMPOUND NAME              | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|----------------------------|-------------------|--------------------|---------|----------------|
|                            |                   | REASON             | ANALYST | DATE           |
| 2-Amino-4,6-dinitrotoluene | 17.35             | Baseline Smoothing | LV5D    | 04/28/23 20:52 |
| 1,3,5-Trinitrobenzene      | 17.41             | Baseline Smoothing | LV5D    | 04/28/23 20:52 |

Lab Sample ID: IC 280-610603/14 Client Sample ID: \_\_\_\_\_Date Analyzed: 04/28/23 20:43 Lab File ID: 04280014.D GC Column: Luna-phenylhe ID: 4.6 (mm)

| COMPOUND NAME              | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|----------------------------|-------------------|--------------------|---------|----------------|
|                            |                   | REASON             | ANALYST | DATE           |
| 2-Amino-4,6-dinitrotoluene | 17.36             | Unspecified        |         |                |
| 1,3,5-Trinitrobenzene      | 17.42             | Baseline Smoothing | LV5D    | 04/29/23 10:22 |

Lab Sample ID: IC 280-610603/15 Client Sample ID: \_\_\_\_\_Date Analyzed: 04/28/23 21:18 Lab File ID: 04280015.D GC Column: Luna-phenylhe ID: 4.6 (mm)

| COMPOUND NAME              | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|----------------------------|-------------------|--------------------|---------|----------------|
|                            |                   | REASON             | ANALYST | DATE           |
| 2-Amino-4,6-dinitrotoluene | 17.37             | Unspecified        |         |                |
| 1,3,5-Trinitrobenzene      | 17.42             | Baseline Smoothing | LV5D    | 04/29/23 10:22 |

Lab Sample ID: IC 280-610603/16 Client Sample ID: \_\_\_\_\_Date Analyzed: 04/28/23 21:53 Lab File ID: 04280016.D GC Column: Luna-phenylhe ID: 4.6 (mm)

| COMPOUND NAME              | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|----------------------------|-------------------|--------------------|---------|----------------|
|                            |                   | REASON             | ANALYST | DATE           |
| 2-Amino-4,6-dinitrotoluene | 17.38             | Unspecified        |         |                |
| 1,3,5-Trinitrobenzene      | 17.43             | Baseline Smoothing | LV5D    | 04/29/23 10:22 |

## HPLC/IC MANUAL INTEGRATION SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Instrument ID: CHHPLC\_X5 Analysis Batch Number: 610603Lab Sample ID: IC 280-610603/17 Client Sample ID: \_\_\_\_\_Date Analyzed: 04/28/23 22:28 Lab File ID: 04280017.D GC Column: Luna-phenylhe ID: 4.6 (mm)

| COMPOUND NAME              | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|----------------------------|-------------------|--------------------|---------|----------------|
|                            |                   | REASON             | ANALYST | DATE           |
| 2-Nitrotoluene             | 15.77             | Baseline Smoothing | LV5D    | 04/29/23 10:23 |
| 4-Nitrotoluene             | 16.00             | Unspecified        |         |                |
| 4-Amino-2,6-dinitrotoluene | 16.56             | Baseline Smoothing | LV5D    | 04/29/23 10:23 |
| 3-Nitrotoluene             | 16.84             | Baseline Smoothing | LV5D    | 04/29/23 10:23 |
| 2-Amino-4,6-dinitrotoluene | 17.41             | Unspecified        |         |                |
| 1,3,5-Trinitrobenzene      | 17.46             | Baseline Smoothing | LV5D    | 04/29/23 10:23 |

Lab Sample ID: IC 280-610603/18 Client Sample ID: \_\_\_\_\_Date Analyzed: 04/28/23 23:03 Lab File ID: 04280018.D GC Column: Luna-phenylhe ID: 4.6 (mm)

| COMPOUND NAME              | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|----------------------------|-------------------|--------------------|---------|----------------|
|                            |                   | REASON             | ANALYST | DATE           |
| 2-Nitrotoluene             | 15.74             | Baseline Smoothing | LV5D    | 04/29/23 10:24 |
| 4-Nitrotoluene             | 15.97             | Unspecified        |         |                |
| 4-Amino-2,6-dinitrotoluene | 16.53             | Baseline Smoothing | LV5D    | 04/29/23 13:01 |
| 3-Nitrotoluene             | 16.84             | Baseline Smoothing | LV5D    | 04/29/23 13:01 |
| 2-Amino-4,6-dinitrotoluene | 17.37             | Unspecified        |         |                |
| 1,3,5-Trinitrobenzene      | 17.44             | Baseline Smoothing | LV5D    | 04/29/23 10:24 |

## HPLC/IC MANUAL INTEGRATION SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Instrument ID: CHHPLC\_X5 Analysis Batch Number: 610603Lab Sample ID: ICV 280-610603/19 Client Sample ID: \_\_\_\_\_Date Analyzed: 04/28/23 23:38 Lab File ID: 04280019.D GC Column: Luna-phenylhe ID: 4.6 (mm)

| COMPOUND NAME              | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|----------------------------|-------------------|--------------------|---------|----------------|
|                            |                   | REASON             | ANALYST | DATE           |
| 3,5-Dinitroaniline         | 14.32             | Baseline Smoothing | LV5D    | 04/29/23 10:27 |
| 1,3-Dinitrobenzene         | 14.59             | Baseline Smoothing | LV5D    | 04/29/23 10:27 |
| 2-Nitrotoluene             | 15.70             | Baseline Smoothing | LV5D    | 04/29/23 10:27 |
| 4-Nitrotoluene             | 15.92             | Baseline Smoothing | LV5D    | 04/29/23 10:27 |
| 4-Amino-2,6-dinitrotoluene | 16.48             | Baseline Smoothing | LV5D    | 04/29/23 10:27 |
| 3-Nitrotoluene             | 16.78             | Baseline Smoothing | LV5D    | 04/29/23 10:27 |
| 2-Amino-4,6-dinitrotoluene | 17.34             | Unspecified        |         |                |
| 1,3,5-Trinitrobenzene      | 17.42             | Baseline Smoothing | LV5D    | 04/29/23 10:27 |



## HPLC/IC MANUAL INTEGRATION SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Instrument ID: CHHPLC\_X5 Analysis Batch Number: 613238Lab Sample ID: CCV 280-613238/7 Client Sample ID: \_\_\_\_\_Date Analyzed: 05/20/23 18:42 Lab File ID: 05200007.D GC Column: Luna-phenylhe ID: 4.6 (mm)

| COMPOUND NAME              | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|----------------------------|-------------------|--------------------|---------|----------------|
|                            |                   | REASON             | ANALYST | DATE           |
| 3-Nitrotoluene             | 16.80             | Split Peak         | K8YG    | 05/22/23 11:10 |
| 2-Amino-4,6-dinitrotoluene | 17.35             | Split Peak         | K8YG    | 05/22/23 11:11 |
| 1,3,5-Trinitrobenzene      | 17.40             | Split Peak         | K8YG    | 05/22/23 11:11 |
| Tetryl                     | 22.33             | Split Peak         | K8YG    | 05/22/23 11:11 |
| 2,4,6-Trinitrotoluene      | 23.13             | Split Peak         | K8YG    | 05/22/23 11:11 |

Lab Sample ID: 280-176674-7 Client Sample ID: FBQmw-173-230401-GW-RDate Analyzed: 05/21/23 00:31 Lab File ID: 05200019.D GC Column: Luna-phenylhe ID: 4.6 (mm)

| COMPOUND NAME      | RETENTION<br>TIME | MANUAL INTEGRATION  |         |                |
|--------------------|-------------------|---------------------|---------|----------------|
|                    |                   | REASON              | ANALYST | DATE           |
| 2,6-Dinitrotoluene |                   | Invalid Compound ID | LV5D    | 05/23/23 13:11 |

Lab Sample ID: CCV 280-613238/20 Client Sample ID: \_\_\_\_\_Date Analyzed: 05/21/23 01:06 Lab File ID: 05200020.D GC Column: Luna-phenylhe ID: 4.6 (mm)

| COMPOUND NAME              | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|----------------------------|-------------------|--------------------|---------|----------------|
|                            |                   | REASON             | ANALYST | DATE           |
| 3,5-Dinitroaniline         | 14.32             | Split Peak         | K8YG    | 05/22/23 11:25 |
| 1,3-Dinitrobenzene         | 14.58             | Split Peak         | K8YG    | 05/22/23 11:25 |
| 2-Nitrotoluene             | 15.65             | Split Peak         | K8YG    | 05/22/23 11:25 |
| 4-Nitrotoluene             | 15.87             | Split Peak         | K8YG    | 05/22/23 11:25 |
| 4-Amino-2,6-dinitrotoluene | 16.43             | Split Peak         | K8YG    | 05/22/23 11:25 |
| 3-Nitrotoluene             | 16.71             | Split Peak         | K8YG    | 05/22/23 11:27 |
| 2-Amino-4,6-dinitrotoluene | 17.25             | Unspecified        | LV5D    | 05/23/23 12:44 |
| 1,3,5-Trinitrobenzene      | 17.37             | Split Peak         | K8YG    | 05/22/23 11:26 |

## HPLC/IC MANUAL INTEGRATION SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Instrument ID: CHHPLC\_X5 Analysis Batch Number: 613683Lab Sample ID: CCV 280-613683/8 Client Sample ID: \_\_\_\_\_Date Analyzed: 05/24/23 19:13 Lab File ID: 05240008.D GC Column: Luna-phenylhe ID: 4.6 (mm)

| COMPOUND NAME              | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|----------------------------|-------------------|--------------------|---------|----------------|
|                            |                   | REASON             | ANALYST | DATE           |
| 2-Amino-4,6-dinitrotoluene | 17.24             | Baseline Smoothing | LV5D    | 05/24/23 19:54 |
| 1,3,5-Trinitrobenzene      | 17.33             | Unspecified        |         |                |

Lab Sample ID: 280-176674-4 Client Sample ID: LL3mw-241-230401-GWDate Analyzed: 05/24/23 19:48 Lab File ID: 05240010.D GC Column: Luna-phenylhe ID: 4.6 (mm)

| COMPOUND NAME | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|---------------|-------------------|--------------------|---------|----------------|
|               |                   | REASON             | ANALYST | DATE           |
| HMX           | 6.73              | Baseline           | K8YG    | 05/25/23 08:40 |
| MNX           | 7.56              | Baseline           | K8YG    | 05/25/23 08:41 |

Lab Sample ID: CCV 280-613683/18 Client Sample ID: \_\_\_\_\_Date Analyzed: 05/25/23 00:28 Lab File ID: 05240018.D GC Column: Luna-phenylhe ID: 4.6 (mm)

| COMPOUND NAME              | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|----------------------------|-------------------|--------------------|---------|----------------|
|                            |                   | REASON             | ANALYST | DATE           |
| 2-Amino-4,6-dinitrotoluene | 17.12             | Baseline           | K8YG    | 05/25/23 09:17 |
| 1,3,5-Trinitrobenzene      | 17.26             | Baseline           | K8YG    | 05/25/23 09:17 |
| Tetryl                     | 22.08             | Baseline           | K8YG    | 05/25/23 09:17 |
| PETN                       | 24.11             | Baseline           | K8YG    | 05/25/23 09:17 |

## GENERAL CHEMISTRY MANUAL INTEGRATION SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Instrument ID: WC\_IonChrom10 Analysis Batch Number: 612961Lab Sample ID: 280-176674-5 Client Sample ID: LL12mw-185-230401-GWDate Analyzed: 05/18/23 16:26 Lab File ID: Info 2\_DENPC179\_Anions\_20 GC Column: Metrosepp A S ID: \_\_\_\_\_

| COMPOUND NAME | RETENTION TIME | MANUAL INTEGRATION |         |                |
|---------------|----------------|--------------------|---------|----------------|
|               |                | REASON             | ANALYST | DATE           |
| Nitrate as N  | 6.91           | Baseline Smoothing | LVW8    | 05/19/23 11:58 |

Lab Sample ID: 280-176674-6 Client Sample ID: LL12mw-185-230402-GWDate Analyzed: 05/18/23 16:41 Lab File ID: Info 2\_DENPC179\_Anions\_20 GC Column: Metrosepp A S ID: \_\_\_\_\_

| COMPOUND NAME | RETENTION TIME | MANUAL INTEGRATION |         |                |
|---------------|----------------|--------------------|---------|----------------|
|               |                | REASON             | ANALYST | DATE           |
| Nitrate as N  | 6.90           | Baseline Smoothing | LVW8    | 05/19/23 11:59 |

Lab Sample ID: CCB 280-612961/25 Client Sample ID: \_\_\_\_\_Date Analyzed: 05/18/23 17:56 Lab File ID: Info 2\_DENPC179\_Anions\_20 GC Column: Metrosepp A S ID: \_\_\_\_\_

| COMPOUND NAME | RETENTION TIME | MANUAL INTEGRATION |         |      |
|---------------|----------------|--------------------|---------|------|
|               |                | REASON             | ANALYST | DATE |
| Nitrite as N  |                | Unspecified        |         |      |

Lab Sample ID: 280-176674-7 Client Sample ID: FBQmw-173-230401-GW-RDate Analyzed: 05/18/23 19:26 Lab File ID: Info 2\_DENPC179\_Anions\_20 GC Column: Metrosepp A S ID: \_\_\_\_\_

| COMPOUND NAME | RETENTION TIME | MANUAL INTEGRATION |         |                |
|---------------|----------------|--------------------|---------|----------------|
|               |                | REASON             | ANALYST | DATE           |
| Nitrate as N  | 6.98           | Baseline Smoothing | LVW8    | 05/19/23 12:01 |

Lab Sample ID: 280-176674-8 Client Sample ID: FBQmw-173-230402-GW-RDate Analyzed: 05/18/23 19:41 Lab File ID: Info 2\_DENPC179\_Anions\_20 GC Column: Metrosepp A S ID: \_\_\_\_\_

| COMPOUND NAME | RETENTION TIME | MANUAL INTEGRATION |         |                |
|---------------|----------------|--------------------|---------|----------------|
|               |                | REASON             | ANALYST | DATE           |
| Nitrate as N  | 6.98           | Baseline Smoothing | LVW8    | 05/19/23 12:02 |

## GENERAL CHEMISTRY MANUAL INTEGRATION SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Instrument ID: WC\_IonChrom10 Analysis Batch Number: 612961Lab Sample ID: CCB 280-612961/37 Client Sample ID: \_\_\_\_\_Date Analyzed: 05/18/23 20:57 Lab File ID: Info 2\_DENPC179\_Anions\_20 GC Column: Metrosepp A S ID: \_\_\_\_\_

| COMPOUND NAME | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|---------------|-------------------|--------------------|---------|----------------|
|               |                   | REASON             | ANALYST | DATE           |
| Nitrite as N  | 5.00              | Baseline Smoothing | LVW8    | 05/19/23 12:03 |
| Nitrate as N  | 6.97              | Baseline Smoothing | LVW8    | 05/19/23 12:03 |

Lab Sample ID: CCB 280-612961/49 Client Sample ID: \_\_\_\_\_Date Analyzed: 05/18/23 23:56 Lab File ID: Info 2\_DENPC179\_Anions\_20 GC Column: Metrosepp A S ID: \_\_\_\_\_

| COMPOUND NAME | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|---------------|-------------------|--------------------|---------|----------------|
|               |                   | REASON             | ANALYST | DATE           |
| Nitrate as N  | 6.98              | Baseline Smoothing | LVW8    | 05/19/23 12:07 |

Lab Sample ID: MB 280-612961/55 Client Sample ID: \_\_\_\_\_Date Analyzed: 05/19/23 01:26 Lab File ID: Info 2\_DENPC179\_Anions\_20 GC Column: Metrosepp A S ID: \_\_\_\_\_

| COMPOUND NAME | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|---------------|-------------------|--------------------|---------|----------------|
|               |                   | REASON             | ANALYST | DATE           |
| Nitrate as N  | 6.98              | Baseline Smoothing | LVW8    | 05/19/23 12:07 |

Lab Sample ID: CCV 280-612961/60 Client Sample ID: \_\_\_\_\_Date Analyzed: 05/19/23 02:42 Lab File ID: Info 2\_DENPC179\_Anions\_20 GC Column: Metrosepp A S ID: \_\_\_\_\_

| COMPOUND NAME | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|---------------|-------------------|--------------------|---------|----------------|
|               |                   | REASON             | ANALYST | DATE           |
| Nitrate as N  | 6.90              | Baseline Smoothing | XAY4    | 05/19/23 09:48 |

Lab Sample ID: CCB 280-612961/61 Client Sample ID: \_\_\_\_\_Date Analyzed: 05/19/23 02:57 Lab File ID: Info 2\_DENPC179\_Anions\_20 GC Column: Metrosepp A S ID: \_\_\_\_\_

| COMPOUND NAME | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|---------------|-------------------|--------------------|---------|----------------|
|               |                   | REASON             | ANALYST | DATE           |
| Nitrate as N  | 7.01              | Baseline Smoothing | LVW8    | 05/19/23 12:09 |

## GENERAL CHEMISTRY MANUAL INTEGRATION SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Instrument ID: WC\_IonChrom10 Analysis Batch Number: 612961Lab Sample ID: CCV 280-612961/72 Client Sample ID: \_\_\_\_\_Date Analyzed: 05/19/23 05:42 Lab File ID: Info 2\_DENPC179\_Anions\_20 GC Column: Metrosepp A S ID: \_\_\_\_\_

| COMPOUND NAME | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|---------------|-------------------|--------------------|---------|----------------|
|               |                   | REASON             | ANALYST | DATE           |
| Nitrite as N  | 5.00              | Baseline Smoothing | XAY4    | 05/19/23 09:48 |
| Nitrate as N  | 6.91              | Baseline Smoothing | LVW8    | 05/19/23 12:11 |

Lab Sample ID: CCB 280-612961/73 Client Sample ID: \_\_\_\_\_Date Analyzed: 05/19/23 05:57 Lab File ID: Info 2\_DENPC179\_Anions\_20 GC Column: Metrosepp A S ID: \_\_\_\_\_

| COMPOUND NAME | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|---------------|-------------------|--------------------|---------|----------------|
|               |                   | REASON             | ANALYST | DATE           |
| Nitrate as N  | 7.00              | Baseline Smoothing | LVW8    | 05/19/23 12:12 |

Lab Sample ID: 280-176674-9 Client Sample ID: LL12mw-245-230401-GWDate Analyzed: 05/19/23 06:42 Lab File ID: Info 2\_DENPC179\_Anions\_20 GC Column: Metrosepp A S ID: \_\_\_\_\_

| COMPOUND NAME | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|---------------|-------------------|--------------------|---------|----------------|
|               |                   | REASON             | ANALYST | DATE           |
| Nitrate as N  | 7.09              | Baseline Smoothing | LVW8    | 05/19/23 12:13 |

Lab Sample ID: CCV 280-612961/84 Client Sample ID: \_\_\_\_\_Date Analyzed: 05/19/23 08:42 Lab File ID: Info 2\_DENPC179\_Anions\_20 GC Column: Metrosepp A S ID: \_\_\_\_\_

| COMPOUND NAME | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|---------------|-------------------|--------------------|---------|----------------|
|               |                   | REASON             | ANALYST | DATE           |
| Nitrite as N  | 5.01              | Baseline Smoothing | XAY4    | 05/19/23 09:48 |

Lab Sample ID: CCB 280-612961/85 Client Sample ID: \_\_\_\_\_Date Analyzed: 05/19/23 08:57 Lab File ID: Info 2\_DENPC179\_Anions\_20 GC Column: Metrosepp A S ID: \_\_\_\_\_

| COMPOUND NAME | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|---------------|-------------------|--------------------|---------|----------------|
|               |                   | REASON             | ANALYST | DATE           |
| Nitrate as N  | 7.00              | Baseline Smoothing | LVW8    | 05/19/23 12:14 |

## GENERAL CHEMISTRY MANUAL INTEGRATION SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Instrument ID: WC\_IonChrom10 Analysis Batch Number: 612962Lab Sample ID: STD 280-612962/2 IC Client Sample ID: \_\_\_\_\_Date Analyzed: 05/18/23 12:11 Lab File ID: Info 2\_DENPC179\_Anions\_20 GC Column: Metrosepp A S ID: \_\_\_\_\_

| COMPOUND NAME | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|---------------|-------------------|--------------------|---------|----------------|
|               |                   | REASON             | ANALYST | DATE           |
| Fluoride      | 3.13              | Baseline Smoothing | LVW8    | 05/18/23 12:58 |

Lab Sample ID: STD 280-612962/3 IC Client Sample ID: \_\_\_\_\_Date Analyzed: 05/18/23 12:26 Lab File ID: Info 2\_DENPC179\_Anions\_20 GC Column: Metrosepp A S ID: \_\_\_\_\_

| COMPOUND NAME | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|---------------|-------------------|--------------------|---------|----------------|
|               |                   | REASON             | ANALYST | DATE           |
| Fluoride      | 3.14              | Baseline Smoothing | LVW8    | 05/18/23 12:58 |

Lab Sample ID: STD 280-612962/4 IC Client Sample ID: \_\_\_\_\_Date Analyzed: 05/18/23 12:40 Lab File ID: Info 2\_DENPC179\_Anions\_20 GC Column: Metrosepp A S ID: \_\_\_\_\_

| COMPOUND NAME | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|---------------|-------------------|--------------------|---------|----------------|
|               |                   | REASON             | ANALYST | DATE           |
| Fluoride      | 3.14              | Baseline Smoothing | LVW8    | 05/18/23 13:29 |

Lab Sample ID: STD 280-612962/5 IC Client Sample ID: \_\_\_\_\_Date Analyzed: 05/18/23 12:55 Lab File ID: Info 2\_DENPC179\_Anions\_20 GC Column: Metrosepp A S ID: \_\_\_\_\_

| COMPOUND NAME | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|---------------|-------------------|--------------------|---------|----------------|
|               |                   | REASON             | ANALYST | DATE           |
| Fluoride      | 3.14              | Baseline Smoothing | LVW8    | 05/18/23 13:29 |

Lab Sample ID: STD 280-612962/6 IC Client Sample ID: \_\_\_\_\_Date Analyzed: 05/18/23 13:10 Lab File ID: Info 2\_DENPC179\_Anions\_20 GC Column: Metrosepp A S ID: \_\_\_\_\_

| COMPOUND NAME | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|---------------|-------------------|--------------------|---------|----------------|
|               |                   | REASON             | ANALYST | DATE           |
| Fluoride      | 3.15              | Baseline Smoothing | LVW8    | 05/18/23 13:30 |

## GENERAL CHEMISTRY MANUAL INTEGRATION SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Instrument ID: WC\_IonChrom10 Analysis Batch Number: 612962Lab Sample ID: ICB 280-612962/8 Client Sample ID: \_\_\_\_\_Date Analyzed: 05/18/23 13:40 Lab File ID: Info 2\_DENPC179\_Anions\_20 GC Column: Metrosepp A S ID: \_\_\_\_\_

| COMPOUND NAME | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|---------------|-------------------|--------------------|---------|----------------|
|               |                   | REASON             | ANALYST | DATE           |
| Sulfate       | 10.29             | Baseline Smoothing | LVW8    | 05/19/23 11:57 |

Lab Sample ID: MB 280-612962/13 Client Sample ID: \_\_\_\_\_Date Analyzed: 05/18/23 14:55 Lab File ID: Info 2\_DENPC179\_Anions\_20 GC Column: Metrosepp A S ID: \_\_\_\_\_

| COMPOUND NAME | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|---------------|-------------------|--------------------|---------|----------------|
|               |                   | REASON             | ANALYST | DATE           |
| Sulfate       | 10.30             | Baseline Smoothing | LVW8    | 05/19/23 11:57 |

Lab Sample ID: CCB 280-612962/25 Client Sample ID: \_\_\_\_\_Date Analyzed: 05/18/23 17:56 Lab File ID: Info 2\_DENPC179\_Anions\_20 GC Column: Metrosepp A S ID: \_\_\_\_\_

| COMPOUND NAME | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|---------------|-------------------|--------------------|---------|----------------|
|               |                   | REASON             | ANALYST | DATE           |
| Sulfate       | 10.31             | Baseline Smoothing | LVW8    | 05/19/23 12:00 |

Lab Sample ID: 280-176674-7 Client Sample ID: FBQmw-173-230401-GW-RDate Analyzed: 05/18/23 19:26 Lab File ID: Info 2\_DENPC179\_Anions\_20 GC Column: Metrosepp A S ID: \_\_\_\_\_

| COMPOUND NAME | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|---------------|-------------------|--------------------|---------|----------------|
|               |                   | REASON             | ANALYST | DATE           |
| Sulfate       | 10.24             | Baseline Smoothing | LVW8    | 05/19/23 12:02 |

Lab Sample ID: CCB 280-612962/37 Client Sample ID: \_\_\_\_\_Date Analyzed: 05/18/23 20:57 Lab File ID: Info 2\_DENPC179\_Anions\_20 GC Column: Metrosepp A S ID: \_\_\_\_\_

| COMPOUND NAME | RETENTION<br>TIME | MANUAL INTEGRATION |         |                |
|---------------|-------------------|--------------------|---------|----------------|
|               |                   | REASON             | ANALYST | DATE           |
| Sulfate       | 10.25             | Baseline Smoothing | LVW8    | 05/19/23 12:03 |

REAGENT TRACEABILITY SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

| Reagent ID               | Exp Date | Prep Date | Dilutant Used                        | Reagent Final Volume | Parent Reagent      |              | Analyte                    | Concentration |
|--------------------------|----------|-----------|--------------------------------------|----------------------|---------------------|--------------|----------------------------|---------------|
|                          |          |           |                                      |                      | Reagent ID          | Volume Added |                            |               |
| <b>3,5-DNA LCS_00042</b> | 01/25/23 | 11/08/22  | Acetonitrile, Lot Acetonitrile 00065 | 100 mL               | 3,5-DNA Stock_00036 | 1 mL         | 3,5-Dinitroaniline         | 10 ug/mL      |
| .3,5-DNA Stock_00036     | 04/25/23 |           | Restek, Lot A0149736                 |                      | (Purchased Reagent) |              | 3,5-Dinitroaniline         | 1000 ug/mL    |
| <b>3,5-DNA LCS_00043</b> | 08/03/23 | 02/03/23  | Acetonitrile, Lot Acetonitrile 00074 | 100 mL               | 3,5-DNA Stock_00038 | 1 mL         | 3,5-Dinitroaniline         | 10 ug/mL      |
| .3,5-DNA Stock_00038     | 02/03/24 |           | Restek, Lot A0164976                 |                      | (Purchased Reagent) |              | 3,5-Dinitroaniline         | 1000 ug/mL    |
| <b>350.1 cal_00583</b>   | 05/24/23 | 05/17/23  | Di Water, Lot na                     | 100 mL               | NH3 CAL STD_00035   | 10 mL        | Ammonia as N               | 100 mg/L      |
| .NH3 CAL STD_00035       | 01/31/24 |           | Ricca, Lot 4208D17                   |                      | (Purchased Reagent) |              | Ammonia as N               | 1000 mg/L     |
| <b>350.1 ICV_00564</b>   | 05/24/23 | 05/17/23  | na, Lot na                           | 100 mL               | NH3 ICV STD_00034   | 10 mL        | Ammonia as N               | 100.2 mg/L    |
| .NH3 ICV STD_00034       | 01/25/25 |           | Inorganic Ventures, Lot S2-NH700817  |                      | (Purchased Reagent) |              | Ammonia as N               | 1002 mg/L     |
| <b>8330 DMT_00013</b>    | 08/24/23 | 02/24/23  | Acetonitrile, Lot ACN_237            | 5 mL                 | MXN,TNX,DNX_00076   | 1 mL         | DNX                        | 20.02 ug/mL   |
|                          |          |           |                                      |                      |                     |              | MXN                        | 23.34 ug/mL   |
|                          |          |           |                                      |                      |                     |              | TNX                        | 20.06 ug/mL   |
| .MXN,TNX,DNX_00076       | 02/29/24 |           | Agilent, Lot 0006725091              |                      | (Purchased Reagent) |              | DNX                        | 100.1 ug/mL   |
|                          |          |           |                                      |                      |                     |              | MXN                        | 116.7 ug/mL   |
|                          |          |           |                                      |                      |                     |              | TNX                        | 100.3 ug/mL   |
| <b>8330 LCS_00121</b>    | 07/21/23 | 01/21/23  | Acetonitrile, Lot Acetonitrile_00073 | 50 mL                | 8330_NG_Stk_00118   | 0.6 mL       | Nitroglycerin              | 100 ug/mL     |
|                          |          |           |                                      |                      | 8330_NG_Stk_00120   | 0.4 mL       | Nitroglycerin              | 100 ug/mL     |
|                          |          |           |                                      |                      | 8330_PETN_Stk_00128 | 1 mL         | PETN                       | 100 ug/mL     |
|                          |          |           |                                      |                      | 8330LCSMxl_00136    | 0.5 mL       | 1,3,5-Trinitrobenzene      | 10 ug/mL      |
|                          |          |           |                                      |                      |                     |              | 1,3-Dinitrobenzene         | 10 ug/mL      |
|                          |          |           |                                      |                      |                     |              | 2,4,6-Trinitrotoluene      | 10 ug/mL      |
|                          |          |           |                                      |                      |                     |              | 2,4-Dinitrotoluene         | 10 ug/mL      |
|                          |          |           |                                      |                      |                     |              | HMX                        | 10 ug/mL      |
|                          |          |           |                                      |                      |                     |              | Nitrobenzene               | 10 ug/mL      |
|                          |          |           |                                      |                      |                     |              | RDX                        | 10 ug/mL      |
|                          |          |           |                                      |                      | 8330LCSmix2_00034   | 0.5 mL       | 2,6-Dinitrotoluene         | 10 ug/mL      |
|                          |          |           |                                      |                      |                     |              | 2-Amino-4,6-dinitrotoluene | 10 ug/mL      |
|                          |          |           |                                      |                      |                     |              | 2-Nitrotoluene             | 10 ug/mL      |
|                          |          |           |                                      |                      |                     |              | 3-Nitrotoluene             | 10 ug/mL      |
|                          |          |           |                                      |                      |                     |              | 4-Amino-2,6-dinitrotoluene | 10 ug/mL      |
|                          |          |           |                                      |                      |                     |              | 4-Nitrotoluene             | 10 ug/mL      |
|                          |          |           |                                      |                      |                     |              | Tetryl                     | 10 ug/mL      |
| .8330_NG_Stk_00118       | 12/24/23 |           | Restek, Lot A0188553                 |                      | (Purchased Reagent) |              | Nitroglycerin              | 5000 ug/mL    |
| .8330_NG_Stk_00120       | 01/21/24 |           | Restek, Lot A0188553                 |                      | (Purchased Reagent) |              | Nitroglycerin              | 5000 ug/mL    |
| .8330_PETN_Stk_00128     | 01/21/24 |           | Restek, Lot A0187506                 |                      | (Purchased Reagent) |              | PETN                       | 5000 ug/mL    |
| .8330LCSMxl_00136        | 10/19/23 |           | Restek, Lot A0171502                 |                      | (Purchased Reagent) |              | 1,3,5-Trinitrobenzene      | 1000 ug/mL    |
|                          |          |           |                                      |                      |                     |              | 1,3-Dinitrobenzene         | 1000 ug/mL    |
|                          |          |           |                                      |                      |                     |              | 2,4,6-Trinitrotoluene      | 1000 ug/mL    |
|                          |          |           |                                      |                      |                     |              | 2,4-Dinitrotoluene         | 1000 ug/mL    |
|                          |          |           |                                      |                      |                     |              | HMX                        | 1000 ug/mL    |
|                          |          |           |                                      |                      |                     |              | Nitrobenzene               | 1000 ug/mL    |
|                          |          |           |                                      |                      |                     |              | RDX                        | 1000 ug/mL    |
| .8330LCSmix2_00034       | 01/21/24 |           | Restek, Lot A0186475                 |                      | (Purchased Reagent) |              | 2,6-Dinitrotoluene         | 1000 ug/mL    |



# REAGENT TRACEABILITY SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

| Reagent ID           | Exp Date | Prep Date | Dilutant Used                        | Reagent Final Volume | Parent Reagent      |              | Analyte                    | Concentration |
|----------------------|----------|-----------|--------------------------------------|----------------------|---------------------|--------------|----------------------------|---------------|
|                      |          |           |                                      |                      | Reagent ID          | Volume Added |                            |               |
|                      |          |           |                                      |                      |                     |              | 2-Amino-4,6-dinitrotoluene | 1000 ug/mL    |
|                      |          |           |                                      |                      |                     |              | 2-Nitrotoluene             | 1000 ug/mL    |
|                      |          |           |                                      |                      |                     |              | 3-Nitrotoluene             | 1000 ug/mL    |
|                      |          |           |                                      |                      |                     |              | 4-Amino-2,6-dinitrotoluene | 1000 ug/mL    |
|                      |          |           |                                      |                      |                     |              | 4-Nitrotoluene             | 1000 ug/mL    |
|                      |          |           |                                      |                      |                     |              | Tetryl                     | 1000 ug/mL    |
| 8330 LCS_00125       | 12/23/23 | 04/03/23  | Acetonitrile, Lot Acetonitrile_00073 | 100 mL               | 3,5-DNA Stock_00040 | 1 mL         | 3,5-Dinitroaniline         | 10 ug/mL      |
|                      |          |           |                                      |                      | 8330 NG Stk 00128   | 1 mL         | Nitroglycerin              | 100 ug/mL     |
|                      |          |           |                                      |                      | 8330 NG Stk 00129   | 1 mL         | Nitroglycerin              | 100 ug/mL     |
|                      |          |           |                                      |                      | 8330 PETN Stk 00133 | 1 mL         | PETN                       | 100 ug/mL     |
|                      |          |           |                                      |                      | 8330 PETN Stk 00134 | 1 mL         | PETN                       | 100 ug/mL     |
|                      |          |           |                                      |                      | 8330LCSMix1_00140   | 1 mL         | 1,3,5-Trinitrobenzene      | 10 ug/mL      |
|                      |          |           |                                      |                      |                     |              | 1,3-Dinitrobenzene         | 10 ug/mL      |
|                      |          |           |                                      |                      |                     |              | 2,4,6-Trinitrotoluene      | 10 ug/mL      |
|                      |          |           |                                      |                      |                     |              | 2,4-Dinitrotoluene         | 10 ug/mL      |
|                      |          |           |                                      |                      |                     |              | HMX                        | 10 ug/mL      |
|                      |          |           |                                      |                      |                     |              | Nitrobenzene               | 10 ug/mL      |
|                      |          |           |                                      |                      |                     |              | RDX                        | 10 ug/mL      |
|                      |          |           |                                      |                      | 8330LCSmix2_00038   | 1 mL         | 2,6-Dinitrotoluene         | 10 ug/mL      |
|                      |          |           |                                      |                      |                     |              | 2-Amino-4,6-dinitrotoluene | 10 ug/mL      |
|                      |          |           |                                      |                      |                     |              | 2-Nitrotoluene             | 10 ug/mL      |
|                      |          |           |                                      |                      |                     |              | 3-Nitrotoluene             | 10 ug/mL      |
|                      |          |           |                                      |                      |                     |              | 4-Amino-2,6-dinitrotoluene | 10 ug/mL      |
|                      |          |           |                                      |                      |                     |              | 4-Nitrotoluene             | 10 ug/mL      |
|                      |          |           |                                      |                      |                     |              | Tetryl                     | 10 ug/mL      |
|                      |          |           |                                      |                      | PicricARestek 00113 | 1 mL         | 2,4,6-Trinitrophenol       | 10 ug/mL      |
| .3,5-DNA Stock 00040 | 04/03/24 |           | Restek, Lot A0172723                 |                      | (Purchased Reagent) |              | 3,5-Dinitroaniline         | 1000 ug/mL    |
| .8330 NG Stk 00128   | 04/03/24 |           | Restek, Lot A0194013                 |                      | (Purchased Reagent) |              | Nitroglycerin              | 5000 ug/mL    |
| .8330 NG Stk 00129   | 04/03/24 |           | Restek, Lot A0194013                 |                      | (Purchased Reagent) |              | Nitroglycerin              | 5000 ug/mL    |
| .8330 PETN Stk 00133 | 04/03/24 |           | Restek, Lot A0188550                 |                      | (Purchased Reagent) |              | PETN                       | 5000 ug/mL    |
| .8330 PETN Stk 00134 | 04/03/24 |           | Restek, Lot A0188550                 |                      | (Purchased Reagent) |              | PETN                       | 5000 ug/mL    |
| .8330LCSMix1_00140   | 04/03/24 |           | Restek, Lot A0183848                 |                      | (Purchased Reagent) |              | 1,3,5-Trinitrobenzene      | 1000 ug/mL    |
|                      |          |           |                                      |                      |                     |              | 1,3-Dinitrobenzene         | 1000 ug/mL    |
|                      |          |           |                                      |                      |                     |              | 2,4,6-Trinitrotoluene      | 1000 ug/mL    |
|                      |          |           |                                      |                      |                     |              | 2,4-Dinitrotoluene         | 1000 ug/mL    |
|                      |          |           |                                      |                      |                     |              | HMX                        | 1000 ug/mL    |
|                      |          |           |                                      |                      |                     |              | Nitrobenzene               | 1000 ug/mL    |
|                      |          |           |                                      |                      |                     |              | RDX                        | 1000 ug/mL    |
| .8330LCSmix2_00038   | 04/03/24 |           | Restek, Lot A0192237                 |                      | (Purchased Reagent) |              | 2,6-Dinitrotoluene         | 1000 ug/mL    |
|                      |          |           |                                      |                      |                     |              | 2-Amino-4,6-dinitrotoluene | 1000 ug/mL    |
|                      |          |           |                                      |                      |                     |              | 2-Nitrotoluene             | 1000 ug/mL    |
|                      |          |           |                                      |                      |                     |              | 3-Nitrotoluene             | 1000 ug/mL    |
|                      |          |           |                                      |                      |                     |              | 4-Amino-2,6-dinitrotoluene | 1000 ug/mL    |
|                      |          |           |                                      |                      |                     |              | 4-Nitrotoluene             | 1000 ug/mL    |
|                      |          |           |                                      |                      |                     |              | Tetryl                     | 1000 ug/mL    |
| .PicricARestek_00113 | 04/03/24 |           | Restek, Lot A0183202                 |                      | (Purchased Reagent) |              | 2,4,6-Trinitrophenol       | 1000 ug/mL    |

## REAGENT TRACEABILITY SUMMARY

Lab Name: Eurofins Denver

Job No.: 280-176674-1

SDG No.:

| Reagent ID           | Exp Date | Prep Date | Dilutant Used                        | Reagent Final Volume | Parent Reagent      |              | Analyte                    | Concentration |
|----------------------|----------|-----------|--------------------------------------|----------------------|---------------------|--------------|----------------------------|---------------|
|                      |          |           |                                      |                      | Reagent ID          | Volume Added |                            |               |
| 8330 LCS_00126       | 04/26/24 | 04/26/23  | Acetonitrile, Lot Acetonitrile_00077 | 100 mL               | 3,5-DNA Stock_00042 | 1 mL         | 3,5-Dinitroaniline         | 10 ug/mL      |
|                      |          |           |                                      |                      | 8330 NG_Stk_00130   | 1 mL         | Nitroglycerin              | 100 ug/mL     |
|                      |          |           |                                      |                      | 8330 NG_Stk_00131   | 1 mL         | Nitroglycerin              | 100 ug/mL     |
|                      |          |           |                                      |                      | 8330 PETN_Stk_00135 | 1 mL         | PETN                       | 100 ug/mL     |
|                      |          |           |                                      |                      | 8330 PETN_Stk_00136 | 1 mL         | PETN                       | 100 ug/mL     |
|                      |          |           |                                      |                      | 8330LCSMix1_00141   | 1 mL         | 1,3,5-Trinitrobenzene      | 10 ug/mL      |
|                      |          |           |                                      |                      |                     |              | 1,3-Dinitrobenzene         | 10 ug/mL      |
|                      |          |           |                                      |                      |                     |              | 2,4,6-Trinitrotoluene      | 10 ug/mL      |
|                      |          |           |                                      |                      |                     |              | 2,4-Dinitrotoluene         | 10 ug/mL      |
|                      |          |           |                                      |                      |                     |              | HMX                        | 10 ug/mL      |
|                      |          |           |                                      |                      |                     |              | Nitrobenzene               | 10 ug/mL      |
|                      |          |           |                                      |                      |                     |              | RDX                        | 10 ug/mL      |
|                      |          |           |                                      |                      | 8330LCSmix2_00039   | 1 mL         | 2,6-Dinitrotoluene         | 10 ug/mL      |
|                      |          |           |                                      |                      |                     |              | 2-Amino-4,6-dinitrotoluene | 10 ug/mL      |
|                      |          |           |                                      |                      |                     |              | 2-Nitrotoluene             | 10 ug/mL      |
|                      |          |           |                                      |                      |                     |              | 3-Nitrotoluene             | 10 ug/mL      |
|                      |          |           |                                      |                      |                     |              | 4-Amino-2,6-dinitrotoluene | 10 ug/mL      |
|                      |          |           |                                      |                      |                     |              | 4-Nitrotoluene             | 10 ug/mL      |
|                      |          |           |                                      |                      |                     |              | Tetryl                     | 10 ug/mL      |
| .3,5-DNA Stock_00042 | 04/26/24 |           | Restek, Lot A0185772                 |                      | (Purchased Reagent) |              | 3,5-Dinitroaniline         | 1000 ug/mL    |
| .8330 NG_Stk_00130   | 04/26/24 |           | Restek, Lot A0194013                 |                      | (Purchased Reagent) |              | Nitroglycerin              | 5000 ug/mL    |
| .8330 NG_Stk_00131   | 04/26/24 |           | Restek, Lot A0194013                 |                      | (Purchased Reagent) |              | Nitroglycerin              | 5000 ug/mL    |
| .8330 PETN_Stk_00135 | 04/26/24 |           | Restek, Lot A0188550                 |                      | (Purchased Reagent) |              | PETN                       | 5000 ug/mL    |
| .8330 PETN_Stk_00136 | 04/26/24 |           | Restek, Lot A0188550                 |                      | (Purchased Reagent) |              | PETN                       | 5000 ug/mL    |
| .8330LCSMix1_00141   | 04/26/24 |           | Restek, Lot A0183848                 |                      | (Purchased Reagent) |              | 1,3,5-Trinitrobenzene      | 1000 ug/mL    |
|                      |          |           |                                      |                      |                     |              | 1,3-Dinitrobenzene         | 1000 ug/mL    |
|                      |          |           |                                      |                      |                     |              | 2,4,6-Trinitrotoluene      | 1000 ug/mL    |
|                      |          |           |                                      |                      |                     |              | 2,4-Dinitrotoluene         | 1000 ug/mL    |
|                      |          |           |                                      |                      |                     |              | HMX                        | 1000 ug/mL    |
|                      |          |           |                                      |                      |                     |              | Nitrobenzene               | 1000 ug/mL    |
| .8330LCSmix2_00039   | 04/26/24 |           | Restek, Lot A0192237                 |                      | (Purchased Reagent) |              | RDX                        | 1000 ug/mL    |
|                      |          |           |                                      |                      |                     |              | 2,6-Dinitrotoluene         | 1000 ug/mL    |
|                      |          |           |                                      |                      |                     |              | 2-Amino-4,6-dinitrotoluene | 1000 ug/mL    |
|                      |          |           |                                      |                      |                     |              | 2-Nitrotoluene             | 1000 ug/mL    |
|                      |          |           |                                      |                      |                     |              | 3-Nitrotoluene             | 1000 ug/mL    |
|                      |          |           |                                      |                      |                     |              | 4-Amino-2,6-dinitrotoluene | 1000 ug/mL    |
| 8330_ADDs_00035      | 07/13/23 | 01/13/23  | Acetonitrile, Lot ACN_237            | 5 mL                 | 2,4Diamino6NT_00094 | 1 mL         | 2,4-diamino-6-nitrotoluene | 20 ug/mL      |
|                      |          |           |                                      |                      | 2,6Diamino4NT_00095 | 1 mL         | 2,6-diamino-4-nitrotoluene | 20 ug/mL      |
|                      |          |           |                                      |                      | 833035DNASTk_00052  | 1 mL         | 3,5-Dinitroaniline         | 20 ug/mL      |
| .2,4Diamino6NT_00094 | 04/10/24 |           | AccuStandard, Lot 220121267-01       |                      | (Purchased Reagent) |              | 2,4-diamino-6-nitrotoluene | 100 ug/mL     |
| .2,6Diamino4NT_00095 | 06/03/24 |           | AccuStandard, Lot 222051002          |                      | (Purchased Reagent) |              | 2,6-diamino-4-nitrotoluene | 100 ug/mL     |
| .833035DNASTk_00052  | 10/06/23 |           | Accustandard, Lot 222011692-01       |                      | (Purchased Reagent) |              | 3,5-Dinitroaniline         | 100 ug/mL     |
| 8330_ADDs_00036      | 07/13/23 | 03/29/23  | Acetonitrile, Lot ACN_237            | 5 mL                 | 2,4Diamino6NT_00094 | 1 mL         | 2,4-diamino-6-nitrotoluene | 20 ug/mL      |

## REAGENT TRACEABILITY SUMMARY

Lab Name: Eurofins Denver

Job No.: 280-176674-1

SDG No.:

| Reagent ID                 | Exp Date | Prep Date | Dilutant Used                        | Reagent Final Volume | Parent Reagent      |              | Analyte                    | Concentration |
|----------------------------|----------|-----------|--------------------------------------|----------------------|---------------------|--------------|----------------------------|---------------|
|                            |          |           |                                      |                      | Reagent ID          | Volume Added |                            |               |
|                            |          |           |                                      |                      | 2,6Diamino4NT 00095 | 1 mL         | 2,6-diamino-4-nitrotoluene | 20 ug/mL      |
| .2,4Diamino6NT 00094       | 04/10/24 |           | AccuStandard, Lot 220121267-01       |                      | (Purchased Reagent) |              | 2,4-diamino-6-nitrotoluene | 100 ug/mL     |
| .2,6Diamino4NT 00095       | 06/03/24 |           | AccuStandard, Lot 222051002          |                      | (Purchased Reagent) |              | 2,6-diamino-4-nitrotoluene | 100 ug/mL     |
| <b>8330_ADDs_00037</b>     | 11/04/23 | 05/04/23  | Acetonitrile, Lot ACN_237            | 5 mL                 | 2,4Diamino6NT_00095 | 1 mL         | 2,4-diamino-6-nitrotoluene | 20 ug/mL      |
|                            |          |           |                                      |                      | 2,6Diamino4NT 00096 | 1 mL         | 2,6-diamino-4-nitrotoluene | 20 ug/mL      |
| .2,4Diamino6NT 00095       | 04/10/24 |           | AccuStandard, Lot 220121267-01       |                      | (Purchased Reagent) |              | 2,4-diamino-6-nitrotoluene | 100 ug/mL     |
| .2,6Diamino4NT 00096       | 05/04/24 |           | AccuStandard, Lot 222051002          |                      | (Purchased Reagent) |              | 2,6-diamino-4-nitrotoluene | 100 ug/mL     |
| <b>8330_OP_DMT_00016</b>   | 02/28/23 | 12/30/22  | Acetonitrile, Lot Acetonitrile_00073 | 10 mL                | MNX,TNX,DNX_00069   | 1 mL         | DNX                        | 10.04 ug/mL   |
|                            |          |           |                                      |                      |                     |              | MNX                        | 11.7 ug/mL    |
|                            |          |           |                                      |                      |                     |              | TNX                        | 10 ug/mL      |
| .MNX,TNX,DNX_00069         | 02/28/23 |           | Agilent, Lot 0006659035              |                      | (Purchased Reagent) |              | DNX                        | 100.4 ug/mL   |
|                            |          |           |                                      |                      |                     |              | MNX                        | 117 ug/mL     |
|                            |          |           |                                      |                      |                     |              | TNX                        | 100 ug/mL     |
| <b>8330_OP_DMT_00017</b>   | 08/28/23 | 02/28/23  | Acetonitrile, Lot Acetonitrile_00073 | 10 mL                | MNX,TNX,DNX_00077   | 1 mL         | DNX                        | 10.01 ug/mL   |
|                            |          |           |                                      |                      |                     |              | MNX                        | 11.67 ug/mL   |
|                            |          |           |                                      |                      |                     |              | TNX                        | 10.03 ug/mL   |
| .MNX,TNX,DNX_00077         | 02/29/24 |           | Agilent, Lot 0006725091              |                      | (Purchased Reagent) |              | DNX                        | 100.1 ug/mL   |
|                            |          |           |                                      |                      |                     |              | MNX                        | 116.7 ug/mL   |
|                            |          |           |                                      |                      |                     |              | TNX                        | 100.3 ug/mL   |
| <b>8330_OP_DMT_00020</b>   | 02/29/24 | 05/09/23  | Acetonitrile, Lot Acetonitrile_00077 | 10 mL                | MNX,TNX,DNX_00079   | 1 mL         | DNX                        | 10.01 ug/mL   |
|                            |          |           |                                      |                      |                     |              | MNX                        | 11.67 ug/mL   |
|                            |          |           |                                      |                      |                     |              | TNX                        | 10.03 ug/mL   |
| .MNX,TNX,DNX_00079         | 02/29/24 |           | Agilent, Lot 0006725091              |                      | (Purchased Reagent) |              | DNX                        | 100.1 ug/mL   |
|                            |          |           |                                      |                      |                     |              | MNX                        | 116.7 ug/mL   |
|                            |          |           |                                      |                      |                     |              | TNX                        | 100.3 ug/mL   |
| <b>8330DiaminLCS_00052</b> | 01/16/23 | 12/19/22  | Acetonitrile, Lot Acetonitrile_00069 | 10 mL                | 2,4Diamino6NT_00092 | 1 mL         | 2,4-diamino-6-nitrotoluene | 10 ug/mL      |
|                            |          |           |                                      |                      | 2,6Diamino4NT 00099 | 1 mL         | 2,6-diamino-4-nitrotoluene | 10 ug/mL      |
| .2,4Diamino6NT 00092       | 01/16/23 |           | AccuStandard, Lot 220121267          |                      | (Purchased Reagent) |              | 2,4-diamino-6-nitrotoluene | 100 ug/mL     |
| .2,6Diamino4NT 00099       | 12/19/23 |           | AccuStandard, Lot 222051002          |                      | (Purchased Reagent) |              | 2,6-diamino-4-nitrotoluene | 100 ug/mL     |
| <b>8330DiaminLCS_00056</b> | 01/19/24 | 01/27/23  | Acetonitrile, Lot Acetonitrile_00074 | 20 mL                | 2,4Diamino6NT_00098 | 1 mL         | 2,4-diamino-6-nitrotoluene | 10 ug/mL      |
|                            |          |           |                                      |                      | 2,4Diamino6NT 00099 | 1 mL         | 2,4-diamino-6-nitrotoluene | 10 ug/mL      |
|                            |          |           |                                      |                      | 2,6Diamino4NT 00100 | 1 mL         | 2,6-diamino-4-nitrotoluene | 10 ug/mL      |
|                            |          |           |                                      |                      | 2,6Diamino4NT 00102 | 1 mL         | 2,6-diamino-4-nitrotoluene | 10 ug/mL      |
| .2,4Diamino6NT 00098       | 04/10/24 |           | AccuStandard, Lot 220121267-01       |                      | (Purchased Reagent) |              | 2,4-diamino-6-nitrotoluene | 100 ug/mL     |
| .2,4Diamino6NT 00099       | 04/10/24 |           | AccuStandard, Lot 220121267-01       |                      | (Purchased Reagent) |              | 2,4-diamino-6-nitrotoluene | 100 ug/mL     |
| .2,6Diamino4NT 00100       | 06/03/24 |           | AccuStandard, Lot 222051002          |                      | (Purchased Reagent) |              | 2,6-diamino-4-nitrotoluene | 100 ug/mL     |
| .2,6Diamino4NT 00102       | 06/03/24 |           | AccuStandard, Lot 222051002          |                      | (Purchased Reagent) |              | 2,6-diamino-4-nitrotoluene | 100 ug/mL     |
| <b>8330DiaminLCS_00057</b> | 04/10/24 | 05/19/23  | Acetonitrile, Lot Acetonitrile_00078 | 10 mL                | 2,4Diamino6NT_00100 | 1 mL         | 2,4-diamino-6-nitrotoluene | 10 ug/mL      |
|                            |          |           |                                      |                      | 2,6Diamino4NT 00105 | 1 mL         | 2,6-diamino-4-nitrotoluene | 10 ug/mL      |

REAGENT TRACEABILITY SUMMARY

Lab Name: Eurofins Denver

Job No.: 280-176674-1

SDG No.:

| Reagent ID            | Exp Date | Prep Date | Dilutant Used                  | Reagent Final Volume | Parent Reagent      |              | Analyte                    | Concentration |
|-----------------------|----------|-----------|--------------------------------|----------------------|---------------------|--------------|----------------------------|---------------|
|                       |          |           |                                |                      | Reagent ID          | Volume Added |                            |               |
| .2,4Diamino6NT_00100  | 04/10/24 |           | AccuStandard, Lot 220121267-01 |                      | (Purchased Reagent) |              | 2,4-diamino-6-nitrotoluene | 100 ug/mL     |
| .2,6Diamino4NT_00105  | 05/19/24 |           | AccuStandard, Lot 222051002    |                      | (Purchased Reagent) |              | 2,6-diamino-4-nitrotoluene | 100 ug/mL     |
| 8330IntermStk_00075   | 08/07/23 | 02/03/23  | Acetonitrile, Lot ACN_238      | 10 mL                | 8330_NG1000_00007   | 1 mL         | Nitroglycerin              | 100 ug/mL     |
|                       |          |           |                                |                      | 8330_PETN1000_00011 | 1 mL         | PETN                       | 100 ug/mL     |
|                       |          |           |                                |                      | 833035DNASTk_00053  | 1 mL         | 3,5-Dinitroaniline         | 10 ug/mL      |
|                       |          |           |                                |                      | 8330ICALStock_00034 | 1 mL         | 1,3,5-Trinitrobenzene      | 10 ug/mL      |
|                       |          |           |                                |                      |                     |              | 1,3-Dinitrobenzene         | 10 ug/mL      |
|                       |          |           |                                |                      |                     |              | 2,4,6-Trinitrotoluene      | 10 ug/mL      |
|                       |          |           |                                |                      |                     |              | 2,4-Dinitrotoluene         | 10 ug/mL      |
|                       |          |           |                                |                      |                     |              | 2,6-Dinitrotoluene         | 10 ug/mL      |
|                       |          |           |                                |                      |                     |              | 2-Amino-4,6-dinitrotoluene | 10 ug/mL      |
|                       |          |           |                                |                      |                     |              | 2-Nitrotoluene             | 10 ug/mL      |
|                       |          |           |                                |                      |                     |              | 3-Nitrotoluene             | 10 ug/mL      |
|                       |          |           |                                |                      |                     |              | 4-Amino-2,6-dinitrotoluene | 10 ug/mL      |
|                       |          |           |                                |                      |                     |              | 4-Nitrotoluene             | 10 ug/mL      |
|                       |          |           |                                |                      |                     |              | HMX                        | 10 ug/mL      |
|                       |          |           |                                |                      |                     |              | Nitrobenzene               | 10 ug/mL      |
|                       |          |           |                                |                      |                     |              | RDX                        | 10 ug/mL      |
|                       |          |           |                                |                      |                     |              | Tetryl                     | 10 ug/mL      |
|                       |          |           |                                |                      |                     |              | 1,2-Dinitrobenzene         | 10 ug/mL      |
|                       |          |           |                                |                      | 8330PASTkPS_00070   | 1 mL         | 2,4,6-Trinitrophenol       | 10 ug/mL      |
| .8330 NG1000_00007    | 02/03/24 |           | Restek, Lot A0175997           |                      | (Purchased Reagent) |              | Nitroglycerin              | 1000 ug/mL    |
| .8330 PETN1000_00011  | 02/03/24 |           | Restek, Lot A0187142           |                      | (Purchased Reagent) |              | PETN                       | 1000 ug/mL    |
| .833035DNASTk_00053   | 10/06/23 |           | Accustandard, Lot 222011692-01 |                      | (Purchased Reagent) |              | 3,5-Dinitroaniline         | 100 ug/mL     |
| .8330ICALStock_00034  | 02/03/24 | 02/03/23  | Acetonitrile, Lot ACN_238      | 10 mL                | 8330 Stock_TS_00023 | 1 mL         | 1,3,5-Trinitrobenzene      | 100 ug/mL     |
|                       |          |           |                                |                      |                     |              | 1,3-Dinitrobenzene         | 100 ug/mL     |
|                       |          |           |                                |                      |                     |              | 2,4,6-Trinitrotoluene      | 100 ug/mL     |
|                       |          |           |                                |                      |                     |              | 2,4-Dinitrotoluene         | 100 ug/mL     |
|                       |          |           |                                |                      |                     |              | 2,6-Dinitrotoluene         | 100 ug/mL     |
|                       |          |           |                                |                      |                     |              | 2-Amino-4,6-dinitrotoluene | 100 ug/mL     |
|                       |          |           |                                |                      |                     |              | 2-Nitrotoluene             | 100 ug/mL     |
|                       |          |           |                                |                      |                     |              | 3-Nitrotoluene             | 100 ug/mL     |
|                       |          |           |                                |                      |                     |              | 4-Amino-2,6-dinitrotoluene | 100 ug/mL     |
|                       |          |           |                                |                      |                     |              | 4-Nitrotoluene             | 100 ug/mL     |
|                       |          |           |                                |                      |                     |              | HMX                        | 100 ug/mL     |
|                       |          |           |                                |                      |                     |              | Nitrobenzene               | 100 ug/mL     |
|                       |          |           |                                |                      |                     |              | RDX                        | 100 ug/mL     |
|                       |          |           |                                |                      |                     |              | Tetryl                     | 100 ug/mL     |
|                       |          |           |                                |                      | 8330SurrStock_00172 | 1 mL         | 1,2-Dinitrobenzene         | 100 ug/mL     |
| ..8330 Stock_TS_00023 | 02/03/24 |           | Agilent, Lot 0006684308        |                      | (Purchased Reagent) |              | 1,3,5-Trinitrobenzene      | 1000 ug/mL    |
|                       |          |           |                                |                      |                     |              | 1,3-Dinitrobenzene         | 1000 ug/mL    |
|                       |          |           |                                |                      |                     |              | 2,4,6-Trinitrotoluene      | 1000 ug/mL    |
|                       |          |           |                                |                      |                     |              | 2,4-Dinitrotoluene         | 1000 ug/mL    |
|                       |          |           |                                |                      |                     |              | 2,6-Dinitrotoluene         | 1000 ug/mL    |
|                       |          |           |                                |                      |                     |              | 2-Amino-4,6-dinitrotoluene | 1000 ug/mL    |
|                       |          |           |                                |                      |                     |              | 2-Nitrotoluene             | 1000 ug/mL    |

# REAGENT TRACEABILITY SUMMARY

Lab Name: Eurofins Denver

Job No.: 280-176674-1

SDG No.:

| Reagent ID                 | Exp Date | Prep Date | Dilutant Used                  | Reagent Final Volume | Parent Reagent      |              | Analyte                    | Concentration |
|----------------------------|----------|-----------|--------------------------------|----------------------|---------------------|--------------|----------------------------|---------------|
|                            |          |           |                                |                      | Reagent ID          | Volume Added |                            |               |
|                            |          |           |                                |                      |                     |              | 3-Nitrotoluene             | 1000 ug/mL    |
|                            |          |           |                                |                      |                     |              | 4-Amino-2,6-dinitrotoluene | 1000 ug/mL    |
|                            |          |           |                                |                      |                     |              | 4-Nitrotoluene             | 1000 ug/mL    |
|                            |          |           |                                |                      |                     |              | HMX                        | 1000 ug/mL    |
|                            |          |           |                                |                      |                     |              | Nitrobenzene               | 1000 ug/mL    |
|                            |          |           |                                |                      |                     |              | RDX                        | 1000 ug/mL    |
|                            |          |           |                                |                      |                     |              | Tetryl                     | 1000 ug/mL    |
| ..8330SurrStock_00172      | 02/03/24 |           | AccuStandard, Lot 219051500    |                      | (Purchased Reagent) |              | 1,2-Dinitrobenzene         | 1000 ug/mL    |
| .8330PASTkPS_00070         | 08/07/23 |           | AccuStandard, Lot 218031154-03 |                      | (Purchased Reagent) |              | 2,4,6-Trinitrophenol       | 100 ug/mL     |
| <b>8330IntermStk_00076</b> | 09/23/23 | 03/23/23  | Acetonitrile, Lot ACN_237      | 10 mL                | 8330_NG1000_00010   | 1 mL         | Nitroglycerin              | 100 ug/mL     |
|                            |          |           |                                |                      | 8330_PETN1000_00010 | 1 mL         | PETN                       | 100 ug/mL     |
|                            |          |           |                                |                      | 833035DNASTk_00054  | 1 mL         | 3,5-Dinitroaniline         | 10 ug/mL      |
|                            |          |           |                                |                      | 8330ICALStock_00034 | 1 mL         | 1,3,5-Trinitrobenzene      | 10 ug/mL      |
|                            |          |           |                                |                      |                     |              | 1,3-Dinitrobenzene         | 10 ug/mL      |
|                            |          |           |                                |                      |                     |              | 2,4,6-Trinitrotoluene      | 10 ug/mL      |
|                            |          |           |                                |                      |                     |              | 2,4-Dinitrotoluene         | 10 ug/mL      |
|                            |          |           |                                |                      |                     |              | 2,6-Dinitrotoluene         | 10 ug/mL      |
|                            |          |           |                                |                      |                     |              | 2-Amino-4,6-dinitrotoluene | 10 ug/mL      |
|                            |          |           |                                |                      |                     |              | 2-Nitrotoluene             | 10 ug/mL      |
|                            |          |           |                                |                      |                     |              | 3-Nitrotoluene             | 10 ug/mL      |
|                            |          |           |                                |                      |                     |              | 4-Amino-2,6-dinitrotoluene | 10 ug/mL      |
|                            |          |           |                                |                      |                     |              | 4-Nitrotoluene             | 10 ug/mL      |
|                            |          |           |                                |                      |                     |              | HMX                        | 10 ug/mL      |
|                            |          |           |                                |                      |                     |              | Nitrobenzene               | 10 ug/mL      |
|                            |          |           |                                |                      |                     |              | RDX                        | 10 ug/mL      |
|                            |          |           |                                |                      |                     |              | Tetryl                     | 10 ug/mL      |
|                            |          |           |                                |                      |                     |              | 1,2-Dinitrobenzene         | 10 ug/mL      |
|                            |          |           |                                |                      | 8330PASTkPS_00071   | 1 mL         | 2,4,6-Trinitrophenol       | 10 ug/mL      |
| .8330_NG1000_00010         | 03/23/24 |           | Restek, Lot A0187280           |                      | (Purchased Reagent) |              | Nitroglycerin              | 1000 ug/mL    |
| .8330_PETN1000_00010       | 03/23/24 |           | Restek, Lot A0187142           |                      | (Purchased Reagent) |              | PETN                       | 1000 ug/mL    |
| .833035DNASTk_00054        | 10/06/23 |           | Accustandard, Lot 222011692-01 |                      | (Purchased Reagent) |              | 3,5-Dinitroaniline         | 100 ug/mL     |
| .8330ICALStock_00034       | 02/03/24 | 02/03/23  | Acetonitrile, Lot ACN_238      | 10 mL                | 8330_Stock_TS_00023 | 1 mL         | 1,3,5-Trinitrobenzene      | 100 ug/mL     |
|                            |          |           |                                |                      |                     |              | 1,3-Dinitrobenzene         | 100 ug/mL     |
|                            |          |           |                                |                      |                     |              | 2,4,6-Trinitrotoluene      | 100 ug/mL     |
|                            |          |           |                                |                      |                     |              | 2,4-Dinitrotoluene         | 100 ug/mL     |
|                            |          |           |                                |                      |                     |              | 2,6-Dinitrotoluene         | 100 ug/mL     |
|                            |          |           |                                |                      |                     |              | 2-Amino-4,6-dinitrotoluene | 100 ug/mL     |
|                            |          |           |                                |                      |                     |              | 2-Nitrotoluene             | 100 ug/mL     |
|                            |          |           |                                |                      |                     |              | 3-Nitrotoluene             | 100 ug/mL     |
|                            |          |           |                                |                      |                     |              | 4-Amino-2,6-dinitrotoluene | 100 ug/mL     |
|                            |          |           |                                |                      |                     |              | 4-Nitrotoluene             | 100 ug/mL     |
|                            |          |           |                                |                      |                     |              | HMX                        | 100 ug/mL     |
|                            |          |           |                                |                      |                     |              | Nitrobenzene               | 100 ug/mL     |
|                            |          |           |                                |                      |                     |              | RDX                        | 100 ug/mL     |
|                            |          |           |                                |                      |                     |              | Tetryl                     | 100 ug/mL     |
|                            |          |           |                                |                      | 8330SurrStock_00172 | 1 mL         | 1,2-Dinitrobenzene         | 100 ug/mL     |

# REAGENT TRACEABILITY SUMMARY

Lab Name: Eurofins Denver

Job No.: 280-176674-1

SDG No.:

| Reagent ID                 | Exp Date | Prep Date | Dilutant Used                        | Reagent Final Volume | Parent Reagent      |              | Analyte                    | Concentration |
|----------------------------|----------|-----------|--------------------------------------|----------------------|---------------------|--------------|----------------------------|---------------|
|                            |          |           |                                      |                      | Reagent ID          | Volume Added |                            |               |
| ..8330 Stock_TS_00023      | 02/03/24 |           | Agilent, Lot 0006684308              |                      | (Purchased Reagent) |              | 1,3,5-Trinitrobenzene      | 1000 ug/mL    |
|                            |          |           |                                      |                      |                     |              | 1,3-Dinitrobenzene         | 1000 ug/mL    |
|                            |          |           |                                      |                      |                     |              | 2,4,6-Trinitrotoluene      | 1000 ug/mL    |
|                            |          |           |                                      |                      |                     |              | 2,4-Dinitrotoluene         | 1000 ug/mL    |
|                            |          |           |                                      |                      |                     |              | 2,6-Dinitrotoluene         | 1000 ug/mL    |
|                            |          |           |                                      |                      |                     |              | 2-Amino-4,6-dinitrotoluene | 1000 ug/mL    |
|                            |          |           |                                      |                      |                     |              | 2-Nitrotoluene             | 1000 ug/mL    |
|                            |          |           |                                      |                      |                     |              | 3-Nitrotoluene             | 1000 ug/mL    |
|                            |          |           |                                      |                      |                     |              | 4-Amino-2,6-dinitrotoluene | 1000 ug/mL    |
|                            |          |           |                                      |                      |                     |              | 4-Nitrotoluene             | 1000 ug/mL    |
|                            |          |           |                                      |                      |                     |              | HMX                        | 1000 ug/mL    |
|                            |          |           |                                      |                      |                     |              | Nitrobenzene               | 1000 ug/mL    |
|                            |          |           |                                      |                      |                     |              | RDX                        | 1000 ug/mL    |
|                            |          |           |                                      |                      |                     |              | Tetryl                     | 1000 ug/mL    |
| ..8330SurrStock_00172      | 02/03/24 |           | AccuStandard, Lot 219051500          |                      | (Purchased Reagent) |              | 1,2-Dinitrobenzene         | 1000 ug/mL    |
| .8330PASTkPS_00071         | 03/23/24 |           | AccuStandard, Lot 218031154-04       |                      | (Purchased Reagent) |              | 2,4,6-Trinitrophenol       | 100 ug/mL     |
| <b>8330Surrogate_00138</b> | 08/02/23 | 02/02/23  | Acetonitrile, Lot Acetonitrile_00074 | 500 mL               | 8330SurrStkSS_00245 | 1 mL         | 1,2-Dinitrobenzene         | 10 ug/mL      |
|                            |          |           |                                      |                      | 8330SurrStkSS_00246 | 1 mL         | 1,2-Dinitrobenzene         | 10 ug/mL      |
|                            |          |           |                                      |                      | 8330SurrStkSS_00247 | 1 mL         | 1,2-Dinitrobenzene         | 10 ug/mL      |
|                            |          |           |                                      |                      | 8330SurrStkSS_00248 | 1 mL         | 1,2-Dinitrobenzene         | 10 ug/mL      |
|                            |          |           |                                      |                      | 8330SurrStkSS_00249 | 1 mL         | 1,2-Dinitrobenzene         | 10 ug/mL      |
| .8330SurrStkSS_00245       | 02/02/24 |           | Restek, Lot A0192220                 |                      | (Purchased Reagent) |              | 1,2-Dinitrobenzene         | 1000 ug/mL    |
| .8330SurrStkSS_00246       | 02/02/24 |           | Restek, Lot A0192220                 |                      | (Purchased Reagent) |              | 1,2-Dinitrobenzene         | 1000 ug/mL    |
| .8330SurrStkSS_00247       | 02/02/24 |           | Restek, Lot A0192220                 |                      | (Purchased Reagent) |              | 1,2-Dinitrobenzene         | 1000 ug/mL    |
| .8330SurrStkSS_00248       | 02/02/24 |           | Restek, Lot A0192220                 |                      | (Purchased Reagent) |              | 1,2-Dinitrobenzene         | 1000 ug/mL    |
| .8330SurrStkSS_00249       | 02/02/24 |           | Restek, Lot A0192220                 |                      | (Purchased Reagent) |              | 1,2-Dinitrobenzene         | 1000 ug/mL    |
| <b>8330Surrogate_00143</b> | 10/12/23 | 04/12/23  | Acetonitrile, Lot Acetonitrile_00077 | 500 mL               | 8330SurrStkSS_00263 | 1 mL         | 1,2-Dinitrobenzene         | 10 ug/mL      |
|                            |          |           |                                      |                      | 8330SurrStkSS_00265 | 1 mL         | 1,2-Dinitrobenzene         | 10 ug/mL      |
|                            |          |           |                                      |                      | 8330SurrStkSS_00266 | 1 mL         | 1,2-Dinitrobenzene         | 10 ug/mL      |
|                            |          |           |                                      |                      | 8330SurrStkSS_00267 | 1 mL         | 1,2-Dinitrobenzene         | 10 ug/mL      |
|                            |          |           |                                      |                      | 8330SurrStkSS_00268 | 1 mL         | 1,2-Dinitrobenzene         | 10 ug/mL      |
| .8330SurrStkSS_00263       | 04/12/24 |           | Restek, Lot A0193444                 |                      | (Purchased Reagent) |              | 1,2-Dinitrobenzene         | 1000 ug/mL    |
| .8330SurrStkSS_00265       | 04/12/24 |           | Restek, Lot A0193444                 |                      | (Purchased Reagent) |              | 1,2-Dinitrobenzene         | 1000 ug/mL    |
| .8330SurrStkSS_00266       | 04/12/24 |           | Restek, Lot A0194831                 |                      | (Purchased Reagent) |              | 1,2-Dinitrobenzene         | 1000 ug/mL    |
| .8330SurrStkSS_00267       | 04/12/24 |           | Restek, Lot A0194831                 |                      | (Purchased Reagent) |              | 1,2-Dinitrobenzene         | 1000 ug/mL    |
| .8330SurrStkSS_00268       | 04/12/24 |           | Restek, Lot A0194831                 |                      | (Purchased Reagent) |              | 1,2-Dinitrobenzene         | 1000 ug/mL    |
| <b>8330Surrogate_00144</b> | 11/17/23 | 05/17/23  | Acetonitrile, Lot Acetonitrile_00078 | 500 mL               | 8330SurrStkSS_00269 | 1 mL         | 1,2-Dinitrobenzene         | 10 ug/mL      |
|                            |          |           |                                      |                      |                     |              | 1,2-Dinitrobenzene (Surr)  | 10 ug/mL      |
|                            |          |           |                                      |                      | 8330SurrStkSS_00270 | 1 mL         | 1,2-Dinitrobenzene         | 10 ug/mL      |
|                            |          |           |                                      |                      |                     |              | 1,2-Dinitrobenzene (Surr)  | 10 ug/mL      |
|                            |          |           |                                      |                      | 8330SurrStkSS_00271 | 1 mL         | 1,2-Dinitrobenzene         | 10 ug/mL      |
|                            |          |           |                                      |                      |                     |              | 1,2-Dinitrobenzene (Surr)  | 10 ug/mL      |
|                            |          |           |                                      |                      | 8330SurrStkSS_00272 | 1 mL         | 1,2-Dinitrobenzene         | 10 ug/mL      |
|                            |          |           |                                      |                      |                     |              | 1,2-Dinitrobenzene (Surr)  | 10 ug/mL      |

REAGENT TRACEABILITY SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

| Reagent ID                 | Exp Date | Prep Date | Dilutant Used                   | Reagent Final Volume | Parent Reagent      |              | Analyte                   | Concentration |
|----------------------------|----------|-----------|---------------------------------|----------------------|---------------------|--------------|---------------------------|---------------|
|                            |          |           |                                 |                      | Reagent ID          | Volume Added |                           |               |
|                            |          |           |                                 |                      | 8330SurrStkSS_00273 | 1 mL         | 1,2-Dinitrobenzene        | 10 ug/mL      |
|                            |          |           |                                 |                      |                     |              | 1,2-Dinitrobenzene (Surr) | 10 ug/mL      |
| .8330SurrStkSS_00269       | 05/17/24 |           | Restek, Lot A0194831            |                      | (Purchased Reagent) |              | 1,2-Dinitrobenzene        | 1000 ug/mL    |
|                            |          |           |                                 |                      |                     |              | 1,2-Dinitrobenzene (Surr) | 1000 ug/mL    |
| .8330SurrStkSS_00270       | 05/17/24 |           | Restek, Lot A0194831            |                      | (Purchased Reagent) |              | 1,2-Dinitrobenzene        | 1000 ug/mL    |
|                            |          |           |                                 |                      |                     |              | 1,2-Dinitrobenzene (Surr) | 1000 ug/mL    |
| .8330SurrStkSS_00271       | 05/17/24 |           | Restek, Lot A0194831            |                      | (Purchased Reagent) |              | 1,2-Dinitrobenzene        | 1000 ug/mL    |
|                            |          |           |                                 |                      |                     |              | 1,2-Dinitrobenzene (Surr) | 1000 ug/mL    |
| .8330SurrStkSS_00272       | 05/17/24 |           | Restek, Lot A0194831            |                      | (Purchased Reagent) |              | 1,2-Dinitrobenzene        | 1000 ug/mL    |
|                            |          |           |                                 |                      |                     |              | 1,2-Dinitrobenzene (Surr) | 1000 ug/mL    |
| .8330SurrStkSS_00273       | 05/17/24 |           | Restek, Lot A0197062            |                      | (Purchased Reagent) |              | 1,2-Dinitrobenzene        | 1000 ug/mL    |
|                            |          |           |                                 |                      |                     |              | 1,2-Dinitrobenzene (Surr) | 1000 ug/mL    |
| <b>Alk daily lcs 01088</b> | 05/26/23 | 05/19/23  | Di Water, Lot na                | 1000 mL              | Alk stk std 00023   | 4 mL         | Alkalinity                | 200 mg/L      |
| .Alk stk std 00023         | 10/30/23 |           | SPEX, Lot 2-106BT               |                      | (Purchased Reagent) |              | Alkalinity                | 50 g/L        |
| <b>IC CAL cl/so4_00480</b> | 05/24/23 | 05/18/23  | Di Water, Lot na                | 100 mL               | IC CL cal 00076     | 12.5 mL      | Chloride                  | 125 mg/L      |
|                            |          |           |                                 |                      | IC SO4 Cal 00008    | 12.5 mL      | Sulfate                   | 125 mg/L      |
| .IC CL cal 00076           | 11/30/23 |           | SPEX CertiPrep, Lot 5-146CL-2X  |                      | (Purchased Reagent) |              | Chloride                  | 1000 mg/L     |
| .IC SO4 Cal 00008          | 07/30/23 |           | SPEX CertiPrep, Lot 5-160SO4-2X |                      | (Purchased Reagent) |              | Sulfate                   | 1000 mg/L     |
| <b>IC Cal low_00709</b>    | 05/24/23 | 05/18/23  | Di Water, Lot NA                | 100 mL               | IC BR ICV 00024     | 2.5 mL       | Bromide                   | 25 mg/L       |
|                            |          |           |                                 |                      | IC FL cal 00023     | 2.5 mL       | Fluoride                  | 25 mg/L       |
|                            |          |           |                                 |                      | IC N03 cal 00028    | 2.5 mL       | Nitrate as N              | 25 mg/L       |
|                            |          |           |                                 |                      | IC NO2 Cal 00007    | 2.5 mL       | Nitrite as N              | 25 mg/L       |
|                            |          |           |                                 |                      | IC P04 cal 00030    | 2.5 mL       | Orthophosphate as P       | 25 mg/L       |
| .IC BR ICV 00024           | 10/31/23 |           | ricca, Lot 1204B86              |                      | (Purchased Reagent) |              | Bromide                   | 1000 mg/L     |
| .IC FL cal 00023           | 01/31/24 |           | Ricca, Lot 2207H13              |                      | (Purchased Reagent) |              | Fluoride                  | 1000 mg/L     |
| .IC N03 cal 00028          | 11/30/23 |           | Ricca, Lot 2205B73              |                      | (Purchased Reagent) |              | Nitrate as N              | 1000 mg/L     |
| .IC NO2 Cal 00007          | 06/30/23 |           | Ricca, Lot 2212826              |                      | (Purchased Reagent) |              | Nitrite as N              | 1000 mg/L     |
| .IC P04 cal 00030          | 07/31/24 |           | RICCA, Lot 4208152              |                      | (Purchased Reagent) |              | Orthophosphate as P       | 1000 mg/L     |
| <b>IC ICV 5_00405</b>      | 05/24/23 | 05/18/23  | Di Water, Lot na                | 100 mL               | IC NO2 ICV 00022    | 2.5 mL       | Nitrite as N              | 25 mg/L       |
|                            |          |           |                                 |                      | IC NO3 ICV 00019    | 2.5 mL       | Nitrate as N              | 25 mg/L       |
| .IC NO2 ICV 00022          | 12/14/23 |           | ERA, Lot 261221m                |                      | (Purchased Reagent) |              | Nitrite as N              | 1000 mg/L     |
| .IC NO3 ICV 00019          | 12/01/23 |           | ERA, Lot 011221m                |                      | (Purchased Reagent) |              | Nitrate as N              | 1000 mg/L     |
| <b>IC LCS_01954</b>        | 05/24/23 | 05/18/23  | Di Water, Lot 27                | 200 mL               | IC Cal low_00709    | 20 mL        | Nitrate as N              | 2.5 mg/L      |
|                            |          |           |                                 |                      |                     |              | Nitrite as N              | 2.5 mg/L      |
|                            |          |           |                                 |                      | IC SO4 Cal 00008    | 10 mL        | Sulfate                   | 50 mg/L       |
| .IC Cal low_00709          | 05/24/23 | 05/18/23  | Di Water, Lot NA                | 100 mL               | IC N03 cal 00028    | 2.5 mL       | Nitrate as N              | 25 mg/L       |
|                            |          |           |                                 |                      | IC NO2 Cal 00007    | 2.5 mL       | Nitrite as N              | 25 mg/L       |
| ..IC N03 cal 00028         | 11/30/23 |           | Ricca, Lot 2205B73              |                      | (Purchased Reagent) |              | Nitrate as N              | 1000 mg/L     |
| ..IC NO2 Cal 00007         | 06/30/23 |           | Ricca, Lot 2212826              |                      | (Purchased Reagent) |              | Nitrite as N              | 1000 mg/L     |
| .IC SO4 Cal 00008          | 07/30/23 |           | SPEX CertiPrep, Lot 5-160SO4-2X |                      | (Purchased Reagent) |              | Sulfate                   | 1000 mg/L     |
| <b>IC SO4 ICV 00024</b>    | 12/01/23 |           | ERA, Lot 021221m                |                      | (Purchased Reagent) |              | Sulfate                   | 1000 mg/L     |
| <b>SFD CAL INT_02403</b>   | 07/26/23 | 04/26/23  | ELGA Water, Lot NA              | 500 mL               | 50% NaOH 00023      | 2 mL         | Sodium Hydroxide          | 2000 mg/L     |
|                            |          |           |                                 |                      | SFD CAL STK 00010   | 3.5656 g     | Sulfide                   | 952.015 mg/L  |
| .50% NaOH 00023            | 01/31/24 |           | RICCA, Lot 1102857              |                      | (Purchased Reagent) |              | Sodium Hydroxide          | 50 %          |
| .SFD CAL STK 00010         | 08/25/25 |           | Thermo Scientific, Lot A0411500 |                      | (Purchased Reagent) |              | Sulfide                   | 0.1335 g/g    |

## REAGENT TRACEABILITY SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

| Reagent ID        | Exp Date | Prep Date                            | Dilutant Used | Reagent Final Volume | Parent Reagent      |              | Analyte                        | Concentration |
|-------------------|----------|--------------------------------------|---------------|----------------------|---------------------|--------------|--------------------------------|---------------|
|                   |          |                                      |               |                      | Reagent ID          | Volume Added |                                |               |
| TOC ICV Std_00051 | 07/01/23 | Ricca, Lot 2208C86                   |               |                      | (Purchased Reagent) |              | TOC Result 1                   | 1000 ppm      |
|                   |          |                                      |               |                      |                     |              | TOC Result 2                   | 1000 ppm      |
|                   |          |                                      |               |                      |                     |              | TOC Result 3                   | 1000 ppm      |
|                   |          |                                      |               |                      |                     |              | TOC Result 4                   | 1000 ppm      |
|                   |          |                                      |               |                      |                     |              | Total Organic Carbon - Average |               |
| TOC LCS Std_00057 | 01/31/25 | Agilent Technologies, Lot 0006719669 |               |                      | (Purchased Reagent) |              | TOC Result 1                   | 1001 ppm      |
|                   |          |                                      |               |                      |                     |              | TOC Result 2                   | 1001 ppm      |
|                   |          |                                      |               |                      |                     |              | TOC Result 3                   | 1001 ppm      |
|                   |          |                                      |               |                      |                     |              | TOC Result 4                   | 1001 ppm      |
|                   |          |                                      |               |                      |                     |              | Total Organic Carbon - Average |               |



Reagent

---

**2,4Diamino6NT\_00092**

# CERTIFICATE OF ANALYSIS

**Catalog No:** M-8330-ADD-12

**Description:** 2,4-Diamino-6-nitrotoluene

**Lot:** 220121267

**Solvent:** Acetonitrile

**Hazards:** Refer to SDS for complete safety information



**Signal Word:** Danger

**Date Certified:** Dec 16, 2020

**Expiration:** Jan 16, 2023

**Sample Size:** 1 mL

**Components:** 1

**Storage Condition:** Refrig (0-5 °C)

**Certified Reference Material**



| Component                  | CAS #     | Purity %<br>(LCMS) | Prepared<br>Concentration <sup>2</sup><br>(µg/mL) | Certified Analyte<br>Concentration <sup>1</sup><br>(µg/mL) |
|----------------------------|-----------|--------------------|---------------------------------------------------|------------------------------------------------------------|
| 2,4-Diamino-6-nitrotoluene | 6629-29-4 | 99.0               | 100.3                                             | 99.3                                                       |

A product with a suffix (-1A, -2B, etc. or -01, -02, etc.) on its lot number has had its expiration date extended and is identical to the same lot number without the suffix.

<sup>2</sup> All weights are traceable through NIST, Test No. 684/289871-17

<sup>1</sup> Certified Analyte Concentration = Purity x Prepared Concentration.

The Uncertainty associated with the certified concentration reported on this certificate is  $\pm 2.4\%$ . This value is the combined expanded uncertainty and represents an estimated standard deviation equal to the positive square root of the total variation of the uncertainty of components. A normal distribution is assumed and a coverage factor of K=2 is chosen using approximately a 95% confidence level.

Labels and certificates follow U.S. Conventions in reporting numerical values: A comma (,) is used to separate units of one-thousand or greater. A period (.) is used as a decimal place marker.

The information on this certificate may not be reproduced without the express permission of the manufacturer. See reverse side for additional information.

Hazard Information: Please refer to the SDS for information regarding the hazards associated with using this material.

This product was prepared according to in-house procedures and is guaranteed to be homogeneous.

Certified By:



Larry Decker, Organic QC Manager

For use in routine laboratory analysis.

Reagent

---

**2,4Diamino6NT\_00094**

# CERTIFICATE OF ANALYSIS

**Catalog No:** M-8330-ADD-12

**Description:** 2,4-Diamino-6-nitrotoluene

**Lot:** 220121267-01

**Solvent:** Acetonitrile

**Hazards:** Refer to SDS for complete safety information

**Date Certified:** Mar 10, 2022

**Expiration:** Apr 10, 2024

**Sample Size:** 1 mL

**Components:** 1

**Storage Condition:** Refrig (0-5 °C)



Signal Word: Danger

**Certified Reference Material**



| Component                  | CAS #     | Purity % | Prepared Concentration <sup>2</sup> | Certified Analyte Concentration <sup>1</sup> |
|----------------------------|-----------|----------|-------------------------------------|----------------------------------------------|
|                            |           | (LCMS)   | (µg/mL)                             | (µg/mL)                                      |
| 2,4-Diamino-6-nitrotoluene | 6629-29-4 | 99.0     | 100.3                               | 99.3                                         |

This Certified Reference Material was verified in accordance with ISO/IEC 17025

A product with a suffix (-1A, -2B, etc. or -01, -02, etc.) on its lot number has had its expiration date extended and is identical to the same lot number without the suffix.

<sup>2</sup> All weights are traceable through NIST, Test No. 684/289871-17

<sup>1</sup> Certified Analyte Concentration = Purity x Prepared Concentration.

The Uncertainty associated with the certified concentration reported on this certificate is  $\pm 2.4\%$ . This value is the combined expanded uncertainty and represents an estimated standard deviation equal to the positive square root of the total variation of the uncertainty of components. A normal distribution is assumed and a coverage factor of K=2 is chosen using approximately a 95% confidence level.

Labels and certificates follow U.S. Conventions in reporting numerical values: A comma (,) is used to separate units of one-thousand or greater. A period (.) is used as a decimal place marker.

The information on this certificate may not be reproduced without the express permission of the manufacturer. See reverse side for additional information

Hazard Information: Please refer to the SDS for information regarding the hazards associated with using this material.

This product was prepared according to in-house procedures and is guaranteed to be homogeneous.

Certified By:



Larry Decker, Organic QC Manager

Reagent

---

**2,4Diamino6NT\_00095**

# CERTIFICATE OF ANALYSIS

**Catalog No:** M-8330-ADD-12

**Description:** 2,4-Diamino-6-nitrotoluene

**Lot:** 220121267-01

**Solvent:** Acetonitrile

**Hazards:** Refer to SDS for complete safety information

**Date Certified:** Mar 10, 2022

**Expiration:** Apr 10, 2024

**Sample Size:** 1 mL

**Components:** 1

**Storage Condition:** Refrig (0-5 °C)



Signal Word: Danger

**Certified Reference Material**



| Component                  | CAS #     | Purity % | Prepared Concentration <sup>2</sup> | Certified Analyte Concentration <sup>1</sup> |
|----------------------------|-----------|----------|-------------------------------------|----------------------------------------------|
|                            |           | (LCMS)   | (µg/mL)                             | (µg/mL)                                      |
| 2,4-Diamino-6-nitrotoluene | 6629-29-4 | 99.0     | 100.3                               | 99.3                                         |

This Certified Reference Material was verified in accordance with ISO/IEC 17025

A product with a suffix (-1A, -2B, etc. or -01, -02, etc.) on its lot number has had its expiration date extended and is identical to the same lot number without the suffix.

<sup>2</sup> All weights are traceable through NIST, Test No. 684/289871-17

<sup>1</sup> Certified Analyte Concentration = Purity x Prepared Concentration.

The Uncertainty associated with the certified concentration reported on this certificate is  $\pm 2.4\%$ . This value is the combined expanded uncertainty and represents an estimated standard deviation equal to the positive square root of the total variation of the uncertainty of components. A normal distribution is assumed and a coverage factor of K=2 is chosen using approximately a 95% confidence level.

Labels and certificates follow U.S. Conventions in reporting numerical values: A comma (,) is used to separate units of one-thousand or greater. A period (.) is used as a decimal place marker.

The information on this certificate may not be reproduced without the express permission of the manufacturer. See reverse side for additional information

Hazard Information: Please refer to the SDS for information regarding the hazards associated with using this material.

This product was prepared according to in-house procedures and is guaranteed to be homogeneous.

Certified By:



Larry Decker, Organic QC Manager

Reagent

---

**2,4Diamino6NT\_00098**

# CERTIFICATE OF ANALYSIS

**Catalog No:** M-8330-ADD-12

**Description:** 2,4-Diamino-6-nitrotoluene

**Lot:** 220121267-01

**Solvent:** Acetonitrile

**Hazards:** Refer to SDS for complete safety information

**Date Certified:** Mar 10, 2022

**Expiration:** Apr 10, 2024

**Sample Size:** 1 mL

**Components:** 1

**Storage Condition:** Refrig (0-5 °C)



**Signal Word:** Danger

**Certified Reference Material**



AR-1463

| Component                  | CAS #     | Purity % | Prepared Concentration <sup>2</sup> | Certified Analyte Concentration <sup>1</sup> |
|----------------------------|-----------|----------|-------------------------------------|----------------------------------------------|
|                            |           | (LCMS)   | (µg/mL)                             | (µg/mL)                                      |
| 2,4-Diamino-6-nitrotoluene | 6629-29-4 | 99.0     | 100.3                               | 99.3                                         |

This Certified Reference Material was verified in accordance with ISO/IEC 17025

A product with a suffix (-1A, -2B, etc. or -01, -02, etc.) on its lot number has had its expiration date extended and is identical to the same lot number without the suffix.

<sup>2</sup> All weights are traceable through NIST, Test No. 684/289871-17

<sup>1</sup> Certified Analyte Concentration = Purity x Prepared Concentration.

The Uncertainty associated with the certified concentration reported on this certificate is  $\pm 2.4\%$ . This value is the combined expanded uncertainty and represents an estimated standard deviation equal to the positive square root of the total variation of the uncertainty of components. A normal distribution is assumed and a coverage factor of K=2 is chosen using approximately a 95% confidence level.

Labels and certificates follow U.S. Conventions in reporting numerical values: A comma (,) is used to separate units of one-thousand or greater. A period (.) is used as a decimal place marker.

The information on this certificate may not be reproduced without the express permission of the manufacturer. See reverse side for additional information

Hazard Information: Please refer to the SDS for information regarding the hazards associated with using this material.

This product was prepared according to in-house procedures and is guaranteed to be homogeneous.

Certified By:

Larry Decker, Organic QC Manager



Reagent

---

**2,4Diamino6NT\_00099**

# CERTIFICATE OF ANALYSIS

**Catalog No:** M-8330-ADD-12

**Description:** 2,4-Diamino-6-nitrotoluene

**Lot:** 220121267-01

**Solvent:** Acetonitrile

**Hazards:** Refer to SDS for complete safety information

**Date Certified:** Mar 10, 2022

**Expiration:** Apr 10, 2024

**Sample Size:** 1 mL

**Components:** 1

**Storage Condition:** Refrig (0-5 °C)



**Signal Word:** Danger

**Certified Reference Material**



AR-1463

| Component                  | CAS #     | Purity % | Prepared Concentration <sup>2</sup> | Certified Analyte Concentration <sup>1</sup> |
|----------------------------|-----------|----------|-------------------------------------|----------------------------------------------|
|                            |           | (LCMS)   | (µg/mL)                             | (µg/mL)                                      |
| 2,4-Diamino-6-nitrotoluene | 6629-29-4 | 99.0     | 100.3                               | 99.3                                         |

This Certified Reference Material was verified in accordance with ISO/IEC 17025

A product with a suffix (-1A, -2B, etc. or -01, -02, etc.) on its lot number has had its expiration date extended and is identical to the same lot number without the suffix.

<sup>2</sup> All weights are traceable through NIST, Test No. 684/289871-17

<sup>1</sup> Certified Analyte Concentration = Purity x Prepared Concentration.

The Uncertainty associated with the certified concentration reported on this certificate is  $\pm 2.4\%$ . This value is the combined expanded uncertainty and represents an estimated standard deviation equal to the positive square root of the total variation of the uncertainty of components. A normal distribution is assumed and a coverage factor of K=2 is chosen using approximately a 95% confidence level.

Labels and certificates follow U.S. Conventions in reporting numerical values: A comma (,) is used to separate units of one-thousand or greater. A period (.) is used as a decimal place marker.

The information on this certificate may not be reproduced without the express permission of the manufacturer. See reverse side for additional information

Hazard Information: Please refer to the SDS for information regarding the hazards associated with using this material.

This product was prepared according to in-house procedures and is guaranteed to be homogeneous.

Certified By:

Larry Decker, Organic QC Manager

Reagent

---

**2,4Diamino6NT\_00100**

# CERTIFICATE OF ANALYSIS

**Catalog No:** M-8330-ADD-12

**Description:** 2,4-Diamino-6-nitrotoluene

**Lot:** 220121267-01

**Solvent:** Acetonitrile

**Hazards:** Refer to SDS for complete safety information

**Date Certified:** Mar 10, 2022

**Expiration:** Apr 10, 2024

**Sample Size:** 1 mL

**Components:** 1

**Storage Condition:** Refrig (0-5 °C)



**Signal Word:** Danger

**Certified Reference Material**



| Component                  | CAS #     | Purity % | Prepared Concentration <sup>2</sup> | Certified Analyte Concentration <sup>1</sup> |
|----------------------------|-----------|----------|-------------------------------------|----------------------------------------------|
|                            |           | (LCMS)   | (µg/mL)                             | (µg/mL)                                      |
| 2,4-Diamino-6-nitrotoluene | 6629-29-4 | 99.0     | 100.3                               | 99.3                                         |

This Certified Reference Material was verified in accordance with ISO/IEC 17025

A product with a suffix (-1A, -2B, etc. or -01, -02, etc.) on its lot number has had its expiration date extended and is identical to the same lot number without the suffix.

<sup>2</sup> All weights are traceable through NIST, Test No. 684/289871-17

<sup>1</sup> Certified Analyte Concentration = Purity x Prepared Concentration.

The Uncertainty associated with the certified concentration reported on this certificate is  $\pm 2.4\%$ . This value is the combined expanded uncertainty and represents an estimated standard deviation equal to the positive square root of the total variation of the uncertainty of components. A normal distribution is assumed and a coverage factor of  $K=2$  is chosen using approximately a 95% confidence level.

Labels and certificates follow U.S. Conventions in reporting numerical values: A comma (,) is used to separate units of one-thousand or greater. A period (.) is used as a decimal place marker.

The information on this certificate may not be reproduced without the express permission of the manufacturer. See reverse side for additional information

Hazard Information: Please refer to the SDS for information regarding the hazards associated with using this material.

This product was prepared according to in-house procedures and is guaranteed to be homogeneous.

Certified By:

Larry Decker, Organic QC Manager

Reagent

---

**2,6Diamino4NT\_00095**

# CERTIFICATE OF ANALYSIS

**Catalog No:** M-8330-ADD-13

**Description:** 2,6-Diamino-4-nitrotoluene

**Lot:** 222051002

**Solvent:** Acetonitrile

**Hazards:** Refer to SDS for complete safety information

**Date Certified:** May 3, 2022

**Expiration:** Jun 3, 2024

**Sample Size:** 1 mL

**Components:** 1

**Storage Condition:** Refrig (0-5 °C)



**Signal Word:** Danger

**Certified Reference Material**



| Component                  | CAS #      | Purity % | Prepared Concentration <sup>2</sup> | Certified Analyte Concentration <sup>1</sup> |
|----------------------------|------------|----------|-------------------------------------|----------------------------------------------|
|                            |            | (GC/FID) | (µg/mL)                             | (µg/mL)                                      |
| 2,6-Diamino-4-nitrotoluene | 59229-75-3 | 99.7     | 101.0                               | 100.7                                        |

This Certified Reference Material was verified in accordance with ISO/IEC 17025

A product with a suffix (-1A, -2B, etc. or -01, -02, etc.) on its lot number has had its expiration date extended and is identical to the same lot number without the suffix.

<sup>2</sup> All weights are traceable through NIST, Test No. 684/289871-17

<sup>1</sup> Certified Analyte Concentration = Purity x Prepared Concentration.

The Uncertainty associated with the certified concentration reported on this certificate is  $\pm 2.4\%$ . This value is the combined expanded uncertainty and represents an estimated standard deviation equal to the positive square root of the total variation of the uncertainty of components. A normal distribution is assumed and a coverage factor of K=2 is chosen using approximately a 95% confidence level.

Labels and certificates follow U.S. Conventions in reporting numerical values: A comma (,) is used to separate units of one-thousand or greater. A period (.) is used as a decimal place marker.

The information on this certificate may not be reproduced without the express permission of the manufacturer. See reverse side for additional information

Hazard Information: Please refer to the SDS for information regarding the hazards associated with using this material.

This product was prepared according to in-house procedures and is guaranteed to be homogeneous.

Certified By:

Larry Decker, Organic QC Manager

Reagent

---

**2,6Diamino4NT\_00096**

# CERTIFICATE OF ANALYSIS

**Catalog No:** M-8330-ADD-13

**Description:** 2,6-Diamino-4-nitrotoluene

**Lot:** 222051002

**Solvent:** Acetonitrile

**Hazards:** Refer to SDS for complete safety information

**Date Certified:** May 3, 2022

**Expiration:** Jun 3, 2024

**Sample Size:** 1 mL

**Components:** 1

**Storage Condition:** Refrig (0-5 °C)



**Signal Word:** Danger

**Certified Reference Material**



| Component                  | CAS #      | Purity % | Prepared Concentration <sup>2</sup> | Certified Analyte Concentration <sup>1</sup> |
|----------------------------|------------|----------|-------------------------------------|----------------------------------------------|
|                            |            | (GC/FID) | (µg/mL)                             | (µg/mL)                                      |
| 2,6-Diamino-4-nitrotoluene | 59229-75-3 | 99.7     | 101.0                               | 100.7                                        |

This Certified Reference Material was verified in accordance with ISO/IEC 17025

A product with a suffix (-1A, -2B, etc. or -01, -02, etc.) on its lot number has had its expiration date extended and is identical to the same lot number without the suffix.

<sup>2</sup> All weights are traceable through NIST, Test No. 684/289871-17

<sup>1</sup> Certified Analyte Concentration = Purity x Prepared Concentration.

The Uncertainty associated with the certified concentration reported on this certificate is  $\pm 2.4\%$ . This value is the combined expanded uncertainty and represents an estimated standard deviation equal to the positive square root of the total variation of the uncertainty of components. A normal distribution is assumed and a coverage factor of K=2 is chosen using approximately a 95% confidence level.

Labels and certificates follow U.S. Conventions in reporting numerical values: A comma (,) is used to separate units of one-thousand or greater. A period (.) is used as a decimal place marker.

The information on this certificate may not be reproduced without the express permission of the manufacturer. See reverse side for additional information

Hazard Information: Please refer to the SDS for information regarding the hazards associated with using this material.

This product was prepared according to in-house procedures and is guaranteed to be homogeneous.

Certified By:

Larry Decker, Organic QC Manager



Reagent

---

**2,6Diamino4NT\_00099**

# CERTIFICATE OF ANALYSIS

**Catalog No:** M-8330-ADD-13

**Description:** 2,6-Diamino-4-nitrotoluene

**Lot:** 222051002

**Solvent:** Acetonitrile

**Hazards:** Refer to SDS for complete safety information

**Date Certified:** May 3, 2022

**Expiration:** Jun 3, 2024

**Sample Size:** 1 mL

**Components:** 1

**Storage Condition:** Refrig (0-5 °C)



Signal Word: Danger

Certified Reference Material



AR-1463

| Component                  | CAS #      | Purity % | Prepared Concentration <sup>2</sup> | Certified Analyte Concentration <sup>1</sup> |
|----------------------------|------------|----------|-------------------------------------|----------------------------------------------|
|                            |            | (GC/FID) | (µg/mL)                             | (µg/mL)                                      |
| 2,6-Diamino-4-nitrotoluene | 59229-75-3 | 99.7     | 101.0                               | 100.7                                        |

This Certified Reference Material was verified in accordance with ISO/IEC 17025

A product with a suffix (-1A, -2B, etc. or -01, -02, etc.) on its lot number has had its expiration date extended and is identical to the same lot number without the suffix.

<sup>2</sup> All weights are traceable through NIST, Test No. 684/289871-17

<sup>1</sup> Certified Analyte Concentration = Purity x Prepared Concentration.

The Uncertainty associated with the certified concentration reported on this certificate is  $\pm 2.4\%$ . This value is the combined expanded uncertainty and represents an estimated standard deviation equal to the positive square root of the total variation of the uncertainty of components. A normal distribution is assumed and a coverage factor of K=2 is chosen using approximately a 95% confidence level.

Labels and certificates follow U.S. Conventions in reporting numerical values: A comma (,) is used to separate units of one-thousand or greater. A period (.) is used as a decimal place marker.

The information on this certificate may not be reproduced without the express permission of the manufacturer. See reverse side for additional information

Hazard Information: Please refer to the SDS for information regarding the hazards associated with using this material.

This product was prepared according to in-house procedures and is guaranteed to be homogeneous.

Certified By:

Larry Decker, Organic QC Manager

Reagent

---

**2,6Diamino4NT\_00100**

# CERTIFICATE OF ANALYSIS

**Catalog No:** M-8330-ADD-13

**Description:** 2,6-Diamino-4-nitrotoluene

**Lot:** 222051002

**Solvent:** Acetonitrile

**Hazards:** Refer to SDS for complete safety information

**Date Certified:** May 3, 2022

**Expiration:** Jun 3, 2024

**Sample Size:** 1 mL

**Components:** 1

**Storage Condition:** Refrig (0-5 °C)



Signal Word: Danger

Certified Reference Material



AR-1463

| Component                  | CAS #      | Purity %<br>(GC/FID) | Prepared<br>Concentration <sup>2</sup><br>(µg/mL) | Certified Analyte<br>Concentration <sup>1</sup><br>(µg/mL) |
|----------------------------|------------|----------------------|---------------------------------------------------|------------------------------------------------------------|
| 2,6-Diamino-4-nitrotoluene | 59229-75-3 | 99.7                 | 101.0                                             | 100.7                                                      |

This Certified Reference Material was verified in accordance with ISO/IEC 17025

A product with a suffix (-1A, -2B, etc. or -01, -02, etc.) on its lot number has had its expiration date extended and is identical to the same lot number without the suffix.

<sup>2</sup> All weights are traceable through NIST, Test No. 684/289871-17

<sup>1</sup> Certified Analyte Concentration = Purity x Prepared Concentration.

The Uncertainty associated with the certified concentration reported on this certificate is  $\pm 2.4\%$ . This value is the combined expanded uncertainty and represents an estimated standard deviation equal to the positive square root of the total variation of the uncertainty of components. A normal distribution is assumed and a coverage factor of K=2 is chosen using approximately a 95% confidence level.

Labels and certificates follow U.S. Conventions in reporting numerical values: A comma (,) is used to separate units of one-thousand or greater. A period (.) is used as a decimal place marker.

The information on this certificate may not be reproduced without the express permission of the manufacturer. See reverse side for additional information

Hazard Information: Please refer to the SDS for information regarding the hazards associated with using this material.

This product was prepared according to in-house procedures and is guaranteed to be homogeneous.

Certified By:

Larry Decker, Organic QC Manager

Reagent

---

**2,6Diamino4NT\_00102**

# CERTIFICATE OF ANALYSIS

**Catalog No:** M-8330-ADD-13

**Description:** 2,6-Diamino-4-nitrotoluene

**Lot:** 222051002

**Solvent:** Acetonitrile

**Hazards:** Refer to SDS for complete safety information

**Date Certified:** May 3, 2022

**Expiration:** Jun 3, 2024

**Sample Size:** 1 mL

**Components:** 1

**Storage Condition:** Refrig (0-5 °C)



Signal Word: Danger

Certified Reference Material



AR-1463

| Component                  | CAS #      | Purity % | Prepared Concentration <sup>2</sup> | Certified Analyte Concentration <sup>1</sup> |
|----------------------------|------------|----------|-------------------------------------|----------------------------------------------|
|                            |            | (GC/FID) | (µg/mL)                             | (µg/mL)                                      |
| 2,6-Diamino-4-nitrotoluene | 59229-75-3 | 99.7     | 101.0                               | 100.7                                        |

This Certified Reference Material was verified in accordance with ISO/IEC 17025

A product with a suffix (-1A, -2B, etc. or -01, -02, etc.) on its lot number has had its expiration date extended and is identical to the same lot number without the suffix.

<sup>2</sup> All weights are traceable through NIST, Test No. 684/289871-17

<sup>1</sup> Certified Analyte Concentration = Purity x Prepared Concentration.

The Uncertainty associated with the certified concentration reported on this certificate is  $\pm 2.4\%$ . This value is the combined expanded uncertainty and represents an estimated standard deviation equal to the positive square root of the total variation of the uncertainty of components. A normal distribution is assumed and a coverage factor of K=2 is chosen using approximately a 95% confidence level.

Labels and certificates follow U.S. Conventions in reporting numerical values: A comma (,) is used to separate units of one-thousand or greater. A period (.) is used as a decimal place marker.

The information on this certificate may not be reproduced without the express permission of the manufacturer. See reverse side for additional information

Hazard Information: Please refer to the SDS for information regarding the hazards associated with using this material.

This product was prepared according to in-house procedures and is guaranteed to be homogeneous.

Certified By:

Larry Decker, Organic QC Manager

Reagent

---

**2,6Diamino4NT\_00105**

# CERTIFICATE OF ANALYSIS

**Catalog No:** M-8330-ADD-13

**Description:** 2,6-Diamino-4-nitrotoluene

**Lot:** 222051002

**Solvent:** Acetonitrile

**Hazards:** Refer to SDS for complete safety information

**Date Certified:** May 3, 2022

**Expiration:** Jun 3, 2024

**Sample Size:** 1 mL

**Components:** 1

**Storage Condition:** Refrig (0-5 °C)



Signal Word: Danger

Certified Reference Material



AR-1463

| Component                  | CAS #      | Purity % | Prepared Concentration <sup>2</sup> | Certified Analyte Concentration <sup>1</sup> |
|----------------------------|------------|----------|-------------------------------------|----------------------------------------------|
|                            |            | (GC/FID) | (µg/mL)                             | (µg/mL)                                      |
| 2,6-Diamino-4-nitrotoluene | 59229-75-3 | 99.7     | 101.0                               | 100.7                                        |

This Certified Reference Material was verified in accordance with ISO/IEC 17025

A product with a suffix (-1A, -2B, etc. or -01, -02, etc.) on its lot number has had its expiration date extended and is identical to the same lot number without the suffix.

<sup>2</sup> All weights are traceable through NIST, Test No. 684/289871-17

<sup>1</sup> Certified Analyte Concentration = Purity x Prepared Concentration.

The Uncertainty associated with the certified concentration reported on this certificate is  $\pm 2.4\%$ . This value is the combined expanded uncertainty and represents an estimated standard deviation equal to the positive square root of the total variation of the uncertainty of components. A normal distribution is assumed and a coverage factor of K=2 is chosen using approximately a 95% confidence level.

Labels and certificates follow U.S. Conventions in reporting numerical values: A comma (,) is used to separate units of one-thousand or greater. A period (.) is used as a decimal place marker.

The information on this certificate may not be reproduced without the express permission of the manufacturer. See reverse side for additional information

Hazard Information: Please refer to the SDS for information regarding the hazards associated with using this material.

This product was prepared according to in-house procedures and is guaranteed to be homogeneous.

Certified By:

Larry Decker, Organic QC Manager



Reagent

---

**3,5-DNA LCS\_00042**

**Preliminary Report**

Eurofins Denver

LCS, Lab Control Sample Report

Sample Path: \\chromfs\Denver\ChromData\CHHPLC\_X\20221108-115900.b\11080002.D

Lims ID: 8330AcidH2O\_00003

Inj. Date: 08-Nov-2022 14:15:10

Worklist ID: 280-0115900-002

Instrument: CHHPLC\_X3

Method: 8330\_X3

| Compound              | Amount Added | Amount Recovered | %Rec | Limits 1<br>3535 | Limits 2<br>3535 | Limits 3<br>3535 |
|-----------------------|--------------|------------------|------|------------------|------------------|------------------|
| 14 3,5-Dinitroaniline | 0.5000       | 0.4857           | 97.1 |                  | 71-117           | 55-119           |

Samples for Limit Group: 1, Lims Prep Method: 3535

570-115659-R-2-A

570-115659-R-3-A

570-115659-R-4-A

570-115659-R-5-A

570-115659-R-6-A

Samples for Limit Group: 2, Lims Prep Method: 3535

280-168716-B-1-A

280-168716-B-2-A

280-168716-B-3-A

280-168716-B-4-A

280-168716-B-5-A

280-168716-B-6-A

280-168716-B-7-A

280-168716-B-8-A

280-168716-B-9-A

680-224687-A-1-A

680-224687-A-2-A

680-224687-D-3-A

680-224687-B-4-A

Samples for Limit Group: 3, Lims Prep Method: 3535

860-36027-C-1-A

Reagent

---

**3,5-DNA Stock\_00036**



110 Benner Circle  
Bellefonte, PA 16823-8812  
Tel: (800)356-1688  
Fax: (814)353-1309

www.restek.com

# CERTIFIED REFERENCE MATERIAL

## Certificate of Analysis



### FOR LABORATORY USE ONLY-READ SDS PRIOR TO USE.

*This Reference Material is intended for Laboratory Use Only as a standard for the qualitative and/or quantitative determination of the analyte(s) listed.*

Catalog No. : 31661 Lot No.: A0149736  
Description : 3,5-Dinitroaniline Standard  
3, 5-Dinitroaniline Std 1000µg/mL, Acetonitrile, 1mL/ampul  
Container Size : 2 mL Pkg Amt: > 1 mL  
Expiration Date : January 31, 2024 Storage: 10°C or colder

### CERTIFIED VALUES

| Elution Order | Compound                                                            | Grav. Conc. (weight/volume) | Expanded Uncertainty (95% C.L.; K=2) |       |             |  |
|---------------|---------------------------------------------------------------------|-----------------------------|--------------------------------------|-------|-------------|--|
| 1             | 3,5-Dinitroaniline<br>CAS # 618-87-1<br>Purity 99%<br>(Lot 10311HS) | 1,000.0 µg/mL               | +/- 10.0737                          | µg/mL | Gravimetric |  |
|               |                                                                     |                             | +/- 31.3469                          | µg/mL | Unstressed  |  |
|               |                                                                     |                             | +/- 31.3469                          | µg/mL | Stressed    |  |

Solvent: Acetonitrile  
CAS # 75-05-8  
Purity 99%

**Column:**  
250mm x 4.6mm  
Ultra C18 (cat.# 9174575)


**Flow Rate:**  
1.0 ml/min.

**Mobile Phase A:**  
water:methanol (44:56 V/V)

**Mobile Phase B:**

**Mobile Phase Composition:**  
100%A

**Det. Type:**  
Wavelength: 210 nm



This chromatogram represents a general set of testing conditions chosen for product acceptance. For optimal results in your lab, conditions should be adjusted for your specific instrument, method, and application.

*Jessica McClenahan*

Jessica McClenahan - Operations Technician I

**Date Mixed:** 03-Jun-2019

**Balance:** 1128360905

*Jennifer L Pollino*

Jennifer Pollino - Operations Tech-ARM QC

**Date Passed:** 05-Jun-2019

Manufactured under Restek's ISO 9001:2015  
Registered Quality System  
Certificate #FM 80397

Reagent

---

**3,5-DNA Stock\_00038**



110 Benner Circle  
Bellefonte, PA 16823-8812  
Tel: (800)356-1688  
Fax: (814)353-1309

www.restek.com

# CERTIFIED REFERENCE MATERIAL

## Certificate of Analysis



### FOR LABORATORY USE ONLY-READ SDS PRIOR TO USE.

*This Reference Material is intended for Laboratory Use Only as a standard for the qualitative and/or quantitative determination of the analyte(s) listed.*

Catalog No. : 31661 Lot No.: A0164976  
Description : 3,5-Dinitroaniline Standard  
3, 5-Dinitroaniline Std 1000µg/mL, Acetonitrile, 1mL/ampul  
Container Size : 2 mL Pkg Amt: > 1 mL  
Expiration Date : May 31, 2025 Storage: 10°C or colder  
Ship: Ambient

### CERTIFIED VALUES

| Elution Order | Compound           | Grav. Conc.<br>(weight/volume) | Expanded Uncertainty<br>(95% C.L.; K=2) |         |                   |
|---------------|--------------------|--------------------------------|-----------------------------------------|---------|-------------------|
| 1             | 3,5-Dinitroaniline | 1,000.8 µg/mL                  | +/-                                     | 10.0818 | µg/mL Gravimetric |
|               | CAS # 618-87-1     |                                | +/-                                     | 31.3720 | µg/mL Unstressed  |
|               | Purity 99%         |                                | +/-                                     | 31.3720 | µg/mL Stressed    |
|               | (Lot 10311HS)      |                                |                                         |         |                   |

Solvent: Acetonitrile  
CAS # 75-05-8  
Purity 99%

**Column:**  
250mm x 4.6mm  
Ultra C18 (cat.# 9174575)


**Flow Rate:**  
1.0 ml/min.

**Mobile Phase A:**  
water:methanol (44:56 V/V)

**Mobile Phase B:**

**Mobile Phase Composition:**  
100%A

**Det. Type:**  
Wavelength: 210nm & 254nm




This chromatogram represents a general set of testing conditions chosen for product acceptance. For optimal results in your lab, conditions should be adjusted for your specific instrument, method, and application.

  
Katelyn McGinn - Operations Tech I

**Date Mixed:** 02-Oct-2020

**Balance:** 1128353505

  
Jennifer Pollino - Operations Tech-ARM QC

**Date Passed:** 06-Oct-2020

Manufactured under Restek's ISO 9001:2015  
Registered Quality System  
Certificate #FM 80397



## General Certified Reference Material Notes

### Expiration Notes:

- Expiration date valid for unopened ampul stored in compliance with the recommended conditions.
- Uncertainty, concentration, and expiration of the CRM are based on the unopened product being stored according to the recommended condition found in the storage field.

### Purity Notes:

- Purity and/or chemical identity are determined by one or more of the following techniques: GC/FID, HPLC, GC/μECD, GC/MS, LC/MS, RI, and/or melting point.
- Compounds with a listed purity of less than 99% have been weight corrected to compensate for impurities and/or salts. A correction factor is used to calculate the amount of compound necessary to achieve the desired concentration of the parent compound in solution.
- Purity of isomeric compounds is reported as the sum of the isomers.
- Purity values are rounded to the nearest whole number.

### Certified Uncertainty Value Notes:

- The uncertainties are determined in accordance with ISO 17034 and Guide 35. The certified combined stressed uncertainty value ( includes gravimetric uncertainty, homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty and were combined using the following formula:

$$U_{combined\ stressed} = k \sqrt{U_{gravimetric}^2 + U_{homogeneity}^2 + U_{storage\ stability}^2 + U_{shipping\ stability}^2}$$

$k$  is a coverage factor of 2, which gives a level of confidence of approximately 95%.

- It is important to note that the shipping stability uncertainty was obtained under temperature extremes for specific time intervals; therefore, the certified combined stressed uncertainty value should only be applied to the product if it was stored at non-standard temperature conditions up to and including 7 days. Contact Restek Technical Service at [www.restek.com/Contact-Us](http://www.restek.com/Contact-Us) for use recommendations if your shipment was in-transit for more than 7 days at non-standard temperature conditions.
- Apply the certified combined unstressed uncertainty value if the product was received under standard shipping conditions. Apply the certified combined stressed uncertainty value if the product was received under non-standard conditions as specified below.

| Label Conditions                                          | Standard Conditions | Non-Standard Conditions |
|-----------------------------------------------------------|---------------------|-------------------------|
| 25°C Nominal (Room Temperature)                           | < 60°C              | ≥ 60°C up to 7 days     |
| 10°C or colder (Refrigerate)                              | < 40°C              | ≥ 40°C up to 7 days     |
| 0°C or colder (Freezer)<br>-20°C or colder (Deep Freezer) | < 25°C              | ≥ 25°C up to 7 days     |

- Separate (not combined) uncertainty values for gravimetric uncertainty are also displayed on the certificate, if needed, separate homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty values are available by contacting Restek Technical Service at [www.restek.com/Contact-Us](http://www.restek.com/Contact-Us).
- The packaged amount is the minimum sample size for which uncertainty is valid. The ampules are over-filled to ensure that the minimum packaged amount can be sufficiently transferred.

### Manufacturing Notes:

- Concentration is based upon gravimetric preparation using either a balance whose calibration has been verified daily using NIST traceable weights, and/or dilutions with Class A glassware.

### Handling Notes:

- Stability of the unopened product, when stored in compliance with the recommended conditions, is guaranteed through the expiration displayed on the product label and certificate. Contact Restek for additional opened product stability information, with the knowledge/understanding that open product stability is subject to the specific handling and environmental conditions to which the product is exposed. For your convenience Restek supplies deactivated vials with most standards packed in 2mL ampules. Larger volume deactivated vials are available through Restek as a custom ordered item. Additionally, Restek sells DMDCS for the purpose of glassware deactivation as catalog number 31861, which includes complete instructions.



Reagent

---

**3,5-DNA Stock\_00040**



# CERTIFIED REFERENCE MATERIAL

110 Benner Circle  
Bellefonte, PA 16823-8812  
Tel: (800)356-1688  
Fax: (814)353-1309

www.restek.com

## Certificate of Analysis



### FOR LABORATORY USE ONLY-READ SDS PRIOR TO USE.

*This Reference Material is intended for Laboratory Use Only as a standard for the qualitative and/or quantitative determination of the analyte(s) listed.*

Catalog No. : 31661 Lot No.: A0172723  
Description : 3,5-Dinitroaniline Standard  
3, 5-Dinitroaniline Std 1000µg/mL, Acetonitrile, 1mL/ampul  
Container Size : 2 mL Pkg Amt: > 1 mL  
Expiration Date : December 31, 2025 Storage: 10°C or colder  
Ship: Ambient

### CERTIFIED VALUES

| Elution Order | Compound                                                            | Grav. Conc. (weight/volume) | Expanded Uncertainty (95% C.L.; K=2)                                                        |
|---------------|---------------------------------------------------------------------|-----------------------------|---------------------------------------------------------------------------------------------|
| 1             | 3,5-Dinitroaniline<br>CAS # 618-87-1<br>Purity 99%<br>(Lot 10311HS) | 1,000.0 µg/mL               | +/- 10.0737 µg/mL Gravimetric<br>+/- 31.3469 µg/mL Unstressed<br>+/- 31.3469 µg/mL Stressed |

Solvent: Acetonitrile  
CAS # 75-05-8  
Purity 99%

**Column:**  
250mm x 4.6mm  
Ultra C18 (cat.# 9174575)


**Flow Rate:**  
1.0 ml/min.

**Mobile Phase A:**  
water:methanol (44:56 V/V)


**Mobile Phase B:**

**Mobile Phase Composition:**  
100%A


**Det. Type:**  
Wavelength: 210nm & 254nm



This chromatogram represents a general set of testing conditions chosen for product acceptance. For optimal results in your lab, conditions should be adjusted for your specific instrument, method, and application.

  
Jeremy Johnson - Mfg. Supervisor

Date Mixed: 25-May-2021 Balance: B251644995

  
Marlene Cowan - Operations Tech I

Date Passed: 27-May-2021

Manufactured under Restek's ISO 9001:2015  
Registered Quality System  
Certificate #FM 80397

## General Certified Reference Material Notes

### Expiration Notes:

- Expiration date valid for unopened ampul stored in compliance with the recommended conditions.
- Uncertainty, concentration, and expiration of the CRM are based on the unopened product being stored according to the recommended condition found in the storage field.

### Purity Notes:

- Purity and/or chemical identity are determined by one or more of the following techniques: GC/FID, HPLC, GC/μECD, GC/MS, LC/MS, RI, and/or melting point.
- Compounds with a listed purity of less than 99% have been weight corrected to compensate for impurities and/or salts. A correction factor is used to calculate the amount of compound necessary to achieve the desired concentration of the parent compound in solution.
- Purity of isomeric compounds is reported as the sum of the isomers.
- Purity values are rounded to the nearest whole number.

### Certified Uncertainty Value Notes:

- The uncertainties are determined in accordance with ISO 17034 and Guide 35. The certified combined stressed uncertainty value ( includes gravimetric uncertainty, homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty and were combined using the following formula:

$$U_{combined\ stressed} = k \sqrt{U_{gravimetric}^2 + U_{homogeneity}^2 + U_{storage\ stability}^2 + U_{shipping\ stability}^2}$$

$k$  is a coverage factor of 2, which gives a level of confidence of approximately 95%.

- It is important to note that the shipping stability uncertainty was obtained under temperature extremes for specific time intervals; therefore, the certified combined stressed uncertainty value should only be applied to the product if it was stored at non-standard temperature conditions up to and including 7 days. Contact Restek Technical Service at [www.restek.com/Contact-Us](http://www.restek.com/Contact-Us) for use recommendations if your shipment was in-transit for more than 7 days at non-standard temperature conditions.
- Apply the certified combined unstressed uncertainty value if the product was received under standard shipping conditions. Apply the certified combined stressed uncertainty value if the product was received under non-standard conditions as specified below.

| Label Conditions                                          | Standard Conditions | Non-Standard Conditions |
|-----------------------------------------------------------|---------------------|-------------------------|
| 25°C Nominal (Room Temperature)                           | < 60°C              | ≥ 60°C up to 7 days     |
| 10°C or colder (Refrigerate)                              | < 40°C              | ≥ 40°C up to 7 days     |
| 0°C or colder (Freezer)<br>-20°C or colder (Deep Freezer) | < 25°C              | ≥ 25°C up to 7 days     |

- Separate (not combined) uncertainty values for gravimetric uncertainty are also displayed on the certificate, if needed, separate homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty values are available by contacting Restek Technical Service at [www.restek.com/Contact-Us](http://www.restek.com/Contact-Us).
- The packaged amount is the minimum sample size for which uncertainty is valid. The ampules are over-filled to ensure that the minimum packaged amount can be sufficiently transferred.

### Manufacturing Notes:

- Concentration is based upon gravimetric preparation using either a balance whose calibration has been verified daily using NIST traceable weights, and/or dilutions with Class A glassware.

### Handling Notes:

- Stability of the unopened product, when stored in compliance with the recommended conditions, is guaranteed through the expiration displayed on the product label and certificate. Contact Restek for additional opened product stability information, with the knowledge/understanding that open product stability is subject to the specific handling and environmental conditions to which the product is exposed. For your convenience Restek supplies deactivated vials with most standards packed in 2mL ampules. Larger volume deactivated vials are available through Restek as a custom ordered item. Additionally, Restek sells DMDCS for the purpose of glassware deactivation as catalog number 31861, which includes complete instructions.



Reagent

---

**3,5-DNA Stock\_00042**





CERTIFIED REFERENCE MATERIAL

110 Benner Circle  
Bellefonte, PA 16823-8812  
Tel: (800)356-1688  
Fax: (814)353-1309

www.restek.com

# Certificate of Analysis



## FOR LABORATORY USE ONLY-READ SDS PRIOR TO USE.

*This Reference Material is intended for Laboratory Use Only as a standard for the qualitative and/or quantitative determination of the analyte(s) listed.*

**Catalog No. :** 31661 **Lot No.:** A0185772  
**Description :** 3,5-Dinitroaniline Standard  
3, 5-Dinitroaniline Std 1000µg/mL, Acetonitrile, 1mL/ampul  
**Container Size :** 2 mL **Pkg Amt:** > 1 mL  
**Expiration Date :** December 31, 2026 **Storage:** 10°C or colder  
**Ship:** Ambient

## CERTIFIED VALUES

| Elution Order | Compound                                                            | Grav. Conc. (weight/volume) | Expanded Uncertainty (95% C.L.; K=2)                                                        |
|---------------|---------------------------------------------------------------------|-----------------------------|---------------------------------------------------------------------------------------------|
| 1             | 3,5-Dinitroaniline<br>CAS # 618-87-1<br>Purity 99%<br>(Lot 10311HS) | 1,008.0 µg/mL               | +/- 10.1543 µg/mL Gravimetric<br>+/- 31.5977 µg/mL Unstressed<br>+/- 31.5977 µg/mL Stressed |

**Solvent:** Acetonitrile  
CAS # 75-05-8  
Purity 99%

**Column:**  
250mm x 4.6mm  
Ultra C18 (cat.# 9174575)


**Flow Rate:**  
1.0 ml/min.

**Mobile Phase A:**  
water:methanol (44:56 V/V)

**Mobile Phase B:**

**Mobile Phase Composition:**  
100%A


**Det. Type:**  
Wavelength: 210nm & 254nm



This chromatogram represents a general set of testing conditions chosen for product acceptance. For optimal results in your lab, conditions should be adjusted for your specific instrument, method, and application.

  
Morgan Craighead - Mix Technician

**Date Mixed:** 27-May-2022      **Balance:** 1128353505

  
Marlina Cowan - Operations Tech I

**Date Passed:** 02-Jun-2022

Manufactured under Restek's ISO 9001:2015  
Registered Quality System  
Certificate #FM 80397

## General Certified Reference Material Notes

### Expiration Notes:

- Expiration date valid for unopened ampul stored in compliance with the recommended conditions.
- Uncertainty, concentration, and expiration of the CRM are based on the unopened product being stored according to the recommended condition found in the storage field.

### Purity Notes:

- Purity and/or chemical identity are determined by one or more of the following techniques: GC/FID, HPLC, GC/μECD, GC/MS, LC/MS, RI, and/or melting point.
- Compounds with a listed purity of less than 99% have been weight corrected to compensate for impurities and/or salts. A correction factor is used to calculate the amount of compound necessary to achieve the desired concentration of the parent compound in solution.
- Purity of isomeric compounds is reported as the sum of the isomers.
- Purity values are rounded to the nearest whole number.

### Certified Uncertainty Value Notes:

- The uncertainties are determined in accordance with ISO 17034 and Guide 35. The certified combined stressed uncertainty value ( includes gravimetric uncertainty, homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty and were combined using the following formula:

$$U_{combined\ stressed} = k \sqrt{U_{gravimetric}^2 + U_{homogeneity}^2 + U_{storage\ stability}^2 + U_{shipping\ stability}^2}$$

$k$  is a coverage factor of 2, which gives a level of confidence of approximately 95%.

- It is important to note that the shipping stability uncertainty was obtained under temperature extremes for specific time intervals; therefore, the certified combined stressed uncertainty value should only be applied to the product if it was stored at non-standard temperature conditions up to and including 7 days. Contact Restek Technical Service at [www.restek.com/Contact-Us](http://www.restek.com/Contact-Us) for use recommendations if your shipment was in-transit for more than 7 days at non-standard temperature conditions.
- Apply the certified combined unstressed uncertainty value if the product was received under standard shipping conditions. Apply the certified combined stressed uncertainty value if the product was received under non-standard conditions as specified below.

| Label Conditions                                          | Standard Conditions | Non-Standard Conditions |
|-----------------------------------------------------------|---------------------|-------------------------|
| 25°C Nominal (Room Temperature)                           | < 60°C              | ≥ 60°C up to 7 days     |
| 10°C or colder (Refrigerate)                              | < 40°C              | ≥ 40°C up to 7 days     |
| 0°C or colder (Freezer)<br>-20°C or colder (Deep Freezer) | < 25°C              | ≥ 25°C up to 7 days     |

- Separate (not combined) uncertainty values for gravimetric uncertainty are also displayed on the certificate, if needed, separate homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty values are available by contacting Restek Technical Service at [www.restek.com/Contact-Us](http://www.restek.com/Contact-Us).
- The packaged amount is the minimum sample size for which uncertainty is valid. The ampules are over-filled to ensure that the minimum packaged amount can be sufficiently transferred.

### Manufacturing Notes:

- Concentration is based upon gravimetric preparation using either a balance whose calibration has been verified daily using NIST traceable weights, and/or dilutions with Class A glassware.

### Handling Notes:

- Stability of the unopened product, when stored in compliance with the recommended conditions, is guaranteed through the expiration displayed on the product label and certificate. Contact Restek for additional opened product stability information, with the knowledge/understanding that open product stability is subject to the specific handling and environmental conditions to which the product is exposed. For your convenience Restek supplies deactivated vials with most standards packed in 2mL ampules. Larger volume deactivated vials are available through Restek as a custom ordered item. Additionally, Restek sells DMDCS for the purpose of glassware deactivation as catalog number 31861, which includes complete instructions.

Reagent

---

**50% NaOH\_00023**

# Certificate of Analysis

## Sodium Hydroxide, 50% (w/w) (760 g/L), Analytical Reagent Grade

**Lot Number:** 1102857

**Product Number:** 7291

**Manufacture Date:** FEB 05, 2021

**Expiration Date:** JAN 2024

The specifications below are equivalent to the ACS Reagent Grade specifications for Sodium Hydroxide pellets, corrected for assay.  
There is no ACS specification for 50% (w/w) Sodium Hydroxide solution.

| Name             | CAS#      | Grade   |
|------------------|-----------|---------|
| Water            | 7732-18-5 | Reagent |
| Sodium Hydroxide | 1310-73-2 | Reagent |

| Test                                                | Specification    | Result     | NIST SRM# |
|-----------------------------------------------------|------------------|------------|-----------|
| Appearance                                          | Colorless liquid | Passed     |           |
| Assay (vs. Sulfuric Acid/Phenolphthalein)           | 50-52 % (w/w)    | 50 % (w/w) | 723       |
| Chloride (Cl)                                       | max 0.002 %      | < 0.002 %  |           |
| Heavy Metals (as Ag)                                | max 0.001 %      | < 0.001 %  |           |
| Iron (Fe)                                           | max 5 ppm        | 0.2 ppm    | 3126      |
| NH <sub>4</sub> OH Precipitate                      | max 0.01 %       | 0.00009 %  |           |
| Nickel (Ni)                                         | max 5 ppm        | < 1 ppm    | 3136      |
| Nitrogen Compounds (as N)                           | max 5 ppm        | < 5 ppm    |           |
| Phosphate (PO <sub>4</sub> )                        | max 5 ppm        | < 5 ppm    |           |
| Potassium (K)                                       | max 0.01 %       | 0.002 %    | 3141      |
| Sodium Carbonate (Na <sub>2</sub> CO <sub>3</sub> ) | max 0.1 %        | 0.02 %     |           |
| Sulfate (SO <sub>4</sub> )                          | max 0.001 %      | < 0.001 %  |           |

| Specification                              | Reference           |
|--------------------------------------------|---------------------|
| Sodium Hydroxide Solution, 50%             | ASTM (D 2187 E)     |
| Sodium Hydroxide Solution, 50%             | ASTM (D 2187 F)     |
| Sodium hydroxide stock solution, about 50% | TAPPI (T 235 cm-85) |
| Sodium Hydroxide Solution                  | ASTM (D 2187 G)     |
| Sodium Hydroxide Solution, 50%             | ASTM (D 4548)       |
| Sodium Hydroxide Solution (50% w/w)        | ASTM (D 2036 A)     |
| Sodium Hydroxide Solution (50% w/w)        | ASTM (D 2036 B)     |
| Sodium Hydroxide Solution (50% w/w)        | ASTM (D 2036 C)     |
| Sodium Hydroxide Solution (50% w/w)        | ASTM (D 2036 D)     |

Volumetric glassware complies with Class A tolerance requirements of ASTM E 288 and NIST Circular 434; it is calibrated before first use and recalibrated regularly in accordance with ASTM E 542 and NIST Procedure NBSIR 74-461. Balances are calibrated regularly with weights certified traceable to the NIST national mass standard. Thermometers and temperature probes are calibrated before first use and recalibrated regularly with a thermometer traceable to NIST standards. All products are prepared according to master documents that assure manufacture according to validated methods. Batch records document raw material traceability and production and testing history for each lot manufactured.

| Part Number | Size / Package Type | Shelf Life (Unopened Container) |
|-------------|---------------------|---------------------------------|
| 7291-1      | 4 L natural poly    | 36 months                       |
| 7291-16     | 500 mL natural poly | 36 months                       |

**Recommended Storage:** 15°C - 30°C (59°F - 86°F)



Israel Alamudun (02/05/2021)

Quality Control Supervisor

This Certificate of Analysis is designed to comply with ISO Guide 31 "Reference Materials -- Contents of Certificates and Labels."

This test report shall not be reproduced, except in full, without the written approval of Ricca Chemical Company.

Reagent

---

**8330 LCS\_00121**

**Preliminary Report**

Eurofins Denver

LCS, Lab Control Sample Report

Sample Path: \\chromfs\Denver\ChromData\CHHPLC\_X\20230121-117963.b\01210012.D

Lims ID: 8330 LCS 121

Inj. Date: 21-Jan-2023 14:04:20

Worklist ID: 280-0117963-012

Instrument: CHHPLC\_X3

Method: 8330\_X3

| Compound                   | Amount Added | Amount Recovered | %Rec  | Limits 1<br>3535 | Limits 2<br>3535 | Limits 3<br>3535 |
|----------------------------|--------------|------------------|-------|------------------|------------------|------------------|
| 4 HMX                      | 0.5000       | 0.4694           | 93.9  | 65-135           | 66-115           |                  |
| 8 RDX                      | 0.5000       | 0.4904           | 98.1  | 68-130           | 69-122           |                  |
| 9 2,4,6-Trinitrophenol     | 0.5000       | 0.5399           | 108.0 | 80-120           | 63-135           |                  |
| 11 1,3,5-Trinitrobenzene   | 0.5000       | 0.5414           | 108.3 | 73-125           | 62-127           |                  |
| 12 1,3-Dinitrobenzene      | 0.5000       | 0.5314           | 106.3 | 78-120           | 59-131           |                  |
| 13 Nitrobenzene            | 0.5000       | 0.5214           | 104.3 | 65-134           | 46-144           |                  |
| 15 Tetryl                  | 0.5000       | 0.5279           | 105.6 | 64-128           | 56-131           |                  |
| 16 Nitroglycerin           | 5.00         | 5.40             | 107.9 | 74-127           | 70-125           |                  |
| 17 2,4,6-Trinitrotoluene   | 0.5000       | 0.4569           | 91.4  | 71-123           | 46-139           |                  |
| 18 4-Amino-2,6-dinitrotolu | 0.5000       | 0.5716           | 114.3 | 76-125           | 43-120           |                  |
| 19 2-Amino-4,6-dinitrotolu | 0.5000       | 0.4927           | 98.5  | 79-120           | 46-124           |                  |
| 20 2,6-Dinitrotoluene      | 0.5000       | 0.5185           | 103.7 | 77-127           | 51-130           |                  |
| 21 2,4-Dinitrotoluene      | 0.5000       | 0.5353           | 107.1 | 78-120           | 53-127           |                  |
| 22 o-Nitrotoluene          | 0.5000       | 0.5106           | 102.1 | 70-127           | 37-138           |                  |
| 23 p-Nitrotoluene          | 0.5000       | 0.4984           | 99.7  | 71-127           | 41-137           |                  |
| 24 m-Nitrotoluene          | 0.5000       | 0.4984           | 99.7  | 73-125           | 31-140           |                  |
| 25 PETN                    | 5.00         | 5.27             | 105.5 | 73-127           | 67-127           |                  |

Samples for Limit Group: 1, Lims Prep Method: 3535

280-171404-A-21-A

280-171267-A-1-A

280-171267-A-2-A

280-171267-A-3-A

280-171267-A-4-A

280-171267-A-5-A

280-171267-A-6-A

280-171267-A-7-A

Samples for Limit Group: 2, Lims Prep Method: 3535

550-196288-AO-1-A

570-123967-A-1-A

Samples for Limit Group: 3, Lims Prep Method: 3535

280-171237-A-1-A

280-171237-A-2-A

280-171237-A-3-A

280-171237-A-4-A



Reagent

---

**8330 LCS\_00126**

**Preliminary Report**

Eurofins Denver

LCS, Lab Control Sample Report

Sample Path: \\chromfs\Denver\ChromData\CHHPLC\_X\20230426-120796.b\04260006.D

Lims ID: C18column:B16162

Inj. Date: 26-Apr-2023 15:04:00

Worklist ID: 280-0120796-006

Instrument: CHHPLC\_X3

Method: 8330\_X3

| Compound                   | Amount Added | Amount Recovered | %Rec  | Limits 1<br>OB_Sonc_ | Limits 2<br>OB_Sonc_ |
|----------------------------|--------------|------------------|-------|----------------------|----------------------|
| 4 HMX                      | 0.5000       | 0.4359           | 87.2  | 80-120               | 74-124               |
| 8 RDX                      | 0.5000       | 0.4718           | 94.4  | 80-124               | 67-129               |
| 9 2,4,6-Trinitrophenol     | 0.5000       | 0.5241           | 104.8 | 38-154               |                      |
| 11 1,3,5-Trinitrobenzene   | 0.5000       | 0.5064           | 101.3 | 80-120               | 80-116               |
| 12 1,3-Dinitrobenzene      | 0.5000       | 0.5004           | 100.1 | 80-120               | 73-119               |
| 13 Nitrobenzene            | 0.5000       | 0.5010           | 100.2 | 76-122               | 67-129               |
| 14 3,5-Dinitroaniline      | 0.5000       | 0.4725           | 94.5  | 80-120               | 86-118               |
| 15 Tetryl                  | 0.5000       | 0.5711           | 114.2 | 80-120               | 68-135               |
| 16 Nitroglycerin           | 5.00         | 5.04             | 100.8 | 75-120               | 73-124               |
| 17 2,4,6-Trinitrotoluene   | 0.5000       | 0.4643           | 92.9  | 80-120               | 71-120               |
| 18 4-Amino-2,6-dinitrotolu | 0.5000       | 0.4973           | 99.5  | 80-120               | 64-127               |
| 19 2-Amino-4,6-dinitrotolu | 0.5000       | 0.4766           | 95.3  | 78-120               | 71-123               |
| 20 2,6-Dinitrotoluene      | 0.5000       | 0.4994           | 99.9  | 80-120               | 79-117               |
| 21 2,4-Dinitrotoluene      | 0.5000       | 0.4894           | 97.9  | 80-120               | 75-121               |
| 22 o-Nitrotoluene          | 0.5000       | 0.4896           | 97.9  | 80-124               | 70-124               |
| 23 p-Nitrotoluene          | 0.5000       | 0.4868           | 97.4  | 80-120               | 71-124               |
| 24 m-Nitrotoluene          | 0.5000       | 0.4890           | 97.8  | 80-122               | 67-129               |
| 25 PETN                    | 5.00         | 5.34             | 106.9 | 80-120               | 72-128               |

Samples for Limit Group: 1, Lims Prep Method: 8330B\_Sonc\_10g

|                   |                   |                   |
|-------------------|-------------------|-------------------|
| 280-174348-A-25-B | 280-174348-A-26-B | 280-174348-A-27-B |
| 280-174348-A-28-B | 280-174348-A-29-B | 280-174348-A-30-D |
| 280-174348-A-31-B | 280-174348-A-16-B | 280-174348-A-18-B |
| 280-174348-A-19-B | 280-174348-A-20-B | 280-174348-A-22-B |
| 280-174348-A-24-B | 280-174485-A-7-B  | 280-174485-A-8-B  |
| 280-174485-A-9-B  | 280-174482-A-1-B  | 280-174482-A-2-B  |
| 280-174482-A-3-B  | 280-174482-A-4-B  | 280-174482-A-5-B  |
| 280-174482-A-6-B  | 280-174482-A-7-B  | 280-174482-A-8-B  |
| 280-174482-A-9-D  | 280-174482-A-10-B | 280-174482-A-11-B |
| 280-174482-A-12-B | 280-174482-A-13-B | 280-174482-A-14-B |
| 280-174482-A-15-B | 280-174482-A-16-B | 280-174482-A-17-B |
| 280-174482-A-18-B | 280-174482-A-19-C |                   |

Samples for Limit Group: 2, Lims Prep Method: 8330B\_Sonc\_10g

|                  |                  |
|------------------|------------------|
| 280-174408-A-1-I | 280-174408-A-2-I |
|------------------|------------------|

Reagent

---

**8330 Stock\_TS\_00023**

## Reference Material Certificate

### Product Information Sheet

**Product Name:** Stock Standard

**Lot Number:** 0006684308

**Product Number:** NAIM-833E-1

**Lot Issue Date:** 01-Jun-2022

**Storage Conditions:** Store at Room Temperature (15° to 30°C).

**Expiration Date:** 30-Jun-2025

| Component Name              | CERTIFIED VALUES |                      |  | CAS#        | Analyte Lot |
|-----------------------------|------------------|----------------------|--|-------------|-------------|
|                             | Concentration    | Expanded Uncertainty |  |             |             |
| HMX                         | 1001             | ± 5 µg/mL            |  | 002691-41-0 | RM06237     |
| RDX                         | 1001             | ± 5 µg/mL            |  | 000121-82-4 | RM10915     |
| 1,3,5-trinitrobenzene       | 1001             | ± 5 µg/mL            |  | 000099-35-4 | RM17843     |
| m-dinitrobenzene            | 1002             | ± 5 µg/mL            |  | 000099-65-0 | RM14290     |
| nitrobenzene                | 1002             | ± 5 µg/mL            |  | 000098-95-3 | RM11472     |
| 2,4,6-trinitrotoluene (TNT) | 1001             | ± 5 µg/mL            |  | 000118-96-7 | RM16204     |
| 2,4-dinitrotoluene          | 1002             | ± 5 µg/mL            |  | 000121-14-2 | RM10279     |
| tetryl                      | 1003             | ± 5 µg/mL            |  | 000479-45-8 | RM14651     |
| 2,6-dinitrotoluene          | 1003             | ± 5 µg/mL            |  | 000606-20-2 | RM16636     |
| 2-nitrotoluene              | 1003             | ± 5 µg/mL            |  | 000088-72-2 | NT01996     |
| 3-nitrotoluene              | 1002             | ± 5 µg/mL            |  | 000099-08-1 | NT02212     |
| 4-nitrotoluene              | 1003             | ± 5 µg/mL            |  | 000099-99-0 | NT02096     |
| 2-amino-4,6-dinitrotoluene  | 1003             | ± 5 µg/mL            |  | 035572-78-2 | RM04232     |
| 4-amino-2,6-dinitrotoluene  | 1004             | ± 5 µg/mL            |  | 019406-51-0 | RM04226     |

**Matrix:** acetonitrile

**Description:**

This document is prepared in accordance with ISO 17034 and Guide 31. This analytical reference material standard was manufactured and verified in accordance with an ISO 9001 registered quality system and analyte concentrations were verified by an ISO 17025 accredited laboratory. The concentration and uncertainty value at the 95% confidence level for each analyte, determined gravimetrically, is listed above.

**Traceability:**

The balances used for these measurements are calibrated with weights traceable to NIST in compliance with ANSI/NCCL Z540.3, ISO 9001, ISO 17025, and ISO 17034. Calibrated Class A glassware is used for volumetric measurements. Thermometers are calibrated against a NIST traceable thermometer in accordance with NIST Special Publication 1088.

**Homogeneity:**

This analytical reference standard was unitized according to an in-house procedure and is guaranteed to be homogeneous. There is no minimum sub-sample size required.

**Instructions for Use:**

Sample aliquots for analysis should be withdrawn at 20°C to 25°C immediately after opening the container and should be processed without delay for the certified values to be valid within the stated uncertainties.

Reagent

---

**8330\_NG\_Stk\_00118**



110 Benner Circle  
Bellefonte, PA 16823-8812  
Tel: (800)356-1688  
Fax: (814)353-1309

www.restek.com

# CERTIFIED REFERENCE MATERIAL

## Certificate of Composition



### FOR LABORATORY USE ONLY-READ SDS PRIOR TO USE.

*This Reference Material is intended for Laboratory Use Only as a standard for the qualitative and/or quantitative determination of the analyte(s) listed.*

Catalog No. : 568871 Lot No.: A0188553  
Description : Custom Nitroglycerin Standard  
Custom Nitroglycerin Standard 5,000µg/mL, Acetonitrile, 1mL/ampul  
Container Size : 2 mL Pkg Amt: > 1 mL  
Expiration Date : August 31, 2025 Storage: 10°C or colder  
Ship: Ambient

### CERTIFIED VALUES

| Elution Order | Compound                                                        | Grav. Conc. (weight/volume) | Expanded Uncertainty (95% C.L.; K=2)                                                          |
|---------------|-----------------------------------------------------------------|-----------------------------|-----------------------------------------------------------------------------------------------|
| 1             | Nitroglycerin<br>CAS # 55-63-0<br>Purity 99%<br>(Lot 200507JLM) | 5,016.0 µg/mL               | +/- 46.6461 µg/mL Gravimetric<br>+/- 277.1256 µg/mL Unstressed<br>+/- 322.4378 µg/mL Stressed |

Solvent: Acetonitrile  
CAS # 75-05-8  
Purity 99%

**Column:**

250mm x 4.6mm  
Ultra C18 (cat.# 9174575)

**Flow Rate:**

1.0 ml/min.

**Mobile Phase A:**


water:methanol (44:56 V/V)

**Mobile Phase B:****Mobile Phase Composition:**

100%A

**Det. Type:**

Wavelength: 210nm & 254nm



This chromatogram represents a general set of testing conditions chosen for product acceptance. For optimal results in your lab, conditions should be adjusted for your specific instrument, method, and application.

  
Matt Fragassi - Mix Technician

Date Mixed: 15-Aug-2022

Balance: 1128353505

  
Jennifer Pollino - Operations Tech III - ARM QC

Date Passed: 16-Aug-2022



Manufactured under Restek's ISO 9001:2015  
Registered Quality System  
Certificate #FM 80397

## General Certified Reference Material Notes

### Expiration Notes:

- Expiration date valid for unopened ampul stored in compliance with the recommended conditions.
- Uncertainty, concentration, and expiration of the CRM are based on the unopened product being stored according to the recommended condition found in the storage field.

### Purity Notes:

- Purity and/or chemical identity are determined by one or more of the following techniques: GC/FID, HPLC, GC/μECD, GC/MS, LC/MS, RI, and/or melting point.
- Compounds with a listed purity of less than 99% have been weight corrected to compensate for impurities and/or salts. A correction factor is used to calculate the amount of compound necessary to achieve the desired concentration of the parent compound in solution.
- Purity of isomeric compounds is reported as the sum of the isomers.
- Purity values are rounded to the nearest whole number.

### Certified Uncertainty Value Notes:

- The uncertainties are determined in accordance with ISO 17034 and Guide 35. The certified combined stressed uncertainty value ( includes gravimetric uncertainty, homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty and were combined using the following formula:

$$U_{combined\ stressed} = k \sqrt{U_{gravimetric}^2 + U_{homogeneity}^2 + U_{storage\ stability}^2 + U_{shipping\ stability}^2}$$

$k$  is a coverage factor of 2, which gives a level of confidence of approximately 95%.

- It is important to note that the shipping stability uncertainty was obtained under temperature extremes for specific time intervals; therefore, the certified combined stressed uncertainty value should only be applied to the product if it was stored at non-standard temperature conditions up to and including 7 days. Contact Restek Technical Service at [www.restek.com/Contact-Us](http://www.restek.com/Contact-Us) for use recommendations if your shipment was in-transit for more than 7 days at non-standard temperature conditions.
- Apply the certified combined unstressed uncertainty value if the product was received under standard shipping conditions. Apply the certified combined stressed uncertainty value if the product was received under non-standard conditions as specified below.

| Label Conditions                                          | Standard Conditions | Non-Standard Conditions |
|-----------------------------------------------------------|---------------------|-------------------------|
| 25°C Nominal (Room Temperature)                           | < 60°C              | ≥ 60°C up to 7 days     |
| 10°C or colder (Refrigerate)                              | < 40°C              | ≥ 40°C up to 7 days     |
| 0°C or colder (Freezer)<br>-20°C or colder (Deep Freezer) | < 25°C              | ≥ 25°C up to 7 days     |

- Separate (not combined) uncertainty values for gravimetric uncertainty are also displayed on the certificate, if needed, separate homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty values are available by contacting Restek Technical Service at [www.restek.com/Contact-Us](http://www.restek.com/Contact-Us).
- The packaged amount is the minimum sample size for which uncertainty is valid. The ampules are over-filled to ensure that the minimum packaged amount can be sufficiently transferred.

### Manufacturing Notes:

- Concentration is based upon gravimetric preparation using either a balance whose calibration has been verified daily using NIST traceable weights, and/or dilutions with Class A glassware.

### Handling Notes:

- Stability of the unopened product, when stored in compliance with the recommended conditions, is guaranteed through the expiration displayed on the product label and certificate. Contact Restek for additional opened product stability information, with the knowledge/understanding that open product stability is subject to the specific handling and environmental conditions to which the product is exposed. For your convenience Restek supplies deactivated vials with most standards packed in 2mL ampules. Larger volume deactivated vials are available through Restek as a custom ordered item. Additionally, Restek sells DMDCS for the purpose of glassware deactivation as catalog number 31861, which includes complete instructions.





Reagent

---

**8330\_NG\_Stk\_00120**



110 Benner Circle  
Bellefonte, PA 16823-8812  
Tel: (800)356-1688  
Fax: (814)353-1309

www.restek.com

# CERTIFIED REFERENCE MATERIAL

## Certificate of Composition



### FOR LABORATORY USE ONLY-READ SDS PRIOR TO USE.

*This Reference Material is intended for Laboratory Use Only as a standard for the qualitative and/or quantitative determination of the analyte(s) listed.*

Catalog No. : 568871 Lot No.: A0188553  
Description : Custom Nitroglycerin Standard  
Custom Nitroglycerin Standard 5,000µg/mL, Acetonitrile, 1mL/ampul  
Container Size : 2 mL Pkg Amt: > 1 mL  
Expiration Date : August 31, 2025 Storage: 10°C or colder  
Ship: Ambient

### CERTIFIED VALUES

| Elution Order | Compound                      | Grav. Conc. (weight/volume) | Expanded Uncertainty (95% C.L.; K=2) |          |       |             |
|---------------|-------------------------------|-----------------------------|--------------------------------------|----------|-------|-------------|
| 1             | Nitroglycerin                 | 5,016.0 µg/mL               | +/-                                  | 46.6461  | µg/mL | Gravimetric |
|               | CAS # 55-63-0 (Lot 200507JLM) |                             | +/-                                  | 277.1256 | µg/mL | Unstressed  |
|               | Purity 99%                    |                             | +/-                                  | 322.4378 | µg/mL | Stressed    |

Solvent: Acetonitrile  
CAS # 75-05-8  
Purity 99%

**Column:**

250mm x 4.6mm

Ultra C18 (cat.# 9174575)

**Flow Rate:**

1.0 ml/min.

**Mobile Phase A:**


water:methanol (44:56 V/V)

**Mobile Phase B:****Mobile Phase Composition:**

100%A

**Det. Type:**

Wavelength: 210nm &amp; 254nm



This chromatogram represents a general set of testing conditions chosen for product acceptance. For optimal results in your lab, conditions should be adjusted for your specific instrument, method, and application.

  
Matt Fragassi - Mix Technician

Date Mixed: 15-Aug-2022

Balance: 1128353505

  
Jennifer Pollino - Operations Tech III - ARM QC

Date Passed: 16-Aug-2022



Manufactured under Restek's ISO 9001:2015  
Registered Quality System  
Certificate #FM 80397

## General Certified Reference Material Notes

### Expiration Notes:

- Expiration date valid for unopened ampul stored in compliance with the recommended conditions.
- Uncertainty, concentration, and expiration of the CRM are based on the unopened product being stored according to the recommended condition found in the storage field.

### Purity Notes:

- Purity and/or chemical identity are determined by one or more of the following techniques: GC/FID, HPLC, GC/μECD, GC/MS, LC/MS, RI, and/or melting point.
- Compounds with a listed purity of less than 99% have been weight corrected to compensate for impurities and/or salts. A correction factor is used to calculate the amount of compound necessary to achieve the desired concentration of the parent compound in solution.
- Purity of isomeric compounds is reported as the sum of the isomers.
- Purity values are rounded to the nearest whole number.

### Certified Uncertainty Value Notes:

- The uncertainties are determined in accordance with ISO 17034 and Guide 35. The certified combined stressed uncertainty value ( includes gravimetric uncertainty, homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty and were combined using the following formula:

$$U_{combined\ stressed} = k \sqrt{U_{gravimetric}^2 + U_{homogeneity}^2 + U_{storage\ stability}^2 + U_{shipping\ stability}^2}$$

$k$  is a coverage factor of 2, which gives a level of confidence of approximately 95%.

- It is important to note that the shipping stability uncertainty was obtained under temperature extremes for specific time intervals; therefore, the certified combined stressed uncertainty value should only be applied to the product if it was stored at non-standard temperature conditions up to and including 7 days. Contact Restek Technical Service at [www.restek.com/Contact-Us](http://www.restek.com/Contact-Us) for use recommendations if your shipment was in-transit for more than 7 days at non-standard temperature conditions.
- Apply the certified combined unstressed uncertainty value if the product was received under standard shipping conditions. Apply the certified combined stressed uncertainty value if the product was received under non-standard conditions as specified below.

| Label Conditions                                          | Standard Conditions | Non-Standard Conditions |
|-----------------------------------------------------------|---------------------|-------------------------|
| 25°C Nominal (Room Temperature)                           | < 60°C              | ≥ 60°C up to 7 days     |
| 10°C or colder (Refrigerate)                              | < 40°C              | ≥ 40°C up to 7 days     |
| 0°C or colder (Freezer)<br>-20°C or colder (Deep Freezer) | < 25°C              | ≥ 25°C up to 7 days     |

- Separate (not combined) uncertainty values for gravimetric uncertainty are also displayed on the certificate, if needed, separate homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty values are available by contacting Restek Technical Service at [www.restek.com/Contact-Us](http://www.restek.com/Contact-Us).
- The packaged amount is the minimum sample size for which uncertainty is valid. The ampules are over-filled to ensure that the minimum packaged amount can be sufficiently transferred.

### Manufacturing Notes:

- Concentration is based upon gravimetric preparation using either a balance whose calibration has been verified daily using NIST traceable weights, and/or dilutions with Class A glassware.

### Handling Notes:

- Stability of the unopened product, when stored in compliance with the recommended conditions, is guaranteed through the expiration displayed on the product label and certificate. Contact Restek for additional opened product stability information, with the knowledge/understanding that open product stability is subject to the specific handling and environmental conditions to which the product is exposed. For your convenience Restek supplies deactivated vials with most standards packed in 2mL ampules. Larger volume deactivated vials are available through Restek as a custom ordered item. Additionally, Restek sells DMDCS for the purpose of glassware deactivation as catalog number 31861, which includes complete instructions.



Reagent

---

**8330\_NG\_Stk\_00128**



110 Benner Circle  
Bellefonte, PA 16823-8812  
Tel: 1-814-353-1300  
Fax: 1-814-353-1309

www.restek.com

## CERTIFIED REFERENCE MATERIAL

# Certificate of Analysis

chromatographic



### FOR LABORATORY USE ONLY-READ SDS PRIOR TO USE.

*This Reference Material is intended for Laboratory Use Only as a standard for the qualitative and/or quantitative determination of the analyte(s) listed.*

**Catalog No. :** 568871 **Lot No.:** A0194013

**Description :** Custom Nitroglycerin Standard  
Custom Nitroglycerin Standard 5,000µg/mL, Acetonitrile, 1mL/ampul

**Container Size :** 2 mL **Pkg Amt:** > 1 mL

**Expiration Date :** January 31, 2026 **Storage:** 10°C or colder

**Ship:** Ambient

### CERTIFIED VALUES

| Elution Order | Compound      | CAS #   | Lot #     | Purity | Grav. Conc. (weight/volume) | Expanded Uncertainty * (95% C.L.; K=2) |
|---------------|---------------|---------|-----------|--------|-----------------------------|----------------------------------------|
| 1             | Nitroglycerin | 55-63-0 | 200507JLM | 99%    | 5,036.0 µg/mL               | +/- 237.6858                           |

\* Expanded Uncertainty displayed in same units as Grav. Conc.

**Solvent:** Acetonitrile  
**CAS #** 75-05-8  
**Purity** 99%



## Quality Confirmation Test

**Column:**

250mm x 4.6mm  
Ultra C18 (cat.# 9174575)

**Flow Rate:**

1.0 ml/min.

**Mobile Phase A:**

water:methanol (44:56 V/V)

**Mobile Phase B:****Mobile Phase Composition:**


100%A

**Det. Type:**

Wavelength: 210nm & 254nm

**Inj. Vol**

5µl



This chromatogram represents a general set of testing conditions chosen for product acceptance. For optimal results in your lab, conditions should be adjusted for your specific instrument, method, and application.

Brandon Reish - Operations Technician II

Date Mixed: 27-Jan-2023

Balance Serial # B251644995

Christie Mills - Operations Tech II - ARM QC

Date Passed: 31-Jan-2023

Manufactured under Restek's ISO 9001:2015  
Registered Quality System  
Certificate #FM 80397

## General Certified Reference Material Notes

### Expiration Notes:

- Expiration date valid for unopened ampul stored in compliance with the recommended conditions.
- Uncertainty, concentration, and expiration of the CRM are based on the unopened product being stored according to the recommended condition found in the storage field.

### Purity Notes:

- Purity and/or chemical identity are determined by one or more of the following techniques: GC/FID, HPLC, GC/μECD, GC/MS, LC/MS, RI, and/or melting point.
- Compounds with a listed purity of less than 99% have been weight corrected to compensate for impurities and/or salts. A correction factor is used to calculate the amount of compound necessary to achieve the desired concentration of the parent compound in solution.
- Purity of isomeric compounds is reported as the sum of the isomers.
- Purity values are rounded to the nearest whole number.

### Certified Uncertainty Value Notes:

- The uncertainties are determined in accordance with ISO 17034 and Guide 35. The certified expanded uncertainty value includes gravimetric uncertainty, homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty and were combined using the following formula:

$$U_{combined\ uncertainty} = k \sqrt{u_{gravimetric}^2 + u_{homogeneity}^2 + u_{storage\ stability}^2 + u_{shipping\ stability}^2}$$

$k$  is a coverage factor of 2, which gives a level of confidence of approximately 95%.

- The packaged amount is the minimum sample size for which uncertainty is valid. The ampuls are over-filled to ensure that the minimum packaged amount can be sufficiently transferred.

### Manufacturing Notes:

- Concentration is based upon gravimetric preparation using either a balance whose calibration has been verified daily using NIST traceable weights, and/or dilutions with Class A glassware.

### Handling Notes:

- Stability of the unopened product, when stored in compliance with the recommended conditions, is guaranteed through the expiration displayed on the product label and certificate. Contact Restek for additional opened product stability information, with the knowledge/understanding that open product stability is subject to the specific handling and environmental conditions to which the product is exposed. For your convenience Restek supplies deactivated vials with most standards packed in 2mL ampuls. Larger volume deactivated vials are available through Restek as a custom ordered item. Additionally, Restek sells DMDCS for the purpose of glassware deactivation as catalog number 31861, which includes complete instructions.
- If any undissolved material is visible inside the ampul, sonicate the unopened ampul until the material is completely dissolved.

Reagent

---

**8330\_NG\_Stk\_00129**



110 Benner Circle  
Bellefonte, PA 16823-8812  
Tel: 1-814-353-1300  
Fax: 1-814-353-1309

www.restek.com

## CERTIFIED REFERENCE MATERIAL

# Certificate of Analysis

chromatographic



### FOR LABORATORY USE ONLY-READ SDS PRIOR TO USE.

*This Reference Material is intended for Laboratory Use Only as a standard for the qualitative and/or quantitative determination of the analyte(s) listed.*

**Catalog No. :** 568871 **Lot No.:** A0194013

**Description :** Custom Nitroglycerin Standard  
Custom Nitroglycerin Standard 5,000µg/mL, Acetonitrile, 1mL/ampul

**Container Size :** 2 mL **Pkg Amt:** > 1 mL

**Expiration Date :** January 31, 2026 **Storage:** 10°C or colder

**Ship:** Ambient

### CERTIFIED VALUES

| Elution Order | Compound      | CAS #   | Lot #     | Purity | Grav. Conc. (weight/volume) | Expanded Uncertainty * (95% C.L.; K=2) |
|---------------|---------------|---------|-----------|--------|-----------------------------|----------------------------------------|
| 1             | Nitroglycerin | 55-63-0 | 200507JLM | 99%    | 5,036.0 µg/mL               | +/- 237.6858                           |

\* Expanded Uncertainty displayed in same units as Grav. Conc.

**Solvent:** Acetonitrile  
**CAS #** 75-05-8  
**Purity** 99%

## Quality Confirmation Test

**Column:**

250mm x 4.6mm  
Ultra C18 (cat.# 9174575)

**Flow Rate:**

1.0 ml/min.

**Mobile Phase A:**

water:methanol (44:56 V/V)

**Mobile Phase B:****Mobile Phase Composition:**


100%A

**Det. Type:**

Wavelength: 210nm & 254nm

**Inj. Vol**

5µl



This chromatogram represents a general set of testing conditions chosen for product acceptance. For optimal results in your lab, conditions should be adjusted for your specific instrument, method, and application.

Brandon Reish - Operations Technician II

Date Mixed: 27-Jan-2023

Balance Serial # B251644995

Christie Mills - Operations Tech II - ARM QC

Date Passed: 31-Jan-2023

Manufactured under Restek's ISO 9001:2015  
Registered Quality System  
Certificate #FM 80397

## General Certified Reference Material Notes

### Expiration Notes:

- Expiration date valid for unopened ampul stored in compliance with the recommended conditions.
- Uncertainty, concentration, and expiration of the CRM are based on the unopened product being stored according to the recommended condition found in the storage field.

### Purity Notes:

- Purity and/or chemical identity are determined by one or more of the following techniques: GC/FID, HPLC, GC/μECD, GC/MS, LC/MS, RI, and/or melting point.
- Compounds with a listed purity of less than 99% have been weight corrected to compensate for impurities and/or salts. A correction factor is used to calculate the amount of compound necessary to achieve the desired concentration of the parent compound in solution.
- Purity of isomeric compounds is reported as the sum of the isomers.
- Purity values are rounded to the nearest whole number.

### Certified Uncertainty Value Notes:

- The uncertainties are determined in accordance with ISO 17034 and Guide 35. The certified expanded uncertainty value includes gravimetric uncertainty, homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty and were combined using the following formula:

$$U_{combined\ uncertainty} = k \sqrt{u_{gravimetric}^2 + u_{homogeneity}^2 + u_{storage\ stability}^2 + u_{shipping\ stability}^2}$$

$k$  is a coverage factor of 2, which gives a level of confidence of approximately 95%.

- The packaged amount is the minimum sample size for which uncertainty is valid. The ampuls are over-filled to ensure that the minimum packaged amount can be sufficiently transferred.

### Manufacturing Notes:

- Concentration is based upon gravimetric preparation using either a balance whose calibration has been verified daily using NIST traceable weights, and/or dilutions with Class A glassware.

### Handling Notes:

- Stability of the unopened product, when stored in compliance with the recommended conditions, is guaranteed through the expiration displayed on the product label and certificate. Contact Restek for additional opened product stability information, with the knowledge/understanding that open product stability is subject to the specific handling and environmental conditions to which the product is exposed. For your convenience Restek supplies deactivated vials with most standards packed in 2mL ampuls. Larger volume deactivated vials are available through Restek as a custom ordered item. Additionally, Restek sells DMDCS for the purpose of glassware deactivation as catalog number 31861, which includes complete instructions.
- If any undissolved material is visible inside the ampul, sonicate the unopened ampul until the material is completely dissolved.

Reagent

---

**8330\_NG\_Stk\_00130**



110 Benner Circle  
Bellefonte, PA 16823-8812  
Tel: 1-814-353-1300  
Fax: 1-814-353-1309

www.restek.com

## CERTIFIED REFERENCE MATERIAL

# Certificate of Analysis

chromatographic



### FOR LABORATORY USE ONLY-READ SDS PRIOR TO USE.

*This Reference Material is intended for Laboratory Use Only as a standard for the qualitative and/or quantitative determination of the analyte(s) listed.*

**Catalog No. :** 568871 **Lot No.:** A0194013

**Description :** Custom Nitroglycerin Standard  
Custom Nitroglycerin Standard 5,000µg/mL, Acetonitrile, 1mL/ampul

**Container Size :** 2 mL **Pkg Amt:** > 1 mL

**Expiration Date :** January 31, 2026 **Storage:** 10°C or colder

**Ship:** Ambient

### CERTIFIED VALUES

| Elution Order | Compound      | CAS #   | Lot #     | Purity | Grav. Conc. (weight/volume) | Expanded Uncertainty * (95% C.L.; K=2) |
|---------------|---------------|---------|-----------|--------|-----------------------------|----------------------------------------|
| 1             | Nitroglycerin | 55-63-0 | 200507JLM | 99%    | 5,036.0 µg/mL               | +/- 237.6858                           |

\* Expanded Uncertainty displayed in same units as Grav. Conc.

**Solvent:** Acetonitrile  
**CAS #** 75-05-8  
**Purity** 99%



## Quality Confirmation Test

**Column:**

250mm x 4.6mm  
Ultra C18 (cat.# 9174575)

**Flow Rate:**

1.0 ml/min.

**Mobile Phase A:**

water:methanol (44:56 V/V)

**Mobile Phase B:****Mobile Phase Composition:**


100%A

**Det. Type:**

Wavelength: 210nm & 254nm

**Inj. Vol**

5µl



This chromatogram represents a general set of testing conditions chosen for product acceptance. For optimal results in your lab, conditions should be adjusted for your specific instrument, method, and application.

Brandon Reish - Operations Technician II

Date Mixed: 27-Jan-2023

Balance Serial # B251644995

Christie Mills - Operations Tech II - ARM QC

Date Passed: 31-Jan-2023

Manufactured under Restek's ISO 9001:2015  
Registered Quality System  
Certificate #FM 80397

## General Certified Reference Material Notes

### Expiration Notes:

- Expiration date valid for unopened ampul stored in compliance with the recommended conditions.
- Uncertainty, concentration, and expiration of the CRM are based on the unopened product being stored according to the recommended condition found in the storage field.

### Purity Notes:

- Purity and/or chemical identity are determined by one or more of the following techniques: GC/FID, HPLC, GC/μECD, GC/MS, LC/MS, RI, and/or melting point.
- Compounds with a listed purity of less than 99% have been weight corrected to compensate for impurities and/or salts. A correction factor is used to calculate the amount of compound necessary to achieve the desired concentration of the parent compound in solution.
- Purity of isomeric compounds is reported as the sum of the isomers.
- Purity values are rounded to the nearest whole number.

### Certified Uncertainty Value Notes:

- The uncertainties are determined in accordance with ISO 17034 and Guide 35. The certified expanded uncertainty value includes gravimetric uncertainty, homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty and were combined using the following formula:

$$U_{combined\ uncertainty} = k \sqrt{u_{gravimetric}^2 + u_{homogeneity}^2 + u_{storage\ stability}^2 + u_{shipping\ stability}^2}$$

$k$  is a coverage factor of 2, which gives a level of confidence of approximately 95%.

- The packaged amount is the minimum sample size for which uncertainty is valid. The ampuls are over-filled to ensure that the minimum packaged amount can be sufficiently transferred.

### Manufacturing Notes:

- Concentration is based upon gravimetric preparation using either a balance whose calibration has been verified daily using NIST traceable weights, and/or dilutions with Class A glassware.

### Handling Notes:

- Stability of the unopened product, when stored in compliance with the recommended conditions, is guaranteed through the expiration displayed on the product label and certificate. Contact Restek for additional opened product stability information, with the knowledge/understanding that open product stability is subject to the specific handling and environmental conditions to which the product is exposed. For your convenience Restek supplies deactivated vials with most standards packed in 2mL ampuls. Larger volume deactivated vials are available through Restek as a custom ordered item. Additionally, Restek sells DMDCS for the purpose of glassware deactivation as catalog number 31861, which includes complete instructions.
- If any undissolved material is visible inside the ampul, sonicate the unopened ampul until the material is completely dissolved.

Reagent

---

**8330\_NG\_Stk\_00131**



110 Benner Circle  
Bellefonte, PA 16823-8812  
Tel: 1-814-353-1300  
Fax: 1-814-353-1309

www.restek.com

## CERTIFIED REFERENCE MATERIAL

# Certificate of Analysis

chromatographic



### FOR LABORATORY USE ONLY-READ SDS PRIOR TO USE.

*This Reference Material is intended for Laboratory Use Only as a standard for the qualitative and/or quantitative determination of the analyte(s) listed.*

**Catalog No. :** 568871 **Lot No.:** A0194013

**Description :** Custom Nitroglycerin Standard  
Custom Nitroglycerin Standard 5,000µg/mL, Acetonitrile, 1mL/ampul

**Container Size :** 2 mL **Pkg Amt:** > 1 mL

**Expiration Date :** January 31, 2026 **Storage:** 10°C or colder

**Ship:** Ambient

### CERTIFIED VALUES

| Elution Order | Compound      | CAS #   | Lot #     | Purity | Grav. Conc. (weight/volume) | Expanded Uncertainty * (95% C.L.; K=2) |
|---------------|---------------|---------|-----------|--------|-----------------------------|----------------------------------------|
| 1             | Nitroglycerin | 55-63-0 | 200507JLM | 99%    | 5,036.0 µg/mL               | +/- 237.6858                           |

\* Expanded Uncertainty displayed in same units as Grav. Conc.

**Solvent:** Acetonitrile  
**CAS #** 75-05-8  
**Purity** 99%

## Quality Confirmation Test

**Column:**

250mm x 4.6mm  
Ultra C18 (cat.# 9174575)

**Flow Rate:**

1.0 ml/min.

**Mobile Phase A:**

water:methanol (44:56 V/V)

**Mobile Phase B:****Mobile Phase Composition:**


100%A

**Det. Type:**

Wavelength: 210nm & 254nm

**Inj. Vol**

5µl



This chromatogram represents a general set of testing conditions chosen for product acceptance. For optimal results in your lab, conditions should be adjusted for your specific instrument, method, and application.

Brandon Reish - Operations Technician II

Date Mixed: 27-Jan-2023

Balance Serial # B251644995

Christie Mills - Operations Tech II - ARM QC

Date Passed: 31-Jan-2023

Manufactured under Restek's ISO 9001:2015  
Registered Quality System  
Certificate #FM 80397

## General Certified Reference Material Notes

### Expiration Notes:

- Expiration date valid for unopened ampul stored in compliance with the recommended conditions.
- Uncertainty, concentration, and expiration of the CRM are based on the unopened product being stored according to the recommended condition found in the storage field.

### Purity Notes:

- Purity and/or chemical identity are determined by one or more of the following techniques: GC/FID, HPLC, GC/μECD, GC/MS, LC/MS, RI, and/or melting point.
- Compounds with a listed purity of less than 99% have been weight corrected to compensate for impurities and/or salts. A correction factor is used to calculate the amount of compound necessary to achieve the desired concentration of the parent compound in solution.
- Purity of isomeric compounds is reported as the sum of the isomers.
- Purity values are rounded to the nearest whole number.

### Certified Uncertainty Value Notes:

- The uncertainties are determined in accordance with ISO 17034 and Guide 35. The certified expanded uncertainty value includes gravimetric uncertainty, homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty and were combined using the following formula:

$$U_{combined\ uncertainty} = k \sqrt{u_{gravimetric}^2 + u_{homogeneity}^2 + u_{storage\ stability}^2 + u_{shipping\ stability}^2}$$

$k$  is a coverage factor of 2, which gives a level of confidence of approximately 95%.

- The packaged amount is the minimum sample size for which uncertainty is valid. The ampuls are over-filled to ensure that the minimum packaged amount can be sufficiently transferred.

### Manufacturing Notes:

- Concentration is based upon gravimetric preparation using either a balance whose calibration has been verified daily using NIST traceable weights, and/or dilutions with Class A glassware.

### Handling Notes:

- Stability of the unopened product, when stored in compliance with the recommended conditions, is guaranteed through the expiration displayed on the product label and certificate. Contact Restek for additional opened product stability information, with the knowledge/understanding that open product stability is subject to the specific handling and environmental conditions to which the product is exposed. For your convenience Restek supplies deactivated vials with most standards packed in 2mL ampuls. Larger volume deactivated vials are available through Restek as a custom ordered item. Additionally, Restek sells DMDCS for the purpose of glassware deactivation as catalog number 31861, which includes complete instructions.
- If any undissolved material is visible inside the ampul, sonicate the unopened ampul until the material is completely dissolved.

Reagent

---

**8330\_NG1000\_00007**



110 Benner Circle  
Bellefonte, PA 16823-8812  
Tel: (800)356-1688  
Fax: (814)353-1309

www.restek.com

## CERTIFIED REFERENCE MATERIAL

### Certificate of Analysis



#### FOR LABORATORY USE ONLY-READ SDS PRIOR TO USE.

*This Reference Material is intended for Laboratory Use Only as a standard for the qualitative and/or quantitative determination of the analyte(s) listed.*

**Catalog No. :** 31498 **Lot No.:** A0175997  
**Description :** Nitroglycerin Standard  
Nitroglycerin Standard 1,000µg/mL, Methanol, 1mL/ampul  
**Container Size :** 2 mL **Pkg Amt:** > 1 mL  
**Expiration Date :** September 30, 2026 **Storage:** 10°C or colder  
**Ship:** Ambient

#### CERTIFIED VALUES

| Elution Order | Compound                                                        | Grav. Conc. (weight/volume) | Expanded Uncertainty (95% C.L.; K=2)                                                       |
|---------------|-----------------------------------------------------------------|-----------------------------|--------------------------------------------------------------------------------------------|
| 1             | Nitroglycerin<br>CAS # 55-63-0<br>Purity 99%<br>(Lot 200507JLM) | 1,000.0 µg/mL               | +/- 5.9397 µg/mL Gravimetric<br>+/- 54.7830 µg/mL Unstressed<br>+/- 63.8824 µg/mL Stressed |

**Solvent:** Methanol  
CAS # 67-56-1  
Purity 99%



## General Certified Reference Material Notes

### Expiration Notes:

- Expiration date valid for unopened ampul stored in compliance with the recommended conditions.
- Uncertainty, concentration, and expiration of the CRM are based on the unopened product being stored according to the recommended condition found in the storage field.

### Purity Notes:

- Purity and/or chemical identity are determined by one or more of the following techniques: GC/FID, HPLC, GC/μECD, GC/MS, LC/MS, RI, and/or melting point.
- Compounds with a listed purity of less than 99% have been weight corrected to compensate for impurities and/or salts. A correction factor is used to calculate the amount of compound necessary to achieve the desired concentration of the parent compound in solution.
- Purity of isomeric compounds is reported as the sum of the isomers.
- Purity values are rounded to the nearest whole number.

### Certified Uncertainty Value Notes:

- The uncertainties are determined in accordance with ISO 17034 and Guide 35. The certified combined stressed uncertainty value ( includes gravimetric uncertainty, homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty and were combined using the following formula:

$$U_{combined\ stressed} = k \sqrt{U_{gravimetric}^2 + U_{homogeneity}^2 + U_{storage\ stability}^2 + U_{shipping\ stability}^2}$$

$k$  is a coverage factor of 2, which gives a level of confidence of approximately 95%.

- It is important to note that the shipping stability uncertainty was obtained under temperature extremes for specific time intervals; therefore, the certified combined stressed uncertainty value should only be applied to the product if it was stored at non-standard temperature conditions up to and including 7 days. Contact Restek Technical Service at [www.restek.com/Contact-Us](http://www.restek.com/Contact-Us) for use recommendations if your shipment was in-transit for more than 7 days at non-standard temperature conditions.
- Apply the certified combined unstressed uncertainty value if the product was received under standard shipping conditions. Apply the certified combined stressed uncertainty value if the product was received under non-standard conditions as specified below.

| Label Conditions                                          | Standard Conditions | Non-Standard Conditions |
|-----------------------------------------------------------|---------------------|-------------------------|
| 25°C Nominal (Room Temperature)                           | < 60°C              | ≥ 60°C up to 7 days     |
| 10°C or colder (Refrigerate)                              | < 40°C              | ≥ 40°C up to 7 days     |
| 0°C or colder (Freezer)<br>-20°C or colder (Deep Freezer) | < 25°C              | ≥ 25°C up to 7 days     |

- Separate (not combined) uncertainty values for gravimetric uncertainty are also displayed on the certificate, if needed, separate homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty values are available by contacting Restek Technical Service at [www.restek.com/Contact-Us](http://www.restek.com/Contact-Us).
- The packaged amount is the minimum sample size for which uncertainty is valid. The ampules are over-filled to ensure that the minimum packaged amount can be sufficiently transferred.

### Manufacturing Notes:

- Concentration is based upon gravimetric preparation using either a balance whose calibration has been verified daily using NIST traceable weights, and/or dilutions with Class A glassware.

### Handling Notes:

- Stability of the unopened product, when stored in compliance with the recommended conditions, is guaranteed through the expiration displayed on the product label and certificate. Contact Restek for additional opened product stability information, with the knowledge/understanding that open product stability is subject to the specific handling and environmental conditions to which the product is exposed. For your convenience Restek supplies deactivated vials with most standards packed in 2mL ampules. Larger volume deactivated vials are available through Restek as a custom ordered item. Additionally, Restek sells DMDCS for the purpose of glassware deactivation as catalog number 31861, which includes complete instructions.

Reagent

---

**8330\_NG1000\_00010**



110 Benner Circle  
Bellefonte, PA 16823-8812  
Tel: (800)356-1688  
Fax: (814)353-1309

www.restek.com

# CERTIFIED REFERENCE MATERIAL

## Certificate of Analysis



### FOR LABORATORY USE ONLY-READ SDS PRIOR TO USE.

*This Reference Material is intended for Laboratory Use Only as a standard for the qualitative and/or quantitative determination of the analyte(s) listed.*

Catalog No. : 31498 Lot No.: A0187280  
Description : Nitroglycerin Standard  
Nitroglycerin Standard 1,000µg/mL, Methanol, 1mL/ampul  
Container Size : 2 mL Pkg Amt: > 1 mL  
Expiration Date : July 31, 2027 Storage: 10°C or colder  
Ship: Ambient

### CERTIFIED VALUES

| Elution Order | Compound                      | Grav. Conc. (weight/volume) | Expanded Uncertainty (95% C.L.; K=2) |       |             |  |
|---------------|-------------------------------|-----------------------------|--------------------------------------|-------|-------------|--|
| 1             | Nitroglycerin                 | 1,004.0 µg/mL               | +/- 5.9635                           | µg/mL | Gravimetric |  |
|               | CAS # 55-63-0 (Lot 200507JLM) |                             | +/- 55.0021                          | µg/mL | Unstressed  |  |
|               | Purity 99%                    |                             | +/- 64.1379                          | µg/mL | Stressed    |  |

Solvent: Methanol  
CAS # 67-56-1  
Purity 99%

Reagent

---

**8330\_OP\_DMT\_00016**

**Preliminary Report**

Eurofins Denver

LCS, Lab Control Sample Report

Sample Path: \\chromfs\Denver\ChromData\CHHPLC\_X\20221230-117448.b\12300069.D

Lims ID: DMT LCS 16

Inj. Date: 30-Dec-2022 14:12:41

Worklist ID: 280-0117448-069

Instrument: CHHPLC\_X3

Method: 8330\_X3

| Compound | Amount Added | Amount Recovered | %Rec  | Limits 1<br>OB_Sonc_ |
|----------|--------------|------------------|-------|----------------------|
| 3 TNX    | 0.5000       | 0.5437           | 108.7 | 50-150               |
| 6 DNX    | 0.5020       | 0.5259           | 104.8 | 50-150               |
| 7 MNX    | 0.5850       | 0.6258           | 107.0 | 50-150               |

Samples for Limit Group: 1, Lims Prep Method: 8330B\_Sonc\_10g

|                   |                   |                   |
|-------------------|-------------------|-------------------|
| 280-170212-A-1-B  | 280-170212-A-2-D  | 280-170212-A-3-B  |
| 280-170212-A-4-B  | 280-170212-A-5-B  | 280-170212-A-6-B  |
| 280-170212-A-7-B  | 280-170212-A-8-B  | 280-170212-A-9-B  |
| 280-170212-A-10-D | 280-170212-A-11-B | 280-170212-A-12-C |
| 280-170212-A-13-B | 280-170212-A-14-B | 280-170212-A-15-B |
| 280-170212-A-16-B | 280-170212-A-17-B | 280-170212-A-18-B |
| 280-170212-A-19-B | 280-170212-A-20-B | 280-170212-A-21-C |
| 280-170212-A-22-B | 280-170212-A-23-B | 280-170212-A-24-F |
| 280-170212-A-25-B | 280-170212-A-26-C | 280-170212-A-27-B |
| 280-170212-A-28-B | 280-170212-A-29-B | 280-170212-A-30-B |
| 280-170212-A-31-B | 280-170212-A-32-B | 280-170212-A-33-B |
| 280-170212-A-34-B | 280-170212-A-35-B | 280-170212-A-36-C |
| 280-170212-A-37-B | 280-170212-A-38-B | 280-170212-A-39-B |
| 280-170212-A-40-B |                   |                   |

Reagent

---

**8330\_OP\_DMT\_00020**

**Preliminary Report**

Eurofins Denver

LCS, Lab Control Sample Report

Sample Path: \\chromfs\Denver\ChromData\CHHPLC\_X\20230509-121258.b\05090011.D

Lims ID: DMT LCS20

Inj. Date: 09-May-2023 16:32:44

Worklist ID: 280-0121258-011

Instrument: CHHPLC\_X3

Method: 8330\_X3

| Compound | Amount Added | Amount Recovered | %Rec  | Limits 1<br>3535 | Limits 2<br>0B_Sonc_ | Limits 3<br>3535 |
|----------|--------------|------------------|-------|------------------|----------------------|------------------|
| 3 TNX    | 0.5015       | 0.5099           | 101.7 | 50-150           |                      |                  |
| 6 DNX    | 0.5005       | 0.5059           | 101.1 | 66-119           | 50-150               |                  |
| 7 MNX    | 0.5835       | 0.6032           | 103.4 | 57-132           | 68-123               |                  |

Samples for Limit Group: 1, Lims Prep Method: 3535

280-175420-A-1-A

280-175420-A-2-A

280-175459-A-1-A

280-175942-A-3-A

280-175942-A-4-A

280-175942-A-5-A

Samples for Limit Group: 2, Lims Prep Method: 8330B\_Sonc\_10g

410-124886-B-1-A

550-201589-N-1-A

550-201670-N-1-A

Samples for Limit Group: 3, Lims Prep Method: 3535

280-175990-A-1-A

280-175990-A-2-A

280-175990-A-3-A

280-175990-A-4-A

860-48196-C-1-A

Reagent

---

**8330\_PETN\_Stk\_00128**





110 Benner Circle  
Bellefonte, PA 16823-8812  
Tel: (800)356-1688  
Fax: (814)353-1309

www.restek.com

## CERTIFIED REFERENCE MATERIAL

# Certificate of Composition



**FOR LABORATORY USE ONLY-READ SDS PRIOR TO USE.**

*This Reference Material is intended for Laboratory Use Only as a standard for the qualitative and/or quantitative determination of the analyte(s) listed.*

**Catalog No. :** 568872 **Lot No.:** A0187506  
**Description :** Custom PETN Standard  
Custom PETN Standard 5,000µg/mL, Acetonitrile, 1mL/ampul  
**Container Size :** 2 mL **Pkg Amt:** > 1 mL  
**Expiration Date :** July 31, 2025 **Storage:** 10°C or colder  
**Ship:** Ambient

### CERTIFIED VALUES

| Elution Order | Compound                                               | Grav. Conc. (weight/volume) | Expanded Uncertainty (95% C.L.; K=2)                                                          |
|---------------|--------------------------------------------------------|-----------------------------|-----------------------------------------------------------------------------------------------|
| 1             | PETN<br>CAS # 78-11-5<br>Purity 99%<br>(Lot 051108JLM) | 5,028.0 µg/mL               | +/- 46.7577 µg/mL Gravimetric<br>+/- 277.7886 µg/mL Unstressed<br>+/- 323.2092 µg/mL Stressed |

**Solvent:** Acetonitrile  
CAS # 75-05-8  
Purity 99%

**Column:**  
250mm x 4.6mm  
Ultra C18 (cat.# 9174575)


**Flow Rate:**  
1.0 ml/min.

**Mobile Phase A:**  
water:methanol (44:56 V/V)

**Mobile Phase B:**

**Mobile Phase Composition:**  
100%A

**Det. Type:**  
Wavelength: 210nm & 254nm



This chromatogram represents a general set of testing conditions chosen for product acceptance. For optimal results in your lab, conditions should be adjusted for your specific instrument, method, and application.

*Miranda Kline*  
Miranda Kline - Operations Technician I

**Date Mixed:** 19-Jul-2022

**Balance:** B345965662

*Fang-Yun Weaver*  
Fang-Yun Weaver - Operations Lead Tech - ARM QC

**Date Passed:** 22-Jul-2022

Manufactured under Restek's ISO 9001:2015  
Registered Quality System  
Certificate #FM 80397

## General Certified Reference Material Notes

### Expiration Notes:

- Expiration date valid for unopened ampul stored in compliance with the recommended conditions.
- Uncertainty, concentration, and expiration of the CRM are based on the unopened product being stored according to the recommended condition found in the storage field.

### Purity Notes:

- Purity and/or chemical identity are determined by one or more of the following techniques: GC/FID, HPLC, GC/μECD, GC/MS, LC/MS, RI, and/or melting point.
- Compounds with a listed purity of less than 99% have been weight corrected to compensate for impurities and/or salts. A correction factor is used to calculate the amount of compound necessary to achieve the desired concentration of the parent compound in solution.
- Purity of isomeric compounds is reported as the sum of the isomers.
- Purity values are rounded to the nearest whole number.

### Certified Uncertainty Value Notes:

- The uncertainties are determined in accordance with ISO 17034 and Guide 35. The certified combined stressed uncertainty value ( includes gravimetric uncertainty, homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty and were combined using the following formula:

$$U_{combined\ stressed} = k \sqrt{U_{gravimetric}^2 + U_{homogeneity}^2 + U_{storage\ stability}^2 + U_{shipping\ stability}^2}$$

$k$  is a coverage factor of 2, which gives a level of confidence of approximately 95%.

- It is important to note that the shipping stability uncertainty was obtained under temperature extremes for specific time intervals; therefore, the certified combined stressed uncertainty value should only be applied to the product if it was stored at non-standard temperature conditions up to and including 7 days. Contact Restek Technical Service at [www.restek.com/Contact-Us](http://www.restek.com/Contact-Us) for use recommendations if your shipment was in-transit for more than 7 days at non-standard temperature conditions.
- Apply the certified combined unstressed uncertainty value if the product was received under standard shipping conditions. Apply the certified combined stressed uncertainty value if the product was received under non-standard conditions as specified below.

| Label Conditions                                          | Standard Conditions | Non-Standard Conditions |
|-----------------------------------------------------------|---------------------|-------------------------|
| 25°C Nominal (Room Temperature)                           | < 60°C              | ≥ 60°C up to 7 days     |
| 10°C or colder (Refrigerate)                              | < 40°C              | ≥ 40°C up to 7 days     |
| 0°C or colder (Freezer)<br>-20°C or colder (Deep Freezer) | < 25°C              | ≥ 25°C up to 7 days     |

- Separate (not combined) uncertainty values for gravimetric uncertainty are also displayed on the certificate, if needed, separate homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty values are available by contacting Restek Technical Service at [www.restek.com/Contact-Us](http://www.restek.com/Contact-Us).
- The packaged amount is the minimum sample size for which uncertainty is valid. The ampules are over-filled to ensure that the minimum packaged amount can be sufficiently transferred.

### Manufacturing Notes:

- Concentration is based upon gravimetric preparation using either a balance whose calibration has been verified daily using NIST traceable weights, and/or dilutions with Class A glassware.

### Handling Notes:

- Stability of the unopened product, when stored in compliance with the recommended conditions, is guaranteed through the expiration displayed on the product label and certificate. Contact Restek for additional opened product stability information, with the knowledge/understanding that open product stability is subject to the specific handling and environmental conditions to which the product is exposed. For your convenience Restek supplies deactivated vials with most standards packed in 2mL ampules. Larger volume deactivated vials are available through Restek as a custom ordered item. Additionally, Restek sells DMDCS for the purpose of glassware deactivation as catalog number 31861, which includes complete instructions.

Reagent

---

**8330\_PETN\_Stk\_00133**



110 Benner Circle  
Bellefonte, PA 16823-8812  
Tel: (800)356-1688  
Fax: (814)353-1309

www.restek.com

## CERTIFIED REFERENCE MATERIAL

# Certificate of Composition



### FOR LABORATORY USE ONLY-READ SDS PRIOR TO USE.

*This Reference Material is intended for Laboratory Use Only as a standard for the qualitative and/or quantitative determination of the analyte(s) listed.*

**Catalog No.:** 568872 **Lot No.:** A0188550  
**Description:** Custom PETN Standard  
Custom PETN Standard 5,000µg/mL, Acetonitrile, 1mL/ampul  
**Container Size:** 2 mL **Pkg Amt:** > 1 mL  
**Expiration Date:** August 31, 2025 **Storage:** 10°C or colder  
**Ship:** Ambient

### CERTIFIED VALUES

| Elution Order | Compound                      | Grav. Conc. (weight/volume) | Expanded Uncertainty (95% C.L.; K=2) |
|---------------|-------------------------------|-----------------------------|--------------------------------------|
| 1             | PETN                          | 5,000.0 µg/mL               | +/- 46.4973 µg/mL Gravimetric        |
|               | CAS # 78-11-5 (Lot 051108JLM) |                             | +/- 276.2417 µg/mL Unstressed        |
|               | Purity 99%                    |                             | +/- 321.4093 µg/mL Stressed          |

**Solvent:** Acetonitrile  
**CAS #** 75-05-8  
**Purity** 99%

**Column:**

250mm x 4.6mm

Ultra C18 (cat.# 9174575)

**Flow Rate:**

1.0 ml/min.

**Mobile Phase A:**


water:methanol (44:56 V/V)

**Mobile Phase B:****Mobile Phase Composition:**

100%A

**Det. Type:**

Wavelength: 210nm &amp; 254nm




This chromatogram represents a general set of testing conditions chosen for product acceptance. For optimal results in your lab, conditions should be adjusted for your specific instrument, method, and application.

  
Matt Fragassi - Mix Technician

Date Mixed: 15-Aug-2022

Balance: 1128353505

  
Jennifer Pollino - Operations Tech III - ARM QC

Date Passed: 16-Aug-2022

Manufactured under Restek's ISO 9001:2015  
Registered Quality System  
Certificate #FM 80397

## General Certified Reference Material Notes

### Expiration Notes:

- Expiration date valid for unopened ampul stored in compliance with the recommended conditions.
- Uncertainty, concentration, and expiration of the CRM are based on the unopened product being stored according to the recommended condition found in the storage field.

### Purity Notes:

- Purity and/or chemical identity are determined by one or more of the following techniques: GC/FID, HPLC, GC/μECD, GC/MS, LC/MS, RI, and/or melting point.
- Compounds with a listed purity of less than 99% have been weight corrected to compensate for impurities and/or salts. A correction factor is used to calculate the amount of compound necessary to achieve the desired concentration of the parent compound in solution.
- Purity of isomeric compounds is reported as the sum of the isomers.
- Purity values are rounded to the nearest whole number.

### Certified Uncertainty Value Notes:

- The uncertainties are determined in accordance with ISO 17034 and Guide 35. The certified combined stressed uncertainty value ( includes gravimetric uncertainty, homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty and were combined using the following formula:

$$U_{combined\ stressed} = k \sqrt{U_{gravimetric}^2 + U_{homogeneity}^2 + U_{storage\ stability}^2 + U_{shipping\ stability}^2}$$

$k$  is a coverage factor of 2, which gives a level of confidence of approximately 95%.

- It is important to note that the shipping stability uncertainty was obtained under temperature extremes for specific time intervals; therefore, the certified combined stressed uncertainty value should only be applied to the product if it was stored at non-standard temperature conditions up to and including 7 days. Contact Restek Technical Service at [www.restek.com/Contact-Us](http://www.restek.com/Contact-Us) for use recommendations if your shipment was in-transit for more than 7 days at non-standard temperature conditions.
- Apply the certified combined unstressed uncertainty value if the product was received under standard shipping conditions. Apply the certified combined stressed uncertainty value if the product was received under non-standard conditions as specified below.

| Label Conditions                                          | Standard Conditions | Non-Standard Conditions |
|-----------------------------------------------------------|---------------------|-------------------------|
| 25°C Nominal (Room Temperature)                           | < 60°C              | ≥ 60°C up to 7 days     |
| 10°C or colder (Refrigerate)                              | < 40°C              | ≥ 40°C up to 7 days     |
| 0°C or colder (Freezer)<br>-20°C or colder (Deep Freezer) | < 25°C              | ≥ 25°C up to 7 days     |

- Separate (not combined) uncertainty values for gravimetric uncertainty are also displayed on the certificate, if needed, separate homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty values are available by contacting Restek Technical Service at [www.restek.com/Contact-Us](http://www.restek.com/Contact-Us).
- The packaged amount is the minimum sample size for which uncertainty is valid. The ampules are over-filled to ensure that the minimum packaged amount can be sufficiently transferred.

### Manufacturing Notes:

- Concentration is based upon gravimetric preparation using either a balance whose calibration has been verified daily using NIST traceable weights, and/or dilutions with Class A glassware.

### Handling Notes:

- Stability of the unopened product, when stored in compliance with the recommended conditions, is guaranteed through the expiration displayed on the product label and certificate. Contact Restek for additional opened product stability information, with the knowledge/understanding that open product stability is subject to the specific handling and environmental conditions to which the product is exposed. For your convenience Restek supplies deactivated vials with most standards packed in 2mL ampules. Larger volume deactivated vials are available through Restek as a custom ordered item. Additionally, Restek sells DMDCS for the purpose of glassware deactivation as catalog number 31861, which includes complete instructions.





Reagent

---

**8330\_PETN\_Stk\_00134**



110 Benner Circle  
Bellefonte, PA 16823-8812  
Tel: (800)356-1688  
Fax: (814)353-1309

www.restek.com

## CERTIFIED REFERENCE MATERIAL

# Certificate of Composition



### FOR LABORATORY USE ONLY-READ SDS PRIOR TO USE.

*This Reference Material is intended for Laboratory Use Only as a standard for the qualitative and/or quantitative determination of the analyte(s) listed.*

**Catalog No.:** 568872 **Lot No.:** A0188550  
**Description:** Custom PETN Standard  
Custom PETN Standard 5,000µg/mL, Acetonitrile, 1mL/ampul  
**Container Size:** 2 mL **Pkg Amt:** > 1 mL  
**Expiration Date:** August 31, 2025 **Storage:** 10°C or colder  
**Ship:** Ambient

### CERTIFIED VALUES

| Elution Order | Compound                      | Grav. Conc. (weight/volume) | Expanded Uncertainty (95% C.L.; K=2) |
|---------------|-------------------------------|-----------------------------|--------------------------------------|
| 1             | PETN                          | 5,000.0 µg/mL               | +/- 46.4973 µg/mL Gravimetric        |
|               | CAS # 78-11-5 (Lot 051108JLM) |                             | +/- 276.2417 µg/mL Unstressed        |
|               | Purity 99%                    |                             | +/- 321.4093 µg/mL Stressed          |

**Solvent:** Acetonitrile  
**CAS #** 75-05-8  
**Purity** 99%

**Column:**

250mm x 4.6mm  
Ultra C18 (cat.# 9174575)

**Flow Rate:**

1.0 ml/min.

**Mobile Phase A:**


water:methanol (44:56 V/V)

**Mobile Phase B:****Mobile Phase Composition:**

100%A

**Det. Type:**

Wavelength: 210nm & 254nm




This chromatogram represents a general set of testing conditions chosen for product acceptance. For optimal results in your lab, conditions should be adjusted for your specific instrument, method, and application.

  
Matt Fragassi - Mix Technician

Date Mixed: 15-Aug-2022

Balance: 1128353505

  
Jennifer Pollino - Operations Tech III - ARM QC

Date Passed: 16-Aug-2022

Manufactured under Restek's ISO 9001:2015  
Registered Quality System  
Certificate #FM 80397

## General Certified Reference Material Notes

### Expiration Notes:

- Expiration date valid for unopened ampul stored in compliance with the recommended conditions.
- Uncertainty, concentration, and expiration of the CRM are based on the unopened product being stored according to the recommended condition found in the storage field.

### Purity Notes:

- Purity and/or chemical identity are determined by one or more of the following techniques: GC/FID, HPLC, GC/μECD, GC/MS, LC/MS, RI, and/or melting point.
- Compounds with a listed purity of less than 99% have been weight corrected to compensate for impurities and/or salts. A correction factor is used to calculate the amount of compound necessary to achieve the desired concentration of the parent compound in solution.
- Purity of isomeric compounds is reported as the sum of the isomers.
- Purity values are rounded to the nearest whole number.

### Certified Uncertainty Value Notes:

- The uncertainties are determined in accordance with ISO 17034 and Guide 35. The certified combined stressed uncertainty value ( includes gravimetric uncertainty, homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty and were combined using the following formula:

$$U_{combined\ stressed} = k \sqrt{U_{gravimetric}^2 + U_{homogeneity}^2 + U_{storage\ stability}^2 + U_{shipping\ stability}^2}$$

$k$  is a coverage factor of 2, which gives a level of confidence of approximately 95%.

- It is important to note that the shipping stability uncertainty was obtained under temperature extremes for specific time intervals; therefore, the certified combined stressed uncertainty value should only be applied to the product if it was stored at non-standard temperature conditions up to and including 7 days. Contact Restek Technical Service at [www.restek.com/Contact-Us](http://www.restek.com/Contact-Us) for use recommendations if your shipment was in-transit for more than 7 days at non-standard temperature conditions.
- Apply the certified combined unstressed uncertainty value if the product was received under standard shipping conditions. Apply the certified combined stressed uncertainty value if the product was received under non-standard conditions as specified below.

| Label Conditions                                          | Standard Conditions | Non-Standard Conditions |
|-----------------------------------------------------------|---------------------|-------------------------|
| 25°C Nominal (Room Temperature)                           | < 60°C              | ≥ 60°C up to 7 days     |
| 10°C or colder (Refrigerate)                              | < 40°C              | ≥ 40°C up to 7 days     |
| 0°C or colder (Freezer)<br>-20°C or colder (Deep Freezer) | < 25°C              | ≥ 25°C up to 7 days     |

- Separate (not combined) uncertainty values for gravimetric uncertainty are also displayed on the certificate, if needed, separate homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty values are available by contacting Restek Technical Service at [www.restek.com/Contact-Us](http://www.restek.com/Contact-Us).
- The packaged amount is the minimum sample size for which uncertainty is valid. The ampules are over-filled to ensure that the minimum packaged amount can be sufficiently transferred.

### Manufacturing Notes:

- Concentration is based upon gravimetric preparation using either a balance whose calibration has been verified daily using NIST traceable weights, and/or dilutions with Class A glassware.

### Handling Notes:

- Stability of the unopened product, when stored in compliance with the recommended conditions, is guaranteed through the expiration displayed on the product label and certificate. Contact Restek for additional opened product stability information, with the knowledge/understanding that open product stability is subject to the specific handling and environmental conditions to which the product is exposed. For your convenience Restek supplies deactivated vials with most standards packed in 2mL ampules. Larger volume deactivated vials are available through Restek as a custom ordered item. Additionally, Restek sells DMDCS for the purpose of glassware deactivation as catalog number 31861, which includes complete instructions.



Reagent

---

**8330\_PETN\_Stk\_00135**



110 Benner Circle  
Bellefonte, PA 16823-8812  
Tel: (800)356-1688  
Fax: (814)353-1309

www.restek.com

## CERTIFIED REFERENCE MATERIAL

# Certificate of Composition



### FOR LABORATORY USE ONLY-READ SDS PRIOR TO USE.

*This Reference Material is intended for Laboratory Use Only as a standard for the qualitative and/or quantitative determination of the analyte(s) listed.*

**Catalog No.:** 568872 **Lot No.:** A0188550  
**Description:** Custom PETN Standard  
Custom PETN Standard 5,000µg/mL, Acetonitrile, 1mL/ampul  
**Container Size:** 2 mL **Pkg Amt:** > 1 mL  
**Expiration Date:** August 31, 2025 **Storage:** 10°C or colder  
**Ship:** Ambient

### CERTIFIED VALUES

| Elution Order | Compound                      | Grav. Conc. (weight/volume) | Expanded Uncertainty (95% C.L.; K=2) |
|---------------|-------------------------------|-----------------------------|--------------------------------------|
| 1             | PETN                          | 5,000.0 µg/mL               | +/- 46.4973 µg/mL Gravimetric        |
|               | CAS # 78-11-5 (Lot 051108JLM) |                             | +/- 276.2417 µg/mL Unstressed        |
|               | Purity 99%                    |                             | +/- 321.4093 µg/mL Stressed          |

**Solvent:** Acetonitrile  
**CAS #** 75-05-8  
**Purity** 99%

**Column:**

250mm x 4.6mm

Ultra C18 (cat.# 9174575)

**Flow Rate:**

1.0 ml/min.

**Mobile Phase A:**


water:methanol (44:56 V/V)

**Mobile Phase B:****Mobile Phase Composition:**

100%A

**Det. Type:**

Wavelength: 210nm &amp; 254nm




This chromatogram represents a general set of testing conditions chosen for product acceptance. For optimal results in your lab, conditions should be adjusted for your specific instrument, method, and application.

  
Matt Fragassi - Mix Technician

Date Mixed: 15-Aug-2022

Balance: 1128353505

  
Jennifer Pollino - Operations Tech III - ARM QC

Date Passed: 16-Aug-2022

Manufactured under Restek's ISO 9001:2015  
Registered Quality System  
Certificate #FM 80397



## General Certified Reference Material Notes

### Expiration Notes:

- Expiration date valid for unopened ampul stored in compliance with the recommended conditions.
- Uncertainty, concentration, and expiration of the CRM are based on the unopened product being stored according to the recommended condition found in the storage field.

### Purity Notes:

- Purity and/or chemical identity are determined by one or more of the following techniques: GC/FID, HPLC, GC/μECD, GC/MS, LC/MS, RI, and/or melting point.
- Compounds with a listed purity of less than 99% have been weight corrected to compensate for impurities and/or salts. A correction factor is used to calculate the amount of compound necessary to achieve the desired concentration of the parent compound in solution.
- Purity of isomeric compounds is reported as the sum of the isomers.
- Purity values are rounded to the nearest whole number.

### Certified Uncertainty Value Notes:

- The uncertainties are determined in accordance with ISO 17034 and Guide 35. The certified combined stressed uncertainty value ( includes gravimetric uncertainty, homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty and were combined using the following formula:

$$U_{combined\ stressed} = k \sqrt{U_{gravimetric}^2 + U_{homogeneity}^2 + U_{storage\ stability}^2 + U_{shipping\ stability}^2}$$

$k$  is a coverage factor of 2, which gives a level of confidence of approximately 95%.

- It is important to note that the shipping stability uncertainty was obtained under temperature extremes for specific time intervals; therefore, the certified combined stressed uncertainty value should only be applied to the product if it was stored at non-standard temperature conditions up to and including 7 days. Contact Restek Technical Service at [www.restek.com/Contact-Us](http://www.restek.com/Contact-Us) for use recommendations if your shipment was in-transit for more than 7 days at non-standard temperature conditions.
- Apply the certified combined unstressed uncertainty value if the product was received under standard shipping conditions. Apply the certified combined stressed uncertainty value if the product was received under non-standard conditions as specified below.

| Label Conditions                                          | Standard Conditions | Non-Standard Conditions |
|-----------------------------------------------------------|---------------------|-------------------------|
| 25°C Nominal (Room Temperature)                           | < 60°C              | ≥ 60°C up to 7 days     |
| 10°C or colder (Refrigerate)                              | < 40°C              | ≥ 40°C up to 7 days     |
| 0°C or colder (Freezer)<br>-20°C or colder (Deep Freezer) | < 25°C              | ≥ 25°C up to 7 days     |

- Separate (not combined) uncertainty values for gravimetric uncertainty are also displayed on the certificate, if needed, separate homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty values are available by contacting Restek Technical Service at [www.restek.com/Contact-Us](http://www.restek.com/Contact-Us).
- The packaged amount is the minimum sample size for which uncertainty is valid. The ampules are over-filled to ensure that the minimum packaged amount can be sufficiently transferred.

### Manufacturing Notes:

- Concentration is based upon gravimetric preparation using either a balance whose calibration has been verified daily using NIST traceable weights, and/or dilutions with Class A glassware.

### Handling Notes:

- Stability of the unopened product, when stored in compliance with the recommended conditions, is guaranteed through the expiration displayed on the product label and certificate. Contact Restek for additional opened product stability information, with the knowledge/understanding that open product stability is subject to the specific handling and environmental conditions to which the product is exposed. For your convenience Restek supplies deactivated vials with most standards packed in 2mL ampules. Larger volume deactivated vials are available through Restek as a custom ordered item. Additionally, Restek sells DMDCS for the purpose of glassware deactivation as catalog number 31861, which includes complete instructions.



Reagent

---

**8330\_PETN\_Stk\_00136**



110 Benner Circle  
Bellefonte, PA 16823-8812  
Tel: (800)356-1688  
Fax: (814)353-1309

www.restek.com

## CERTIFIED REFERENCE MATERIAL

# Certificate of Composition



### FOR LABORATORY USE ONLY-READ SDS PRIOR TO USE.

*This Reference Material is intended for Laboratory Use Only as a standard for the qualitative and/or quantitative determination of the analyte(s) listed.*

**Catalog No.:** 568872 **Lot No.:** A0188550  
**Description:** Custom PETN Standard  
Custom PETN Standard 5,000µg/mL, Acetonitrile, 1mL/ampul  
**Container Size:** 2 mL **Pkg Amt:** > 1 mL  
**Expiration Date:** August 31, 2025 **Storage:** 10°C or colder  
**Ship:** Ambient

### CERTIFIED VALUES

| Elution Order | Compound                      | Grav. Conc. (weight/volume) | Expanded Uncertainty (95% C.L.; K=2) |
|---------------|-------------------------------|-----------------------------|--------------------------------------|
| 1             | PETN                          | 5,000.0 µg/mL               | +/- 46.4973 µg/mL Gravimetric        |
|               | CAS # 78-11-5 (Lot 051108JLM) |                             | +/- 276.2417 µg/mL Unstressed        |
|               | Purity 99%                    |                             | +/- 321.4093 µg/mL Stressed          |

**Solvent:** Acetonitrile  
**CAS #** 75-05-8  
**Purity** 99%

**Column:**

250mm x 4.6mm  
Ultra C18 (cat.# 9174575)

**Flow Rate:**

1.0 ml/min.

**Mobile Phase A:**


water:methanol (44:56 V/V)

**Mobile Phase B:****Mobile Phase Composition:**

100%A

**Det. Type:**

Wavelength: 210nm & 254nm




This chromatogram represents a general set of testing conditions chosen for product acceptance. For optimal results in your lab, conditions should be adjusted for your specific instrument, method, and application.

  
Matt Fragassi - Mix Technician

Date Mixed: 15-Aug-2022

Balance: 1128353505

  
Jennifer Pollino - Operations Tech III - ARM QC

Date Passed: 16-Aug-2022



Manufactured under Restek's ISO 9001:2015  
Registered Quality System  
Certificate #FM 80397

## General Certified Reference Material Notes

### Expiration Notes:

- Expiration date valid for unopened ampul stored in compliance with the recommended conditions.
- Uncertainty, concentration, and expiration of the CRM are based on the unopened product being stored according to the recommended condition found in the storage field.

### Purity Notes:

- Purity and/or chemical identity are determined by one or more of the following techniques: GC/FID, HPLC, GC/μECD, GC/MS, LC/MS, RI, and/or melting point.
- Compounds with a listed purity of less than 99% have been weight corrected to compensate for impurities and/or salts. A correction factor is used to calculate the amount of compound necessary to achieve the desired concentration of the parent compound in solution.
- Purity of isomeric compounds is reported as the sum of the isomers.
- Purity values are rounded to the nearest whole number.

### Certified Uncertainty Value Notes:

- The uncertainties are determined in accordance with ISO 17034 and Guide 35. The certified combined stressed uncertainty value ( includes gravimetric uncertainty, homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty and were combined using the following formula:

$$U_{combined\ stressed} = k \sqrt{U_{gravimetric}^2 + U_{homogeneity}^2 + U_{storage\ stability}^2 + U_{shipping\ stability}^2}$$

$k$  is a coverage factor of 2, which gives a level of confidence of approximately 95%.

- It is important to note that the shipping stability uncertainty was obtained under temperature extremes for specific time intervals; therefore, the certified combined stressed uncertainty value should only be applied to the product if it was stored at non-standard temperature conditions up to and including 7 days. Contact Restek Technical Service at [www.restek.com/Contact-Us](http://www.restek.com/Contact-Us) for use recommendations if your shipment was in-transit for more than 7 days at non-standard temperature conditions.
- Apply the certified combined unstressed uncertainty value if the product was received under standard shipping conditions. Apply the certified combined stressed uncertainty value if the product was received under non-standard conditions as specified below.

| Label Conditions                                          | Standard Conditions | Non-Standard Conditions |
|-----------------------------------------------------------|---------------------|-------------------------|
| 25°C Nominal (Room Temperature)                           | < 60°C              | ≥ 60°C up to 7 days     |
| 10°C or colder (Refrigerate)                              | < 40°C              | ≥ 40°C up to 7 days     |
| 0°C or colder (Freezer)<br>-20°C or colder (Deep Freezer) | < 25°C              | ≥ 25°C up to 7 days     |

- Separate (not combined) uncertainty values for gravimetric uncertainty are also displayed on the certificate, if needed, separate homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty values are available by contacting Restek Technical Service at [www.restek.com/Contact-Us](http://www.restek.com/Contact-Us).
- The packaged amount is the minimum sample size for which uncertainty is valid. The ampules are over-filled to ensure that the minimum packaged amount can be sufficiently transferred.

### Manufacturing Notes:

- Concentration is based upon gravimetric preparation using either a balance whose calibration has been verified daily using NIST traceable weights, and/or dilutions with Class A glassware.

### Handling Notes:

- Stability of the unopened product, when stored in compliance with the recommended conditions, is guaranteed through the expiration displayed on the product label and certificate. Contact Restek for additional opened product stability information, with the knowledge/understanding that open product stability is subject to the specific handling and environmental conditions to which the product is exposed. For your convenience Restek supplies deactivated vials with most standards packed in 2mL ampules. Larger volume deactivated vials are available through Restek as a custom ordered item. Additionally, Restek sells DMDCS for the purpose of glassware deactivation as catalog number 31861, which includes complete instructions.



Reagent

---

**8330\_PETN1000\_00010**





110 Benner Circle  
Bellefonte, PA 16823-8812  
Tel: (800)356-1688  
Fax: (814)353-1309

www.restek.com

# CERTIFIED REFERENCE MATERIAL

## Certificate of Analysis



### FOR LABORATORY USE ONLY-READ SDS PRIOR TO USE.

*This Reference Material is intended for Laboratory Use Only as a standard for the qualitative and/or quantitative determination of the analyte(s) listed.*

Catalog No. : 31600 Lot No.: A0187142  
Description : PETN Standard  
PETN Standard 1000µg/mL, Methanol, 1mL/ampul  
Container Size : 2 mL Pkg Amt: > 1 mL  
Expiration Date : July 31, 2027 Storage: 10°C or colder  
Handling: Sonicate prior to use. Ship: Ambient

### CERTIFIED VALUES

| Elution Order | Compound                                               | Grav. Conc. (weight/volume) | Expanded Uncertainty (95% C.L.; K=2)                                                       |
|---------------|--------------------------------------------------------|-----------------------------|--------------------------------------------------------------------------------------------|
| I             | PETN<br>CAS # 78-11-5<br>Purity 99%<br>(Lot 051108JLM) | 1,001.0 µg/mL               | +/- 5.9456 µg/mL Gravimetric<br>+/- 54.8378 µg/mL Unstressed<br>+/- 63.9463 µg/mL Stressed |

Solvent: Methanol  
CAS # 67-56-1  
Purity 99%

Reagent

---

**8330\_PETN1000\_00011**



110 Benner Circle  
Bellefonte, PA 16823-8812  
Tel: (800)356-1688  
Fax: (814)353-1309

www.restek.com

# CERTIFIED REFERENCE MATERIAL

## Certificate of Analysis



### FOR LABORATORY USE ONLY-READ SDS PRIOR TO USE.

*This Reference Material is intended for Laboratory Use Only as a standard for the qualitative and/or quantitative determination of the analyte(s) listed.*

**Catalog No. :** 31600 **Lot No.:** A0187142  
**Description :** PETN Standard  
PETN Standard 1000µg/mL, Methanol, 1mL/ampul  
**Container Size :** 2 mL **Pkg Amt:** > 1 mL  
**Expiration Date :** July 31, 2027 **Storage:** 10°C or colder  
**Handling:** Sonicate prior to use. **Ship:** Ambient

### CERTIFIED VALUES

| Elution Order | Compound                                               | Grav. Conc. (weight/volume) | Expanded Uncertainty (95% C.L.; K=2)                                                       |
|---------------|--------------------------------------------------------|-----------------------------|--------------------------------------------------------------------------------------------|
| I             | PETN<br>CAS # 78-11-5<br>Purity 99%<br>(Lot 051108JLM) | 1,001.0 µg/mL               | +/- 5.9456 µg/mL Gravimetric<br>+/- 54.8378 µg/mL Unstressed<br>+/- 63.9463 µg/mL Stressed |

**Solvent:** Methanol  
CAS # 67-56-1  
Purity 99%

Reagent

---

**833035DNASTk\_00052**

# CERTIFICATE OF ANALYSIS

**Catalog No:** M-8330-ADD-4  
**Description:** 3,5-Dinitroaniline  
**Lot:** 222011692-01

**Solvent:** Methanol (50%)  
Acetonitrile (50%)

**Hazards:** Refer to SDS for complete safety information



Signal Word: Danger

**Date Certified:** Sep 6, 2022

**Expiration:** Oct 6, 2023

**Sample Size:** 1 mL

**Components:** 1

**Storage Condition:** Ambient (>5 °C)

**Certified Reference Material**



| Component          | CAS #    | Purity % | Prepared Concentration <sup>2</sup> | Certified Analyte Concentration <sup>1</sup> |
|--------------------|----------|----------|-------------------------------------|----------------------------------------------|
|                    |          | (GC/MS)  | (µg/mL)                             | (µg/mL)                                      |
| 3,5-Dinitroaniline | 618-87-1 | 100.0    | 100.6                               | 100.6                                        |

This Certified Reference Material was verified in accordance with ISO/IEC 17025

A product with a suffix (-1A, -2B, etc. or -01, -02, etc.) on its lot number has had its expiration date extended and is identical to the same lot number without the suffix.

<sup>2</sup> All weights are traceable through NIST, Test No. 684/289871-17

<sup>1</sup> Certified Analyte Concentration = Purity x Prepared Concentration.

The Uncertainty associated with the certified concentration reported on this certificate is  $\pm 2.4\%$ . This value is the combined expanded uncertainty and represents an estimated standard deviation equal to the positive square root of the total variation of the uncertainty of components. A normal distribution is assumed and a coverage factor of K=2 is chosen using approximately a 95% confidence level.

Labels and certificates follow U.S. Conventions in reporting numerical values: A comma (,) is used to separate units of one-thousand or greater. A period (.) is used as a decimal place marker.

The information on this certificate may not be reproduced without the express permission of the manufacturer. See reverse side for additional information

Hazard Information: Please refer to the SDS for information regarding the hazards associated with using this material.

This product was prepared according to in-house procedures and is guaranteed to be homogeneous.

Certified By:



Larry Decker, Organic QC Manager

Reagent

---

**833035DNASTk\_00053**

# CERTIFICATE OF ANALYSIS

**Catalog No:** M-8330-ADD-4  
**Description:** 3,5-Dinitroaniline  
**Lot:** 222011692-01

**Solvent:** Methanol (50%)  
Acetonitrile (50%)

**Hazards:** Refer to SDS for complete safety information



Signal Word: Danger

**Date Certified:** Sep 6, 2022

**Expiration:** Oct 6, 2023

**Sample Size:** 1 mL

**Components:** 1

**Storage Condition:** Ambient (>5 °C)

**Certified Reference Material**



| Component          | CAS #    | Purity % | Prepared Concentration <sup>2</sup> | Certified Analyte Concentration <sup>1</sup> |
|--------------------|----------|----------|-------------------------------------|----------------------------------------------|
|                    |          | (GC/MS)  | (µg/mL)                             | (µg/mL)                                      |
| 3,5-Dinitroaniline | 618-87-1 | 100.0    | 100.6                               | 100.6                                        |

This Certified Reference Material was verified in accordance with ISO/IEC 17025

A product with a suffix (-1A, -2B, etc. or -01, -02, etc.) on its lot number has had its expiration date extended and is identical to the same lot number without the suffix.

<sup>2</sup> All weights are traceable through NIST, Test No. 684/289871-17

<sup>1</sup> Certified Analyte Concentration = Purity x Prepared Concentration.

The Uncertainty associated with the certified concentration reported on this certificate is  $\pm 2.4\%$ . This value is the combined expanded uncertainty and represents an estimated standard deviation equal to the positive square root of the total variation of the uncertainty of components. A normal distribution is assumed and a coverage factor of K=2 is chosen using approximately a 95% confidence level.

Labels and certificates follow U.S. Conventions in reporting numerical values: A comma (,) is used to separate units of one-thousand or greater. A period (.) is used as a decimal place marker.

The information on this certificate may not be reproduced without the express permission of the manufacturer. See reverse side for additional information

Hazard Information: Please refer to the SDS for information regarding the hazards associated with using this material.

This product was prepared according to in-house procedures and is guaranteed to be homogeneous.

Certified By:



Larry Decker, Organic QC Manager

Reagent

---

**833035DNASTk\_00054**



# CERTIFICATE OF ANALYSIS

**Catalog No:** M-8330-ADD-4  
**Description:** 3,5-Dinitroaniline  
**Lot:** 222011692-01

**Solvent:** Methanol (50%)  
Acetonitrile (50%)

**Hazards:** Refer to SDS for complete safety information



Signal Word: Danger

**Date Certified:** Sep 6, 2022

**Expiration:** Oct 6, 2023

**Sample Size:** 1 mL

**Components:** 1

**Storage Condition:** Ambient (>5 °C)

**Certified Reference Material**



| Component          | CAS #    | Purity % | Prepared Concentration <sup>2</sup> | Certified Analyte Concentration <sup>1</sup> |
|--------------------|----------|----------|-------------------------------------|----------------------------------------------|
|                    |          | (GC/MS)  | (µg/mL)                             | (µg/mL)                                      |
| 3,5-Dinitroaniline | 618-87-1 | 100.0    | 100.6                               | 100.6                                        |

This Certified Reference Material was verified in accordance with ISO/IEC 17025

A product with a suffix (-1A, -2B, etc. or -01, -02, etc.) on its lot number has had its expiration date extended and is identical to the same lot number without the suffix.

<sup>2</sup> All weights are traceable through NIST, Test No. 684/289871-17

<sup>1</sup> Certified Analyte Concentration = Purity x Prepared Concentration.

The Uncertainty associated with the certified concentration reported on this certificate is  $\pm 2.4\%$ . This value is the combined expanded uncertainty and represents an estimated standard deviation equal to the positive square root of the total variation of the uncertainty of components. A normal distribution is assumed and a coverage factor of K=2 is chosen using approximately a 95% confidence level.

Labels and certificates follow U.S. Conventions in reporting numerical values: A comma (,) is used to separate units of one-thousand or greater. A period (.) is used as a decimal place marker.

The information on this certificate may not be reproduced without the express permission of the manufacturer. See reverse side for additional information

Hazard Information: Please refer to the SDS for information regarding the hazards associated with using this material.

This product was prepared according to in-house procedures and is guaranteed to be homogeneous.

Certified By:

Larry Decker, Organic QC Manager

Reagent

---

**8330DiaminLCS\_00052**

**Preliminary Report**

Eurofins Denver

LCS, Lab Control Sample Report

Sample Path: \\chromfs\Denver\ChromData\CHHPLC\_X\20221221-117211.b\12210014.D

Lims ID: 8330Dia LCS52

Inj. Date: 21-Dec-2022 13:28:20

Worklist ID: 280-0117211-014

Instrument: CHHPLC\_X3

Method: 8330\_X3

| Compound                  | Amount Added | Amount Recovered | %Rec  | Limits 1<br>0B_Sonc_ | Limits 2<br>3535 |
|---------------------------|--------------|------------------|-------|----------------------|------------------|
| 2 2,6-diamino-4-nitrotolu | 0.5000       | 0.5110           | 102.2 | 10-150               |                  |
| 5 2,4-diamino-6-nitrotolu | 0.5000       | 0.5040           | 100.8 | 10-150               |                  |

Samples for Limit Group: 1, Lims Prep Method: 8330B\_Sonc\_10g

280-170073-A-48-B

280-170073-A-43-B

Samples for Limit Group: 2, Lims Prep Method: 3535

280-170078-A-20-C

280-170078-A-21-C

280-170078-A-22-C

280-170078-A-23-C

280-170078-A-24-C

280-170078-A-25-B

280-170078-A-26-B

280-170078-A-27-F

280-170078-A-28-B

280-170078-A-29-B

280-170078-A-30-B

280-170078-A-31-B

280-170078-A-32-B

280-170078-A-33-B

280-170078-E-34-B

Reagent

---

**8330LCsMix1\_00136**



CERTIFIED REFERENCE MATERIAL

110 Benner Circle  
Bellefonte, PA 16823-8812  
Tel: (800)356-1688  
Fax: (814)353-1309

www.restek.com

# Certificate of Analysis



FOR LABORATORY USE ONLY-READ SDS PRIOR TO USE.

*This Reference Material is intended for Laboratory Use Only as a standard for the qualitative and/or quantitative determination of the analyte(s) listed.*

**Catalog No. :** 31450 **Lot No.:** A0183848

**Description :** 8330 Calibration Mix #1  
8330 Calibration Std #1 1000µg/mL, Acetonitrile, 1mL/ampul

**Container Size :** 2 mL **Pkg Amt:** > 1 mL

**Expiration Date :** April 30, 2027 **Storage:** 10°C or colder  
**Ship:** Ambient

## CERTIFIED VALUES

| Elution Order | Compound              |                           | Grav. Conc.<br>(weight/volume) | Expanded Uncertainty<br>(95% C.L.; K=2) |         |                   |
|---------------|-----------------------|---------------------------|--------------------------------|-----------------------------------------|---------|-------------------|
| 1             | HMX                   |                           | 1,008.0 µg/mL                  | +/-                                     | 5.9872  | µg/mL Gravimetric |
|               | CAS #                 | 2691-41-0 (Lot 210324JLM) |                                | +/-                                     | 55.2213 | µg/mL Unstressed  |
|               | Purity                | 99%                       |                                | +/-                                     | 64.3934 | µg/mL Stressed    |
|               |                       |                           |                                |                                         |         |                   |
| 2             | RDX                   |                           | 1,007.0 µg/mL                  | +/-                                     | 5.9813  | µg/mL Gravimetric |
|               | CAS #                 | 121-82-4 (Lot 080228JLM)  |                                | +/-                                     | 55.1665 | µg/mL Unstressed  |
|               | Purity                | 99%                       |                                | +/-                                     | 64.3295 | µg/mL Stressed    |
|               |                       |                           |                                |                                         |         |                   |
| 3             | 1,3,5-Trinitrobenzene |                           | 1,008.0 µg/mL                  | +/-                                     | 5.9872  | µg/mL Gravimetric |
|               | CAS #                 | 99-35-4 (Lot A6TDK)       |                                | +/-                                     | 55.2213 | µg/mL Unstressed  |
|               | Purity                | 99%                       |                                | +/-                                     | 64.3934 | µg/mL Stressed    |
|               |                       |                           |                                |                                         |         |                   |
| 4             | 1,3-Dinitrobenzene    |                           | 1,006.0 µg/mL                  | +/-                                     | 5.9753  | µg/mL Gravimetric |
|               | CAS #                 | 99-65-0 (Lot 1-DXX-24-1)  |                                | +/-                                     | 55.1117 | µg/mL Unstressed  |
|               | Purity                | 99%                       |                                | +/-                                     | 64.2657 | µg/mL Stressed    |
|               |                       |                           |                                |                                         |         |                   |
| 5             | Nitrobenzene          |                           | 1,007.0 µg/mL                  | +/-                                     | 5.9813  | µg/mL Gravimetric |
|               | CAS #                 | 98-95-3 (Lot 10224044)    |                                | +/-                                     | 55.1665 | µg/mL Unstressed  |
|               | Purity                | 99%                       |                                | +/-                                     | 64.3295 | µg/mL Stressed    |
|               |                       |                           |                                |                                         |         |                   |
| 6             | 2,4,6-Trinitrotoluene |                           | 1,002.0 µg/mL                  | +/-                                     | 5.9516  | µg/mL Gravimetric |
|               | CAS #                 | 118-96-7 (Lot D11836200)  |                                | +/-                                     | 54.8926 | µg/mL Unstressed  |
|               | Purity                | 99%                       |                                | +/-                                     | 64.0101 | µg/mL Stressed    |
|               |                       |                           |                                |                                         |         |                   |
| 7             | 2,4-Dinitrotoluene    |                           | 1,005.0 µg/mL                  | +/-                                     | 5.9694  | µg/mL Gravimetric |
|               | CAS #                 | 121-14-2 (Lot MKAA0690V)  |                                | +/-                                     | 55.0569 | µg/mL Unstressed  |
|               | Purity                | 99%                       |                                | +/-                                     | 64.2018 | µg/mL Stressed    |
|               |                       |                           |                                |                                         |         |                   |

**Solvent:** Acetonitrile  
**CAS #** 75-05-8  
**Purity** 99%

**Column:**  
250mm x 4.6mm  
Ultra C18 (cat.# 9174575)


**Flow Rate:**  
1.0 ml/min.

**Mobile Phase A:**  
water:methanol (44:56 V/V)


**Mobile Phase B:**

**Mobile Phase Composition:**  
100%A

**Det. Type:**  
Wavelength: 210nm & 254nm



This chromatogram represents a general set of testing conditions chosen for product acceptance. For optimal results in your lab, conditions should be adjusted for your specific instrument, method, and application.

  
Tom Suckar - Mix Technician

**Date Mixed:** 08-Apr-2022

**Balance:** B251644995

  
Fang-Yan Lo - QC Analyst

**Date Passed:** 13-Apr-2022

Manufactured under Restek's ISO 9001:2015  
Registered Quality System  
Certificate #FM 80397

## General Certified Reference Material Notes

### Expiration Notes:

- Expiration date valid for unopened ampul stored in compliance with the recommended conditions.
- Uncertainty, concentration, and expiration of the CRM are based on the unopened product being stored according to the recommended condition found in the storage field.

### Purity Notes:

- Purity and/or chemical identity are determined by one or more of the following techniques: GC/FID, HPLC, GC/μECD, GC/MS, LC/MS, RI, and/or melting point.
- Compounds with a listed purity of less than 99% have been weight corrected to compensate for impurities and/or salts. A correction factor is used to calculate the amount of compound necessary to achieve the desired concentration of the parent compound in solution.
- Purity of isomeric compounds is reported as the sum of the isomers.
- Purity values are rounded to the nearest whole number.

### Certified Uncertainty Value Notes:

- The uncertainties are determined in accordance with ISO 17034 and Guide 35. The certified combined stressed uncertainty value ( includes gravimetric uncertainty, homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty and were combined using the following formula:

$$U_{combined\ stressed} = k \sqrt{U_{gravimetric}^2 + U_{homogeneity}^2 + U_{storage\ stability}^2 + U_{shipping\ stability}^2}$$

$k$  is a coverage factor of 2, which gives a level of confidence of approximately 95%.

- It is important to note that the shipping stability uncertainty was obtained under temperature extremes for specific time intervals; therefore, the certified combined stressed uncertainty value should only be applied to the product if it was stored at non-standard temperature conditions up to and including 7 days. Contact Restek Technical Service at [www.restek.com/Contact-Us](http://www.restek.com/Contact-Us) for use recommendations if your shipment was in-transit for more than 7 days at non-standard temperature conditions.
- Apply the certified combined unstressed uncertainty value if the product was received under standard shipping conditions. Apply the certified combined stressed uncertainty value if the product was received under non-standard conditions as specified below.

| Label Conditions                                          | Standard Conditions | Non-Standard Conditions |
|-----------------------------------------------------------|---------------------|-------------------------|
| 25°C Nominal (Room Temperature)                           | < 60°C              | ≥ 60°C up to 7 days     |
| 10°C or colder (Refrigerate)                              | < 40°C              | ≥ 40°C up to 7 days     |
| 0°C or colder (Freezer)<br>-20°C or colder (Deep Freezer) | < 25°C              | ≥ 25°C up to 7 days     |

- Separate (not combined) uncertainty values for gravimetric uncertainty are also displayed on the certificate, if needed, separate homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty values are available by contacting Restek Technical Service at [www.restek.com/Contact-Us](http://www.restek.com/Contact-Us).
- The packaged amount is the minimum sample size for which uncertainty is valid. The ampules are over-filled to ensure that the minimum packaged amount can be sufficiently transferred.

### Manufacturing Notes:

- Concentration is based upon gravimetric preparation using either a balance whose calibration has been verified daily using NIST traceable weights, and/or dilutions with Class A glassware.

### Handling Notes:

- Stability of the unopened product, when stored in compliance with the recommended conditions, is guaranteed through the expiration displayed on the product label and certificate. Contact Restek for additional opened product stability information, with the knowledge/understanding that open product stability is subject to the specific handling and environmental conditions to which the product is exposed. For your convenience Restek supplies deactivated vials with most standards packed in 2mL ampules. Larger volume deactivated vials are available through Restek as a custom ordered item. Additionally, Restek sells DMDCS for the purpose of glassware deactivation as catalog number 31861, which includes complete instructions.





Reagent

---

**8330LCSMix1\_00140**



CERTIFIED REFERENCE MATERIAL

110 Benner Circle  
Bellefonte, PA 16823-8812  
Tel: (800)356-1688  
Fax: (814)353-1309

www.restek.com

# Certificate of Analysis



FOR LABORATORY USE ONLY-READ SDS PRIOR TO USE.

*This Reference Material is intended for Laboratory Use Only as a standard for the qualitative and/or quantitative determination of the analyte(s) listed.*

**Catalog No. :** 31450 **Lot No.:** A0183848

**Description :** 8330 Calibration Mix #1  
8330 Calibration Std #1 1000µg/mL, Acetonitrile, 1mL/ampul

**Container Size :** 2 mL **Pkg Amt:** > 1 mL

**Expiration Date :** April 30, 2027 **Storage:** 10°C or colder

**Ship:** Ambient

## CERTIFIED VALUES

| Elution Order | Compound              |           |                  | Grav. Conc. (weight/volume) | Expanded Uncertainty (95% C.L.; K=2) |       |             |
|---------------|-----------------------|-----------|------------------|-----------------------------|--------------------------------------|-------|-------------|
| 1             | HMX                   |           |                  | 1,008.0 µg/mL               | +/- 5.9872                           | µg/mL | Gravimetric |
|               | CAS #                 | 2691-41-0 | (Lot 210324JLM)  |                             | +/- 55.2213                          | µg/mL | Unstressed  |
|               | Purity                | 99%       |                  |                             | +/- 64.3934                          | µg/mL | Stressed    |
| 2             | RDX                   |           |                  | 1,007.0 µg/mL               | +/- 5.9813                           | µg/mL | Gravimetric |
|               | CAS #                 | 121-82-4  | (Lot 080228JLM)  |                             | +/- 55.1665                          | µg/mL | Unstressed  |
|               | Purity                | 99%       |                  |                             | +/- 64.3295                          | µg/mL | Stressed    |
| 3             | 1,3,5-Trinitrobenzene |           |                  | 1,008.0 µg/mL               | +/- 5.9872                           | µg/mL | Gravimetric |
|               | CAS #                 | 99-35-4   | (Lot A6TDK)      |                             | +/- 55.2213                          | µg/mL | Unstressed  |
|               | Purity                | 99%       |                  |                             | +/- 64.3934                          | µg/mL | Stressed    |
| 4             | 1,3-Dinitrobenzene    |           |                  | 1,006.0 µg/mL               | +/- 5.9753                           | µg/mL | Gravimetric |
|               | CAS #                 | 99-65-0   | (Lot 1-DXX-24-1) |                             | +/- 55.1117                          | µg/mL | Unstressed  |
|               | Purity                | 99%       |                  |                             | +/- 64.2657                          | µg/mL | Stressed    |
| 5             | Nitrobenzene          |           |                  | 1,007.0 µg/mL               | +/- 5.9813                           | µg/mL | Gravimetric |
|               | CAS #                 | 98-95-3   | (Lot 10224044)   |                             | +/- 55.1665                          | µg/mL | Unstressed  |
|               | Purity                | 99%       |                  |                             | +/- 64.3295                          | µg/mL | Stressed    |
| 6             | 2,4,6-Trinitrotoluene |           |                  | 1,002.0 µg/mL               | +/- 5.9516                           | µg/mL | Gravimetric |
|               | CAS #                 | 118-96-7  | (Lot D11836200)  |                             | +/- 54.8926                          | µg/mL | Unstressed  |
|               | Purity                | 99%       |                  |                             | +/- 64.0101                          | µg/mL | Stressed    |
| 7             | 2,4-Dinitrotoluene    |           |                  | 1,005.0 µg/mL               | +/- 5.9694                           | µg/mL | Gravimetric |
|               | CAS #                 | 121-14-2  | (Lot MKAA0690V)  |                             | +/- 55.0569                          | µg/mL | Unstressed  |
|               | Purity                | 99%       |                  |                             | +/- 64.2018                          | µg/mL | Stressed    |

Solvent: Acetonitrile  
CAS # 75-05-8  
Purity 99%

Column:  
250mm x 4.6mm  
Ultra C18 (cat.# 9174575)


Flow Rate:  
1.0 ml/min.

Mobile Phase A:  
water:methanol (44:56 V/V)

Mobile Phase B:

Mobile Phase Composition:  
100%A

Det. Type:  
Wavelength: 210nm & 254nm



This chromatogram represents a general set of testing conditions chosen for product acceptance. For optimal results in your lab, conditions should be adjusted for your specific instrument, method, and application.

  
Tom Suckar - Mix Technician

Date Mixed: 08-Apr-2022      Balance: B251644995

  
Fang-Yun Lo - QC Analyst

Date Passed: 13-Apr-2022

Manufactured under Restek's ISO 9001:2015  
Registered Quality System  
Certificate #FM 80397

## General Certified Reference Material Notes

### Expiration Notes:

- Expiration date valid for unopened ampul stored in compliance with the recommended conditions.
- Uncertainty, concentration, and expiration of the CRM are based on the unopened product being stored according to the recommended condition found in the storage field.

### Purity Notes:

- Purity and/or chemical identity are determined by one or more of the following techniques: GC/FID, HPLC, GC/μECD, GC/MS, LC/MS, RI, and/or melting point.
- Compounds with a listed purity of less than 99% have been weight corrected to compensate for impurities and/or salts. A correction factor is used to calculate the amount of compound necessary to achieve the desired concentration of the parent compound in solution.
- Purity of isomeric compounds is reported as the sum of the isomers.
- Purity values are rounded to the nearest whole number.

### Certified Uncertainty Value Notes:

- The uncertainties are determined in accordance with ISO 17034 and Guide 35. The certified combined stressed uncertainty value ( includes gravimetric uncertainty, homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty and were combined using the following formula:

$$U_{combined\ stressed} = k \sqrt{U_{gravimetric}^2 + U_{homogeneity}^2 + U_{storage\ stability}^2 + U_{shipping\ stability}^2}$$

$k$  is a coverage factor of 2, which gives a level of confidence of approximately 95%.

- It is important to note that the shipping stability uncertainty was obtained under temperature extremes for specific time intervals; therefore, the certified combined stressed uncertainty value should only be applied to the product if it was stored at non-standard temperature conditions up to and including 7 days. Contact Restek Technical Service at [www.restek.com/Contact-Us](http://www.restek.com/Contact-Us) for use recommendations if your shipment was in-transit for more than 7 days at non-standard temperature conditions.
- Apply the certified combined unstressed uncertainty value if the product was received under standard shipping conditions. Apply the certified combined stressed uncertainty value if the product was received under non-standard conditions as specified below.

| Label Conditions                                          | Standard Conditions | Non-Standard Conditions |
|-----------------------------------------------------------|---------------------|-------------------------|
| 25°C Nominal (Room Temperature)                           | < 60°C              | ≥ 60°C up to 7 days     |
| 10°C or colder (Refrigerate)                              | < 40°C              | ≥ 40°C up to 7 days     |
| 0°C or colder (Freezer)<br>-20°C or colder (Deep Freezer) | < 25°C              | ≥ 25°C up to 7 days     |

- Separate (not combined) uncertainty values for gravimetric uncertainty are also displayed on the certificate, if needed, separate homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty values are available by contacting Restek Technical Service at [www.restek.com/Contact-Us](http://www.restek.com/Contact-Us).
- The packaged amount is the minimum sample size for which uncertainty is valid. The ampules are over-filled to ensure that the minimum packaged amount can be sufficiently transferred.

### Manufacturing Notes:

- Concentration is based upon gravimetric preparation using either a balance whose calibration has been verified daily using NIST traceable weights, and/or dilutions with Class A glassware.

### Handling Notes:

- Stability of the unopened product, when stored in compliance with the recommended conditions, is guaranteed through the expiration displayed on the product label and certificate. Contact Restek for additional opened product stability information, with the knowledge/understanding that open product stability is subject to the specific handling and environmental conditions to which the product is exposed. For your convenience Restek supplies deactivated vials with most standards packed in 2mL ampules. Larger volume deactivated vials are available through Restek as a custom ordered item. Additionally, Restek sells DMDCS for the purpose of glassware deactivation as catalog number 31861, which includes complete instructions.

Reagent

---

**8330LCSMix1\_00141**



CERTIFIED REFERENCE MATERIAL

110 Benner Circle  
Bellefonte, PA 16823-8812  
Tel: (800)356-1688  
Fax: (814)353-1309

www.restek.com

# Certificate of Analysis



FOR LABORATORY USE ONLY-READ SDS PRIOR TO USE.

*This Reference Material is intended for Laboratory Use Only as a standard for the qualitative and/or quantitative determination of the analyte(s) listed.*

**Catalog No. :** 31450 **Lot No.:** A0183848

**Description :** 8330 Calibration Mix #1  
8330 Calibration Std #1 1000µg/mL, Acetonitrile, 1mL/ampul

**Container Size :** 2 mL **Pkg Amt:** > 1 mL

**Expiration Date :** April 30, 2027 **Storage:** 10°C or colder

**Ship:** Ambient

## CERTIFIED VALUES

| Elution Order | Compound              |           |                  | Grav. Conc. (weight/volume) | Expanded Uncertainty (95% C.L.; K=2) |       |             |
|---------------|-----------------------|-----------|------------------|-----------------------------|--------------------------------------|-------|-------------|
| 1             | HMX                   |           |                  | 1,008.0 µg/mL               | +/- 5.9872                           | µg/mL | Gravimetric |
|               | CAS #                 | 2691-41-0 | (Lot 210324JLM)  |                             | +/- 55.2213                          | µg/mL | Unstressed  |
|               | Purity                | 99%       |                  |                             | +/- 64.3934                          | µg/mL | Stressed    |
| 2             | RDX                   |           |                  | 1,007.0 µg/mL               | +/- 5.9813                           | µg/mL | Gravimetric |
|               | CAS #                 | 121-82-4  | (Lot 080228JLM)  |                             | +/- 55.1665                          | µg/mL | Unstressed  |
|               | Purity                | 99%       |                  |                             | +/- 64.3295                          | µg/mL | Stressed    |
| 3             | 1,3,5-Trinitrobenzene |           |                  | 1,008.0 µg/mL               | +/- 5.9872                           | µg/mL | Gravimetric |
|               | CAS #                 | 99-35-4   | (Lot A6TDK)      |                             | +/- 55.2213                          | µg/mL | Unstressed  |
|               | Purity                | 99%       |                  |                             | +/- 64.3934                          | µg/mL | Stressed    |
| 4             | 1,3-Dinitrobenzene    |           |                  | 1,006.0 µg/mL               | +/- 5.9753                           | µg/mL | Gravimetric |
|               | CAS #                 | 99-65-0   | (Lot 1-DXX-24-1) |                             | +/- 55.1117                          | µg/mL | Unstressed  |
|               | Purity                | 99%       |                  |                             | +/- 64.2657                          | µg/mL | Stressed    |
| 5             | Nitrobenzene          |           |                  | 1,007.0 µg/mL               | +/- 5.9813                           | µg/mL | Gravimetric |
|               | CAS #                 | 98-95-3   | (Lot 10224044)   |                             | +/- 55.1665                          | µg/mL | Unstressed  |
|               | Purity                | 99%       |                  |                             | +/- 64.3295                          | µg/mL | Stressed    |
| 6             | 2,4,6-Trinitrotoluene |           |                  | 1,002.0 µg/mL               | +/- 5.9516                           | µg/mL | Gravimetric |
|               | CAS #                 | 118-96-7  | (Lot D11836200)  |                             | +/- 54.8926                          | µg/mL | Unstressed  |
|               | Purity                | 99%       |                  |                             | +/- 64.0101                          | µg/mL | Stressed    |
| 7             | 2,4-Dinitrotoluene    |           |                  | 1,005.0 µg/mL               | +/- 5.9694                           | µg/mL | Gravimetric |
|               | CAS #                 | 121-14-2  | (Lot MKAA0690V)  |                             | +/- 55.0569                          | µg/mL | Unstressed  |
|               | Purity                | 99%       |                  |                             | +/- 64.2018                          | µg/mL | Stressed    |

Solvent: Acetonitrile  
CAS # 75-05-8  
Purity 99%

Column:  
250mm x 4.6mm  
Ultra C18 (cat.# 9174575)


Flow Rate:  
1.0 ml/min.

Mobile Phase A:  
water:methanol (44:56 V/V)

Mobile Phase B:

Mobile Phase Composition:  
100%A

Det. Type:  
Wavelength: 210nm & 254nm



This chromatogram represents a general set of testing conditions chosen for product acceptance. For optimal results in your lab, conditions should be adjusted for your specific instrument, method, and application.

  
Tom Suckar - Mix Technician

Date Mixed: 08-Apr-2022      Balance: B251644995

  
Fang-Yun Lo - QC Analyst

Date Passed: 13-Apr-2022

Manufactured under Restek's ISO 9001:2015  
Registered Quality System  
Certificate #FM 80397

## General Certified Reference Material Notes

### Expiration Notes:

- Expiration date valid for unopened ampul stored in compliance with the recommended conditions.
- Uncertainty, concentration, and expiration of the CRM are based on the unopened product being stored according to the recommended condition found in the storage field.

### Purity Notes:

- Purity and/or chemical identity are determined by one or more of the following techniques: GC/FID, HPLC, GC/μECD, GC/MS, LC/MS, RI, and/or melting point.
- Compounds with a listed purity of less than 99% have been weight corrected to compensate for impurities and/or salts. A correction factor is used to calculate the amount of compound necessary to achieve the desired concentration of the parent compound in solution.
- Purity of isomeric compounds is reported as the sum of the isomers.
- Purity values are rounded to the nearest whole number.

### Certified Uncertainty Value Notes:

- The uncertainties are determined in accordance with ISO 17034 and Guide 35. The certified combined stressed uncertainty value ( includes gravimetric uncertainty, homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty and were combined using the following formula:

$$U_{combined\ stressed} = k \sqrt{U_{gravimetric}^2 + U_{homogeneity}^2 + U_{storage\ stability}^2 + U_{shipping\ stability}^2}$$

$k$  is a coverage factor of 2, which gives a level of confidence of approximately 95%.

- It is important to note that the shipping stability uncertainty was obtained under temperature extremes for specific time intervals; therefore, the certified combined stressed uncertainty value should only be applied to the product if it was stored at non-standard temperature conditions up to and including 7 days. Contact Restek Technical Service at [www.restek.com/Contact-Us](http://www.restek.com/Contact-Us) for use recommendations if your shipment was in-transit for more than 7 days at non-standard temperature conditions.
- Apply the certified combined unstressed uncertainty value if the product was received under standard shipping conditions. Apply the certified combined stressed uncertainty value if the product was received under non-standard conditions as specified below.

| Label Conditions                                          | Standard Conditions | Non-Standard Conditions |
|-----------------------------------------------------------|---------------------|-------------------------|
| 25°C Nominal (Room Temperature)                           | < 60°C              | ≥ 60°C up to 7 days     |
| 10°C or colder (Refrigerate)                              | < 40°C              | ≥ 40°C up to 7 days     |
| 0°C or colder (Freezer)<br>-20°C or colder (Deep Freezer) | < 25°C              | ≥ 25°C up to 7 days     |

- Separate (not combined) uncertainty values for gravimetric uncertainty are also displayed on the certificate, if needed, separate homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty values are available by contacting Restek Technical Service at [www.restek.com/Contact-Us](http://www.restek.com/Contact-Us).
- The packaged amount is the minimum sample size for which uncertainty is valid. The ampules are over-filled to ensure that the minimum packaged amount can be sufficiently transferred.

### Manufacturing Notes:

- Concentration is based upon gravimetric preparation using either a balance whose calibration has been verified daily using NIST traceable weights, and/or dilutions with Class A glassware.

### Handling Notes:

- Stability of the unopened product, when stored in compliance with the recommended conditions, is guaranteed through the expiration displayed on the product label and certificate. Contact Restek for additional opened product stability information, with the knowledge/understanding that open product stability is subject to the specific handling and environmental conditions to which the product is exposed. For your convenience Restek supplies deactivated vials with most standards packed in 2mL ampules. Larger volume deactivated vials are available through Restek as a custom ordered item. Additionally, Restek sells DMDCS for the purpose of glassware deactivation as catalog number 31861, which includes complete instructions.



Reagent

---

**8330LCSmix2\_00034**



110 Benner Circle  
Bellefonte, PA 16823-8812  
Tel: (800)356-1688  
Fax: (814)353-1309

www.restek.com

# CERTIFIED REFERENCE MATERIAL

## Certificate of Analysis



### FOR LABORATORY USE ONLY-READ SDS PRIOR TO USE.

*This Reference Material is intended for Laboratory Use Only as a standard for the qualitative and/or quantitative determination of the analyte(s) listed.*

**Catalog No.:** 31451 **Lot No.:** A0186475

**Description:** 8330 Calibration Mix #2  
8330 Calibration Std #2 1000µg/mL, Acetonitrile, 1mL/ampul

**Container Size:** 2 mL **Pkg Amt:** > 1 mL

**Expiration Date:** June 30, 2027 **Storage:** 10°C or colder

**Ship:** Ambient

### CERTIFIED VALUES

| Elution Order | Compound                                                                       | Grav. Conc. (weight/volume) | Expanded Uncertainty (95% C.L.; K=2)                       |                                       |
|---------------|--------------------------------------------------------------------------------|-----------------------------|------------------------------------------------------------|---------------------------------------|
| 1             | Tetryl<br>CAS # 479-45-8 (Lot 211028JLM)<br>Purity 99%                         | 1,004.0 µg/mL               | +/- 5.9635 µg/mL<br>+/- 55.0021 µg/mL<br>+/- 64.1379 µg/mL | Gravimetric<br>Unstressed<br>Stressed |
| 2             | 4-Amino-2,6-dinitrotoluene<br>CAS # 19406-51-0 (Lot ER070908-01)<br>Purity 99% | 1,002.0 µg/mL               | +/- 5.9516 µg/mL<br>+/- 54.8926 µg/mL<br>+/- 64.0101 µg/mL | Gravimetric<br>Unstressed<br>Stressed |
| 3             | 2-Amino-4,6-dinitrotoluene<br>CAS # 35572-78-2 (Lot 29550-55)<br>Purity 99%    | 1,002.0 µg/mL               | +/- 5.9516 µg/mL<br>+/- 54.8926 µg/mL<br>+/- 64.0101 µg/mL | Gravimetric<br>Unstressed<br>Stressed |
| 4             | 2,6-Dinitrotoluene<br>CAS # 606-20-2 (Lot BCBB8606)<br>Purity 99%              | 1,004.0 µg/mL               | +/- 5.9635 µg/mL<br>+/- 55.0021 µg/mL<br>+/- 64.1379 µg/mL | Gravimetric<br>Unstressed<br>Stressed |
| 5             | 2-Nitrotoluene<br>CAS # 88-72-2 (Lot BCBZ7826)<br>Purity 99%                   | 1,002.0 µg/mL               | +/- 5.9516 µg/mL<br>+/- 54.8926 µg/mL<br>+/- 64.0101 µg/mL | Gravimetric<br>Unstressed<br>Stressed |
| 6             | 4-Nitrotoluene<br>CAS # 99-99-0 (Lot BCCB0171)<br>Purity 99%                   | 1,002.0 µg/mL               | +/- 5.9516 µg/mL<br>+/- 54.8926 µg/mL<br>+/- 64.0101 µg/mL | Gravimetric<br>Unstressed<br>Stressed |
| 7             | 3-Nitrotoluene<br>CAS # 99-08-1 (Lot FBO01)<br>Purity 99%                      | 1,000.0 µg/mL               | +/- 5.9397 µg/mL<br>+/- 54.7830 µg/mL<br>+/- 63.8824 µg/mL | Gravimetric<br>Unstressed<br>Stressed |

**Solvent:** Acetonitrile  
**CAS #** 75-05-8  
**Purity** 99%

**Column:**  
250mm x 4.6mm  
Ultra C18 (cat.# 9174575)


**Flow Rate:**  
1.0 ml/min.

**Mobile Phase A:**  
water:methanol (44:56 V/V)

**Mobile Phase B:**

**Mobile Phase Composition:**  
100%A

**Det. Type:**  
Wavelength: 210nm & 254nm



This chromatogram represents a general set of testing conditions chosen for product acceptance. For optimal results in your lab, conditions should be adjusted for your specific instrument, method, and application.

  
Ashley Frantz - Quoting Technician

**Date Mixed:** 21-Jun-2022 **Balance:** 1128360905

  
Christie Mills - Operations Tech II - ARM QC

**Date Passed:** 29-Jun-2022

Manufactured under Restek's ISO 9001:2015  
Registered Quality System  
Certificate #FM 80397

## General Certified Reference Material Notes

### Expiration Notes:

- Expiration date valid for unopened ampul stored in compliance with the recommended conditions.
- Uncertainty, concentration, and expiration of the CRM are based on the unopened product being stored according to the recommended condition found in the storage field.

### Purity Notes:

- Purity and/or chemical identity are determined by one or more of the following techniques: GC/FID, HPLC, GC/μECD, GC/MS, LC/MS, RI, and/or melting point.
- Compounds with a listed purity of less than 99% have been weight corrected to compensate for impurities and/or salts. A correction factor is used to calculate the amount of compound necessary to achieve the desired concentration of the parent compound in solution.
- Purity of isomeric compounds is reported as the sum of the isomers.
- Purity values are rounded to the nearest whole number.

### Certified Uncertainty Value Notes:

- The uncertainties are determined in accordance with ISO 17034 and Guide 35. The certified combined stressed uncertainty value ( includes gravimetric uncertainty, homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty and were combined using the following formula:

$$U_{combined\ stressed} = k \sqrt{U_{gravimetric}^2 + U_{homogeneity}^2 + U_{storage\ stability}^2 + U_{shipping\ stability}^2}$$

$k$  is a coverage factor of 2, which gives a level of confidence of approximately 95%.

- It is important to note that the shipping stability uncertainty was obtained under temperature extremes for specific time intervals; therefore, the certified combined stressed uncertainty value should only be applied to the product if it was stored at non-standard temperature conditions up to and including 7 days. Contact Restek Technical Service at [www.restek.com/Contact-Us](http://www.restek.com/Contact-Us) for use recommendations if your shipment was in-transit for more than 7 days at non-standard temperature conditions.
- Apply the certified combined unstressed uncertainty value if the product was received under standard shipping conditions. Apply the certified combined stressed uncertainty value if the product was received under non-standard conditions as specified below.

| Label Conditions                                          | Standard Conditions | Non-Standard Conditions |
|-----------------------------------------------------------|---------------------|-------------------------|
| 25°C Nominal (Room Temperature)                           | < 60°C              | ≥ 60°C up to 7 days     |
| 10°C or colder (Refrigerate)                              | < 40°C              | ≥ 40°C up to 7 days     |
| 0°C or colder (Freezer)<br>-20°C or colder (Deep Freezer) | < 25°C              | ≥ 25°C up to 7 days     |

- Separate (not combined) uncertainty values for gravimetric uncertainty are also displayed on the certificate, if needed, separate homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty values are available by contacting Restek Technical Service at [www.restek.com/Contact-Us](http://www.restek.com/Contact-Us).
- The packaged amount is the minimum sample size for which uncertainty is valid. The ampules are over-filled to ensure that the minimum packaged amount can be sufficiently transferred.

### Manufacturing Notes:

- Concentration is based upon gravimetric preparation using either a balance whose calibration has been verified daily using NIST traceable weights, and/or dilutions with Class A glassware.

### Handling Notes:

- Stability of the unopened product, when stored in compliance with the recommended conditions, is guaranteed through the expiration displayed on the product label and certificate. Contact Restek for additional opened product stability information, with the knowledge/understanding that open product stability is subject to the specific handling and environmental conditions to which the product is exposed. For your convenience Restek supplies deactivated vials with most standards packed in 2mL ampules. Larger volume deactivated vials are available through Restek as a custom ordered item. Additionally, Restek sells DMDCS for the purpose of glassware deactivation as catalog number 31861, which includes complete instructions.



Reagent

---

**8330LCSmix2\_00038**



110 Benner Circle  
Bellefonte, PA 16823-8812  
Tel: 1-814-353-1300  
Fax: 1-814-353-1309

www.restek.com

## CERTIFIED REFERENCE MATERIAL

# Certificate of Analysis

chromatographic plus



### FOR LABORATORY USE ONLY-READ SDS PRIOR TO USE.

*This Reference Material is intended for Laboratory Use Only as a standard for the qualitative and/or quantitative determination of the analyte(s) listed.*

**Catalog No. :** 31451 **Lot No.:** A0192237

**Description :** 8330 Calibration Mix #2  
8330 Calibration Std #2 1000µg/mL, Acetonitrile, 1mL/ampul

**Container Size :** 2 mL **Pkg Amt:** > 1 mL

**Expiration Date :** December 31, 2027 **Storage:** 10°C or colder

**Ship:** Ambient

### CERTIFIED VALUES

| Elution Order | Compound                   | CAS #      | Lot #       | Purity | Grav. Conc. (weight/volume) | Expanded Uncertainty * (95% C.L.; K=2) |
|---------------|----------------------------|------------|-------------|--------|-----------------------------|----------------------------------------|
| 1             | Tetryl                     | 479-45-8   | 211028JLM   | 99%    | 1,006.0 µg/mL               | +/- 46.9317                            |
| 2             | 4-Amino-2,6-dinitrotoluene | 19406-51-0 | ER070908-01 | 99%    | 1,004.0 µg/mL               | +/- 46.8384                            |
| 3             | 2-Amino-4,6-dinitrotoluene | 35572-78-2 | A210503-001 | 99%    | 1,006.0 µg/mL               | +/- 46.9317                            |
| 4             | 2,6-Dinitrotoluene         | 606-20-2   | BCCG1833    | 99%    | 1,004.0 µg/mL               | +/- 46.8384                            |
| 5             | 2-Nitrotoluene             | 88-72-2    | BCBZ7826    | 99%    | 1,008.0 µg/mL               | +/- 47.0250                            |
| 6             | 4-Nitrotoluene             | 99-99-0    | BCCB0171    | 99%    | 1,002.0 µg/mL               | +/- 46.7451                            |
| 7             | 3-Nitrotoluene             | 99-08-1    | 07329LG     | 99%    | 1,008.0 µg/mL               | +/- 47.0250                            |

\* Expanded Uncertainty displayed in same units as Grav. Conc.

**Solvent:** Acetonitrile  
**CAS #** 75-05-8  
**Purity** 99%

# Quality Confirmation Test

**Column:**

250mm x 4.6mm  
Ultra C18 (cat.# 9174575)

**Flow Rate:**

1.0 ml/min.

**Mobile Phase A:**

water:methanol (44:56 V/V)

**Mobile Phase B:****Mobile Phase Composition:**


100%A

**Det. Type:**

Wavelength: 210nm & 254nm

**Inj. Vol**

5µl



This chromatogram represents a general set of testing conditions chosen for product acceptance. For optimal results in your lab, conditions should be adjusted for your specific instrument, method, and application.

John Friedline - Operations Technician I

Date Mixed: 05-Dec-2022

Balance Serial # 1128342314

Christie Mills - Operations Tech II - ARM QC

Date Passed: 07-Dec-2022

Manufactured under Restek's ISO 9001:2015  
Registered Quality System  
Certificate #FM 80397



## General Certified Reference Material Notes

### Expiration Notes:

- Expiration date valid for unopened ampul stored in compliance with the recommended conditions.
- Uncertainty, concentration, and expiration of the CRM are based on the unopened product being stored according to the recommended condition found in the storage field.

### Purity Notes:

- Purity and/or chemical identity are determined by one or more of the following techniques: GC/FID, HPLC, GC/μECD, GC/MS, LC/MS, RI, and/or melting point.
- Compounds with a listed purity of less than 99% have been weight corrected to compensate for impurities and/or salts. A correction factor is used to calculate the amount of compound necessary to achieve the desired concentration of the parent compound in solution.
- Purity of isomeric compounds is reported as the sum of the isomers.
- Purity values are rounded to the nearest whole number.

### Certified Uncertainty Value Notes:

- The uncertainties are determined in accordance with ISO 17034 and Guide 35. The certified expanded uncertainty value includes gravimetric uncertainty, homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty and were combined using the following formula:

$$U_{combined\ uncertainty} = k \sqrt{u_{gravimetric}^2 + u_{homogeneity}^2 + u_{storage\ stability}^2 + u_{shipping\ stability}^2}$$

$k$  is a coverage factor of 2, which gives a level of confidence of approximately 95%.

- The packaged amount is the minimum sample size for which uncertainty is valid. The ampuls are over-filled to ensure that the minimum packaged amount can be sufficiently transferred.

### Manufacturing Notes:

- Concentration is based upon gravimetric preparation using either a balance whose calibration has been verified daily using NIST traceable weights, and/or dilutions with Class A glassware.

### Handling Notes:

- Stability of the unopened product, when stored in compliance with the recommended conditions, is guaranteed through the expiration displayed on the product label and certificate. Contact Restek for additional opened product stability information, with the knowledge/understanding that open product stability is subject to the specific handling and environmental conditions to which the product is exposed. For your convenience Restek supplies deactivated vials with most standards packed in 2mL ampuls. Larger volume deactivated vials are available through Restek as a custom ordered item. Additionally, Restek sells DMDCS for the purpose of glassware deactivation as catalog number 31861, which includes complete instructions.
- If any undissolved material is visible inside the ampul, sonicate the unopened ampul until the material is completely dissolved.

Reagent

---

**8330LCSmix2\_00039**



110 Benner Circle  
Bellefonte, PA 16823-8812  
Tel: 1-814-353-1300  
Fax: 1-814-353-1309

www.restek.com

## CERTIFIED REFERENCE MATERIAL

# Certificate of Analysis

chromatographic plus



### FOR LABORATORY USE ONLY-READ SDS PRIOR TO USE.

*This Reference Material is intended for Laboratory Use Only as a standard for the qualitative and/or quantitative determination of the analyte(s) listed.*

**Catalog No. :** 31451 **Lot No.:** A0192237

**Description :** 8330 Calibration Mix #2  
8330 Calibration Std #2 1000µg/mL, Acetonitrile, 1mL/ampul

**Container Size :** 2 mL **Pkg Amt:** > 1 mL

**Expiration Date :** December 31, 2027 **Storage:** 10°C or colder

**Ship:** Ambient

### CERTIFIED VALUES

| Elution Order | Compound                   | CAS #      | Lot #       | Purity | Grav. Conc. (weight/volume) | Expanded Uncertainty * (95% C.L.; K=2) |
|---------------|----------------------------|------------|-------------|--------|-----------------------------|----------------------------------------|
| 1             | Tetryl                     | 479-45-8   | 211028JLM   | 99%    | 1,006.0 µg/mL               | +/- 46.9317                            |
| 2             | 4-Amino-2,6-dinitrotoluene | 19406-51-0 | ER070908-01 | 99%    | 1,004.0 µg/mL               | +/- 46.8384                            |
| 3             | 2-Amino-4,6-dinitrotoluene | 35572-78-2 | A210503-001 | 99%    | 1,006.0 µg/mL               | +/- 46.9317                            |
| 4             | 2,6-Dinitrotoluene         | 606-20-2   | BCCG1833    | 99%    | 1,004.0 µg/mL               | +/- 46.8384                            |
| 5             | 2-Nitrotoluene             | 88-72-2    | BCBZ7826    | 99%    | 1,008.0 µg/mL               | +/- 47.0250                            |
| 6             | 4-Nitrotoluene             | 99-99-0    | BCCB0171    | 99%    | 1,002.0 µg/mL               | +/- 46.7451                            |
| 7             | 3-Nitrotoluene             | 99-08-1    | 07329LG     | 99%    | 1,008.0 µg/mL               | +/- 47.0250                            |

\* Expanded Uncertainty displayed in same units as Grav. Conc.

**Solvent:** Acetonitrile  
**CAS #** 75-05-8  
**Purity** 99%

# Quality Confirmation Test

**Column:**  
250mm x 4.6mm  
Ultra C18 (cat.# 9174575)

**Flow Rate:**  
1.0 ml/min.


**Mobile Phase A:**  
water:methanol (44:56 V/V)

**Mobile Phase B:**

**Mobile Phase Composition:**  
100%A

**Det. Type:**  
Wavelength: 210nm & 254nm

**Inj. Vol**  
5µl



This chromatogram represents a general set of testing conditions chosen for product acceptance. For optimal results in your lab, conditions should be adjusted for your specific instrument, method, and application.

John Friedline - Operations Technician I

**Date Mixed:** 05-Dec-2022

**Balance Serial #** 1128342314

Christie Mills - Operations Tech II - ARM QC

**Date Passed:** 07-Dec-2022

Manufactured under Restek's ISO 9001:2015  
Registered Quality System  
Certificate #FM 80397

## General Certified Reference Material Notes

### Expiration Notes:

- Expiration date valid for unopened ampul stored in compliance with the recommended conditions.
- Uncertainty, concentration, and expiration of the CRM are based on the unopened product being stored according to the recommended condition found in the storage field.

### Purity Notes:

- Purity and/or chemical identity are determined by one or more of the following techniques: GC/FID, HPLC, GC/μECD, GC/MS, LC/MS, RI, and/or melting point.
- Compounds with a listed purity of less than 99% have been weight corrected to compensate for impurities and/or salts. A correction factor is used to calculate the amount of compound necessary to achieve the desired concentration of the parent compound in solution.
- Purity of isomeric compounds is reported as the sum of the isomers.
- Purity values are rounded to the nearest whole number.

### Certified Uncertainty Value Notes:

- The uncertainties are determined in accordance with ISO 17034 and Guide 35. The certified expanded uncertainty value includes gravimetric uncertainty, homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty and were combined using the following formula:

$$U_{combined\ uncertainty} = k \sqrt{u_{gravimetric}^2 + u_{homogeneity}^2 + u_{storage\ stability}^2 + u_{shipping\ stability}^2}$$

$k$  is a coverage factor of 2, which gives a level of confidence of approximately 95%.

- The packaged amount is the minimum sample size for which uncertainty is valid. The ampuls are over-filled to ensure that the minimum packaged amount can be sufficiently transferred.

### Manufacturing Notes:

- Concentration is based upon gravimetric preparation using either a balance whose calibration has been verified daily using NIST traceable weights, and/or dilutions with Class A glassware.

### Handling Notes:

- Stability of the unopened product, when stored in compliance with the recommended conditions, is guaranteed through the expiration displayed on the product label and certificate. Contact Restek for additional opened product stability information, with the knowledge/understanding that open product stability is subject to the specific handling and environmental conditions to which the product is exposed. For your convenience Restek supplies deactivated vials with most standards packed in 2mL ampuls. Larger volume deactivated vials are available through Restek as a custom ordered item. Additionally, Restek sells DMDCS for the purpose of glassware deactivation as catalog number 31861, which includes complete instructions.
- If any undissolved material is visible inside the ampul, sonicate the unopened ampul until the material is completely dissolved.

Reagent

---

**8330PAsTkPS\_00070**

# CERTIFICATE OF ANALYSIS

Catalog No: M-8330-ADD-3

Description: Picric acid

Lot: 218031154-03

Solvent: Acetonitrile (50%)

Methanol (50%)

Hazards: Refer to SDS for complete safety information

Date Certified: Jul 7, 2021

Expiration: Aug 7, 2023

Sample Size: 1 mL

Components: 1

Storage Condition: Ambient (>5 °C)

Certified Reference Material



Signal Word: Danger



| Component   | CAS #   | Purity % | Prepared<br>Concentration <sup>2</sup> | Certified Analyte<br>Concentration <sup>1</sup> |
|-------------|---------|----------|----------------------------------------|-------------------------------------------------|
|             |         | (HPLC)   | (µg/mL)                                | (µg/mL)                                         |
| Picric acid | 88-89-1 | 99.1     | 100.1                                  | 99.2                                            |

31499

This Certified Reference Material was verified in accordance with ISO/IEC 17025

A product with a suffix (-1A, -2B, etc. or -01, -02, etc.) on its lot number has had its expiration date extended and is identical to the same lot number without the suffix.

<sup>2</sup> All weights are traceable through NIST, Test No. 822-275872-11

<sup>1</sup> Certified Analyte Concentration = Purity x Prepared Concentration.

The Uncertainty associated with the certified concentration reported on this certificate is  $\pm 2.4\%$ . This value is the combined expanded uncertainty and represents an estimated standard deviation equal to the positive square root of the total variation of the uncertainty of components. A normal distribution is assumed and a coverage factor of K=2 is chosen using approximately a 95% confidence level.


Labels and certificates follow U.S. Conventions in reporting numerical values: A comma (,) is used to separate units of one-thousand or greater. A period (.) is used as a decimal place marker.

The information on this certificate may not be reproduced without the express permission of the manufacturer. See reverse side for additional information

Hazard Information: Please refer to the SDS for information regarding the hazards associated with using this material.

This product was prepared according to in-house procedures and is guaranteed to be homogeneous.

Certified By:

  
Larry Decker, Organic QC Manager

Reagent

---

**8330PASTkPS\_00071**



# CERTIFICATE OF ANALYSIS

**Catalog No:** M-8330-ADD-3

**Description:** Picric acid

**Lot:** 218031154-04

**Solvent:** Acetonitrile (50%)

Methanol (50%)

**Hazards:** Refer to SDS for complete safety information

**Date Certified:** Aug 30, 2022

**Expiration:** Sep 30, 2024

**Sample Size:** 1 mL

**Components:** 1

**Storage Condition:** Ambient (>5 °C)

**Certified Reference Material**



**Signal Word:** Danger



| Component   | CAS #   | Purity % | Prepared<br>Concentration <sup>2</sup> | Certified Analyte<br>Concentration <sup>1</sup> |
|-------------|---------|----------|----------------------------------------|-------------------------------------------------|
|             |         | (HPLC)   | (µg/mL)                                | (µg/mL)                                         |
| Picric acid | 88-89-1 | 99.1     | 100.1                                  | 99.2                                            |

This Certified Reference Material was verified in accordance with ISO/IEC 17025

A product with a suffix (-1A, -2B, etc. or -01, -02, etc.) on its lot number has had its expiration date extended and is identical to the same lot number without the suffix.

<sup>2</sup> All weights are traceable through NIST, Test No. 822-275872-11

<sup>1</sup> Certified Analyte Concentration = Purity x Prepared Concentration.

The Uncertainty associated with the certified concentration reported on this certificate is  $\pm 2.4\%$ . This value is the combined expanded uncertainty and represents an estimated standard deviation equal to the positive square root of the total variation of the uncertainty of components. A normal distribution is assumed and a coverage factor of  $k=2$  is chosen using approximately a 95% confidence level.

Labels and certificates follow U.S. Conventions in reporting numerical values: A comma (,) is used to separate units of one-thousand or greater. A period (.) is used as a decimal place marker.

The information on this certificate may not be reproduced without the express permission of the manufacturer. See reverse side for additional information

Hazard Information: Please refer to the SDS for information regarding the hazards associated with using this material.

This product was prepared according to in-house procedures and is guaranteed to be homogeneous.

Certified By:



Larry Decker, Organic QC Manager

Reagent

---

**8330Surrogate\_00138**

**Preliminary Report**

Eurofins Denver

LCS, Lab Control Sample Report

Sample Path: \\chromfs\Denver\ChromData\CHHPLC\_X\20230202-118305.b\02020011.D

Lims ID: Surr138

Inj. Date: 02-Feb-2023 19:11:17

Worklist ID: 280-0118305-011

Instrument: CHHPLC\_X3

Method: 8330\_X3

| Compound                 | Amount Added | Amount Recovered | %Rec  | Limits 1<br>OB_Sonc_ |
|--------------------------|--------------|------------------|-------|----------------------|
| \$ 10 1,2-Dinitrobenzene | 0.5000       | 0.5112           | 102.2 | 78-119               |

Samples for Limit Group: 1, Lims Prep Method: 8330B\_Sonc\_10g

|                   |                   |                   |
|-------------------|-------------------|-------------------|
| 280-171573-C-4-B  | 280-171573-C-5-B  | 280-171573-C-6-B  |
| 280-171573-C-7-B  | 280-171573-C-8-B  | 280-171573-C-9-D  |
| 280-171573-C-10-B | 280-171573-C-11-B | 280-171573-B-12-B |
| 280-171573-B-13-B | 280-171573-B-14-B | 280-171573-B-15-B |
| 280-171573-B-16-B | 280-171573-C-42-B | 280-171573-B-43-B |
| 280-171586-A-1-C  | 280-171586-A-2-C  | 280-171586-A-3-C  |
| 280-171586-A-4-C  | 280-171586-A-5-C  | 280-171586-A-6-J  |
| 280-171586-A-7-C  |                   |                   |

Reagent

---

**8330Surrogate\_00143**

**Preliminary Report**

Eurofins Denver

LCS, Lab Control Sample Report

Sample Path: \\chromfs\Denver\ChromData\CHHPLC\_X\20230413-120394.b\04130011.D

Lims ID: 8330Surr143

Inj. Date: 13-Apr-2023 17:28:02

Worklist ID: 280-0120394-011

Instrument: CHHPLC\_X3

Method: 8330\_X3

| Compound                 | Amount Added | Amount Recovered | %Rec  | Limits 1<br>OB_Sonc_ | Limits 2<br>OB_Sonc_ |
|--------------------------|--------------|------------------|-------|----------------------|----------------------|
| \$ 10 1,2-Dinitrobenzene | 0.5000       | 0.5398           | 108.0 | 83-122               | 78-119               |

Samples for Limit Group: 1, Lims Prep Method: 8330B\_Sonc\_10g

|                   |                   |                   |
|-------------------|-------------------|-------------------|
| 280-174145-A-1-E  | 280-174145-A-2-C  | 280-174145-A-3-B  |
| 280-174145-A-4-B  | 280-174145-A-5-B  | 280-174145-A-6-B  |
| 280-174145-A-7-C  | 280-174145-A-8-B  | 280-174145-A-9-B  |
| 280-174145-A-10-B | 280-174145-A-11-B | 280-174145-A-12-B |
| 280-174145-A-13-B | 280-174145-A-14-B | 280-174145-A-15-D |
| 280-174145-A-16-B | 280-174145-A-17-B | 280-174145-A-18-B |
| 280-174145-A-19-B | 280-174145-A-20-B |                   |

Samples for Limit Group: 2, Lims Prep Method: 8330B\_Sonc\_10g

|                   |                   |                   |
|-------------------|-------------------|-------------------|
| 280-174388-A-1-C  | 280-174388-A-2-C  | 280-174388-A-3-C  |
| 280-174388-A-4-C  | 280-174388-A-5-C  | 280-174388-A-6-C  |
| 280-174388-A-7-G  | 280-174388-A-8-G  | 280-174388-A-9-C  |
| 280-174388-A-10-C | 280-174388-A-11-C | 280-174388-A-12-C |
| 280-174388-A-13-C | 280-174388-A-14-C | 280-174388-A-15-C |
| 280-174388-A-16-C | 280-174388-A-17-C | 280-174388-A-18-C |
| 280-174388-A-19-C | 280-174388-A-20-C |                   |

Reagent

---

**8330Surrogate\_00144**

**Preliminary Report**

Eurofins Denver

LCS, Lab Control Sample Report

Sample Path: \\chromfs\Denver\ChromData\CHHPLC\_X\20230517-121556.b\05170006.D

Lims ID: C18column:B16162

Inj. Date: 17-May-2023 17:57:37

Worklist ID: 280-0121556-006

Instrument: CHHPLC\_X3

Method: 8330\_X3

| Compound                 | Amount Added | Amount Recovered | %Rec  | Limits 1<br>3535 | Limits 2<br>3535 | Limits 3<br>3535 |
|--------------------------|--------------|------------------|-------|------------------|------------------|------------------|
| \$ 10 1,2-Dinitrobenzene | 0.5000       | 0.5038           | 100.8 | 83-119           | 63-127           |                  |

Samples for Limit Group: 1, Lims Prep Method: 3535

|                  |                  |                  |
|------------------|------------------|------------------|
| 280-176194-A-9-A | 280-176284-C-2-A | 280-176284-B-3-A |
| 280-176284-C-5-A | 280-176284-C-7-A | 280-176434-I-1-A |
| 280-176434-E-4-A | 280-176434-B-7-A | 280-176284-B-4-C |
| 280-176434-H-2-A | 280-176524-B-1-A | 280-176524-A-2-A |
| 280-176524-A-3-A | 280-176524-A-4-A | 280-176335-A-1-A |
| 280-176335-A-2-A | 280-176335-B-3-A | 280-176335-A-4-A |

Samples for Limit Group: 2, Lims Prep Method: 3535

|                  |                  |                  |
|------------------|------------------|------------------|
| 570-137880-R-1-A | 570-137880-S-2-A | 570-137880-S-3-A |
| 570-137880-R-4-A | 570-138095-R-2-A | 570-138095-S-3-A |
| 570-138095-R-4-A | 570-138095-S-5-A | 550-201994-O-1-A |
| 550-201994-O-3-A |                  |                  |

Samples for Limit Group: 3, Lims Prep Method: 3535

|                  |                 |                  |
|------------------|-----------------|------------------|
| 380-47016-AV-1-A | 380-47022-Q-1-A | 280-176543-B-1-A |
| 280-176543-B-2-A |                 |                  |

Reagent

---

**8330SurrStkSS\_00245**





110 Benner Circle  
Bellefonte, PA 16823-8812  
Tel: 1-814-353-1300  
Fax: 1-814-353-1309

www.restek.com

## CERTIFIED REFERENCE MATERIAL

# Certificate of Analysis

*chromatographic plus*



### FOR LABORATORY USE ONLY-READ SDS PRIOR TO USE.

*This Reference Material is intended for Laboratory Use Only as a standard for the qualitative and/or quantitative determination of the analyte(s) listed.*

**Catalog No. :** 31453 **Lot No.:** A0192220  
**Description :** 8330 Surrogate Mix  
8330 Surrogate Mix 1000 µg/mL, Methanol, 1mL/ampul  
**Container Size :** 2 mL **Pkg Amt:** > 1 mL  
**Expiration Date :** December 31, 2027 **Storage:** 10°C or colder  
**Ship:** Ambient

### CERTIFIED VALUES

| Elution Order | Compound           | CAS #    | Lot #    | Purity | Grav. Conc. (weight/volume) | Expanded Uncertainty* (95% C.L.; K=2) |
|---------------|--------------------|----------|----------|--------|-----------------------------|---------------------------------------|
| 1             | 1,2-Dinitrobenzene | 528-29-0 | MKCH6067 | 99%    | 1,000.0 µg/mL               | +/- 56.1888                           |

\* Expanded Uncertainty displayed in same units as Grav. Conc.

**Solvent:** Methanol  
**CAS #** 67-56-1  
**Purity** 99%

## Quality Confirmation Test

**Column:**

250mm x 4.6mm

Ultra C18 (cat.# 9174575)

**Flow Rate:**

1.0 ml/min.

**Mobile Phase A:**

water:methanol (44:56 V/V)

**Mobile Phase B:****Mobile Phase Composition:**


100%A

**Det. Type:**

Wavelength: 210nm &amp; 254nm

**Inj. Vol**


2µl



This chromatogram represents a general set of testing conditions chosen for product acceptance. For optimal results in your lab, conditions should be adjusted for your specific instrument, method, and application.

  
Alicia Leathers - Operation Technician I

**Date Mixed:** 04-Dec-2022**Balance Serial #** 1128353505

  
Jennifer Pollino - Operations Tech III - ARM QC

**Date Passed:** 06-Dec-2022

Manufactured under Restek's ISO 9001:2015  
Registered Quality System  
Certificate #FM 80397

## General Certified Reference Material Notes

### Expiration Notes:

- Expiration date valid for unopened ampul stored in compliance with the recommended conditions.
- Uncertainty, concentration, and expiration of the CRM are based on the unopened product being stored according to the recommended condition found in the storage field.

### Purity Notes:

- Purity and/or chemical identity are determined by one or more of the following techniques: GC/FID, HPLC, GC/μECD, GC/MS, LC/MS, RI, and/or melting point.
- Compounds with a listed purity of less than 99% have been weight corrected to compensate for impurities and/or salts. A correction factor is used to calculate the amount of compound necessary to achieve the desired concentration of the parent compound in solution.
- Purity of isomeric compounds is reported as the sum of the isomers.
- Purity values are rounded to the nearest whole number.

### Certified Uncertainty Value Notes:

- The uncertainties are determined in accordance with ISO 17034 and Guide 35. The certified expanded uncertainty value includes gravimetric uncertainty, homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty and were combined using the following formula:

$$U_{combined\ uncertainty} = k \sqrt{u_{gravimetric}^2 + u_{homogeneity}^2 + u_{storage\ stability}^2 + u_{shipping\ stability}^2}$$

$k$  is a coverage factor of 2, which gives a level of confidence of approximately 95%.

- The packaged amount is the minimum sample size for which uncertainty is valid. The ampuls are over-filled to ensure that the minimum packaged amount can be sufficiently transferred.

### Manufacturing Notes:

- Concentration is based upon gravimetric preparation using either a balance whose calibration has been verified daily using NIST traceable weights, and/or dilutions with Class A glassware.

### Handling Notes:

- Stability of the unopened product, when stored in compliance with the recommended conditions, is guaranteed through the expiration displayed on the product label and certificate. Contact Restek for additional opened product stability information, with the knowledge/understanding that open product stability is subject to the specific handling and environmental conditions to which the product is exposed. For your convenience Restek supplies deactivated vials with most standards packed in 2mL ampuls. Larger volume deactivated vials are available through Restek as a custom ordered item. Additionally, Restek sells DMDCS for the purpose of glassware deactivation as catalog number 31861, which includes complete instructions.
- If any undissolved material is visible inside the ampul, sonicate the unopened ampul until the material is completely dissolved.



Reagent

---

**8330SurrStkSS\_00246**



110 Benner Circle  
Bellefonte, PA 16823-8812  
Tel: 1-814-353-1300  
Fax: 1-814-353-1309

www.restek.com

## CERTIFIED REFERENCE MATERIAL

# Certificate of Analysis

*chromatographic plus*



### FOR LABORATORY USE ONLY-READ SDS PRIOR TO USE.

*This Reference Material is intended for Laboratory Use Only as a standard for the qualitative and/or quantitative determination of the analyte(s) listed.*

**Catalog No. :** 31453 **Lot No.:** A0192220

**Description :** 8330 Surrogate Mix  
8330 Surrogate Mix 1000 µg/mL, Methanol, 1mL/ampul

**Container Size :** 2 mL **Pkg Amt:** > 1 mL

**Expiration Date :** December 31, 2027 **Storage:** 10°C or colder

**Ship:** Ambient

### CERTIFIED VALUES

| Elution Order | Compound           | CAS #    | Lot #    | Purity | Grav. Conc. (weight/volume) | Expanded Uncertainty* (95% C.L.; K=2) |
|---------------|--------------------|----------|----------|--------|-----------------------------|---------------------------------------|
| 1             | 1,2-Dinitrobenzene | 528-29-0 | MKCH6067 | 99%    | 1,000.0 µg/mL               | +/- 56.1888                           |

\* Expanded Uncertainty displayed in same units as Grav. Conc.

**Solvent:** Methanol  
**CAS #** 67-56-1  
**Purity** 99%

## Quality Confirmation Test

**Column:**

250mm x 4.6mm

Ultra C18 (cat.# 9174575)

**Flow Rate:**

1.0 ml/min.

**Mobile Phase A:**

water:methanol (44:56 V/V)

**Mobile Phase B:****Mobile Phase Composition:**


100%A

**Det. Type:**

Wavelength: 210nm &amp; 254nm

**Inj. Vol**

2µl




This chromatogram represents a general set of testing conditions chosen for product acceptance. For optimal results in your lab, conditions should be adjusted for your specific instrument, method, and application.

  
Alicia Leathers - Operation Technician I

Date Mixed: 04-Dec-2022

Balance Serial # 1128353505

  
Jennifer Pollino - Operations Tech III - ARM QC

Date Passed: 06-Dec-2022

Manufactured under Restek's ISO 9001:2015  
Registered Quality System  
Certificate #FM 80397

## General Certified Reference Material Notes

### Expiration Notes:

- Expiration date valid for unopened ampul stored in compliance with the recommended conditions.
- Uncertainty, concentration, and expiration of the CRM are based on the unopened product being stored according to the recommended condition found in the storage field.

### Purity Notes:

- Purity and/or chemical identity are determined by one or more of the following techniques: GC/FID, HPLC, GC/μECD, GC/MS, LC/MS, RI, and/or melting point.
- Compounds with a listed purity of less than 99% have been weight corrected to compensate for impurities and/or salts. A correction factor is used to calculate the amount of compound necessary to achieve the desired concentration of the parent compound in solution.
- Purity of isomeric compounds is reported as the sum of the isomers.
- Purity values are rounded to the nearest whole number.

### Certified Uncertainty Value Notes:

- The uncertainties are determined in accordance with ISO 17034 and Guide 35. The certified expanded uncertainty value includes gravimetric uncertainty, homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty and were combined using the following formula:

$$U_{combined\ uncertainty} = k \sqrt{u_{gravimetric}^2 + u_{homogeneity}^2 + u_{storage\ stability}^2 + u_{shipping\ stability}^2}$$

$k$  is a coverage factor of 2, which gives a level of confidence of approximately 95%.

- The packaged amount is the minimum sample size for which uncertainty is valid. The ampuls are over-filled to ensure that the minimum packaged amount can be sufficiently transferred.

### Manufacturing Notes:

- Concentration is based upon gravimetric preparation using either a balance whose calibration has been verified daily using NIST traceable weights, and/or dilutions with Class A glassware.

### Handling Notes:

- Stability of the unopened product, when stored in compliance with the recommended conditions, is guaranteed through the expiration displayed on the product label and certificate. Contact Restek for additional opened product stability information, with the knowledge/understanding that open product stability is subject to the specific handling and environmental conditions to which the product is exposed. For your convenience Restek supplies deactivated vials with most standards packed in 2mL ampuls. Larger volume deactivated vials are available through Restek as a custom ordered item. Additionally, Restek sells DMDCS for the purpose of glassware deactivation as catalog number 31861, which includes complete instructions.
- If any undissolved material is visible inside the ampul, sonicate the unopened ampul until the material is completely dissolved.





Reagent

---

**8330SurrStkSS\_00247**



110 Benner Circle  
Bellefonte, PA 16823-8812  
Tel: 1-814-353-1300  
Fax: 1-814-353-1309

www.restek.com

## CERTIFIED REFERENCE MATERIAL

# Certificate of Analysis

*chromatographic plus*



### FOR LABORATORY USE ONLY-READ SDS PRIOR TO USE.

*This Reference Material is intended for Laboratory Use Only as a standard for the qualitative and/or quantitative determination of the analyte(s) listed.*

**Catalog No. :** 31453 **Lot No.:** A0192220

**Description :** 8330 Surrogate Mix  
8330 Surrogate Mix 1000 µg/mL, Methanol, 1mL/ampul

**Container Size :** 2 mL **Pkg Amt:** > 1 mL

**Expiration Date :** December 31, 2027 **Storage:** 10°C or colder

**Ship:** Ambient

### CERTIFIED VALUES

| Elution Order | Compound           | CAS #    | Lot #    | Purity | Grav. Conc. (weight/volume) | Expanded Uncertainty* (95% C.L.; K=2) |
|---------------|--------------------|----------|----------|--------|-----------------------------|---------------------------------------|
| 1             | 1,2-Dinitrobenzene | 528-29-0 | MKCH6067 | 99%    | 1,000.0 µg/mL               | +/- 56.1888                           |

\* Expanded Uncertainty displayed in same units as Grav. Conc.

**Solvent:** Methanol  
**CAS #** 67-56-1  
**Purity** 99%

## Quality Confirmation Test

**Column:**

250mm x 4.6mm  
Ultra C18 (cat.# 9174575)

**Flow Rate:**

1.0 ml/min.

**Mobile Phase A:**

water:methanol (44:56 V/V)

**Mobile Phase B:****Mobile Phase Composition:**


100%A

**Det. Type:**

Wavelength: 210nm & 254nm

**Inj. Vol**

2µl




This chromatogram represents a general set of testing conditions chosen for product acceptance. For optimal results in your lab, conditions should be adjusted for your specific instrument, method, and application.

  
Alicia Leathers - Operation Technician I

Date Mixed: 04-Dec-2022

Balance Serial # 1128353505

  
Jennifer Pollino - Operations Tech III - ARM QC

Date Passed: 06-Dec-2022

Manufactured under Restek's ISO 9001:2015  
Registered Quality System  
Certificate #FM 80397

## General Certified Reference Material Notes

### Expiration Notes:

- Expiration date valid for unopened ampul stored in compliance with the recommended conditions.
- Uncertainty, concentration, and expiration of the CRM are based on the unopened product being stored according to the recommended condition found in the storage field.

### Purity Notes:

- Purity and/or chemical identity are determined by one or more of the following techniques: GC/FID, HPLC, GC/μECD, GC/MS, LC/MS, RI, and/or melting point.
- Compounds with a listed purity of less than 99% have been weight corrected to compensate for impurities and/or salts. A correction factor is used to calculate the amount of compound necessary to achieve the desired concentration of the parent compound in solution.
- Purity of isomeric compounds is reported as the sum of the isomers.
- Purity values are rounded to the nearest whole number.

### Certified Uncertainty Value Notes:

- The uncertainties are determined in accordance with ISO 17034 and Guide 35. The certified expanded uncertainty value includes gravimetric uncertainty, homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty and were combined using the following formula:

$$U_{combined\ uncertainty} = k \sqrt{u_{gravimetric}^2 + u_{homogeneity}^2 + u_{storage\ stability}^2 + u_{shipping\ stability}^2}$$

$k$  is a coverage factor of 2, which gives a level of confidence of approximately 95%.

- The packaged amount is the minimum sample size for which uncertainty is valid. The ampuls are over-filled to ensure that the minimum packaged amount can be sufficiently transferred.

### Manufacturing Notes:

- Concentration is based upon gravimetric preparation using either a balance whose calibration has been verified daily using NIST traceable weights, and/or dilutions with Class A glassware.

### Handling Notes:

- Stability of the unopened product, when stored in compliance with the recommended conditions, is guaranteed through the expiration displayed on the product label and certificate. Contact Restek for additional opened product stability information, with the knowledge/understanding that open product stability is subject to the specific handling and environmental conditions to which the product is exposed. For your convenience Restek supplies deactivated vials with most standards packed in 2mL ampuls. Larger volume deactivated vials are available through Restek as a custom ordered item. Additionally, Restek sells DMDCS for the purpose of glassware deactivation as catalog number 31861, which includes complete instructions.
- If any undissolved material is visible inside the ampul, sonicate the unopened ampul until the material is completely dissolved.



Reagent

---

**8330SurrStkSS\_00248**



110 Benner Circle  
Bellefonte, PA 16823-8812  
Tel: 1-814-353-1300  
Fax: 1-814-353-1309

www.restek.com

## CERTIFIED REFERENCE MATERIAL

# Certificate of Analysis

*chromatographic plus*



### FOR LABORATORY USE ONLY-READ SDS PRIOR TO USE.

*This Reference Material is intended for Laboratory Use Only as a standard for the qualitative and/or quantitative determination of the analyte(s) listed.*

**Catalog No. :** 31453 **Lot No.:** A0192220

**Description :** 8330 Surrogate Mix  
8330 Surrogate Mix 1000 µg/mL, Methanol, 1mL/ampul

**Container Size :** 2 mL **Pkg Amt:** > 1 mL

**Expiration Date :** December 31, 2027 **Storage:** 10°C or colder

**Ship:** Ambient

### CERTIFIED VALUES

| Elution Order | Compound           | CAS #    | Lot #    | Purity | Grav. Conc. (weight/volume) | Expanded Uncertainty* (95% C.L.; K=2) |
|---------------|--------------------|----------|----------|--------|-----------------------------|---------------------------------------|
| 1             | 1,2-Dinitrobenzene | 528-29-0 | MKCH6067 | 99%    | 1,000.0 µg/mL               | +/- 56.1888                           |

\* Expanded Uncertainty displayed in same units as Grav. Conc.

**Solvent:** Methanol  
**CAS #** 67-56-1  
**Purity** 99%



## Quality Confirmation Test

**Column:**

250mm x 4.6mm

Ultra C18 (cat.# 9174575)

**Flow Rate:**

1.0 ml/min.

**Mobile Phase A:**

water:methanol (44:56 V/V)

**Mobile Phase B:****Mobile Phase Composition:**


100%A

**Det. Type:**

Wavelength: 210nm &amp; 254nm

**Inj. Vol**

2µl




This chromatogram represents a general set of testing conditions chosen for product acceptance. For optimal results in your lab, conditions should be adjusted for your specific instrument, method, and application.

  
Alicia Leathers - Operation Technician I

Date Mixed: 04-Dec-2022

Balance Serial # 1128353505

  
Jennifer Pollino - Operations Tech III - ARM QC

Date Passed: 06-Dec-2022

Manufactured under Restek's ISO 9001:2015  
Registered Quality System  
Certificate #FM 80397

## General Certified Reference Material Notes

### Expiration Notes:

- Expiration date valid for unopened ampul stored in compliance with the recommended conditions.
- Uncertainty, concentration, and expiration of the CRM are based on the unopened product being stored according to the recommended condition found in the storage field.

### Purity Notes:

- Purity and/or chemical identity are determined by one or more of the following techniques: GC/FID, HPLC, GC/μECD, GC/MS, LC/MS, RI, and/or melting point.
- Compounds with a listed purity of less than 99% have been weight corrected to compensate for impurities and/or salts. A correction factor is used to calculate the amount of compound necessary to achieve the desired concentration of the parent compound in solution.
- Purity of isomeric compounds is reported as the sum of the isomers.
- Purity values are rounded to the nearest whole number.

### Certified Uncertainty Value Notes:

- The uncertainties are determined in accordance with ISO 17034 and Guide 35. The certified expanded uncertainty value includes gravimetric uncertainty, homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty and were combined using the following formula:

$$U_{combined\ uncertainty} = k \sqrt{u_{gravimetric}^2 + u_{homogeneity}^2 + u_{storage\ stability}^2 + u_{shipping\ stability}^2}$$

$k$  is a coverage factor of 2, which gives a level of confidence of approximately 95%.

- The packaged amount is the minimum sample size for which uncertainty is valid. The ampuls are over-filled to ensure that the minimum packaged amount can be sufficiently transferred.

### Manufacturing Notes:

- Concentration is based upon gravimetric preparation using either a balance whose calibration has been verified daily using NIST traceable weights, and/or dilutions with Class A glassware.

### Handling Notes:

- Stability of the unopened product, when stored in compliance with the recommended conditions, is guaranteed through the expiration displayed on the product label and certificate. Contact Restek for additional opened product stability information, with the knowledge/understanding that open product stability is subject to the specific handling and environmental conditions to which the product is exposed. For your convenience Restek supplies deactivated vials with most standards packed in 2mL ampuls. Larger volume deactivated vials are available through Restek as a custom ordered item. Additionally, Restek sells DMDCS for the purpose of glassware deactivation as catalog number 31861, which includes complete instructions.
- If any undissolved material is visible inside the ampul, sonicate the unopened ampul until the material is completely dissolved.



Reagent

---

**8330SurrStkSS\_00249**



110 Benner Circle  
Bellefonte, PA 16823-8812  
Tel: 1-814-353-1300  
Fax: 1-814-353-1309

www.restek.com

## CERTIFIED REFERENCE MATERIAL

# Certificate of Analysis

*chromatographic plus*



### FOR LABORATORY USE ONLY-READ SDS PRIOR TO USE.

*This Reference Material is intended for Laboratory Use Only as a standard for the qualitative and/or quantitative determination of the analyte(s) listed.*

**Catalog No. :** 31453 **Lot No.:** A0192220  
**Description :** 8330 Surrogate Mix  
8330 Surrogate Mix 1000 µg/mL, Methanol, 1mL/ampul  
**Container Size :** 2 mL **Pkg Amt:** > 1 mL  
**Expiration Date :** December 31, 2027 **Storage:** 10°C or colder  
**Ship:** Ambient

### CERTIFIED VALUES

| Elution Order | Compound           | CAS #    | Lot #    | Purity | Grav. Conc. (weight/volume) | Expanded Uncertainty* (95% C.L.; K=2) |
|---------------|--------------------|----------|----------|--------|-----------------------------|---------------------------------------|
| 1             | 1,2-Dinitrobenzene | 528-29-0 | MKCH6067 | 99%    | 1,000.0 µg/mL               | +/- 56.1888                           |

\* Expanded Uncertainty displayed in same units as Grav. Conc.

**Solvent:** Methanol  
**CAS #** 67-56-1  
**Purity** 99%

## Quality Confirmation Test

**Column:**

250mm x 4.6mm  
Ultra C18 (cat.# 9174575)

**Flow Rate:**

1.0 ml/min.

**Mobile Phase A:**

water:methanol (44:56 V/V)

**Mobile Phase B:****Mobile Phase Composition:**


100%A

**Det. Type:**

Wavelength: 210nm & 254nm

**Inj. Vol**

2µl



This chromatogram represents a general set of testing conditions chosen for product acceptance. For optimal results in your lab, conditions should be adjusted for your specific instrument, method, and application.

  
Alicia Leathers - Operation Technician I

Date Mixed: 04-Dec-2022

Balance Serial # 1128353505

  
Jennifer Pollino - Operations Tech III - ARM QC

Date Passed: 06-Dec-2022

Manufactured under Restek's ISO 9001:2015  
Registered Quality System  
Certificate #FM 80397

## General Certified Reference Material Notes

### Expiration Notes:

- Expiration date valid for unopened ampul stored in compliance with the recommended conditions.
- Uncertainty, concentration, and expiration of the CRM are based on the unopened product being stored according to the recommended condition found in the storage field.

### Purity Notes:

- Purity and/or chemical identity are determined by one or more of the following techniques: GC/FID, HPLC, GC/ $\mu$ ECD, GC/MS, LC/MS, RI, and/or melting point.
- Compounds with a listed purity of less than 99% have been weight corrected to compensate for impurities and/or salts. A correction factor is used to calculate the amount of compound necessary to achieve the desired concentration of the parent compound in solution.
- Purity of isomeric compounds is reported as the sum of the isomers.
- Purity values are rounded to the nearest whole number.

### Certified Uncertainty Value Notes:

- The uncertainties are determined in accordance with ISO 17034 and Guide 35. The certified expanded uncertainty value includes gravimetric uncertainty, homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty and were combined using the following formula:

$$U_{combined\ uncertainty} = k \sqrt{u_{gravimetric}^2 + u_{homogeneity}^2 + u_{storage\ stability}^2 + u_{shipping\ stability}^2}$$

$k$  is a coverage factor of 2, which gives a level of confidence of approximately 95%.

- The packaged amount is the minimum sample size for which uncertainty is valid. The ampuls are over-filled to ensure that the minimum packaged amount can be sufficiently transferred.

### Manufacturing Notes:

- Concentration is based upon gravimetric preparation using either a balance whose calibration has been verified daily using NIST traceable weights, and/or dilutions with Class A glassware.

### Handling Notes:

- Stability of the unopened product, when stored in compliance with the recommended conditions, is guaranteed through the expiration displayed on the product label and certificate. Contact Restek for additional opened product stability information, with the knowledge/understanding that open product stability is subject to the specific handling and environmental conditions to which the product is exposed. For your convenience Restek supplies deactivated vials with most standards packed in 2mL ampuls. Larger volume deactivated vials are available through Restek as a custom ordered item. Additionally, Restek sells DMDCS for the purpose of glassware deactivation as catalog number 31861, which includes complete instructions.
- If any undissolved material is visible inside the ampul, sonicate the unopened ampul until the material is completely dissolved.





Reagent

---

**8330SurrStkSS\_00263**



110 Benner Circle  
Bellefonte, PA 16823-8812  
Tel: 1-814-353-1300  
Fax: 1-814-353-1309

www.restek.com

## CERTIFIED REFERENCE MATERIAL

# Certificate of Analysis

*chromatographic plus*



### FOR LABORATORY USE ONLY-READ SDS PRIOR TO USE.

*This Reference Material is intended for Laboratory Use Only as a standard for the qualitative and/or quantitative determination of the analyte(s) listed.*

**Catalog No. :** 31453 **Lot No.:** A0193444

**Description :** 8330 Surrogate Mix  
8330 Surrogate Mix 1000 µg/mL, Methanol, 1mL/ampul

**Container Size :** 2 mL **Pkg Amt:** > 1 mL

**Expiration Date :** January 31, 2028 **Storage:** 10°C or colder

**Ship:** Ambient

### CERTIFIED VALUES

| Elution Order | Compound           | CAS #    | Lot #    | Purity | Grav. Conc. (weight/volume) | Expanded Uncertainty * (95% C.L.; K=2) |
|---------------|--------------------|----------|----------|--------|-----------------------------|----------------------------------------|
| 1             | 1,2-Dinitrobenzene | 528-29-0 | MKCH6067 | 99%    | 1,004.0 µg/mL               | +/- 56.4136                            |

\* Expanded Uncertainty displayed in same units as Grav. Conc.

**Solvent:** Methanol  
**CAS #** 67-56-1  
**Purity** 99%

## Quality Confirmation Test

**Column:**

250mm x 4.6mm  
Ultra C18 (cat.# 9174575)

**Flow Rate:**

1.0 ml/min.

**Mobile Phase A:**

water:methanol (44:56 V/V)

**Mobile Phase B:****Mobile Phase Composition:**


100%A

**Det. Type:**

Wavelength: 210nm & 254nm

**Inj. Vol**

2µl



This chromatogram represents a general set of testing conditions chosen for product acceptance. For optimal results in your lab, conditions should be adjusted for your specific instrument, method, and application.

  
Bryan Snyder - Operations Tech I

Date Mixed: 11-Jan-2023

Balance Serial # 1127510105

  
Christie Mills - Operations Tech II - ARM QC

Date Passed: 16-Jan-2023

Manufactured under Restek's ISO 9001:2015  
Registered Quality System  
Certificate #FM 80397

## General Certified Reference Material Notes

### Expiration Notes:

- Expiration date valid for unopened ampul stored in compliance with the recommended conditions.
- Uncertainty, concentration, and expiration of the CRM are based on the unopened product being stored according to the recommended condition found in the storage field.

### Purity Notes:

- Purity and/or chemical identity are determined by one or more of the following techniques: GC/FID, HPLC, GC/μECD, GC/MS, LC/MS, RI, and/or melting point.
- Compounds with a listed purity of less than 99% have been weight corrected to compensate for impurities and/or salts. A correction factor is used to calculate the amount of compound necessary to achieve the desired concentration of the parent compound in solution.
- Purity of isomeric compounds is reported as the sum of the isomers.
- Purity values are rounded to the nearest whole number.

### Certified Uncertainty Value Notes:

- The uncertainties are determined in accordance with ISO 17034 and Guide 35. The certified expanded uncertainty value includes gravimetric uncertainty, homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty and were combined using the following formula:

$$U_{combined\ uncertainty} = k \sqrt{u_{gravimetric}^2 + u_{homogeneity}^2 + u_{storage\ stability}^2 + u_{shipping\ stability}^2}$$

$k$  is a coverage factor of 2, which gives a level of confidence of approximately 95%.

- The packaged amount is the minimum sample size for which uncertainty is valid. The ampuls are over-filled to ensure that the minimum packaged amount can be sufficiently transferred.

### Manufacturing Notes:

- Concentration is based upon gravimetric preparation using either a balance whose calibration has been verified daily using NIST traceable weights, and/or dilutions with Class A glassware.

### Handling Notes:

- Stability of the unopened product, when stored in compliance with the recommended conditions, is guaranteed through the expiration displayed on the product label and certificate. Contact Restek for additional opened product stability information, with the knowledge/understanding that open product stability is subject to the specific handling and environmental conditions to which the product is exposed. For your convenience Restek supplies deactivated vials with most standards packed in 2mL ampuls. Larger volume deactivated vials are available through Restek as a custom ordered item. Additionally, Restek sells DMDCS for the purpose of glassware deactivation as catalog number 31861, which includes complete instructions.
- If any undissolved material is visible inside the ampul, sonicate the unopened ampul until the material is completely dissolved.

Reagent

---

**8330SurrStkSS\_00265**



110 Benner Circle  
Bellefonte, PA 16823-8812  
Tel: 1-814-353-1300  
Fax: 1-814-353-1309

www.restek.com

## CERTIFIED REFERENCE MATERIAL

# Certificate of Analysis

*chromatographic plus*



### FOR LABORATORY USE ONLY-READ SDS PRIOR TO USE.

*This Reference Material is intended for Laboratory Use Only as a standard for the qualitative and/or quantitative determination of the analyte(s) listed.*

**Catalog No. :** 31453 **Lot No.:** A0193444  
**Description :** 8330 Surrogate Mix  
8330 Surrogate Mix 1000 µg/mL, Methanol, 1mL/ampul  
**Container Size :** 2 mL **Pkg Amt:** > 1 mL  
**Expiration Date :** January 31, 2028 **Storage:** 10°C or colder  
**Ship:** Ambient

### CERTIFIED VALUES

| Elution Order | Compound           | CAS #    | Lot #    | Purity | Grav. Conc. (weight/volume) | Expanded Uncertainty * (95% C.L.; K=2) |
|---------------|--------------------|----------|----------|--------|-----------------------------|----------------------------------------|
| 1             | 1,2-Dinitrobenzene | 528-29-0 | MKCH6067 | 99%    | 1,004.0 µg/mL               | +/- 56.4136                            |

\* Expanded Uncertainty displayed in same units as Grav. Conc.

**Solvent:** Methanol  
**CAS #** 67-56-1  
**Purity** 99%

# Quality Confirmation Test

**Column:**

250mm x 4.6mm  
Ultra C18 (cat.# 9174575)

**Flow Rate:**

1.0 ml/min.

**Mobile Phase A:**

water:methanol (44:56 V/V)

**Mobile Phase B:****Mobile Phase Composition:**


100%A

**Det. Type:**

Wavelength: 210nm & 254nm

**Inj. Vol**

2µl



This chromatogram represents a general set of testing conditions chosen for product acceptance. For optimal results in your lab, conditions should be adjusted for your specific instrument, method, and application.

  
Bryan Snyder - Operations Tech I

Date Mixed: 11-Jan-2023

Balance Serial # 1127510105

  
Christie Mills - Operations Tech II - ARM QC

Date Passed: 16-Jan-2023

Manufactured under Restek's ISO 9001:2015  
Registered Quality System  
Certificate #FM 80397

## General Certified Reference Material Notes

### Expiration Notes:

- Expiration date valid for unopened ampul stored in compliance with the recommended conditions.
- Uncertainty, concentration, and expiration of the CRM are based on the unopened product being stored according to the recommended condition found in the storage field.

### Purity Notes:

- Purity and/or chemical identity are determined by one or more of the following techniques: GC/FID, HPLC, GC/μECD, GC/MS, LC/MS, RI, and/or melting point.
- Compounds with a listed purity of less than 99% have been weight corrected to compensate for impurities and/or salts. A correction factor is used to calculate the amount of compound necessary to achieve the desired concentration of the parent compound in solution.
- Purity of isomeric compounds is reported as the sum of the isomers.
- Purity values are rounded to the nearest whole number.

### Certified Uncertainty Value Notes:

- The uncertainties are determined in accordance with ISO 17034 and Guide 35. The certified expanded uncertainty value includes gravimetric uncertainty, homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty and were combined using the following formula:

$$U_{combined\ uncertainty} = k \sqrt{u_{gravimetric}^2 + u_{homogeneity}^2 + u_{storage\ stability}^2 + u_{shipping\ stability}^2}$$

$k$  is a coverage factor of 2, which gives a level of confidence of approximately 95%.

- The packaged amount is the minimum sample size for which uncertainty is valid. The ampuls are over-filled to ensure that the minimum packaged amount can be sufficiently transferred.

### Manufacturing Notes:

- Concentration is based upon gravimetric preparation using either a balance whose calibration has been verified daily using NIST traceable weights, and/or dilutions with Class A glassware.

### Handling Notes:

- Stability of the unopened product, when stored in compliance with the recommended conditions, is guaranteed through the expiration displayed on the product label and certificate. Contact Restek for additional opened product stability information, with the knowledge/understanding that open product stability is subject to the specific handling and environmental conditions to which the product is exposed. For your convenience Restek supplies deactivated vials with most standards packed in 2mL ampuls. Larger volume deactivated vials are available through Restek as a custom ordered item. Additionally, Restek sells DMDCS for the purpose of glassware deactivation as catalog number 31861, which includes complete instructions.
- If any undissolved material is visible inside the ampul, sonicate the unopened ampul until the material is completely dissolved.



Reagent

---

**8330SurrStkSS\_00266**



110 Benner Circle  
Bellefonte, PA 16823-8812  
Tel: 1-814-353-1300  
Fax: 1-814-353-1309

www.restek.com

## CERTIFIED REFERENCE MATERIAL

# Certificate of Analysis

chromatographic plus



### FOR LABORATORY USE ONLY-READ SDS PRIOR TO USE.

*This Reference Material is intended for Laboratory Use Only as a standard for the qualitative and/or quantitative determination of the analyte(s) listed.*

Catalog No. : 31453 Lot No.: A0194831  
Description : 8330 Surrogate Mix  
8330 Surrogate Mix 1000 µg/mL, Methanol, 1mL/ampul  
Container Size : 2 mL Pkg Amt: > 1 mL  
Expiration Date : February 29, 2028 Storage: 10°C or colder  
Ship: Ambient

### CERTIFIED VALUES

| Elution Order | Compound           | CAS #    | Lot #    | Purity | Grav. Conc. (weight/volume) | Expanded Uncertainty * (95% C.L.; K=2) |
|---------------|--------------------|----------|----------|--------|-----------------------------|----------------------------------------|
| 1             | 1,2-Dinitrobenzene | 528-29-0 | MKCH6067 | 99%    | 1,004.8 µg/mL               | +/- 56.4585                            |

\* Expanded Uncertainty displayed in same units as Grav. Conc.

Solvent: Methanol  
CAS # 67-56-1  
Purity 99%

## Quality Confirmation Test

**Column:**

250mm x 4.6mm  
Ultra C18 (cat.# 9174575)

**Flow Rate:**

1.0 ml/min.

**Mobile Phase A:**

water:methanol (44:56 V/V)

**Mobile Phase B:****Mobile Phase Composition:**


100%A

**Det. Type:**

Wavelength: 210nm & 254nm

**Inj. Vol**

2µl



This chromatogram represents a general set of testing conditions chosen for product acceptance. For optimal results in your lab, conditions should be adjusted for your specific instrument, method, and application.

Stacey Wanner - Operations Technician I

Date Mixed: 20-Feb-2023

Balance Serial # B442140311

Jennifer Pollino - Operations Tech III - ARM QC

Date Passed: 21-Feb-2023

Manufactured under Restek's ISO 9001:2015  
Registered Quality System  
Certificate #FM 80397

## General Certified Reference Material Notes

### Expiration Notes:

- Expiration date valid for unopened ampul stored in compliance with the recommended conditions.
- Uncertainty, concentration, and expiration of the CRM are based on the unopened product being stored according to the recommended condition found in the storage field.

### Purity Notes:

- Purity and/or chemical identity are determined by one or more of the following techniques: GC/FID, HPLC, GC/μECD, GC/MS, LC/MS, RI, and/or melting point.
- Compounds with a listed purity of less than 99% have been weight corrected to compensate for impurities and/or salts. A correction factor is used to calculate the amount of compound necessary to achieve the desired concentration of the parent compound in solution.
- Purity of isomeric compounds is reported as the sum of the isomers.
- Purity values are rounded to the nearest whole number.

### Certified Uncertainty Value Notes:

- The uncertainties are determined in accordance with ISO 17034 and Guide 35. The certified expanded uncertainty value includes gravimetric uncertainty, homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty and were combined using the following formula:

$$U_{combined\ uncertainty} = k \sqrt{u_{gravimetric}^2 + u_{homogeneity}^2 + u_{storage\ stability}^2 + u_{shipping\ stability}^2}$$

$k$  is a coverage factor of 2, which gives a level of confidence of approximately 95%.

- The packaged amount is the minimum sample size for which uncertainty is valid. The ampuls are over-filled to ensure that the minimum packaged amount can be sufficiently transferred.

### Manufacturing Notes:

- Concentration is based upon gravimetric preparation using either a balance whose calibration has been verified daily using NIST traceable weights, and/or dilutions with Class A glassware.

### Handling Notes:

- Stability of the unopened product, when stored in compliance with the recommended conditions, is guaranteed through the expiration displayed on the product label and certificate. Contact Restek for additional opened product stability information, with the knowledge/understanding that open product stability is subject to the specific handling and environmental conditions to which the product is exposed. For your convenience Restek supplies deactivated vials with most standards packed in 2mL ampuls. Larger volume deactivated vials are available through Restek as a custom ordered item. Additionally, Restek sells DMDCS for the purpose of glassware deactivation as catalog number 31861, which includes complete instructions.
- If any undissolved material is visible inside the ampul, sonicate the unopened ampul until the material is completely dissolved.



Reagent

---

**8330SurrStkSS\_00267**



110 Benner Circle  
Bellefonte, PA 16823-8812  
Tel: 1-814-353-1300  
Fax: 1-814-353-1309

www.restek.com

## CERTIFIED REFERENCE MATERIAL

# Certificate of Analysis

chromatographic plus



### FOR LABORATORY USE ONLY-READ SDS PRIOR TO USE.

*This Reference Material is intended for Laboratory Use Only as a standard for the qualitative and/or quantitative determination of the analyte(s) listed.*

Catalog No. : 31453 Lot No.: A0194831  
Description : 8330 Surrogate Mix  
8330 Surrogate Mix 1000 µg/mL, Methanol, 1mL/ampul  
Container Size : 2 mL Pkg Amt: > 1 mL  
Expiration Date : February 29, 2028 Storage: 10°C or colder  
Ship: Ambient

### CERTIFIED VALUES

| Elution Order | Compound           | CAS #    | Lot #    | Purity | Grav. Conc. (weight/volume) | Expanded Uncertainty * (95% C.L.; K=2) |
|---------------|--------------------|----------|----------|--------|-----------------------------|----------------------------------------|
| 1             | 1,2-Dinitrobenzene | 528-29-0 | MKCH6067 | 99%    | 1,004.8 µg/mL               | +/- 56.4585                            |

\* Expanded Uncertainty displayed in same units as Grav. Conc.

Solvent: Methanol  
CAS # 67-56-1  
Purity 99%

## Quality Confirmation Test

**Column:**

250mm x 4.6mm  
Ultra C18 (cat.# 9174575)

**Flow Rate:**

1.0 ml/min.

**Mobile Phase A:**

water:methanol (44:56 V/V)

**Mobile Phase B:****Mobile Phase Composition:**


100%A

**Det. Type:**

Wavelength: 210nm & 254nm

**Inj. Vol**

2µl



This chromatogram represents a general set of testing conditions chosen for product acceptance. For optimal results in your lab, conditions should be adjusted for your specific instrument, method, and application.

Stacey Wanner - Operations Technician I

Date Mixed: 20-Feb-2023

Balance Serial # B442140311

Jennifer Pollino - Operations Tech III - ARM QC

Date Passed: 21-Feb-2023

Manufactured under Restek's ISO 9001:2015  
Registered Quality System  
Certificate #FM 80397



## General Certified Reference Material Notes

### Expiration Notes:

- Expiration date valid for unopened ampul stored in compliance with the recommended conditions.
- Uncertainty, concentration, and expiration of the CRM are based on the unopened product being stored according to the recommended condition found in the storage field.

### Purity Notes:

- Purity and/or chemical identity are determined by one or more of the following techniques: GC/FID, HPLC, GC/ $\mu$ ECD, GC/MS, LC/MS, RI, and/or melting point.
- Compounds with a listed purity of less than 99% have been weight corrected to compensate for impurities and/or salts. A correction factor is used to calculate the amount of compound necessary to achieve the desired concentration of the parent compound in solution.
- Purity of isomeric compounds is reported as the sum of the isomers.
- Purity values are rounded to the nearest whole number.

### Certified Uncertainty Value Notes:

- The uncertainties are determined in accordance with ISO 17034 and Guide 35. The certified expanded uncertainty value includes gravimetric uncertainty, homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty and were combined using the following formula:

$$U_{combined\ uncertainty} = k \sqrt{u_{gravimetric}^2 + u_{homogeneity}^2 + u_{storage\ stability}^2 + u_{shipping\ stability}^2}$$

$k$  is a coverage factor of 2, which gives a level of confidence of approximately 95%.

- The packaged amount is the minimum sample size for which uncertainty is valid. The ampuls are over-filled to ensure that the minimum packaged amount can be sufficiently transferred.

### Manufacturing Notes:

- Concentration is based upon gravimetric preparation using either a balance whose calibration has been verified daily using NIST traceable weights, and/or dilutions with Class A glassware.

### Handling Notes:

- Stability of the unopened product, when stored in compliance with the recommended conditions, is guaranteed through the expiration displayed on the product label and certificate. Contact Restek for additional opened product stability information, with the knowledge/understanding that open product stability is subject to the specific handling and environmental conditions to which the product is exposed. For your convenience Restek supplies deactivated vials with most standards packed in 2mL ampuls. Larger volume deactivated vials are available through Restek as a custom ordered item. Additionally, Restek sells DMDCS for the purpose of glassware deactivation as catalog number 31861, which includes complete instructions.
- If any undissolved material is visible inside the ampul, sonicate the unopened ampul until the material is completely dissolved.



Reagent

---

**8330SurrStkSS\_00268**



110 Benner Circle  
Bellefonte, PA 16823-8812  
Tel: 1-814-353-1300  
Fax: 1-814-353-1309

www.restek.com

## CERTIFIED REFERENCE MATERIAL

# Certificate of Analysis

*chromatographic plus*



### FOR LABORATORY USE ONLY-READ SDS PRIOR TO USE.

*This Reference Material is intended for Laboratory Use Only as a standard for the qualitative and/or quantitative determination of the analyte(s) listed.*

**Catalog No. :** 31453 **Lot No.:** A0194831

**Description :** 8330 Surrogate Mix  
8330 Surrogate Mix 1000 µg/mL, Methanol, 1mL/ampul

**Container Size :** 2 mL **Pkg Amt:** > 1 mL

**Expiration Date :** February 29, 2028 **Storage:** 10°C or colder  
**Ship:** Ambient

### CERTIFIED VALUES

| Elution Order | Compound           | CAS #    | Lot #    | Purity | Grav. Conc. (weight/volume) | Expanded Uncertainty * (95% C.L.; K=2) |
|---------------|--------------------|----------|----------|--------|-----------------------------|----------------------------------------|
| 1             | 1,2-Dinitrobenzene | 528-29-0 | MKCH6067 | 99%    | 1,004.8 µg/mL               | +/- 56.4585                            |

\* Expanded Uncertainty displayed in same units as Grav. Conc.

**Solvent:** Methanol  
**CAS #** 67-56-1  
**Purity** 99%

## Quality Confirmation Test

**Column:**

250mm x 4.6mm  
Ultra C18 (cat.# 9174575)

**Flow Rate:**

1.0 ml/min.

**Mobile Phase A:**

water:methanol (44:56 V/V)

**Mobile Phase B:****Mobile Phase Composition:**


100%A

**Det. Type:**

Wavelength: 210nm & 254nm

**Inj. Vol**

2µl



This chromatogram represents a general set of testing conditions chosen for product acceptance. For optimal results in your lab, conditions should be adjusted for your specific instrument, method, and application.

Stacey Wanner - Operations Technician I

Date Mixed: 20-Feb-2023

Balance Serial # B442140311

Jennifer Pollino - Operations Tech III - ARM QC

Date Passed: 21-Feb-2023

Manufactured under Restek's ISO 9001:2015  
Registered Quality System  
Certificate #FM 80397

## General Certified Reference Material Notes

### Expiration Notes:

- Expiration date valid for unopened ampul stored in compliance with the recommended conditions.
- Uncertainty, concentration, and expiration of the CRM are based on the unopened product being stored according to the recommended condition found in the storage field.

### Purity Notes:

- Purity and/or chemical identity are determined by one or more of the following techniques: GC/FID, HPLC, GC/ $\mu$ ECD, GC/MS, LC/MS, RI, and/or melting point.
- Compounds with a listed purity of less than 99% have been weight corrected to compensate for impurities and/or salts. A correction factor is used to calculate the amount of compound necessary to achieve the desired concentration of the parent compound in solution.
- Purity of isomeric compounds is reported as the sum of the isomers.
- Purity values are rounded to the nearest whole number.

### Certified Uncertainty Value Notes:

- The uncertainties are determined in accordance with ISO 17034 and Guide 35. The certified expanded uncertainty value includes gravimetric uncertainty, homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty and were combined using the following formula:

$$U_{combined\ uncertainty} = k \sqrt{u_{gravimetric}^2 + u_{homogeneity}^2 + u_{storage\ stability}^2 + u_{shipping\ stability}^2}$$

$k$  is a coverage factor of 2, which gives a level of confidence of approximately 95%.

- The packaged amount is the minimum sample size for which uncertainty is valid. The ampuls are over-filled to ensure that the minimum packaged amount can be sufficiently transferred.

### Manufacturing Notes:

- Concentration is based upon gravimetric preparation using either a balance whose calibration has been verified daily using NIST traceable weights, and/or dilutions with Class A glassware.

### Handling Notes:

- Stability of the unopened product, when stored in compliance with the recommended conditions, is guaranteed through the expiration displayed on the product label and certificate. Contact Restek for additional opened product stability information, with the knowledge/understanding that open product stability is subject to the specific handling and environmental conditions to which the product is exposed. For your convenience Restek supplies deactivated vials with most standards packed in 2mL ampuls. Larger volume deactivated vials are available through Restek as a custom ordered item. Additionally, Restek sells DMDCS for the purpose of glassware deactivation as catalog number 31861, which includes complete instructions.
- If any undissolved material is visible inside the ampul, sonicate the unopened ampul until the material is completely dissolved.



Reagent

---

**8330SurrStkSS\_00269**





110 Benner Circle  
Bellefonte, PA 16823-8812  
Tel: 1-814-353-1300  
Fax: 1-814-353-1309

www.restek.com

## CERTIFIED REFERENCE MATERIAL

# Certificate of Analysis

*chromatographic plus*



### FOR LABORATORY USE ONLY-READ SDS PRIOR TO USE.

*This Reference Material is intended for Laboratory Use Only as a standard for the qualitative and/or quantitative determination of the analyte(s) listed.*

**Catalog No. :** 31453 **Lot No.:** A0194831

**Description :** 8330 Surrogate Mix  
8330 Surrogate Mix 1000 µg/mL, Methanol, 1mL/ampul

**Container Size :** 2 mL **Pkg Amt:** > 1 mL

**Expiration Date :** February 29, 2028 **Storage:** 10°C or colder  
**Ship:** Ambient

### CERTIFIED VALUES

| Elution Order | Compound           | CAS #    | Lot #    | Purity | Grav. Conc. (weight/volume) | Expanded Uncertainty * (95% C.L.; K=2) |
|---------------|--------------------|----------|----------|--------|-----------------------------|----------------------------------------|
| 1             | 1,2-Dinitrobenzene | 528-29-0 | MKCH6067 | 99%    | 1,004.8 µg/mL               | +/- 56.4585                            |

\* Expanded Uncertainty displayed in same units as Grav. Conc.

**Solvent:** Methanol  
**CAS #** 67-56-1  
**Purity** 99%

## Quality Confirmation Test

**Column:**

250mm x 4.6mm  
Ultra C18 (cat.# 9174575)

**Flow Rate:**

1.0 ml/min.

**Mobile Phase A:**

water:methanol (44:56 V/V)

**Mobile Phase B:****Mobile Phase Composition:**


100%A

**Det. Type:**

Wavelength: 210nm & 254nm

**Inj. Vol**

2µl



This chromatogram represents a general set of testing conditions chosen for product acceptance. For optimal results in your lab, conditions should be adjusted for your specific instrument, method, and application.

Stacey Wanner - Operations Technician I

Date Mixed: 20-Feb-2023

Balance Serial # B442140311

Jennifer Pollino - Operations Tech III - ARM QC

Date Passed: 21-Feb-2023

Manufactured under Restek's ISO 9001:2015  
Registered Quality System  
Certificate #FM 80397

## General Certified Reference Material Notes

### Expiration Notes:

- Expiration date valid for unopened ampul stored in compliance with the recommended conditions.
- Uncertainty, concentration, and expiration of the CRM are based on the unopened product being stored according to the recommended condition found in the storage field.

### Purity Notes:

- Purity and/or chemical identity are determined by one or more of the following techniques: GC/FID, HPLC, GC/μECD, GC/MS, LC/MS, RI, and/or melting point.
- Compounds with a listed purity of less than 99% have been weight corrected to compensate for impurities and/or salts. A correction factor is used to calculate the amount of compound necessary to achieve the desired concentration of the parent compound in solution.
- Purity of isomeric compounds is reported as the sum of the isomers.
- Purity values are rounded to the nearest whole number.

### Certified Uncertainty Value Notes:

- The uncertainties are determined in accordance with ISO 17034 and Guide 35. The certified expanded uncertainty value includes gravimetric uncertainty, homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty and were combined using the following formula:

$$U_{combined\ uncertainty} = k \sqrt{u_{gravimetric}^2 + u_{homogeneity}^2 + u_{storage\ stability}^2 + u_{shipping\ stability}^2}$$

$k$  is a coverage factor of 2, which gives a level of confidence of approximately 95%.

- The packaged amount is the minimum sample size for which uncertainty is valid. The ampuls are over-filled to ensure that the minimum packaged amount can be sufficiently transferred.

### Manufacturing Notes:

- Concentration is based upon gravimetric preparation using either a balance whose calibration has been verified daily using NIST traceable weights, and/or dilutions with Class A glassware.

### Handling Notes:

- Stability of the unopened product, when stored in compliance with the recommended conditions, is guaranteed through the expiration displayed on the product label and certificate. Contact Restek for additional opened product stability information, with the knowledge/understanding that open product stability is subject to the specific handling and environmental conditions to which the product is exposed. For your convenience Restek supplies deactivated vials with most standards packed in 2mL ampuls. Larger volume deactivated vials are available through Restek as a custom ordered item. Additionally, Restek sells DMDCS for the purpose of glassware deactivation as catalog number 31861, which includes complete instructions.
- If any undissolved material is visible inside the ampul, sonicate the unopened ampul until the material is completely dissolved.



Reagent

---

**8330SurrStkSS\_00270**



110 Benner Circle  
Bellefonte, PA 16823-8812  
Tel: 1-814-353-1300  
Fax: 1-814-353-1309

www.restek.com

## CERTIFIED REFERENCE MATERIAL

# Certificate of Analysis

chromatographic plus



### FOR LABORATORY USE ONLY-READ SDS PRIOR TO USE.

*This Reference Material is intended for Laboratory Use Only as a standard for the qualitative and/or quantitative determination of the analyte(s) listed.*

Catalog No. : 31453 Lot No.: A0194831  
Description : 8330 Surrogate Mix  
8330 Surrogate Mix 1000 µg/mL, Methanol, 1mL/ampul  
Container Size : 2 mL Pkg Amt: > 1 mL  
Expiration Date : February 29, 2028 Storage: 10°C or colder  
Ship: Ambient

### CERTIFIED VALUES

| Elution Order | Compound           | CAS #    | Lot #    | Purity | Grav. Conc. (weight/volume) | Expanded Uncertainty * (95% C.L.; K=2) |
|---------------|--------------------|----------|----------|--------|-----------------------------|----------------------------------------|
| 1             | 1,2-Dinitrobenzene | 528-29-0 | MKCH6067 | 99%    | 1,004.8 µg/mL               | +/- 56.4585                            |

\* Expanded Uncertainty displayed in same units as Grav. Conc.

Solvent: Methanol  
CAS # 67-56-1  
Purity 99%

## Quality Confirmation Test

**Column:**

250mm x 4.6mm  
Ultra C18 (cat.# 9174575)

**Flow Rate:**

1.0 ml/min.

**Mobile Phase A:**

water:methanol (44:56 V/V)

**Mobile Phase B:****Mobile Phase Composition:**


100%A

**Det. Type:**

Wavelength: 210nm & 254nm

**Inj. Vol**

2µl



This chromatogram represents a general set of testing conditions chosen for product acceptance. For optimal results in your lab, conditions should be adjusted for your specific instrument, method, and application.

Stacey Wanner - Operations Technician I

Date Mixed: 20-Feb-2023

Balance Serial # B442140311

Jennifer Pollino - Operations Tech III - ARM QC

Date Passed: 21-Feb-2023

Manufactured under Restek's ISO 9001:2015  
Registered Quality System  
Certificate #FM 80397

## General Certified Reference Material Notes

### Expiration Notes:

- Expiration date valid for unopened ampul stored in compliance with the recommended conditions.
- Uncertainty, concentration, and expiration of the CRM are based on the unopened product being stored according to the recommended condition found in the storage field.

### Purity Notes:

- Purity and/or chemical identity are determined by one or more of the following techniques: GC/FID, HPLC, GC/μECD, GC/MS, LC/MS, RI, and/or melting point.
- Compounds with a listed purity of less than 99% have been weight corrected to compensate for impurities and/or salts. A correction factor is used to calculate the amount of compound necessary to achieve the desired concentration of the parent compound in solution.
- Purity of isomeric compounds is reported as the sum of the isomers.
- Purity values are rounded to the nearest whole number.

### Certified Uncertainty Value Notes:

- The uncertainties are determined in accordance with ISO 17034 and Guide 35. The certified expanded uncertainty value includes gravimetric uncertainty, homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty and were combined using the following formula:

$$U_{combined\ uncertainty} = k \sqrt{u_{gravimetric}^2 + u_{homogeneity}^2 + u_{storage\ stability}^2 + u_{shipping\ stability}^2}$$

$k$  is a coverage factor of 2, which gives a level of confidence of approximately 95%.

- The packaged amount is the minimum sample size for which uncertainty is valid. The ampuls are over-filled to ensure that the minimum packaged amount can be sufficiently transferred.

### Manufacturing Notes:

- Concentration is based upon gravimetric preparation using either a balance whose calibration has been verified daily using NIST traceable weights, and/or dilutions with Class A glassware.

### Handling Notes:

- Stability of the unopened product, when stored in compliance with the recommended conditions, is guaranteed through the expiration displayed on the product label and certificate. Contact Restek for additional opened product stability information, with the knowledge/understanding that open product stability is subject to the specific handling and environmental conditions to which the product is exposed. For your convenience Restek supplies deactivated vials with most standards packed in 2mL ampuls. Larger volume deactivated vials are available through Restek as a custom ordered item. Additionally, Restek sells DMDCS for the purpose of glassware deactivation as catalog number 31861, which includes complete instructions.
- If any undissolved material is visible inside the ampul, sonicate the unopened ampul until the material is completely dissolved.





Reagent

---

**8330SurrStkSS\_00271**



110 Benner Circle  
Bellefonte, PA 16823-8812  
Tel: 1-814-353-1300  
Fax: 1-814-353-1309

www.restek.com

## CERTIFIED REFERENCE MATERIAL

# Certificate of Analysis

chromatographic plus



### FOR LABORATORY USE ONLY-READ SDS PRIOR TO USE.

*This Reference Material is intended for Laboratory Use Only as a standard for the qualitative and/or quantitative determination of the analyte(s) listed.*

**Catalog No. :** 31453 **Lot No.:** A0194831

**Description :** 8330 Surrogate Mix  
8330 Surrogate Mix 1000 µg/mL, Methanol, 1mL/ampul

**Container Size :** 2 mL **Pkg Amt:** > 1 mL

**Expiration Date :** February 29, 2028 **Storage:** 10°C or colder

**Ship:** Ambient

### CERTIFIED VALUES

| Elution Order | Compound           | CAS #    | Lot #    | Purity | Grav. Conc. (weight/volume) | Expanded Uncertainty * (95% C.L.; K=2) |
|---------------|--------------------|----------|----------|--------|-----------------------------|----------------------------------------|
| 1             | 1,2-Dinitrobenzene | 528-29-0 | MKCH6067 | 99%    | 1,004.8 µg/mL               | +/- 56.4585                            |

\* Expanded Uncertainty displayed in same units as Grav. Conc.

**Solvent:** Methanol  
**CAS #** 67-56-1  
**Purity** 99%

## Quality Confirmation Test

**Column:**

250mm x 4.6mm  
Ultra C18 (cat.# 9174575)

**Flow Rate:**

1.0 ml/min.

**Mobile Phase A:**

water:methanol (44:56 V/V)

**Mobile Phase B:****Mobile Phase Composition:**


100%A

**Det. Type:**

Wavelength: 210nm & 254nm

**Inj. Vol**

2µl



This chromatogram represents a general set of testing conditions chosen for product acceptance. For optimal results in your lab, conditions should be adjusted for your specific instrument, method, and application.

Stacey Wanner - Operations Technician I

Date Mixed: 20-Feb-2023

Balance Serial # B442140311

Jennifer Pollino - Operations Tech III - ARM QC

Date Passed: 21-Feb-2023

Manufactured under Restek's ISO 9001:2015  
Registered Quality System  
Certificate #FM 80397

## General Certified Reference Material Notes

### Expiration Notes:

- Expiration date valid for unopened ampul stored in compliance with the recommended conditions.
- Uncertainty, concentration, and expiration of the CRM are based on the unopened product being stored according to the recommended condition found in the storage field.

### Purity Notes:

- Purity and/or chemical identity are determined by one or more of the following techniques: GC/FID, HPLC, GC/ $\mu$ ECD, GC/MS, LC/MS, RI, and/or melting point.
- Compounds with a listed purity of less than 99% have been weight corrected to compensate for impurities and/or salts. A correction factor is used to calculate the amount of compound necessary to achieve the desired concentration of the parent compound in solution.
- Purity of isomeric compounds is reported as the sum of the isomers.
- Purity values are rounded to the nearest whole number.

### Certified Uncertainty Value Notes:

- The uncertainties are determined in accordance with ISO 17034 and Guide 35. The certified expanded uncertainty value includes gravimetric uncertainty, homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty and were combined using the following formula:

$$U_{combined\ uncertainty} = k \sqrt{u_{gravimetric}^2 + u_{homogeneity}^2 + u_{storage\ stability}^2 + u_{shipping\ stability}^2}$$

$k$  is a coverage factor of 2, which gives a level of confidence of approximately 95%.

- The packaged amount is the minimum sample size for which uncertainty is valid. The ampuls are over-filled to ensure that the minimum packaged amount can be sufficiently transferred.

### Manufacturing Notes:

- Concentration is based upon gravimetric preparation using either a balance whose calibration has been verified daily using NIST traceable weights, and/or dilutions with Class A glassware.

### Handling Notes:

- Stability of the unopened product, when stored in compliance with the recommended conditions, is guaranteed through the expiration displayed on the product label and certificate. Contact Restek for additional opened product stability information, with the knowledge/understanding that open product stability is subject to the specific handling and environmental conditions to which the product is exposed. For your convenience Restek supplies deactivated vials with most standards packed in 2mL ampuls. Larger volume deactivated vials are available through Restek as a custom ordered item. Additionally, Restek sells DMDCS for the purpose of glassware deactivation as catalog number 31861, which includes complete instructions.
- If any undissolved material is visible inside the ampul, sonicate the unopened ampul until the material is completely dissolved.



Reagent

---

**8330SurrStkSS\_00272**



110 Benner Circle  
Bellefonte, PA 16823-8812  
Tel: 1-814-353-1300  
Fax: 1-814-353-1309

www.restek.com

## CERTIFIED REFERENCE MATERIAL

# Certificate of Analysis

chromatographic plus



### FOR LABORATORY USE ONLY-READ SDS PRIOR TO USE.

*This Reference Material is intended for Laboratory Use Only as a standard for the qualitative and/or quantitative determination of the analyte(s) listed.*

Catalog No. : 31453 Lot No.: A0194831  
Description : 8330 Surrogate Mix  
8330 Surrogate Mix 1000 µg/mL, Methanol, 1mL/ampul  
Container Size : 2 mL Pkg Amt: > 1 mL  
Expiration Date : February 29, 2028 Storage: 10°C or colder  
Ship: Ambient

### CERTIFIED VALUES

| Elution Order | Compound           | CAS #    | Lot #    | Purity | Grav. Conc. (weight/volume) | Expanded Uncertainty * (95% C.L.; K=2) |
|---------------|--------------------|----------|----------|--------|-----------------------------|----------------------------------------|
| 1             | 1,2-Dinitrobenzene | 528-29-0 | MKCH6067 | 99%    | 1,004.8 µg/mL               | +/- 56.4585                            |

\* Expanded Uncertainty displayed in same units as Grav. Conc.

Solvent: Methanol  
CAS # 67-56-1  
Purity 99%



## Quality Confirmation Test

**Column:**

250mm x 4.6mm  
Ultra C18 (cat.# 9174575)

**Flow Rate:**

1.0 ml/min.

**Mobile Phase A:**

water:methanol (44:56 V/V)

**Mobile Phase B:****Mobile Phase Composition:**


100%A

**Det. Type:**

Wavelength: 210nm & 254nm

**Inj. Vol**

2µl



This chromatogram represents a general set of testing conditions chosen for product acceptance. For optimal results in your lab, conditions should be adjusted for your specific instrument, method, and application.

Stacey Wanner - Operations Technician I

Date Mixed: 20-Feb-2023

Balance Serial # B442140311

Jennifer Pollino - Operations Tech III - ARM QC

Date Passed: 21-Feb-2023

Manufactured under Restek's ISO 9001:2015  
Registered Quality System  
Certificate #FM 80397

## General Certified Reference Material Notes

### Expiration Notes:

- Expiration date valid for unopened ampul stored in compliance with the recommended conditions.
- Uncertainty, concentration, and expiration of the CRM are based on the unopened product being stored according to the recommended condition found in the storage field.

### Purity Notes:

- Purity and/or chemical identity are determined by one or more of the following techniques: GC/FID, HPLC, GC/μECD, GC/MS, LC/MS, RI, and/or melting point.
- Compounds with a listed purity of less than 99% have been weight corrected to compensate for impurities and/or salts. A correction factor is used to calculate the amount of compound necessary to achieve the desired concentration of the parent compound in solution.
- Purity of isomeric compounds is reported as the sum of the isomers.
- Purity values are rounded to the nearest whole number.

### Certified Uncertainty Value Notes:

- The uncertainties are determined in accordance with ISO 17034 and Guide 35. The certified expanded uncertainty value includes gravimetric uncertainty, homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty and were combined using the following formula:

$$U_{combined\ uncertainty} = k \sqrt{u_{gravimetric}^2 + u_{homogeneity}^2 + u_{storage\ stability}^2 + u_{shipping\ stability}^2}$$

$k$  is a coverage factor of 2, which gives a level of confidence of approximately 95%.

- The packaged amount is the minimum sample size for which uncertainty is valid. The ampuls are over-filled to ensure that the minimum packaged amount can be sufficiently transferred.

### Manufacturing Notes:

- Concentration is based upon gravimetric preparation using either a balance whose calibration has been verified daily using NIST traceable weights, and/or dilutions with Class A glassware.

### Handling Notes:

- Stability of the unopened product, when stored in compliance with the recommended conditions, is guaranteed through the expiration displayed on the product label and certificate. Contact Restek for additional opened product stability information, with the knowledge/understanding that open product stability is subject to the specific handling and environmental conditions to which the product is exposed. For your convenience Restek supplies deactivated vials with most standards packed in 2mL ampuls. Larger volume deactivated vials are available through Restek as a custom ordered item. Additionally, Restek sells DMDCS for the purpose of glassware deactivation as catalog number 31861, which includes complete instructions.
- If any undissolved material is visible inside the ampul, sonicate the unopened ampul until the material is completely dissolved.



Reagent

---

**8330SurrStkSS\_00273**



110 Benner Circle  
Bellefonte, PA 16823-8812  
Tel: 1-814-353-1300  
Fax: 1-814-353-1309

www.restek.com

## CERTIFIED REFERENCE MATERIAL

# Certificate of Analysis

*chromatographic plus*



### FOR LABORATORY USE ONLY-READ SDS PRIOR TO USE.

*This Reference Material is intended for Laboratory Use Only as a standard for the qualitative and/or quantitative determination of the analyte(s) listed.*

**Catalog No. :** 31453 **Lot No.:** A0197062

**Description :** 8330 Surrogate Mix  
8330 Surrogate Mix 1000 µg/mL, Methanol, 1mL/ampul

**Container Size :** 2 mL **Pkg Amt:** > 1 mL

**Expiration Date :** April 30, 2028 **Storage:** 10°C or colder

**Ship:** Ambient

### CERTIFIED VALUES

| Elution Order | Compound           | CAS #    | Lot #    | Purity | Grav. Conc. (weight/volume) | Expanded Uncertainty * (95% C.L.; K=2) |
|---------------|--------------------|----------|----------|--------|-----------------------------|----------------------------------------|
| 1             | 1,2-Dinitrobenzene | 528-29-0 | STBK9056 | 99%    | 1,001.6 µg/mL               | +/- 56.2801                            |

\* Expanded Uncertainty displayed in same units as Grav. Conc.

**Solvent:** Methanol  
**CAS #** 67-56-1  
**Purity** 99%

## Quality Confirmation Test

**Column:**

250mm x 4.6mm  
Ultra C18 (cat.# 9174575)

**Flow Rate:**

1.0 ml/min.

**Mobile Phase A:**

water:methanol (44:56 V/V)

**Mobile Phase B:****Mobile Phase Composition:**


100%A

**Det. Type:**

Wavelength: 210nm & 254nm

**Inj. Vol**

2.0µl



This chromatogram represents a general set of testing conditions chosen for product acceptance. For optimal results in your lab, conditions should be adjusted for your specific instrument, method, and application.

Penelope Riglin - Operations Tech I

Date Mixed: 17-Apr-2023

Balance Serial # 1128360905

Jennifer Pollino - Operations Tech III - ARM QC

Date Passed: 19-Apr-2023

Manufactured under Restek's ISO 9001:2015  
Registered Quality System  
Certificate #FM 80397

## General Certified Reference Material Notes

### Expiration Notes:

- Expiration date valid for unopened ampul stored in compliance with the recommended conditions.
- Uncertainty, concentration, and expiration of the CRM are based on the unopened product being stored according to the recommended condition found in the storage field.

### Purity Notes:

- Purity and/or chemical identity are determined by one or more of the following techniques: GC/FID, HPLC, GC/μECD, GC/MS, LC/MS, RI, and/or melting point.
- Compounds with a listed purity of less than 99% have been weight corrected to compensate for impurities and/or salts. A correction factor is used to calculate the amount of compound necessary to achieve the desired concentration of the parent compound in solution.
- Purity of isomeric compounds is reported as the sum of the isomers.
- Purity values are rounded to the nearest whole number.

### Certified Uncertainty Value Notes:

- The uncertainties are determined in accordance with ISO 17034 and Guide 35. The certified expanded uncertainty value includes gravimetric uncertainty, homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty and were combined using the following formula:

$$U_{combined\ uncertainty} = k \sqrt{u_{gravimetric}^2 + u_{homogeneity}^2 + u_{storage\ stability}^2 + u_{shipping\ stability}^2}$$

$k$  is a coverage factor of 2, which gives a level of confidence of approximately 95%.

- The packaged amount is the minimum sample size for which uncertainty is valid. The ampuls are over-filled to ensure that the minimum packaged amount can be sufficiently transferred.

### Manufacturing Notes:

- Concentration is based upon gravimetric preparation using either a balance whose calibration has been verified daily using NIST traceable weights, and/or dilutions with Class A glassware.

### Handling Notes:

- Stability of the unopened product, when stored in compliance with the recommended conditions, is guaranteed through the expiration displayed on the product label and certificate. Contact Restek for additional opened product stability information, with the knowledge/understanding that open product stability is subject to the specific handling and environmental conditions to which the product is exposed. For your convenience Restek supplies deactivated vials with most standards packed in 2mL ampuls. Larger volume deactivated vials are available through Restek as a custom ordered item. Additionally, Restek sells DMDCS for the purpose of glassware deactivation as catalog number 31861, which includes complete instructions.
- If any undissolved material is visible inside the ampul, sonicate the unopened ampul until the material is completely dissolved.

Reagent

---

**8330SurrStock\_00172**



# CERTIFICATE OF ANALYSIS

**Catalog No:** M-8330-SS

**Description:** 1,2-Dinitrobenzene

**Lot:** 219051500

**Solvent:** Methanol

**Hazards:** Refer to SDS for complete safety information

**Date Certified:** May 22, 2019

**Expiration:** May 22, 2029

**Sample Size:** 1 mL

**Components:** 1

**Storage Condition:** Ambient (>5 °C)



**Signal Word:** Danger

**Certified Reference Material**



| Component          | CAS #    | Purity % | Prepared Concentration <sup>2</sup> | Certified Analyte Concentration <sup>1</sup> |
|--------------------|----------|----------|-------------------------------------|----------------------------------------------|
|                    |          | (GC/FID) | (µg/mL)                             | (µg/mL)                                      |
| 1,2-Dinitrobenzene | 528-29-0 | 100.0    | 1002                                | 1002                                         |

This Certified Reference Material was verified in accordance with ISO/IEC 17025

A product with a suffix (-1A, -2B, etc. or -01, -02, etc.) on its lot number has had its expiration date extended and is identical to the same lot number without the suffix.

<sup>2</sup> All weights are traceable through NIST, Test No. 684/289871-17

<sup>1</sup> Certified Analyte Concentration = Purity x Prepared Concentration.

The Uncertainty associated with the certified concentration reported on this certificate is  $\pm 2.4\%$ . This value is the combined expanded uncertainty and represents an estimated standard deviation equal to the positive square root of the total variation of the uncertainty of components. A normal distribution is assumed and a coverage factor of K=2 is chosen using approximately a 95% confidence level.

Labels and certificates follow U.S. Conventions in reporting numerical values: A comma (,) is used to separate units of one-thousand or greater. A period (.) is used as a decimal place marker.

The information on this certificate may not be reproduced without the express permission of the manufacturer. See reverse side for additional information

Hazard Information: Please refer to the SDS for information regarding the hazards associated with using this material.

This product was prepared according to in-house procedures and is guaranteed to be homogeneous.

Certified By:

Larry Decker, Organic QC Manager

Reagent

---

**Alk stk std\_00023**



**SPEXertificate**<sup>®</sup>  
*Certificate of Reference Material*



**Catalog Number:** VTACO-3 **Lot No.:** 2-106BT  
**Description:** Custom Standard  
**Matrix:** H<sub>2</sub>O

This Certified Reference Material, CRM, is intended primarily for use as a calibration standard or quality control standard for conductance meters and conductivity cells used in water purity measurements. It can be employed in USEPA, ASTM and other methods relevant to the certified properties listed below.

The CRM is prepared gravimetrically using high purity sodium carbonate (Na<sub>2</sub>CO<sub>3</sub>) Lot# H04155 and ASTM Type I water, using Class A laboratory ware to give precise concentration.

**Certified Value:**

Analyte: Na<sub>2</sub>CO<sub>3</sub>  
Alkalinity: 50,000 µg/mL as CaCO<sub>3</sub>  
Uncertainty: ± 500 µg/mL  
Traceable to: Standard made from NIST SRM # 351a

Balances are calibrated regularly with weight sets traceable to NIST#s 32856, 32867. This CRM is guaranteed stable to ±1% of the certified concentration inclusive of uncertainty of measurements and other components such as homogeneity, long and short-term stability for a period of one year from the date of certification. This guarantee is valid only when the material is kept tightly capped and stored under ambient laboratory conditions.

Date of Certification: OCT - - 2022

Certifying Officer: Katherine Cullinan  
Katherine Cullinan, QC Manager

Page 1 of 2  
Rev. 0

# Report of Certification

This Certified Reference Material (CRM) has been prepared and certified under an ISO 9001 (certified by DQS), ISO/IEC 17025 (accredited by A2LA) and ISO 17034 (accredited by A2LA) quality system consistent with the following guides:

- ISO 9001: Quality management systems – Requirements
- ISO/IEC 17025: General requirements for the competence of testing and calibration laboratories
- ISO 17034: General requirements for the competence of reference material producers
- ISO Guide 30: Reference Materials – Selected terms and definitions
- ISO Guide 31: Reference Materials – Contents of certificates, labels, and accompanying documentation
- ISO Guide 35: Reference Materials – Guidance for characterization and assessment of homogeneity and stability
- Guide to the Expression of Uncertainty in Measurement, 2008
- EURACHEM/CITAC Guide: Qualifying Uncertainty in Analytical Measurement – Third Edition
- NIST Technical Note 1297

## Material Source:

All analytes and matrix materials are obtained and verified by Spex CertiPrep from pre-qualified vendors as per ISO 9001, ISO/IEC 17025 and ISO 17034 guidelines. Vendor identifications are proprietary; however, sources of all materials used in the preparation and testing of Spex CertiPrep CRMs are tracked and documented. For further assistance, please contact Sales Support at USMet-CRMSales@antylia.com.

## Instructions for Use:

Primary usage of this CRM is in neat form or diluted serially with matrix of a purity at or greater than the purity of the original matrix solution. If dilution is required, the diluent must be compatible with all certified analytes and contain stabilizers appropriate for the period of intended use. The CRM can also be used as a spike or with a spike, again with appropriate compatibility considerations. All solutions should be thoroughly mixed, by shaking, prior to use and never pipetted directly from the bottle. Do not return excess solution to the bottle. All surfaces that come in contact with the solution must be thoroughly cleaned and leached prior to use. Dilutions should be performed only with Class A volumetric glassware. See SDS for health and safety information.

## Method of Preparation:

Clean laboratory procedures and techniques have been used throughout the preparation. All materials, equipment, analytical instrumentation and personnel have been qualified prior to use. The highest purity acids applicable, 18 megohm, double deionized water, acid-leached triple-rinsed bottles (where appropriate), and Class A/calibrated volumetrics have been used in all preparations.

## Homogeneity:

The homogeneity of the CRM has been confirmed by procedures consistent with ISO/IEC 17025, ISO 17034, and ASTM D6362-98 Appendix X2. Random, replicate samples of the final, packaged material have been analyzed to prove homogeneity in accordance with our internal procedure 4600-HOMOGEN-1A. Since the product is highly homogeneous, any sample size taken for analysis would be within the uncertainty budget. This is consistent with the intended use of the CRM.

## Statistical Estimator and Confidence Limits:

The certified value 'X' listed on the reverse of this document is at the 95% level of confidence and can be expressed as:

- $X = x \pm U$  where  $X$  = certified value,  $U$  = expanded uncertainty,  $x$  = property value
- $U = k u_c$  where  $k = 2$  is the coverage factor at the 95% confidence level
- $u_c$  = combined standard uncertainty obtained by combining the individual element standard uncertainty components  $u_i$  and  $u_c = \sqrt{\sum u_i^2}$

## Certification Report:

All certified values reported were derived from the Certification Report, Spex CertiPrep's traceability documentation, identified by the lot number of this CRM. During the stated period of validity, the purchaser will be notified if this product is recalled due to any significant changes in the stability of the solution. For further assistance, please contact Sales Support at USMet-CRMSales@antylia.com.

## Legal Notice:

Spex CertiPrep reference materials are not for any cosmetic, drug or household application and are to be used only by qualified individuals who are trained in appropriate procedures. No claims against Spex CertiPrep, LLC. of any kind whatsoever, whether based on breach of warranty, alleged negligence, or otherwise, with respect to this Reference Material shall be greater than the purchase price. In no event shall Spex CertiPrep, LLC. be liable for any loss of profits or any incidental, special, or consequential damages.

Reagent

---

**IC BR ICV\_00024**

# Certificate of Analysis

## Bromide Standard, 1000 ppm Br<sup>-</sup>

**Lot Number:** 1204B86

**Product Number:** 1180

**Manufacture Date:** APR 26, 2022

**Expiration Date:** OCT 2023

The certified value reported is the prepared value based upon the method of preparation of the material. The uncertainty in the prepared value is based upon the volumetric method of preparation.

| Name           | CAS#      | Grade           |
|----------------|-----------|-----------------|
| Water          | 7732-18-5 | ACS/ASTM/USP/EP |
| Sodium Bromide | 7647-15-6 | High Purity     |

| Test         | Specification    | Result   |
|--------------|------------------|----------|
| Appearance   | Colorless liquid | Passed   |
| Bromide (Br) | 995-1005 ppm     | 1000 ppm |

| Specification                                               | Reference           |
|-------------------------------------------------------------|---------------------|
| Bromide Solution, Standard (1 mL = 1 mg Br <sup>-</sup> )   | ASTM (D 3869 D)     |
| Standard Bromide Solution, 1000 mg/L                        | APHA (4110 B)       |
| Bromide Stock Solution (1.00 mL = 1.00 mg Br <sup>-</sup> ) | EPA (SW-846) (9056) |
| Sodium Bromide Standard Solution, 1000 mg/L                 | ASTM (D 1246)       |
| Bromide Stock Solution (1.00 mL = 1.00 mg Br <sup>-</sup> ) | ASTM (D 4327)       |

Volumetric glassware complies with Class A tolerance requirements of ASTM E 288 and NIST Circular 434; it is calibrated before first use and recalibrated regularly in accordance with ASTM E 542 and NIST Procedure NBSIR 74-461. Balances are calibrated regularly with weights certified traceable to the NIST national mass standard. Thermometers and temperature probes are calibrated before first use and recalibrated regularly with a thermometer traceable to NIST standards. All products are prepared according to master documents that assure manufacture according to validated methods. Batch records document raw material traceability and production and testing history for each lot manufactured.

| Part Number | Size / Package Type | Shelf Life (Unopened Container) |
|-------------|---------------------|---------------------------------|
| 1180-16     | 500 mL natural poly | 18 months                       |

**Recommended Storage:** 15°C - 30°C (59°F - 86°F)



Jason Thomas (04/26/2022)

Operations Manager

This Certificate of Analysis is designed to comply with ISO Guide 31 "Reference Materials -- Contents of Certificates and Labels."

This test report shall not be reproduced, except in full, without the written approval of Ricca Chemical Company.

Reagent

---

**IC FL cal\_00023**

# Certificate of Analysis

**Fluoride Standard, 1000 ppm F<sup>-</sup>**

**Lot Number:** 2207H13

**Product Number:** 3173

**Manufacture Date:** JUL 29, 2022

**Expiration Date:** JAN 2024

The certified value reported is the prepared value based upon the method of preparation of the material. The uncertainty in the prepared value is based upon the volumetric method of preparation.

The concentration is confirmed by Fluoride ISE and is certified traceable to NIST SRM 2203.

| Name            | CAS#      | Grade           |
|-----------------|-----------|-----------------|
| Water           | 7732-18-5 | ACS/ASTM/USP/EP |
| Sodium Fluoride | 7681-49-4 | High Purity     |

| Test         | Specification    | Result   |
|--------------|------------------|----------|
| Appearance   | Colorless liquid | Passed   |
| Fluoride (F) | 995-1005 ppm     | 1005 ppm |

| Specification                                                    | Reference           |
|------------------------------------------------------------------|---------------------|
| Fluoride Solution, Stock (1.00 mL = 1.00 mg F)                   | ASTM (D 5542)       |
| Fluoride Stock Solution (1.00 mL = 1.00 mg F)                    | EPA (SW-846) (9056) |
| Fluoride Calibration Stock Solution (1,000 mg/L F <sup>-</sup> ) | EPA (SW-846) (9214) |
| Stock Solution, 1.0 mL = 1.0 mg F                                | EPA (340.3)         |
| Fluoride Solution, Stock (1.00 mL = 1.00 mg F)                   | ASTM (D 5996)       |
| Fluoride Stock Solution (1.00 mL = 1.00 mg F)                    | ASTM (D 4327)       |
| Fluoride Stock Standard Solution (1 mg of F in 1 mL)             | ACS (N/A)           |

Volumetric glassware complies with Class A tolerance requirements of ASTM E 288 and NIST Circular 434; it is calibrated before first use and recalibrated regularly in accordance with ASTM E 542 and NIST Procedure NBSIR 74-461. Balances are calibrated regularly with weights certified traceable to the NIST national mass standard. Thermometers and temperature probes are calibrated before first use and recalibrated regularly with a thermometer traceable to NIST standards. All products are prepared according to master documents that assure manufacture according to validated methods. Batch records document raw material traceability and production and testing history for each lot manufactured.

| Part Number | Size / Package Type | Shelf Life (Unopened Container) |
|-------------|---------------------|---------------------------------|
| 3173-16     | 500 mL natural poly | 18 months                       |
| 3173-32     | 1 L natural poly    | 18 months                       |

**Recommended Storage:** 15°C - 30°C (59°F - 86°F)



Myrlande Gilles (07/29/2022)

Quality Control

This Certificate of Analysis is designed to comply with ISO Guide 31 "Reference Materials -- Contents of Certificates and Labels."

This test report shall not be reproduced, except in full, without the written approval of Ricca Chemical Company.



Reagent

---

**IC N03 cal\_00028**

# Certificate of Analysis

**Nitrate Nitrogen Standard, 1000 ppm N (4427 ppm NO<sub>3</sub>)**

**Lot Number:** 2205B73

**Product Number:** 5459

**Manufacture Date:** MAY 11, 2022

**Expiration Date:** NOV 2023

The certified value reported is the prepared value based upon the method of preparation of the material. The uncertainty in the prepared value is based upon the volumetric method of preparation.

| Name              | CAS#      | Grade           |
|-------------------|-----------|-----------------|
| Water             | 7732-18-5 | ACS/ASTM/USP/EP |
| Potassium Nitrate | 7757-79-1 | High Purity     |
| Chloroform        | 67-66-3   |                 |

| Test         | Specification    | Result   |
|--------------|------------------|----------|
| Appearance   | Colorless liquid | Passed   |
| Nitrogen (N) | 995-1005 ppm     | 1000 ppm |

| Specification                                                | Reference       |
|--------------------------------------------------------------|-----------------|
| Nitrate Solution, Stock (1.0 mL = 1.0 mg NO <sub>3</sub> -N) | ASTM (D 3867 A) |
| Nitrate Solution, Stock (1.0 mL = 1.0 mg NO <sub>3</sub> -N) | ASTM (D 3867 B) |
| Stock Nitrate Solution: 1 mL = 1.0 mg NO <sub>3</sub> -N     | EPA (353.2)     |
| Stock Nitrate Solution: 1.0 mL = 1.00 mg NO <sub>3</sub> -N  | EPA (353.3)     |

Volumetric glassware complies with Class A tolerance requirements of ASTM E 288 and NIST Circular 434; it is calibrated before first use and recalibrated regularly in accordance with ASTM E 542 and NIST Procedure NBSIR 74-461. Balances are calibrated regularly with weights certified traceable to the NIST national mass standard. Thermometers and temperature probes are calibrated before first use and recalibrated regularly with a thermometer traceable to NIST standards. All products are prepared according to master documents that assure manufacture according to validated methods. Batch records document raw material traceability and production and testing history for each lot manufactured.

| Part Number | Size / Package Type | Shelf Life (Unopened Container) |
|-------------|---------------------|---------------------------------|
| 5459-16     | 500 mL natural poly | 18 months                       |

**Recommended Storage:** 15°C - 30°C (59°F - 86°F)



Myrlande Gilles (05/11/2022)

Quality Control

This Certificate of Analysis is designed to comply with ISO Guide 31 "Reference Materials -- Contents of Certificates and Labels."

This test report shall not be reproduced, except in full, without the written approval of Ricca Chemical Company.

Reagent

---

**IC NO2 Cal\_00007**

# Certificate of Analysis

**Nitrite Nitrogen Standard, 1000 ppm N (3285 ppm NO<sub>2</sub>)**

**Lot Number:** 2212826

**Product Number:** 5461

**Manufacture Date:** DEC 08, 2022

**Expiration Date:** JUN 2023

| Name              | CAS#      | Grade           |
|-------------------|-----------|-----------------|
| Water             | 7732-18-5 | ACS/ASTM/USP/EP |
| Potassium Nitrite | 7758-09-0 | ACS             |
| Chloroform        | 67-66-3   |                 |

| Test                               | Specification    | Result     | NIST SRM# |
|------------------------------------|------------------|------------|-----------|
| Appearance                         | Colorless liquid | Passed     |           |
| Assay (vs. Potassium Permanganate) | 995-1005 ppm N   | 1005 ppm N | 8040      |

| Specification                                                | Reference       |
|--------------------------------------------------------------|-----------------|
| Nitrite Solution, Stock (1.0 mL = 1.0 mg NO <sub>2</sub> -N) | ASTM (D 3867 A) |
| Nitrite Solution, Stock (1.0 mL = 1.0 mg NO <sub>2</sub> -N) | ASTM (D 3867 B) |
| Stock Nitrite Solution: 1 mL = 1.0 mg NO <sub>2</sub> -N     | EPA (353.2)     |
| Stock Nitrite Solution: 1.0 mL = 1.00 mg NO <sub>2</sub> -N  | EPA (353.3)     |

Volumetric glassware complies with Class A tolerance requirements of ASTM E 288 and NIST Circular 434; it is calibrated before first use and recalibrated regularly in accordance with ASTM E 542 and NIST Procedure NBSIR 74-461. Balances are calibrated regularly with weights certified traceable to the NIST national mass standard. Thermometers and temperature probes are calibrated before first use and recalibrated regularly with a thermometer traceable to NIST standards. All products are prepared according to master documents that assure manufacture according to validated methods. Batch records document raw material traceability and production and testing history for each lot manufactured.

| Part Number | Size / Package Type | Shelf Life (Unopened Container) |
|-------------|---------------------|---------------------------------|
| 5461-16     | 500 mL natural poly | 6 months                        |
| 5461-4      | 120 mL natural poly | 6 months                        |

**Recommended Storage:** 2°C - 8°C (36°F - 46°F)



Myrlande Gilles (12/08/2022)

Quality Control

This Certificate of Analysis is designed to comply with ISO Guide 31 "Reference Materials -- Contents of Certificates and Labels."

This test report shall not be reproduced, except in full, without the written approval of Ricca Chemical Company.

Reagent

---

**IC NO2 ICV\_00022**



A Waters Company

Certified Reference Material

▪ **Certificate of Analysis** ▪

**Product:** 1000 mg/L Nitrite as N (NO<sub>2</sub>-N)  
**Catalog Number:** 053-125mL, 990-500mL  
**Lot No.** 261221m  
**Starting Material:** Sodium Nitrite (NaNO<sub>2</sub>)  
**Matrix:** 18 megohm deionized water  
**Density:** 1.0022 ± 0.0005 g/mL . 21.3 °C and 740 mm Hg  
**Verification Method:** Ion Chromatography  
**Certificate Issue Date:** December 21, 2021  
**Expiration Date:** December 14, 2023  
**Revision Number:** Original

**CERTIFICATION**

| Parameter                         | Certified Value <sup>1</sup> | Uncertainty <sup>2</sup> | NIST Traceability       |            |
|-----------------------------------|------------------------------|--------------------------|-------------------------|------------|
|                                   | mg/L                         | %                        | SRM Number <sup>3</sup> | Recovery % |
| Nitrite as N (NO <sub>2</sub> -N) | 1000                         | 0.948                    | -                       | -          |

## Certified Reference Material

# • Certificate of Analysis •

1. The **Certified Value** is the actual gravimetric/volumetric "made-to" concentration confirmed by ERA analytical verification. The certified value is monitored and the purchasers will be notified of any significant changes resulting in recertification or withdrawal of this certified reference material during the period of validity of this certificate.

2. The **Uncertainty** represents an expanded uncertainty and approximates a 95% confidence interval. The uncertainty is based on the characterization, homogeneity and stability characteristics of the product, multiplied by a coverage factor ( $k=2$ ). The uncertainty applies to the product as supplied and does not take into account any required or optional dilution and/or preparations the laboratory may perform while using this product. The formula used to calculate the expanded uncertainty is:

$$U_{\text{expanded}} = k * \text{SQRT}((U_{\text{char}})^2 + (U_{\text{homogen}})^2 + (U_{\text{LTS}})^2 + (U_{\text{STS}})^2 + (U_{\text{RSS}})^2)$$

Where:

$U_{\text{expanded}}$  = Expanded uncertainty.

$k$  = Coverage factor.

$U_{\text{char}}$  = Combined standard uncertainty of the manufacturing and/or analytical verification assessment.

$U_{\text{homogen}}$  = Standard uncertainty of the homogeneity assessment.

$U_{\text{LTS}}$  = Standard uncertainty associated with long-term stability.

$U_{\text{STS}}$  = Standard uncertainty associated with short-term (transport) stability.

$U_{\text{RSS}}$  = Standard uncertainty associated with repeated sampling of the product (where permitted by product use instructions).

3. Where NIST Standard Reference Materials (SRMs) are available, each analyte has been analytically traced to the NIST SRM listed. **Analytical Traceability Recovery (%)** =  $[(\% \text{ recovery ERA certified reference material})/(\% \text{ recovery NIST SRM})]*100$

The traceability data shown were compiled by analyzing this ERA certified reference material and/or its associated stock solution(s) against the applicable NIST SRMs.

4. **Metrological Traceability.** This certified reference material is metrologically traceable to NIST mass reference materials through an unbroken chain of comparisons.

5. **Storage:** 20-25°C

6. **Intended Use:** This standard is intended to be used to calibrate your analytical process and/or as a quality control check of the entire process for the analytes/matrix included in the standard.

7. **Minimum Sample Size:** ERA suggests that when subsampling this product prior to analysis, you use a minimum sample size of at least 1 mL. Using a smaller sample size may invalidate the assigned value and/or uncertainty shown.

8. **Repeat Sampling:** Repeated Sampling of this product is permitted, provided minimum sample sizes and storage instructions are adhered to.

9. **Safety:** ERA products may be hazardous and are intended for use by professional laboratory personnel trained in the competent handling of such materials. Responsibility for the safe use of these products rests entirely with the buyer and/or user. Safety Data Sheets (SDS) for all ERA products are available through our website: [www.eraqc.com](http://www.eraqc.com)

If you have any questions or need technical assistance, please call ERA technical assistance at 1-800-372-0122 or send an email to [info@eraqc.com](mailto:info@eraqc.com).

Certifying Officer

Brian Miller



Quality Officer

Matthew Seebeck



ISO 17034:2016



ISO/IEC 17025:2017



Reagent

---

**IC NO3 ICV\_00019**





A Waters Company

Certified Reference Material

## ▪ Certificate of Analysis ▪

**Product:** 1000 mg/L Nitrate as N ( $\text{NO}_3\text{-N}$ )  
**Catalog Number:** 052-125mL, 991-500mL  
**Lot No.** 011221m  
**Starting Material:** Potassium Nitrate ( $\text{KNO}_3$ )  
**Matrix:** 18 megohm deionized water  
**Density:**  $1.0029 \pm 0.0005$  g/mL 20.2 °C and 749 mm Hg  
**Verification Method:** Ion Chromatography  
**Certificate Issue Date:** December 21, 2021  
**Expiration Date:** December 1, 2023  
**Revision Number:** Original

### CERTIFICATION

| Parameter                               | Certified Value <sup>1</sup> | Uncertainty <sup>2</sup> | NIST Traceability       |            |
|-----------------------------------------|------------------------------|--------------------------|-------------------------|------------|
|                                         | mg/L                         | %                        | SRM Number <sup>3</sup> | Recovery % |
| Nitrate as N ( $\text{NO}_3\text{-N}$ ) | 1000                         | 1.57                     | 3185                    | 101        |

## Certified Reference Material

# • Certificate of Analysis •

1. The **Certified Value** is the actual gravimetric/volumetric "made-to" concentration confirmed by ERA analytical verification. The certified value is monitored and the purchasers will be notified of any significant changes resulting in recertification or withdrawal of this certified reference material during the period of validity of this certificate.

2. The **Uncertainty** represents an expanded uncertainty and approximates a 95% confidence interval. The uncertainty is based on the characterization, homogeneity and stability characteristics of the product, multiplied by a coverage factor ( $k=2$ ). The uncertainty applies to the product as supplied and does not take into account any required or optional dilution and/or preparations the laboratory may perform while using this product. The formula used to calculate the expanded uncertainty is:

$$U_{\text{expanded}} = k * \text{SQRT}((U_{\text{char}})^2 + (U_{\text{homogen}})^2 + (U_{\text{LTS}})^2 + (U_{\text{STS}})^2 + (U_{\text{RSS}})^2)$$

Where:

$U_{\text{expanded}}$  = Expanded uncertainty.

$k$  = Coverage factor.

$U_{\text{char}}$  = Combined standard uncertainty of the manufacturing and/or analytical verification assessment.

$U_{\text{homogen}}$  = Standard uncertainty of the homogeneity assessment.

$U_{\text{LTS}}$  = Standard uncertainty associated with long-term stability.

$U_{\text{STS}}$  = Standard uncertainty associated with short-term (transport) stability.

$U_{\text{RSS}}$  = Standard uncertainty associated with repeated sampling of the product (where permitted by product use instructions).

3. Where NIST Standard Reference Materials (SRMs) are available, each analyte has been analytically traced to the NIST SRM listed. **Analytical Traceability Recovery (%)** =  $\{(\% \text{ recovery ERA certified reference material})/(\% \text{ recovery NIST SRM})\} * 100$

The traceability data shown were compiled by analyzing this ERA certified reference material and/or its associated stock solution(s) against the applicable NIST SRMs.

4. **Metrological Traceability.** This certified reference material is metrologically traceable to NIST mass reference materials through an unbroken chain of comparisons.

5. **Storage:** 20-25°C

6. **Intended Use:** This standard is intended to be used to calibrate your analytical process and/or as a quality control check of the entire process for the analytes/matrix included in the standard.

7. **Minimum Sample Size:** ERA suggests that when subsampling this product prior to analysis, you use a minimum sample size of at least 1 mL. Using a smaller sample size may invalidate the assigned value and/or uncertainty shown.

8. **Repeat Sampling:** Repeated Sampling of this product is permitted, provided minimum sample sizes and storage instructions are adhered to.

9. **Safety:** ERA products may be hazardous and are intended for use by professional laboratory personnel trained in the competent handling of such materials. Responsibility for the safe use of these products rests entirely with the buyer and/or user. Safety Data Sheets (SDS) for all ERA products are available through our website: [www.eraqc.com](http://www.eraqc.com)

If you have any questions or need technical assistance, please call ERA technical assistance at 1-800-372-0122 or send an email to [info@eraqc.com](mailto:info@eraqc.com).

Certifying Officer  
Brian Miller

Quality Officer  
Matthew Seebeck



Reagent

---

**IC P04 cal\_00030**

# Certificate of Analysis

## Phosphorus AA Standard, 1000 ppm P in H<sub>2</sub>O

**Lot Number:** 4208152**Product Number:** AP1KW**Manufacture Date:** AUG 02, 2022**Expiration Date:** JUL 2024

This is a single element solution that was prepared volumetrically to contain the certified value reported. The uncertainty associated with the certified value is the sum of the estimated errors due to the purity of the raw material, the volumetric preparation of the solution, and transpiration of the solution through the container wall.

The final solution concentration is confirmed by AA, ICP, or ICP-MS, and is traceable to NIST Standard Reference Material 3139.

| Name                          | CAS#      | Grade           |
|-------------------------------|-----------|-----------------|
| Water                         | 7732-18-5 | ACS/ASTM/USP/EP |
| Ammonium Dihydrogen Phosphate | 7722-76-1 | High Purity     |

| Test           | Specification    | Result   | NIST SRM# |
|----------------|------------------|----------|-----------|
| Appearance     | Colorless liquid | Passed   |           |
| Phosphorus (P) | 995-1005 ppm     | 1000 ppm | 3139      |

Volumetric glassware complies with Class A tolerance requirements of ASTM E 288 and NIST Circular 434; it is calibrated before first use and recalibrated regularly in accordance with ASTM E 542 and NIST Procedure NBSIR 74-461. Balances are calibrated regularly with weights certified traceable to the NIST national mass standard. Thermometers and temperature probes are calibrated before first use and recalibrated regularly with a thermometer traceable to NIST standards. All products are prepared according to master documents that assure manufacture according to validated methods. Batch records document raw material traceability and production and testing history for each lot manufactured.

| Part Number | Size / Package Type | Shelf Life (Unopened Container) |
|-------------|---------------------|---------------------------------|
| AP1KW-500   | 500 mL natural poly | 24 months                       |

**Recommended Storage:** 15°C - 30°C (59°F - 86°F)

Paul Brandon (08/02/2022)

Production Manager

This Certificate of Analysis is designed to comply with ISO Guide 31 "Reference Materials -- Contents of Certificates and Labels."

This test report shall not be reproduced, except in full, without the written approval of Ricca Chemical Company.

Reagent

---

**IC SO4 Cal\_00008**

# SPEXertificate<sup>®</sup>

## Certificate of Reference Material

**Catalog Number:** AS-SO49-2X

**Lot No.** 5-160SO4-2X

**Description:** 1000 µg/mL Sulfate

**Matrix:** H<sub>2</sub>O

This **Ion Chromatography** Certified Reference Material, CRM, is intended primarily for use as a calibration standard or quality control standard for ion chromatography instrumentation. It can be employed in USEPA, ASTM and other methods relevant to the certified properties listed below.

**Certified Value:** 1000 µg/mL ±5 µg/mL

**Certified Value is Traceable to:** 3181\*

\* - indicates NIST SRM    † - indicates SPEX CertiPrep CRM (when NIST SRM is not available)    ‡ - prepared gravimetrically

The CRM is prepared gravimetrically using high purity Potassium Sulfate, Lot# AN07201A. The certified value listed is the average of values obtained by classical wet assay and ion chromatography analysis.

Refer to side 2 for details of measurement uncertainties.

**Classical Wet Assay:** 999 µg/mL

**Method:** Precipitation using Barium Chloride. Filter, ignite, and weigh as BaSO<sub>4</sub>.

**Instrumental Analysis by Ion Chromatography:** 1000 µg/mL

**Uncertified Properties**

**Trace Ionic Impurities in the Actual Solution via IC Analysis:**

| Element         | µg/mL  | Element                       | µg/mL |
|-----------------|--------|-------------------------------|-------|
| Br <sup>-</sup> | <0.01  | NO <sub>2</sub> <sup>-</sup>  | <0.01 |
| Cl <sup>-</sup> | <0.03  | NO <sub>3</sub> <sup>-</sup>  | <0.02 |
| F <sup>-</sup>  | <0.003 | PO <sub>4</sub> <sup>-3</sup> | <0.08 |

Balances are calibrated regularly with weight sets traceable to NIST #32856, #32867 and others. This CRM is guaranteed stable and accurate to +/- 0.5% of the certified value. This includes uncertainty components due to preparation, homogeneity by the most precise method, and short-term and long-term stability. This guarantee is valid for a period of one year from the date of certification only when the material is kept tightly capped and stored under ambient laboratory conditions.

Date of Certification:

**JUL - - 2022**

Certifying Officer:

  
Katherine Cullinan, QC Manager

Page 1 of 2  
Rev. 0

# Report of Certification

This Certified Reference Material (CRM) has been prepared and certified under an ISO 9001 (certified by DQS), ISO/IEC 17025 (accredited by A2LA) and ISO 17034 (accredited by A2LA) quality system consistent with the following guides:

- ISO 9001: Quality management systems – Requirements
- ISO/IEC 17025: General requirements for the competence of testing and calibration laboratories
- ISO 17034: General requirements for the competence of reference material producers
- ISO Guide 30: Reference Materials – Selected terms and definitions
- ISO Guide 31: Reference Materials – Contents of certificates, labels, and accompanying documentation
- ISO Guide 35: Reference Materials – Guidance for characterization and assessment of homogeneity and stability
- Guide to the Expression of Uncertainty in Measurement, 2008
- EURACHEM/CITAC Guide: Qualifying Uncertainty in Analytical Measurement – Third Edition
- NIST Technical Note 1297

## Material Source:

All analytes and matrix materials are obtained and verified by Spex CertiPrep from pre-qualified vendors as per ISO 9001, ISO/IEC 17025 and ISO 17034 guidelines. Vendor identifications are proprietary; however, sources of all materials used in the preparation and testing of Spex CertiPrep CRMs are tracked and documented. For further assistance, please contact Sales Support at USMet-CRMSales@antylia.com.

## Instructions for Use:

Primary usage of this CRM is in neat form or diluted serially with matrix of a purity at or greater than the purity of the original matrix solution. If dilution is required, the diluent must be compatible with all certified analytes and contain stabilizers appropriate for the period of intended use. The CRM can also be used as a spike or with a spike, again with appropriate compatibility considerations. All solutions should be thoroughly mixed, by shaking, prior to use and never pipetted directly from the bottle. Do not return excess solution to the bottle. All surfaces that come in contact with the solution must be thoroughly cleaned and leached prior to use. Dilutions should be performed only with Class A volumetric glassware. See SDS for health and safety information.

## Method of Preparation:

Clean laboratory procedures and techniques have been used throughout the preparation. All materials, equipment, analytical instrumentation and personnel have been qualified prior to use. The highest purity acids applicable, 18 megohm, double deionized water, acid-leached triple-rinsed bottles (where appropriate), and Class A/calibrated volumetrics have been used in all preparations.

## Homogeneity:

The homogeneity of the CRM has been confirmed by procedures consistent with ISO/IEC 17025, ISO 17034, and ASTM D6362-98 Appendix X2. Random, replicate samples of the final, packaged material have been analyzed to prove homogeneity in accordance with our internal procedure 4600-HOMOGEN-1A. Since the product is highly homogeneous, any sample size taken for analysis would be within the uncertainty budget. This is consistent with the intended use of the CRM.

## Statistical Estimator and Confidence Limits:

The certified value 'X' listed on the reverse of this document is at the 95% level of confidence and can be expressed as:

- $X = x \pm U$  where  $\bar{X}$  = certified value,  $U$  = expanded uncertainty,  $x$  = property value
- $U = k u_c$  where  $k = 2$  is the coverage factor at the 95% confidence level
- $u_c$  = combined standard uncertainty obtained by combining the individual element standard uncertainty components  $u_i$  and  $u_c = \sqrt{\sum u_i^2}$

## Certification Report:

All certified values reported were derived from the Certification Report, Spex CertiPrep's traceability documentation, identified by the lot number of this CRM. During the stated period of validity, the purchaser will be notified if this product is recalled due to any significant changes in the stability of the solution. For further assistance, please contact Sales Support at USMet-CRMSales@antylia.com.

## Legal Notice:

Spex CertiPrep reference materials are not for any cosmetic, drug or household application and are to be used only by qualified individuals who are trained in appropriate procedures. No claims against Spex CertiPrep, LLC. of any kind whatsoever, whether based on breach of warranty, alleged negligence, or otherwise, with respect to this Reference Material shall be greater than the purchase price. In no event shall Spex CertiPrep, LLC. be liable for any loss of profits or any incidental, special, or consequential damages.

Reagent

---

**IC SO4 ICV\_00024**





A Waters Company

Certified Reference Material

▪ **Certificate of Analysis** ▪

**Product:** 1000 mg/L Sulfate  
**Catalog Number:** 062-125mL, 995-500mL  
**Lot No.** 021221m  
**Starting Material:** Potassium Sulfate ( $K_2SO_4$ )  
**Matrix:** 18 megohm deionized water  
**Density:**  $1.0002 \pm 0.0005$  g/mL 20.4 °C and 748 mm Hg  
**Verification Method:** Ion Chromatography  
**Certificate Issue Date:** December 22, 2021  
**Expiration Date:** December 1, 2023  
**Revision Number:** Original

**CERTIFICATION**

| Parameter | Certified Value <sup>1</sup> | Uncertainty <sup>2</sup> | NIST Traceability       |            |
|-----------|------------------------------|--------------------------|-------------------------|------------|
|           | mg/L                         | %                        | SRM Number <sup>3</sup> | Recovery % |
| Sulfate   | 1000                         | 0.994                    | 3181                    | 97.5       |

# ▪ Certificate of Analysis ▪

1. The **Certified Value** is the actual gravimetric/volumetric "made-to" concentration confirmed by ERA analytical verification. The certified value is monitored and the purchasers will be notified of any significant changes resulting in recertification or withdrawal of this certified reference material during the period of validity of this certificate.

2. The **Uncertainty** represents an expanded uncertainty and approximates a 95% confidence interval. The uncertainty is based on the characterization, homogeneity and stability characteristics of the product, multiplied by a coverage factor ( $k=2$ ). The uncertainty applies to the product as supplied and does not take into account any required or optional dilution and/or preparations the laboratory may perform while using this product. The formula used to calculate the expanded uncertainty is:  

$$U_{\text{expanded}} = k \cdot \text{SQRT}((U_{\text{char}})^2 + (U_{\text{homogen}})^2 + (U_{\text{LTS}})^2 + (U_{\text{STS}})^2 + (U_{\text{RSS}})^2)$$

Where:

$U_{\text{expanded}}$  = Expanded uncertainty.

$k$  = Coverage factor.

$U_{\text{char}}$  = Combined standard uncertainty of the manufacturing and/or analytical verification assessment.

$U_{\text{homogen}}$  = Standard uncertainty of the homogeneity assessment.

$U_{\text{LTS}}$  = Standard uncertainty associated with long-term stability.

$U_{\text{STS}}$  = Standard uncertainty associated with short-term (transport) stability.

$U_{\text{RSS}}$  = Standard uncertainty associated with repeated sampling of the product (where permitted by product use instructions).

3. Where NIST Standard Reference Materials (SRMs) are available, each analyte has been analytically traced to the NIST SRM listed. **Analytical Traceability Recovery (%)** =  $[(\% \text{ recovery ERA certified reference material})/(\% \text{ recovery NIST SRM})] \cdot 100$

The traceability data shown were compiled by analyzing this ERA certified reference material and/or its associated stock solution(s) against the applicable NIST SRMs.

4. **Metrological Traceability.** This certified reference material is metrologically traceable to NIST mass reference materials through an unbroken chain of comparisons.

5. **Storage:** 20-25°C

6. **Intended Use:** This standard is intended to be used to calibrate your analytical process and/or as a quality control check of the entire process for the analytes/matrix included in the standard.

7. **Minimum Sample Size:** ERA suggests that when subsampling this product prior to analysis, you use a minimum sample size of at least 1 mL. Using a smaller sample size may invalidate the assigned value and/or uncertainty shown.

8. **Repeat Sampling:** Repeated Sampling of this product is permitted, provided minimum sample sizes and storage instructions are adhered to.

9. **Safety:** ERA products may be hazardous and are intended for use by professional laboratory personnel trained in the competent handling of such materials. Responsibility for the safe use of these products rests entirely with the buyer and/or user. Safety Data Sheets (SDS) for all ERA products are available through our website: [www.eraqc.com](http://www.eraqc.com)

If you have any questions or need technical assistance, please call ERA technical assistance at 1-800-372-0122 or send an email to [info@eraqc.com](mailto:info@eraqc.com).

Certifying Officer

Brian Miller



Quality Officer

Matthew Seebeck



ISO 17034:2016




ISO/IEC 17025:2017



Reagent

---

**MNX , TNX , DNX \_ 00069**



## Reference Material Certificate

**Product Name:** Custom Standard**Lot Number:** 0006659035**Product Number:** CUS-23984**Lot Issue Date:** 11-Jan-2022**Storage Conditions:** Store at Room Temperature (15° to 30°C).**Expiration Date:** 28-Feb-2023

| Component Name                                      | CERTIFIED VALUES |                      | CAS# | Analyte Lot |
|-----------------------------------------------------|------------------|----------------------|------|-------------|
|                                                     | Concentration    | Expanded Uncertainty |      |             |
| 1,3,5-trinitroso-1,3,5-triazacyclohexane (TNX)      | 100.4            | ± 0.5 µg/mL          | N/A  | RM12426     |
| 1-nitro-3,5-dinitroso-1,3,5-triazacyclohexane (DNX) | 100.4            | ± 0.5 µg/mL          | N/A  | RM12428     |
| 1-nitroso-3,5-dinitro-1,3,5-triazacyclohexane (MNX) | 117.0            | ± 0.6 µg/mL          | N/A  | RM12428     |

**Matrix:** acetonitrile**Description:**

This document is prepared in accordance with ISO 17034 and Guide 31. This analytical reference material standard was manufactured and verified in accordance with an ISO 9001 registered quality system and analyte concentrations were verified by an ISO 17025 accredited laboratory. The concentration and uncertainty value at the 95% confidence level for each analyte, determined gravimetrically, is listed above.

**Traceability:**

The balances used for these measurements are calibrated with weights traceable to NIST in compliance with ANSI/NCCL Z540.3, ISO 9001, ISO 17025, and ISO 17034. Calibrated Class A glassware is used for volumetric measurements. Thermometers are calibrated against a NIST traceable thermometer in accordance with NIST Special Publication 1088.

**Homogeneity:**

This analytical reference standard was unitized according to an in-house procedure and is guaranteed to be homogeneous. There is no minimum sub-sample size required.

**Instructions for Use:**

Sample aliquots for analysis should be withdrawn at 20°C to 25°C immediately after opening the container and should be processed without delay for the certified values to be valid within the stated uncertainties.

**Safety:**

Refer to the Safety Data Sheet on [www.agilent.com](http://www.agilent.com) for information regarding this analytical reference material.

**Intended Use:**

This analytical reference standard is intended for the preparation of working reference samples for use in routine laboratory analyses, calibration of instruments, validation of analytical methods, assessments of measurement methods, and continuing calibration verification.

**Expiration of Certification:**

Reagent

---

**MNX , TNX , DNX \_ 00076**

**Reference Material Certificate**  
**Product Information Sheet**

**Product Name:** Custom Standard

**Lot Number:** 0006725091

**Product Number:** CUS-23984

**Lot Issue Date:** 31-Jan-2023

**Storage Conditions:** Store at Room Temperature (15° to 30°C).

**Expiration Date:** 29-Feb-2024

| Component Name                                      | Concentration | Uncertainty | CAS# | Analyte Lot |
|-----------------------------------------------------|---------------|-------------|------|-------------|
| 1,3,5-trinitroso-1,3,5-triazacyclohexane (TNX)      | 100.3 ±       | 0.5 µg/mL   | N/A  | RM12426     |
| 1-nitro-3,5-dinitroso-1,3,5-triazacyclohexane (DNX) | 100.1 ±       | 0.5 µg/mL   | N/A  | RM12428     |
| 1-nitroso-3,5-dinitro-1,3,5-triazacyclohexane (MNX) | 116.7 ±       | 0.6 µg/mL   | N/A  | RM12428     |

**Matrix:** acetonitrile

**Description:**

This document is prepared in accordance with ISO 17034 and Guide 31. This analytical reference material (RM) standard was manufactured and verified in accordance with an ISO 9001 registered quality system and analyte concentrations were verified by an ISO 17025 accredited laboratory. The concentration and uncertainty value at the 95% confidence level for each analyte, determined gravimetrically, is listed above. Purity values are taken from approved vendor raw material certificates.

**Traceability:**

The balances used for these measurements are calibrated with weights traceable to NIST in compliance with ANSI/NC SL Z540.3, ISO 9001, ISO 17025, and ISO 17034. Calibrated Class A glassware is used for volumetric measurements. Thermometers are calibrated against a NIST traceable thermometer in accordance with NIST Special Publication 1088.

**Homogeneity:**

This analytical reference (RM) standard was unitized according to an in-house procedure and is guaranteed to be homogeneous. There is no minimum sub-sample size required.

**Instructions for Use:**

Sample aliquots for analysis should be withdrawn at 20°C to 25°C immediately after opening the container and should be processed without delay for the certified values to be valid within the stated uncertainties.

**Safety:**

Refer to the Safety Data Sheet on [www.agilent.com](http://www.agilent.com) for information regarding this analytical reference material.

**Intended Use:**

This analytical reference (RM) standard is intended for the preparation of working reference samples for use in routine laboratory analyses, calibration of instruments, validation of analytical methods, assessments of measurement methods, and continuing calibration verification.

**Expiration of Certification:**

Reagent

---

**MNX , TNX , DNX \_ 00077**

**Reference Material Certificate**  
**Product Information Sheet**

**Product Name:** Custom Standard

**Lot Number:** 0006725091

**Product Number:** CUS-23984

**Lot Issue Date:** 31-Jan-2023

**Storage Conditions:** Store at Room Temperature (15° to 30°C).

**Expiration Date:** 29-Feb-2024

| Component Name                                      | Concentration | Uncertainty | CAS# | Analyte Lot |
|-----------------------------------------------------|---------------|-------------|------|-------------|
| 1,3,5-trinitroso-1,3,5-triazacyclohexane (TNX)      | 100.3 ±       | 0.5 µg/mL   | N/A  | RM12426     |
| 1-nitro-3,5-dinitroso-1,3,5-triazacyclohexane (DNX) | 100.1 ±       | 0.5 µg/mL   | N/A  | RM12428     |
| 1-nitroso-3,5-dinitro-1,3,5-triazacyclohexane (MNX) | 116.7 ±       | 0.6 µg/mL   | N/A  | RM12428     |

**Matrix:** acetonitrile

**Description:**

This document is prepared in accordance with ISO 17034 and Guide 31. This analytical reference material (RM) standard was manufactured and verified in accordance with an ISO 9001 registered quality system and analyte concentrations were verified by an ISO 17025 accredited laboratory. The concentration and uncertainty value at the 95% confidence level for each analyte, determined gravimetrically, is listed above. Purity values are taken from approved vendor raw material certificates.

**Traceability:**

The balances used for these measurements are calibrated with weights traceable to NIST in compliance with ANSI/NC SL Z540.3, ISO 9001, ISO 17025, and ISO 17034. Calibrated Class A glassware is used for volumetric measurements. Thermometers are calibrated against a NIST traceable thermometer in accordance with NIST Special Publication 1088.

**Homogeneity:**

This analytical reference (RM) standard was unitized according to an in-house procedure and is guaranteed to be homogeneous. There is no minimum sub-sample size required.

**Instructions for Use:**

Sample aliquots for analysis should be withdrawn at 20°C to 25°C immediately after opening the container and should be processed without delay for the certified values to be valid within the stated uncertainties.

**Safety:**

Refer to the Safety Data Sheet on [www.agilent.com](http://www.agilent.com) for information regarding this analytical reference material.

**Intended Use:**

This analytical reference (RM) standard is intended for the preparation of working reference samples for use in routine laboratory analyses, calibration of instruments, validation of analytical methods, assessments of measurement methods, and continuing calibration verification.

**Expiration of Certification:**



Reagent

---

**MNX , TNX , DNX \_ 00079**

**Reference Material Certificate**  
**Product Information Sheet**

**Product Name:** Custom Standard

**Lot Number:** 0006725091

**Product Number:** CUS-23984

**Lot Issue Date:** 31-Jan-2023

**Storage Conditions:** Store at Room Temperature (15° to 30°C).

**Expiration Date:** 29-Feb-2024

| Component Name                                      | Concentration | Uncertainty | CAS# | Analyte Lot |
|-----------------------------------------------------|---------------|-------------|------|-------------|
| 1,3,5-trinitroso-1,3,5-triazacyclohexane (TNX)      | 100.3 ±       | 0.5 µg/mL   | N/A  | RM12426     |
| 1-nitro-3,5-dinitroso-1,3,5-triazacyclohexane (DNX) | 100.1 ±       | 0.5 µg/mL   | N/A  | RM12428     |
| 1-nitroso-3,5-dinitro-1,3,5-triazacyclohexane (MNX) | 116.7 ±       | 0.6 µg/mL   | N/A  | RM12428     |

**Matrix:** acetonitrile

**Description:**

This document is prepared in accordance with ISO 17034 and Guide 31. This analytical reference material (RM) standard was manufactured and verified in accordance with an ISO 9001 registered quality system and analyte concentrations were verified by an ISO 17025 accredited laboratory. The concentration and uncertainty value at the 95% confidence level for each analyte, determined gravimetrically, is listed above. Purity values are taken from approved vendor raw material certificates.

**Traceability:**

The balances used for these measurements are calibrated with weights traceable to NIST in compliance with ANSI/NC SL Z540.3, ISO 9001, ISO 17025, and ISO 17034. Calibrated Class A glassware is used for volumetric measurements. Thermometers are calibrated against a NIST traceable thermometer in accordance with NIST Special Publication 1088.

**Homogeneity:**

This analytical reference (RM) standard was unitized according to an in-house procedure and is guaranteed to be homogeneous. There is no minimum sub-sample size required.

**Instructions for Use:**

Sample aliquots for analysis should be withdrawn at 20°C to 25°C immediately after opening the container and should be processed without delay for the certified values to be valid within the stated uncertainties.

**Safety:**

Refer to the Safety Data Sheet on [www.agilent.com](http://www.agilent.com) for information regarding this analytical reference material.

**Intended Use:**

This analytical reference (RM) standard is intended for the preparation of working reference samples for use in routine laboratory analyses, calibration of instruments, validation of analytical methods, assessments of measurement methods, and continuing calibration verification.

**Expiration of Certification:**

Reagent

---

**NH3 CAL STD\_00035**

# Certificate of Analysis

**Ammonia Nitrogen Standard, 1000 ppm N (1216 ppm NH<sub>3</sub>)**

**Lot Number:** 4208D17

**Product Number:** 5455

**Manufacture Date:** AUG 05, 2022

**Expiration Date:** JAN 2024

The certified value reported is the prepared value based upon the method of preparation of the material. The uncertainty in the prepared value is based upon the volumetric method of preparation.

| Name              | CAS#       | Grade           |
|-------------------|------------|-----------------|
| Water             | 7732-18-5  | ACS/ASTM/USP/EP |
| Ammonium Chloride | 12125-02-9 | High Purity     |

| Test         | Specification    | Result   |
|--------------|------------------|----------|
| Appearance   | Colorless liquid | Passed   |
| Nitrogen (N) | 995-1005 ppm     | 1000 ppm |

| Specification                                                         | Reference         |
|-----------------------------------------------------------------------|-------------------|
| Ammonia Solution, Stock (1.0 mL = 1.0 mg ammonia nitrogen)            | ASTM (D 3590 A)   |
| Ammonia Solution, Stock (1.0 mL = 1.0 mg ammonium nitrogen)           | ASTM (D 3590 B)   |
| Stock Ammonium Chloride Solution                                      | APHA (4500-CN- L) |
| Stock Ammonium Solution                                               | APHA (4500-NH3 C) |
| Stock Ammonium chloride Solution                                      | APHA (4500-NH3 D) |
| Stock Ammonium Solution                                               | APHA (4500-NH3 F) |
| Ammonium Chloride, Stock Solution, 1.0 mL = 1.0 mg NH <sub>3</sub> -N | EPA (351.2)       |
| Ammonium Chloride, Stock Solution, 1.0 mL = 1.0 mg NH <sub>3</sub> -N | EPA (350.2)       |
| Ammonium Chloride, Stock Solution, 1.0 mL = 1.0 mg NH <sub>3</sub> -N | EPA (350.3)       |
| Ammonium Chloride, Stock Solution, 1.0 mL = 1.0 mg NH <sub>3</sub> -N | EPA (351.4)       |
| Stock Solution, 1.0 mL = 1.0 mg NH <sub>3</sub> -N                    | EPA (350.1)       |
| Ammonium Chloride, Stock Solution, 1.0 mL = 1.0 mg NH <sub>3</sub> -N | EPA (351.3)       |

Volumetric glassware complies with Class A tolerance requirements of ASTM E 288 and NIST Circular 434; it is calibrated before first use and recalibrated regularly in accordance with ASTM E 542 and NIST Procedure NBSIR 74-461. Balances are calibrated regularly with weights certified traceable to the NIST national mass standard. Thermometers and temperature probes are calibrated before first use and recalibrated regularly with a thermometer traceable to NIST standards. All products are prepared according to master documents that assure manufacture according to validated methods. Batch records document raw material traceability and production and testing history for each lot manufactured.

| Part Number | Size / Package Type | Shelf Life (Unopened Container) |
|-------------|---------------------|---------------------------------|
| 5455-16     | 500 mL natural poly | 18 months                       |

**Recommended Storage:** 15°C - 30°C (59°F - 86°F)



Paul Brandon (08/05/2022)

Production Manager

This Certificate of Analysis is designed to comply with ISO Guide 31 "Reference Materials -- Contents of Certificates and Labels."

This test report shall not be reproduced, except in full, without the written approval of Ricca Chemical Company.

Reagent

---

**NH3 ICV STD\_00034**

## 1.0 ACCREDITATION / REGISTRATION

**INORGANIC VENTURES** is accredited to ISO 17034, "General Requirements for the Competence of Reference Material Producers" and ISO/IEC 17025, "General Requirements for the Competence of Testing and Calibration Laboratories". Inorganic Ventures is also an ISO 9001 registered manufacturer (QSR Certificate Number QSR-1034).



## 2.0 PRODUCT DESCRIPTION

Product Code: Single Analyte Ion Chromatography Solution  
Catalog Number: ICNNH41  
Lot Number: S2-NH700817  
Matrix: H<sub>2</sub>O  
Value / Analyte(s): 1 000 µg/mL ea:  
Ammonium as N  
Starting Material: Ammonium chloride  
Starting Material Lot#: 1736  
Starting Material Purity: 99.8500%

## 3.0 CERTIFIED VALUES AND UNCERTAINTIES

**Certified Value:** 999 ± 4 µg/mL  
**Density:** 0.999 g/mL (measured at 20 ± 4 °C)

### Assay Information:

|                        |                                                                  |
|------------------------|------------------------------------------------------------------|
| <b>Assay Method #1</b> | <b>995 ± 4 µg/mL</b><br>Fajans NIST SRM 999c Lot Number: 999c    |
| <b>Assay Method #2</b> | <b>1002 ± 4 µg/mL</b><br>IC Assay NIST SRM 194a Lot Number: 194a |

- The Calculated Value is a value calculated from the weight of a starting material that has been certified directly vs. a National Institute of Standards and Technology (NIST) SRM/RM. See Sec 4.2 for balance traceability.

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of  $k = 2$ .

#### Characterization of CRM/RM by Two or More Methods

Certified Value,  $X_{\text{CRM/RM}}$ , where two or more methods of characterization are used is the weighted mean of the results:

$$X_{\text{CRM/RM}} = \sum (w_i) (X_i)$$

$X_i$  = mean of Assay Method  $i$  with standard uncertainty  $u_{\text{char } i}$

$w_i$  = the weighting factors for each method calculated using the inverse square of the variance:

$$w_i = (1/u_{\text{char } i}^2) / (\sum (1/u_{\text{char } i}^2))$$

$$\text{CRM/RM Expanded Uncertainty } (\pm) = U_{\text{CRM/RM}} = k (u_{\text{char}}^2 + u_{\text{bb}}^2 + u_{\text{Its}}^2 + u_{\text{ts}}^2)^{1/2}$$

$k$  = coverage factor = 2

$u_{\text{char}} = [\sum (w_i)^2 (u_{\text{char } i}^2)]^{1/2}$  where  $u_{\text{char } i}$  are the errors from each characterization method

$u_{\text{bb}}$  = bottle to bottle homogeneity standard uncertainty

$u_{\text{Its}}$  = long term stability standard uncertainty (storage)

$u_{\text{ts}}$  = transport stability standard uncertainty

#### Characterization of CRM/RM by One Method

Certified Value,  $X_{\text{CRM/RM}}$ , where one method of characterization is used is the mean of individual results:

$$X_{\text{CRM/RM}} = (X_a) (u_{\text{char } a})$$

$X_a$  = mean of Assay Method A with

$u_{\text{char } a}$  = the standard uncertainty of characterization Method A

$$\text{CRM/RM Expanded Uncertainty } (\pm) = U_{\text{CRM/RM}} = k (u_{\text{char } a}^2 + u_{\text{bb}}^2 + u_{\text{Its}}^2 + u_{\text{ts}}^2)^{1/2}$$

$k$  = coverage factor = 2

$u_{\text{char } a}$  = the errors from characterization

$u_{\text{bb}}$  = bottle to bottle homogeneity standard uncertainty

$u_{\text{Its}}$  = long term stability standard uncertainty (storage)

$u_{\text{ts}}$  = transport stability standard uncertainty

## 4.0 TRACEABILITY TO NIST

- This product is traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM/RM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRM/RM are available, the term 'in-house std.' is specified.

### 4.1 Thermometer Calibration

- All thermometers are NIST traceable through thermometers that are calibrated by an accredited calibration laboratory.

### 4.2 Balance Calibration


- All analytical balances are calibrated by an accredited calibration laboratory and procedure. The weights used for testing are annually compared to master weights and are traceable to NIST.

### 4.3 Glassware Calibration

- An in-house procedure is used to calibrate all Class A glassware used in the manufacturing and quality control of CRM/RMs.

## 5.0 CHROMATOGRAM





#### Dionex ICS-2000 Ion Chromatograph

**Analytical Column:** IonPac CS18 2 x 250 mm

**Eluent:** 10 mM MSA

**Guard Column:** IonPac CG18 2 x 50 mm

**Eluent Flow Rate:** 0.25 mL/min

**Anion Self Regen**

**Column Temp:** 30°C

**Suppressor/** N/A

**Cell Temp:** 35°C

**Chemical**

**Suppression:**

**Scale X-Axis:** minutes

**Cation Self Regen**

**Scale Y-Axis:** 8 µS/cm

**Suppressor/** CERS 500 2mm

**Chemical**

**Concentration:** 10 µg/g

**Suppression:**

**Suppressor** 8 mA

**Current/ Chemical**

**Suppressant:**

#### 6.0 INTENDED USE

- For the calibration of analytical instruments and validation of analytical methods as appropriate.

#### 7.0 INSTRUCTIONS FOR THE CORRECT USE OF THIS REFERENCE MATERIAL

##### 7.1 Storage and Handling Recommendations

- Store between approximately 4° - 30° C while in sealed TCT bag.

- While stored in the sealed TCT bag, transpiration of this CRM/RM is negligible. After opening the sealed TCT bag transpiration of the CRM/RM will occur, resulting in a gradual increase in the analyte concentration(s). It is the responsibility of the user to account for this effect. When the bottle is weighed both before and after being placed in storage, the mass difference observed will be a measure of transpiration mass loss.

- After opening the sealed TCT bag, keep cap tightly sealed when not in use and store between 4° - 24° C to minimize the effects of transpiration. Use at 20° ± 4° C to minimize volumetric dilution error when using the reported density. Do not pipette from the container. Do not return removed aliquots to container.

- For more information, visit [www.inorganicventures.com/TCT](http://www.inorganicventures.com/TCT)

## **8.0 HAZARDOUS INFORMATION**

- Please refer to the Safety Data Sheet for information regarding this CRM/RM.

## **9.0 HOMOGENEITY**

- This solution was mixed according to an in-house procedure and is guaranteed to be homogeneous. Homogeneity data indicate that the end user should take a minimum sample size of 0.2 mL to assure homogeneity.

## **10.0 QUALITY STANDARD DOCUMENTATION**

### **10.1 ISO 9001 Quality Management System Registration**

- QSR Certificate Number QSR-1034

### **10.2 ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories"**

- Chemical Testing - Accredited / A2LA Certificate Number 883.01

### **10.3 ISO 17034 "General Requirements for the Competence of Reference Material Producers"**

- Reference Material Producer - Accredited / A2LA Certificate Number 883.02

Inorganic Ventures, 300 Technology Drive, Christiansburg, Va. 24073, USA; Telephone: 800.669.6799; 540.585.3030, Fax: 540.585.3012; [inorganicventures.com](http://inorganicventures.com); [info@inorganicventures.com](mailto:info@inorganicventures.com)

## **11.0 CERTIFICATION, LOT EXPIRATION AND PERIOD OF VALIDITY**

### **11.1 Certification Issue Date**

January 25, 2021

- The certification is valid within the measurement uncertainty specified provided the CRM/RM is stored and handled in accordance with instructions given in Sec 7.1. This certification is nullified if instructions in Sec 7.1 are not followed or if the CRM/RM is damaged, contaminated, or otherwise modified.

### **11.2 Lot Expiration Date**

- **January 25, 2025**

- The date after which this CRM/RM should not be used.

- The lot expiration date reflects the period of time that the stability of a CRM/RM can be supported by long term stability studies conducted on properly stored and handled CRM/RMs. Lot expiration is limited primarily by transpiration (loss of water from the solution) and infrequently by chemical stability.

### 11.3 Period of Validity

- Sealed TCT Bag Open Date: \_\_\_\_\_

- This CRM/RM should not be used longer than one year (or six months in the case of a 30 mL bottle) from the date of opening the aluminized bag or after the date given in Sec. 11.2, whichever comes first. This is contingent upon the CRM/RM being stored and handled in accordance with the instructions given in Sec. 7.1.

## 12.0 NAMES AND SIGNATURES OF CERTIFYING OFFICERS

### Certificate Prepared By:

Uyen Truong  
Supervisor, Product Documentation



### Certificate Approved By:

Michael Booth  
Director, Quality Control



### Certifying Officer:

Paul Gaines  
Chairman / Senior Technical Director



Reagent

---

**PicricARestek\_00113**



CERTIFIED REFERENCE MATERIAL

110 Benner Circle  
Bellefonte, PA 16823-8812  
Tel: (800)356-1688  
Fax: (814)353-1309

www.restek.com

## Certificate of Analysis



FOR LABORATORY USE ONLY-READ SDS PRIOR TO USE.

*This Reference Material is intended for Laboratory Use Only as a standard for the qualitative and/or quantitative determination of the analyte(s) listed.*

**Catalog No.:** 31499 **Lot No.:** A0183202

**Description:** Picric Acid Standard

Picric Acid Standard 1000µg/mL, Methanol, 1mL/1000µg/mL \*PGI BOX REQUIRED\* SHIP FED EX GROUND ONLY

**Container Size:** 2 mL **Pkg Amt:** > 1 mL

**Expiration Date:** March 31, 2027 **Storage:** 10°C or colder

**Ship:** Ambient

### CERTIFIED VALUES

| Elution Order | Compound                    | Grav. Conc. (weight/volume) | Expanded Uncertainty (95% C.L.; K=2) |
|---------------|-----------------------------|-----------------------------|--------------------------------------|
| 1             | Picric Acid                 | 1,004.0 µg/mL               | +/- 5.9635 µg/mL Gravimetric         |
|               | CAS # 88-89-1 (Lot 06130CU) |                             | +/- 55.0021 µg/mL Unstressed         |
|               | Purity 99%                  |                             | +/- 64.1379 µg/mL Stressed           |

**Solvent:** Methanol

**CAS #** 67-56-1

**Purity** 99%

**Column:**  
250mm x 4.6mm  
Ultra C18 (cat.# 9174575)


**Flow Rate:**  
1.0 ml/min.

**Mobile Phase A:**  
water:methanol (44:56 V/V)


**Mobile Phase B:**

**Mobile Phase Composition:**  
100%A

**Det. Type:**  
Wavelength: 210nm & 254nm



This chromatogram represents a general set of testing conditions chosen for product acceptance. For optimal results in your lab, conditions should be adjusted for your specific instrument, method, and application.

  
Brandon Reish - Mix Technician

Date Mixed: 23-Mar-2022      Balance: 1128360905

  
Amanda Miller - Operations Tech-ARM QC

Date Passed: 28-Mar-2022

Manufactured under Restek's ISO 9001:2015  
Registered Quality System  
Certificate #FM 80397

## General Certified Reference Material Notes

### Expiration Notes:

- Expiration date valid for unopened ampul stored in compliance with the recommended conditions.
- Uncertainty, concentration, and expiration of the CRM are based on the unopened product being stored according to the recommended condition found in the storage field.

### Purity Notes:

- Purity and/or chemical identity are determined by one or more of the following techniques: GC/FID, HPLC, GC/μECD, GC/MS, LC/MS, RI, and/or melting point.
- Compounds with a listed purity of less than 99% have been weight corrected to compensate for impurities and/or salts. A correction factor is used to calculate the amount of compound necessary to achieve the desired concentration of the parent compound in solution.
- Purity of isomeric compounds is reported as the sum of the isomers.
- Purity values are rounded to the nearest whole number.

### Certified Uncertainty Value Notes:

- The uncertainties are determined in accordance with ISO 17034 and Guide 35. The certified combined stressed uncertainty value ( includes gravimetric uncertainty, homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty and were combined using the following formula:

$$U_{combined\ stressed} = k \sqrt{U_{gravimetric}^2 + U_{homogeneity}^2 + U_{storage\ stability}^2 + U_{shipping\ stability}^2}$$

$k$  is a coverage factor of 2, which gives a level of confidence of approximately 95%.

- It is important to note that the shipping stability uncertainty was obtained under temperature extremes for specific time intervals; therefore, the certified combined stressed uncertainty value should only be applied to the product if it was stored at non-standard temperature conditions up to and including 7 days. Contact Restek Technical Service at [www.restek.com/Contact-Us](http://www.restek.com/Contact-Us) for use recommendations if your shipment was in-transit for more than 7 days at non-standard temperature conditions.
- Apply the certified combined unstressed uncertainty value if the product was received under standard shipping conditions. Apply the certified combined stressed uncertainty value if the product was received under non-standard conditions as specified below.

| Label Conditions                                          | Standard Conditions | Non-Standard Conditions |
|-----------------------------------------------------------|---------------------|-------------------------|
| 25°C Nominal (Room Temperature)                           | < 60°C              | ≥ 60°C up to 7 days     |
| 10°C or colder (Refrigerate)                              | < 40°C              | ≥ 40°C up to 7 days     |
| 0°C or colder (Freezer)<br>-20°C or colder (Deep Freezer) | < 25°C              | ≥ 25°C up to 7 days     |

- Separate (not combined) uncertainty values for gravimetric uncertainty are also displayed on the certificate, if needed, separate homogeneity between-ampul uncertainty, storage stability uncertainty and shipping stability uncertainty values are available by contacting Restek Technical Service at [www.restek.com/Contact-Us](http://www.restek.com/Contact-Us).
- The packaged amount is the minimum sample size for which uncertainty is valid. The ampules are over-filled to ensure that the minimum packaged amount can be sufficiently transferred.

### Manufacturing Notes:

- Concentration is based upon gravimetric preparation using either a balance whose calibration has been verified daily using NIST traceable weights, and/or dilutions with Class A glassware.

### Handling Notes:

- Stability of the unopened product, when stored in compliance with the recommended conditions, is guaranteed through the expiration displayed on the product label and certificate. Contact Restek for additional opened product stability information, with the knowledge/understanding that open product stability is subject to the specific handling and environmental conditions to which the product is exposed. For your convenience Restek supplies deactivated vials with most standards packed in 2mL ampules. Larger volume deactivated vials are available through Restek as a custom ordered item. Additionally, Restek sells DMDCS for the purpose of glassware deactivation as catalog number 31861, which includes complete instructions.





Reagent

---

**SFD CAL STK\_00010**



|                             |                                        |
|-----------------------------|----------------------------------------|
| Version                     | 02                                     |
| Molecular weight            | 240.18                                 |
| Quality Test / Release Date | 08/28/2022                             |
| Molecular Formula           | Na <sub>2</sub> S . 9 H <sub>2</sub> O |
| CAS No                      | 1313-84-4                              |
| Linear Formula              | Na <sub>2</sub> S.9H <sub>2</sub> O    |
| Flash Point (°C)            |                                        |

## Certificate of Analysis

Products are processed under ISO 9001:2015 quality management systems and samples are tested for conformance to the noted specifications. Certain data may have been supplied by third parties. We disclaim the implied warranties of merchantability and fitness for a particular purpose, and the accuracy of third party data or information associated with the product. Products are for research and development use only. Products are not for direct administration to humans or animals. It is the responsibility of the final formulator or end user to determine suitability, and to qualify and/or validate each product for its intended use.

|                       |                                            |                             |            |
|-----------------------|--------------------------------------------|-----------------------------|------------|
| Catalog Number        | 38706                                      | Quality Test / Release Date | 08/28/2022 |
| Lot Number            | A0411500                                   | Suggested retest date       | 08/28/2025 |
| Description           | Sodium sulfide nonahydrate,98+%,extra pure |                             |            |
| Country of Origin     | INDIA                                      |                             |            |
| Declaration of Origin | synthetic                                  |                             |            |

|          |  |
|----------|--|
| BSE/TSE  |  |
| Chemical |  |

| Result name                                   | Specifications                     | Test Value                   |
|-----------------------------------------------|------------------------------------|------------------------------|
| Appearance (Color)                            | Colorless to light yellow          | off white                    |
| Appearance (Form)                             | Adhering crystals and/or chunks    | Adhering crystals and chunks |
| Titration Iodimetric                          | 32.0 to 38.0 % (Na <sub>2</sub> S) | 32.0 % (Na <sub>2</sub> S)   |
| Total nitrogen (as N)                         | =<0.005 %                          | =<0.005 %                    |
| Sulfite (as SO <sub>2</sub> )                 | =<2000 ppm                         | =<1000 ppm                   |
| Thiosulphate (S <sub>2</sub> O <sub>3</sub> ) | =<5000 ppm                         | =<1000 ppm                   |

C. Wygaerts, QA Manager

Issued: 08-29-2022

Acros Organics BV  
ENA23, zone1, nr 1350, Janssen Pharmaceuticaaan 3a, B-2440 Geel, Belgium  
Tel +32 14/57.52.11 - Fax+32 14/59.34.34 Internet: <https://www.thermofisher.com>  
1 Reagent Lane, Fair Lawn, NJ 07410, USA Fax 201-796-1329

Reagent

---

**TOC ICV Std\_00051**

# Certificate of Analysis

## Organic Carbon Standard, 1000 ppm C

**Lot Number:** 2208C86

**Product Number:** 1847

**Manufacture Date:** AUG 08, 2022

**Expiration Date:** JUL 2023

The certified value reported is the prepared value based upon the method of preparation of the material. The uncertainty in the prepared value is based upon the volumetric method of preparation.

| Name                     | CAS#      | Grade           |
|--------------------------|-----------|-----------------|
| Water                    | 7732-18-5 | ACS/ASTM/USP/EP |
| Phosphoric Acid          | 7664-38-2 | ACS             |
| Potassium Acid Phthalate | 877-24-7  | ACS Acidimetric |

| Test       | Specification    | Result   |
|------------|------------------|----------|
| Appearance | Colorless liquid | Passed   |
| Carbon (C) | 995-1005 ppm     | 1000 ppm |

| Specification                                                      | Reference           |
|--------------------------------------------------------------------|---------------------|
| Organic Carbon Stock Solution                                      | APHA (5310 B)       |
| Potassium Hydrogen Phthalate, Stock Solution                       | EPA (SW-846) (9060) |
| Potassium Hydrogen Phthalate, Stock Solution, 1000 mg Carbon/liter | EPA (415.1)         |
| Organic Carbon Solution, Standard (1 mL = 1 mg C)                  | ASTM (D 2579)       |

Volumetric glassware complies with Class A tolerance requirements of ASTM E 288 and NIST Circular 434; it is calibrated before first use and recalibrated regularly in accordance with ASTM E 542 and NIST Procedure NBSIR 74-461. Balances are calibrated regularly with weights certified traceable to the NIST national mass standard. Thermometers and temperature probes are calibrated before first use and recalibrated regularly with a thermometer traceable to NIST standards. All products are prepared according to master documents that assure manufacture according to validated methods. Batch records document raw material traceability and production and testing history for each lot manufactured.

| Part Number | Size / Package Type | Shelf Life (Unopened Container) |
|-------------|---------------------|---------------------------------|
| 1847-16     | 500 mL amber glass  | 12 months                       |
| 1847-32     | 1 L amber glass     | 12 months                       |

**Recommended Storage:** 15°C - 30°C (59°F - 86°F)



Myrlande Gilles (08/08/2022)

### Quality Control

This Certificate of Analysis is designed to comply with ISO Guide 31 "Reference Materials -- Contents of Certificates and Labels."

This test report shall not be reproduced, except in full, without the written approval of Ricca Chemical Company.

# Safety Data Sheet

Classified According to OSHA Hazard Communication Standard (HCS)

## SECTION 1: Identification

### 1.1. Product Identifier

**Trade Name or Designation:** Organic Carbon Standard, 1000 ppm C

**Product Number:** 1847

**Other Identifying Product Numbers:** 1847-1, 1847-16, 1847-2.5, 1847-32, 1847-4, 1847-5, 1847-8

### 1.2. Recommended Use and Restrictions on Use

General Laboratory Reagent

### 1.3. Details of the Supplier of the Safety Data Sheet

**Company:** Ricca Chemical Company

**Address:** 448 West Fork Drive

Arlington, TX 76012 USA

**Telephone:** 888-467-4222

### 1.4. Emergency Telephone Number (24 hours)

CHEMTREC (USA) 800-424-9300

CHEMTREC (International) 1+ 703-527-3887

## SECTION 2: Hazard(s) Identification

### 2.1. Classification of the Substance or Mixture

For the full text of the Hazard and Precautionary Statements listed below, see Section 16.

| Hazard Class                                 | Category   | Hazard Statements | Precautionary Statements: |
|----------------------------------------------|------------|-------------------|---------------------------|
| Hazardous to the Aquatic Environment (Acute) | Category 3 | H402              | P273, P501                |

### 2.2. GHS Label Elements

**Pictograms:**

# Safety Data Sheet

**Signal Word:**

**Hazard Statements:**

| Hazard Number | Hazard Statement         |
|---------------|--------------------------|
| H402          | Harmful to aquatic life. |

**Precautionary Statements:**

| Precautionary Number | Precautionary Statement                                                                     |
|----------------------|---------------------------------------------------------------------------------------------|
| P273                 | Avoid release to the environment.                                                           |
| P501                 | Dispose of contents in accordance with local, state, federal and international regulations. |

## 2.4. Hazards not Otherwise Classified or Covered by GHS

Data not available.

## SECTION 3: Composition / Information on Ingredients

### 3.1. Components of Substance or Mixture

| Chemical Name            | Formula                                | Molecular Weight | CAS Number | Weight% |
|--------------------------|----------------------------------------|------------------|------------|---------|
| Water                    | H <sub>2</sub> O                       | 18.01 g/mol      | 7732-18-5  | 99.39   |
| Phosphoric Acid          | H <sub>3</sub> PO <sub>4</sub>         | 97.99 g/mol      | 7664-38-2  | 0.40    |
| Potassium Acid Phthalate | HOCOC <sub>6</sub> H <sub>4</sub> COOK | 204.22 g/mol     | 877-24-7   | 0.21    |

## SECTION 4: First-Aid Measures

### 4.1. General First Aid Information

**Eye Contact:** May cause irritation, redness, pain, and tearing.

**Inhalation:** Not expected to require first aid. If necessary, remove to fresh air.

**Skin Contact:** May cause slight irritation.

**Ingestion:** Dilute with water or milk. Call a physician if necessary.

### 4.2. Most Important Symptoms and Effects, Acute and Delayed

Does not present any significant health hazards. Wash areas of contact with water. EYE CONTACT: May cause irritation, redness, pain, and tearing. SKIN CONTACT: May cause slight irritation.

### 4.3. Medical Attention or Special Treatment Needed

Irrigate immediately with large quantity of water for at least 15 minutes. Call a physician if irritation develops. Remove to fresh air. Give artificial respiration if necessary. If breathing is difficult, give oxygen. Flush with plenty of water for at least 15 minutes. Call a physician if irritation develops. Dilute with water or milk. Call a physician if necessary.

## Safety Data Sheet

### SECTION 5: Fire-Fighting Measures

#### 5.1. Extinguishing Media

Use any means suitable for extinguishing surrounding fire.

#### 5.2. Specific Hazards Arising from the Substance or Mixture

Not considered to be a fire or explosion hazard.

#### 5.3. Special Protective Equipment for Firefighters

Use protective clothing and breathing equipment appropriate for the surrounding fire.

### SECTION 6: Accidental Release Measures

#### 6.1. Personal Precautions, Protective Equipment and Emergency Procedures

Wear appropriate PPE for the size and nature of the spill. As a general rule, wear safety glasses and gloves.

#### 6.2. Cleanup and Containment Methods and Materials

Collect liquid and dilute with water. Neutralize with dilute base solutions. Resulting solution may be flushed to the sewer with excess water. For larger spills, absorb with suitable material and dispose of in accordance with local regulations.

### SECTION 7: Handling and Storage

#### 7.1. Precautions for Safe Handling and Storage Conditions

As with all chemicals, wash hands thoroughly after handling. Avoid contact with eyes and skin. Protect from freezing and physical damage. Keep tightly capped. Store in the refrigerator.

### SECTION 8: Exposure Controls / Personal Protection

#### 8.1 Control Parameters

| Chemical Name               | Limit Type | Country | Exposure Limit           | Information Source                                                     |
|-----------------------------|------------|---------|--------------------------|------------------------------------------------------------------------|
| Phosphoric Acid (7664-38-2) | TWA        | USA     | 1 mg/m <sup>3</sup> TWA  | U.S. - OSHA - Final PELs - Time Weighted Averages (TWAs)               |
| Phosphoric Acid (7664-38-2) | TLV-STEL   | USA     | 3 mg/m <sup>3</sup> STEL | ACGIH - Threshold Limit Values - Short Term Exposure Limits (TLV-STEL) |
| Phosphoric Acid (7664-38-2) | TLV-TWA    | USA     | 1 mg/m <sup>3</sup> TWA  | ACGIH - Threshold Limit Values - Time Weighted Averages (TLV-TWA)      |

## Safety Data Sheet

### 8.2. Exposure Controls

**Engineering Controls:** No specific controls are needed. Normal room ventilation is adequate.

**Respiratory Protection:** Normal room ventilation is adequate.

**Skin Protection:** Chemical resistant gloves.

**Eye Protection:** Safety glasses or goggles.

### 8.3. Personal Protective Equipment

Normal room ventilation is adequate. Chemical resistant gloves. Safety glasses or goggles.

## SECTION 9: Physical and Chemical Properties

### 9.1. Basic Physical and Chemical Properties

**Appearance:** Colorless liquid

**Physical State:** Liquid

**Odor:** Data not available.

**Odor Threshold:** Data not available.

**pH:** < 2

**Melting/Freezing Point:** 0.0°C

**Initial Boiling Point/Range:** 100°C - 100°C

**Flash Point:** Data not available.

**Evaporation Rate:** Data not available.

**Flammability:** Data not available.

**Flammability/Explosive Limits:** Data not available.

**Vapor Pressure:** Data not available.

**Vapor Density:** Data not available.

**Relative Density:** 1.00

**Solubility:** Miscible

**Partition Coefficient:** Data not available.

**Auto-Ignition Temperature:** Data not available.

**Decomposition Temperature:** Data not available.

**Viscosity:** Data not available.

**Explosive Properties:** Data not available.

**Oxidizing Properties:** Data not available.



# Safety Data Sheet

## SECTION 10: Stability and Reactivity

### 10.1. Reactivity and Chemical Stability

Stable under normal conditions of use and storage.

### 10.2. Possibility of Hazardous Reactions

Data not available.

### 10.3. Conditions to Avoid and Incompatible Materials

Nitric Acid.

### 10.4. Hazardous Decomposition Products

Will not occur.

## SECTION 11: Toxicological Information

### 11.1. Information on Toxicological Effects

#### Acute Toxicity - Oral Exposure:

Not applicable.

#### Acute Toxicity - Dermal Exposure:

Not applicable.

#### Acute Toxicity - Inhalation Exposure:

Not applicable.

#### Acute Toxicity - Other Information:

LD50, Oral, Rat: >3200 mg/kg (Potassium Acid Phthalate), details of toxic effects not reported other than lethal dose value. LD50, Oral, Rat: (Phosphoric Acid) 1530 mg/kg, behavioral, kidney, bladder, hair effects noted. LD50, Dermal, Rat: 2740 mg/kg (Phosphoric Acid), behavioral effects noted.

#### Skin Corrosion and Irritation:

Not applicable.

#### Serious Eye Damage and Irritation:

Not applicable.

#### Respiratory Sensitization:

Not applicable.

#### Skin Sensitization:

Not applicable.

#### Germ Cell Mutagenicity:

Not applicable.

#### Carcinogenicity:

Not applicable.

## Safety Data Sheet

**Reproductive Toxicity:**

Not applicable.

**Specific Target Organ Toxicity from Single Exposure:**

Not applicable.

**Specific Target Organ Toxicity from Repeated Exposure:**

Not applicable.

**Aspiration Hazard:**

Not applicable.

**Additional Toxicology Information:**

Data not available.

### SECTION 12: Ecological Information

**12.1. Ecotoxicity**

Harmful to aquatic life. Avoid release to the environment. Dispose of contents in accordance with local, state, federal and international regulations.

**12.2. Persistence and Degradability**

Data not available.

**12.3. Bioaccumulative Potential**

Data not available.

**12.4. Mobility in Soil**

Data not available.

**12.5. Other Adverse Ecological Effects**

Data not available.

### SECTION 13: Disposal Considerations

**13.1. Waste Treatment Methods**

Data not available.

## Safety Data Sheet

### SECTION 14: Transportation Information

#### 14.1. Transportation by Land-Department of Transportation (DOT, United States of America)

Not regulated according to DOT Regulations.

#### 14.2. Transportation by Air - International Air Transport Association (IATA)

Not regulated according to IATA Dangerous Goods Regulations.

#### 14.3 Transportation of Dangerous Goods (TDG, Canada)

Not regulated according to TDG Regulations.

# Safety Data Sheet

## SECTION 15: Regulatory Information

### 15.1. Occupational Safety and Health Administration (OSHA) Hazards

Not listed.

### 15.2. Superfund Amendments and Reauthorization Act (SARA) 302 Extremely Hazardous Substances

Not listed.

### 15.3. Superfund Amendments and Reauthorization Act (SARA) 311/312 Hazardous Chemicals

Phosphoric Acid (CAS # 7664-38-2): 5000 lb final RQ; 2270 kg final RQ

### 15.4. Superfund Amendments and Reauthorization Act (SARA) 313 Toxic Release Inventory (TRI)

Not listed.

### 15.5. Massachusetts Right-to-Know Substance List

Phosphoric Acid (CAS # 7664-38-2): Present

### 15.6. Pennsylvania Right-to-Know Hazardous Substances

Phosphoric Acid (CAS # 7664-38-2): Environmental hazard

Phosphoric Acid (CAS # 7664-38-2): Present

Water (CAS # 7732-18-5): Present

### 15.7. New Jersey Worker and Community Right-to-Know Components

Phosphoric Acid (CAS # 7664-38-2): corrosive

Phosphoric Acid (CAS # 7664-38-2): sn 1516

### 15.8. California Proposition 65

Not listed.

### 15.9. Canada Domestic Substances List / Non-Domestic Substances List (DSL/NDSL)

Phosphoric Acid (CAS # 7664-38-2): Present (DSL)

Water (CAS # 7732-18-5): Present (DSL)

Potassium Acid Phthalate (CAS # 877-24-7): Present (DSL)

### 15.10. United States of America Toxic Substances Control Act (TSCA) List

**All components of this solution are listed as active on the TSCA Inventory or are mixtures (hydrates) of active items listed on the TSCA Inventory.**

Phosphoric Acid (CAS # 7664-38-2): Present

Water (CAS # 7732-18-5): Present

Potassium Acid Phthalate (CAS # 877-24-7): Present

## Safety Data Sheet

### 15.11. European Inventory of Existing Commercial Chemical Substances (EINECS), European List of Notified Chemical Substances (ELINCS), and No Longer Polymers (NLP)

Phosphoric Acid (CAS # 7664-38-2): 231-633-2

Water (CAS # 7732-18-5): 231-791-2

Water (CAS # 7732-18-5): 232-148-9

Potassium Acid Phthalate (CAS # 877-24-7): 212-889-4

### SECTION 16: Other Information

#### 16.1. Full Text of Hazard Statements and Precautionary Statements

Harmful to aquatic life.

Avoid release to the environment.

Dispose of contents in accordance with local, state, federal and international regulations.

#### 16.2. Miscellaneous Hazard Classes

**Canadian Carcinogenicity Hazard Class:** Not Applicable.


**Physical Hazards Not Otherwise Classified (PHNOC):** Not Applicable.

**Health Hazards Not Otherwise Classified (HHNOC):** Not Applicable.

**Biohazardous Infectious Materials Hazard Class:** Not Applicable.

#### 16.3. National Fire Protection Association (NFPA) Rating

Health: 0  
 Flammability: 0  
 Reactivity: 0  
 Special Hazard:



#### 16.4. Document Revision

**Last Revision Date:** 2023-02-14

## DISCLAIMER

When handled properly by qualified personnel, the product described herein does not present a significant health or safety hazard. Alteration of its characteristics by concentration, evaporation, addition of other substances, or other means may present hazards not specifically addressed herein and which must be evaluated by the user. The information furnished herein is believed to be accurate and represents the best data currently available to us. No warranty, expressed or implied, is made and RICCA CHEMICAL COMPANY assumes no legal responsibility or liability whatsoever resulting from its use.

Reagent

---

**TOC LCS Std\_00057**

## Reference Material Certificate

### Product Information Sheet

**Product Name:** Total Organic Carbon (TOC) Standard  
**Product Number:** IQC-106  
**Storage Conditions:** Store at Room Temperature (15° to 30°C).

**Lot Number:** 0006719669  
**Lot Issue Date:** 28-Dec-2022  
**Expiration Date:** 31-Jan-2025

| Component Name                              | CERTIFIED VALUES |                      |  | CAS#        | Analyte Lot |
|---------------------------------------------|------------------|----------------------|--|-------------|-------------|
|                                             | Concentration    | Expanded Uncertainty |  |             |             |
| potassium hydrogen phthalate (KHP) (as TOC) | 1001             | ± 5 mg/L             |  | 000877-24-7 | RM19275     |

**Matrix:** water (low TOC, < 50 ppb)

#### Description:

This document is prepared in accordance with ISO 17034 and Guide 31. This analytical reference material (RM) was manufactured and verified in accordance with an ISO 9001 registered quality system. The analyte concentration(s) were prepared and verified by an ISO 17034 / ISO 17025 accredited laboratory and compared to calibration standards independently prepared using NIST SRM(s) when available. The certified value and uncertainty value at the 95% confidence level for each analyte is determined gravimetrically.

#### Traceability:

The balances used for these measurements are calibrated with weights traceable to NIST in compliance with ANSI/NCSL Z540.3, ISO 9001, ISO 17025, and ISO 17034. Calibrated Class A glassware is used for volumetric measurements. Thermometers are calibrated against a NIST traceable thermometer in accordance with NIST Special Publication 1088.

#### Homogeneity:

This analytical reference standard was unitized according to an in-house procedure and is guaranteed to be homogeneous. There is no minimum sub-sample size required.

#### Instructions for Use:

Sample aliquots for analysis should be withdrawn at 20°C to 25°C immediately after opening the container and should be processed without delay for the certified values to be valid within the stated uncertainties.

#### Safety:

Refer to the Safety Data Sheet on [www.agilent.com](http://www.agilent.com) for information regarding this analytical reference material.

#### Intended Use:

This analytical reference standard is intended for the preparation of working reference samples for use in routine laboratory analyses, calibration of instruments, validation of analytical methods, assessments of measurement methods, and continuing calibration verification.

#### Expiration of Certification:

The certification of this analytical reference standard is valid until the expiration date specified above, provided the material is handled and stored in accordance with the instructions given in this certificate. This certification is nullified if the material is damaged, contaminated, or otherwise modified.





**Maintenance of Certification:**

If substantive changes are noted that affect the certification before the expiration of this certificate, Agilent will notify the purchaser.

**Sample lot approver:**

Monica Bourgeois  
QMS Representative



ISO 17034 Cert  
No. AR-1936

RM was produced in accordance with the TUV/SUD registered ISO  
9001:2015 Quality Management System. Cert# 951215321

Page: 2 of 2

[www.agilent.com/quality/](http://www.agilent.com/quality/)  
CSD-QA-040.1

ISO 17025

# 8330B\_DOD5

---

Nitroaromatics and Nitramines (HPLC)

FORM II  
HPLC/IC SURROGATE RECOVERY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
 SDG No.: \_\_\_\_\_  
 Matrix: Water Level: Low  
 GC Column (1): UltraCarb5u ID: 4.6 (mm) GC Column (2): Luna-phenyl 4.6 (mm)

| Client Sample ID              | Lab Sample ID          | 12DNB1 # | 12DNB2 # |
|-------------------------------|------------------------|----------|----------|
| FWGmw-015-230401-G<br>W       | 280-176674-2           | 85 M     |          |
| LL3mw-241-230401-G<br>W       | 280-176674-4           | 93       |          |
| LL3mw-241-230401-G<br>W       | 280-176674-4           |          | 88       |
| FBQmw-173-230401-G<br>W-R     | 280-176674-7           | 91 M     |          |
| FBQmw-173-230401-G<br>W-R     | 280-176674-7           |          | 85       |
| FBQmw-175-230401-G<br>W-R     | 280-176674-10          | 95 M     |          |
| FBQmw-175-230402-G<br>W-R     | 280-176674-11          | 93 M     |          |
|                               | MB<br>280-613095/1-A   | 102      |          |
|                               | LCS<br>280-613095/2-A  | 93       |          |
|                               | LCS<br>280-613095/4-A  | 98       |          |
|                               | LCSD<br>280-613095/3-A | 94       |          |
|                               | LCSD<br>280-613095/5-A | 98       |          |
| FBQmw-175-230401-G<br>W-R MS  | 280-176674-10 MS       | 82 Q     |          |
| FBQmw-175-230401-G<br>W-R MS  | 280-176674-10 MS       | 95 M     |          |
| FBQmw-175-230401-G<br>W-R MSD | 280-176674-10<br>MSD   | 98       |          |
| FBQmw-175-230401-G<br>W-R MSD | 280-176674-10<br>MSD   | 109      |          |

12DNB = 1,2-Dinitrobenzene

QC LIMITS  
83-119

# Column to be used to flag recovery values

FORM II 8330B

FORM III  
HPLC/IC LAB CONTROL SAMPLE RECOVERY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
 SDG No.: \_\_\_\_\_  
 Matrix: Water Level: Low Lab File ID: 05190035.D  
 Lab ID: LCS 280-613095/2-A Client ID: \_\_\_\_\_

| COMPOUND                   | SPIKE<br>ADDED<br>(ug/L) | LCS<br>CONCENTRATION<br>(ug/L) | LCS<br>%<br>REC | QC<br>LIMITS<br>REC | # |
|----------------------------|--------------------------|--------------------------------|-----------------|---------------------|---|
| 1,3,5-Trinitrobenzene      | 2.00                     | 1.98                           | 99              | 73-125              |   |
| 1,3-Dinitrobenzene         | 2.00                     | 1.90                           | 95              | 78-120              |   |
| 2,4,6-Trinitrotoluene      | 2.00                     | 1.83                           | 91              | 71-123              |   |
| 2,4-Dinitrotoluene         | 2.00                     | 1.82                           | 91              | 78-120              |   |
| 2,6-Dinitrotoluene         | 2.00                     | 1.84                           | 92              | 77-127              |   |
| 2-Amino-4,6-dinitrotoluene | 2.00                     | 1.80                           | 90              | 79-120              |   |
| 2-Nitrotoluene             | 2.00                     | 1.56                           | 78              | 70-127              |   |
| 3,5-Dinitroaniline         | 2.00                     | 1.73                           | 87              | 71-117              |   |
| 3-Nitrotoluene             | 2.00                     | 1.54                           | 77              | 73-125              |   |
| 4-Amino-2,6-dinitrotoluene | 2.00                     | 1.78                           | 89              | 76-125              |   |
| 4-Nitrotoluene             | 2.00                     | 1.56                           | 78              | 71-127              |   |
| DNX                        | 2.00                     | 1.93                           | 96              | 66-119              | M |
| HMX                        | 2.00                     | 1.62                           | 81              | 65-135              | M |
| MNX                        | 2.33                     | 2.45                           | 105             | 57-132              |   |
| Nitrobenzene               | 2.00                     | 1.74                           | 87              | 65-134              |   |
| Nitroglycerin              | 20.0                     | 18.9                           | 95              | 74-127              |   |
| PETN                       | 20.0                     | 20.5                           | 102             | 73-127              |   |
| RDX                        | 2.00                     | 1.79                           | 89              | 68-130              |   |
| Tetryl                     | 2.00                     | 2.07                           | 103             | 64-128              |   |
| TNX                        | 2.01                     | 1.88                           | 94              | 50-150              | M |

# Column to be used to flag recovery and RPD values

FORM III  
HPLC/IC LAB CONTROL SAMPLE RECOVERY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
SDG No.: \_\_\_\_\_  
Matrix: Water Level: Low Lab File ID: 05190037.D  
Lab ID: LCS 280-613095/4-A Client ID: \_\_\_\_\_

| COMPOUND                   | SPIKE<br>ADDED<br>(ug/L) | LCS<br>CONCENTRATION<br>(ug/L) | LCS<br>%<br>REC | QC<br>LIMITS<br>REC | # |
|----------------------------|--------------------------|--------------------------------|-----------------|---------------------|---|
| 2,4-diamino-6-nitrotoluene | 2.00                     | 1.61                           | 81              | 68-122              | M |
| 2,6-diamino-4-nitrotoluene | 2.00                     | 1.68                           | 84              | 72-122              | M |

# Column to be used to flag recovery and RPD values

FORM III  
HPLC/IC LAB CONTROL SAMPLE DUPLICATE RECOVERY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
 SDG No.: \_\_\_\_\_  
 Matrix: Water Level: Low Lab File ID: 05190036.D  
 Lab ID: LCSD 280-613095/3-A Client ID: \_\_\_\_\_

| COMPOUND                   | SPIKE<br>ADDED<br>(ug/L) | LCSD<br>CONCENTRATION<br>(ug/L) | LCSD<br>%<br>REC | %<br>RPD | QC LIMITS |        | # |
|----------------------------|--------------------------|---------------------------------|------------------|----------|-----------|--------|---|
|                            |                          |                                 |                  |          | RPD       | REC    |   |
| 1,3,5-Trinitrobenzene      | 2.00                     | 2.04                            | 102              | 3        | 20        | 73-125 |   |
| 1,3-Dinitrobenzene         | 2.00                     | 1.96                            | 98               | 3        | 20        | 78-120 |   |
| 2,4,6-Trinitrotoluene      | 2.00                     | 1.88                            | 94               | 3        | 20        | 71-123 |   |
| 2,4-Dinitrotoluene         | 2.00                     | 1.87                            | 93               | 3        | 20        | 78-120 |   |
| 2,6-Dinitrotoluene         | 2.00                     | 1.95                            | 98               | 6        | 20        | 77-127 |   |
| 2-Amino-4,6-dinitrotoluene | 2.00                     | 1.84                            | 92               | 2        | 20        | 79-120 |   |
| 2-Nitrotoluene             | 2.00                     | 1.64                            | 82               | 5        | 20        | 70-127 |   |
| 3,5-Dinitroaniline         | 2.00                     | 1.83                            | 91               | 5        | 20        | 71-117 |   |
| 3-Nitrotoluene             | 2.00                     | 1.62                            | 81               | 5        | 20        | 73-125 |   |
| 4-Amino-2,6-dinitrotoluene | 2.00                     | 1.82                            | 91               | 2        | 20        | 76-125 |   |
| 4-Nitrotoluene             | 2.00                     | 1.61                            | 80               | 3        | 20        | 71-127 |   |
| DNX                        | 2.00                     | 1.95                            | 97               | 1        | 20        | 66-119 | M |
| HMX                        | 2.00                     | 1.71                            | 86               | 6        | 20        | 65-135 | M |
| MNX                        | 2.33                     | 2.47                            | 106              | 1        | 20        | 57-132 |   |
| Nitrobenzene               | 2.00                     | 1.81                            | 91               | 4        | 20        | 65-134 |   |
| Nitroglycerin              | 20.0                     | 19.5                            | 97               | 3        | 20        | 74-127 |   |
| PETN                       | 20.0                     | 21.1                            | 106              | 3        | 20        | 73-127 |   |
| RDX                        | 2.00                     | 1.84                            | 92               | 3        | 20        | 68-130 |   |
| Tetryl                     | 2.00                     | 2.06                            | 103              | 0        | 20        | 64-128 |   |
| TNX                        | 2.01                     | 1.91                            | 95               | 2        | 20        | 50-150 | M |

# Column to be used to flag recovery and RPD values

FORM III  
HPLC/IC LAB CONTROL SAMPLE DUPLICATE RECOVERY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
SDG No.: \_\_\_\_\_  
Matrix: Water Level: Low Lab File ID: 05190038.D  
Lab ID: LCSD 280-613095/5-A Client ID: \_\_\_\_\_

| COMPOUND                   | SPIKE<br>ADDED<br>(ug/L) | LCSD<br>CONCENTRATION<br>(ug/L) | LCSD<br>%<br>REC | %<br>RPD | QC LIMITS |        | # |
|----------------------------|--------------------------|---------------------------------|------------------|----------|-----------|--------|---|
|                            |                          |                                 |                  |          | RPD       | REC    |   |
| 2,4-diamino-6-nitrotoluene | 2.00                     | 1.47                            | 74               | 9        | 20        | 68-122 | M |
| 2,6-diamino-4-nitrotoluene | 2.00                     | 1.56                            | 78               | 8        | 20        | 72-122 | M |

# Column to be used to flag recovery and RPD values

FORM III  
HPLC/IC MATRIX SPIKE RECOVERY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
 SDG No.: \_\_\_\_\_  
 Matrix: Water Level: Low Lab File ID: 05190049.D  
 Lab ID: 280-176674-10 MS Client ID: FBQmw-175-230401-GW-R MS

| COMPOUND                   | SPIKE<br>ADDED<br>(ug/L) | SAMPLE<br>CONCENTRATION<br>(ug/L) | MS<br>CONCENTRATION<br>(ug/L) | MS<br>%<br>REC | QC<br>LIMITS<br>REC | #    |
|----------------------------|--------------------------|-----------------------------------|-------------------------------|----------------|---------------------|------|
| 1,3,5-Trinitrobenzene      | 2.00                     | 0.20 U                            | 1.75                          | 87             | 73-125              | Q    |
| 1,3-Dinitrobenzene         | 2.00                     | 0.10 U                            | 1.66                          | 83             | 78-120              | Q    |
| 2,4,6-Trinitrotoluene      | 2.00                     | 0.10 U                            | 1.56                          | 78             | 71-123              | Q    |
| 2,4-Dinitrotoluene         | 2.00                     | 0.080 U                           | 1.50                          | 75             | 78-120              | Q J1 |
| 2,6-Dinitrotoluene         | 2.00                     | 0.080 U                           | 1.58                          | 79             | 77-127              | Q    |
| 2-Amino-4,6-dinitrotoluene | 2.00                     | 0.10 U                            | 1.46                          | 73             | 79-120              | Q J1 |
| 2-Nitrotoluene             | 2.00                     | 0.20 U                            | 1.34                          | 67             | 70-127              | Q J1 |
| 3,5-Dinitroaniline         | 2.00                     | 0.30 U                            | 1.46                          | 73             | 71-117              | Q    |
| 3-Nitrotoluene             | 2.00                     | 0.35 U                            | 1.23                          | 61             | 73-125              | Q J1 |
| 4-Amino-2,6-dinitrotoluene | 2.00                     | 0.12 U                            | 1.44                          | 72             | 76-125              | Q J1 |
| 4-Nitrotoluene             | 2.00                     | 0.40 U                            | 1.30                          | 65             | 71-127              | Q J1 |
| DNX                        | 2.00                     | 0.25 U                            | 1.56                          | 78             | 66-119              | M Q  |
| HMX                        | 2.00                     | 0.20 U                            | 1.46                          | 73             | 65-135              | M Q  |
| MNX                        | 2.34                     | 0.29 U                            | 1.98                          | 85             | 57-132              | M Q  |
| Nitrobenzene               | 2.00                     | 0.20 U                            | 1.47                          | 73             | 65-134              | Q    |
| Nitroglycerin              | 20.0                     | 2.0 U                             | 17.5                          | 87             | 74-127              | Q    |
| PETN                       | 20.0                     | 1.0 U                             | 18.1                          | 90             | 73-127              | Q    |
| RDX                        | 2.00                     | 0.20 U                            | 1.49                          | 74             | 68-130              | M Q  |
| Tetryl                     | 2.00                     | 0.10 U                            | 1.72                          | 86             | 64-128              | Q    |
| TNX                        | 2.01                     | 0.25 U                            | 1.70                          | 85             | 50-150              | M Q  |

# Column to be used to flag recovery and RPD values



FORM III  
HPLC/IC MATRIX SPIKE RECOVERY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
SDG No.: \_\_\_\_\_  
Matrix: Water Level: Low Lab File ID: 05190051.D  
Lab ID: 280-176674-10 MS Client ID: FBQmw-175-230401-GW-R MS

| COMPOUND                   | SPIKE<br>ADDED<br>(ug/L) | SAMPLE<br>CONCENTRATION<br>(ug/L) | MS<br>CONCENTRATION<br>(ug/L) | MS<br>%<br>REC | QC<br>LIMITS<br>REC | #    |
|----------------------------|--------------------------|-----------------------------------|-------------------------------|----------------|---------------------|------|
| 2,4-diamino-6-nitrotoluene | 2.01                     | 0.90 U                            | 1.25                          | 62             | 68-122              | M J1 |
| 2,6-diamino-4-nitrotoluene | 2.01                     | 0.90 U                            | 1.60                          | 80             | 72-122              | M    |

# Column to be used to flag recovery and RPD values

FORM III  
HPLC/IC MATRIX SPIKE DUPLICATE RECOVERY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
 SDG No.: \_\_\_\_\_  
 Matrix: Water Level: Low Lab File ID: 05190050.D  
 Lab ID: 280-176674-10 MSD Client ID: FBQmw-175-230401-GW-R MSD

| COMPOUND                   | SPIKE<br>ADDED<br>(ug/L) | MSD<br>CONCENTRATION<br>(ug/L) | MSD<br>%<br>REC | %<br>RPD | QC LIMITS |        | #    |
|----------------------------|--------------------------|--------------------------------|-----------------|----------|-----------|--------|------|
|                            |                          |                                |                 |          | RPD       | REC    |      |
| 1,3,5-Trinitrobenzene      | 2.19                     | 2.21                           | 101             | 23       | 20        | 73-125 | J1   |
| 1,3-Dinitrobenzene         | 2.19                     | 2.13                           | 97              | 25       | 20        | 78-120 | J1   |
| 2,4,6-Trinitrotoluene      | 2.19                     | 2.00                           | 92              | 25       | 20        | 71-123 | J1   |
| 2,4-Dinitrotoluene         | 2.19                     | 2.00                           | 92              | 29       | 20        | 78-120 | J1   |
| 2,6-Dinitrotoluene         | 2.19                     | 2.02                           | 93              | 25       | 20        | 77-127 | J1   |
| 2-Amino-4,6-dinitrotoluene | 2.19                     | 1.96                           | 90              | 29       | 20        | 79-120 | J1   |
| 2-Nitrotoluene             | 2.19                     | 1.78                           | 82              | 28       | 20        | 70-127 | J1   |
| 3,5-Dinitroaniline         | 2.19                     | 1.86                           | 85              | 24       | 20        | 71-117 | J1   |
| 3-Nitrotoluene             | 2.19                     | 1.67                           | 76              | 30       | 20        | 73-125 | J1   |
| 4-Amino-2,6-dinitrotoluene | 2.19                     | 1.87                           | 86              | 26       | 20        | 76-125 | J1   |
| 4-Nitrotoluene             | 2.19                     | 1.75                           | 80              | 29       | 20        | 71-127 | J1   |
| DNX                        | 2.19                     | 1.93                           | 88              | 21       | 20        | 66-119 | M J1 |
| HMX                        | 2.19                     | 1.77                           | 81              | 19       | 20        | 65-135 | M    |
| MNX                        | 2.55                     | 2.37                           | 93              | 18       | 20        | 57-132 | M    |
| Nitrobenzene               | 2.19                     | 1.94                           | 89              | 28       | 20        | 65-134 | J1   |
| Nitroglycerin              | 21.9                     | 21.5                           | 98              | 20       | 20        | 74-127 |      |
| PETN                       | 21.9                     | 22.4                           | 102             | 21       | 20        | 73-127 | J1   |
| RDX                        | 2.19                     | 1.83                           | 84              | 21       | 20        | 68-130 | M J1 |
| Tetryl                     | 2.19                     | 2.23                           | 102             | 26       | 20        | 64-128 | J1   |
| TNX                        | 2.19                     | 1.98                           | 90              | 15       | 20        | 50-150 | M    |

# Column to be used to flag recovery and RPD values

FORM III  
HPLC/IC MATRIX SPIKE DUPLICATE RECOVERY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
SDG No.: \_\_\_\_\_  
Matrix: Water Level: Low Lab File ID: 05190052.D  
Lab ID: 280-176674-10 MSD Client ID: FBQmw-175-230401-GW-R MSD

| COMPOUND                   | SPIKE<br>ADDED<br>(ug/L) | MSD<br>CONCENTRATION<br>(ug/L) | MSD<br>%<br>REC | %<br>RPD | QC LIMITS |        | #    |
|----------------------------|--------------------------|--------------------------------|-----------------|----------|-----------|--------|------|
|                            |                          |                                |                 |          | RPD       | REC    |      |
| 2,4-diamino-6-nitrotoluene | 2.00                     | 1.03                           | 51              | 19       | 20        | 68-122 | M J1 |
| 2,6-diamino-4-nitrotoluene | 2.00                     | 1.32                           | 66              | 19       | 20        | 72-122 | M J1 |

# Column to be used to flag recovery and RPD values

FORM IV  
HPLC/IC METHOD BLANK SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
 SDG No.: \_\_\_\_\_  
 Lab Sample ID: MB 280-613095/1-A  
 Matrix: Water Date Extracted: 05/19/2023 13:39  
 Lab File ID: (1) 05190034.D Lab File ID: (2) \_\_\_\_\_  
 Date Analyzed: (1) 05/19/2023 23:27 Date Analyzed: (2) \_\_\_\_\_  
 Instrument ID: (1) CHHPLC\_X3 Instrument ID: (2) CHHPLC\_X5  
 GC Column: (1) UltraCarb5uO ID: 4.6 (mm) GC Column: (2) Luna-phenylh ID: 4.6 (mm)

THIS METHOD BLANK APPLIES TO THE FOLLOWING SAMPLES:

| CLIENT SAMPLE ID          | LAB SAMPLE ID       | DATE ANALYZED 1  | DATE ANALYZED 2  |
|---------------------------|---------------------|------------------|------------------|
|                           | LCS 280-613095/2-A  | 05/19/2023 23:50 |                  |
|                           | LCSD 280-613095/3-A | 05/20/2023 00:13 |                  |
|                           | LCS 280-613095/4-A  | 05/20/2023 00:36 |                  |
|                           | LCSD 280-613095/5-A | 05/20/2023 00:59 |                  |
| FWGmw-015-230401-GW       | 280-176674-2        | 05/20/2023 02:53 |                  |
| LL3mw-241-230401-GW       | 280-176674-4        | 05/20/2023 04:02 | 05/24/2023 19:48 |
| FBQmw-173-230401-GW-R     | 280-176674-7        | 05/20/2023 04:25 | 05/21/2023 00:31 |
| FBQmw-175-230401-GW-R     | 280-176674-10       | 05/20/2023 04:48 |                  |
| FBQmw-175-230401-GW-R MS  | 280-176674-10 MS    | 05/20/2023 05:11 |                  |
| FBQmw-175-230401-GW-R MSD | 280-176674-10 MSD   | 05/20/2023 05:34 |                  |
| FBQmw-175-230401-GW-R MS  | 280-176674-10 MS    | 05/20/2023 05:57 |                  |
| FBQmw-175-230401-GW-R MSD | 280-176674-10 MSD   | 05/20/2023 06:20 |                  |
| FBQmw-175-230402-GW-R     | 280-176674-11       | 05/20/2023 06:43 |                  |

FORM X  
IDENTIFICATION SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
 SDG No.: \_\_\_\_\_  
 Client Sample ID: LL3mw-241-230401-GW Lab Sample ID: 280-176674-4  
 Instrument ID (1): CHHPLC\_X3 Instrument ID (2): CHHPLC\_X5  
 Date Analyzed (1): 05/20/2023 04:02 Date Analyzed (2): 05/24/2023 19:48  
 GC Column (1): UltraCarb5uOD ID: 4.6(mm) GC Column (2): Luna-phenylhe ID: 4.6(mm)

| ANALYTE                    | COL | PEAK | RT    | RT WINDOW |       | CONCENTRATION |      | RPD  |
|----------------------------|-----|------|-------|-----------|-------|---------------|------|------|
|                            |     |      |       | FROM      | TO    | PEAK          | MEAN |      |
| HMX                        | 1   |      | 6.55  | 6.40      | 6.70  | 0.14          |      | 6.6  |
|                            | 2   |      | 6.73  | 6.54      | 6.84  | 0.15          |      |      |
| RDX                        | 1   |      | 7.57  | 7.42      | 7.72  | 1.0           |      | 72.2 |
|                            | 2   |      | 8.97  | 8.80      | 9.10  | 0.49          |      |      |
| 1,3,5-Trinitrobenzene      | 1   |      | 8.65  | 8.50      | 8.80  | 4.6           |      | 29.0 |
|                            | 2   |      | 17.34 | 17.18     | 17.48 | 6.1           |      |      |
| 3,5-Dinitroaniline         | 1   |      | 9.88  | 9.73      | 10.03 | 0.36          |      | 3.0  |
|                            | 2   |      | 14.28 | 14.14     | 14.44 | 0.37          |      |      |
| 2,4,6-Trinitrotoluene      | 1   |      | 10.91 | 10.81     | 11.01 | 2.0           |      | 8.5  |
|                            | 2   |      | 23.01 | 22.89     | 23.19 | 1.8           |      |      |
| 4-Amino-2,6-dinitrotoluene | 1   |      | 11.11 | 11.02     | 11.22 | 1.7           |      | 0.2  |
|                            | 2   |      | 16.37 | 16.26     | 16.56 | 1.7           |      |      |
| 2,4-Dinitrotoluene         | 1   |      | 11.67 | 11.58     | 11.78 | 0.12          |      | 89.4 |
|                            | 2   |      | 18.91 | 18.80     | 19.10 | 0.045         |      |      |

FORM X  
IDENTIFICATION SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
 SDG No.: \_\_\_\_\_  
 Client Sample ID: FBQmw-173-230401-GW-R Lab Sample ID: 280-176674-7  
 Instrument ID (1): CHHPLC\_X3 Instrument ID (2): CHHPLC\_X5  
 Date Analyzed (1): 05/20/2023 04:25 Date Analyzed (2): 05/21/2023 00:31  
 GC Column (1): UltraCarb5uOD ID: 4.6(mm) GC Column (2): Luna-phenylhe ID: 4.6(mm)

| ANALYTE                    | COL | PEAK | RT    | RT WINDOW |       | CONCENTRATION |      | RPD  |
|----------------------------|-----|------|-------|-----------|-------|---------------|------|------|
|                            |     |      |       | FROM      | TO    | PEAK          | MEAN |      |
| 3,5-Dinitroaniline         | 1   |      | 9.89  | 9.73      | 10.03 | 0.15          |      | 0.4  |
|                            | 2   |      | 14.34 | 14.23     | 14.53 | 0.15          |      |      |
| 2,4,6-Trinitrotoluene      | 1   |      | 10.91 | 10.81     | 11.01 | 0.11          |      | 0.1  |
|                            | 2   |      | 23.08 | 22.98     | 23.28 | 0.11          |      |      |
| 4-Amino-2,6-dinitrotoluene | 1   |      | 11.11 | 11.02     | 11.22 | 1.5           |      | 14.0 |
|                            | 2   |      | 16.43 | 16.36     | 16.66 | 1.8           |      |      |
| 2-Amino-4,6-dinitrotoluene | 1   |      | 11.37 | 11.27     | 11.47 | 1.2           |      | 9.0  |
|                            | 2   |      | 17.24 | 17.20     | 17.50 | 1.3           |      |      |

FORM I  
HPLC/IC ORGANICS ANALYSIS DATA SHEET

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Client Sample ID: FWGmw-015-230401-GW Lab Sample ID: 280-176674-2

Matrix: Water Lab File ID: 05190043.D

Analysis Method: 8330B Date Collected: 05/17/2023 14:06

Extraction Method: 3535 Date Extracted: 05/19/2023 13:39

Sample wt/vol: 481.7 (mL) Date Analyzed: 05/20/2023 02:53

Con. Extract Vol.: 5 (mL) Dilution Factor: 1

Injection Volume: 100 (uL) GC Column: UltraCarb5uODS ID: 4.6 (mm)

% Moisture: \_\_\_\_\_ % Solids: \_\_\_\_\_ GPC Cleanup: (Y/N) N

Cleanup Factor: \_\_\_\_\_

Analysis Batch No.: 613168 Units: ug/L

| CAS NO.    | COMPOUND NAME              | RESULT | Q   | LOQ  | LOD   | DL    |
|------------|----------------------------|--------|-----|------|-------|-------|
| 99-35-4    | 1,3,5-Trinitrobenzene      | 0.21   | U   | 0.22 | 0.21  | 0.087 |
| 99-65-0    | 1,3-Dinitrobenzene         | 0.10   | U   | 0.11 | 0.10  | 0.038 |
| 118-96-7   | 2,4,6-Trinitrotoluene      | 0.10   | U   | 0.11 | 0.10  | 0.047 |
| 121-14-2   | 2,4-Dinitrotoluene         | 0.083  | U   | 0.10 | 0.083 | 0.028 |
| 606-20-2   | 2,6-Dinitrotoluene         | 0.083  | U   | 0.10 | 0.083 | 0.042 |
| 35572-78-2 | 2-Amino-4,6-dinitrotoluene | 0.10   | U   | 0.11 | 0.10  | 0.053 |
| 88-72-2    | 2-Nitrotoluene             | 0.21   | U   | 0.22 | 0.21  | 0.089 |
| 99-08-1    | 3-Nitrotoluene             | 0.36   | U   | 0.42 | 0.36  | 0.20  |
| 19406-51-0 | 4-Amino-2,6-dinitrotoluene | 0.12   | U   | 0.16 | 0.12  | 0.060 |
| 99-99-0    | 4-Nitrotoluene             | 0.42   | U   | 0.43 | 0.42  | 0.10  |
| 2691-41-0  | HMX                        | 0.21   | U M | 0.22 | 0.21  | 0.091 |
| 98-95-3    | Nitrobenzene               | 0.21   | U   | 0.22 | 0.21  | 0.094 |
| 55-63-0    | Nitroglycerin              | 2.1    | U   | 2.2  | 2.1   | 0.96  |
| 78-11-5    | PETN                       | 1.0    | U   | 1.1  | 1.0   | 0.46  |
| 121-82-4   | RDX                        | 0.21   | U   | 0.22 | 0.21  | 0.053 |
| 479-45-8   | Tetryl                     | 0.10   | U   | 0.11 | 0.10  | 0.033 |

| CAS NO.  | SURROGATE          | %REC | Q | LIMITS |
|----------|--------------------|------|---|--------|
| 528-29-0 | 1,2-Dinitrobenzene | 85   | M | 83-119 |

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\05190043.D  
 Lims ID: 280-176674-B-2-A  
 Client ID: FWGmw-015-230401-GW  
 Sample Type: Client  
 Inject. Date: 20-May-2023 02:53:48 ALS Bottle#: 43 Worklist Smp#: 43  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-176674-B-2-A  
 Operator ID: JZ/JG Instrument ID: CHHPLC\_X3  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 20-May-2023 11:22:00 Calib Date: 24-Feb-2023 17:41:20  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1663

First Level Reviewer: LV5D

Date: 20-May-2023 11:05:10

| Compound                      | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | OnCol Amt<br>ug/mL | Flags |
|-------------------------------|-----|--------------|------------------|------------------|----------|--------------------|-------|
| 4 HMX                         | 1   |              | 6.550            |                  |          | ND                 | U     |
| 8 RDX                         | 1   |              | 7.570            |                  |          | ND                 |       |
| \$ 10 1,2-Dinitrobenzene      | 1   | 8.529        | 8.523            | 0.006            | 21540    | 0.1705             | M     |
| 11 1,3,5-Trinitrobenzene      | 1   |              | 8.650            |                  |          | ND                 |       |
| 12 1,3-Dinitrobenzene         | 1   |              | 9.277            |                  |          | ND                 |       |
| 13 Nitrobenzene               | 1   |              | 9.643            |                  |          | ND                 |       |
| 15 Tetryl                     | 1   |              | 10.017           |                  |          | ND                 |       |
| 16 Nitroglycerin              | 2   |              | 10.470           |                  |          | ND                 |       |
| 17 2,4,6-Trinitrotoluene      | 1   |              | 10.910           |                  |          | ND                 |       |
| 18 4-Amino-2,6-dinitrotoluene | 1   |              | 11.117           |                  |          | ND                 |       |
| 19 2-Amino-4,6-dinitrotoluene | 1   |              | 11.370           |                  |          | ND                 |       |
| 20 2,6-Dinitrotoluene         | 1   |              | 11.517           |                  |          | ND                 |       |
| 21 2,4-Dinitrotoluene         | 1   |              | 11.683           |                  |          | ND                 |       |
| 22 o-Nitrotoluene             | 1   |              | 12.517           |                  |          | ND                 |       |
| 23 p-Nitrotoluene             | 1   |              | 12.937           |                  |          | ND                 |       |
| 24 m-Nitrotoluene             | 1   |              | 13.517           |                  |          | ND                 |       |
| 25 PETN                       | 2   |              | 14.697           |                  |          | ND                 |       |

**QC Flag Legend**

Processing Flags

Review Flags

M - Manually Integrated

U - Marked Undetected



Report Date: 20-May-2023 11:22:12

Chrom Revision: 2.3 16-May-2023 16:37:30

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190043.d

Injection Date: 20-May-2023 02:53:48

Instrument ID: CHHPLC\_X3

Operator ID: JZ/JG

Lims ID: 280-176674-B-2-A

Lab Sample ID: 280-176674-2

Worklist Smp#: 43

Client ID: FWGmw-015-230401-GW

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 43

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



Eurofins Denver  
Recovery Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\05190043.D  
Lims ID: 280-176674-B-2-A  
Client ID: FWGmw-015-230401-GW  
Sample Type: Client  
Inject. Date: 20-May-2023 02:53:48 ALS Bottle#: 43 Worklist Smp#: 43  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Sample Info: 280-176674-B-2-A  
Operator ID: JZ/JG Instrument ID: CHHPLC\_X3  
Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\8330\_X3.m  
Limit Group: GCSV - 8330  
Last Update: 20-May-2023 11:22:00 Calib Date: 24-Feb-2023 17:41:20  
Integrator: Falcon  
Quant Method: External Standard Quant By: Initial Calibration  
Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
Process Host: CTX1663

First Level Reviewer: LV5D

Date: 20-May-2023 11:05:10

| Compound                 | Amount Added | Amount Recovered | % Rec. |
|--------------------------|--------------|------------------|--------|
| \$ 10 1,2-Dinitrobenzene | 0.2000       | 0.1705           | 85.27  |

## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190043.d  
Injection Date: 20-May-2023 02:53:48 Instrument ID: CHHPLC\_X3  
Lims ID: 280-176674-B-2-A Lab Sample ID: 280-176674-2  
Client ID: FWGmw-015-230401-GW  
Operator ID: JZ/JG ALS Bottle#: 43 Worklist Smp#: 43  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector: LC DAD1B, 254 nm

## \$ 10 1,2-Dinitrobenzene, CAS: 528-29-0

Signal: 1


RT: 8.53  
Area: 108785  
Amount: 0.861263  
Amount Units: ug/mL

## Processing Integration Results



RT: 8.53  
Area: 21540  
Amount: 0.170535  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 11:05:09 -06:00:00 (UTC)

Audit Action: Manually Integrated


Audit Reason: Baseline

## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190043.d  
Injection Date: 20-May-2023 02:53:48 Instrument ID: CHHPLC\_X3  
Lims ID: 280-176674-B-2-A Lab Sample ID: 280-176674-2  
Client ID: FWGmw-015-230401-GW  
Operator ID: JZ/JG ALS Bottle#: 43 Worklist Smp#: 43  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

4 HMX, CAS: 2691-41-0, Signal: 1

RT: 6.52  
Response: 1924031  
Amount: 20.571678



Reviewer: LV5D, 20-May-2023 11:05:10

Audit Action: Marked Compound Undetected

Audit Reason: Invalid Compound ID

FORM I  
HPLC/IC ORGANICS ANALYSIS DATA SHEET

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Client Sample ID: LL3mw-241-230401-GW Lab Sample ID: 280-176674-4

Matrix: Water Lab File ID: 05190046.D

Analysis Method: 8330B Date Collected: 05/17/2023 15:00

Extraction Method: 3535 Date Extracted: 05/19/2023 13:39

Sample wt/vol: 488.4 (mL) Date Analyzed: 05/20/2023 04:02

Con. Extract Vol.: 5 (mL) Dilution Factor: 1

Injection Volume: 100 (uL) GC Column: UltraCarb5uODS ID: 4.6 (mm)

% Moisture: \_\_\_\_\_ % Solids: \_\_\_\_\_ GPC Cleanup: (Y/N) N

Cleanup Factor: \_\_\_\_\_

Analysis Batch No.: 613168 Units: ug/L

| CAS NO.    | COMPOUND NAME              | RESULT | Q   | LOQ  | LOD   | DL    |
|------------|----------------------------|--------|-----|------|-------|-------|
| 99-35-4    | 1,3,5-Trinitrobenzene      | 4.6    |     | 0.21 | 0.20  | 0.086 |
| 118-96-7   | 2,4,6-Trinitrotoluene      | 2.0    |     | 0.11 | 0.10  | 0.046 |
| 6629-29-4  | 2,4-diamino-6-nitrotoluene | 0.92   | U   | 1.0  | 0.92  | 0.44  |
| 121-14-2   | 2,4-Dinitrotoluene         | 0.12   | J1  | 0.10 | 0.082 | 0.028 |
| 59229-75-3 | 2,6-diamino-4-nitrotoluene | 0.92   | U M | 1.0  | 0.92  | 0.23  |
| 606-20-2   | 2,6-Dinitrotoluene         | 0.082  | U   | 0.10 | 0.082 | 0.041 |
| 88-72-2    | 2-Nitrotoluene             | 0.20   | U   | 0.21 | 0.20  | 0.088 |
| 618-87-1   | 3,5-Dinitroaniline         | 0.36   | J   | 0.41 | 0.31  | 0.13  |
| 99-08-1    | 3-Nitrotoluene             | 0.36   | U   | 0.41 | 0.36  | 0.20  |
| 19406-51-0 | 4-Amino-2,6-dinitrotoluene | 1.7    |     | 0.15 | 0.12  | 0.059 |
| 99-99-0    | 4-Nitrotoluene             | 0.41   | U M | 0.42 | 0.41  | 0.10  |
| 80251-29-2 | DNX                        | 0.26   | U M | 0.51 | 0.26  | 0.10  |
| 2691-41-0  | HMX                        | 0.14   | J M | 0.21 | 0.20  | 0.090 |
| 98-95-3    | Nitrobenzene               | 0.20   | U M | 0.21 | 0.20  | 0.093 |
| 55-63-0    | Nitroglycerin              | 2.0    | U   | 2.1  | 2.0   | 0.94  |
| 78-11-5    | PETN                       | 1.0    | U   | 1.1  | 1.0   | 0.46  |
| 121-82-4   | RDX                        | 1.0    | J1  | 0.21 | 0.20  | 0.053 |
| 479-45-8   | Tetryl                     | 0.10   | U M | 0.11 | 0.10  | 0.033 |
| 13980-04-6 | TNX                        | 0.26   | U   | 0.51 | 0.26  | 0.081 |

| CAS NO.  | SURROGATE          | %REC | Q | LIMITS |
|----------|--------------------|------|---|--------|
| 528-29-0 | 1,2-Dinitrobenzene | 93   |   | 83-119 |

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\05190046.D  
 Lims ID: 280-176674-A-4-A  
 Client ID: LL3mw-241-230401-GW  
 Sample Type: Client  
 Inject. Date: 20-May-2023 04:02:40 ALS Bottle#: 46 Worklist Smp#: 46  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-176674-A-4-A  
 Operator ID: JZ/JG Instrument ID: CHHPLC\_X3  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 20-May-2023 11:22:13 Calib Date: 24-Feb-2023 17:41:20  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1663

First Level Reviewer: LV5D

Date: 20-May-2023 11:06:03

| Compound                      | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | OnCol Amt<br>ug/mL | Flags |
|-------------------------------|-----|--------------|------------------|------------------|----------|--------------------|-------|
| 3 TNX                         | 1   |              | 6.430            |                  |          | ND                 |       |
| 2 2,6-diamino-4-nitrotoluene  | 1   |              | 6.431            |                  |          | ND                 | U     |
| 4 HMX                         | 1   | 6.549        | 6.550            | -0.001           | 1277     | 0.0137             | M     |
| 5 2,4-diamino-6-nitrotoluene  | 1   |              | 6.618            |                  |          | ND                 |       |
| 6 DNx                         | 1   |              | 6.757            |                  |          | ND                 | U     |
| 7 MNx                         | 1   | 7.162        | 7.190            | -0.028           | 3472     | 0.0265             |       |
| 8 RDX                         | 1   | 7.569        | 7.570            | -0.001           | 10823    | 0.1017             |       |
| \$ 10 1,2-Dinitrobenzene      | 1   | 8.522        | 8.523            | -0.001           | 23545    | 0.1864             |       |
| 11 1,3,5-Trinitrobenzene      | 1   | 8.649        | 8.650            | -0.001           | 96646    | 0.4451             |       |
| 12 1,3-Dinitrobenzene         | 1   | 9.262        | 9.277            | -0.015           | 7703     | 0.0262             |       |
| 13 Nitrobenzene               | 1   |              | 9.643            |                  |          | ND                 | U     |
| 14 3,5-Dinitroaniline         | 1   | 9.882        | 9.883            | -0.001           | 7908     | 0.0352             |       |
| 15 Tetryl                     | 1   |              | 10.017           |                  |          | ND                 | U     |
| 16 Nitroglycerin              | 2   |              | 10.470           |                  |          | ND                 |       |
| 17 2,4,6-Trinitrotoluene      | 1   | 10.909       | 10.910           | -0.001           | 41127    | 0.1949             |       |
| 18 4-Amino-2,6-dinitrotoluene | 1   | 11.109       | 11.117           | -0.008           | 26325    | 0.1699             |       |
| 19 2-Amino-4,6-dinitrotoluene | 1   | 11.369       | 11.370           | -0.001           | 34224    | 0.1699             |       |
| 20 2,6-Dinitrotoluene         | 1   |              | 11.517           |                  |          | ND                 |       |
| 21 2,4-Dinitrotoluene         | 1   | 11.669       | 11.683           | -0.014           | 3431     | 0.0116             |       |
| 22 o-Nitrotoluene             | 1   |              | 12.517           |                  |          | ND                 |       |
| 23 p-Nitrotoluene             | 1   |              | 12.937           |                  |          | ND                 | U     |
| 24 m-Nitrotoluene             | 1   |              | 13.517           |                  |          | ND                 |       |
| 25 PETN                       | 2   |              | 14.697           |                  |          | ND                 |       |

### QC Flag Legend

Processing Flags

Review Flags

M - Manually Integrated

U - Marked Undetected

Report Date: 20-May-2023 11:22:13

Chrom Revision: 2.3 16-May-2023 16:37:30

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190046.d

Injection Date: 20-May-2023 04:02:40

Instrument ID: CHHPLC\_X3

Operator ID: JZ/JG

Lims ID: 280-176674-A-4-A

Lab Sample ID: 280-176674-4

Worklist Smp#: 46

Client ID: LL3mw-241-230401-GW

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 46

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1





Eurofins Denver  
Recovery Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\05190046.D  
Lims ID: 280-176674-A-4-A  
Client ID: LL3mw-241-230401-GW  
Sample Type: Client  
Inject. Date: 20-May-2023 04:02:40 ALS Bottle#: 46 Worklist Smp#: 46  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Sample Info: 280-176674-A-4-A  
Operator ID: JZ/JG Instrument ID: CHHPLC\_X3  
Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\8330\_X3.m  
Limit Group: GCSV - 8330  
Last Update: 20-May-2023 11:22:13 Calib Date: 24-Feb-2023 17:41:20  
Integrator: Falcon  
Quant Method: External Standard Quant By: Initial Calibration  
Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
Process Host: CTX1663

First Level Reviewer: LV5D

Date: 20-May-2023 11:06:03

| Compound                 | Amount Added | Amount Recovered | % Rec. |
|--------------------------|--------------|------------------|--------|
| \$ 10 1,2-Dinitrobenzene | 0.2000       | 0.1864           | 93.20  |


## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190046.d  
Injection Date: 20-May-2023 04:02:40 Instrument ID: CHHPLC\_X3  
Lims ID: 280-176674-A-4-A Lab Sample ID: 280-176674-4  
Client ID: LL3mw-241-230401-GW  
Operator ID: JZ/JG ALS Bottle#: 46 Worklist Smp#: 46  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**2,2,6-diamino-4-nitrotoluene, CAS: 59229-75-3, Signal: 1**

RT: 6.44  
Response: 190824  
Amount: 0.827016



Reviewer: LV5D, 20-May-2023 11:06:03

Audit Action: Marked Compound Undetected

Audit Reason: Invalid Compound ID

## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190046.d  
Injection Date: 20-May-2023 04:02:40 Instrument ID: CHHPLC\_X3  
Lims ID: 280-176674-A-4-A Lab Sample ID: 280-176674-4  
Client ID: LL3mw-241-230401-GW  
Operator ID: JZ/JG  
Injection Vol: 100.0 ul  
Method: 8330\_X3

ALS Bottle#: 46 Worklist Smp#: 46  
Dil. Factor: 1.0000  
Limit Group: GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**23 p-Nitrotoluene, CAS: 99-99-0, Signal: 1**

RT: 13.03  
Response: 1671  
Amount: 0.014936



Reviewer: LV5D, 20-May-2023 11:06:03

Audit Action: Marked Compound Undetected

Audit Reason: Invalid Compound ID


## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190046.d  
Injection Date: 20-May-2023 04:02:40 Instrument ID: CHHPLC\_X3  
Lims ID: 280-176674-A-4-A Lab Sample ID: 280-176674-4  
Client ID: LL3mw-241-230401-GW  
Operator ID: JZ/JG ALS Bottle#: 46 Worklist Smp#: 46  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

6 DNX, CAS: 80251-29-2, Signal: 1

RT: 6.69  
Response: 36624  
Amount: 0.253519



Reviewer: LV5D, 20-May-2023 11:06:03

Audit Action: Marked Compound Undetected

Audit Reason: Invalid Compound ID

## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190046.d  
Injection Date: 20-May-2023 04:02:40 Instrument ID: CHHPLC\_X3  
Lims ID: 280-176674-A-4-A Lab Sample ID: 280-176674-4  
Client ID: LL3mw-241-230401-GW  
Operator ID: JZ/JG ALS Bottle#: 46 Worklist Smp#: 46  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector: LC DAD1B, 254 nm

**4 HMX, CAS: 2691-41-0**

Signal: 1


RT: 6.55  
Area: 48805  
Amount: 0.521822  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.55  
Area: 1277  
Amount: 0.013654  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 11:05:44 -06:00:00 (UTC)

Audit Action: Manually Integrated

Audit Reason: Baseline


## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190046.d  
Injection Date: 20-May-2023 04:02:40 Instrument ID: CHHPLC\_X3  
Lims ID: 280-176674-A-4-A Lab Sample ID: 280-176674-4  
Client ID: LL3mw-241-230401-GW  
Operator ID: JZ/JG ALS Bottle#: 46 Worklist Smp#: 46  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**13 Nitrobenzene, CAS: 98-95-3, Signal: 1**

RT: 9.58  
Response: 3658  
Amount: 0.019127



Reviewer: LV5D, 20-May-2023 11:06:03

Audit Action: Marked Compound Undetected

Audit Reason: Invalid Compound ID


## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190046.d  
Injection Date: 20-May-2023 04:02:40 Instrument ID: CHHPLC\_X3  
Lims ID: 280-176674-A-4-A Lab Sample ID: 280-176674-4  
Client ID: LL3mw-241-230401-GW  
Operator ID: JZ/JG ALS Bottle#: 46 Worklist Smp#: 46  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**15 Tetryl, CAS: 479-45-8, Signal: 1**

RT: 10.08  
Response: 2214  
Amount: 0.013490



Reviewer: LV5D, 20-May-2023 11:06:03

Audit Action: Marked Compound Undetected

Audit Reason: Invalid Compound ID

FORM I  
HPLC/IC ORGANICS ANALYSIS DATA SHEET

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Client Sample ID: LL3mw-241-230401-GW Lab Sample ID: 280-176674-4

Matrix: Water Lab File ID: 05240010.D

Analysis Method: 8330B Date Collected: 05/17/2023 15:00

Extraction Method: 3535 Date Extracted: 05/19/2023 13:39

Sample wt/vol: 488.4 (mL) Date Analyzed: 05/24/2023 19:48

Con. Extract Vol.: 5 (mL) Dilution Factor: 1

Injection Volume: 100 (uL) GC Column: Luna-phenylhex ID: 4.6 (mm)

% Moisture: \_\_\_\_\_ % Solids: \_\_\_\_\_ GPC Cleanup: (Y/N) N

Cleanup Factor: \_\_\_\_\_

Analysis Batch No.: 613683 Units: ug/L

| CAS NO.    | COMPOUND NAME              | RESULT | Q    | LOQ  | LOD   | DL    |
|------------|----------------------------|--------|------|------|-------|-------|
| 99-65-0    | 1,3-Dinitrobenzene         | 0.10   | U    | 0.11 | 0.10  | 0.038 |
| 121-14-2   | 2,4-Dinitrotoluene         | 0.045  | J J1 | 0.10 | 0.082 | 0.028 |
| 35572-78-2 | 2-Amino-4,6-dinitrotoluene | 0.10   | U    | 0.11 | 0.10  | 0.052 |
| 5755-27-1  | MNX                        | 0.30   | U M  | 0.51 | 0.30  | 0.095 |
| 121-82-4   | RDX                        | 0.49   | J1   | 0.21 | 0.20  | 0.053 |

| CAS NO.  | SURROGATE          | %REC | Q | LIMITS |
|----------|--------------------|------|---|--------|
| 528-29-0 | 1,2-Dinitrobenzene | 88   |   | 83-119 |



Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230524-121801.b\05240010.D  
 Lims ID: 280-176674-A-4-A  
 Client ID: LL3mw-241-230401-GW  
 Sample Type: Client  
 Inject. Date: 24-May-2023 19:48:22 ALS Bottle#: 10 Worklist Smp#: 10  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-176674-A-4-A  
 Operator ID: JZ/JG Instrument ID: CHHPLC\_X5  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230524-121801.b\8330\_X5\_Luna.m  
 Limit Group: GCSV - 8330  
 Last Update: 25-May-2023 13:06:09 Calib Date: 29-Apr-2023 04:18:01  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280027.D  
 Column 1 : Luna-Phenyl hexyl ( 4.60 mm) Det: LC mwd1A, 254 nm  
 Process Host: CTX1654

First Level Reviewer: LV5D

Date: 24-May-2023 20:24:45

| Compound                      | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | OnCol Amt<br>ug/ml | Flags |
|-------------------------------|-----|--------------|------------------|------------------|----------|--------------------|-------|
| 1 2,6-diamino-4-nitrotoluene  | 1   |              | 4.214            |                  |          | ND                 |       |
| 2 2,4-diamino-6-nitrotoluene  | 1   |              | 4.767            |                  |          | ND                 |       |
| 3 TNX                         | 1   |              | 5.206            |                  |          | ND                 |       |
| 4 DNX                         | 1   |              | 6.046            |                  |          | ND                 |       |
| 5 HMX                         | 1   | 6.731        | 6.686            | 0.045            | 3536     | 0.0146             | M     |
| 6 MNX                         | 1   | 7.557        | 7.546            | 0.011            | 2363     | 0.009010           | M     |
| 8 RDX                         | 1   | 8.971        | 8.946            | 0.025            | 10969    | 0.0477             |       |
| 9 Nitrobenzene                | 1   |              | 11.459           |                  |          | ND                 |       |
| \$ 10 1,2-Dinitrobenzene      | 1   | 12.437       | 12.426           | 0.011            | 49455    | 0.1762             |       |
| 11 3,5-Dinitroaniline         | 1   | 14.277       | 14.292           | -0.015           | 16722    | 0.0362             |       |
| 12 1,3-Dinitrobenzene         | 1   |              | 14.552           |                  |          | ND                 |       |
| 13 Nitroglycerin              | 2   |              | 15.066           |                  |          | ND                 |       |
| 14 o-Nitrotoluene             | 1   | 15.511       | 15.619           | -0.108           | 3979     | 0.0161             |       |
| 16 p-Nitrotoluene             | 1   |              | 15.846           |                  |          | ND                 |       |
| 17 4-Amino-2,6-dinitrotoluene | 1   | 16.371       | 16.406           | -0.035           | 50092    | 0.1703             |       |
| 18 m-Nitrotoluene             | 1   | 16.764       | 16.692           | 0.072            | 34294    | 0.1164             |       |
| 19 2-Amino-4,6-dinitrotoluene | 1   |              | 17.239           |                  |          | ND                 |       |
| 20 1,3,5-Trinitrobenzene      | 1   | 17.337       | 17.332           | 0.005            | 241842   | 0.5961             |       |
| 21 2,6-Dinitrotoluene         | 1   |              | 18.492           |                  |          | ND                 |       |
| 22 2,4-Dinitrotoluene         | 1   | 18.911       | 18.952           | -0.041           | 2584     | 0.004421           |       |
| 23 Tetryl                     | 1   |              | 22.239           |                  |          | ND                 | 7     |
| 24 2,4,6-Trinitrotoluene      | 1   | 23.011       | 23.039           | -0.028           | 78556    | 0.1790             |       |
| 25 PETN                       | 2   |              | 24.279           |                  |          | ND                 |       |

**QC Flag Legend**

Processing Flags

7 - Failed Limit of Detection

Review Flags

M - Manually Integrated

Report Date: 25-May-2023 13:06:10

Chrom Revision: 2.3 23-May-2023 13:55:56

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230524-121801.b\05240010.D

Injection Date: 24-May-2023 19:48:22

Instrument ID: CHHPLC\_X5

Operator ID: JZ/JG

Lims ID: 280-176674-A-4-A

Lab Sample ID: 280-176674-4

Worklist Smp#: 10

Client ID: LL3mw-241-230401-GW

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 10

Method: 8330\_X5\_Luna

Limit Group: GCSV - 8330


Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



Eurofins Denver  
Recovery Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230524-121801.b\05240010.D  
Lims ID: 280-176674-A-4-A  
Client ID: LL3mw-241-230401-GW  
Sample Type: Client  
Inject. Date: 24-May-2023 19:48:22 ALS Bottle#: 10 Worklist Smp#: 10  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Sample Info: 280-176674-A-4-A  
Operator ID: JZ/JG Instrument ID: CHHPLC\_X5  
Method: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230524-121801.b\8330\_X5\_Luna.m  
Limit Group: GCSV - 8330  
Last Update: 25-May-2023 13:06:09 Calib Date: 29-Apr-2023 04:18:01  
Integrator: Falcon  
Quant Method: External Standard Quant By: Initial Calibration  
Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280027.D  
Column 1 : Luna-Phenyl hexyl ( 4.60 mm) Det: LC mwd1A, 254 nm  
Process Host: CTX1654

First Level Reviewer: LV5D

Date: 24-May-2023 20:24:45

| Compound                 | Amount Added | Amount Recovered | % Rec. |
|--------------------------|--------------|------------------|--------|
| \$ 10 1,2-Dinitrobenzene | 0.2000       | 0.1762           | 88.09  |

## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230524-121801.b\05240010.D  
Injection Date: 24-May-2023 19:48:22 Instrument ID: CHHPLC\_X5  
Lims ID: 280-176674-A-4-A Lab Sample ID: 280-176674-4  
Client ID: LL3mw-241-230401-GW  
Operator ID: JZ/JG ALS Bottle#: 10 Worklist Smp#: 10  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X5\_Luna Limit Group: GCSV - 8330  
Column: Luna-Phenyl hexyl ( 4.60 mm) Detector LC mwd1A, 254 nm

**5 HMX, CAS: 2691-41-0**

Signal: 1


RT: 6.73  
Area: 8111  
Amount: 0.039256  
Amount Units: ug/ml

## Processing Integration Results



RT: 6.73  
Area: 3536  
Amount: 0.014583  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: K8YG, 25-May-2023 08:40:59 -06:00:00 (UTC)

Audit Action: Manually Integrated

Audit Reason: Baseline

## Eurofins Denver


|                 |                                                                   |                |                  |
|-----------------|-------------------------------------------------------------------|----------------|------------------|
| Data File:      | \\chromfs\Denver\ChromData\CHHPLC_X5\20230524-121801.b\05240010.D |                |                  |
| Injection Date: | 24-May-2023 19:48:22                                              | Instrument ID: | CHHPLC_X5        |
| Lims ID:        | 280-176674-A-4-A                                                  | Lab Sample ID: | 280-176674-4     |
| Client ID:      | LL3mw-241-230401-GW                                               |                |                  |
| Operator ID:    | JZ/JG                                                             | ALS Bottle#:   | 10               |
| Injection Vol:  | 100.0 ul                                                          | Dil. Factor:   | 1.0000           |
| Method:         | 8330_X5_Luna                                                      | Limit Group:   | GCSV - 8330      |
| Column:         | Luna-Phenyl hexyl ( 4.60 mm)                                      | Detector       | LC mwd1A, 254 nm |
|                 |                                                                   | Worklist Smp#: | 10               |

**6 MNX, CAS: 5755-27-1**

Signal: 1


RT: 7.56  
Area: 7992  
Amount: 0.030473  
Amount Units: ug/ml

## Processing Integration Results



RT: 7.56  
Area: 2363  
Amount: 0.009010  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: K8YG, 25-May-2023 08:41:04 -06:00:00 (UTC)

Audit Action: Manually Integrated

Audit Reason: Baseline

FORM I  
HPLC/IC ORGANICS ANALYSIS DATA SHEET

|                                                |                                                      |
|------------------------------------------------|------------------------------------------------------|
| Lab Name: <u>Eurofins Denver</u>               | Job No.: <u>280-176674-1</u>                         |
| SDG No.: _____                                 |                                                      |
| Client Sample ID: <u>FBQmw-173-230401-GW-R</u> | Lab Sample ID: <u>280-176674-7</u>                   |
| Matrix: <u>Water</u>                           | Lab File ID: <u>05190047.D</u>                       |
| Analysis Method: <u>8330B</u>                  | Date Collected: <u>05/17/2023 11:47</u>              |
| Extraction Method: <u>3535</u>                 | Date Extracted: <u>05/19/2023 13:39</u>              |
| Sample wt/vol: <u>483.2 (mL)</u>               | Date Analyzed: <u>05/20/2023 04:25</u>               |
| Con. Extract Vol.: <u>5 (mL)</u>               | Dilution Factor: <u>1</u>                            |
| Injection Volume: <u>100 (uL)</u>              | GC Column: <u>UltraCarb5uODS</u> ID: <u>4.6 (mm)</u> |
| % Moisture: _____ % Solids: _____              | GPC Cleanup: (Y/N) <u>N</u>                          |
| Cleanup Factor: _____                          |                                                      |
| Analysis Batch No.: <u>613168</u>              | Units: <u>ug/L</u>                                   |

| CAS NO.    | COMPOUND NAME              | RESULT | Q   | LOQ  | LOD   | DL    |
|------------|----------------------------|--------|-----|------|-------|-------|
| 99-35-4    | 1,3,5-Trinitrobenzene      | 0.21   | U   | 0.22 | 0.21  | 0.087 |
| 118-96-7   | 2,4,6-Trinitrotoluene      | 0.11   |     | 0.11 | 0.10  | 0.047 |
| 6629-29-4  | 2,4-diamino-6-nitrotoluene | 0.93   | U   | 1.0  | 0.93  | 0.45  |
| 121-14-2   | 2,4-Dinitrotoluene         | 0.083  | U   | 0.10 | 0.083 | 0.028 |
| 59229-75-3 | 2,6-diamino-4-nitrotoluene | 0.93   | U   | 1.0  | 0.93  | 0.23  |
| 606-20-2   | 2,6-Dinitrotoluene         | 0.083  | U   | 0.10 | 0.083 | 0.041 |
| 35572-78-2 | 2-Amino-4,6-dinitrotoluene | 1.2    |     | 0.11 | 0.10  | 0.052 |
| 88-72-2    | 2-Nitrotoluene             | 0.21   | U   | 0.22 | 0.21  | 0.088 |
| 618-87-1   | 3,5-Dinitroaniline         | 0.15   | J   | 0.41 | 0.31  | 0.14  |
| 99-08-1    | 3-Nitrotoluene             | 0.36   | U   | 0.41 | 0.36  | 0.20  |
| 19406-51-0 | 4-Amino-2,6-dinitrotoluene | 1.5    |     | 0.16 | 0.12  | 0.060 |
| 99-99-0    | 4-Nitrotoluene             | 0.41   | U   | 0.42 | 0.41  | 0.10  |
| 80251-29-2 | DNX                        | 0.26   | U   | 0.52 | 0.26  | 0.10  |
| 2691-41-0  | HMX                        | 0.21   | U   | 0.22 | 0.21  | 0.091 |
| 5755-27-1  | MNX                        | 0.30   | U   | 0.52 | 0.30  | 0.096 |
| 98-95-3    | Nitrobenzene               | 0.21   | U   | 0.22 | 0.21  | 0.094 |
| 55-63-0    | Nitroglycerin              | 2.1    | U   | 2.2  | 2.1   | 0.95  |
| 78-11-5    | PETN                       | 1.0    | U   | 1.1  | 1.0   | 0.46  |
| 121-82-4   | RDX                        | 0.21   | U M | 0.22 | 0.21  | 0.053 |
| 479-45-8   | Tetryl                     | 0.10   | U   | 0.11 | 0.10  | 0.033 |
| 13980-04-6 | TNX                        | 0.26   | U M | 0.52 | 0.26  | 0.082 |

| CAS NO.  | SURROGATE          | %REC | Q | LIMITS |
|----------|--------------------|------|---|--------|
| 528-29-0 | 1,2-Dinitrobenzene | 91   | M | 83-119 |

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\05190047.D  
 Lims ID: 280-176674-C-7-A  
 Client ID: FBQmw-173-230401-GW-R  
 Sample Type: Client  
 Inject. Date: 20-May-2023 04:25:37 ALS Bottle#: 47 Worklist Smp#: 47  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-176674-C-7-A  
 Operator ID: JZ/JG Instrument ID: CHHPLC\_X3  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 20-May-2023 11:22:13 Calib Date: 24-Feb-2023 17:41:20  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1663

First Level Reviewer: LV5D

Date: 20-May-2023 11:06:24

| Compound                      | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | OnCol Amt<br>ug/mL | Flags |
|-------------------------------|-----|--------------|------------------|------------------|----------|--------------------|-------|
| 3 TNX                         | 1   |              | 6.430            |                  |          | ND                 | U     |
| 2 2,6-diamino-4-nitrotoluene  | 1   |              | 6.431            |                  |          | ND                 |       |
| 4 HMX                         | 1   |              | 6.550            |                  |          | ND                 |       |
| 5 2,4-diamino-6-nitrotoluene  | 1   |              | 6.618            |                  |          | ND                 |       |
| 6 DNX                         | 1   |              | 6.757            |                  |          | ND                 |       |
| 7 MNX                         | 1   |              | 7.190            |                  |          | ND                 |       |
| 8 RDX                         | 1   |              | 7.570            |                  |          | ND                 | U     |
| \$ 10 1,2-Dinitrobenzene      | 1   | 8.527        | 8.523            | 0.004            | 23003    | 0.1821             | M     |
| 11 1,3,5-Trinitrobenzene      | 1   |              | 8.650            |                  |          | ND                 |       |
| 12 1,3-Dinitrobenzene         | 1   | 9.273        | 9.277            | -0.004           | 1111     | 0.003774           | 7M    |
| 13 Nitrobenzene               | 1   |              | 9.643            |                  |          | ND                 | 7     |
| 14 3,5-Dinitroaniline         | 1   | 9.887        | 9.883            | 0.004            | 3120     | 0.0143             |       |
| 15 Tetryl                     | 1   |              | 10.017           |                  |          | ND                 |       |
| 16 Nitroglycerin              | 2   |              | 10.470           |                  |          | ND                 |       |
| 17 2,4,6-Trinitrotoluene      | 1   | 10.907       | 10.910           | -0.003           | 2156     | 0.0102             |       |
| 18 4-Amino-2,6-dinitrotoluene | 1   | 11.113       | 11.117           | -0.004           | 22844    | 0.1474             |       |
| 19 2-Amino-4,6-dinitrotoluene | 1   | 11.367       | 11.370           | -0.003           | 23618    | 0.1173             |       |
| 20 2,6-Dinitrotoluene         | 1   |              | 11.517           |                  |          | ND                 |       |
| 21 2,4-Dinitrotoluene         | 1   |              | 11.683           |                  |          | ND                 | 7     |
| 22 o-Nitrotoluene             | 1   |              | 12.517           |                  |          | ND                 |       |
| 23 p-Nitrotoluene             | 1   |              | 12.937           |                  |          | ND                 |       |
| 24 m-Nitrotoluene             | 1   |              | 13.517           |                  |          | ND                 |       |
| 25 PETN                       | 2   |              | 14.697           |                  |          | ND                 |       |

**QC Flag Legend**

Processing Flags

7 - Failed Limit of Detection

Review Flags

M - Manually Integrated

U - Marked Undetected



Report Date: 20-May-2023 11:22:14

Chrom Revision: 2.3 16-May-2023 16:37:30

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190047.d

Injection Date: 20-May-2023 04:25:37

Instrument ID: CHHPLC\_X3

Operator ID: JZ/JG

Lims ID: 280-176674-C-7-A

Lab Sample ID: 280-176674-7

Worklist Smp#: 47

Client ID: FBQmw-173-230401-GW-R

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 47

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



Eurofins Denver  
Recovery Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\05190047.D  
Lims ID: 280-176674-C-7-A  
Client ID: FBQmw-173-230401-GW-R  
Sample Type: Client  
Inject. Date: 20-May-2023 04:25:37 ALS Bottle#: 47 Worklist Smp#: 47  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Sample Info: 280-176674-C-7-A  
Operator ID: JZ/JG Instrument ID: CHHPLC\_X3  
Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\8330\_X3.m  
Limit Group: GCSV - 8330  
Last Update: 20-May-2023 11:22:13 Calib Date: 24-Feb-2023 17:41:20  
Integrator: Falcon  
Quant Method: External Standard Quant By: Initial Calibration  
Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
Process Host: CTX1663

First Level Reviewer: LV5D

Date: 20-May-2023 11:06:24

| Compound                 | Amount Added | Amount Recovered | % Rec. |
|--------------------------|--------------|------------------|--------|
| \$ 10 1,2-Dinitrobenzene | 0.2000       | 0.1821           | 91.06  |

## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190047.d  
Injection Date: 20-May-2023 04:25:37 Instrument ID: CHHPLC\_X3  
Lims ID: 280-176674-C-7-A Lab Sample ID: 280-176674-7  
Client ID: FBQmw-173-230401-GW-R  
Operator ID: JZ/JG ALS Bottle#: 47 Worklist Smp#: 47  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector: LC DAD1B, 254 nm

## \$ 10 1,2-Dinitrobenzene, CAS: 528-29-0

Signal: 1


RT: 8.53  
Area: 25403  
Amount: 0.201118  
Amount Units: ug/mL

## Processing Integration Results



RT: 8.53  
Area: 23003  
Amount: 0.182117  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 11:06:13 -06:00:00 (UTC)

Audit Action: Manually Integrated

Audit Reason: Baseline

## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190047.d  
Injection Date: 20-May-2023 04:25:37 Instrument ID: CHHPLC\_X3  
Lims ID: 280-176674-C-7-A Lab Sample ID: 280-176674-7  
Client ID: FBQmw-173-230401-GW-R  
Operator ID: JZ/JG ALS Bottle#: 47 Worklist Smp#: 47  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector: LC DAD1B, 254 nm

**12 1,3-Dinitrobenzene, CAS: 99-65-0**

Signal: 1


RT: 9.27  
Area: 1740  
Amount: 0.005910  
Amount Units: ug/mL

## Processing Integration Results



RT: 9.27  
Area: 1111  
Amount: 0.003774  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 11:06:19 -06:00:00 (UTC)

Audit Action: Split an Integrated Peak


Audit Reason: Baseline

## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190047.d  
Injection Date: 20-May-2023 04:25:37 Instrument ID: CHHPLC\_X3  
Lims ID: 280-176674-C-7-A Lab Sample ID: 280-176674-7  
Client ID: FBQmw-173-230401-GW-R  
Operator ID: JZ/JG ALS Bottle#: 47 Worklist Smp#: 47  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**8 RDX, CAS: 121-82-4, Signal: 1**

RT: 7.69  
Response: 4608  
Amount: 0.043316



Reviewer: LV5D, 20-May-2023 11:06:24

Audit Action: Marked Compound Undetected


Audit Reason: Invalid Compound ID

## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190047.d  
Injection Date: 20-May-2023 04:25:37 Instrument ID: CHHPLC\_X3  
Lims ID: 280-176674-C-7-A Lab Sample ID: 280-176674-7  
Client ID: FBQmw-173-230401-GW-R  
Operator ID: JZ/JG ALS Bottle#: 47 Worklist Smp#: 47  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

3 TNX, CAS: 13980-04-6, Signal: 1

RT: 6.43  
Response: 140600  
Amount: 0.711224



Reviewer: LV5D, 20-May-2023 11:06:24

Audit Action: Marked Compound Undetected

Audit Reason: Invalid Compound ID

FORM I  
HPLC/IC ORGANICS ANALYSIS DATA SHEET

Lab Name: Eurofins Denver Job No.: 280-176674-1  
SDG No.: \_\_\_\_\_  
Client Sample ID: FBQmw-173-230401-GW-R Lab Sample ID: 280-176674-7  
Matrix: Water Lab File ID: 05200019.D  
Analysis Method: 8330B Date Collected: 05/17/2023 11:47  
Extraction Method: 3535 Date Extracted: 05/19/2023 13:39  
Sample wt/vol: 483.2 (mL) Date Analyzed: 05/21/2023 00:31  
Con. Extract Vol.: 5 (mL) Dilution Factor: 1  
Injection Volume: 100 (uL) GC Column: Luna-phenylhex ID: 4.6 (mm)  
% Moisture: \_\_\_\_\_ % Solids: \_\_\_\_\_ GPC Cleanup: (Y/N) N  
Cleanup Factor: \_\_\_\_\_  
Analysis Batch No.: 613238 Units: ug/L

| CAS NO. | COMPOUND NAME      | RESULT | Q | LOQ  | LOD  | DL    |
|---------|--------------------|--------|---|------|------|-------|
| 99-65-0 | 1,3-Dinitrobenzene | 0.10   | U | 0.11 | 0.10 | 0.038 |

| CAS NO.  | SURROGATE          | %REC | Q | LIMITS |
|----------|--------------------|------|---|--------|
| 528-29-0 | 1,2-Dinitrobenzene | 85   |   | 83-119 |

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230520-121670.b\05200019.D  
 Lims ID: 280-176674-C-7-A  
 Client ID: FBQmw-173-230401-GW-R  
 Sample Type: Client  
 Inject. Date: 21-May-2023 00:31:28 ALS Bottle#: 19 Worklist Smp#: 19  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-176674-C-7-A  
 Operator ID: JZ/JG Instrument ID: CHHPLC\_X5  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230520-121670.b\8330\_X5\_Luna.m  
 Limit Group: GCSV - 8330  
 Last Update: 23-May-2023 13:12:40 Calib Date: 29-Apr-2023 04:18:01  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280027.D  
 Column 1 : Luna-Phenyl hexyl ( 4.60 mm) Det: LC mwd1A, 254 nm  
 Process Host: CTX1626

First Level Reviewer: K8YG

Date: 22-May-2023 11:24:46

| Compound                      | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | OnCol Amt<br>ug/ml | Flags |
|-------------------------------|-----|--------------|------------------|------------------|----------|--------------------|-------|
| 1 2,6-diamino-4-nitrotoluene  | 1   |              | 4.214            |                  |          | ND                 | 7     |
| 2 2,4-diamino-6-nitrotoluene  | 1   |              | 4.767            |                  |          | ND                 |       |
| 3 TNX                         | 1   |              | 5.221            |                  |          | ND                 |       |
| 4 DNX                         | 1   |              | 6.067            |                  |          | ND                 |       |
| 5 HMX                         | 1   |              | 6.758            |                  |          | ND                 |       |
| 6 MNX                         | 1   |              | 7.567            |                  |          | ND                 |       |
| 8 RDX                         | 1   |              | 9.005            |                  |          | ND                 |       |
| 9 Nitrobenzene                | 1   |              | 11.531           |                  |          | ND                 |       |
| \$ 10 1,2-Dinitrobenzene      | 1   | 12.484       | 12.511           | -0.027           | 47511    | 0.1693             |       |
| 11 3,5-Dinitroaniline         | 1   | 14.344       | 14.378           | -0.034           | 6984     | 0.0143             |       |
| 12 1,3-Dinitrobenzene         | 1   |              | 14.638           |                  |          | ND                 |       |
| 13 Nitroglycerin              | 2   |              | 15.158           |                  |          | ND                 |       |
| 14 o-Nitrotoluene             | 1   | 15.711       | 15.718           | -0.007           | 2921     | 0.0118             |       |
| 16 p-Nitrotoluene             | 1   |              | 15.945           |                  |          | ND                 |       |
| 17 4-Amino-2,6-dinitrotoluene | 1   | 16.431       | 16.505           | -0.074           | 49909    | 0.1696             |       |
| 18 m-Nitrotoluene             | 1   |              | 16.798           |                  |          | ND                 |       |
| 19 2-Amino-4,6-dinitrotoluene | 1   | 17.244       | 17.351           | -0.107           | 54929    | 0.1284             |       |
| 20 1,3,5-Trinitrobenzene      | 1   |              | 17.398           |                  |          | ND                 |       |
| 21 2,6-Dinitrotoluene         | 1   |              | 18.605           |                  |          | ND                 | U     |
| 22 2,4-Dinitrotoluene         | 1   |              | 19.051           |                  |          | ND                 | 7     |
| 23 Tetryl                     | 1   |              | 22.331           |                  |          | ND                 |       |
| 24 2,4,6-Trinitrotoluene      | 1   | 23.077       | 23.125           | -0.048           | 4488     | 0.0102             |       |
| 25 PETN                       | 2   |              | 24.365           |                  |          | ND                 |       |

**QC Flag Legend**

Processing Flags

7 - Failed Limit of Detection

Review Flags

U - Marked Undetected



Report Date: 23-May-2023 13:12:53

Chrom Revision: 2.3 16-May-2023 16:37:30

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230520-121670.b\05200019.D

Injection Date: 21-May-2023 00:31:28

Instrument ID: CHHPLC\_X5

Operator ID: JZ/JG

Lims ID: 280-176674-C-7-A

Lab Sample ID: 280-176674-7

Worklist Smp#: 19

Client ID: FBQmw-173-230401-GW-R

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 19

Method: 8330\_X5\_Luna


Limit Group: GCSV - 8330

Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



Eurofins Denver  
Recovery Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230520-121670.b\05200019.D  
Lims ID: 280-176674-C-7-A  
Client ID: FBQmw-173-230401-GW-R  
Sample Type: Client  
Inject. Date: 21-May-2023 00:31:28 ALS Bottle#: 19 Worklist Smp#: 19  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Sample Info: 280-176674-C-7-A  
Operator ID: JZ/JG Instrument ID: CHHPLC\_X5  
Method: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230520-121670.b\8330\_X5\_Luna.m  
Limit Group: GCSV - 8330  
Last Update: 23-May-2023 13:12:40 Calib Date: 29-Apr-2023 04:18:01  
Integrator: Falcon  
Quant Method: External Standard Quant By: Initial Calibration  
Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280027.D  
Column 1 : Luna-Phenyl hexyl ( 4.60 mm) Det: LC mwd1A, 254 nm  
Process Host: CTX1626

First Level Reviewer: K8YG

Date: 22-May-2023 11:24:46


| Compound                 | Amount Added | Amount Recovered | % Rec. |
|--------------------------|--------------|------------------|--------|
| \$ 10 1,2-Dinitrobenzene | 0.2000       | 0.1693           | 84.63  |

## Eurofins Denver

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230520-121670.b\05200019.D  
Injection Date: 21-May-2023 00:31:28 Instrument ID: CHHPLC\_X5  
Lims ID: 280-176674-C-7-A Lab Sample ID: 280-176674-7  
Client ID: FBQmw-173-230401-GW-R  
Operator ID: JZ/JG ALS Bottle#: 19 Worklist Smp#: 19  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X5\_Luna Limit Group: GCSV - 8330  
Column: Luna-Phenyl hexyl ( 4.60 mm) Detector LC mwd1A, 254 nm

**21 2,6-Dinitrotoluene, CAS: 606-20-2, Signal: 1**

RT: 18.48  
Response: 1888  
Amount: 0.006473



Reviewer: LV5D, 23-May-2023 13:11:55

Audit Action: Marked Compound Undetected

Audit Reason: Invalid Compound ID

FORM I  
HPLC/IC ORGANICS ANALYSIS DATA SHEET

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Client Sample ID: FBQmw-175-230401-GW-R Lab Sample ID: 280-176674-10

Matrix: Water Lab File ID: 05190048.D

Analysis Method: 8330B Date Collected: 05/17/2023 08:58

Extraction Method: 3535 Date Extracted: 05/19/2023 13:39

Sample wt/vol: 497.5 (mL) Date Analyzed: 05/20/2023 04:48

Con. Extract Vol.: 5 (mL) Dilution Factor: 1

Injection Volume: 100 (uL) GC Column: UltraCarb5uODS ID: 4.6 (mm)

% Moisture: \_\_\_\_\_ % Solids: \_\_\_\_\_ GPC Cleanup: (Y/N) N

Cleanup Factor: \_\_\_\_\_

Analysis Batch No.: 613168 Units: ug/L

| CAS NO.    | COMPOUND NAME              | RESULT | Q         | LOQ  | LOD   | DL    |
|------------|----------------------------|--------|-----------|------|-------|-------|
| 99-35-4    | 1,3,5-Trinitrobenzene      | 0.20   | U J1      | 0.21 | 0.20  | 0.085 |
| 99-65-0    | 1,3-Dinitrobenzene         | 0.10   | U J1      | 0.11 | 0.10  | 0.037 |
| 118-96-7   | 2,4,6-Trinitrotoluene      | 0.10   | U J1      | 0.11 | 0.10  | 0.045 |
| 6629-29-4  | 2,4-diamino-6-nitrotoluene | 0.90   | U J1      | 1.0  | 0.90  | 0.44  |
| 121-14-2   | 2,4-Dinitrotoluene         | 0.080  | U J1      | 0.10 | 0.080 | 0.028 |
| 59229-75-3 | 2,6-diamino-4-nitrotoluene | 0.90   | U J1      | 1.0  | 0.90  | 0.22  |
| 606-20-2   | 2,6-Dinitrotoluene         | 0.080  | U J1      | 0.10 | 0.080 | 0.040 |
| 35572-78-2 | 2-Amino-4,6-dinitrotoluene | 0.10   | U J1      | 0.11 | 0.10  | 0.051 |
| 88-72-2    | 2-Nitrotoluene             | 0.20   | U J1      | 0.21 | 0.20  | 0.086 |
| 618-87-1   | 3,5-Dinitroaniline         | 0.30   | U J1      | 0.40 | 0.30  | 0.13  |
| 99-08-1    | 3-Nitrotoluene             | 0.35   | U J1      | 0.40 | 0.35  | 0.20  |
| 19406-51-0 | 4-Amino-2,6-dinitrotoluene | 0.12   | U J1      | 0.15 | 0.12  | 0.058 |
| 99-99-0    | 4-Nitrotoluene             | 0.40   | U J1      | 0.41 | 0.40  | 0.10  |
| 80251-29-2 | DNX                        | 0.25   | U J1      | 0.50 | 0.25  | 0.098 |
| 2691-41-0  | HMX                        | 0.20   | U         | 0.21 | 0.20  | 0.088 |
| 5755-27-1  | MNX                        | 0.29   | U         | 0.50 | 0.29  | 0.093 |
| 98-95-3    | Nitrobenzene               | 0.20   | U J1      | 0.21 | 0.20  | 0.091 |
| 55-63-0    | Nitroglycerin              | 2.0    | U         | 2.1  | 2.0   | 0.93  |
| 78-11-5    | PETN                       | 1.0    | U J1      | 1.1  | 1.0   | 0.45  |
| 121-82-4   | RDX                        | 0.20   | U M<br>J1 | 0.21 | 0.20  | 0.052 |
| 479-45-8   | Tetryl                     | 0.10   | U J1      | 0.11 | 0.10  | 0.032 |
| 13980-04-6 | TNX                        | 0.25   | U M       | 0.50 | 0.25  | 0.080 |

| CAS NO.  | SURROGATE          | %REC | Q | LIMITS |
|----------|--------------------|------|---|--------|
| 528-29-0 | 1,2-Dinitrobenzene | 95   | M | 83-119 |

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\05190048.D  
 Lims ID: 280-176674-A-10-A  
 Client ID: FBQmw-175-230401-GW-R  
 Sample Type: Client  
 Inject. Date: 20-May-2023 04:48:36 ALS Bottle#: 48 Worklist Smp#: 48  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-176674-A-10-A  
 Operator ID: JZ/JG Instrument ID: CHHPLC\_X3  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 20-May-2023 11:22:13 Calib Date: 24-Feb-2023 17:41:20  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1663

First Level Reviewer: LV5D

Date: 20-May-2023 11:14:43

| Compound                      | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | OnCol Amt<br>ug/mL | Flags |
|-------------------------------|-----|--------------|------------------|------------------|----------|--------------------|-------|
| 1 Triamine Trinitrobenzene    | 1   |              | 2.444            |                  |          | ND                 |       |
| 3 TNX                         | 1   |              | 6.430            |                  |          | ND                 | U     |
| 2 2,6-diamino-4-nitrotoluene  | 1   |              | 6.431            |                  |          | ND                 |       |
| 4 HMX                         | 1   |              | 6.550            |                  |          | ND                 |       |
| 5 2,4-diamino-6-nitrotoluene  | 1   |              | 6.618            |                  |          | ND                 |       |
| 6 DNX                         | 1   |              | 6.757            |                  |          | ND                 |       |
| 7 MNX                         | 1   |              | 7.190            |                  |          | ND                 |       |
| 8 RDX                         | 1   |              | 7.570            |                  |          | ND                 | U     |
| 9 2,4,6-Trinitrophenol        | 1   |              | 7.963            |                  |          | ND                 |       |
| \$ 10 1,2-Dinitrobenzene      | 1   | 8.527        | 8.523            | 0.004            | 23925    | 0.1894             | M     |
| 11 1,3,5-Trinitrobenzene      | 1   |              | 8.650            |                  |          | ND                 |       |
| 12 1,3-Dinitrobenzene         | 1   |              | 9.277            |                  |          | ND                 |       |
| 13 Nitrobenzene               | 1   |              | 9.643            |                  |          | ND                 |       |
| 14 3,5-Dinitroaniline         | 1   |              | 9.883            |                  |          | ND                 |       |
| 15 Tetryl                     | 1   |              | 10.017           |                  |          | ND                 |       |
| 16 Nitroglycerin              | 2   |              | 10.470           |                  |          | ND                 |       |
| 17 2,4,6-Trinitrotoluene      | 1   |              | 10.910           |                  |          | ND                 |       |
| 18 4-Amino-2,6-dinitrotoluene | 1   |              | 11.117           |                  |          | ND                 |       |
| 19 2-Amino-4,6-dinitrotoluene | 1   |              | 11.370           |                  |          | ND                 |       |
| 20 2,6-Dinitrotoluene         | 1   |              | 11.517           |                  |          | ND                 |       |
| 21 2,4-Dinitrotoluene         | 1   |              | 11.683           |                  |          | ND                 |       |
| 22 o-Nitrotoluene             | 1   |              | 12.517           |                  |          | ND                 |       |
| 23 p-Nitrotoluene             | 1   |              | 12.937           |                  |          | ND                 |       |
| 24 m-Nitrotoluene             | 1   |              | 13.517           |                  |          | ND                 |       |
| 25 PETN                       | 2   |              | 14.697           |                  |          | ND                 |       |
| 26 Ammonium Picrate           | 1   |              | 0.000            |                  |          | ND                 |       |

### QC Flag Legend

Processing Flags

Review Flags

M - Manually Integrated

U - Marked Undetected

Report Date: 20-May-2023 11:22:14

Chrom Revision: 2.3 16-May-2023 16:37:30

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190048.d

Injection Date: 20-May-2023 04:48:36

Instrument ID: CHHPLC\_X3

Operator ID: JZ/JG

Lims ID: 280-176674-A-10-A

Lab Sample ID: 280-176674-10

Worklist Smp#: 48

Client ID: FBQmw-175-230401-GW-R

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 48

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



Eurofins Denver  
Recovery Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\05190048.D  
Lims ID: 280-176674-A-10-A  
Client ID: FBQmw-175-230401-GW-R  
Sample Type: Client  
Inject. Date: 20-May-2023 04:48:36 ALS Bottle#: 48 Worklist Smp#: 48  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Sample Info: 280-176674-A-10-A  
Operator ID: JZ/JG Instrument ID: CHHPLC\_X3  
Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\8330\_X3.m  
Limit Group: GCSV - 8330  
Last Update: 20-May-2023 11:22:13 Calib Date: 24-Feb-2023 17:41:20  
Integrator: Falcon  
Quant Method: External Standard Quant By: Initial Calibration  
Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
Process Host: CTX1663

First Level Reviewer: LV5D

Date: 20-May-2023 11:14:43

| Compound                 | Amount Added | Amount Recovered | % Rec. |
|--------------------------|--------------|------------------|--------|
| \$ 10 1,2-Dinitrobenzene | 0.2000       | 0.1894           | 94.71  |



## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190048.d  
Injection Date: 20-May-2023 04:48:36 Instrument ID: CHHPLC\_X3  
Lims ID: 280-176674-A-10-A Lab Sample ID: 280-176674-10  
Client ID: FBQmw-175-230401-GW-R  
Operator ID: JZ/JG ALS Bottle#: 48 Worklist Smp#: 48  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector: LC DAD1B, 254 nm

## \$ 10 1,2-Dinitrobenzene, CAS: 528-29-0

Signal: 1


RT: 8.53  
Area: 28479  
Amount: 0.225471  
Amount Units: ug/mL

## Processing Integration Results



RT: 8.53  
Area: 23925  
Amount: 0.189417  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 11:14:42 -06:00:00 (UTC)

Audit Action: Manually Integrated


Audit Reason: Baseline

## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190048.d  
Injection Date: 20-May-2023 04:48:36 Instrument ID: CHHPLC\_X3  
Lims ID: 280-176674-A-10-A Lab Sample ID: 280-176674-10  
Client ID: FBQmw-175-230401-GW-R  
Operator ID: JZ/JG ALS Bottle#: 48 Worklist Smp#: 48  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**8 RDX, CAS: 121-82-4, Signal: 1**

RT: 7.69  
Response: 101948  
Amount: 0.958340



Reviewer: LV5D, 20-May-2023 11:14:43

Audit Action: Marked Compound Undetected


Audit Reason: Invalid Compound ID

## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190048.d  
Injection Date: 20-May-2023 04:48:36 Instrument ID: CHHPLC\_X3  
Lims ID: 280-176674-A-10-A Lab Sample ID: 280-176674-10  
Client ID: FBQmw-175-230401-GW-R  
Operator ID: JZ/JG ALS Bottle#: 48 Worklist Smp#: 48  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

3 TNX, CAS: 13980-04-6, Signal: 1

RT: 6.37  
Response: 1691772  
Amount: 8.557817



Reviewer: LV5D, 20-May-2023 11:14:43

Audit Action: Marked Compound Undetected

Audit Reason: Invalid Compound ID

FORM I  
HPLC/IC ORGANICS ANALYSIS DATA SHEET

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Client Sample ID: FBQmw-175-230402-GW-R Lab Sample ID: 280-176674-11

Matrix: Water Lab File ID: 05190053.D

Analysis Method: 8330B Date Collected: 05/17/2023 08:58

Extraction Method: 3535 Date Extracted: 05/19/2023 13:39

Sample wt/vol: 497 (mL) Date Analyzed: 05/20/2023 06:43

Con. Extract Vol.: 5 (mL) Dilution Factor: 1

Injection Volume: 100 (uL) GC Column: UltraCarb5uODS ID: 4.6 (mm)

% Moisture: \_\_\_\_\_ % Solids: \_\_\_\_\_ GPC Cleanup: (Y/N) N

Cleanup Factor: \_\_\_\_\_

Analysis Batch No.: 613168 Units: ug/L

| CAS NO.    | COMPOUND NAME              | RESULT | Q   | LOQ  | LOD   | DL    |
|------------|----------------------------|--------|-----|------|-------|-------|
| 99-35-4    | 1,3,5-Trinitrobenzene      | 0.20   | U   | 0.21 | 0.20  | 0.085 |
| 99-65-0    | 1,3-Dinitrobenzene         | 0.10   | U   | 0.11 | 0.10  | 0.037 |
| 118-96-7   | 2,4,6-Trinitrotoluene      | 0.10   | U   | 0.11 | 0.10  | 0.045 |
| 6629-29-4  | 2,4-diamino-6-nitrotoluene | 0.91   | U   | 1.0  | 0.91  | 0.44  |
| 121-14-2   | 2,4-Dinitrotoluene         | 0.080  | U   | 0.10 | 0.080 | 0.028 |
| 59229-75-3 | 2,6-diamino-4-nitrotoluene | 0.91   | U   | 1.0  | 0.91  | 0.22  |
| 606-20-2   | 2,6-Dinitrotoluene         | 0.080  | U   | 0.10 | 0.080 | 0.040 |
| 35572-78-2 | 2-Amino-4,6-dinitrotoluene | 0.10   | U   | 0.11 | 0.10  | 0.051 |
| 88-72-2    | 2-Nitrotoluene             | 0.20   | U M | 0.21 | 0.20  | 0.086 |
| 618-87-1   | 3,5-Dinitroaniline         | 0.30   | U   | 0.40 | 0.30  | 0.13  |
| 99-08-1    | 3-Nitrotoluene             | 0.35   | U   | 0.40 | 0.35  | 0.20  |
| 19406-51-0 | 4-Amino-2,6-dinitrotoluene | 0.12   | U   | 0.15 | 0.12  | 0.058 |
| 99-99-0    | 4-Nitrotoluene             | 0.40   | U   | 0.41 | 0.40  | 0.10  |
| 80251-29-2 | DNX                        | 0.25   | U   | 0.50 | 0.25  | 0.098 |
| 2691-41-0  | HMX                        | 0.20   | U   | 0.21 | 0.20  | 0.088 |
| 5755-27-1  | MNX                        | 0.29   | U   | 0.50 | 0.29  | 0.093 |
| 98-95-3    | Nitrobenzene               | 0.20   | U   | 0.21 | 0.20  | 0.092 |
| 55-63-0    | Nitroglycerin              | 2.0    | U   | 2.1  | 2.0   | 0.93  |
| 78-11-5    | PETN                       | 1.0    | U   | 1.1  | 1.0   | 0.45  |
| 121-82-4   | RDX                        | 0.20   | U M | 0.21 | 0.20  | 0.052 |
| 479-45-8   | Tetryl                     | 0.10   | U   | 0.11 | 0.10  | 0.032 |
| 13980-04-6 | TNX                        | 0.25   | U M | 0.50 | 0.25  | 0.080 |

| CAS NO.  | SURROGATE          | %REC | Q | LIMITS |
|----------|--------------------|------|---|--------|
| 528-29-0 | 1,2-Dinitrobenzene | 93   | M | 83-119 |

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\05190053.D  
 Lims ID: 280-176674-A-11-A  
 Client ID: FBQmw-175-230402-GW-R  
 Sample Type: Client  
 Inject. Date: 20-May-2023 06:43:26 ALS Bottle#: 53 Worklist Smp#: 53  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-176674-A-11-A  
 Operator ID: JZ/JG Instrument ID: CHHPLC\_X3  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 20-May-2023 11:22:13 Calib Date: 24-Feb-2023 17:41:20  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1663

First Level Reviewer: LV5D

Date: 20-May-2023 11:16:13

| Compound                      | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | OnCol Amt<br>ug/mL | Flags |
|-------------------------------|-----|--------------|------------------|------------------|----------|--------------------|-------|
| 3 TNX                         | 1   |              | 6.430            |                  |          | ND                 | U     |
| 2 2,6-diamino-4-nitrotoluene  | 1   |              | 6.431            |                  |          | ND                 |       |
| 4 HMX                         | 1   |              | 6.550            |                  |          | ND                 |       |
| 5 2,4-diamino-6-nitrotoluene  | 1   |              | 6.618            |                  |          | ND                 |       |
| 6 DNX                         | 1   |              | 6.757            |                  |          | ND                 |       |
| 7 MNX                         | 1   |              | 7.190            |                  |          | ND                 |       |
| 8 RDX                         | 1   |              | 7.570            |                  |          | ND                 | U     |
| \$ 10 1,2-Dinitrobenzene      | 1   | 8.530        | 8.523            | 0.007            | 23513    | 0.1862             | M     |
| 11 1,3,5-Trinitrobenzene      | 1   |              | 8.650            |                  |          | ND                 |       |
| 12 1,3-Dinitrobenzene         | 1   |              | 9.277            |                  |          | ND                 |       |
| 13 Nitrobenzene               | 1   |              | 9.643            |                  |          | ND                 |       |
| 14 3,5-Dinitroaniline         | 1   |              | 9.883            |                  |          | ND                 |       |
| 15 Tetryl                     | 1   |              | 10.017           |                  |          | ND                 |       |
| 16 Nitroglycerin              | 2   |              | 10.470           |                  |          | ND                 |       |
| 17 2,4,6-Trinitrotoluene      | 1   |              | 10.910           |                  |          | ND                 |       |
| 18 4-Amino-2,6-dinitrotoluene | 1   |              | 11.117           |                  |          | ND                 |       |
| 19 2-Amino-4,6-dinitrotoluene | 1   |              | 11.370           |                  |          | ND                 |       |
| 20 2,6-Dinitrotoluene         | 1   |              | 11.517           |                  |          | ND                 |       |
| 21 2,4-Dinitrotoluene         | 1   |              | 11.683           |                  |          | ND                 |       |
| 22 o-Nitrotoluene             | 1   | 12.503       | 12.517           | -0.014           | 881      | 0.006888           | 7M    |
| 23 p-Nitrotoluene             | 1   |              | 12.937           |                  |          | ND                 |       |
| 24 m-Nitrotoluene             | 1   |              | 13.517           |                  |          | ND                 |       |
| 25 PETN                       | 2   |              | 14.697           |                  |          | ND                 |       |

**QC Flag Legend**

Processing Flags

7 - Failed Limit of Detection

Review Flags

M - Manually Integrated

U - Marked Undetected

Report Date: 20-May-2023 11:22:17

Chrom Revision: 2.3 16-May-2023 16:37:30

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190053.d

Injection Date: 20-May-2023 06:43:26

Instrument ID: CHHPLC\_X3

Operator ID: JZ/JG

Lims ID: 280-176674-A-11-A

Lab Sample ID: 280-176674-11

Worklist Smp#: 53

Client ID: FBQmw-175-230402-GW-R

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 53

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



Eurofins Denver  
Recovery Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\05190053.D  
Lims ID: 280-176674-A-11-A  
Client ID: FBQmw-175-230402-GW-R  
Sample Type: Client  
Inject. Date: 20-May-2023 06:43:26 ALS Bottle#: 53 Worklist Smp#: 53  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Sample Info: 280-176674-A-11-A  
Operator ID: JZ/JG Instrument ID: CHHPLC\_X3  
Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\8330\_X3.m  
Limit Group: GCSV - 8330  
Last Update: 20-May-2023 11:22:13 Calib Date: 24-Feb-2023 17:41:20  
Integrator: Falcon  
Quant Method: External Standard Quant By: Initial Calibration  
Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
Process Host: CTX1663

First Level Reviewer: LV5D

Date: 20-May-2023 11:16:13

| Compound                 | Amount Added | Amount Recovered | % Rec. |
|--------------------------|--------------|------------------|--------|
| \$ 10 1,2-Dinitrobenzene | 0.2000       | 0.1862           | 93.08  |



## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190053.d  
Injection Date: 20-May-2023 06:43:26 Instrument ID: CHHPLC\_X3  
Lims ID: 280-176674-A-11-A Lab Sample ID: 280-176674-11  
Client ID: FBQmw-175-230402-GW-R  
Operator ID: JZ/JG ALS Bottle#: 53 Worklist Smp#: 53  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector: LC DAD1B, 254 nm

## \$ 10 1,2-Dinitrobenzene, CAS: 528-29-0

Signal: 1


RT: 8.53  
Area: 28365  
Amount: 0.224569  
Amount Units: ug/mL

## Processing Integration Results



RT: 8.53  
Area: 23513  
Amount: 0.186155  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 11:16:08 -06:00:00 (UTC)

Audit Action: Manually Integrated

Audit Reason: Baseline

## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190053.d  
Injection Date: 20-May-2023 06:43:26 Instrument ID: CHHPLC\_X3  
Lims ID: 280-176674-A-11-A Lab Sample ID: 280-176674-11  
Client ID: FBQmw-175-230402-GW-R  
Operator ID: JZ/JG ALS Bottle#: 53 Worklist Smp#: 53  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**22 o-Nitrotoluene, CAS: 88-72-2**

Signal: 1


RT: 12.51  
Area: 1718  
Amount: 0.013433  
Amount Units: ug/mL

## Processing Integration Results



RT: 12.50  
Area: 881  
Amount: 0.006888  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 11:16:12 -06:00:00 (UTC)

Audit Action: Manually Integrated


Audit Reason: Baseline

## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190053.d  
Injection Date: 20-May-2023 06:43:26 Instrument ID: CHHPLC\_X3  
Lims ID: 280-176674-A-11-A Lab Sample ID: 280-176674-11  
Client ID: FBQmw-175-230402-GW-R  
Operator ID: JZ/JG ALS Bottle#: 53 Worklist Smp#: 53  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**8 RDX, CAS: 121-82-4, Signal: 1**

RT: 7.70  
Response: 127518  
Amount: 1.198705



Reviewer: LV5D, 20-May-2023 11:16:13

Audit Action: Marked Compound Undetected


Audit Reason: Invalid Compound ID

## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190053.d  
Injection Date: 20-May-2023 06:43:26 Instrument ID: CHHPLC\_X3  
Lims ID: 280-176674-A-11-A Lab Sample ID: 280-176674-11  
Client ID: FBQmw-175-230402-GW-R  
Operator ID: JZ/JG ALS Bottle#: 53 Worklist Smp#: 53  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

3 TNX, CAS: 13980-04-6, Signal: 1

RT: 6.40  
Response: 1704684  
Amount: 8.623132



Reviewer: LV5D, 20-May-2023 11:16:13

Audit Action: Marked Compound Undetected

Audit Reason: Invalid Compound ID

FORM VI  
HPLC/IC BY EXTERNAL STANDARD - INITIAL CALIBRATION DATA  
RETENTION TIME SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1 Analy Batch No.: 599427  
SDG No.: \_\_\_\_\_  
Instrument ID: CHHPLC\_X3 GC Column: UltraCarb5u ID: 4.6 (mm) Heated Purge: (Y/N) N  
Calibration Start Date: 01/17/2023 00:33 Calibration End Date: 01/17/2023 03:16 Calibration ID: 75939

Calibration Files

|         |                  |              |
|---------|------------------|--------------|
| LEVEL:  | LAB SAMPLE ID:   | LAB FILE ID: |
| Level 1 | IC 280-599427/18 | 01160046.D   |
| Level 2 | IC 280-599427/17 | 01160045.D   |
| Level 3 | IC 280-599427/16 | 01160044.D   |
| Level 4 | IC 280-599427/15 | 01160043.D   |
| Level 5 | IC 280-599427/14 | 01160042.D   |
| Level 6 | IC 280-599427/13 | 01160041.D   |
| Level 7 | IC 280-599427/12 | 01160040.D   |
| Level 8 | IC 280-599427/11 | 01160039.D   |

| ANALYTE                    | LVL 1 | LVL 2 | LVL 3 | LVL 4 | LVL 5 | LVL 6 | LVL 7 | LVL 8 |  |  | RT WINDOW     | AVG RT |
|----------------------------|-------|-------|-------|-------|-------|-------|-------|-------|--|--|---------------|--------|
| 2,6-diamino-4-nitrotoluene | 6.459 | 6.459 | 6.459 | 6.460 | 6.459 | 6.460 | 6.457 | 6.432 |  |  | 6.305 - 6.605 | 6.456  |
| 2,4-diamino-6-nitrotoluene | 6.639 | 6.639 | 6.632 | 6.640 | 6.639 | 6.640 | 6.637 | 6.618 |  |  | 6.619 - 6.919 | 6.636  |
| 3,5-Dinitroaniline         | 9.805 | 9.798 | 9.812 | 9.800 | 9.806 | 9.800 | 9.804 | 9.798 |  |  | 9.650 - 9.950 | 9.803  |

FORM VI  
HPLC/IC BY EXTERNAL STANDARD - INITIAL CALIBRATION DATA  
CURVE EVALUATION

Lab Name: Eurofins Denver Job No.: 280-176674-1 Analy Batch No.: 599427  
SDG No.: \_\_\_\_\_  
Instrument ID: CHHPLC\_X3 GC Column: UltraCarb5u ID: 4.6(mm) Heated Purge: (Y/N) N  
Calibration Start Date: 01/17/2023 00:33 Calibration End Date: 01/17/2023 03:16 Calibration ID: 75939

Calibration Files

|         |                  |              |
|---------|------------------|--------------|
| LEVEL:  | LAB SAMPLE ID:   | LAB FILE ID: |
| Level 1 | IC 280-599427/18 | 01160046.D   |
| Level 2 | IC 280-599427/17 | 01160045.D   |
| Level 3 | IC 280-599427/16 | 01160044.D   |
| Level 4 | IC 280-599427/15 | 01160043.D   |
| Level 5 | IC 280-599427/14 | 01160042.D   |
| Level 6 | IC 280-599427/13 | 01160041.D   |
| Level 7 | IC 280-599427/12 | 01160040.D   |
| Level 8 | IC 280-599427/11 | 01160039.D   |

| ANALYTE                    | CF               |                  |                  |                  | CURVE<br>TYPE | COEFFICIENT    |                |    | # | MIN CF | %RSD<br>/RSE | # | MAX<br>%RSD<br>/RSE | R^2<br>OR COD | # | MIN R^2<br>OR COD |
|----------------------------|------------------|------------------|------------------|------------------|---------------|----------------|----------------|----|---|--------|--------------|---|---------------------|---------------|---|-------------------|
|                            | LVL 1<br>LVL 5   | LVL 2<br>LVL 6   | LVL 3<br>LVL 7   | LVL 4<br>LVL 8   |               | B              | M1             | M2 |   |        |              |   |                     |               |   |                   |
| 2,6-diamino-4-nitrotoluene | 266300<br>235998 | 260780<br>221767 | 247020<br>234891 | 225236<br>222744 | Lin2          | 864.13225<br>4 | 229692.99<br>3 |    |   |        |              |   |                     | 0.9980        |   | 0.9900            |
| 2,4-diamino-6-nitrotoluene | 152100<br>147503 | 163500<br>136319 | 149080<br>148050 | 138680<br>136546 | Ave           |                | 146472.13<br>4 |    |   |        | 6.3          |   | 20.0                |               |   |                   |
| 3,5-Dinitroaniline         | 219850<br>241580 | 252000<br>241484 | 246040<br>245164 | 237836<br>244381 | Lin2          | -392.4275<br>7 | 245423.81<br>0 |    |   |        |              |   |                     | 0.9990        |   | 0.9900            |

Note: The M1 coefficient is the same as Ave CF for an Ave curve type. RSD is calculated for Ave curve types. RSE is used for all other types.

FORM VI  
HPLC/IC BY EXTERNAL STANDARD - INITIAL CALIBRATION DATA  
RESPONSE AND CONCENTRATION

Lab Name: Eurofins Denver Job No.: 280-176674-1 Analy Batch No.: 599427

SDG No.: \_\_\_\_\_

Instrument ID: CHHPLC\_X3 GC Column: UltraCarb5u ID: 4.6(mm) Heated Purge: (Y/N) N

Calibration Start Date: 01/17/2023 00:33 Calibration End Date: 01/17/2023 03:16 Calibration ID: 75939

Calibration Files

|         |                  |              |
|---------|------------------|--------------|
| LEVEL:  | LAB SAMPLE ID:   | LAB FILE ID: |
| Level 1 | IC 280-599427/18 | 01160046.D   |
| Level 2 | IC 280-599427/17 | 01160045.D   |
| Level 3 | IC 280-599427/16 | 01160044.D   |
| Level 4 | IC 280-599427/15 | 01160043.D   |
| Level 5 | IC 280-599427/14 | 01160042.D   |
| Level 6 | IC 280-599427/13 | 01160041.D   |
| Level 7 | IC 280-599427/12 | 01160040.D   |
| Level 8 | IC 280-599427/11 | 01160039.D   |

| ANALYTE                    | CURVE TYPE | RESPONSE       |                |                |       |       | CONCENTRATION (UG/ML) |                |                |       |       |
|----------------------------|------------|----------------|----------------|----------------|-------|-------|-----------------------|----------------|----------------|-------|-------|
|                            |            | LVL 1<br>LVL 6 | LVL 2<br>LVL 7 | LVL 3<br>LVL 8 | LVL 4 | LVL 5 | LVL 1<br>LVL 6        | LVL 2<br>LVL 7 | LVL 3<br>LVL 8 | LVL 4 | LVL 5 |
| 2,6-diamino-4-nitrotoluene | Lin2       | 5326           | 13039          | 24702          | 56309 | 94399 | 0.0200                | 0.0500         | 0.100          | 0.250 | 0.400 |
|                            |            | 155237         | 234891         | 556860         |       |       | 0.700                 | 1.00           | 2.50           |       |       |
| 2,4-diamino-6-nitrotoluene | Ave        | 3042           | 8175           | 14908          | 34670 | 59001 | 0.0200                | 0.0500         | 0.100          | 0.250 | 0.400 |
|                            |            | 95423          | 148050         | 341365         |       |       | 0.700                 | 1.00           | 2.50           |       |       |
| 3,5-Dinitroaniline         | Lin2       | 4397           | 12600          | 24604          | 59459 | 96632 | 0.0200                | 0.0500         | 0.100          | 0.250 | 0.400 |
|                            |            | 169039         | 245164         | 610953         |       |       | 0.700                 | 1.00           | 2.50           |       |       |

Curve Type Legend

Ave = Average  
Lin2 = Linear 1/conc^2

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230114-117775.b\01160039.D  
 Lims ID: IC ADD 8  
 Client ID:  
 Sample Type: IC Calib Level: 8  
 Inject. Date: 17-Jan-2023 00:33:01 ALS Bottle#: 42 Worklist Smp#: 11  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: ICADD8  
 Operator ID: JZ/MAR Instrument ID: CHHPLC\_X3  
 Sublist: chrom-8330\_X3\*sub10  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230114-117775.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 17-Jan-2023 12:50:38 Calib Date: 17-Jan-2023 06:47:07  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230114-117775.b\01160055.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1672

First Level Reviewer: NQH7

Date: 17-Jan-2023 08:20:12

| Compound                     | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/mL | OnCol Amt<br>ug/mL | Flags |
|------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 2 2,6-diamino-4-nitrotoluene | 1   | 6.432        | 6.455            | -0.023           | 556860   | 2.50             | 2.42               | M     |
| 5 2,4-diamino-6-nitrotoluene | 1   | 6.618        | 6.769            | -0.151           | 341365   | 2.50             | 2.33               | M     |
| 14 3,5-Dinitroaniline        | 1   | 9.798        | 9.800            | -0.002           | 610953   | 2.50             | 2.49               |       |

**QC Flag Legend**

Processing Flags

Review Flags

M - Manually Integrated

**Reagents:**

8330\_ADDs\_00035

Amount Added: 125.00

Units: uL



Report Date: 17-Jan-2023 12:50:38

Chrom Revision: 2.3 20-Dec-2022 14:14:06

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230114-117775.b\01160039.d

Injection Date: 17-Jan-2023 00:33:01

Instrument ID: CHHPLC\_X3

Operator ID: JZ/MAR

Lims ID: IC ADD 8

Worklist Smp#: 11

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 42

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230114-117775.b\01160039.d  
Injection Date: 17-Jan-2023 00:33:01 Instrument ID: CHHPLC\_X3  
Lims ID: IC ADD 8  
Client ID:  
Operator ID: JZ/MAR ALS Bottle#: 42 Worklist Smp#: 11  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**2,2,6-diamino-4-nitrotoluene, CAS: 59229-75-3**

Signal: 1


RT: 6.43  
Area: 581780  
Amount: 2.370843  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.43  
Area: 556860  
Amount: 2.420604  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 17-Jan-2023 12:22:46  
Audit Action: Split an Integrated Peak

Audit Reason: Baseline

## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230114-117775.b\01160039.d

Injection Date: 17-Jan-2023 00:33:01

Instrument ID: CHHPLC\_X3

Lims ID: IC ADD 8

Client ID:

Operator ID: JZ/MAR

ALS Bottle#:

42

Worklist Smp#: 11

Injection Vol: 100.0 ul

Dil. Factor:

1.0000

Method: 8330\_X3

Limit Group:

GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm)

Detector


LC DAD1B, 254 nm

**5 2,4-diamino-6-nitrotoluene, CAS: 6629-29-4**

Signal: 1


RT: 6.62  
Area: 364706  
Amount: 2.146312  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.62  
Area: 341365  
Amount: 2.330580  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 17-Jan-2023 12:22:47

Audit Action: Split an Integrated Peak

Audit Reason: Baseline

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230114-117775.b\01160040.D  
 Lims ID: IC ADD 7  
 Client ID:  
 Sample Type: IC Calib Level: 7  
 Inject. Date: 17-Jan-2023 00:56:21 ALS Bottle#: 43 Worklist Smp#: 12  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: ICADD7  
 Operator ID: JZ/MAR Instrument ID: CHHPLC\_X3  
 Sublist: chrom-8330\_X3\*sub10  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230114-117775.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 17-Jan-2023 12:50:39 Calib Date: 17-Jan-2023 06:47:07  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230114-117775.b\01160055.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1672

First Level Reviewer: LV5D

Date: 17-Jan-2023 12:23:10

| Compound                     | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/mL | OnCol Amt<br>ug/mL | Flags |
|------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 2 2,6-diamino-4-nitrotoluene | 1   | 6.457        | 6.455            | 0.002            | 234891   | 1.00             | 1.02               | M     |
| 5 2,4-diamino-6-nitrotoluene | 1   | 6.637        | 6.769            | -0.132           | 148050   | 1.00             | 1.01               | M     |
| 14 3,5-Dinitroaniline        | 1   | 9.804        | 9.800            | 0.004            | 245164   | 1.00             | 1.00               |       |

**QC Flag Legend**

Processing Flags

Review Flags

M - Manually Integrated

**Reagents:**

8330\_ADDs\_00035

Amount Added: 50.00

Units: uL

Report Date: 17-Jan-2023 12:50:39

Chrom Revision: 2.3 20-Dec-2022 14:14:06

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230114-117775.b\01160040.d

Injection Date: 17-Jan-2023 00:56:21

Instrument ID: CHHPLC\_X3

Operator ID: JZ/MAR

Lims ID: IC ADD 7

Worklist Smp#: 12

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 43

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230114-117775.b\01160040.d

Injection Date: 17-Jan-2023 00:56:21

Instrument ID: CHHPLC\_X3

Lims ID: IC ADD 7

Client ID:

Operator ID: JZ/MAR

ALS Bottle#:

43

Worklist Smp#: 12

Injection Vol: 100.0 ul

Dil. Factor:

1.0000

Method: 8330\_X3

Limit Group:

GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm)

Detector


LC DAD1B, 254 nm

**2,2,6-diamino-4-nitrotoluene, CAS: 59229-75-3**

Signal: 1


RT: 6.46  
Area: 243731  
Amount: 0.963471  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.46  
Area: 234891  
Amount: 1.018868  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 17-Jan-2023 12:23:08

Audit Action: Split an Integrated Peak

Audit Reason: Baseline

## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230114-117775.b\01160040.d

Injection Date: 17-Jan-2023 00:56:21

Instrument ID: CHHPLC\_X3

Lims ID: IC ADD 7

Client ID:

Operator ID: JZ/MAR

ALS Bottle#:

43

Worklist Smp#: 12

Injection Vol: 100.0 ul

Dil. Factor:

1.0000

Method: 8330\_X3

Limit Group:

GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm)

Detector


LC DAD1B, 254 nm

**5 2,4-diamino-6-nitrotoluene, CAS: 6629-29-4**

Signal: 1


RT: 6.64  
Area: 156243  
Amount: 0.928768  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.64  
Area: 148050  
Amount: 1.010772  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 17-Jan-2023 12:23:09

Audit Action: Split an Integrated Peak

Audit Reason: Baseline

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230114-117775.b\01160041.D  
 Lims ID: IC ADD 6  
 Client ID:  
 Sample Type: IC Calib Level: 6  
 Inject. Date: 17-Jan-2023 01:19:43 ALS Bottle#: 44 Worklist Smp#: 13  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: ICADD6  
 Operator ID: JZ/MAR Instrument ID: CHHPLC\_X3  
 Sublist: chrom-8330\_X3\*sub10  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230114-117775.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 17-Jan-2023 12:50:39 Calib Date: 17-Jan-2023 06:47:07  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230114-117775.b\01160055.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1672

First Level Reviewer: LV5D

Date: 17-Jan-2023 12:23:22

| Compound                     | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/mL | OnCol Amt<br>ug/mL | Flags |
|------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 2 2,6-diamino-4-nitrotoluene | 1   | 6.460        | 6.455            | 0.005            | 155237   | 0.7000           | 0.6721             | M     |
| 5 2,4-diamino-6-nitrotoluene | 1   | 6.640        | 6.769            | -0.129           | 95423    | 0.7000           | 0.6515             | M     |
| 14 3,5-Dinitroaniline        | 1   | 9.800        | 9.800            | 0.000            | 169039   | 0.7000           | 0.6904             |       |

**QC Flag Legend**

Review Flags

M - Manually Integrated

**Reagents:**

8330\_ADDs\_00035

Amount Added: 35.00

Units: uL



Report Date: 17-Jan-2023 12:50:39

Chrom Revision: 2.3 20-Dec-2022 14:14:06

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230114-117775.b\01160041.d

Injection Date: 17-Jan-2023 01:19:43

Instrument ID: CHHPLC\_X3

Operator ID: JZ/MAR

Lims ID: IC ADD 6

Worklist Smp#: 13

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 44

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230114-117775.b\01160041.d  
Injection Date: 17-Jan-2023 01:19:43 Instrument ID: CHHPLC\_X3  
Lims ID: IC ADD 6  
Client ID:  
Operator ID: JZ/MAR ALS Bottle#: 44 Worklist Smp#: 13  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**2,2,6-diamino-4-nitrotoluene, CAS: 59229-75-3**

Signal: 1


RT: 6.46  
Area: 174539  
Amount: 0.677131  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.46  
Area: 155237  
Amount: 0.672083  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 17-Jan-2023 12:23:20  
Audit Action: Split an Integrated Peak

Audit Reason: Baseline

## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230114-117775.b\01160041.d  
Injection Date: 17-Jan-2023 01:19:43 Instrument ID: CHHPLC\_X3  
Lims ID: IC ADD 6  
Client ID:  
Operator ID: JZ/MAR ALS Bottle#: 44 Worklist Smp#: 13  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**5 2,4-diamino-6-nitrotoluene, CAS: 6629-29-4**

Signal: 1


RT: 6.64  
Area: 114251  
Amount: 0.682909  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.64  
Area: 95423  
Amount: 0.651475  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 17-Jan-2023 12:23:21  
Audit Action: Split an Integrated Peak

Audit Reason: Baseline

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230114-117775.b\01160042.D  
 Lims ID: IC ADD 5  
 Client ID:  
 Sample Type: IC Calib Level: 5  
 Inject. Date: 17-Jan-2023 01:43:11 ALS Bottle#: 45 Worklist Smp#: 14  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: ICADD5  
 Operator ID: JZ/MAR Instrument ID: CHHPLC\_X3  
 Sublist: chrom-8330\_X3\*sub10  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230114-117775.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 17-Jan-2023 12:50:40 Calib Date: 17-Jan-2023 06:47:07  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230114-117775.b\01160055.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1672

First Level Reviewer: LV5D

Date: 17-Jan-2023 12:30:57

| Compound                     | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/mL | OnCol Amt<br>ug/mL | Flags |
|------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 2 2,6-diamino-4-nitrotoluene | 1   | 6.459        | 6.455            | 0.004            | 94399    | 0.4000           | 0.4072             | M     |
| 5 2,4-diamino-6-nitrotoluene | 1   | 6.639        | 6.769            | -0.130           | 59001    | 0.4000           | 0.4028             | M     |
| 14 3,5-Dinitroaniline        | 1   | 9.806        | 9.800            | 0.006            | 96632    | 0.4000           | 0.3953             |       |

**QC Flag Legend**

Review Flags

M - Manually Integrated

**Reagents:**

8330\_ADDs\_00035

Amount Added: 20.00

Units: uL

Report Date: 17-Jan-2023 12:50:40

Chrom Revision: 2.3 20-Dec-2022 14:14:06

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230114-117775.b\01160042.d

Injection Date: 17-Jan-2023 01:43:11

Instrument ID: CHHPLC\_X3

Operator ID: JZ/MAR

Lims ID: IC ADD 5

Worklist Smp#: 14

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 45

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230114-117775.b\01160042.d

Injection Date: 17-Jan-2023 01:43:11

Instrument ID: CHHPLC\_X3

Lims ID: IC ADD 5

Client ID:

Operator ID: JZ/MAR

ALS Bottle#:

45

Worklist Smp#: 14

Injection Vol: 100.0 ul

Dil. Factor:

1.0000

Method: 8330\_X3

Limit Group:

GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm)

Detector


LC DAD1B, 254 nm

**2,2,6-diamino-4-nitrotoluene, CAS: 59229-75-3**

Signal: 1


RT: 6.46  
Area: 108476  
Amount: 0.410648  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.46  
Area: 94399  
Amount: 0.407217  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 17-Jan-2023 12:30:54

Audit Action: Split an Integrated Peak

Audit Reason: Baseline

## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230114-117775.b\01160042.d

Injection Date: 17-Jan-2023 01:43:11

Instrument ID: CHHPLC\_X3

Lims ID: IC ADD 5

Client ID:

Operator ID: JZ/MAR

ALS Bottle#:

45

Worklist Smp#: 14

Injection Vol: 100.0 ul

Dil. Factor:

1.0000

Method: 8330\_X3

Limit Group:

GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm)

Detector


LC DAD1B, 254 nm

**5 2,4-diamino-6-nitrotoluene, CAS: 6629-29-4**

Signal: 1


RT: 6.64  
Area: 71991  
Amount: 0.437489  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.64  
Area: 59001  
Amount: 0.402814  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 17-Jan-2023 12:30:55

Audit Action: Split an Integrated Peak

Audit Reason: Baseline

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230114-117775.b\01160043.D  
 Lims ID: IC ADD 4  
 Client ID:  
 Sample Type: IC Calib Level: 4  
 Inject. Date: 17-Jan-2023 02:06:34 ALS Bottle#: 46 Worklist Smp#: 15  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: ICADD4  
 Operator ID: JZ/MAR Instrument ID: CHHPLC\_X3  
 Sublist: chrom-8330\_X3\*sub10  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230114-117775.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 17-Jan-2023 12:50:41 Calib Date: 17-Jan-2023 06:47:07  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230114-117775.b\01160055.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1672

First Level Reviewer: LV5D

Date: 17-Jan-2023 12:31:13

| Compound                     | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/mL | OnCol Amt<br>ug/mL | Flags |
|------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 2 2,6-diamino-4-nitrotoluene | 1   | 6.460        | 6.455            | 0.005            | 56309    | 0.2500           | 0.2414             | M     |
| 5 2,4-diamino-6-nitrotoluene | 1   | 6.640        | 6.769            | -0.129           | 34670    | 0.2500           | 0.2367             | M     |
| 14 3,5-Dinitroaniline        | 1   | 9.800        | 9.800            | 0.000            | 59459    | 0.2500           | 0.2439             |       |

**QC Flag Legend**

Review Flags

M - Manually Integrated

**Reagents:**

8330\_ADDs\_00035

Amount Added: 12.50

Units: uL



Report Date: 17-Jan-2023 12:50:41

Chrom Revision: 2.3 20-Dec-2022 14:14:06

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230114-117775.b\01160043.d

Injection Date: 17-Jan-2023 02:06:34

Instrument ID: CHHPLC\_X3

Operator ID: JZ/MAR

Lims ID: IC ADD 4

Worklist Smp#: 15

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 46

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230114-117775.b\01160043.d

Injection Date: 17-Jan-2023 02:06:34

Instrument ID: CHHPLC\_X3

Lims ID: IC ADD 4

Client ID:

Operator ID: JZ/MAR

ALS Bottle#:

46

Worklist Smp#: 15

Injection Vol: 100.0 ul

Dil. Factor:

1.0000

Method: 8330\_X3

Limit Group:

GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm)

Detector


LC DAD1B, 254 nm

**2,2,6-diamino-4-nitrotoluene, CAS: 59229-75-3**

Signal: 1


RT: 6.46  
Area: 74424  
Amount: 0.275328  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.46  
Area: 56309  
Amount: 0.241387  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 17-Jan-2023 12:31:11

Audit Action: Assigned New Baseline

Audit Reason: Baseline

## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230114-117775.b\01160043.d

Injection Date: 17-Jan-2023 02:06:34

Instrument ID: CHHPLC\_X3

Lims ID: IC ADD 4

Client ID:

Operator ID: JZ/MAR

ALS Bottle#:

46

Worklist Smp#: 15

Injection Vol: 100.0 ul

Dil. Factor:

1.0000

Method: 8330\_X3

Limit Group:

GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm)

Detector


LC DAD1B, 254 nm

**5 2,4-diamino-6-nitrotoluene, CAS: 6629-29-4**

Signal: 1


RT: 6.64  
Area: 50478  
Amount: 0.310727  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.64  
Area: 34670  
Amount: 0.236700  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 17-Jan-2023 12:31:11

Audit Action: Assigned New Baseline

Audit Reason: Baseline

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230114-117775.b\01160044.D  
 Lims ID: IC ADD 3  
 Client ID:  
 Sample Type: IC Calib Level: 3  
 Inject. Date: 17-Jan-2023 02:30:00 ALS Bottle#: 47 Worklist Smp#: 16  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: ICADD3  
 Operator ID: JZ/MAR Instrument ID: CHHPLC\_X3  
 Sublist: chrom-8330\_X3\*sub10  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230114-117775.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 17-Jan-2023 12:50:41 Calib Date: 17-Jan-2023 06:47:07  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230114-117775.b\01160055.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1672

First Level Reviewer: LV5D

Date: 17-Jan-2023 12:31:24

| Compound                     | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/mL | OnCol Amt<br>ug/mL | Flags |
|------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 2 2,6-diamino-4-nitrotoluene | 1   | 6.459        | 6.455            | 0.004            | 24702    | 0.1000           | 0.1038             | M     |
| 5 2,4-diamino-6-nitrotoluene | 1   | 6.632        | 6.769            | -0.137           | 14908    | 0.1000           | 0.1018             | M     |
| 14 3,5-Dinitroaniline        | 1   | 9.812        | 9.800            | 0.012            | 24604    | 0.1000           | 0.1019             |       |

**QC Flag Legend**

Review Flags

M - Manually Integrated

**Reagents:**

8330\_ADDs\_00035

Amount Added: 5.00

Units: uL

Report Date: 17-Jan-2023 12:50:41

Chrom Revision: 2.3 20-Dec-2022 14:14:06

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230114-117775.b\01160044.d

Injection Date: 17-Jan-2023 02:30:00

Instrument ID: CHHPLC\_X3

Operator ID: JZ/MAR

Lims ID: IC ADD 3

Worklist Smp#: 16

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 47

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230114-117775.b\01160044.d

Injection Date: 17-Jan-2023 02:30:00

Instrument ID: CHHPLC\_X3

Lims ID: IC ADD 3

Client ID:

Operator ID: JZ/MAR

ALS Bottle#:

47

Worklist Smp#: 16

Injection Vol: 100.0 ul

Dil. Factor:

1.0000

Method: 8330\_X3

Limit Group:

GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm)

Detector


LC DAD1B, 254 nm

**2,2,6-diamino-4-nitrotoluene, CAS: 59229-75-3**

Signal: 1


RT: 6.46  
Area: 43811  
Amount: 0.142153  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.46  
Area: 24702  
Amount: 0.103781  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 17-Jan-2023 12:31:22

Audit Action: Split an Integrated Peak

Audit Reason: Baseline

## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230114-117775.b\01160044.d

Injection Date: 17-Jan-2023 02:30:00

Instrument ID: CHHPLC\_X3

Lims ID: IC ADD 3

Client ID:

Operator ID: JZ/MAR

ALS Bottle#:

47

Worklist Smp#: 16

Injection Vol: 100.0 ul

Dil. Factor:

1.0000

Method: 8330\_X3

Limit Group:

GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm)

Detector


LC DAD1B, 254 nm

**5 2,4-diamino-6-nitrotoluene, CAS: 6629-29-4**

Signal: 1


RT: 6.63  
Area: 28852  
Amount: 0.169920  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.63  
Area: 14908  
Amount: 0.101780  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 17-Jan-2023 12:31:22

Audit Action: Split an Integrated Peak

Audit Reason: Baseline

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230114-117775.b\01160045.D  
 Lims ID: IC ADD 2  
 Client ID:  
 Sample Type: IC Calib Level: 2  
 Inject. Date: 17-Jan-2023 02:53:21 ALS Bottle#: 48 Worklist Smp#: 17  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: ICADD2  
 Operator ID: JZ/MAR Instrument ID: CHHPLC\_X3  
 Sublist: chrom-8330\_X3\*sub10  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230114-117775.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 17-Jan-2023 12:50:42 Calib Date: 17-Jan-2023 06:47:07  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230114-117775.b\01160055.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1672

First Level Reviewer: LV5D

Date: 17-Jan-2023 12:31:36

| Compound                     | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/mL | OnCol Amt<br>ug/mL | Flags |
|------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 2 2,6-diamino-4-nitrotoluene | 1   | 6.459        | 6.455            | 0.004            | 13039    | 0.0500           | 0.0530             | M     |
| 5 2,4-diamino-6-nitrotoluene | 1   | 6.639        | 6.769            | -0.130           | 8175     | 0.0500           | 0.0558             | M     |
| 14 3,5-Dinitroaniline        | 1   | 9.798        | 9.800            | -0.002           | 12600    | 0.0500           | 0.0529             |       |

**QC Flag Legend**

Review Flags

M - Manually Integrated

**Reagents:**

8330\_ADDs\_00035

Amount Added: 2.50

Units: uL



Report Date: 17-Jan-2023 12:50:42

Chrom Revision: 2.3 20-Dec-2022 14:14:06

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230114-117775.b\01160045.d

Injection Date: 17-Jan-2023 02:53:21

Instrument ID: CHHPLC\_X3

Operator ID: JZ/MAR

Lims ID: IC ADD 2

Worklist Smp#: 17

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 48

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230114-117775.b\01160045.d

Injection Date: 17-Jan-2023 02:53:21

Instrument ID: CHHPLC\_X3

Lims ID: IC ADD 2

Client ID:

Operator ID: JZ/MAR

ALS Bottle#:

48

Worklist Smp#: 17

Injection Vol: 100.0 ul

Dil. Factor:

1.0000

Method: 8330\_X3

Limit Group:

GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm)

Detector


LC DAD1B, 254 nm

**2,2,6-diamino-4-nitrotoluene, CAS: 59229-75-3**

Signal: 1


RT: 6.46  
Area: 31801  
Amount: 0.089109  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.46  
Area: 13039  
Amount: 0.053005  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 17-Jan-2023 12:31:33

Audit Action: Split an Integrated Peak

Audit Reason: Baseline

## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230114-117775.b\01160045.d

Injection Date: 17-Jan-2023 02:53:21

Instrument ID: CHHPLC\_X3

Lims ID: IC ADD 2

Client ID:

Operator ID: JZ/MAR

ALS Bottle#:

48

Worklist Smp#: 17

Injection Vol: 100.0 ul

Dil. Factor:

1.0000

Method: 8330\_X3

Limit Group:

GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm)

Detector


LC DAD1B, 254 nm

**5 2,4-diamino-6-nitrotoluene, CAS: 6629-29-4**

Signal: 1


RT: 6.64  
Area: 23838  
Amount: 0.125460  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.64  
Area: 8175  
Amount: 0.055813  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 17-Jan-2023 12:31:34

Audit Action: Split an Integrated Peak

Audit Reason: Baseline

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230114-117775.b\01160046.D  
 Lims ID: IC ADD 1  
 Client ID:  
 Sample Type: IC Calib Level: 1  
 Inject. Date: 17-Jan-2023 03:16:44 ALS Bottle#: 49 Worklist Smp#: 18  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: ICADD1  
 Operator ID: JZ/MAR Instrument ID: CHHPLC\_X3  
 Sublist: chrom-8330\_X3\*sub10  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230114-117775.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 17-Jan-2023 12:50:42 Calib Date: 17-Jan-2023 06:47:07  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230114-117775.b\01160055.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1672

First Level Reviewer: LV5D

Date: 17-Jan-2023 12:31:50

| Compound                     | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/mL | OnCol Amt<br>ug/mL | Flags |
|------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 2 2,6-diamino-4-nitrotoluene | 1   | 6.459        | 6.455            | 0.004            | 5326     | 0.0200           | 0.0194             | M     |
| 5 2,4-diamino-6-nitrotoluene | 1   | 6.639        | 6.769            | -0.130           | 3042     | 0.0200           | 0.0208             | M     |
| 14 3,5-Dinitroaniline        | 1   | 9.805        | 9.800            | 0.005            | 4397     | 0.0200           | 0.0195             |       |

**QC Flag Legend**

Processing Flags

Review Flags

M - Manually Integrated

**Reagents:**

8330\_ADDs\_00035

Amount Added: 1.00

Units: uL

Report Date: 17-Jan-2023 12:50:42

Chrom Revision: 2.3 20-Dec-2022 14:14:06

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230114-117775.b\01160046.d

Injection Date: 17-Jan-2023 03:16:44

Instrument ID: CHHPLC\_X3

Operator ID: JZ/MAR

Lims ID: IC ADD 1

Worklist Smp#: 18

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 49

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230114-117775.b\01160046.d  
Injection Date: 17-Jan-2023 03:16:44 Instrument ID: CHHPLC\_X3  
Lims ID: IC ADD 1  
Client ID:  
Operator ID: JZ/MAR ALS Bottle#: 49 Worklist Smp#: 18  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**2,2,6-diamino-4-nitrotoluene, CAS: 59229-75-3**

Signal: 1


RT: 6.46  
Area: 24501  
Amount: 0.068213  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.46  
Area: 5326  
Amount: 0.019425  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 17-Jan-2023 12:31:46  
Audit Action: Split an Integrated Peak

Audit Reason: Baseline

## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230114-117775.b\01160046.d

Injection Date: 17-Jan-2023 03:16:44

Instrument ID: CHHPLC\_X3

Lims ID: IC ADD 1

Client ID:

Operator ID: JZ/MAR

ALS Bottle#:

49

Worklist Smp#: 18

Injection Vol: 100.0 ul

Dil. Factor:

1.0000

Method: 8330\_X3

Limit Group:

GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm)

Detector


LC DAD1B, 254 nm

**5 2,4-diamino-6-nitrotoluene, CAS: 6629-29-4**

Signal: 1


RT: 6.64  
Area: 15283  
Amount: 0.068540  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.64  
Area: 3042  
Amount: 0.020768  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 17-Jan-2023 12:31:48

Audit Action: Split an Integrated Peak

Audit Reason: Baseline

# Calibration

/ 2,6-diamino-4-nitrotoluene

Curve Type: Linear  
 Weighting: Conc\_Sq  
 Origin: None  
 Dependency: Response  
 Calib Mode: ESTD  
 Response Base: AREA  
 RF Rounding: 0


## Curve Coefficients

Intercept: 864.1  
 Slope: 2.297E+05

## Error Coefficients

Standard Error: 8170  
 Relative Standard Error: 4.2  
 Correlation Coefficient: 0.999  
 Coefficient of Determination (Adjusted): 0.998

| ID | Level            | Concentration | Response | IS Amount | IS Response | RF            | Used |
|----|------------------|---------------|----------|-----------|-------------|---------------|------|
| 1  | IC 280-599427/18 | 0.02          | 5326.0   |           |             | 266300.0      | Y    |
| 2  | IC 280-599427/17 | 0.05          | 13039.0  |           |             | 260780.0      | Y    |
| 3  | IC 280-599427/16 | 0.1           | 24702.0  |           |             | 247020.0      | Y    |
| 4  | IC 280-599427/15 | 0.25          | 56309.0  |           |             | 225236.0      | Y    |
| 5  | IC 280-599427/14 | 0.4           | 94399.0  |           |             | 235997.5      | Y    |
| 6  | IC 280-599427/13 | 0.7           | 155237.0 |           |             | 221767.142857 | Y    |
| 7  | IC 280-599427/12 | 1.0           | 234891.0 |           |             | 234891.0      | Y    |
| 8  | IC 280-599427/11 | 2.5           | 556860.0 |           |             | 222744.0      | Y    |





## Calibration

/ 2,4-diamino-6-nitrotoluene

Curve Type: Average  
Weighting: Conc\_Sq  
Origin: Force  
Dependency: Response  
Calib Mode: ESTD  
Response Base: AREA  
RF Rounding: 0


## Curve Coefficients

Intercept: 0  
Slope: 1.465E+05

## Error Coefficients

Standard Error: 9810  
Relative Standard Error: 6.3  
Correlation Coefficient: 0.999  
Coefficient of Determination (Adjusted): 0.995

| ID | Level            | Concentration | Response | IS Amount | IS Response | RF            | Used |
|----|------------------|---------------|----------|-----------|-------------|---------------|------|
| 1  | IC 280-599427/18 | 0.02          | 3042.0   |           |             | 152100.0      | Y    |
| 2  | IC 280-599427/17 | 0.05          | 8175.0   |           |             | 163500.0      | Y    |
| 3  | IC 280-599427/16 | 0.1           | 14908.0  |           |             | 149080.0      | Y    |
| 4  | IC 280-599427/15 | 0.25          | 34670.0  |           |             | 138680.0      | Y    |
| 5  | IC 280-599427/14 | 0.4           | 59001.0  |           |             | 147502.5      | Y    |
| 6  | IC 280-599427/13 | 0.7           | 95423.0  |           |             | 136318.571429 | Y    |
| 7  | IC 280-599427/12 | 1.0           | 148050.0 |           |             | 148050.0      | Y    |
| 8  | IC 280-599427/11 | 2.5           | 341365.0 |           |             | 136546.0      | Y    |



## Calibration

/ 3,5-Dinitroaniline

Curve Type: Linear  
Weighting: Conc\_Sq  
Origin: None  
Dependency: Response  
Calib Mode: ESTD  
Response Base: AREA  
RF Rounding: 0


## Curve Coefficients

Intercept: -392.4  
Slope: 2.454E+05

## Error Coefficients

Standard Error: 1570  
Relative Standard Error: 3.0  
Correlation Coefficient: 1.000  
Coefficient of Determination (Adjusted): 0.999

| ID | Level            | Concentration | Response | IS Amount | IS Response | RF            | Used |
|----|------------------|---------------|----------|-----------|-------------|---------------|------|
| 1  | IC 280-599427/18 | 0.02          | 4397.0   |           |             | 219850.0      | Y    |
| 2  | IC 280-599427/17 | 0.05          | 12600.0  |           |             | 252000.0      | Y    |
| 3  | IC 280-599427/16 | 0.1           | 24604.0  |           |             | 246040.0      | Y    |
| 4  | IC 280-599427/15 | 0.25          | 59459.0  |           |             | 237836.0      | Y    |
| 5  | IC 280-599427/14 | 0.4           | 96632.0  |           |             | 241580.0      | Y    |
| 6  | IC 280-599427/13 | 0.7           | 169039.0 |           |             | 241484.285714 | Y    |
| 7  | IC 280-599427/12 | 1.0           | 245164.0 |           |             | 245164.0      | Y    |
| 8  | IC 280-599427/11 | 2.5           | 610953.0 |           |             | 244381.2      | Y    |



FORM VI  
HPLC/IC BY EXTERNAL STANDARD - INITIAL CALIBRATION DATA  
RETENTION TIME SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1 Analy Batch No.: 601664

SDG No.: \_\_\_\_\_

Instrument ID: CHHPLC\_X3 GC Column: UltraCarb5u ID: 4.6(mm) Heated Purge: (Y/N) N

Calibration Start Date: 02/08/2023 15:38 Calibration End Date: 02/08/2023 18:42 Calibration ID: 76877

Calibration Files

|         |                  |              |
|---------|------------------|--------------|
| LEVEL:  | LAB SAMPLE ID:   | LAB FILE ID: |
| Level 1 | IC 280-601664/19 | 02080019.D   |
| Level 2 | IC 280-601664/18 | 02080018.D   |
| Level 3 | IC 280-601664/17 | 02080017.D   |
| Level 4 | IC 280-601664/16 | 02080016.D   |
| Level 5 | IC 280-601664/15 | 02080015.D   |
| Level 6 | IC 280-601664/14 | 02080014.D   |
| Level 7 | IC 280-601664/13 | 02080013.D   |
| Level 8 | IC 280-601664/12 | 02080012.D   |
| Level 9 | IC 280-601664/11 | 02080011.D   |

| ANALYTE                    | LVL 1  | LVL 2  | LVL 3  | LVL 4  | LVL 5  | LVL 6  | LVL 7  | LVL 8  | LVL 9  |  | RT WINDOW       | AVG RT |
|----------------------------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--|-----------------|--------|
| HMX                        | 6.550  | 6.549  | 6.547  | 6.544  | 6.548  | 6.548  | 6.544  | 6.547  | 6.546  |  | 6.394 - 6.694   | 6.547  |
| RDX                        | 7.563  | 7.563  | 7.567  | 7.564  | 7.568  | 7.562  | 7.564  | 7.567  | 7.559  |  | 7.414 - 7.714   | 7.564  |
| Picric acid                | 8.050  | 8.049  | 8.047  | 8.044  | 8.041  | 8.028  | 8.018  | 8.014  | 7.966  |  | 7.894 - 8.194   | 8.029  |
| 1,3,5-Trinitrobenzene      | 8.643  | 8.643  | 8.647  | 8.644  | 8.641  | 8.642  | 8.638  | 8.647  | 8.639  |  | 8.494 - 8.794   | 8.643  |
| 1,3-Dinitrobenzene         | 9.263  | 9.256  | 9.260  | 9.257  | 9.255  | 9.255  | 9.251  | 9.260  | 9.246  |  | 9.107 - 9.407   | 9.256  |
| Nitrobenzene               | 9.636  | 9.636  | 9.633  | 9.631  | 9.628  | 9.628  | 9.624  | 9.634  | 9.619  |  | 9.481 - 9.781   | 9.630  |
| 3,5-Dinitroaniline         | 9.836  | 9.836  | 9.840  | 9.831  | 9.828  | 9.828  | 9.824  | 9.840  | 9.819  |  | 9.681 - 9.981   | 9.831  |
| Tetryl                     | 9.983  | 9.982  | 9.987  | 9.977  | 9.975  | 9.968  | 9.971  | 9.987  | 9.966  |  | 9.827 - 10.127  | 9.977  |
| Nitroglycerin              | 10.450 | 10.449 | 10.453 | 10.437 | 10.441 | 10.435 | 10.437 | 10.454 | 10.426 |  | 10.287 - 10.587 | 10.442 |
| 2,4,6-Trinitrotoluene      | 10.876 | 10.876 | 10.880 | 10.864 | 10.868 | 10.868 | 10.864 | 10.880 | 10.859 |  | 10.764 - 10.964 | 10.871 |
| 4-Amino-2,6-dinitrotoluene | 11.043 | 11.042 | 11.040 | 11.031 | 11.028 | 11.022 | 11.024 | 11.040 | 11.013 |  | 10.931 - 11.131 | 11.031 |
| 2-Amino-4,6-dinitrotoluene | 11.296 | 11.289 | 11.293 | 11.277 | 11.275 | 11.275 | 11.271 | 11.287 | 11.259 |  | 11.177 - 11.377 | 11.280 |
| 2,6-Dinitrotoluene         | 11.470 | 11.469 | 11.467 | 11.451 | 11.455 | 11.448 | 11.451 | 11.467 | 11.439 |  | 11.351 - 11.551 | 11.457 |
| 2,4-Dinitrotoluene         | 11.636 | 11.636 | 11.633 | 11.617 | 11.621 | 11.615 | 11.617 | 11.634 | 11.606 |  | 11.517 - 11.717 | 11.624 |
| 2-Nitrotoluene             | 12.456 | 12.456 | 12.453 | 12.437 | 12.441 | 12.435 | 12.431 | 12.447 | 12.426 |  | 12.287 - 12.587 | 12.442 |
| 4-Nitrotoluene             | 12.863 | 12.869 | 12.867 | 12.844 | 12.848 | 12.842 | 12.844 | 12.860 | 12.833 |  | 12.694 - 12.994 | 12.852 |
| 3-Nitrotoluene             | 13.430 | 13.442 | 13.433 | 13.411 | 13.415 | 13.408 | 13.404 | 13.420 | 13.399 |  | 13.261 - 13.561 | 13.418 |
| PETN                       | 14.583 | 14.589 | 14.580 | 14.551 | 14.555 | 14.548 | 14.544 | 14.567 | 14.539 |  | 14.401 - 14.701 | 14.562 |
| 1,2-Dinitrobenzene         | 8.517  | 8.516  | 8.520  | 8.517  | 8.515  | 8.515  | 8.511  | 8.521  | 8.506  |  | 8.367 - 8.667   | 8.515  |

FORM VI  
HPLC/IC BY EXTERNAL STANDARD - INITIAL CALIBRATION DATA  
CURVE EVALUATION

Lab Name: Eurofins Denver Job No.: 280-176674-1 Analy Batch No.: 601664  
SDG No.: \_\_\_\_\_  
Instrument ID: CHHPLC\_X3 GC Column: UltraCarb5u ID: 4.6(mm) Heated Purge: (Y/N) N  
Calibration Start Date: 02/08/2023 15:38 Calibration End Date: 02/08/2023 18:42 Calibration ID: 76877

Calibration Files

|         |                  |              |
|---------|------------------|--------------|
| LEVEL:  | LAB SAMPLE ID:   | LAB FILE ID: |
| Level 1 | IC 280-601664/19 | 02080019.D   |
| Level 2 | IC 280-601664/18 | 02080018.D   |
| Level 3 | IC 280-601664/17 | 02080017.D   |
| Level 4 | IC 280-601664/16 | 02080016.D   |
| Level 5 | IC 280-601664/15 | 02080015.D   |
| Level 6 | IC 280-601664/14 | 02080014.D   |
| Level 7 | IC 280-601664/13 | 02080013.D   |
| Level 8 | IC 280-601664/12 | 02080012.D   |
| Level 9 | IC 280-601664/11 | 02080011.D   |

| ANALYTE               | CF                         |                  |                  |                  | CURVE TYPE | COEFFICIENT    |                |    | # | MIN CF | %RSD /RSE | # | MAX %RSD /RSE | R^2 OR COD | # | MIN R^2 OR COD |
|-----------------------|----------------------------|------------------|------------------|------------------|------------|----------------|----------------|----|---|--------|-----------|---|---------------|------------|---|----------------|
|                       | LVL 1<br>LVL 5<br>LVL 9    | LVL 2<br>LVL 6   | LVL 3<br>LVL 7   | LVL 4<br>LVL 8   |            | B              | M1             | M2 |   |        |           |   |               |            |   |                |
| HMX                   | 96200<br>93960<br>91737    | 98050<br>92130   | 92960<br>91317   | 93870<br>91529   | Ave        |                | 93528.149<br>2 |    |   |        | 2.5       |   | 20.0          |            |   |                |
| RDX                   | 120300<br>102168<br>106491 | 109800<br>103248 | 104720<br>102589 | 101120<br>106983 | Ave        |                | 106379.80<br>8 |    |   |        | 5.5       |   | 20.0          |            |   |                |
| Picric acid           | 67300<br>77436<br>78734    | 74600<br>77483   | 75360<br>77481   | 76400<br>77676   | Ave        |                | 75829.947<br>6 |    |   |        | 4.5       |   | 20.0          |            |   |                |
| 1,3,5-Trinitrobenzene | 219600<br>217528<br>217838 | 215650<br>218445 | 213580<br>218654 | 213990<br>219038 | Ave        |                | 217147.03<br>2 |    |   |        | 1.0       |   | 20.0          |            |   |                |
| 1,3-Dinitrobenzene    | 290200<br>296196<br>297242 | 292450<br>296278 | 293280<br>294821 | 292490<br>296615 | Ave        |                | 294396.92<br>5 |    |   |        | 0.8       |   | 20.0          |            |   |                |
| Nitrobenzene          | 190500<br>190360<br>192956 | 190800<br>192660 | 192100<br>189984 | 190190<br>191651 | Ave        |                | 191244.54<br>3 |    |   |        | 0.6       |   | 20.0          |            |   |                |
| 3,5-Dinitroaniline    | 214700<br>228212<br>230510 | 218850<br>227498 | 228520<br>228160 | 227230<br>229927 | Lin2       | -154.2826<br>1 | 229201.89<br>0 |    |   |        |           |   |               | 1.0000     |   | 0.9900         |
| Tetryl                | 161100<br>165680<br>166480 | 153650<br>171658 | 163060<br>165663 | 164380<br>165416 | Ave        |                | 164120.70<br>6 |    |   |        | 3.0       |   | 20.0          |            |   |                |

Note: The M1 coefficient is the same as Ave CF for an Ave curve type. RSD is calculated for Ave curve types. RSE is used for all other types.

FORM VI  
HPLC/IC BY EXTERNAL STANDARD - INITIAL CALIBRATION DATA  
CURVE EVALUATION

Lab Name: Eurofins Denver Job No.: 280-176674-1 Analy Batch No.: 601664

SDG No.: \_\_\_\_\_

Instrument ID: CHHPLC\_X3 GC Column: UltraCarb5u ID: 4.6(mm) Heated Purge: (Y/N) N

Calibration Start Date: 02/08/2023 15:38 Calibration End Date: 02/08/2023 18:42 Calibration ID: 76877

| ANALYTE                    | CF                         |                  |                  |                  | CURVE<br>TYPE | COEFFICIENT |                |    | # | MIN CF | %RSD<br>/RSE | # | MAX<br>%RSD<br>/RSE | R^2<br>OR COD | # | MIN R^2<br>OR COD |
|----------------------------|----------------------------|------------------|------------------|------------------|---------------|-------------|----------------|----|---|--------|--------------|---|---------------------|---------------|---|-------------------|
|                            | LVL 1<br>LVL 5<br>LVL 9    | LVL 2<br>LVL 6   | LVL 3<br>LVL 7   | LVL 4<br>LVL 8   |               | B           | M1             | M2 |   |        |              |   |                     |               |   |                   |
| Nitroglycerin              | 61590<br>64985<br>65842    | 59395<br>66398   | 63298<br>66047   | 63170<br>65903   | Ave           |             | 64069.924<br>3 |    |   |        | 3.7          |   | 20.0                |               |   |                   |
| 2,4,6-Trinitrotoluene      | 216700<br>210476<br>209648 | 213100<br>208335 | 214160<br>208960 | 207870<br>210110 | Ave           |             | 211039.88<br>9 |    |   |        | 1.4          |   | 20.0                |               |   |                   |
| 4-Amino-2,6-dinitrotoluene | 167000<br>152884<br>152878 | 155300<br>152868 | 157500<br>150916 | 154070<br>150980 | Ave           |             | 154932.80<br>2 |    |   |        | 3.2          |   | 20.0                |               |   |                   |
| 2-Amino-4,6-dinitrotoluene | 221400<br>198056<br>200205 | 201050<br>198963 | 201440<br>196574 | 197650<br>197356 | Ave           |             | 201410.39<br>8 |    |   |        | 3.8          |   | 20.0                |               |   |                   |
| 2,6-Dinitrotoluene         | 148600<br>143452<br>140114 | 145550<br>140035 | 142180<br>142239 | 139670<br>142869 | Ave           |             | 142745.39<br>7 |    |   |        | 2.0          |   | 20.0                |               |   |                   |
| 2,4-Dinitrotoluene         | 314100<br>293052<br>294340 | 299250<br>294023 | 297440<br>290757 | 295390<br>291656 | Ave           |             | 296667.47<br>1 |    |   |        | 2.4          |   | 20.0                |               |   |                   |
| 2-Nitrotoluene             | 137700<br>124752<br>126235 | 130550<br>125108 | 132460<br>124116 | 125480<br>124668 | Ave           |             | 127896.49<br>0 |    |   |        | 3.7          |   | 20.0                |               |   |                   |
| 4-Nitrotoluene             | 126300<br>108612<br>108849 | 119650<br>108228 | 111680<br>107359 | 108570<br>107671 | Ave           |             | 111879.80<br>8 |    |   |        | 5.9          |   | 20.0                |               |   |                   |
| 3-Nitrotoluene             | 147300<br>137292<br>138526 | 153100<br>136790 | 141920<br>136640 | 136660<br>136203 | Ave           |             | 140492.28<br>9 |    |   |        | 4.2          |   | 20.0                |               |   |                   |
| PETN                       | 72730<br>68498<br>68809    | 67310<br>68653   | 68114<br>68586   | 68161<br>68748   | Ave           |             | 68845.425<br>1 |    |   |        | 2.2          |   | 20.0                |               |   |                   |
| 1,2-Dinitrobenzene         | 118800<br>127432<br>129675 | 124550<br>128775 | 125560<br>126837 | 127270<br>127879 | Ave           |             | 126308.70<br>5 |    |   |        | 2.5          |   | 20.0                |               |   |                   |

Note: The M1 coefficient is the same as Ave CF for an Ave curve type. RSD is calculated for Ave curve types. RSE is used for all other types.

FORM VI  
HPLC/IC BY EXTERNAL STANDARD - INITIAL CALIBRATION DATA  
RESPONSE AND CONCENTRATION

Lab Name: Eurofins Denver Job No.: 280-176674-1 Analy Batch No.: 601664

SDG No.: \_\_\_\_\_

Instrument ID: CHHPLC\_X3 GC Column: UltraCarb5u ID: 4.6(mm) Heated Purge: (Y/N) N

Calibration Start Date: 02/08/2023 15:38 Calibration End Date: 02/08/2023 18:42 Calibration ID: 76877

Calibration Files

|         |                  |              |
|---------|------------------|--------------|
| LEVEL:  | LAB SAMPLE ID:   | LAB FILE ID: |
| Level 1 | IC 280-601664/19 | 02080019.D   |
| Level 2 | IC 280-601664/18 | 02080018.D   |
| Level 3 | IC 280-601664/17 | 02080017.D   |
| Level 4 | IC 280-601664/16 | 02080016.D   |
| Level 5 | IC 280-601664/15 | 02080015.D   |
| Level 6 | IC 280-601664/14 | 02080014.D   |
| Level 7 | IC 280-601664/13 | 02080013.D   |
| Level 8 | IC 280-601664/12 | 02080012.D   |
| Level 9 | IC 280-601664/11 | 02080011.D   |

| ANALYTE                    | CURVE TYPE | RESPONSE       |                |                |                |        | CONCENTRATION (UG/ML) |                |                |                |       |
|----------------------------|------------|----------------|----------------|----------------|----------------|--------|-----------------------|----------------|----------------|----------------|-------|
|                            |            | LVL 1<br>LVL 6 | LVL 2<br>LVL 7 | LVL 3<br>LVL 8 | LVL 4<br>LVL 9 | LVL 5  | LVL 1<br>LVL 6        | LVL 2<br>LVL 7 | LVL 3<br>LVL 8 | LVL 4<br>LVL 9 | LVL 5 |
| HMX                        | Ave        | 962            | 1961           | 4648           | 9387           | 23490  | 0.0100                | 0.0200         | 0.0500         | 0.100          | 0.250 |
|                            |            | 36852          | 63922          | 91529          | 229343         |        | 0.400                 | 0.700          | 1.00           | 2.50           |       |
| RDX                        | Ave        | 1203           | 2196           | 5236           | 10112          | 25542  | 0.0100                | 0.0200         | 0.0500         | 0.100          | 0.250 |
|                            |            | 41299          | 71812          | 106983         | 266228         |        | 0.400                 | 0.700          | 1.00           | 2.50           |       |
| Picric acid                | Ave        | 673            | 1492           | 3768           | 7640           | 19359  | 0.0100                | 0.0200         | 0.0500         | 0.100          | 0.250 |
|                            |            | 30993          | 54237          | 77676          | 196834         |        | 0.400                 | 0.700          | 1.00           | 2.50           |       |
| 1,3,5-Trinitrobenzene      | Ave        | 2196           | 4313           | 10679          | 21399          | 54382  | 0.0100                | 0.0200         | 0.0500         | 0.100          | 0.250 |
|                            |            | 87378          | 153058         | 219038         | 544595         |        | 0.400                 | 0.700          | 1.00           | 2.50           |       |
| 1,3-Dinitrobenzene         | Ave        | 2902           | 5849           | 14664          | 29249          | 74049  | 0.0100                | 0.0200         | 0.0500         | 0.100          | 0.250 |
|                            |            | 118511         | 206375         | 296615         | 743106         |        | 0.400                 | 0.700          | 1.00           | 2.50           |       |
| Nitrobenzene               | Ave        | 1905           | 3816           | 9605           | 19019          | 47590  | 0.0100                | 0.0200         | 0.0500         | 0.100          | 0.250 |
|                            |            | 77064          | 132989         | 191651         | 482389         |        | 0.400                 | 0.700          | 1.00           | 2.50           |       |
| 3,5-Dinitroaniline         | Lin2       | 2147           | 4377           | 11426          | 22723          | 57053  | 0.0100                | 0.0200         | 0.0500         | 0.100          | 0.250 |
|                            |            | 90999          | 159712         | 229927         | 576276         |        | 0.400                 | 0.700          | 1.00           | 2.50           |       |
| Tetryl                     | Ave        | 1611           | 3073           | 8153           | 16438          | 41420  | 0.0100                | 0.0200         | 0.0500         | 0.100          | 0.250 |
|                            |            | 68663          | 115964         | 165416         | 416200         |        | 0.400                 | 0.700          | 1.00           | 2.50           |       |
| Nitroglycerin              | Ave        | 6159           | 11879          | 31649          | 63170          | 162463 | 0.100                 | 0.200          | 0.500          | 1.00           | 2.50  |
|                            |            | 265593         | 462332         | 659032         | 1646056        |        | 4.00                  | 7.00           | 10.0           | 25.0           |       |
| 2,4,6-Trinitrotoluene      | Ave        | 2167           | 4262           | 10708          | 20787          | 52619  | 0.0100                | 0.0200         | 0.0500         | 0.100          | 0.250 |
|                            |            | 83334          | 146272         | 210110         | 524120         |        | 0.400                 | 0.700          | 1.00           | 2.50           |       |
| 4-Amino-2,6-dinitrotoluene | Ave        | 1670           | 3106           | 7875           | 15407          | 38221  | 0.0100                | 0.0200         | 0.0500         | 0.100          | 0.250 |
|                            |            | 61147          | 105641         | 150980         | 382195         |        | 0.400                 | 0.700          | 1.00           | 2.50           |       |
| 2-Amino-4,6-dinitrotoluene | Ave        | 2214           | 4021           | 10072          | 19765          | 49514  | 0.0100                | 0.0200         | 0.0500         | 0.100          | 0.250 |
|                            |            | 79585          | 137602         | 197356         | 500512         |        | 0.400                 | 0.700          | 1.00           | 2.50           |       |
| 2,6-Dinitrotoluene         | Ave        | 1486           | 2911           | 7109           | 13967          | 35863  | 0.0100                | 0.0200         | 0.0500         | 0.100          | 0.250 |

FORM VI  
HPLC/IC BY EXTERNAL STANDARD - INITIAL CALIBRATION DATA  
RESPONSE AND CONCENTRATION

Lab Name: Eurofins Denver Job No.: 280-176674-1 Analy Batch No.: 601664

SDG No.: \_\_\_\_\_

Instrument ID: CHHPLC\_X3 GC Column: UltraCarb5u ID: 4.6(mm) Heated Purge: (Y/N) N

Calibration Start Date: 02/08/2023 15:38 Calibration End Date: 02/08/2023 18:42 Calibration ID: 76877

| ANALYTE            | CURVE<br>TYPE | RESPONSE       |                 |                 |                  |        | CONCENTRATION (UG/ML) |                 |                |                |       |
|--------------------|---------------|----------------|-----------------|-----------------|------------------|--------|-----------------------|-----------------|----------------|----------------|-------|
|                    |               | LVL 1<br>LVL 6 | LVL 2<br>LVL 7  | LVL 3<br>LVL 8  | LVL 4<br>LVL 9   | LVL 5  | LVL 1<br>LVL 6        | LVL 2<br>LVL 7  | LVL 3<br>LVL 8 | LVL 4<br>LVL 9 | LVL 5 |
|                    |               | 56014          | 99567           | 142869          | 350285           |        | 0.400                 | 0.700           | 1.00           | 2.50           |       |
| 2,4-Dinitrotoluene | Ave           | 3141<br>117609 | 5985<br>203530  | 14872<br>291656 | 29539<br>735849  | 73263  | 0.0100<br>0.400       | 0.0200<br>0.700 | 0.0500<br>1.00 | 0.100<br>2.50  | 0.250 |
| 2-Nitrotoluene     | Ave           | 1377<br>50043  | 2611<br>86881   | 6623<br>124668  | 12548<br>315588  | 31188  | 0.0100<br>0.400       | 0.0200<br>0.700 | 0.0500<br>1.00 | 0.100<br>2.50  | 0.250 |
| 4-Nitrotoluene     | Ave           | 1263<br>43291  | 2393<br>75151   | 5584<br>107671  | 10857<br>272123  | 27153  | 0.0100<br>0.400       | 0.0200<br>0.700 | 0.0500<br>1.00 | 0.100<br>2.50  | 0.250 |
| 3-Nitrotoluene     | Ave           | 1473<br>54716  | 3062<br>95648   | 7096<br>136203  | 13666<br>346314  | 34323  | 0.0100<br>0.400       | 0.0200<br>0.700 | 0.0500<br>1.00 | 0.100<br>2.50  | 0.250 |
| PETN               | Ave           | 7273<br>274612 | 13462<br>480104 | 34057<br>687477 | 68161<br>1720221 | 171245 | 0.100<br>4.00         | 0.200<br>7.00   | 0.500<br>10.0  | 1.00<br>25.0   | 2.50  |
| 1,2-Dinitrobenzene | Ave           | 1188<br>51510  | 2491<br>88786   | 6278<br>127879  | 12727<br>324188  | 31858  | 0.0100<br>0.400       | 0.0200<br>0.700 | 0.0500<br>1.00 | 0.100<br>2.50  | 0.250 |

Curve Type Legend

Ave = Average  
Lin2 = Linear 1/conc^2

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230208-118465.b\02080011.D  
 Lims ID: IC INT 9  
 Client ID:  
 Sample Type: IC Calib Level: 9  
 Inject. Date: 08-Feb-2023 15:38:53 ALS Bottle#: 11 Worklist Smp#: 11  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: IC INT 9  
 Operator ID: JZ/MAR Instrument ID: CHHPLC\_X3  
 Sublist: chrom-8330\_X3\*sub26  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230208-118465.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 09-Feb-2023 12:48:02 Calib Date: 08-Feb-2023 18:42:13  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230208-118465.b\02080019.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1604

First Level Reviewer: LV5D

Date: 08-Feb-2023 16:28:55

| Compound                      | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/mL | OnCol Amt<br>ug/mL | Flags |
|-------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 4 HMX                         | 1   | 6.546        | 6.544            | 0.002            | 229343   | 2.50             | 2.45               | M     |
| 8 RDX                         | 1   | 7.559        | 7.564            | -0.005           | 266228   | 2.50             | 2.50               |       |
| 9 2,4,6-Trinitrophenol        | 1   | 7.966        | 8.044            | -0.078           | 196834   | 2.50             | 2.60               |       |
| \$ 10 1,2-Dinitrobenzene      | 1   | 8.506        | 8.517            | -0.011           | 324188   | 2.50             | 2.57               |       |
| 11 1,3,5-Trinitrobenzene      | 1   | 8.639        | 8.644            | -0.005           | 544595   | 2.50             | 2.51               |       |
| 12 1,3-Dinitrobenzene         | 1   | 9.246        | 9.257            | -0.011           | 743106   | 2.50             | 2.52               |       |
| 13 Nitrobenzene               | 1   | 9.619        | 9.631            | -0.012           | 482389   | 2.50             | 2.52               |       |
| 14 3,5-Dinitroaniline         | 1   | 9.819        | 9.831            | -0.012           | 576276   | 2.50             | 2.51               |       |
| 15 Tetryl                     | 1   | 9.966        | 9.977            | -0.011           | 416200   | 2.50             | 2.54               |       |
| 16 Nitroglycerin              | 2   | 10.426       | 10.437           | -0.011           | 1646056  | 25.0             | 25.7               |       |
| 17 2,4,6-Trinitrotoluene      | 1   | 10.859       | 10.864           | -0.005           | 524120   | 2.50             | 2.48               |       |
| 18 4-Amino-2,6-dinitrotoluene | 1   | 11.013       | 11.031           | -0.018           | 382195   | 2.50             | 2.47               |       |
| 19 2-Amino-4,6-dinitrotoluene | 1   | 11.259       | 11.277           | -0.018           | 500512   | 2.50             | 2.49               |       |
| 20 2,6-Dinitrotoluene         | 1   | 11.439       | 11.451           | -0.012           | 350285   | 2.50             | 2.45               |       |
| 21 2,4-Dinitrotoluene         | 1   | 11.606       | 11.617           | -0.011           | 735849   | 2.50             | 2.48               |       |
| 22 o-Nitrotoluene             | 1   | 12.426       | 12.437           | -0.011           | 315588   | 2.50             | 2.47               |       |
| 23 p-Nitrotoluene             | 1   | 12.833       | 12.844           | -0.011           | 272123   | 2.50             | 2.43               |       |
| 24 m-Nitrotoluene             | 1   | 13.399       | 13.411           | -0.012           | 346314   | 2.50             | 2.47               |       |
| 25 PETN                       | 2   | 14.539       | 14.551           | -0.012           | 1720221  | 25.0             | 25.0               |       |

**QC Flag Legend**

Processing Flags

Review Flags

M - Manually Integrated

**Reagents:**

8330IntermStk\_00075

Amount Added: 250.00

Units: uL



Report Date: 09-Feb-2023 12:48:02

Chrom Revision: 2.3 01-Feb-2023 13:23:06

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230208-118465.b\02080011.d

Injection Date: 08-Feb-2023 15:38:53

Instrument ID: CHHPLC\_X3

Operator ID: JZ/MAR

Lims ID: IC INT 9

Worklist Smp#: 11

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 11

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230208-118465.b\02080011.d  
Injection Date: 08-Feb-2023 15:38:53 Instrument ID: CHHPLC\_X3  
Lims ID: IC INT 9  
Client ID:  
Operator ID: JZ/MAR ALS Bottle#: 11 Worklist Smp#: 11  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**4 HMX, CAS: 2691-41-0**

Signal: 1


RT: 6.55  
Area: 245170  
Amount: 2.586222  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.55  
Area: 229343  
Amount: 2.452128  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 08-Feb-2023 16:28:37

Audit Action: Manually Integrated

Audit Reason: Baseline Smoothing

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230208-118465.b\02080012.D  
 Lims ID: IC INT 8  
 Client ID:  
 Sample Type: IC Calib Level: 8  
 Inject. Date: 08-Feb-2023 16:01:51 ALS Bottle#: 12 Worklist Smp#: 12  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: IC INT 8  
 Operator ID: JZ/MAR Instrument ID: CHHPLC\_X3  
 Sublist: chrom-8330\_X3\*sub26  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230208-118465.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 09-Feb-2023 12:48:03 Calib Date: 08-Feb-2023 18:42:13  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230208-118465.b\02080019.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1604

First Level Reviewer: LV5D

Date: 08-Feb-2023 16:28:33

| Compound                      | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/mL | OnCol Amt<br>ug/mL | Flags |
|-------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 4 HMX                         | 1   | 6.547        | 6.544            | 0.003            | 91529    | 1.00             | 0.9786             | M     |
| 8 RDX                         | 1   | 7.567        | 7.564            | 0.003            | 106983   | 1.00             | 1.01               |       |
| 9 2,4,6-Trinitrophenol        | 1   | 8.014        | 8.044            | -0.030           | 77676    | 1.00             | 1.02               |       |
| \$ 10 1,2-Dinitrobenzene      | 1   | 8.521        | 8.517            | 0.004            | 127879   | 1.00             | 1.01               |       |
| 11 1,3,5-Trinitrobenzene      | 1   | 8.647        | 8.644            | 0.003            | 219038   | 1.00             | 1.01               |       |
| 12 1,3-Dinitrobenzene         | 1   | 9.260        | 9.257            | 0.003            | 296615   | 1.00             | 1.01               |       |
| 13 Nitrobenzene               | 1   | 9.634        | 9.631            | 0.003            | 191651   | 1.00             | 1.00               |       |
| 14 3,5-Dinitroaniline         | 1   | 9.840        | 9.831            | 0.009            | 229927   | 1.00             | 1.00               |       |
| 15 Tetryl                     | 1   | 9.987        | 9.977            | 0.010            | 165416   | 1.00             | 1.01               |       |
| 16 Nitroglycerin              | 2   | 10.454       | 10.437           | 0.017            | 659032   | 10.0             | 10.3               |       |
| 17 2,4,6-Trinitrotoluene      | 1   | 10.880       | 10.864           | 0.016            | 210110   | 1.00             | 1.00               |       |
| 18 4-Amino-2,6-dinitrotoluene | 1   | 11.040       | 11.031           | 0.009            | 150980   | 1.00             | 0.9745             |       |
| 19 2-Amino-4,6-dinitrotoluene | 1   | 11.287       | 11.277           | 0.010            | 197356   | 1.00             | 0.9799             |       |
| 20 2,6-Dinitrotoluene         | 1   | 11.467       | 11.451           | 0.016            | 142869   | 1.00             | 1.00               |       |
| 21 2,4-Dinitrotoluene         | 1   | 11.634       | 11.617           | 0.017            | 291656   | 1.00             | 0.9831             |       |
| 22 o-Nitrotoluene             | 1   | 12.447       | 12.437           | 0.010            | 124668   | 1.00             | 0.9748             |       |
| 23 p-Nitrotoluene             | 1   | 12.860       | 12.844           | 0.016            | 107671   | 1.00             | 0.9624             |       |
| 24 m-Nitrotoluene             | 1   | 13.420       | 13.411           | 0.009            | 136203   | 1.00             | 0.9695             |       |
| 25 PETN                       | 2   | 14.567       | 14.551           | 0.016            | 687477   | 10.0             | 9.99               |       |

**QC Flag Legend**

Processing Flags

Review Flags

M - Manually Integrated

**Reagents:**

8330IntermStk\_00075

Amount Added: 100.00

Units: uL

Report Date: 09-Feb-2023 12:48:03

Chrom Revision: 2.3 01-Feb-2023 13:23:06

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230208-118465.b\02080012.d

Injection Date: 08-Feb-2023 16:01:51

Instrument ID: CHHPLC\_X3

Operator ID: JZ/MAR

Lims ID: IC INT 8

Worklist Smp#: 12

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 12

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230208-118465.b\02080012.d  
Injection Date: 08-Feb-2023 16:01:51 Instrument ID: CHHPLC\_X3  
Lims ID: IC INT 8  
Client ID:  
Operator ID: JZ/MAR ALS Bottle#: 12 Worklist Smp#: 12  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**4 HMX, CAS: 2691-41-0**

Signal: 1


RT: 6.55  
Area: 107005  
Amount: 1.043580  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.55  
Area: 91529  
Amount: 0.978625  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 08-Feb-2023 16:28:16  
Audit Action: Manually Integrated

Audit Reason: Baseline Smoothing

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230208-118465.b\02080013.D  
 Lims ID: IC INT 7  
 Client ID:  
 Sample Type: IC Calib Level: 7  
 Inject. Date: 08-Feb-2023 16:24:45 ALS Bottle#: 13 Worklist Smp#: 13  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: IC INT 7  
 Operator ID: JZ/MAR Instrument ID: CHHPLC\_X3  
 Sublist: chrom-8330\_X3\*sub26  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230208-118465.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 09-Feb-2023 12:48:04 Calib Date: 08-Feb-2023 18:42:13  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230208-118465.b\02080019.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1604

First Level Reviewer: LV5D

Date: 08-Feb-2023 16:55:38

| Compound                      | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/mL | OnCol Amt<br>ug/mL | Flags |
|-------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 4 HMX                         | 1   | 6.544        | 6.544            | 0.000            | 63922    | 0.7000           | 0.6835             | M     |
| 8 RDX                         | 1   | 7.564        | 7.564            | 0.000            | 71812    | 0.7000           | 0.6751             | M     |
| 9 2,4,6-Trinitrophenol        | 1   | 8.018        | 8.044            | -0.026           | 54237    | 0.7000           | 0.7152             |       |
| \$ 10 1,2-Dinitrobenzene      | 1   | 8.511        | 8.517            | -0.006           | 88786    | 0.7000           | 0.7029             |       |
| 11 1,3,5-Trinitrobenzene      | 1   | 8.638        | 8.644            | -0.006           | 153058   | 0.7000           | 0.7049             |       |
| 12 1,3-Dinitrobenzene         | 1   | 9.251        | 9.257            | -0.006           | 206375   | 0.7000           | 0.7010             |       |
| 13 Nitrobenzene               | 1   | 9.624        | 9.631            | -0.007           | 132989   | 0.7000           | 0.6954             |       |
| 14 3,5-Dinitroaniline         | 1   | 9.824        | 9.831            | -0.007           | 159712   | 0.7000           | 0.6975             |       |
| 15 Tetryl                     | 1   | 9.971        | 9.977            | -0.006           | 115964   | 0.7000           | 0.7066             |       |
| 16 Nitroglycerin              | 2   | 10.437       | 10.437           | 0.000            | 462332   | 7.00             | 7.22               |       |
| 17 2,4,6-Trinitrotoluene      | 1   | 10.864       | 10.864           | 0.000            | 146272   | 0.7000           | 0.6931             |       |
| 18 4-Amino-2,6-dinitrotoluene | 1   | 11.024       | 11.031           | -0.007           | 105641   | 0.7000           | 0.6819             |       |
| 19 2-Amino-4,6-dinitrotoluene | 1   | 11.271       | 11.277           | -0.006           | 137602   | 0.7000           | 0.6832             |       |
| 20 2,6-Dinitrotoluene         | 1   | 11.451       | 11.451           | 0.000            | 99567    | 0.7000           | 0.6975             |       |
| 21 2,4-Dinitrotoluene         | 1   | 11.617       | 11.617           | 0.000            | 203530   | 0.7000           | 0.6861             |       |
| 22 o-Nitrotoluene             | 1   | 12.431       | 12.437           | -0.006           | 86881    | 0.7000           | 0.6793             |       |
| 23 p-Nitrotoluene             | 1   | 12.844       | 12.844           | 0.000            | 75151    | 0.7000           | 0.6717             |       |
| 24 m-Nitrotoluene             | 1   | 13.404       | 13.411           | -0.007           | 95648    | 0.7000           | 0.6808             |       |
| 25 PETN                       | 2   | 14.544       | 14.551           | -0.007           | 480104   | 7.00             | 6.97               |       |

**QC Flag Legend**

Processing Flags

Review Flags

M - Manually Integrated

**Reagents:**

8330IntermStk\_00075

Amount Added: 70.00

Units: uL

Report Date: 09-Feb-2023 12:48:04

Chrom Revision: 2.3 01-Feb-2023 13:23:06

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230208-118465.b\02080013.d

Injection Date: 08-Feb-2023 16:24:45

Instrument ID: CHHPLC\_X3

Operator ID: JZ/MAR

Lims ID: IC INT 7

Worklist Smp#: 13

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 13

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230208-118465.b\02080013.d

Injection Date: 08-Feb-2023 16:24:45

Instrument ID: CHHPLC\_X3

Lims ID: IC INT 7

Client ID:

Operator ID: JZ/MAR

ALS Bottle#:

13

Worklist Smp#: 13

Injection Vol: 100.0 ul

Dil. Factor: 1.0000

Method: 8330\_X3

Limit Group: GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm)


Detector: LC DAD1B, 254 nm

**4 HMX, CAS: 2691-41-0**

Signal: 1


RT: 6.54  
Area: 70424  
Amount: 0.744251  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.54  
Area: 63922  
Amount: 0.683452  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 08-Feb-2023 16:55:19

Audit Action: Manually Integrated

Audit Reason: Baseline Smoothing



## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230208-118465.b\02080013.d  
Injection Date: 08-Feb-2023 16:24:45 Instrument ID: CHHPLC\_X3  
Lims ID: IC INT 7  
Client ID:  
Operator ID: JZ/MAR ALS Bottle#: 13 Worklist Smp#: 13  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**8 RDX, CAS: 121-82-4**

Signal: 1


RT: 7.56  
Area: 75013  
Amount: 0.701791  
Amount Units: ug/mL

## Processing Integration Results



RT: 7.56  
Area: 71812  
Amount: 0.675053  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 08-Feb-2023 19:09:56  
Audit Action: Manually Integrated

Audit Reason: Baseline Smoothing

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230208-118465.b\02080014.D  
 Lims ID: IC INT 6  
 Client ID:  
 Sample Type: IC Calib Level: 6  
 Inject. Date: 08-Feb-2023 16:47:39 ALS Bottle#: 14 Worklist Smp#: 14  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: IC INT 6  
 Operator ID: JZ/MAR Instrument ID: CHHPLC\_X3  
 Sublist: chrom-8330\_X3\*sub26  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230208-118465.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 09-Feb-2023 12:48:05 Calib Date: 08-Feb-2023 18:42:13  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230208-118465.b\02080019.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1604

First Level Reviewer: LV5D

Date: 08-Feb-2023 17:11:22

| Compound                      | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/mL | OnCol Amt<br>ug/mL | Flags |
|-------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 4 HMX                         | 1   | 6.548        | 6.544            | 0.004            | 36852    | 0.4000           | 0.3940             | M     |
| 8 RDX                         | 1   | 7.562        | 7.564            | -0.002           | 41299    | 0.4000           | 0.3882             | M     |
| 9 2,4,6-Trinitrophenol        | 1   | 8.028        | 8.044            | -0.016           | 30993    | 0.4000           | 0.4087             |       |
| \$ 10 1,2-Dinitrobenzene      | 1   | 8.515        | 8.517            | -0.002           | 51510    | 0.4000           | 0.4078             |       |
| 11 1,3,5-Trinitrobenzene      | 1   | 8.642        | 8.644            | -0.002           | 87378    | 0.4000           | 0.4024             |       |
| 12 1,3-Dinitrobenzene         | 1   | 9.255        | 9.257            | -0.002           | 118511   | 0.4000           | 0.4026             |       |
| 13 Nitrobenzene               | 1   | 9.628        | 9.631            | -0.003           | 77064    | 0.4000           | 0.4030             |       |
| 14 3,5-Dinitroaniline         | 1   | 9.828        | 9.831            | -0.003           | 90999    | 0.4000           | 0.3977             |       |
| 15 Tetryl                     | 1   | 9.968        | 9.977            | -0.009           | 68663    | 0.4000           | 0.4184             |       |
| 16 Nitroglycerin              | 2   | 10.435       | 10.437           | -0.002           | 265593   | 4.00             | 4.15               |       |
| 17 2,4,6-Trinitrotoluene      | 1   | 10.868       | 10.864           | 0.004            | 83334    | 0.4000           | 0.3949             |       |
| 18 4-Amino-2,6-dinitrotoluene | 1   | 11.022       | 11.031           | -0.009           | 61147    | 0.4000           | 0.3947             |       |
| 19 2-Amino-4,6-dinitrotoluene | 1   | 11.275       | 11.277           | -0.002           | 79585    | 0.4000           | 0.3951             |       |
| 20 2,6-Dinitrotoluene         | 1   | 11.448       | 11.451           | -0.003           | 56014    | 0.4000           | 0.3924             |       |
| 21 2,4-Dinitrotoluene         | 1   | 11.615       | 11.617           | -0.002           | 117609   | 0.4000           | 0.3964             |       |
| 22 o-Nitrotoluene             | 1   | 12.435       | 12.437           | -0.002           | 50043    | 0.4000           | 0.3913             |       |
| 23 p-Nitrotoluene             | 1   | 12.842       | 12.844           | -0.002           | 43291    | 0.4000           | 0.3869             |       |
| 24 m-Nitrotoluene             | 1   | 13.408       | 13.411           | -0.003           | 54716    | 0.4000           | 0.3895             |       |
| 25 PETN                       | 2   | 14.548       | 14.551           | -0.003           | 274612   | 4.00             | 3.99               |       |

**QC Flag Legend**

Processing Flags

Review Flags

M - Manually Integrated

**Reagents:**

8330IntermStk\_00075

Amount Added: 40.00

Units: uL

Report Date: 09-Feb-2023 12:48:05

Chrom Revision: 2.3 01-Feb-2023 13:23:06

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230208-118465.b\02080014.d

Injection Date: 08-Feb-2023 16:47:39

Instrument ID: CHHPLC\_X3

Operator ID: JZ/MAR

Lims ID: IC INT 6

Worklist Smp#: 14

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 14

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230208-118465.b\02080014.d

Injection Date: 08-Feb-2023 16:47:39

Instrument ID: CHHPLC\_X3

Lims ID: IC INT 6

Client ID:

Operator ID: JZ/MAR

ALS Bottle#: 14

Worklist Smp#: 14

Injection Vol: 100.0 ul

Dil. Factor: 1.0000

Method: 8330\_X3

Limit Group: GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm)


Detector: LC DAD1B, 254 nm

**4 HMX, CAS: 2691-41-0**

Signal: 1


RT: 6.55  
Area: 37332  
Amount: 0.405878  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.55  
Area: 36852  
Amount: 0.394020  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 08-Feb-2023 17:11:10

Audit Action: Manually Integrated

Audit Reason: Baseline Smoothing

## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230208-118465.b\02080014.d  
Injection Date: 08-Feb-2023 16:47:39 Instrument ID: CHHPLC\_X3  
Lims ID: IC INT 6  
Client ID:  
Operator ID: JZ/MAR ALS Bottle#: 14 Worklist Smp#: 14  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**8 RDX, CAS: 121-82-4**

Signal: 1


RT: 7.56  
Area: 43581  
Amount: 0.405322  
Amount Units: ug/mL

## Processing Integration Results



RT: 7.56  
Area: 41299  
Amount: 0.388222  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 08-Feb-2023 19:09:51  
Audit Action: Manually Integrated

Audit Reason: Baseline Smoothing

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230208-118465.b\02080015.D  
 Lims ID: IC INT 5  
 Client ID:  
 Sample Type: IC Calib Level: 5  
 Inject. Date: 08-Feb-2023 17:10:30 ALS Bottle#: 15 Worklist Smp#: 15  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: IC INT 5  
 Operator ID: JZ/MAR Instrument ID: CHHPLC\_X3  
 Sublist: chrom-8330\_X3\*sub26  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230208-118465.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 09-Feb-2023 12:48:05 Calib Date: 08-Feb-2023 18:42:13  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230208-118465.b\02080019.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1604

First Level Reviewer: LV5D

Date: 08-Feb-2023 17:44:54

| Compound                      | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/mL | OnCol Amt<br>ug/mL | Flags |
|-------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 4 HMX                         | 1   | 6.548        | 6.544            | 0.004            | 23490    | 0.2500           | 0.2512             |       |
| 8 RDX                         | 1   | 7.568        | 7.564            | 0.004            | 25542    | 0.2500           | 0.2401             | M     |
| 9 2,4,6-Trinitrophenol        | 1   | 8.041        | 8.044            | -0.003           | 19359    | 0.2500           | 0.2553             |       |
| \$ 10 1,2-Dinitrobenzene      | 1   | 8.515        | 8.517            | -0.002           | 31858    | 0.2500           | 0.2522             |       |
| 11 1,3,5-Trinitrobenzene      | 1   | 8.641        | 8.644            | -0.003           | 54382    | 0.2500           | 0.2504             |       |
| 12 1,3-Dinitrobenzene         | 1   | 9.255        | 9.257            | -0.002           | 74049    | 0.2500           | 0.2515             |       |
| 13 Nitrobenzene               | 1   | 9.628        | 9.631            | -0.003           | 47590    | 0.2500           | 0.2488             |       |
| 14 3,5-Dinitroaniline         | 1   | 9.828        | 9.831            | -0.003           | 57053    | 0.2500           | 0.2496             |       |
| 15 Tetryl                     | 1   | 9.975        | 9.977            | -0.002           | 41420    | 0.2500           | 0.2524             |       |
| 16 Nitroglycerin              | 2   | 10.441       | 10.437           | 0.004            | 162463   | 2.50             | 2.54               |       |
| 17 2,4,6-Trinitrotoluene      | 1   | 10.868       | 10.864           | 0.004            | 52619    | 0.2500           | 0.2493             |       |
| 18 4-Amino-2,6-dinitrotoluene | 1   | 11.028       | 11.031           | -0.003           | 38221    | 0.2500           | 0.2467             |       |
| 19 2-Amino-4,6-dinitrotoluene | 1   | 11.275       | 11.277           | -0.002           | 49514    | 0.2500           | 0.2458             |       |
| 20 2,6-Dinitrotoluene         | 1   | 11.455       | 11.451           | 0.004            | 35863    | 0.2500           | 0.2512             |       |
| 21 2,4-Dinitrotoluene         | 1   | 11.621       | 11.617           | 0.004            | 73263    | 0.2500           | 0.2470             |       |
| 22 o-Nitrotoluene             | 1   | 12.441       | 12.437           | 0.004            | 31188    | 0.2500           | 0.2439             |       |
| 23 p-Nitrotoluene             | 1   | 12.848       | 12.844           | 0.004            | 27153    | 0.2500           | 0.2427             |       |
| 24 m-Nitrotoluene             | 1   | 13.415       | 13.411           | 0.004            | 34323    | 0.2500           | 0.2443             |       |
| 25 PETN                       | 2   | 14.555       | 14.551           | 0.004            | 171245   | 2.50             | 2.49               |       |

**QC Flag Legend**

Processing Flags

Review Flags

M - Manually Integrated

**Reagents:**

8330IntermStk\_00075

Amount Added: 25.00

Units: uL

Report Date: 09-Feb-2023 12:48:06

Chrom Revision: 2.3 01-Feb-2023 13:23:06

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230208-118465.b\02080015.d

Injection Date: 08-Feb-2023 17:10:30

Instrument ID: CHHPLC\_X3

Operator ID: JZ/MAR

Lims ID: IC INT 5

Worklist Smp#: 15

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 15

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230208-118465.b\02080015.d  
Injection Date: 08-Feb-2023 17:10:30 Instrument ID: CHHPLC\_X3  
Lims ID: IC INT 5  
Client ID:  
Operator ID: JZ/MAR ALS Bottle#: 15 Worklist Smp#: 15  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**8 RDX, CAS: 121-82-4**

Signal: 1


RT: 7.57  
Area: 27559  
Amount: 0.254192  
Amount Units: ug/mL

## Processing Integration Results



RT: 7.57  
Area: 25542  
Amount: 0.240102  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 08-Feb-2023 19:09:46  
Audit Action: Manually Integrated

Audit Reason: Baseline Smoothing



Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230208-118465.b\02080016.D  
 Lims ID: IC INT 4  
 Client ID:  
 Sample Type: IC Calib Level: 4  
 Inject. Date: 08-Feb-2023 17:33:25 ALS Bottle#: 16 Worklist Smp#: 16  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: IC INT 4  
 Operator ID: JZ/MAR Instrument ID: CHHPLC\_X3  
 Sublist: chrom-8330\_X3\*sub26  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230208-118465.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 09-Feb-2023 12:48:06 Calib Date: 08-Feb-2023 18:42:13  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230208-118465.b\02080019.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1604

First Level Reviewer: LV5D

Date: 08-Feb-2023 18:05:31

| Compound                      | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/mL | OnCol Amt<br>ug/mL | Flags |
|-------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 4 HMX                         | 1   | 6.544        | 6.544            | 0.000            | 9387     | 0.1000           | 0.1004             |       |
| 8 RDX                         | 1   | 7.564        | 7.564            | 0.000            | 10112    | 0.1000           | 0.0951             | M     |
| 9 2,4,6-Trinitrophenol        | 1   | 8.044        | 8.044            | 0.000            | 7640     | 0.1000           | 0.1008             |       |
| \$ 10 1,2-Dinitrobenzene      | 1   | 8.517        | 8.517            | 0.000            | 12727    | 0.1000           | 0.1008             |       |
| 11 1,3,5-Trinitrobenzene      | 1   | 8.644        | 8.644            | 0.000            | 21399    | 0.1000           | 0.0985             |       |
| 12 1,3-Dinitrobenzene         | 1   | 9.257        | 9.257            | 0.000            | 29249    | 0.1000           | 0.0994             |       |
| 13 Nitrobenzene               | 1   | 9.631        | 9.631            | 0.000            | 19019    | 0.1000           | 0.0994             |       |
| 14 3,5-Dinitroaniline         | 1   | 9.831        | 9.831            | 0.000            | 22723    | 0.1000           | 0.0998             |       |
| 15 Tetryl                     | 1   | 9.977        | 9.977            | 0.000            | 16438    | 0.1000           | 0.1002             |       |
| 16 Nitroglycerin              | 2   | 10.437       | 10.437           | 0.000            | 63170    | 1.00             | 0.9860             |       |
| 17 2,4,6-Trinitrotoluene      | 1   | 10.864       | 10.864           | 0.000            | 20787    | 0.1000           | 0.0985             |       |
| 18 4-Amino-2,6-dinitrotoluene | 1   | 11.031       | 11.031           | 0.000            | 15407    | 0.1000           | 0.0994             |       |
| 19 2-Amino-4,6-dinitrotoluene | 1   | 11.277       | 11.277           | 0.000            | 19765    | 0.1000           | 0.0981             |       |
| 20 2,6-Dinitrotoluene         | 1   | 11.451       | 11.451           | 0.000            | 13967    | 0.1000           | 0.0978             |       |
| 21 2,4-Dinitrotoluene         | 1   | 11.617       | 11.617           | 0.000            | 29539    | 0.1000           | 0.0996             |       |
| 22 o-Nitrotoluene             | 1   | 12.437       | 12.437           | 0.000            | 12548    | 0.1000           | 0.0981             |       |
| 23 p-Nitrotoluene             | 1   | 12.844       | 12.844           | 0.000            | 10857    | 0.1000           | 0.0970             |       |
| 24 m-Nitrotoluene             | 1   | 13.411       | 13.411           | 0.000            | 13666    | 0.1000           | 0.0973             |       |
| 25 PETN                       | 2   | 14.551       | 14.551           | 0.000            | 68161    | 1.00             | 0.99               | M     |

**QC Flag Legend**

Processing Flags

Review Flags

M - Manually Integrated

**Reagents:**

8330IntermStk\_00075

Amount Added: 10.00

Units: uL

Report Date: 09-Feb-2023 12:48:06

Chrom Revision: 2.3 01-Feb-2023 13:23:06

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230208-118465.b\02080016.d

Injection Date: 08-Feb-2023 17:33:25

Instrument ID: CHHPLC\_X3

Operator ID: JZ/MAR

Lims ID: IC INT 4

Worklist Smp#: 16

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 16

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230208-118465.b\02080016.d  
Injection Date: 08-Feb-2023 17:33:25 Instrument ID: CHHPLC\_X3  
Lims ID: IC INT 4  
Client ID:  
Operator ID: JZ/MAR ALS Bottle#: 16 Worklist Smp#: 16  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**8 RDX, CAS: 121-82-4**

Signal: 1


RT: 7.56  
Area: 11152  
Amount: 0.101776  
Amount Units: ug/mL

## Processing Integration Results



RT: 7.56  
Area: 10112  
Amount: 0.095056  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 08-Feb-2023 19:09:38  
Audit Action: Manually Integrated

Audit Reason: Baseline Smoothing

## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230208-118465.b\02080016.d  
Injection Date: 08-Feb-2023 17:33:25 Instrument ID: CHHPLC\_X3  
Lims ID: IC INT 4  
Client ID:  
Operator ID: JZ/MAR ALS Bottle#: 16 Worklist Smp#: 16  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1C, 215 nm

**25 PETN, CAS: 78-11-5**

Signal: 1


RT: 14.55  
Area: 66752  
Amount: 0.971802  
Amount Units: ug/mL

## Processing Integration Results



RT: 14.55  
Area: 68161  
Amount: 0.990059  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 08-Feb-2023 19:08:31

Audit Action: Manually Integrated

Audit Reason: Baseline Smoothing

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230208-118465.b\02080017.D  
 Lims ID: IC INT 3  
 Client ID:  
 Sample Type: IC Calib Level: 3  
 Inject. Date: 08-Feb-2023 17:56:21 ALS Bottle#: 17 Worklist Smp#: 17  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: IC INT 3  
 Operator ID: JZ/MAR Instrument ID: CHHPLC\_X3  
 Sublist: chrom-8330\_X3\*sub26  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230208-118465.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 09-Feb-2023 12:48:07 Calib Date: 08-Feb-2023 18:42:13  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230208-118465.b\02080019.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1604

First Level Reviewer: LV5D

Date: 08-Feb-2023 19:07:50

| Compound                      | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/mL | OnCol Amt<br>ug/mL | Flags |
|-------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 4 HMX                         | 1   | 6.547        | 6.544            | 0.003            | 4648     | 0.0500           | 0.0497             | M     |
| 8 RDX                         | 1   | 7.567        | 7.564            | 0.003            | 5236     | 0.0500           | 0.0492             | M     |
| 9 2,4,6-Trinitrophenol        | 1   | 8.047        | 8.044            | 0.003            | 3768     | 0.0500           | 0.0497             |       |
| \$ 10 1,2-Dinitrobenzene      | 1   | 8.520        | 8.517            | 0.003            | 6278     | 0.0500           | 0.0497             |       |
| 11 1,3,5-Trinitrobenzene      | 1   | 8.647        | 8.644            | 0.003            | 10679    | 0.0500           | 0.0492             |       |
| 12 1,3-Dinitrobenzene         | 1   | 9.260        | 9.257            | 0.003            | 14664    | 0.0500           | 0.0498             |       |
| 13 Nitrobenzene               | 1   | 9.633        | 9.631            | 0.002            | 9605     | 0.0500           | 0.0502             |       |
| 14 3,5-Dinitroaniline         | 1   | 9.840        | 9.831            | 0.009            | 11426    | 0.0500           | 0.0505             |       |
| 15 Tetryl                     | 1   | 9.987        | 9.977            | 0.010            | 8153     | 0.0500           | 0.0497             |       |
| 16 Nitroglycerin              | 2   | 10.453       | 10.437           | 0.016            | 31649    | 0.5000           | 0.4940             |       |
| 17 2,4,6-Trinitrotoluene      | 1   | 10.880       | 10.864           | 0.016            | 10708    | 0.0500           | 0.0507             |       |
| 18 4-Amino-2,6-dinitrotoluene | 1   | 11.040       | 11.031           | 0.009            | 7875     | 0.0500           | 0.0508             |       |
| 19 2-Amino-4,6-dinitrotoluene | 1   | 11.293       | 11.277           | 0.016            | 10072    | 0.0500           | 0.0500             |       |
| 20 2,6-Dinitrotoluene         | 1   | 11.467       | 11.451           | 0.016            | 7109     | 0.0500           | 0.0498             |       |
| 21 2,4-Dinitrotoluene         | 1   | 11.633       | 11.617           | 0.016            | 14872    | 0.0500           | 0.0501             |       |
| 22 o-Nitrotoluene             | 1   | 12.453       | 12.437           | 0.016            | 6623     | 0.0500           | 0.0518             |       |
| 23 p-Nitrotoluene             | 1   | 12.867       | 12.844           | 0.023            | 5584     | 0.0500           | 0.0499             |       |
| 24 m-Nitrotoluene             | 1   | 13.433       | 13.411           | 0.022            | 7096     | 0.0500           | 0.0505             |       |
| 25 PETN                       | 2   | 14.580       | 14.551           | 0.029            | 34057    | 0.5000           | 0.4947             | M     |

**QC Flag Legend**

Processing Flags

Review Flags

M - Manually Integrated

**Reagents:**

8330IntermStk\_00075

Amount Added: 5.00

Units: uL

Report Date: 09-Feb-2023 12:48:07

Chrom Revision: 2.3 01-Feb-2023 13:23:06

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230208-118465.b\02080017.d

Injection Date: 08-Feb-2023 17:56:21

Instrument ID: CHHPLC\_X3

Operator ID: JZ/MAR

Lims ID: IC INT 3

Worklist Smp#: 17

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 17

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230208-118465.b\02080017.d

Injection Date: 08-Feb-2023 17:56:21

Instrument ID: CHHPLC\_X3

Lims ID: IC INT 3

Client ID:

Operator ID: JZ/MAR

ALS Bottle#:

17

Worklist Smp#: 17

Injection Vol: 100.0 ul

Dil. Factor:

1.0000

Method: 8330\_X3

Limit Group:

GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm)

Detector


LC DAD1B, 254 nm

**4 HMX, CAS: 2691-41-0**

Signal: 1


RT: 6.55  
Area: 4942  
Amount: 0.050700  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.55  
Area: 4648  
Amount: 0.049696  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 08-Feb-2023 19:07:46

Audit Action: Manually Integrated

Audit Reason: Baseline Smoothing

## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230208-118465.b\02080017.d  
Injection Date: 08-Feb-2023 17:56:21 Instrument ID: CHHPLC\_X3  
Lims ID: IC INT 3  
Client ID:  
Operator ID: JZ/MAR ALS Bottle#: 17 Worklist Smp#: 17  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**8 RDX, CAS: 121-82-4**

Signal: 1


RT: 7.57  
Area: 6041  
Amount: 0.054246  
Amount Units: ug/mL

## Processing Integration Results



RT: 7.57  
Area: 5236  
Amount: 0.049220  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 08-Feb-2023 19:09:26  
Audit Action: Manually Integrated

Audit Reason: Baseline Smoothing



## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230208-118465.b\02080017.d

Injection Date: 08-Feb-2023 17:56:21

Instrument ID: CHHPLC\_X3

Lims ID: IC INT 3

Client ID:

Operator ID: JZ/MAR

ALS Bottle#:

17

Worklist Smp#: 17

Injection Vol: 100.0 ul

Dil. Factor:

1.0000

Method: 8330\_X3

Limit Group:

GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm)

Detector


LC DAD1C, 215 nm

**25 PETN, CAS: 78-11-5**

Signal: 1


RT: 14.58  
Area: 32403  
Amount: 0.474274  
Amount Units: ug/mL

## Processing Integration Results



RT: 14.58  
Area: 34057  
Amount: 0.494688  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 08-Feb-2023 19:08:23

Audit Action: Manually Integrated

Audit Reason: Baseline Smoothing

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230208-118465.b\02080018.D  
 Lims ID: IC INT 2  
 Client ID:  
 Sample Type: IC Calib Level: 2  
 Inject. Date: 08-Feb-2023 18:19:14 ALS Bottle#: 18 Worklist Smp#: 18  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: IC INT 2  
 Operator ID: JZ/MAR Instrument ID: CHHPLC\_X3  
 Sublist: chrom-8330\_X3\*sub26  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230208-118465.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 09-Feb-2023 12:48:08 Calib Date: 08-Feb-2023 18:42:13  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230208-118465.b\02080019.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1604

First Level Reviewer: LV5D

Date: 08-Feb-2023 19:08:08

| Compound                      | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/mL | OnCol Amt<br>ug/mL | Flags |
|-------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 4 HMX                         | 1   | 6.549        | 6.544            | 0.005            | 1961     | 0.0200           | 0.0210             | M     |
| 8 RDX                         | 1   | 7.563        | 7.564            | -0.001           | 2196     | 0.0200           | 0.0206             | M     |
| 9 2,4,6-Trinitrophenol        | 1   | 8.049        | 8.044            | 0.005            | 1492     | 0.0200           | 0.0197             |       |
| \$ 10 1,2-Dinitrobenzene      | 1   | 8.516        | 8.517            | -0.001           | 2491     | 0.0200           | 0.0197             |       |
| 11 1,3,5-Trinitrobenzene      | 1   | 8.643        | 8.644            | -0.001           | 4313     | 0.0200           | 0.0199             |       |
| 12 1,3-Dinitrobenzene         | 1   | 9.256        | 9.257            | -0.001           | 5849     | 0.0200           | 0.0199             |       |
| 13 Nitrobenzene               | 1   | 9.636        | 9.631            | 0.005            | 3816     | 0.0200           | 0.0200             |       |
| 14 3,5-Dinitroaniline         | 1   | 9.836        | 9.831            | 0.005            | 4377     | 0.0200           | 0.0198             |       |
| 15 Tetryl                     | 1   | 9.982        | 9.977            | 0.005            | 3073     | 0.0200           | 0.0187             |       |
| 16 Nitroglycerin              | 2   | 10.449       | 10.437           | 0.012            | 11879    | 0.2000           | 0.1854             |       |
| 17 2,4,6-Trinitrotoluene      | 1   | 10.876       | 10.864           | 0.012            | 4262     | 0.0200           | 0.0202             |       |
| 18 4-Amino-2,6-dinitrotoluene | 1   | 11.042       | 11.031           | 0.011            | 3106     | 0.0200           | 0.0200             |       |
| 19 2-Amino-4,6-dinitrotoluene | 1   | 11.289       | 11.277           | 0.012            | 4021     | 0.0200           | 0.0200             |       |
| 20 2,6-Dinitrotoluene         | 1   | 11.469       | 11.451           | 0.018            | 2911     | 0.0200           | 0.0204             |       |
| 21 2,4-Dinitrotoluene         | 1   | 11.636       | 11.617           | 0.019            | 5985     | 0.0200           | 0.0202             |       |
| 22 o-Nitrotoluene             | 1   | 12.456       | 12.437           | 0.019            | 2611     | 0.0200           | 0.0204             |       |
| 23 p-Nitrotoluene             | 1   | 12.869       | 12.844           | 0.025            | 2393     | 0.0200           | 0.0214             |       |
| 24 m-Nitrotoluene             | 1   | 13.442       | 13.411           | 0.031            | 3062     | 0.0200           | 0.0218             |       |
| 25 PETN                       | 2   | 14.589       | 14.551           | 0.038            | 13462    | 0.2000           | 0.1955             | M     |

**QC Flag Legend**

Processing Flags

Review Flags

M - Manually Integrated

**Reagents:**

8330IntermStk\_00075

Amount Added: 2.00

Units: uL

Report Date: 09-Feb-2023 12:48:08

Chrom Revision: 2.3 01-Feb-2023 13:23:06

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230208-118465.b\02080018.d

Injection Date: 08-Feb-2023 18:19:14

Instrument ID: CHHPLC\_X3

Operator ID: JZ/MAR

Lims ID: IC INT 2

Worklist Smp#: 18

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 18

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230208-118465.b\02080018.d

Injection Date: 08-Feb-2023 18:19:14

Instrument ID: CHHPLC\_X3

Lims ID: IC INT 2

Client ID:

Operator ID: JZ/MAR

ALS Bottle#:

18

Worklist Smp#: 18

Injection Vol: 100.0 ul

Dil. Factor:

1.0000

Method: 8330\_X3

Limit Group:

GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm)

Detector


LC DAD1B, 254 nm

**4 HMX, CAS: 2691-41-0**

Signal: 1


RT: 6.55  
Area: 2102  
Amount: 0.022288  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.55  
Area: 1961  
Amount: 0.020967  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 08-Feb-2023 19:09:07

Audit Action: Manually Integrated

Audit Reason: Baseline Smoothing

## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230208-118465.b\02080018.d  
Injection Date: 08-Feb-2023 18:19:14 Instrument ID: CHHPLC\_X3  
Lims ID: IC INT 2  
Client ID:  
Operator ID: JZ/MAR ALS Bottle#: 18 Worklist Smp#: 18  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**8 RDX, CAS: 121-82-4**

Signal: 1


RT: 7.56  
Area: 2668  
Amount: 0.023407  
Amount Units: ug/mL

## Processing Integration Results



RT: 7.56  
Area: 2196  
Amount: 0.020643  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 08-Feb-2023 19:09:19  
Audit Action: Manually Integrated

Audit Reason: Baseline Smoothing

## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230208-118465.b\02080018.d

Injection Date: 08-Feb-2023 18:19:14

Instrument ID: CHHPLC\_X3

Lims ID: IC INT 2

Client ID:

Operator ID: JZ/MAR

ALS Bottle#:

18

Worklist Smp#: 18

Injection Vol: 100.0 ul

Dil. Factor: 1.0000

Method: 8330\_X3

Limit Group: GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm)

Detector


LC DAD1C, 215 nm

**25 PETN, CAS: 78-11-5**

Signal: 1


RT: 14.59  
Area: 11186  
Amount: 0.166814  
Amount Units: ug/mL

## Processing Integration Results



RT: 14.59  
Area: 13462  
Amount: 0.195540  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 08-Feb-2023 19:08:19

Audit Action: Manually Integrated

Audit Reason: Baseline Smoothing

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230208-118465.b\02080019.D  
 Lims ID: IC INT 1  
 Client ID:  
 Sample Type: IC Calib Level: 1  
 Inject. Date: 08-Feb-2023 18:42:13 ALS Bottle#: 19 Worklist Smp#: 19  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: IC INT 1  
 Operator ID: JZ/MAR Instrument ID: CHHPLC\_X3  
 Sublist: chrom-8330\_X3\*sub26  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230208-118465.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 09-Feb-2023 12:48:08 Calib Date: 08-Feb-2023 18:42:13  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230208-118465.b\02080019.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1604

First Level Reviewer: LV5D

Date: 08-Feb-2023 19:08:17

| Compound                      | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/mL | OnCol Amt<br>ug/mL | Flags |
|-------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 4 HMX                         | 1   | 6.550        | 6.544            | 0.006            | 962      | 0.0100           | 0.0103             | M     |
| 8 RDX                         | 1   | 7.563        | 7.564            | -0.001           | 1203     | 0.0100           | 0.0113             | M     |
| 9 2,4,6-Trinitrophenol        | 1   | 8.050        | 8.044            | 0.006            | 673      | 0.0100           | 0.008875           |       |
| \$ 10 1,2-Dinitrobenzene      | 1   | 8.517        | 8.517            | 0.000            | 1188     | 0.0100           | 0.009406           |       |
| 11 1,3,5-Trinitrobenzene      | 1   | 8.643        | 8.644            | -0.001           | 2196     | 0.0100           | 0.0101             |       |
| 12 1,3-Dinitrobenzene         | 1   | 9.263        | 9.257            | 0.006            | 2902     | 0.0100           | 0.009857           |       |
| 13 Nitrobenzene               | 1   | 9.636        | 9.631            | 0.005            | 1905     | 0.0100           | 0.0100             |       |
| 14 3,5-Dinitroaniline         | 1   | 9.836        | 9.831            | 0.005            | 2147     | 0.0100           | 0.0100             |       |
| 15 Tetryl                     | 1   | 9.983        | 9.977            | 0.006            | 1611     | 0.0100           | 0.009816           |       |
| 16 Nitroglycerin              | 2   | 10.450       | 10.437           | 0.013            | 6159     | 0.1000           | 0.0961             |       |
| 17 2,4,6-Trinitrotoluene      | 1   | 10.876       | 10.864           | 0.012            | 2167     | 0.0100           | 0.0103             |       |
| 18 4-Amino-2,6-dinitrotoluene | 1   | 11.043       | 11.031           | 0.012            | 1670     | 0.0100           | 0.0108             |       |
| 19 2-Amino-4,6-dinitrotoluene | 1   | 11.296       | 11.277           | 0.019            | 2214     | 0.0100           | 0.0110             |       |
| 20 2,6-Dinitrotoluene         | 1   | 11.470       | 11.451           | 0.019            | 1486     | 0.0100           | 0.0104             |       |
| 21 2,4-Dinitrotoluene         | 1   | 11.636       | 11.617           | 0.019            | 3141     | 0.0100           | 0.0106             |       |
| 22 o-Nitrotoluene             | 1   | 12.456       | 12.437           | 0.019            | 1377     | 0.0100           | 0.0108             |       |
| 23 p-Nitrotoluene             | 1   | 12.863       | 12.844           | 0.019            | 1263     | 0.0100           | 0.0113             |       |
| 24 m-Nitrotoluene             | 1   | 13.430       | 13.411           | 0.019            | 1473     | 0.0100           | 0.0105             |       |
| 25 PETN                       | 2   | 14.583       | 14.551           | 0.032            | 7273     | 0.1000           | 0.1056             | M     |

**QC Flag Legend**

Processing Flags

Review Flags

M - Manually Integrated

**Reagents:**

8330IntermStk\_00075

Amount Added: 1.00

Units: uL

Report Date: 09-Feb-2023 12:48:08

Chrom Revision: 2.3 01-Feb-2023 13:23:06

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230208-118465.b\02080019.d

Injection Date: 08-Feb-2023 18:42:13

Instrument ID: CHHPLC\_X3

Operator ID: JZ/MAR

Lims ID: IC INT 1

Worklist Smp#: 19

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 19

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1





## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230208-118465.b\02080019.d

Injection Date: 08-Feb-2023 18:42:13

Instrument ID: CHHPLC\_X3

Lims ID: IC INT 1

Client ID:

Operator ID: JZ/MAR

ALS Bottle#:

19

Worklist Smp#: 19

Injection Vol: 100.0 ul

Dil. Factor:

1.0000

Method: 8330\_X3

Limit Group:

GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm)

Detector


LC DAD1B, 254 nm

**4 HMX, CAS: 2691-41-0**

Signal: 1


RT: 6.55  
Area: 1188  
Amount: 0.012270  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.55  
Area: 962  
Amount: 0.010286  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 08-Feb-2023 19:09:00

Audit Action: Manually Integrated

Audit Reason: Baseline Smoothing

## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230208-118465.b\02080019.d

Injection Date: 08-Feb-2023 18:42:13

Instrument ID: CHHPLC\_X3

Lims ID: IC INT 1

Client ID:

Operator ID: JZ/MAR

ALS Bottle#:

19

Worklist Smp#: 19

Injection Vol: 100.0 ul

Dil. Factor:

1.0000

Method: 8330\_X3

Limit Group:

GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm)

Detector


LC DAD1B, 254 nm

**8 RDX, CAS: 121-82-4**

Signal: 1


RT: 7.56  
Area: 1576  
Amount: 0.013341  
Amount Units: ug/mL

## Processing Integration Results



RT: 7.56  
Area: 1203  
Amount: 0.011309  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 08-Feb-2023 19:08:52

Audit Action: Manually Integrated

Audit Reason: Baseline Smoothing

FORM VI  
HPLC/IC BY EXTERNAL STANDARD - INITIAL CALIBRATION DATA  
RETENTION TIME SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1 Analy Batch No.: 603284  
SDG No.: \_\_\_\_\_  
Instrument ID: CHHPLC\_X3 GC Column: UltraCarb5u ID: 4.6(mm) Heated Purge: (Y/N) N  
Calibration Start Date: 02/24/2023 15:00 Calibration End Date: 02/24/2023 17:41 Calibration ID: 77558

Calibration Files

| LEVEL:  | LAB SAMPLE ID:   | LAB FILE ID: |
|---------|------------------|--------------|
| Level 1 | IC 280-603284/18 | 02240018.D   |
| Level 2 | IC 280-603284/17 | 02240017.D   |
| Level 3 | IC 280-603284/16 | 02240016.D   |
| Level 4 | IC 280-603284/15 | 02240015.D   |
| Level 5 | IC 280-603284/14 | 02240014.D   |
| Level 6 | IC 280-603284/13 | 02240013.D   |
| Level 7 | IC 280-603284/12 | 02240012.D   |
| Level 8 | IC 280-603284/11 | 02240011.D   |

| ANALYTE | LVL 1 | LVL 2 | LVL 3 | LVL 4 | LVL 5 | LVL 6 | LVL 7 | LVL 8 |  |  | RT WINDOW     | AVG RT |
|---------|-------|-------|-------|-------|-------|-------|-------|-------|--|--|---------------|--------|
| TNX     | 6.444 | 6.448 | 6.446 | 6.443 | 6.442 | 6.441 | 6.439 | 6.442 |  |  | 6.343 - 6.543 | 6.443  |
| DNX     | 6.771 | 6.768 | 6.766 | 6.770 | 6.769 | 6.767 | 6.759 | 6.762 |  |  | 6.670 - 6.870 | 6.767  |
| MNX     | 7.197 | 7.195 | 7.199 | 7.196 | 7.195 | 7.194 | 7.192 | 7.189 |  |  | 7.046 - 7.346 | 7.195  |

FORM VI  
HPLC/IC BY EXTERNAL STANDARD - INITIAL CALIBRATION DATA  
CURVE EVALUATION

Lab Name: Eurofins Denver Job No.: 280-176674-1 Analy Batch No.: 603284  
SDG No.: \_\_\_\_\_  
Instrument ID: CHHPLC\_X3 GC Column: UltraCarb5u ID: 4.6(mm) Heated Purge: (Y/N) N  
Calibration Start Date: 02/24/2023 15:00 Calibration End Date: 02/24/2023 17:41 Calibration ID: 77558

Calibration Files

|         |                  |              |
|---------|------------------|--------------|
| LEVEL:  | LAB SAMPLE ID:   | LAB FILE ID: |
| Level 1 | IC 280-603284/18 | 02240018.D   |
| Level 2 | IC 280-603284/17 | 02240017.D   |
| Level 3 | IC 280-603284/16 | 02240016.D   |
| Level 4 | IC 280-603284/15 | 02240015.D   |
| Level 5 | IC 280-603284/14 | 02240014.D   |
| Level 6 | IC 280-603284/13 | 02240013.D   |
| Level 7 | IC 280-603284/12 | 02240012.D   |
| Level 8 | IC 280-603284/11 | 02240011.D   |

| ANALYTE | CF               |                  |                  |                  | CURVE<br>TYPE | COEFFICIENT |                |    | # | MIN CF | %RSD<br>/RSE | # | MAX<br>%RSD<br>/RSE | R^2<br>OR COD | # | MIN R^2<br>OR COD |
|---------|------------------|------------------|------------------|------------------|---------------|-------------|----------------|----|---|--------|--------------|---|---------------------|---------------|---|-------------------|
|         | LVL 1<br>LVL 5   | LVL 2<br>LVL 6   | LVL 3<br>LVL 7   | LVL 4<br>LVL 8   |               | B           | M1             | M2 |   |        |              |   |                     |               |   |                   |
| TNX     | 187089<br>192298 | 197747<br>202669 | 198325<br>203463 | 196782<br>203127 | Ave           |             | 197687.33<br>1 |    |   |        | 2.9          |   | 20.0                |               |   |                   |
| DNX     | 139760<br>138674 | 146613<br>146977 | 145554<br>146290 | 145335<br>146496 | Ave           |             | 144462.48<br>6 |    |   |        | 2.3          |   | 20.0                |               |   |                   |
| MNX     | 127763<br>125139 | 130334<br>133592 | 131602<br>134086 | 134629<br>132227 | Ave           |             | 131171.56<br>5 |    |   |        | 2.5          |   | 20.0                |               |   |                   |

Note: The M1 coefficient is the same as Ave CF for an Ave curve type. RSD is calculated for Ave curve types. RSE is used for all other types.

FORM VI  
HPLC/IC BY EXTERNAL STANDARD - INITIAL CALIBRATION DATA  
RESPONSE AND CONCENTRATION

Lab Name: Eurofins Denver Job No.: 280-176674-1 Analy Batch No.: 603284

SDG No.: \_\_\_\_\_

Instrument ID: CHHPLC\_X3 GC Column: UltraCarb5u ID: 4.6(mm) Heated Purge: (Y/N) N

Calibration Start Date: 02/24/2023 15:00 Calibration End Date: 02/24/2023 17:41 Calibration ID: 77558

Calibration Files

|         |                  |              |
|---------|------------------|--------------|
| LEVEL:  | LAB SAMPLE ID:   | LAB FILE ID: |
| Level 1 | IC 280-603284/18 | 02240018.D   |
| Level 2 | IC 280-603284/17 | 02240017.D   |
| Level 3 | IC 280-603284/16 | 02240016.D   |
| Level 4 | IC 280-603284/15 | 02240015.D   |
| Level 5 | IC 280-603284/14 | 02240014.D   |
| Level 6 | IC 280-603284/13 | 02240013.D   |
| Level 7 | IC 280-603284/12 | 02240012.D   |
| Level 8 | IC 280-603284/11 | 02240011.D   |

| ANALYTE | CURVE TYPE | RESPONSE       |                |                |       |       | CONCENTRATION (UG/ML) |                |                |       |       |
|---------|------------|----------------|----------------|----------------|-------|-------|-----------------------|----------------|----------------|-------|-------|
|         |            | LVL 1<br>LVL 6 | LVL 2<br>LVL 7 | LVL 3<br>LVL 8 | LVL 4 | LVL 5 | LVL 1<br>LVL 6        | LVL 2<br>LVL 7 | LVL 3<br>LVL 8 | LVL 4 | LVL 5 |
| TNX     | Ave        | 3753           | 9917           | 19892          | 49343 | 77150 | 0.0201                | 0.0502         | 0.100          | 0.251 | 0.401 |
|         |            | 142294         | 204073         | 509340         |       |       | 0.702                 | 1.00           | 2.51           |       |       |
| DNX     | Ave        | 2798           | 7338           | 14570          | 36370 | 55525 | 0.0200                | 0.0501         | 0.100          | 0.250 | 0.400 |
|         |            | 102987         | 146436         | 366607         |       |       | 0.701                 | 1.00           | 2.50           |       |       |
| MNX     | Ave        | 2982           | 7605           | 15358          | 39278 | 58415 | 0.0233                | 0.0584         | 0.117          | 0.292 | 0.467 |
|         |            | 109131         | 156478         | 385772         |       |       | 0.817                 | 1.17           | 2.92           |       |       |

Curve Type Legend

Ave = Average

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240011.D  
 Lims ID: IC DMT 8  
 Client ID:  
 Sample Type: IC Calib Level: 8  
 Inject. Date: 24-Feb-2023 15:00:10 ALS Bottle#: 11 Worklist Smp#: 11  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: IC DMT 8  
 Operator ID: JZ/MAR Instrument ID: CHHPLC\_X3  
 Sublist: chrom-8330\_X3\*sub17  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 28-Feb-2023 13:17:25 Calib Date: 24-Feb-2023 17:41:20  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1602

First Level Reviewer: LV5D

Date: 24-Feb-2023 15:28:35

| Compound | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/mL | OnCol Amt<br>ug/mL | Flags |
|----------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 3 TNX    | 1   | 6.442        | 6.443            | -0.001           | 509340   | 2.51             | 2.58               | M     |
| 6 DNX    | 1   | 6.762        | 6.770            | -0.008           | 366607   | 2.50             | 2.54               | M     |
| 7 MNX    | 1   | 7.189        | 7.196            | -0.007           | 385772   | 2.92             | 2.94               | M     |

**QC Flag Legend**

Processing Flags

Review Flags

M - Manually Integrated

**Reagents:**

8330 DMT\_00013

Amount Added: 125.00

Units: uL

Report Date: 28-Feb-2023 13:17:26

Chrom Revision: 2.3 15-Feb-2023 20:44:50

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230224-118952.b\02240011.d

Injection Date: 24-Feb-2023 15:00:10

Instrument ID: CHHPLC\_X3

Operator ID: JZ/MAR

Lims ID: IC DMT 8

Worklist Smp#: 11

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 11

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230224-118952.b\02240011.d  
Injection Date: 24-Feb-2023 15:00:10 Instrument ID: CHHPLC\_X3  
Lims ID: IC DMT 8  
Client ID:  
Operator ID: JZ/MAR ALS Bottle#: 11 Worklist Smp#: 11  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**3 TNX, CAS: 13980-04-6**

Signal: 1


RT: 6.44  
Area: 548148  
Amount: 2.507500  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.44  
Area: 509340  
Amount: 2.576493  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 24-Feb-2023 16:41:21  
Audit Action: Manually Integrated

Audit Reason: Baseline



## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230224-118952.b\02240011.d  
Injection Date: 24-Feb-2023 15:00:10 Instrument ID: CHHPLC\_X3  
Lims ID: IC DMT 8  
Client ID:  
Operator ID: JZ/MAR ALS Bottle#: 11 Worklist Smp#: 11  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**6 DNX, CAS: 80251-29-2**

Signal: 1


RT: 6.76  
Area: 399339  
Amount: 2.502500  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.76  
Area: 366607  
Amount: 2.537731  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 24-Feb-2023 16:41:24  
Audit Action: Manually Integrated

Audit Reason: Baseline

## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230224-118952.b\02240011.d  
Injection Date: 24-Feb-2023 15:00:10 Instrument ID: CHHPLC\_X3  
Lims ID: IC DMT 8  
Client ID:  
Operator ID: JZ/MAR ALS Bottle#: 11 Worklist Smp#: 11  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**7 MNX, CAS: 5755-27-1**

Signal: 1


RT: 7.19  
Area: 413424  
Amount: 2.917500  
Amount Units: ug/mL

## Processing Integration Results



RT: 7.19  
Area: 385772  
Amount: 2.940973  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 24-Feb-2023 16:41:14  
Audit Action: Manually Integrated

Audit Reason: Baseline

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240012.D  
 Lims ID: IC DMT 7  
 Client ID:  
 Sample Type: IC Calib Level: 7  
 Inject. Date: 24-Feb-2023 15:23:08 ALS Bottle#: 12 Worklist Smp#: 12  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: IC DMT 7  
 Operator ID: JZ/MAR Instrument ID: CHHPLC\_X3  
 Sublist: chrom-8330\_X3\*sub17  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 28-Feb-2023 13:17:26 Calib Date: 24-Feb-2023 17:41:20  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1602

First Level Reviewer: LV5D

Date: 24-Feb-2023 16:07:37

| Compound | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/mL | OnCol Amt<br>ug/mL | Flags |
|----------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 3 TNX    | 1   | 6.439        | 6.443            | -0.004           | 204073   | 1.00             | 1.03               | M     |
| 6 DNX    | 1   | 6.759        | 6.770            | -0.011           | 146436   | 1.00             | 1.01               | M     |
| 7 MNX    | 1   | 7.192        | 7.196            | -0.004           | 156478   | 1.17             | 1.19               | M     |

**QC Flag Legend**

Processing Flags

Review Flags

M - Manually Integrated

**Reagents:**

8330 DMT\_00013

Amount Added: 50.00

Units: uL

Report Date: 28-Feb-2023 13:17:26

Chrom Revision: 2.3 15-Feb-2023 20:44:50

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230224-118952.b\02240012.d

Injection Date: 24-Feb-2023 15:23:08

Instrument ID: CHHPLC\_X3

Operator ID: JZ/MAR

Lims ID: IC DMT 7

Worklist Smp#: 12

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 12

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230224-118952.b\02240012.d  
Injection Date: 24-Feb-2023 15:23:08 Instrument ID: CHHPLC\_X3  
Lims ID: IC DMT 7  
Client ID:  
Operator ID: JZ/MAR ALS Bottle#: 12 Worklist Smp#: 12  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**3 TNX, CAS: 13980-04-6**

Signal: 1


RT: 6.44  
Area: 232230  
Amount: 1.069897  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.44  
Area: 204073  
Amount: 1.032302  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 24-Feb-2023 16:41:30  
Audit Action: Manually Integrated

Audit Reason: Baseline

## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230224-118952.b\02240012.d  
Injection Date: 24-Feb-2023 15:23:08 Instrument ID: CHHPLC\_X3  
Lims ID: IC DMT 7  
Client ID:  
Operator ID: JZ/MAR ALS Bottle#: 12 Worklist Smp#: 12  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**6 DNX, CAS: 80251-29-2**

Signal: 1


RT: 6.76  
Area: 164380  
Amount: 1.059050  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.76  
Area: 146436  
Amount: 1.013661  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 24-Feb-2023 16:41:34  
Audit Action: Manually Integrated

Audit Reason: Baseline

## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230224-118952.b\02240012.d  
Injection Date: 24-Feb-2023 15:23:08 Instrument ID: CHHPLC\_X3  
Lims ID: IC DMT 7  
Client ID:  
Operator ID: JZ/MAR ALS Bottle#: 12 Worklist Smp#: 12  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**7 MNX, CAS: 5755-27-1**

Signal: 1


RT: 7.19  
Area: 166389  
Amount: 1.221740  
Amount Units: ug/mL

## Processing Integration Results



RT: 7.19  
Area: 156478  
Amount: 1.192926  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 24-Feb-2023 16:41:38  
Audit Action: Manually Integrated

Audit Reason: Baseline

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240013.D  
 Lims ID: IC DMT 6  
 Client ID:  
 Sample Type: IC Calib Level: 6  
 Inject. Date: 24-Feb-2023 15:46:11 ALS Bottle#: 13 Worklist Smp#: 13  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: IC DMT 6  
 Operator ID: JZ/MAR Instrument ID: CHHPLC\_X3  
 Sublist: chrom-8330\_X3\*sub17  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 28-Feb-2023 13:17:27 Calib Date: 24-Feb-2023 17:41:20  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1602

First Level Reviewer: LV5D

Date: 24-Feb-2023 16:40:45

| Compound | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/mL | OnCol Amt<br>ug/mL | Flags |
|----------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 3 TNX    | 1   | 6.441        | 6.443            | -0.002           | 142294   | 0.7021           | 0.7198             | M     |
| 6 DNX    | 1   | 6.767        | 6.770            | -0.003           | 102987   | 0.7007           | 0.7129             | M     |
| 7 MNX    | 1   | 7.194        | 7.196            | -0.002           | 109131   | 0.8169           | 0.8320             | M     |

**QC Flag Legend**

Processing Flags

Review Flags

M - Manually Integrated

**Reagents:**

8330 DMT\_00013

Amount Added: 35.00

Units: uL



Report Date: 28-Feb-2023 13:17:27

Chrom Revision: 2.3 15-Feb-2023 20:44:50

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230224-118952.b\02240013.d

Injection Date: 24-Feb-2023 15:46:11

Instrument ID: CHHPLC\_X3

Operator ID: JZ/MAR

Lims ID: IC DMT 6

Worklist Smp#: 13

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 13

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230224-118952.b\02240013.d  
Injection Date: 24-Feb-2023 15:46:11 Instrument ID: CHHPLC\_X3  
Lims ID: IC DMT 6  
Client ID:  
Operator ID: JZ/MAR ALS Bottle#: 13 Worklist Smp#: 13  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**3 TNX, CAS: 13980-04-6**

Signal: 1


RT: 6.44  
Area: 149855  
Amount: 0.731638  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.44  
Area: 142294  
Amount: 0.719793  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 24-Feb-2023 16:41:49  
Audit Action: Manually Integrated

Audit Reason: Baseline

## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230224-118952.b\02240013.d  
Injection Date: 24-Feb-2023 15:46:11 Instrument ID: CHHPLC\_X3  
Lims ID: IC DMT 6  
Client ID:  
Operator ID: JZ/MAR ALS Bottle#: 13 Worklist Smp#: 13  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**6 DNX, CAS: 80251-29-2**

Signal: 1


RT: 6.77  
Area: 106819  
Amount: 0.729190  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.77  
Area: 102987  
Amount: 0.712898  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 24-Feb-2023 16:41:44  
Audit Action: Manually Integrated

Audit Reason: Baseline

## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230224-118952.b\02240013.d  
Injection Date: 24-Feb-2023 15:46:11 Instrument ID: CHHPLC\_X3  
Lims ID: IC DMT 6  
Client ID:  
Operator ID: JZ/MAR ALS Bottle#: 13 Worklist Smp#: 13  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**7 MNX, CAS: 5755-27-1**

Signal: 1


RT: 7.19  
Area: 110998  
Amount: 0.844687  
Amount Units: ug/mL

## Processing Integration Results



RT: 7.19  
Area: 109131  
Amount: 0.831971  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 24-Feb-2023 16:41:46  
Audit Action: Manually Integrated

Audit Reason: Baseline

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240014.D  
 Lims ID: IC DMT 5  
 Client ID:  
 Sample Type: IC Calib Level: 5  
 Inject. Date: 24-Feb-2023 16:09:13 ALS Bottle#: 14 Worklist Smp#: 14  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: IC DMT 5  
 Operator ID: JZ/MAR Instrument ID: CHHPLC\_X3  
 Sublist: chrom-8330\_X3\*sub17  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 28-Feb-2023 13:17:27 Calib Date: 24-Feb-2023 17:41:20  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1602

First Level Reviewer: LV5D

Date: 24-Feb-2023 16:41:09

| Compound | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/mL | OnCol Amt<br>ug/mL | Flags |
|----------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 3 TNX    | 1   | 6.442        | 6.443            | -0.001           | 77150    | 0.4012           | 0.3903             | M     |
| 6 DNX    | 1   | 6.769        | 6.770            | -0.001           | 55525    | 0.4004           | 0.3844             | M     |
| 7 MNX    | 1   | 7.195        | 7.196            | -0.001           | 58415    | 0.4668           | 0.4453             | M     |

**QC Flag Legend**

Processing Flags

Review Flags

M - Manually Integrated

**Reagents:**

8330 DMT\_00013

Amount Added: 20.00

Units: uL

Report Date: 28-Feb-2023 13:17:27

Chrom Revision: 2.3 15-Feb-2023 20:44:50

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230224-118952.b\02240014.d

Injection Date: 24-Feb-2023 16:09:13

Instrument ID: CHHPLC\_X3

Operator ID: JZ/MAR

Lims ID: IC DMT 5

Worklist Smp#: 14

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 14

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230224-118952.b\02240014.d

Injection Date: 24-Feb-2023 16:09:13

Instrument ID: CHHPLC\_X3

Lims ID: IC DMT 5

Client ID:

Operator ID: JZ/MAR

ALS Bottle#:

14

Worklist Smp#: 14

Injection Vol: 100.0 ul

Dil. Factor: 1.0000

Method: 8330\_X3

Limit Group: GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm)

Detector


LC DAD1B, 254 nm

**3 TNX, CAS: 13980-04-6**

Signal: 1


RT: 6.44  
Area: 81283  
Amount: 0.403064  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.44  
Area: 77150  
Amount: 0.390263  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 24-Feb-2023 16:41:55

Audit Action: Manually Integrated

Audit Reason: Baseline

## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230224-118952.b\02240014.d  
Injection Date: 24-Feb-2023 16:09:13 Instrument ID: CHHPLC\_X3  
Lims ID: IC DMT 5  
Client ID:  
Operator ID: JZ/MAR ALS Bottle#: 14 Worklist Smp#: 14  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**6 DNX, CAS: 80251-29-2**

Signal: 1


RT: 6.77  
Area: 57199  
Amount: 0.395371  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.77  
Area: 55525  
Amount: 0.384356  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 24-Feb-2023 16:41:57

Audit Action: Manually Integrated

Audit Reason: Baseline



## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230224-118952.b\02240014.d

Injection Date: 24-Feb-2023 16:09:13

Instrument ID: CHHPLC\_X3

Lims ID: IC DMT 5

Client ID:

Operator ID: JZ/MAR

ALS Bottle#:

14

Worklist Smp#: 14

Injection Vol: 100.0 ul

Dil. Factor: 1.0000

Method: 8330\_X3

Limit Group: GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm)

Detector


LC DAD1B, 254 nm

**7 MNX, CAS: 5755-27-1**

Signal: 1


RT: 7.20  
Area: 60440  
Amount: 0.462241  
Amount Units: ug/mL

## Processing Integration Results



RT: 7.20  
Area: 58415  
Amount: 0.445333  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 24-Feb-2023 16:40:59

Audit Action: Manually Integrated

Audit Reason: Baseline

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240015.D  
 Lims ID: IC DMT 4  
 Client ID:  
 Sample Type: IC Calib Level: 4  
 Inject. Date: 24-Feb-2023 16:32:12 ALS Bottle#: 15 Worklist Smp#: 15  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: IC DMT 4  
 Operator ID: JZ/MAR Instrument ID: CHHPLC\_X3  
 Sublist: chrom-8330\_X3\*sub17  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 28-Feb-2023 13:17:28 Calib Date: 24-Feb-2023 17:41:20  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1602

First Level Reviewer: LV5D

Date: 24-Feb-2023 17:15:42

| Compound | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/mL | OnCol Amt<br>ug/mL | Flags |
|----------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 3 TNX    | 1   | 6.443        | 6.443            | 0.000            | 49343    | 0.2508           | 0.2496             | M     |
| 6 DNX    | 1   | 6.770        | 6.770            | 0.000            | 36370    | 0.2503           | 0.2518             | M     |
| 7 MNX    | 1   | 7.196        | 7.196            | 0.000            | 39278    | 0.2918           | 0.2994             |       |

**QC Flag Legend**

Processing Flags

Review Flags

M - Manually Integrated

**Reagents:**

8330 DMT\_00013

Amount Added: 12.50

Units: uL

Report Date: 28-Feb-2023 13:17:28

Chrom Revision: 2.3 15-Feb-2023 20:44:50

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230224-118952.b\02240015.d

Injection Date: 24-Feb-2023 16:32:12

Instrument ID: CHHPLC\_X3

Operator ID: JZ/MAR

Lims ID: IC DMT 4

Worklist Smp#: 15

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 15

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230224-118952.b\02240015.d

Injection Date: 24-Feb-2023 16:32:12

Instrument ID: CHHPLC\_X3

Lims ID: IC DMT 4

Client ID:

Operator ID: JZ/MAR

ALS Bottle#:

15

Worklist Smp#: 15

Injection Vol: 100.0 ul

Dil. Factor: 1.0000

Method: 8330\_X3

Limit Group: GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm)

Detector


LC DAD1B, 254 nm

**3 TNX, CAS: 13980-04-6**

Signal: 1


RT: 6.44  
Area: 73619  
Amount: 0.336113  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.44  
Area: 49343  
Amount: 0.249601  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 24-Feb-2023 17:15:24

Audit Action: Manually Integrated

Audit Reason: Baseline

## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230224-118952.b\02240015.d  
Injection Date: 24-Feb-2023 16:32:12 Instrument ID: CHHPLC\_X3  
Lims ID: IC DMT 4  
Client ID:  
Operator ID: JZ/MAR ALS Bottle#: 15 Worklist Smp#: 15  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**6 DNX, CAS: 80251-29-2**

Signal: 1


RT: 6.77  
Area: 45251  
Amount: 0.297994  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.77  
Area: 36370  
Amount: 0.251761  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 24-Feb-2023 17:15:28  
Audit Action: Manually Integrated

Audit Reason: Baseline

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240016.D  
 Lims ID: IC DMT 3  
 Client ID:  
 Sample Type: IC Calib Level: 3  
 Inject. Date: 24-Feb-2023 16:55:16 ALS Bottle#: 16 Worklist Smp#: 16  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: IC DMT 3  
 Operator ID: JZ/MAR Instrument ID: CHHPLC\_X3  
 Sublist: chrom-8330\_X3\*sub17  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 28-Feb-2023 13:17:28 Calib Date: 24-Feb-2023 17:41:20  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1602

First Level Reviewer: LV5D

Date: 24-Feb-2023 17:30:15

| Compound | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/mL | OnCol Amt<br>ug/mL | Flags |
|----------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 3 TNX    | 1   | 6.446        | 6.443            | 0.003            | 19892    | 0.1003           | 0.1006             | M     |
| 6 DNX    | 1   | 6.766        | 6.770            | -0.004           | 14570    | 0.1001           | 0.1009             | M     |
| 7 MNX    | 1   | 7.199        | 7.196            | 0.003            | 15358    | 0.1167           | 0.1171             | M     |

**QC Flag Legend**

Processing Flags

Review Flags

M - Manually Integrated

**Reagents:**

8330 DMT\_00013

Amount Added: 5.00

Units: uL

Report Date: 28-Feb-2023 13:17:28

Chrom Revision: 2.3 15-Feb-2023 20:44:50

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230224-118952.b\02240016.d

Injection Date: 24-Feb-2023 16:55:16

Instrument ID: CHHPLC\_X3

Operator ID: JZ/MAR

Lims ID: IC DMT 3

Worklist Smp#: 16

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 16

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230224-118952.b\02240016.d  
Injection Date: 24-Feb-2023 16:55:16 Instrument ID: CHHPLC\_X3  
Lims ID: IC DMT 3  
Client ID:  
Operator ID: JZ/MAR ALS Bottle#: 16 Worklist Smp#: 16  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**3 TNX, CAS: 13980-04-6**

Signal: 1


RT: 6.45  
Area: 29096  
Amount: 0.135495  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.45  
Area: 19892  
Amount: 0.100624  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 24-Feb-2023 17:30:03  
Audit Action: Manually Integrated

Audit Reason: Baseline



## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230224-118952.b\02240016.d  
Injection Date: 24-Feb-2023 16:55:16 Instrument ID: CHHPLC\_X3  
Lims ID: IC DMT 3  
Client ID:  
Operator ID: JZ/MAR ALS Bottle#: 16 Worklist Smp#: 16  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**6 DNX, CAS: 80251-29-2**

Signal: 1


RT: 6.77  
Area: 15665  
Amount: 0.106775  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.77  
Area: 14570  
Amount: 0.100857  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 24-Feb-2023 17:30:07  
Audit Action: Manually Integrated

Audit Reason: Baseline

## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230224-118952.b\02240016.d

Injection Date: 24-Feb-2023 16:55:16

Instrument ID: CHHPLC\_X3

Lims ID: IC DMT 3

Client ID:

Operator ID: JZ/MAR

ALS Bottle#:

16

Worklist Smp#: 16

Injection Vol: 100.0 ul

Dil. Factor:

1.0000

Method: 8330\_X3

Limit Group:

GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm)

Detector


LC DAD1B, 254 nm

**7 MNX, CAS: 5755-27-1**

Signal: 1


RT: 7.20  
Area: 15598  
Amount: 0.117968  
Amount Units: ug/mL

## Processing Integration Results



RT: 7.20  
Area: 15358  
Amount: 0.117083  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 24-Feb-2023 17:30:11

Audit Action: Manually Integrated

Audit Reason: Baseline

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240017.D  
 Lims ID: IC DMT 2  
 Client ID:  
 Sample Type: IC Calib Level: 2  
 Inject. Date: 24-Feb-2023 17:18:19 ALS Bottle#: 17 Worklist Smp#: 17  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: IC DMT 2  
 Operator ID: JZ/MAR Instrument ID: CHHPLC\_X3  
 Sublist: chrom-8330\_X3\*sub17  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 28-Feb-2023 13:17:29 Calib Date: 24-Feb-2023 17:41:20  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1602

First Level Reviewer: LV5D

Date: 24-Feb-2023 17:51:38

| Compound | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/mL | OnCol Amt<br>ug/mL | Flags |
|----------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 3 TNX    | 1   | 6.448        | 6.443            | 0.005            | 9917     | 0.0502           | 0.0502             | M     |
| 6 DNX    | 1   | 6.768        | 6.770            | -0.002           | 7338     | 0.0501           | 0.0508             | M     |
| 7 MNX    | 1   | 7.195        | 7.196            | -0.001           | 7605     | 0.0584           | 0.0580             |       |

**QC Flag Legend**

Processing Flags

Review Flags

M - Manually Integrated

**Reagents:**

8330 DMT\_00013

Amount Added: 2.50

Units: uL

Report Date: 28-Feb-2023 13:17:29

Chrom Revision: 2.3 15-Feb-2023 20:44:50

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230224-118952.b\02240017.d

Injection Date: 24-Feb-2023 17:18:19

Instrument ID: CHHPLC\_X3

Operator ID: JZ/MAR

Lims ID: IC DMT 2

Worklist Smp#: 17

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 17

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230224-118952.b\02240017.d  
Injection Date: 24-Feb-2023 17:18:19 Instrument ID: CHHPLC\_X3  
Lims ID: IC DMT 2  
Client ID:  
Operator ID: JZ/MAR ALS Bottle#: 17 Worklist Smp#: 17  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**3 TNX, CAS: 13980-04-6**

Signal: 1


RT: 6.45  
Area: 21070  
Amount: 0.091223  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.45  
Area: 9917  
Amount: 0.050165  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 24-Feb-2023 17:51:28  
Audit Action: Manually Integrated

Audit Reason: Baseline

## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230224-118952.b\02240017.d

Injection Date: 24-Feb-2023 17:18:19

Instrument ID: CHHPLC\_X3

Lims ID: IC DMT 2

Client ID:

Operator ID: JZ/MAR

ALS Bottle#:

17

Worklist Smp#: 17

Injection Vol: 100.0 ul

Dil. Factor: 1.0000

Method: 8330\_X3

Limit Group: GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm)

Detector


LC DAD1B, 254 nm

**6 DNX, CAS: 80251-29-2**

Signal: 1


RT: 6.77  
Area: 8177  
Amount: 0.055426  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.77  
Area: 7338  
Amount: 0.050795  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 24-Feb-2023 17:51:32

Audit Action: Manually Integrated

Audit Reason: Baseline

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
 Lims ID: IC DMT 1  
 Client ID:  
 Sample Type: IC Calib Level: 1  
 Inject. Date: 24-Feb-2023 17:41:20 ALS Bottle#: 18 Worklist Smp#: 18  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: IC DMT 1  
 Operator ID: JZ/MAR Instrument ID: CHHPLC\_X3  
 Sublist: chrom-8330\_X3\*sub17  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 28-Feb-2023 13:17:30 Calib Date: 24-Feb-2023 17:41:20  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1602

First Level Reviewer: LV5D

Date: 24-Feb-2023 18:10:19

| Compound | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/mL | OnCol Amt<br>ug/mL | Flags |
|----------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 3 TNX    | 1   | 6.444        | 6.443            | 0.001            | 3753     | 0.0201           | 0.0190             | M     |
| 6 DNX    | 1   | 6.771        | 6.770            | 0.001            | 2798     | 0.0200           | 0.0194             | M     |
| 7 MNX    | 1   | 7.197        | 7.196            | 0.001            | 2982     | 0.0233           | 0.0227             |       |

**QC Flag Legend**

Processing Flags

Review Flags

M - Manually Integrated

**Reagents:**

8330 DMT\_00013

Amount Added: 1.00

Units: uL

Report Date: 28-Feb-2023 13:17:30

Chrom Revision: 2.3 15-Feb-2023 20:44:50

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230224-118952.b\02240018.d

Injection Date: 24-Feb-2023 17:41:20

Instrument ID: CHHPLC\_X3

Operator ID: JZ/MAR

Lims ID: IC DMT 1

Worklist Smp#: 18

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 18

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1





## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230224-118952.b\02240018.d  
Injection Date: 24-Feb-2023 17:41:20 Instrument ID: CHHPLC\_X3  
Lims ID: IC DMT 1  
Client ID:  
Operator ID: JZ/MAR ALS Bottle#: 18 Worklist Smp#: 18  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**3 TNX, CAS: 13980-04-6**

Signal: 1


RT: 6.44  
Area: 12281  
Amount: 0.048962  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.44  
Area: 3753  
Amount: 0.018985  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 24-Feb-2023 18:10:04  
Audit Action: Manually Integrated

Audit Reason: Baseline

## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230224-118952.b\02240018.d  
Injection Date: 24-Feb-2023 17:41:20 Instrument ID: CHHPLC\_X3  
Lims ID: IC DMT 1  
Client ID:  
Operator ID: JZ/MAR ALS Bottle#: 18 Worklist Smp#: 18  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**6 DNX, CAS: 80251-29-2**

Signal: 1


RT: 6.77  
Area: 3465  
Amount: 0.023313  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.77  
Area: 2798  
Amount: 0.019368  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 24-Feb-2023 18:10:11  
Audit Action: Manually Integrated

Audit Reason: Baseline

# Calibration

/ TNX

Curve Type: Average  
Weighting: Conc\_Sq  
Origin: Force  
Dependency: Response  
Calib Mode: ESTD  
Response Base: AREA  
RF Rounding: 0


## Curve Coefficients

Intercept: 0  
Slope: 1.977E+05

## Error Coefficients

Standard Error: 5810  
Relative Standard Error: 2.9  
Correlation Coefficient: 1.000  
Coefficient of Determination (Adjusted): 0.999

| ID | Level            | Concentration | Response | IS Amount | IS Response | RF            | Used |
|----|------------------|---------------|----------|-----------|-------------|---------------|------|
| 1  | IC 280-603284/18 | 0.02006       | 3753.0   |           |             | 187088.733799 | Y    |
| 2  | IC 280-603284/17 | 0.05015       | 9917.0   |           |             | 197746.759721 | Y    |
| 3  | IC 280-603284/16 | 0.1003        | 19892.0  |           |             | 198325.024925 | Y    |
| 4  | IC 280-603284/15 | 0.25075       | 49343.0  |           |             | 196781.655035 | Y    |
| 5  | IC 280-603284/14 | 0.4012        | 77150.0  |           |             | 192298.105683 | Y    |
| 6  | IC 280-603284/13 | 0.7021        | 142294.0 |           |             | 202669.135451 | Y    |
| 7  | IC 280-603284/12 | 1.003         | 204073.0 |           |             | 203462.612164 | Y    |
| 8  | IC 280-603284/11 | 2.5075        | 509340.0 |           |             | 203126.62014  | Y    |



## Calibration

/ DNX

Curve Type: Average  
 Weighting: Conc\_Sq  
 Origin: Force  
 Dependency: Response  
 Calib Mode: ESTD  
 Response Base: AREA  
 RF Rounding: 0


### Curve Coefficients

Intercept: 0  
 Slope: 1.445E+05

### Error Coefficients

Standard Error: 2320  
 Relative Standard Error: 2.3  
 Correlation Coefficient: 1.000  
 Coefficient of Determination (Adjusted): 0.999

| ID | Level            | Concentration | Response | IS Amount | IS Response | RF            | Used |
|----|------------------|---------------|----------|-----------|-------------|---------------|------|
| 1  | IC 280-603284/18 | 0.02002       | 2798.0   |           |             | 139760.23976  | Y    |
| 2  | IC 280-603284/17 | 0.05005       | 7338.0   |           |             | 146613.386613 | Y    |
| 3  | IC 280-603284/16 | 0.1001        | 14570.0  |           |             | 145554.445554 | Y    |
| 4  | IC 280-603284/15 | 0.25025       | 36370.0  |           |             | 145334.665335 | Y    |
| 5  | IC 280-603284/14 | 0.4004        | 55525.0  |           |             | 138673.826174 | Y    |
| 6  | IC 280-603284/13 | 0.7007        | 102987.0 |           |             | 146977.308406 | Y    |
| 7  | IC 280-603284/12 | 1.001         | 146436.0 |           |             | 146289.71029  | Y    |
| 8  | IC 280-603284/11 | 2.5025        | 366607.0 |           |             | 146496.303696 | Y    |



# Calibration

/ MNX

Curve Type: Average  
Weighting: Conc\_Sq  
Origin: Force  
Dependency: Response  
Calib Mode: ESTD  
Response Base: AREA  
RF Rounding: 0


## Curve Coefficients

Intercept: 0  
Slope: 1.312E+05

## Error Coefficients

Standard Error: 2200  
Relative Standard Error: 2.5  
Correlation Coefficient: 1.000  
Coefficient of Determination (Adjusted): 0.999

| ID | Level            | Concentration | Response | IS Amount | IS Response | RF            | Used |
|----|------------------|---------------|----------|-----------|-------------|---------------|------|
| 1  | IC 280-603284/18 | 0.02334       | 2982.0   |           |             | 127763.496144 | Y    |
| 2  | IC 280-603284/17 | 0.05835       | 7605.0   |           |             | 130334.190231 | Y    |
| 3  | IC 280-603284/16 | 0.1167        | 15358.0  |           |             | 131602.399314 | Y    |
| 4  | IC 280-603284/15 | 0.29175       | 39278.0  |           |             | 134628.963153 | Y    |
| 5  | IC 280-603284/14 | 0.4668        | 58415.0  |           |             | 125139.24593  | Y    |
| 6  | IC 280-603284/13 | 0.8169        | 109131.0 |           |             | 133591.626882 | Y    |
| 7  | IC 280-603284/12 | 1.167         | 156478.0 |           |             | 134085.689803 | Y    |
| 8  | IC 280-603284/11 | 2.9175        | 385772.0 |           |             | 132226.906598 | Y    |



FORM VI  
HPLC/IC BY EXTERNAL STANDARD - INITIAL CALIBRATION DATA  
RETENTION TIME SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1 Analy Batch No.: 610603

SDG No.: \_\_\_\_\_

Instrument ID: CHHPLC\_X5 GC Column: Luna-phenyl ID: 4.6(mm) Heated Purge: (Y/N) N

Calibration Start Date: 04/28/2023 18:24 Calibration End Date: 04/28/2023 23:03 Calibration ID: 79725

Calibration Files

|         |                  |              |
|---------|------------------|--------------|
| LEVEL:  | LAB SAMPLE ID:   | LAB FILE ID: |
| Level 1 | IC 280-610603/18 | 04280018.D   |
| Level 2 | IC 280-610603/17 | 04280017.D   |
| Level 3 | IC 280-610603/16 | 04280016.D   |
| Level 4 | IC 280-610603/15 | 04280015.D   |
| Level 5 | IC 280-610603/14 | 04280014.D   |
| Level 6 | IC 280-610603/13 | 04280013.D   |
| Level 7 | IC 280-610603/12 | 04280012.D   |
| Level 8 | IC 280-610603/11 | 04280011.D   |
| Level 9 | IC 280-610603/10 | 04280010.D   |

| ANALYTE                    | LVL 1  | LVL 2  | LVL 3  | LVL 4  | LVL 5  | LVL 6  | LVL 7  | LVL 8  | LVL 9  |  | RT WINDOW       | AVG RT |
|----------------------------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--|-----------------|--------|
| TNX                        | 5.116  | 5.136  | 5.116  | 5.116  | 5.115  | 5.108  | 5.096  | 5.095  | 5.014  |  | 4.965 - 5.265   | 5.101  |
| DNX                        | 5.963  | 5.983  | 5.969  | 5.963  | 5.955  | 5.954  | 5.943  | 5.942  | 5.854  |  | 5.805 - 6.105   | 5.947  |
| HMX                        | ++++   | 6.630  | 6.603  | 6.603  | 6.595  | 6.594  | 6.589  | 6.588  | 6.521  |  | 6.445 - 6.745   | 6.590  |
| MNX                        | 7.496  | 7.523  | 7.503  | 7.496  | 7.489  | 7.488  | 7.476  | 7.475  | 7.387  |  | 7.339 - 7.639   | 7.481  |
| Picric acid                | 8.063  | 8.116  | 8.076  | 8.056  | 8.029  | 8.008  | 7.956  | 7.935  | ++++   |  | 7.879 - 8.179   | 8.030  |
| RDX                        | 8.923  | 8.956  | 8.929  | 8.923  | 8.915  | 8.908  | 8.903  | 8.902  | 8.814  |  | 8.765 - 9.065   | 8.908  |
| Nitrobenzene               | 11.489 | 11.523 | 11.489 | 11.483 | 11.482 | 11.461 | 11.462 | 11.462 | 11.374 |  | 11.332 - 11.632 | 11.469 |
| 3,5-Dinitroaniline         | 14.376 | 14.410 | 14.369 | 14.363 | 14.355 | 14.334 | 14.329 | 14.355 | 14.367 |  | 14.205 - 14.505 | 14.362 |
| 1,3-Dinitrobenzene         | 14.636 | 14.663 | 14.629 | 14.623 | 14.622 | 14.601 | 14.596 | 14.615 | 14.527 |  | 14.472 - 14.772 | 14.612 |
| Nitroglycerin              | 15.196 | 15.216 | 15.176 | 15.176 | 15.182 | 15.154 | 15.149 | 15.175 | 15.121 |  | 15.032 - 15.332 | 15.172 |
| 2-Nitrotoluene             | 15.736 | 15.770 | 15.723 | 15.730 | 15.729 | 15.708 | 15.696 | 15.722 | 15.661 |  | 15.579 - 15.879 | 15.719 |
| 4-Nitrotoluene             | 15.969 | 15.996 | 15.956 | 15.956 | 15.955 | 15.934 | 15.922 | 15.948 | 15.834 |  | 15.805 - 16.105 | 15.941 |
| 4-Amino-2,6-dinitrotoluene | 16.529 | 16.556 | 16.523 | 16.523 | 16.515 | 16.501 | 16.482 | 16.528 | 16.441 |  | 16.365 - 16.665 | 16.511 |
| 3-Nitrotoluene             | 16.843 | 16.843 | 16.816 | 16.816 | 16.815 | 16.801 | 16.782 | 16.828 | 16.734 |  | 16.665 - 16.965 | 16.809 |
| 2-Amino-4,6-dinitrotoluene | 17.369 | 17.410 | 17.376 | 17.370 | 17.362 | 17.354 | 17.329 | 17.402 | 17.367 |  | 17.212 - 17.512 | 17.371 |
| 1,3,5-Trinitrobenzene      | 17.443 | 17.456 | 17.429 | 17.423 | 17.422 | 17.414 | 17.389 | 17.442 | 17.367 |  | 17.272 - 17.572 | 17.421 |
| 2,6-Dinitrotoluene         | 18.643 | 18.670 | 18.643 | 18.636 | 18.629 | 18.628 | 18.622 | 18.675 | 18.561 |  | 18.479 - 18.779 | 18.634 |
| 2,4-Dinitrotoluene         | 19.096 | 19.123 | 19.096 | 19.090 | 19.082 | 19.081 | 19.076 | 19.128 | 19.021 |  | 18.932 - 19.232 | 19.088 |
| Tetryl                     | 22.363 | 22.396 | 22.369 | 22.363 | 22.389 | 22.361 | 22.389 | 22.442 | 22.341 |  | 22.239 - 22.539 | 22.379 |
| 2,4,6-Trinitrotoluene      | 23.176 | 23.196 | 23.163 | 23.156 | 23.182 | 23.154 | 23.183 | 23.228 | 23.134 |  | 23.032 - 23.332 | 23.175 |
| PETN                       | 24.430 | 24.450 | 24.423 | 24.416 | 24.475 | 24.428 | 24.469 | 24.508 | 24.421 |  | 24.325 - 24.625 | 24.447 |
| 1,2-Dinitrobenzene         | 12.483 | 12.516 | 12.476 | 12.470 | 12.469 | 12.448 | 12.442 | 12.455 | 12.374 |  | 12.319 - 12.619 | 12.459 |

FORM VI  
HPLC/IC BY EXTERNAL STANDARD - INITIAL CALIBRATION DATA  
CURVE EVALUATION

Lab Name: Eurofins Denver Job No.: 280-176674-1 Analy Batch No.: 610603  
SDG No.: \_\_\_\_\_  
Instrument ID: CHHPLC\_X5 GC Column: Luna-phenyl ID: 4.6(mm) Heated Purge: (Y/N) N  
Calibration Start Date: 04/28/2023 18:24 Calibration End Date: 04/28/2023 23:03 Calibration ID: 79725

Calibration Files

|         |                  |              |
|---------|------------------|--------------|
| LEVEL:  | LAB SAMPLE ID:   | LAB FILE ID: |
| Level 1 | IC 280-610603/18 | 04280018.D   |
| Level 2 | IC 280-610603/17 | 04280017.D   |
| Level 3 | IC 280-610603/16 | 04280016.D   |
| Level 4 | IC 280-610603/15 | 04280015.D   |
| Level 5 | IC 280-610603/14 | 04280014.D   |
| Level 6 | IC 280-610603/13 | 04280013.D   |
| Level 7 | IC 280-610603/12 | 04280012.D   |
| Level 8 | IC 280-610603/11 | 04280011.D   |
| Level 9 | IC 280-610603/10 | 04280010.D   |

| ANALYTE            | CF                         |                  |                  |                  | CURVE TYPE | COEFFICIENT    |                |    | # | MIN CF | %RSD /RSE | # | MAX %RSD /RSE | R^2 OR COD | # | MIN R^2 OR COD |
|--------------------|----------------------------|------------------|------------------|------------------|------------|----------------|----------------|----|---|--------|-----------|---|---------------|------------|---|----------------|
|                    | LVL 1<br>LVL 5<br>LVL 9    | LVL 2<br>LVL 6   | LVL 3<br>LVL 7   | LVL 4<br>LVL 8   |            | B              | M1             | M2 |   |        |           |   |               |            |   |                |
| TNX                | 453838<br>343462<br>370331 | 407727<br>368303 | 404407<br>394444 | 416291<br>390186 | Ave        |                | 394332.07<br>5 |    |   |        | 8.1       |   | 20.0          |            |   |                |
| DNX                | 280120<br>247976<br>265367 | 286114<br>268282 | 294366<br>283929 | 286124<br>281374 | Ave        |                | 277072.29<br>0 |    |   |        | 5.1       |   | 20.0          |            |   |                |
| HMX                | ++++<br>184440<br>193632   | 222450<br>169983 | 215220<br>187240 | 195730<br>189027 | Lin2       | 831.97563<br>4 | 185425.33<br>8 |    |   |        |           |   |               | 0.9970     |   | 0.9900         |
| MNX                | 312425<br>228103<br>265422 | 263710<br>245259 | 268963<br>258928 | 261611<br>255982 | Ave        |                | 262267.07<br>3 |    |   |        | 8.6       |   | 20.0          |            |   |                |
| Picric acid        | 219700<br>160808<br>++++   | 172650<br>149730 | 179000<br>163827 | 170600<br>165319 | Lin2       | 546.11107<br>9 | 159807.27<br>0 |    |   |        |           |   |               | 0.9960     |   | 0.9900         |
| RDX                | 285200<br>214572<br>225981 | 243150<br>197783 | 232420<br>219114 | 228270<br>221588 | Lin2       | 672.50349<br>1 | 215639.31<br>8 |    |   |        |           |   |               | 0.9980     |   | 0.9900         |
| Nitrobenzene       | 426900<br>386724<br>401056 | 393550<br>354198 | 417660<br>395110 | 415850<br>398454 | Ave        |                | 398833.54<br>4 |    |   |        | 5.3       |   | 20.0          |            |   |                |
| 3,5-Dinitroaniline | 505300<br>452216<br>407228 | 461050<br>414810 | 483520<br>461404 | 469800<br>460771 | Lin2       | 607.26133<br>1 | 444569.70<br>7 |    |   |        |           |   |               | 0.9970     |   | 0.9900         |

Note: The M1 coefficient is the same as Ave CF for an Ave curve type. RSD is calculated for Ave curve types. RSE is used for all other types.

FORM VI  
HPLC/IC BY EXTERNAL STANDARD - INITIAL CALIBRATION DATA  
CURVE EVALUATION

Lab Name: Eurofins Denver Job No.: 280-176674-1 Analy Batch No.: 610603

SDG No.: \_\_\_\_\_

Instrument ID: CHHPLC\_X5 GC Column: Luna-phenyl ID: 4.6(mm) Heated Purge: (Y/N) N

Calibration Start Date: 04/28/2023 18:24 Calibration End Date: 04/28/2023 23:03 Calibration ID: 79725

| ANALYTE                    | CF                         |                  |                  |                  | CURVE<br>TYPE | COEFFICIENT |                |    | # | MIN CF | %RSD<br>/RSE | # | MAX<br>%RSD<br>/RSE | R^2<br>OR COD | # | MIN R^2<br>OR COD |
|----------------------------|----------------------------|------------------|------------------|------------------|---------------|-------------|----------------|----|---|--------|--------------|---|---------------------|---------------|---|-------------------|
|                            | LVL 1<br>LVL 5<br>LVL 9    | LVL 2<br>LVL 6   | LVL 3<br>LVL 7   | LVL 4<br>LVL 8   |               | B           | M1             | M2 |   |        |              |   |                     |               |   |                   |
| 1,3-Dinitrobenzene         | 694600<br>610632<br>687042 | 615350<br>554495 | 629180<br>616069 | 625960<br>630835 | Ave           |             | 629351.39<br>7 |    |   |        | 6.6          |   | 20.0                |               |   |                   |
| Nitroglycerin              | 146030<br>130878<br>134334 | 137045<br>120030 | 142168<br>132805 | 134661<br>134666 | Ave           |             | 134735.15<br>4 |    |   |        | 5.4          |   | 20.0                |               |   |                   |
| 2-Nitrotoluene             | 233600<br>242696<br>244772 | 264050<br>225665 | 259020<br>249076 | 252900<br>250892 | Ave           |             | 246963.36<br>8 |    |   |        | 4.8          |   | 20.0                |               |   |                   |
| 4-Nitrotoluene             | 220600<br>226220<br>241286 | 223300<br>203570 | 243580<br>227886 | 230020<br>230951 | Ave           |             | 227490.30<br>2 |    |   |        | 5.2          |   | 20.0                |               |   |                   |
| 4-Amino-2,6-dinitrotoluene | 327500<br>286712<br>289176 | 303200<br>260730 | 305300<br>288184 | 294720<br>292369 | Ave           |             | 294210.14<br>3 |    |   |        | 6.1          |   | 20.0                |               |   |                   |
| 3-Nitrotoluene             | 326900<br>284164<br>298841 | 291150<br>260845 | 310920<br>290920 | 295690<br>292277 | Ave           |             | 294634.08<br>9 |    |   |        | 6.1          |   | 20.0                |               |   |                   |
| 2-Amino-4,6-dinitrotoluene | 391200<br>434020<br>467007 | 412850<br>397298 | 463240<br>420357 | 439480<br>426004 | Ave           |             | 427939.49<br>4 |    |   |        | 6.1          |   | 20.0                |               |   |                   |
| 1,3,5-Trinitrobenzene      | 443800<br>392672<br>378606 | 389950<br>354488 | 441180<br>413963 | 419460<br>417113 | Ave           |             | 405692.41<br>7 |    |   |        | 7.2          |   | 20.0                |               |   |                   |
| 2,6-Dinitrotoluene         | 309300<br>283472<br>290190 | 285800<br>257870 | 317520<br>287934 | 302600<br>290390 | Ave           |             | 291675.18<br>7 |    |   |        | 5.9          |   | 20.0                |               |   |                   |
| 2,4-Dinitrotoluene         | 598700<br>572392<br>588866 | 612050<br>519843 | 608140<br>580639 | 594330<br>585209 | Ave           |             | 584463.07<br>5 |    |   |        | 4.7          |   | 20.0                |               |   |                   |
| Tetryl                     | 328200<br>273308<br>283850 | 317400<br>246538 | 284660<br>274789 | 275340<br>278494 | Ave           |             | 284730.94<br>1 |    |   |        | 8.6          |   | 20.0                |               |   |                   |
| 2,4,6-Trinitrotoluene      | 491900<br>418244<br>436729 | 496250<br>380470 | 443680<br>421144 | 431360<br>430103 | Ave           |             | 438875.56<br>5 |    |   |        | 8.2          |   | 20.0                |               |   |                   |

Note: The M1 coefficient is the same as Ave CF for an Ave curve type. RSD is calculated for Ave curve types. RSE is used for all other types.



FORM VI  
HPLC/IC BY EXTERNAL STANDARD - INITIAL CALIBRATION DATA  
CURVE EVALUATION

Lab Name: Eurofins Denver Job No.: 280-176674-1 Analy Batch No.: 610603

SDG No.: \_\_\_\_\_

Instrument ID: CHHPLC\_X5 GC Column: Luna-phenyl ID: 4.6 (mm) Heated Purge: (Y/N) N

Calibration Start Date: 04/28/2023 18:24 Calibration End Date: 04/28/2023 23:03 Calibration ID: 79725

| ANALYTE            | CF                         |                  |                  |                  | CURVE<br>TYPE | COEFFICIENT |                |    | # | MIN CF | %RSD<br>/RSE | # | MAX<br>%RSD<br>/RSE | R^2<br>OR COD | # | MIN R^2<br>OR COD |
|--------------------|----------------------------|------------------|------------------|------------------|---------------|-------------|----------------|----|---|--------|--------------|---|---------------------|---------------|---|-------------------|
|                    | LVL 1<br>LVL 5<br>LVL 9    | LVL 2<br>LVL 6   | LVL 3<br>LVL 7   | LVL 4<br>LVL 8   |               | B           | M1             | M2 |   |        |              |   |                     |               |   |                   |
| PETN               | 163440<br>140765<br>142965 | 144025<br>128148 | 145886<br>141788 | 146450<br>141973 | Ave           |             | 143937.77<br>0 |    |   |        | 6.3          |   | 20.0                |               |   |                   |
| 1,2-Dinitrobenzene | 316900<br>269384<br>276674 | 301900<br>245698 | 287580<br>273104 | 278930<br>276103 | Ave           |             | 280697.02<br>1 |    |   |        | 7.2          |   | 20.0                |               |   |                   |

Note: The M1 coefficient is the same as Ave CF for an Ave curve type. RSD is calculated for Ave curve types. RSE is used for all other types.

FORM VI  
HPLC/IC BY EXTERNAL STANDARD - INITIAL CALIBRATION DATA  
RESPONSE AND CONCENTRATION

Lab Name: Eurofins Denver Job No.: 280-176674-1 Analy Batch No.: 610603

SDG No.: \_\_\_\_\_

Instrument ID: CHHPLC\_X5 GC Column: Luna-phenyl ID: 4.6(mm) Heated Purge: (Y/N) N

Calibration Start Date: 04/28/2023 18:24 Calibration End Date: 04/28/2023 23:03 Calibration ID: 79725

Calibration Files

|         |                  |              |
|---------|------------------|--------------|
| LEVEL:  | LAB SAMPLE ID:   | LAB FILE ID: |
| Level 1 | IC 280-610603/18 | 04280018.D   |
| Level 2 | IC 280-610603/17 | 04280017.D   |
| Level 3 | IC 280-610603/16 | 04280016.D   |
| Level 4 | IC 280-610603/15 | 04280015.D   |
| Level 5 | IC 280-610603/14 | 04280014.D   |
| Level 6 | IC 280-610603/13 | 04280013.D   |
| Level 7 | IC 280-610603/12 | 04280012.D   |
| Level 8 | IC 280-610603/11 | 04280011.D   |
| Level 9 | IC 280-610603/10 | 04280010.D   |

| ANALYTE                    | CURVE TYPE | RESPONSE        |                 |                  |                   |        | CONCENTRATION (UG/ML) |                 |                |                |       |
|----------------------------|------------|-----------------|-----------------|------------------|-------------------|--------|-----------------------|-----------------|----------------|----------------|-------|
|                            |            | LVL 1<br>LVL 6  | LVL 2<br>LVL 7  | LVL 3<br>LVL 8   | LVL 4<br>LVL 9    | LVL 5  | LVL 1<br>LVL 6        | LVL 2<br>LVL 7  | LVL 3<br>LVL 8 | LVL 4<br>LVL 9 | LVL 5 |
| TNX                        | Ave        | 4552<br>147763  | 8179<br>276939  | 20281<br>391357  | 41754<br>928605   | 86123  | 0.0100<br>0.401       | 0.0201<br>0.702 | 0.0502<br>1.00 | 0.100<br>2.51  | 0.251 |
| DNX                        | Ave        | 2804<br>107420  | 5728<br>198949  | 14733<br>281655  | 28641<br>664081   | 62056  | 0.0100<br>0.400       | 0.0200<br>0.701 | 0.0501<br>1.00 | 0.100<br>2.50  | 0.250 |
| HMX                        | Lin2       | ++++<br>67993   | 4449<br>131068  | 10761<br>189027  | 19573<br>484081   | 46110  | ++++<br>0.400         | 0.0200<br>0.700 | 0.0500<br>1.00 | 0.100<br>2.50  | 0.250 |
| MNX                        | Ave        | 3646<br>114487  | 6155<br>211518  | 15694<br>298731  | 30530<br>774370   | 66549  | 0.0117<br>0.467       | 0.0233<br>0.817 | 0.0584<br>1.17 | 0.117<br>2.92  | 0.292 |
| Picric acid                | Lin2       | 2197<br>59892   | 3453<br>114679  | 8950<br>165319   | 17060<br>++++     | 40202  | 0.0100<br>0.400       | 0.0200<br>0.700 | 0.0500<br>1.00 | 0.100<br>++++  | 0.250 |
| RDX                        | Lin2       | 2852<br>79113   | 4863<br>153380  | 11621<br>221588  | 22827<br>564953   | 53643  | 0.0100<br>0.400       | 0.0200<br>0.700 | 0.0500<br>1.00 | 0.100<br>2.50  | 0.250 |
| Nitrobenzene               | Ave        | 4269<br>141679  | 7871<br>276577  | 20883<br>398454  | 41585<br>1002641  | 96681  | 0.0100<br>0.400       | 0.0200<br>0.700 | 0.0500<br>1.00 | 0.100<br>2.50  | 0.250 |
| 3,5-Dinitroaniline         | Lin2       | 5053<br>165924  | 9221<br>322983  | 24176<br>460771  | 46980<br>1018070  | 113054 | 0.0100<br>0.400       | 0.0200<br>0.700 | 0.0500<br>1.00 | 0.100<br>2.50  | 0.250 |
| 1,3-Dinitrobenzene         | Ave        | 6946<br>221798  | 12307<br>431248 | 31459<br>630835  | 62596<br>1717605  | 152658 | 0.0100<br>0.400       | 0.0200<br>0.700 | 0.0500<br>1.00 | 0.100<br>2.50  | 0.250 |
| Nitroglycerin              | Ave        | 14603<br>480119 | 27409<br>929634 | 71084<br>1346659 | 134661<br>3358347 | 327195 | 0.100<br>4.00         | 0.200<br>7.00   | 0.500<br>10.0  | 1.00<br>25.0   | 2.50  |
| 2-Nitrotoluene             | Ave        | 2336<br>90266   | 5281<br>174353  | 12951<br>250892  | 25290<br>611929   | 60674  | 0.0100<br>0.400       | 0.0200<br>0.700 | 0.0500<br>1.00 | 0.100<br>2.50  | 0.250 |
| 4-Nitrotoluene             | Ave        | 2206<br>81428   | 4466<br>159520  | 12179<br>230951  | 23002<br>603215   | 56555  | 0.0100<br>0.400       | 0.0200<br>0.700 | 0.0500<br>1.00 | 0.100<br>2.50  | 0.250 |
| 4-Amino-2,6-dinitrotoluene | Ave        | 3275            | 6064            | 15265            | 29472             | 71678  | 0.0100                | 0.0200          | 0.0500         | 0.100          | 0.250 |

FORM VI  
HPLC/IC BY EXTERNAL STANDARD - INITIAL CALIBRATION DATA  
RESPONSE AND CONCENTRATION

Lab Name: Eurofins Denver Job No.: 280-176674-1 Analy Batch No.: 610603

SDG No.: \_\_\_\_\_

Instrument ID: CHHPLC\_X5 GC Column: Luna-phenyl ID: 4.6(mm) Heated Purge: (Y/N) N

Calibration Start Date: 04/28/2023 18:24 Calibration End Date: 04/28/2023 23:03 Calibration ID: 79725

| ANALYTE                    | CURVE<br>TYPE | RESPONSE       |                |                |                |        | CONCENTRATION (UG/ML) |                |                |                |       |
|----------------------------|---------------|----------------|----------------|----------------|----------------|--------|-----------------------|----------------|----------------|----------------|-------|
|                            |               | LVL 1<br>LVL 6 | LVL 2<br>LVL 7 | LVL 3<br>LVL 8 | LVL 4<br>LVL 9 | LVL 5  | LVL 1<br>LVL 6        | LVL 2<br>LVL 7 | LVL 3<br>LVL 8 | LVL 4<br>LVL 9 | LVL 5 |
|                            |               | 104292         | 201729         | 292369         | 722940         |        | 0.400                 | 0.700          | 1.00           | 2.50           |       |
| 3-Nitrotoluene             | Ave           | 3269           | 5823           | 15546          | 29569          | 71041  | 0.0100                | 0.0200         | 0.0500         | 0.100          | 0.250 |
|                            |               | 104338         | 203644         | 292277         | 747102         |        | 0.400                 | 0.700          | 1.00           | 2.50           |       |
| 2-Amino-4,6-dinitrotoluene | Ave           | 3912           | 8257           | 23162          | 43948          | 108505 | 0.0100                | 0.0200         | 0.0500         | 0.100          | 0.250 |
|                            |               | 158919         | 294250         | 426004         | 1167517        |        | 0.400                 | 0.700          | 1.00           | 2.50           |       |
| 1,3,5-Trinitrobenzene      | Ave           | 4438           | 7799           | 22059          | 41946          | 98168  | 0.0100                | 0.0200         | 0.0500         | 0.100          | 0.250 |
|                            |               | 141795         | 289774         | 417113         | 946516         |        | 0.400                 | 0.700          | 1.00           | 2.50           |       |
| 2,6-Dinitrotoluene         | Ave           | 3093           | 5716           | 15876          | 30260          | 70868  | 0.0100                | 0.0200         | 0.0500         | 0.100          | 0.250 |
|                            |               | 103148         | 201554         | 290390         | 725476         |        | 0.400                 | 0.700          | 1.00           | 2.50           |       |
| 2,4-Dinitrotoluene         | Ave           | 5987           | 12241          | 30407          | 59433          | 143098 | 0.0100                | 0.0200         | 0.0500         | 0.100          | 0.250 |
|                            |               | 207937         | 406447         | 585209         | 1472164        |        | 0.400                 | 0.700          | 1.00           | 2.50           |       |
| Tetryl                     | Ave           | 3282           | 6348           | 14233          | 27534          | 68327  | 0.0100                | 0.0200         | 0.0500         | 0.100          | 0.250 |
|                            |               | 98615          | 192352         | 278494         | 709626         |        | 0.400                 | 0.700          | 1.00           | 2.50           |       |
| 2,4,6-Trinitrotoluene      | Ave           | 4919           | 9925           | 22184          | 43136          | 104561 | 0.0100                | 0.0200         | 0.0500         | 0.100          | 0.250 |
|                            |               | 152188         | 294801         | 430103         | 1091822        |        | 0.400                 | 0.700          | 1.00           | 2.50           |       |
| PETN                       | Ave           | 16344          | 28805          | 72943          | 146450         | 351913 | 0.100                 | 0.200          | 0.500          | 1.00           | 2.50  |
|                            |               | 512590         | 992519         | 1419730        | 3574120        |        | 4.00                  | 7.00           | 10.0           | 25.0           |       |
| 1,2-Dinitrobenzene         | Ave           | 3169           | 6038           | 14379          | 27893          | 67346  | 0.0100                | 0.0200         | 0.0500         | 0.100          | 0.250 |
|                            |               | 98279          | 191173         | 276103         | 691686         |        | 0.400                 | 0.700          | 1.00           | 2.50           |       |

Curve Type Legend

Ave = Average  
Lin2 = Linear 1/conc^2

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280010.D  
 Lims ID: IC INT 9  
 Client ID:  
 Sample Type: IC Calib Level: 9  
 Inject. Date: 28-Apr-2023 18:24:11 ALS Bottle#: 10 Worklist Smp#: 10  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: IC INT 9  
 Operator ID: JZ/JG Instrument ID: CHHPLC\_X5  
 Sublist: chrom-8330\_X5\_Luna\*sub6  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\8330\_X5\_Luna.m  
 Limit Group: GCSV - 8330  
 Last Update: 29-Apr-2023 13:04:16 Calib Date: 29-Apr-2023 04:18:01  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280027.D  
 Column 1 : Luna-Phenyl hexyl ( 4.60 mm) Det: LC mwd1A, 254 nm  
 Process Host: CTX1621

First Level Reviewer: LV5D

Date: 29-Apr-2023 13:01:46

| Compound                      | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 3 TNX                         | 1   | 5.014        | 5.115            | -0.101           | 928605   | 2.51             | 2.35               |       |
| 4 DNx                         | 1   | 5.854        | 5.955            | -0.101           | 664081   | 2.50             | 2.40               |       |
| 5 HMX                         | 1   | 6.521        | 6.595            | -0.074           | 484081   | 2.50             | 2.61               |       |
| 6 MNX                         | 1   | 7.387        | 7.489            | -0.102           | 774370   | 2.92             | 2.95               |       |
| 7 2,4,6-Trinitrophenol        | 1   | 7.867        | 8.029            | -0.162           | 285146   | 2.50             | 1.78               |       |
| 8 RDX                         | 1   | 8.814        | 8.915            | -0.101           | 564953   | 2.50             | 2.62               |       |
| 9 Nitrobenzene                | 1   | 11.374       | 11.482           | -0.108           | 1002641  | 2.50             | 2.51               |       |
| \$ 10 1,2-Dinitrobenzene      | 1   | 12.374       | 12.469           | -0.095           | 691686   | 2.50             | 2.46               |       |
| 11 3,5-Dinitroaniline         | 1   | 14.367       | 14.355           | 0.012            | 1018070  | 2.50             | 2.29               | Ma    |
| 12 1,3-Dinitrobenzene         | 1   | 14.527       | 14.622           | -0.095           | 1717605  | 2.50             | 2.73               | M     |
| 13 Nitroglycerin              | 2   | 15.121       | 15.182           | -0.061           | 3358347  | 25.0             | 24.9               |       |
| 14 o-Nitrotoluene             | 1   | 15.661       | 15.729           | -0.068           | 611929   | 2.50             | 2.48               | M     |
| 16 p-Nitrotoluene             | 1   | 15.834       | 15.955           | -0.121           | 603215   | 2.50             | 2.65               | M     |
| 17 4-Amino-2,6-dinitrotoluene | 1   | 16.441       | 16.515           | -0.074           | 722940   | 2.50             | 2.46               | M     |
| 18 m-Nitrotoluene             | 1   | 16.734       | 16.815           | -0.081           | 747102   | 2.50             | 2.54               | M     |
| 19 2-Amino-4,6-dinitrotoluene | 1   | 17.367       | 17.362           | 0.005            | 1167517  | 2.50             | 2.73               | M     |
| 20 1,3,5-Trinitrobenzene      | 1   | 17.367       | 17.422           | -0.055           | 946516   | 2.50             | 2.33               | M     |
| 21 2,6-Dinitrotoluene         | 1   | 18.561       | 18.629           | -0.068           | 725476   | 2.50             | 2.49               | M     |
| 22 2,4-Dinitrotoluene         | 1   | 19.021       | 19.082           | -0.061           | 1472164  | 2.50             | 2.52               | M     |
| 23 Tetryl                     | 1   | 22.341       | 22.389           | -0.048           | 709626   | 2.50             | 2.49               |       |
| 24 2,4,6-Trinitrotoluene      | 1   | 23.134       | 23.182           | -0.048           | 1091822  | 2.50             | 2.49               |       |
| 25 PETN                       | 2   | 24.421       | 24.475           | -0.054           | 3574120  | 25.0             | 24.8               |       |

## QC Flag Legend

Processing Flags

Review Flags

M - Manually Integrated

a - User Assigned ID

Reagents:

8330IntermStk\_00076

Amount Added: 250.00

Units: uL

8330 DMT\_00013

Amount Added: 125.00

Units: uL

Report Date: 29-Apr-2023 13:04:16

Chrom Revision: 2.3 29-Mar-2023 18:39:10

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280010.D

Injection Date: 28-Apr-2023 18:24:11

Instrument ID: CHHPLC\_X5

Operator ID: JZ/JG

Lims ID: IC INT 9

Worklist Smp#: 10

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 10

Method: 8330\_X5\_Luna

Limit Group: GCSV - 8330


Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280010.D  
Injection Date: 28-Apr-2023 18:24:11 Instrument ID: CHHPLC\_X5  
Lims ID: IC INT 9  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 10 Worklist Smp#: 10  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X5\_Luna Limit Group: GCSV - 8330  
Column: Luna-Phenyl hexyl ( 4.60 mm) Detector LC mwd1A, 254 nm

**11 3,5-Dinitroaniline, CAS: 618-87-1**

Signal: 1


RT: 14.37  
Area: 1013469  
Amount: 2.724745  
Amount Units: ug/ml

## Processing Integration Results



RT: 14.37  
Area: 1018070  
Amount: 2.288646  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LV5D, 28-Apr-2023 19:01:40  
Audit Action: Assigned New Baseline

Audit Reason: Baseline Smoothing

## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280010.D  
Injection Date: 28-Apr-2023 18:24:11 Instrument ID: CHHPLC\_X5  
Lims ID: IC INT 9  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 10 Worklist Smp#: 10  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X5\_Luna Limit Group: GCSV - 8330  
Column: Luna-Phenyl hexyl ( 4.60 mm) Detector LC mwd1A, 254 nm

**12 1,3-Dinitrobenzene, CAS: 99-65-0**

Signal: 1


RT: 14.53  
Area: 2715993  
Amount: 3.616571  
Amount Units: ug/ml

## Processing Integration Results



RT: 14.53  
Area: 1717605  
Amount: 2.729167  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LV5D, 28-Apr-2023 19:01:40  
Audit Action: Assigned New Baseline

Audit Reason: Baseline Smoothing



## Eurofins Denver

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280010.D

Injection Date: 28-Apr-2023 18:24:11

Instrument ID: CHHPLC\_X5

Lims ID: IC INT 9

Client ID:

Operator ID: JZ/JG

ALS Bottle#: 10 Worklist Smp#: 10

Injection Vol: 100.0 ul

Dil. Factor: 1.0000

Method: 8330\_X5\_Luna

Limit Group: GCSV - 8330

Column: Luna-Phenyl hexyl ( 4.60 mm)


Detector LC mwd1A, 254 nm

**14 o-Nitrotoluene, CAS: 88-72-2**

Signal: 1


RT: 15.66  
Area: 811716  
Amount: 3.605542  
Amount Units: ug/ml

## Processing Integration Results



RT: 15.66  
Area: 611929  
Amount: 2.477813  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LV5D, 28-Apr-2023 19:01:40

Audit Action: Assigned New Baseline

Audit Reason: Baseline Smoothing

## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280010.D  
Injection Date: 28-Apr-2023 18:24:11 Instrument ID: CHHPLC\_X5  
Lims ID: IC INT 9  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 10 Worklist Smp#: 10  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X5\_Luna Limit Group: GCSV - 8330  
Column: Luna-Phenyl hexyl ( 4.60 mm) Detector LC mwd1A, 254 nm

**16 p-Nitrotoluene, CAS: 99-99-0**

Signal: 1


RT: 15.83  
Area: 577007  
Amount: 2.272934  
Amount Units: ug/ml

## Processing Integration Results



RT: 15.83  
Area: 603215  
Amount: 2.651608  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LV5D, 28-Apr-2023 19:01:40  
Audit Action: Assigned New Baseline

Audit Reason: Baseline Smoothing

## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280010.D  
Injection Date: 28-Apr-2023 18:24:11 Instrument ID: CHHPLC\_X5  
Lims ID: IC INT 9  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 10 Worklist Smp#: 10  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X5\_Luna Limit Group: GCSV - 8330  
Column: Luna-Phenyl hexyl ( 4.60 mm) Detector LC mwd1A, 254 nm

**17 4-Amino-2,6-dinitrotoluene, CAS: 19406-51-0**

Signal: 1


RT: 16.44  
Area: 691163  
Amount: 2.645631  
Amount Units: ug/ml

## Processing Integration Results



RT: 16.44  
Area: 722940  
Amount: 2.457223  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LV5D, 28-Apr-2023 19:01:40  
Audit Action: Assigned New Baseline

Audit Reason: Baseline Smoothing

## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280010.D  
Injection Date: 28-Apr-2023 18:24:11 Instrument ID: CHHPLC\_X5  
Lims ID: IC INT 9  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 10 Worklist Smp#: 10  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X5\_Luna Limit Group: GCSV - 8330  
Column: Luna-Phenyl hexyl ( 4.60 mm) Detector LC mwd1A, 254 nm

**18 m-Nitrotoluene, CAS: 99-08-1**

Signal: 1


RT: 16.73  
Area: 719908  
Amount: 2.277668  
Amount Units: ug/ml

## Processing Integration Results



RT: 16.73  
Area: 747102  
Amount: 2.535694  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LV5D, 28-Apr-2023 19:01:40  
Audit Action: Assigned New Baseline

Audit Reason: Baseline Smoothing

## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280010.D  
Injection Date: 28-Apr-2023 18:24:11 Instrument ID: CHHPLC\_X5  
Lims ID: IC INT 9  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 10 Worklist Smp#: 10  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X5\_Luna Limit Group: GCSV - 8330  
Column: Luna-Phenyl hexyl ( 4.60 mm) Detector LC mwd1A, 254 nm

**20 1,3,5-Trinitrobenzene, CAS: 99-35-4**

Signal: 1


RT: 17.37  
Area: 2081660  
Amount: 3.645357  
Amount Units: ug/ml

## Processing Integration Results



RT: 17.37  
Area: 946516  
Amount: 2.333088  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LV5D, 28-Apr-2023 19:01:40  
Audit Action: Assigned New Baseline

Audit Reason: Baseline Smoothing

## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280010.D  
Injection Date: 28-Apr-2023 18:24:11 Instrument ID: CHHPLC\_X5  
Lims ID: IC INT 9  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 10 Worklist Smp#: 10  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X5\_Luna Limit Group: GCSV - 8330  
Column: Luna-Phenyl hexyl ( 4.60 mm) Detector LC mwd1A, 254 nm

**19 2-Amino-4,6-dinitrotoluene, CAS: 35572-78-2**

Signal: 1


RT: 17.37  
Area: 1147839  
Amount: 3.846971  
Amount Units: ug/ml

## Processing Integration Results



RT: 17.37  
Area: 1167517  
Amount: 2.728229  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LV5D, 28-Apr-2023 19:01:40  
Audit Action: Assigned New Baseline

Audit Reason: Baseline Smoothing

## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280010.D  
Injection Date: 28-Apr-2023 18:24:11 Instrument ID: CHHPLC\_X5  
Lims ID: IC INT 9  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 10 Worklist Smp#: 10  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X5\_Luna Limit Group: GCSV - 8330  
Column: Luna-Phenyl hexyl ( 4.60 mm) Detector: LC mwd1A, 254 nm

**21 2,6-Dinitrotoluene, CAS: 606-20-2**

Signal: 1


RT: 18.56  
Area: 716752  
Amount: 2.903314  
Amount Units: ug/ml

## Processing Integration Results



RT: 18.56  
Area: 725476  
Amount: 2.487274  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LV5D, 28-Apr-2023 19:01:40  
Audit Action: Assigned New Baseline

Audit Reason: Baseline Smoothing

## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280010.D  
Injection Date: 28-Apr-2023 18:24:11 Instrument ID: CHHPLC\_X5  
Lims ID: IC INT 9  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 10 Worklist Smp#: 10  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X5\_Luna Limit Group: GCSV - 8330  
Column: Luna-Phenyl hexyl ( 4.60 mm) Detector LC mwd1A, 254 nm

**22 2,4-Dinitrotoluene, CAS: 121-14-2**

Signal: 1


RT: 19.02  
Area: 1458826  
Amount: 2.601930  
Amount Units: ug/ml

## Processing Integration Results



RT: 19.02  
Area: 1472164  
Amount: 2.518831  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LV5D, 28-Apr-2023 19:01:40  
Audit Action: Assigned New Baseline

Audit Reason: Baseline Smoothing



Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280011.D  
 Lims ID: IC INT 8  
 Client ID:  
 Sample Type: IC Calib Level: 8  
 Inject. Date: 28-Apr-2023 18:59:07 ALS Bottle#: 11 Worklist Smp#: 11  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: IC INT 8  
 Operator ID: JZ/JG Instrument ID: CHHPLC\_X5  
 Sublist: chrom-8330\_X5\_Luna\*sub6  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\8330\_X5\_Luna.m  
 Limit Group: GCSV - 8330  
 Last Update: 29-Apr-2023 13:04:17 Calib Date: 29-Apr-2023 04:18:01  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280027.D  
 Column 1 : Luna-Phenyl hexyl ( 4.60 mm) Det: LC mwd1A, 254 nm  
 Process Host: CTX1621

First Level Reviewer: LV5D

Date: 28-Apr-2023 20:43:46

| Compound                      | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 3 TNX                         | 1   | 5.095        | 5.115            | -0.020           | 391357   | 1.00             | 0.99               |       |
| 4 DNX                         | 1   | 5.942        | 5.955            | -0.013           | 281655   | 1.00             | 1.02               |       |
| 5 HMX                         | 1   | 6.588        | 6.595            | -0.007           | 189027   | 1.00             | 1.01               |       |
| 6 MNX                         | 1   | 7.475        | 7.489            | -0.014           | 298731   | 1.17             | 1.14               |       |
| 7 2,4,6-Trinitrophenol        | 1   | 7.935        | 8.029            | -0.094           | 165319   | 1.00             | 1.03               |       |
| 8 RDX                         | 1   | 8.902        | 8.915            | -0.013           | 221588   | 1.00             | 1.02               |       |
| 9 Nitrobenzene                | 1   | 11.462       | 11.482           | -0.020           | 398454   | 1.00             | 1.00               |       |
| \$ 10 1,2-Dinitrobenzene      | 1   | 12.455       | 12.469           | -0.014           | 276103   | 1.00             | 0.9836             |       |
| 11 3,5-Dinitroaniline         | 1   | 14.355       | 14.355           | 0.000            | 460771   | 1.00             | 1.04               |       |
| 12 1,3-Dinitrobenzene         | 1   | 14.615       | 14.622           | -0.007           | 630835   | 1.00             | 1.00               |       |
| 13 Nitroglycerin              | 2   | 15.175       | 15.182           | -0.007           | 1346659  | 10.0             | 10.0               |       |
| 14 o-Nitrotoluene             | 1   | 15.722       | 15.729           | -0.007           | 250892   | 1.00             | 1.02               |       |
| 16 p-Nitrotoluene             | 1   | 15.948       | 15.955           | -0.007           | 230951   | 1.00             | 1.02               |       |
| 17 4-Amino-2,6-dinitrotoluene | 1   | 16.528       | 16.515           | 0.013            | 292369   | 1.00             | 0.99               |       |
| 18 m-Nitrotoluene             | 1   | 16.828       | 16.815           | 0.013            | 292277   | 1.00             | 0.99               |       |
| 19 2-Amino-4,6-dinitrotoluene | 1   | 17.402       | 17.362           | 0.040            | 426004   | 1.00             | 1.00               | Ma    |
| 20 1,3,5-Trinitrobenzene      | 1   | 17.442       | 17.422           | 0.020            | 417113   | 1.00             | 1.03               | M     |
| 21 2,6-Dinitrotoluene         | 1   | 18.675       | 18.629           | 0.046            | 290390   | 1.00             | 1.00               |       |
| 22 2,4-Dinitrotoluene         | 1   | 19.128       | 19.082           | 0.046            | 585209   | 1.00             | 1.00               |       |
| 23 Tetryl                     | 1   | 22.442       | 22.389           | 0.053            | 278494   | 1.00             | 0.9781             |       |
| 24 2,4,6-Trinitrotoluene      | 1   | 23.228       | 23.182           | 0.046            | 430103   | 1.00             | 0.9800             |       |
| 25 PETN                       | 2   | 24.508       | 24.475           | 0.033            | 1419730  | 10.0             | 9.86               |       |

### QC Flag Legend

Processing Flags

Review Flags

M - Manually Integrated

a - User Assigned ID

Reagents:

8330IntermStk\_00076

Amount Added: 100.00

Units: uL

8330 DMT\_00013

Amount Added: 50.00

Units: uL

Report Date: 29-Apr-2023 13:04:17

Chrom Revision: 2.3 29-Mar-2023 18:39:10

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280011.D

Injection Date: 28-Apr-2023 18:59:07

Instrument ID: CHHPLC\_X5

Operator ID: JZ/JG

Lims ID: IC INT 8

Worklist Smp#: 11

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 11

Method: 8330\_X5\_Luna

Limit Group: GCSV - 8330


Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280011.D  
Injection Date: 28-Apr-2023 18:59:07 Instrument ID: CHHPLC\_X5  
Lims ID: IC INT 8  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 11 Worklist Smp#: 11  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X5\_Luna Limit Group: GCSV - 8330  
Column: Luna-Phenyl hexyl ( 4.60 mm) Detector LC mwd1A, 254 nm

**19 2-Amino-4,6-dinitrotoluene, CAS: 35572-78-2**

Signal: 1


Not Detected  
Expected RT: 17.36

## Processing Integration Results



RT: 17.40  
Area: 426004  
Amount: 0.995477  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LV5D, 28-Apr-2023 20:43:43

Audit Action: Manually Integrated/Assigned Compound ID Audit Reason: Baseline Smoothing

## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280011.D  
Injection Date: 28-Apr-2023 18:59:07 Instrument ID: CHHPLC\_X5  
Lims ID: IC INT 8  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 11 Worklist Smp#: 11  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X5\_Luna Limit Group: GCSV - 8330  
Column: Luna-Phenyl hexyl ( 4.60 mm) Detector LC mwd1A, 254 nm

**20 1,3,5-Trinitrobenzene, CAS: 99-35-4**

Signal: 1


RT: 17.40  
Area: 843108  
Amount: 1.537916  
Amount Units: ug/ml

## Processing Integration Results



RT: 17.44  
Area: 417113  
Amount: 1.028151  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LV5D, 28-Apr-2023 20:43:40  
Audit Action: Split an Integrated Peak

Audit Reason: Baseline Smoothing

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280012.D  
 Lims ID: IC INT 7  
 Client ID:  
 Sample Type: IC Calib Level: 7  
 Inject. Date: 28-Apr-2023 19:34:01 ALS Bottle#: 12 Worklist Smp#: 12  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: IC INT 7  
 Operator ID: JZ/JG Instrument ID: CHHPLC\_X5  
 Sublist: chrom-8330\_X5\_Luna\*sub6  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\8330\_X5\_Luna.m  
 Limit Group: GCSV - 8330  
 Last Update: 29-Apr-2023 13:04:17 Calib Date: 29-Apr-2023 04:18:01  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280027.D  
 Column 1 : Luna-Phenyl hexyl ( 4.60 mm) Det: LC mwd1A, 254 nm  
 Process Host: CTX1621

First Level Reviewer: LV5D

Date: 28-Apr-2023 20:44:11

| Compound                      | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 3 TNX                         | 1   | 5.096        | 5.115            | -0.019           | 276939   | 0.7021           | 0.7023             |       |
| 4 DNx                         | 1   | 5.943        | 5.955            | -0.012           | 198949   | 0.7007           | 0.7180             |       |
| 5 HMX                         | 1   | 6.589        | 6.595            | -0.006           | 131068   | 0.7000           | 0.7024             |       |
| 6 MNX                         | 1   | 7.476        | 7.489            | -0.013           | 211518   | 0.8169           | 0.8065             |       |
| 7 2,4,6-Trinitrophenol        | 1   | 7.956        | 8.029            | -0.073           | 114679   | 0.7000           | 0.7142             |       |
| 8 RDX                         | 1   | 8.903        | 8.915            | -0.012           | 153380   | 0.7000           | 0.7082             |       |
| 9 Nitrobenzene                | 1   | 11.462       | 11.482           | -0.020           | 276577   | 0.7000           | 0.6935             |       |
| \$ 10 1,2-Dinitrobenzene      | 1   | 12.442       | 12.469           | -0.027           | 191173   | 0.7000           | 0.6811             |       |
| 11 3,5-Dinitroaniline         | 1   | 14.329       | 14.355           | -0.026           | 322983   | 0.7000           | 0.7251             |       |
| 12 1,3-Dinitrobenzene         | 1   | 14.596       | 14.622           | -0.026           | 431248   | 0.7000           | 0.6852             |       |
| 13 Nitroglycerin              | 2   | 15.149       | 15.182           | -0.033           | 929634   | 7.00             | 6.90               |       |
| 14 o-Nitrotoluene             | 1   | 15.696       | 15.729           | -0.033           | 174353   | 0.7000           | 0.7060             |       |
| 16 p-Nitrotoluene             | 1   | 15.922       | 15.955           | -0.033           | 159520   | 0.7000           | 0.7012             |       |
| 17 4-Amino-2,6-dinitrotoluene | 1   | 16.482       | 16.515           | -0.033           | 201729   | 0.7000           | 0.6857             |       |
| 18 m-Nitrotoluene             | 1   | 16.782       | 16.815           | -0.033           | 203644   | 0.7000           | 0.6912             |       |
| 19 2-Amino-4,6-dinitrotoluene | 1   | 17.329       | 17.362           | -0.033           | 294250   | 0.7000           | 0.6876             | M     |
| 20 1,3,5-Trinitrobenzene      | 1   | 17.389       | 17.422           | -0.033           | 289774   | 0.7000           | 0.7143             | Ma    |
| 21 2,6-Dinitrotoluene         | 1   | 18.622       | 18.629           | -0.007           | 201554   | 0.7000           | 0.6910             |       |
| 22 2,4-Dinitrotoluene         | 1   | 19.076       | 19.082           | -0.006           | 406447   | 0.7000           | 0.6954             |       |
| 23 Tetryl                     | 1   | 22.389       | 22.389           | 0.000            | 192352   | 0.7000           | 0.6756             |       |
| 24 2,4,6-Trinitrotoluene      | 1   | 23.183       | 23.182           | 0.001            | 294801   | 0.7000           | 0.6717             |       |
| 25 PETN                       | 2   | 24.469       | 24.475           | -0.006           | 992519   | 7.00             | 6.90               |       |

## QC Flag Legend

Processing Flags

Review Flags

M - Manually Integrated

a - User Assigned ID

**Reagents:**

8330IntermStk\_00076

Amount Added: 70.00

Units: uL

8330 DMT\_00013

Amount Added: 35.00

Units: uL

Report Date: 29-Apr-2023 13:04:17

Chrom Revision: 2.3 29-Mar-2023 18:39:10

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280012.D

Injection Date: 28-Apr-2023 19:34:01

Instrument ID: CHHPLC\_X5

Operator ID: JZ/JG

Lims ID: IC INT 7

Worklist Smp#: 12

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 12

Method: 8330\_X5\_Luna

Limit Group: GCSV - 8330


Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1





## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280012.D  
Injection Date: 28-Apr-2023 19:34:01 Instrument ID: CHHPLC\_X5  
Lims ID: IC INT 7  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 12 Worklist Smp#: 12  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X5\_Luna Limit Group: GCSV - 8330  
Column: Luna-Phenyl hexyl ( 4.60 mm) Detector LC mwd1A, 254 nm

**19 2-Amino-4,6-dinitrotoluene, CAS: 35572-78-2**

Signal: 1


RT: 17.33  
Area: 584005  
Amount: 1.540308  
Amount Units: ug/ml

## Processing Integration Results



RT: 17.33  
Area: 294250  
Amount: 0.687597  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LV5D, 28-Apr-2023 20:43:51  
Audit Action: Split an Integrated Peak

Audit Reason: Baseline Smoothing

## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280012.D  
Injection Date: 28-Apr-2023 19:34:01 Instrument ID: CHHPLC\_X5  
Lims ID: IC INT 7  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 12 Worklist Smp#: 12  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X5\_Luna Limit Group: GCSV - 8330  
Column: Luna-Phenyl hexyl ( 4.60 mm) Detector LC mwd1A, 254 nm

**20 1,3,5-Trinitrobenzene, CAS: 99-35-4**

Signal: 1


Not Detected  
Expected RT: 17.42

## Processing Integration Results



RT: 17.39  
Area: 289774  
Amount: 0.714270  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LV5D, 28-Apr-2023 20:43:53

Audit Action: Manually Integrated/Assigned Compound ID Audit Reason: Baseline Smoothing

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280013.D  
 Lims ID: IC INT 6  
 Client ID:  
 Sample Type: IC Calib Level: 6  
 Inject. Date: 28-Apr-2023 20:08:59 ALS Bottle#: 13 Worklist Smp#: 13  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: IC INT 6  
 Operator ID: JZ/JG Instrument ID: CHHPLC\_X5  
 Sublist: chrom-8330\_X5\_Luna\*sub6  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\8330\_X5\_Luna.m  
 Limit Group: GCSV - 8330  
 Last Update: 29-Apr-2023 13:04:18 Calib Date: 29-Apr-2023 04:18:01  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280027.D  
 Column 1 : Luna-Phenyl hexyl ( 4.60 mm) Det: LC mwd1A, 254 nm  
 Process Host: CTX1621

First Level Reviewer: LV5D

Date: 28-Apr-2023 20:52:44

| Compound                      | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 3 TNX                         | 1   | 5.108        | 5.115            | -0.007           | 147763   | 0.4012           | 0.3747             |       |
| 4 DNx                         | 1   | 5.954        | 5.955            | -0.001           | 107420   | 0.4004           | 0.3877             |       |
| 5 HMx                         | 1   | 6.594        | 6.595            | -0.001           | 67993    | 0.4000           | 0.3622             |       |
| 6 MNx                         | 1   | 7.488        | 7.489            | -0.001           | 114487   | 0.4668           | 0.4365             |       |
| 7 2,4,6-Trinitrophenol        | 1   | 8.008        | 8.029            | -0.021           | 59892    | 0.4000           | 0.3714             |       |
| 8 RDX                         | 1   | 8.908        | 8.915            | -0.007           | 79113    | 0.4000           | 0.3638             |       |
| 9 Nitrobenzene                | 1   | 11.461       | 11.482           | -0.021           | 141679   | 0.4000           | 0.3552             |       |
| \$ 10 1,2-Dinitrobenzene      | 1   | 12.448       | 12.469           | -0.021           | 98279    | 0.4000           | 0.3501             |       |
| 11 3,5-Dinitroaniline         | 1   | 14.334       | 14.355           | -0.021           | 165924   | 0.4000           | 0.3719             |       |
| 12 1,3-Dinitrobenzene         | 1   | 14.601       | 14.622           | -0.021           | 221798   | 0.4000           | 0.3524             |       |
| 13 Nitroglycerin              | 2   | 15.154       | 15.182           | -0.028           | 480119   | 4.00             | 3.56               |       |
| 14 o-Nitrotoluene             | 1   | 15.708       | 15.729           | -0.021           | 90266    | 0.4000           | 0.3655             |       |
| 16 p-Nitrotoluene             | 1   | 15.934       | 15.955           | -0.021           | 81428    | 0.4000           | 0.3579             |       |
| 17 4-Amino-2,6-dinitrotoluene | 1   | 16.501       | 16.515           | -0.014           | 104292   | 0.4000           | 0.3545             |       |
| 18 m-Nitrotoluene             | 1   | 16.801       | 16.815           | -0.014           | 104338   | 0.4000           | 0.3541             |       |
| 19 2-Amino-4,6-dinitrotoluene | 1   | 17.354       | 17.362           | -0.008           | 158919   | 0.4000           | 0.3714             | Ma    |
| 20 1,3,5-Trinitrobenzene      | 1   | 17.414       | 17.422           | -0.008           | 141795   | 0.4000           | 0.3495             | Ma    |
| 21 2,6-Dinitrotoluene         | 1   | 18.628       | 18.629           | -0.001           | 103148   | 0.4000           | 0.3536             |       |
| 22 2,4-Dinitrotoluene         | 1   | 19.081       | 19.082           | -0.001           | 207937   | 0.4000           | 0.3558             |       |
| 23 Tetryl                     | 1   | 22.361       | 22.389           | -0.028           | 98615    | 0.4000           | 0.3463             |       |
| 24 2,4,6-Trinitrotoluene      | 1   | 23.154       | 23.182           | -0.028           | 152188   | 0.4000           | 0.3468             |       |
| 25 PETN                       | 2   | 24.428       | 24.475           | -0.047           | 512590   | 4.00             | 3.56               |       |

### QC Flag Legend

Processing Flags

Review Flags

M - Manually Integrated

a - User Assigned ID

**Reagents:**

8330IntermStk\_00076

Amount Added: 40.00

Units: uL

8330 DMT\_00013

Amount Added: 20.00

Units: uL

Report Date: 29-Apr-2023 13:04:18

Chrom Revision: 2.3 29-Mar-2023 18:39:10

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280013.D

Injection Date: 28-Apr-2023 20:08:59

Instrument ID: CHHPLC\_X5

Operator ID: JZ/JG

Lims ID: IC INT 6

Worklist Smp#: 13

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 13

Method: 8330\_X5\_Luna

Limit Group: GCSV - 8330


Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280013.D  
Injection Date: 28-Apr-2023 20:08:59 Instrument ID: CHHPLC\_X5  
Lims ID: IC INT 6  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 13 Worklist Smp#: 13  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X5\_Luna Limit Group: GCSV - 8330  
Column: Luna-Phenyl hexyl ( 4.60 mm) Detector LC mwd1A, 254 nm

**19 2-Amino-4,6-dinitrotoluene, CAS: 35572-78-2**

Signal: 1


Not Detected  
Expected RT: 17.36

## Processing Integration Results



RT: 17.35  
Area: 158919  
Amount: 0.371359  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LV5D, 28-Apr-2023 20:52:42

Audit Action: Manually Integrated/Assigned Compound ID Audit Reason: Baseline Smoothing

## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280013.D  
Injection Date: 28-Apr-2023 20:08:59 Instrument ID: CHHPLC\_X5  
Lims ID: IC INT 6  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 13 Worklist Smp#: 13  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X5\_Luna Limit Group: GCSV - 8330  
Column: Luna-Phenyl hexyl ( 4.60 mm) Detector LC mwd1A, 254 nm

**20 1,3,5-Trinitrobenzene, CAS: 99-35-4**

Signal: 1


RT: 17.35  
Area: 300702  
Amount: 0.613228  
Amount Units: ug/ml

## Processing Integration Results



RT: 17.41  
Area: 141795  
Amount: 0.349514  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LV5D, 28-Apr-2023 20:52:40

Audit Action: Manually Integrated/Assigned Compound ID Audit Reason: Baseline Smoothing

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280014.D  
 Lims ID: IC INT 5  
 Client ID:  
 Sample Type: IC Calib Level: 5  
 Inject. Date: 28-Apr-2023 20:43:54 ALS Bottle#: 14 Worklist Smp#: 14  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: IC INT 5  
 Operator ID: JZ/JG Instrument ID: CHHPLC\_X5  
 Sublist: chrom-8330\_X5\_Luna\*sub6  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\8330\_X5\_Luna.m  
 Limit Group: GCSV - 8330  
 Last Update: 29-Apr-2023 13:04:18 Calib Date: 29-Apr-2023 04:18:01  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280027.D  
 Column 1 : Luna-Phenyl hexyl ( 4.60 mm) Det: LC mwd1A, 254 nm  
 Process Host: CTX1621

First Level Reviewer: LV5D

Date: 29-Apr-2023 10:22:21

| Compound                      | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 3 TNX                         | 1   | 5.115        | 5.115            | 0.000            | 86123    | 0.2508           | 0.2184             |       |
| 4 DNx                         | 1   | 5.955        | 5.955            | 0.000            | 62056    | 0.2503           | 0.2240             |       |
| 5 HMX                         | 1   | 6.595        | 6.595            | 0.000            | 46110    | 0.2500           | 0.2442             |       |
| 6 MNX                         | 1   | 7.489        | 7.489            | 0.000            | 66549    | 0.2918           | 0.2537             |       |
| 7 2,4,6-Trinitrophenol        | 1   | 8.029        | 8.029            | 0.000            | 40202    | 0.2500           | 0.2481             |       |
| 8 RDX                         | 1   | 8.915        | 8.915            | 0.000            | 53643    | 0.2500           | 0.2456             |       |
| 9 Nitrobenzene                | 1   | 11.482       | 11.482           | 0.000            | 96681    | 0.2500           | 0.2424             |       |
| \$ 10 1,2-Dinitrobenzene      | 1   | 12.469       | 12.469           | 0.000            | 67346    | 0.2500           | 0.2399             |       |
| 11 3,5-Dinitroaniline         | 1   | 14.355       | 14.355           | 0.000            | 113054   | 0.2500           | 0.2529             |       |
| 12 1,3-Dinitrobenzene         | 1   | 14.622       | 14.622           | 0.000            | 152658   | 0.2500           | 0.2426             |       |
| 13 Nitroglycerin              | 2   | 15.182       | 15.182           | 0.000            | 327195   | 2.50             | 2.43               |       |
| 14 o-Nitrotoluene             | 1   | 15.729       | 15.729           | 0.000            | 60674    | 0.2500           | 0.2457             |       |
| 16 p-Nitrotoluene             | 1   | 15.955       | 15.955           | 0.000            | 56555    | 0.2500           | 0.2486             |       |
| 17 4-Amino-2,6-dinitrotoluene | 1   | 16.515       | 16.515           | 0.000            | 71678    | 0.2500           | 0.2436             |       |
| 18 m-Nitrotoluene             | 1   | 16.815       | 16.815           | 0.000            | 71041    | 0.2500           | 0.2411             |       |
| 19 2-Amino-4,6-dinitrotoluene | 1   | 17.362       | 17.362           | 0.000            | 108505   | 0.2500           | 0.2536             | M     |
| 20 1,3,5-Trinitrobenzene      | 1   | 17.422       | 17.422           | 0.000            | 98168    | 0.2500           | 0.2420             | M     |
| 21 2,6-Dinitrotoluene         | 1   | 18.629       | 18.629           | 0.000            | 70868    | 0.2500           | 0.2430             |       |
| 22 2,4-Dinitrotoluene         | 1   | 19.082       | 19.082           | 0.000            | 143098   | 0.2500           | 0.2448             |       |
| 23 Tetryl                     | 1   | 22.389       | 22.389           | 0.000            | 68327    | 0.2500           | 0.2400             |       |
| 24 2,4,6-Trinitrotoluene      | 1   | 23.182       | 23.182           | 0.000            | 104561   | 0.2500           | 0.2382             |       |
| 25 PETN                       | 2   | 24.475       | 24.475           | 0.000            | 351913   | 2.50             | 2.44               |       |

**QC Flag Legend**

Processing Flags

Review Flags

M - Manually Integrated



**Reagents:**

8330IntermStk\_00076

Amount Added: 25.00

Units: uL

8330 DMT\_00013

Amount Added: 12.50

Units: uL

Report Date: 29-Apr-2023 13:04:18

Chrom Revision: 2.3 29-Mar-2023 18:39:10

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280014.D

Injection Date: 28-Apr-2023 20:43:54

Instrument ID: CHHPLC\_X5

Operator ID: JZ/JG

Lims ID: IC INT 5

Worklist Smp#: 14

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 14

Method: 8330\_X5\_Luna

Limit Group: GCSV - 8330


Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280014.D  
Injection Date: 28-Apr-2023 20:43:54 Instrument ID: CHHPLC\_X5  
Lims ID: IC INT 5  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 14 Worklist Smp#: 14  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X5\_Luna Limit Group: GCSV - 8330  
Column: Luna-Phenyl hexyl ( 4.60 mm) Detector LC mwd1A, 254 nm

**20 1,3,5-Trinitrobenzene, CAS: 99-35-4**

Signal: 1


RT: 17.36  
Area: 206659  
Amount: 0.382644  
Amount Units: ug/ml

## Processing Integration Results



RT: 17.42  
Area: 98168  
Amount: 0.241976  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LV5D, 29-Apr-2023 10:22:17

Audit Action: Split an Integrated Peak

Audit Reason: Baseline Smoothing

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280015.D  
 Lims ID: IC INT 4  
 Client ID:  
 Sample Type: IC Calib Level: 4  
 Inject. Date: 28-Apr-2023 21:18:48 ALS Bottle#: 15 Worklist Smp#: 15  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: IC INT 4  
 Operator ID: JZ/JG Instrument ID: CHHPLC\_X5  
 Sublist: chrom-8330\_X5\_Luna\*sub6  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\8330\_X5\_Luna.m  
 Limit Group: GCSV - 8330  
 Last Update: 29-Apr-2023 13:04:19 Calib Date: 29-Apr-2023 04:18:01  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280027.D  
 Column 1 : Luna-Phenyl hexyl ( 4.60 mm) Det: LC mwd1A, 254 nm  
 Process Host: CTX1621

First Level Reviewer: LV5D

Date: 29-Apr-2023 10:22:33

| Compound                      | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 3 TNX                         | 1   | 5.116        | 5.115            | 0.001            | 41754    | 0.1003           | 0.1059             |       |
| 4 DNx                         | 1   | 5.963        | 5.955            | 0.008            | 28641    | 0.1001           | 0.1034             |       |
| 5 HMX                         | 1   | 6.603        | 6.595            | 0.008            | 19573    | 0.1000           | 0.1011             |       |
| 6 MNX                         | 1   | 7.496        | 7.489            | 0.007            | 30530    | 0.1167           | 0.1164             |       |
| 7 2,4,6-Trinitrophenol        | 1   | 8.056        | 8.029            | 0.027            | 17060    | 0.1000           | 0.1033             |       |
| 8 RDX                         | 1   | 8.923        | 8.915            | 0.008            | 22827    | 0.1000           | 0.1027             |       |
| 9 Nitrobenzene                | 1   | 11.483       | 11.482           | 0.001            | 41585    | 0.1000           | 0.1043             |       |
| \$ 10 1,2-Dinitrobenzene      | 1   | 12.470       | 12.469           | 0.001            | 27893    | 0.1000           | 0.0994             |       |
| 11 3,5-Dinitroaniline         | 1   | 14.363       | 14.355           | 0.008            | 46980    | 0.1000           | 0.1043             |       |
| 12 1,3-Dinitrobenzene         | 1   | 14.623       | 14.622           | 0.001            | 62596    | 0.1000           | 0.0995             |       |
| 13 Nitroglycerin              | 2   | 15.176       | 15.182           | -0.006           | 134661   | 1.00             | 1.00               |       |
| 14 o-Nitrotoluene             | 1   | 15.730       | 15.729           | 0.001            | 25290    | 0.1000           | 0.1024             |       |
| 16 p-Nitrotoluene             | 1   | 15.956       | 15.955           | 0.001            | 23002    | 0.1000           | 0.1011             |       |
| 17 4-Amino-2,6-dinitrotoluene | 1   | 16.523       | 16.515           | 0.008            | 29472    | 0.1000           | 0.1002             |       |
| 18 m-Nitrotoluene             | 1   | 16.816       | 16.815           | 0.001            | 29569    | 0.1000           | 0.1004             |       |
| 19 2-Amino-4,6-dinitrotoluene | 1   | 17.370       | 17.362           | 0.008            | 43948    | 0.1000           | 0.1027             | M     |
| 20 1,3,5-Trinitrobenzene      | 1   | 17.423       | 17.422           | 0.001            | 41946    | 0.1000           | 0.1034             | M     |
| 21 2,6-Dinitrotoluene         | 1   | 18.636       | 18.629           | 0.007            | 30260    | 0.1000           | 0.1037             |       |
| 22 2,4-Dinitrotoluene         | 1   | 19.090       | 19.082           | 0.008            | 59433    | 0.1000           | 0.1017             |       |
| 23 Tetryl                     | 1   | 22.363       | 22.389           | -0.026           | 27534    | 0.1000           | 0.0967             |       |
| 24 2,4,6-Trinitrotoluene      | 1   | 23.156       | 23.182           | -0.026           | 43136    | 0.1000           | 0.0983             |       |
| 25 PETN                       | 2   | 24.416       | 24.475           | -0.059           | 146450   | 1.00             | 1.02               |       |

**QC Flag Legend**

Processing Flags

Review Flags

M - Manually Integrated

**Reagents:**

8330IntermStk\_00076

Amount Added: 10.00

Units: uL

8330 DMT\_00013

Amount Added: 5.00

Units: uL

Report Date: 29-Apr-2023 13:04:19

Chrom Revision: 2.3 29-Mar-2023 18:39:10

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280015.D

Injection Date: 28-Apr-2023 21:18:48

Instrument ID: CHHPLC\_X5

Operator ID: JZ/JG

Lims ID: IC INT 4

Worklist Smp#: 15

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 15

Method: 8330\_X5\_Luna

Limit Group: GCSV - 8330


Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280015.D  
Injection Date: 28-Apr-2023 21:18:48 Instrument ID: CHHPLC\_X5  
Lims ID: IC INT 4  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 15 Worklist Smp#: 15  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X5\_Luna Limit Group: GCSV - 8330  
Column: Luna-Phenyl hexyl ( 4.60 mm) Detector LC mwd1A, 254 nm

**20 1,3,5-Trinitrobenzene, CAS: 99-35-4**

Signal: 1


RT: 17.37  
Area: 85885  
Amount: 0.174611  
Amount Units: ug/ml

## Processing Integration Results



RT: 17.42  
Area: 41946  
Amount: 0.103394  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LV5D, 29-Apr-2023 10:22:30  
Audit Action: Split an Integrated Peak

Audit Reason: Baseline Smoothing

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280016.D  
 Lims ID: IC INT 3  
 Client ID:  
 Sample Type: IC Calib Level: 3  
 Inject. Date: 28-Apr-2023 21:53:46 ALS Bottle#: 16 Worklist Smp#: 16  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: IC INT 3  
 Operator ID: JZ/JG Instrument ID: CHHPLC\_X5  
 Sublist: chrom-8330\_X5\_Luna\*sub6  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\8330\_X5\_Luna.m  
 Limit Group: GCSV - 8330  
 Last Update: 29-Apr-2023 13:04:19 Calib Date: 29-Apr-2023 04:18:01  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280027.D  
 Column 1 : Luna-Phenyl hexyl ( 4.60 mm) Det: LC mwd1A, 254 nm  
 Process Host: CTX1621

First Level Reviewer: LV5D

Date: 29-Apr-2023 10:22:43

| Compound                      | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 3 TNX                         | 1   | 5.116        | 5.115            | 0.001            | 20281    | 0.0502           | 0.0514             |       |
| 4 DNX                         | 1   | 5.969        | 5.955            | 0.014            | 14733    | 0.0501           | 0.0532             |       |
| 5 HMX                         | 1   | 6.603        | 6.595            | 0.008            | 10761    | 0.0500           | 0.0535             |       |
| 6 MNX                         | 1   | 7.503        | 7.489            | 0.014            | 15694    | 0.0584           | 0.0598             |       |
| 7 2,4,6-Trinitrophenol        | 1   | 8.076        | 8.029            | 0.047            | 8950     | 0.0500           | 0.0526             |       |
| 8 RDX                         | 1   | 8.929        | 8.915            | 0.014            | 11621    | 0.0500           | 0.0508             |       |
| 9 Nitrobenzene                | 1   | 11.489       | 11.482           | 0.007            | 20883    | 0.0500           | 0.0524             |       |
| \$ 10 1,2-Dinitrobenzene      | 1   | 12.476       | 12.469           | 0.007            | 14379    | 0.0500           | 0.0512             |       |
| 11 3,5-Dinitroaniline         | 1   | 14.369       | 14.355           | 0.014            | 24176    | 0.0500           | 0.0530             |       |
| 12 1,3-Dinitrobenzene         | 1   | 14.629       | 14.622           | 0.007            | 31459    | 0.0500           | 0.0500             |       |
| 13 Nitroglycerin              | 2   | 15.176       | 15.182           | -0.006           | 71084    | 0.5000           | 0.5276             |       |
| 14 o-Nitrotoluene             | 1   | 15.723       | 15.729           | -0.006           | 12951    | 0.0500           | 0.0524             |       |
| 16 p-Nitrotoluene             | 1   | 15.956       | 15.955           | 0.001            | 12179    | 0.0500           | 0.0535             |       |
| 17 4-Amino-2,6-dinitrotoluene | 1   | 16.523       | 16.515           | 0.008            | 15265    | 0.0500           | 0.0519             |       |
| 18 m-Nitrotoluene             | 1   | 16.816       | 16.815           | 0.001            | 15546    | 0.0500           | 0.0528             |       |
| 19 2-Amino-4,6-dinitrotoluene | 1   | 17.376       | 17.362           | 0.014            | 23162    | 0.0500           | 0.0541             | M     |
| 20 1,3,5-Trinitrobenzene      | 1   | 17.429       | 17.422           | 0.007            | 22059    | 0.0500           | 0.0544             | M     |
| 21 2,6-Dinitrotoluene         | 1   | 18.643       | 18.629           | 0.014            | 15876    | 0.0500           | 0.0544             |       |
| 22 2,4-Dinitrotoluene         | 1   | 19.096       | 19.082           | 0.014            | 30407    | 0.0500           | 0.0520             |       |
| 23 Tetryl                     | 1   | 22.369       | 22.389           | -0.020           | 14233    | 0.0500           | 0.0500             |       |
| 24 2,4,6-Trinitrotoluene      | 1   | 23.163       | 23.182           | -0.019           | 22184    | 0.0500           | 0.0505             |       |
| 25 PETN                       | 2   | 24.423       | 24.475           | -0.052           | 72943    | 0.5000           | 0.5068             |       |

**QC Flag Legend**

Processing Flags

Review Flags

M - Manually Integrated



**Reagents:**

8330IntermStk\_00076

Amount Added: 5.00

Units: uL

8330 DMT\_00013

Amount Added: 2.50

Units: uL

Report Date: 29-Apr-2023 13:04:19

Chrom Revision: 2.3 29-Mar-2023 18:39:10

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280016.D

Injection Date: 28-Apr-2023 21:53:46

Instrument ID: CHHPLC\_X5

Operator ID: JZ/JG

Lims ID: IC INT 3

Worklist Smp#: 16

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 16

Method: 8330\_X5\_Luna

Limit Group: GCSV - 8330


Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280016.D  
Injection Date: 28-Apr-2023 21:53:46 Instrument ID: CHHPLC\_X5  
Lims ID: IC INT 3  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 16 Worklist Smp#: 16  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X5\_Luna Limit Group: GCSV - 8330  
Column: Luna-Phenyl hexyl ( 4.60 mm) Detector LC mwd1A, 254 nm

**20 1,3,5-Trinitrobenzene, CAS: 99-35-4**

Signal: 1


RT: 17.38  
Area: 45204  
Amount: 0.102031  
Amount Units: ug/ml

## Processing Integration Results



RT: 17.43  
Area: 22059  
Amount: 0.054374  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LV5D, 29-Apr-2023 10:22:40

Audit Action: Split an Integrated Peak

Audit Reason: Baseline Smoothing

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280017.D  
 Lims ID: IC INT 2  
 Client ID:  
 Sample Type: IC Calib Level: 2  
 Inject. Date: 28-Apr-2023 22:28:39 ALS Bottle#: 17 Worklist Smp#: 17  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: IC INT 2  
 Operator ID: JZ/JG Instrument ID: CHHPLC\_X5  
 Sublist: chrom-8330\_X5\_Luna\*sub6  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\8330\_X5\_Luna.m  
 Limit Group: GCSV - 8330  
 Last Update: 29-Apr-2023 13:04:20 Calib Date: 29-Apr-2023 04:18:01  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280027.D  
 Column 1 : Luna-Phenyl hexyl ( 4.60 mm) Det: LC mwd1A, 254 nm  
 Process Host: CTX1621

First Level Reviewer: LV5D

Date: 29-Apr-2023 10:23:35

| Compound                      | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 3 TNX                         | 1   | 5.136        | 5.115            | 0.021            | 8179     | 0.0201           | 0.0207             |       |
| 4 DNx                         | 1   | 5.983        | 5.955            | 0.028            | 5728     | 0.0200           | 0.0207             |       |
| 5 HMX                         | 1   | 6.630        | 6.595            | 0.035            | 4449     | 0.0200           | 0.0195             |       |
| 6 MNX                         | 1   | 7.523        | 7.489            | 0.034            | 6155     | 0.0233           | 0.0235             |       |
| 7 2,4,6-Trinitrophenol        | 1   | 8.116        | 8.029            | 0.087            | 3453     | 0.0200           | 0.0182             |       |
| 8 RDX                         | 1   | 8.956        | 8.915            | 0.041            | 4863     | 0.0200           | 0.0194             |       |
| 9 Nitrobenzene                | 1   | 11.523       | 11.482           | 0.041            | 7871     | 0.0200           | 0.0197             |       |
| \$ 10 1,2-Dinitrobenzene      | 1   | 12.516       | 12.469           | 0.047            | 6038     | 0.0200           | 0.0215             |       |
| 11 3,5-Dinitroaniline         | 1   | 14.410       | 14.355           | 0.055            | 9221     | 0.0200           | 0.0194             |       |
| 12 1,3-Dinitrobenzene         | 1   | 14.663       | 14.622           | 0.041            | 12307    | 0.0200           | 0.0196             |       |
| 13 Nitroglycerin              | 2   | 15.216       | 15.182           | 0.034            | 27409    | 0.2000           | 0.2034             |       |
| 14 o-Nitrotoluene             | 1   | 15.770       | 15.729           | 0.041            | 5281     | 0.0200           | 0.0214             | M     |
| 16 p-Nitrotoluene             | 1   | 15.996       | 15.955           | 0.041            | 4466     | 0.0200           | 0.0196             | M     |
| 17 4-Amino-2,6-dinitrotoluene | 1   | 16.556       | 16.515           | 0.041            | 6064     | 0.0200           | 0.0206             | M     |
| 18 m-Nitrotoluene             | 1   | 16.843       | 16.815           | 0.028            | 5823     | 0.0200           | 0.0198             | M     |
| 19 2-Amino-4,6-dinitrotoluene | 1   | 17.410       | 17.362           | 0.048            | 8257     | 0.0200           | 0.0193             | M     |
| 20 1,3,5-Trinitrobenzene      | 1   | 17.456       | 17.422           | 0.034            | 7799     | 0.0200           | 0.0192             | M     |
| 21 2,6-Dinitrotoluene         | 1   | 18.670       | 18.629           | 0.041            | 5716     | 0.0200           | 0.0196             |       |
| 22 2,4-Dinitrotoluene         | 1   | 19.123       | 19.082           | 0.041            | 12241    | 0.0200           | 0.0209             |       |
| 23 Tetryl                     | 1   | 22.396       | 22.389           | 0.007            | 6348     | 0.0200           | 0.0223             |       |
| 24 2,4,6-Trinitrotoluene      | 1   | 23.196       | 23.182           | 0.014            | 9925     | 0.0200           | 0.0226             |       |
| 25 PETN                       | 2   | 24.450       | 24.475           | -0.025           | 28805    | 0.2000           | 0.2001             |       |

**QC Flag Legend**

Processing Flags

Review Flags

M - Manually Integrated

**Reagents:**

8330IntermStk\_00076

Amount Added: 2.00

Units: uL

8330 DMT\_00013

Amount Added: 1.00

Units: uL

Report Date: 29-Apr-2023 13:04:20

Chrom Revision: 2.3 29-Mar-2023 18:39:10

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280017.D

Injection Date: 28-Apr-2023 22:28:39

Instrument ID: CHHPLC\_X5

Operator ID: JZ/JG

Lims ID: IC INT 2

Worklist Smp#: 17

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 17

Method: 8330\_X5\_Luna

Limit Group: GCSV - 8330


Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280017.D  
Injection Date: 28-Apr-2023 22:28:39 Instrument ID: CHHPLC\_X5  
Lims ID: IC INT 2  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 17 Worklist Smp#: 17  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X5\_Luna Limit Group: GCSV - 8330  
Column: Luna-Phenyl hexyl ( 4.60 mm) Detector LC mwd1A, 254 nm

**14 o-Nitrotoluene, CAS: 88-72-2**

Signal: 1


RT: 15.77  
Area: 8594  
Amount: 0.026168  
Amount Units: ug/ml

## Processing Integration Results



RT: 15.77  
Area: 5281  
Amount: 0.021384  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LV5D, 29-Apr-2023 10:23:05  
Audit Action: Split an Integrated Peak

Audit Reason: Baseline Smoothing

## Eurofins Denver

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280017.D

Injection Date: 28-Apr-2023 22:28:39

Instrument ID: CHHPLC\_X5

Lims ID: IC INT 2

Client ID:

Operator ID: JZ/JG

ALS Bottle#: 17 Worklist Smp#: 17

Injection Vol: 100.0 ul

Dil. Factor: 1.0000

Method: 8330\_X5\_Luna

Limit Group: GCSV - 8330

Column: Luna-Phenyl hexyl ( 4.60 mm)


Detector LC mwd1A, 254 nm

**17 4-Amino-2,6-dinitrotoluene, CAS: 19406-51-0**

Signal: 1


RT: 16.56  
Area: 4498  
Amount: 0.016304  
Amount Units: ug/ml

## Processing Integration Results



RT: 16.56  
Area: 6064  
Amount: 0.020611  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LV5D, 29-Apr-2023 10:23:33

Audit Action: Assigned New Baseline

Audit Reason: Baseline Smoothing



## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280017.D  
Injection Date: 28-Apr-2023 22:28:39 Instrument ID: CHHPLC\_X5  
Lims ID: IC INT 2  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 17 Worklist Smp#: 17  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X5\_Luna Limit Group: GCSV - 8330  
Column: Luna-Phenyl hexyl ( 4.60 mm) Detector LC mwd1A, 254 nm

**18 m-Nitrotoluene, CAS: 99-08-1**

Signal: 1


RT: 16.84  
Area: 4234  
Amount: 0.015427  
Amount Units: ug/ml

## Processing Integration Results



RT: 16.84  
Area: 5823  
Amount: 0.019763  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LV5D, 29-Apr-2023 10:23:33  
Audit Action: Assigned New Baseline

Audit Reason: Baseline Smoothing

## Eurofins Denver

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280017.D

Injection Date: 28-Apr-2023 22:28:39

Instrument ID: CHHPLC\_X5

Lims ID: IC INT 2

Client ID:

Operator ID: JZ/JG

ALS Bottle#:

17

Worklist Smp#: 17

Injection Vol: 100.0 ul

Dil. Factor: 1.0000

Method: 8330\_X5\_Luna

Limit Group: GCSV - 8330

Column: Luna-Phenyl hexyl ( 4.60 mm)

Detector


LC mwd1A, 254 nm

**20 1,3,5-Trinitrobenzene, CAS: 99-35-4**

Signal: 1


RT: 17.41  
Area: 6616  
Amount: 0.016894  
Amount Units: ug/ml

## Processing Integration Results



RT: 17.46  
Area: 7799  
Amount: 0.019224  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LV5D, 29-Apr-2023 10:23:23

Audit Action: Split an Integrated Peak

Audit Reason: Baseline Smoothing

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280018.D  
 Lims ID: IC INT 1  
 Client ID:  
 Sample Type: IC Calib Level: 1  
 Inject. Date: 28-Apr-2023 23:03:35 ALS Bottle#: 18 Worklist Smp#: 18  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: IC INT 1  
 Operator ID: JZ/JG Instrument ID: CHHPLC\_X5  
 Sublist: chrom-8330\_X5\_Luna\*sub6  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\8330\_X5\_Luna.m  
 Limit Group: GCSV - 8330  
 Last Update: 29-Apr-2023 13:04:20 Calib Date: 29-Apr-2023 04:18:01  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280027.D  
 Column 1 : Luna-Phenyl hexyl ( 4.60 mm) Det: LC mwd1A, 254 nm  
 Process Host: CTX1621

First Level Reviewer: LV5D

Date: 29-Apr-2023 10:24:45

| Compound                      | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 3 TNX                         | 1   | 5.116        | 5.115            | 0.001            | 4552     | 0.0100           | 0.0115             |       |
| 4 DNx                         | 1   | 5.963        | 5.955            | 0.008            | 2804     | 0.0100           | 0.0101             |       |
| 5 HMX                         | 1   | 6.596        | 6.595            | 0.001            | 2178     | 0.0100           | 0.007259           |       |
| 6 MNX                         | 1   | 7.496        | 7.489            | 0.007            | 3646     | 0.0117           | 0.0139             |       |
| 7 2,4,6-Trinitrophenol        | 1   | 8.063        | 8.029            | 0.034            | 2197     | 0.0100           | 0.0103             |       |
| 8 RDX                         | 1   | 8.923        | 8.915            | 0.008            | 2852     | 0.0100           | 0.0101             |       |
| 9 Nitrobenzene                | 1   | 11.489       | 11.482           | 0.007            | 4269     | 0.0100           | 0.0107             |       |
| \$ 10 1,2-Dinitrobenzene      | 1   | 12.483       | 12.469           | 0.014            | 3169     | 0.0100           | 0.0113             |       |
| 11 3,5-Dinitroaniline         | 1   | 14.376       | 14.355           | 0.021            | 5053     | 0.0100           | 0.0100             |       |
| 12 1,3-Dinitrobenzene         | 1   | 14.636       | 14.622           | 0.014            | 6946     | 0.0100           | 0.0110             |       |
| 13 Nitroglycerin              | 2   | 15.196       | 15.182           | 0.014            | 14603    | 0.1000           | 0.1084             |       |
| 14 o-Nitrotoluene             | 1   | 15.736       | 15.729           | 0.007            | 2336     | 0.0100           | 0.009459           | M     |
| 16 p-Nitrotoluene             | 1   | 15.969       | 15.955           | 0.014            | 2206     | 0.0100           | 0.009697           | M     |
| 17 4-Amino-2,6-dinitrotoluene | 1   | 16.529       | 16.515           | 0.014            | 3275     | 0.0100           | 0.0111             | M     |
| 18 m-Nitrotoluene             | 1   | 16.843       | 16.815           | 0.028            | 3269     | 0.0100           | 0.0111             | M     |
| 19 2-Amino-4,6-dinitrotoluene | 1   | 17.369       | 17.362           | 0.007            | 3912     | 0.0100           | 0.009141           | M     |
| 20 1,3,5-Trinitrobenzene      | 1   | 17.443       | 17.422           | 0.021            | 4438     | 0.0100           | 0.0109             | M     |
| 21 2,6-Dinitrotoluene         | 1   | 18.643       | 18.629           | 0.014            | 3093     | 0.0100           | 0.0106             |       |
| 22 2,4-Dinitrotoluene         | 1   | 19.096       | 19.082           | 0.014            | 5987     | 0.0100           | 0.0102             |       |
| 23 Tetryl                     | 1   | 22.363       | 22.389           | -0.026           | 3282     | 0.0100           | 0.0115             |       |
| 24 2,4,6-Trinitrotoluene      | 1   | 23.176       | 23.182           | -0.006           | 4919     | 0.0100           | 0.0112             |       |
| 25 PETN                       | 2   | 24.430       | 24.475           | -0.045           | 16344    | 0.1000           | 0.1135             |       |

## QC Flag Legend

Processing Flags

Review Flags

M - Manually Integrated

**Reagents:**

8330IntermStk\_00076

Amount Added: 1.00

Units: uL

8330 DMT\_00013

Amount Added: 0.50

Units: uL

Report Date: 29-Apr-2023 13:04:20

Chrom Revision: 2.3 29-Mar-2023 18:39:10

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280018.D

Injection Date: 28-Apr-2023 23:03:35

Instrument ID: CHHPLC\_X5

Operator ID: JZ/JG

Lims ID: IC INT 1

Worklist Smp#: 18

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 18

Method: 8330\_X5\_Luna

Limit Group: GCSV - 8330


Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280018.D  
Injection Date: 28-Apr-2023 23:03:35 Instrument ID: CHHPLC\_X5  
Lims ID: IC INT 1  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 18 Worklist Smp#: 18  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X5\_Luna Limit Group: GCSV - 8330  
Column: Luna-Phenyl hexyl ( 4.60 mm) Detector LC mwd1A, 254 nm

**14 o-Nitrotoluene, CAS: 88-72-2**

Signal: 1


RT: 15.74  
Area: 4530  
Amount: 0.013152  
Amount Units: ug/ml

## Processing Integration Results



RT: 15.74  
Area: 2336  
Amount: 0.009459  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LV5D, 29-Apr-2023 10:24:06  
Audit Action: Split an Integrated Peak

Audit Reason: Baseline Smoothing

## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280018.D  
Injection Date: 28-Apr-2023 23:03:35 Instrument ID: CHHPLC\_X5  
Lims ID: IC INT 1  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 18 Worklist Smp#: 18  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X5\_Luna Limit Group: GCSV - 8330  
Column: Luna-Phenyl hexyl ( 4.60 mm) Detector LC mwd1A, 254 nm

**17 4-Amino-2,6-dinitrotoluene, CAS: 19406-51-0**

Signal: 1


RT: 16.53  
Area: 2409  
Amount: 0.008465  
Amount Units: ug/ml

## Processing Integration Results



RT: 16.53  
Area: 3275  
Amount: 0.011131  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LV5D, 29-Apr-2023 13:01:58  
Audit Action: Assigned New Baseline

Audit Reason: Baseline Smoothing

## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280018.D  
Injection Date: 28-Apr-2023 23:03:35 Instrument ID: CHHPLC\_X5  
Lims ID: IC INT 1  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 18 Worklist Smp#: 18  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X5\_Luna Limit Group: GCSV - 8330  
Column: Luna-Phenyl hexyl ( 4.60 mm) Detector LC mwd1A, 254 nm

**18 m-Nitrotoluene, CAS: 99-08-1**

Signal: 1


RT: 16.84  
Area: 2248  
Amount: 0.007935  
Amount Units: ug/ml

## Processing Integration Results



RT: 16.84  
Area: 3269  
Amount: 0.011095  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LV5D, 29-Apr-2023 13:01:58  
Audit Action: Assigned New Baseline

Audit Reason: Baseline Smoothing



## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280018.D  
Injection Date: 28-Apr-2023 23:03:35 Instrument ID: CHHPLC\_X5  
Lims ID: IC INT 1  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 18 Worklist Smp#: 18  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X5\_Luna Limit Group: GCSV - 8330  
Column: Luna-Phenyl hexyl ( 4.60 mm) Detector LC mwd1A, 254 nm

**20 1,3,5-Trinitrobenzene, CAS: 99-35-4**

Signal: 1


RT: 17.38  
Area: 3762  
Amount: 0.009448  
Amount Units: ug/ml

## Processing Integration Results



RT: 17.44  
Area: 4438  
Amount: 0.010939  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LV5D, 29-Apr-2023 10:24:15  
Audit Action: Split an Integrated Peak

Audit Reason: Baseline Smoothing

# Calibration

/ TNX

Curve Type: Average  
Weighting: Conc\_Sq  
Origin: Force  
Dependency: Response  
Calib Mode: ESTD  
Response Base: AREA  
RF Rounding: 0


## Curve Coefficients

Intercept: 0  
Slope: 3.943E+05

## Error Coefficients

Standard Error: 22100  
Relative Standard Error: 8.1  
Correlation Coefficient: 0.999  
Coefficient of Determination (Adjusted): 0.991

| ID | Level            | Concentration | Response | IS Amount | IS Response | RF            | Used |
|----|------------------|---------------|----------|-----------|-------------|---------------|------|
| 1  | IC 280-610603/18 | 0.01003       | 4552.0   |           |             | 453838.484546 | Y    |
| 2  | IC 280-610603/17 | 0.02006       | 8179.0   |           |             | 407726.819541 | Y    |
| 3  | IC 280-610603/16 | 0.05015       | 20281.0  |           |             | 404406.779661 | Y    |
| 4  | IC 280-610603/15 | 0.1003        | 41754.0  |           |             | 416291.12662  | Y    |
| 5  | IC 280-610603/14 | 0.25075       | 86123.0  |           |             | 343461.615155 | Y    |
| 6  | IC 280-610603/13 | 0.4012        | 147763.0 |           |             | 368302.592223 | Y    |
| 7  | IC 280-610603/12 | 0.7021        | 276939.0 |           |             | 394443.811423 | Y    |
| 8  | IC 280-610603/11 | 1.003         | 391357.0 |           |             | 390186.440678 | Y    |
| 9  | IC 280-610603/10 | 2.5075        | 928605.0 |           |             | 370331.006979 | Y    |



# Calibration

/ DNX

Curve Type: Average  
Weighting: Conc\_Sq  
Origin: Force  
Dependency: Response  
Calib Mode: ESTD  
Response Base: AREA  
RF Rounding: 0


## Curve Coefficients

Intercept: 0  
Slope: 2.771E+05

## Error Coefficients

Standard Error: 11000  
Relative Standard Error: 5.1  
Correlation Coefficient: 0.999  
Coefficient of Determination (Adjusted): 0.997

| ID | Level            | Concentration | Response | IS Amount | IS Response | RF            | Used |
|----|------------------|---------------|----------|-----------|-------------|---------------|------|
| 1  | IC 280-610603/18 | 0.01001       | 2804.0   |           |             | 280119.88012  | Y    |
| 2  | IC 280-610603/17 | 0.02002       | 5728.0   |           |             | 286113.886114 | Y    |
| 3  | IC 280-610603/16 | 0.05005       | 14733.0  |           |             | 294365.634366 | Y    |
| 4  | IC 280-610603/15 | 0.1001        | 28641.0  |           |             | 286123.876124 | Y    |
| 5  | IC 280-610603/14 | 0.25025       | 62056.0  |           |             | 247976.023976 | Y    |
| 6  | IC 280-610603/13 | 0.4004        | 107420.0 |           |             | 268281.718282 | Y    |
| 7  | IC 280-610603/12 | 0.7007        | 198949.0 |           |             | 283928.928215 | Y    |
| 8  | IC 280-610603/11 | 1.001         | 281655.0 |           |             | 281373.626374 | Y    |
| 9  | IC 280-610603/10 | 2.5025        | 664081.0 |           |             | 265367.032967 | Y    |



# Calibration

/ HMX

Curve Type: Linear  
Weighting: Conc\_Sq  
Origin: None  
Dependency: Response  
Calib Mode: ESTD  
Response Base: AREA  
RF Rounding: 0


## Curve Coefficients

Intercept: 832  
Slope: 1.854E+05

## Error Coefficients

Standard Error: 8620  
Relative Standard Error: 5.4  
Correlation Coefficient: 0.999  
Coefficient of Determination (Adjusted): 0.997

| ID | Level            | Concentration | Response | IS Amount | IS Response | RF       | Used |
|----|------------------|---------------|----------|-----------|-------------|----------|------|
| 1  | IC 280-610603/18 | 0.01          | 2178.0   |           |             | 217800.0 | N    |
| 2  | IC 280-610603/17 | 0.02          | 4449.0   |           |             | 222450.0 | Y    |
| 3  | IC 280-610603/16 | 0.05          | 10761.0  |           |             | 215220.0 | Y    |
| 4  | IC 280-610603/15 | 0.1           | 19573.0  |           |             | 195730.0 | Y    |
| 5  | IC 280-610603/14 | 0.25          | 46110.0  |           |             | 184440.0 | Y    |
| 6  | IC 280-610603/13 | 0.4           | 67993.0  |           |             | 169982.5 | Y    |
| 7  | IC 280-610603/12 | 0.7           | 131068.0 |           |             | 187240.0 | Y    |
| 8  | IC 280-610603/11 | 1.0           | 189027.0 |           |             | 189027.0 | Y    |
| 9  | IC 280-610603/10 | 2.5           | 484081.0 |           |             | 193632.4 | Y    |



# Calibration

/ MNX

Curve Type: Average  
 Weighting: Conc\_Sq  
 Origin: Force  
 Dependency: Response  
 Calib Mode: ESTD  
 Response Base: AREA  
 RF Rounding: 0


## Curve Coefficients

Intercept: 0  
 Slope: 2.623E+05

## Error Coefficients

Standard Error: 6210  
 Relative Standard Error: 8.6  
 Correlation Coefficient: 1.000  
 Coefficient of Determination (Adjusted): 0.990

| ID | Level            | Concentration | Response | IS Amount | IS Response | RF            | Used |
|----|------------------|---------------|----------|-----------|-------------|---------------|------|
| 1  | IC 280-610603/18 | 0.01167       | 3646.0   |           |             | 312425.021422 | Y    |
| 2  | IC 280-610603/17 | 0.02334       | 6155.0   |           |             | 263710.368466 | Y    |
| 3  | IC 280-610603/16 | 0.05835       | 15694.0  |           |             | 268963.153385 | Y    |
| 4  | IC 280-610603/15 | 0.1167        | 30530.0  |           |             | 261610.968295 | Y    |
| 5  | IC 280-610603/14 | 0.29175       | 66549.0  |           |             | 228102.827763 | Y    |
| 6  | IC 280-610603/13 | 0.4668        | 114487.0 |           |             | 245259.211654 | Y    |
| 7  | IC 280-610603/12 | 0.8169        | 211518.0 |           |             | 258927.653324 | Y    |
| 8  | IC 280-610603/11 | 1.167         | 298731.0 |           |             | 255982.005141 | Y    |
| 9  | IC 280-610603/10 | 2.9175        | 774370.0 |           |             | 265422.450728 | Y    |



# Calibration

/ 2,4,6-Trinitrophenol

Curve Type: Linear  
 Weighting: Conc\_Sq  
 Origin: None  
 Dependency: Response  
 Calib Mode: ESTD  
 Response Base: AREA  
 RF Rounding: 0


## Curve Coefficients

Intercept: 546.1  
 Slope: 1.598E+05

## Error Coefficients

Standard Error: 2930  
 Relative Standard Error: 5.7  
 Correlation Coefficient: 0.999  
 Coefficient of Determination (Adjusted): 0.996

| ID | Level            | Concentration | Response | IS Amount | IS Response | RF            | Used |
|----|------------------|---------------|----------|-----------|-------------|---------------|------|
| 1  | IC 280-610603/18 | 0.01          | 2197.0   |           |             | 219700.0      | Y    |
| 2  | IC 280-610603/17 | 0.02          | 3453.0   |           |             | 172650.0      | Y    |
| 3  | IC 280-610603/16 | 0.05          | 8950.0   |           |             | 179000.0      | Y    |
| 4  | IC 280-610603/15 | 0.1           | 17060.0  |           |             | 170600.0      | Y    |
| 5  | IC 280-610603/14 | 0.25          | 40202.0  |           |             | 160808.0      | Y    |
| 6  | IC 280-610603/13 | 0.4           | 59892.0  |           |             | 149730.0      | Y    |
| 7  | IC 280-610603/12 | 0.7           | 114679.0 |           |             | 163827.142857 | Y    |
| 8  | IC 280-610603/11 | 1.0           | 165319.0 |           |             | 165319.0      | Y    |
| 9  | IC 280-610603/10 | 2.5           | 285146.0 |           |             | 114058.4      | N    |



## Calibration

/ RDX

**Curve Type:** Linear  
**Weighting:** Conc\_Sq  
**Origin:** None  
**Dependency:** Response  
**Calib Mode:** ESTD  
**Response Base:** AREA  
**RF Rounding:** 0


### Curve Coefficients

**Intercept:** 672.5  
**Slope:** 2.156E+05

### Error Coefficients

**Standard Error:** 10200  
**Relative Standard Error:** 4.4  
**Correlation Coefficient:** 1.000  
**Coefficient of Determination (Adjusted):** 0.998

| ID | Level            | Concentration | Response | IS Amount | IS Response | RF            | Used |
|----|------------------|---------------|----------|-----------|-------------|---------------|------|
| 1  | IC 280-610603/18 | 0.01          | 2852.0   |           |             | 285200.0      | Y    |
| 2  | IC 280-610603/17 | 0.02          | 4863.0   |           |             | 243150.0      | Y    |
| 3  | IC 280-610603/16 | 0.05          | 11621.0  |           |             | 232420.0      | Y    |
| 4  | IC 280-610603/15 | 0.1           | 22827.0  |           |             | 228270.0      | Y    |
| 5  | IC 280-610603/14 | 0.25          | 53643.0  |           |             | 214572.0      | Y    |
| 6  | IC 280-610603/13 | 0.4           | 79113.0  |           |             | 197782.5      | Y    |
| 7  | IC 280-610603/12 | 0.7           | 153380.0 |           |             | 219114.285714 | Y    |
| 8  | IC 280-610603/11 | 1.0           | 221588.0 |           |             | 221588.0      | Y    |
| 9  | IC 280-610603/10 | 2.5           | 564953.0 |           |             | 225981.2      | Y    |



# Calibration

/ Nitrobenzene

Curve Type: Average  
Weighting: Conc\_Sq  
Origin: Force  
Dependency: Response  
Calib Mode: ESTD  
Response Base: AREA  
RF Rounding: 0


## Curve Coefficients

Intercept: 0  
Slope: 3.988E+05

## Error Coefficients

Standard Error: 6800  
Relative Standard Error: 5.3  
Correlation Coefficient: 1.000  
Coefficient of Determination (Adjusted): 0.996

| ID | Level            | Concentration | Response  | IS Amount | IS Response | RF       | Used |
|----|------------------|---------------|-----------|-----------|-------------|----------|------|
| 1  | IC 280-610603/18 | 0.01          | 4269.0    |           |             | 426900.0 | Y    |
| 2  | IC 280-610603/17 | 0.02          | 7871.0    |           |             | 393550.0 | Y    |
| 3  | IC 280-610603/16 | 0.05          | 20883.0   |           |             | 417660.0 | Y    |
| 4  | IC 280-610603/15 | 0.1           | 41585.0   |           |             | 415850.0 | Y    |
| 5  | IC 280-610603/14 | 0.25          | 96681.0   |           |             | 386724.0 | Y    |
| 6  | IC 280-610603/13 | 0.4           | 141679.0  |           |             | 354197.5 | Y    |
| 7  | IC 280-610603/12 | 0.7           | 276577.0  |           |             | 395110.0 | Y    |
| 8  | IC 280-610603/11 | 1.0           | 398454.0  |           |             | 398454.0 | Y    |
| 9  | IC 280-610603/10 | 2.5           | 1002641.0 |           |             | 401056.4 | Y    |





# Calibration

/ 1,2-Dinitrobenzene

Curve Type: Average  
 Weighting: Conc\_Sq  
 Origin: Force  
 Dependency: Response  
 Calib Mode: ESTD  
 Response Base: AREA  
 RF Rounding: 0


## Curve Coefficients

Intercept: 0  
 Slope: 2.807E+05

## Error Coefficients

Standard Error: 6660  
 Relative Standard Error: 7.2  
 Correlation Coefficient: 1.000  
 Coefficient of Determination (Adjusted): 0.993

| ID | Level            | Concentration | Response | IS Amount | IS Response | RF            | Used |
|----|------------------|---------------|----------|-----------|-------------|---------------|------|
| 1  | IC 280-610603/18 | 0.01          | 3169.0   |           |             | 316900.0      | Y    |
| 2  | IC 280-610603/17 | 0.02          | 6038.0   |           |             | 301900.0      | Y    |
| 3  | IC 280-610603/16 | 0.05          | 14379.0  |           |             | 287580.0      | Y    |
| 4  | IC 280-610603/15 | 0.1           | 27893.0  |           |             | 278930.0      | Y    |
| 5  | IC 280-610603/14 | 0.25          | 67346.0  |           |             | 269384.0      | Y    |
| 6  | IC 280-610603/13 | 0.4           | 98279.0  |           |             | 245697.5      | Y    |
| 7  | IC 280-610603/12 | 0.7           | 191173.0 |           |             | 273104.285714 | Y    |
| 8  | IC 280-610603/11 | 1.0           | 276103.0 |           |             | 276103.0      | Y    |
| 9  | IC 280-610603/10 | 2.5           | 691686.0 |           |             | 276674.4      | Y    |



# Calibration

/ 3,5-Dinitroaniline

Curve Type: Linear  
 Weighting: Conc\_Sq  
 Origin: None  
 Dependency: Response  
 Calib Mode: ESTD  
 Response Base: AREA  
 RF Rounding: 0


## Curve Coefficients

Intercept: 607.3  
 Slope: 4.446E+05

## Error Coefficients

Standard Error: 36600  
 Relative Standard Error: 5.5  
 Correlation Coefficient: 0.997  
 Coefficient of Determination (Adjusted): 0.997

| ID | Level            | Concentration | Response  | IS Amount | IS Response | RF            | Used |
|----|------------------|---------------|-----------|-----------|-------------|---------------|------|
| 1  | IC 280-610603/18 | 0.01          | 5053.0    |           |             | 505300.0      | Y    |
| 2  | IC 280-610603/17 | 0.02          | 9221.0    |           |             | 461050.0      | Y    |
| 3  | IC 280-610603/16 | 0.05          | 24176.0   |           |             | 483520.0      | Y    |
| 4  | IC 280-610603/15 | 0.1           | 46980.0   |           |             | 469800.0      | Y    |
| 5  | IC 280-610603/14 | 0.25          | 113054.0  |           |             | 452216.0      | Y    |
| 6  | IC 280-610603/13 | 0.4           | 165924.0  |           |             | 414810.0      | Y    |
| 7  | IC 280-610603/12 | 0.7           | 322983.0  |           |             | 461404.285714 | Y    |
| 8  | IC 280-610603/11 | 1.0           | 460771.0  |           |             | 460771.0      | Y    |
| 9  | IC 280-610603/10 | 2.5           | 1018070.0 |           |             | 407228.0      | Y    |



## Calibration

/ 1,3-Dinitrobenzene

Curve Type: Average  
Weighting: Conc\_Sq  
Origin: Force  
Dependency: Response  
Calib Mode: ESTD  
Response Base: AREA  
RF Rounding: 0


## Curve Coefficients

Intercept: 0  
Slope: 6.294E+05

## Error Coefficients

Standard Error: 52200  
Relative Standard Error: 6.6  
Correlation Coefficient: 0.998  
Coefficient of Determination (Adjusted): 0.994

| ID | Level            | Concentration | Response  | IS Amount | IS Response | RF            | Used |
|----|------------------|---------------|-----------|-----------|-------------|---------------|------|
| 1  | IC 280-610603/18 | 0.01          | 6946.0    |           |             | 694600.0      | Y    |
| 2  | IC 280-610603/17 | 0.02          | 12307.0   |           |             | 615350.0      | Y    |
| 3  | IC 280-610603/16 | 0.05          | 31459.0   |           |             | 629180.0      | Y    |
| 4  | IC 280-610603/15 | 0.1           | 62596.0   |           |             | 625960.0      | Y    |
| 5  | IC 280-610603/14 | 0.25          | 152658.0  |           |             | 610632.0      | Y    |
| 6  | IC 280-610603/13 | 0.4           | 221798.0  |           |             | 554495.0      | Y    |
| 7  | IC 280-610603/12 | 0.7           | 431248.0  |           |             | 616068.571429 | Y    |
| 8  | IC 280-610603/11 | 1.0           | 630835.0  |           |             | 630835.0      | Y    |
| 9  | IC 280-610603/10 | 2.5           | 1717605.0 |           |             | 687042.0      | Y    |



## Calibration

/ Nitroglycerin

Curve Type: Average  
Weighting: Conc\_Sq  
Origin: Force  
Dependency: Response  
Calib Mode: ESTD  
Response Base: AREA  
RF Rounding: 0


## Curve Coefficients

Intercept: 0  
Slope: 1.347E+05

## Error Coefficients

Standard Error: 21900  
Relative Standard Error: 5.4  
Correlation Coefficient: 1.000  
Coefficient of Determination (Adjusted): 0.996

| ID | Level            | Concentration | Response  | IS Amount | IS Response | RF            | Used |
|----|------------------|---------------|-----------|-----------|-------------|---------------|------|
| 1  | IC 280-610603/18 | 0.1           | 14603.0   |           |             | 146030.0      | Y    |
| 2  | IC 280-610603/17 | 0.2           | 27409.0   |           |             | 137045.0      | Y    |
| 3  | IC 280-610603/16 | 0.5           | 71084.0   |           |             | 142168.0      | Y    |
| 4  | IC 280-610603/15 | 1.0           | 134661.0  |           |             | 134661.0      | Y    |
| 5  | IC 280-610603/14 | 2.5           | 327195.0  |           |             | 130878.0      | Y    |
| 6  | IC 280-610603/13 | 4.0           | 480119.0  |           |             | 120029.75     | Y    |
| 7  | IC 280-610603/12 | 7.0           | 929634.0  |           |             | 132804.857143 | Y    |
| 8  | IC 280-610603/11 | 10.0          | 1346659.0 |           |             | 134665.9      | Y    |
| 9  | IC 280-610603/10 | 25.0          | 3358347.0 |           |             | 134333.88     | Y    |



# Calibration

/ o-Nitrotoluene

Curve Type: Average  
 Weighting: Conc\_Sq  
 Origin: Force  
 Dependency: Response  
 Calib Mode: ESTD  
 Response Base: AREA  
 RF Rounding: 0


## Curve Coefficients

Intercept: 0  
 Slope: 2.47E+05

## Error Coefficients

Standard Error: 3910  
 Relative Standard Error: 4.8  
 Correlation Coefficient: 1.000  
 Coefficient of Determination (Adjusted): 0.997

| ID | Level            | Concentration | Response | IS Amount | IS Response | RF            | Used |
|----|------------------|---------------|----------|-----------|-------------|---------------|------|
| 1  | IC 280-610603/18 | 0.01          | 2336.0   |           |             | 233600.0      | Y    |
| 2  | IC 280-610603/17 | 0.02          | 5281.0   |           |             | 264050.0      | Y    |
| 3  | IC 280-610603/16 | 0.05          | 12951.0  |           |             | 259020.0      | Y    |
| 4  | IC 280-610603/15 | 0.1           | 25290.0  |           |             | 252900.0      | Y    |
| 5  | IC 280-610603/14 | 0.25          | 60674.0  |           |             | 242696.0      | Y    |
| 6  | IC 280-610603/13 | 0.4           | 90266.0  |           |             | 225665.0      | Y    |
| 7  | IC 280-610603/12 | 0.7           | 174353.0 |           |             | 249075.714286 | Y    |
| 8  | IC 280-610603/11 | 1.0           | 250892.0 |           |             | 250892.0      | Y    |
| 9  | IC 280-610603/10 | 2.5           | 611929.0 |           |             | 244771.6      | Y    |



## Calibration

/ p-Nitrotoluene

Curve Type: Average  
Weighting: Conc\_Sq  
Origin: Force  
Dependency: Response  
Calib Mode: ESTD  
Response Base: AREA  
RF Rounding: 0


## Curve Coefficients

Intercept: 0  
Slope: 2.275E+05

## Error Coefficients

Standard Error: 12700  
Relative Standard Error: 5.2  
Correlation Coefficient: 0.999  
Coefficient of Determination (Adjusted): 0.997

| ID | Level            | Concentration | Response | IS Amount | IS Response | RF            | Used |
|----|------------------|---------------|----------|-----------|-------------|---------------|------|
| 1  | IC 280-610603/18 | 0.01          | 2206.0   |           |             | 220600.0      | Y    |
| 2  | IC 280-610603/17 | 0.02          | 4466.0   |           |             | 223300.0      | Y    |
| 3  | IC 280-610603/16 | 0.05          | 12179.0  |           |             | 243580.0      | Y    |
| 4  | IC 280-610603/15 | 0.1           | 23002.0  |           |             | 230020.0      | Y    |
| 5  | IC 280-610603/14 | 0.25          | 56555.0  |           |             | 226220.0      | Y    |
| 6  | IC 280-610603/13 | 0.4           | 81428.0  |           |             | 203570.0      | Y    |
| 7  | IC 280-610603/12 | 0.7           | 159520.0 |           |             | 227885.714286 | Y    |
| 8  | IC 280-610603/11 | 1.0           | 230951.0 |           |             | 230951.0      | Y    |
| 9  | IC 280-610603/10 | 2.5           | 603215.0 |           |             | 241286.0      | Y    |



## Calibration

/ 4-Amino-2,6-dinitrotoluene

Curve Type: Average  
Weighting: Conc\_Sq  
Origin: Force  
Dependency: Response  
Calib Mode: ESTD  
Response Base: AREA  
RF Rounding: 0


## Curve Coefficients

Intercept: 0  
Slope: 2.942E+05

## Error Coefficients

Standard Error: 6740  
Relative Standard Error: 6.1  
Correlation Coefficient: 1.000  
Coefficient of Determination (Adjusted): 0.995

| ID | Level            | Concentration | Response | IS Amount | IS Response | RF            | Used |
|----|------------------|---------------|----------|-----------|-------------|---------------|------|
| 1  | IC 280-610603/18 | 0.01          | 3275.0   |           |             | 327500.0      | Y    |
| 2  | IC 280-610603/17 | 0.02          | 6064.0   |           |             | 303200.0      | Y    |
| 3  | IC 280-610603/16 | 0.05          | 15265.0  |           |             | 305300.0      | Y    |
| 4  | IC 280-610603/15 | 0.1           | 29472.0  |           |             | 294720.0      | Y    |
| 5  | IC 280-610603/14 | 0.25          | 71678.0  |           |             | 286712.0      | Y    |
| 6  | IC 280-610603/13 | 0.4           | 104292.0 |           |             | 260730.0      | Y    |
| 7  | IC 280-610603/12 | 0.7           | 201729.0 |           |             | 288184.285714 | Y    |
| 8  | IC 280-610603/11 | 1.0           | 292369.0 |           |             | 292369.0      | Y    |
| 9  | IC 280-610603/10 | 2.5           | 722940.0 |           |             | 289176.0      | Y    |



# Calibration

/ m-Nitrotoluene

Curve Type: Average  
Weighting: Conc\_Sq  
Origin: Force  
Dependency: Response  
Calib Mode: ESTD  
Response Base: AREA  
RF Rounding: 0


## Curve Coefficients

Intercept: 0  
Slope: 2.946E+05

## Error Coefficients

Standard Error: 6260  
Relative Standard Error: 6.1  
Correlation Coefficient: 1.000  
Coefficient of Determination (Adjusted): 0.995

| ID | Level            | Concentration | Response | IS Amount | IS Response | RF       | Used |
|----|------------------|---------------|----------|-----------|-------------|----------|------|
| 1  | IC 280-610603/18 | 0.01          | 3269.0   |           |             | 326900.0 | Y    |
| 2  | IC 280-610603/17 | 0.02          | 5823.0   |           |             | 291150.0 | Y    |
| 3  | IC 280-610603/16 | 0.05          | 15546.0  |           |             | 310920.0 | Y    |
| 4  | IC 280-610603/15 | 0.1           | 29569.0  |           |             | 295690.0 | Y    |
| 5  | IC 280-610603/14 | 0.25          | 71041.0  |           |             | 284164.0 | Y    |
| 6  | IC 280-610603/13 | 0.4           | 104338.0 |           |             | 260845.0 | Y    |
| 7  | IC 280-610603/12 | 0.7           | 203644.0 |           |             | 290920.0 | Y    |
| 8  | IC 280-610603/11 | 1.0           | 292277.0 |           |             | 292277.0 | Y    |
| 9  | IC 280-610603/10 | 2.5           | 747102.0 |           |             | 298840.8 | Y    |





## Calibration

/ 2-Amino-4,6-dinitrotoluene

Curve Type: Average  
Weighting: Conc\_Sq  
Origin: Force  
Dependency: Response  
Calib Mode: ESTD  
Response Base: AREA  
RF Rounding: 0


## Curve Coefficients

Intercept: 0  
Slope: 4.279E+05

## Error Coefficients

Standard Error: 34900  
Relative Standard Error: 6.1  
Correlation Coefficient: 0.998  
Coefficient of Determination (Adjusted): 0.995

| ID | Level            | Concentration | Response  | IS Amount | IS Response | RF            | Used |
|----|------------------|---------------|-----------|-----------|-------------|---------------|------|
| 1  | IC 280-610603/18 | 0.01          | 3912.0    |           |             | 391200.0      | Y    |
| 2  | IC 280-610603/17 | 0.02          | 8257.0    |           |             | 412850.0      | Y    |
| 3  | IC 280-610603/16 | 0.05          | 23162.0   |           |             | 463240.0      | Y    |
| 4  | IC 280-610603/15 | 0.1           | 43948.0   |           |             | 439480.0      | Y    |
| 5  | IC 280-610603/14 | 0.25          | 108505.0  |           |             | 434020.0      | Y    |
| 6  | IC 280-610603/13 | 0.4           | 158919.0  |           |             | 397297.5      | Y    |
| 7  | IC 280-610603/12 | 0.7           | 294250.0  |           |             | 420357.142857 | Y    |
| 8  | IC 280-610603/11 | 1.0           | 426004.0  |           |             | 426004.0      | Y    |
| 9  | IC 280-610603/10 | 2.5           | 1167517.0 |           |             | 467006.8      | Y    |



# Calibration

/ 1,3,5-Trinitrobenzene

Curve Type: Average  
Weighting: Conc\_Sq  
Origin: Force  
Dependency: Response  
Calib Mode: ESTD  
Response Base: AREA  
RF Rounding: 0


## Curve Coefficients

Intercept: 0  
Slope: 4.057E+05

## Error Coefficients

Standard Error: 25500  
Relative Standard Error: 7.2  
Correlation Coefficient: 0.998  
Coefficient of Determination (Adjusted): 0.993

| ID | Level            | Concentration | Response | IS Amount | IS Response | RF            | Used |
|----|------------------|---------------|----------|-----------|-------------|---------------|------|
| 1  | IC 280-610603/18 | 0.01          | 4438.0   |           |             | 443800.0      | Y    |
| 2  | IC 280-610603/17 | 0.02          | 7799.0   |           |             | 389950.0      | Y    |
| 3  | IC 280-610603/16 | 0.05          | 22059.0  |           |             | 441180.0      | Y    |
| 4  | IC 280-610603/15 | 0.1           | 41946.0  |           |             | 419460.0      | Y    |
| 5  | IC 280-610603/14 | 0.25          | 98168.0  |           |             | 392672.0      | Y    |
| 6  | IC 280-610603/13 | 0.4           | 141795.0 |           |             | 354487.5      | Y    |
| 7  | IC 280-610603/12 | 0.7           | 289774.0 |           |             | 413962.857143 | Y    |
| 8  | IC 280-610603/11 | 1.0           | 417113.0 |           |             | 417113.0      | Y    |
| 9  | IC 280-610603/10 | 2.5           | 946516.0 |           |             | 378606.4      | Y    |



## Calibration

/ 2,6-Dinitrotoluene

Curve Type: Average  
Weighting: Conc\_Sq  
Origin: Force  
Dependency: Response  
Calib Mode: ESTD  
Response Base: AREA  
RF Rounding: 0


## Curve Coefficients

Intercept: 0  
Slope: 2.917E+05

## Error Coefficients

Standard Error: 5150  
Relative Standard Error: 5.9  
Correlation Coefficient: 1.000  
Coefficient of Determination (Adjusted): 0.995

| ID | Level            | Concentration | Response | IS Amount | IS Response | RF            | Used |
|----|------------------|---------------|----------|-----------|-------------|---------------|------|
| 1  | IC 280-610603/18 | 0.01          | 3093.0   |           |             | 309300.0      | Y    |
| 2  | IC 280-610603/17 | 0.02          | 5716.0   |           |             | 285800.0      | Y    |
| 3  | IC 280-610603/16 | 0.05          | 15876.0  |           |             | 317520.0      | Y    |
| 4  | IC 280-610603/15 | 0.1           | 30260.0  |           |             | 302600.0      | Y    |
| 5  | IC 280-610603/14 | 0.25          | 70868.0  |           |             | 283472.0      | Y    |
| 6  | IC 280-610603/13 | 0.4           | 103148.0 |           |             | 257870.0      | Y    |
| 7  | IC 280-610603/12 | 0.7           | 201554.0 |           |             | 287934.285714 | Y    |
| 8  | IC 280-610603/11 | 1.0           | 290390.0 |           |             | 290390.0      | Y    |
| 9  | IC 280-610603/10 | 2.5           | 725476.0 |           |             | 290190.4      | Y    |



# Calibration

/ 2,4-Dinitrotoluene

Curve Type: Average  
 Weighting: Conc\_Sq  
 Origin: Force  
 Dependency: Response  
 Calib Mode: ESTD  
 Response Base: AREA  
 RF Rounding: 0


## Curve Coefficients

Intercept: 0  
 Slope: 5.845E+05

## Error Coefficients

Standard Error: 10100  
 Relative Standard Error: 4.7  
 Correlation Coefficient: 1.000  
 Coefficient of Determination (Adjusted): 0.997

| ID | Level            | Concentration | Response  | IS Amount | IS Response | RF            | Used |
|----|------------------|---------------|-----------|-----------|-------------|---------------|------|
| 1  | IC 280-610603/18 | 0.01          | 5987.0    |           |             | 598700.0      | Y    |
| 2  | IC 280-610603/17 | 0.02          | 12241.0   |           |             | 612050.0      | Y    |
| 3  | IC 280-610603/16 | 0.05          | 30407.0   |           |             | 608140.0      | Y    |
| 4  | IC 280-610603/15 | 0.1           | 59433.0   |           |             | 594330.0      | Y    |
| 5  | IC 280-610603/14 | 0.25          | 143098.0  |           |             | 572392.0      | Y    |
| 6  | IC 280-610603/13 | 0.4           | 207937.0  |           |             | 519842.5      | Y    |
| 7  | IC 280-610603/12 | 0.7           | 406447.0  |           |             | 580638.571429 | Y    |
| 8  | IC 280-610603/11 | 1.0           | 585209.0  |           |             | 585209.0      | Y    |
| 9  | IC 280-610603/10 | 2.5           | 1472164.0 |           |             | 588865.6      | Y    |



## Calibration

/ Tetryl

Curve Type: Average  
 Weighting: Conc\_Sq  
 Origin: Force  
 Dependency: Response  
 Calib Mode: ESTD  
 Response Base: AREA  
 RF Rounding: 0


### Curve Coefficients

Intercept: 0  
 Slope: 2.847E+05

### Error Coefficients

Standard Error: 6470  
 Relative Standard Error: 8.6  
 Correlation Coefficient: 1.000  
 Coefficient of Determination (Adjusted): 0.990

| ID | Level            | Concentration | Response | IS Amount | IS Response | RF            | Used |
|----|------------------|---------------|----------|-----------|-------------|---------------|------|
| 1  | IC 280-610603/18 | 0.01          | 3282.0   |           |             | 328200.0      | Y    |
| 2  | IC 280-610603/17 | 0.02          | 6348.0   |           |             | 317400.0      | Y    |
| 3  | IC 280-610603/16 | 0.05          | 14233.0  |           |             | 284660.0      | Y    |
| 4  | IC 280-610603/15 | 0.1           | 27534.0  |           |             | 275340.0      | Y    |
| 5  | IC 280-610603/14 | 0.25          | 68327.0  |           |             | 273308.0      | Y    |
| 6  | IC 280-610603/13 | 0.4           | 98615.0  |           |             | 246537.5      | Y    |
| 7  | IC 280-610603/12 | 0.7           | 192352.0 |           |             | 274788.571429 | Y    |
| 8  | IC 280-610603/11 | 1.0           | 278494.0 |           |             | 278494.0      | Y    |
| 9  | IC 280-610603/10 | 2.5           | 709626.0 |           |             | 283850.4      | Y    |



# Calibration

/ 2,4,6-Trinitrotoluene

Curve Type: Average  
 Weighting: Conc\_Sq  
 Origin: Force  
 Dependency: Response  
 Calib Mode: ESTD  
 Response Base: AREA  
 RF Rounding: 0


## Curve Coefficients

Intercept: 0  
 Slope: 4.389E+05

## Error Coefficients

Standard Error: 10200  
 Relative Standard Error: 8.2  
 Correlation Coefficient: 1.000  
 Coefficient of Determination (Adjusted): 0.991

| ID | Level            | Concentration | Response  | IS Amount | IS Response | RF            | Used |
|----|------------------|---------------|-----------|-----------|-------------|---------------|------|
| 1  | IC 280-610603/18 | 0.01          | 4919.0    |           |             | 491900.0      | Y    |
| 2  | IC 280-610603/17 | 0.02          | 9925.0    |           |             | 496250.0      | Y    |
| 3  | IC 280-610603/16 | 0.05          | 22184.0   |           |             | 443680.0      | Y    |
| 4  | IC 280-610603/15 | 0.1           | 43136.0   |           |             | 431360.0      | Y    |
| 5  | IC 280-610603/14 | 0.25          | 104561.0  |           |             | 418244.0      | Y    |
| 6  | IC 280-610603/13 | 0.4           | 152188.0  |           |             | 380470.0      | Y    |
| 7  | IC 280-610603/12 | 0.7           | 294801.0  |           |             | 421144.285714 | Y    |
| 8  | IC 280-610603/11 | 1.0           | 430103.0  |           |             | 430103.0      | Y    |
| 9  | IC 280-610603/10 | 2.5           | 1091822.0 |           |             | 436728.8      | Y    |



# Calibration

/ PETN

Curve Type: Average  
Weighting: Conc\_Sq  
Origin: Force  
Dependency: Response  
Calib Mode: ESTD  
Response Base: AREA  
RF Rounding: 0


## Curve Coefficients

Intercept: 0  
Slope: 1.439E+05

## Error Coefficients

Standard Error: 25700  
Relative Standard Error: 6.3  
Correlation Coefficient: 1.000  
Coefficient of Determination (Adjusted): 0.995

| ID | Level            | Concentration | Response  | IS Amount | IS Response | RF            | Used |
|----|------------------|---------------|-----------|-----------|-------------|---------------|------|
| 1  | IC 280-610603/18 | 0.1           | 16344.0   |           |             | 163440.0      | Y    |
| 2  | IC 280-610603/17 | 0.2           | 28805.0   |           |             | 144025.0      | Y    |
| 3  | IC 280-610603/16 | 0.5           | 72943.0   |           |             | 145886.0      | Y    |
| 4  | IC 280-610603/15 | 1.0           | 146450.0  |           |             | 146450.0      | Y    |
| 5  | IC 280-610603/14 | 2.5           | 351913.0  |           |             | 140765.2      | Y    |
| 6  | IC 280-610603/13 | 4.0           | 512590.0  |           |             | 128147.5      | Y    |
| 7  | IC 280-610603/12 | 7.0           | 992519.0  |           |             | 141788.428571 | Y    |
| 8  | IC 280-610603/11 | 10.0          | 1419730.0 |           |             | 141973.0      | Y    |
| 9  | IC 280-610603/10 | 25.0          | 3574120.0 |           |             | 142964.8      | Y    |



FORM VI  
HPLC/IC BY EXTERNAL STANDARD - INITIAL CALIBRATION DATA  
RETENTION TIME SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1 Analy Batch No.: 610603

SDG No.: \_\_\_\_\_

Instrument ID: CHHPLC\_X5 GC Column: Luna-phenyl ID: 4.6 (mm) Heated Purge: (Y/N) N

Calibration Start Date: 04/29/2023 00:13 Calibration End Date: 04/29/2023 04:18 Calibration ID: 79726

Calibration Files

|         |                  |              |
|---------|------------------|--------------|
| LEVEL:  | LAB SAMPLE ID:   | LAB FILE ID: |
| Level 1 | IC 280-610603/27 | 04280027.D   |
| Level 2 | IC 280-610603/26 | 04280026.D   |
| Level 3 | IC 280-610603/25 | 04280025.D   |
| Level 4 | IC 280-610603/24 | 04280024.D   |
| Level 5 | IC 280-610603/23 | 04280023.D   |
| Level 6 | IC 280-610603/22 | 04280022.D   |
| Level 7 | IC 280-610603/21 | 04280021.D   |
| Level 8 | IC 280-610603/20 | 04280020.D   |

| ANALYTE                    | LVL 1 | LVL 2 | LVL 3 | LVL 4 | LVL 5 | LVL 6 | LVL 7 | LVL 8 |  |  | RT WINDOW     | AVG RT |
|----------------------------|-------|-------|-------|-------|-------|-------|-------|-------|--|--|---------------|--------|
| 2,6-diamino-4-nitrotoluene | +++++ | +++++ | 4.112 | 4.109 | 4.113 | 4.111 | 4.105 | 4.080 |  |  | 3.963 - 4.263 | 4.105  |
| 2,4-diamino-6-nitrotoluene | +++++ | +++++ | 4.652 | 4.643 | 4.653 | 4.651 | 4.638 | 4.620 |  |  | 4.503 - 4.803 | 4.643  |



FORM VI  
HPLC/IC BY EXTERNAL STANDARD - INITIAL CALIBRATION DATA  
CURVE EVALUATION

Lab Name: Eurofins Denver Job No.: 280-176674-1 Analy Batch No.: 610603  
 SDG No.: \_\_\_\_\_  
 Instrument ID: CHHPLC\_X5 GC Column: Luna-phenyl ID: 4.6(mm) Heated Purge: (Y/N) N  
 Calibration Start Date: 04/29/2023 00:13 Calibration End Date: 04/29/2023 04:18 Calibration ID: 79726

Calibration Files

|         |                  |              |
|---------|------------------|--------------|
| LEVEL:  | LAB SAMPLE ID:   | LAB FILE ID: |
| Level 1 | IC 280-610603/27 | 04280027.D   |
| Level 2 | IC 280-610603/26 | 04280026.D   |
| Level 3 | IC 280-610603/25 | 04280025.D   |
| Level 4 | IC 280-610603/24 | 04280024.D   |
| Level 5 | IC 280-610603/23 | 04280023.D   |
| Level 6 | IC 280-610603/22 | 04280022.D   |
| Level 7 | IC 280-610603/21 | 04280021.D   |
| Level 8 | IC 280-610603/20 | 04280020.D   |

| ANALYTE                    | CF              |                 |                  |                  | CURVE<br>TYPE | COEFFICIENT |                |    | # | MIN CF | %RSD<br>/RSE | # | MAX<br>%RSD<br>/RSE | R^2<br>OR COD | # | MIN R^2<br>OR COD |
|----------------------------|-----------------|-----------------|------------------|------------------|---------------|-------------|----------------|----|---|--------|--------------|---|---------------------|---------------|---|-------------------|
|                            | LVL 1<br>LVL 5  | LVL 2<br>LVL 6  | LVL 3<br>LVL 7   | LVL 4<br>LVL 8   |               | B           | M1             | M2 |   |        |              |   |                     |               |   |                   |
| 2,6-diamino-4-nitrotoluene | +++++<br>404133 | +++++<br>405346 | 414510<br>402965 | 401736<br>391111 | Ave           |             | 403300.00<br>2 |    |   |        | 1.9          |   | 20.0                |               |   |                   |
| 2,4-diamino-6-nitrotoluene | +++++<br>219290 | +++++<br>221550 | 230940<br>222678 | 231676<br>225549 | Ave           |             | 225280.53<br>3 |    |   |        | 2.3          |   | 20.0                |               |   |                   |

Note: The M1 coefficient is the same as Ave CF for an Ave curve type. RSD is calculated for Ave curve types. RSE is used for all other types.

FORM VI  
HPLC/IC BY EXTERNAL STANDARD - INITIAL CALIBRATION DATA  
RESPONSE AND CONCENTRATION

Lab Name: Eurofins Denver Job No.: 280-176674-1 Analy Batch No.: 610603

SDG No.: \_\_\_\_\_

Instrument ID: CHHPLC\_X5 GC Column: Luna-phenyl ID: 4.6 (mm) Heated Purge: (Y/N) N

Calibration Start Date: 04/29/2023 00:13 Calibration End Date: 04/29/2023 04:18 Calibration ID: 79726

Calibration Files

|         |                  |              |
|---------|------------------|--------------|
| LEVEL:  | LAB SAMPLE ID:   | LAB FILE ID: |
| Level 1 | IC 280-610603/27 | 04280027.D   |
| Level 2 | IC 280-610603/26 | 04280026.D   |
| Level 3 | IC 280-610603/25 | 04280025.D   |
| Level 4 | IC 280-610603/24 | 04280024.D   |
| Level 5 | IC 280-610603/23 | 04280023.D   |
| Level 6 | IC 280-610603/22 | 04280022.D   |
| Level 7 | IC 280-610603/21 | 04280021.D   |
| Level 8 | IC 280-610603/20 | 04280020.D   |

| ANALYTE                    | CURVE TYPE | RESPONSE        |                 |                 |        |        | CONCENTRATION (UG/ML) |                |                |       |       |
|----------------------------|------------|-----------------|-----------------|-----------------|--------|--------|-----------------------|----------------|----------------|-------|-------|
|                            |            | LVL 1<br>LVL 6  | LVL 2<br>LVL 7  | LVL 3<br>LVL 8  | LVL 4  | LVL 5  | LVL 1<br>LVL 6        | LVL 2<br>LVL 7 | LVL 3<br>LVL 8 | LVL 4 | LVL 5 |
| 2,6-diamino-4-nitrotoluene | Ave        | +++++<br>283742 | +++++<br>402965 | 41451<br>977777 | 100434 | 161653 | +++++<br>0.700        | +++++<br>1.00  | 0.100<br>2.50  | 0.250 | 0.400 |
| 2,4-diamino-6-nitrotoluene | Ave        | +++++<br>155085 | +++++<br>222678 | 23094<br>563873 | 57919  | 87716  | +++++<br>0.700        | +++++<br>1.00  | 0.100<br>2.50  | 0.250 | 0.400 |

Curve Type Legend

Ave = Average

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280020.D  
 Lims ID: IC ADD 8  
 Client ID:  
 Sample Type: IC Calib Level: 8  
 Inject. Date: 29-Apr-2023 00:13:27 ALS Bottle#: 20 Worklist Smp#: 20  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: IC ADD 8  
 Operator ID: JZ/JG Instrument ID: CHHPLC\_X5  
 Sublist: chrom-8330\_X5\_Luna\*sub10  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\8330\_X5\_Luna.m  
 Limit Group: GCSV - 8330  
 Last Update: 29-Apr-2023 13:04:21 Calib Date: 29-Apr-2023 04:18:01  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280027.D  
 Column 1 : Luna-Phenyl hexyl ( 4.60 mm) Det: LC mwd1A, 254 nm  
 Process Host: CTX1621

First Level Reviewer: LV5D

Date: 29-Apr-2023 10:25:00

| Compound                     | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 1 2,6-diamino-4-nitrotoluene | 1   | 4.080        | 4.113            | -0.033           | 977777   | 2.50             | 2.42               |       |
| 2 2,4-diamino-6-nitrotoluene | 1   | 4.620        | 4.653            | -0.033           | 563873   | 2.50             | 2.50               |       |

**QC Flag Legend**

Processing Flags

**Reagents:**

8330\_ADDs\_00036

Amount Added: 125.00

Units: uL

Report Date: 29-Apr-2023 13:04:21

Chrom Revision: 2.3 29-Mar-2023 18:39:10

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280020.D

Injection Date: 29-Apr-2023 00:13:27

Instrument ID: CHHPLC\_X5

Operator ID: JZ/JG

Lims ID: IC ADD 8

Worklist Smp#: 20

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 20

Method: 8330\_X5\_Luna

Limit Group: GCSV - 8330


Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280021.D  
Lims ID: IC ADD 7  
Client ID:  
Sample Type: IC Calib Level: 7  
Inject. Date: 29-Apr-2023 00:48:26 ALS Bottle#: 21 Worklist Smp#: 21  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Sample Info: IC ADD 7  
Operator ID: JZ/JG Instrument ID: CHHPLC\_X5  
Sublist: chrom-8330\_X5\_Luna\*sub10  
Method: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\8330\_X5\_Luna.m  
Limit Group: GCSV - 8330  
Last Update: 29-Apr-2023 13:04:21 Calib Date: 29-Apr-2023 04:18:01  
Integrator: Falcon  
Quant Method: External Standard Quant By: Initial Calibration  
Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280027.D  
Column 1 : Luna-Phenyl hexyl ( 4.60 mm) Det: LC mwd1A, 254 nm  
Process Host: CTX1621

| Compound                     | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 1 2,6-diamino-4-nitrotoluene | 1   | 4.105        | 4.113            | -0.008           | 402965   | 1.00             | 1.00               |       |
| 2 2,4-diamino-6-nitrotoluene | 1   | 4.638        | 4.653            | -0.015           | 222678   | 1.00             | 0.9884             |       |

## Reagents:

8330\_ADDs\_00036

Amount Added: 50.00

Units: uL

Report Date: 29-Apr-2023 13:04:21

Chrom Revision: 2.3 29-Mar-2023 18:39:10

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280021.D

Injection Date: 29-Apr-2023 00:48:26

Instrument ID: CHHPLC\_X5

Operator ID: JZ/JG

Lims ID: IC ADD 7

Worklist Smp#: 21

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 21

Method: 8330\_X5\_Luna

Limit Group: GCSV - 8330


Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280022.D  
 Lims ID: IC ADD 6  
 Client ID:  
 Sample Type: IC Calib Level: 6  
 Inject. Date: 29-Apr-2023 01:23:20 ALS Bottle#: 22 Worklist Smp#: 22  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: IC ADD 6  
 Operator ID: JZ/JG Instrument ID: CHHPLC\_X5  
 Sublist: chrom-8330\_X5\_Luna\*sub10  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\8330\_X5\_Luna.m  
 Limit Group: GCSV - 8330  
 Last Update: 29-Apr-2023 13:04:22 Calib Date: 29-Apr-2023 04:18:01  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280027.D  
 Column 1 : Luna-Phenyl hexyl ( 4.60 mm) Det: LC mwd1A, 254 nm  
 Process Host: CTX1621

| Compound                     | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 1 2,6-diamino-4-nitrotoluene | 1   | 4.111        | 4.113            | -0.002           | 283742   | 0.7000           | 0.7036             |       |
| 2 2,4-diamino-6-nitrotoluene | 1   | 4.651        | 4.653            | -0.002           | 155085   | 0.7000           | 0.6884             |       |

**Reagents:**

8330\_ADDs\_00036

Amount Added: 35.00

Units: uL

Report Date: 29-Apr-2023 13:04:22

Chrom Revision: 2.3 29-Mar-2023 18:39:10

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280022.D

Injection Date: 29-Apr-2023 01:23:20

Instrument ID: CHHPLC\_X5

Operator ID: JZ/JG

Lims ID: IC ADD 6

Worklist Smp#: 22

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 22

Method: 8330\_X5\_Luna

Limit Group: GCSV - 8330


Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1





Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280023.D  
 Lims ID: IC ADD 5  
 Client ID:  
 Sample Type: IC Calib Level: 5  
 Inject. Date: 29-Apr-2023 01:58:17 ALS Bottle#: 23 Worklist Smp#: 23  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: IC ADD 5  
 Operator ID: JZ/JG Instrument ID: CHHPLC\_X5  
 Sublist: chrom-8330\_X5\_Luna\*sub10  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\8330\_X5\_Luna.m  
 Limit Group: GCSV - 8330  
 Last Update: 29-Apr-2023 13:04:22 Calib Date: 29-Apr-2023 04:18:01  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280027.D  
 Column 1 : Luna-Phenyl hexyl ( 4.60 mm) Det: LC mwd1A, 254 nm  
 Process Host: CTX1621

| Compound                     | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 1 2,6-diamino-4-nitrotoluene | 1   | 4.113        | 4.113            | 0.000            | 161653   | 0.4000           | 0.4008             |       |
| 2 2,4-diamino-6-nitrotoluene | 1   | 4.653        | 4.653            | 0.000            | 87716    | 0.4000           | 0.3894             |       |

**Reagents:**

8330\_ADDs\_00036

Amount Added: 20.00

Units: uL

Report Date: 29-Apr-2023 13:04:22

Chrom Revision: 2.3 29-Mar-2023 18:39:10

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280023.D

Injection Date: 29-Apr-2023 01:58:17

Instrument ID: CHHPLC\_X5

Operator ID: JZ/JG

Lims ID: IC ADD 5

Worklist Smp#: 23

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 23

Method: 8330\_X5\_Luna

Limit Group: GCSV - 8330


Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280024.D  
Lims ID: IC ADD 4  
Client ID:  
Sample Type: IC Calib Level: 4  
Inject. Date: 29-Apr-2023 02:33:12 ALS Bottle#: 24 Worklist Smp#: 24  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Sample Info: IC ADD 4  
Operator ID: JZ/JG Instrument ID: CHHPLC\_X5  
Sublist: chrom-8330\_X5\_Luna\*sub10  
Method: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\8330\_X5\_Luna.m  
Limit Group: GCSV - 8330  
Last Update: 29-Apr-2023 13:04:22 Calib Date: 29-Apr-2023 04:18:01  
Integrator: Falcon  
Quant Method: External Standard Quant By: Initial Calibration  
Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280027.D  
Column 1 : Luna-Phenyl hexyl ( 4.60 mm) Det: LC mwd1A, 254 nm  
Process Host: CTX1621

| Compound                     | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 1 2,6-diamino-4-nitrotoluene | 1   | 4.109        | 4.113            | -0.004           | 100434   | 0.2500           | 0.2490             |       |
| 2 2,4-diamino-6-nitrotoluene | 1   | 4.643        | 4.653            | -0.010           | 57919    | 0.2500           | 0.2571             |       |

## Reagents:

8330\_ADDs\_00036

Amount Added: 12.50

Units: uL

Report Date: 29-Apr-2023 13:04:22

Chrom Revision: 2.3 29-Mar-2023 18:39:10

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280024.D

Injection Date: 29-Apr-2023 02:33:12

Instrument ID: CHHPLC\_X5

Operator ID: JZ/JG

Lims ID: IC ADD 4

Worklist Smp#: 24

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 24

Method: 8330\_X5\_Luna

Limit Group: GCSV - 8330


Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280025.D  
Lims ID: IC ADD 3  
Client ID:  
Sample Type: IC Calib Level: 3  
Inject. Date: 29-Apr-2023 03:08:08 ALS Bottle#: 25 Worklist Smp#: 25  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Sample Info: IC ADD 3  
Operator ID: JZ/JG Instrument ID: CHHPLC\_X5  
Sublist: chrom-8330\_X5\_Luna\*sub10  
Method: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\8330\_X5\_Luna.m  
Limit Group: GCSV - 8330  
Last Update: 29-Apr-2023 13:04:23 Calib Date: 29-Apr-2023 04:18:01  
Integrator: Falcon  
Quant Method: External Standard Quant By: Initial Calibration  
Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280027.D  
Column 1 : Luna-Phenyl hexyl ( 4.60 mm) Det: LC mwd1A, 254 nm  
Process Host: CTX1621

| Compound                     | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 1 2,6-diamino-4-nitrotoluene | 1   | 4.112        | 4.113            | -0.001           | 41451    | 0.1000           | 0.1028             |       |
| 2 2,4-diamino-6-nitrotoluene | 1   | 4.652        | 4.653            | -0.001           | 23094    | 0.1000           | 0.1025             |       |

## Reagents:

8330\_ADDs\_00036

Amount Added: 5.00

Units: uL

Report Date: 29-Apr-2023 13:04:23

Chrom Revision: 2.3 29-Mar-2023 18:39:10

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280025.D

Injection Date: 29-Apr-2023 03:08:08

Instrument ID: CHHPLC\_X5

Operator ID: JZ/JG

Lims ID: IC ADD 3

Worklist Smp#: 25

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 25

Method: 8330\_X5\_Luna

Limit Group: GCSV - 8330


Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280026.D  
Lims ID: IC ADD 2  
Client ID:  
Sample Type: IC Calib Level: 2  
Inject. Date: 29-Apr-2023 03:43:02 ALS Bottle#: 26 Worklist Smp#: 26  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Sample Info: IC ADD 2  
Operator ID: JZ/JG Instrument ID: CHHPLC\_X5  
Sublist: chrom-8330\_X5\_Luna\*sub10  
Method: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\8330\_X5\_Luna.m  
Limit Group: GCSV - 8330  
Last Update: 29-Apr-2023 13:04:23 Calib Date: 29-Apr-2023 04:18:01  
Integrator: Falcon  
Quant Method: External Standard Quant By: Initial Calibration  
Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280027.D  
Column 1 : Luna-Phenyl hexyl ( 4.60 mm) Det: LC mwd1A, 254 nm  
Process Host: CTX1621

| Compound                     | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 1 2,6-diamino-4-nitrotoluene | 1   | 4.113        | 4.113            | 0.000            | 21691    | 0.0400           | 0.0538             |       |
| 2 2,4-diamino-6-nitrotoluene | 1   | 4.653        | 4.653            | 0.000            | 12626    | 0.0400           | 0.0560             |       |

## Reagents:

8330\_ADDs\_00036

Amount Added: 2.00

Units: uL

Report Date: 29-Apr-2023 13:04:23

Chrom Revision: 2.3 29-Mar-2023 18:39:10

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280026.D

Injection Date: 29-Apr-2023 03:43:02

Instrument ID: CHHPLC\_X5

Operator ID: JZ/JG

Lims ID: IC ADD 2

Worklist Smp#: 26

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 26

Method: 8330\_X5\_Luna

Limit Group: GCSV - 8330


Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1





Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280027.D  
 Lims ID: IC ADD 1  
 Client ID:  
 Sample Type: IC Calib Level: 1  
 Inject. Date: 29-Apr-2023 04:18:01 ALS Bottle#: 27 Worklist Smp#: 27  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: IC ADD 1  
 Operator ID: JZ/JG Instrument ID: CHHPLC\_X5  
 Sublist: chrom-8330\_X5\_Luna\*sub10  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\8330\_X5\_Luna.m  
 Limit Group: GCSV - 8330  
 Last Update: 29-Apr-2023 13:04:23 Calib Date: 29-Apr-2023 04:18:01  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280027.D  
 Column 1 : Luna-Phenyl hexyl ( 4.60 mm) Det: LC mwd1A, 254 nm  
 Process Host: CTX1621

| Compound                     | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 1 2,6-diamino-4-nitrotoluene | 1   | 4.113        | 4.113            | 0.000            | 8084     | 0.0200           | 0.0200             |       |
| 2 2,4-diamino-6-nitrotoluene | 1   | 4.653        | 4.653            | 0.000            | 5388     | 0.0200           | 0.0239             |       |

**Reagents:**

8330\_ADDs\_00036

Amount Added: 1.00

Units: uL

Report Date: 29-Apr-2023 13:04:24

Chrom Revision: 2.3 29-Mar-2023 18:39:10

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280027.D

Injection Date: 29-Apr-2023 04:18:01

Instrument ID: CHHPLC\_X5

Operator ID: JZ/JG

Lims ID: IC ADD 1

Worklist Smp#: 27

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 27

Method: 8330\_X5\_Luna

Limit Group: GCSV - 8330


Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



# Calibration

/ 2,6-diamino-4-nitrotoluene

Curve Type: Average  
Weighting: Conc\_Sq  
Origin: Force  
Dependency: Response  
Calib Mode: ESTD  
Response Base: AREA  
RF Rounding: 0


## Curve Coefficients

Intercept: 0  
Slope: 4.033E+05

## Error Coefficients

Standard Error: 13700  
Relative Standard Error: 1.9  
Correlation Coefficient: 1.000  
Coefficient of Determination (Adjusted): 0.999

| ID | Level            | Concentration | Response | IS Amount | IS Response | RF            | Used |
|----|------------------|---------------|----------|-----------|-------------|---------------|------|
| 1  | IC 280-610603/27 | 0.02          | 8084.0   |           |             | 404200.0      | N    |
| 2  | IC 280-610603/26 | 0.04          | 21691.0  |           |             | 542275.0      | N    |
| 3  | IC 280-610603/25 | 0.1           | 41451.0  |           |             | 414510.0      | Y    |
| 4  | IC 280-610603/24 | 0.25          | 100434.0 |           |             | 401736.0      | Y    |
| 5  | IC 280-610603/23 | 0.4           | 161653.0 |           |             | 404132.5      | Y    |
| 6  | IC 280-610603/22 | 0.7           | 283742.0 |           |             | 405345.714286 | Y    |
| 7  | IC 280-610603/21 | 1.0           | 402965.0 |           |             | 402965.0      | Y    |
| 8  | IC 280-610603/20 | 2.5           | 977777.0 |           |             | 391110.8      | Y    |



## Calibration

/ 2,4-diamino-6-nitrotoluene

Curve Type: Average  
Weighting: Conc\_Sq  
Origin: Force  
Dependency: Response  
Calib Mode: ESTD  
Response Base: AREA  
RF Rounding: 0


## Curve Coefficients

Intercept: 0  
Slope: 2.253E+05

## Error Coefficients

Standard Error: 2130  
Relative Standard Error: 2.3  
Correlation Coefficient: 1.000  
Coefficient of Determination (Adjusted): 0.999

| ID | Level            | Concentration | Response | IS Amount | IS Response | RF       | Used |
|----|------------------|---------------|----------|-----------|-------------|----------|------|
| 1  | IC 280-610603/27 | 0.02          | 5388.0   |           |             | 269400.0 | N    |
| 2  | IC 280-610603/26 | 0.04          | 12626.0  |           |             | 315650.0 | N    |
| 3  | IC 280-610603/25 | 0.1           | 23094.0  |           |             | 230940.0 | Y    |
| 4  | IC 280-610603/24 | 0.25          | 57919.0  |           |             | 231676.0 | Y    |
| 5  | IC 280-610603/23 | 0.4           | 87716.0  |           |             | 219290.0 | Y    |
| 6  | IC 280-610603/22 | 0.7           | 155085.0 |           |             | 221550.0 | Y    |
| 7  | IC 280-610603/21 | 1.0           | 222678.0 |           |             | 222678.0 | Y    |
| 8  | IC 280-610603/20 | 2.5           | 563873.0 |           |             | 225549.2 | Y    |



FORM VII  
HPLC/IC CONTINUING CALIBRATION DATA

Lab Name: Eurofins Denver Job No.: 280-176674-1  
SDG No.: \_\_\_\_\_  
Lab Sample ID: ICV 280-599427/19 Calibration Date: 01/17/2023 03:40  
Instrument ID: CHHPLC\_X3 Calib Start Date: 01/17/2023 00:33  
GC Column: UltraCarb5uODS ID: 4.60 (mm) Calib End Date: 01/17/2023 03:16  
Lab File ID: 01160047.D Conc. Units: ug/L

| ANALYTE                    | CURVE<br>TYPE | AVE CF | CF     | MIN CF | CALC<br>AMOUNT | SPIKE<br>AMOUNT | %D   | MAX<br>%D |
|----------------------------|---------------|--------|--------|--------|----------------|-----------------|------|-----------|
| 2,6-diamino-4-nitrotoluene | Lin2          |        | 224434 |        | 485            | 500             | -3.0 | 20.0      |
| 2,4-diamino-6-nitrotoluene | Ave           | 146472 | 141938 |        | 485            | 500             | -3.1 | 20.0      |
| 3,5-Dinitroaniline         | Lin2          |        | 235360 |        | 481            | 500             | -3.8 | 20.0      |

FORM VII  
HPLC/IC CONTINUING CALIBRATION RETENTION TIME SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
SDG No.: \_\_\_\_\_  
Lab Sample ID: ICV 280-599427/19 Calibration Date: 01/17/2023 03:40  
Instrument ID: CHHPLC\_X3 Calib Start Date: 01/17/2023 00:33  
GC Column: UltraCarb5uODS ID: 4.60 (mm) Calib End Date: 01/17/2023 03:16  
Lab File ID: 01160047.D

| Analyte                    | RT   | RT WINDOW |      |
|----------------------------|------|-----------|------|
|                            |      | FROM      | TO   |
| 2,6-diamino-4-nitrotoluene | 6.45 | 6.31      | 6.61 |
| 2,4-diamino-6-nitrotoluene | 6.63 | 6.62      | 6.92 |
| 3,5-Dinitroaniline         | 9.80 | 9.65      | 9.95 |

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230114-117775.b\01160047.D  
 Lims ID: ICV ADD  
 Client ID:  
 Sample Type: ICV  
 Inject. Date: 17-Jan-2023 03:40:08 ALS Bottle#: 50 Worklist Smp#: 19  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: ICVADD  
 Operator ID: JZ/MAR Instrument ID: CHHPLC\_X3  
 Sublist:  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230114-117775.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 17-Jan-2023 12:54:36 Calib Date: 17-Jan-2023 03:16:44  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230114-117775.b\01160046.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1672

First Level Reviewer: LV5D

Date: 17-Jan-2023 12:32:06

| Compound                     | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/mL | OnCol Amt<br>ug/mL | Flags |
|------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 2 2,6-diamino-4-nitrotoluene | 1   | 6.454        | 6.455            | -0.001           | 112217   | 0.5000           | 0.4848             | M     |
| 5 2,4-diamino-6-nitrotoluene | 1   | 6.634        | 6.769            | -0.135           | 70969    | 0.5000           | 0.4845             | M     |
| 14 3,5-Dinitroaniline        | 1   | 9.801        | 9.800            | 0.001            | 117680   | 0.5000           | 0.4811             |       |

**QC Flag Legend**

Processing Flags

Review Flags

M - Manually Integrated

**Reagents:**

8330DiaminLCS\_00052

Amount Added: 50.00

Units: uL

3,5-DNA LCS\_00042

Amount Added: 50.00

Units: uL

Report Date: 17-Jan-2023 12:54:42

Chrom Revision: 2.3 20-Dec-2022 14:14:06

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230114-117775.b\01160047.d

Injection Date: 17-Jan-2023 03:40:08

Instrument ID: CHHPLC\_X3

Operator ID: JZ/MAR

Lims ID: ICV ADD

Worklist Smp#: 19

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 50

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1





## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230114-117775.b\01160047.d  
Injection Date: 17-Jan-2023 03:40:08 Instrument ID: CHHPLC\_X3  
Lims ID: ICV ADD  
Client ID:  
Operator ID: JZ/MAR ALS Bottle#: 50 Worklist Smp#: 19  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**2,2,6-diamino-4-nitrotoluene, CAS: 59229-75-3**

Signal: 1


RT: 6.45  
Area: 126714  
Amount: 0.529427  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.45  
Area: 112217  
Amount: 0.484790  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 17-Jan-2023 12:32:03  
Audit Action: Split an Integrated Peak

Audit Reason: Baseline

## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230114-117775.b\01160047.d  
Injection Date: 17-Jan-2023 03:40:08 Instrument ID: CHHPLC\_X3  
Lims ID: ICV ADD  
Client ID:  
Operator ID: JZ/MAR ALS Bottle#: 50 Worklist Smp#: 19  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**5 2,4-diamino-6-nitrotoluene, CAS: 6629-29-4**

Signal: 1


RT: 6.63  
Area: 84952  
Amount: 0.579987  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.63  
Area: 70969  
Amount: 0.484522  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 17-Jan-2023 12:32:04  
Audit Action: Split an Integrated Peak

Audit Reason: Baseline

FORM VII  
HPLC/IC CONTINUING CALIBRATION DATA

Lab Name: Eurofins Denver Job No.: 280-176674-1  
 SDG No.: \_\_\_\_\_  
 Lab Sample ID: ICV 280-601664/20 Calibration Date: 02/08/2023 19:05  
 Instrument ID: CHHPLC\_X3 Calib Start Date: 02/08/2023 15:38  
 GC Column: UltraCarb5uODS ID: 4.60 (mm) Calib End Date: 02/08/2023 18:42  
 Lab File ID: 02080020.D Conc. Units: ug/L

| ANALYTE                    | CURVE<br>TYPE | AVE CF | CF     | MIN CF | CALC<br>AMOUNT | SPIKE<br>AMOUNT | %D   | MAX<br>%D |
|----------------------------|---------------|--------|--------|--------|----------------|-----------------|------|-----------|
| HMX                        | Ave           | 93528  | 88148  |        | 471            | 500             | -5.8 | 20.0      |
| RDX                        | Ave           | 106380 | 105642 |        | 497            | 500             | -0.7 | 20.0      |
| Picric acid                | Ave           | 75830  | 83408  |        | 550            | 500             | 10.0 | 20.0      |
| 1,3,5-Trinitrobenzene      | Ave           | 217147 | 243558 |        | 561            | 500             | 12.2 | 20.0      |
| 1,3-Dinitrobenzene         | Ave           | 294397 | 317476 |        | 539            | 500             | 7.8  | 20.0      |
| Nitrobenzene               | Ave           | 191245 | 205682 |        | 538            | 500             | 7.5  | 20.0      |
| 3,5-Dinitroaniline         | Lin2          |        | 223846 |        | 489            | 500             | -2.2 | 20.0      |
| Tetryl                     | Ave           | 164121 | 182414 |        | 556            | 500             | 11.1 | 20.0      |
| Nitroglycerin              | Ave           | 64070  | 67923  |        | 5300           | 5000            | 6.0  | 20.0      |
| 2,4,6-Trinitrotoluene      | Ave           | 211040 | 213428 |        | 506            | 500             | 1.1  | 20.0      |
| 4-Amino-2,6-dinitrotoluene | Ave           | 154933 | 158824 |        | 513            | 500             | 2.5  | 20.0      |
| 2-Amino-4,6-dinitrotoluene | Ave           | 201410 | 204888 |        | 509            | 500             | 1.7  | 20.0      |
| 2,6-Dinitrotoluene         | Ave           | 142745 | 147024 |        | 515            | 500             | 3.0  | 20.0      |
| 2,4-Dinitrotoluene         | Ave           | 296667 | 301684 |        | 508            | 500             | 1.7  | 20.0      |
| 2-Nitrotoluene             | Ave           | 127896 | 128960 |        | 504            | 500             | 0.8  | 20.0      |
| 4-Nitrotoluene             | Ave           | 111880 | 111646 |        | 499            | 500             | -0.2 | 20.0      |
| 3-Nitrotoluene             | Ave           | 140492 | 138998 |        | 495            | 500             | -1.1 | 20.0      |
| PETN                       | Ave           | 68845  | 76520  |        | 5560           | 5000            | 11.1 | 20.0      |
| 1,2-Dinitrobenzene         | Ave           | 126309 | 125412 |        | 496            | 500             | -0.7 | 20.0      |

FORM VII  
HPLC/IC CONTINUING CALIBRATION RETENTION TIME SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
 SDG No.: \_\_\_\_\_  
 Lab Sample ID: ICV 280-601664/20 Calibration Date: 02/08/2023 19:05  
 Instrument ID: CHHPLC\_X3 Calib Start Date: 02/08/2023 15:38  
 GC Column: UltraCarb5uODS ID: 4.60 (mm) Calib End Date: 02/08/2023 18:42  
 Lab File ID: 02080020.D

| Analyte                    | RT    | RT WINDOW |       |
|----------------------------|-------|-----------|-------|
|                            |       | FROM      | TO    |
| HMX                        | 6.55  | 6.39      | 6.69  |
| RDX                        | 7.56  | 7.41      | 7.71  |
| Picric acid                | 8.02  | 7.89      | 8.19  |
| 1,3,5-Trinitrobenzene      | 8.65  | 8.49      | 8.79  |
| 1,3-Dinitrobenzene         | 9.26  | 9.11      | 9.41  |
| Nitrobenzene               | 9.64  | 9.48      | 9.78  |
| 3,5-Dinitroaniline         | 9.84  | 9.68      | 9.98  |
| Tetryl                     | 9.99  | 9.83      | 10.13 |
| Nitroglycerin              | 10.46 | 10.29     | 10.59 |
| 2,4,6-Trinitrotoluene      | 10.88 | 10.76     | 10.96 |
| 4-Amino-2,6-dinitrotoluene | 11.05 | 10.93     | 11.13 |
| 2-Amino-4,6-dinitrotoluene | 11.30 | 11.18     | 11.38 |
| 2,6-Dinitrotoluene         | 11.48 | 11.35     | 11.55 |
| 2,4-Dinitrotoluene         | 11.64 | 11.52     | 11.72 |
| 2-Nitrotoluene             | 12.46 | 12.29     | 12.59 |
| 4-Nitrotoluene             | 12.88 | 12.69     | 12.99 |
| 3-Nitrotoluene             | 13.44 | 13.26     | 13.56 |
| PETN                       | 14.60 | 14.40     | 14.70 |
| 1,2-Dinitrobenzene         | 8.52  | 8.37      | 8.67  |

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230208-118465.b\02080020.D  
 Lims ID: ICV INT  
 Client ID:  
 Sample Type: ICV  
 Inject. Date: 08-Feb-2023 19:05:10 ALS Bottle#: 20 Worklist Smp#: 20  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: ICV INT  
 Operator ID: JZ/MAR Instrument ID: CHHPLC\_X3  
 Sublist:  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230208-118465.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 09-Feb-2023 13:03:49 Calib Date: 08-Feb-2023 18:42:13  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230208-118465.b\02080019.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1604

First Level Reviewer: LV5D

Date: 08-Feb-2023 19:31:36

| Compound                      | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/mL | OnCol Amt<br>ug/mL | Flags |
|-------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 4 HMX                         | 1   | 6.549        | 6.544            | 0.005            | 44074    | 0.5000           | 0.4712             | M     |
| 8 RDX                         | 1   | 7.563        | 7.564            | -0.001           | 52821    | 0.5000           | 0.4965             | M     |
| 9 2,4,6-Trinitrophenol        | 1   | 8.023        | 8.044            | -0.021           | 41704    | 0.5000           | 0.5500             |       |
| \$ 10 1,2-Dinitrobenzene      | 1   | 8.523        | 8.517            | 0.006            | 62706    | 0.5000           | 0.4965             |       |
| 11 1,3,5-Trinitrobenzene      | 1   | 8.649        | 8.644            | 0.005            | 121779   | 0.5000           | 0.5608             |       |
| 12 1,3-Dinitrobenzene         | 1   | 9.263        | 9.257            | 0.006            | 158738   | 0.5000           | 0.5392             |       |
| 13 Nitrobenzene               | 1   | 9.636        | 9.631            | 0.005            | 102841   | 0.5000           | 0.5377             |       |
| 14 3,5-Dinitroaniline         | 1   | 9.843        | 9.831            | 0.012            | 111923   | 0.5000           | 0.4890             |       |
| 15 Tetryl                     | 1   | 9.989        | 9.977            | 0.012            | 91207    | 0.5000           | 0.5557             |       |
| 16 Nitroglycerin              | 2   | 10.456       | 10.437           | 0.019            | 339614   | 5.00             | 5.30               |       |
| 17 2,4,6-Trinitrotoluene      | 1   | 10.883       | 10.864           | 0.019            | 106714   | 0.5000           | 0.5057             |       |
| 18 4-Amino-2,6-dinitrotoluene | 1   | 11.049       | 11.031           | 0.018            | 79412    | 0.5000           | 0.5126             |       |
| 19 2-Amino-4,6-dinitrotoluene | 1   | 11.296       | 11.277           | 0.019            | 102444   | 0.5000           | 0.5086             |       |
| 20 2,6-Dinitrotoluene         | 1   | 11.476       | 11.451           | 0.025            | 73512    | 0.5000           | 0.5150             |       |
| 21 2,4-Dinitrotoluene         | 1   | 11.643       | 11.617           | 0.026            | 150842   | 0.5000           | 0.5085             |       |
| 22 o-Nitrotoluene             | 1   | 12.463       | 12.437           | 0.026            | 64480    | 0.5000           | 0.5042             |       |
| 23 p-Nitrotoluene             | 1   | 12.876       | 12.844           | 0.032            | 55823    | 0.5000           | 0.4990             |       |
| 24 m-Nitrotoluene             | 1   | 13.443       | 13.411           | 0.032            | 69499    | 0.5000           | 0.4947             |       |
| 25 PETN                       | 2   | 14.596       | 14.551           | 0.045            | 382602   | 5.00             | 5.56               |       |

**QC Flag Legend**

Processing Flags

Review Flags

M - Manually Integrated

**Reagents:**

8330Surrogate\_00138

Amount Added: 50.00

Units: uL

8330 LCS\_00121

Amount Added: 50.00

Units: uL

3,5-DNA LCS\_00043

Amount Added: 50.00

Units: uL

Report Date: 09-Feb-2023 13:04:01

Chrom Revision: 2.3 01-Feb-2023 13:23:06

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230208-118465.b\02080020.d

Injection Date: 08-Feb-2023 19:05:10

Instrument ID: CHHPLC\_X3

Operator ID: JZ/MAR

Lims ID: ICV INT

Worklist Smp#: 20

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 20

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230208-118465.b\02080020.d  
Injection Date: 08-Feb-2023 19:05:10 Instrument ID: CHHPLC\_X3  
Lims ID: ICV INT  
Client ID:  
Operator ID: JZ/MAR ALS Bottle#: 20 Worklist Smp#: 20  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**4 HMX, CAS: 2691-41-0**

Signal: 1


RT: 6.55  
Area: 49010  
Amount: 0.524013  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.55  
Area: 44074  
Amount: 0.471238  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 08-Feb-2023 19:31:27  
Audit Action: Manually Integrated

Audit Reason: Baseline Smoothing



## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230208-118465.b\02080020.d  
Injection Date: 08-Feb-2023 19:05:10 Instrument ID: CHHPLC\_X3  
Lims ID: ICV INT  
Client ID:  
Operator ID: JZ/MAR ALS Bottle#: 20 Worklist Smp#: 20  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**8 RDX, CAS: 121-82-4**

Signal: 1


RT: 7.56  
Area: 54996  
Amount: 0.516978  
Amount Units: ug/mL

## Processing Integration Results



RT: 7.56  
Area: 52821  
Amount: 0.496532  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 08-Feb-2023 19:32:23  
Audit Action: Manually Integrated

Audit Reason: Baseline Smoothing

FORM VII  
HPLC/IC CONTINUING CALIBRATION DATA

Lab Name: Eurofins Denver Job No.: 280-176674-1  
SDG No.: \_\_\_\_\_  
Lab Sample ID: ICV 280-603284/19 Calibration Date: 02/24/2023 18:04  
Instrument ID: CHHPLC\_X3 Calib Start Date: 02/24/2023 15:00  
GC Column: UltraCarb5uODS ID: 4.60 (mm) Calib End Date: 02/24/2023 17:41  
Lab File ID: 02240019.D Conc. Units: ug/L

| ANALYTE | CURVE<br>TYPE | AVE CF | CF     | MIN CF | CALC<br>AMOUNT | SPIKE<br>AMOUNT | %D   | MAX<br>%D |
|---------|---------------|--------|--------|--------|----------------|-----------------|------|-----------|
| TNX     | Ave           | 197687 | 197288 |        | 499            | 500             | -0.2 | 20.0      |
| DNX     | Ave           | 144462 | 141211 |        | 491            | 502             | -2.3 | 20.0      |
| MNX     | Ave           | 131172 | 134903 |        | 602            | 585             | 2.8  | 20.0      |

FORM VII  
HPLC/IC CONTINUING CALIBRATION RETENTION TIME SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
SDG No.: \_\_\_\_\_  
Lab Sample ID: ICV 280-603284/19 Calibration Date: 02/24/2023 18:04  
Instrument ID: CHHPLC\_X3 Calib Start Date: 02/24/2023 15:00  
GC Column: UltraCarb5uODS ID: 4.60 (mm) Calib End Date: 02/24/2023 17:41  
Lab File ID: 02240019.D

| Analyte | RT   | RT WINDOW |      |
|---------|------|-----------|------|
|         |      | FROM      | TO   |
| TNX     | 6.45 | 6.34      | 6.54 |
| DNX     | 6.77 | 6.67      | 6.87 |
| MNX     | 7.20 | 7.05      | 7.35 |

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240019.D  
 Lims ID: ICV DMT  
 Client ID:  
 Sample Type: ICV  
 Inject. Date: 24-Feb-2023 18:04:21 ALS Bottle#: 19 Worklist Smp#: 19  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: ICV DMT  
 Operator ID: JZ/MAR Instrument ID: CHHPLC\_X3  
 Sublist:  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 28-Feb-2023 13:33:34 Calib Date: 24-Feb-2023 17:41:20  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1602

First Level Reviewer: LV5D

Date: 28-Feb-2023 11:45:18

| Compound | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/mL | OnCol Amt<br>ug/mL | Flags |
|----------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 3 TNX    | 1   | 6.446        | 6.443            | 0.003            | 98644    | 0.5000           | 0.4990             | M     |
| 6 DNX    | 1   | 6.772        | 6.770            | 0.002            | 70888    | 0.5020           | 0.4907             | M     |
| 7 MNX    | 1   | 7.199        | 7.196            | 0.003            | 78918    | 0.5850           | 0.6016             |       |

**QC Flag Legend**

Processing Flags

Review Flags

M - Manually Integrated

**Reagents:**

8330\_OP\_DMT\_00016

Amount Added: 50.00

Units: uL

Report Date: 28-Feb-2023 13:33:59

Chrom Revision: 2.3 15-Feb-2023 20:44:50

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230224-118952.b\02240019.d

Injection Date: 24-Feb-2023 18:04:21

Instrument ID: CHHPLC\_X3

Operator ID: JZ/MAR

Lims ID: ICV DMT

Worklist Smp#: 19

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 19

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230224-118952.b\02240019.d  
Injection Date: 24-Feb-2023 18:04:21 Instrument ID: CHHPLC\_X3  
Lims ID: ICV DMT  
Client ID:  
Operator ID: JZ/MAR ALS Bottle#: 19 Worklist Smp#: 19  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**3 TNX, CAS: 13980-04-6**

Signal: 1


RT: 6.45  
Area: 120613  
Amount: 0.610120  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.45  
Area: 98644  
Amount: 0.498990  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 28-Feb-2023 11:45:12  
Audit Action: Manually Integrated

Audit Reason: Baseline

## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230224-118952.b\02240019.d  
Injection Date: 24-Feb-2023 18:04:21 Instrument ID: CHHPLC\_X3  
Lims ID: ICV DMT  
Client ID:  
Operator ID: JZ/MAR ALS Bottle#: 19 Worklist Smp#: 19  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**6 DNX, CAS: 80251-29-2**

Signal: 1


RT: 6.77  
Area: 83407  
Amount: 0.577361  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.77  
Area: 70888  
Amount: 0.490702  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 28-Feb-2023 11:45:16  
Audit Action: Manually Integrated

Audit Reason: Baseline

FORM VII  
HPLC/IC CONTINUING CALIBRATION DATA

Lab Name: Eurofins Denver Job No.: 280-176674-1  
 SDG No.: \_\_\_\_\_  
 Lab Sample ID: CCV 280-613168/32 Calibration Date: 05/19/2023 22:41  
 Instrument ID: CHHPLC\_X3 Calib Start Date: 02/08/2023 15:38  
 GC Column: UltraCarb5uODS ID: 4.60 (mm) Calib End Date: 02/08/2023 18:42  
 Lab File ID: 05190032.D Conc. Units: ug/L

| ANALYTE                    | CURVE<br>TYPE | AVE CF | CF     | MIN CF | CALC<br>AMOUNT | SPIKE<br>AMOUNT | %D    | MAX<br>%D |
|----------------------------|---------------|--------|--------|--------|----------------|-----------------|-------|-----------|
| HMX                        | Ave           | 93528  | 97276  |        | 260            | 250             | 4.0   | 20.0      |
| RDX                        | Ave           | 106380 | 114836 |        | 270            | 250             | 7.9   | 20.0      |
| Picric acid                | Ave           | 75830  | 83564  |        | 275            | 250             | 10.2  | 20.0      |
| 1,3,5-Trinitrobenzene      | Ave           | 217147 | 230412 |        | 265            | 250             | 6.1   | 20.0      |
| 1,3-Dinitrobenzene         | Ave           | 294397 | 314180 |        | 267            | 250             | 6.7   | 20.0      |
| Nitrobenzene               | Ave           | 191245 | 195144 |        | 255            | 250             | 2.0   | 20.0      |
| 3,5-Dinitroaniline         | Lin2          |        | 242216 |        | 265            | 250             | 5.9   | 20.0      |
| Tetryl                     | Ave           | 164121 | 147104 |        | 224            | 250             | -10.4 | 20.0      |
| Nitroglycerin              | Ave           | 64070  | 68084  |        | 2660           | 2500            | 6.3   | 20.0      |
| 2,4,6-Trinitrotoluene      | Ave           | 211040 | 215868 |        | 256            | 250             | 2.3   | 20.0      |
| 4-Amino-2,6-dinitrotoluene | Ave           | 154933 | 161716 |        | 261            | 250             | 4.4   | 20.0      |
| 2-Amino-4,6-dinitrotoluene | Ave           | 201410 | 206292 |        | 256            | 250             | 2.4   | 20.0      |
| 2,6-Dinitrotoluene         | Ave           | 142745 | 155028 |        | 272            | 250             | 8.6   | 20.0      |
| 2,4-Dinitrotoluene         | Ave           | 296667 | 312288 |        | 263            | 250             | 5.3   | 20.0      |
| 2-Nitrotoluene             | Ave           | 127896 | 126744 |        | 248            | 250             | -0.9  | 20.0      |
| 4-Nitrotoluene             | Ave           | 111880 | 108304 |        | 242            | 250             | -3.2  | 20.0      |
| 3-Nitrotoluene             | Ave           | 140492 | 138916 |        | 247            | 250             | -1.1  | 20.0      |
| PETN                       | Ave           | 68845  | 74332  |        | 2700           | 2500            | 8.0   | 20.0      |
| 1,2-Dinitrobenzene         | Ave           | 126309 | 141520 |        | 280            | 250             | 12.0  | 20.0      |



FORM VII  
HPLC/IC CONTINUING CALIBRATION RETENTION TIME SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
 SDG No.: \_\_\_\_\_  
 Lab Sample ID: CCV 280-613168/32 Calibration Date: 05/19/2023 22:41  
 Instrument ID: CHHPLC\_X3 Calib Start Date: 02/08/2023 15:38  
 GC Column: UltraCarb5uODS ID: 4.60 (mm) Calib End Date: 02/08/2023 18:42  
 Lab File ID: 05190032.D

| Analyte                    | RT    | RT WINDOW |       |
|----------------------------|-------|-----------|-------|
|                            |       | FROM      | TO    |
| HMX                        | 6.55  | 6.40      | 6.70  |
| RDX                        | 7.57  | 7.42      | 7.72  |
| Picric acid                | 7.97  | 7.81      | 8.11  |
| 1,3,5-Trinitrobenzene      | 8.64  | 8.50      | 8.80  |
| 1,3-Dinitrobenzene         | 9.27  | 9.13      | 9.43  |
| Nitrobenzene               | 9.64  | 9.49      | 9.79  |
| 3,5-Dinitroaniline         | 9.88  | 9.73      | 10.03 |
| Tetryl                     | 10.00 | 9.87      | 10.17 |
| Nitroglycerin              | 10.46 | 10.32     | 10.62 |
| 2,4,6-Trinitrotoluene      | 10.89 | 10.81     | 11.01 |
| 4-Amino-2,6-dinitrotoluene | 11.09 | 11.02     | 11.22 |
| 2-Amino-4,6-dinitrotoluene | 11.34 | 11.27     | 11.47 |
| 2,6-Dinitrotoluene         | 11.49 | 11.42     | 11.62 |
| 2,4-Dinitrotoluene         | 11.66 | 11.58     | 11.78 |
| 2-Nitrotoluene             | 12.48 | 12.37     | 12.67 |
| 4-Nitrotoluene             | 12.90 | 12.79     | 13.09 |
| 3-Nitrotoluene             | 13.47 | 13.37     | 13.67 |
| PETN                       | 14.63 | 14.55     | 14.85 |
| 1,2-Dinitrobenzene         | 8.52  | 8.37      | 8.67  |

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\05190032.D  
 Lims ID: CCV DMT  
 Client ID:  
 Sample Type: CCV  
 Inject. Date: 19-May-2023 22:41:21 ALS Bottle#: 7 Worklist Smp#: 32  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: CCV DMT  
 Operator ID: JZ/JG Instrument ID: CHHPLC\_X3  
 Sublist: chrom-8330\_X3\*sub27  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 20-May-2023 11:22:31 Calib Date: 24-Feb-2023 17:41:20  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1663

First Level Reviewer: LV5D

Date: 20-May-2023 10:26:09

| Compound                      | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/mL | OnCol Amt<br>ug/mL | Flags |
|-------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 3 TNX                         | 1   | 6.429        | 6.430            | -0.001           | 49929    | 0.2508           | 0.2526             | M     |
| 4 HMX                         | 1   | 6.549        | 6.550            | -0.001           | 24319    | 0.2500           | 0.2600             | M     |
| 6 DNx                         | 1   | 6.756        | 6.757            | -0.001           | 37017    | 0.2503           | 0.2562             | M     |
| 7 MNx                         | 1   | 7.183        | 7.190            | -0.007           | 40240    | 0.2918           | 0.3068             |       |
| 8 RDX                         | 1   | 7.569        | 7.570            | -0.001           | 28709    | 0.2500           | 0.2699             |       |
| 9 2,4,6-Trinitrophenol        | 1   | 7.969        | 7.963            | 0.006            | 20891    | 0.2500           | 0.2755             |       |
| \$ 10 1,2-Dinitrobenzene      | 1   | 8.523        | 8.523            | 0.000            | 35380    | 0.2500           | 0.2801             |       |
| 11 1,3,5-Trinitrobenzene      | 1   | 8.643        | 8.650            | -0.007           | 57603    | 0.2500           | 0.2653             |       |
| 12 1,3-Dinitrobenzene         | 1   | 9.269        | 9.277            | -0.008           | 78545    | 0.2500           | 0.2668             |       |
| 13 Nitrobenzene               | 1   | 9.636        | 9.643            | -0.007           | 48786    | 0.2500           | 0.2551             |       |
| 14 3,5-Dinitroaniline         | 1   | 9.876        | 9.883            | -0.007           | 60554    | 0.2500           | 0.2649             |       |
| 15 Tetryl                     | 1   | 10.003       | 10.017           | -0.014           | 36776    | 0.2500           | 0.2241             |       |
| 16 Nitroglycerin              | 2   | 10.456       | 10.470           | -0.014           | 170210   | 2.50             | 2.66               |       |
| 17 2,4,6-Trinitrotoluene      | 1   | 10.889       | 10.910           | -0.021           | 53967    | 0.2500           | 0.2557             |       |
| 18 4-Amino-2,6-dinitrotoluene | 1   | 11.089       | 11.117           | -0.028           | 40429    | 0.2500           | 0.2609             |       |
| 19 2-Amino-4,6-dinitrotoluene | 1   | 11.343       | 11.370           | -0.027           | 51573    | 0.2500           | 0.2561             |       |
| 20 2,6-Dinitrotoluene         | 1   | 11.489       | 11.517           | -0.028           | 38757    | 0.2500           | 0.2715             |       |
| 21 2,4-Dinitrotoluene         | 1   | 11.663       | 11.683           | -0.020           | 78072    | 0.2500           | 0.2632             |       |
| 22 o-Nitrotoluene             | 1   | 12.483       | 12.517           | -0.034           | 31686    | 0.2500           | 0.2477             |       |
| 23 p-Nitrotoluene             | 1   | 12.896       | 12.937           | -0.041           | 27076    | 0.2500           | 0.2420             |       |
| 24 m-Nitrotoluene             | 1   | 13.469       | 13.517           | -0.048           | 34729    | 0.2500           | 0.2472             |       |
| 25 PETN                       | 2   | 14.629       | 14.697           | -0.068           | 185830   | 2.50             | 2.70               |       |

### QC Flag Legend

Processing Flags

Review Flags

M - Manually Integrated

**Reagents:**

8330 DMT\_00013

Amount Added: 12.50

Units: uL

8330IntermStk\_00075

Amount Added: 25.00

Units: uL

Report Date: 20-May-2023 11:22:31

Chrom Revision: 2.3 16-May-2023 16:37:30

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190032.d

Injection Date: 19-May-2023 22:41:21

Instrument ID: CHHPLC\_X3

Operator ID: JZ/JG

Lims ID: CCV DMT

Worklist Smp#: 32

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 7

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190032.d

Injection Date: 19-May-2023 22:41:21

Instrument ID: CHHPLC\_X3

Lims ID: CCV DMT

Client ID:

Operator ID: JZ/JG

ALS Bottle#:

7

Worklist Smp#: 32

Injection Vol: 100.0 ul

Dil. Factor: 1.0000

Method: 8330\_X3

Limit Group: GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm)

Detector


LC DAD1B, 254 nm

**4 HMX, CAS: 2691-41-0**

Signal: 1


RT: 6.55  
Area: 28256  
Amount: 0.302112  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.55  
Area: 24319  
Amount: 0.260018  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 10:26:03 -06:00:00 (UTC)

Audit Action: Split an Integrated Peak

Audit Reason: Baseline

FORM VII  
HPLC/IC CONTINUING CALIBRATION DATA

Lab Name: Eurofins Denver Job No.: 280-176674-1  
SDG No.: \_\_\_\_\_  
Lab Sample ID: CCV 280-613168/32 Calibration Date: 05/19/2023 22:41  
Instrument ID: CHHPLC\_X3 Calib Start Date: 02/24/2023 15:00  
GC Column: UltraCarb5uODS ID: 4.60 (mm) Calib End Date: 02/24/2023 17:41  
Lab File ID: 05190032.D Conc. Units: ug/L

| ANALYTE | CURVE<br>TYPE | AVE CF | CF     | MIN CF | CALC<br>AMOUNT | SPIKE<br>AMOUNT | %D  | MAX<br>%D |
|---------|---------------|--------|--------|--------|----------------|-----------------|-----|-----------|
| TNX     | Ave           | 197687 | 199119 |        | 253            | 251             | 0.7 | 20.0      |
| DNX     | Ave           | 144462 | 147920 |        | 256            | 250             | 2.4 | 20.0      |
| MNX     | Ave           | 131172 | 137926 |        | 307            | 292             | 5.1 | 20.0      |

FORM VII  
HPLC/IC CONTINUING CALIBRATION RETENTION TIME SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
SDG No.: \_\_\_\_\_  
Lab Sample ID: CCV 280-613168/32 Calibration Date: 05/19/2023 22:41  
Instrument ID: CHHPLC\_X3 Calib Start Date: 02/24/2023 15:00  
GC Column: UltraCarb5uODS ID: 4.60 (mm) Calib End Date: 02/24/2023 17:41  
Lab File ID: 05190032.D

| Analyte | RT   | RT WINDOW |      |
|---------|------|-----------|------|
|         |      | FROM      | TO   |
| TNX     | 6.43 | 6.33      | 6.53 |
| DNX     | 6.76 | 6.66      | 6.86 |
| MNX     | 7.18 | 7.04      | 7.34 |

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\05190032.D  
 Lims ID: CCV DMT  
 Client ID:  
 Sample Type: CCV  
 Inject. Date: 19-May-2023 22:41:21 ALS Bottle#: 7 Worklist Smp#: 32  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: CCV DMT  
 Operator ID: JZ/JG Instrument ID: CHHPLC\_X3  
 Sublist: chrom-8330\_X3\*sub27  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 20-May-2023 11:22:31 Calib Date: 24-Feb-2023 17:41:20  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1663

First Level Reviewer: LV5D

Date: 20-May-2023 10:26:09

| Compound                      | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/mL | OnCol Amt<br>ug/mL | Flags |
|-------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 3 TNX                         | 1   | 6.429        | 6.430            | -0.001           | 49929    | 0.2508           | 0.2526             | M     |
| 4 HMX                         | 1   | 6.549        | 6.550            | -0.001           | 24319    | 0.2500           | 0.2600             | M     |
| 6 DNx                         | 1   | 6.756        | 6.757            | -0.001           | 37017    | 0.2503           | 0.2562             | M     |
| 7 MNx                         | 1   | 7.183        | 7.190            | -0.007           | 40240    | 0.2918           | 0.3068             |       |
| 8 RDX                         | 1   | 7.569        | 7.570            | -0.001           | 28709    | 0.2500           | 0.2699             |       |
| 9 2,4,6-Trinitrophenol        | 1   | 7.969        | 7.963            | 0.006            | 20891    | 0.2500           | 0.2755             |       |
| \$ 10 1,2-Dinitrobenzene      | 1   | 8.523        | 8.523            | 0.000            | 35380    | 0.2500           | 0.2801             |       |
| 11 1,3,5-Trinitrobenzene      | 1   | 8.643        | 8.650            | -0.007           | 57603    | 0.2500           | 0.2653             |       |
| 12 1,3-Dinitrobenzene         | 1   | 9.269        | 9.277            | -0.008           | 78545    | 0.2500           | 0.2668             |       |
| 13 Nitrobenzene               | 1   | 9.636        | 9.643            | -0.007           | 48786    | 0.2500           | 0.2551             |       |
| 14 3,5-Dinitroaniline         | 1   | 9.876        | 9.883            | -0.007           | 60554    | 0.2500           | 0.2649             |       |
| 15 Tetryl                     | 1   | 10.003       | 10.017           | -0.014           | 36776    | 0.2500           | 0.2241             |       |
| 16 Nitroglycerin              | 2   | 10.456       | 10.470           | -0.014           | 170210   | 2.50             | 2.66               |       |
| 17 2,4,6-Trinitrotoluene      | 1   | 10.889       | 10.910           | -0.021           | 53967    | 0.2500           | 0.2557             |       |
| 18 4-Amino-2,6-dinitrotoluene | 1   | 11.089       | 11.117           | -0.028           | 40429    | 0.2500           | 0.2609             |       |
| 19 2-Amino-4,6-dinitrotoluene | 1   | 11.343       | 11.370           | -0.027           | 51573    | 0.2500           | 0.2561             |       |
| 20 2,6-Dinitrotoluene         | 1   | 11.489       | 11.517           | -0.028           | 38757    | 0.2500           | 0.2715             |       |
| 21 2,4-Dinitrotoluene         | 1   | 11.663       | 11.683           | -0.020           | 78072    | 0.2500           | 0.2632             |       |
| 22 o-Nitrotoluene             | 1   | 12.483       | 12.517           | -0.034           | 31686    | 0.2500           | 0.2477             |       |
| 23 p-Nitrotoluene             | 1   | 12.896       | 12.937           | -0.041           | 27076    | 0.2500           | 0.2420             |       |
| 24 m-Nitrotoluene             | 1   | 13.469       | 13.517           | -0.048           | 34729    | 0.2500           | 0.2472             |       |
| 25 PETN                       | 2   | 14.629       | 14.697           | -0.068           | 185830   | 2.50             | 2.70               |       |

**QC Flag Legend**

Processing Flags

Review Flags

M - Manually Integrated



**Reagents:**

8330 DMT\_00013

Amount Added: 12.50

Units: uL

8330IntermStk\_00075

Amount Added: 25.00

Units: uL

Report Date: 20-May-2023 11:22:31

Chrom Revision: 2.3 16-May-2023 16:37:30

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190032.d

Injection Date: 19-May-2023 22:41:21

Instrument ID: CHHPLC\_X3

Operator ID: JZ/JG

Lims ID: CCV DMT

Worklist Smp#: 32

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 7

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190032.d  
Injection Date: 19-May-2023 22:41:21 Instrument ID: CHHPLC\_X3  
Lims ID: CCV DMT  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 7 Worklist Smp#: 32  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**3 TNX, CAS: 13980-04-6**

Signal: 1


RT: 6.43  
Area: 53608  
Amount: 0.271176  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.43  
Area: 49929  
Amount: 0.252566  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 10:26:02 -06:00:00 (UTC)

Audit Action: Split an Integrated Peak

Audit Reason: Baseline

## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190032.d

Injection Date: 19-May-2023 22:41:21

Instrument ID: CHHPLC\_X3

Lims ID: CCV DMT

Client ID:

Operator ID: JZ/JG

ALS Bottle#:

7

Worklist Smp#: 32

Injection Vol: 100.0 ul

Dil. Factor: 1.0000

Method: 8330\_X3

Limit Group: GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm)

Detector


LC DAD1B, 254 nm

**6 DNX, CAS: 80251-29-2**

Signal: 1


RT: 6.76  
Area: 40382  
Amount: 0.279533  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.76  
Area: 37017  
Amount: 0.256240  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 10:26:08 -06:00:00 (UTC)

Audit Action: Manually Integrated

Audit Reason: Baseline

FORM VII  
HPLC/IC CONTINUING CALIBRATION DATA

Lab Name: Eurofins Denver Job No.: 280-176674-1  
SDG No.: \_\_\_\_\_  
Lab Sample ID: CCV 280-613168/33 Calibration Date: 05/19/2023 23:04  
Instrument ID: CHHPLC\_X3 Calib Start Date: 01/17/2023 00:33  
GC Column: UltraCarb5uODS ID: 4.60 (mm) Calib End Date: 01/17/2023 03:16  
Lab File ID: 05190033.D Conc. Units: ug/L

| ANALYTE                    | CURVE<br>TYPE | AVE CF | CF     | MIN CF | CALC<br>AMOUNT | SPIKE<br>AMOUNT | %D   | MAX<br>%D |
|----------------------------|---------------|--------|--------|--------|----------------|-----------------|------|-----------|
| 2,6-diamino-4-nitrotoluene | Lin2          |        | 221920 |        | 238            | 250             | -4.9 | 20.0      |
| 2,4-diamino-6-nitrotoluene | Ave           | 146472 | 136996 |        | 234            | 250             | -6.5 | 20.0      |

FORM VII  
HPLC/IC CONTINUING CALIBRATION RETENTION TIME SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
SDG No.: \_\_\_\_\_  
Lab Sample ID: CCV 280-613168/33 Calibration Date: 05/19/2023 23:04  
Instrument ID: CHHPLC\_X3 Calib Start Date: 01/17/2023 00:33  
GC Column: UltraCarb5uODS ID: 4.60 (mm) Calib End Date: 01/17/2023 03:16  
Lab File ID: 05190033.D

| Analyte                    | RT   | RT WINDOW |      |
|----------------------------|------|-----------|------|
|                            |      | FROM      | TO   |
| 2,6-diamino-4-nitrotoluene | 6.43 | 6.28      | 6.58 |
| 2,4-diamino-6-nitrotoluene | 6.62 | 6.47      | 6.77 |

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\05190033.D  
 Lims ID: CCV ADD  
 Client ID:  
 Sample Type: CCV  
 Inject. Date: 19-May-2023 23:04:19 ALS Bottle#: 8 Worklist Smp#: 33  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: CCV ADD  
 Operator ID: JZ/JG Instrument ID: CHHPLC\_X3  
 Sublist: chrom-8330\_X3\*sub28  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 20-May-2023 11:22:31 Calib Date: 24-Feb-2023 17:41:20  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1663

First Level Reviewer: LV5D

Date: 20-May-2023 10:26:15

| Compound                     | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/mL | OnCol Amt<br>ug/mL | Flags |
|------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 2 2,6-diamino-4-nitrotoluene | 1   | 6.425        | 6.431            | -0.006           | 55480    | 0.2500           | 0.2378             | M     |
| 5 2,4-diamino-6-nitrotoluene | 1   | 6.618        | 6.618            | 0.000            | 34249    | 0.2500           | 0.2338             | M     |

**QC Flag Legend**

Review Flags

M - Manually Integrated

**Reagents:**

8330\_ADDs\_00037

Amount Added: 12.50

Units: uL

Report Date: 20-May-2023 11:22:31

Chrom Revision: 2.3 16-May-2023 16:37:30

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190033.d

Injection Date: 19-May-2023 23:04:19

Instrument ID: CHHPLC\_X3

Operator ID: JZ/JG

Lims ID: CCV ADD

Worklist Smp#: 33

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 8

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1





## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190033.d  
Injection Date: 19-May-2023 23:04:19 Instrument ID: CHHPLC\_X3  
Lims ID: CCV ADD  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 8 Worklist Smp#: 33  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**2,2,6-diamino-4-nitrotoluene, CAS: 59229-75-3**

Signal: 1


RT: 6.42  
Area: 56120  
Amount: 0.240564  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.42  
Area: 55480  
Amount: 0.237778  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 10:26:13 -06:00:00 (UTC)

Audit Action: Split an Integrated Peak

Audit Reason: Baseline

## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190033.d

Injection Date: 19-May-2023 23:04:19

Instrument ID: CHHPLC\_X3

Lims ID: CCV ADD

Client ID:

Operator ID: JZ/JG

ALS Bottle#:

8

Worklist Smp#: 33

Injection Vol: 100.0 ul

Dil. Factor:

1.0000

Method: 8330\_X3

Limit Group:

GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm)

Detector


LC DAD1B, 254 nm

**5 2,4-diamino-6-nitrotoluene, CAS: 6629-29-4**

Signal: 1


RT: 6.62  
Area: 36354  
Amount: 0.248197  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.62  
Area: 34249  
Amount: 0.233826  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 10:26:14 -06:00:00 (UTC)

Audit Action: Split an Integrated Peak

Audit Reason: Baseline

FORM VII  
HPLC/IC CONTINUING CALIBRATION DATA

Lab Name: Eurofins Denver Job No.: 280-176674-1  
 SDG No.: \_\_\_\_\_  
 Lab Sample ID: CCV 280-613168/44 Calibration Date: 05/20/2023 03:16  
 Instrument ID: CHHPLC\_X3 Calib Start Date: 02/08/2023 15:38  
 GC Column: UltraCarb5uODS ID: 4.60 (mm) Calib End Date: 02/08/2023 18:42  
 Lab File ID: 05190044.D Conc. Units: ug/L

| ANALYTE                    | CURVE<br>TYPE | AVE CF | CF     | MIN CF | CALC<br>AMOUNT | SPIKE<br>AMOUNT | %D   | MAX<br>%D |
|----------------------------|---------------|--------|--------|--------|----------------|-----------------|------|-----------|
| HMX                        | Ave           | 93528  | 94944  |        | 254            | 250             | 1.5  | 20.0      |
| RDX                        | Ave           | 106380 | 115160 |        | 271            | 250             | 8.3  | 20.0      |
| Picric acid                | Ave           | 75830  | 84904  |        | 280            | 250             | 12.0 | 20.0      |
| 1,3,5-Trinitrobenzene      | Ave           | 217147 | 225564 |        | 260            | 250             | 3.9  | 20.0      |
| 1,3-Dinitrobenzene         | Ave           | 294397 | 315568 |        | 268            | 250             | 7.2  | 20.0      |
| Nitrobenzene               | Ave           | 191245 | 191912 |        | 251            | 250             | 0.3  | 20.0      |
| 3,5-Dinitroaniline         | Lin2          |        | 235100 |        | 257            | 250             | 2.8  | 20.0      |
| Tetryl                     | Ave           | 164121 | 149828 |        | 228            | 250             | -8.7 | 20.0      |
| Nitroglycerin              | Ave           | 64070  | 68442  |        | 2670           | 2500            | 6.8  | 20.0      |
| 2,4,6-Trinitrotoluene      | Ave           | 211040 | 217880 |        | 258            | 250             | 3.2  | 20.0      |
| 4-Amino-2,6-dinitrotoluene | Ave           | 154933 | 161584 |        | 261            | 250             | 4.3  | 20.0      |
| 2-Amino-4,6-dinitrotoluene | Ave           | 201410 | 205692 |        | 255            | 250             | 2.1  | 20.0      |
| 2,6-Dinitrotoluene         | Ave           | 142745 | 156192 |        | 274            | 250             | 9.4  | 20.0      |
| 2,4-Dinitrotoluene         | Ave           | 296667 | 314000 |        | 265            | 250             | 5.8  | 20.0      |
| 2-Nitrotoluene             | Ave           | 127896 | 125248 |        | 245            | 250             | -2.1 | 20.0      |
| 4-Nitrotoluene             | Ave           | 111880 | 107640 |        | 241            | 250             | -3.8 | 20.0      |
| 3-Nitrotoluene             | Ave           | 140492 | 136876 |        | 244            | 250             | -2.6 | 20.0      |
| PETN                       | Ave           | 68845  | 74703  |        | 2710           | 2500            | 8.5  | 20.0      |
| 1,2-Dinitrobenzene         | Ave           | 126309 | 140776 |        | 279            | 250             | 11.5 | 20.0      |

FORM VII  
HPLC/IC CONTINUING CALIBRATION RETENTION TIME SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
 SDG No.: \_\_\_\_\_  
 Lab Sample ID: CCV 280-613168/44 Calibration Date: 05/20/2023 03:16  
 Instrument ID: CHHPLC\_X3 Calib Start Date: 02/08/2023 15:38  
 GC Column: UltraCarb5uODS ID: 4.60 (mm) Calib End Date: 02/08/2023 18:42  
 Lab File ID: 05190044.D

| Analyte                    | RT    | RT WINDOW |       |
|----------------------------|-------|-----------|-------|
|                            |       | FROM      | TO    |
| HMX                        | 6.55  | 6.40      | 6.70  |
| RDX                        | 7.57  | 7.42      | 7.72  |
| Picric acid                | 7.97  | 7.81      | 8.11  |
| 1,3,5-Trinitrobenzene      | 8.65  | 8.50      | 8.80  |
| 1,3-Dinitrobenzene         | 9.28  | 9.13      | 9.43  |
| Nitrobenzene               | 9.65  | 9.49      | 9.79  |
| 3,5-Dinitroaniline         | 9.89  | 9.73      | 10.03 |
| Tetryl                     | 10.03 | 9.87      | 10.17 |
| Nitroglycerin              | 10.48 | 10.32     | 10.62 |
| 2,4,6-Trinitrotoluene      | 10.92 | 10.81     | 11.01 |
| 4-Amino-2,6-dinitrotoluene | 11.12 | 11.02     | 11.22 |
| 2-Amino-4,6-dinitrotoluene | 11.37 | 11.27     | 11.47 |
| 2,6-Dinitrotoluene         | 11.52 | 11.42     | 11.62 |
| 2,4-Dinitrotoluene         | 11.69 | 11.58     | 11.78 |
| 2-Nitrotoluene             | 12.51 | 12.37     | 12.67 |
| 4-Nitrotoluene             | 12.93 | 12.79     | 13.09 |
| 3-Nitrotoluene             | 13.51 | 13.37     | 13.67 |
| PETN                       | 14.69 | 14.55     | 14.85 |
| 1,2-Dinitrobenzene         | 8.53  | 8.37      | 8.67  |

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\05190044.D  
 Lims ID: CCV DMT  
 Client ID:  
 Sample Type: CCV  
 Inject. Date: 20-May-2023 03:16:42 ALS Bottle#: 7 Worklist Smp#: 44  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: CCV DMT  
 Operator ID: JZ/JG Instrument ID: CHHPLC\_X3  
 Sublist: chrom-8330\_X3\*sub27  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 20-May-2023 11:22:12 Calib Date: 24-Feb-2023 17:41:20  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1663

First Level Reviewer: LV5D

Date: 20-May-2023 11:05:27

| Compound                      | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/mL | OnCol Amt<br>ug/mL | Flags |
|-------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 3 TNX                         | 1   | 6.432        | 6.430            | 0.002            | 49301    | 0.2508           | 0.2494             | M     |
| 4 HMX                         | 1   | 6.552        | 6.550            | 0.002            | 23736    | 0.2500           | 0.2538             | M     |
| 6 DNx                         | 1   | 6.758        | 6.757            | 0.001            | 36943    | 0.2503           | 0.2557             | M     |
| 7 MNx                         | 1   | 7.185        | 7.190            | -0.005           | 39943    | 0.2918           | 0.3045             |       |
| 8 RDX                         | 1   | 7.572        | 7.570            | 0.002            | 28790    | 0.2500           | 0.2706             |       |
| 9 2,4,6-Trinitrophenol        | 1   | 7.972        | 7.963            | 0.009            | 21226    | 0.2500           | 0.2799             |       |
| \$ 10 1,2-Dinitrobenzene      | 1   | 8.532        | 8.523            | 0.009            | 35194    | 0.2500           | 0.2786             |       |
| 11 1,3,5-Trinitrobenzene      | 1   | 8.652        | 8.650            | 0.002            | 56391    | 0.2500           | 0.2597             |       |
| 12 1,3-Dinitrobenzene         | 1   | 9.278        | 9.277            | 0.001            | 78892    | 0.2500           | 0.2680             |       |
| 13 Nitrobenzene               | 1   | 9.651        | 9.643            | 0.008            | 47978    | 0.2500           | 0.2509             |       |
| 14 3,5-Dinitroaniline         | 1   | 9.891        | 9.883            | 0.008            | 58775    | 0.2500           | 0.2571             |       |
| 15 Tetryl                     | 1   | 10.025       | 10.017           | 0.008            | 37457    | 0.2500           | 0.2282             |       |
| 16 Nitroglycerin              | 2   | 10.478       | 10.470           | 0.008            | 171104   | 2.50             | 2.67               |       |
| 17 2,4,6-Trinitrotoluene      | 1   | 10.918       | 10.910           | 0.008            | 54470    | 0.2500           | 0.2581             |       |
| 18 4-Amino-2,6-dinitrotoluene | 1   | 11.118       | 11.117           | 0.001            | 40396    | 0.2500           | 0.2607             |       |
| 19 2-Amino-4,6-dinitrotoluene | 1   | 11.371       | 11.370           | 0.001            | 51423    | 0.2500           | 0.2553             |       |
| 20 2,6-Dinitrotoluene         | 1   | 11.518       | 11.517           | 0.001            | 39048    | 0.2500           | 0.2735             |       |
| 21 2,4-Dinitrotoluene         | 1   | 11.691       | 11.683           | 0.008            | 78500    | 0.2500           | 0.2646             |       |
| 22 o-Nitrotoluene             | 1   | 12.511       | 12.517           | -0.006           | 31312    | 0.2500           | 0.2448             |       |
| 23 p-Nitrotoluene             | 1   | 12.931       | 12.937           | -0.006           | 26910    | 0.2500           | 0.2405             |       |
| 24 m-Nitrotoluene             | 1   | 13.511       | 13.517           | -0.006           | 34219    | 0.2500           | 0.2436             |       |
| 25 PETN                       | 2   | 14.685       | 14.697           | -0.012           | 186758   | 2.50             | 2.71               |       |

**QC Flag Legend**

Processing Flags

Review Flags

M - Manually Integrated

**Reagents:**

8330 DMT\_00013

Amount Added: 12.50

Units: uL

8330IntermStk\_00075

Amount Added: 25.00

Units: uL

Report Date: 20-May-2023 11:22:13

Chrom Revision: 2.3 16-May-2023 16:37:30

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190044.d

Injection Date: 20-May-2023 03:16:42

Instrument ID: CHHPLC\_X3

Operator ID: JZ/JG

Lims ID: CCV DMT

Worklist Smp#: 44

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 7

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190044.d

Injection Date: 20-May-2023 03:16:42

Instrument ID: CHHPLC\_X3

Lims ID: CCV DMT

Client ID:

Operator ID: JZ/JG

ALS Bottle#:

7

Worklist Smp#: 44

Injection Vol: 100.0 ul

Dil. Factor: 1.0000

Method: 8330\_X3

Limit Group: GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm)

Detector


LC DAD1B, 254 nm

**4 HMX, CAS: 2691-41-0**

Signal: 1


RT: 6.55  
Area: 28479  
Amount: 0.304497  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.55  
Area: 23736  
Amount: 0.253785  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 11:05:20 -06:00:00 (UTC)

Audit Action: Split an Integrated Peak

Audit Reason: Baseline



FORM VII  
HPLC/IC CONTINUING CALIBRATION DATA

Lab Name: Eurofins Denver Job No.: 280-176674-1  
SDG No.: \_\_\_\_\_  
Lab Sample ID: CCV 280-613168/44 Calibration Date: 05/20/2023 03:16  
Instrument ID: CHHPLC\_X3 Calib Start Date: 02/24/2023 15:00  
GC Column: UltraCarb5uODS ID: 4.60 (mm) Calib End Date: 02/24/2023 17:41  
Lab File ID: 05190044.D Conc. Units: ug/L

| ANALYTE | CURVE<br>TYPE | AVE CF | CF     | MIN CF | CALC<br>AMOUNT | SPIKE<br>AMOUNT | %D   | MAX<br>%D |
|---------|---------------|--------|--------|--------|----------------|-----------------|------|-----------|
| TNX     | Ave           | 197687 | 196614 |        | 249            | 251             | -0.5 | 20.0      |
| DNX     | Ave           | 144462 | 147624 |        | 256            | 250             | 2.2  | 20.0      |
| MNX     | Ave           | 131172 | 136908 |        | 305            | 292             | 4.4  | 20.0      |

FORM VII  
HPLC/IC CONTINUING CALIBRATION RETENTION TIME SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
SDG No.: \_\_\_\_\_  
Lab Sample ID: CCV 280-613168/44 Calibration Date: 05/20/2023 03:16  
Instrument ID: CHHPLC\_X3 Calib Start Date: 02/24/2023 15:00  
GC Column: UltraCarb5uODS ID: 4.60 (mm) Calib End Date: 02/24/2023 17:41  
Lab File ID: 05190044.D

| Analyte | RT   | RT WINDOW |      |
|---------|------|-----------|------|
|         |      | FROM      | TO   |
| TNX     | 6.43 | 6.33      | 6.53 |
| DNX     | 6.76 | 6.66      | 6.86 |
| MNX     | 7.19 | 7.04      | 7.34 |

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\05190044.D  
 Lims ID: CCV DMT  
 Client ID:  
 Sample Type: CCV  
 Inject. Date: 20-May-2023 03:16:42 ALS Bottle#: 7 Worklist Smp#: 44  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: CCV DMT  
 Operator ID: JZ/JG Instrument ID: CHHPLC\_X3  
 Sublist: chrom-8330\_X3\*sub27  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 20-May-2023 11:22:12 Calib Date: 24-Feb-2023 17:41:20  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1663

First Level Reviewer: LV5D

Date: 20-May-2023 11:05:27

| Compound                      | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/mL | OnCol Amt<br>ug/mL | Flags |
|-------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 3 TNX                         | 1   | 6.432        | 6.430            | 0.002            | 49301    | 0.2508           | 0.2494             | M     |
| 4 HMX                         | 1   | 6.552        | 6.550            | 0.002            | 23736    | 0.2500           | 0.2538             | M     |
| 6 DNx                         | 1   | 6.758        | 6.757            | 0.001            | 36943    | 0.2503           | 0.2557             | M     |
| 7 MNx                         | 1   | 7.185        | 7.190            | -0.005           | 39943    | 0.2918           | 0.3045             |       |
| 8 RDX                         | 1   | 7.572        | 7.570            | 0.002            | 28790    | 0.2500           | 0.2706             |       |
| 9 2,4,6-Trinitrophenol        | 1   | 7.972        | 7.963            | 0.009            | 21226    | 0.2500           | 0.2799             |       |
| \$ 10 1,2-Dinitrobenzene      | 1   | 8.532        | 8.523            | 0.009            | 35194    | 0.2500           | 0.2786             |       |
| 11 1,3,5-Trinitrobenzene      | 1   | 8.652        | 8.650            | 0.002            | 56391    | 0.2500           | 0.2597             |       |
| 12 1,3-Dinitrobenzene         | 1   | 9.278        | 9.277            | 0.001            | 78892    | 0.2500           | 0.2680             |       |
| 13 Nitrobenzene               | 1   | 9.651        | 9.643            | 0.008            | 47978    | 0.2500           | 0.2509             |       |
| 14 3,5-Dinitroaniline         | 1   | 9.891        | 9.883            | 0.008            | 58775    | 0.2500           | 0.2571             |       |
| 15 Tetryl                     | 1   | 10.025       | 10.017           | 0.008            | 37457    | 0.2500           | 0.2282             |       |
| 16 Nitroglycerin              | 2   | 10.478       | 10.470           | 0.008            | 171104   | 2.50             | 2.67               |       |
| 17 2,4,6-Trinitrotoluene      | 1   | 10.918       | 10.910           | 0.008            | 54470    | 0.2500           | 0.2581             |       |
| 18 4-Amino-2,6-dinitrotoluene | 1   | 11.118       | 11.117           | 0.001            | 40396    | 0.2500           | 0.2607             |       |
| 19 2-Amino-4,6-dinitrotoluene | 1   | 11.371       | 11.370           | 0.001            | 51423    | 0.2500           | 0.2553             |       |
| 20 2,6-Dinitrotoluene         | 1   | 11.518       | 11.517           | 0.001            | 39048    | 0.2500           | 0.2735             |       |
| 21 2,4-Dinitrotoluene         | 1   | 11.691       | 11.683           | 0.008            | 78500    | 0.2500           | 0.2646             |       |
| 22 o-Nitrotoluene             | 1   | 12.511       | 12.517           | -0.006           | 31312    | 0.2500           | 0.2448             |       |
| 23 p-Nitrotoluene             | 1   | 12.931       | 12.937           | -0.006           | 26910    | 0.2500           | 0.2405             |       |
| 24 m-Nitrotoluene             | 1   | 13.511       | 13.517           | -0.006           | 34219    | 0.2500           | 0.2436             |       |
| 25 PETN                       | 2   | 14.685       | 14.697           | -0.012           | 186758   | 2.50             | 2.71               |       |

## QC Flag Legend

Processing Flags

Review Flags

M - Manually Integrated

**Reagents:**

8330 DMT\_00013

Amount Added: 12.50

Units: uL

8330IntermStk\_00075

Amount Added: 25.00

Units: uL

Report Date: 20-May-2023 11:22:13

Chrom Revision: 2.3 16-May-2023 16:37:30

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190044.d

Injection Date: 20-May-2023 03:16:42

Instrument ID: CHHPLC\_X3

Operator ID: JZ/JG

Lims ID: CCV DMT

Worklist Smp#: 44

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 7

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190044.d  
Injection Date: 20-May-2023 03:16:42 Instrument ID: CHHPLC\_X3  
Lims ID: CCV DMT  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 7 Worklist Smp#: 44  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**3 TNX, CAS: 13980-04-6**

Signal: 1


RT: 6.43  
Area: 53439  
Amount: 0.270321  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.43  
Area: 49301  
Amount: 0.249389  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 11:05:19 -06:00:00 (UTC)

Audit Action: Split an Integrated Peak

Audit Reason: Baseline

## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190044.d

Injection Date: 20-May-2023 03:16:42

Instrument ID: CHHPLC\_X3

Lims ID: CCV DMT

Client ID:

Operator ID: JZ/JG

ALS Bottle#:

7

Worklist Smp#: 44

Injection Vol: 100.0 ul

Dil. Factor: 1.0000

Method: 8330\_X3

Limit Group: GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm)

Detector


LC DAD1B, 254 nm

**6 DNX, CAS: 80251-29-2**

Signal: 1


RT: 6.76  
Area: 39950  
Amount: 0.276542  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.76  
Area: 36943  
Amount: 0.255727  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 11:05:26 -06:00:00 (UTC)

Audit Action: Manually Integrated

Audit Reason: Baseline

FORM VII  
HPLC/IC CONTINUING CALIBRATION DATA

Lab Name: Eurofins Denver Job No.: 280-176674-1  
SDG No.: \_\_\_\_\_  
Lab Sample ID: CCV 280-613168/45 Calibration Date: 05/20/2023 03:39  
Instrument ID: CHHPLC\_X3 Calib Start Date: 01/17/2023 00:33  
GC Column: UltraCarb5uODS ID: 4.60 (mm) Calib End Date: 01/17/2023 03:16  
Lab File ID: 05190045.D Conc. Units: ug/L

| ANALYTE                    | CURVE<br>TYPE | AVE CF | CF     | MIN CF | CALC<br>AMOUNT | SPIKE<br>AMOUNT | %D   | MAX<br>%D |
|----------------------------|---------------|--------|--------|--------|----------------|-----------------|------|-----------|
| 2,6-diamino-4-nitrotoluene | Lin2          |        | 220240 |        | 236            | 250             | -5.6 | 20.0      |
| 2,4-diamino-6-nitrotoluene | Ave           | 146472 | 137824 |        | 235            | 250             | -5.9 | 20.0      |



FORM VII  
HPLC/IC CONTINUING CALIBRATION RETENTION TIME SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
SDG No.: \_\_\_\_\_  
Lab Sample ID: CCV 280-613168/45 Calibration Date: 05/20/2023 03:39  
Instrument ID: CHHPLC\_X3 Calib Start Date: 01/17/2023 00:33  
GC Column: UltraCarb5uODS ID: 4.60 (mm) Calib End Date: 01/17/2023 03:16  
Lab File ID: 05190045.D

| Analyte                    | RT   | RT WINDOW |      |
|----------------------------|------|-----------|------|
|                            |      | FROM      | TO   |
| 2,6-diamino-4-nitrotoluene | 6.43 | 6.28      | 6.58 |
| 2,4-diamino-6-nitrotoluene | 6.62 | 6.47      | 6.77 |

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\05190045.D  
 Lims ID: CCV ADD  
 Client ID:  
 Sample Type: CCV  
 Inject. Date: 20-May-2023 03:39:42 ALS Bottle#: 8 Worklist Smp#: 45  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: CCV ADD  
 Operator ID: JZ/JG Instrument ID: CHHPLC\_X3  
 Sublist: chrom-8330\_X3\*sub28  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 20-May-2023 11:22:13 Calib Date: 24-Feb-2023 17:41:20  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1663

First Level Reviewer: LV5D

Date: 20-May-2023 11:05:33

| Compound                     | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/mL | OnCol Amt<br>ug/mL | Flags |
|------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 2 2,6-diamino-4-nitrotoluene | 1   | 6.426        | 6.431            | -0.005           | 55060    | 0.2500           | 0.2359             | M     |
| 5 2,4-diamino-6-nitrotoluene | 1   | 6.619        | 6.618            | 0.001            | 34456    | 0.2500           | 0.2352             | M     |

**QC Flag Legend**

Review Flags

M - Manually Integrated

**Reagents:**

8330\_ADDs\_00037

Amount Added: 12.50

Units: uL

Report Date: 20-May-2023 11:22:13

Chrom Revision: 2.3 16-May-2023 16:37:30

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190045.d

Injection Date: 20-May-2023 03:39:42

Instrument ID: CHHPLC\_X3

Operator ID: JZ/JG

Lims ID: CCV ADD

Worklist Smp#: 45

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 8

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190045.d  
Injection Date: 20-May-2023 03:39:42 Instrument ID: CHHPLC\_X3  
Lims ID: CCV ADD  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 8 Worklist Smp#: 45  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**2,2,6-diamino-4-nitrotoluene, CAS: 59229-75-3**

Signal: 1


RT: 6.43  
Area: 55550  
Amount: 0.238082  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.43  
Area: 55060  
Amount: 0.235949  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 11:05:31 -06:00:00 (UTC)

Audit Action: Split an Integrated Peak

Audit Reason: Baseline

## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190045.d

Injection Date: 20-May-2023 03:39:42

Instrument ID: CHHPLC\_X3

Lims ID: CCV ADD

Client ID:

Operator ID: JZ/JG

ALS Bottle#:

8

Worklist Smp#: 45

Injection Vol: 100.0 ul

Dil. Factor:

1.0000

Method: 8330\_X3

Limit Group:

GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm)

Detector


LC DAD1B, 254 nm

**5 2,4-diamino-6-nitrotoluene, CAS: 6629-29-4**

Signal: 1


RT: 6.62  
Area: 35768  
Amount: 0.244197  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.62  
Area: 34456  
Amount: 0.235239  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 11:05:32 -06:00:00 (UTC)

Audit Action: Split an Integrated Peak

Audit Reason: Baseline

FORM VII  
HPLC/IC CONTINUING CALIBRATION DATA

Lab Name: Eurofins Denver Job No.: 280-176674-1  
 SDG No.: \_\_\_\_\_  
 Lab Sample ID: CCV 280-613168/56 Calibration Date: 05/20/2023 07:52  
 Instrument ID: CHHPLC\_X3 Calib Start Date: 02/08/2023 15:38  
 GC Column: UltraCarb5uODS ID: 4.60 (mm) Calib End Date: 02/08/2023 18:42  
 Lab File ID: 05190056.D Conc. Units: ug/L

| ANALYTE                    | CURVE<br>TYPE | AVE CF | CF     | MIN CF | CALC<br>AMOUNT | SPIKE<br>AMOUNT | %D   | MAX<br>%D |
|----------------------------|---------------|--------|--------|--------|----------------|-----------------|------|-----------|
| HMX                        | Ave           | 93528  | 102472 |        | 274            | 250             | 9.6  | 20.0      |
| RDX                        | Ave           | 106380 | 116316 |        | 273            | 250             | 9.3  | 20.0      |
| Picric acid                | Ave           | 75830  | 84504  |        | 279            | 250             | 11.4 | 20.0      |
| 1,3,5-Trinitrobenzene      | Ave           | 217147 | 232972 |        | 268            | 250             | 7.3  | 20.0      |
| 1,3-Dinitrobenzene         | Ave           | 294397 | 316060 |        | 268            | 250             | 7.4  | 20.0      |
| Nitrobenzene               | Ave           | 191245 | 192044 |        | 251            | 250             | 0.4  | 20.0      |
| 3,5-Dinitroaniline         | Lin2          |        | 243092 |        | 266            | 250             | 6.3  | 20.0      |
| Tetryl                     | Ave           | 164121 | 149672 |        | 228            | 250             | -8.8 | 20.0      |
| Nitroglycerin              | Ave           | 64070  | 68638  |        | 2680           | 2500            | 7.1  | 20.0      |
| 2,4,6-Trinitrotoluene      | Ave           | 211040 | 218544 |        | 259            | 250             | 3.6  | 20.0      |
| 4-Amino-2,6-dinitrotoluene | Ave           | 154933 | 161800 |        | 261            | 250             | 4.4  | 20.0      |
| 2-Amino-4,6-dinitrotoluene | Ave           | 201410 | 211024 |        | 262            | 250             | 4.8  | 20.0      |
| 2,6-Dinitrotoluene         | Ave           | 142745 | 151340 |        | 265            | 250             | 6.0  | 20.0      |
| 2,4-Dinitrotoluene         | Ave           | 296667 | 318628 |        | 269            | 250             | 7.4  | 20.0      |
| 2-Nitrotoluene             | Ave           | 127896 | 124148 |        | 243            | 250             | -2.9 | 20.0      |
| 4-Nitrotoluene             | Ave           | 111880 | 106836 |        | 239            | 250             | -4.5 | 20.0      |
| 3-Nitrotoluene             | Ave           | 140492 | 136984 |        | 244            | 250             | -2.5 | 20.0      |
| PETN                       | Ave           | 68845  | 74686  |        | 2710           | 2500            | 8.5  | 20.0      |
| 1,2-Dinitrobenzene         | Ave           | 126309 | 140368 |        | 278            | 250             | 11.1 | 20.0      |

FORM VII  
HPLC/IC CONTINUING CALIBRATION RETENTION TIME SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
 SDG No.: \_\_\_\_\_  
 Lab Sample ID: CCV 280-613168/56 Calibration Date: 05/20/2023 07:52  
 Instrument ID: CHHPLC\_X3 Calib Start Date: 02/08/2023 15:38  
 GC Column: UltraCarb5uODS ID: 4.60 (mm) Calib End Date: 02/08/2023 18:42  
 Lab File ID: 05190056.D

| Analyte                    | RT    | RT WINDOW |       |
|----------------------------|-------|-----------|-------|
|                            |       | FROM      | TO    |
| HMX                        | 6.55  | 6.40      | 6.70  |
| RDX                        | 7.57  | 7.42      | 7.72  |
| Picric acid                | 7.97  | 7.81      | 8.11  |
| 1,3,5-Trinitrobenzene      | 8.65  | 8.50      | 8.80  |
| 1,3-Dinitrobenzene         | 9.27  | 9.13      | 9.43  |
| Nitrobenzene               | 9.64  | 9.49      | 9.79  |
| 3,5-Dinitroaniline         | 9.87  | 9.73      | 10.03 |
| Tetryl                     | 10.00 | 9.87      | 10.17 |
| Nitroglycerin              | 10.45 | 10.32     | 10.62 |
| 2,4,6-Trinitrotoluene      | 10.89 | 10.81     | 11.01 |
| 4-Amino-2,6-dinitrotoluene | 11.09 | 11.02     | 11.22 |
| 2-Amino-4,6-dinitrotoluene | 11.34 | 11.27     | 11.47 |
| 2,6-Dinitrotoluene         | 11.49 | 11.42     | 11.62 |
| 2,4-Dinitrotoluene         | 11.66 | 11.58     | 11.78 |
| 2-Nitrotoluene             | 12.49 | 12.37     | 12.67 |
| 4-Nitrotoluene             | 12.91 | 12.79     | 13.09 |
| 3-Nitrotoluene             | 13.48 | 13.37     | 13.67 |
| PETN                       | 14.65 | 14.55     | 14.85 |
| 1,2-Dinitrobenzene         | 8.52  | 8.37      | 8.67  |

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\05190056.D  
 Lims ID: CCV DMT  
 Client ID:  
 Sample Type: CCV  
 Inject. Date: 20-May-2023 07:52:17 ALS Bottle#: 7 Worklist Smp#: 56  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: CCV DMT  
 Operator ID: JZ/JG Instrument ID: CHHPLC\_X3  
 Sublist: chrom-8330\_X3\*sub27  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 20-May-2023 11:22:18 Calib Date: 24-Feb-2023 17:41:20  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1663

First Level Reviewer: LV5D

Date: 20-May-2023 11:16:53

| Compound                      | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/mL | OnCol Amt<br>ug/mL | Flags |
|-------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 3 TNX                         | 1   | 6.433        | 6.430            | 0.003            | 51212    | 0.2508           | 0.2591             | M     |
| 4 HMX                         | 1   | 6.553        | 6.550            | 0.003            | 25618    | 0.2500           | 0.2739             | M     |
| 6 DNx                         | 1   | 6.760        | 6.757            | 0.003            | 37116    | 0.2503           | 0.2569             | M     |
| 7 MNx                         | 1   | 7.187        | 7.190            | -0.003           | 40248    | 0.2918           | 0.3068             |       |
| 8 RDX                         | 1   | 7.573        | 7.570            | 0.003            | 29079    | 0.2500           | 0.2734             |       |
| 9 2,4,6-Trinitrophenol        | 1   | 7.967        | 7.963            | 0.004            | 21126    | 0.2500           | 0.2786             |       |
| \$ 10 1,2-Dinitrobenzene      | 1   | 8.520        | 8.523            | -0.003           | 35092    | 0.2500           | 0.2778             |       |
| 11 1,3,5-Trinitrobenzene      | 1   | 8.647        | 8.650            | -0.003           | 58243    | 0.2500           | 0.2682             |       |
| 12 1,3-Dinitrobenzene         | 1   | 9.267        | 9.277            | -0.010           | 79015    | 0.2500           | 0.2684             |       |
| 13 Nitrobenzene               | 1   | 9.640        | 9.643            | -0.003           | 48011    | 0.2500           | 0.2510             |       |
| 14 3,5-Dinitroaniline         | 1   | 9.873        | 9.883            | -0.010           | 60773    | 0.2500           | 0.2658             |       |
| 15 Tetryl                     | 1   | 10.000       | 10.017           | -0.017           | 37418    | 0.2500           | 0.2280             |       |
| 16 Nitroglycerin              | 2   | 10.453       | 10.470           | -0.017           | 171596   | 2.50             | 2.68               |       |
| 17 2,4,6-Trinitrotoluene      | 1   | 10.893       | 10.910           | -0.017           | 54636    | 0.2500           | 0.2589             |       |
| 18 4-Amino-2,6-dinitrotoluene | 1   | 11.087       | 11.117           | -0.030           | 40450    | 0.2500           | 0.2611             |       |
| 19 2-Amino-4,6-dinitrotoluene | 1   | 11.340       | 11.370           | -0.030           | 52756    | 0.2500           | 0.2619             |       |
| 20 2,6-Dinitrotoluene         | 1   | 11.493       | 11.517           | -0.024           | 37835    | 0.2500           | 0.2651             |       |
| 21 2,4-Dinitrotoluene         | 1   | 11.660       | 11.683           | -0.023           | 79657    | 0.2500           | 0.2685             |       |
| 22 o-Nitrotoluene             | 1   | 12.487       | 12.517           | -0.030           | 31037    | 0.2500           | 0.2427             |       |
| 23 p-Nitrotoluene             | 1   | 12.907       | 12.937           | -0.030           | 26709    | 0.2500           | 0.2387             |       |
| 24 m-Nitrotoluene             | 1   | 13.480       | 13.517           | -0.037           | 34246    | 0.2500           | 0.2438             |       |
| 25 PETN                       | 2   | 14.653       | 14.697           | -0.044           | 186714   | 2.50             | 2.71               |       |

**QC Flag Legend**

Processing Flags

Review Flags

M - Manually Integrated



**Reagents:**

8330 DMT\_00013

Amount Added: 12.50

Units: uL

8330IntermStk\_00075

Amount Added: 25.00

Units: uL

Report Date: 20-May-2023 11:22:18

Chrom Revision: 2.3 16-May-2023 16:37:30

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190056.d

Injection Date: 20-May-2023 07:52:17

Instrument ID: CHHPLC\_X3

Operator ID: JZ/JG

Lims ID: CCV DMT

Worklist Smp#: 56

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 7

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190056.d

Injection Date: 20-May-2023 07:52:17

Instrument ID: CHHPLC\_X3

Lims ID: CCV DMT

Client ID:

Operator ID: JZ/JG

ALS Bottle#:

7

Worklist Smp#: 56

Injection Vol: 100.0 ul

Dil. Factor: 1.0000

Method: 8330\_X3

Limit Group: GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm)


Detector: LC DAD1B, 254 nm

**4 HMX, CAS: 2691-41-0**

Signal: 1


RT: 6.55  
Area: 28233  
Amount: 0.301866  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.55  
Area: 25618  
Amount: 0.273907  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 11:16:49 -06:00:00 (UTC)

Audit Action: Split an Integrated Peak

Audit Reason: Baseline

FORM VII  
HPLC/IC CONTINUING CALIBRATION DATA

Lab Name: Eurofins Denver Job No.: 280-176674-1  
SDG No.: \_\_\_\_\_  
Lab Sample ID: CCV 280-613168/56 Calibration Date: 05/20/2023 07:52  
Instrument ID: CHHPLC\_X3 Calib Start Date: 02/24/2023 15:00  
GC Column: UltraCarb5uODS ID: 4.60 (mm) Calib End Date: 02/24/2023 17:41  
Lab File ID: 05190056.D Conc. Units: ug/L

| ANALYTE | CURVE<br>TYPE | AVE CF | CF     | MIN CF | CALC<br>AMOUNT | SPIKE<br>AMOUNT | %D  | MAX<br>%D |
|---------|---------------|--------|--------|--------|----------------|-----------------|-----|-----------|
| TNX     | Ave           | 197687 | 204235 |        | 259            | 251             | 3.3 | 20.0      |
| DNX     | Ave           | 144462 | 148316 |        | 257            | 250             | 2.7 | 20.0      |
| MNX     | Ave           | 131172 | 137954 |        | 307            | 292             | 5.2 | 20.0      |

FORM VII  
HPLC/IC CONTINUING CALIBRATION RETENTION TIME SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
SDG No.: \_\_\_\_\_  
Lab Sample ID: CCV 280-613168/56 Calibration Date: 05/20/2023 07:52  
Instrument ID: CHHPLC\_X3 Calib Start Date: 02/24/2023 15:00  
GC Column: UltraCarb5uODS ID: 4.60 (mm) Calib End Date: 02/24/2023 17:41  
Lab File ID: 05190056.D

| Analyte | RT   | RT WINDOW |      |
|---------|------|-----------|------|
|         |      | FROM      | TO   |
| TNX     | 6.43 | 6.33      | 6.53 |
| DNX     | 6.76 | 6.66      | 6.86 |
| MNX     | 7.19 | 7.04      | 7.34 |

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\05190056.D  
 Lims ID: CCV DMT  
 Client ID:  
 Sample Type: CCV  
 Inject. Date: 20-May-2023 07:52:17 ALS Bottle#: 7 Worklist Smp#: 56  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: CCV DMT  
 Operator ID: JZ/JG Instrument ID: CHHPLC\_X3  
 Sublist: chrom-8330\_X3\*sub27  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 20-May-2023 11:22:18 Calib Date: 24-Feb-2023 17:41:20  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1663

First Level Reviewer: LV5D

Date: 20-May-2023 11:16:53

| Compound                      | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/mL | OnCol Amt<br>ug/mL | Flags |
|-------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 3 TNX                         | 1   | 6.433        | 6.430            | 0.003            | 51212    | 0.2508           | 0.2591             | M     |
| 4 HMX                         | 1   | 6.553        | 6.550            | 0.003            | 25618    | 0.2500           | 0.2739             | M     |
| 6 DNx                         | 1   | 6.760        | 6.757            | 0.003            | 37116    | 0.2503           | 0.2569             | M     |
| 7 MNX                         | 1   | 7.187        | 7.190            | -0.003           | 40248    | 0.2918           | 0.3068             |       |
| 8 RDX                         | 1   | 7.573        | 7.570            | 0.003            | 29079    | 0.2500           | 0.2734             |       |
| 9 2,4,6-Trinitrophenol        | 1   | 7.967        | 7.963            | 0.004            | 21126    | 0.2500           | 0.2786             |       |
| \$ 10 1,2-Dinitrobenzene      | 1   | 8.520        | 8.523            | -0.003           | 35092    | 0.2500           | 0.2778             |       |
| 11 1,3,5-Trinitrobenzene      | 1   | 8.647        | 8.650            | -0.003           | 58243    | 0.2500           | 0.2682             |       |
| 12 1,3-Dinitrobenzene         | 1   | 9.267        | 9.277            | -0.010           | 79015    | 0.2500           | 0.2684             |       |
| 13 Nitrobenzene               | 1   | 9.640        | 9.643            | -0.003           | 48011    | 0.2500           | 0.2510             |       |
| 14 3,5-Dinitroaniline         | 1   | 9.873        | 9.883            | -0.010           | 60773    | 0.2500           | 0.2658             |       |
| 15 Tetryl                     | 1   | 10.000       | 10.017           | -0.017           | 37418    | 0.2500           | 0.2280             |       |
| 16 Nitroglycerin              | 2   | 10.453       | 10.470           | -0.017           | 171596   | 2.50             | 2.68               |       |
| 17 2,4,6-Trinitrotoluene      | 1   | 10.893       | 10.910           | -0.017           | 54636    | 0.2500           | 0.2589             |       |
| 18 4-Amino-2,6-dinitrotoluene | 1   | 11.087       | 11.117           | -0.030           | 40450    | 0.2500           | 0.2611             |       |
| 19 2-Amino-4,6-dinitrotoluene | 1   | 11.340       | 11.370           | -0.030           | 52756    | 0.2500           | 0.2619             |       |
| 20 2,6-Dinitrotoluene         | 1   | 11.493       | 11.517           | -0.024           | 37835    | 0.2500           | 0.2651             |       |
| 21 2,4-Dinitrotoluene         | 1   | 11.660       | 11.683           | -0.023           | 79657    | 0.2500           | 0.2685             |       |
| 22 o-Nitrotoluene             | 1   | 12.487       | 12.517           | -0.030           | 31037    | 0.2500           | 0.2427             |       |
| 23 p-Nitrotoluene             | 1   | 12.907       | 12.937           | -0.030           | 26709    | 0.2500           | 0.2387             |       |
| 24 m-Nitrotoluene             | 1   | 13.480       | 13.517           | -0.037           | 34246    | 0.2500           | 0.2438             |       |
| 25 PETN                       | 2   | 14.653       | 14.697           | -0.044           | 186714   | 2.50             | 2.71               |       |

**QC Flag Legend**

Processing Flags

Review Flags

M - Manually Integrated

**Reagents:**

8330 DMT\_00013

Amount Added: 12.50

Units: uL

8330IntermStk\_00075

Amount Added: 25.00

Units: uL

Report Date: 20-May-2023 11:22:18

Chrom Revision: 2.3 16-May-2023 16:37:30

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190056.d

Injection Date: 20-May-2023 07:52:17

Instrument ID: CHHPLC\_X3

Operator ID: JZ/JG

Lims ID: CCV DMT

Worklist Smp#: 56

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 7

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1





## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190056.d  
Injection Date: 20-May-2023 07:52:17 Instrument ID: CHHPLC\_X3  
Lims ID: CCV DMT  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 7 Worklist Smp#: 56  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**3 TNX, CAS: 13980-04-6**

Signal: 1


RT: 6.43  
Area: 53726  
Amount: 0.271773  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.43  
Area: 51212  
Amount: 0.259056  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 11:16:48 -06:00:00 (UTC)

Audit Action: Split an Integrated Peak

Audit Reason: Baseline

## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190056.d  
Injection Date: 20-May-2023 07:52:17 Instrument ID: CHHPLC\_X3  
Lims ID: CCV DMT  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 7 Worklist Smp#: 56  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector: LC DAD1B, 254 nm

**6 DNX, CAS: 80251-29-2**

Signal: 1


RT: 6.76  
Area: 41085  
Amount: 0.284399  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.76  
Area: 37116  
Amount: 0.256925  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 11:16:52 -06:00:00 (UTC)

Audit Action: Manually Integrated

Audit Reason: Baseline

FORM VII  
HPLC/IC CONTINUING CALIBRATION DATA

Lab Name: Eurofins Denver Job No.: 280-176674-1  
SDG No.: \_\_\_\_\_  
Lab Sample ID: CCV 280-613168/57 Calibration Date: 05/20/2023 08:15  
Instrument ID: CHHPLC\_X3 Calib Start Date: 01/17/2023 00:33  
GC Column: UltraCarb5uODS ID: 4.60 (mm) Calib End Date: 01/17/2023 03:16  
Lab File ID: 05190057.D Conc. Units: ug/L

| ANALYTE                    | CURVE<br>TYPE | AVE CF | CF     | MIN CF | CALC<br>AMOUNT | SPIKE<br>AMOUNT | %D   | MAX<br>%D |
|----------------------------|---------------|--------|--------|--------|----------------|-----------------|------|-----------|
| 2,6-diamino-4-nitrotoluene | Lin2          |        | 221428 |        | 237            | 250             | -5.1 | 20.0      |
| 2,4-diamino-6-nitrotoluene | Ave           | 146472 | 135328 |        | 231            | 250             | -7.6 | 20.0      |

FORM VII  
HPLC/IC CONTINUING CALIBRATION RETENTION TIME SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
SDG No.: \_\_\_\_\_  
Lab Sample ID: CCV 280-613168/57 Calibration Date: 05/20/2023 08:15  
Instrument ID: CHHPLC\_X3 Calib Start Date: 01/17/2023 00:33  
GC Column: UltraCarb5uODS ID: 4.60 (mm) Calib End Date: 01/17/2023 03:16  
Lab File ID: 05190057.D

| Analyte                    | RT   | RT WINDOW |      |
|----------------------------|------|-----------|------|
|                            |      | FROM      | TO   |
| 2,6-diamino-4-nitrotoluene | 6.43 | 6.28      | 6.58 |
| 2,4-diamino-6-nitrotoluene | 6.62 | 6.47      | 6.77 |

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\05190057.D  
 Lims ID: CCV ADD  
 Client ID:  
 Sample Type: CCV  
 Inject. Date: 20-May-2023 08:15:18 ALS Bottle#: 8 Worklist Smp#: 57  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: CCV ADD  
 Operator ID: JZ/JG Instrument ID: CHHPLC\_X3  
 Sublist: chrom-8330\_X3\*sub28  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 20-May-2023 11:22:18 Calib Date: 24-Feb-2023 17:41:20  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1663

First Level Reviewer: LV5D

Date: 20-May-2023 11:17:00

| Compound                     | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/mL | OnCol Amt<br>ug/mL | Flags |
|------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 2 2,6-diamino-4-nitrotoluene | 1   | 6.429        | 6.431            | -0.002           | 55357    | 0.2500           | 0.2372             | M     |
| 5 2,4-diamino-6-nitrotoluene | 1   | 6.622        | 6.618            | 0.004            | 33832    | 0.2500           | 0.2310             | M     |

**QC Flag Legend**

Review Flags

M - Manually Integrated

**Reagents:**

8330\_ADDs\_00037

Amount Added: 12.50

Units: uL

Report Date: 20-May-2023 11:22:18

Chrom Revision: 2.3 16-May-2023 16:37:30

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190057.d

Injection Date: 20-May-2023 08:15:18

Instrument ID: CHHPLC\_X3

Operator ID: JZ/JG

Lims ID: CCV ADD

Worklist Smp#: 57

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 8

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190057.d  
Injection Date: 20-May-2023 08:15:18 Instrument ID: CHHPLC\_X3  
Lims ID: CCV ADD  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 8 Worklist Smp#: 57  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**2,2,6-diamino-4-nitrotoluene, CAS: 59229-75-3**

Signal: 1


RT: 6.43  
Area: 55871  
Amount: 0.239480  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.43  
Area: 55357  
Amount: 0.237242  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 11:16:57 -06:00:00 (UTC)

Audit Action: Split an Integrated Peak

Audit Reason: Baseline

## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190057.d

Injection Date: 20-May-2023 08:15:18

Instrument ID: CHHPLC\_X3

Lims ID: CCV ADD

Client ID:

Operator ID: JZ/JG

ALS Bottle#:

8

Worklist Smp#: 57

Injection Vol: 100.0 ul

Dil. Factor:

1.0000

Method: 8330\_X3

Limit Group:

GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm)

Detector


LC DAD1B, 254 nm

**5 2,4-diamino-6-nitrotoluene, CAS: 6629-29-4**

Signal: 1


RT: 6.62  
Area: 35329  
Amount: 0.241199  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.62  
Area: 33832  
Amount: 0.230979  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 11:16:59 -06:00:00 (UTC)

Audit Action: Split an Integrated Peak

Audit Reason: Baseline



FORM VII  
HPLC/IC CONTINUING CALIBRATION DATA

Lab Name: Eurofins Denver Job No.: 280-176674-1  
 SDG No.: \_\_\_\_\_  
 Lab Sample ID: ICV 280-610603/19 Calibration Date: 04/28/2023 23:38  
 Instrument ID: CHHPLC\_X5 Calib Start Date: 04/28/2023 18:24  
 GC Column: Luna-phenylhex ID: 4.60 (mm) Calib End Date: 04/28/2023 23:03  
 Lab File ID: 04280019.D Conc. Units: ug/L

| ANALYTE                    | CURVE<br>TYPE | AVE CF | CF     | MIN CF | CALC<br>AMOUNT | SPIKE<br>AMOUNT | %D    | MAX<br>%D |
|----------------------------|---------------|--------|--------|--------|----------------|-----------------|-------|-----------|
| TNX                        | Ave           | 394332 | 392929 |        | 500            | 502             | -0.4  | 20.0      |
| DNX                        | Ave           | 277072 | 289930 |        | 524            | 501             | 4.6   | 20.0      |
| HMX                        | Lin2          |        | 156404 |        | 417            | 500             | -16.5 | 20.0      |
| MNX                        | Ave           | 262267 | 264334 |        | 588            | 584             | 0.8   | 20.0      |
| Picric acid                | Lin2          |        | 156686 |        | 487            | 500             | -2.6  | 20.0      |
| RDX                        | Lin2          |        | 195166 |        | 449            | 500             | -10.1 | 20.0      |
| Nitrobenzene               | Ave           | 398834 | 373736 |        | 469            | 500             | -6.3  | 20.0      |
| 3,5-Dinitroaniline         | Lin2          |        | 425372 |        | 477            | 500             | -4.6  | 20.0      |
| 1,3-Dinitrobenzene         | Ave           | 629351 | 577162 |        | 459            | 500             | -8.3  | 20.0      |
| Nitroglycerin              | Ave           | 134735 | 124537 |        | 4620           | 5000            | -7.6  | 20.0      |
| 2-Nitrotoluene             | Ave           | 246963 | 231686 |        | 469            | 500             | -6.2  | 20.0      |
| 4-Nitrotoluene             | Ave           | 227490 | 217326 |        | 478            | 500             | -4.5  | 20.0      |
| 4-Amino-2,6-dinitrotoluene | Ave           | 294210 | 272862 |        | 464            | 500             | -7.3  | 20.0      |
| 3-Nitrotoluene             | Ave           | 294634 | 265188 |        | 450            | 500             | -10.0 | 20.0      |
| 2-Amino-4,6-dinitrotoluene | Ave           | 427939 | 400432 |        | 468            | 500             | -6.4  | 20.0      |
| 1,3,5-Trinitrobenzene      | Ave           | 405692 | 392558 |        | 484            | 500             | -3.2  | 20.0      |
| 2,6-Dinitrotoluene         | Ave           | 291675 | 264938 |        | 454            | 500             | -9.2  | 20.0      |
| 2,4-Dinitrotoluene         | Ave           | 584463 | 526892 |        | 451            | 500             | -9.9  | 20.0      |
| Tetryl                     | Ave           | 284731 | 322214 |        | 566            | 500             | 13.2  | 20.0      |
| 2,4,6-Trinitrotoluene      | Ave           | 438876 | 392624 |        | 447            | 500             | -10.5 | 20.0      |
| PETN                       | Ave           | 143938 | 138644 |        | 4820           | 5000            | -3.7  | 20.0      |
| 1,2-Dinitrobenzene         | Ave           | 280697 | 261152 |        | 465            | 500             | -7.0  | 20.0      |

FORM VII  
HPLC/IC CONTINUING CALIBRATION RETENTION TIME SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
 SDG No.: \_\_\_\_\_  
 Lab Sample ID: ICV 280-610603/19 Calibration Date: 04/28/2023 23:38  
 Instrument ID: CHHPLC\_X5 Calib Start Date: 04/28/2023 18:24  
 GC Column: Luna-phenylhex ID: 4.60 (mm) Calib End Date: 04/28/2023 23:03  
 Lab File ID: 04280019.D

| Analyte                    | RT    | RT WINDOW |       |
|----------------------------|-------|-----------|-------|
|                            |       | FROM      | TO    |
| TNX                        | 5.08  | 4.97      | 5.27  |
| DNX                        | 5.92  | 5.81      | 6.11  |
| HMX                        | 6.58  | 6.45      | 6.75  |
| MNX                        | 7.46  | 7.34      | 7.64  |
| Picric acid                | 7.94  | 7.88      | 8.18  |
| RDX                        | 8.88  | 8.77      | 9.07  |
| Nitrobenzene               | 11.45 | 11.33     | 11.63 |
| 3,5-Dinitroaniline         | 14.32 | 14.21     | 14.51 |
| 1,3-Dinitrobenzene         | 14.59 | 14.47     | 14.77 |
| Nitroglycerin              | 15.16 | 15.03     | 15.33 |
| 2-Nitrotoluene             | 15.70 | 15.58     | 15.88 |
| 4-Nitrotoluene             | 15.92 | 15.81     | 16.11 |
| 4-Amino-2,6-dinitrotoluene | 16.48 | 16.37     | 16.67 |
| 3-Nitrotoluene             | 16.78 | 16.67     | 16.97 |
| 2-Amino-4,6-dinitrotoluene | 17.34 | 17.21     | 17.51 |
| 1,3,5-Trinitrobenzene      | 17.42 | 17.27     | 17.57 |
| 2,6-Dinitrotoluene         | 18.60 | 18.48     | 18.78 |
| 2,4-Dinitrotoluene         | 19.06 | 18.93     | 19.23 |
| Tetryl                     | 22.33 | 22.24     | 22.54 |
| 2,4,6-Trinitrotoluene      | 23.12 | 23.03     | 23.33 |
| PETN                       | 24.39 | 24.33     | 24.63 |
| 1,2-Dinitrobenzene         | 12.44 | 12.32     | 12.62 |

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280019.D  
 Lims ID: ICV INT  
 Client ID:  
 Sample Type: ICV  
 Inject. Date: 28-Apr-2023 23:38:32 ALS Bottle#: 19 Worklist Smp#: 19  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: ICV INT  
 Operator ID: JZ/JG Instrument ID: CHHPLC\_X5  
 Sublist:  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\8330\_X5\_Luna.m  
 Limit Group: GCSV - 8330  
 Last Update: 29-Apr-2023 13:04:23 Calib Date: 28-Apr-2023 23:03:35  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280018.D  
 Column 1 : Luna-Phenyl hexyl ( 4.60 mm) Det: LC mwd1A, 254 nm  
 Process Host: CTX1621

First Level Reviewer: LV5D

Date: 29-Apr-2023 10:26:41

| Compound                      | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 3 TNX                         | 1   | 5.082        | 5.115            | -0.033           | 197054   | 0.5015           | 0.4997             |       |
| 4 DNx                         | 1   | 5.922        | 5.955            | -0.033           | 145110   | 0.5005           | 0.5237             |       |
| 5 HMX                         | 1   | 6.576        | 6.595            | -0.019           | 78202    | 0.5000           | 0.4173             |       |
| 6 MNX                         | 1   | 7.456        | 7.489            | -0.033           | 154239   | 0.5835           | 0.5881             |       |
| 7 2,4,6-Trinitrophenol        | 1   | 7.936        | 8.029            | -0.093           | 78343    | 0.5000           | 0.4868             |       |
| 8 RDX                         | 1   | 8.882        | 8.915            | -0.033           | 97583    | 0.5000           | 0.4494             |       |
| 9 Nitrobenzene                | 1   | 11.449       | 11.482           | -0.033           | 186868   | 0.5000           | 0.4685             |       |
| \$ 10 1,2-Dinitrobenzene      | 1   | 12.436       | 12.469           | -0.033           | 130576   | 0.5000           | 0.4652             |       |
| 11 3,5-Dinitroaniline         | 1   | 14.322       | 14.355           | -0.033           | 212686   | 0.5000           | 0.4770             | M     |
| 12 1,3-Dinitrobenzene         | 1   | 14.589       | 14.622           | -0.033           | 288581   | 0.5000           | 0.4585             | M     |
| 13 Nitroglycerin              | 2   | 15.156       | 15.182           | -0.026           | 622685   | 5.00             | 4.62               |       |
| 14 o-Nitrotoluene             | 1   | 15.696       | 15.729           | -0.033           | 115843   | 0.5000           | 0.4691             | M     |
| 16 p-Nitrotoluene             | 1   | 15.922       | 15.955           | -0.033           | 108663   | 0.5000           | 0.4777             | M     |
| 17 4-Amino-2,6-dinitrotoluene | 1   | 16.482       | 16.515           | -0.033           | 136431   | 0.5000           | 0.4637             | M     |
| 18 m-Nitrotoluene             | 1   | 16.782       | 16.815           | -0.033           | 132594   | 0.5000           | 0.4500             | M     |
| 19 2-Amino-4,6-dinitrotoluene | 1   | 17.342       | 17.362           | -0.020           | 200216   | 0.5000           | 0.4679             | M     |
| 20 1,3,5-Trinitrobenzene      | 1   | 17.422       | 17.422           | 0.000            | 196279   | 0.5000           | 0.4838             | M     |
| 21 2,6-Dinitrotoluene         | 1   | 18.602       | 18.629           | -0.027           | 132469   | 0.5000           | 0.4542             |       |
| 22 2,4-Dinitrotoluene         | 1   | 19.056       | 19.082           | -0.026           | 263446   | 0.5000           | 0.4507             |       |
| 23 Tetryl                     | 1   | 22.329       | 22.389           | -0.060           | 161107   | 0.5000           | 0.5658             |       |
| 24 2,4,6-Trinitrotoluene      | 1   | 23.122       | 23.182           | -0.060           | 196312   | 0.5000           | 0.4473             |       |
| 25 PETN                       | 2   | 24.389       | 24.475           | -0.086           | 693219   | 5.00             | 4.82               |       |

**QC Flag Legend**

Processing Flags

Review Flags

M - Manually Integrated

**Reagents:**

|                     |                     |           |
|---------------------|---------------------|-----------|
| 8330Surrogate_00143 | Amount Added: 50.00 | Units: uL |
| 8330_OP_DMT_00017   | Amount Added: 50.00 | Units: uL |
| 8330 LCS_00126      | Amount Added: 50.00 | Units: uL |

Report Date: 29-Apr-2023 13:06:04

Chrom Revision: 2.3 29-Mar-2023 18:39:10

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280019.D

Injection Date: 28-Apr-2023 23:38:32

Instrument ID: CHHPLC\_X5

Operator ID: JZ/JG

Lims ID: ICV INT

Worklist Smp#: 19

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 19

Method: 8330\_X5\_Luna

Limit Group: GCSV - 8330


Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280019.D  
Injection Date: 28-Apr-2023 23:38:32 Instrument ID: CHHPLC\_X5  
Lims ID: ICV INT  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 19 Worklist Smp#: 19  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X5\_Luna Limit Group: GCSV - 8330  
Column: Luna-Phenyl hexyl ( 4.60 mm) Detector: LC mwd1A, 254 nm

**11 3,5-Dinitroaniline, CAS: 618-87-1**

Signal: 1


RT: 14.32  
Area: 204903  
Amount: 0.459536  
Amount Units: ug/ml

## Processing Integration Results



RT: 14.32  
Area: 212686  
Amount: 0.477043  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LV5D, 29-Apr-2023 10:27:10  
Audit Action: Assigned New Baseline

Audit Reason: Baseline Smoothing

## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280019.D  
Injection Date: 28-Apr-2023 23:38:32 Instrument ID: CHHPLC\_X5  
Lims ID: ICV INT  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 19 Worklist Smp#: 19  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X5\_Luna Limit Group: GCSV - 8330  
Column: Luna-Phenyl hexyl ( 4.60 mm) Detector LC mwd1A, 254 nm

**12 1,3-Dinitrobenzene, CAS: 99-65-0**

Signal: 1


RT: 14.59  
Area: 278067  
Amount: 0.441831  
Amount Units: ug/ml

## Processing Integration Results



RT: 14.59  
Area: 288581  
Amount: 0.458537  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LV5D, 29-Apr-2023 10:27:10  
Audit Action: Assigned New Baseline

Audit Reason: Baseline Smoothing

## Eurofins Denver

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280019.D

Injection Date: 28-Apr-2023 23:38:32

Instrument ID: CHHPLC\_X5

Lims ID: ICV INT

Client ID:

Operator ID: JZ/JG

ALS Bottle#:

19

Worklist Smp#: 19

Injection Vol: 100.0 ul

Dil. Factor: 1.0000

Method: 8330\_X5\_Luna

Limit Group: GCSV - 8330

Column: Luna-Phenyl hexyl ( 4.60 mm)


Detector: LC mwd1A, 254 nm

**14 o-Nitrotoluene, CAS: 88-72-2**

Signal: 1


RT: 15.70  
Area: 94233  
Amount: 0.381567  
Amount Units: ug/ml

## Processing Integration Results



RT: 15.70  
Area: 115843  
Amount: 0.469070  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LV5D, 29-Apr-2023 10:27:13

Audit Action: Assigned New Baseline

Audit Reason: Baseline Smoothing



## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280019.D  
Injection Date: 28-Apr-2023 23:38:32 Instrument ID: CHHPLC\_X5  
Lims ID: ICV INT  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 19 Worklist Smp#: 19  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X5\_Luna Limit Group: GCSV - 8330  
Column: Luna-Phenyl hexyl ( 4.60 mm) Detector LC mwd1A, 254 nm

**16 p-Nitrotoluene, CAS: 99-99-0**

Signal: 1


RT: 15.92  
Area: 83952  
Amount: 0.369036  
Amount Units: ug/ml

## Processing Integration Results



RT: 15.92  
Area: 108663  
Amount: 0.477660  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LV5D, 29-Apr-2023 10:27:13  
Audit Action: Assigned New Baseline

Audit Reason: Baseline Smoothing

## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280019.D  
Injection Date: 28-Apr-2023 23:38:32 Instrument ID: CHHPLC\_X5  
Lims ID: ICV INT  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 19 Worklist Smp#: 19  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X5\_Luna Limit Group: GCSV - 8330  
Column: Luna-Phenyl hexyl ( 4.60 mm) Detector LC mwd1A, 254 nm

**17 4-Amino-2,6-dinitrotoluene, CAS: 19406-51-0**

Signal: 1


RT: 16.48  
Area: 103981  
Amount: 0.365374  
Amount Units: ug/ml

## Processing Integration Results



RT: 16.48  
Area: 136431  
Amount: 0.463720  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LV5D, 29-Apr-2023 10:27:13  
Audit Action: Assigned New Baseline

Audit Reason: Baseline Smoothing

## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280019.D  
Injection Date: 28-Apr-2023 23:38:32 Instrument ID: CHHPLC\_X5  
Lims ID: ICV INT  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 19 Worklist Smp#: 19  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X5\_Luna Limit Group: GCSV - 8330  
Column: Luna-Phenyl hexyl ( 4.60 mm) Detector LC mwd1A, 254 nm

**18 m-Nitrotoluene, CAS: 99-08-1**

Signal: 1


RT: 16.78  
Area: 97090  
Amount: 0.342723  
Amount Units: ug/ml

## Processing Integration Results



RT: 16.78  
Area: 132594  
Amount: 0.450029  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LV5D, 29-Apr-2023 10:27:13  
Audit Action: Assigned New Baseline

Audit Reason: Baseline Smoothing

## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280019.D  
Injection Date: 28-Apr-2023 23:38:32 Instrument ID: CHHPLC\_X5  
Lims ID: ICV INT  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 19 Worklist Smp#: 19  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X5\_Luna Limit Group: GCSV - 8330  
Column: Luna-Phenyl hexyl ( 4.60 mm) Detector LC mwd1A, 254 nm

**20 1,3,5-Trinitrobenzene, CAS: 99-35-4**

Signal: 1


RT: 17.40  
Area: 328214  
Amount: 0.809022  
Amount Units: ug/ml

## Processing Integration Results



RT: 17.42  
Area: 196279  
Amount: 0.483812  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LV5D, 29-Apr-2023 10:27:01  
Audit Action: Split an Integrated Peak

Audit Reason: Baseline Smoothing

FORM VII  
HPLC/IC CONTINUING CALIBRATION DATA

Lab Name: Eurofins Denver Job No.: 280-176674-1  
SDG No.: \_\_\_\_\_  
Lab Sample ID: ICV 280-610603/28 Calibration Date: 04/29/2023 04:52  
Instrument ID: CHHPLC\_X5 Calib Start Date: 04/29/2023 00:13  
GC Column: Luna-phenylhex ID: 4.60 (mm) Calib End Date: 04/29/2023 04:18  
Lab File ID: 04280028.D Conc. Units: ug/L

| ANALYTE                    | CURVE<br>TYPE | AVE CF | CF     | MIN CF | CALC<br>AMOUNT | SPIKE<br>AMOUNT | %D   | MAX<br>%D |
|----------------------------|---------------|--------|--------|--------|----------------|-----------------|------|-----------|
| 2,6-diamino-4-nitrotoluene | Ave           | 403300 | 412518 |        | 511            | 500             | 2.3  | 20.0      |
| 2,4-diamino-6-nitrotoluene | Ave           | 225281 | 220158 |        | 489            | 500             | -2.3 | 20.0      |

FORM VII  
HPLC/IC CONTINUING CALIBRATION RETENTION TIME SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
SDG No.: \_\_\_\_\_  
Lab Sample ID: ICV 280-610603/28 Calibration Date: 04/29/2023 04:52  
Instrument ID: CHHPLC\_X5 Calib Start Date: 04/29/2023 00:13  
GC Column: Luna-phenylhex ID: 4.60 (mm) Calib End Date: 04/29/2023 04:18  
Lab File ID: 04280028.D

| Analyte                    | RT   | RT WINDOW |      |
|----------------------------|------|-----------|------|
|                            |      | FROM      | TO   |
| 2,6-diamino-4-nitrotoluene | 4.10 | 3.96      | 4.26 |
| 2,4-diamino-6-nitrotoluene | 4.64 | 4.50      | 4.80 |

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280028.D  
Lims ID: ICV ADD  
Client ID:  
Sample Type: ICV  
Inject. Date: 29-Apr-2023 04:52:55 ALS Bottle#: 28 Worklist Smp#: 28  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Sample Info: ICV ADD  
Operator ID: JZ/JG Instrument ID: CHHPLC\_X5  
Sublist:  
Method: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\8330\_X5\_Luna.m  
Limit Group: GCSV - 8330  
Last Update: 29-Apr-2023 13:04:23 Calib Date: 29-Apr-2023 04:18:01  
Integrator: Falcon  
Quant Method: External Standard Quant By: Initial Calibration  
Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280027.D  
Column 1 : Luna-Phenyl hexyl ( 4.60 mm) Det: LC mwd1A, 254 nm  
Process Host: CTX1621

First Level Reviewer: LV5D

Date: 29-Apr-2023 10:31:01

| Compound                     | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 1 2,6-diamino-4-nitrotoluene | 1   | 4.102        | 4.113            | -0.011           | 206259   | 0.5000           | 0.5114             |       |
| 2 2,4-diamino-6-nitrotoluene | 1   | 4.635        | 4.653            | -0.018           | 110079   | 0.5000           | 0.4886             |       |

## QC Flag Legend

Processing Flags

## Reagents:

8330DiaminLCS\_00056

Amount Added: 50.00

Units: uL

Report Date: 29-Apr-2023 13:06:04

Chrom Revision: 2.3 29-Mar-2023 18:39:10

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280028.D

Injection Date: 29-Apr-2023 04:52:55

Instrument ID: CHHPLC\_X5

Operator ID: JZ/JG

Lims ID: ICV ADD

Worklist Smp#: 28

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 28

Method: 8330\_X5\_Luna

Limit Group: GCSV - 8330


Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1





FORM VII  
HPLC/IC CONTINUING CALIBRATION DATA

Lab Name: Eurofins Denver Job No.: 280-176674-1  
 SDG No.: \_\_\_\_\_  
 Lab Sample ID: CCV 280-613238/7 Calibration Date: 05/20/2023 18:42  
 Instrument ID: CHHPLC\_X5 Calib Start Date: 04/28/2023 18:24  
 GC Column: Luna-phenylhex ID: 4.60 (mm) Calib End Date: 04/28/2023 23:03  
 Lab File ID: 05200007.D Conc. Units: ug/L

| ANALYTE                    | CURVE TYPE | AVE CF | CF     | MIN CF | CALC AMOUNT | SPIKE AMOUNT | %D   | MAX %D |
|----------------------------|------------|--------|--------|--------|-------------|--------------|------|--------|
| HMX                        | Lin2       |        | 193344 |        | 256         | 250          | 2.5  | 20.0   |
| Picric acid                | Lin2       |        | 165940 |        | 256         | 250          | 2.5  | 20.0   |
| RDX                        | Lin2       |        | 225160 |        | 258         | 250          | 3.2  | 20.0   |
| Nitrobenzene               | Ave        | 398834 | 399416 |        | 250         | 250          | 0.1  | 20.0   |
| 3,5-Dinitroaniline         | Lin2       |        | 474016 |        | 265         | 250          | 6.1  | 20.0   |
| 1,3-Dinitrobenzene         | Ave        | 629351 | 630392 |        | 250         | 250          | 0.2  | 20.0   |
| Nitroglycerin              | Ave        | 134735 | 136176 |        | 2530        | 2500         | 1.1  | 20.0   |
| 2-Nitrotoluene             | Ave        | 246963 | 255896 |        | 259         | 250          | 3.6  | 20.0   |
| 4-Nitrotoluene             | Ave        | 227490 | 235428 |        | 259         | 250          | 3.5  | 20.0   |
| 4-Amino-2,6-dinitrotoluene | Ave        | 294210 | 299492 |        | 254         | 250          | 1.8  | 20.0   |
| 3-Nitrotoluene             | Ave        | 294634 | 302560 |        | 257         | 250          | 2.7  | 20.0   |
| 2-Amino-4,6-dinitrotoluene | Ave        | 427939 | 426644 |        | 249         | 250          | -0.3 | 20.0   |
| 1,3,5-Trinitrobenzene      | Ave        | 405692 | 444348 |        | 274         | 250          | 9.5  | 20.0   |
| 2,6-Dinitrotoluene         | Ave        | 291675 | 295332 |        | 253         | 250          | 1.3  | 20.0   |
| 2,4-Dinitrotoluene         | Ave        | 584463 | 597064 |        | 255         | 250          | 2.2  | 20.0   |
| Tetryl                     | Ave        | 284731 | 290312 |        | 255         | 250          | 2.0  | 20.0   |
| 2,4,6-Trinitrotoluene      | Ave        | 438876 | 444388 |        | 253         | 250          | 1.3  | 20.0   |
| PETN                       | Ave        | 143938 | 135777 |        | 2360        | 2500         | -5.7 | 20.0   |
| 1,2-Dinitrobenzene         | Ave        | 280697 | 280972 |        | 250         | 250          | 0.1  | 20.0   |

FORM VII  
HPLC/IC CONTINUING CALIBRATION RETENTION TIME SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
 SDG No.: \_\_\_\_\_  
 Lab Sample ID: CCV 280-613238/7 Calibration Date: 05/20/2023 18:42  
 Instrument ID: CHHPLC\_X5 Calib Start Date: 04/28/2023 18:24  
 GC Column: Luna-phenylhex ID: 4.60 (mm) Calib End Date: 04/28/2023 23:03  
 Lab File ID: 05200007.D

| Analyte                    | RT    | RT WINDOW |       |
|----------------------------|-------|-----------|-------|
|                            |       | FROM      | TO    |
| HMX                        | 6.76  | 6.61      | 6.91  |
| Picric acid                | 8.49  | 8.34      | 8.64  |
| RDX                        | 9.01  | 8.86      | 9.16  |
| Nitrobenzene               | 11.53 | 11.38     | 11.68 |
| 3,5-Dinitroaniline         | 14.38 | 14.23     | 14.53 |
| 1,3-Dinitrobenzene         | 14.64 | 14.49     | 14.79 |
| Nitroglycerin              | 15.16 | 15.01     | 15.31 |
| 2-Nitrotoluene             | 15.72 | 15.57     | 15.87 |
| 4-Nitrotoluene             | 15.95 | 15.80     | 16.10 |
| 4-Amino-2,6-dinitrotoluene | 16.51 | 16.36     | 16.66 |
| 3-Nitrotoluene             | 16.80 | 16.65     | 16.95 |
| 2-Amino-4,6-dinitrotoluene | 17.35 | 17.20     | 17.50 |
| 1,3,5-Trinitrobenzene      | 17.40 | 17.25     | 17.55 |
| 2,6-Dinitrotoluene         | 18.61 | 18.46     | 18.76 |
| 2,4-Dinitrotoluene         | 19.05 | 18.90     | 19.20 |
| Tetryl                     | 22.33 | 22.18     | 22.48 |
| 2,4,6-Trinitrotoluene      | 23.13 | 22.98     | 23.28 |
| PETN                       | 24.37 | 24.22     | 24.52 |
| 1,2-Dinitrobenzene         | 12.51 | 12.36     | 12.66 |

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230520-121670.b\05200007.D  
 Lims ID: CCV  
 Client ID:  
 Sample Type: CCV  
 Inject. Date: 20-May-2023 18:42:13 ALS Bottle#: 7 Worklist Smp#: 7  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: CCV  
 Operator ID: JZ/JG Instrument ID: CHHPLC\_X5  
 Sublist: chrom-8330\_X5\_Luna\*sub7  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230520-121670.b\8330\_X5\_Luna.m  
 Limit Group: GCSV - 8330  
 Last Update: 23-May-2023 13:12:40 Calib Date: 29-Apr-2023 04:18:01  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280027.D  
 Column 1 : Luna-Phenyl hexyl ( 4.60 mm) Det: LC mwd1A, 254 nm  
 Process Host: CTX1626

First Level Reviewer: LV5D

Date: 23-May-2023 12:44:07

| Compound                      | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 5 HMX                         | 1   | 6.758        | 6.758            | 0.000            | 48336    | 0.2500           | 0.2562             |       |
| 7 2,4,6-Trinitrophenol        | 1   | 8.491        | 8.491            | 0.000            | 41485    | 0.2500           | 0.2562             |       |
| 8 RDX                         | 1   | 9.005        | 9.005            | 0.000            | 56290    | 0.2500           | 0.2579             |       |
| 9 Nitrobenzene                | 1   | 11.531       | 11.531           | 0.000            | 99854    | 0.2500           | 0.2504             |       |
| \$ 10 1,2-Dinitrobenzene      | 1   | 12.511       | 12.511           | 0.000            | 70243    | 0.2500           | 0.2502             |       |
| 11 3,5-Dinitroaniline         | 1   | 14.378       | 14.378           | 0.000            | 118504   | 0.2500           | 0.2652             |       |
| 12 1,3-Dinitrobenzene         | 1   | 14.638       | 14.638           | 0.000            | 157598   | 0.2500           | 0.2504             |       |
| 13 Nitroglycerin              | 2   | 15.158       | 15.158           | 0.000            | 340440   | 2.50             | 2.53               |       |
| 14 o-Nitrotoluene             | 1   | 15.718       | 15.718           | 0.000            | 63974    | 0.2500           | 0.2590             |       |
| 16 p-Nitrotoluene             | 1   | 15.945       | 15.945           | 0.000            | 58857    | 0.2500           | 0.2587             |       |
| 17 4-Amino-2,6-dinitrotoluene | 1   | 16.505       | 16.505           | 0.000            | 74873    | 0.2500           | 0.2545             |       |
| 18 m-Nitrotoluene             | 1   | 16.798       | 16.798           | 0.000            | 75640    | 0.2500           | 0.2567             | a     |
| 19 2-Amino-4,6-dinitrotoluene | 1   | 17.351       | 17.351           | 0.000            | 106661   | 0.2500           | 0.2492             | Ma    |
| 20 1,3,5-Trinitrobenzene      | 1   | 17.398       | 17.398           | 0.000            | 111087   | 0.2500           | 0.2738             | Ma    |
| 21 2,6-Dinitrotoluene         | 1   | 18.605       | 18.605           | 0.000            | 73833    | 0.2500           | 0.2531             |       |
| 22 2,4-Dinitrotoluene         | 1   | 19.051       | 19.051           | 0.000            | 149266   | 0.2500           | 0.2554             |       |
| 23 Tetryl                     | 1   | 22.331       | 22.331           | 0.000            | 72578    | 0.2500           | 0.2549             | M     |
| 24 2,4,6-Trinitrotoluene      | 1   | 23.125       | 23.125           | 0.000            | 111097   | 0.2500           | 0.2531             | M     |
| 25 PETN                       | 2   | 24.365       | 24.365           | 0.000            | 339443   | 2.50             | 2.36               |       |

**QC Flag Legend**

Processing Flags

Review Flags

M - Manually Integrated

a - User Assigned ID

**Reagents:**

8330IntermStk\_00076

Amount Added: 25.00

Units: uL

## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230520-121670.b\05200007.D  
Injection Date: 20-May-2023 18:42:13 Instrument ID: CHHPLC\_X5  
Lims ID: CCV  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 7 Worklist Smp#: 7  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X5\_Luna Limit Group: GCSV - 8330  
Column: Luna-Phenyl hexyl ( 4.60 mm) Detector LC mwd1A, 254 nm

**18 m-Nitrotoluene, CAS: 99-08-1**

Signal: 1


RT: 16.80  
Area: 75640  
Amount: 0.256725  
Amount Units: ug/ml

## Processing Integration Results



RT: 16.80  
Area: 75640  
Amount: 0.256725  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: K8YG, 22-May-2023 11:10:41 -06:00:00 (UTC)

Audit Action: Assigned Compound ID

Audit Reason: Split Peak

## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230520-121670.b\05200007.D  
Injection Date: 20-May-2023 18:42:13 Instrument ID: CHHPLC\_X5  
Lims ID: CCV  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 7 Worklist Smp#: 7  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X5\_Luna Limit Group: GCSV - 8330  
Column: Luna-Phenyl hexyl ( 4.60 mm) Detector LC mwd1A, 254 nm

**19 2-Amino-4,6-dinitrotoluene, CAS: 35572-78-2**

Signal: 1


Not Detected  
Expected RT: 17.35

## Processing Integration Results



RT: 17.35  
Area: 106661  
Amount: 0.249243  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: K8YG, 22-May-2023 11:11:04 -06:00:00 (UTC)

Audit Action: Split an Integrated Peak

Audit Reason: Split Peak

## Eurofins Denver

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230520-121670.b\05200007.D

Injection Date: 20-May-2023 18:42:13

Instrument ID: CHHPLC\_X5

Lims ID: CCV

Client ID:

Operator ID: JZ/JG

ALS Bottle#:

7

Worklist Smp#: 7

Injection Vol: 100.0 ul

Dil. Factor: 1.0000

Method: 8330\_X5\_Luna

Limit Group: GCSV - 8330

Column: Luna-Phenyl hexyl ( 4.60 mm)


Detector LC mwd1A, 254 nm

**20 1,3,5-Trinitrobenzene, CAS: 99-35-4**

Signal: 1


RT: 17.35  
Area: 217742  
Amount: 0.536717  
Amount Units: ug/ml

## Processing Integration Results



RT: 17.40  
Area: 111087  
Amount: 0.273821  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: K8YG, 22-May-2023 11:11:10 -06:00:00 (UTC)

Audit Action: Manually Integrated/Assigned Compound ID Audit Reason: Split Peak

## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230520-121670.b\05200007.D  
Injection Date: 20-May-2023 18:42:13 Instrument ID: CHHPLC\_X5  
Lims ID: CCV  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 7 Worklist Smp#: 7  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X5\_Luna Limit Group: GCSV - 8330  
Column: Luna-Phenyl hexyl ( 4.60 mm) Detector LC mwd1A, 254 nm

**23 Tetryl, CAS: 479-45-8**

Signal: 1


RT: 22.33  
Area: 69631  
Amount: 0.244550  
Amount Units: ug/ml

## Processing Integration Results



RT: 22.33  
Area: 72578  
Amount: 0.254900  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: K8YG, 22-May-2023 11:11:23 -06:00:00 (UTC)

Audit Action: Assigned New Baseline

Audit Reason: Split Peak

## Eurofins Denver

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230520-121670.b\05200007.D

Injection Date: 20-May-2023 18:42:13

Instrument ID: CHHPLC\_X5

Lims ID: CCV

Client ID:

Operator ID: JZ/JG

ALS Bottle#:

7

Worklist Smp#: 7

Injection Vol: 100.0 ul

Dil. Factor: 1.0000

Method: 8330\_X5\_Luna

Limit Group: GCSV - 8330

Column: Luna-Phenyl hexyl ( 4.60 mm)


Detector LC mwd1A, 254 nm

**24 2,4,6-Trinitrotoluene, CAS: 118-96-7**

Signal: 1


RT: 23.12  
Area: 107588  
Amount: 0.245145  
Amount Units: ug/ml

## Processing Integration Results



RT: 23.12  
Area: 111097  
Amount: 0.253140  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: K8YG, 22-May-2023 11:11:23 -06:00:00 (UTC)

Audit Action: Assigned New Baseline

Audit Reason: Split Peak



FORM VII  
HPLC/IC CONTINUING CALIBRATION DATA

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Lab Sample ID: CCV 280-613238/20 Calibration Date: 05/21/2023 01:06

Instrument ID: CHHPLC\_X5 Calib Start Date: 04/28/2023 18:24

GC Column: Luna-phenylhex ID: 4.60 (mm) Calib End Date: 04/28/2023 23:03

Lab File ID: 05200020.D Conc. Units: ug/L

| ANALYTE                    | CURVE<br>TYPE | AVE CF | CF     | MIN CF | CALC<br>AMOUNT | SPIKE<br>AMOUNT | %D   | MAX<br>%D |
|----------------------------|---------------|--------|--------|--------|----------------|-----------------|------|-----------|
| HMX                        | Lin2          |        | 195372 |        | 259            | 250             | 3.6  | 20.0      |
| Picric acid                | Lin2          |        | 167800 |        | 259            | 250             | 3.6  | 20.0      |
| RDX                        | Lin2          |        | 229832 |        | 263            | 250             | 5.3  | 20.0      |
| Nitrobenzene               | Ave           | 398834 | 395272 |        | 248            | 250             | -0.9 | 20.0      |
| 3,5-Dinitroaniline         | Lin2          |        | 473648 |        | 265            | 250             | 6.0  | 20.0      |
| 1,3-Dinitrobenzene         | Ave           | 629351 | 624272 |        | 248            | 250             | -0.8 | 20.0      |
| Nitroglycerin              | Ave           | 134735 | 140513 |        | 2610           | 2500            | 4.3  | 20.0      |
| 2-Nitrotoluene             | Ave           | 246963 | 250940 |        | 254            | 250             | 1.6  | 20.0      |
| 4-Nitrotoluene             | Ave           | 227490 | 229120 |        | 252            | 250             | 0.7  | 20.0      |
| 4-Amino-2,6-dinitrotoluene | Ave           | 294210 | 296208 |        | 252            | 250             | 0.7  | 20.0      |
| 3-Nitrotoluene             | Ave           | 294634 | 291452 |        | 247            | 250             | -1.1 | 20.0      |
| 2-Amino-4,6-dinitrotoluene | Ave           | 427939 | 423056 |        | 247            | 250             | -1.1 | 20.0      |
| 1,3,5-Trinitrobenzene      | Ave           | 405692 | 416692 |        | 257            | 250             | 2.7  | 20.0      |
| 2,6-Dinitrotoluene         | Ave           | 291675 | 294404 |        | 252            | 250             | 0.9  | 20.0      |
| 2,4-Dinitrotoluene         | Ave           | 584463 | 597332 |        | 256            | 250             | 2.2  | 20.0      |
| Tetryl                     | Ave           | 284731 | 290192 |        | 255            | 250             | 1.9  | 20.0      |
| 2,4,6-Trinitrotoluene      | Ave           | 438876 | 442008 |        | 252            | 250             | 0.7  | 20.0      |
| PETN                       | Ave           | 143938 | 143362 |        | 2490           | 2500            | -0.4 | 20.0      |
| 1,2-Dinitrobenzene         | Ave           | 280697 | 281760 |        | 251            | 250             | 0.4  | 20.0      |

FORM VII  
HPLC/IC CONTINUING CALIBRATION RETENTION TIME SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
 SDG No.: \_\_\_\_\_  
 Lab Sample ID: CCV 280-613238/20 Calibration Date: 05/21/2023 01:06  
 Instrument ID: CHHPLC\_X5 Calib Start Date: 04/28/2023 18:24  
 GC Column: Luna-phenylhex ID: 4.60 (mm) Calib End Date: 04/28/2023 23:03  
 Lab File ID: 05200020.D

| Analyte                    | RT    | RT WINDOW |       |
|----------------------------|-------|-----------|-------|
|                            |       | FROM      | TO    |
| HMX                        | 6.75  | 6.61      | 6.91  |
| Picric acid                | 8.48  | 8.34      | 8.64  |
| RDX                        | 8.99  | 8.86      | 9.16  |
| Nitrobenzene               | 11.49 | 11.38     | 11.68 |
| 3,5-Dinitroaniline         | 14.32 | 14.23     | 14.53 |
| 1,3-Dinitrobenzene         | 14.58 | 14.49     | 14.79 |
| Nitroglycerin              | 15.09 | 15.01     | 15.31 |
| 2-Nitrotoluene             | 15.65 | 15.57     | 15.87 |
| 4-Nitrotoluene             | 15.87 | 15.80     | 16.10 |
| 4-Amino-2,6-dinitrotoluene | 16.43 | 16.36     | 16.66 |
| 3-Nitrotoluene             | 16.71 | 16.65     | 16.95 |
| 2-Amino-4,6-dinitrotoluene | 17.25 | 17.20     | 17.50 |
| 1,3,5-Trinitrobenzene      | 17.37 | 17.25     | 17.55 |
| 2,6-Dinitrotoluene         | 18.52 | 18.46     | 18.76 |
| 2,4-Dinitrotoluene         | 18.97 | 18.90     | 19.20 |
| Tetryl                     | 22.25 | 22.18     | 22.48 |
| 2,4,6-Trinitrotoluene      | 23.05 | 22.98     | 23.28 |
| PETN                       | 24.28 | 24.22     | 24.52 |
| 1,2-Dinitrobenzene         | 12.45 | 12.36     | 12.66 |

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230520-121670.b\05200020.D  
 Lims ID: CCV  
 Client ID:  
 Sample Type: CCV  
 Inject. Date: 21-May-2023 01:06:23 ALS Bottle#: 7 Worklist Smp#: 20  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: CCV  
 Operator ID: JZ/JG Instrument ID: CHHPLC\_X5  
 Sublist: chrom-8330\_X5\_Luna\*sub7  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230520-121670.b\8330\_X5\_Luna.m  
 Limit Group: GCSV - 8330  
 Last Update: 23-May-2023 13:12:54 Calib Date: 29-Apr-2023 04:18:01  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280027.D  
 Column 1 : Luna-Phenyl hexyl ( 4.60 mm) Det: LC mwd1A, 254 nm  
 Process Host: CTX1626

First Level Reviewer: K8YG

Date: 22-May-2023 11:27:42

| Compound                      | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 5 HMX                         | 1   | 6.745        | 6.758            | -0.013           | 48843    | 0.2500           | 0.2589             |       |
| 7 2,4,6-Trinitrophenol        | 1   | 8.478        | 8.491            | -0.013           | 41950    | 0.2500           | 0.2591             |       |
| 8 RDX                         | 1   | 8.992        | 9.005            | -0.013           | 57458    | 0.2500           | 0.2633             |       |
| 9 Nitrobenzene                | 1   | 11.485       | 11.531           | -0.046           | 98818    | 0.2500           | 0.2478             |       |
| \$ 10 1,2-Dinitrobenzene      | 1   | 12.452       | 12.511           | -0.059           | 70440    | 0.2500           | 0.2509             |       |
| 11 3,5-Dinitroaniline         | 1   | 14.318       | 14.378           | -0.060           | 118412   | 0.2500           | 0.2650             | M     |
| 12 1,3-Dinitrobenzene         | 1   | 14.578       | 14.638           | -0.060           | 156068   | 0.2500           | 0.2480             | M     |
| 13 Nitroglycerin              | 2   | 15.085       | 15.158           | -0.073           | 351283   | 2.50             | 2.61               |       |
| 14 o-Nitrotoluene             | 1   | 15.645       | 15.718           | -0.073           | 62735    | 0.2500           | 0.2540             | M     |
| 16 p-Nitrotoluene             | 1   | 15.872       | 15.945           | -0.073           | 57280    | 0.2500           | 0.2518             | M     |
| 17 4-Amino-2,6-dinitrotoluene | 1   | 16.425       | 16.505           | -0.080           | 74052    | 0.2500           | 0.2517             | M     |
| 18 m-Nitrotoluene             | 1   | 16.712       | 16.798           | -0.086           | 72863    | 0.2500           | 0.2473             | Ma    |
| 19 2-Amino-4,6-dinitrotoluene | 1   | 17.252       | 17.351           | -0.099           | 105764   | 0.2500           | 0.2471             | Ma    |
| 20 1,3,5-Trinitrobenzene      | 1   | 17.365       | 17.398           | -0.033           | 104173   | 0.2500           | 0.2568             | Ma    |
| 21 2,6-Dinitrotoluene         | 1   | 18.518       | 18.605           | -0.087           | 73601    | 0.2500           | 0.2523             |       |
| 22 2,4-Dinitrotoluene         | 1   | 18.972       | 19.051           | -0.079           | 149333   | 0.2500           | 0.2555             |       |
| 23 Tetryl                     | 1   | 22.252       | 22.331           | -0.079           | 72548    | 0.2500           | 0.2548             |       |
| 24 2,4,6-Trinitrotoluene      | 1   | 23.052       | 23.125           | -0.073           | 110502   | 0.2500           | 0.2518             |       |
| 25 PETN                       | 2   | 24.278       | 24.365           | -0.087           | 358404   | 2.50             | 2.49               |       |

**QC Flag Legend**

Processing Flags

Review Flags

M - Manually Integrated

a - User Assigned ID

**Reagents:**

8330IntermStk\_00076

Amount Added: 25.00

Units: uL

Report Date: 23-May-2023 13:12:54

Chrom Revision: 2.3 16-May-2023 16:37:30

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230520-121670.b\05200020.D

Injection Date: 21-May-2023 01:06:23

Instrument ID: CHHPLC\_X5

Operator ID: JZ/JG

Lims ID: CCV

Worklist Smp#: 20

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 7

Method: 8330\_X5\_Luna

Limit Group: GCSV - 8330


Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230520-121670.b\05200020.D  
Injection Date: 21-May-2023 01:06:23 Instrument ID: CHHPLC\_X5  
Lims ID: CCV  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 7 Worklist Smp#: 20  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X5\_Luna Limit Group: GCSV - 8330  
Column: Luna-Phenyl hexyl ( 4.60 mm) Detector: LC mwd1A, 254 nm

**11 3,5-Dinitroaniline, CAS: 618-87-1**

Signal: 1


RT: 14.32  
Area: 116911  
Amount: 0.261610  
Amount Units: ug/ml

## Processing Integration Results



RT: 14.32  
Area: 118412  
Amount: 0.264986  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: K8YG, 22-May-2023 11:25:27 -06:00:00 (UTC)

Audit Action: Assigned New Baseline

Audit Reason: Split Peak

## Eurofins Denver

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230520-121670.b\05200020.D

Injection Date: 21-May-2023 01:06:23

Instrument ID: CHHPLC\_X5

Lims ID: CCV

Client ID:

Operator ID: JZ/JG

ALS Bottle#:

7

Worklist Smp#: 20

Injection Vol: 100.0 ul

Dil. Factor: 1.0000

Method: 8330\_X5\_Luna

Limit Group: GCSV - 8330

Column: Luna-Phenyl hexyl ( 4.60 mm)


Detector LC mwd1A, 254 nm

**12 1,3-Dinitrobenzene, CAS: 99-65-0**

Signal: 1


RT: 14.58  
Area: 153969  
Amount: 0.244647  
Amount Units: ug/ml

## Processing Integration Results



RT: 14.58  
Area: 156068  
Amount: 0.247982  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: K8YG, 22-May-2023 11:25:27 -06:00:00 (UTC)

Audit Action: Assigned New Baseline

Audit Reason: Split Peak

## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230520-121670.b\05200020.D  
Injection Date: 21-May-2023 01:06:23 Instrument ID: CHHPLC\_X5  
Lims ID: CCV  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 7 Worklist Smp#: 20  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X5\_Luna Limit Group: GCSV - 8330  
Column: Luna-Phenyl hexyl ( 4.60 mm) Detector LC mwd1A, 254 nm

**14 o-Nitrotoluene, CAS: 88-72-2**

Signal: 1


RT: 15.64  
Area: 55718  
Amount: 0.225612  
Amount Units: ug/ml

## Processing Integration Results



RT: 15.64  
Area: 62735  
Amount: 0.254026  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: K8YG, 22-May-2023 11:25:22 -06:00:00 (UTC)

Audit Action: Assigned New Baseline

Audit Reason: Split Peak

## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230520-121670.b\05200020.D  
Injection Date: 21-May-2023 01:06:23 Instrument ID: CHHPLC\_X5  
Lims ID: CCV  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 7 Worklist Smp#: 20  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X5\_Luna Limit Group: GCSV - 8330  
Column: Luna-Phenyl hexyl ( 4.60 mm) Detector LC mwd1A, 254 nm

**16 p-Nitrotoluene, CAS: 99-99-0**

Signal: 1


RT: 15.87  
Area: 49119  
Amount: 0.215917  
Amount Units: ug/ml

## Processing Integration Results



RT: 15.87  
Area: 57280  
Amount: 0.251791  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: K8YG, 22-May-2023 11:25:22 -06:00:00 (UTC)

Audit Action: Assigned New Baseline

Audit Reason: Split Peak



## Eurofins Denver

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230520-121670.b\05200020.D

Injection Date: 21-May-2023 01:06:23

Instrument ID: CHHPLC\_X5

Lims ID: CCV

Client ID:

Operator ID: JZ/JG

ALS Bottle#:

7

Worklist Smp#: 20

Injection Vol: 100.0 ul

Dil. Factor:

1.0000

Method: 8330\_X5\_Luna

Limit Group:

GCSV - 8330

Column: Luna-Phenyl hexyl ( 4.60 mm)

Detector


LC mwd1A, 254 nm

**17 4-Amino-2,6-dinitrotoluene, CAS: 19406-51-0**

Signal: 1


RT: 16.42  
Area: 63437  
Amount: 0.215618  
Amount Units: ug/ml

## Processing Integration Results



RT: 16.42  
Area: 74052  
Amount: 0.251698  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: K8YG, 22-May-2023 11:25:22 -06:00:00 (UTC)

Audit Action: Assigned New Baseline

Audit Reason: Split Peak

## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230520-121670.b\05200020.D  
Injection Date: 21-May-2023 01:06:23 Instrument ID: CHHPLC\_X5  
Lims ID: CCV  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 7 Worklist Smp#: 20  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X5\_Luna Limit Group: GCSV - 8330  
Column: Luna-Phenyl hexyl ( 4.60 mm) Detector LC mwd1A, 254 nm

**18 m-Nitrotoluene, CAS: 99-08-1**

Signal: 1


RT: 16.71  
Area: 60099  
Amount: 0.203978  
Amount Units: ug/ml

## Processing Integration Results



RT: 16.71  
Area: 72863  
Amount: 0.247300  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: K8YG, 22-May-2023 11:27:33 -06:00:00 (UTC)

Audit Action: Manually Integrated/Assigned Compound ID Audit Reason: Split Peak

## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230520-121670.b\05200020.D  
Injection Date: 21-May-2023 01:06:23 Instrument ID: CHHPLC\_X5  
Lims ID: CCV  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 7 Worklist Smp#: 20  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X5\_Luna Limit Group: GCSV - 8330  
Column: Luna-Phenyl hexyl ( 4.60 mm) Detector LC mwd1A, 254 nm

**19 2-Amino-4,6-dinitrotoluene, CAS: 35572-78-2**

Signal: 1


RT: 17.25  
Area: 190797  
Amount: 0.445850  
Amount Units: ug/ml

## Processing Integration Results



RT: 17.25  
Area: 105764  
Amount: 0.247147  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LV5D, 23-May-2023 12:44:14 -06:00:00 (UTC)

Audit Action: Manually Integrated/Assigned Compound ID Audit Reason:

## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230520-121670.b\05200020.D  
Injection Date: 21-May-2023 01:06:23 Instrument ID: CHHPLC\_X5  
Lims ID: CCV  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 7 Worklist Smp#: 20  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X5\_Luna Limit Group: GCSV - 8330  
Column: Luna-Phenyl hexyl ( 4.60 mm) Detector LC mwd1A, 254 nm

20 1,3,5-Trinitrobenzene, CAS: 99-35-4

Signal: 1


Not Detected  
Expected RT: 17.40

## Processing Integration Results



RT: 17.36  
Area: 104173  
Amount: 0.256778  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: K8YG, 22-May-2023 11:26:20 -06:00:00 (UTC)

Audit Action: Split an Integrated Peak

Audit Reason: Split Peak

FORM VII  
HPLC/IC CONTINUING CALIBRATION DATA

Lab Name: Eurofins Denver Job No.: 280-176674-1  
 SDG No.: \_\_\_\_\_  
 Lab Sample ID: CCV 280-613683/8 Calibration Date: 05/24/2023 19:13  
 Instrument ID: CHHPLC\_X5 Calib Start Date: 04/28/2023 18:24  
 GC Column: Luna-phenylhex ID: 4.60 (mm) Calib End Date: 04/28/2023 23:03  
 Lab File ID: 05240008.D Conc. Units: ug/L

| ANALYTE                    | CURVE<br>TYPE | AVE CF | CF     | MIN CF | CALC<br>AMOUNT | SPIKE<br>AMOUNT | %D   | MAX<br>%D |
|----------------------------|---------------|--------|--------|--------|----------------|-----------------|------|-----------|
| TNX                        | Ave           | 394332 | 368052 |        | 234            | 251             | -6.7 | 20.0      |
| DNX                        | Ave           | 277072 | 275600 |        | 249            | 250             | -0.5 | 20.0      |
| HMX                        | Lin2          |        | 196824 |        | 261            | 250             | 4.4  | 20.0      |
| MNX                        | Ave           | 262267 | 253656 |        | 282            | 292             | -3.3 | 20.0      |
| Picric acid                | Lin2          |        | 165684 |        | 256            | 250             | 2.3  | 20.0      |
| RDX                        | Lin2          |        | 228684 |        | 262            | 250             | 4.8  | 20.0      |
| Nitrobenzene               | Ave           | 398834 | 398688 |        | 250            | 250             | -0.0 | 20.0      |
| 3,5-Dinitroaniline         | Lin2          |        | 469884 |        | 263            | 250             | 5.1  | 20.0      |
| 1,3-Dinitrobenzene         | Ave           | 629351 | 624260 |        | 248            | 250             | -0.8 | 20.0      |
| Nitroglycerin              | Ave           | 134735 | 137371 |        | 2550           | 2500            | 2.0  | 20.0      |
| 2-Nitrotoluene             | Ave           | 246963 | 253404 |        | 257            | 250             | 2.6  | 20.0      |
| 4-Nitrotoluene             | Ave           | 227490 | 233124 |        | 256            | 250             | 2.5  | 20.0      |
| 4-Amino-2,6-dinitrotoluene | Ave           | 294210 | 296312 |        | 252            | 250             | 0.7  | 20.0      |
| 3-Nitrotoluene             | Ave           | 294634 | 299556 |        | 254            | 250             | 1.7  | 20.0      |
| 2-Amino-4,6-dinitrotoluene | Ave           | 427939 | 430188 |        | 251            | 250             | 0.5  | 20.0      |
| 1,3,5-Trinitrobenzene      | Ave           | 405692 | 423372 |        | 261            | 250             | 4.4  | 20.0      |
| 2,6-Dinitrotoluene         | Ave           | 291675 | 312336 |        | 268            | 250             | 7.1  | 20.0      |
| 2,4-Dinitrotoluene         | Ave           | 584463 | 590292 |        | 252            | 250             | 1.0  | 20.0      |
| Tetryl                     | Ave           | 284731 | 285060 |        | 250            | 250             | 0.1  | 20.0      |
| 2,4,6-Trinitrotoluene      | Ave           | 438876 | 438628 |        | 250            | 250             | -0.0 | 20.0      |
| PETN                       | Ave           | 143938 | 141864 |        | 2460           | 2500            | -1.4 | 20.0      |
| 1,2-Dinitrobenzene         | Ave           | 280697 | 279280 |        | 249            | 250             | -0.5 | 20.0      |

FORM VII  
HPLC/IC CONTINUING CALIBRATION RETENTION TIME SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
 SDG No.: \_\_\_\_\_  
 Lab Sample ID: CCV 280-613683/8 Calibration Date: 05/24/2023 19:13  
 Instrument ID: CHHPLC\_X5 Calib Start Date: 04/28/2023 18:24  
 GC Column: Luna-phenylhex ID: 4.60 (mm) Calib End Date: 04/28/2023 23:03  
 Lab File ID: 05240008.D

| Analyte                    | RT    | RT WINDOW |       |
|----------------------------|-------|-----------|-------|
|                            |       | FROM      | TO    |
| TNX                        | 5.21  | 5.06      | 5.36  |
| DNX                        | 6.05  | 5.90      | 6.20  |
| HMX                        | 6.69  | 6.54      | 6.84  |
| MNX                        | 7.55  | 7.40      | 7.70  |
| Picric acid                | 8.54  | 8.39      | 8.69  |
| RDX                        | 8.95  | 8.80      | 9.10  |
| Nitrobenzene               | 11.46 | 11.31     | 11.61 |
| 3,5-Dinitroaniline         | 14.29 | 14.14     | 14.44 |
| 1,3-Dinitrobenzene         | 14.55 | 14.40     | 14.70 |
| Nitroglycerin              | 15.07 | 14.92     | 15.22 |
| 2-Nitrotoluene             | 15.62 | 15.47     | 15.77 |
| 4-Nitrotoluene             | 15.85 | 15.70     | 16.00 |
| 4-Amino-2,6-dinitrotoluene | 16.41 | 16.26     | 16.56 |
| 3-Nitrotoluene             | 16.69 | 16.54     | 16.84 |
| 2-Amino-4,6-dinitrotoluene | 17.24 | 17.09     | 17.39 |
| 1,3,5-Trinitrobenzene      | 17.33 | 17.18     | 17.48 |
| 2,6-Dinitrotoluene         | 18.49 | 18.34     | 18.64 |
| 2,4-Dinitrotoluene         | 18.95 | 18.80     | 19.10 |
| Tetryl                     | 22.24 | 22.09     | 22.39 |
| 2,4,6-Trinitrotoluene      | 23.04 | 22.89     | 23.19 |
| PETN                       | 24.28 | 24.13     | 24.43 |
| 1,2-Dinitrobenzene         | 12.43 | 12.28     | 12.58 |

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230524-121801.b\05240008.D  
 Lims ID: CCV DMT  
 Client ID:  
 Sample Type: CCV  
 Inject. Date: 24-May-2023 19:13:25 ALS Bottle#: 7 Worklist Smp#: 8  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: CCV DMT  
 Operator ID: JZ/JG Instrument ID: CHHPLC\_X5  
 Sublist: chrom-8330\_X5\_Luna\*sub6  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230524-121801.b\8330\_X5\_Luna.m  
 Limit Group: GCSV - 8330  
 Last Update: 25-May-2023 13:06:09 Calib Date: 29-Apr-2023 04:18:01  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280027.D  
 Column 1 : Luna-Phenyl hexyl ( 4.60 mm) Det: LC mwd1A, 254 nm  
 Process Host: CTX1654

First Level Reviewer: LV5D

Date: 24-May-2023 19:54:18

| Compound                      | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 3 TNX                         | 1   | 5.206        | 5.206            | 0.000            | 92289    | 0.2508           | 0.2340             |       |
| 4 DNx                         | 1   | 6.046        | 6.046            | 0.000            | 68969    | 0.2503           | 0.2489             |       |
| 5 HMX                         | 1   | 6.686        | 6.686            | 0.000            | 49206    | 0.2500           | 0.2609             |       |
| 6 MNX                         | 1   | 7.546        | 7.546            | 0.000            | 74004    | 0.2918           | 0.2822             |       |
| 7 2,4,6-Trinitrophenol        | 1   | 8.539        | 8.539            | 0.000            | 41421    | 0.2500           | 0.2558             |       |
| 8 RDX                         | 1   | 8.946        | 8.946            | 0.000            | 57171    | 0.2500           | 0.2620             |       |
| 9 Nitrobenzene                | 1   | 11.459       | 11.459           | 0.000            | 99672    | 0.2500           | 0.2499             |       |
| \$ 10 1,2-Dinitrobenzene      | 1   | 12.426       | 12.426           | 0.000            | 69820    | 0.2500           | 0.2487             |       |
| 11 3,5-Dinitroaniline         | 1   | 14.292       | 14.292           | 0.000            | 117471   | 0.2500           | 0.2629             |       |
| 12 1,3-Dinitrobenzene         | 1   | 14.552       | 14.552           | 0.000            | 156065   | 0.2500           | 0.2480             |       |
| 13 Nitroglycerin              | 2   | 15.066       | 15.066           | 0.000            | 343427   | 2.50             | 2.55               |       |
| 14 o-Nitrotoluene             | 1   | 15.619       | 15.619           | 0.000            | 63351    | 0.2500           | 0.2565             |       |
| 16 p-Nitrotoluene             | 1   | 15.846       | 15.846           | 0.000            | 58281    | 0.2500           | 0.2562             |       |
| 17 4-Amino-2,6-dinitrotoluene | 1   | 16.406       | 16.406           | 0.000            | 74078    | 0.2500           | 0.2518             |       |
| 18 m-Nitrotoluene             | 1   | 16.692       | 16.692           | 0.000            | 74889    | 0.2500           | 0.2542             |       |
| 19 2-Amino-4,6-dinitrotoluene | 1   | 17.239       | 17.239           | 0.000            | 107547   | 0.2500           | 0.2513             | M     |
| 20 1,3,5-Trinitrobenzene      | 1   | 17.332       | 17.332           | 0.000            | 105843   | 0.2500           | 0.2609             | M     |
| 21 2,6-Dinitrotoluene         | 1   | 18.492       | 18.492           | 0.000            | 78084    | 0.2500           | 0.2677             |       |
| 22 2,4-Dinitrotoluene         | 1   | 18.952       | 18.952           | 0.000            | 147573   | 0.2500           | 0.2525             |       |
| 23 Tetryl                     | 1   | 22.239       | 22.239           | 0.000            | 71265    | 0.2500           | 0.2503             |       |
| 24 2,4,6-Trinitrotoluene      | 1   | 23.039       | 23.039           | 0.000            | 109657   | 0.2500           | 0.2499             |       |
| 25 PETN                       | 2   | 24.279       | 24.279           | 0.000            | 354659   | 2.50             | 2.46               |       |

**QC Flag Legend**

Processing Flags

Review Flags

M - Manually Integrated

**Reagents:**

8330 DMT\_00013

Amount Added: 12.50

Units: uL

8330IntermStk\_00076

Amount Added: 25.00

Units: uL



Report Date: 25-May-2023 13:06:09

Chrom Revision: 2.3 23-May-2023 13:55:56

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230524-121801.b\05240008.D

Injection Date: 24-May-2023 19:13:25

Instrument ID: CHHPLC\_X5

Operator ID: JZ/JG

Lims ID: CCV DMT

Worklist Smp#: 8

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 7

Method: 8330\_X5\_Luna

Limit Group: GCSV - 8330


Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230524-121801.b\05240008.D  
Injection Date: 24-May-2023 19:13:25 Instrument ID: CHHPLC\_X5  
Lims ID: CCV DMT  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 7 Worklist Smp#: 8  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X5\_Luna Limit Group: GCSV - 8330  
Column: Luna-Phenyl hexyl ( 4.60 mm) Detector LC mwd1A, 254 nm

**19 2-Amino-4,6-dinitrotoluene, CAS: 35572-78-2**

Signal: 1


RT: 17.24  
Area: 213371  
Amount: 0.498601  
Amount Units: ug/ml

## Processing Integration Results



RT: 17.24  
Area: 107547  
Amount: 0.251314  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LV5D, 24-May-2023 19:54:12 -06:00:00 (UTC)

Audit Action: Split an Integrated Peak

Audit Reason: Baseline Smoothing

FORM VII  
HPLC/IC CONTINUING CALIBRATION DATA

Lab Name: Eurofins Denver Job No.: 280-176674-1  
 SDG No.: \_\_\_\_\_  
 Lab Sample ID: CCV 280-613683/18 Calibration Date: 05/25/2023 00:28  
 Instrument ID: CHHPLC\_X5 Calib Start Date: 04/28/2023 18:24  
 GC Column: Luna-phenylhex ID: 4.60 (mm) Calib End Date: 04/28/2023 23:03  
 Lab File ID: 05240018.D Conc. Units: ug/L

| ANALYTE                    | CURVE<br>TYPE | AVE CF | CF     | MIN CF | CALC<br>AMOUNT | SPIKE<br>AMOUNT | %D   | MAX<br>%D |
|----------------------------|---------------|--------|--------|--------|----------------|-----------------|------|-----------|
| TNX                        | Ave           | 394332 | 372790 |        | 237            | 251             | -5.5 | 20.0      |
| DNX                        | Ave           | 277072 | 281950 |        | 255            | 250             | 1.8  | 20.0      |
| HMX                        | Lin2          |        | 201336 |        | 267            | 250             | 6.8  | 20.0      |
| MNX                        | Ave           | 262267 | 254687 |        | 283            | 292             | -2.9 | 20.0      |
| Picric acid                | Lin2          |        | 168120 |        | 260            | 250             | 3.8  | 20.0      |
| RDX                        | Lin2          |        | 231928 |        | 266            | 250             | 6.3  | 20.0      |
| Nitrobenzene               | Ave           | 398834 | 396548 |        | 249            | 250             | -0.6 | 20.0      |
| 3,5-Dinitroaniline         | Lin2          |        | 468764 |        | 262            | 250             | 4.9  | 20.0      |
| 1,3-Dinitrobenzene         | Ave           | 629351 | 626840 |        | 249            | 250             | -0.4 | 20.0      |
| Nitroglycerin              | Ave           | 134735 | 135780 |        | 2520           | 2500            | 0.8  | 20.0      |
| 2-Nitrotoluene             | Ave           | 246963 | 248760 |        | 252            | 250             | 0.7  | 20.0      |
| 4-Nitrotoluene             | Ave           | 227490 | 227292 |        | 250            | 250             | -0.0 | 20.0      |
| 4-Amino-2,6-dinitrotoluene | Ave           | 294210 | 295572 |        | 251            | 250             | 0.5  | 20.0      |
| 3-Nitrotoluene             | Ave           | 294634 | 291356 |        | 247            | 250             | -1.1 | 20.0      |
| 2-Amino-4,6-dinitrotoluene | Ave           | 427939 | 415388 |        | 243            | 250             | -2.9 | 20.0      |
| 1,3,5-Trinitrobenzene      | Ave           | 405692 | 433976 |        | 267            | 250             | 7.0  | 20.0      |
| 2,6-Dinitrotoluene         | Ave           | 291675 | 294584 |        | 252            | 250             | 1.0  | 20.0      |
| 2,4-Dinitrotoluene         | Ave           | 584463 | 588900 |        | 252            | 250             | 0.8  | 20.0      |
| Tetryl                     | Ave           | 284731 | 284520 |        | 250            | 250             | -0.0 | 20.0      |
| 2,4,6-Trinitrotoluene      | Ave           | 438876 | 439980 |        | 251            | 250             | 0.3  | 20.0      |
| PETN                       | Ave           | 143938 | 141745 |        | 2460           | 2500            | -1.5 | 20.0      |
| 1,2-Dinitrobenzene         | Ave           | 280697 | 279344 |        | 249            | 250             | -0.5 | 20.0      |

FORM VII  
HPLC/IC CONTINUING CALIBRATION RETENTION TIME SUMMARY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
 SDG No.: \_\_\_\_\_  
 Lab Sample ID: CCV 280-613683/18 Calibration Date: 05/25/2023 00:28  
 Instrument ID: CHHPLC\_X5 Calib Start Date: 04/28/2023 18:24  
 GC Column: Luna-phenylhex ID: 4.60 (mm) Calib End Date: 04/28/2023 23:03  
 Lab File ID: 05240018.D

| Analyte                    | RT    | RT WINDOW |       |
|----------------------------|-------|-----------|-------|
|                            |       | FROM      | TO    |
| TNX                        | 5.21  | 5.06      | 5.36  |
| DNX                        | 6.05  | 5.90      | 6.20  |
| HMX                        | 6.68  | 6.54      | 6.84  |
| MNX                        | 7.54  | 7.40      | 7.70  |
| Picric acid                | 8.48  | 8.39      | 8.69  |
| RDX                        | 8.92  | 8.80      | 9.10  |
| Nitrobenzene               | 11.41 | 11.31     | 11.61 |
| 3,5-Dinitroaniline         | 14.21 | 14.14     | 14.44 |
| 1,3-Dinitrobenzene         | 14.48 | 14.40     | 14.70 |
| Nitroglycerin              | 14.97 | 14.92     | 15.22 |
| 2-Nitrotoluene             | 15.52 | 15.47     | 15.77 |
| 4-Nitrotoluene             | 15.75 | 15.70     | 16.00 |
| 4-Amino-2,6-dinitrotoluene | 16.30 | 16.26     | 16.56 |
| 3-Nitrotoluene             | 16.59 | 16.54     | 16.84 |
| 2-Amino-4,6-dinitrotoluene | 17.12 | 17.09     | 17.39 |
| 1,3,5-Trinitrobenzene      | 17.26 | 17.18     | 17.48 |
| 2,6-Dinitrotoluene         | 18.38 | 18.34     | 18.64 |
| 2,4-Dinitrotoluene         | 18.84 | 18.80     | 19.10 |
| Tetryl                     | 22.08 | 22.09     | 22.39 |
| 2,4,6-Trinitrotoluene      | 22.90 | 22.89     | 23.19 |
| PETN                       | 24.11 | 24.13     | 24.43 |
| 1,2-Dinitrobenzene         | 12.36 | 12.28     | 12.58 |

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230524-121801.b\05240018.D  
 Lims ID: CCV DMT  
 Client ID:  
 Sample Type: CCV  
 Inject. Date: 25-May-2023 00:28:05 ALS Bottle#: 7 Worklist Smp#: 18  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: CCV DMT  
 Operator ID: JZ/JG Instrument ID: CHHPLC\_X5  
 Sublist: chrom-8330\_X5\_Luna\*sub6  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230524-121801.b\8330\_X5\_Luna.m  
 Limit Group: GCSV - 8330  
 Last Update: 25-May-2023 13:06:15 Calib Date: 29-Apr-2023 04:18:01  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230428-120896.b\04280027.D  
 Column 1 : Luna-Phenyl hexyl ( 4.60 mm) Det: LC mwd1A, 254 nm  
 Process Host: CTX1654

First Level Reviewer: K8YG

Date: 25-May-2023 09:19:36

| Compound                      | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 3 TNX                         | 1   | 5.209        | 5.206            | 0.003            | 93477    | 0.2508           | 0.2371             |       |
| 4 DNx                         | 1   | 6.049        | 6.046            | 0.003            | 70558    | 0.2503           | 0.2547             |       |
| 5 HMX                         | 1   | 6.683        | 6.686            | -0.003           | 50334    | 0.2500           | 0.2670             |       |
| 6 MNX                         | 1   | 7.536        | 7.546            | -0.010           | 74305    | 0.2918           | 0.2833             |       |
| 7 2,4,6-Trinitrophenol        | 1   | 8.483        | 8.539            | -0.056           | 42030    | 0.2500           | 0.2596             |       |
| 8 RDX                         | 1   | 8.923        | 8.946            | -0.023           | 57982    | 0.2500           | 0.2658             |       |
| 9 Nitrobenzene                | 1   | 11.409       | 11.459           | -0.050           | 99137    | 0.2500           | 0.2486             |       |
| \$ 10 1,2-Dinitrobenzene      | 1   | 12.363       | 12.426           | -0.063           | 69836    | 0.2500           | 0.2488             |       |
| 11 3,5-Dinitroaniline         | 1   | 14.209       | 14.292           | -0.083           | 117191   | 0.2500           | 0.2622             |       |
| 12 1,3-Dinitrobenzene         | 1   | 14.476       | 14.552           | -0.076           | 156710   | 0.2500           | 0.2490             |       |
| 13 Nitroglycerin              | 2   | 14.969       | 15.066           | -0.097           | 339450   | 2.50             | 2.52               |       |
| 14 o-Nitrotoluene             | 1   | 15.523       | 15.619           | -0.096           | 62190    | 0.2500           | 0.2518             |       |
| 16 p-Nitrotoluene             | 1   | 15.749       | 15.846           | -0.097           | 56823    | 0.2500           | 0.2498             |       |
| 17 4-Amino-2,6-dinitrotoluene | 1   | 16.296       | 16.406           | -0.110           | 73893    | 0.2500           | 0.2512             |       |
| 18 m-Nitrotoluene             | 1   | 16.589       | 16.692           | -0.103           | 72839    | 0.2500           | 0.2472             |       |
| 19 2-Amino-4,6-dinitrotoluene | 1   | 17.116       | 17.239           | -0.123           | 103847   | 0.2500           | 0.2427             | a     |
| 20 1,3,5-Trinitrobenzene      | 1   | 17.256       | 17.332           | -0.076           | 108494   | 0.2500           | 0.2674             | a     |
| 21 2,6-Dinitrotoluene         | 1   | 18.383       | 18.492           | -0.109           | 73646    | 0.2500           | 0.2525             |       |
| 22 2,4-Dinitrotoluene         | 1   | 18.836       | 18.952           | -0.116           | 147225   | 0.2500           | 0.2519             |       |
| 23 Tetryl                     | 1   | 22.083       | 22.239           | -0.156           | 71130    | 0.2500           | 0.2498             | a     |
| 24 2,4,6-Trinitrotoluene      | 1   | 22.896       | 23.039           | -0.143           | 109995   | 0.2500           | 0.2506             |       |
| 25 PETN                       | 2   | 24.109       | 24.279           | -0.170           | 354363   | 2.50             | 2.46               | a     |

**QC Flag Legend**

Processing Flags

Review Flags

a - User Assigned ID

**Reagents:**

8330 DMT\_00013

Amount Added: 12.50

Units: uL

8330IntermStk\_00076

Amount Added: 25.00

Units: uL

Report Date: 25-May-2023 13:06:16

Chrom Revision: 2.3 23-May-2023 13:55:56

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230524-121801.b\05240018.D

Injection Date: 25-May-2023 00:28:05

Instrument ID: CHHPLC\_X5

Operator ID: JZ/JG

Lims ID: CCV DMT

Worklist Smp#: 18

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 7

Method: 8330\_X5\_Luna

Limit Group: GCSV - 8330


Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: Luna-Phenyl hexyl ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230524-121801.b\05240018.D  
Injection Date: 25-May-2023 00:28:05 Instrument ID: CHHPLC\_X5  
Lims ID: CCV DMT  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 7 Worklist Smp#: 18  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X5\_Luna Limit Group: GCSV - 8330  
Column: Luna-Phenyl hexyl ( 4.60 mm) Detector LC mwd1A, 254 nm

**19 2-Amino-4,6-dinitrotoluene, CAS: 35572-78-2**

Signal: 1


RT: 17.26  
Area: 108494  
Amount: 0.253526  
Amount Units: ug/ml

## Processing Integration Results



RT: 17.12  
Area: 103847  
Amount: 0.242667  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: K8YG, 25-May-2023 09:17:32 -06:00:00 (UTC)

Audit Action: Assigned Compound ID

Audit Reason: Baseline



## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230524-121801.b\05240018.D  
Injection Date: 25-May-2023 00:28:05 Instrument ID: CHHPLC\_X5  
Lims ID: CCV DMT  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 7 Worklist Smp#: 18  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X5\_Luna Limit Group: GCSV - 8330  
Column: Luna-Phenyl hexyl ( 4.60 mm) Detector LC mwd1A, 254 nm

20 1,3,5-Trinitrobenzene, CAS: 99-35-4

Signal: 1


Not Detected  
Expected RT: 17.33

## Processing Integration Results



RT: 17.26  
Area: 108494  
Amount: 0.267429  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: K8YG, 25-May-2023 09:17:34 -06:00:00 (UTC)

Audit Action: Assigned Compound ID

Audit Reason: Baseline

## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230524-121801.b\05240018.D  
Injection Date: 25-May-2023 00:28:05 Instrument ID: CHHPLC\_X5  
Lims ID: CCV DMT  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 7 Worklist Smp#: 18  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X5\_Luna Limit Group: GCSV - 8330  
Column: Luna-Phenyl hexyl ( 4.60 mm) Detector LC mwd1A, 254 nm

**23 Tetryl, CAS: 479-45-8**

Signal: 1


Not Detected  
Expected RT: 22.24

## Processing Integration Results



RT: 22.08  
Area: 71130  
Amount: 0.249815  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: K8YG, 25-May-2023 09:17:39 -06:00:00 (UTC)

Audit Action: Assigned Compound ID

Audit Reason: Baseline

## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\CHHPLC\_X5\20230524-121801.b\05240018.D  
Injection Date: 25-May-2023 00:28:05 Instrument ID: CHHPLC\_X5  
Lims ID: CCV DMT  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 7 Worklist Smp#: 18  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X5\_Luna Limit Group: GCSV - 8330  
Column: Luna-Phenyl hexyl ( 4.60 mm) Detector LC mwd1B, 215 nm

**25 PETN, CAS: 78-11-5**

Signal: 1


Not Detected  
Expected RT: 24.28

## Processing Integration Results



RT: 24.11  
Area: 354363  
Amount: 2.461918  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: K8YG, 25-May-2023 09:17:43 -06:00:00 (UTC)

Audit Action: Assigned Compound ID

Audit Reason: Baseline

FORM I  
HPLC/IC ORGANICS ANALYSIS DATA SHEET

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Client Sample ID: \_\_\_\_\_ Lab Sample ID: MB 280-613095/1-A

Matrix: Water Lab File ID: 05190034.D

Analysis Method: 8330B Date Collected: \_\_\_\_\_

Extraction Method: 3535 Date Extracted: 05/19/2023 13:39

Sample wt/vol: 500 (mL) Date Analyzed: 05/19/2023 23:27

Con. Extract Vol.: 5 (mL) Dilution Factor: 1

Injection Volume: 100 (uL) GC Column: UltraCarb5uODS ID: 4.6 (mm)

% Moisture: \_\_\_\_\_ % Solids: \_\_\_\_\_ GPC Cleanup: (Y/N) N

Cleanup Factor: \_\_\_\_\_

Analysis Batch No.: 613168 Units: ug/L

| CAS NO.    | COMPOUND NAME              | RESULT | Q   | LOQ  | LOD   | DL    |
|------------|----------------------------|--------|-----|------|-------|-------|
| 99-35-4    | 1,3,5-Trinitrobenzene      | 0.20   | U M | 0.21 | 0.20  | 0.084 |
| 99-65-0    | 1,3-Dinitrobenzene         | 0.10   | U   | 0.11 | 0.10  | 0.037 |
| 118-96-7   | 2,4,6-Trinitrotoluene      | 0.10   | U   | 0.11 | 0.10  | 0.045 |
| 6629-29-4  | 2,4-diamino-6-nitrotoluene | 0.90   | U   | 1.0  | 0.90  | 0.43  |
| 121-14-2   | 2,4-Dinitrotoluene         | 0.080  | U   | 0.10 | 0.080 | 0.027 |
| 59229-75-3 | 2,6-diamino-4-nitrotoluene | 0.90   | U M | 1.0  | 0.90  | 0.22  |
| 606-20-2   | 2,6-Dinitrotoluene         | 0.080  | U   | 0.10 | 0.080 | 0.040 |
| 35572-78-2 | 2-Amino-4,6-dinitrotoluene | 0.10   | U   | 0.11 | 0.10  | 0.051 |
| 88-72-2    | 2-Nitrotoluene             | 0.20   | U   | 0.21 | 0.20  | 0.086 |
| 618-87-1   | 3,5-Dinitroaniline         | 0.30   | U   | 0.40 | 0.30  | 0.13  |
| 99-08-1    | 3-Nitrotoluene             | 0.35   | U   | 0.40 | 0.35  | 0.20  |
| 19406-51-0 | 4-Amino-2,6-dinitrotoluene | 0.12   | U   | 0.15 | 0.12  | 0.058 |
| 99-99-0    | 4-Nitrotoluene             | 0.40   | U   | 0.41 | 0.40  | 0.10  |
| 80251-29-2 | DNX                        | 0.25   | U   | 0.50 | 0.25  | 0.097 |
| 2691-41-0  | HMX                        | 0.20   | U   | 0.21 | 0.20  | 0.088 |
| 5755-27-1  | MNX                        | 0.29   | U   | 0.50 | 0.29  | 0.093 |
| 98-95-3    | Nitrobenzene               | 0.20   | U   | 0.21 | 0.20  | 0.091 |
| 55-63-0    | Nitroglycerin              | 2.0    | U   | 2.1  | 2.0   | 0.92  |
| 78-11-5    | PETN                       | 1.0    | U   | 1.1  | 1.0   | 0.45  |
| 121-82-4   | RDX                        | 0.20   | U   | 0.21 | 0.20  | 0.052 |
| 479-45-8   | Tetryl                     | 0.10   | U   | 0.11 | 0.10  | 0.032 |
| 13980-04-6 | TNX                        | 0.25   | U   | 0.50 | 0.25  | 0.079 |

| CAS NO.  | SURROGATE          | %REC | Q | LIMITS |
|----------|--------------------|------|---|--------|
| 528-29-0 | 1,2-Dinitrobenzene | 102  |   | 83-119 |

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\05190034.D  
 Lims ID: MB 280-613095/1-A  
 Client ID:  
 Sample Type: MB  
 Inject. Date: 19-May-2023 23:27:21 ALS Bottle#: 34 Worklist Smp#: 34  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: MB 280-613095/1-A  
 Operator ID: JZ/JG Instrument ID: CHHPLC\_X3  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 20-May-2023 11:22:31 Calib Date: 24-Feb-2023 17:41:20  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1663

First Level Reviewer: LV5D

Date: 20-May-2023 10:26:22

| Compound                      | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/mL | OnCol Amt<br>ug/mL | Flags |
|-------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 1 Triamine Trinitrobenzene    | 1   |              | 2.444            |                  |          |                  | ND                 |       |
| 3 TNX                         | 1   |              | 6.430            |                  |          |                  | ND                 |       |
| 2 2,6-diamino-4-nitrotoluene  | 1   |              | 6.431            |                  |          |                  | ND                 | U     |
| 4 HMX                         | 1   |              | 6.550            |                  |          |                  | ND                 |       |
| 5 2,4-diamino-6-nitrotoluene  | 1   |              | 6.618            |                  |          |                  | ND                 |       |
| 6 DNX                         | 1   |              | 6.757            |                  |          |                  | ND                 |       |
| 7 MNX                         | 1   |              | 7.190            |                  |          |                  | ND                 |       |
| 8 RDX                         | 1   |              | 7.570            |                  |          |                  | ND                 |       |
| 9 2,4,6-Trinitrophenol        | 1   |              | 7.963            |                  |          |                  | ND                 |       |
| \$ 10 1,2-Dinitrobenzene      | 1   | 8.526        | 8.523            | 0.003            | 25862    | 0.2000           | 0.2048             | U     |
| 11 1,3,5-Trinitrobenzene      | 1   |              | 8.650            |                  |          |                  | ND                 |       |
| 12 1,3-Dinitrobenzene         | 1   |              | 9.277            |                  |          |                  | ND                 |       |
| 13 Nitrobenzene               | 1   |              | 9.643            |                  |          |                  | ND                 |       |
| 14 3,5-Dinitroaniline         | 1   |              | 9.883            |                  |          |                  | ND                 |       |
| 15 Tetryl                     | 1   |              | 10.017           |                  |          |                  | ND                 |       |
| 16 Nitroglycerin              | 2   |              | 10.470           |                  |          |                  | ND                 |       |
| 17 2,4,6-Trinitrotoluene      | 1   |              | 10.910           |                  |          |                  | ND                 |       |
| 18 4-Amino-2,6-dinitrotoluene | 1   |              | 11.117           |                  |          |                  | ND                 |       |
| 19 2-Amino-4,6-dinitrotoluene | 1   |              | 11.370           |                  |          |                  | ND                 |       |
| 20 2,6-Dinitrotoluene         | 1   |              | 11.517           |                  |          |                  | ND                 |       |
| 21 2,4-Dinitrotoluene         | 1   |              | 11.683           |                  |          |                  | ND                 |       |
| 22 o-Nitrotoluene             | 1   |              | 12.517           |                  |          |                  | ND                 |       |
| 23 p-Nitrotoluene             | 1   |              | 12.937           |                  |          |                  | ND                 |       |
| 24 m-Nitrotoluene             | 1   |              | 13.517           |                  |          |                  | ND                 |       |
| 25 PETN                       | 2   |              | 14.697           |                  |          |                  | ND                 |       |
| 26 Ammonium Picrate           | 1   |              | 0.000            |                  |          |                  | ND                 |       |

## QC Flag Legend

Processing Flags

Review Flags

U - Marked Undetected

Report Date: 20-May-2023 11:22:32

Chrom Revision: 2.3 16-May-2023 16:37:30

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190034.d

Injection Date: 19-May-2023 23:27:21

Instrument ID: CHHPLC\_X3

Operator ID: JZ/JG

Lims ID: MB 280-613095/1-A

Worklist Smp#: 34

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 34

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



Eurofins Denver  
Recovery Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\05190034.D  
Lims ID: MB 280-613095/1-A  
Client ID:  
Sample Type: MB  
Inject. Date: 19-May-2023 23:27:21 ALS Bottle#: 34 Worklist Smp#: 34  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Sample Info: MB 280-613095/1-A  
Operator ID: JZ/JG Instrument ID: CHHPLC\_X3  
Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\8330\_X3.m  
Limit Group: GCSV - 8330  
Last Update: 20-May-2023 11:22:31 Calib Date: 24-Feb-2023 17:41:20  
Integrator: Falcon  
Quant Method: External Standard Quant By: Initial Calibration  
Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
Process Host: CTX1663

First Level Reviewer: LV5D

Date: 20-May-2023 10:26:22

| Compound                 | Amount Added | Amount Recovered | % Rec. |
|--------------------------|--------------|------------------|--------|
| \$ 10 1,2-Dinitrobenzene | 0.2000       | 0.2048           | 102.38 |



## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190034.d

Injection Date: 19-May-2023 23:27:21

Instrument ID: CHHPLC\_X3

Lims ID: MB 280-613095/1-A

Client ID:

Operator ID: JZ/JG

ALS Bottle#:

34

Worklist Smp#:

34

Injection Vol: 100.0 ul

Dil. Factor:

1.0000

Method: 8330\_X3

Limit Group:

GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm)

Detector


LC DAD1B, 254 nm

**11 1,3,5-Trinitrobenzene, CAS: 99-35-4, Signal: 1**

RT: 8.80

Response: 1046

Amount: 0.004817



Reviewer: LV5D, 20-May-2023 10:26:22

Audit Action: Marked Compound Undetected

Audit Reason: Invalid Compound ID

## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190034.d

Injection Date: 19-May-2023 23:27:21

Instrument ID: CHHPLC\_X3

Lims ID: MB 280-613095/1-A

Client ID:

Operator ID: JZ/JG

ALS Bottle#:

34

Worklist Smp#:

34

Injection Vol: 100.0 ul

Dil. Factor:

1.0000

Method: 8330\_X3

Limit Group:

GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm)

Detector


LC DAD1B, 254 nm

**2,2,6-diamino-4-nitrotoluene, CAS: 59229-75-3, Signal: 1**

RT: 6.44

Response: 9407

Amount: 0.037193



Reviewer: LV5D, 20-May-2023 10:26:22

Audit Action: Marked Compound Undetected

Audit Reason: Invalid Compound ID

FORM I  
HPLC/IC ORGANICS ANALYSIS DATA SHEET

|                                   |                                                      |
|-----------------------------------|------------------------------------------------------|
| Lab Name: <u>Eurofins Denver</u>  | Job No.: <u>280-176674-1</u>                         |
| SDG No.: _____                    |                                                      |
| Client Sample ID: _____           | Lab Sample ID: <u>LCS 280-613095/2-A</u>             |
| Matrix: <u>Water</u>              | Lab File ID: <u>05190035.D</u>                       |
| Analysis Method: <u>8330B</u>     | Date Collected: _____                                |
| Extraction Method: <u>3535</u>    | Date Extracted: <u>05/19/2023 13:39</u>              |
| Sample wt/vol: <u>500 (mL)</u>    | Date Analyzed: <u>05/19/2023 23:50</u>               |
| Con. Extract Vol.: <u>5 (mL)</u>  | Dilution Factor: <u>1</u>                            |
| Injection Volume: <u>100 (uL)</u> | GC Column: <u>UltraCarb5uODS</u> ID: <u>4.6 (mm)</u> |
| % Moisture: _____ % Solids: _____ | GPC Cleanup: (Y/N) <u>N</u>                          |
| Cleanup Factor: _____             |                                                      |
| Analysis Batch No.: <u>613168</u> | Units: <u>ug/L</u>                                   |

| CAS NO.    | COMPOUND NAME              | RESULT | Q | LOQ  | LOD   | DL    |
|------------|----------------------------|--------|---|------|-------|-------|
| 99-35-4    | 1,3,5-Trinitrobenzene      | 1.98   |   | 0.21 | 0.20  | 0.084 |
| 99-65-0    | 1,3-Dinitrobenzene         | 1.90   |   | 0.11 | 0.10  | 0.037 |
| 118-96-7   | 2,4,6-Trinitrotoluene      | 1.83   |   | 0.11 | 0.10  | 0.045 |
| 121-14-2   | 2,4-Dinitrotoluene         | 1.82   |   | 0.10 | 0.080 | 0.027 |
| 606-20-2   | 2,6-Dinitrotoluene         | 1.84   |   | 0.10 | 0.080 | 0.040 |
| 35572-78-2 | 2-Amino-4,6-dinitrotoluene | 1.80   |   | 0.11 | 0.10  | 0.051 |
| 88-72-2    | 2-Nitrotoluene             | 1.56   |   | 0.21 | 0.20  | 0.086 |
| 618-87-1   | 3,5-Dinitroaniline         | 1.73   |   | 0.40 | 0.30  | 0.13  |
| 99-08-1    | 3-Nitrotoluene             | 1.54   |   | 0.40 | 0.35  | 0.20  |
| 19406-51-0 | 4-Amino-2,6-dinitrotoluene | 1.78   |   | 0.15 | 0.12  | 0.058 |
| 99-99-0    | 4-Nitrotoluene             | 1.56   |   | 0.41 | 0.40  | 0.10  |
| 80251-29-2 | DNX                        | 1.93   | M | 0.50 | 0.25  | 0.097 |
| 2691-41-0  | HMX                        | 1.62   | M | 0.21 | 0.20  | 0.088 |
| 5755-27-1  | MNX                        | 2.45   |   | 0.50 | 0.29  | 0.093 |
| 98-95-3    | Nitrobenzene               | 1.74   |   | 0.21 | 0.20  | 0.091 |
| 55-63-0    | Nitroglycerin              | 18.9   |   | 2.1  | 2.0   | 0.92  |
| 78-11-5    | PETN                       | 20.5   |   | 1.1  | 1.0   | 0.45  |
| 121-82-4   | RDX                        | 1.79   |   | 0.21 | 0.20  | 0.052 |
| 479-45-8   | Tetryl                     | 2.07   |   | 0.11 | 0.10  | 0.032 |
| 13980-04-6 | TNX                        | 1.88   | M | 0.50 | 0.25  | 0.079 |

| CAS NO.  | SURROGATE          | %REC | Q | LIMITS |
|----------|--------------------|------|---|--------|
| 528-29-0 | 1,2-Dinitrobenzene | 93   |   | 83-119 |

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\05190035.D  
 Lims ID: LCS 280-613095/2-A  
 Client ID:  
 Sample Type: LCS  
 Inject. Date: 19-May-2023 23:50:14 ALS Bottle#: 35 Worklist Smp#: 35  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: LCS 280-613095/2-A  
 Operator ID: JZ/JG Instrument ID: CHHPLC\_X3  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 20-May-2023 11:22:31 Calib Date: 24-Feb-2023 17:41:20  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1663

First Level Reviewer: LV5D

Date: 20-May-2023 10:26:38

| Compound                      | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/mL | OnCol Amt<br>ug/mL | Flags |
|-------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 3 TNX                         | 1   | 6.429        | 6.430            | -0.001           | 37089    | 0.2006           | 0.1876             | M     |
| 4 HMX                         | 1   | 6.549        | 6.550            | -0.001           | 15150    | 0.2000           | 0.1620             | M     |
| 6 DNX                         | 1   | 6.755        | 6.757            | -0.002           | 27816    | 0.2002           | 0.1925             | M     |
| 7 MNX                         | 1   | 7.182        | 7.190            | -0.008           | 32117    | 0.2334           | 0.2448             |       |
| 8 RDX                         | 1   | 7.569        | 7.570            | -0.001           | 19040    | 0.2000           | 0.1790             |       |
| 9 2,4,6-Trinitrophenol        | 1   | 7.962        | 7.963            | -0.001           | 16851    | 0.2000           | 0.2222             |       |
| \$ 10 1,2-Dinitrobenzene      | 1   | 8.522        | 8.523            | -0.001           | 23565    | 0.2000           | 0.1866             |       |
| 11 1,3,5-Trinitrobenzene      | 1   | 8.649        | 8.650            | -0.001           | 42888    | 0.2000           | 0.1975             |       |
| 12 1,3-Dinitrobenzene         | 1   | 9.275        | 9.277            | -0.002           | 55963    | 0.2000           | 0.1901             |       |
| 13 Nitrobenzene               | 1   | 9.642        | 9.643            | -0.001           | 33367    | 0.2000           | 0.1745             |       |
| 14 3,5-Dinitroaniline         | 1   | 9.882        | 9.883            | -0.001           | 39507    | 0.2000           | 0.1730             |       |
| 15 Tetryl                     | 1   | 10.015       | 10.017           | -0.002           | 33961    | 0.2000           | 0.2069             |       |
| 16 Nitroglycerin              | 2   | 10.462       | 10.470           | -0.008           | 121383   | 2.00             | 1.89               |       |
| 17 2,4,6-Trinitrotoluene      | 1   | 10.902       | 10.910           | -0.008           | 38574    | 0.2000           | 0.1828             |       |
| 18 4-Amino-2,6-dinitrotoluene | 1   | 11.102       | 11.117           | -0.015           | 27601    | 0.2000           | 0.1781             |       |
| 19 2-Amino-4,6-dinitrotoluene | 1   | 11.355       | 11.370           | -0.015           | 36267    | 0.2000           | 0.1801             |       |
| 20 2,6-Dinitrotoluene         | 1   | 11.502       | 11.517           | -0.015           | 26241    | 0.2000           | 0.1838             |       |
| 21 2,4-Dinitrotoluene         | 1   | 11.668       | 11.683           | -0.015           | 53943    | 0.2000           | 0.1818             |       |
| 22 o-Nitrotoluene             | 1   | 12.488       | 12.517           | -0.029           | 19964    | 0.2000           | 0.1561             |       |
| 23 p-Nitrotoluene             | 1   | 12.908       | 12.937           | -0.029           | 17433    | 0.2000           | 0.1558             |       |
| 24 m-Nitrotoluene             | 1   | 13.482       | 13.517           | -0.035           | 21634    | 0.2000           | 0.1540             |       |
| 25 PETN                       | 2   | 14.655       | 14.697           | -0.042           | 140997   | 2.00             | 2.05               |       |

**QC Flag Legend**

Processing Flags

Review Flags

M - Manually Integrated

Report Date: 20-May-2023 11:22:32

Chrom Revision: 2.3 16-May-2023 16:37:30

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190035.d

Injection Date: 19-May-2023 23:50:14

Instrument ID: CHHPLC\_X3

Operator ID: JZ/JG

Lims ID: LCS 280-613095/2-A

Worklist Smp#: 35

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 35

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



Eurofins Denver  
Recovery Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\05190035.D  
Lims ID: LCS 280-613095/2-A  
Client ID:  
Sample Type: LCS  
Inject. Date: 19-May-2023 23:50:14 ALS Bottle#: 35 Worklist Smp#: 35  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Sample Info: LCS 280-613095/2-A  
Operator ID: JZ/JG Instrument ID: CHHPLC\_X3  
Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\8330\_X3.m  
Limit Group: GCSV - 8330  
Last Update: 20-May-2023 11:22:31 Calib Date: 24-Feb-2023 17:41:20  
Integrator: Falcon  
Quant Method: External Standard Quant By: Initial Calibration  
Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
Process Host: CTX1663

First Level Reviewer: LV5D

Date: 20-May-2023 10:26:38

| Compound                 | Amount Added | Amount Recovered | % Rec. |
|--------------------------|--------------|------------------|--------|
| \$ 10 1,2-Dinitrobenzene | 0.2000       | 0.1866           | 93.28  |

## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190035.d  
Injection Date: 19-May-2023 23:50:14 Instrument ID: CHHPLC\_X3  
Lims ID: LCS 280-613095/2-A  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 35 Worklist Smp#: 35  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**6 DNX, CAS: 80251-29-2**

Signal: 1


RT: 6.76  
Area: 39586  
Amount: 0.274023  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.76  
Area: 27816  
Amount: 0.192548  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 10:26:37 -06:00:00 (UTC)

Audit Action: Manually Integrated

Audit Reason: Baseline

## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190035.d

Injection Date: 19-May-2023 23:50:14

Instrument ID: CHHPLC\_X3

Lims ID: LCS 280-613095/2-A

Client ID:

Operator ID: JZ/JG

ALS Bottle#:

35

Worklist Smp#: 35

Injection Vol: 100.0 ul

Dil. Factor: 1.0000

Method: 8330\_X3

Limit Group: GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm)


Detector: LC DAD1B, 254 nm

**4 HMX, CAS: 2691-41-0**

Signal: 1


RT: 6.55  
Area: 25726  
Amount: 0.275062  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.55  
Area: 15150  
Amount: 0.161983  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 10:26:34 -06:00:00 (UTC)

Audit Action: Split an Integrated Peak

Audit Reason: Baseline



## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190035.d

Injection Date: 19-May-2023 23:50:14

Instrument ID: CHHPLC\_X3

Lims ID: LCS 280-613095/2-A

Client ID:

Operator ID: JZ/JG

ALS Bottle#:

35

Worklist Smp#: 35

Injection Vol: 100.0 ul

Dil. Factor: 1.0000

Method: 8330\_X3

Limit Group: GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm)

Detector


LC DAD1B, 254 nm

**3 TNX, CAS: 13980-04-6**

Signal: 1


RT: 6.43  
Area: 54654  
Amount: 0.276467  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.43  
Area: 37089  
Amount: 0.187614  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 10:26:32 -06:00:00 (UTC)

Audit Action: Split an Integrated Peak

Audit Reason: Baseline

FORM I  
HPLC/IC ORGANICS ANALYSIS DATA SHEET

Lab Name: Eurofins Denver Job No.: 280-176674-1  
 SDG No.: \_\_\_\_\_  
 Client Sample ID: \_\_\_\_\_ Lab Sample ID: LCS 280-613095/4-A  
 Matrix: Water Lab File ID: 05190037.D  
 Analysis Method: 8330B Date Collected: \_\_\_\_\_  
 Extraction Method: 3535 Date Extracted: 05/19/2023 13:39  
 Sample wt/vol: 500 (mL) Date Analyzed: 05/20/2023 00:36  
 Con. Extract Vol.: 5 (mL) Dilution Factor: 1  
 Injection Volume: 100 (uL) GC Column: UltraCarb5uODS ID: 4.6 (mm)  
 % Moisture: \_\_\_\_\_ % Solids: \_\_\_\_\_ GPC Cleanup: (Y/N) N  
 Cleanup Factor: \_\_\_\_\_  
 Analysis Batch No.: 613168 Units: ug/L

| CAS NO.    | COMPOUND NAME              | RESULT | Q | LOQ | LOD  | DL   |
|------------|----------------------------|--------|---|-----|------|------|
| 6629-29-4  | 2,4-diamino-6-nitrotoluene | 1.61   | M | 1.0 | 0.90 | 0.43 |
| 59229-75-3 | 2,6-diamino-4-nitrotoluene | 1.68   | M | 1.0 | 0.90 | 0.22 |

| CAS NO.  | SURROGATE          | %REC | Q | LIMITS |
|----------|--------------------|------|---|--------|
| 528-29-0 | 1,2-Dinitrobenzene | 98   |   | 83-119 |

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\05190037.D  
 Lims ID: LCS 280-613095/4-A  
 Client ID:  
 Sample Type: LCS  
 Inject. Date: 20-May-2023 00:36:13 ALS Bottle#: 37 Worklist Smp#: 37  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: LCS 280-613095/4-A  
 Operator ID: JZ/JG Instrument ID: CHHPLC\_X3  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 20-May-2023 11:22:31 Calib Date: 24-Feb-2023 17:41:20  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1663

First Level Reviewer: LV5D

Date: 20-May-2023 10:27:04

| Compound                     | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/mL | OnCol Amt<br>ug/mL | Flags |
|------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 2 2,6-diamino-4-nitrotoluene | 1   | 6.418        | 6.431            | -0.013           | 39464    | 0.2000           | 0.1680             | M     |
| 5 2,4-diamino-6-nitrotoluene | 1   | 6.605        | 6.618            | -0.013           | 23586    | 0.2000           | 0.1610             | M     |
| \$ 10 1,2-Dinitrobenzene     | 1   | 8.531        | 8.523            | 0.008            | 24832    | 0.2000           | 0.1966             |       |

**QC Flag Legend**

Review Flags

M - Manually Integrated

Report Date: 20-May-2023 11:22:33

Chrom Revision: 2.3 16-May-2023 16:37:30

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190037.d

Injection Date: 20-May-2023 00:36:13

Instrument ID: CHHPLC\_X3

Operator ID: JZ/JG

Lims ID: LCS 280-613095/4-A

Worklist Smp#: 37

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 37

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



Eurofins Denver  
Recovery Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\05190037.D  
Lims ID: LCS 280-613095/4-A  
Client ID:  
Sample Type: LCS  
Inject. Date: 20-May-2023 00:36:13 ALS Bottle#: 37 Worklist Smp#: 37  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Sample Info: LCS 280-613095/4-A  
Operator ID: JZ/JG Instrument ID: CHHPLC\_X3  
Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\8330\_X3.m  
Limit Group: GCSV - 8330  
Last Update: 20-May-2023 11:22:31 Calib Date: 24-Feb-2023 17:41:20  
Integrator: Falcon  
Quant Method: External Standard Quant By: Initial Calibration  
Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
Process Host: CTX1663

First Level Reviewer: LV5D

Date: 20-May-2023 10:27:04

| Compound                 | Amount Added | Amount Recovered | % Rec. |
|--------------------------|--------------|------------------|--------|
| \$ 10 1,2-Dinitrobenzene | 0.2000       | 0.1966           | 98.30  |

## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190037.d

Injection Date: 20-May-2023 00:36:13

Instrument ID: CHHPLC\_X3

Lims ID: LCS 280-613095/4-A

Client ID:

Operator ID: JZ/JG

ALS Bottle#:

37

Worklist Smp#: 37

Injection Vol: 100.0 ul

Dil. Factor:

1.0000

Method: 8330\_X3

Limit Group:

GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm)

Detector


LC DAD1B, 254 nm

**5 2,4-diamino-6-nitrotoluene, CAS: 6629-29-4**

Signal: 1


RT: 6.60  
Area: 34389  
Amount: 0.234782  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.60  
Area: 23586  
Amount: 0.161027  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 10:27:04 -06:00:00 (UTC)

Audit Action: Split an Integrated Peak

Audit Reason: Baseline

## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190037.d

Injection Date: 20-May-2023 00:36:13

Instrument ID: CHHPLC\_X3

Lims ID: LCS 280-613095/4-A

Client ID:

Operator ID: JZ/JG

ALS Bottle#:

37

Worklist Smp#: 37

Injection Vol: 100.0 ul

Dil. Factor:

1.0000

Method: 8330\_X3

Limit Group:

GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm)

Detector


LC DAD1B, 254 nm

**2,2,6-diamino-4-nitrotoluene, CAS: 59229-75-3**

Signal: 1


RT: 6.42  
Area: 53667  
Amount: 0.229885  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.42  
Area: 39464  
Amount: 0.168050  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 10:27:02 -06:00:00 (UTC)

Audit Action: Split an Integrated Peak

Audit Reason: Baseline

FORM I  
HPLC/IC ORGANICS ANALYSIS DATA SHEET

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Client Sample ID: \_\_\_\_\_ Lab Sample ID: LCSD 280-613095/3-A

Matrix: Water Lab File ID: 05190036.D

Analysis Method: 8330B Date Collected: \_\_\_\_\_

Extraction Method: 3535 Date Extracted: 05/19/2023 13:39

Sample wt/vol: 500 (mL) Date Analyzed: 05/20/2023 00:13

Con. Extract Vol.: 5 (mL) Dilution Factor: 1

Injection Volume: 100 (uL) GC Column: UltraCarb5uODS ID: 4.6 (mm)

% Moisture: \_\_\_\_\_ % Solids: \_\_\_\_\_ GPC Cleanup: (Y/N) N

Cleanup Factor: \_\_\_\_\_

Analysis Batch No.: 613168 Units: ug/L

| CAS NO.    | COMPOUND NAME              | RESULT | Q | LOQ  | LOD   | DL    |
|------------|----------------------------|--------|---|------|-------|-------|
| 99-35-4    | 1,3,5-Trinitrobenzene      | 2.04   |   | 0.21 | 0.20  | 0.084 |
| 99-65-0    | 1,3-Dinitrobenzene         | 1.96   |   | 0.11 | 0.10  | 0.037 |
| 118-96-7   | 2,4,6-Trinitrotoluene      | 1.88   |   | 0.11 | 0.10  | 0.045 |
| 121-14-2   | 2,4-Dinitrotoluene         | 1.87   |   | 0.10 | 0.080 | 0.027 |
| 606-20-2   | 2,6-Dinitrotoluene         | 1.95   |   | 0.10 | 0.080 | 0.040 |
| 35572-78-2 | 2-Amino-4,6-dinitrotoluene | 1.84   |   | 0.11 | 0.10  | 0.051 |
| 88-72-2    | 2-Nitrotoluene             | 1.64   |   | 0.21 | 0.20  | 0.086 |
| 618-87-1   | 3,5-Dinitroaniline         | 1.83   |   | 0.40 | 0.30  | 0.13  |
| 99-08-1    | 3-Nitrotoluene             | 1.62   |   | 0.40 | 0.35  | 0.20  |
| 19406-51-0 | 4-Amino-2,6-dinitrotoluene | 1.82   |   | 0.15 | 0.12  | 0.058 |
| 99-99-0    | 4-Nitrotoluene             | 1.61   |   | 0.41 | 0.40  | 0.10  |
| 80251-29-2 | DNX                        | 1.95   | M | 0.50 | 0.25  | 0.097 |
| 2691-41-0  | HMX                        | 1.71   | M | 0.21 | 0.20  | 0.088 |
| 5755-27-1  | MNX                        | 2.47   |   | 0.50 | 0.29  | 0.093 |
| 98-95-3    | Nitrobenzene               | 1.81   |   | 0.21 | 0.20  | 0.091 |
| 55-63-0    | Nitroglycerin              | 19.5   |   | 2.1  | 2.0   | 0.92  |
| 78-11-5    | PETN                       | 21.1   |   | 1.1  | 1.0   | 0.45  |
| 121-82-4   | RDX                        | 1.84   |   | 0.21 | 0.20  | 0.052 |
| 479-45-8   | Tetryl                     | 2.06   |   | 0.11 | 0.10  | 0.032 |
| 13980-04-6 | TNX                        | 1.91   | M | 0.50 | 0.25  | 0.079 |

| CAS NO.  | SURROGATE          | %REC | Q | LIMITS |
|----------|--------------------|------|---|--------|
| 528-29-0 | 1,2-Dinitrobenzene | 94   |   | 83-119 |



Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\05190036.D  
 Lims ID: LCSD 280-613095/3-A  
 Client ID:  
 Sample Type: LCSD  
 Inject. Date: 20-May-2023 00:13:13 ALS Bottle#: 36 Worklist Smp#: 36  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: LCSD 280-613095/3-A  
 Operator ID: JZ/JG Instrument ID: CHHPLC\_X3  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 20-May-2023 11:22:31 Calib Date: 24-Feb-2023 17:41:20  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1663

First Level Reviewer: LV5D

Date: 20-May-2023 10:26:51

| Compound                      | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/mL | OnCol Amt<br>ug/mL | Flags |
|-------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 3 TNX                         | 1   | 6.430        | 6.430            | 0.000            | 37692    | 0.2006           | 0.1907             | M     |
| 4 HMX                         | 1   | 6.550        | 6.550            | 0.000            | 16032    | 0.2000           | 0.1714             | M     |
| 6 DNx                         | 1   | 6.757        | 6.757            | 0.000            | 28144    | 0.2002           | 0.1948             | M     |
| 7 MNx                         | 1   | 7.184        | 7.190            | -0.006           | 32355    | 0.2334           | 0.2467             |       |
| 8 RDX                         | 1   | 7.570        | 7.570            | 0.000            | 19542    | 0.2000           | 0.1837             |       |
| 9 2,4,6-Trinitrophenol        | 1   | 7.957        | 7.963            | -0.006           | 17381    | 0.2000           | 0.2292             |       |
| \$ 10 1,2-Dinitrobenzene      | 1   | 8.517        | 8.523            | -0.006           | 23733    | 0.2000           | 0.1879             |       |
| 11 1,3,5-Trinitrobenzene      | 1   | 8.644        | 8.650            | -0.006           | 44269    | 0.2000           | 0.2039             |       |
| 12 1,3-Dinitrobenzene         | 1   | 9.270        | 9.277            | -0.007           | 57656    | 0.2000           | 0.1958             |       |
| 13 Nitrobenzene               | 1   | 9.637        | 9.643            | -0.006           | 34660    | 0.2000           | 0.1812             |       |
| 14 3,5-Dinitroaniline         | 1   | 9.877        | 9.883            | -0.006           | 41720    | 0.2000           | 0.1827             |       |
| 15 Tetryl                     | 1   | 10.010       | 10.017           | -0.007           | 33876    | 0.2000           | 0.2064             |       |
| 16 Nitroglycerin              | 2   | 10.464       | 10.470           | -0.006           | 124924   | 2.00             | 1.95               |       |
| 17 2,4,6-Trinitrotoluene      | 1   | 10.904       | 10.910           | -0.006           | 39578    | 0.2000           | 0.1875             |       |
| 18 4-Amino-2,6-dinitrotoluene | 1   | 11.104       | 11.117           | -0.013           | 28204    | 0.2000           | 0.1820             |       |
| 19 2-Amino-4,6-dinitrotoluene | 1   | 11.357       | 11.370           | -0.013           | 37107    | 0.2000           | 0.1842             |       |
| 20 2,6-Dinitrotoluene         | 1   | 11.504       | 11.517           | -0.013           | 27882    | 0.2000           | 0.1953             |       |
| 21 2,4-Dinitrotoluene         | 1   | 11.670       | 11.683           | -0.013           | 55373    | 0.2000           | 0.1867             |       |
| 22 o-Nitrotoluene             | 1   | 12.497       | 12.517           | -0.020           | 21009    | 0.2000           | 0.1643             |       |
| 23 p-Nitrotoluene             | 1   | 12.924       | 12.937           | -0.013           | 17980    | 0.2000           | 0.1607             |       |
| 24 m-Nitrotoluene             | 1   | 13.497       | 13.517           | -0.020           | 22802    | 0.2000           | 0.1623             |       |
| 25 PETN                       | 2   | 14.670       | 14.697           | -0.027           | 145328   | 2.00             | 2.11               |       |

### QC Flag Legend

Processing Flags

Review Flags

M - Manually Integrated

Report Date: 20-May-2023 11:22:33

Chrom Revision: 2.3 16-May-2023 16:37:30

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190036.d

Injection Date: 20-May-2023 00:13:13

Instrument ID: CHHPLC\_X3

Operator ID: JZ/JG

Lims ID: LCSD 280-613095/3-A

Worklist Smp#: 36

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 36

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



Eurofins Denver  
Recovery Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\05190036.D  
Lims ID: LCSD 280-613095/3-A  
Client ID:  
Sample Type: LCSD  
Inject. Date: 20-May-2023 00:13:13 ALS Bottle#: 36 Worklist Smp#: 36  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Sample Info: LCSD 280-613095/3-A  
Operator ID: JZ/JG Instrument ID: CHHPLC\_X3  
Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\8330\_X3.m  
Limit Group: GCSV - 8330  
Last Update: 20-May-2023 11:22:31 Calib Date: 24-Feb-2023 17:41:20  
Integrator: Falcon  
Quant Method: External Standard Quant By: Initial Calibration  
Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
Process Host: CTX1663

First Level Reviewer: LV5D

Date: 20-May-2023 10:26:51

| Compound                 | Amount Added | Amount Recovered | % Rec. |
|--------------------------|--------------|------------------|--------|
| \$ 10 1,2-Dinitrobenzene | 0.2000       | 0.1879           | 93.95  |

## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190036.d  
Injection Date: 20-May-2023 00:13:13 Instrument ID: CHHPLC\_X3  
Lims ID: LCSD 280-613095/3-A  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 36 Worklist Smp#: 36  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**6 DNX, CAS: 80251-29-2**

Signal: 1


RT: 6.76  
Area: 39615  
Amount: 0.274223  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.76  
Area: 28144  
Amount: 0.194819  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 10:26:49 -06:00:00 (UTC)

Audit Action: Manually Integrated

Audit Reason: Baseline

## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190036.d

Injection Date: 20-May-2023 00:13:13

Instrument ID: CHHPLC\_X3

Lims ID: LCSD 280-613095/3-A

Client ID:

Operator ID: JZ/JG

ALS Bottle#:

36

Worklist Smp#: 36

Injection Vol: 100.0 ul

Dil. Factor: 1.0000

Method: 8330\_X3

Limit Group: GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm)


Detector: LC DAD1B, 254 nm

**4 HMX, CAS: 2691-41-0**

Signal: 1


RT: 6.55  
Area: 27366  
Amount: 0.292596  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.55  
Area: 16032  
Amount: 0.171414  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 10:26:46 -06:00:00 (UTC)

Audit Action: Split an Integrated Peak

Audit Reason: Baseline

## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190036.d

Injection Date: 20-May-2023 00:13:13

Instrument ID: CHHPLC\_X3

Lims ID: LCSD 280-613095/3-A

Client ID:

Operator ID: JZ/JG

ALS Bottle#:

36

Worklist Smp#: 36

Injection Vol: 100.0 ul

Dil. Factor: 1.0000

Method: 8330\_X3

Limit Group: GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm)


Detector: LC DAD1B, 254 nm

**3 TNX, CAS: 13980-04-6**

Signal: 1


RT: 6.43  
Area: 56112  
Amount: 0.283842  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.43  
Area: 37692  
Amount: 0.190665  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 10:26:45 -06:00:00 (UTC)

Audit Action: Split an Integrated Peak

Audit Reason: Baseline

FORM I  
HPLC/IC ORGANICS ANALYSIS DATA SHEET

Lab Name: Eurofins Denver Job No.: 280-176674-1  
SDG No.: \_\_\_\_\_  
Client Sample ID: \_\_\_\_\_ Lab Sample ID: LCSD 280-613095/5-A  
Matrix: Water Lab File ID: 05190038.D  
Analysis Method: 8330B Date Collected: \_\_\_\_\_  
Extraction Method: 3535 Date Extracted: 05/19/2023 13:39  
Sample wt/vol: 500 (mL) Date Analyzed: 05/20/2023 00:59  
Con. Extract Vol.: 5 (mL) Dilution Factor: 1  
Injection Volume: 100 (uL) GC Column: UltraCarb5uODS ID: 4.6 (mm)  
% Moisture: \_\_\_\_\_ % Solids: \_\_\_\_\_ GPC Cleanup: (Y/N) N  
Cleanup Factor: \_\_\_\_\_  
Analysis Batch No.: 613168 Units: ug/L

| CAS NO.    | COMPOUND NAME              | RESULT | Q | LOQ | LOD  | DL   |
|------------|----------------------------|--------|---|-----|------|------|
| 6629-29-4  | 2,4-diamino-6-nitrotoluene | 1.47   | M | 1.0 | 0.90 | 0.43 |
| 59229-75-3 | 2,6-diamino-4-nitrotoluene | 1.56   | M | 1.0 | 0.90 | 0.22 |

| CAS NO.  | SURROGATE          | %REC | Q | LIMITS |
|----------|--------------------|------|---|--------|
| 528-29-0 | 1,2-Dinitrobenzene | 98   |   | 83-119 |

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\05190038.D  
 Lims ID: LCSD 280-613095/5-A  
 Client ID:  
 Sample Type: LCSD  
 Inject. Date: 20-May-2023 00:59:08 ALS Bottle#: 38 Worklist Smp#: 38  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: LCSD 280-613095/5-A  
 Operator ID: JZ/JG Instrument ID: CHHPLC\_X3  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 20-May-2023 11:22:31 Calib Date: 24-Feb-2023 17:41:20  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1663

First Level Reviewer: LV5D

Date: 20-May-2023 10:27:14

| Compound                     | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/mL | OnCol Amt<br>ug/mL | Flags |
|------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 2 2,6-diamino-4-nitrotoluene | 1   | 6.419        | 6.431            | -0.012           | 36645    | 0.2000           | 0.1558             | M     |
| 5 2,4-diamino-6-nitrotoluene | 1   | 6.605        | 6.618            | -0.013           | 21570    | 0.2000           | 0.1473             | M     |
| \$ 10 1,2-Dinitrobenzene     | 1   | 8.525        | 8.523            | 0.002            | 24672    | 0.2000           | 0.1953             |       |

**QC Flag Legend**

Review Flags

M - Manually Integrated



Report Date: 20-May-2023 11:22:34

Chrom Revision: 2.3 16-May-2023 16:37:30

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190038.d

Injection Date: 20-May-2023 00:59:08

Instrument ID: CHHPLC\_X3

Operator ID: JZ/JG

Lims ID: LCSD 280-613095/5-A

Worklist Smp#: 38

Client ID:

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 38

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



Eurofins Denver  
Recovery Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\05190038.D  
Lims ID: LCSD 280-613095/5-A  
Client ID:  
Sample Type: LCSD  
Inject. Date: 20-May-2023 00:59:08 ALS Bottle#: 38 Worklist Smp#: 38  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Sample Info: LCSD 280-613095/5-A  
Operator ID: JZ/JG Instrument ID: CHHPLC\_X3  
Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\8330\_X3.m  
Limit Group: GCSV - 8330  
Last Update: 20-May-2023 11:22:31 Calib Date: 24-Feb-2023 17:41:20  
Integrator: Falcon  
Quant Method: External Standard Quant By: Initial Calibration  
Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
Process Host: CTX1663

First Level Reviewer: LV5D

Date: 20-May-2023 10:27:14

| Compound                 | Amount Added | Amount Recovered | % Rec. |
|--------------------------|--------------|------------------|--------|
| \$ 10 1,2-Dinitrobenzene | 0.2000       | 0.1953           | 97.67  |

## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190038.d  
Injection Date: 20-May-2023 00:59:08 Instrument ID: CHHPLC\_X3  
Lims ID: LCSD 280-613095/5-A  
Client ID:  
Operator ID: JZ/JG ALS Bottle#: 38 Worklist Smp#: 38  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**5 2,4-diamino-6-nitrotoluene, CAS: 6629-29-4**

Signal: 1


RT: 6.61  
Area: 32318  
Amount: 0.220643  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.61  
Area: 21570  
Amount: 0.147264  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 10:27:13 -06:00:00 (UTC)

Audit Action: Split an Integrated Peak

Audit Reason: Baseline

## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190038.d

Injection Date: 20-May-2023 00:59:08

Instrument ID: CHHPLC\_X3

Lims ID: LCSD 280-613095/5-A

Client ID:

Operator ID: JZ/JG

ALS Bottle#:

38

Worklist Smp#: 38

Injection Vol: 100.0 ul

Dil. Factor:

1.0000

Method: 8330\_X3

Limit Group:

GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm)

Detector


LC DAD1B, 254 nm

**2,2,6-diamino-4-nitrotoluene, CAS: 59229-75-3**

Signal: 1


RT: 6.42  
Area: 51579  
Amount: 0.220794  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.42  
Area: 36645  
Amount: 0.155777  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 10:27:12 -06:00:00 (UTC)

Audit Action: Split an Integrated Peak

Audit Reason: Baseline

FORM I  
HPLC/IC ORGANICS ANALYSIS DATA SHEET

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Client Sample ID: FBQmw-175-230401-GW-R MS Lab Sample ID: 280-176674-10 MS

Matrix: Water Lab File ID: 05190049.D

Analysis Method: 8330B Date Collected: 05/17/2023 08:58

Extraction Method: 3535 Date Extracted: 05/19/2023 13:39

Sample wt/vol: 499.4 (mL) Date Analyzed: 05/20/2023 05:11

Con. Extract Vol.: 5 (mL) Dilution Factor: 1

Injection Volume: 100 (uL) GC Column: UltraCarb5uODS ID: 4.6 (mm)

% Moisture: \_\_\_\_\_ % Solids: \_\_\_\_\_ GPC Cleanup: (Y/N) N

Cleanup Factor: \_\_\_\_\_

Analysis Batch No.: 613168 Units: ug/L

| CAS NO.    | COMPOUND NAME              | RESULT | Q    | LOQ  | LOD   | DL    |
|------------|----------------------------|--------|------|------|-------|-------|
| 99-35-4    | 1,3,5-Trinitrobenzene      | 1.75   | Q    | 0.21 | 0.20  | 0.084 |
| 99-65-0    | 1,3-Dinitrobenzene         | 1.66   | Q    | 0.11 | 0.10  | 0.037 |
| 118-96-7   | 2,4,6-Trinitrotoluene      | 1.56   | Q    | 0.11 | 0.10  | 0.045 |
| 121-14-2   | 2,4-Dinitrotoluene         | 1.50   | Q J1 | 0.10 | 0.080 | 0.027 |
| 606-20-2   | 2,6-Dinitrotoluene         | 1.58   | Q    | 0.10 | 0.080 | 0.040 |
| 35572-78-2 | 2-Amino-4,6-dinitrotoluene | 1.46   | Q J1 | 0.11 | 0.10  | 0.051 |
| 88-72-2    | 2-Nitrotoluene             | 1.34   | Q J1 | 0.21 | 0.20  | 0.086 |
| 618-87-1   | 3,5-Dinitroaniline         | 1.46   | Q    | 0.40 | 0.30  | 0.13  |
| 99-08-1    | 3-Nitrotoluene             | 1.23   | Q J1 | 0.40 | 0.35  | 0.20  |
| 19406-51-0 | 4-Amino-2,6-dinitrotoluene | 1.44   | Q J1 | 0.15 | 0.12  | 0.058 |
| 99-99-0    | 4-Nitrotoluene             | 1.30   | Q J1 | 0.41 | 0.40  | 0.10  |
| 80251-29-2 | DNX                        | 1.56   | M Q  | 0.50 | 0.25  | 0.097 |
| 2691-41-0  | HMX                        | 1.46   | M Q  | 0.21 | 0.20  | 0.088 |
| 5755-27-1  | MNX                        | 1.98   | M Q  | 0.50 | 0.29  | 0.093 |
| 98-95-3    | Nitrobenzene               | 1.47   | Q    | 0.21 | 0.20  | 0.091 |
| 55-63-0    | Nitroglycerin              | 17.5   | Q    | 2.1  | 2.0   | 0.92  |
| 78-11-5    | PETN                       | 18.1   | Q    | 1.1  | 1.0   | 0.45  |
| 121-82-4   | RDX                        | 1.49   | M Q  | 0.21 | 0.20  | 0.052 |
| 479-45-8   | Tetryl                     | 1.72   | Q    | 0.11 | 0.10  | 0.032 |
| 13980-04-6 | TNX                        | 1.70   | M Q  | 0.50 | 0.25  | 0.079 |

| CAS NO.  | SURROGATE          | %REC | Q | LIMITS |
|----------|--------------------|------|---|--------|
| 528-29-0 | 1,2-Dinitrobenzene | 82   | Q | 83-119 |

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\05190049.D  
 Lims ID: 280-176674-B-10-A MS  
 Client ID: FBQmw-175-230401-GW-R  
 Sample Type: MS  
 Inject. Date: 20-May-2023 05:11:34 ALS Bottle#: 49 Worklist Smp#: 49  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-176674-B-10-A MS  
 Operator ID: JZ/JG Instrument ID: CHHPLC\_X3  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 20-May-2023 11:22:13 Calib Date: 24-Feb-2023 17:41:20  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1663

First Level Reviewer: LV5D

Date: 20-May-2023 11:15:10

| Compound                      | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/mL | OnCol Amt<br>ug/mL | Flags |
|-------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 3 TNX                         | 1   | 6.428        | 6.430            | -0.002           | 33570    | 0.2006           | 0.1698             | M     |
| 4 HMX                         | 1   | 6.548        | 6.550            | -0.002           | 13659    | 0.2000           | 0.1460             | M     |
| 6 DNx                         | 1   | 6.755        | 6.757            | -0.002           | 22572    | 0.2002           | 0.1562             | M     |
| 7 MNx                         | 1   | 7.181        | 7.190            | -0.009           | 25990    | 0.2334           | 0.1981             | M     |
| 8 RDX                         | 1   | 7.568        | 7.570            | -0.002           | 15813    | 0.2000           | 0.1486             | M     |
| 9 2,4,6-Trinitrophenol        | 1   | 7.948        | 7.963            | -0.015           | 13879    | 0.2000           | 0.1830             | M     |
| \$ 10 1,2-Dinitrobenzene      | 1   | 8.521        | 8.523            | -0.002           | 20772    | 0.2000           | 0.1645             |       |
| 11 1,3,5-Trinitrobenzene      | 1   | 8.648        | 8.650            | -0.002           | 37864    | 0.2000           | 0.1744             |       |
| 12 1,3-Dinitrobenzene         | 1   | 9.268        | 9.277            | -0.009           | 48740    | 0.2000           | 0.1656             |       |
| 13 Nitrobenzene               | 1   | 9.634        | 9.643            | -0.009           | 28032    | 0.2000           | 0.1466             |       |
| 14 3,5-Dinitroaniline         | 1   | 9.874        | 9.883            | -0.009           | 33202    | 0.2000           | 0.1455             |       |
| 15 Tetryl                     | 1   | 10.008       | 10.017           | -0.009           | 28215    | 0.2000           | 0.1719             |       |
| 16 Nitroglycerin              | 2   | 10.454       | 10.470           | -0.016           | 111863   | 2.00             | 1.75               |       |
| 17 2,4,6-Trinitrotoluene      | 1   | 10.894       | 10.910           | -0.016           | 32823    | 0.2000           | 0.1555             |       |
| 18 4-Amino-2,6-dinitrotoluene | 1   | 11.094       | 11.117           | -0.023           | 22268    | 0.2000           | 0.1437             |       |
| 19 2-Amino-4,6-dinitrotoluene | 1   | 11.354       | 11.370           | -0.016           | 29397    | 0.2000           | 0.1460             |       |
| 20 2,6-Dinitrotoluene         | 1   | 11.494       | 11.517           | -0.023           | 22479    | 0.2000           | 0.1575             |       |
| 21 2,4-Dinitrotoluene         | 1   | 11.668       | 11.683           | -0.015           | 44520    | 0.2000           | 0.1501             |       |
| 22 o-Nitrotoluene             | 1   | 12.488       | 12.517           | -0.029           | 17175    | 0.2000           | 0.1343             |       |
| 23 p-Nitrotoluene             | 1   | 12.901       | 12.937           | -0.036           | 14552    | 0.2000           | 0.1301             |       |
| 24 m-Nitrotoluene             | 1   | 13.474       | 13.517           | -0.043           | 17266    | 0.2000           | 0.1229             |       |
| 25 PETN                       | 2   | 14.641       | 14.697           | -0.056           | 124233   | 2.00             | 1.80               |       |

**QC Flag Legend**

Processing Flags

Review Flags

M - Manually Integrated

Report Date: 20-May-2023 11:22:15

Chrom Revision: 2.3 16-May-2023 16:37:30

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190049.d

Injection Date: 20-May-2023 05:11:34

Instrument ID: CHHPLC\_X3

Operator ID: JZ/JG

Lims ID: 280-176674-B-10-A MS

Worklist Smp#: 49

Client ID: FBQmw-175-230401-GW-R

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 49

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



Eurofins Denver  
Recovery Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\05190049.D  
Lims ID: 280-176674-B-10-A MS  
Client ID: FBQmw-175-230401-GW-R  
Sample Type: MS  
Inject. Date: 20-May-2023 05:11:34 ALS Bottle#: 49 Worklist Smp#: 49  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Sample Info: 280-176674-B-10-A MS  
Operator ID: JZ/JG Instrument ID: CHHPLC\_X3  
Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\8330\_X3.m  
Limit Group: GCSV - 8330  
Last Update: 20-May-2023 11:22:13 Calib Date: 24-Feb-2023 17:41:20  
Integrator: Falcon  
Quant Method: External Standard Quant By: Initial Calibration  
Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
Process Host: CTX1663

First Level Reviewer: LV5D

Date: 20-May-2023 11:15:10

| Compound                 | Amount Added | Amount Recovered | % Rec. |
|--------------------------|--------------|------------------|--------|
| \$ 10 1,2-Dinitrobenzene | 0.2000       | 0.1645           | 82.23  |



## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190049.d  
Injection Date: 20-May-2023 05:11:34 Instrument ID: CHHPLC\_X3  
Lims ID: 280-176674-B-10-A MS  
Client ID: FBQmw-175-230401-GW-R  
Operator ID: JZ/JG ALS Bottle#: 49 Worklist Smp#: 49  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**6 DNX, CAS: 80251-29-2**

Signal: 1


RT: 6.75  
Area: 379225  
Amount: 2.625076  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.75  
Area: 22572  
Amount: 0.156248  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 11:14:54 -06:00:00 (UTC)

Audit Action: Manually Integrated

Audit Reason: Baseline

## Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190049.d

Injection Date: 20-May-2023 05:11:34

Instrument ID: CHHPLC\_X3

Lims ID: 280-176674-B-10-A MS

Client ID: FBQmw-175-230401-GW-R

Operator ID: JZ/JG

ALS Bottle#:

49

Worklist Smp#:

49

Injection Vol: 100.0 ul

Dil. Factor: 1.0000

Method: 8330\_X3

Limit Group: GCSV - 8330

Column: UltraCarb5uODS (20) ( 4.60 mm)

Detector


LC DAD1B, 254 nm

## 4 HMX, CAS: 2691-41-0

Signal: 1


RT: 6.55  
Area: 265382  
Amount: 2.837456  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.55  
Area: 13659  
Amount: 0.146042  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 11:14:49 -06:00:00 (UTC)

Audit Action: Manually Integrated

Audit Reason: Baseline

## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190049.d  
Injection Date: 20-May-2023 05:11:34 Instrument ID: CHHPLC\_X3  
Lims ID: 280-176674-B-10-A MS  
Client ID: FBQmw-175-230401-GW-R  
Operator ID: JZ/JG ALS Bottle#: 49 Worklist Smp#: 49  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**7 MNX, CAS: 5755-27-1**

Signal: 1


RT: 7.18  
Area: 196895  
Amount: 1.501049  
Amount Units: ug/mL

## Processing Integration Results



RT: 7.18  
Area: 25990  
Amount: 0.198137  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 11:14:58 -06:00:00 (UTC)

Audit Action: Manually Integrated

Audit Reason: Baseline

## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190049.d  
Injection Date: 20-May-2023 05:11:34 Instrument ID: CHHPLC\_X3  
Lims ID: 280-176674-B-10-A MS  
Client ID: FBQmw-175-230401-GW-R  
Operator ID: JZ/JG ALS Bottle#: 49 Worklist Smp#: 49  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**8 RDX, CAS: 121-82-4**

Signal: 1


RT: 7.57  
Area: 54141  
Amount: 0.508941  
Amount Units: ug/mL

## Processing Integration Results



RT: 7.57  
Area: 15813  
Amount: 0.148647  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 11:15:03 -06:00:00 (UTC)

Audit Action: Manually Integrated

Audit Reason: Baseline

## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190049.d  
Injection Date: 20-May-2023 05:11:34 Instrument ID: CHHPLC\_X3  
Lims ID: 280-176674-B-10-A MS  
Client ID: FBQmw-175-230401-GW-R  
Operator ID: JZ/JG ALS Bottle#: 49 Worklist Smp#: 49  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector: LC DAD1B, 254 nm

**3 TNX, CAS: 13980-04-6**

Signal: 1


RT: 6.43  
Area: 686004  
Amount: 3.470147  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.43  
Area: 33570  
Amount: 0.169814  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 11:14:47 -06:00:00 (UTC)

Audit Action: Manually Integrated

Audit Reason: Baseline

FORM I  
HPLC/IC ORGANICS ANALYSIS DATA SHEET

Lab Name: Eurofins Denver Job No.: 280-176674-1  
SDG No.: \_\_\_\_\_  
Client Sample ID: FBQmw-175-230401-GW-R MS Lab Sample ID: 280-176674-10 MS  
Matrix: Water Lab File ID: 05190051.D  
Analysis Method: 8330B Date Collected: 05/17/2023 08:58  
Extraction Method: 3535 Date Extracted: 05/19/2023 13:39  
Sample wt/vol: 498.7 (mL) Date Analyzed: 05/20/2023 05:57  
Con. Extract Vol.: 5 (mL) Dilution Factor: 1  
Injection Volume: 100 (uL) GC Column: UltraCarb5uODS ID: 4.6 (mm)  
% Moisture: \_\_\_\_\_ % Solids: \_\_\_\_\_ GPC Cleanup: (Y/N) N  
Cleanup Factor: \_\_\_\_\_  
Analysis Batch No.: 613168 Units: ug/L

| CAS NO.    | COMPOUND NAME              | RESULT | Q    | LOQ | LOD  | DL   |
|------------|----------------------------|--------|------|-----|------|------|
| 6629-29-4  | 2,4-diamino-6-nitrotoluene | 1.25   | M J1 | 1.0 | 0.90 | 0.44 |
| 59229-75-3 | 2,6-diamino-4-nitrotoluene | 1.60   | M    | 1.0 | 0.90 | 0.22 |

| CAS NO.  | SURROGATE          | %REC | Q | LIMITS |
|----------|--------------------|------|---|--------|
| 528-29-0 | 1,2-Dinitrobenzene | 95   | M | 83-119 |

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\05190051.D  
 Lims ID: 280-176674-C-10-A MS  
 Client ID: FBQmw-175-230401-GW-R  
 Sample Type: MS  
 Inject. Date: 20-May-2023 05:57:36 ALS Bottle#: 51 Worklist Smp#: 51  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-176674-C-10-A MS  
 Operator ID: JZ/JG Instrument ID: CHHPLC\_X3  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 20-May-2023 11:22:13 Calib Date: 24-Feb-2023 17:41:20  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1663

First Level Reviewer: LV5D

Date: 20-May-2023 11:15:52

| Compound                     | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/mL | OnCol Amt<br>ug/mL | Flags |
|------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 2 2,6-diamino-4-nitrotoluene | 1   | 6.418        | 6.431            | -0.013           | 37412    | 0.2000           | 0.1591             | M     |
| 5 2,4-diamino-6-nitrotoluene | 1   | 6.604        | 6.618            | -0.014           | 18219    | 0.2000           | 0.1244             | M     |
| \$ 10 1,2-Dinitrobenzene     | 1   | 8.524        | 8.523            | 0.001            | 24080    | 0.2000           | 0.1906             | M     |

**QC Flag Legend**

Processing Flags

Review Flags

M - Manually Integrated

Report Date: 20-May-2023 11:22:16

Chrom Revision: 2.3 16-May-2023 16:37:30

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190051.d

Injection Date: 20-May-2023 05:57:36

Instrument ID: CHHPLC\_X3

Operator ID: JZ/JG

Lims ID: 280-176674-C-10-A MS

Worklist Smp#: 51

Client ID: FBQmw-175-230401-GW-R

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 51

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1





Eurofins Denver  
Recovery Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\05190051.D  
Lims ID: 280-176674-C-10-A MS  
Client ID: FBQmw-175-230401-GW-R  
Sample Type: MS  
Inject. Date: 20-May-2023 05:57:36 ALS Bottle#: 51 Worklist Smp#: 51  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Sample Info: 280-176674-C-10-A MS  
Operator ID: JZ/JG Instrument ID: CHHPLC\_X3  
Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\8330\_X3.m  
Limit Group: GCSV - 8330  
Last Update: 20-May-2023 11:22:13 Calib Date: 24-Feb-2023 17:41:20  
Integrator: Falcon  
Quant Method: External Standard Quant By: Initial Calibration  
Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
Process Host: CTX1663

First Level Reviewer: LV5D

Date: 20-May-2023 11:15:52

| Compound                 | Amount Added | Amount Recovered | % Rec. |
|--------------------------|--------------|------------------|--------|
| \$ 10 1,2-Dinitrobenzene | 0.2000       | 0.1906           | 95.32  |

## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190051.d  
Injection Date: 20-May-2023 05:57:36 Instrument ID: CHHPLC\_X3  
Lims ID: 280-176674-C-10-A MS  
Client ID: FBQmw-175-230401-GW-R  
Operator ID: JZ/JG ALS Bottle#: 51 Worklist Smp#: 51  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

## \$ 10 1,2-Dinitrobenzene, CAS: 528-29-0

Signal: 1


RT: 8.52  
Area: 28021  
Amount: 0.221845  
Amount Units: ug/mL

## Processing Integration Results



RT: 8.52  
Area: 24080  
Amount: 0.190644  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 11:15:46 -06:00:00 (UTC)

Audit Action: Manually Integrated

Audit Reason: Baseline

## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190051.d  
Injection Date: 20-May-2023 05:57:36 Instrument ID: CHHPLC\_X3  
Lims ID: 280-176674-C-10-A MS  
Client ID: FBQmw-175-230401-GW-R  
Operator ID: JZ/JG ALS Bottle#: 51 Worklist Smp#: 51  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**5 2,4-diamino-6-nitrotoluene, CAS: 6629-29-4**

Signal: 1


RT: 6.60  
Area: 825706  
Amount: 5.637291  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.60  
Area: 18219  
Amount: 0.124385  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 11:15:51 -06:00:00 (UTC)

Audit Action: Manually Integrated

Audit Reason: Baseline

## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190051.d  
Injection Date: 20-May-2023 05:57:36 Instrument ID: CHHPLC\_X3  
Lims ID: 280-176674-C-10-A MS  
Client ID: FBQmw-175-230401-GW-R  
Operator ID: JZ/JG ALS Bottle#: 51 Worklist Smp#: 51  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**2,2,6-diamino-4-nitrotoluene, CAS: 59229-75-3**

Signal: 1


RT: 6.42  
Area: 743677  
Amount: 3.233938  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.42  
Area: 37412  
Amount: 0.159116  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 11:15:49 -06:00:00 (UTC)

Audit Action: Manually Integrated

Audit Reason: Baseline

FORM I  
HPLC/IC ORGANICS ANALYSIS DATA SHEET

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Client Sample ID: FBQmw-175-230401-GW-R MSD Lab Sample ID: 280-176674-10 MSD

Matrix: Water Lab File ID: 05190050.D

Analysis Method: 8330B Date Collected: 05/17/2023 08:58

Extraction Method: 3535 Date Extracted: 05/19/2023 13:39

Sample wt/vol: 457.4 (mL) Date Analyzed: 05/20/2023 05:34

Con. Extract Vol.: 5 (mL) Dilution Factor: 1

Injection Volume: 100 (uL) GC Column: UltraCarb5uODS ID: 4.6 (mm)

% Moisture: \_\_\_\_\_ % Solids: \_\_\_\_\_ GPC Cleanup: (Y/N) N

Cleanup Factor: \_\_\_\_\_

Analysis Batch No.: 613168 Units: ug/L

| CAS NO.    | COMPOUND NAME              | RESULT | Q    | LOQ  | LOD   | DL    |
|------------|----------------------------|--------|------|------|-------|-------|
| 99-35-4    | 1,3,5-Trinitrobenzene      | 2.21   | J1   | 0.23 | 0.22  | 0.092 |
| 99-65-0    | 1,3-Dinitrobenzene         | 2.13   | J1   | 0.12 | 0.11  | 0.040 |
| 118-96-7   | 2,4,6-Trinitrotoluene      | 2.00   | J1   | 0.12 | 0.11  | 0.049 |
| 121-14-2   | 2,4-Dinitrotoluene         | 2.00   | J1   | 0.11 | 0.087 | 0.030 |
| 606-20-2   | 2,6-Dinitrotoluene         | 2.02   | J1   | 0.11 | 0.087 | 0.044 |
| 35572-78-2 | 2-Amino-4,6-dinitrotoluene | 1.96   | J1   | 0.12 | 0.11  | 0.055 |
| 88-72-2    | 2-Nitrotoluene             | 1.78   | J1   | 0.23 | 0.22  | 0.093 |
| 618-87-1   | 3,5-Dinitroaniline         | 1.86   | J1   | 0.44 | 0.33  | 0.14  |
| 99-08-1    | 3-Nitrotoluene             | 1.67   | J1   | 0.44 | 0.38  | 0.21  |
| 19406-51-0 | 4-Amino-2,6-dinitrotoluene | 1.87   | J1   | 0.16 | 0.13  | 0.063 |
| 99-99-0    | 4-Nitrotoluene             | 1.75   | J1   | 0.45 | 0.44  | 0.11  |
| 80251-29-2 | DNX                        | 1.93   | M J1 | 0.55 | 0.27  | 0.11  |
| 2691-41-0  | HMX                        | 1.77   | M    | 0.23 | 0.22  | 0.096 |
| 5755-27-1  | MNX                        | 2.37   | M    | 0.55 | 0.32  | 0.10  |
| 98-95-3    | Nitrobenzene               | 1.94   | J1   | 0.23 | 0.22  | 0.099 |
| 55-63-0    | Nitroglycerin              | 21.5   |      | 2.3  | 2.2   | 1.0   |
| 78-11-5    | PETN                       | 22.4   | J1   | 1.2  | 1.1   | 0.49  |
| 121-82-4   | RDX                        | 1.83   | M J1 | 0.23 | 0.22  | 0.056 |
| 479-45-8   | Tetryl                     | 2.23   | J1   | 0.12 | 0.11  | 0.035 |
| 13980-04-6 | TNX                        | 1.98   | M    | 0.55 | 0.27  | 0.087 |

| CAS NO.  | SURROGATE          | %REC | Q | LIMITS |
|----------|--------------------|------|---|--------|
| 528-29-0 | 1,2-Dinitrobenzene | 98   |   | 83-119 |

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\05190050.D  
 Lims ID: 280-176674-B-10-B MSD  
 Client ID: FBQmw-175-230401-GW-R  
 Sample Type: MSD  
 Inject. Date: 20-May-2023 05:34:33 ALS Bottle#: 50 Worklist Smp#: 50  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-176674-B-10-B MSD  
 Operator ID: JZ/JG Instrument ID: CHHPLC\_X3  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 20-May-2023 11:22:13 Calib Date: 24-Feb-2023 17:41:20  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1663

First Level Reviewer: LV5D

Date: 20-May-2023 11:15:40

| Compound                      | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/mL | OnCol Amt<br>ug/mL | Flags |
|-------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 3 TNX                         | 1   | 6.433        | 6.430            | 0.003            | 35886    | 0.2006           | 0.1815             | M     |
| 4 HMX                         | 1   | 6.546        | 6.550            | -0.004           | 15136    | 0.2000           | 0.1618             | M     |
| 6 DNx                         | 1   | 6.753        | 6.757            | -0.004           | 25551    | 0.2002           | 0.1769             | M     |
| 7 MNx                         | 1   | 7.179        | 7.190            | -0.011           | 28473    | 0.2334           | 0.2171             | M     |
| 8 RDX                         | 1   | 7.566        | 7.570            | -0.004           | 17855    | 0.2000           | 0.1678             | M     |
| 9 2,4,6-Trinitrophenol        | 1   | 7.946        | 7.963            | -0.017           | 15361    | 0.2000           | 0.2026             | M     |
| \$ 10 1,2-Dinitrobenzene      | 1   | 8.519        | 8.523            | -0.004           | 24659    | 0.2000           | 0.1952             |       |
| 11 1,3,5-Trinitrobenzene      | 1   | 8.646        | 8.650            | -0.004           | 43820    | 0.2000           | 0.2018             |       |
| 12 1,3-Dinitrobenzene         | 1   | 9.266        | 9.277            | -0.011           | 57363    | 0.2000           | 0.1948             |       |
| 13 Nitrobenzene               | 1   | 9.639        | 9.643            | -0.004           | 33936    | 0.2000           | 0.1774             |       |
| 14 3,5-Dinitroaniline         | 1   | 9.879        | 9.883            | -0.004           | 38777    | 0.2000           | 0.1699             |       |
| 15 Tetryl                     | 1   | 10.006       | 10.017           | -0.011           | 33421    | 0.2000           | 0.2036             |       |
| 16 Nitroglycerin              | 2   | 10.466       | 10.470           | -0.004           | 125847   | 2.00             | 1.96               |       |
| 17 2,4,6-Trinitrotoluene      | 1   | 10.899       | 10.910           | -0.011           | 38698    | 0.2000           | 0.1834             |       |
| 18 4-Amino-2,6-dinitrotoluene | 1   | 11.099       | 11.117           | -0.018           | 26512    | 0.2000           | 0.1711             |       |
| 19 2-Amino-4,6-dinitrotoluene | 1   | 11.359       | 11.370           | -0.011           | 36151    | 0.2000           | 0.1795             |       |
| 20 2,6-Dinitrotoluene         | 1   | 11.506       | 11.517           | -0.011           | 26414    | 0.2000           | 0.1850             |       |
| 21 2,4-Dinitrotoluene         | 1   | 11.673       | 11.683           | -0.010           | 54377    | 0.2000           | 0.1833             |       |
| 22 o-Nitrotoluene             | 1   | 12.499       | 12.517           | -0.018           | 20853    | 0.2000           | 0.1630             |       |
| 23 p-Nitrotoluene             | 1   | 12.919       | 12.937           | -0.018           | 17936    | 0.2000           | 0.1603             |       |
| 24 m-Nitrotoluene             | 1   | 13.499       | 13.517           | -0.018           | 21452    | 0.2000           | 0.1527             |       |
| 25 PETN                       | 2   | 14.673       | 14.697           | -0.024           | 141016   | 2.00             | 2.05               |       |

**QC Flag Legend**

Processing Flags

Review Flags

M - Manually Integrated

Report Date: 20-May-2023 11:22:15

Chrom Revision: 2.3 16-May-2023 16:37:30

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190050.d

Injection Date: 20-May-2023 05:34:33

Instrument ID: CHHPLC\_X3

Operator ID: JZ/JG

Lims ID: 280-176674-B-10-B MSD

Worklist Smp#: 50

Client ID: FBQmw-175-230401-GW-R

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 50

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1



Eurofins Denver  
Recovery Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\05190050.D  
Lims ID: 280-176674-B-10-B MSD  
Client ID: FBQmw-175-230401-GW-R  
Sample Type: MSD  
Inject. Date: 20-May-2023 05:34:33 ALS Bottle#: 50 Worklist Smp#: 50  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Sample Info: 280-176674-B-10-B MSD  
Operator ID: JZ/JG Instrument ID: CHHPLC\_X3  
Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\8330\_X3.m  
Limit Group: GCSV - 8330  
Last Update: 20-May-2023 11:22:13 Calib Date: 24-Feb-2023 17:41:20  
Integrator: Falcon  
Quant Method: External Standard Quant By: Initial Calibration  
Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
Process Host: CTX1663

First Level Reviewer: LV5D

Date: 20-May-2023 11:15:40

| Compound                 | Amount Added | Amount Recovered | % Rec. |
|--------------------------|--------------|------------------|--------|
| \$ 10 1,2-Dinitrobenzene | 0.2000       | 0.1952           | 97.61  |



## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190050.d  
Injection Date: 20-May-2023 05:34:33 Instrument ID: CHHPLC\_X3  
Lims ID: 280-176674-B-10-B MSD  
Client ID: FBQmw-175-230401-GW-R  
Operator ID: JZ/JG ALS Bottle#: 50 Worklist Smp#: 50  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**6 DNX, CAS: 80251-29-2**

Signal: 1


RT: 6.75  
Area: 316028  
Amount: 2.187613  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.75  
Area: 25551  
Amount: 0.176869  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 11:15:24 -06:00:00 (UTC)

Audit Action: Manually Integrated

Audit Reason: Baseline

## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190050.d  
Injection Date: 20-May-2023 05:34:33 Instrument ID: CHHPLC\_X3  
Lims ID: 280-176674-B-10-B MSD  
Client ID: FBQmw-175-230401-GW-R  
Operator ID: JZ/JG ALS Bottle#: 50 Worklist Smp#: 50  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

## 4 HMX, CAS: 2691-41-0

Signal: 1


RT: 6.55  
Area: 240862  
Amount: 2.575289  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.55  
Area: 15136  
Amount: 0.161834  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 11:15:18 -06:00:00 (UTC)

Audit Action: Manually Integrated

Audit Reason: Baseline

## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190050.d  
Injection Date: 20-May-2023 05:34:33 Instrument ID: CHHPLC\_X3  
Lims ID: 280-176674-B-10-B MSD  
Client ID: FBQmw-175-230401-GW-R  
Operator ID: JZ/JG ALS Bottle#: 50 Worklist Smp#: 50  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**7 MNX, CAS: 5755-27-1**

Signal: 1


RT: 7.18  
Area: 138248  
Amount: 1.053948  
Amount Units: ug/mL

## Processing Integration Results



RT: 7.18  
Area: 28473  
Amount: 0.217067  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 11:15:29 -06:00:00 (UTC)

Audit Action: Manually Integrated

Audit Reason: Baseline

## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190050.d  
Injection Date: 20-May-2023 05:34:33 Instrument ID: CHHPLC\_X3  
Lims ID: 280-176674-B-10-B MSD  
Client ID: FBQmw-175-230401-GW-R  
Operator ID: JZ/JG ALS Bottle#: 50 Worklist Smp#: 50  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**8 RDX, CAS: 121-82-4**

Signal: 1


RT: 7.57  
Area: 31817  
Amount: 0.299089  
Amount Units: ug/mL

## Processing Integration Results



RT: 7.57  
Area: 17855  
Amount: 0.167842  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 11:15:33 -06:00:00 (UTC)

Audit Action: Manually Integrated

Audit Reason: Baseline

## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190050.d  
Injection Date: 20-May-2023 05:34:33 Instrument ID: CHHPLC\_X3  
Lims ID: 280-176674-B-10-B MSD  
Client ID: FBQmw-175-230401-GW-R  
Operator ID: JZ/JG ALS Bottle#: 50 Worklist Smp#: 50  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector: LC DAD1B, 254 nm

**3 TNX, CAS: 13980-04-6**

Signal: 1


RT: 6.43  
Area: 622576  
Amount: 3.149296  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.43  
Area: 35886  
Amount: 0.181529  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 11:15:16 -06:00:00 (UTC)

Audit Action: Manually Integrated

Audit Reason: Baseline

FORM I  
HPLC/IC ORGANICS ANALYSIS DATA SHEET

Lab Name: Eurofins Denver Job No.: 280-176674-1  
SDG No.: \_\_\_\_\_  
Client Sample ID: FBQmw-175-230401-GW-R MSD Lab Sample ID: 280-176674-10 MSD  
Matrix: Water Lab File ID: 05190052.D  
Analysis Method: 8330B Date Collected: 05/17/2023 08:58  
Extraction Method: 3535 Date Extracted: 05/19/2023 13:39  
Sample wt/vol: 498.8 (mL) Date Analyzed: 05/20/2023 06:20  
Con. Extract Vol.: 5 (mL) Dilution Factor: 1  
Injection Volume: 100 (uL) GC Column: UltraCarb5uODS ID: 4.6 (mm)  
% Moisture: \_\_\_\_\_ % Solids: \_\_\_\_\_ GPC Cleanup: (Y/N) N  
Cleanup Factor: \_\_\_\_\_  
Analysis Batch No.: 613168 Units: ug/L

| CAS NO.    | COMPOUND NAME              | RESULT | Q    | LOQ | LOD  | DL   |
|------------|----------------------------|--------|------|-----|------|------|
| 6629-29-4  | 2,4-diamino-6-nitrotoluene | 1.03   | M J1 | 1.0 | 0.90 | 0.44 |
| 59229-75-3 | 2,6-diamino-4-nitrotoluene | 1.32   | M J1 | 1.0 | 0.90 | 0.22 |

| CAS NO.  | SURROGATE          | %REC | Q | LIMITS |
|----------|--------------------|------|---|--------|
| 528-29-0 | 1,2-Dinitrobenzene | 109  |   | 83-119 |

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\05190052.D  
 Lims ID: 280-176674-D-10-A MSD  
 Client ID: FBQmw-175-230401-GW-R  
 Sample Type: MSD  
 Inject. Date: 20-May-2023 06:20:33 ALS Bottle#: 52 Worklist Smp#: 52  
 Injection Vol: 100.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-176674-D-10-A MSD  
 Operator ID: JZ/JG Instrument ID: CHHPLC\_X3  
 Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\8330\_X3.m  
 Limit Group: GCSV - 8330  
 Last Update: 20-May-2023 11:22:13 Calib Date: 24-Feb-2023 17:41:20  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
 Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
 Process Host: CTX1663

First Level Reviewer: LV5D

Date: 20-May-2023 11:15:58

| Compound                     | Det | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/mL | OnCol Amt<br>ug/mL | Flags |
|------------------------------|-----|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 2 2,6-diamino-4-nitrotoluene | 1   | 6.417        | 6.431            | -0.014           | 31161    | 0.2000           | 0.1319             | M     |
| 5 2,4-diamino-6-nitrotoluene | 1   | 6.604        | 6.618            | -0.014           | 15008    | 0.2000           | 0.1025             | M     |
| \$ 10 1,2-Dinitrobenzene     | 1   | 8.524        | 8.523            | 0.001            | 27619    | 0.2000           | 0.2187             |       |

**QC Flag Legend**

Processing Flags

Review Flags

M - Manually Integrated

Report Date: 20-May-2023 11:22:16

Chrom Revision: 2.3 16-May-2023 16:37:30

Eurofins Denver

Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190052.d

Injection Date: 20-May-2023 06:20:33

Instrument ID: CHHPLC\_X3

Operator ID: JZ/JG

Lims ID: 280-176674-D-10-A MSD

Worklist Smp#: 52

Client ID: FBQmw-175-230401-GW-R

Injection Vol: 100.0 ul

Dil. Factor: 1.0000


ALS Bottle#: 52

Method: 8330\_X3

Limit Group: GCSV - 8330


Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Target: 1



Column: UltraCarb5uODS (20) ( 4.60 mm)

Y Scaling: Method Defined: Scale to the Nth Largest Peak: 1





Eurofins Denver  
Recovery Report

Data File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\05190052.D  
Lims ID: 280-176674-D-10-A MSD  
Client ID: FBQmw-175-230401-GW-R  
Sample Type: MSD  
Inject. Date: 20-May-2023 06:20:33 ALS Bottle#: 52 Worklist Smp#: 52  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Sample Info: 280-176674-D-10-A MSD  
Operator ID: JZ/JG Instrument ID: CHHPLC\_X3  
Method: \\chromfs\Denver\ChromData\CHHPLC\_X\20230519-121647.b\8330\_X3.m  
Limit Group: GCSV - 8330  
Last Update: 20-May-2023 11:22:13 Calib Date: 24-Feb-2023 17:41:20  
Integrator: Falcon  
Quant Method: External Standard Quant By: Initial Calibration  
Last ICal File: \\chromfs\Denver\ChromData\CHHPLC\_X\20230224-118952.b\02240018.D  
Column 1 : UltraCarb5uODS (20) ( 4.60 mm) Det: LC DAD1B, 254 nm  
Process Host: CTX1663

First Level Reviewer: LV5D

Date: 20-May-2023 11:15:58

| Compound                 | Amount Added | Amount Recovered | % Rec. |
|--------------------------|--------------|------------------|--------|
| \$ 10 1,2-Dinitrobenzene | 0.2000       | 0.2187           | 109.33 |

## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190052.d  
Injection Date: 20-May-2023 06:20:33 Instrument ID: CHHPLC\_X3  
Lims ID: 280-176674-D-10-A MSD  
Client ID: FBQmw-175-230401-GW-R  
Operator ID: JZ/JG ALS Bottle#: 52 Worklist Smp#: 52  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**5 2,4-diamino-6-nitrotoluene, CAS: 6629-29-4**

Signal: 1


RT: 6.60  
Area: 293570  
Amount: 2.004272  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.60  
Area: 15008  
Amount: 0.102463  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 11:15:57 -06:00:00 (UTC)

Audit Action: Manually Integrated

Audit Reason: Baseline

## Eurofins Denver


Data File: \\chromfs\denver\chromdata\chhplc\_x\20230519-121647.b\05190052.d  
Injection Date: 20-May-2023 06:20:33 Instrument ID: CHHPLC\_X3  
Lims ID: 280-176674-D-10-A MSD  
Client ID: FBQmw-175-230401-GW-R  
Operator ID: JZ/JG ALS Bottle#: 52 Worklist Smp#: 52  
Injection Vol: 100.0 ul Dil. Factor: 1.0000  
Method: 8330\_X3 Limit Group: GCSV - 8330  
Column: UltraCarb5uODS (20) ( 4.60 mm) Detector LC DAD1B, 254 nm

**2,2,6-diamino-4-nitrotoluene, CAS: 59229-75-3**

Signal: 1


RT: 6.42  
Area: 644343  
Amount: 2.801474  
Amount Units: ug/mL

## Processing Integration Results



RT: 6.42  
Area: 31161  
Amount: 0.131902  
Amount Units: ug/mL

## Manual Integration Results



Reviewer: LV5D, 20-May-2023 11:15:56 -06:00:00 (UTC)

Audit Action: Manually Integrated

Audit Reason: Baseline

## HPLC/IC ANALYSIS RUN LOG

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Instrument ID: CHHPLC\_X3 Start Date: 01/17/2023 00:33Analysis Batch Number: 599427 End Date: 01/17/2023 07:10

| LAB SAMPLE ID     | CLIENT SAMPLE ID | DATE ANALYZED    | DILUTION<br>FACTOR | LAB FILE ID | COLUMN ID               |
|-------------------|------------------|------------------|--------------------|-------------|-------------------------|
| IC 280-599427/11  |                  | 01/17/2023 00:33 | 1                  | 01160039.D  | UltraCarb5uODS 4.6 (mm) |
| IC 280-599427/12  |                  | 01/17/2023 00:56 | 1                  | 01160040.D  | UltraCarb5uODS 4.6 (mm) |
| IC 280-599427/13  |                  | 01/17/2023 01:19 | 1                  | 01160041.D  | UltraCarb5uODS 4.6 (mm) |
| IC 280-599427/14  |                  | 01/17/2023 01:43 | 1                  | 01160042.D  | UltraCarb5uODS 4.6 (mm) |
| IC 280-599427/15  |                  | 01/17/2023 02:06 | 1                  | 01160043.D  | UltraCarb5uODS 4.6 (mm) |
| IC 280-599427/16  |                  | 01/17/2023 02:30 | 1                  | 01160044.D  | UltraCarb5uODS 4.6 (mm) |
| IC 280-599427/17  |                  | 01/17/2023 02:53 | 1                  | 01160045.D  | UltraCarb5uODS 4.6 (mm) |
| IC 280-599427/18  |                  | 01/17/2023 03:16 | 1                  | 01160046.D  | UltraCarb5uODS 4.6 (mm) |
| ICV 280-599427/19 |                  | 01/17/2023 03:40 | 1                  | 01160047.D  | UltraCarb5uODS 4.6 (mm) |
| IC 280-599427/20  |                  | 01/17/2023 04:03 | 1                  |             | UltraCarb5uODS 4.6 (mm) |
| IC 280-599427/21  |                  | 01/17/2023 04:27 | 1                  |             | UltraCarb5uODS 4.6 (mm) |
| IC 280-599427/22  |                  | 01/17/2023 04:50 | 1                  |             | UltraCarb5uODS 4.6 (mm) |
| IC 280-599427/23  |                  | 01/17/2023 05:13 | 1                  |             | UltraCarb5uODS 4.6 (mm) |
| IC 280-599427/24  |                  | 01/17/2023 05:37 | 1                  |             | UltraCarb5uODS 4.6 (mm) |
| IC 280-599427/25  |                  | 01/17/2023 06:00 | 1                  |             | UltraCarb5uODS 4.6 (mm) |
| IC 280-599427/26  |                  | 01/17/2023 06:24 | 1                  |             | UltraCarb5uODS 4.6 (mm) |
| IC 280-599427/27  |                  | 01/17/2023 06:47 | 1                  |             | UltraCarb5uODS 4.6 (mm) |
| ICV 280-599427/28 |                  | 01/17/2023 07:10 | 1                  |             | UltraCarb5uODS 4.6 (mm) |

## HPLC/IC ANALYSIS RUN LOG

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Instrument ID: CHHPLC\_X3 Start Date: 02/08/2023 15:38Analysis Batch Number: 601664 End Date: 02/08/2023 19:05

| LAB SAMPLE ID     | CLIENT SAMPLE ID | DATE ANALYZED    | DILUTION<br>FACTOR | LAB FILE ID | COLUMN ID               |
|-------------------|------------------|------------------|--------------------|-------------|-------------------------|
| IC 280-601664/11  |                  | 02/08/2023 15:38 | 1                  | 02080011.D  | UltraCarb5uODS 4.6 (mm) |
| IC 280-601664/12  |                  | 02/08/2023 16:01 | 1                  | 02080012.D  | UltraCarb5uODS 4.6 (mm) |
| IC 280-601664/13  |                  | 02/08/2023 16:24 | 1                  | 02080013.D  | UltraCarb5uODS 4.6 (mm) |
| IC 280-601664/14  |                  | 02/08/2023 16:47 | 1                  | 02080014.D  | UltraCarb5uODS 4.6 (mm) |
| IC 280-601664/15  |                  | 02/08/2023 17:10 | 1                  | 02080015.D  | UltraCarb5uODS 4.6 (mm) |
| IC 280-601664/16  |                  | 02/08/2023 17:33 | 1                  | 02080016.D  | UltraCarb5uODS 4.6 (mm) |
| IC 280-601664/17  |                  | 02/08/2023 17:56 | 1                  | 02080017.D  | UltraCarb5uODS 4.6 (mm) |
| IC 280-601664/18  |                  | 02/08/2023 18:19 | 1                  | 02080018.D  | UltraCarb5uODS 4.6 (mm) |
| IC 280-601664/19  |                  | 02/08/2023 18:42 | 1                  | 02080019.D  | UltraCarb5uODS 4.6 (mm) |
| ICV 280-601664/20 |                  | 02/08/2023 19:05 | 1                  | 02080020.D  | UltraCarb5uODS 4.6 (mm) |

## HPLC/IC ANALYSIS RUN LOG

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Instrument ID: CHHPLC\_X3 Start Date: 02/24/2023 15:00Analysis Batch Number: 603284 End Date: 02/24/2023 18:04

| LAB SAMPLE ID     | CLIENT SAMPLE ID | DATE ANALYZED    | DILUTION<br>FACTOR | LAB FILE ID | COLUMN ID               |
|-------------------|------------------|------------------|--------------------|-------------|-------------------------|
| IC 280-603284/11  |                  | 02/24/2023 15:00 | 1                  | 02240011.D  | UltraCarb5uODS 4.6 (mm) |
| IC 280-603284/12  |                  | 02/24/2023 15:23 | 1                  | 02240012.D  | UltraCarb5uODS 4.6 (mm) |
| IC 280-603284/13  |                  | 02/24/2023 15:46 | 1                  | 02240013.D  | UltraCarb5uODS 4.6 (mm) |
| IC 280-603284/14  |                  | 02/24/2023 16:09 | 1                  | 02240014.D  | UltraCarb5uODS 4.6 (mm) |
| IC 280-603284/15  |                  | 02/24/2023 16:32 | 1                  | 02240015.D  | UltraCarb5uODS 4.6 (mm) |
| IC 280-603284/16  |                  | 02/24/2023 16:55 | 1                  | 02240016.D  | UltraCarb5uODS 4.6 (mm) |
| IC 280-603284/17  |                  | 02/24/2023 17:18 | 1                  | 02240017.D  | UltraCarb5uODS 4.6 (mm) |
| IC 280-603284/18  |                  | 02/24/2023 17:41 | 1                  | 02240018.D  | UltraCarb5uODS 4.6 (mm) |
| ICV 280-603284/19 |                  | 02/24/2023 18:04 | 1                  | 02240019.D  | UltraCarb5uODS 4.6 (mm) |

## HPLC/IC ANALYSIS RUN LOG

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Instrument ID: CHHPLC\_X5 Start Date: 04/28/2023 18:24Analysis Batch Number: 610603 End Date: 04/29/2023 04:52

| LAB SAMPLE ID     | CLIENT SAMPLE ID | DATE ANALYZED    | DILUTION<br>FACTOR | LAB FILE ID | COLUMN ID               |
|-------------------|------------------|------------------|--------------------|-------------|-------------------------|
| IC 280-610603/10  |                  | 04/28/2023 18:24 | 1                  | 04280010.D  | Luna-phenylhex 4.6 (mm) |
| IC 280-610603/11  |                  | 04/28/2023 18:59 | 1                  | 04280011.D  | Luna-phenylhex 4.6 (mm) |
| IC 280-610603/12  |                  | 04/28/2023 19:34 | 1                  | 04280012.D  | Luna-phenylhex 4.6 (mm) |
| IC 280-610603/13  |                  | 04/28/2023 20:08 | 1                  | 04280013.D  | Luna-phenylhex 4.6 (mm) |
| IC 280-610603/14  |                  | 04/28/2023 20:43 | 1                  | 04280014.D  | Luna-phenylhex 4.6 (mm) |
| IC 280-610603/15  |                  | 04/28/2023 21:18 | 1                  | 04280015.D  | Luna-phenylhex 4.6 (mm) |
| IC 280-610603/16  |                  | 04/28/2023 21:53 | 1                  | 04280016.D  | Luna-phenylhex 4.6 (mm) |
| IC 280-610603/17  |                  | 04/28/2023 22:28 | 1                  | 04280017.D  | Luna-phenylhex 4.6 (mm) |
| IC 280-610603/18  |                  | 04/28/2023 23:03 | 1                  | 04280018.D  | Luna-phenylhex 4.6 (mm) |
| ICV 280-610603/19 |                  | 04/28/2023 23:38 | 1                  | 04280019.D  | Luna-phenylhex 4.6 (mm) |
| IC 280-610603/20  |                  | 04/29/2023 00:13 | 1                  | 04280020.D  | Luna-phenylhex 4.6 (mm) |
| IC 280-610603/21  |                  | 04/29/2023 00:48 | 1                  | 04280021.D  | Luna-phenylhex 4.6 (mm) |
| IC 280-610603/22  |                  | 04/29/2023 01:23 | 1                  | 04280022.D  | Luna-phenylhex 4.6 (mm) |
| IC 280-610603/23  |                  | 04/29/2023 01:58 | 1                  | 04280023.D  | Luna-phenylhex 4.6 (mm) |
| IC 280-610603/24  |                  | 04/29/2023 02:33 | 1                  | 04280024.D  | Luna-phenylhex 4.6 (mm) |
| IC 280-610603/25  |                  | 04/29/2023 03:08 | 1                  | 04280025.D  | Luna-phenylhex 4.6 (mm) |
| IC 280-610603/26  |                  | 04/29/2023 03:43 | 1                  | 04280026.D  | Luna-phenylhex 4.6 (mm) |
| IC 280-610603/27  |                  | 04/29/2023 04:18 | 1                  | 04280027.D  | Luna-phenylhex 4.6 (mm) |
| ICV 280-610603/28 |                  | 04/29/2023 04:52 | 1                  | 04280028.D  | Luna-phenylhex 4.6 (mm) |

## HPLC/IC ANALYSIS RUN LOG

Lab Name: Eurofins Denver

Job No.: 280-176674-1

SDG No.:

Instrument ID: CHHPLC\_X3

Start Date: 05/19/2023 22:41

Analysis Batch Number: 613168

End Date: 05/20/2023 12:04

| LAB SAMPLE ID       | CLIENT SAMPLE ID             | DATE ANALYZED    | DILUTION<br>FACTOR | LAB FILE ID | COLUMN ID               |
|---------------------|------------------------------|------------------|--------------------|-------------|-------------------------|
| CCV 280-613168/32   |                              | 05/19/2023 22:41 | 1                  | 05190032.D  | UltraCarb5uODS 4.6 (mm) |
| CCV 280-613168/33   |                              | 05/19/2023 23:04 | 1                  | 05190033.D  | UltraCarb5uODS 4.6 (mm) |
| MB 280-613095/1-A   |                              | 05/19/2023 23:27 | 1                  | 05190034.D  | UltraCarb5uODS 4.6 (mm) |
| LCS 280-613095/2-A  |                              | 05/19/2023 23:50 | 1                  | 05190035.D  | UltraCarb5uODS 4.6 (mm) |
| LCSD 280-613095/3-A |                              | 05/20/2023 00:13 | 1                  | 05190036.D  | UltraCarb5uODS 4.6 (mm) |
| LCS 280-613095/4-A  |                              | 05/20/2023 00:36 | 1                  | 05190037.D  | UltraCarb5uODS 4.6 (mm) |
| LCSD 280-613095/5-A |                              | 05/20/2023 00:59 | 1                  | 05190038.D  | UltraCarb5uODS 4.6 (mm) |
| ZZZZZ               |                              | 05/20/2023 01:22 | 1                  |             | UltraCarb5uODS 4.6 (mm) |
| ZZZZZ               |                              | 05/20/2023 02:07 | 1                  |             | UltraCarb5uODS 4.6 (mm) |
| ZZZZZ               |                              | 05/20/2023 02:30 | 1                  |             | UltraCarb5uODS 4.6 (mm) |
| 280-176674-2        | FWGmw-015-230401-GW          | 05/20/2023 02:53 | 1                  | 05190043.D  | UltraCarb5uODS 4.6 (mm) |
| CCV 280-613168/44   |                              | 05/20/2023 03:16 | 1                  | 05190044.D  | UltraCarb5uODS 4.6 (mm) |
| CCV 280-613168/45   |                              | 05/20/2023 03:39 | 1                  | 05190045.D  | UltraCarb5uODS 4.6 (mm) |
| 280-176674-4        | LL3mw-241-230401-GW          | 05/20/2023 04:02 | 1                  | 05190046.D  | UltraCarb5uODS 4.6 (mm) |
| 280-176674-7        | FBQmw-173-230401-GW-R        | 05/20/2023 04:25 | 1                  | 05190047.D  | UltraCarb5uODS 4.6 (mm) |
| 280-176674-10       | FBQmw-175-230401-GW-R        | 05/20/2023 04:48 | 1                  | 05190048.D  | UltraCarb5uODS 4.6 (mm) |
| 280-176674-10 MS    | FBQmw-175-230401-GW-R<br>MS  | 05/20/2023 05:11 | 1                  | 05190049.D  | UltraCarb5uODS 4.6 (mm) |
| 280-176674-10 MSD   | FBQmw-175-230401-GW-R<br>MSD | 05/20/2023 05:34 | 1                  | 05190050.D  | UltraCarb5uODS 4.6 (mm) |
| 280-176674-10 MS    | FBQmw-175-230401-GW-R<br>MS  | 05/20/2023 05:57 | 1                  | 05190051.D  | UltraCarb5uODS 4.6 (mm) |
| 280-176674-10 MSD   | FBQmw-175-230401-GW-R<br>MSD | 05/20/2023 06:20 | 1                  | 05190052.D  | UltraCarb5uODS 4.6 (mm) |
| 280-176674-11       | FBQmw-175-230402-GW-R        | 05/20/2023 06:43 | 1                  | 05190053.D  | UltraCarb5uODS 4.6 (mm) |
| ZZZZZ               |                              | 05/20/2023 07:06 | 1                  |             | UltraCarb5uODS 4.6 (mm) |
| ZZZZZ               |                              | 05/20/2023 07:29 | 1                  |             | UltraCarb5uODS 4.6 (mm) |
| CCV 280-613168/56   |                              | 05/20/2023 07:52 | 1                  | 05190056.D  | UltraCarb5uODS 4.6 (mm) |
| CCV 280-613168/57   |                              | 05/20/2023 08:15 | 1                  | 05190057.D  | UltraCarb5uODS 4.6 (mm) |
| ZZZZZ               |                              | 05/20/2023 08:38 | 1                  |             | UltraCarb5uODS 4.6 (mm) |
| ZZZZZ               |                              | 05/20/2023 09:01 | 1                  |             | UltraCarb5uODS 4.6 (mm) |
| ZZZZZ               |                              | 05/20/2023 09:24 | 1                  |             | UltraCarb5uODS 4.6 (mm) |
| ZZZZZ               |                              | 05/20/2023 09:47 | 1                  |             | UltraCarb5uODS 4.6 (mm) |
| ZZZZZ               |                              | 05/20/2023 10:09 | 1                  |             | UltraCarb5uODS 4.6 (mm) |
| ZZZZZ               |                              | 05/20/2023 10:32 | 1                  |             | UltraCarb5uODS 4.6 (mm) |
| ZZZZZ               |                              | 05/20/2023 10:55 | 1                  |             | UltraCarb5uODS 4.6 (mm) |
| ZZZZZ               |                              | 05/20/2023 11:18 | 1                  |             | UltraCarb5uODS 4.6 (mm) |
| CCV 280-613168/66   |                              | 05/20/2023 11:41 | 1                  |             | UltraCarb5uODS 4.6 (mm) |
| CCV 280-613168/67   |                              | 05/20/2023 12:04 | 1                  |             | UltraCarb5uODS 4.6 (mm) |



## HPLC/IC ANALYSIS RUN LOG

Lab Name: Eurofins Denver Job No.: 280-176674-1  
SDG No.: \_\_\_\_\_  
Instrument ID: CHHPLC\_X5 Start Date: 05/20/2023 18:42  
Analysis Batch Number: 613238 End Date: 05/21/2023 03:26

| LAB SAMPLE ID     | CLIENT SAMPLE ID      | DATE ANALYZED    | DILUTION<br>FACTOR | LAB FILE ID | COLUMN ID               |
|-------------------|-----------------------|------------------|--------------------|-------------|-------------------------|
| CCV 280-613238/7  |                       | 05/20/2023 18:42 | 1                  | 05200007.D  | Luna-phenylhex 4.6 (mm) |
| ZZZZZ             |                       | 05/20/2023 19:17 | 1                  |             | Luna-phenylhex 4.6 (mm) |
| ZZZZZ             |                       | 05/20/2023 19:52 | 1                  |             | Luna-phenylhex 4.6 (mm) |
| ZZZZZ             |                       | 05/20/2023 20:27 | 1                  |             | Luna-phenylhex 4.6 (mm) |
| ZZZZZ             |                       | 05/20/2023 21:02 | 1                  |             | Luna-phenylhex 4.6 (mm) |
| ZZZZZ             |                       | 05/20/2023 21:36 | 10000              |             | Luna-phenylhex 4.6 (mm) |
| ZZZZZ             |                       | 05/20/2023 22:11 | 1                  |             | Luna-phenylhex 4.6 (mm) |
| ZZZZZ             |                       | 05/20/2023 22:46 | 5                  |             | Luna-phenylhex 4.6 (mm) |
| ZZZZZ             |                       | 05/20/2023 23:21 | 1                  |             | Luna-phenylhex 4.6 (mm) |
| 280-176674-7      | FBQmw-173-230401-GW-R | 05/21/2023 00:31 | 1                  | 05200019.D  | Luna-phenylhex 4.6 (mm) |
| CCV 280-613238/20 |                       | 05/21/2023 01:06 | 1                  | 05200020.D  | Luna-phenylhex 4.6 (mm) |
| ZZZZZ             |                       | 05/21/2023 01:41 | 1                  |             | Luna-phenylhex 4.6 (mm) |
| ZZZZZ             |                       | 05/21/2023 02:16 | 5                  |             | Luna-phenylhex 4.6 (mm) |
| ZZZZZ             |                       | 05/21/2023 02:51 | 1                  |             | Luna-phenylhex 4.6 (mm) |
| CCV 280-613238/24 |                       | 05/21/2023 03:26 | 1                  |             | Luna-phenylhex 4.6 (mm) |

## HPLC/IC ANALYSIS RUN LOG

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Instrument ID: CHHPLC\_X5 Start Date: 05/24/2023 19:13Analysis Batch Number: 613683 End Date: 05/25/2023 00:28

| LAB SAMPLE ID     | CLIENT SAMPLE ID    | DATE ANALYZED    | DILUTION<br>FACTOR | LAB FILE ID | COLUMN ID               |
|-------------------|---------------------|------------------|--------------------|-------------|-------------------------|
| CCV 280-613683/8  |                     | 05/24/2023 19:13 | 1                  | 05240008.D  | Luna-phenylhex 4.6 (mm) |
| 280-176674-4      | LL3mw-241-230401-GW | 05/24/2023 19:48 | 1                  | 05240010.D  | Luna-phenylhex 4.6 (mm) |
| ZZZZZ             |                     | 05/24/2023 20:23 | 1                  |             | Luna-phenylhex 4.6 (mm) |
| ZZZZZ             |                     | 05/24/2023 20:58 | 1                  |             | Luna-phenylhex 4.6 (mm) |
| ZZZZZ             |                     | 05/24/2023 21:33 | 1                  |             | Luna-phenylhex 4.6 (mm) |
| ZZZZZ             |                     | 05/24/2023 22:08 | 1                  |             | Luna-phenylhex 4.6 (mm) |
| ZZZZZ             |                     | 05/24/2023 22:43 | 1                  |             | Luna-phenylhex 4.6 (mm) |
| ZZZZZ             |                     | 05/24/2023 23:18 | 1                  |             | Luna-phenylhex 4.6 (mm) |
| ZZZZZ             |                     | 05/24/2023 23:53 | 1                  |             | Luna-phenylhex 4.6 (mm) |
| CCV 280-613683/18 |                     | 05/25/2023 00:28 | 1                  | 05240018.D  | Luna-phenylhex 4.6 (mm) |

## HPLC/IC BATCH WORKSHEET

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Batch Number: 599427 Batch Start Date: 01/17/23 00:33 Batch Analyst: Zhang, JianBatch Method: 8330B Batch End Date: \_\_\_\_\_

| Lab Sample ID        | Client Sample ID | Method Chain | Basis | FinalAmount | 3,5-DNA LCS<br>00042 | 8330_ADDs 00035 | 8330DiaminLCS<br>00052 |  |  |
|----------------------|------------------|--------------|-------|-------------|----------------------|-----------------|------------------------|--|--|
| IC<br>280-599427/11  |                  | 8330B        |       | 1 mL        |                      | 125 uL          |                        |  |  |
| IC<br>280-599427/12  |                  | 8330B        |       | 1 mL        |                      | 50 uL           |                        |  |  |
| IC<br>280-599427/13  |                  | 8330B        |       | 1 mL        |                      | 35 uL           |                        |  |  |
| IC<br>280-599427/14  |                  | 8330B        |       | 1 mL        |                      | 20 uL           |                        |  |  |
| IC<br>280-599427/15  |                  | 8330B        |       | 1 mL        |                      | 12.5 uL         |                        |  |  |
| IC<br>280-599427/16  |                  | 8330B        |       | 1 mL        |                      | 5 uL            |                        |  |  |
| IC<br>280-599427/17  |                  | 8330B        |       | 1 mL        |                      | 2.5 uL          |                        |  |  |
| IC<br>280-599427/18  |                  | 8330B        |       | 1 mL        |                      | 1 uL            |                        |  |  |
| ICV<br>280-599427/19 |                  | 8330B        |       | 1 mL        | 50 uL                |                 | 50 uL                  |  |  |

| Batch Notes |        |
|-------------|--------|
| Methanol ID | 221054 |

| Basis | Basis Description |
|-------|-------------------|
|       |                   |

The pound sign (#) in the amount added field denotes that the reagent was used undiluted. All calculations are performed using the stated concentration for this reagent.

8330B

Page 1 of 1

## HPLC/IC BATCH WORKSHEET

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Batch Number: 601664 Batch Start Date: 02/08/23 15:38 Batch Analyst: Zhang, JianBatch Method: 8330B Batch End Date: \_\_\_\_\_

| Lab Sample ID        | Client Sample ID | Method Chain | Basis | FinalAmount | 3,5-DNA LCS<br>00043 | 8330 LCS 00121 | 8330IntermStk<br>00075 | 8330Surrogate<br>00138 |  |
|----------------------|------------------|--------------|-------|-------------|----------------------|----------------|------------------------|------------------------|--|
| IC<br>280-601664/11  |                  | 8330B        |       | 1 mL        |                      |                | 250 uL                 |                        |  |
| IC<br>280-601664/12  |                  | 8330B        |       | 1 mL        |                      |                | 100 uL                 |                        |  |
| IC<br>280-601664/13  |                  | 8330B        |       | 1 mL        |                      |                | 70 uL                  |                        |  |
| IC<br>280-601664/14  |                  | 8330B        |       | 1 mL        |                      |                | 40 uL                  |                        |  |
| IC<br>280-601664/15  |                  | 8330B        |       | 1 mL        |                      |                | 25 uL                  |                        |  |
| IC<br>280-601664/16  |                  | 8330B        |       | 1 mL        |                      |                | 10 uL                  |                        |  |
| IC<br>280-601664/17  |                  | 8330B        |       | 1 mL        |                      |                | 5 uL                   |                        |  |
| IC<br>280-601664/18  |                  | 8330B        |       | 1 mL        |                      |                | 2 uL                   |                        |  |
| IC<br>280-601664/19  |                  | 8330B        |       | 1 mL        |                      |                | 1 uL                   |                        |  |
| ICV<br>280-601664/20 |                  | 8330B        |       | 1 mL        | 50 uL                | 50 uL          |                        | 50 uL                  |  |

| Batch Notes |        |
|-------------|--------|
| Methanol ID | 221054 |

| Basis | Basis Description |
|-------|-------------------|
|       |                   |

The pound sign (#) in the amount added field denotes that the reagent was used undiluted. All calculations are performed using the stated concentration for this reagent.

## HPLC/IC BATCH WORKSHEET

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Batch Number: 603284 Batch Start Date: 02/24/23 15:00 Batch Analyst: Zhang, JianBatch Method: 8330B Batch End Date: \_\_\_\_\_

| Lab Sample ID        | Client Sample ID | Method Chain | Basis | FinalAmount | 8330 DMT 00013 | 8330_OP_DMT<br>00016 |  |  |  |
|----------------------|------------------|--------------|-------|-------------|----------------|----------------------|--|--|--|
| IC<br>280-603284/11  |                  | 8330B        |       | 1 mL        | 125 uL         |                      |  |  |  |
| IC<br>280-603284/12  |                  | 8330B        |       | 1 mL        | 50 uL          |                      |  |  |  |
| IC<br>280-603284/13  |                  | 8330B        |       | 1 mL        | 35 uL          |                      |  |  |  |
| IC<br>280-603284/14  |                  | 8330B        |       | 1 mL        | 20 uL          |                      |  |  |  |
| IC<br>280-603284/15  |                  | 8330B        |       | 1 mL        | 12.5 uL        |                      |  |  |  |
| IC<br>280-603284/16  |                  | 8330B        |       | 1 mL        | 5 uL           |                      |  |  |  |
| IC<br>280-603284/17  |                  | 8330B        |       | 1 mL        | 2.5 uL         |                      |  |  |  |
| IC<br>280-603284/18  |                  | 8330B        |       | 1 mL        | 1 uL           |                      |  |  |  |
| ICV<br>280-603284/19 |                  | 8330B        |       | 1 mL        |                | 50 uL                |  |  |  |

| Batch Notes |        |
|-------------|--------|
| Methanol ID | 221054 |

| Basis | Basis Description |
|-------|-------------------|
|       |                   |

The pound sign (#) in the amount added field denotes that the reagent was used undiluted. All calculations are performed using the stated concentration for this reagent.

## HPLC/IC BATCH WORKSHEET

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Batch Number: 610603 Batch Start Date: 04/28/23 18:24 Batch Analyst: Zhang, JianBatch Method: 8330B Batch End Date: \_\_\_\_\_

| Lab Sample ID        | Client Sample ID | Method Chain | Basis | FinalAmount | 8330 DMT 00013 | 8330 LCS 00126 | 8330_ADDs 00036 | 8330_OP_DMT<br>00017 | 8330DiaminLCS<br>00056 |
|----------------------|------------------|--------------|-------|-------------|----------------|----------------|-----------------|----------------------|------------------------|
| IC<br>280-610603/10  |                  | 8330B        |       | 1 mL        | 125 uL         |                |                 |                      |                        |
| IC<br>280-610603/11  |                  | 8330B        |       | 1 mL        | 50 uL          |                |                 |                      |                        |
| IC<br>280-610603/12  |                  | 8330B        |       | 1 mL        | 35 uL          |                |                 |                      |                        |
| IC<br>280-610603/13  |                  | 8330B        |       | 1 mL        | 20 uL          |                |                 |                      |                        |
| IC<br>280-610603/14  |                  | 8330B        |       | 1 mL        | 12.5 uL        |                |                 |                      |                        |
| IC<br>280-610603/15  |                  | 8330B        |       | 1 mL        | 5 uL           |                |                 |                      |                        |
| IC<br>280-610603/16  |                  | 8330B        |       | 1 mL        | 2.5 uL         |                |                 |                      |                        |
| IC<br>280-610603/17  |                  | 8330B        |       | 1 mL        | 1 uL           |                |                 |                      |                        |
| IC<br>280-610603/18  |                  | 8330B        |       | 1 mL        | 0.5 uL         |                |                 |                      |                        |
| ICV<br>280-610603/19 |                  | 8330B        |       | 1 mL        |                | 50 uL          |                 | 50 uL                |                        |
| IC<br>280-610603/20  |                  | 8330B        |       | 1 mL        |                |                | 125 uL          |                      |                        |
| IC<br>280-610603/21  |                  | 8330B        |       | 1 mL        |                |                | 50 uL           |                      |                        |
| IC<br>280-610603/22  |                  | 8330B        |       | 1 mL        |                |                | 35 uL           |                      |                        |
| IC<br>280-610603/23  |                  | 8330B        |       | 1 mL        |                |                | 20 uL           |                      |                        |
| IC<br>280-610603/24  |                  | 8330B        |       | 1 mL        |                |                | 12.5 uL         |                      |                        |
| IC<br>280-610603/25  |                  | 8330B        |       | 1 mL        |                |                | 5 uL            |                      |                        |
| IC<br>280-610603/26  |                  | 8330B        |       | 1 mL        |                |                | 2 uL            |                      |                        |
| IC<br>280-610603/27  |                  | 8330B        |       | 1 mL        |                |                | 1 uL            |                      |                        |
| ICV<br>280-610603/28 |                  | 8330B        |       | 1 mL        |                |                |                 |                      | 50 uL                  |

| Lab Sample ID | Client Sample ID | Method Chain | Basis | 8330IntermStk<br>00076 | 8330Surrogate<br>00143 |  |  |  |  |
|---------------|------------------|--------------|-------|------------------------|------------------------|--|--|--|--|
|---------------|------------------|--------------|-------|------------------------|------------------------|--|--|--|--|

The pound sign (#) in the amount added field denotes that the reagent was used undiluted. All calculations are performed using the stated concentration for this reagent.

## HPLC/IC BATCH WORKSHEET

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Batch Number: 610603 Batch Start Date: 04/28/23 18:24 Batch Analyst: Zhang, JianBatch Method: 8330B Batch End Date: \_\_\_\_\_

| Lab Sample ID        | Client Sample ID | Method Chain | Basis | 8330IntermStk<br>00076 | 8330Surrogate<br>00143 |  |  |  |  |
|----------------------|------------------|--------------|-------|------------------------|------------------------|--|--|--|--|
| IC<br>280-610603/10  |                  | 8330B        |       | 250 uL                 |                        |  |  |  |  |
| IC<br>280-610603/11  |                  | 8330B        |       | 100 uL                 |                        |  |  |  |  |
| IC<br>280-610603/12  |                  | 8330B        |       | 70 uL                  |                        |  |  |  |  |
| IC<br>280-610603/13  |                  | 8330B        |       | 40 uL                  |                        |  |  |  |  |
| IC<br>280-610603/14  |                  | 8330B        |       | 25 uL                  |                        |  |  |  |  |
| IC<br>280-610603/15  |                  | 8330B        |       | 10 uL                  |                        |  |  |  |  |
| IC<br>280-610603/16  |                  | 8330B        |       | 5 uL                   |                        |  |  |  |  |
| IC<br>280-610603/17  |                  | 8330B        |       | 2 uL                   |                        |  |  |  |  |
| IC<br>280-610603/18  |                  | 8330B        |       | 1 uL                   |                        |  |  |  |  |
| ICV<br>280-610603/19 |                  | 8330B        |       |                        | 50 uL                  |  |  |  |  |
| IC<br>280-610603/20  |                  | 8330B        |       |                        |                        |  |  |  |  |
| IC<br>280-610603/21  |                  | 8330B        |       |                        |                        |  |  |  |  |
| IC<br>280-610603/22  |                  | 8330B        |       |                        |                        |  |  |  |  |
| IC<br>280-610603/23  |                  | 8330B        |       |                        |                        |  |  |  |  |
| IC<br>280-610603/24  |                  | 8330B        |       |                        |                        |  |  |  |  |
| IC<br>280-610603/25  |                  | 8330B        |       |                        |                        |  |  |  |  |
| IC<br>280-610603/26  |                  | 8330B        |       |                        |                        |  |  |  |  |
| IC<br>280-610603/27  |                  | 8330B        |       |                        |                        |  |  |  |  |
| ICV<br>280-610603/28 |                  | 8330B        |       |                        |                        |  |  |  |  |

| Batch Notes |        |
|-------------|--------|
| Methanol ID | 221056 |

The pound sign (#) in the amount added field denotes that the reagent was used undiluted. All calculations are performed using the stated concentration for this reagent.

8330B

Page 2 of 3

## HPLC/IC BATCH WORKSHEET

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Batch Number: 610603 Batch Start Date: 04/28/23 18:24 Batch Analyst: Zhang, JianBatch Method: 8330B Batch End Date: \_\_\_\_\_

| Basis | Basis Description |
|-------|-------------------|
|       |                   |

The pound sign (#) in the amount added field denotes that the reagent was used undiluted. All calculations are performed using the stated concentration for this reagent.



## HPLC/IC BATCH WORKSHEET

Lab Name: Eurofins DenverJob No.: 280-176674-1

SDG No.: \_\_\_\_\_

Batch Number: 613095Batch Start Date: 05/19/23 13:39Batch Analyst: Hermanova, EvaBatch Method: 3535Batch End Date: 05/19/23 16:39

| Lab Sample ID          | Client Sample ID          | Method Chain | Basis | GrossWeight | TareWeight | InitialAmount | FinalAmount | 8330 LCS 00125 | 8330_OP_DMT<br>00020 |
|------------------------|---------------------------|--------------|-------|-------------|------------|---------------|-------------|----------------|----------------------|
| MB 280-613095/1        |                           | 3535, 8330B  |       |             |            | 500 mL        | 5 mL        |                |                      |
| LCS<br>280-613095/2    |                           | 3535, 8330B  |       |             |            | 500 mL        | 5 mL        | 0.1 mL         | 0.1 mL               |
| LCSD<br>280-613095/3   |                           | 3535, 8330B  |       |             |            | 500 mL        | 5 mL        | 0.1 mL         | 0.1 mL               |
| LCS<br>280-613095/4    |                           | 3535, 8330B  |       |             |            | 500 mL        | 5 mL        |                |                      |
| LCSD<br>280-613095/5   |                           | 3535, 8330B  |       |             |            | 500 mL        | 5 mL        |                |                      |
| 280-176674-B-2         | FWGmw-015-230401<br>-GW   | 3535, 8330B  | T     | 772.1 g     | 290.4 g    | 481.7 mL      | 5 mL        |                |                      |
| 280-176674-A-4         | LL3mw-241-230401<br>-GW   | 3535, 8330B  | T     | 774.5 g     | 286.1 g    | 488.4 mL      | 5 mL        |                |                      |
| 280-176674-C-7         | FBQmw-173-230401<br>-GW-R | 3535, 8330B  | T     | 769.7 g     | 286.5 g    | 483.2 mL      | 5 mL        |                |                      |
| 280-176674-A-10        | FBQmw-175-230401<br>-GW-R | 3535, 8330B  | T     | 783.4 g     | 285.9 g    | 497.5 mL      | 5 mL        |                |                      |
| 280-176674-B-10<br>MS  | FBQmw-175-230401<br>-GW-R | 3535, 8330B  | T     | 784.4 g     | 285.0 g    | 499.4 mL      | 5 mL        | 0.1 mL         | 0.1 mL               |
| 280-176674-B-10<br>MSD | FBQmw-175-230401<br>-GW-R | 3535, 8330B  | T     | 743.4 g     | 286.0 g    | 457.4 mL      | 5 mL        | 0.1 mL         | 0.1 mL               |
| 280-176674-C-10<br>MS  | FBQmw-175-230401<br>-GW-R | 3535, 8330B  | T     | 784.2 g     | 285.5 g    | 498.7 mL      | 5 mL        |                |                      |
| 280-176674-D-10<br>MSD | FBQmw-175-230401<br>-GW-R | 3535, 8330B  | T     | 785.1 g     | 286.3 g    | 498.8 mL      | 5 mL        |                |                      |
| 280-176674-A-11        | FBQmw-175-230402<br>-GW-R | 3535, 8330B  | T     | 783.9 g     | 286.9 g    | 497 mL        | 5 mL        |                |                      |

| Lab Sample ID        | Client Sample ID        | Method Chain | Basis | 8330DiaminLCS<br>00057 | 8330Surrogate<br>00144 | AnalysisComment |  |  |  |
|----------------------|-------------------------|--------------|-------|------------------------|------------------------|-----------------|--|--|--|
| MB 280-613095/1      |                         | 3535, 8330B  |       |                        | 0.1 mL                 |                 |  |  |  |
| LCS<br>280-613095/2  |                         | 3535, 8330B  |       |                        | 0.1 mL                 | DMT             |  |  |  |
| LCSD<br>280-613095/3 |                         | 3535, 8330B  |       |                        | 0.1 mL                 | DMT             |  |  |  |
| LCS<br>280-613095/4  |                         | 3535, 8330B  |       | 0.1 mL                 | 0.1 mL                 | Diamino         |  |  |  |
| LCSD<br>280-613095/5 |                         | 3535, 8330B  |       | 0.1 mL                 | 0.1 mL                 | Diamino         |  |  |  |
| 280-176674-B-2       | FWGmw-015-230401<br>-GW | 3535, 8330B  | T     |                        | 0.1 mL                 |                 |  |  |  |

The pound sign (#) in the amount added field denotes that the reagent was used undiluted. All calculations are performed using the stated concentration for this reagent.

## HPLC/IC BATCH WORKSHEET

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Batch Number: 613095 Batch Start Date: 05/19/23 13:39 Batch Analyst: Hermanova, EvaBatch Method: 3535 Batch End Date: 05/19/23 16:39

| Lab Sample ID          | Client Sample ID          | Method Chain | Basis | 8330DiaminLCS<br>00057 | 8330Surrogate<br>00144 | AnalysisComment |  |  |  |
|------------------------|---------------------------|--------------|-------|------------------------|------------------------|-----------------|--|--|--|
| 280-176674-A-4         | LL3mw-241-230401<br>-GW   | 3535, 8330B  | T     |                        | 0.1 mL                 |                 |  |  |  |
| 280-176674-C-7         | FBQmw-173-230401<br>-GW-R | 3535, 8330B  | T     |                        | 0.1 mL                 |                 |  |  |  |
| 280-176674-A-10        | FBQmw-175-230401<br>-GW-R | 3535, 8330B  | T     |                        | 0.1 mL                 |                 |  |  |  |
| 280-176674-B-10<br>MS  | FBQmw-175-230401<br>-GW-R | 3535, 8330B  | T     |                        | 0.1 mL                 | DMT             |  |  |  |
| 280-176674-B-10<br>MSD | FBQmw-175-230401<br>-GW-R | 3535, 8330B  | T     |                        | 0.1 mL                 | DMT             |  |  |  |
| 280-176674-C-10<br>MS  | FBQmw-175-230401<br>-GW-R | 3535, 8330B  | T     | 0.1 mL                 | 0.1 mL                 | Diamino         |  |  |  |
| 280-176674-D-10<br>MSD | FBQmw-175-230401<br>-GW-R | 3535, 8330B  | T     | 0.1 mL                 | 0.1 mL                 | Diamino         |  |  |  |
| 280-176674-A-11        | FBQmw-175-230402<br>-GW-R | 3535, 8330B  | T     |                        | 0.1 mL                 |                 |  |  |  |

The pound sign (#) in the amount added field denotes that the reagent was used undiluted. All calculations are performed using the stated concentration for this reagent.

## HPLC/IC BATCH WORKSHEET

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Batch Number: 613095 Batch Start Date: 05/19/23 13:39 Batch Analyst: Hermanova, EvaBatch Method: 3535 Batch End Date: 05/19/23 16:39

| Batch Notes                        |                            |
|------------------------------------|----------------------------|
| First Start time                   | 05/19/2023 13:59           |
| First End time                     | 05/19/2023 16:17           |
| SPE Cartridge Type                 | Sep-Pak Porapak Rdx        |
| SPE Cartridge Lot ID               | 005332286A                 |
| Balance ID                         | 24350888                   |
| Manifold ID                        | Manifold: A, B             |
| QC Bottle Lot ID                   | 0420301F, 02024011         |
| Pipette/Syringe/Dispenser ID       | Dobby/DOD/Pugsley          |
| Solvent Name                       | CaCl2                      |
| Solvent Lot #                      | CaCl2_Sol_00085            |
| Rinse Solvent Name                 | Acetonitrile               |
| Rinse Solvent Lot                  | Acetonitrile_00078         |
| Acid Name                          | 0.1% AAinACN               |
| Acid ID                            | 0.1% AAinACN_00207         |
| Analyst ID - Spike Analyst         | EH (Trainee), DL (Trainer) |
| Analyst ID - Spike Witness Analyst | Reviewer: NC               |
| Batch Comment                      | DV-OP-0017                 |

| Basis | Basis Description |
|-------|-------------------|
| T     | Total/NA          |

The pound sign (#) in the amount added field denotes that the reagent was used undiluted. All calculations are performed using the stated concentration for this reagent.

# GENERAL CHEMISTRY

COVER PAGE  
GENERAL CHEMISTRY

Lab Name: Eurofins Denver Job Number: 280-176674-1

SDG No.:

Project: RVAAP FWGW

| Client Sample ID      | Lab Sample ID |
|-----------------------|---------------|
| LL12mw-187-230401-GW  | 280-176674-3  |
| LL12mw-185-230401-GW  | 280-176674-5  |
| LL12mw-185-230402-GW  | 280-176674-6  |
| FBQmw-173-230401-GW-R | 280-176674-7  |
| FBQmw-173-230402-GW-R | 280-176674-8  |
| LL12mw-245-230401-GW  | 280-176674-9  |

Comments:

1B-IN  
INORGANIC ANALYSIS DATA SHEET  
GENERAL CHEMISTRY

Client Sample ID: LL12mw-187-230401-GW

Lab Sample ID: 280-176674-3

Lab Name: Eurofins Denver

Job No.: 280-176674-1

SDG ID.:

Matrix: Water

Date Sampled: 05/17/2023 13:40

Reporting Basis: WET

Date Received: 05/18/2023 09:45

| Analyte      | Result | LOQ  | LOD | DL  | Units | C | Q   | DIL  | Method |
|--------------|--------|------|-----|-----|-------|---|-----|------|--------|
| Ammonia as N | 730    | 20   | 10  | 5.8 | mg/L  |   |     | 200  | 350.1  |
| Nitrate as N | 1500   | 1000 | 400 | 180 | mg/L  |   | H D | 2000 | 9056   |

1B-IN  
INORGANIC ANALYSIS DATA SHEET  
GENERAL CHEMISTRY

Client Sample ID: LL12mw-185-230401-GW

Lab Sample ID: 280-176674-5

Lab Name: Eurofins Denver

Job No.: 280-176674-1

SDG ID.:

Matrix: Water

Date Sampled: 05/17/2023 09:45

Reporting Basis: WET

Date Received: 05/18/2023 09:45

| Analyte      | Result | LOQ  | LOD   | DL    | Units | C | Q   | DIL | Method |
|--------------|--------|------|-------|-------|-------|---|-----|-----|--------|
| Ammonia as N | 0.050  | 0.10 | 0.050 | 0.029 | mg/L  | U |     | 1   | 350.1  |
| Nitrate as N | 67     | 10   | 4.0   | 1.8   | mg/L  |   | D M | 20  | 9056   |

1B-IN  
INORGANIC ANALYSIS DATA SHEET  
GENERAL CHEMISTRY

Client Sample ID: LL12mw-185-230402-GW

Lab Sample ID: 280-176674-6

Lab Name: Eurofins Denver

Job No.: 280-176674-1

SDG ID.:

Matrix: Water

Date Sampled: 05/17/2023 09:45

Reporting Basis: WET

Date Received: 05/18/2023 09:45

| Analyte      | Result | LOQ  | LOD   | DL    | Units | C | Q   | DIL | Method |
|--------------|--------|------|-------|-------|-------|---|-----|-----|--------|
| Ammonia as N | 0.050  | 0.10 | 0.050 | 0.029 | mg/L  | U |     | 1   | 350.1  |
| Nitrate as N | 74     | 10   | 4.0   | 1.8   | mg/L  |   | D M | 20  | 9056   |



1B-IN  
INORGANIC ANALYSIS DATA SHEET  
GENERAL CHEMISTRY

Client Sample ID: FBQmw-173-230401-GW-R

Lab Sample ID: 280-176674-7

Lab Name: Eurofins Denver

Job No.: 280-176674-1

SDG ID.:

Matrix: Water

Date Sampled: 05/17/2023 11:47

Reporting Basis: WET

Date Received: 05/18/2023 09:45

| Analyte                        | Result | LOQ  | LOD  | DL    | Units | C | Q | DIL | Method   |
|--------------------------------|--------|------|------|-------|-------|---|---|-----|----------|
| Sulfide                        | 2.0    | 4.0  | 2.0  | 1.6   | mg/L  | U |   | 1   | 9034     |
| Nitrate as N                   | 0.20   | 0.50 | 0.20 | 0.090 | mg/L  | J | M | 1   | 9056     |
| Nitrite as N                   | 0.10   | 0.50 | 0.10 | 0.049 | mg/L  | U |   | 1   | 9056     |
| Sulfate                        | 42     | 5.0  | 2.5  | 1.0   | mg/L  |   | M | 1   | 9056A    |
| Total Organic Carbon - Average | 0.85   | 1.0  | 0.80 | 0.35  | mg/L  | J |   | 1   | 9060     |
| Alkalinity                     | 18     | 10   | 6.4  | 3.1   | mg/L  |   |   | 1   | SM 2320B |

1B-IN  
INORGANIC ANALYSIS DATA SHEET  
GENERAL CHEMISTRY

Client Sample ID: FBQmw-173-230402-GW-R

Lab Sample ID: 280-176674-8

Lab Name: Eurofins Denver

Job No.: 280-176674-1

SDG ID.:

Matrix: Water

Date Sampled: 05/17/2023 11:47

Reporting Basis: WET

Date Received: 05/18/2023 09:45

| Analyte      | Result | LOQ  | LOD  | DL    | Units | C | Q | DIL | Method   |
|--------------|--------|------|------|-------|-------|---|---|-----|----------|
| Sulfide      | 2.0    | 4.0  | 2.0  | 1.6   | mg/L  | U |   | 1   | 9034     |
| Nitrate as N | 0.21   | 0.50 | 0.20 | 0.090 | mg/L  | J | M | 1   | 9056     |
| Nitrite as N | 0.10   | 0.50 | 0.10 | 0.049 | mg/L  | U |   | 1   | 9056     |
| Sulfate      | 45     | 5.0  | 2.5  | 1.0   | mg/L  |   |   | 1   | 9056A    |
| Alkalinity   | 18     | 10   | 6.4  | 3.1   | mg/L  |   |   | 1   | SM 2320B |

1B-IN  
INORGANIC ANALYSIS DATA SHEET  
GENERAL CHEMISTRY

Client Sample ID: LL12mw-245-230401-GW

Lab Sample ID: 280-176674-9

Lab Name: Eurofins Denver

Job No.: 280-176674-1

SDG ID.:

Matrix: Water

Date Sampled: 05/17/2023 14:26

Reporting Basis: WET

Date Received: 05/18/2023 09:45

| Analyte      | Result | LOQ  | LOD   | DL    | Units | C | Q | DIL | Method |
|--------------|--------|------|-------|-------|-------|---|---|-----|--------|
| Ammonia as N | 0.050  | 0.10 | 0.050 | 0.029 | mg/L  | U |   | 1   | 350.1  |
| Nitrate as N | 0.20   | 0.50 | 0.20  | 0.090 | mg/L  | U | M | 1   | 9056   |

2-IN  
CALIBRATION QUALITY CONTROL  
GENERAL CHEMISTRY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
 SDG No.: \_\_\_\_\_  
 Analyst: MMP Batch Start Date: 05/23/2023  
 Reporting Units: mg/L Analytical Batch No.: 613561

| Sample Number | QC Type | Time  | Analyte      | Result | Spike Amount | (%) Recovery | Limits | Qual | Reagent         |
|---------------|---------|-------|--------------|--------|--------------|--------------|--------|------|-----------------|
| 14            | ICV     | 08:51 | Ammonia as N | 2.49   | 2.51         | 99           | 90-110 |      | 350.1 ICV_00564 |
| 15            | ICVL    | 08:53 | Ammonia as N | 0.491  | 0.501        | 98           | 90-110 |      | 350.1 ICV_00564 |
| 16            | ICB     | 08:56 | Ammonia as N | 0.050  |              |              |        | U    |                 |
| 83            | CCV     | 11:54 | Ammonia as N | 2.58   | 2.50         | 103          | 90-110 |      | 350.1 cal_00583 |
| 84            | CCVL    | 11:57 | Ammonia as N | 0.520  | 0.500        | 104          | 90-110 |      | 350.1 cal_00583 |
| 85            | CCB     | 12:00 | Ammonia as N | 0.050  |              |              |        | U    |                 |
| 96            | CCV     | 12:29 | Ammonia as N | 2.47   | 2.50         | 99           | 90-110 |      | 350.1 cal_00583 |
| 97            | CCVL    | 12:32 | Ammonia as N | 0.512  | 0.500        | 102          | 90-110 |      | 350.1 cal_00583 |
| 98            | CCB     | 12:34 | Ammonia as N | 0.050  |              |              |        | U    |                 |
| 111           | CCV     | 13:09 | Ammonia as N | 2.55   | 2.50         | 102          | 90-110 |      | 350.1 cal_00583 |
| 112           | CCVL    | 13:12 | Ammonia as N | 0.494  | 0.500        | 99           | 90-110 |      | 350.1 cal_00583 |
| 113           | CCB     | 13:14 | Ammonia as N | 0.050  |              |              |        | U    |                 |
| 124           | CCV     | 13:46 | Ammonia as N | 2.74   | 2.50         | 109          | 90-110 |      | 350.1 cal_00583 |
| 125           | CCVL    | 13:49 | Ammonia as N | 0.551  | 0.500        | 110          | 90-110 |      | 350.1 cal_00583 |
| 195           | CCV     | 16:56 | Ammonia as N | 2.56   | 2.50         | 102          | 90-110 |      | 350.1 cal_00583 |
| 196           | CCVL    | 16:58 | Ammonia as N | 0.479  | 0.500        | 96           | 90-110 |      | 350.1 cal_00583 |
| 197           | CCB     | 17:01 | Ammonia as N | 0.050  |              |              |        | U    |                 |
| 209           | CCV     | 17:33 | Ammonia as N | 2.54   | 2.50         | 102          | 90-110 |      | 350.1 cal_00583 |
| 210           | CCVL    | 17:36 | Ammonia as N | 0.477  | 0.500        | 95           | 90-110 |      | 350.1 cal_00583 |
| 211           | CCB     | 17:38 | Ammonia as N | 0.050  |              |              |        | U    |                 |

Note! Calculations are performed before rounding to avoid round-off errors in calculated results.

FORM II-IN

2-IN  
CALIBRATION QUALITY CONTROL  
GENERAL CHEMISTRY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
 SDG No.: \_\_\_\_\_  
 Analyst: MEC Batch Start Date: 05/18/2023  
 Reporting Units: mg/L Analytical Batch No.: 612961

| Sample Number | QC Type | Time  | Analyte      | Result | Spike Amount | (%) Recovery | Limits | Qual | Reagent        |
|---------------|---------|-------|--------------|--------|--------------|--------------|--------|------|----------------|
| 7             | ICV     | 13:25 | Nitrate as N | 1.95   | 2.00         | 98           | 90-110 |      | IC ICV 5_00405 |
|               |         |       | Nitrite as N | 1.99   | 2.00         | 100          | 90-110 |      | IC ICV 5_00405 |
| 8             | ICB     | 13:40 | Nitrate as N | 0.20   |              |              |        | U    |                |
|               |         |       | Nitrite as N | 0.10   |              |              |        | U    |                |
| 24            | CCV     | 17:41 | Nitrate as N | 2.43   | 2.50         | 97           | 90-110 |      | IC LCS_01954   |
|               |         |       | Nitrite as N | 2.54   | 2.50         | 101          | 90-110 |      | IC LCS_01954   |
| 25            | CCB     | 17:56 | Nitrate as N | 0.20   |              |              |        | U    |                |
|               |         |       | Nitrite as N | 0.10   |              |              |        | U    |                |
| 36            | CCV     | 20:42 | Nitrate as N | 2.43   | 2.50         | 97           | 90-110 |      | IC LCS_01954   |
|               |         |       | Nitrite as N | 2.53   | 2.50         | 101          | 90-110 |      | IC LCS_01954   |
| 37            | CCB     | 20:57 | Nitrate as N | 0.20   |              |              |        | U    |                |
|               |         |       | Nitrite as N | 0.10   |              |              |        | U    |                |
| 48            | CCV     | 23:41 | Nitrate as N | 2.52   | 2.50         | 101          | 90-110 |      | IC LCS_01954   |
| 49            | CCB     | 23:56 | Nitrate as N | 0.20   |              |              |        | U    |                |
| 60            | CCV     | 02:42 | Nitrate as N | 2.68   | 2.50         | 107          | 90-110 |      | IC LCS_01954   |
| 61            | CCB     | 02:57 | Nitrate as N | 0.20   |              |              |        | U    |                |
| 72            | CCV     | 05:42 | Nitrate as N | 2.67   | 2.50         | 107          | 90-110 |      | IC LCS_01954   |
| 73            | CCB     | 05:57 | Nitrate as N | 0.20   |              |              |        | U    |                |
| 84            | CCV     | 08:42 | Nitrate as N | 2.55   | 2.50         | 102          | 90-110 |      | IC LCS_01954   |
| 85            | CCB     | 08:57 | Nitrate as N | 0.20   |              |              |        | U    |                |

Note! Calculations are performed before rounding to avoid round-off errors in calculated results.

2-IN  
CALIBRATION QUALITY CONTROL  
GENERAL CHEMISTRY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
SDG No.: \_\_\_\_\_  
Analyst: EJS Batch Start Date: 05/22/2023  
Reporting Units: mg/L Analytical Batch No.: 613351

| Sample Number | QC Type | Time  | Analyte      | Result | Spike Amount | (%) Recovery | Limits | Qual | Reagent      |
|---------------|---------|-------|--------------|--------|--------------|--------------|--------|------|--------------|
| 1             | CCV     | 16:00 | Nitrate as N | 2.43   | 2.50         | 97           | 90-110 |      | IC LCS_01954 |
| 2             | CCB     | 16:14 | Nitrate as N | 0.20   |              |              |        | U    |              |
| 17            | CCV     | 19:59 | Nitrate as N | 2.42   | 2.50         | 97           | 90-110 |      | IC LCS_01954 |
| 18            | CCB     | 20:15 | Nitrate as N | 0.20   |              |              |        | U    |              |

Note! Calculations are performed before rounding to avoid round-off errors in calculated results.

FORM II-IN

2-IN  
CALIBRATION QUALITY CONTROL  
GENERAL CHEMISTRY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
SDG No.: \_\_\_\_\_  
Analyst: MEC Batch Start Date: 05/18/2023  
Reporting Units: mg/L Analytical Batch No.: 612962

| Sample Number | QC Type | Time  | Analyte | Result | Spike Amount | (%) Recovery | Limits | Qual | Reagent          |
|---------------|---------|-------|---------|--------|--------------|--------------|--------|------|------------------|
| 7             | ICV     | 13:25 | Sulfate | 81.8   | 80.0         | 102          | 90-110 |      | IC SO4 ICV_00024 |
| 8             | ICB     | 13:40 | Sulfate | 2.5    |              |              |        | U    |                  |
| 24            | CCV     | 17:41 | Sulfate | 48.9   | 50.0         | 98           | 90-110 |      | IC LCS_01954     |
| 25            | CCB     | 17:56 | Sulfate | 2.5    |              |              |        | U    |                  |
| 36            | CCV     | 20:42 | Sulfate | 49.4   | 50.0         | 99           | 90-110 |      | IC LCS_01954     |
| 37            | CCB     | 20:57 | Sulfate | 2.5    |              |              |        | U    |                  |

Note! Calculations are performed before rounding to avoid round-off errors in calculated results.

FORM II-IN

2-IN  
CALIBRATION QUALITY CONTROL  
GENERAL CHEMISTRY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
SDG No.: \_\_\_\_\_  
Analyst: ABW Batch Start Date: 05/19/2023  
Reporting Units: mg/L Analytical Batch No.: 613268

| Sample Number | QC Type | Time  | Analyte                        | Result | Spike Amount | (%) Recovery | Limits | Qual | Reagent           |
|---------------|---------|-------|--------------------------------|--------|--------------|--------------|--------|------|-------------------|
| 1             | ICV     | 18:03 | Total Organic Carbon - Average | 19.7   | 20.0         | 99           | 90-110 |      | TOC ICV Std_00051 |
| 2             | ICB     | 18:17 | Total Organic Carbon - Average | 0.80   |              |              |        | U    |                   |
| 15            | CCV     | 21:30 | Total Organic Carbon - Average | 24.5   | 25.0         | 98           | 90-110 |      | TOC LCS Std_00057 |
| 16            | CCB     | 21:45 | Total Organic Carbon - Average | 0.80   |              |              |        | U    |                   |

Note! Calculations are performed before rounding to avoid round-off errors in calculated results.

FORM II-IN



2-IN  
CALIBRATION QUALITY CONTROL  
GENERAL CHEMISTRY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
SDG No.: \_\_\_\_\_  
Analyst: KEG Batch Start Date: 05/19/2023  
Reporting Units: mg/L Analytical Batch No.: 613263

| Sample Number | QC Type | Time  | Analyte    | Result | Spike Amount | (%) Recovery | Limits | Qual | Reagent             |
|---------------|---------|-------|------------|--------|--------------|--------------|--------|------|---------------------|
| 29            | CCV     | 15:42 | Alkalinity | 182    | 200          | 91           | 90-110 |      | Alk daily lcs 01088 |
| 30            | CCB     | 15:49 | Alkalinity | 3.43   |              |              |        | J    |                     |
| 44            | CCV     | 17:30 | Alkalinity | 182    | 200          | 91           | 90-110 |      | Alk daily lcs 01088 |
| 45            | CCB     | 17:36 | Alkalinity | 3.58   |              |              |        | J    |                     |

Note! Calculations are performed before rounding to avoid round-off errors in calculated results.

FORM II-IN

3-IN  
METHOD BLANK  
GENERAL CHEMISTRY

Lab Name: Eurofins Denver

Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

| Method                                                                            | Lab Sample ID     | Analyte                        | Result | Qual | Units | LOQ  | Dil |
|-----------------------------------------------------------------------------------|-------------------|--------------------------------|--------|------|-------|------|-----|
| Batch ID: 613561 Date: 05/23/2023 12:13                                           |                   |                                |        |      |       |      |     |
| 350.1                                                                             | MB 280-613561/90  | Ammonia as N                   | 0.050  | U    | mg/L  | 0.10 | 1   |
| Batch ID: 613208 Date: 05/19/2023 17:27 Prep Batch: 613197 Date: 05/19/2023 15:58 |                   |                                |        |      |       |      |     |
| 9034                                                                              | MB 280-613197/2-A | Sulfide                        | 2.0    | U    | mg/L  | 4.0  | 1   |
| Batch ID: 612961 Date: 05/18/2023 14:55                                           |                   |                                |        |      |       |      |     |
| 9056                                                                              | MB 280-612961/13  | Nitrate as N                   | 0.20   | U    | mg/L  | 0.50 | 1   |
| 9056                                                                              | MB 280-612961/13  | Nitrite as N                   | 0.10   | U    | mg/L  | 0.50 | 1   |
| Batch ID: 612961 Date: 05/19/2023 01:26                                           |                   |                                |        |      |       |      |     |
| 9056                                                                              | MB 280-612961/55  | Nitrate as N                   | 0.20   | U M  | mg/L  | 0.50 | 1   |
| Batch ID: 613351 Date: 05/22/2023 17:14                                           |                   |                                |        |      |       |      |     |
| 9056                                                                              | MB 280-613351/6   | Nitrate as N                   | 0.20   | U    | mg/L  | 0.50 | 1   |
| Batch ID: 612962 Date: 05/18/2023 14:55                                           |                   |                                |        |      |       |      |     |
| 9056A                                                                             | MB 280-612962/13  | Sulfate                        | 2.5    | U M  | mg/L  | 5.0  | 1   |
| Batch ID: 613268 Date: 05/19/2023 18:46                                           |                   |                                |        |      |       |      |     |
| 9060                                                                              | MB 280-613268/4   | Total Organic Carbon - Average | 0.80   | U    | mg/L  | 1.0  | 1   |
| Batch ID: 613263 Date: 05/19/2023 16:11                                           |                   |                                |        |      |       |      |     |
| SM 2320B                                                                          | MB 280-613263/33  | Alkalinity                     | 3.58   | J    | mg/L  | 10   | 1   |

5-IN  
MATRIX SPIKE SAMPLE RECOVERY  
GENERAL CHEMISTRY

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Matrix: Water

| Method                                  | Lab Sample ID | Analyte                        | Result | C | Unit | Spike Amount | Pct. Rec. | Limits | RPD | RPD Limit | Q |
|-----------------------------------------|---------------|--------------------------------|--------|---|------|--------------|-----------|--------|-----|-----------|---|
| Batch ID: 613561 Date: 05/23/2023 13:04 |               |                                |        |   |      |              |           |        |     |           |   |
| 350.1                                   | 280-176674-5  | Ammonia as N                   | 0.050  | U | mg/L |              |           |        |     |           |   |
| 350.1                                   | 280-176674-5  | Ammonia as N                   | 1.05   |   | mg/L | 1.00         | 105       | 90-110 |     |           |   |
| MS                                      |               |                                |        |   |      |              |           |        |     |           |   |
| Batch ID: 613268 Date: 05/19/2023 19:32 |               |                                |        |   |      |              |           |        |     |           |   |
| 9060                                    | 280-176674-7  | Total Organic Carbon - Average | 0.85   | J | mg/L |              |           |        |     |           |   |
| 9060                                    | 280-176674-7  | Total Organic Carbon - Average | 25.0   |   | mg/L | 25.0         | 97        | 88-112 |     |           |   |
| MS                                      |               |                                |        |   |      |              |           |        |     |           |   |

Calculations are performed before rounding to avoid round-off errors in calculated results.  
Note - Results and Reporting Limits have been adjusted for dry weight.

5-IN  
MATRIX SPIKE DUPLICATE SAMPLE RECOVERY  
GENERAL CHEMISTRY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
SDG No.: \_\_\_\_\_  
Matrix: Water

| Method                                  | Lab Sample ID | Analyte                | Result | C | Unit | Spike Amount | Pct. Rec. | Limits | RPD | RPD Limit | Q |
|-----------------------------------------|---------------|------------------------|--------|---|------|--------------|-----------|--------|-----|-----------|---|
| Batch ID: 613561 Date: 05/23/2023 13:06 |               |                        |        |   |      |              |           |        |     |           |   |
| 350.1                                   | 280-176674-5  | Ammonia as N           | 1.03   |   | mg/L | 1.00         | 103       | 90-110 | 1   | 10        |   |
| MSD                                     |               |                        |        |   |      |              |           |        |     |           |   |
| Batch ID: 613268 Date: 05/19/2023 19:46 |               |                        |        |   |      |              |           |        |     |           |   |
| 9060                                    | 280-176674-7  | Total Organic Carbon - | 25.1   |   | mg/L | 25.0         | 97        | 88-112 | 0   | 15        |   |
| MSD Average                             |               |                        |        |   |      |              |           |        |     |           |   |

Calculations are performed before rounding to avoid round-off errors in calculated results.  
Note - Results and Reporting Limits have been adjusted for dry weight.

7A-IN  
LAB CONTROL SAMPLE  
GENERAL CHEMISTRY

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Matrix: Water

| Method                                                                            | Lab Sample ID         | Analyte                        | Result | C | Unit | Spike Amount | Pct. Rec. | Limits | RPD | RPD Limit | Q |
|-----------------------------------------------------------------------------------|-----------------------|--------------------------------|--------|---|------|--------------|-----------|--------|-----|-----------|---|
| Batch ID: 613561 Date: 05/23/2023 12:10                                           |                       |                                |        |   |      |              |           |        |     |           |   |
| 350.1                                                                             | LCS<br>280-613561/89  | Ammonia as N                   | 2.58   |   | mg/L | 2.50         | 103       | 90-110 |     |           |   |
| LCS Source: 350.1 cal_00583                                                       |                       |                                |        |   |      |              |           |        |     |           |   |
| Batch ID: 613208 Date: 05/19/2023 17:27 Prep Batch: 613197 Date: 05/19/2023 15:58 |                       |                                |        |   |      |              |           |        |     |           |   |
| 9034                                                                              | LCS<br>280-613197/1-A | Sulfide                        | 16.0   |   | mg/L | 19.0         | 84        | 44-110 |     |           |   |
| LCS Source: SFD CAL INT_02403                                                     |                       |                                |        |   |      |              |           |        |     |           |   |
| Batch ID: 612961 Date: 05/18/2023 14:25                                           |                       |                                |        |   |      |              |           |        |     |           |   |
| 9056                                                                              | LCS<br>280-612961/11  | Nitrate as N                   | 2.46   |   | mg/L | 2.50         | 98        | 88-111 | 0   | 10        |   |
| 9056                                                                              | LCS<br>280-612961/11  | Nitrite as N                   | 2.42   |   | mg/L | 2.50         | 97        | 87-111 | 1   | 10        |   |
| LCS Source: IC LCS_01954                                                          |                       |                                |        |   |      |              |           |        |     |           |   |
| Batch ID: 612961 Date: 05/19/2023 00:56                                           |                       |                                |        |   |      |              |           |        |     |           |   |
| 9056                                                                              | LCS<br>280-612961/53  | Nitrate as N                   | 2.50   |   | mg/L | 2.50         | 100       | 88-111 | 2   | 10        |   |
| LCS Source: IC LCS_01954                                                          |                       |                                |        |   |      |              |           |        |     |           |   |
| Batch ID: 613351 Date: 05/22/2023 16:44                                           |                       |                                |        |   |      |              |           |        |     |           |   |
| 9056                                                                              | LCS<br>280-613351/4   | Nitrate as N                   | 2.41   |   | mg/L | 2.50         | 96        | 88-111 | 0   | 10        |   |
| LCS Source: IC LCS_01954                                                          |                       |                                |        |   |      |              |           |        |     |           |   |
| Batch ID: 612962 Date: 05/18/2023 14:25                                           |                       |                                |        |   |      |              |           |        |     |           |   |
| 9056A                                                                             | LCS<br>280-612962/11  | Sulfate                        | 50.5   |   | mg/L | 50.0         | 101       | 87-112 | 1   | 10        |   |
| LCS Source: IC LCS_01954                                                          |                       |                                |        |   |      |              |           |        |     |           |   |
| Batch ID: 613268 Date: 05/19/2023 18:31                                           |                       |                                |        |   |      |              |           |        |     |           |   |
| 9060                                                                              | LCS<br>280-613268/3   | Total Organic Carbon - Average | 24.8   |   | mg/L | 25.0         | 99        | 88-112 |     |           |   |
| LCS Source: TOC LCS Std_00057                                                     |                       |                                |        |   |      |              |           |        |     |           |   |
| Batch ID: 613263 Date: 05/19/2023 15:57                                           |                       |                                |        |   |      |              |           |        |     |           |   |
| SM<br>2320B                                                                       | LCS<br>280-613263/31  | Alkalinity                     | 181    |   | mg/L | 200          | 91        | 89-110 | 1   | 10        |   |
| LCS Source: Alk daily lcs_01088                                                   |                       |                                |        |   |      |              |           |        |     |           |   |

Calculations are performed before rounding to avoid round-off errors in calculated results.

FORM VIIA-IN

7A-IN  
LAB CONTROL SAMPLE DUPLICATE  
GENERAL CHEMISTRY

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Matrix: Water

| Method                                  | Lab Sample ID         | Analyte      | Result | C | Unit | Spike Amount                     | Pct. Rec. | Limits | RPD | RPD Limit | Q |
|-----------------------------------------|-----------------------|--------------|--------|---|------|----------------------------------|-----------|--------|-----|-----------|---|
| Batch ID: 612961 Date: 05/18/2023 14:40 |                       |              |        |   |      |                                  |           |        |     |           |   |
|                                         |                       |              |        |   |      | LCSD Source: IC LCS_01954        |           |        |     |           |   |
| 9056                                    | LCSD<br>280-612961/12 | Nitrate as N | 2.45   |   | mg/L | 2.50                             | 98        | 88-111 | 0   | 10        |   |
| 9056                                    | LCSD<br>280-612961/12 | Nitrite as N | 2.41   |   | mg/L | 2.50                             | 96        | 87-111 | 1   | 10        |   |
| Batch ID: 612961 Date: 05/19/2023 01:11 |                       |              |        |   |      |                                  |           |        |     |           |   |
|                                         |                       |              |        |   |      | LCSD Source: IC LCS_01954        |           |        |     |           |   |
| 9056                                    | LCSD<br>280-612961/54 | Nitrate as N | 2.45   |   | mg/L | 2.50                             | 98        | 88-111 | 2   | 10        |   |
| Batch ID: 613351 Date: 05/22/2023 16:59 |                       |              |        |   |      |                                  |           |        |     |           |   |
|                                         |                       |              |        |   |      | LCSD Source: IC LCS_01954        |           |        |     |           |   |
| 9056                                    | LCSD<br>280-613351/5  | Nitrate as N | 2.42   |   | mg/L | 2.50                             | 97        | 88-111 | 0   | 10        |   |
| Batch ID: 612962 Date: 05/18/2023 14:40 |                       |              |        |   |      |                                  |           |        |     |           |   |
|                                         |                       |              |        |   |      | LCSD Source: IC LCS_01954        |           |        |     |           |   |
| 9056A                                   | LCSD<br>280-612962/12 | Sulfate      | 50.0   |   | mg/L | 50.0                             | 100       | 87-112 | 1   | 10        |   |
| Batch ID: 613263 Date: 05/19/2023 16:04 |                       |              |        |   |      |                                  |           |        |     |           |   |
|                                         |                       |              |        |   |      | LCSD Source: Alk daily lcs_01088 |           |        |     |           |   |
| SM<br>2320B                             | LCSD<br>280-613263/32 | Alkalinity   | 183    |   | mg/L | 200                              | 92        | 89-110 | 1   | 10        |   |

Calculations are performed before rounding to avoid round-off errors in calculated results.

FORM VIIA-IN

7A-IN  
METHOD REPORTING LIMIT CHECK  
GENERAL CHEMISTRY

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Matrix: Water

| Method                                  | Lab Sample ID        | Analyte      | Result | C | Unit | Spike Amount                    | Pct. Rec. | Limits | RPD | RPD Limit | Q |
|-----------------------------------------|----------------------|--------------|--------|---|------|---------------------------------|-----------|--------|-----|-----------|---|
| Batch ID: 612961 Date: 05/18/2023 14:10 |                      |              |        |   |      |                                 |           |        |     |           |   |
|                                         |                      |              |        |   |      | LCS Source: IC Cal low_00709    |           |        |     |           |   |
| 9056                                    | MRL<br>280-612961/10 | Nitrate as N | 0.248  | J | mg/L | 0.250                           | 99        | 50-150 |     |           |   |
| 9056                                    | MRL<br>280-612961/10 | Nitrite as N | 0.221  | J | mg/L | 0.250                           | 89        | 50-150 |     |           |   |
| Batch ID: 613351 Date: 05/22/2023 16:29 |                      |              |        |   |      |                                 |           |        |     |           |   |
|                                         |                      |              |        |   |      | LCS Source: IC Cal low_00709    |           |        |     |           |   |
| 9056                                    | MRL<br>280-613351/3  | Nitrate as N | 0.239  | J | mg/L | 0.250                           | 95        | 50-150 |     |           |   |
| Batch ID: 612962 Date: 05/18/2023 14:10 |                      |              |        |   |      |                                 |           |        |     |           |   |
|                                         |                      |              |        |   |      | LCS Source: IC CAL cl/so4_00480 |           |        |     |           |   |
| 9056A                                   | MRL<br>280-612962/10 | Sulfate      | 2.98   | J | mg/L | 2.50                            | 119       | 50-150 |     |           |   |

Calculations are performed before rounding to avoid round-off errors in calculated results.

FORM VIIA-IN

16A-IN

Lab Name: Eurofins Denver Job No: 280-176674-1

SDG No.: \_\_\_\_\_ Analysis Batch No.: 612961

Instrument ID: WC\_IonChrom10 Calibration ID: 80543

Start Date: 05/18/2023 11:56 End Date: 05/18/2023 13:10

Analytical Method: 9056

| Analyte      | Corr. Coeff. | Slope    | Intercept | Calib. Type | Weighting    |
|--------------|--------------|----------|-----------|-------------|--------------|
| Nitrate as N | 0.9995       | 45000000 | -1160000  | WLR         | Inverse Conc |
| Nitrite as N | 1.0000       | 42900000 | -968000   | WLR         | Inverse Conc |



16B-IN

Lab Name: Eurofins Denver                      Job No:     280-176674-1

SDG No.: Analysis Batch No.: 612961

Instrument ID: WC\_IonChrom10                      Calibration ID: 80543

Start Date: 05/18/2023 11:56                      End Date: 05/18/2023 13:10

Analytical Method: 9056                      Concentration Units: ug/mL

| Analyte      | True | Found | %D | True | Found | %D | True | Found | %D |
|--------------|------|-------|----|------|-------|----|------|-------|----|
| Nitrate as N | 0.20 | 0.19  | -5 | 0.50 | 0.54  | 7  | 1.0  | 0.98  | -2 |
| Nitrite as N | 0.20 | 0.19  | -6 | 0.50 | 0.52  | 4  | 1.0  | 1.0   | 3  |

16B-IN

Lab Name: Eurofins Denver

Job No: 280-176674-1

SDG No.:

Analysis Batch No.: 612961

Instrument ID: WC\_IonChrom10

Calibration ID: 80543

Start Date: 05/18/2023 11:56

End Date: 05/18/2023 13:10

Analytical Method: 9056

Concentration Units: ug/mL

| Analyte      | True | Found | %D | True | Found | %D | True | Found | %D |
|--------------|------|-------|----|------|-------|----|------|-------|----|
| Nitrate as N | 2.5  | 2.5   | -1 | 5.0  | 5.0   | 0  |      |       |    |
| Nitrite as N | 2.5  | 2.5   | 0  | 5.0  | 5.0   | -1 |      |       |    |

16A-IN

|                              |                            |
|------------------------------|----------------------------|
| Lab Name: Eurofins Denver    | Job No: 280-176674-1       |
| SDG No.:                     | Analysis Batch No.: 612962 |
| Instrument ID: WC_IonChrom10 | Calibration ID: 80544      |
| Start Date: 05/18/2023 11:56 | End Date: 05/18/2023 13:10 |
| Analytical Method: 9056A     |                            |

| Analyte | Corr. Coeff. | Slope    | Intercept | Calib. Type | Weighting    |
|---------|--------------|----------|-----------|-------------|--------------|
| Sulfate | 0.9995       | 13800000 | 4650000   | WLR         | Inverse Conc |

16B-IN

Lab Name: Eurofins Denver                      Job No:     280-176674-1

SDG No.: Analysis Batch No.: 612962

Instrument ID: WC\_IonChrom10                      Calibration ID: 80544

Start Date: 05/18/2023 11:56                      End Date: 05/18/2023 13:10

Analytical Method: 9056A                      Concentration Units: ug/mL

| Analyte | True | Found | %D | True | Found | %D  | True | Found | %D |
|---------|------|-------|----|------|-------|-----|------|-------|----|
| Sulfate | 0.50 | 0.48  | -4 | 1.3  | 1.1   | -11 | 2.5  | 2.8   | 12 |

16B-IN

Lab Name: Eurofins Denver                      Job No:     280-176674-1

SDG No.: Analysis Batch No.: 612962

Instrument ID: WC\_IonChrom10                      Calibration ID: 80544

Start Date: 05/18/2023 11:56                      End Date: 05/18/2023 13:10

Analytical Method: 9056A                      Concentration Units: ug/mL

| Analyte | True | Found | %D | True | Found | %D | True | Found | %D |
|---------|------|-------|----|------|-------|----|------|-------|----|
| Sulfate | 30.0 | 31.0  | 3  | 60.0 | 59.9  | 0  | 100  | 99.0  | -1 |

9-IN  
DETECTION LIMITS  
GENERAL CHEMISTRY

Lab Name: Eurofins Denver Job Number: 280-176674-1  
SDG Number: \_\_\_\_\_  
Matrix: Water Instrument ID: WC\_SKALAR\_01  
Method: 350.1 DL Date: 04/29/2022 13:46

| Analyte      | Wavelength/<br>Mass | LOQ<br>(mg/L) | DL<br>(mg/L) |
|--------------|---------------------|---------------|--------------|
| Ammonia as N |                     | 0.1           | 0.029        |

9-IN  
CALIBRATION BLANK DETECTION LIMITS  
GENERAL CHEMISTRY

Lab Name: Eurofins Denver Job Number: 280-176674-1  
SDG Number: \_\_\_\_\_  
Matrix: Water Instrument ID: WC\_SKALAR\_01  
Method: 350.1 XMDL Date: 03/28/2011 13:26

| Analyte      | Wavelength/<br>Mass | XRL<br>(mg/L) | XMDL<br>(mg/L) |
|--------------|---------------------|---------------|----------------|
| Ammonia as N |                     | 0.1           | 0.0225         |

9-IN  
DETECTION LIMITS  
GENERAL CHEMISTRY

Lab Name: Eurofins Denver Job Number: 280-176674-1  
SDG Number: \_\_\_\_\_  
Matrix: Water Instrument ID: NOEQUIP  
Method: 9034 DL Date: 04/15/2022 13:08  
Prep Method: 9030B

| Analyte | Wavelength/<br>Mass | LOQ<br>(mg/L) | DL<br>(mg/L) |
|---------|---------------------|---------------|--------------|
| Sulfide |                     | 4             | 1.6          |



9-IN  
CALIBRATION BLANK DETECTION LIMITS  
GENERAL CHEMISTRY

Lab Name: Eurofins Denver Job Number: 280-176674-1  
SDG Number: \_\_\_\_\_  
Matrix: Water Instrument ID: NOEQUIP  
Method: 9034 XMDL Date: 04/15/2022 13:09

| Analyte | Wavelength/<br>Mass | XRL<br>(mg/L) | XMDL<br>(mg/L) |
|---------|---------------------|---------------|----------------|
| Sulfide |                     | 4             | 1.6            |

9-IN  
DETECTION LIMITS  
GENERAL CHEMISTRY

Lab Name: Eurofins Denver Job Number: 280-176674-1  
SDG Number: \_\_\_\_\_  
Matrix: Water Instrument ID: WC\_IonChrom10  
Method: 9056 DL Date: 06/21/2019 00:00

| Analyte      | Wavelength/<br>Mass | LOQ<br>(mg/L) | DL<br>(mg/L) |
|--------------|---------------------|---------------|--------------|
| Nitrate as N |                     | 0.5           | 0.0901       |
| Nitrite as N |                     | 0.5           | 0.049        |

9-IN  
CALIBRATION BLANK DETECTION LIMITS  
GENERAL CHEMISTRY

Lab Name: Eurofins Denver

Job Number: 280-176674-1

SDG Number: \_\_\_\_\_

Matrix: Water

Instrument ID: WC\_IonChrom10

Method: 9056

XMDL Date: 06/21/2019 00:00

| Analyte      | Wavelength/<br>Mass | XRL<br>(mg/L) | XMDL<br>(mg/L) |
|--------------|---------------------|---------------|----------------|
| Nitrate as N |                     | 0.5           | 0.0901         |
| Nitrite as N |                     | 0.5           | 0.049          |

9-IN  
DETECTION LIMITS  
GENERAL CHEMISTRY

Lab Name: Eurofins Denver Job Number: 280-176674-1  
SDG Number: \_\_\_\_\_  
Matrix: Water Instrument ID: WC\_IonChrom10  
Method: 9056A DL Date: 06/21/2019 00:00

| Analyte | Wavelength/<br>Mass | LOQ<br>(mg/L) | DL<br>(mg/L) |
|---------|---------------------|---------------|--------------|
| Sulfate |                     | 5             | 1.03         |

9-IN  
CALIBRATION BLANK DETECTION LIMITS  
GENERAL CHEMISTRY

Lab Name: Eurofins Denver Job Number: 280-176674-1  
SDG Number: \_\_\_\_\_  
Matrix: Water Instrument ID: WC\_IonChrom10  
Method: 9056A XMDL Date: 06/21/2019 00:00

| Analyte | Wavelength/<br>Mass | XRL<br>(mg/L) | XMDL<br>(mg/L) |
|---------|---------------------|---------------|----------------|
| Sulfate |                     | 5             | 1.03           |

9-IN  
DETECTION LIMITS  
GENERAL CHEMISTRY

Lab Name: Eurofins Denver Job Number: 280-176674-1  
SDG Number: \_\_\_\_\_  
Matrix: Water Instrument ID: WC\_SHI5  
Method: 9060 DL Date: 06/21/2019 00:00

| Analyte                           | Wavelength/<br>Mass | LOQ<br>(mg/L) | DL<br>(mg/L) |
|-----------------------------------|---------------------|---------------|--------------|
| Total Organic Carbon -<br>Average |                     | 1             | 0.345        |

9-IN  
CALIBRATION BLANK DETECTION LIMITS  
GENERAL CHEMISTRY

Lab Name: Eurofins Denver Job Number: 280-176674-1  
SDG Number: \_\_\_\_\_  
Matrix: Water Instrument ID: WC\_SHI5  
Method: 9060 XMDL Date: 06/21/2019 00:00

| Analyte                           | Wavelength/<br>Mass | XRL<br>(mg/L) | XMDL<br>(mg/L) |
|-----------------------------------|---------------------|---------------|----------------|
| Total Organic Carbon -<br>Average |                     | 1             | 0.345          |

9-IN  
DETECTION LIMITS  
GENERAL CHEMISTRY

Lab Name: Eurofins Denver Job Number: 280-176674-1  
SDG Number: \_\_\_\_\_  
Matrix: Water Instrument ID: WC\_AT4  
Method: SM 2320B DL Date: 02/03/2019 00:00

| Analyte    | Wavelength/<br>Mass | LOQ<br>(mg/L) | DL<br>(mg/L) |
|------------|---------------------|---------------|--------------|
| Alkalinity |                     | 10            | 3.08         |



9-IN  
CALIBRATION BLANK DETECTION LIMITS  
GENERAL CHEMISTRY

Lab Name: Eurofins Denver Job Number: 280-176674-1  
SDG Number: \_\_\_\_\_  
Matrix: Water Instrument ID: WC\_AT4  
Method: SM 2320B XMDL Date: 02/03/2019 00:00

| Analyte    | Wavelength/<br>Mass | XRL<br>(mg/L) | XMDL<br>(mg/L) |
|------------|---------------------|---------------|----------------|
| Alkalinity |                     | 10            | 3.08           |

12-IN  
PREPARATION LOG  
GENERAL CHEMISTRY

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Prep Method: 9030B

| Lab<br>Sample<br>ID | Preparation<br>Date | Prep<br>Batch | Initial<br>Weight | Initial<br>Volume<br>(mL) | Final<br>Volume<br>(mL) |
|---------------------|---------------------|---------------|-------------------|---------------------------|-------------------------|
| LCS 280-613197/1-A  | 05/19/2023 15:58    | 613197        |                   | 50                        | 50                      |
| MB 280-613197/2-A   | 05/19/2023 15:58    | 613197        |                   | 50                        | 50                      |
| 280-176674-7        | 05/19/2023 15:58    | 613197        |                   | 50                        | 50                      |
| 280-176674-8        | 05/19/2023 15:58    | 613197        |                   | 50                        | 50                      |

13-IN  
ANALYSIS RUN LOG  
GENERAL CHEMISTRY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
 SDG No.: \_\_\_\_\_  
 Instrument ID: WC\_SKALAR\_01 Analysis Method: 350.1  
 Start Date: 05/23/2023 08:15 End Date: 05/23/2023 17:47

| Lab Sample Id      | D/F | T<br>y<br>p<br>e | Time  | Analytes    |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
|--------------------|-----|------------------|-------|-------------|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|
|                    |     |                  |       | N<br>H<br>3 |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 08:15 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 08:17 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 08:21 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 08:24 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 08:27 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 08:29 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 08:32 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 08:34 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 08:37 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 08:40 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 08:42 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 08:45 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 08:48 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ICV 280-613561/14  | 1   |                  | 08:51 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ICVL 280-613561/15 | 1   |                  | 08:53 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ICB 280-613561/16  | 1   |                  | 08:56 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 08:58 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 09:01 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 09:04 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 09:06 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 09:09 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 09:12 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 09:14 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 09:17 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 09:20 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 09:22 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-613561/27  |     |                  | 09:25 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCVL 280-613561/28 |     |                  | 09:28 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-613561/29  |     |                  | 09:30 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 09:33 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 09:36 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 09:38 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 09:41 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 09:43 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 09:46 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 09:49 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 09:52 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 09:54 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 09:57 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-613561/40  |     |                  | 09:59 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCVL 280-613561/41 |     |                  | 10:02 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |

13-IN  
ANALYSIS RUN LOG  
GENERAL CHEMISTRY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
 SDG No.: \_\_\_\_\_  
 Instrument ID: WC\_SKALAR\_01 Analysis Method: 350.1  
 Start Date: 05/23/2023 08:15 End Date: 05/23/2023 17:47

| Lab Sample Id      | D/F | T<br>y<br>p<br>e | Time  | Analytes    |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
|--------------------|-----|------------------|-------|-------------|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|
|                    |     |                  |       | N<br>H<br>3 |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-613561/42  |     |                  | 10:05 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 10:08 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 10:10 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 10:14 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 10:16 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 10:18 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 10:21 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 10:24 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 10:26 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 10:29 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 10:32 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 10:34 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 10:37 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-613561/55  |     |                  | 10:40 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCVL 280-613561/56 |     |                  | 10:42 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-613561/57  |     |                  | 10:45 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 10:48 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 10:50 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 10:53 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 10:56 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 10:58 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 11:01 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 11:04 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 11:06 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 11:09 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 11:12 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-613561/68  |     |                  | 11:14 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCVL 280-613561/69 |     |                  | 11:17 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-613561/70  |     |                  | 11:20 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 11:22 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 11:25 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 11:28 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 11:30 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 11:33 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 11:36 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 11:38 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 11:41 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 11:44 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 11:46 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 11:49 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 11:52 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |

13-IN  
ANALYSIS RUN LOG  
GENERAL CHEMISTRY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
SDG No.: \_\_\_\_\_  
Instrument ID: WC\_SKALAR\_01 Analysis Method: 350.1  
Start Date: 05/23/2023 08:15 End Date: 05/23/2023 17:47

| Lab Sample Id       | D/F | T<br>y<br>p<br>e | Time  | Analytes    |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
|---------------------|-----|------------------|-------|-------------|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|
|                     |     |                  |       | N<br>H<br>3 |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-613561/83   | 1   |                  | 11:54 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCVL 280-613561/84  | 1   |                  | 11:57 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-613561/85   | 1   |                  | 12:00 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 12:02 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 12:05 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 12:08 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| LCS 280-613561/89   | 1   | T                | 12:10 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| MB 280-613561/90    | 1   | T                | 12:13 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 12:16 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 12:18 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 12:21 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 12:24 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 12:26 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-613561/96   | 1   |                  | 12:29 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCVL 280-613561/97  | 1   |                  | 12:32 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-613561/98   | 1   |                  | 12:34 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 12:37 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 12:40 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 12:43 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 12:45 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 12:48 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 12:50 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 12:53 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 12:56 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 12:58 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| 280-176674-5        | 1   | T                | 13:01 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| 280-176674-5 MS     | 1   | T                | 13:04 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| 280-176674-5 MSD    | 1   | T                | 13:06 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-613561/111  | 1   |                  | 13:09 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCVL 280-613561/112 | 1   |                  | 13:12 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-613561/113  | 1   |                  | 13:14 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 13:17 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 13:20 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 13:22 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| 280-176674-6        | 1   | T                | 13:25 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| 280-176674-9        | 1   | T                | 13:28 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 13:30 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 13:33 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 13:37 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 13:39 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 13:43 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |

13-IN  
ANALYSIS RUN LOG  
GENERAL CHEMISTRY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
 SDG No.: \_\_\_\_\_  
 Instrument ID: WC\_SKALAR\_01 Analysis Method: 350.1  
 Start Date: 05/23/2023 08:15 End Date: 05/23/2023 17:47

| Lab Sample Id       | D/F | T<br>y<br>p<br>e | Time  | Analytes    |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
|---------------------|-----|------------------|-------|-------------|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|
|                     |     |                  |       | N<br>H<br>3 |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-613561/124  | 1   |                  | 13:46 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCVL 280-613561/125 | 1   |                  | 13:49 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-613561/126  |     |                  | 13:51 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 13:54 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 13:57 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 13:59 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 14:02 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 14:05 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 14:08 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 14:10 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 14:13 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 14:15 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 14:18 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 14:21 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 14:24 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-613561/139  |     |                  | 14:26 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCVL 280-613561/140 |     |                  | 14:29 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-613561/141  |     |                  | 14:32 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 14:34 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 14:37 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 14:40 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 14:42 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 14:45 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 14:48 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 14:50 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 14:53 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 14:56 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 14:58 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-613561/152  |     |                  | 15:01 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCVL 280-613561/153 |     |                  | 15:04 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-613561/154  |     |                  | 15:06 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 15:09 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 15:12 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 15:15 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 15:17 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 15:20 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 15:22 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 15:25 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 15:28 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 15:30 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 15:33 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |

13-IN  
ANALYSIS RUN LOG  
GENERAL CHEMISTRY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
 SDG No.: \_\_\_\_\_  
 Instrument ID: WC\_SKALAR\_01 Analysis Method: 350.1  
 Start Date: 05/23/2023 08:15 End Date: 05/23/2023 17:47

| Lab Sample Id       | D/F | T<br>y<br>p<br>e | Time  | Analytes    |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
|---------------------|-----|------------------|-------|-------------|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|
|                     |     |                  |       | N<br>H<br>3 |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 15:35 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 15:38 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-613561/167  |     |                  | 15:41 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCVL 280-613561/168 |     |                  | 15:44 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-613561/169  |     |                  | 15:46 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 15:49 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 15:51 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 15:54 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 15:57 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 16:00 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 16:02 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 16:05 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 16:08 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 16:10 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 16:13 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-613561/180  |     |                  | 16:16 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCVL 280-613561/181 |     |                  | 16:18 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-613561/182  |     |                  | 16:21 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 16:24 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 16:26 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 16:30 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 16:32 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 16:34 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 16:37 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 16:39 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 16:42 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 16:45 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 16:48 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 16:50 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 16:53 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-613561/195  | 1   |                  | 16:56 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCVL 280-613561/196 | 1   |                  | 16:58 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-613561/197  | 1   |                  | 17:01 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| 280-176674-3        | 200 | T                | 17:04 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 17:06 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 17:09 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 17:12 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 17:14 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 17:17 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 17:19 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 17:22 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |

13-IN  
ANALYSIS RUN LOG  
GENERAL CHEMISTRY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
 SDG No.: \_\_\_\_\_  
 Instrument ID: WC\_SKALAR\_01 Analysis Method: 350.1  
 Start Date: 05/23/2023 08:15 End Date: 05/23/2023 17:47

| Lab Sample Id       | D/F | T<br>y<br>p<br>e | Time  | Analytes    |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
|---------------------|-----|------------------|-------|-------------|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|
|                     |     |                  |       | N<br>H<br>3 |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 17:25 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 17:28 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 17:30 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-613561/209  | 1   |                  | 17:33 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCVL 280-613561/210 | 1   |                  | 17:36 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-613561/211  | 1   |                  | 17:38 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 17:41 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 17:43 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 17:47 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |

Prep Types: \_\_\_\_\_  
 T = Total/NA



13-IN  
ANALYSIS RUN LOG  
GENERAL CHEMISTRY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
 SDG No.: \_\_\_\_\_  
 Instrument ID: NOEQUIP Analysis Method: 9034  
 Start Date: 05/19/2023 17:27 End Date: 05/19/2023 17:27

| Lab Sample Id      | D/F | T<br>y<br>p<br>e | Time  | Analytes |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
|--------------------|-----|------------------|-------|----------|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|
|                    |     |                  |       | S<br>2   |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| LCS 280-613197/1-A | 1   | T                | 17:27 | X        |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| MB 280-613197/2-A  | 1   | T                | 17:27 | X        |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 17:27 |          |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 17:27 |          |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 17:27 |          |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 17:27 |          |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 17:27 |          |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 17:27 |          |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 17:27 |          |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 17:27 |          |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| 280-176674-7       | 1   | T                | 17:27 | X        |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| 280-176674-8       | 1   | T                | 17:27 | X        |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 17:27 |          |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |

Prep Types: \_\_\_\_\_  
 T = Total/NA

13-IN  
ANALYSIS RUN LOG  
GENERAL CHEMISTRY

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Instrument ID: WC\_IonChrom10 Analysis Method: 9056

Start Date: 05/18/2023 11:56 End Date: 05/19/2023 11:12

| Lab Sample Id       | D/F | T<br>y<br>p<br>e | Time  | Analytes              |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
|---------------------|-----|------------------|-------|-----------------------|-------------|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|
|                     |     |                  |       | N<br>O<br>2<br>-<br>N | N<br>O<br>3 |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| STD 280-612961/1 IC |     |                  | 11:56 | X                     | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| STD 280-612961/2 IC | 1   |                  | 12:11 | X                     | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| STD 280-612961/3 IC | 1   |                  | 12:26 | X                     | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| STD 280-612961/4 IC | 1   |                  | 12:40 | X                     | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| STD 280-612961/5 IC | 1   |                  | 12:55 | X                     | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| STD 280-612961/6 IC | 1   |                  | 13:10 | X                     | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ICV 280-612961/7    | 1   |                  | 13:25 | X                     | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ICB 280-612961/8    | 1   |                  | 13:40 | X                     | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| RTC 280-612961/9    |     |                  | 13:55 |                       |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| MRL 280-612961/10   | 1   | T                | 14:10 | X                     | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| LCS 280-612961/11   | 1   | T                | 14:25 | X                     | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| LCSD 280-612961/12  | 1   | T                | 14:40 | X                     | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| MB 280-612961/13    | 1   | T                | 14:55 | X                     | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 15:10 |                       |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 15:25 |                       |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 15:40 |                       |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 15:56 |                       |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 16:11 |                       |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| 280-176674-5        | 20  | T                | 16:26 |                       | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| 280-176674-6        | 20  | T                | 16:41 |                       | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 16:56 |                       |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 17:11 |                       |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 17:26 |                       |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-612961/24   | 1   |                  | 17:41 | X                     | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-612961/25   | 1   |                  | 17:56 | X                     | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 18:11 |                       |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 18:26 |                       |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 18:41 |                       |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 18:56 |                       |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 19:11 |                       |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| 280-176674-7        | 1   | T                | 19:26 | X                     | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| 280-176674-8        | 1   | T                | 19:41 | X                     | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 19:56 |                       |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 20:11 |                       |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 20:27 |                       |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-612961/36   | 1   |                  | 20:42 | X                     | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-612961/37   | 1   |                  | 20:57 | X                     | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 21:12 |                       |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 21:27 |                       |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |

|                              |                            |
|------------------------------|----------------------------|
| Lab Name: Eurofins Denver    | Job No.: 280-176674-1      |
| SDG No.:                     |                            |
| Instrument ID: WC_IonChrom10 | Analysis Method: 9056      |
| Start Date: 05/18/2023 11:56 | End Date: 05/19/2023 11:12 |

5/31/2023 9:08 PM

13-IN  
ANALYSIS RUN LOG  
GENERAL CHEMISTRY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
 SDG No.: \_\_\_\_\_  
 Instrument ID: WC\_IonChrom10 Analysis Method: 9056  
 Start Date: 05/18/2023 11:56 End Date: 05/19/2023 11:12

| Lab Sample Id      | D/F | T<br>y<br>p<br>e | Time  | Analytes              |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
|--------------------|-----|------------------|-------|-----------------------|-------------|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|
|                    |     |                  |       | N<br>O<br>2<br>-<br>N | N<br>O<br>3 |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 07:27 |                       |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 07:42 |                       |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 07:57 |                       |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 08:12 |                       |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 08:27 |                       |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-612961/84  | 1   |                  | 08:42 | X                     |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-612961/85  | 1   |                  | 08:57 | X                     |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 09:12 |                       |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 09:27 |                       |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 09:42 |                       |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 09:57 |                       |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 10:12 |                       |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 10:27 |                       |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 10:42 |                       |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-612961/93  |     |                  | 10:57 |                       |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-612961/107 |     |                  | 11:12 |                       |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |

Prep Types: \_\_\_\_\_  
 T = Total/NA

13-IN  
ANALYSIS RUN LOG  
GENERAL CHEMISTRY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
 SDG No.: \_\_\_\_\_  
 Instrument ID: WC\_IonChrom10 Analysis Method: 9056  
 Start Date: 05/22/2023 16:00 End Date: 05/23/2023 03:30

| Lab Sample Id     | D/F  | T<br>y<br>p<br>e | Time  | Analytes    |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
|-------------------|------|------------------|-------|-------------|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|
|                   |      |                  |       | N<br>O<br>3 |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-613351/1  | 1    |                  | 16:00 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-613351/2  | 1    |                  | 16:14 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| MRL 280-613351/3  | 1    | T                | 16:29 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| LCS 280-613351/4  | 1    | T                | 16:44 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| LCSD 280-613351/5 | 1    | T                | 16:59 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| MB 280-613351/6   | 1    | T                | 17:14 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |      |                  | 17:29 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |      |                  | 17:44 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |      |                  | 17:59 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |      |                  | 18:14 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |      |                  | 18:29 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |      |                  | 18:44 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| 280-176674-3      | 2000 | T                | 18:59 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |      |                  | 19:14 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |      |                  | 19:29 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |      |                  | 19:44 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-613351/17 | 1    |                  | 19:59 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-613351/18 | 1    |                  | 20:15 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |      |                  | 20:30 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |      |                  | 20:45 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |      |                  | 21:00 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-613351/29 |      |                  | 23:00 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-613351/30 |      |                  | 23:15 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |      |                  | 01:31 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |      |                  | 01:46 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-613351/41 |      |                  | 02:01 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-613351/42 |      |                  | 02:16 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |      |                  | 02:31 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |      |                  | 02:46 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |      |                  | 03:01 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-613351/46 |      |                  | 03:16 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-613351/47 |      |                  | 03:30 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |

Prep Types: \_\_\_\_\_  
 T = Total/NA

13-IN  
ANALYSIS RUN LOG  
GENERAL CHEMISTRY

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Instrument ID: WC\_IonChrom10 Analysis Method: 9056A

Start Date: 05/18/2023 11:56 End Date: 05/19/2023 11:12

| Lab Sample Id       | D/F | T<br>y<br>p<br>e | Time  | Analytes    |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
|---------------------|-----|------------------|-------|-------------|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|
|                     |     |                  |       | S<br>O<br>4 |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| STD 280-612962/1 IC | 1   |                  | 11:56 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| STD 280-612962/2 IC | 1   |                  | 12:11 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| STD 280-612962/3 IC | 1   |                  | 12:26 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| STD 280-612962/4 IC | 1   |                  | 12:40 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| STD 280-612962/5 IC | 1   |                  | 12:55 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| STD 280-612962/6 IC | 1   |                  | 13:10 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ICV 280-612962/7    | 1   |                  | 13:25 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ICB 280-612962/8    | 1   |                  | 13:40 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| RTC 280-612962/9    |     |                  | 13:55 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| MRL 280-612962/10   | 1   | T                | 14:10 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| LCS 280-612962/11   | 1   | T                | 14:25 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| LCSD 280-612962/12  | 1   | T                | 14:40 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| MB 280-612962/13    | 1   | T                | 14:55 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 15:25 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 15:40 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 15:56 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 16:11 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 16:26 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 16:41 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 16:56 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 17:11 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 17:26 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-612962/24   | 1   |                  | 17:41 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-612962/25   | 1   |                  | 17:56 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 18:11 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 18:26 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 18:41 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 18:56 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 19:11 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| 280-176674-7        | 1   | T                | 19:26 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| 280-176674-8        | 1   | T                | 19:41 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 19:56 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 20:11 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 20:27 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-612962/36   | 1   |                  | 20:42 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-612962/37   | 1   |                  | 20:57 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 21:12 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 21:27 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 21:42 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 21:57 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ              |     |                  | 22:12 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |

13-IN  
ANALYSIS RUN LOG  
GENERAL CHEMISTRY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
 SDG No.: \_\_\_\_\_  
 Instrument ID: WC\_IonChrom10 Analysis Method: 9056A  
 Start Date: 05/18/2023 11:56 End Date: 05/19/2023 11:12

| Lab Sample Id     | D/F | T<br>y<br>p<br>e | Time  | Analytes    |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
|-------------------|-----|------------------|-------|-------------|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|
|                   |     |                  |       | S<br>O<br>4 |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 22:27 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 22:42 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 22:57 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 23:12 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 23:26 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-612962/48 |     |                  | 23:41 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-612962/49 |     |                  | 23:56 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 00:11 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 00:26 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 00:41 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 00:56 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 01:11 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 01:26 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 01:41 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 01:56 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 02:11 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 02:26 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-612962/60 |     |                  | 02:42 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-612962/61 |     |                  | 02:57 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 03:12 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 03:27 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 03:42 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 03:57 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 04:12 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 04:27 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 04:57 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 05:12 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 05:27 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-612962/72 |     |                  | 05:42 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-612962/73 |     |                  | 05:57 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 06:12 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 06:27 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 06:42 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 06:57 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 07:12 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 07:27 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 07:42 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 07:57 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 08:12 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 08:27 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-612962/84 |     |                  | 08:42 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |

13-IN  
ANALYSIS RUN LOG  
GENERAL CHEMISTRY

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Instrument ID: WC\_IonChrom10 Analysis Method: 9056A

Start Date: 05/18/2023 11:56 End Date: 05/19/2023 11:12

| Lab Sample Id      | D/F | T<br>y<br>p<br>e | Time  | Analytes    |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
|--------------------|-----|------------------|-------|-------------|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|
|                    |     |                  |       | S<br>O<br>4 |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-612962/85  |     |                  | 08:57 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 09:12 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 09:27 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 09:42 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 09:57 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 10:12 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 10:27 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 10:42 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-612962/93  |     |                  | 10:57 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-612962/107 |     |                  | 11:12 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |

Prep Types: \_\_\_\_\_

T = Total/NA



13-IN  
ANALYSIS RUN LOG  
GENERAL CHEMISTRY

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Instrument ID: WC\_SHI5 Analysis Method: 9060

Start Date: 05/19/2023 18:03 End Date: 05/20/2023 17:05

| Lab Sample Id     | D/F | T<br>y<br>p<br>e | Time  | Analytes         |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
|-------------------|-----|------------------|-------|------------------|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|
|                   |     |                  |       | T<br>O<br>C<br>Q |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ICV 280-613268/1  | 1   |                  | 18:03 | X                |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ICB 280-613268/2  | 1   |                  | 18:17 | X                |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| LCS 280-613268/3  | 1   | T                | 18:31 | X                |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| MB 280-613268/4   | 1   | T                | 18:46 | X                |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 19:00 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| 280-176674-7      | 1   | T                | 19:17 | X                |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| 280-176674-7 MS   | 1   | T                | 19:32 | X                |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| 280-176674-7 MSD  | 1   | T                | 19:46 | X                |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 20:02 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 20:17 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 20:33 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 20:47 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 21:01 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 21:16 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-613268/15 | 1   |                  | 21:30 | X                |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-613268/16 | 1   |                  | 21:45 | X                |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 22:00 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 22:16 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 22:31 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 22:45 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 22:59 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 23:14 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 23:30 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 23:44 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 23:59 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 00:13 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-613268/27 |     |                  | 00:28 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-613268/28 |     |                  | 00:42 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 00:57 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 01:13 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 01:27 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 01:45 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 02:02 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 02:16 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 02:31 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 02:45 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 02:59 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 03:13 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-613268/39 |     |                  | 03:30 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-613268/40 |     |                  | 03:45 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |

13-IN  
ANALYSIS RUN LOG  
GENERAL CHEMISTRY

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Instrument ID: WC\_SHI5 Analysis Method: 9060

Start Date: 05/19/2023 18:03 End Date: 05/20/2023 17:05

| Lab Sample Id     | D/F | T<br>y<br>p<br>e | Time  | Analytes         |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
|-------------------|-----|------------------|-------|------------------|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|
|                   |     |                  |       | T<br>O<br>C<br>Q |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 03:59 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 04:14 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 04:30 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 04:44 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 05:02 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 05:19 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 05:33 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 05:47 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 06:02 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 06:16 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-613268/51 |     |                  | 06:32 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-613268/52 |     |                  | 06:47 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 07:01 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 07:17 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 07:33 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 07:48 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 08:04 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 08:21 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 08:37 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 08:51 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 09:07 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 09:22 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-613268/63 |     |                  | 09:36 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-613268/64 |     |                  | 09:51 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 10:05 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 10:19 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 10:33 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 10:48 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 11:02 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 11:18 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 11:34 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 11:49 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 12:03 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 12:19 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-613268/75 |     |                  | 12:36 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-613268/76 |     |                  | 12:50 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 13:05 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 13:19 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 13:35 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 13:49 |                  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |

13-IN  
ANALYSIS RUN LOG  
GENERAL CHEMISTRY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
 SDG No.: \_\_\_\_\_  
 Instrument ID: WC\_SHI5 Analysis Method: 9060  
 Start Date: 05/19/2023 18:03 End Date: 05/20/2023 17:05

| Lab Sample Id     | D/F | T<br>y<br>p<br>e | Time  | Analytes    |   |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
|-------------------|-----|------------------|-------|-------------|---|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|
|                   |     |                  |       | T<br>O<br>C | Q |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 14:06 |             |   |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 14:21 |             |   |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 14:35 |             |   |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 14:49 |             |   |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 15:06 |             |   |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 15:20 |             |   |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-613268/87 |     |                  | 15:35 |             |   |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-613268/88 |     |                  | 15:50 |             |   |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 16:04 |             |   |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 16:20 |             |   |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 16:34 |             |   |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-613268/92 |     |                  | 16:51 |             |   |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-613268/93 |     |                  | 17:05 |             |   |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |

Prep Types: \_\_\_\_\_  
 T = Total/NA

13-IN  
ANALYSIS RUN LOG  
GENERAL CHEMISTRY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
 SDG No.: \_\_\_\_\_  
 Instrument ID: WC\_AT4 Analysis Method: SM 2320B  
 Start Date: 05/19/2023 12:32 End Date: 05/20/2023 07:16

| Lab Sample Id      | D/F | T<br>y<br>p<br>e | Time  | Analytes    |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
|--------------------|-----|------------------|-------|-------------|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|
|                    |     |                  |       | A<br>l<br>k |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| RINSE 280-613263/1 |     |                  | 12:32 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 12:40 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 12:43 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 12:51 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 12:58 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 13:04 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 13:12 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 13:21 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 13:26 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 13:33 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 13:39 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 13:47 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 13:53 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 14:00 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 14:04 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 14:11 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-613263/17  |     |                  | 14:19 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-613263/18  |     |                  | 14:25 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 14:31 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 14:37 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 14:42 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 14:47 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 14:52 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 14:59 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 15:06 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 15:15 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 15:26 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 15:35 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-613263/29  |     | 1                | 15:42 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-613263/30  |     | 1                | 15:49 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| LCS 280-613263/31  |     | 1 T              | 15:57 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| LCSD 280-613263/32 |     | 1 T              | 16:04 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| MB 280-613263/33   |     | 1 T              | 16:11 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 16:22 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 16:34 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 16:40 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 16:48 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 16:55 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 17:01 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 17:06 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 17:11 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |

13-IN  
ANALYSIS RUN LOG  
GENERAL CHEMISTRY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
 SDG No.: \_\_\_\_\_  
 Instrument ID: WC\_AT4 Analysis Method: SM 2320B  
 Start Date: 05/19/2023 12:32 End Date: 05/20/2023 07:16

| Lab Sample Id     | D/F | T<br>y<br>p<br>e | Time  | Analytes    |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
|-------------------|-----|------------------|-------|-------------|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|
|                   |     |                  |       | A<br>l<br>k |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| 280-176674-7      | 1   | T                | 17:16 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| 280-176674-8      | 1   | T                | 17:22 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-613263/44 | 1   |                  | 17:30 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-613263/45 | 1   |                  | 17:36 | X           |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 17:43 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 17:50 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 17:55 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 18:01 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 18:14 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 18:21 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 18:28 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 18:38 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 18:50 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 19:01 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-613263/56 |     |                  | 19:09 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-613263/57 |     |                  | 19:14 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 19:22 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 19:28 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 19:40 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 19:51 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 20:05 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 20:21 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 20:34 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 20:49 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 20:55 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 21:00 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 21:07 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 21:15 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-613263/70 |     |                  | 21:23 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-613263/71 |     |                  | 21:29 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 21:36 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 21:45 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 21:53 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 22:02 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 22:13 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 22:18 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 22:25 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 22:32 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 22:52 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ            |     |                  | 23:12 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-613263/82 |     |                  | 23:21 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |

13-IN  
ANALYSIS RUN LOG  
GENERAL CHEMISTRY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
 SDG No.: \_\_\_\_\_  
 Instrument ID: WC\_AT4 Analysis Method: SM 2320B  
 Start Date: 05/19/2023 12:32 End Date: 05/20/2023 07:16

| Lab Sample Id      | D/F | T<br>y<br>p<br>e | Time  | Analytes    |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
|--------------------|-----|------------------|-------|-------------|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|
|                    |     |                  |       | A<br>l<br>k |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-613263/83  |     |                  | 23:27 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 23:35 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 23:41 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 23:47 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 23:53 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 23:59 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 00:05 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 00:11 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 00:17 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 00:23 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 00:29 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 00:37 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 00:43 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-613263/96  |     |                  | 00:51 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-613263/97  |     |                  | 00:57 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 01:12 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 01:20 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 01:41 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 02:03 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 02:10 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 02:16 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 02:24 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 02:34 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 02:42 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 02:51 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-613263/108 |     |                  | 02:59 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-613263/109 |     |                  | 03:05 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 03:12 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 03:18 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 03:26 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 03:33 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 03:39 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 03:45 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 03:51 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 03:57 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 04:03 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 04:09 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 04:15 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ             |     |                  | 04:22 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-613263/122 |     |                  | 04:30 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-613263/123 |     |                  | 04:36 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |

13-IN  
ANALYSIS RUN LOG  
GENERAL CHEMISTRY

Lab Name: Eurofins Denver Job No.: 280-176674-1  
 SDG No.: \_\_\_\_\_  
 Instrument ID: WC\_AT4 Analysis Method: SM 2320B  
 Start Date: 05/19/2023 12:32 End Date: 05/20/2023 07:16

| Lab Sample Id        | D/F | T<br>y<br>p<br>e | Time  | Analytes    |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
|----------------------|-----|------------------|-------|-------------|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|
|                      |     |                  |       | A<br>l<br>k |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ               |     |                  | 04:44 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ               |     |                  | 04:51 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ               |     |                  | 04:58 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ               |     |                  | 05:04 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ               |     |                  | 05:15 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ               |     |                  | 05:35 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ               |     |                  | 05:56 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ               |     |                  | 06:15 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ               |     |                  | 06:35 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| ZZZZZZ               |     |                  | 06:55 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| RINSE 280-613263/134 |     |                  | 07:02 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCV 280-613263/135   |     |                  | 07:10 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| CCB 280-613263/136   |     |                  | 07:16 |             |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |

Prep Types: \_\_\_\_\_  
 T = Total/NA

## GENERAL CHEMISTRY BATCH WORKSHEET

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Batch Number: 613561 Batch Start Date: 05/23/23 08:15 Batch Analyst: Peterson, McKenzie MBatch Method: 350.1 Batch End Date: \_\_\_\_\_

| Lab Sample ID          | Client Sample ID         | Method Chain | Basis | ClResPres | InitialAmount | FinalAmount | 350.1 cal 00583 | 350.1 ICV 00564 |  |
|------------------------|--------------------------|--------------|-------|-----------|---------------|-------------|-----------------|-----------------|--|
| ICV<br>280-613561/14   |                          | 350.1        |       | no        | 100 mL        | 100 mL      |                 | 2.5 mL          |  |
| ICVL<br>280-613561/15  |                          | 350.1        |       | no        | 100 mL        | 100 mL      |                 | 0.5 mL          |  |
| ICB<br>280-613561/16   |                          | 350.1        |       | no        | 30 mL         | 30 mL       |                 |                 |  |
| CCV<br>280-613561/83   |                          | 350.1        |       | no        | 100 mL        | 100 mL      | 2.5 mL          |                 |  |
| CCVL<br>280-613561/84  |                          | 350.1        |       | no        | 100 mL        | 100 mL      | 0.5 mL          |                 |  |
| CCB<br>280-613561/85   |                          | 350.1        |       | no        | 10 mL         | 10 mL       |                 |                 |  |
| LCS<br>280-613561/89   |                          | 350.1        |       | no        | 100 mL        | 100 mL      | 2.5 mL          |                 |  |
| MB<br>280-613561/90    |                          | 350.1        |       | no        | 30 mL         | 30 mL       |                 |                 |  |
| CCV<br>280-613561/96   |                          | 350.1        |       | no        | 100 mL        | 100 mL      | 2.5 mL          |                 |  |
| CCVL<br>280-613561/97  |                          | 350.1        |       | no        | 100 mL        | 100 mL      | 0.5 mL          |                 |  |
| CCB<br>280-613561/98   |                          | 350.1        |       | no        | 10 mL         | 10 mL       |                 |                 |  |
| 280-176674-B-5         | LL12mw-185-23040<br>1-GW | 350.1        | T     | no        | 10 mL         | 10 mL       |                 |                 |  |
| 280-176674-B-5<br>MS   | LL12mw-185-23040<br>1-GW | 350.1        | T     | no        | 10 mL         | 10 mL       | 0.1 mL          |                 |  |
| 280-176674-B-5<br>MSD  | LL12mw-185-23040<br>1-GW | 350.1        | T     | no        | 10 mL         | 10 mL       | 0.1 mL          |                 |  |
| CCV<br>280-613561/111  |                          | 350.1        |       | no        | 100 mL        | 100 mL      | 2.5 mL          |                 |  |
| CCVL<br>280-613561/112 |                          | 350.1        |       | no        | 100 mL        | 100 mL      | 0.5 mL          |                 |  |
| CCB<br>280-613561/113  |                          | 350.1        |       | no        | 10 mL         | 10 mL       |                 |                 |  |
| 280-176674-B-6         | LL12mw-185-23040<br>2-GW | 350.1        | T     | no        | 10 mL         | 10 mL       |                 |                 |  |
| 280-176674-B-9         | LL12mw-245-23040<br>1-GW | 350.1        | T     | no        | 10 mL         | 10 mL       |                 |                 |  |
| CCV<br>280-613561/124  |                          | 350.1        |       | no        | 100 mL        | 100 mL      | 2.5 mL          |                 |  |
| CCVL<br>280-613561/125 |                          | 350.1        |       | no        | 100 mL        | 100 mL      | 0.5 mL          |                 |  |

The pound sign (#) in the amount added field denotes that the reagent was used undiluted. All calculations are performed using the stated concentration for this reagent.

350.1

Page 1 of 2



## GENERAL CHEMISTRY BATCH WORKSHEET

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Batch Number: 613561 Batch Start Date: 05/23/23 08:15 Batch Analyst: Peterson, McKenzie MBatch Method: 350.1 Batch End Date: \_\_\_\_\_

| Lab Sample ID          | Client Sample ID         | Method Chain | Basis | ClResPres | InitialAmount | FinalAmount | 350.1 cal 00583 | 350.1 ICV 00564 |  |
|------------------------|--------------------------|--------------|-------|-----------|---------------|-------------|-----------------|-----------------|--|
| CCV<br>280-613561/195  |                          | 350.1        |       | no        | 100 mL        | 100 mL      | 2.5 mL          |                 |  |
| CCVL<br>280-613561/196 |                          | 350.1        |       | no        | 100 mL        | 100 mL      | 0.5 mL          |                 |  |
| CCB<br>280-613561/197  |                          | 350.1        |       | no        | 10 mL         | 10 mL       |                 |                 |  |
| 280-176674-B-3         | LL12mw-187-23040<br>1-GW | 350.1        | T     | no        | 10 mL         | 10 mL       |                 |                 |  |
| CCV<br>280-613561/209  |                          | 350.1        |       | no        | 100 mL        | 100 mL      | 2.5 mL          |                 |  |
| CCVL<br>280-613561/210 |                          | 350.1        |       | no        | 100 mL        | 100 mL      | 0.5 mL          |                 |  |
| CCB<br>280-613561/211  |                          | 350.1        |       | no        | 10 mL         | 10 mL       |                 |                 |  |

| Batch Notes                  |                                   |
|------------------------------|-----------------------------------|
| Dechlorination Agent ID      | Na2S2O3_00002                     |
| pH Indicator ID              | HC293086                          |
| Acid used for pH adjustment  | SulfuricAcid_00288                |
| Sodium Nitroprusside ID      | Na Nitro_00050                    |
| Hypochlorite ID              | Na Hypo_00050                     |
| EDTA Buffer ID               | Buffer A_00039                    |
| Potassium Sodium Tartrate ID | Buffer B_00058                    |
| Carrier Identification       | Ammonia Rinse_00056               |
| Sodium Salicylate ID         | Sodium Sal_00028                  |
| Pipette/Syringe/Dispenser ID | SAH500, 1000 INTERCESSOR, BWH 200 |

| Basis | Basis Description |
|-------|-------------------|
| T     | Total/NA          |

The pound sign (#) in the amount added field denotes that the reagent was used undiluted. All calculations are performed using the stated concentration for this reagent.

350.1

Page 2 of 2

## GENERAL CHEMISTRY BATCH WORKSHEET

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Batch Number: 613197 Batch Start Date: 05/19/23 15:58 Batch Analyst: Ince, Connor ABatch Method: 9030B Batch End Date: \_\_\_\_\_

| Lab Sample ID       | Client Sample ID          | Method Chain | Basis | InitialAmount | FinalAmount | DistillUnitPort | Initial pH | Final pH | AcidVolAdded |
|---------------------|---------------------------|--------------|-------|---------------|-------------|-----------------|------------|----------|--------------|
| LCS<br>280-613197/1 |                           | 9030B, 9034  |       | 50 mL         | 50 mL       | 1               | 7 SU       | <2 SU    | 5 mL         |
| MB 280-613197/2     |                           | 9030B, 9034  |       | 50 mL         | 50 mL       | 2               | 7 SU       | <2 SU    | 5 mL         |
| 280-176674-E-7      | FBQmw-173-230401<br>-GW-R | 9030B, 9034  | T     | 50 mL         | 50 mL       | 11              | 12+ SU     | <2 SU    | 5 mL         |
| 280-176674-B-8      | FBQmw-173-230402<br>-GW-R | 9030B, 9034  | T     | 50 mL         | 50 mL       | 12              | 12+ SU     | <2 SU    | 5 mL         |

| Lab Sample ID       | Client Sample ID          | Method Chain | Basis | SFD CAL INT<br>02403 |  |  |  |  |  |
|---------------------|---------------------------|--------------|-------|----------------------|--|--|--|--|--|
| LCS<br>280-613197/1 |                           | 9030B, 9034  |       | 1 mL                 |  |  |  |  |  |
| MB 280-613197/2     |                           | 9030B, 9034  |       |                      |  |  |  |  |  |
| 280-176674-E-7      | FBQmw-173-230401<br>-GW-R | 9030B, 9034  | T     |                      |  |  |  |  |  |
| 280-176674-B-8      | FBQmw-173-230402<br>-GW-R | 9030B, 9034  | T     |                      |  |  |  |  |  |

| Batch Notes                       |                    |
|-----------------------------------|--------------------|
| Formaldehyde ID                   | Form_00125         |
| Zinc Acetate Buffer ID            | znac_00127         |
| Pipette/Syringe/Dispenser ID      | WC-5000BOD1        |
| Sulfuric Acid Reagent ID Number   | SulfuricAcid_00288 |
| Balance ID                        | M19170             |
| pH Indicator ID                   | HC157843           |
| Oven, Bath or Block Temperature 1 | 80 Degrees C       |
| Distillation Start Time           | 05/19/2023 12:30   |
| Uncorrected Temperature 2         | 80 Degrees C       |
| Oven, Bath or Block Temperature 2 | 80 Degrees C       |
| Distillation End Time             | 05/19/2023 14:00   |
| Batch Comment                     | CAI                |

| Basis | Basis Description |
|-------|-------------------|
| T     | Total/NA          |

The pound sign (#) in the amount added field denotes that the reagent was used undiluted. All calculations are performed using the stated concentration for this reagent.

## GENERAL CHEMISTRY BATCH WORKSHEET

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Batch Number: 613208 Batch Start Date: 05/19/23 17:27 Batch Analyst: Ince, Connor ABatch Method: 9034 Batch End Date: \_\_\_\_\_

| Lab Sample ID         | Client Sample ID          | Method Chain | Basis | BuretStart1 | BuretStop1 | IodineAmount | TitrantVolume1            | FinalAmount |  |
|-----------------------|---------------------------|--------------|-------|-------------|------------|--------------|---------------------------|-------------|--|
| LCS<br>280-613197/1-A |                           | 9034         |       | 0.0 mL      | 3.0 mL     | 1.25 mL      | 3 mL                      | 50 mL       |  |
| MB<br>280-613197/2-A  |                           | 9034         |       | 3.0 mL      | 4.0 mL     | 0.25 mL      | 1 mL                      | 50 mL       |  |
| 280-176674-E-7-A      | FBQmw-173-230401<br>-GW-R | 9034         | T     | 16.1 mL     | 17.1 mL    | 0.25 mL      | 1 mL                      | 50 mL       |  |
| 280-176674-B-8-A      | FBQmw-173-230402<br>-GW-R | 9034         | T     | 17.1 mL     | 18.0 mL    | 0.25 mL      | 0.89999999999999<br>99 mL | 50 mL       |  |

| Batch Notes                  |                  |
|------------------------------|------------------|
| Normality of First Titrant   | 0.025 N          |
| Normality of Iodine Solution | 0.1 N            |
| Zinc Acetate Buffer ID       | znac_00127       |
| Iodine ID                    | Iod_00279        |
| Hydrochloric Acid ID         | HCL Sol_00207    |
| Sodium Thiosulfate ID        | Na Thio_00181    |
| Starch Reagent ID            | Starch Ind_00071 |
| Pipette/Syringe/Dispenser ID | WC-5000BOD1      |
| Batch Comment                | CAI              |

| Basis | Basis Description |
|-------|-------------------|
| T     | Total/NA          |

The pound sign (#) in the amount added field denotes that the reagent was used undiluted. All calculations are performed using the stated concentration for this reagent.

## GENERAL CHEMISTRY BATCH WORKSHEET

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Batch Number: 612961 Batch Start Date: 05/18/23 11:56 Batch Analyst: Collins, Michael EBatch Method: 9056 Batch End Date: \_\_\_\_\_

| Lab Sample ID          | Client Sample ID          | Method Chain | Basis | InitialAmount | FinalAmount | Cl ICV Std<br>00006 | IC CAL cl/so4<br>00480 | IC Cal low<br>00709 | IC ICV 5 00405 |
|------------------------|---------------------------|--------------|-------|---------------|-------------|---------------------|------------------------|---------------------|----------------|
| STD<br>280-612961/2 IC |                           | 9056         |       | 10 mL         | 10 mL       |                     | 0.1 mL                 | 0.08 mL             |                |
| STD<br>280-612961/3 IC |                           | 9056         |       | 10 mL         | 10 mL       |                     | 0.2 mL                 | 0.2 mL              |                |
| STD<br>280-612961/4 IC |                           | 9056         |       | 10 mL         | 10 mL       |                     | 2.4 mL                 | 0.4 mL              |                |
| STD<br>280-612961/5 IC |                           | 9056         |       | 10 mL         | 10 mL       |                     | 4.8 mL                 | 1 mL                |                |
| STD<br>280-612961/6 IC |                           | 9056         |       | 10 mL         | 10 mL       |                     | 8 mL                   | 2 mL                |                |
| ICV<br>280-612961/7    |                           | 9056         |       | 10 mL         | 10 mL       | 0.8 mL              |                        |                     | 0.8 mL         |
| ICB<br>280-612961/8    |                           | 9056         |       | 10 mL         | 10 mL       |                     |                        |                     |                |
| MRL<br>280-612961/10   |                           | 9056         |       | 10 mL         | 10 mL       |                     | 0.1 mL                 | 0.05 mL             |                |
| LCS<br>280-612961/11   |                           | 9056         |       | 10 mL         | 10 mL       |                     |                        |                     |                |
| LCSD<br>280-612961/12  |                           | 9056         |       | 10 mL         | 10 mL       |                     |                        |                     |                |
| MB<br>280-612961/13    |                           | 9056         |       | 10 mL         | 10 mL       |                     |                        |                     |                |
| 280-176674-A-5         | LL12mw-185-23040<br>1-GW  | 9056         | T     | 10 mL         | 10 mL       |                     |                        |                     |                |
| 280-176674-A-6         | LL12mw-185-23040<br>2-GW  | 9056         | T     | 10 mL         | 10 mL       |                     |                        |                     |                |
| CCV<br>280-612961/24   |                           | 9056         |       | 10 mL         | 10 mL       |                     |                        |                     |                |
| CCB<br>280-612961/25   |                           | 9056         |       | 10 mL         | 10 mL       |                     |                        |                     |                |
| 280-176674-A-7         | FBQmw-173-230401<br>-GW-R | 9056         | T     | 10 mL         | 10 mL       |                     |                        |                     |                |
| 280-176674-A-8         | FBQmw-173-230402<br>-GW-R | 9056         | T     | 10 mL         | 10 mL       |                     |                        |                     |                |
| CCV<br>280-612961/36   |                           | 9056         |       | 10 mL         | 10 mL       |                     |                        |                     |                |
| CCB<br>280-612961/37   |                           | 9056         |       | 10 mL         | 10 mL       |                     |                        |                     |                |
| CCV<br>280-612961/48   |                           | 9056         |       | 10 mL         | 10 mL       |                     |                        |                     |                |
| CCB<br>280-612961/49   |                           | 9056         |       | 10 mL         | 10 mL       |                     |                        |                     |                |

The pound sign (#) in the amount added field denotes that the reagent was used undiluted. All calculations are performed using the stated concentration for this reagent.

## GENERAL CHEMISTRY BATCH WORKSHEET

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Batch Number: 612961 Batch Start Date: 05/18/23 11:56 Batch Analyst: Collins, Michael EBatch Method: 9056 Batch End Date: \_\_\_\_\_

| Lab Sample ID         | Client Sample ID         | Method Chain | Basis | InitialAmount | FinalAmount | Cl ICV Std<br>00006 | IC CAL cl/so4<br>00480 | IC Cal low<br>00709 | IC ICV 5 00405 |
|-----------------------|--------------------------|--------------|-------|---------------|-------------|---------------------|------------------------|---------------------|----------------|
| LCS<br>280-612961/53  |                          | 9056         |       | 10 mL         | 10 mL       |                     |                        |                     |                |
| LCSD<br>280-612961/54 |                          | 9056         |       | 10 mL         | 10 mL       |                     |                        |                     |                |
| MB<br>280-612961/55   |                          | 9056         |       | 10 mL         | 10 mL       |                     |                        |                     |                |
| CCV<br>280-612961/60  |                          | 9056         |       | 10 mL         | 10 mL       |                     |                        |                     |                |
| CCB<br>280-612961/61  |                          | 9056         |       | 10 mL         | 10 mL       |                     |                        |                     |                |
| CCV<br>280-612961/72  |                          | 9056         |       | 10 mL         | 10 mL       |                     |                        |                     |                |
| CCB<br>280-612961/73  |                          | 9056         |       | 10 mL         | 10 mL       |                     |                        |                     |                |
| 280-176674-A-9        | LL12mw-245-23040<br>1-GW | 9056         | T     | 10 mL         | 10 mL       |                     |                        |                     |                |
| CCV<br>280-612961/84  |                          | 9056         |       | 10 mL         | 10 mL       |                     |                        |                     |                |
| CCB<br>280-612961/85  |                          | 9056         |       | 10 mL         | 10 mL       |                     |                        |                     |                |

| Lab Sample ID          | Client Sample ID | Method Chain | Basis | IC LCS 01954 | IC SO4 ICV<br>00024 |  |  |  |  |
|------------------------|------------------|--------------|-------|--------------|---------------------|--|--|--|--|
| STD<br>280-612961/2 IC |                  | 9056         |       |              |                     |  |  |  |  |
| STD<br>280-612961/3 IC |                  | 9056         |       |              |                     |  |  |  |  |
| STD<br>280-612961/4 IC |                  | 9056         |       |              |                     |  |  |  |  |
| STD<br>280-612961/5 IC |                  | 9056         |       |              |                     |  |  |  |  |
| STD<br>280-612961/6 IC |                  | 9056         |       |              |                     |  |  |  |  |
| ICV<br>280-612961/7    |                  | 9056         |       |              | 0.8 mL              |  |  |  |  |
| ICB<br>280-612961/8    |                  | 9056         |       |              |                     |  |  |  |  |
| MRL<br>280-612961/10   |                  | 9056         |       |              |                     |  |  |  |  |
| LCS<br>280-612961/11   |                  | 9056         |       | 5 mL         |                     |  |  |  |  |

The pound sign (#) in the amount added field denotes that the reagent was used undiluted. All calculations are performed using the stated concentration for this reagent.

## GENERAL CHEMISTRY BATCH WORKSHEET

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Batch Number: 612961 Batch Start Date: 05/18/23 11:56 Batch Analyst: Collins, Michael EBatch Method: 9056 Batch End Date: \_\_\_\_\_

| Lab Sample ID         | Client Sample ID          | Method Chain | Basis | IC LCS 01954 | IC S04 ICV<br>00024 |  |  |  |  |
|-----------------------|---------------------------|--------------|-------|--------------|---------------------|--|--|--|--|
| LCSD<br>280-612961/12 |                           | 9056         |       | 5 mL         |                     |  |  |  |  |
| MB<br>280-612961/13   |                           | 9056         |       |              |                     |  |  |  |  |
| 280-176674-A-5        | LL12mw-185-23040<br>1-GW  | 9056         | T     |              |                     |  |  |  |  |
| 280-176674-A-6        | LL12mw-185-23040<br>2-GW  | 9056         | T     |              |                     |  |  |  |  |
| CCV<br>280-612961/24  |                           | 9056         |       | 10 mL        |                     |  |  |  |  |
| CCB<br>280-612961/25  |                           | 9056         |       |              |                     |  |  |  |  |
| 280-176674-A-7        | FBQmw-173-230401<br>-GW-R | 9056         | T     |              |                     |  |  |  |  |
| 280-176674-A-8        | FBQmw-173-230402<br>-GW-R | 9056         | T     |              |                     |  |  |  |  |
| CCV<br>280-612961/36  |                           | 9056         |       | 10 mL        |                     |  |  |  |  |
| CCB<br>280-612961/37  |                           | 9056         |       |              |                     |  |  |  |  |
| CCV<br>280-612961/48  |                           | 9056         |       | 10 mL        |                     |  |  |  |  |
| CCB<br>280-612961/49  |                           | 9056         |       |              |                     |  |  |  |  |
| LCS<br>280-612961/53  |                           | 9056         |       | 10 mL        |                     |  |  |  |  |
| LCSD<br>280-612961/54 |                           | 9056         |       | 10 mL        |                     |  |  |  |  |
| MB<br>280-612961/55   |                           | 9056         |       |              |                     |  |  |  |  |
| CCV<br>280-612961/60  |                           | 9056         |       | 10 mL        |                     |  |  |  |  |
| CCB<br>280-612961/61  |                           | 9056         |       |              |                     |  |  |  |  |
| CCV<br>280-612961/72  |                           | 9056         |       | 10 mL        |                     |  |  |  |  |
| CCB<br>280-612961/73  |                           | 9056         |       |              |                     |  |  |  |  |
| 280-176674-A-9        | LL12mw-245-23040<br>1-GW  | 9056         | T     |              |                     |  |  |  |  |
| CCV<br>280-612961/84  |                           | 9056         |       | 10 mL        |                     |  |  |  |  |

The pound sign (#) in the amount added field denotes that the reagent was used undiluted. All calculations are performed using the stated concentration for this reagent.

## GENERAL CHEMISTRY BATCH WORKSHEET

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Batch Number: 612961 Batch Start Date: 05/18/23 11:56 Batch Analyst: Collins, Michael EBatch Method: 9056 Batch End Date: \_\_\_\_\_

|                      |                  |              |       |              |                     |  |  |  |  |
|----------------------|------------------|--------------|-------|--------------|---------------------|--|--|--|--|
| Lab Sample ID        | Client Sample ID | Method Chain | Basis | IC LCS 01954 | IC S04 ICV<br>00024 |  |  |  |  |
| CCB<br>280-612961/85 |                  | 9056         |       |              |                     |  |  |  |  |

| Batch Notes                    |                               |
|--------------------------------|-------------------------------|
| Filter ID                      | SF020E                        |
| Pipette/Syringe/Dispenser ID   | ARM5000, 100HEX, 200CJ, 100IC |
| Sufficient Volume for Batch QC | yes                           |
| Eluent 1 ID                    | IC10 Eluent_00006             |
| Batch Comment                  | MEC                           |

| Basis | Basis Description |
|-------|-------------------|
| T     | Total/NA          |

The pound sign (#) in the amount added field denotes that the reagent was used undiluted. All calculations are performed using the stated concentration for this reagent.

## GENERAL CHEMISTRY BATCH WORKSHEET

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Batch Number: 613351 Batch Start Date: 05/22/23 16:00 Batch Analyst: Sherman, Erik JBatch Method: 9056 Batch End Date: \_\_\_\_\_

| Lab Sample ID        | Client Sample ID         | Method Chain | Basis | InitialAmount | FinalAmount | IC CAL cl/so4<br>00480 | IC Cal low<br>00709 | IC LCS 01954 |  |
|----------------------|--------------------------|--------------|-------|---------------|-------------|------------------------|---------------------|--------------|--|
| CCV<br>280-613351/1  |                          | 9056         |       | 10 mL         | 10 mL       |                        |                     | 10 mL        |  |
| CCB<br>280-613351/2  |                          | 9056         |       | 10 mL         | 10 mL       |                        |                     |              |  |
| MRL<br>280-613351/3  |                          | 9056         |       | 10 mL         | 10 mL       | 0.2 mL                 | 0.1 mL              |              |  |
| LCS<br>280-613351/4  |                          | 9056         |       | 10 mL         | 10 mL       |                        |                     | 10 mL        |  |
| LCSD<br>280-613351/5 |                          | 9056         |       | 10 mL         | 10 mL       |                        |                     | 10 mL        |  |
| MB 280-613351/6      |                          | 9056         |       | 10 mL         | 10 mL       |                        |                     |              |  |
| 280-176674-A-3       | LL12mw-187-23040<br>1-GW | 9056         | T     | 10 mL         | 10 mL       |                        |                     |              |  |
| CCV<br>280-613351/17 |                          | 9056         |       | 10 mL         | 10 mL       |                        |                     | 10 mL        |  |
| CCB<br>280-613351/18 |                          | 9056         |       | 10 mL         | 10 mL       |                        |                     |              |  |

| Batch Notes                    |                                                 |
|--------------------------------|-------------------------------------------------|
| Filter ID                      | SF020E                                          |
| Pipette/Syringe/Dispenser ID   | ARM5000, 1000HEX, 200CJ, IC100, AB8A100, BMF100 |
| Sufficient Volume for Batch QC | Y                                               |
| Eluent 1 ID                    | IC10 ELUENT_00006                               |
| Batch Comment                  | EJS                                             |

| Basis | Basis Description |
|-------|-------------------|
| T     | Total/NA          |

The pound sign (#) in the amount added field denotes that the reagent was used undiluted. All calculations are performed using the stated concentration for this reagent.



## GENERAL CHEMISTRY BATCH WORKSHEET

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Batch Number: 612962 Batch Start Date: 05/18/23 11:56 Batch Analyst: Collins, Michael EBatch Method: 9056A Batch End Date: \_\_\_\_\_

| Lab Sample ID          | Client Sample ID          | Method Chain | Basis | InitialAmount | FinalAmount | Cl ICV Std<br>00006 | IC CAL cl/so4<br>00480 | IC Cal low<br>00709 | IC ICV 5 00405 |
|------------------------|---------------------------|--------------|-------|---------------|-------------|---------------------|------------------------|---------------------|----------------|
| STD<br>280-612962/1 IC |                           | 9056A        |       | 10 mL         | 10 mL       |                     | 0.04 mL                |                     |                |
| STD<br>280-612962/2 IC |                           | 9056A        |       | 10 mL         | 10 mL       |                     | 0.1 mL                 | 0.08 mL             |                |
| STD<br>280-612962/3 IC |                           | 9056A        |       | 10 mL         | 10 mL       |                     | 0.2 mL                 | 0.2 mL              |                |
| STD<br>280-612962/4 IC |                           | 9056A        |       | 10 mL         | 10 mL       |                     | 2.4 mL                 | 0.4 mL              |                |
| STD<br>280-612962/5 IC |                           | 9056A        |       | 10 mL         | 10 mL       |                     | 4.8 mL                 | 1 mL                |                |
| STD<br>280-612962/6 IC |                           | 9056A        |       | 10 mL         | 10 mL       |                     | 8 mL                   | 2 mL                |                |
| ICV<br>280-612962/7    |                           | 9056A        |       | 10 mL         | 10 mL       | 0.8 mL              |                        |                     | 0.8 mL         |
| ICB<br>280-612962/8    |                           | 9056A        |       | 10 mL         | 10 mL       |                     |                        |                     |                |
| MRL<br>280-612962/10   |                           | 9056A        |       | 10 mL         | 10 mL       |                     | 0.1 mL                 | 0.05 mL             |                |
| LCS<br>280-612962/11   |                           | 9056A        |       | 10 mL         | 10 mL       |                     |                        |                     |                |
| LCSD<br>280-612962/12  |                           | 9056A        |       | 10 mL         | 10 mL       |                     |                        |                     |                |
| MB<br>280-612962/13    |                           | 9056A        |       | 10 mL         | 10 mL       |                     |                        |                     |                |
| CCV<br>280-612962/24   |                           | 9056A        |       | 10 mL         | 10 mL       |                     |                        |                     |                |
| CCB<br>280-612962/25   |                           | 9056A        |       | 10 mL         | 10 mL       |                     |                        |                     |                |
| 280-176674-A-7         | FBQmw-173-230401<br>-GW-R | 9056A        | T     | 10 mL         | 10 mL       |                     |                        |                     |                |
| 280-176674-A-8         | FBQmw-173-230402<br>-GW-R | 9056A        | T     | 10 mL         | 10 mL       |                     |                        |                     |                |
| CCV<br>280-612962/36   |                           | 9056A        |       | 10 mL         | 10 mL       |                     |                        |                     |                |
| CCB<br>280-612962/37   |                           | 9056A        |       | 10 mL         | 10 mL       |                     |                        |                     |                |

| Lab Sample ID          | Client Sample ID | Method Chain | Basis | IC LCS 01954 | IC SO4 ICV<br>00024 |  |  |  |  |
|------------------------|------------------|--------------|-------|--------------|---------------------|--|--|--|--|
| STD<br>280-612962/1 IC |                  | 9056A        |       |              |                     |  |  |  |  |

The pound sign (#) in the amount added field denotes that the reagent was used undiluted. All calculations are performed using the stated concentration for this reagent.

9056A

Page 1 of 3

## GENERAL CHEMISTRY BATCH WORKSHEET

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Batch Number: 612962 Batch Start Date: 05/18/23 11:56 Batch Analyst: Collins, Michael EBatch Method: 9056A Batch End Date: \_\_\_\_\_

| Lab Sample ID          | Client Sample ID          | Method Chain | Basis | IC LCS 01954 | IC S04 ICV<br>00024 |  |  |  |  |
|------------------------|---------------------------|--------------|-------|--------------|---------------------|--|--|--|--|
| STD<br>280-612962/2 IC |                           | 9056A        |       |              |                     |  |  |  |  |
| STD<br>280-612962/3 IC |                           | 9056A        |       |              |                     |  |  |  |  |
| STD<br>280-612962/4 IC |                           | 9056A        |       |              |                     |  |  |  |  |
| STD<br>280-612962/5 IC |                           | 9056A        |       |              |                     |  |  |  |  |
| STD<br>280-612962/6 IC |                           | 9056A        |       |              |                     |  |  |  |  |
| ICV<br>280-612962/7    |                           | 9056A        |       |              | 0.8 mL              |  |  |  |  |
| ICB<br>280-612962/8    |                           | 9056A        |       |              |                     |  |  |  |  |
| MRL<br>280-612962/10   |                           | 9056A        |       |              |                     |  |  |  |  |
| LCS<br>280-612962/11   |                           | 9056A        |       | 5 mL         |                     |  |  |  |  |
| LCS<br>280-612962/12   |                           | 9056A        |       | 5 mL         |                     |  |  |  |  |
| MB<br>280-612962/13    |                           | 9056A        |       |              |                     |  |  |  |  |
| CCV<br>280-612962/24   |                           | 9056A        |       | 10 mL        |                     |  |  |  |  |
| CCB<br>280-612962/25   |                           | 9056A        |       |              |                     |  |  |  |  |
| 280-176674-A-7         | FBQmw-173-230401<br>-GW-R | 9056A        | T     |              |                     |  |  |  |  |
| 280-176674-A-8         | FBQmw-173-230402<br>-GW-R | 9056A        | T     |              |                     |  |  |  |  |
| CCV<br>280-612962/36   |                           | 9056A        |       | 10 mL        |                     |  |  |  |  |
| CCB<br>280-612962/37   |                           | 9056A        |       |              |                     |  |  |  |  |

The pound sign (#) in the amount added field denotes that the reagent was used undiluted. All calculations are performed using the stated concentration for this reagent.

9056A

Page 2 of 3

## GENERAL CHEMISTRY BATCH WORKSHEET

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Batch Number: 612962 Batch Start Date: 05/18/23 11:56 Batch Analyst: Collins, Michael EBatch Method: 9056A Batch End Date: \_\_\_\_\_

| Batch Notes                    |                                |
|--------------------------------|--------------------------------|
| Filter ID                      | SF020E                         |
| Pipette/Syringe/Dispenser ID   | ARM5000, 1000HEX, 200CJ, 100IC |
| Sufficient Volume for Batch QC | yes                            |
| Eluent 1 ID                    | IC10 Eluent_00006              |
| Batch Comment                  | MEC                            |

| Basis | Basis Description |
|-------|-------------------|
| T     | Total/NA          |

The pound sign (#) in the amount added field denotes that the reagent was used undiluted. All calculations are performed using the stated concentration for this reagent.

9056A

Page 3 of 3

## GENERAL CHEMISTRY BATCH WORKSHEET

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Batch Number: 613268 Batch Start Date: 05/19/23 18:03 Batch Analyst: Wondimu, Akewak BBatch Method: 9060 Batch End Date: \_\_\_\_\_

| Lab Sample ID         | Client Sample ID          | Method Chain | Basis | InitialAmount | FinalAmount | Initial pH | TOC ICV Std<br>00051 | TOC LCS Std<br>00057 | AnalysisComment                                                                    |
|-----------------------|---------------------------|--------------|-------|---------------|-------------|------------|----------------------|----------------------|------------------------------------------------------------------------------------|
| ICV<br>280-613268/1   |                           | 9060         |       | 50 mL         | 50 mL       | <2 SU      | 1 mL                 |                      | The final volume of these samples was fixed at 50 by the TOCControlV preprocessor. |
| ICB<br>280-613268/2   |                           | 9060         |       | 20 mL         | 20 mL       | <2 SU      |                      |                      | The final volume of these samples was fixed at 50 by the TOCControlV preprocessor. |
| LCS<br>280-613268/3   |                           | 9060         |       | 200 mL        | 200 mL      | <2 SU      |                      | 5 mL                 | The final volume of these samples was fixed at 50 by the TOCControlV preprocessor. |
| MB 280-613268/4       |                           | 9060         |       | 20 mL         | 20 mL       | <2 SU      |                      |                      | The final volume of these samples was fixed at 50 by the TOCControlV preprocessor. |
| 280-176674-D-7        | FBQmw-173-230401<br>-GW-R | 9060         | T     | 20 mL         | 20 mL       | <2 SU      |                      |                      | The final volume of these samples was fixed at 50 by the TOCControlV preprocessor. |
| 280-176674-D-7<br>MS  | FBQmw-173-230401<br>-GW-R | 9060         | T     | 50 mL         | 50 mL       | <2 SU      |                      | 1.25 mL              | The final volume of these samples was fixed at 50 by the TOCControlV preprocessor. |
| 280-176674-D-7<br>MSD | FBQmw-173-230401<br>-GW-R | 9060         | T     | 50 mL         | 50 mL       | <2 SU      |                      | 1.25 mL              | The final volume of these samples was fixed at 50 by the TOCControlV preprocessor. |

The pound sign (#) in the amount added field denotes that the reagent was used undiluted. All calculations are performed using the stated concentration for this reagent.

9060

Page 1 of 2

## GENERAL CHEMISTRY BATCH WORKSHEET

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Batch Number: 613268 Batch Start Date: 05/19/23 18:03 Batch Analyst: Wondimu, Akewak BBatch Method: 9060 Batch End Date: \_\_\_\_\_

| Lab Sample ID        | Client Sample ID | Method Chain | Basis | InitialAmount | FinalAmount | Initial pH | TOC ICV Std<br>00051 | TOC LCS Std<br>00057 | AnalysisComment                                                                                   |
|----------------------|------------------|--------------|-------|---------------|-------------|------------|----------------------|----------------------|---------------------------------------------------------------------------------------------------|
| CCV<br>280-613268/15 |                  | 9060         |       | 200 mL        | 200 mL      | <2 SU      |                      | 5 mL                 | The final<br>volume of these<br>samples was<br>fixed at 50 by<br>the TOCControlV<br>preprocessor. |
| CCB<br>280-613268/16 |                  | 9060         |       | 20 mL         | 20 mL       | <2 SU      |                      |                      | The final<br>volume of these<br>samples was<br>fixed at 50 by<br>the TOCControlV<br>preprocessor. |

| Batch Notes                  |                               |
|------------------------------|-------------------------------|
| Acid ID                      | 0.2% H2SO4_00462, H2SO4_00214 |
| Combustion Catalyst ID       | S638-92069-02                 |
| Pipette/Syringe/Dispenser ID | BWH 5000,1000 MS              |
| Vial Lot Number              | 705212201761                  |
| Batch Comment                | Run and reviewed by ABW       |

| Basis | Basis Description |
|-------|-------------------|
| T     | Total/NA          |

The pound sign (#) in the amount added field denotes that the reagent was used undiluted. All calculations are performed using the stated concentration for this reagent.

## GENERAL CHEMISTRY BATCH WORKSHEET

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Batch Number: 613270 Batch Start Date: 05/19/23 18:03 Batch Analyst: Wondimu, Akewak BBatch Method: 9060 Batch End Date: \_\_\_\_\_

| Lab Sample ID         | Client Sample ID          | Method Chain | Basis | InitialAmount | FinalAmount | Initial pH | Final pH | TOC ICV Std<br>00051 | TOC LCS Std<br>00057 |
|-----------------------|---------------------------|--------------|-------|---------------|-------------|------------|----------|----------------------|----------------------|
| LCS<br>280-613270/3   |                           | 9060         |       | 200 mL        | 200 mL      | <2 SU      | <2 SU    |                      | 5 mL                 |
| MB 280-613270/4       |                           | 9060         |       | 20 mL         | 20 mL       | <2 SU      | <2 SU    |                      |                      |
| 280-176674-D-7        | FBQmw-173-230401<br>-GW-R | 9060         | T     | 20 mL         | 20 mL       | <2 SU      | <2 SU    |                      |                      |
| 280-176674-D-7<br>MS  | FBQmw-173-230401<br>-GW-R | 9060         | T     | 50 mL         | 50 mL       | <2 SU      | <2 SU    |                      | 1.25 mL              |
| 280-176674-D-7<br>MSD | FBQmw-173-230401<br>-GW-R | 9060         | T     | 50 mL         | 50 mL       | <2 SU      | <2 SU    |                      | 1.25 mL              |
| LCSD<br>280-613270/35 |                           | 9060         |       | 100 mL        | 100 mL      | <2 SU      | <2 SU    | 2.5 mL               |                      |

| Lab Sample ID        | Client Sample ID          | Method Chain | Basis | AnalysisComment                                                                    |  |  |  |  |  |
|----------------------|---------------------------|--------------|-------|------------------------------------------------------------------------------------|--|--|--|--|--|
| LCS<br>280-613270/3  |                           | 9060         |       | The final volume of these samples was fixed at 50 by the TOCControlV preprocessor. |  |  |  |  |  |
| MB 280-613270/4      |                           | 9060         |       | The final volume of these samples was fixed at 50 by the TOCControlV preprocessor. |  |  |  |  |  |
| 280-176674-D-7       | FBQmw-173-230401<br>-GW-R | 9060         | T     | The final volume of these samples was fixed at 50 by the TOCControlV preprocessor. |  |  |  |  |  |
| 280-176674-D-7<br>MS | FBQmw-173-230401<br>-GW-R | 9060         | T     | The final volume of these samples was fixed at 50 by the TOCControlV preprocessor. |  |  |  |  |  |

The pound sign (#) in the amount added field denotes that the reagent was used undiluted. All calculations are performed using the stated concentration for this reagent.

## GENERAL CHEMISTRY BATCH WORKSHEET

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Batch Number: 613270 Batch Start Date: 05/19/23 18:03 Batch Analyst: Wondimu, Akewak BBatch Method: 9060 Batch End Date: \_\_\_\_\_

| Lab Sample ID         | Client Sample ID          | Method Chain | Basis | AnalysisComment                                                                                   |  |  |  |  |  |
|-----------------------|---------------------------|--------------|-------|---------------------------------------------------------------------------------------------------|--|--|--|--|--|
| 280-176674-D-7<br>MSD | FBQmw-173-230401<br>-GW-R | 9060         | T     | The final<br>volume of these<br>samples was<br>fixed at 50 by<br>the TOCControlV<br>preprocessor. |  |  |  |  |  |
| LCSD<br>280-613270/35 |                           | 9060         |       | The final<br>volume of these<br>samples was<br>fixed at 50 by<br>the TOCControlV<br>preprocessor. |  |  |  |  |  |

| Batch Notes                  |                               |
|------------------------------|-------------------------------|
| Acid ID                      | 0.2% H2SO4_00462, H2SO4_00214 |
| Combustion Catylyst ID       | S638-92069-02                 |
| Pipette/Syringe/Dispenser ID | BWH 5000,1000 MS              |
| pH Paper ID                  | HC205981                      |
| Batch Comment                | Run and reviewed by ABW       |

| Basis | Basis Description |
|-------|-------------------|
| T     | Total/NA          |

The pound sign (#) in the amount added field denotes that the reagent was used undiluted. All calculations are performed using the stated concentration for this reagent.

9060

Page 2 of 2

## GENERAL CHEMISTRY BATCH WORKSHEET

Lab Name: Eurofins Denver Job No.: 280-176674-1

SDG No.: \_\_\_\_\_

Batch Number: 613263 Batch Start Date: 05/19/23 12:32 Batch Analyst: Guillen, Kristela EBatch Method: SM 2320B Batch End Date: \_\_\_\_\_

| Lab Sample ID          | Client Sample ID          | Method Chain | Basis | CalcMsg                   | Alk daily lcs<br>01088 |  |  |  |  |
|------------------------|---------------------------|--------------|-------|---------------------------|------------------------|--|--|--|--|
| CCV<br>280-613263/29   |                           | SM 2320B     |       | InitialAmount<br>is blank | 10 mL                  |  |  |  |  |
| CCB<br>280-613263/30   |                           | SM 2320B     |       | InitialAmount<br>is blank |                        |  |  |  |  |
| LCS<br>280-613263/31   |                           | SM 2320B     |       | InitialAmount<br>is blank | 10 mL                  |  |  |  |  |
| LCS D<br>280-613263/32 |                           | SM 2320B     |       | InitialAmount<br>is blank | 10 mL                  |  |  |  |  |
| MB<br>280-613263/33    |                           | SM 2320B     |       | InitialAmount<br>is blank |                        |  |  |  |  |
| 280-176674-A-7         | FBQmw-173-230401<br>-GW-R | SM 2320B     | T     | InitialAmount<br>is blank |                        |  |  |  |  |
| 280-176674-A-8         | FBQmw-173-230402<br>-GW-R | SM 2320B     | T     | InitialAmount<br>is blank |                        |  |  |  |  |
| CCV<br>280-613263/44   |                           | SM 2320B     |       | InitialAmount<br>is blank | 10 mL                  |  |  |  |  |
| CCB<br>280-613263/45   |                           | SM 2320B     |       | InitialAmount<br>is blank |                        |  |  |  |  |

| Batch Notes                |                      |
|----------------------------|----------------------|
| Nominal Amount Used        | 10 mL                |
| Acid ID                    | 0.02 H2SO4_00303     |
| Normality of First Titrant | 0.02 N               |
| Sodium Carbonate ID        | Alk Stk Std_00023    |
| pH Buffer 1 ID             | pH 2.0 buffer_00100  |
| pH Buffer 2 ID             | pH 4.0 buffer_00211  |
| pH Buffer 3 ID             | pH 7.0 buffer_00323  |
| pH Buffer 4 ID             | pH 7.0 buffer_00322  |
| pH Buffer 5 ID             | pH 10 buffer_00163   |
| pH Buffer 6 ID             | pH 12 buffer_00182   |
| Probe ID                   | PCE 80 pH 1200D_8918 |
| Batch Comment              | KG                   |

| Basis | Basis Description |
|-------|-------------------|
| T     | Total/NA          |

The pound sign (#) in the amount added field denotes that the reagent was used undiluted. All calculations are performed using the stated concentration for this reagent.

SM 2320B

Page 1 of 1




|   | SerialNumber | Cup position | Sample Type | Identity | Concentration | Corrected Height | Result | Use this                            | Y Residuals | Relative Error(%) |
|---|--------------|--------------|-------------|----------|---------------|------------------|--------|-------------------------------------|-------------|-------------------|
| 1 | 4            | ST1          | S1          | 0 mg/L   | 0.000         | -0.001           | -0.009 | <input checked="" type="checkbox"/> |             |                   |
| 2 | 5            | ST2          | S2          | 0.05     | 0.050         | 0.021            | 0.055  | <input checked="" type="checkbox"/> | 8.6%        | 10.5%             |
| 3 | 6            | ST3          | S3          | 0.1      | 0.100         | 0.038            | 0.104  | <input checked="" type="checkbox"/> | 3.91%       | 4.32%             |
| 4 | 7            | ST4          | S4          | 0.5      | 0.500         | 0.172            | 0.496  | <input checked="" type="checkbox"/> | -0.75%      | -0.76%            |
| 5 | 8            | ST5          | S5          | 1        | 1.000         | 0.342            | 1.005  | <input checked="" type="checkbox"/> | 0.45%       | 0.47%             |
| 6 | 9            | ST6          | S6          | 2.5      | 2.500         | 0.816            | 2.498  | <input checked="" type="checkbox"/> | -0.09%      | -0.09%            |
| 7 | 10           | ST7          | S7          | 5        | 5.000         | 1.524            | 5.000  | <input checked="" type="checkbox"/> | 0.01%       | 0.01%             |

|                         |                   |                         |             |              |
|-------------------------|-------------------|-------------------------|-------------|--------------|
| Method Name             | Ammonia           | a                       | 0.00199894  | R Squared    |
| Module Name             | NH3 and TKN       | b                       | 0.34728059  | Constant Sx0 |
| Calibration Order       | II Order ISO 8466 | c                       | -0.00856484 | Constant Vx0 |
| Residual Std. Dev. (Sy) | 0.00217970902050  | d                       |             |              |
|                         |                   | Correlation Coefficient | 0.99933221  |              |

Calibration Order


II Order ISO 8466

▼



# FlowAccessV3

Date: May 23 2023 5:29PM



$a = 0.0019989397283142$   $b = 0.3472805936030370$   $c = -0.0085648353783778$   $RSD = 0.00217970902050$   
 $r = 0.99933221456885$   $R^2 = 0.99999013246634$

Run Name : Skalar20230523A0

User Name : Administrator Operator Name : Administrator



## FlowAccessV3 Results Report

Run Name : Skalar20230523A0

DateTime :May 23 2023 8:05AM

User Name : Administrator Operator Name : Administrator

|    | Position | SampleType | SampleIdentity      | Comments | ExternalDilution |
|----|----------|------------|---------------------|----------|------------------|
| 1  | IW       | IW         | InitialWash         |          | 1.0000           |
| 2  | ST7      | T          | Tracer              |          | 1.0000           |
| 3  | ST5      | D          | Drift               |          | 1.0000           |
| 4  | WT       | W          | Wash                |          | 1.0000           |
| 5  | ST1      | S1         | 0 mg/L              |          | 1.0000           |
| 6  | ST2      | S2         | 0.05                |          | 1.0000           |
| 7  | ST3      | S3         | 0.1                 |          | 1.0000           |
| 8  | ST4      | S4         | 0.5                 |          | 1.0000           |
| 9  | ST5      | S5         | 1.0                 |          | 1.0000           |
| 10 | ST6      | S6         | 2.5                 |          | 1.0000           |
| 11 | ST7      | S7         | 5.0                 |          | 1.0000           |
| 12 | ST5      | D          | Drift               |          | 1.0000           |
| 13 | WT       | W          | Wash                |          | 1.0000           |
| 14 | ST8      | U          | ICV                 |          | 1.0000           |
| 15 | ST9      | U          | ICVL                |          | 1.0000           |
| 16 | ST1      | U          | ICB                 |          | 1.0000           |
| 17 | D29      | U          | LCS                 |          | 1.0000           |
| 18 | D30      | U          | LCSD                |          | 1.0000           |
| 19 | ST1      | U          | MB                  |          | 1.0000           |
| 20 | A1       | U          | 280-176003-a-23     |          | 1.0000           |
| 21 | A2       | U          | 280-176003-a-23 ms  |          | 1.0000           |
| 22 | A3       | U          | 280-176003-a-23 msd |          | 1.0000           |
| 23 | A4       | U          | 280-176003-a-5      |          | 1.0000           |
| 24 | A5       | U          | 280-176003-a-14     |          | 1.0000           |
| 25 | A6       | U          | 280-176114-a-12     |          | 1.0000           |
| 26 | A7       | U          | 280-176114-a-19     |          | 1.0000           |
| 27 | ST6      | U          | CCV                 |          | 1.0000           |
| 28 | ST4      | U          | CCVL                |          | 1.0000           |
| 29 | WT       | U          | CCB                 |          | 1.0000           |
| 30 | A8       | U          | 280-176117-a-9      |          | 1.0000           |
| 31 | A9       | U          | 280-176119-a-3      |          | 1.0000           |
| 32 | A10      | U          | 280-176185-g-3      |          | 5.0000           |
| 33 | A11      | U          | 280-176185-g-4      |          | 5.0000           |
| 34 | A12      | U          | 280-176185-g-5      |          | 1.0000           |
| 35 | A13      | U          | 280-176114-a-4      |          | 1.0000           |
| 36 | A14      | U          | 280-176114-a-4 ms   |          | 1.0000           |
| 37 | A15      | U          | 280-176114-a-4 msd  |          | 1.0000           |
| 38 | A16      | U          | 280-176185-g-6      |          | 1.0000           |
| 39 | A17      | U          | 280-176185-g-7      |          | 1.0000           |
| 40 | ST6      | U          | CCV                 |          | 1.0000           |

Page:1 /18

DateTime :May 23 2023 5:52PM

## FlowAccessV3 Results Report

Run Name : Skalar20230523A0

DateTime :May 23 2023 8:05AM

User Name : Administrator Operator Name : Administrator

|    | Position | SampleType | SampleIdentity      | Comments | ExternalDilution |
|----|----------|------------|---------------------|----------|------------------|
| 41 | ST4      | U          | CCVL                |          | 1.0000           |
| 42 | WT       | U          | CCB                 |          | 1.0000           |
| 43 | ST5      | D          | Drift               |          | 1.0000           |
| 44 | WT       | W          | Wash                |          | 1.0000           |
| 45 | A18      | U          | 280-176239-g-3      |          | 1.0000           |
| 46 | A19      | U          | 280-176239-g-4      |          | 1.0000           |
| 47 | A20      | U          | 280-176239-g-5      |          | 1.0000           |
| 48 | A21      | U          | 280-176239-g-6      |          | 1.0000           |
| 49 | A22      | U          | 280-176239-g-7      |          | 1.0000           |
| 50 | A23      | U          | 280-176208-b-2      |          | 1.0000           |
| 51 | A24      | U          | 280-176208-b-4      |          | 1.0000           |
| 52 | D29      | U          | LCS                 |          | 1.0000           |
| 53 | ST1      | U          | MB                  |          | 1.0000           |
| 54 | A25      | U          | 280-176123-a-36     |          | 1.0000           |
| 55 | ST6      | U          | CCV                 |          | 1.0000           |
| 56 | ST4      | U          | CCVL                |          | 1.0000           |
| 57 | WT       | U          | CCB                 |          | 1.0000           |
| 58 | A26      | U          | 280-176123-a-36 MS  |          | 1.0000           |
| 59 | A27      | U          | 280-176123-a-36 MSD |          | 1.0000           |
| 60 | A28      | U          | 280-176123-a-1      |          | 1.0000           |
| 61 | A29      | U          | 280-176123-a-2      |          | 1.0000           |
| 62 | A30      | U          | 280-176123-a-3      |          | 1.0000           |
| 63 | A31      | U          | 280-176123-a-4      |          | 1.0000           |
| 64 | A32      | U          | 280-176123-a-5      |          | 1.0000           |
| 65 | A33      | U          | 280-176229-b-2      |          | 1.0000           |
| 66 | A34      | U          | 280-176229-b-3      |          | 1.0000           |
| 67 | A35      | U          | 280-176229-b-5      |          | 1.0000           |
| 68 | ST6      | U          | CCV                 |          | 1.0000           |
| 69 | ST4      | U          | CCVL                |          | 1.0000           |
| 70 | WT       | U          | CCB                 |          | 1.0000           |
| 71 | ST5      | D          | Drift               |          | 1.0000           |
| 72 | WT       | W          | Wash                |          | 1.0000           |
| 73 | B1       | U          | 280-176229-b-6      |          | 1.0000           |
| 74 | B2       | U          | 280-176229-b-4      |          | 1.0000           |
| 75 | B3       | U          | 280-176229-b-4 MS   |          | 1.0000           |
| 76 | B4       | U          | 280-176229-b-4 MSD  |          | 1.0000           |
| 77 | B5       | U          | 280-176208-a-3      |          | 1.0000           |
| 78 | B6       | U          | 280-176229-b-9      |          | 1.0000           |
| 79 | B7       | U          | 280-176229-b-7      |          | 1.0000           |
| 80 | B8       | U          | 280-176229-b-8      |          | 1.0000           |

## FlowAccessV3 Results Report

Run Name : Skalar20230523A0

DateTime :May 23 2023 8:05AM

User Name : Administrator Operator Name : Administrator

|     | Position | SampleType | SampleIdentity     | Comments | ExternalDilution |
|-----|----------|------------|--------------------|----------|------------------|
| 81  | B9       | U          | 280-176213-j-1     |          | 20.0000          |
| 82  | B10      | U          | 280-176228-q-1     |          | 1.0000           |
| 83  | ST6      | U          | CCV                |          | 1.0000           |
| 84  | ST4      | U          | CCVL               |          | 1.0000           |
| 85  | WT       | U          | CCB                |          | 1.0000           |
| 86  | B11      | U          | 280-176193-a-1     |          | 1.0000           |
| 87  | B12      | U          | 280-176193-a-3     |          | 1.0000           |
| 88  | B13      | U          | 280-176193-a-4     |          | 1.0000           |
| 89  | D29      | U          | LCS                |          | 1.0000           |
| 90  | ST1      | U          | MB                 |          | 1.0000           |
| 91  | B14      | U          | 280-176193-a-2     |          | 1.0000           |
| 92  | B15      | U          | 280-176193-a-2 MS  |          | 1.0000           |
| 93  | B16      | U          | 280-176193-a-2 MSD |          | 1.0000           |
| 94  | B17      | U          | 280-176193-b-5     |          | 1.0000           |
| 95  | B18      | U          | 280-176271-b-1     |          | 1.0000           |
| 96  | ST6      | U          | CCV                |          | 1.0000           |
| 97  | ST4      | U          | CCVL               |          | 1.0000           |
| 98  | WT       | U          | CCB                |          | 1.0000           |
| 99  | ST5      | D          | Drift              |          | 1.0000           |
| 100 | WT       | W          | Wash               |          | 1.0000           |
| 101 | B19      | U          | 280-176271-b-2     |          | 1.0000           |
| 102 | B20      | U          | 280-176436-i-1     |          | 200.0000         |
| 103 | B21      | U          | 280-176274-d-2     |          | 100.0000         |
| 104 | B22      | U          | 590-20546-c-1      |          | 2.0000           |
| 105 | B23      | U          | 280-176302-c-2     |          | 1.0000           |
| 106 | B24      | U          | 280-176302-c-3     |          | 1.0000           |
| 107 | B25      | U          | 280-176302-c-4     |          | 1.0000           |
| 108 | B26      | U          | 280-176674-b-5     |          | 1.0000           |
| 109 | B27      | U          | 280-176674-b-5 MS  |          | 1.0000           |
| 110 | B28      | U          | 280-176674-b-5 MSD |          | 1.0000           |
| 111 | ST6      | U          | CCV                |          | 1.0000           |
| 112 | ST4      | U          | CCVL               |          | 1.0000           |
| 113 | WT       | U          | CCB                |          | 1.0000           |
| 114 | B29      | U          | 280-176302-c-5     |          | 1.0000           |
| 115 | B30      | U          | 280-176302-c-6     |          | 1.0000           |
| 116 | B31      | U          | 280-175902-q-1     |          | 5.0000           |
| 117 | B32      | U          | 280-176674-b-6     |          | 1.0000           |
| 118 | B33      | U          | 280-176674-b-9     |          | 1.0000           |
| 119 | B34      | U          | 590-20396-b-2      |          | 1.0000           |
| 120 | B35      | U          | 280-176674-b-3     |          | 1.0000           |

## FlowAccessV3 Results Report

Run Name : Skalar20230523A0

DateTime : May 23 2023 8:05AM

User Name : Administrator Operator Name : Administrator

|     | Position | SampleType | SampleIdentity     | Comments                       | ExternalDilution |
|-----|----------|------------|--------------------|--------------------------------|------------------|
| 121 | C1       | U          | 280-176472-a-1     |                                | 1.0000           |
| 122 | C2       | U          | 590-20500-d-2      |                                | 2.0000           |
| 123 | D29      | U          | LCS                |                                | 1.0000           |
| 124 | ST6      | U          | CCV                |                                | 1.0000           |
| 125 | ST4      | U          | CCVL               |                                | 1.0000           |
| 126 | WT       | U          | CCB                |                                | 1.0000           |
| 127 | ST5      | D          | Drift              |                                | 1.0000           |
| 128 | WT       | W          | Wash               |                                | 1.0000           |
| 129 | ST1      | U          | MB                 |                                | 1.0000           |
| 130 | C3       | U          | 280-176515-c-2     |                                | 1.0000           |
| 131 | C4       | U          | 280-176515-c-2 MS  |                                | 1.0000           |
| 132 | C5       | U          | 280-176515-c-2 MSD |                                | 1.0000           |
| 133 | C6       | U          | 590-20500-d-1      |                                | 2.0000           |
| 134 | C7       | U          | 280-176210-b-2     |                                | 20.0000          |
| 135 | C8       | U          | 280-176210-b-3     |                                | 20.0000          |
| 136 | C9       | U          | 280-176326-b-3     |                                | 1.0000           |
| 137 | C10      | U          | 280-176326-b-4     |                                | 1.0000           |
| 138 | B31      | U          | 280-175902-q-1     | rerun @ no dilution            | 1.0000           |
| 139 | ST6      | U          | CCV                |                                | 1.0000           |
| 140 | ST4      | U          | CCVL               |                                | 1.0000           |
| 141 | WT       | U          | CCB                |                                | 1.0000           |
| 142 | C2       | U          | 590-20500-d-2      | rerun @ no dilution            | 1.0000           |
| 143 | B35      | U          | 280-176674-b-3     | rerun @ a dilution             | 5.0000           |
| 144 | C1       | U          | 280-176472-a-1     | rerun @ a dilution             | 10.0000          |
| 145 | C7       | U          | 280-176210-b-2     | rerun @ a higher dilution      | 40.0000          |
| 146 | C11      | U          | 280-176326-a-5     |                                | 1.0000           |
| 147 | C12      | U          | 280-176363-c-1     |                                | 1.0000           |
| 148 | C13      | U          | 280-176363-a-10    |                                | 1.0000           |
| 149 | C14      | U          | 280-176363-C-11    |                                | 200.0000         |
| 150 | C15      | U          | 280-176365-c-4     |                                | 1.0000           |
| 151 | C16      | U          | 280-176365-c-4 MS  |                                | 1.0000           |
| 152 | ST6      | U          | CCV                |                                | 1.0000           |
| 153 | ST4      | U          | CCVL               |                                | 1.0000           |
| 154 | WT       | U          | CCB                |                                | 1.0000           |
| 155 | ST5      | D          | Drift              |                                | 1.0000           |
| 156 | WT       | W          | Wash               |                                | 1.0000           |
| 157 | C17      | U          | 280-176365-c-4 MSD |                                | 1.0000           |
| 158 | C1       | U          | 280-176472-a-1     | rerun @ 1x for elevated result | 1.0000           |
| 159 | B35      | U          | 280-176674-b-3     | rerun @ a higher dilution      | 100.0000         |
| 160 | C18      | U          | 280-176363-c-12    |                                | 1.0000           |

Page:4 /18

DateTime : May 23 2023 5:52PM



## FlowAccessV3 Results Report

Run Name : Skalar20230523A0

DateTime :May 23 2023 8:05AM

User Name : Administrator Operator Name : Administrator

|     | Position | SampleType | SampleIdentity     | Comments                  | ExternalDilution |
|-----|----------|------------|--------------------|---------------------------|------------------|
| 161 | C19      | U          | 280-176363-c-13    |                           | 100.0000         |
| 162 | C20      | U          | 280-176363-d-14    |                           | 1.0000           |
| 163 | C21      | U          | 280-176363-d-15    |                           | 50.0000          |
| 164 | C22      | U          | 280-176363-c-16    |                           | 1.0000           |
| 165 | C23      | U          | 280-176634-D-1     |                           | 1.0000           |
| 166 | C24      | U          | 280-176634-D-4     |                           | 1.0000           |
| 167 | ST6      | U          | CCV                |                           | 1.0000           |
| 168 | ST4      | U          | CCVL               |                           | 1.0000           |
| 169 | WT       | U          | CCB                |                           | 1.0000           |
| 170 | C25      | U          | 280-176634-d-5     |                           | 1.0000           |
| 171 | C26      | U          | 280-176634-d-7     |                           | 1.0000           |
| 172 | D29      | U          | LCS                |                           | 1.0000           |
| 173 | D30      | U          | LCSD               |                           | 1.0000           |
| 174 | ST1      | U          | MB                 |                           | 1.0000           |
| 175 | C27      | U          | 280-176296-g-3     |                           | 1.0000           |
| 176 | C28      | U          | 280-176296-g-3 MS  |                           | 1.0000           |
| 177 | C29      | U          | 280-176296-g-3 MSD |                           | 1.0000           |
| 178 | C30      | U          | 280-176296-g-4     |                           | 2.0000           |
| 179 | C31      | U          | 280-176296-g-5     |                           | 1.0000           |
| 180 | ST6      | U          | CCV                |                           | 1.0000           |
| 181 | ST4      | U          | CCVL               |                           | 1.0000           |
| 182 | WT       | U          | CCB                |                           | 1.0000           |
| 183 | ST5      | D          | Drift              |                           | 1.0000           |
| 184 | WT       | W          | Wash               |                           | 1.0000           |
| 185 | C32      | U          | 280-176296-g-7     |                           | 1.0000           |
| 186 | C33      | U          | 280-176296-g-8     |                           | 1.0000           |
| 187 | C34      | U          | 280-176296-g-9     |                           | 1.0000           |
| 188 | C35      | U          | 280-176572-f-13    |                           | 1.0000           |
| 189 | D1       | U          | 280-176634-d-8     |                           | 1.0000           |
| 190 | D2       | U          | 280-176634-d-10    |                           | 1.0000           |
| 191 | D3       | U          | 280-176634-d-11    |                           | 1.0000           |
| 192 | D4       | U          | 280-176853-a-7     |                           | 1.0000           |
| 193 | D5       | U          | 280-176853-a-7 MS  |                           | 1.0000           |
| 194 | D6       | U          | 280-176853-a-7 MSD |                           | 1.0000           |
| 195 | ST6      | U          | CCV                |                           | 1.0000           |
| 196 | ST4      | U          | CCVL               |                           | 1.0000           |
| 197 | WT       | U          | CCB                |                           | 1.0000           |
| 198 | B35      | U          | 280-176674-b-3     | rerun @ a higher dilution | 200.0000         |
| 199 | D29      | U          | LCS                | FILTRATION                | 1.0000           |
| 200 | ST1      | U          | MB                 | FILTRATION                | 1.0000           |

Page:5 /18

DateTime :May 23 2023 5:52PM

## FlowAccessV3 Results Report

Run Name : Skalar20230523A0

DateTime :May 23 2023 8:05AM

User Name : Administrator Operator Name : Administrator

|     | Position | SampleType | SampleIdentity     | Comments                           | ExternalDilution |
|-----|----------|------------|--------------------|------------------------------------|------------------|
| 201 | D8       | U          | 280-175902-s-1     | FILTRATION                         | 1.0000           |
| 202 | D9       | U          | 280-175902-s-1 MS  | FILTRATION                         | 1.0000           |
| 203 | D10      | U          | 280-175902-s-1 MSD | FILTRATION                         | 1.0000           |
| 204 | D11      | U          | 280-176589-e-4     |                                    | 1.0000           |
| 205 | D12      | U          | 280-176589-g-3     | upreserved                         | 1.0000           |
| 206 | D13      | U          | 280-176589-F-3     | TREATED chlorine/iodine/peroxide+  | 1.0000           |
| 207 | D14      | U          | 280-176589-f-3     | UNTREATED chlorine/iodine/peroxide | 1.0000           |
| 208 | WT       | W          | Wash               |                                    | 1.0000           |
| 209 | ST6      | U          | CCV                |                                    | 1.0000           |
| 210 | ST4      | U          | CCVL               |                                    | 1.0000           |
| 211 | WT       | U          | CCB                |                                    | 1.0000           |
| 212 | ST5      | D          | Drift              |                                    | 1.0000           |
| 213 | WT       | W          | Wash               |                                    | 1.0000           |
| 214 | E        | E          | EndRun             |                                    | 1.0000           |

## FlowAccessV3 Results Report

Run Name : Skalar20230523A0

DateTime :May 23 2023 8:05AM

User Name : Administrator Operator Name : Administrator

|    | NH3 and TKN- Results | NH3 and TKN- | NH3 and TKN-       |
|----|----------------------|--------------|--------------------|
| 1  | -0.006               | 0.0000       | May 23 2023 8:15AM |
| 2  | 5.107                | 1.5522       | May 23 2023 8:17AM |
| 3  | 1.025                | 0.3491       | May 23 2023 8:21AM |
| 4  | -0.006               | 0.0000       | May 23 2023 8:24AM |
| 5  | -0.009               | -0.0010      | May 23 2023 8:27AM |
| 6  | 0.055                | 0.0212       | May 23 2023 8:29AM |
| 7  | 0.104                | 0.0381       | May 23 2023 8:32AM |
| 8  | 0.496                | 0.1722       | May 23 2023 8:34AM |
| 9  | 1.005                | 0.3423       | May 23 2023 8:37AM |
| 10 | 2.498                | 0.8160       | May 23 2023 8:40AM |
| 11 | 5.000                | 1.5244       | May 23 2023 8:42AM |
| 12 | 1.043                | 0.3547       | May 23 2023 8:45AM |
| 13 | -0.006               | 0.0000       | May 23 2023 8:48AM |
| 14 | 2.487                | 0.8128       | May 23 2023 8:51AM |
| 15 | 0.491                | 0.1706       | May 23 2023 8:53AM |
| 16 | -0.011               | -0.0018      | May 23 2023 8:56AM |
| 17 | 2.524                | 0.8238       | May 23 2023 8:58AM |
| 18 | 2.546                | 0.8306       | May 23 2023 9:01AM |
| 19 | -0.001               | 0.0015       | May 23 2023 9:04AM |
| 20 | 0.056                | 0.0213       | May 23 2023 9:06AM |
| 21 | 1.079                | 0.3667       | May 23 2023 9:09AM |
| 22 | 1.123                | 0.3811       | May 23 2023 9:12AM |
| 23 | 0.003                | 0.0031       | May 23 2023 9:14AM |
| 24 | -0.012               | -0.0021      | May 23 2023 9:17AM |
| 25 | -0.009               | -0.0012      | May 23 2023 9:20AM |
| 26 | -0.008               | -0.0007      | May 23 2023 9:22AM |
| 27 | 2.570                | 0.8378       | May 23 2023 9:25AM |
| 28 | 0.512                | 0.1774       | May 23 2023 9:28AM |
| 29 | 0.000                | 0.0019       | May 23 2023 9:30AM |
| 30 | -0.010               | -0.0014      | May 23 2023 9:33AM |
| 31 | 0.047                | 0.0184       | May 23 2023 9:36AM |
| 32 | 6.198                | 0.4193       | May 23 2023 9:38AM |
| 33 | 9.719                | 0.6447       | May 23 2023 9:41AM |
| 34 | 2.588                | 0.8435       | May 23 2023 9:43AM |
| 35 | -0.003               | 0.0009       | May 23 2023 9:46AM |
| 36 | 1.041                | 0.3542       | May 23 2023 9:49AM |
| 37 | 1.008                | 0.3434       | May 23 2023 9:52AM |
| 38 | 0.005                | 0.0038       | May 23 2023 9:54AM |
| 39 | 0.020                | 0.0089       | May 23 2023 9:57AM |
| 40 | 2.597                | 0.8461       | May 23 2023 9:59AM |

## FlowAccessV3 Results Report

Run Name : Skalar20230523A0

DateTime :May 23 2023 8:05AM

User Name : Administrator Operator Name : Administrator

|    | NH3 and TKN- Results | NH3 and TKN- | NH3 and TKN-        |
|----|----------------------|--------------|---------------------|
| 41 | 0.517                | 0.1794       | May 23 2023 10:02AM |
| 42 | 0.000                | 0.0020       | May 23 2023 10:05AM |
| 43 | 1.041                | 0.3544       | May 23 2023 10:08AM |
| 44 | -0.006               | 0.0000       | May 23 2023 10:10AM |
| 45 | 0.003                | 0.0029       | May 23 2023 10:14AM |
| 46 | -0.004               | 0.0006       | May 23 2023 10:16AM |
| 47 | -0.004               | 0.0007       | May 23 2023 10:18AM |
| 48 | 0.144                | 0.0519       | May 23 2023 10:21AM |
| 49 | 0.003                | 0.0029       | May 23 2023 10:24AM |
| 50 | 0.404                | 0.1408       | May 23 2023 10:26AM |
| 51 | 0.568                | 0.1964       | May 23 2023 10:29AM |
| 52 | 2.568                | 0.8373       | May 23 2023 10:32AM |
| 53 | -0.001               | 0.0017       | May 23 2023 10:34AM |
| 54 | 0.072                | 0.0270       | May 23 2023 10:37AM |
| 55 | 2.575                | 0.8396       | May 23 2023 10:40AM |
| 56 | 0.514                | 0.1784       | May 23 2023 10:42AM |
| 57 | 0.004                | 0.0033       | May 23 2023 10:45AM |
| 58 | 1.128                | 0.3830       | May 23 2023 10:48AM |
| 59 | 1.158                | 0.3928       | May 23 2023 10:50AM |
| 60 | 0.007                | 0.0045       | May 23 2023 10:53AM |
| 61 | 0.011                | 0.0057       | May 23 2023 10:56AM |
| 62 | 0.012                | 0.0061       | May 23 2023 10:58AM |
| 63 | 0.455                | 0.1582       | May 23 2023 11:01AM |
| 64 | 0.005                | 0.0038       | May 23 2023 11:04AM |
| 65 | 0.003                | 0.0031       | May 23 2023 11:06AM |
| 66 | 1.655                | 0.5533       | May 23 2023 11:09AM |
| 67 | 0.175                | 0.0625       | May 23 2023 11:12AM |
| 68 | 2.554                | 0.8332       | May 23 2023 11:14AM |
| 69 | 0.525                | 0.1820       | May 23 2023 11:17AM |
| 70 | 0.007                | 0.0046       | May 23 2023 11:20AM |
| 71 | 1.063                | 0.3615       | May 23 2023 11:22AM |
| 72 | -0.006               | 0.0000       | May 23 2023 11:25AM |
| 73 | 0.119                | 0.0432       | May 23 2023 11:28AM |
| 74 | 0.095                | 0.0348       | May 23 2023 11:30AM |
| 75 | 0.996                | 0.3395       | May 23 2023 11:33AM |
| 76 | 1.051                | 0.3575       | May 23 2023 11:36AM |
| 77 | 0.338                | 0.1184       | May 23 2023 11:38AM |
| 78 | 0.009                | 0.0051       | May 23 2023 11:41AM |
| 79 | 1.769                | 0.5897       | May 23 2023 11:44AM |
| 80 | 1.096                | 0.3724       | May 23 2023 11:46AM |

## FlowAccessV3 Results Report

Run Name : Skalar20230523A0

DateTime :May 23 2023 8:05AM

User Name : Administrator Operator Name : Administrator

|     | NH3 and TKN- Results | NH3 and TKN- | NH3 and TKN-        |
|-----|----------------------|--------------|---------------------|
| 81  | 56.952               | 0.9215       | May 23 2023 11:49AM |
| 82  | 0.048                | 0.0187       | May 23 2023 11:52AM |
| 83  | 2.575                | 0.8395       | May 23 2023 11:54AM |
| 84  | 0.520                | 0.1802       | May 23 2023 11:57AM |
| 85  | 0.007                | 0.0043       | May 23 2023 12:00PM |
| 86  | 0.034                | 0.0137       | May 23 2023 12:02PM |
| 87  | 0.182                | 0.0650       | May 23 2023 12:05PM |
| 88  | 0.030                | 0.0125       | May 23 2023 12:08PM |
| 89  | 2.575                | 0.8394       | May 23 2023 12:10PM |
| 90  | 0.004                | 0.0034       | May 23 2023 12:13PM |
| 91  | 0.100                | 0.0366       | May 23 2023 12:16PM |
| 92  | 1.140                | 0.3869       | May 23 2023 12:18PM |
| 93  | 1.173                | 0.3975       | May 23 2023 12:21PM |
| 94  | 0.028                | 0.0118       | May 23 2023 12:24PM |
| 95  | 0.138                | 0.0499       | May 23 2023 12:26PM |
| 96  | 2.471                | 0.8078       | May 23 2023 12:29PM |
| 97  | 0.512                | 0.1777       | May 23 2023 12:32PM |
| 98  | -0.002               | 0.0013       | May 23 2023 12:34PM |
| 99  | 1.073                | 0.3648       | May 23 2023 12:37PM |
| 100 | -0.006               | 0.0000       | May 23 2023 12:40PM |
| 101 | 0.029                | 0.0122       | May 23 2023 12:43PM |
| 102 | 289.772              | 0.4872       | May 23 2023 12:45PM |
| 103 | 270.798              | 0.8796       | May 23 2023 12:48PM |
| 104 | 2.930                | 0.4923       | May 23 2023 12:50PM |
| 105 | -0.004               | 0.0005       | May 23 2023 12:53PM |
| 106 | -0.010               | -0.0013      | May 23 2023 12:56PM |
| 107 | -0.001               | 0.0016       | May 23 2023 12:58PM |
| 108 | 0.000                | 0.0020       | May 23 2023 1:01PM  |
| 109 | 1.045                | 0.3555       | May 23 2023 1:04PM  |
| 110 | 1.034                | 0.3520       | May 23 2023 1:06PM  |
| 111 | 2.553                | 0.8329       | May 23 2023 1:09PM  |
| 112 | 0.494                | 0.1716       | May 23 2023 1:12PM  |
| 113 | -0.028               | -0.0076      | May 23 2023 1:14PM  |
| 114 | 0.046                | 0.0180       | May 23 2023 1:17PM  |
| 115 | -0.023               | -0.0060      | May 23 2023 1:20PM  |
| 116 | 1.257                | 0.0888       | May 23 2023 1:22PM  |
| 117 | -0.014               | -0.0029      | May 23 2023 1:25PM  |
| 118 | 0.003                | 0.0031       | May 23 2023 1:28PM  |
| 119 | -0.044               | -0.0134      | May 23 2023 1:30PM  |
| 120 | 13.946               | 3.1794       | May 23 2023 1:33PM  |

## FlowAccessV3 Results Report

Run Name : Skalar20230523A0

DateTime :May 23 2023 8:05AM

User Name : Administrator Operator Name : Administrator

|     | NH3 and TKN- Results | NH3 and TKN- | NH3 and TKN-       |
|-----|----------------------|--------------|--------------------|
| 121 | 0.971                | 0.3312       | May 23 2023 1:37PM |
| 122 | 0.752                | 0.1313       | May 23 2023 1:39PM |
| 123 | 2.728                | 0.8856       | May 23 2023 1:43PM |
| 124 | 2.736                | 0.8881       | May 23 2023 1:46PM |
| 125 | 0.551                | 0.1906       | May 23 2023 1:49PM |
| 126 | 0.011                | 0.0059       | May 23 2023 1:51PM |
| 127 | 1.075                | 0.3655       | May 23 2023 1:54PM |
| 128 | -0.006               | 0.0000       | May 23 2023 1:57PM |
| 129 | -0.019               | -0.0046      | May 23 2023 1:59PM |
| 130 | 0.046                | 0.0178       | May 23 2023 2:02PM |
| 131 | 1.037                | 0.3528       | May 23 2023 2:05PM |
| 132 | 1.054                | 0.3586       | May 23 2023 2:08PM |
| 133 | 0.028                | 0.0069       | May 23 2023 2:10PM |
| 134 | 113.581              | 1.6980       | May 23 2023 2:13PM |
| 135 | 11.457               | 0.1981       | May 23 2023 2:15PM |
| 136 | -0.003               | 0.0010       | May 23 2023 2:18PM |
| 137 | -0.015               | -0.0032      | May 23 2023 2:21PM |
| 138 | 1.624                | 0.5433       | May 23 2023 2:24PM |
| 139 | 2.573                | 0.8387       | May 23 2023 2:26PM |
| 140 | 0.506                | 0.1754       | May 23 2023 2:29PM |
| 141 | -0.016               | -0.0036      | May 23 2023 2:32PM |
| 142 | 0.055                | 0.0209       | May 23 2023 2:34PM |
| 143 | 66.259               | 3.1000       | May 23 2023 2:37PM |
| 144 | 3.572                | 0.1250       | May 23 2023 2:40PM |
| 145 | 124.186              | 0.9976       | May 23 2023 2:42PM |
| 146 | 0.089                | 0.0329       | May 23 2023 2:45PM |
| 147 | 1.000                | 0.3408       | May 23 2023 2:48PM |
| 148 | 2.337                | 0.7667       | May 23 2023 2:50PM |
| 149 | 170.694              | 0.2922       | May 23 2023 2:53PM |
| 150 | 0.124                | 0.0448       | May 23 2023 2:56PM |
| 151 | 1.181                | 0.4002       | May 23 2023 2:58PM |
| 152 | 2.524                | 0.8241       | May 23 2023 3:01PM |
| 153 | 0.513                | 0.1779       | May 23 2023 3:04PM |
| 154 | -0.002               | 0.0013       | May 23 2023 3:06PM |
| 155 | 1.106                | 0.3755       | May 23 2023 3:09PM |
| 156 | -0.006               | 0.0000       | May 23 2023 3:12PM |
| 157 | 1.192                | 0.4037       | May 23 2023 3:15PM |
| 158 | 0.787                | 0.2701       | May 23 2023 3:17PM |
| 159 | 656.181              | 1.9120       | May 23 2023 3:20PM |
| 160 | 1.073                | 0.3647       | May 23 2023 3:22PM |

## FlowAccessV3 Results Report

Run Name : Skalar20230523A0

DateTime :May 23 2023 8:05AM

User Name : Administrator Operator Name : Administrator

|     | NH3 and TKN- Results | NH3 and TKN- | NH3 and TKN-       |
|-----|----------------------|--------------|--------------------|
| 161 | 238.901              | 0.7828       | May 23 2023 3:25PM |
| 162 | 0.829                | 0.2840       | May 23 2023 3:28PM |
| 163 | 68.651               | 0.4627       | May 23 2023 3:30PM |
| 164 | 0.908                | 0.3102       | May 23 2023 3:33PM |
| 165 | 0.150                | 0.0539       | May 23 2023 3:35PM |
| 166 | 1.326                | 0.4475       | May 23 2023 3:38PM |
| 167 | 2.595                | 0.8455       | May 23 2023 3:41PM |
| 168 | 0.525                | 0.1818       | May 23 2023 3:44PM |
| 169 | 0.000                | 0.0018       | May 23 2023 3:46PM |
| 170 | 0.367                | 0.1283       | May 23 2023 3:49PM |
| 171 | 0.325                | 0.1141       | May 23 2023 3:51PM |
| 172 | 2.632                | 0.8569       | May 23 2023 3:54PM |
| 173 | 2.618                | 0.8525       | May 23 2023 3:57PM |
| 174 | 0.004                | 0.0034       | May 23 2023 4:00PM |
| 175 | 0.621                | 0.2144       | May 23 2023 4:02PM |
| 176 | 1.766                | 0.5886       | May 23 2023 4:05PM |
| 177 | 1.663                | 0.5558       | May 23 2023 4:08PM |
| 178 | 8.030                | 1.2583       | May 23 2023 4:10PM |
| 179 | 2.189                | 0.7212       | May 23 2023 4:13PM |
| 180 | 2.643                | 0.8599       | May 23 2023 4:16PM |
| 181 | 0.527                | 0.1826       | May 23 2023 4:18PM |
| 182 | -0.007               | -0.0005      | May 23 2023 4:21PM |
| 183 | 1.070                | 0.3639       | May 23 2023 4:24PM |
| 184 | -0.006               | 0.0000       | May 23 2023 4:26PM |
| 185 | 1.697                | 0.5666       | May 23 2023 4:30PM |
| 186 | -0.012               | -0.0021      | May 23 2023 4:32PM |
| 187 | -0.004               | 0.0007       | May 23 2023 4:34PM |
| 188 | 0.634                | 0.2188       | May 23 2023 4:37PM |
| 189 | 1.131                | 0.3839       | May 23 2023 4:39PM |
| 190 | 0.275                | 0.0968       | May 23 2023 4:42PM |
| 191 | 0.604                | 0.2087       | May 23 2023 4:45PM |
| 192 | -0.028               | -0.0077      | May 23 2023 4:48PM |
| 193 | 1.073                | 0.3648       | May 23 2023 4:50PM |
| 194 | 1.037                | 0.3529       | May 23 2023 4:53PM |
| 195 | 2.555                | 0.8335       | May 23 2023 4:56PM |
| 196 | 0.479                | 0.1665       | May 23 2023 4:58PM |
| 197 | -0.051               | -0.0156      | May 23 2023 5:01PM |
| 198 | 729.998              | 1.1555       | May 23 2023 5:04PM |
| 199 | 2.877                | 0.9302       | May 23 2023 5:06PM |
| 200 | -0.057               | -0.0178      | May 23 2023 5:09PM |

Page:11 /18

DateTime :May 23 2023 5:52PM

## FlowAccessV3 Results Report

Run Name : Skalar20230523A0

DateTime :May 23 2023 8:05AM

User Name : Administrator Operator Name : Administrator

|     | NH3 and TKN- Results | NH3 and TKN- | NH3 and TKN-       |
|-----|----------------------|--------------|--------------------|
| 201 | 1.168                | 0.3959       | May 23 2023 5:12PM |
| 202 | 2.192                | 0.7221       | May 23 2023 5:14PM |
| 203 | 2.252                | 0.7407       | May 23 2023 5:17PM |
| 204 | -0.062               | -0.0196      | May 23 2023 5:19PM |
| 205 | 12.812               | 3.0455       | May 23 2023 5:22PM |
| 206 | 12.798               | 3.0437       | May 23 2023 5:25PM |
| 207 | 12.983               | 3.0671       | May 23 2023 5:28PM |
| 208 | -0.006               | 0.0000       | May 23 2023 5:30PM |
| 209 | 2.543                | 0.8297       | May 23 2023 5:33PM |
| 210 | 0.477                | 0.1657       | May 23 2023 5:36PM |
| 211 | -0.035               | -0.0101      | May 23 2023 5:38PM |
| 212 | 1.076                | 0.3659       | May 23 2023 5:41PM |
| 213 | -0.006               | 0.0000       | May 23 2023 5:43PM |
| 214 | -0.006               | 0.0000       | May 23 2023 5:47PM |



## FlowAccessV3 Results Report

Run Name : Skalar20230523A0

DateTime :May 23 2023 8:05AM

User Name : Administrator Operator Name : Administrator

|    | NH3 and TKN- |
|----|--------------|
| 1  | -0.006       |
| 2  | 5.107        |
| 3  | 1.025        |
| 4  | -0.006       |
| 5  | -0.009       |
| 6  | 0.055        |
| 7  | 0.104        |
| 8  | 0.496        |
| 9  | 1.005        |
| 10 | 2.498        |
| 11 | 5.000        |
| 12 | 1.043        |
| 13 | -0.006       |
| 14 | 2.487        |
| 15 | 0.491        |
| 16 | -0.011       |
| 17 | 2.524        |
| 18 | 2.546        |
| 19 | -0.001       |
| 20 | 0.056        |
| 21 | 1.079        |
| 22 | 1.123        |
| 23 | 0.003        |
| 24 | -0.012       |
| 25 | -0.009       |
| 26 | -0.008       |
| 27 | 2.570        |
| 28 | 0.512        |
| 29 | 0.000        |
| 30 | -0.010       |
| 31 | 0.047        |
| 32 | 1.240        |
| 33 | 1.944        |
| 34 | 2.588        |
| 35 | -0.003       |
| 36 | 1.041        |
| 37 | 1.008        |
| 38 | 0.005        |
| 39 | 0.020        |
| 40 | 2.597        |

## FlowAccessV3 Results Report

Run Name : Skalar20230523A0

DateTime :May 23 2023 8:05AM

User Name : Administrator Operator Name : Administrator

|    | NH3 and TKN- |
|----|--------------|
| 41 | 0.517        |
| 42 | 0.000        |
| 43 | 1.041        |
| 44 | -0.006       |
| 45 | 0.003        |
| 46 | -0.004       |
| 47 | -0.004       |
| 48 | 0.144        |
| 49 | 0.003        |
| 50 | 0.404        |
| 51 | 0.568        |
| 52 | 2.568        |
| 53 | -0.001       |
| 54 | 0.072        |
| 55 | 2.575        |
| 56 | 0.514        |
| 57 | 0.004        |
| 58 | 1.128        |
| 59 | 1.158        |
| 60 | 0.007        |
| 61 | 0.011        |
| 62 | 0.012        |
| 63 | 0.455        |
| 64 | 0.005        |
| 65 | 0.003        |
| 66 | 1.655        |
| 67 | 0.175        |
| 68 | 2.554        |
| 69 | 0.525        |
| 70 | 0.007        |
| 71 | 1.063        |
| 72 | -0.006       |
| 73 | 0.119        |
| 74 | 0.095        |
| 75 | 0.996        |
| 76 | 1.051        |
| 77 | 0.338        |
| 78 | 0.009        |
| 79 | 1.769        |
| 80 | 1.096        |

## FlowAccessV3 Results Report

Run Name : Skalar20230523A0

DateTime :May 23 2023 8:05AM

User Name : Administrator Operator Name : Administrator

|     | NH3 and TKN- |
|-----|--------------|
| 81  | 2.848        |
| 82  | 0.048        |
| 83  | 2.575        |
| 84  | 0.520        |
| 85  | 0.007        |
| 86  | 0.034        |
| 87  | 0.182        |
| 88  | 0.030        |
| 89  | 2.575        |
| 90  | 0.004        |
| 91  | 0.100        |
| 92  | 1.140        |
| 93  | 1.173        |
| 94  | 0.028        |
| 95  | 0.138        |
| 96  | 2.471        |
| 97  | 0.512        |
| 98  | -0.002       |
| 99  | 1.073        |
| 100 | -0.006       |
| 101 | 0.029        |
| 102 | 1.449        |
| 103 | 2.708        |
| 104 | 1.465        |
| 105 | -0.004       |
| 106 | -0.010       |
| 107 | -0.001       |
| 108 | 0.000        |
| 109 | 1.045        |
| 110 | 1.034        |
| 111 | 2.553        |
| 112 | 0.494        |
| 113 | -0.028       |
| 114 | 0.046        |
| 115 | -0.023       |
| 116 | 0.251        |
| 117 | -0.014       |
| 118 | 0.003        |
| 119 | -0.044       |
| 120 | 13.946       |

## FlowAccessV3 Results Report

Run Name : Skalar20230523A0

DateTime :May 23 2023 8:05AM

User Name : Administrator Operator Name : Administrator

|     | NH3 and TKN- |
|-----|--------------|
| 121 | 0.971        |
| 122 | 0.376        |
| 123 | 2.728        |
| 124 | 2.736        |
| 125 | 0.551        |
| 126 | 0.011        |
| 127 | 1.075        |
| 128 | -0.006       |
| 129 | -0.019       |
| 130 | 0.046        |
| 131 | 1.037        |
| 132 | 1.054        |
| 133 | 0.014        |
| 134 | 5.679        |
| 135 | 0.573        |
| 136 | -0.003       |
| 137 | -0.015       |
| 138 | 1.624        |
| 139 | 2.573        |
| 140 | 0.506        |
| 141 | -0.016       |
| 142 | 0.055        |
| 143 | 13.252       |
| 144 | 0.357        |
| 145 | 3.105        |
| 146 | 0.089        |
| 147 | 1.000        |
| 148 | 2.337        |
| 149 | 0.853        |
| 150 | 0.124        |
| 151 | 1.181        |
| 152 | 2.524        |
| 153 | 0.513        |
| 154 | -0.002       |
| 155 | 1.106        |
| 156 | -0.006       |
| 157 | 1.192        |
| 158 | 0.787        |
| 159 | 6.562        |
| 160 | 1.073        |

## FlowAccessV3 Results Report

Run Name : Skalar20230523A0

DateTime :May 23 2023 8:05AM

User Name : Administrator Operator Name : Administrator

|     | NH3 and TKN- |
|-----|--------------|
| 161 | 2.389        |
| 162 | 0.829        |
| 163 | 1.373        |
| 164 | 0.908        |
| 165 | 0.150        |
| 166 | 1.326        |
| 167 | 2.595        |
| 168 | 0.525        |
| 169 | 0.000        |
| 170 | 0.367        |
| 171 | 0.325        |
| 172 | 2.632        |
| 173 | 2.618        |
| 174 | 0.004        |
| 175 | 0.621        |
| 176 | 1.766        |
| 177 | 1.663        |
| 178 | 4.015        |
| 179 | 2.189        |
| 180 | 2.643        |
| 181 | 0.527        |
| 182 | -0.007       |
| 183 | 1.070        |
| 184 | -0.006       |
| 185 | 1.697        |
| 186 | -0.012       |
| 187 | -0.004       |
| 188 | 0.634        |
| 189 | 1.131        |
| 190 | 0.275        |
| 191 | 0.604        |
| 192 | -0.028       |
| 193 | 1.073        |
| 194 | 1.037        |
| 195 | 2.555        |
| 196 | 0.479        |
| 197 | -0.051       |
| 198 | 3.650        |
| 199 | 2.877        |
| 200 | -0.057       |

# FlowAccessV3 Results Report

Run Name : Skalar20230523A0

DateTime :May 23 2023 8:05AM

User Name : Administrator Operator Name : Administrator

|     | NH3 and TKN- |
|-----|--------------|
| 201 | 1.168        |
| 202 | 2.192        |
| 203 | 2.252        |
| 204 | -0.062       |
| 205 | 12.812       |
| 206 | 12.798       |
| 207 | 12.983       |
| 208 | -0.006       |
| 209 | 2.543        |
| 210 | 0.477        |
| 211 | -0.035       |
| 212 | 1.076        |
| 213 | -0.006       |
| 214 | -0.006       |

# General Chemistry Raw Data Report

Job ID: 280-176674-1

Batch: 613208

Method: 9034

Analyst Initials: CAI

Instrument: No Equipment Used for this Test

Lab Sample ID: LCS 280-613197/1-A

Analysis Date: May 19, 2023 17:27

| Analyte | Detector | Dilution | Raw Result | Unit | Initial Amount | Final Amount |
|---------|----------|----------|------------|------|----------------|--------------|
| Sulfide | None     | 1        | 16         | mg/L | 50 mL          | 50 mL        |

Lab Sample ID: MB 280-613197/2-A

Analysis Date: May 19, 2023 17:27

| Analyte | Detector | Dilution | Raw Result | Unit | Initial Amount | Final Amount |
|---------|----------|----------|------------|------|----------------|--------------|
| Sulfide | None     | 1        | 0          | mg/L | 50 mL          | 50 mL        |

Lab Sample ID: 280-176674-E-7-A

Analysis Date: May 19, 2023 17:27

| Analyte | Detector | Dilution | Raw Result | Unit | Initial Amount | Final Amount |
|---------|----------|----------|------------|------|----------------|--------------|
| Sulfide | None     | 1        | 0          | mg/L | 50 mL          | 50 mL        |

Lab Sample ID: 280-176674-B-8-A

Analysis Date: May 19, 2023 17:27

| Analyte | Detector | Dilution | Raw Result | Unit | Initial Amount | Final Amount |
|---------|----------|----------|------------|------|----------------|--------------|
| Sulfide | None     | 1        | 0.8        | mg/L | 50 mL          | 50 mL        |

# Sulfide by Titration

|                     |                    |                         |           |
|---------------------|--------------------|-------------------------|-----------|
| Analyst:            | CAI                | SOP Information:        |           |
| Date:               | 5/19/2023          | Number:                 | 91        |
| Titration Solutions |                    | Calibration Information |           |
| Solution 1:         | Iodine             | Source/Ver-Lot#:        | INT_02386 |
| TALS ID             | Iod_00279          | Prep Date:              | 4/26/2023 |
| Normality:          | 0.1                | Made By:                | CAI       |
| Solution 2:         | sodium thiosulfate | Concentration:          | 952       |
| TALS ID             | Na Thio_00181      | Expiration Date:        | 7/26/2023 |
| Normality:          | 0.025              |                         |           |
|                     | Starch Indicator   |                         |           |
| TALS ID             | Starch IND_00071   |                         |           |

| ICV Information  |           |
|------------------|-----------|
| Source/Ver-Lot#: | INT_02291 |
| Prep Date:       | 4/26/2023 |
| Made By:         | CAI       |
| Concentration:   | 920       |
| Expiration Date: | 7/26/2023 |
|                  |           |
|                  |           |
|                  |           |

|     | CAL    | Buret | Buret | mL     | Final | Conc    |
|-----|--------|-------|-------|--------|-------|---------|
|     | Volume | Start | End   | Iodine | mL    | mg/L    |
| CAL | 5      | 0.00  | 8.20  | 5      | 8.20  | 944.000 |
| CAL | 5      | 8.20  | 16.20 | 5      | 8.00  | 960.000 |

|     | CAL    | Buret | Buret | mL     | Final | Conc    |
|-----|--------|-------|-------|--------|-------|---------|
|     | Volume | Start | End   | Iodine | mL    | mg/L    |
| ICV | 5      | 0.00  | 8.50  | 5      | 8.50  | 920.000 |
| ICV | 5      | 8.50  | 17.00 | 5      | 8.50  | 920.000 |



# TALS Raw Data Report

Job Number: 280-176524-1  
LIMS Batch: 613208  
Equipment: NOEQUIP

Laboratory: Eurofins Denver

|       |                              |                               |               |                      |       |         |       |         |
|-------|------------------------------|-------------------------------|---------------|----------------------|-------|---------|-------|---------|
| RS# 1 | Lab ID: LCS 280-613197/1-A   | Inj Date: 5/19/2023 5:27:00PM | Dil: 1.0      | Meth: 9034           |       |         |       |         |
|       | Analyte                      | Rspnse                        | Raw Res/Units | Final Res/Qual/Units | % Rec | Rec Lmt | % RPD | RPD Lmt |
|       | Sulfide                      |                               | 16 mg/L       | mg/L                 | 84    | 44 110  |       |         |
|       | Sulfide as H2S               |                               | 17 mg/L       | mg/L                 |       |         |       |         |
| RS# 2 | Lab ID: MB 280-613197/2-A    | Inj Date: 5/19/2023 5:27:00PM | Dil: 1.0      | Meth: 9034           |       |         |       |         |
|       | Analyte                      | Rspnse                        | Raw Res/Units | Final Res/Qual/Units | % Rec | Rec Lmt | % RPD | RPD Lmt |
|       | Sulfide                      |                               | 0 mg/L        | 2.0 U mg/L           |       |         |       |         |
|       | Sulfide as H2S               |                               | 0 mg/L        | 2.1 U mg/L           |       |         |       |         |
| RS# 3 | Lab ID: 280-176524-E-1-A     | Inj Date: 5/19/2023 5:27:00PM | Dil: 1.0      | Meth: 9034           |       |         |       |         |
|       | Analyte                      | Rspnse                        | Raw Res/Units | Final Res/Qual/Units | % Rec | Rec Lmt | % RPD | RPD Lmt |
|       | Sulfide                      |                               | 0 mg/L        | 2.0 U mg/L           |       |         |       |         |
|       | Sulfide as H2S               |                               | 0 mg/L        | 2.1 U mg/L           |       |         |       |         |
| RS# 4 | Lab ID: 280-176524-E-1-B MS  | Inj Date: 5/19/2023 5:27:00PM | Dil: 1.0      | Meth: 9034           |       |         |       |         |
|       | Analyte                      | Rspnse                        | Raw Res/Units | Final Res/Qual/Units | % Rec | Rec Lmt | % RPD | RPD Lmt |
|       | Sulfide                      |                               | 15.2 mg/L     | mg/L                 | 80    | 44 110  |       |         |
|       | Sulfide as H2S               |                               | 16.15 mg/L    | mg/L                 |       |         |       |         |
| RS# 5 | Lab ID: 280-176524-E-1-C MSD | Inj Date: 5/19/2023 5:27:00PM | Dil: 1.0      | Meth: 9034           |       |         |       |         |
|       | Analyte                      | Rspnse                        | Raw Res/Units | Final Res/Qual/Units | % Rec | Rec Lmt | % RPD | RPD Lmt |
|       | Sulfide                      |                               | 15.2 mg/L     | mg/L                 | 80    | 44 110  | 0     | 20      |
|       | Sulfide as H2S               |                               | 16.15 mg/L    | mg/L                 |       |         |       |         |
| RS# 6 | Lab ID: 280-176524-E-4-A     | Inj Date: 5/19/2023 5:27:00PM | Dil: 1.0      | Meth: 9034           |       |         |       |         |
|       | Analyte                      | Rspnse                        | Raw Res/Units | Final Res/Qual/Units | % Rec | Rec Lmt | % RPD | RPD Lmt |
|       | Sulfide                      |                               | 0 mg/L        | 2.0 U mg/L           |       |         |       |         |
|       | Sulfide as H2S               |                               | 0 mg/L        | 2.1 U mg/L           |       |         |       |         |

# TALS Raw Data Report

Job Number: 280-176526-1  
LIMS Batch: 613208  
Equipment: NOEQUIP

Laboratory: Eurofins Denver

|        |                              |                               |               |                      |       |          |            |         |
|--------|------------------------------|-------------------------------|---------------|----------------------|-------|----------|------------|---------|
| RS# 1  | Lab ID: LCS 280-613197/1-A   | Inj Date: 5/19/2023 5:27:00PM |               |                      |       | Dil: 1.0 | Meth: 9034 |         |
|        | Analyte                      | Rspnse                        | Raw Res/Units | Final Res/Qual/Units | % Rec | Rec Lmt  | % RPD      | RPD Lmt |
|        | Sulfide                      |                               | 16 mg/L       | mg/L                 | 84    | 44 110   |            |         |
|        | Sulfide as H2S               |                               | 17 mg/L       | mg/L                 |       |          |            |         |
| RS# 2  | Lab ID: MB 280-613197/2-A    | Inj Date: 5/19/2023 5:27:00PM |               |                      |       | Dil: 1.0 | Meth: 9034 |         |
|        | Analyte                      | Rspnse                        | Raw Res/Units | Final Res/Qual/Units | % Rec | Rec Lmt  | % RPD      | RPD Lmt |
|        | Sulfide                      |                               | 0 mg/L        | mg/L                 |       |          |            |         |
|        | Sulfide as H2S               |                               | 0 mg/L        | mg/L                 |       |          |            |         |
| RS# 3  | Lab ID: 280-176524-E-1-A     | Inj Date: 5/19/2023 5:27:00PM |               |                      |       | Dil: 1.0 | Meth: 9034 |         |
|        | Analyte                      | Rspnse                        | Raw Res/Units | Final Res/Qual/Units | % Rec | Rec Lmt  | % RPD      | RPD Lmt |
|        | Sulfide                      |                               | 0 mg/L        | mg/L                 |       |          |            |         |
|        | Sulfide as H2S               |                               | 0 mg/L        | mg/L                 |       |          |            |         |
| RS# 4  | Lab ID: 280-176524-E-1-B MS  | Inj Date: 5/19/2023 5:27:00PM |               |                      |       | Dil: 1.0 | Meth: 9034 |         |
|        | Analyte                      | Rspnse                        | Raw Res/Units | Final Res/Qual/Units | % Rec | Rec Lmt  | % RPD      | RPD Lmt |
|        | Sulfide                      |                               | 15.2 mg/L     | mg/L                 | 80    | 44 110   |            |         |
|        | Sulfide as H2S               |                               | 16.15 mg/L    | mg/L                 |       |          |            |         |
| RS# 5  | Lab ID: 280-176524-E-1-C MSD | Inj Date: 5/19/2023 5:27:00PM |               |                      |       | Dil: 1.0 | Meth: 9034 |         |
|        | Analyte                      | Rspnse                        | Raw Res/Units | Final Res/Qual/Units | % Rec | Rec Lmt  | % RPD      | RPD Lmt |
|        | Sulfide                      |                               | 15.2 mg/L     | mg/L                 | 80    | 44 110   | 0          | 20      |
|        | Sulfide as H2S               |                               | 16.15 mg/L    | mg/L                 |       |          |            |         |
| RS# 13 | Lab ID: 280-176526-O-1-A     | Inj Date: 5/19/2023 5:27:00PM |               |                      |       | Dil: 1.0 | Meth: 9034 |         |
|        | Analyte                      | Rspnse                        | Raw Res/Units | Final Res/Qual/Units | % Rec | Rec Lmt  | % RPD      | RPD Lmt |
|        | Sulfide                      |                               | 0 mg/L        | mg/L                 |       |          |            |         |
|        | Sulfide as H2S               |                               | 0 mg/L        | mg/L                 |       |          |            |         |

# TALS Raw Data Report

Job Number: 280-176591-1  
LIMS Batch: 613208  
Equipment: NOEQUIP

Laboratory: Eurofins Denver

|        |                              |                               |               |                      |       |         |       |         |
|--------|------------------------------|-------------------------------|---------------|----------------------|-------|---------|-------|---------|
| RS# 1  | Lab ID: LCS 280-613197/1-A   | Inj Date: 5/19/2023 5:27:00PM | Dil: 1.0      | Meth: 9034           |       |         |       |         |
|        | Analyte                      | Rspnse                        | Raw Res/Units | Final Res/Qual/Units | % Rec | Rec Lmt | % RPD | RPD Lmt |
|        | Sulfide                      |                               | 16 mg/L       | mg/L                 | 84    | 44      | 110   |         |
|        | Sulfide as H2S               |                               | 17 mg/L       | mg/L                 |       |         |       |         |
| RS# 2  | Lab ID: MB 280-613197/2-A    | Inj Date: 5/19/2023 5:27:00PM | Dil: 1.0      | Meth: 9034           |       |         |       |         |
|        | Analyte                      | Rspnse                        | Raw Res/Units | Final Res/Qual/Units | % Rec | Rec Lmt | % RPD | RPD Lmt |
|        | Sulfide                      |                               | 0 mg/L        | 2.0 U mg/L           |       |         |       |         |
|        | Sulfide as H2S               |                               | 0 mg/L        | 2.1 U mg/L           |       |         |       |         |
| RS# 3  | Lab ID: 280-176524-E-1-A     | Inj Date: 5/19/2023 5:27:00PM | Dil: 1.0      | Meth: 9034           |       |         |       |         |
|        | Analyte                      | Rspnse                        | Raw Res/Units | Final Res/Qual/Units | % Rec | Rec Lmt | % RPD | RPD Lmt |
|        | Sulfide                      |                               | 0 mg/L        | 2.0 U mg/L           |       |         |       |         |
|        | Sulfide as H2S               |                               | 0 mg/L        | 2.1 U mg/L           |       |         |       |         |
| RS# 4  | Lab ID: 280-176524-E-1-B MS  | Inj Date: 5/19/2023 5:27:00PM | Dil: 1.0      | Meth: 9034           |       |         |       |         |
|        | Analyte                      | Rspnse                        | Raw Res/Units | Final Res/Qual/Units | % Rec | Rec Lmt | % RPD | RPD Lmt |
|        | Sulfide                      |                               | 15.2 mg/L     | mg/L                 | 80    | 44      | 110   |         |
|        | Sulfide as H2S               |                               | 16.15 mg/L    | mg/L                 |       |         |       |         |
| RS# 5  | Lab ID: 280-176524-E-1-C MSD | Inj Date: 5/19/2023 5:27:00PM | Dil: 1.0      | Meth: 9034           |       |         |       |         |
|        | Analyte                      | Rspnse                        | Raw Res/Units | Final Res/Qual/Units | % Rec | Rec Lmt | % RPD | RPD Lmt |
|        | Sulfide                      |                               | 15.2 mg/L     | mg/L                 | 80    | 44      | 110   | 0 20    |
|        | Sulfide as H2S               |                               | 16.15 mg/L    | mg/L                 |       |         |       |         |
| RS# 7  | Lab ID: 280-176591-E-2-A     | Inj Date: 5/19/2023 5:27:00PM | Dil: 1.0      | Meth: 9034           |       |         |       |         |
|        | Analyte                      | Rspnse                        | Raw Res/Units | Final Res/Qual/Units | % Rec | Rec Lmt | % RPD | RPD Lmt |
|        | Sulfide                      |                               | 0 mg/L        | 2.0 U mg/L           |       |         |       |         |
|        | Sulfide as H2S               |                               | 0 mg/L        | 2.1 U mg/L           |       |         |       |         |
| RS# 8  | Lab ID: 280-176591-E-4-A     | Inj Date: 5/19/2023 5:27:00PM | Dil: 1.0      | Meth: 9034           |       |         |       |         |
|        | Analyte                      | Rspnse                        | Raw Res/Units | Final Res/Qual/Units | % Rec | Rec Lmt | % RPD | RPD Lmt |
|        | Sulfide                      |                               | 0.8 mg/L      | 2.0 U mg/L           |       |         |       |         |
|        | Sulfide as H2S               |                               | 0.85 mg/L     | 2.1 U mg/L           |       |         |       |         |
| RS# 9  | Lab ID: 280-176591-E-6-A     | Inj Date: 5/19/2023 5:27:00PM | Dil: 1.0      | Meth: 9034           |       |         |       |         |
|        | Analyte                      | Rspnse                        | Raw Res/Units | Final Res/Qual/Units | % Rec | Rec Lmt | % RPD | RPD Lmt |
|        | Sulfide                      |                               | 0 mg/L        | 2.0 U mg/L           |       |         |       |         |
|        | Sulfide as H2S               |                               | 0 mg/L        | 2.1 U mg/L           |       |         |       |         |
| RS# 10 | Lab ID: 280-176591-E-7-A     | Inj Date: 5/19/2023 5:27:00PM | Dil: 1.0      | Meth: 9034           |       |         |       |         |
|        | Analyte                      | Rspnse                        | Raw Res/Units | Final Res/Qual/Units | % Rec | Rec Lmt | % RPD | RPD Lmt |
|        | Sulfide                      |                               | 0 mg/L        | 2.0 U mg/L           |       |         |       |         |
|        | Sulfide as H2S               |                               | 0 mg/L        | 2.1 U mg/L           |       |         |       |         |

# TALS Raw Data Report

Job Number: 280-176674-1  
LIMS Batch: 613208  
Equipment: NOEQUIP

Laboratory: Eurofins Denver

|        |                              |                               |               |                      |       |          |       |            |  |
|--------|------------------------------|-------------------------------|---------------|----------------------|-------|----------|-------|------------|--|
| RS# 1  | Lab ID: LCS 280-613197/1-A   | Inj Date: 5/19/2023 5:27:00PM |               |                      |       | Dil: 1.0 |       | Meth: 9034 |  |
|        | Analyte                      | Rspnse                        | Raw Res/Units | Final Res/Qual/Units | % Rec | Rec Lmt  | % RPD | RPD Lmt    |  |
|        | Sulfide                      |                               | 16 mg/L       | mg/L                 | 84    | 44 110   |       |            |  |
|        | Sulfide as H2S               |                               | 17 mg/L       | mg/L                 |       |          |       |            |  |
| RS# 2  | Lab ID: MB 280-613197/2-A    | Inj Date: 5/19/2023 5:27:00PM |               |                      |       | Dil: 1.0 |       | Meth: 9034 |  |
|        | Analyte                      | Rspnse                        | Raw Res/Units | Final Res/Qual/Units | % Rec | Rec Lmt  | % RPD | RPD Lmt    |  |
|        | Sulfide                      |                               | 0 mg/L        | 2.0 U mg/L           |       |          |       |            |  |
|        | Sulfide as H2S               |                               | 0 mg/L        | 2.1 U mg/L           |       |          |       |            |  |
| RS# 3  | Lab ID: 280-176524-E-1-A     | Inj Date: 5/19/2023 5:27:00PM |               |                      |       | Dil: 1.0 |       | Meth: 9034 |  |
|        | Analyte                      | Rspnse                        | Raw Res/Units | Final Res/Qual/Units | % Rec | Rec Lmt  | % RPD | RPD Lmt    |  |
|        | Sulfide                      |                               | 0 mg/L        | 2.0 U mg/L           |       |          |       |            |  |
|        | Sulfide as H2S               |                               | 0 mg/L        | 2.1 U mg/L           |       |          |       |            |  |
| RS# 4  | Lab ID: 280-176524-E-1-B MS  | Inj Date: 5/19/2023 5:27:00PM |               |                      |       | Dil: 1.0 |       | Meth: 9034 |  |
|        | Analyte                      | Rspnse                        | Raw Res/Units | Final Res/Qual/Units | % Rec | Rec Lmt  | % RPD | RPD Lmt    |  |
|        | Sulfide                      |                               | 15.2 mg/L     | mg/L                 | 80    | 44 110   |       |            |  |
|        | Sulfide as H2S               |                               | 16.15 mg/L    | mg/L                 |       |          |       |            |  |
| RS# 5  | Lab ID: 280-176524-E-1-C MSD | Inj Date: 5/19/2023 5:27:00PM |               |                      |       | Dil: 1.0 |       | Meth: 9034 |  |
|        | Analyte                      | Rspnse                        | Raw Res/Units | Final Res/Qual/Units | % Rec | Rec Lmt  | % RPD | RPD Lmt    |  |
|        | Sulfide                      |                               | 15.2 mg/L     | mg/L                 | 80    | 44 110   | 0     | 20         |  |
|        | Sulfide as H2S               |                               | 16.15 mg/L    | mg/L                 |       |          |       |            |  |
| RS# 11 | Lab ID: 280-176674-E-7-A     | Inj Date: 5/19/2023 5:27:00PM |               |                      |       | Dil: 1.0 |       | Meth: 9034 |  |
|        | Analyte                      | Rspnse                        | Raw Res/Units | Final Res/Qual/Units | % Rec | Rec Lmt  | % RPD | RPD Lmt    |  |
|        | Sulfide                      |                               | 0 mg/L        | 2.0 U mg/L           |       |          |       |            |  |
|        | Sulfide as H2S               |                               | 0 mg/L        | 2.1 U mg/L           |       |          |       |            |  |
| RS# 12 | Lab ID: 280-176674-B-8-A     | Inj Date: 5/19/2023 5:27:00PM |               |                      |       | Dil: 1.0 |       | Meth: 9034 |  |
|        | Analyte                      | Rspnse                        | Raw Res/Units | Final Res/Qual/Units | % Rec | Rec Lmt  | % RPD | RPD Lmt    |  |
|        | Sulfide                      |                               | 0.8 mg/L      | 2.0 U mg/L           |       |          |       |            |  |
|        | Sulfide as H2S               |                               | 0.85 mg/L     | 2.1 U mg/L           |       |          |       |            |  |



Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-121594.b  
 Lims ID: STD L2  
 Client ID:  
 Sample Type: IC Calib Level: 2  
 Inject. Date: 18-May-2023 12:11:00 ALS Bottle#: 0 Worklist Smp#: 2  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121594-002  
 Misc. Info.: 280-0121594-002  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Sublist: chrom-Anions\_IC10\*sub5  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
 Limit Group: Wet - Anions  
 Last Update: 19-May-2023 11:53:46 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-121594.b

Column 1 : Det: Info 2\_091554\_1

Process Host: CTX1603

First Level Reviewer: LVW8

Date: 18-May-2023 12:58:14

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 1 Fluoride            | 3.133        | 3.133            | 0.000            | 5835197  | NC               | NC                 | Ma    |
| 2 Chloride            | 4.300        | 4.300            | 0.000            | 20218551 | NC               | NC                 |       |
| 3 Nitrite as N        | 5.002        | 5.002            | 0.000            | 7062392  | 0.2000           | 0.1873             |       |
| 4 Bromide             | 6.165        | 6.165            | 0.000            | 2255477  | NC               | NC                 |       |
| 5 Nitrate as N        | 6.953        | 6.953            | 0.000            | 7411099  | 0.2000           | 0.1905             |       |
| 6 Orthophosphate as P | 9.093        | 9.093            | 0.000            | 8660779  | 0.2000           | 0.1157             |       |
| 7 Sulfate             | 10.245       | 10.245           | 0.000            | 20026015 | NC               | NC                 |       |

### QC Flag Legend

Processing Flags

NC - Not Calibrated

Review Flags

M - Manually Integrated

a - User Assigned ID

### Reagents:

IC CAL c/so4\_00480 Amount Added: 0.10 Units: mL  
 IC Cal low\_00709 Amount Added: 0.08 Units: mL

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-122558.d

Injection Date: 18-May-2023 12:11:00

Instrument ID: WC\_IonChrom10

Operator ID: wetchemd

Lims ID: STD L2

Worklist Smp#: 2

Client ID:


Injection Vol: 5.0 ul

Dil. Factor: 1.0000

ALS Bottle#: 0

Method: Anions\_IC10

Limit Group: Wet - Anions



Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-121594.b  
 Lims ID: STD L3  
 Client ID:  
 Sample Type: IC Calib Level: 3  
 Inject. Date: 18-May-2023 12:26:00 ALS Bottle#: 0 Worklist Smp#: 3  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121594-003  
 Misc. Info.: 280-0121594-003  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Sublist: chrom-Anions\_IC10\*sub5  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
 Limit Group: Wet - Anions  
 Last Update: 19-May-2023 11:53:47 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-121594.b  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1603

First Level Reviewer: LVW8

Date: 18-May-2023 12:58:34

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 1 Fluoride            | 3.137        | 3.137            | 0.000            | 18179755 | NC               | NC                 | Ma    |
| 2 Chloride            | 4.300        | 4.300            | 0.000            | 48348641 | NC               | NC                 |       |
| 3 Nitrite as N        | 5.000        | 5.000            | 0.000            | 21399863 | 0.5000           | 0.5218             |       |
| 4 Bromide             | 6.162        | 6.162            | 0.000            | 4768597  | NC               | NC                 |       |
| 5 Nitrate as N        | 6.938        | 6.938            | 0.000            | 22979043 | 0.5000           | 0.5365             |       |
| 6 Orthophosphate as P | 9.075        | 9.075            | 0.000            | 20921591 | 0.5000           | 0.6419             |       |
| 7 Sulfate             | 10.247       | 10.247           | 0.000            | 43450426 | NC               | NC                 |       |

**QC Flag Legend**

Processing Flags

NC - Not Calibrated

Review Flags

M - Manually Integrated

a - User Assigned ID

**Reagents:**

IC CAL c/so4\_00480

Amount Added: 0.20

Units: mL

IC Cal low\_00709

Amount Added: 0.20

Units: mL



Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-124056.d

Injection Date: 18-May-2023 12:26:00

Instrument ID: WC\_IonChrom10

Operator ID: wetchemd

Lims ID: STD L3

Worklist Smp#: 3

Client ID:


Injection Vol: 5.0 ul

Dil. Factor: 1.0000

ALS Bottle#: 0

Method: Anions\_IC10

Limit Group: Wet - Anions



Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-121594.b  
 Lims ID: STD L4  
 Client ID:  
 Sample Type: IC Calib Level: 4  
 Inject. Date: 18-May-2023 12:40:00 ALS Bottle#: 0 Worklist Smp#: 4  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121594-004  
 Misc. Info.: 280-0121594-004  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Sublist: chrom-Anions\_IC10\*sub5  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
 Limit Group: Wet - Anions  
 Last Update: 19-May-2023 11:53:47 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-121594.b

Column 1 : Det: Info 2\_091554\_1

Process Host: CTX1603

First Level Reviewer: LVW8

Date: 18-May-2023 13:29:45

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response  | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|-----------|------------------|--------------------|-------|
| 1 Fluoride            | 3.138        | 3.138            | 0.000            | 39638142  | NC               | NC                 | Ma    |
| 2 Chloride            | 4.305        | 4.305            | 0.000            | 570289992 | NC               | NC                 |       |
| 3 Nitrite as N        | 5.000        | 5.000            | 0.000            | 43146865  | 1.00             | 1.03               |       |
| 4 Bromide             | 6.162        | 6.162            | 0.000            | 7900278   | NC               | NC                 |       |
| 5 Nitrate as N        | 6.928        | 6.928            | 0.000            | 42854846  | 1.00             | 0.9783             |       |
| 6 Orthophosphate as P | 9.068        | 9.068            | 0.000            | 33022328  | 1.00             | 1.16               |       |
| 7 Sulfate             | 10.242       | 10.242           | 0.000            | 432139991 | NC               | NC                 |       |

### QC Flag Legend

Processing Flags

NC - Not Calibrated

Review Flags

M - Manually Integrated

a - User Assigned ID

### Reagents:

IC CAL c/so4\_00480

Amount Added: 2.40

Units: mL


IC Cal low\_00709

Amount Added: 0.40

Units: mL

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-125554.d  
Injection Date: 18-May-2023 12:40:00 Instrument ID: WC\_IonChrom10 Operator ID: wetchemd  
Lims ID: STD L4 Worklist Smp#: 4  
Client ID:  
Injection Vol: 5.0 ul Dil. Factor: 1.0000 ALS Bottle#: 0  
Method: Anions\_IC10 Limit Group: Wet - Anions



Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13  
 Lims ID: STD L5  
 Client ID:  
 Sample Type: IC Calib Level: 5  
 Inject. Date: 18-May-2023 12:55:00 ALS Bottle#: 0 Worklist Smp#: 5  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121594-005  
 Misc. Info.: 280-0121594-005  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Sublist: chrom-Anions\_IC10\*sub5  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
 Limit Group: Wet - Anions  
 Last Update: 19-May-2023 11:53:48 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1603

First Level Reviewer: LVW8

Date: 18-May-2023 13:30:00

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response   | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|------------|------------------|--------------------|-------|
| 1 Fluoride            | 3.140        | 3.140            | 0.000            | 102461401  | NC               | NC                 | Ma    |
| 2 Chloride            | 4.313        | 4.313            | 0.000            | 1123523133 | NC               | NC                 |       |
| 3 Nitrite as N        | 4.995        | 4.995            | 0.000            | 105836223  | 2.50             | 2.49               |       |
| 4 Bromide             | 6.152        | 6.152            | 0.000            | 18980516   | NC               | NC                 |       |
| 5 Nitrate as N        | 6.903        | 6.903            | 0.000            | 110713949  | 2.50             | 2.49               |       |
| 6 Orthophosphate as P | 9.043        | 9.043            | 0.000            | 66560323   | 2.50             | 2.60               |       |
| 7 Sulfate             | 10.227       | 10.227           | 0.000            | 831159688  | NC               | NC                 |       |

**QC Flag Legend**

Processing Flags

NC - Not Calibrated

Review Flags

M - Manually Integrated

a - User Assigned ID

**Reagents:**

IC CAL c/so4\_00480

Amount Added: 4.80

Units: mL


IC Cal low\_00709

Amount Added: 1.00

Units: mL

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-131053.d  
Injection Date: 18-May-2023 12:55:00 Instrument ID: WC\_IonChrom10 Operator ID: wetchemd  
Lims ID: STD L5 Worklist Smp#: 5  
Client ID:  
Injection Vol: 5.0 ul Dil. Factor: 1.0000 ALS Bottle#: 0  
Method: Anions\_IC10 Limit Group: Wet - Anions



Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13  
 Lims ID: STD L6  
 Client ID:  
 Sample Type: IC Calib Level: 6  
 Inject. Date: 18-May-2023 13:10:00 ALS Bottle#: 0 Worklist Smp#: 6  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121594-006  
 Misc. Info.: 280-0121594-006  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Sublist: chrom-Anions\_IC10\*sub5  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
 Limit Group: Wet - Anions  
 Last Update: 19-May-2023 11:53:49 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1603

First Level Reviewer: LVW8

Date: 18-May-2023 13:42:47

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response   | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|------------|------------------|--------------------|-------|
| 1 Fluoride            | 3.145        | 3.145            | 0.000            | 200510278  | NC               | NC                 | Ma    |
| 2 Chloride            | 4.333        | 4.333            | 0.000            | 1864000362 | NC               | NC                 |       |
| 3 Nitrite as N        | 4.997        | 4.997            | 0.000            | 212086804  | 5.00             | 4.97               |       |
| 4 Bromide             | 6.150        | 6.150            | 0.000            | 38244681   | NC               | NC                 |       |
| 5 Nitrate as N        | 6.890        | 6.890            | 0.000            | 224161262  | 5.00             | 5.01               |       |
| 6 Orthophosphate as P | 9.035        | 9.035            | 0.000            | 115018456  | 5.00             | 4.68               |       |
| 7 Sulfate             | 10.220       | 10.220           | 0.000            | 1371227441 | NC               | NC                 |       |

**QC Flag Legend**

Processing Flags

NC - Not Calibrated

Review Flags

M - Manually Integrated

a - User Assigned ID

**Reagents:**

IC CAL c/so4\_00480

Amount Added: 8.00

Units: mL

IC Cal low\_00709

Amount Added: 2.00

Units: mL

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-132551.d

Injection Date: 18-May-2023 13:10:00

Instrument ID: WC\_IonChrom10

Operator ID: wetchemd

Lims ID: STD L6

Worklist Smp#: 6

Client ID:


Injection Vol: 5.0 ul

Dil. Factor: 1.0000

ALS Bottle#: 0

Method: Anions\_IC10

Limit Group: Wet - Anions



# IC Instrument Information

WL: 121594 Inst ID: 1010 Analysis Date: 5/18/23 Analyst: MEC

| Rush                     | Job No.       | Samples   | Anions                                                     | QC Req | HT Exp |
|--------------------------|---------------|-----------|------------------------------------------------------------|--------|--------|
| <input type="checkbox"/> | <u>176637</u> | <u>1</u>  | F Cl NO2 Br <u>NO3</u> PO4 SO4                             | MS/D   | _____  |
| <input type="checkbox"/> | <u>176589</u> | <u>1</u>  | F Cl <u>NO2</u> Br NO3 PO4 SO4                             | MS/D   | _____  |
| <input type="checkbox"/> | <u>176670</u> | <u>7</u>  | F <u>Cl</u> NO2 Br <u>NO3</u> PO4 <u>SO4</u>               | MS/D   | _____  |
| <input type="checkbox"/> | <u>176676</u> | <u>15</u> | <u>F</u> <u>Cl</u> <u>NO2</u> Br <u>NO3</u> PO4 <u>SO4</u> | MS/D   | _____  |
| <input type="checkbox"/> | <u>176674</u> | <u>4</u>  | F Cl <u>NO2</u> Br <u>NO3</u> PO4 <u>SO4</u>               | MS/D   | _____  |
| <input type="checkbox"/> | <u>176678</u> | <u>5</u>  | F <u>Cl</u> <u>NO2</u> Br <u>NO3</u> PO4 <u>SO4</u>        | MS/D   | _____  |
| <input type="checkbox"/> | <u>176669</u> | <u>4</u>  | F <u>Cl</u> NO2 Br <u>NO3</u> PO4 <u>SO4</u>               | MS/D   | _____  |
| <input type="checkbox"/> | <u>176683</u> | <u>3</u>  | F <u>Cl</u> NO2 Br <u>NO3</u> PO4 <u>SO4</u>               | MS/D   | _____  |
| <input type="checkbox"/> | <u>176668</u> | <u>1</u>  | F Cl NO2 Br <u>NO3</u> PO4 SO4                             | MS/D   | _____  |
| <input type="checkbox"/> | _____         | _____     | F Cl NO2 Br NO3 PO4 SO4                                    | MS/D   | _____  |
| <input type="checkbox"/> | _____         | _____     | F Cl NO2 Br NO3 PO4 SO4                                    | MS/D   | _____  |
| <input type="checkbox"/> | _____         | _____     | F Cl NO2 Br NO3 PO4 SO4                                    | MS/D   | _____  |
| <input type="checkbox"/> | _____         | _____     | F Cl NO2 Br NO3 PO4 SO4                                    | MS/D   | _____  |
| <input type="checkbox"/> | _____         | _____     | F Cl NO2 Br NO3 PO4 SO4                                    | MS/D   | _____  |

## Dilutions

| Job No. | Samples | Anions                  | Dilution | Reason |
|---------|---------|-------------------------|----------|--------|
| _____   | _____   | F Cl NO2 Br NO3 PO4 SO4 | _____    | _____  |
| _____   | _____   | F Cl NO2 Br NO3 PO4 SO4 | _____    | _____  |
| _____   | _____   | F Cl NO2 Br NO3 PO4 SO4 | _____    | _____  |
| _____   | _____   | F Cl NO2 Br NO3 PO4 SO4 | _____    | _____  |
| _____   | _____   | F Cl NO2 Br NO3 PO4 SO4 | _____    | _____  |
| _____   | _____   | F Cl NO2 Br NO3 PO4 SO4 | _____    | _____  |
| _____   | _____   | F Cl NO2 Br NO3 PO4 SO4 | _____    | _____  |
| _____   | _____   | F Cl NO2 Br NO3 PO4 SO4 | _____    | _____  |
| _____   | _____   | F Cl NO2 Br NO3 PO4 SO4 | _____    | _____  |



Eurofins Environment Testing America  
Initial Calibration Report

Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230519-121638.b\Anions\_IC10.m

Instrument: WC\_IonChrom10

Lims Location: 280

Lock State: Initial Calib Locked

Cpnd Order: Retention Time

Integrator: Falcon

Last Modified: 19-May-2023 12:31:06

No. Compounds: 7

Sublist: chrom-Anions\_IC10\*sub5

Limit Group: Wet - Anions

## Detectors

Detector: 1, Info 2\_091554\_1

Data Type: ic Spec Type: none

Supports Extracted Chromatograms: False

Run Time: 0.000-14.900

No. Points: 7196

## Calibration File Names

|          |                                                                                                  |
|----------|--------------------------------------------------------------------------------------------------|
| Level: 1 | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230518-12110 |
|          | Inj Date: 18-May-2023 11:56:00 Worklist: 121594 Sample#: 1                                       |
| Level: 2 | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230518-12255 |
|          | Inj Date: 18-May-2023 12:11:00 Worklist: 121594 Sample#: 2                                       |
| Level: 3 | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230518-12405 |
|          | Inj Date: 18-May-2023 12:26:00 Worklist: 121594 Sample#: 3                                       |
| Level: 4 | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230518-12555 |
|          | Inj Date: 18-May-2023 12:40:00 Worklist: 121594 Sample#: 4                                       |
| Level: 5 | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230518-13105 |
|          | Inj Date: 18-May-2023 12:55:00 Worklist: 121594 Sample#: 5                                       |
| Level: 6 | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230518-13255 |
|          | Inj Date: 18-May-2023 13:10:00 Worklist: 121594 Sample#: 6                                       |

Start Cal Date: 18-May-2023 11:56:00

End Cal Date: 19-May-2023 12:01:00

## Individual Compound Calibration Parameters

Quant Method: ESTD

RF Calibration: Replace

Rule Name: Linear1

Curve: Linear

Weighting: Conc

Origin: None

Error: raw\_COD

Error Limit: 1.00

RF %Dif: 0.0

SPCC Limit: 0.0

CCC Limit: 0.0

Dependent Variable: Resp

Number of Compounds: 3

| RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | b | m1 | m2 | Error |
|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|---|----|----|-------|
|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|---|----|----|-------|

| RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | b | m1 | m2 | Error |
|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|---|----|----|-------|
|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|---|----|----|-------|

## 3 Nitrite as N

Signal: 1

|             |             |            |            |            |  |         |          |  |       |
|-------------|-------------|------------|------------|------------|--|---------|----------|--|-------|
| 35311960    | 42799726    | 43146865   | 42334489   | 42417361   |  | -967573 |          |  | 1.000 |
| 0.200000(2) | 0.500000(3) | 1.0000 (4) | 2.5000 (5) | 5.0000 (6) |  |         | 42866306 |  |       |
| 7062392     | 21399863    | 43146865   | 105836223  | 212086804  |  |         |          |  |       |
| 121594(2)   | 121594(3)   | 121594(4)  | 121594(5)  | 121594(6)  |  |         |          |  |       |

## 5 Nitrate as N

Signal: 1

|             |             |            |            |            |  |          |          |  |       |
|-------------|-------------|------------|------------|------------|--|----------|----------|--|-------|
| 37055495    | 45958086    | 42854846   | 44285580   | 44832252   |  | -1159043 |          |  | 0.999 |
| 0.200000(2) | 0.500000(3) | 1.0000 (4) | 2.5000 (5) | 5.0000 (6) |  |          | 44990806 |  |       |
| 7411099     | 22979043    | 42854846   | 110713949  | 224161262  |  |          |          |  |       |
| 121594(2)   | 121594(3)   | 121594(4)  | 121594(5)  | 121594(6)  |  |          |          |  |       |

## 6 Orthophosphate as P

Signal: 1

|             |             |            |            |            |  |         |          |  |        |
|-------------|-------------|------------|------------|------------|--|---------|----------|--|--------|
| 43303895    | 41843182    | 33022328   | 26624129   | 23003691   |  | 5964839 | R1       |  | *0.980 |
| 0.200000(2) | 0.500000(3) | 1.0000 (4) | 2.5000 (5) | 5.0000 (6) |  |         | 23299922 |  |        |
| 8660779     | 20921591    | 33022328   | 66560323   | 115018456  |  |         |          |  |        |
| 121594(2)   | 121594(3)   | 121594(4)  | 121594(5)  | 121594(6)  |  |         |          |  |        |

ICalib Error Legend

R1, Curve Fit Fail Error Limit Test

Eurofins Environment Testing America  
Initial Calibration Report

Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230519-121638.b\Anions\_IC10.m

Instrument: WC\_IonChrom10

Lims Location: 280

Lock State: Initial Calib Locked

Cpnd Order: Retention Time

Integrator: Falcon

Last Modified: 19-May-2023 12:31:06

No.Compounds:7

Sublist: chrom-Anions\_IC10\*sub5

Limit Group: Wet - Anions 28D

## Detectors

Detector: 1, Info 2\_091554\_1

Data Type: ic Spec Type: none

Supports Extracted Chromatograms: False

Run Time: 0.000-14.900

No. Points: 7196

## Calibration File Names

|          |                                                                                                  |
|----------|--------------------------------------------------------------------------------------------------|
| Level: 1 | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230518-12110 |
|          | Inj Date: 18-May-2023 11:56:00 Worklist: 121594 Sample#: 1                                       |
| Level: 2 | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230518-12255 |
|          | Inj Date: 18-May-2023 12:11:00 Worklist: 121594 Sample#: 2                                       |
| Level: 3 | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230518-12405 |
|          | Inj Date: 18-May-2023 12:26:00 Worklist: 121594 Sample#: 3                                       |
| Level: 4 | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230518-12555 |
|          | Inj Date: 18-May-2023 12:40:00 Worklist: 121594 Sample#: 4                                       |
| Level: 5 | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230518-13105 |
|          | Inj Date: 18-May-2023 12:55:00 Worklist: 121594 Sample#: 5                                       |
| Level: 6 | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230518-13255 |
|          | Inj Date: 18-May-2023 13:10:00 Worklist: 121594 Sample#: 6                                       |

Start Cal Date: 18-May-2023 11:56:00

End Cal Date: 19-May-2023 12:01:00

## Individual Compound Calibration Parameters

Quant Method: ESTD

RF Calibration: Replace

Rule Name: Linear1

Curve: Linear

Weighting: Conc

Origin: None

Error: raw\_COD

Error Limit: 1.00

RF %Dif: 0.0

SPCC Limit: 0.0

CCC Limit: 0.0

Dependent Variable: Resp

Number of Compounds: 4

| RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | b | m1 | m2 | Error |
|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|---|----|----|-------|
|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|---|----|----|-------|

| RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | b         | m1       | m2 | Error |
|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|-----------|----------|----|-------|
| 1 Fluoride                         |                                    |                                    |                                    |                                    |                                    | Signal: 1 |          |    |       |
| 29175985                           | 36359510                           | 39638142                           | 40984560                           | 40102056                           |                                    | -2208602  |          |    | 1.000 |
| 0.200000(2)                        | 0.500000(3)                        | 1.0000 (4)                         | 2.5000 (5)                         | 5.0000 (6)                         |                                    |           | 41050846 |    |       |
| M5835197                           | M18179755                          | M39638142                          | M102461401                         | M200510278                         |                                    |           |          |    |       |
| 121594(2)                          | 121594(3)                          | 121594(4)                          | 121594(5)                          | 121594(6)                          |                                    |           |          |    |       |
| 2 Chloride                         |                                    |                                    |                                    |                                    |                                    | Signal: 1 |          |    |       |
| 16780722                           | 16174841                           | 19339456                           | 19009666                           | 18725386                           | 18640004                           | -1171132  |          |    | 1.000 |
| 0.500000(1)                        | 1.2500 (2)                         | 2.5000 (3)                         | 30.0 (4)                           | 60.0 (5)                           | 100.0 (6)                          |           | 18747994 |    |       |
| 8390361                            | 20218551                           | 48348641                           | 570289992                          | 1123523133                         | 1864000362                         |           |          |    |       |
| 121594(1)                          | 121594(2)                          | 121594(3)                          | 121594(4)                          | 121594(5)                          | 121594(6)                          |           |          |    |       |
| 4 Bromide                          |                                    |                                    |                                    |                                    |                                    | Signal: 1 |          |    |       |
| 11277385                           | 9537194                            | 7900278                            | 7592206                            | 7648936                            |                                    | 806494    |          |    | 0.999 |
| 0.200000(2)                        | 0.500000(3)                        | 1.0000 (4)                         | 2.5000 (5)                         | 5.0000 (6)                         |                                    |           | 7404030  |    |       |
| 2255477                            | 4768597                            | 7900278                            | 18980516                           | 38244681                           |                                    |           |          |    |       |
| 121594(2)                          | 121594(3)                          | 121594(4)                          | 121594(5)                          | 121594(6)                          |                                    |           |          |    |       |
| 7 Sulfate                          |                                    |                                    |                                    |                                    |                                    | Signal: 1 |          |    |       |
| 22625290                           | 16020812                           | 17380170                           | 14404666                           | 13852661                           | 13712274                           | 4654296   |          |    | 0.999 |
| 0.500000(1)                        | 1.2500 (2)                         | 2.5000 (3)                         | 30.0 (4)                           | 60.0 (5)                           | 100.0 (6)                          |           | 13803812 |    |       |
| 11312645                           | 20026015                           | 43450426                           | 432139991                          | 831159688                          | 1371227441                         |           |          |    |       |
| 121594(1)                          | 121594(2)                          | 121594(3)                          | 121594(4)                          | 121594(5)                          | 121594(6)                          |           |          |    |       |

## Preliminary Report

Eurofins Denver  
ICV, ICal Verification Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13  
Lims ID: ICV  
Client ID:  
Sample Type: ICV  
Inject. Date: 18-May-2023 13:25:00 ALS Bottle#: 0 Worklist Smp#: 7  
Injection Vol: 5.0 ul Dil. Factor: 1.0000  
Sample Info: 280-0121594-007  
Misc. Info.: 280-0121594-007  
Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
Sublist: chrom-Anions\_IC7\*sub4  
Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
Limit Group: Wet - Anions 28D  
Last Update: 19-May-2023 12:32:45 Calib Date: 18-May-2023 13:10:00  
Integrator: Falcon  
Quant Method: External Standard Quant By: Initial Calibration  
Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13  
Column 1 : Det: Info 2\_091554\_1  
Process Host: CTX1603  
First Level Reviewer: LVW8 Date: 19-May-2023 12:38:30  
Start Cal Date: 18-May-2023 11:56:00  
End Cal Date: 18-May-2023 13:10:00

| Compound   | Amount Added | Amount Detected | %Drift | Max. %Drift | %Rec  | %Rec Limits |
|------------|--------------|-----------------|--------|-------------|-------|-------------|
| 1 Fluoride | 2.00         | 2.34            | * 16.8 | 10.0        | 116.8 |             |
| 2 Chloride | 80.0         | 83.0            | 3.8    | 10.0        | 103.8 |             |
| 4 Bromide  | 2.00         | 1.91            | -4.7   | 10.0        | 95.3  |             |
| 7 Sulfate  | 80.0         | 81.8            | 2.3    | 10.0        | 102.3 |             |

## Preliminary Report

Eurofins Denver  
ICV, ICal Verification Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13-  
Lims ID: ICV  
Client ID:  
Sample Type: ICV  
Inject. Date: 18-May-2023 13:25:00 ALS Bottle#: 0 Worklist Smp#: 7  
Injection Vol: 5.0 ul Dil. Factor: 1.0000  
Sample Info: 280-0121594-007  
Misc. Info.: 280-0121594-007  
Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
Sublist: chrom-Anions\_IC7\*sub4  
Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
Limit Group: Wet - Anions  
Last Update: 19-May-2023 12:38:46 Calib Date: 18-May-2023 13:10:00  
Integrator: Falcon  
Quant Method: External Standard Quant By: Initial Calibration  
Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13-  
Column 1 : Det: Info 2\_091554\_1  
Process Host: CTX1603  
First Level Reviewer: LVW8 Date: 19-May-2023 12:38:46  
Start Cal Date: 18-May-2023 11:56:00  
End Cal Date: 18-May-2023 13:10:00

| Compound            | Amount Added | Amount Detected | %Drift | Max. %Drift | %Rec  | %Rec Limits |
|---------------------|--------------|-----------------|--------|-------------|-------|-------------|
| 3 Nitrite as N      | 2.00         | 1.99            | -0.3   | 10.0        | 99.7  |             |
| 5 Nitrate as N      | 2.00         | 1.95            | -2.3   | 10.0        | 97.7  |             |
| 6 Orthophosphate as | 2.00         | 2.30            | * 14.8 | 10.0        | 114.8 |             |

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13  
 Lims ID: ICV  
 Client ID:  
 Sample Type: ICV  
 Inject. Date: 18-May-2023 13:25:00 ALS Bottle#: 0 Worklist Smp#: 7  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121594-007  
 Misc. Info.: 280-0121594-007  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Sublist: chrom-Anions\_IC7\*sub4  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
 Limit Group: Wet - Anions  
 Last Update: 19-May-2023 12:38:46 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1603

First Level Reviewer: LVW8

Date: 19-May-2023 12:38:46

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response   | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|------------|------------------|--------------------|-------|
| 1 Fluoride            | 3.145        | 3.145            | 0.000            | 93724806   | NC               | NC                 | Ma    |
| 2 Chloride            | 4.328        | 4.328            | 0.000            | 1555403505 | NC               | NC                 |       |
| 3 Nitrite as N        | 5.000        | 5.000            | 0.000            | 84485309   | 2.00             | 1.99               |       |
| 4 Bromide             | 6.162        | 6.162            | 0.000            | 14920856   | NC               | NC                 |       |
| 5 Nitrate as N        | 6.917        | 6.917            | 0.000            | 86743393   | 2.00             | 1.95               |       |
| 6 Orthophosphate as P | 9.058        | 9.058            | 0.000            | 59448017   | 2.00             | 2.30               |       |
| 7 Sulfate             | 10.230       | 10.230           | 0.000            | 1134012406 | NC               | NC                 |       |

**QC Flag Legend**

Processing Flags

NC - Not Calibrated

Review Flags

M - Manually Integrated

a - User Assigned ID

**Reagents:**

|                  |                    |           |
|------------------|--------------------|-----------|
| IC SO4 ICV_00024 | Amount Added: 0.80 | Units: mL |
| CI ICV Std_00006 | Amount Added: 0.80 | Units: mL |
| IC ICV 5_00405   | Amount Added: 0.80 | Units: mL |

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-134050.d

Injection Date: 18-May-2023 13:25:00

Instrument ID: WC\_IonChrom10

Operator ID: wetchemd

Lims ID: ICV

Worklist Smp#: 7

Client ID:


Injection Vol: 5.0 ul

Dil. Factor: 1.0000

ALS Bottle#: 0

Method: Anions\_IC10

Limit Group: Wet - Anions





Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13  
 Lims ID: ICB  
 Client ID:  
 Sample Type: ICB  
 Inject. Date: 18-May-2023 13:40:00 ALS Bottle#: 0 Worklist Smp#: 8  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121594-008  
 Misc. Info.: 280-0121594-008  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
 Limit Group: Wet - Anions  
 Last Update: 19-May-2023 12:31:30 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1603

First Level Reviewer: LVW8

Date: 19-May-2023 11:57:30

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 1 Fluoride            |              | 3.145            |                  |          |                  | ND                 |       |
| 2 Chloride            |              | 4.328            |                  |          |                  | ND                 |       |
| 3 Nitrite as N        |              | 5.000            |                  |          |                  | ND                 |       |
| 4 Bromide             |              | 6.162            |                  |          |                  | ND                 | U     |
| 5 Nitrate as N        |              | 6.917            |                  |          |                  | ND                 |       |
| 6 Orthophosphate as P | 9.123        | 9.058            | 0.065            | 12987322 |                  | 0.3014             |       |
| 7 Sulfate             | 10.285       | 10.230           | 0.055            | 1071622  |                  | NC                 | M     |

**QC Flag Legend**

Processing Flags

NC - Not Calibrated

Review Flags

M - Manually Integrated

U - Marked Undetected

Report Date: 19-May-2023 12:31:31

Chrom Revision: 2.3 16-May-2023 16:37:30

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-135549.d

Injection Date: 18-May-2023 13:40:00

Instrument ID: WC\_IonChrom10

Operator ID: wetchemd

Lims ID: ICB

Worklist Smp#: 8

Client ID:


Injection Vol: 5.0 ul

Dil. Factor: 1.0000

ALS Bottle#: 0

Method: Anions\_IC10

Limit Group: Wet - Anions



Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-142  
 Lims ID: MRL  
 Client ID:  
 Sample Type: MRL  
 Inject. Date: 18-May-2023 14:10:00 ALS Bottle#: 0 Worklist Smp#: 10  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121594-010  
 Misc. Info.: 280-0121594-010  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
 Limit Group: Wet - Anions  
 Last Update: 19-May-2023 12:31:32 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-132  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1603

First Level Reviewer: LVW8

Date: 18-May-2023 14:58:41

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 1 Fluoride            |              | 3.145            |                  |          | ND               | ND                 |       |
| 2 Chloride            | 4.302        | 4.313            | -0.011           | 43926233 | NC               | NC                 |       |
| 3 Nitrite as N        | 5.003        | 4.998            | 0.005            | 8516946  | 0.2500           | 0.2213             |       |
| 4 Bromide             | 6.178        | 6.163            | 0.015            | 2736457  | NC               | NC                 |       |
| 5 Nitrate as N        | 6.958        | 6.913            | 0.045            | 9978227  | 0.2500           | 0.2475             |       |
| 6 Orthophosphate as P | 9.105        | 9.060            | 0.045            | 20685657 | 0.2500           | 0.6318             |       |
| 7 Sulfate             | 10.265       | 10.247           | 0.018            | 45774744 | NC               | NC                 |       |

**QC Flag Legend**

Processing Flags

NC - Not Calibrated

ND - Not Detected or Marked ND

**Reagents:**

IC CAL cl/so4\_00480

Amount Added: 0.10

Units: mL

IC Cal low\_00709

Amount Added: 0.05

Units: mL

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-142550.d

Injection Date: 18-May-2023 14:10:00

Instrument ID: WC\_IonChrom10

Operator ID: wetchemd

Lims ID: MRL

Worklist Smp#: 10

Client ID:


Injection Vol: 5.0 ul

Dil. Factor: 1.0000

ALS Bottle#: 0

Method: Anions\_IC10

Limit Group: Wet - Anions



Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-14  
 Lims ID: LCS  
 Client ID:  
 Sample Type: LCS  
 Inject. Date: 18-May-2023 14:25:00 ALS Bottle#: 0 Worklist Smp#: 11  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121594-011  
 Misc. Info.: 280-0121594-011  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
 Limit Group: Wet - Anions  
 Last Update: 19-May-2023 12:31:32 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13

Column 1 : Det: Info 2\_091554\_1

Process Host: CTX1603

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response  | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|-----------|------------------|--------------------|-------|
| 1 Fluoride            |              | 3.145            |                  |           | ND               | ND                 |       |
| 2 Chloride            | 4.313        | 4.313            | 0.000            | 924355606 | NC               | NC                 |       |
| 3 Nitrite as N        | 4.998        | 4.998            | 0.000            | 102973098 | 2.50             | 2.42               |       |
| 4 Bromide             | 6.163        | 6.163            | 0.000            | 18068153  | NC               | NC                 |       |
| 5 Nitrate as N        | 6.912        | 6.913            | -0.001           | 109333247 | 2.50             | 2.46               |       |
| 6 Orthophosphate as P | 9.062        | 9.060            | 0.002            | 71396675  | 2.50             | 2.81               |       |
| 7 Sulfate             | 10.248       | 10.247           | 0.001            | 701445617 | NC               | NC                 |       |

### QC Flag Legend

Processing Flags

NC - Not Calibrated

ND - Not Detected or Marked ND

### Reagents:


IC LCS\_01954

Amount Added: 5.00

Units: mL

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-144050.d  
Injection Date: 18-May-2023 14:25:00 Instrument ID: WC\_IonChrom10 Operator ID: wetchemd  
Lims ID: LCS Worklist Smp#: 11  
Client ID:  
Injection Vol: 5.0 ul Dil. Factor: 1.0000 ALS Bottle#: 0  
Method: Anions\_IC10 Limit Group: Wet - Anions



Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-14  
 Lims ID: LCSD  
 Client ID:  
 Sample Type: LCSD  
 Inject. Date: 18-May-2023 14:40:00 ALS Bottle#: 0 Worklist Smp#: 12  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121594-012  
 Misc. Info.: 280-0121594-012  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
 Limit Group: Wet - Anions  
 Last Update: 19-May-2023 12:31:32 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1603

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response  | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|-----------|------------------|--------------------|-------|
| 1 Fluoride            |              | 3.145            |                  |           | ND               | ND                 |       |
| 2 Chloride            | 4.313        | 4.313            | 0.000            | 921265063 | NC               | NC                 |       |
| 3 Nitrite as N        | 4.998        | 4.998            | 0.000            | 102379982 | 2.50             | 2.41               |       |
| 4 Bromide             | 6.163        | 6.163            | 0.000            | 17988727  | NC               | NC                 |       |
| 5 Nitrate as N        | 6.912        | 6.913            | -0.001           | 108961704 | 2.50             | 2.45               |       |
| 6 Orthophosphate as P | 9.062        | 9.060            | 0.002            | 71915979  | 2.50             | 2.83               |       |
| 7 Sulfate             | 10.250       | 10.247           | 0.003            | 695470892 | NC               | NC                 |       |

### QC Flag Legend

Processing Flags

NC - Not Calibrated

ND - Not Detected or Marked ND

### Reagents:


IC LCS\_01954

Amount Added: 5.00

Units: mL

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-145552.d  
Injection Date: 18-May-2023 14:40:00 Instrument ID: WC\_IonChrom10 Operator ID: wetchemd  
Lims ID: LCSD Worklist Smp#: 12  
Client ID:  
Injection Vol: 5.0 ul Dil. Factor: 1.0000 ALS Bottle#: 0  
Method: Anions\_IC10 Limit Group: Wet - Anions





Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-15  
 Lims ID: MB  
 Client ID:  
 Sample Type: MB  
 Inject. Date: 18-May-2023 14:55:00 ALS Bottle#: 0 Worklist Smp#: 13  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121594-013  
 Misc. Info.: 280-0121594-013  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
 Limit Group: Wet - Anions  
 Last Update: 19-May-2023 12:31:35 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1603

First Level Reviewer: LVW8

Date: 19-May-2023 11:57:53

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 1 Fluoride            |              | 3.145            |                  |          |                  | ND                 |       |
| 2 Chloride            |              | 4.313            |                  |          |                  | ND                 |       |
| 3 Nitrite as N        |              | 4.998            |                  |          |                  | ND                 |       |
| 4 Bromide             |              | 6.163            |                  |          |                  | ND                 | U     |
| 5 Nitrate as N        |              | 6.913            |                  |          |                  | ND                 |       |
| 6 Orthophosphate as P | 9.127        | 9.127            | 0.000            | 13698318 |                  | 0.3319             |       |
| 7 Sulfate             | 10.300       | 10.300           | 0.000            | 1160087  |                  | NC                 | M     |

**QC Flag Legend**

Processing Flags

NC - Not Calibrated

Review Flags

M - Manually Integrated

U - Marked Undetected

Report Date: 19-May-2023 12:31:35

Chrom Revision: 2.3 16-May-2023 16:37:30

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-151052.d

Injection Date: 18-May-2023 14:55:00

Instrument ID: WC\_IonChrom10

Operator ID: wetchemd

Lims ID: MB

Worklist Smp#: 13

Client ID:


Injection Vol: 5.0 ul

Dil. Factor: 1.0000

ALS Bottle#: 0

Method: Anions\_IC10

Limit Group: Wet - Anions



Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-16  
 Lims ID: 280-176674-A-5  
 Client ID: LL12mw-185-230401-GW  
 Sample Type: Client  
 Inject. Date: 18-May-2023 16:26:00 ALS Bottle#: 0 Worklist Smp#: 19  
 Injection Vol: 5.0 ul Dil. Factor: 20.0000  
 Sample Info: 280-0121594-019  
 Misc. Info.: 280-0121594-019  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
 Limit Group: Wet - Anions  
 Last Update: 19-May-2023 12:31:35 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1603

First Level Reviewer: LVW8

Date: 19-May-2023 11:59:05

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response  | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|-----------|--------------------|-------|
| 1 Fluoride            |              | 3.145            |                  |           | ND                 |       |
| 2 Chloride            | 4.312        | 4.313            | -0.001           | 766614887 | NC                 | M     |
| 3 Nitrite as N        |              | 4.998            |                  |           | ND                 |       |
| 4 Bromide             | 6.152        | 6.163            | -0.011           | 1354596   | NC                 |       |
| 5 Nitrate as N        | 6.905        | 6.913            | -0.008           | 149556743 | 3.35               | M     |
| 6 Orthophosphate as P | 9.153        | 9.127            | 0.026            | 1762122   | -0.1804            |       |
| 7 Sulfate             | 10.247       | 10.300           | -0.053           | 224738382 | NC                 | M     |

**QC Flag Legend**

Processing Flags


NC - Not Calibrated

Review Flags

M - Manually Integrated

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-164102.d  
Injection Date: 18-May-2023 16:26:00 Instrument ID: WC\_IonChrom10 Operator ID: wetchemd  
Lims ID: 280-176674-A-5 Lab Sample ID: 280-176674-5 Worklist Smp#: 19  
Client ID: LL12mw-185-230401-GW  
Injection Vol: 5.0 ul Dil. Factor: 20.0000 ALS Bottle#: 0  
Method: Anions\_IC10 Limit Group: Wet - Anions



## Eurofins Denver


|                 |                                                                                                     |                |                 |                |    |
|-----------------|-----------------------------------------------------------------------------------------------------|----------------|-----------------|----------------|----|
| Data File:      | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230518-164102.d |                |                 |                |    |
| Injection Date: | 18-May-2023 16:26:00                                                                                | Instrument ID: | WC_IonChrom10   |                |    |
| Lims ID:        | 280-176674-A-5                                                                                      | Lab Sample ID: | 280-176674-5    |                |    |
| Client ID:      | LL12mw-185-230401-GW                                                                                |                |                 |                |    |
| Operator ID:    | wetchemd                                                                                            | ALS Bottle#:   | 0               | Worklist Smp#: | 19 |
| Injection Vol:  | 5.0 ul                                                                                              | Dil. Factor:   | 20.0000         |                |    |
| Method:         | Anions_IC10                                                                                         | Limit Group:   | Wet - Anions    |                |    |
| Column:         |                                                                                                     | Detector:      | Info 2_091554_1 |                |    |

**5 Nitrate as N, CAS: 14797-55-8**

Signal: 1


RT: 6.91  
Area: 150795607  
Amount: 3.377460  
Amount Units: ug/ml

## Processing Integration Results



RT: 6.91  
Area: 149556743  
Amount: 3.349924  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LVW8, 19-May-2023 11:58:59 -06:00:00 (UTC)

Audit Action: Manually Integrated

Audit Reason: Baseline Smoothing

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-16  
 Lims ID: 280-176674-A-6  
 Client ID: LL12mw-185-230402-GW  
 Sample Type: Client  
 Inject. Date: 18-May-2023 16:41:00 ALS Bottle#: 0 Worklist Smp#: 20  
 Injection Vol: 5.0 ul Dil. Factor: 20.0000  
 Sample Info: 280-0121594-020  
 Misc. Info.: 280-0121594-020  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
 Limit Group: Wet - Anions  
 Last Update: 19-May-2023 12:31:35 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1603

First Level Reviewer: LVW8

Date: 19-May-2023 11:59:17

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response  | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|-----------|--------------------|-------|
| 1 Fluoride            |              | 3.145            |                  |           | ND                 |       |
| 2 Chloride            | 4.310        | 4.313            | -0.003           | 841073588 | NC                 |       |
| 3 Nitrite as N        |              | 4.998            |                  |           | ND                 |       |
| 4 Bromide             | 6.182        | 6.163            | 0.019            | 1031418   | NC                 |       |
| 5 Nitrate as N        | 6.900        | 6.913            | -0.013           | 165199089 | 3.70               | M     |
| 6 Orthophosphate as P | 9.192        | 9.127            | 0.065            | 2106909   | -0.1656            |       |
| 7 Sulfate             | 10.243       | 10.300           | -0.057           | 248437087 | NC                 | M     |

**QC Flag Legend**

Processing Flags

NC - Not Calibrated

Review Flags


M - Manually Integrated

Report Date: 19-May-2023 12:31:39

Chrom Revision: 2.3 16-May-2023 16:37:30

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-165605.d  
Injection Date: 18-May-2023 16:41:00 Instrument ID: WC\_IonChrom10 Operator ID: wetchemd  
Lims ID: 280-176674-A-6 Lab Sample ID: 280-176674-6 Worklist Smp#: 20  
Client ID: LL12mw-185-230402-GW  
Injection Vol: 5.0 ul Dil. Factor: 20.0000 ALS Bottle#: 0  
Method: Anions\_IC10 Limit Group: Wet - Anions



## Eurofins Denver


|                 |                                                                                                     |                |                 |                |    |
|-----------------|-----------------------------------------------------------------------------------------------------|----------------|-----------------|----------------|----|
| Data File:      | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230518-165605.d |                |                 |                |    |
| Injection Date: | 18-May-2023 16:41:00                                                                                | Instrument ID: | WC_IonChrom10   |                |    |
| Lims ID:        | 280-176674-A-6                                                                                      | Lab Sample ID: | 280-176674-6    |                |    |
| Client ID:      | LL12mw-185-230402-GW                                                                                |                |                 |                |    |
| Operator ID:    | wetchemd                                                                                            | ALS Bottle#:   | 0               | Worklist Smp#: | 20 |
| Injection Vol:  | 5.0 ul                                                                                              | Dil. Factor:   | 20.0000         |                |    |
| Method:         | Anions_IC10                                                                                         | Limit Group:   | Wet - Anions    |                |    |
| Column:         |                                                                                                     | Detector       | Info 2_091554_1 |                |    |

**5 Nitrate as N, CAS: 14797-55-8**

Signal: 1


RT: 6.90  
Area: 167240152  
Amount: 3.742969  
Amount Units: ug/ml

## Processing Integration Results



RT: 6.90  
Area: 165199089  
Amount: 3.697603  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LVW8, 19-May-2023 11:59:11 -06:00:00 (UTC)

Audit Action: Manually Integrated

Audit Reason: Baseline Smoothing



Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-179  
 Lims ID: ccv  
 Client ID:  
 Sample Type: CCV  
 Inject. Date: 18-May-2023 17:41:00 ALS Bottle#: 0 Worklist Smp#: 24  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121594-024  
 Misc. Info.: 280-0121594-024  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Sublist: chrom-Anions\_IC10\*sub5  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
 Limit Group: Wet - Anions  
 Last Update: 19-May-2023 12:31:42 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-132  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1603

First Level Reviewer: XAY4

Date: 19-May-2023 09:47:24

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response  | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|-----------|------------------|--------------------|-------|
| 1 Fluoride            |              | 3.145            |                  |           | ND               | ND                 |       |
| 2 Chloride            | 4.312        | 4.313            | -0.001           | 924998506 | NC               | NC                 |       |
| 3 Nitrite as N        | 5.000        | 4.998            | 0.002            | 107797836 | 2.50             | 2.54               |       |
| 4 Bromide             | 6.140        | 6.163            | -0.023           | 17602305  | NC               | NC                 |       |
| 5 Nitrate as N        | 6.903        | 6.913            | -0.010           | 108124722 | 2.50             | 2.43               |       |
| 6 Orthophosphate as P | 9.045        | 9.127            | -0.082           | 53503464  | 2.50             | 2.04               |       |
| 7 Sulfate             | 10.223       | 10.300           | -0.077           | 680174954 | NC               | NC                 |       |

**QC Flag Legend**

Processing Flags

NC - Not Calibrated

ND - Not Detected or Marked ND

**Reagents:**

IC LCS\_01954

Amount Added: 10.00

Units: mL

Report Date: 19-May-2023 12:31:42

Chrom Revision: 2.3 16-May-2023 16:37:30

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-175619.d

Injection Date: 18-May-2023 17:41:00

Instrument ID: WC\_IonChrom10

Operator ID: wetchemd

Lims ID: ccv

Worklist Smp#: 24

Client ID:


Injection Vol: 5.0 ul

Dil. Factor: 1.0000

ALS Bottle#: 0

Method: Anions\_IC10

Limit Group: Wet - Anions



Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-18  
 Lims ID: ccb  
 Client ID:  
 Sample Type: CCB  
 Inject. Date: 18-May-2023 17:56:00 ALS Bottle#: 0 Worklist Smp#: 25  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121594-025  
 Misc. Info.: 280-0121594-025  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
 Limit Group: Wet - Anions  
 Last Update: 19-May-2023 12:31:42 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1603

First Level Reviewer: LVW8

Date: 19-May-2023 12:00:26

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 1 Fluoride            |              | 3.145            |                  |          |                  | ND                 |       |
| 2 Chloride            |              | 4.313            |                  |          |                  | ND                 |       |
| 3 Nitrite as N        |              | 4.998            |                  |          |                  | ND                 | U     |
| 4 Bromide             |              | 6.163            |                  |          |                  | ND                 | U     |
| 5 Nitrate as N        |              | 6.913            |                  |          |                  | ND                 |       |
| 6 Orthophosphate as P | 9.133        | 9.127            | 0.006            | 5344408  |                  | -0.0266            |       |
| 7 Sulfate             | 10.305       | 10.300           | 0.005            | 810551   |                  | NC                 | M     |

**QC Flag Legend**

Processing Flags

NC - Not Calibrated

Review Flags

M - Manually Integrated

U - Marked Undetected

Report Date: 19-May-2023 12:31:43

Chrom Revision: 2.3 16-May-2023 16:37:30

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-181124.d

Injection Date: 18-May-2023 17:56:00

Instrument ID: WC\_IonChrom10

Operator ID: wetchemd

Lims ID: ccb

Worklist Smp#: 25

Client ID:


Injection Vol: 5.0 ul

Dil. Factor: 1.0000

ALS Bottle#: 0

Method: Anions\_IC10

Limit Group: Wet - Anions




## Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-181124.d  
Injection Date: 18-May-2023 17:56:00 Instrument ID: WC\_IonChrom10  
Lims ID: ccb  
Client ID:  
Operator ID: wetchemd ALS Bottle#: 0 Worklist Smp#: 25  
Injection Vol: 5.0 ul Dil. Factor: 1.0000  
Method: Anions\_IC10 Limit Group: Wet - Anions  
Column: Detector Info 2\_091554\_1

**3 Nitrite as N, CAS: 14797-65-0, Signal: 1**

RT: 4.78  
Response: 44417  
Amount: 0



Reviewer: LVW8, 19-May-2023 12:00:26

Audit Action: Manually Integrated

Audit Reason:

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-19  
 Lims ID: 280-176674-A-7  
 Client ID: FBQmw-173-230401-GW-R  
 Sample Type: Client  
 Inject. Date: 18-May-2023 19:26:00 ALS Bottle#: 0 Worklist Smp#: 31  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121594-031  
 Misc. Info.: 280-0121594-031  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
 Limit Group: Wet - Anions  
 Last Update: 19-May-2023 12:31:42 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1603

First Level Reviewer: LVW8

Date: 19-May-2023 12:02:08

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response  | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|-----------|--------------------|-------|
| 1 Fluoride            |              | 3.145            |                  |           | ND                 |       |
| 2 Chloride            | 4.303        | 4.313            | -0.010           | 19574865  | NC                 | M     |
| 3 Nitrite as N        |              | 4.998            |                  |           | ND                 |       |
| 4 Bromide             | 6.212        | 6.163            | 0.049            | 1197331   | NC                 |       |
| 5 Nitrate as N        | 6.975        | 6.913            | 0.062            | 7707821   | 0.1971             | M     |
| 6 Orthophosphate as P | 9.138        | 9.127            | 0.011            | 8994927   | 0.1300             |       |
| 7 Sulfate             | 10.243       | 10.300           | -0.057           | 580169865 | NC                 | M     |

**QC Flag Legend**

Processing Flags


NC - Not Calibrated

Review Flags

M - Manually Integrated

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-194147.d  
Injection Date: 18-May-2023 19:26:00 Instrument ID: WC\_IonChrom10 Operator ID: wetchemd  
Lims ID: 280-176674-A-7 Lab Sample ID: 280-176674-7 Worklist Smp#: 31  
Client ID: FBQmw-173-230401-GW-R  
Injection Vol: 5.0 ul Dil. Factor: 1.0000 ALS Bottle#: 0  
Method: Anions\_IC10 Limit Group: Wet - Anions



## Eurofins Denver


|                 |                                                                                                     |                |                 |                |    |
|-----------------|-----------------------------------------------------------------------------------------------------|----------------|-----------------|----------------|----|
| Data File:      | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230518-194147.d |                |                 |                |    |
| Injection Date: | 18-May-2023 19:26:00                                                                                | Instrument ID: | WC_IonChrom10   |                |    |
| Lims ID:        | 280-176674-A-7                                                                                      | Lab Sample ID: | 280-176674-7    |                |    |
| Client ID:      | FBQmw-173-230401-GW-R                                                                               |                |                 |                |    |
| Operator ID:    | wetchemd                                                                                            | ALS Bottle#:   | 0               | Worklist Smp#: | 31 |
| Injection Vol:  | 5.0 ul                                                                                              | Dil. Factor:   | 1.0000          |                |    |
| Method:         | Anions_IC10                                                                                         | Limit Group:   | Wet - Anions    |                |    |
| Column:         |                                                                                                     | Detector:      | Info 2_091554_1 |                |    |

**5 Nitrate as N, CAS: 14797-55-8**

Signal: 1


RT: 6.98  
Area: 9538418  
Amount: 0.237770  
Amount Units: ug/ml

## Processing Integration Results



RT: 6.98  
Area: 7707821  
Amount: 0.197082  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LVW8, 19-May-2023 12:01:59 -06:00:00 (UTC)

Audit Action: Manually Integrated

Audit Reason: Baseline Smoothing



Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-19  
 Lims ID: 280-176674-A-8  
 Client ID: FBQmw-173-230402-GW-R  
 Sample Type: Client  
 Inject. Date: 18-May-2023 19:41:00 ALS Bottle#: 0 Worklist Smp#: 32  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121594-032  
 Misc. Info.: 280-0121594-032  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
 Limit Group: Wet - Anions  
 Last Update: 19-May-2023 12:31:42 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1603

First Level Reviewer: LVW8

Date: 19-May-2023 12:02:17

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response  | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|-----------|--------------------|-------|
| 1 Fluoride            |              | 3.145            |                  |           | ND                 |       |
| 2 Chloride            | 4.302        | 4.313            | -0.011           | 20099453  | NC                 | M     |
| 3 Nitrite as N        |              | 4.998            |                  |           | ND                 |       |
| 4 Bromide             | 6.210        | 6.163            | 0.047            | 1338653   | NC                 |       |
| 5 Nitrate as N        | 6.977        | 6.913            | 0.064            | 8138449   | 0.2067             | M     |
| 6 Orthophosphate as P | 9.148        | 9.127            | 0.021            | 5690064   | -0.0118            |       |
| 7 Sulfate             | 10.245       | 10.300           | -0.055           | 627177033 | NC                 |       |

**QC Flag Legend**

Processing Flags


NC - Not Calibrated

Review Flags

M - Manually Integrated

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-195651.d  
Injection Date: 18-May-2023 19:41:00 Instrument ID: WC\_IonChrom10 Operator ID: wetchemd  
Lims ID: 280-176674-A-8 Lab Sample ID: 280-176674-8 Worklist Smp#: 32  
Client ID: FBQmw-173-230402-GW-R  
Injection Vol: 5.0 ul Dil. Factor: 1.0000 ALS Bottle#: 0  
Method: Anions\_IC10 Limit Group: Wet - Anions



## Eurofins Denver


|                 |                                                                                                     |                |                 |                |    |
|-----------------|-----------------------------------------------------------------------------------------------------|----------------|-----------------|----------------|----|
| Data File:      | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230518-195651.d |                |                 |                |    |
| Injection Date: | 18-May-2023 19:41:00                                                                                | Instrument ID: | WC_IonChrom10   |                |    |
| Lims ID:        | 280-176674-A-8                                                                                      | Lab Sample ID: | 280-176674-8    |                |    |
| Client ID:      | FBQmw-173-230402-GW-R                                                                               |                |                 |                |    |
| Operator ID:    | wetchemd                                                                                            | ALS Bottle#:   | 0               | Worklist Smp#: | 32 |
| Injection Vol:  | 5.0 ul                                                                                              | Dil. Factor:   | 1.0000          |                |    |
| Method:         | Anions_IC10                                                                                         | Limit Group:   | Wet - Anions    |                |    |
| Column:         |                                                                                                     | Detector       | Info 2_091554_1 |                |    |

**5 Nitrate as N, CAS: 14797-55-8**

Signal: 1


RT: 6.98  
Area: 10229136  
Amount: 0.253122  
Amount Units: ug/ml

## Processing Integration Results



RT: 6.98  
Area: 8138449  
Amount: 0.206653  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LVW8, 19-May-2023 12:02:12 -06:00:00 (UTC)

Audit Action: Manually Integrated

Audit Reason: Baseline Smoothing

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-20230518-131594.b  
 Lims ID: ccv  
 Client ID:  
 Sample Type: CCV  
 Inject. Date: 18-May-2023 20:42:00 ALS Bottle#: 0 Worklist Smp#: 36  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121594-036  
 Misc. Info.: 280-0121594-036  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Sublist: chrom-Anions\_IC10\*sub5  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
 Limit Group: Wet - Anions  
 Last Update: 19-May-2023 12:31:51 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-131594.b

Column 1 : Det: Info 2\_091554\_1

Process Host: CTX1603

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response  | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|-----------|------------------|--------------------|-------|
| 1 Fluoride            |              | 3.145            |                  |           | ND               | ND                 |       |
| 2 Chloride            | 4.315        | 4.313            | 0.002            | 926185158 | NC               | NC                 |       |
| 3 Nitrite as N        | 5.003        | 4.998            | 0.005            | 107572491 | 2.50             | 2.53               |       |
| 4 Bromide             | 6.145        | 6.163            | -0.018           | 17667283  | NC               | NC                 |       |
| 5 Nitrate as N        | 6.910        | 6.913            | -0.003           | 108155910 | 2.50             | 2.43               |       |
| 6 Orthophosphate as P | 9.048        | 9.127            | -0.079           | 53833009  | 2.50             | 2.05               |       |
| 7 Sulfate             | 10.232       | 10.300           | -0.068           | 686565829 | NC               | NC                 |       |

### QC Flag Legend

Processing Flags

NC - Not Calibrated

ND - Not Detected or Marked ND

### Reagents:


IC LCS\_01954

Amount Added: 10.00

Units: mL

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-205701.d  
Injection Date: 18-May-2023 20:42:00 Instrument ID: WC\_IonChrom10 Operator ID: wetchemd  
Lims ID: ccv Worklist Smp#: 36  
Client ID:  
Injection Vol: 5.0 ul Dil. Factor: 1.0000 ALS Bottle#: 0  
Method: Anions\_IC10 Limit Group: Wet - Anions



Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-21.  
 Lims ID: ccb  
 Client ID:  
 Sample Type: CCB  
 Inject. Date: 18-May-2023 20:57:00 ALS Bottle#: 0 Worklist Smp#: 37  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121594-037  
 Misc. Info.: 280-0121594-037  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
 Limit Group: Wet - Anions  
 Last Update: 19-May-2023 12:31:51 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13.  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1603

First Level Reviewer: LVW8

Date: 19-May-2023 12:03:41

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 1 Fluoride            |              | 3.145            |                  |          |                  | ND                 |       |
| 2 Chloride            | 4.305        | 4.313            | -0.008           | 7127336  |                  | NC                 | M     |
| 3 Nitrite as N        | 5.003        | 4.998            | 0.005            | 116037   |                  | 0.0253             | M     |
| 4 Bromide             | 6.167        | 6.163            | 0.004            | 189030   |                  | NC                 | M     |
| 5 Nitrate as N        | 6.973        | 6.913            | 0.060            | 76369    |                  | 0.0275             | M     |
| 6 Orthophosphate as P | 9.052        | 9.127            | -0.075           | 6191477  |                  | 0.009727           |       |
| 7 Sulfate             | 10.250       | 10.300           | -0.050           | 11085723 |                  | NC                 | M     |

**QC Flag Legend**

Processing Flags

NC - Not Calibrated

Review Flags

M - Manually Integrated

Report Date: 19-May-2023 12:31:51

Chrom Revision: 2.3 16-May-2023 16:37:30

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-211202.d

Injection Date: 18-May-2023 20:57:00

Instrument ID: WC\_IonChrom10

Operator ID: wetchemd

Lims ID: ccb

Worklist Smp#: 37

Client ID:


Injection Vol: 5.0 ul

Dil. Factor: 1.0000

ALS Bottle#: 0

Method: Anions\_IC10

Limit Group: Wet - Anions



## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-211202.d  
Injection Date: 18-May-2023 20:57:00 Instrument ID: WC\_IonChrom10  
Lims ID: ccb  
Client ID:  
Operator ID: wetchemd ALS Bottle#: 0 Worklist Smp#: 37  
Injection Vol: 5.0 ul Dil. Factor: 1.0000  
Method: Anions\_IC10 Limit Group: Wet - Anions  
Column: Detector Info 2\_091554\_1

**5 Nitrate as N, CAS: 14797-55-8**

Signal: 1


RT: 6.98  
Area: 265918  
Amount: 0.031672  
Amount Units: ug/ml

## Processing Integration Results



RT: 6.97  
Area: 76369  
Amount: 0.027459  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LVW8, 19-May-2023 12:03:29 -06:00:00 (UTC)

Audit Action: Manually Integrated

Audit Reason: Baseline Smoothing



## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-211202.d  
Injection Date: 18-May-2023 20:57:00 Instrument ID: WC\_IonChrom10  
Lims ID: ccb  
Client ID:  
Operator ID: wetchemd ALS Bottle#: 0 Worklist Smp#: 37  
Injection Vol: 5.0 ul Dil. Factor: 1.0000  
Method: Anions\_IC10 Limit Group: Wet - Anions  
Column: Detector Info 2\_091554\_1

**3 Nitrite as N, CAS: 14797-65-0**

Signal: 1


RT: 5.01  
Area: 191185  
Amount: 0.027032  
Amount Units: ug/ml

## Processing Integration Results



RT: 5.00  
Area: 116037  
Amount: 0.025279  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LVW8, 19-May-2023 12:03:26 -06:00:00 (UTC)

Audit Action: Manually Integrated

Audit Reason: Baseline Smoothing

## Eurofins Denver

|                 |                                                                                                |                |                      |                |    |
|-----------------|------------------------------------------------------------------------------------------------|----------------|----------------------|----------------|----|
| Data File:      | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230518-230 |                |                      |                |    |
| Lims ID:        | ccv                                                                                            |                |                      |                |    |
| Client ID:      |                                                                                                |                |                      |                |    |
| Sample Type:    | CCV                                                                                            |                |                      |                |    |
| Inject. Date:   | 18-May-2023 23:41:00                                                                           | ALS Bottle#:   | 0                    | Worklist Smp#: | 48 |
| Injection Vol:  | 5.0 ul                                                                                         | Dil. Factor:   | 1.0000               |                |    |
| Sample Info:    | 280-0121594-048                                                                                |                |                      |                |    |
| Misc. Info.:    | 280-0121594-048                                                                                |                |                      |                |    |
| Operator ID:    | wetchemd                                                                                       | Instrument ID: | WC_IonChrom10        |                |    |
| Sublist:        | chrom-Anions_IC10*sub5                                                                         |                |                      |                |    |
| Method:         | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Anions_IC10.m                       |                |                      |                |    |
| Limit Group:    | Wet - Anions                                                                                   |                |                      |                |    |
| Last Update:    | 19-May-2023 12:32:00                                                                           | Calib Date:    | 18-May-2023 13:10:00 |                |    |
| Integrator:     | Falcon                                                                                         |                |                      |                |    |
| Quant Method:   | External Standard                                                                              | Quant By:      | Initial Calibration  |                |    |
| Last ICal File: | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230518-130 |                |                      |                |    |
| Column 1 :      | Det: Info 2_091554_1                                                                           |                |                      |                |    |
| Process Host:   | CTX1603                                                                                        |                |                      |                |    |

Date: 19-May-2023 10:54:12

| Compound              | RT (min.) | Exp RT (min.) | Dlt RT (min.) | Response   | Cal Amt ug/ml | OnCol Amt ug/ml | Flags |
|-----------------------|-----------|---------------|---------------|------------|---------------|-----------------|-------|
| 1 Fluoride            |           | 3.145         |               |            | NC            | ND              |       |
| 2 Chloride            | 4.313     | 4.313         | 0.000         | 1006185849 | NC            | NC              |       |
| 3 Nitrite as N        | 5.000     | 4.998         | 0.002         | 113584540  | 2.50          | 2.67            |       |
| 4 Bromide             | 6.127     | 6.163         | -0.036        | 17948809   | NC            | NC              |       |
| 5 Nitrate as N        | 6.897     | 6.913         | -0.016        | 112164302  | 2.50          | 2.52            |       |
| 6 Orthophosphate as P | 9.023     | 9.127         | -0.104        | 59801890   | 2.50          | 2.31            |       |
| 7 Sulfate             | 10.212    | 10.300        | -0.088        | 780948377  | NC            | NC              | M     |

## QC Flag Legend

## Processing Flags

NC - Not Calibrated

## Review Flags

M - Manually Integrated

Reagents:

IC LCS\_01954                      Amount Added: 10.00                      Units: mL

Report Date: 19-May-2023 12:32:00

Chrom Revision: 2.3 16-May-2023 16:37:30

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-235651.d

Injection Date: 18-May-2023 23:41:00

Instrument ID: WC\_IonChrom10

Operator ID: wetchemd

Lims ID: ccv

Worklist Smp#: 48

Client ID:


Injection Vol: 5.0 ul

Dil. Factor: 1.0000

ALS Bottle#: 0

Method: Anions\_IC10

Limit Group: Wet - Anions



## Eurofins Denver

|                 |                                                                                                |                |                      |                |    |
|-----------------|------------------------------------------------------------------------------------------------|----------------|----------------------|----------------|----|
| Data File:      | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230519-003 |                |                      |                |    |
| Lims ID:        | ccb                                                                                            |                |                      |                |    |
| Client ID:      |                                                                                                |                |                      |                |    |
| Sample Type:    | CCB                                                                                            |                |                      |                |    |
| Inject. Date:   | 18-May-2023 23:56:00                                                                           | ALS Bottle#:   | 0                    | Worklist Smp#: | 49 |
| Injection Vol:  | 5.0 ul                                                                                         | Dil. Factor:   | 1.0000               |                |    |
| Sample Info:    | 280-0121594-049                                                                                |                |                      |                |    |
| Misc. Info.:    | 280-0121594-049                                                                                |                |                      |                |    |
| Operator ID:    | wetchemd                                                                                       | Instrument ID: | WC_IonChrom10        |                |    |
| Method:         | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Anions_IC10.m                       |                |                      |                |    |
| Limit Group:    | Wet - Anions                                                                                   |                |                      |                |    |
| Last Update:    | 19-May-2023 12:32:00                                                                           | Calib Date:    | 18-May-2023 13:10:00 |                |    |
| Integrator:     | Falcon                                                                                         |                |                      |                |    |
| Quant Method:   | External Standard                                                                              | Quant By:      | Initial Calibration  |                |    |
| Last ICal File: | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230518-132 |                |                      |                |    |
| Column 1 :      | Det: Info 2_091554_1                                                                           |                |                      |                |    |
| Process Host:   | CTX1603                                                                                        |                |                      |                |    |

First Level Reviewer: XAY4 Date: 19-May-2023 10:59:02

| Compound              | RT (min.) | Exp RT (min.) | Dlt RT (min.) | Response | Cal Amt ug/ml | OnCol Amt ug/ml | Flags |
|-----------------------|-----------|---------------|---------------|----------|---------------|-----------------|-------|
| 1 Fluoride            |           | 3.145         |               |          |               | ND              |       |
| 2 Chloride            | 4.303     | 4.313         | -0.010        | 6683950  |               | NC              | M     |
| 3 Nitrite as N        | 5.003     | 4.998         | 0.005         | 102978   |               | 0.0250          | M     |
| 4 Bromide             | 6.148     | 6.163         | -0.015        | 157855   |               | NC              | M     |
| 5 Nitrate as N        | 6.975     | 6.913         | 0.062         | 54025    |               | 0.0270          | M     |
| 6 Orthophosphate as P | 9.042     | 9.127         | -0.085        | 11087615 |               | 0.2199          |       |
| 7 Sulfate             | 10.235    | 10.300        | -0.065        | 10128974 |               | NC              | M     |

## QC Flag Legend

## Processing Flags

NC - Not Calibrated

## Review Flags

M - Manually Integrated

Report Date: 19-May-2023 12:32:00

Chrom Revision: 2.3 16-May-2023 16:37:30

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230519-001149.d

Injection Date: 18-May-2023 23:56:00

Instrument ID: WC\_IonChrom10

Operator ID: wetchemd

Lims ID: ccb

Worklist Smp#: 49

Client ID:


Injection Vol: 5.0 ul

Dil. Factor: 1.0000

ALS Bottle#: 0

Method: Anions\_IC10

Limit Group: Wet - Anions



## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230519-001149.d  
Injection Date: 18-May-2023 23:56:00 Instrument ID: WC\_IonChrom10  
Lims ID: ccb  
Client ID:  
Operator ID: wetchemd ALS Bottle#: 0 Worklist Smp#: 49  
Injection Vol: 5.0 ul Dil. Factor: 1.0000  
Method: Anions\_IC10 Limit Group: Wet - Anions  
Column: Detector Info 2\_091554\_1

**5 Nitrate as N, CAS: 14797-55-8**

Signal: 1


RT: 6.98  
Area: 107227  
Amount: 0  
Amount Units: ug/ml

## Processing Integration Results



RT: 6.98  
Area: 54025  
Amount: 0.026963  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LVW8, 19-May-2023 12:07:00 -06:00:00 (UTC)

Audit Action: Manually Integrated

Audit Reason: Baseline Smoothing

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230519-01.  
 Lims ID: lcs  
 Client ID:  
 Sample Type: LCS  
 Inject. Date: 19-May-2023 00:56:00 ALS Bottle#: 0 Worklist Smp#: 53  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121594-053  
 Misc. Info.: 280-0121594-053  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
 Limit Group: Wet - Anions  
 Last Update: 19-May-2023 12:32:00 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13.  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1603

First Level Reviewer: XAY4

Date: 19-May-2023 10:54:45

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response  | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|-----------|------------------|--------------------|-------|
| 1 Fluoride            |              | 3.145            |                  |           | NC               | ND                 |       |
| 2 Chloride            | 4.315        | 4.313            | 0.002            | 985100117 | NC               | NC                 |       |
| 3 Nitrite as N        | 5.003        | 4.998            | 0.005            | 112986038 | 2.50             | 2.66               |       |
| 4 Bromide             | 6.132        | 6.163            | -0.031           | 18390855  | NC               | NC                 |       |
| 5 Nitrate as N        | 6.905        | 6.913            | -0.008           | 111290713 | 2.50             | 2.50               |       |
| 6 Orthophosphate as P | 9.028        | 9.127            | -0.099           | 58043149  | 2.50             | 2.24               |       |
| 7 Sulfate             | 10.212       | 10.300           | -0.088           | 876224058 | NC               | NC                 | M     |

**QC Flag Legend**

Processing Flags

NC - Not Calibrated

Review Flags

M - Manually Integrated

**Reagents:**

IC LCS\_01954

Amount Added: 10.00

Units: mL

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230519-011150.d

Injection Date: 19-May-2023 00:56:00

Instrument ID: WC\_IonChrom10

Operator ID: wetchemd

Lims ID: lcs

Worklist Smp#: 53

Client ID:


Injection Vol: 5.0 ul

Dil. Factor: 1.0000

ALS Bottle#: 0

Method: Anions\_IC10

Limit Group: Wet - Anions





Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230519-012  
 Lims ID: lcsd  
 Client ID:  
 Sample Type: LCSD  
 Inject. Date: 19-May-2023 01:11:00 ALS Bottle#: 0 Worklist Smp#: 54  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121594-054  
 Misc. Info.: 280-0121594-054  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
 Limit Group: Wet - Anions  
 Last Update: 19-May-2023 12:32:00 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-132  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1603

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response  | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|-----------|------------------|--------------------|-------|
| 1 Fluoride            |              | 3.145            |                  |           | ND               | ND                 |       |
| 2 Chloride            | 4.313        | 4.313            | 0.000            | 961779421 | NC               | NC                 |       |
| 3 Nitrite as N        | 5.002        | 4.998            | 0.004            | 110662881 | 2.50             | 2.60               |       |
| 4 Bromide             | 6.130        | 6.163            | -0.033           | 18424761  | NC               | NC                 |       |
| 5 Nitrate as N        | 6.898        | 6.913            | -0.015           | 109150694 | 2.50             | 2.45               |       |
| 6 Orthophosphate as P | 9.020        | 9.127            | -0.107           | 57102319  | 2.50             | 2.19               |       |
| 7 Sulfate             | 10.215       | 10.300           | -0.085           | 717016076 | NC               | NC                 |       |

### QC Flag Legend

Processing Flags

NC - Not Calibrated

ND - Not Detected or Marked ND

### Reagents:

IC LCS\_01954

Amount Added: 10.00

Units: mL

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230519-012652.d

Injection Date: 19-May-2023 01:11:00

Instrument ID: WC\_IonChrom10

Operator ID: wetchemd

Lims ID: lcsd

Worklist Smp#: 54

Client ID:


Injection Vol: 5.0 ul

Dil. Factor: 1.0000

ALS Bottle#: 0

Method: Anions\_IC10

Limit Group: Wet - Anions



Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230519-01  
 Lims ID: mb  
 Client ID:  
 Sample Type: MB  
 Inject. Date: 19-May-2023 01:26:00 ALS Bottle#: 0 Worklist Smp#: 55  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121594-055  
 Misc. Info.: 280-0121594-055  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
 Limit Group: Wet - Anions  
 Last Update: 19-May-2023 12:32:00 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1603

First Level Reviewer: XAY4

Date: 19-May-2023 10:59:34

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 1 Fluoride            |              | 3.145            |                  |          |                  | ND                 |       |
| 2 Chloride            | 4.303        | 4.313            | -0.010           | 6687549  |                  | NC                 | M     |
| 3 Nitrite as N        | 5.003        | 4.998            | 0.005            | 132102   |                  | 0.0257             | M     |
| 4 Bromide             | 6.158        | 6.163            | -0.005           | 137170   |                  | NC                 | M     |
| 5 Nitrate as N        | 6.982        | 6.913            | 0.069            | 60360    |                  | 0.0271             | M     |
| 6 Orthophosphate as P | 9.037        | 9.127            | -0.090           | 10239561 |                  | 0.1835             |       |
| 7 Sulfate             | 10.237       | 10.300           | -0.063           | 9929845  |                  | NC                 | M     |

**QC Flag Legend**

Processing Flags

NC - Not Calibrated

Review Flags

M - Manually Integrated

Report Date: 19-May-2023 12:32:04

Chrom Revision: 2.3 16-May-2023 16:37:30

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230519-014153.d

Injection Date: 19-May-2023 01:26:00

Instrument ID: WC\_IonChrom10

Operator ID: wetchemd

Lims ID: mb

Worklist Smp#: 55

Client ID:


Injection Vol: 5.0 ul

Dil. Factor: 1.0000

ALS Bottle#: 0

Method: Anions\_IC10

Limit Group: Wet - Anions



## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230519-014153.d  
Injection Date: 19-May-2023 01:26:00 Instrument ID: WC\_IonChrom10  
Lims ID: mb  
Client ID:  
Operator ID: wetchemd ALS Bottle#: 0 Worklist Smp#: 55  
Injection Vol: 5.0 ul Dil. Factor: 1.0000  
Method: Anions\_IC10 Limit Group: Wet - Anions  
Column: Detector Info 2\_091554\_1

**5 Nitrate as N, CAS: 14797-55-8**

Signal: 1


RT: 6.98  
Area: 170622  
Amount: 0  
Amount Units: ug/ml

## Processing Integration Results



RT: 6.98  
Area: 60360  
Amount: 0.027103  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LVW8, 19-May-2023 12:07:59 -06:00:00 (UTC)

Audit Action: Manually Integrated

Audit Reason: Baseline Smoothing

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230519-02  
 Lims ID: ccv  
 Client ID:  
 Sample Type: CCV  
 Inject. Date: 19-May-2023 02:42:00 ALS Bottle#: 0 Worklist Smp#: 60  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121594-060  
 Misc. Info.: 280-0121594-060  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Sublist: chrom-Anions\_IC10\*sub5  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
 Limit Group: Wet - Anions  
 Last Update: 19-May-2023 12:32:07 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1603

First Level Reviewer: XAY4

Date: 19-May-2023 09:48:09

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response  | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|-----------|------------------|--------------------|-------|
| 1 Fluoride            |              | 3.145            |                  |           | NC               | ND                 |       |
| 2 Chloride            | 4.315        | 4.313            | 0.002            | 951474195 | NC               | NC                 |       |
| 3 Nitrite as N        | 5.003        | 4.998            | 0.005            | 111742013 | 2.50             | 2.63               |       |
| 4 Bromide             | 6.138        | 6.163            | -0.025           | 17981985  | NC               | NC                 |       |
| 5 Nitrate as N        | 6.903        | 6.913            | -0.010           | 119448879 | 2.50             | 2.68               | M     |
| 6 Orthophosphate as P | 9.030        | 9.127            | -0.097           | 57271017  | 2.50             | 2.20               |       |
| 7 Sulfate             | 10.220       | 10.300           | -0.080           | 712406832 | NC               | NC                 |       |

**QC Flag Legend**

Processing Flags

NC - Not Calibrated

Review Flags

M - Manually Integrated

**Reagents:**

IC LCS\_01954

Amount Added: 10.00

Units: mL

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230519-025700.d

Injection Date: 19-May-2023 02:42:00

Instrument ID: WC\_IonChrom10

Operator ID: wetchemd

Lims ID: ccv

Worklist Smp#: 60

Client ID:


Injection Vol: 5.0 ul

Dil. Factor: 1.0000

ALS Bottle#: 0

Method: Anions\_IC10

Limit Group: Wet - Anions



## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230519-025700.d  
Injection Date: 19-May-2023 02:42:00 Instrument ID: WC\_IonChrom10  
Lims ID: ccv  
Client ID:  
Operator ID: wetchemd ALS Bottle#: 0 Worklist Smp#: 60  
Injection Vol: 5.0 ul Dil. Factor: 1.0000  
Method: Anions\_IC10 Limit Group: Wet - Anions  
Column: Detector Info 2\_091554\_1

**5 Nitrate as N, CAS: 14797-55-8**

Signal: 1


RT: 6.90  
Area: 122651005  
Amount: 2.751897  
Amount Units: ug/ml

## Processing Integration Results



RT: 6.90  
Area: 119448879  
Amount: 2.680724  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: XAY4, 19-May-2023 09:48:05 -06:00:00 (UTC)

Audit Action: Manually Integrated

Audit Reason: Baseline Smoothing



Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230519-03  
 Lims ID: ccb  
 Client ID:  
 Sample Type: CCB  
 Inject. Date: 19-May-2023 02:57:00 ALS Bottle#: 0 Worklist Smp#: 61  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121594-061  
 Misc. Info.: 280-0121594-061  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
 Limit Group: Wet - Anions  
 Last Update: 19-May-2023 12:32:07 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1603

First Level Reviewer: LVW8

Date: 19-May-2023 12:09:33

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 1 Fluoride            |              | 3.145            |                  |          |                  | ND                 |       |
| 2 Chloride            | 4.307        | 4.313            | -0.006           | 7010633  |                  | NC                 | M     |
| 3 Nitrite as N        | 5.005        | 4.998            | 0.007            | 181021   |                  | 0.0268             | M     |
| 4 Bromide             | 6.160        | 6.163            | -0.003           | 189023   |                  | NC                 | M     |
| 5 Nitrate as N        | 7.007        | 6.913            | 0.094            | 431220   |                  | 0.0353             | M     |
| 6 Orthophosphate as P | 9.048        | 9.127            | -0.079           | 9019381  |                  | 0.1311             |       |
| 7 Sulfate             | 10.240       | 10.300           | -0.060           | 10195287 |                  | NC                 | M     |

**QC Flag Legend**

Processing Flags

NC - Not Calibrated

Review Flags

M - Manually Integrated

Report Date: 19-May-2023 12:32:08

Chrom Revision: 2.3 16-May-2023 16:37:30

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230519-031202.d

Injection Date: 19-May-2023 02:57:00

Instrument ID: WC\_IonChrom10

Operator ID: wetchemd

Lims ID: ccb

Worklist Smp#: 61

Client ID:


Injection Vol: 5.0 ul

Dil. Factor: 1.0000

ALS Bottle#: 0

Method: Anions\_IC10

Limit Group: Wet - Anions



## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230519-031202.d  
Injection Date: 19-May-2023 02:57:00 Instrument ID: WC\_IonChrom10  
Lims ID: ccb  
Client ID:  
Operator ID: wetchemd ALS Bottle#: 0 Worklist Smp#: 61  
Injection Vol: 5.0 ul Dil. Factor: 1.0000  
Method: Anions\_IC10 Limit Group: Wet - Anions  
Column: Detector Info 2\_091554\_1

**5 Nitrate as N, CAS: 14797-55-8**

Signal: 1


RT: 7.01  
Area: 633844  
Amount: 0.039850  
Amount Units: ug/ml

## Processing Integration Results



RT: 7.01  
Area: 431220  
Amount: 0.035346  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LVW8, 19-May-2023 12:09:23 -06:00:00 (UTC)

Audit Action: Manually Integrated

Audit Reason: Baseline Smoothing

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230519-05  
 Lims ID: ccv  
 Client ID:  
 Sample Type: CCV  
 Inject. Date: 19-May-2023 05:42:00 ALS Bottle#: 0 Worklist Smp#: 72  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121594-072  
 Misc. Info.: 280-0121594-072  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Sublist: chrom-Anions\_IC10\*sub5  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
 Limit Group: Wet - Anions  
 Last Update: 19-May-2023 12:32:15 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1603

First Level Reviewer: XAY4

Date: 19-May-2023 09:48:31

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response   | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|------------|------------------|--------------------|-------|
| 1 Fluoride            |              | 3.145            |                  |            | NC               | ND                 |       |
| 2 Chloride            | 4.313        | 4.313            | 0.000            | 1048693012 | NC               | NC                 | M     |
| 3 Nitrite as N        | 5.003        | 4.998            | 0.005            | 114016951  | 2.50             | 2.68               | M     |
| 4 Bromide             | 6.137        | 6.163            | -0.026           | 18725437   | NC               | NC                 |       |
| 5 Nitrate as N        | 6.910        | 6.913            | -0.003           | 119161316  | 2.50             | 2.67               | M     |
| 6 Orthophosphate as P | 9.032        | 9.127            | -0.095           | 69679815   | 2.50             | 2.73               |       |
| 7 Sulfate             | 10.187       | 10.300           | -0.113           | 1768620183 | NC               | NC                 | M     |

**QC Flag Legend**

Processing Flags

NC - Not Calibrated

Review Flags

M - Manually Integrated

**Reagents:**


IC LCS\_01954

Amount Added: 10.00

Units: mL

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230519-055741.d  
Injection Date: 19-May-2023 05:42:00 Instrument ID: WC\_IonChrom10 Operator ID: wetchemd  
Lims ID: ccv Worklist Smp#: 72  
Client ID:  
Injection Vol: 5.0 ul Dil. Factor: 1.0000 ALS Bottle#: 0  
Method: Anions\_IC10 Limit Group: Wet - Anions



## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230519-055741.d  
Injection Date: 19-May-2023 05:42:00 Instrument ID: WC\_IonChrom10  
Lims ID: ccv  
Client ID:  
Operator ID: wetchemd ALS Bottle#: 0 Worklist Smp#: 72  
Injection Vol: 5.0 ul Dil. Factor: 1.0000  
Method: Anions\_IC10 Limit Group: Wet - Anions  
Column: Detector Info 2\_091554\_1

**5 Nitrate as N, CAS: 14797-55-8**

Signal: 1


RT: 6.91  
Area: 119310357  
Amount: 2.677645  
Amount Units: ug/ml

## Processing Integration Results



RT: 6.91  
Area: 119161316  
Amount: 2.674332  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LVW8, 19-May-2023 12:11:59 -06:00:00 (UTC)

Audit Action: Manually Integrated

Audit Reason: Baseline Smoothing

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230519-063  
 Lims ID: ccb  
 Client ID:  
 Sample Type: CCB  
 Inject. Date: 19-May-2023 05:57:00 ALS Bottle#: 0 Worklist Smp#: 73  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121594-073  
 Misc. Info.: 280-0121594-073  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
 Limit Group: Wet - Anions  
 Last Update: 19-May-2023 12:32:15 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-132  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1603

First Level Reviewer: LVW8

Date: 19-May-2023 11:14:15

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 1 Fluoride            |              | 3.145            |                  |          |                  | ND                 |       |
| 2 Chloride            | 4.307        | 4.313            | -0.006           | 7258099  |                  | NC                 | M     |
| 3 Nitrite as N        | 5.010        | 4.998            | 0.012            | 183242   |                  | 0.0268             | M     |
| 4 Bromide             | 6.160        | 6.163            | -0.003           | 166879   |                  | NC                 | M     |
| 5 Nitrate as N        | 7.003        | 6.913            | 0.090            | 191074   |                  | 0.0300             | M     |
| 6 Orthophosphate as P | 9.042        | 9.127            | -0.085           | 10457998 |                  | 0.1928             |       |
| 7 Sulfate             | 10.243       | 10.300           | -0.057           | 11421829 |                  | NC                 | M     |

**QC Flag Legend**

Processing Flags

NC - Not Calibrated

Review Flags

M - Manually Integrated

Report Date: 19-May-2023 12:32:16

Chrom Revision: 2.3 16-May-2023 16:37:30

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230519-061244.d

Injection Date: 19-May-2023 05:57:00

Instrument ID: WC\_IonChrom10

Operator ID: wetchemd

Lims ID: ccb

Worklist Smp#: 73

Client ID:


Injection Vol: 5.0 ul

Dil. Factor: 1.0000

ALS Bottle#: 0

Method: Anions\_IC10

Limit Group: Wet - Anions





## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230519-061244.d  
Injection Date: 19-May-2023 05:57:00 Instrument ID: WC\_IonChrom10  
Lims ID: ccb  
Client ID:  
Operator ID: wetchemd ALS Bottle#: 0 Worklist Smp#: 73  
Injection Vol: 5.0 ul Dil. Factor: 1.0000  
Method: Anions\_IC10 Limit Group: Wet - Anions  
Column: Detector Info 2\_091554\_1

**5 Nitrate as N, CAS: 14797-55-8**

Signal: 1


RT: 7.01  
Area: 293258  
Amount: 0.032280  
Amount Units: ug/ml

## Processing Integration Results



RT: 7.00  
Area: 191074  
Amount: 0.030009  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LVW8, 19-May-2023 12:12:21 -06:00:00 (UTC)

Audit Action: Manually Integrated

Audit Reason: Baseline Smoothing

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230519-065  
 Lims ID: 280-176674-A-9  
 Client ID: LL12mw-245-230401-GW  
 Sample Type: Client  
 Inject. Date: 19-May-2023 06:42:00 ALS Bottle#: 0 Worklist Smp#: 76  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121594-076  
 Misc. Info.: 280-0121594-076  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
 Limit Group: Wet - Anions  
 Last Update: 19-May-2023 12:32:15 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-132  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1603

First Level Reviewer: LVW8

Date: 19-May-2023 12:13:21

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response   | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|------------|--------------------|-------|
| 1 Fluoride            |              | 3.145            |                  |            | ND                 |       |
| 2 Chloride            | 4.342        | 4.313            | 0.029            | 185834867  | NC                 | M     |
| 3 Nitrite as N        |              | 4.998            |                  |            | ND                 |       |
| 4 Bromide             | 6.292        | 6.163            | 0.129            | 1063857    | NC                 | M     |
| 5 Nitrate as N        | 7.090        | 6.913            | 0.177            | 938335     | 0.0466             | M     |
| 6 Orthophosphate as P | 9.183        | 9.127            | 0.056            | 4723047    | -0.0533            |       |
| 7 Sulfate             | 10.220       | 10.300           | -0.080           | 3278605455 | NC                 | M     |

**QC Flag Legend**

Processing Flags

NC - Not Calibrated

Review Flags

M - Manually Integrated

Report Date: 19-May-2023 12:32:18

Chrom Revision: 2.3 16-May-2023 16:37:30

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230519-065750.d

Injection Date: 19-May-2023 06:42:00

Instrument ID: WC\_IonChrom10

Operator ID: wetchemd

Lims ID: 280-176674-A-9

Lab Sample ID: 280-176674-9

Worklist Smp#: 76

Client ID: LL12mw-245-230401-GW


Injection Vol: 5.0 ul

Dil. Factor: 1.0000

ALS Bottle#: 0

Method: Anions\_IC10

Limit Group: Wet - Anions



## Eurofins Denver


|                 |                                                                                                     |                |                 |                |    |
|-----------------|-----------------------------------------------------------------------------------------------------|----------------|-----------------|----------------|----|
| Data File:      | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230519-065750.d |                |                 |                |    |
| Injection Date: | 19-May-2023 06:42:00                                                                                | Instrument ID: | WC_IonChrom10   |                |    |
| Lims ID:        | 280-176674-A-9                                                                                      | Lab Sample ID: | 280-176674-9    |                |    |
| Client ID:      | LL12mw-245-230401-GW                                                                                |                |                 |                |    |
| Operator ID:    | wetchemd                                                                                            | ALS Bottle#:   | 0               | Worklist Smp#: | 76 |
| Injection Vol:  | 5.0 ul                                                                                              | Dil. Factor:   | 1.0000          |                |    |
| Method:         | Anions_IC10                                                                                         | Limit Group:   | Wet - Anions    |                |    |
| Column:         |                                                                                                     | Detector       | Info 2_091554_1 |                |    |

**5 Nitrate as N, CAS: 14797-55-8**

Signal: 1


RT: 7.09  
Area: 2335704  
Amount: 0.077677  
Amount Units: ug/ml

## Processing Integration Results



RT: 7.09  
Area: 938335  
Amount: 0.046618  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LVW8, 19-May-2023 12:13:13 -06:00:00 (UTC)

Audit Action: Manually Integrated

Audit Reason: Baseline Smoothing

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230519-085  
 Lims ID: ccv  
 Client ID:  
 Sample Type: CCV  
 Inject. Date: 19-May-2023 08:42:00 ALS Bottle#: 0 Worklist Smp#: 84  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121594-084  
 Misc. Info.: 280-0121594-084  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Sublist: chrom-Anions\_IC10\*sub5  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
 Limit Group: Wet - Anions  
 Last Update: 19-May-2023 12:32:23 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13

Column 1 : Det: Info 2\_091554\_1

Process Host: CTX1603

First Level Reviewer: XAY4

Date: 19-May-2023 09:48:49

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response  | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|-----------|------------------|--------------------|-------|
| 1 Fluoride            |              | 3.145            |                  |           | NC               | ND                 |       |
| 2 Chloride            | 4.315        | 4.313            | 0.002            | 990692167 | NC               | NC                 |       |
| 3 Nitrite as N        | 5.005        | 4.998            | 0.007            | 111229675 | 2.50             | 2.62               | M     |
| 4 Bromide             | 6.135        | 6.163            | -0.028           | 19447213  | NC               | NC                 |       |
| 5 Nitrate as N        | 6.910        | 6.913            | -0.003           | 113790663 | 2.50             | 2.55               |       |
| 6 Orthophosphate as P | 9.018        | 9.127            | -0.109           | 54727773  | 2.50             | 2.09               |       |
| 7 Sulfate             | 10.210       | 10.300           | -0.090           | 953187267 | NC               | NC                 |       |

### QC Flag Legend

Processing Flags

NC - Not Calibrated

Review Flags

M - Manually Integrated

### Reagents:


IC LCS\_01954

Amount Added: 10.00

Units: mL

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230519-085748.d  
Injection Date: 19-May-2023 08:42:00 Instrument ID: WC\_IonChrom10 Operator ID: wetchemd  
Lims ID: ccv Worklist Smp#: 84  
Client ID:  
Injection Vol: 5.0 ul Dil. Factor: 1.0000 ALS Bottle#: 0  
Method: Anions\_IC10 Limit Group: Wet - Anions



Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230519-09  
Lims ID: ccb  
Client ID:  
Sample Type: CCB  
Inject. Date: 19-May-2023 08:57:00 ALS Bottle#: 0 Worklist Smp#: 85  
Injection Vol: 5.0 ul Dil. Factor: 1.0000  
Sample Info: 280-0121594-085  
Misc. Info.: 280-0121594-085  
Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
Limit Group: Wet - Anions  
Last Update: 19-May-2023 12:32:23 Calib Date: 18-May-2023 13:10:00  
Integrator: Falcon  
Quant Method: External Standard Quant By: Initial Calibration  
Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13  
Column 1 : Det: Info 2\_091554\_1  
Process Host: CTX1603

First Level Reviewer: LVW8

Date: 19-May-2023 11:14:37

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 1 Fluoride            |              | 3.145            |                  |          |                  | ND                 |       |
| 2 Chloride            | 4.305        | 4.313            | -0.008           | 7181062  |                  | NC                 | M     |
| 3 Nitrite as N        | 5.008        | 4.998            | 0.010            | 167665   |                  | 0.0265             | M     |
| 4 Bromide             | 6.162        | 6.163            | -0.001           | 168568   |                  | NC                 | M     |
| 5 Nitrate as N        | 6.997        | 6.913            | 0.084            | 103094   |                  | 0.0281             | M     |
| 6 Orthophosphate as P | 9.037        | 9.127            | -0.090           | 7455528  |                  | 0.0640             |       |
| 7 Sulfate             | 10.245       | 10.300           | -0.055           | 10601160 |                  | NC                 | M     |

## QC Flag Legend

Processing Flags

NC - Not Calibrated

Review Flags

M - Manually Integrated

Report Date: 19-May-2023 12:32:24

Chrom Revision: 2.3 16-May-2023 16:37:30

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230519-091245.d

Injection Date: 19-May-2023 08:57:00

Instrument ID: WC\_IonChrom10

Operator ID: wetchemd

Lims ID: ccb

Worklist Smp#: 85

Client ID:


Injection Vol: 5.0 ul

Dil. Factor: 1.0000

ALS Bottle#: 0

Method: Anions\_IC10

Limit Group: Wet - Anions





## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230519-091245.d  
Injection Date: 19-May-2023 08:57:00 Instrument ID: WC\_IonChrom10  
Lims ID: ccb  
Client ID:  
Operator ID: wetchemd ALS Bottle#: 0 Worklist Smp#: 85  
Injection Vol: 5.0 ul Dil. Factor: 1.0000  
Method: Anions\_IC10 Limit Group: Wet - Anions  
Column: Detector Info 2\_091554\_1

**5 Nitrate as N, CAS: 14797-55-8**

Signal: 1


RT: 7.00  
Area: 187510  
Amount: 0.029930  
Amount Units: ug/ml

## Processing Integration Results



RT: 7.00  
Area: 103094  
Amount: 0.028053  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LVW8, 19-May-2023 12:14:03 -06:00:00 (UTC)

Audit Action: Manually Integrated

Audit Reason: Baseline Smoothing

# IC Instrument Information

 WL: 121704 Inst ID: IL10 Analysis Date: 5/22/23 Analyst: MEC

| Rush                     | Job No.       | Samples  | Anions                         | QC Req | HT Exp |
|--------------------------|---------------|----------|--------------------------------|--------|--------|
| <input type="checkbox"/> | <u>176839</u> | <u>1</u> | F <u>Cl</u> NO2 Br NO3 PO4 SO4 | MS/D   |        |
| <input type="checkbox"/> | <u>176725</u> | <u>1</u> | F Cl NO2 Br NO3 PO4 SO4        | MS/D   |        |
| <input type="checkbox"/> | <u>176674</u> | <u>1</u> | F Cl NO2 Br NO3 PO4 SO4        | MS/D   |        |
| <input type="checkbox"/> | <u>176676</u> | <u>1</u> | F Cl NO2 Br NO3 PO4 SO4        | MS/D   |        |
| <input type="checkbox"/> | <u>176678</u> | <u>2</u> | F Cl NO2 Br NO3 PO4 SO4        | MS/D   |        |
| <input type="checkbox"/> | <u>176715</u> | <u>2</u> | F Cl NO2 Br NO3 PO4 SO4        | MS/D   |        |
| <input type="checkbox"/> | <u>176737</u> | <u>1</u> | F Cl NO2 Br NO3 PO4 SO4        | MS/D   |        |
| <input type="checkbox"/> | <u>176111</u> | <u>5</u> | F Cl NO2 Br NO3 PO4 SO4        | MS/D   |        |
| <input type="checkbox"/> | <u>176725</u> |          | F Cl NO2 Br NO3 PO4 SO4        | MS/D   |        |
| <input type="checkbox"/> | <u>176859</u> | <u>1</u> | F Cl NO2 Br NO3 PO4 SO4        | MS/D   |        |
| <input type="checkbox"/> |               |          | F Cl NO2 Br NO3 PO4 SO4        | MS/D   |        |
| <input type="checkbox"/> |               |          | F Cl NO2 Br NO3 PO4 SO4        | MS/D   |        |
| <input type="checkbox"/> |               |          | F Cl NO2 Br NO3 PO4 SO4        | MS/D   |        |
| <input type="checkbox"/> |               |          | F Cl NO2 Br NO3 PO4 SO4        | MS/D   |        |

## Dilutions

| Job No.    | Samples  | Anions                  | Dilution                            | Reason        |
|------------|----------|-------------------------|-------------------------------------|---------------|
| <u>839</u> | <u>1</u> | F Cl NO2 Br NO3 PO4 SO4 | <u>100x</u>                         | <u>test</u>   |
| <u>674</u> | <u>1</u> | F Cl NO2 Br NO3 PO4 SO4 | <u>2000x</u>                        | <u>test</u>   |
| <u>676</u> | <u>1</u> | F Cl NO2 Br NO3 PO4 SO4 | <u>20x</u>                          | <u>test</u>   |
| <u>678</u> | <u>2</u> | F Cl NO2 Br NO3 PO4 SO4 | <u>2x</u>                           | <u>test</u>   |
| <u>715</u> | <u>2</u> | F Cl NO2 Br NO3 PO4 SO4 | <u>2.5x</u>                         | <u>test</u>   |
| <u>737</u> | <u>1</u> | F Cl NO2 Br NO3 PO4 SO4 | <u>100x</u>                         | <u>test</u>   |
| <u>111</u> | <u>5</u> | F Cl NO2 Br NO3 PO4 SO4 | <u>20, 50, 100, 200, 500, 1000x</u> | <u>- test</u> |
|            |          | F Cl NO2 Br NO3 PO4 SO4 |                                     |               |
|            |          | F Cl NO2 Br NO3 PO4 SO4 |                                     |               |

Eurofins Environment Testing America  
Initial Calibration Report

Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230519-121638.b\Anions\_IC10.m

Instrument: WC\_IonChrom10

Lims Location: 280

Lock State: Initial Calib Locked

Cpnd Order: Retention Time

Integrator: Falcon

Last Modified: 19-May-2023 12:31:06

No. Compounds: 7

Sublist: chrom-Anions\_IC10\*sub5

Limit Group: Wet - Anions

## Detectors

Detector: 1, Info 2\_091554\_1

Data Type: ic Spec Type: none

Supports Extracted Chromatograms: False

Run Time: 0.000-14.900

No. Points: 7196

## Calibration File Names

|          |                                                                                                  |
|----------|--------------------------------------------------------------------------------------------------|
| Level: 1 | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230518-12110 |
|          | Inj Date: 18-May-2023 11:56:00 Worklist: 121594 Sample#: 1                                       |
| Level: 2 | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230518-12255 |
|          | Inj Date: 18-May-2023 12:11:00 Worklist: 121594 Sample#: 2                                       |
| Level: 3 | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230518-12405 |
|          | Inj Date: 18-May-2023 12:26:00 Worklist: 121594 Sample#: 3                                       |
| Level: 4 | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230518-12555 |
|          | Inj Date: 18-May-2023 12:40:00 Worklist: 121594 Sample#: 4                                       |
| Level: 5 | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230518-13105 |
|          | Inj Date: 18-May-2023 12:55:00 Worklist: 121594 Sample#: 5                                       |
| Level: 6 | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230518-13255 |
|          | Inj Date: 18-May-2023 13:10:00 Worklist: 121594 Sample#: 6                                       |

Start Cal Date: 18-May-2023 11:56:00

End Cal Date: 19-May-2023 12:01:00

## Individual Compound Calibration Parameters

Quant Method: ESTD

RF Calibration: Replace

Rule Name: Linear1

Curve: Linear

Weighting: Conc

Origin: None

Error: raw\_COD

Error Limit: 1.00

RF %Dif: 0.0

SPCC Limit: 0.0

CCC Limit: 0.0

Dependent Variable: Resp

Number of Compounds: 3

| RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | b | m1 | m2 | Error |
|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|---|----|----|-------|
|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|---|----|----|-------|

| RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | b | m1 | m2 | Error |
|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|---|----|----|-------|
|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|---|----|----|-------|

## 3 Nitrite as N

Signal: 1

|             |             |            |            |            |  |         |          |  |       |
|-------------|-------------|------------|------------|------------|--|---------|----------|--|-------|
| 35311960    | 42799726    | 43146865   | 42334489   | 42417361   |  | -967573 |          |  | 1.000 |
| 0.200000(2) | 0.500000(3) | 1.0000 (4) | 2.5000 (5) | 5.0000 (6) |  |         | 42866306 |  |       |
| 7062392     | 21399863    | 43146865   | 105836223  | 212086804  |  |         |          |  |       |
| 121594(2)   | 121594(3)   | 121594(4)  | 121594(5)  | 121594(6)  |  |         |          |  |       |

## 5 Nitrate as N

Signal: 1

|             |             |            |            |            |  |          |          |  |       |
|-------------|-------------|------------|------------|------------|--|----------|----------|--|-------|
| 37055495    | 45958086    | 42854846   | 44285580   | 44832252   |  | -1159043 |          |  | 0.999 |
| 0.200000(2) | 0.500000(3) | 1.0000 (4) | 2.5000 (5) | 5.0000 (6) |  |          | 44990806 |  |       |
| 7411099     | 22979043    | 42854846   | 110713949  | 224161262  |  |          |          |  |       |
| 121594(2)   | 121594(3)   | 121594(4)  | 121594(5)  | 121594(6)  |  |          |          |  |       |

## 6 Orthophosphate as P

Signal: 1

|             |             |            |            |            |  |         |          |  |        |
|-------------|-------------|------------|------------|------------|--|---------|----------|--|--------|
| 43303895    | 41843182    | 33022328   | 26624129   | 23003691   |  | 5964839 | R1       |  | *0.980 |
| 0.200000(2) | 0.500000(3) | 1.0000 (4) | 2.5000 (5) | 5.0000 (6) |  |         | 23299922 |  |        |
| 8660779     | 20921591    | 33022328   | 66560323   | 115018456  |  |         |          |  |        |
| 121594(2)   | 121594(3)   | 121594(4)  | 121594(5)  | 121594(6)  |  |         |          |  |        |

ICalib Error Legend

R1, Curve Fit Fail Error Limit Test

Eurofins Environment Testing America  
Initial Calibration Report

Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230519-121638.b\Anions\_IC10.m

Instrument: WC\_IonChrom10

Lims Location: 280

Lock State: Initial Calib Locked

Cpnd Order: Retention Time

Integrator: Falcon

Last Modified: 19-May-2023 12:31:06

No.Compounds:7

Sublist: chrom-Anions\_IC10\*sub5

Limit Group: Wet - Anions 28D

## Detectors

Detector: 1, Info 2\_091554\_1

Data Type: ic Spec Type: none

Supports Extracted Chromatograms: False

Run Time: 0.000-14.900

No. Points: 7196

## Calibration File Names

|          |                                                                                                  |
|----------|--------------------------------------------------------------------------------------------------|
| Level: 1 | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230518-12110 |
|          | Inj Date: 18-May-2023 11:56:00 Worklist: 121594 Sample#: 1                                       |
| Level: 2 | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230518-12255 |
|          | Inj Date: 18-May-2023 12:11:00 Worklist: 121594 Sample#: 2                                       |
| Level: 3 | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230518-12405 |
|          | Inj Date: 18-May-2023 12:26:00 Worklist: 121594 Sample#: 3                                       |
| Level: 4 | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230518-12555 |
|          | Inj Date: 18-May-2023 12:40:00 Worklist: 121594 Sample#: 4                                       |
| Level: 5 | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230518-13105 |
|          | Inj Date: 18-May-2023 12:55:00 Worklist: 121594 Sample#: 5                                       |
| Level: 6 | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230518-13255 |
|          | Inj Date: 18-May-2023 13:10:00 Worklist: 121594 Sample#: 6                                       |

Start Cal Date: 18-May-2023 11:56:00

End Cal Date: 19-May-2023 12:01:00

## Individual Compound Calibration Parameters

Quant Method: ESTD

RF Calibration: Replace

Rule Name: Linear1

Curve: Linear

Weighting: Conc

Origin: None

Error: raw\_COD

Error Limit: 1.00

RF %Dif: 0.0

SPCC Limit: 0.0

CCC Limit: 0.0

Dependent Variable: Resp

Number of Compounds: 4

| RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | b | m1 | m2 | Error |
|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|---|----|----|-------|
|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|---|----|----|-------|

| RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | b         | m1       | m2 | Error |
|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|-----------|----------|----|-------|
| 1 Fluoride                         |                                    |                                    |                                    |                                    |                                    | Signal: 1 |          |    |       |
| 29175985                           | 36359510                           | 39638142                           | 40984560                           | 40102056                           |                                    | -2208602  |          |    | 1.000 |
| 0.200000(2)                        | 0.500000(3)                        | 1.0000 (4)                         | 2.5000 (5)                         | 5.0000 (6)                         |                                    |           | 41050846 |    |       |
| M5835197                           | M18179755                          | M39638142                          | M102461401                         | M200510278                         |                                    |           |          |    |       |
| 121594(2)                          | 121594(3)                          | 121594(4)                          | 121594(5)                          | 121594(6)                          |                                    |           |          |    |       |
| 2 Chloride                         |                                    |                                    |                                    |                                    |                                    | Signal: 1 |          |    |       |
| 16780722                           | 16174841                           | 19339456                           | 19009666                           | 18725386                           | 18640004                           | -1171132  |          |    | 1.000 |
| 0.500000(1)                        | 1.2500 (2)                         | 2.5000 (3)                         | 30.0 (4)                           | 60.0 (5)                           | 100.0 (6)                          |           | 18747994 |    |       |
| 8390361                            | 20218551                           | 48348641                           | 570289992                          | 1123523133                         | 1864000362                         |           |          |    |       |
| 121594(1)                          | 121594(2)                          | 121594(3)                          | 121594(4)                          | 121594(5)                          | 121594(6)                          |           |          |    |       |
| 4 Bromide                          |                                    |                                    |                                    |                                    |                                    | Signal: 1 |          |    |       |
| 11277385                           | 9537194                            | 7900278                            | 7592206                            | 7648936                            |                                    | 806494    |          |    | 0.999 |
| 0.200000(2)                        | 0.500000(3)                        | 1.0000 (4)                         | 2.5000 (5)                         | 5.0000 (6)                         |                                    |           | 7404030  |    |       |
| 2255477                            | 4768597                            | 7900278                            | 18980516                           | 38244681                           |                                    |           |          |    |       |
| 121594(2)                          | 121594(3)                          | 121594(4)                          | 121594(5)                          | 121594(6)                          |                                    |           |          |    |       |
| 7 Sulfate                          |                                    |                                    |                                    |                                    |                                    | Signal: 1 |          |    |       |
| 22625290                           | 16020812                           | 17380170                           | 14404666                           | 13852661                           | 13712274                           | 4654296   |          |    | 0.999 |
| 0.500000(1)                        | 1.2500 (2)                         | 2.5000 (3)                         | 30.0 (4)                           | 60.0 (5)                           | 100.0 (6)                          |           | 13803812 |    |       |
| 11312645                           | 20026015                           | 43450426                           | 432139991                          | 831159688                          | 1371227441                         |           |          |    |       |
| 121594(1)                          | 121594(2)                          | 121594(3)                          | 121594(4)                          | 121594(5)                          | 121594(6)                          |           |          |    |       |

## Preliminary Report

Eurofins Denver  
ICV, ICal Verification Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13  
Lims ID: ICV  
Client ID:  
Sample Type: ICV  
Inject. Date: 18-May-2023 13:25:00 ALS Bottle#: 0 Worklist Smp#: 7  
Injection Vol: 5.0 ul Dil. Factor: 1.0000  
Sample Info: 280-0121594-007  
Misc. Info.: 280-0121594-007  
Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
Sublist: chrom-Anions\_IC7\*sub4

Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
Limit Group: Wet - Anions 28D  
Last Update: 19-May-2023 12:32:45 Calib Date: 18-May-2023 13:10:00  
Integrator: Falcon  
Quant Method: External Standard Quant By: Initial Calibration  
Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13

Column 1 : Det: Info 2\_091554\_1

Process Host: CTX1603

First Level Reviewer: LVW8

Date: 19-May-2023 12:38:30

Start Cal Date: 18-May-2023 11:56:00

End Cal Date: 18-May-2023 13:10:00

| Compound   | Amount Added | Amount Detected | %Drift | Max. %Drift | %Rec  | %Rec Limits |
|------------|--------------|-----------------|--------|-------------|-------|-------------|
| 1 Fluoride | 2.00         | 2.34            | * 16.8 | 10.0        | 116.8 |             |
| 2 Chloride | 80.0         | 83.0            | 3.8    | 10.0        | 103.8 |             |
| 4 Bromide  | 2.00         | 1.91            | -4.7   | 10.0        | 95.3  |             |
| 7 Sulfate  | 80.0         | 81.8            | 2.3    | 10.0        | 102.3 |             |

## Preliminary Report

Eurofins Denver  
ICV, ICal Verification Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13  
Lims ID: ICV  
Client ID:  
Sample Type: ICV  
Inject. Date: 18-May-2023 13:25:00 ALS Bottle#: 0 Worklist Smp#: 7  
Injection Vol: 5.0 ul Dil. Factor: 1.0000  
Sample Info: 280-0121594-007  
Misc. Info.: 280-0121594-007  
Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
Sublist: chrom-Anions\_IC7\*sub4  
Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
Limit Group: Wet - Anions  
Last Update: 19-May-2023 12:38:46 Calib Date: 18-May-2023 13:10:00  
Integrator: Falcon  
Quant Method: External Standard Quant By: Initial Calibration  
Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13  
Column 1 : Det: Info 2\_091554\_1  
Process Host: CTX1603  
First Level Reviewer: LVW8 Date: 19-May-2023 12:38:46  
Start Cal Date: 18-May-2023 11:56:00  
End Cal Date: 18-May-2023 13:10:00

| Compound            | Amount Added | Amount Detected | %Drift | Max. %Drift | %Rec  | %Rec Limits |
|---------------------|--------------|-----------------|--------|-------------|-------|-------------|
| 3 Nitrite as N      | 2.00         | 1.99            | -0.3   | 10.0        | 99.7  |             |
| 5 Nitrate as N      | 2.00         | 1.95            | -2.3   | 10.0        | 97.7  |             |
| 6 Orthophosphate as | 2.00         | 2.30            | * 14.8 | 10.0        | 114.8 |             |



Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230522-121704.b\Info 2\_DENPC179\_Anions\_20230522-163  
 Lims ID: ccv  
 Client ID:  
 Sample Type: CCV  
 Inject. Date: 22-May-2023 16:00:00 ALS Bottle#: 0 Worklist Smp#: 1  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121704-001  
 Misc. Info.: 280-0121704-001  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Sublist: chrom-Anions\_IC10\*sub5  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230522-121704.b\Anions\_IC10.m  
 Limit Group: Wet - Anions  
 Last Update: 23-May-2023 10:59:14 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1623

First Level Reviewer: LVW8

Date: 22-May-2023 16:19:27

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response  | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|-----------|------------------|--------------------|-------|
| 1 Fluoride            |              | 3.140            |                  |           | ND               | ND                 |       |
| 2 Chloride            | 4.315        | 4.315            | 0.000            | 926518595 | NC               | NC                 |       |
| 3 Nitrite as N        | 4.998        | 4.998            | 0.000            | 103103875 | 2.50             | 2.43               |       |
| 4 Bromide             | 6.152        | 6.152            | 0.000            | 18468888  | NC               | NC                 |       |
| 5 Nitrate as N        | 6.905        | 6.905            | 0.000            | 108362993 | 2.50             | 2.43               |       |
| 6 Orthophosphate as P | 9.155        | 9.155            | 0.000            | 2988274   | 2.50             | -0.1278            |       |
| 7 Sulfate             | 10.285       | 10.285           | 0.000            | 703902223 | NC               | NC                 |       |

**QC Flag Legend**

Processing Flags

NC - Not Calibrated

ND - Not Detected or Marked ND

**Reagents:**

IC LCS\_01954

Amount Added: 10.00

Units: mL

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230522-121704.b\Info 2\_DENPC179\_Anions\_20230522-162  
 Lims ID: ccb  
 Client ID:  
 Sample Type: CCB  
 Inject. Date: 22-May-2023 16:14:00 ALS Bottle#: 0 Worklist Smp#: 2  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121704-002  
 Misc. Info.: 280-0121704-002  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230522-121704.b\Anions\_IC10.m  
 Limit Group: Wet - Anions  
 Last Update: 23-May-2023 10:59:15 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-132  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1623

First Level Reviewer: XAY4

Date: 23-May-2023 10:20:49

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 1 Fluoride            |              | 3.140            |                  |          |                  | ND                 |       |
| 2 Chloride            |              | 4.315            |                  |          |                  | ND                 |       |
| 3 Nitrite as N        | 4.762        | 4.762            | 0.000            | 67101    |                  | 0.0241             |       |
| 4 Bromide             |              | 6.152            |                  |          |                  | ND                 | U     |
| 5 Nitrate as N        |              | 6.905            |                  |          |                  | ND                 |       |
| 6 Orthophosphate as P | 9.140        | 9.140            | 0.000            | 4188666  |                  | -0.0762            |       |
| 7 Sulfate             | 10.360       | 10.360           | 0.000            | 3894449  |                  | NC                 |       |

**QC Flag Legend**

Processing Flags

NC - Not Calibrated

Review Flags

U - Marked Undetected

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230522-121704.b\Info 2\_DENPC179\_Anions\_20230522-16  
 Lims ID: mrl  
 Client ID:  
 Sample Type: MRL  
 Inject. Date: 22-May-2023 16:29:00 ALS Bottle#: 0 Worklist Smp#: 3  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121704-003  
 Misc. Info.: 280-0121704-003  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230522-121704.b\Anions\_IC10.m  
 Limit Group: Wet - Anions  
 Last Update: 23-May-2023 10:59:16 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1623

First Level Reviewer: LVW8

Date: 22-May-2023 17:51:11

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 1 Fluoride            |              | 3.140            |                  |          | ND               | ND                 |       |
| 2 Chloride            | 4.302        | 4.302            | 0.000            | 43970753 | NC               | NC                 |       |
| 3 Nitrite as N        | 5.002        | 5.002            | 0.000            | 8614732  | 0.2500           | 0.2235             |       |
| 4 Bromide             | 6.172        | 6.172            | 0.000            | 2650290  | NC               | NC                 |       |
| 5 Nitrate as N        | 6.953        | 6.953            | 0.000            | 9581243  | 0.2500           | 0.2387             |       |
| 6 Orthophosphate as P | 9.157        | 9.157            | 0.000            | 2494989  | 0.2500           | -0.1489            |       |
| 7 Sulfate             | 10.305       | 10.305           | 0.000            | 43075457 | NC               | NC                 |       |

**QC Flag Legend**

Processing Flags

NC - Not Calibrated

ND - Not Detected or Marked ND

**Reagents:**

IC CAL cl/so4\_00480

Amount Added: 0.20

Units: mL

IC Cal low\_00709

Amount Added: 0.10

Units: mL

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230522-121704.b\Info 2\_DENPC179\_Anions\_20230522-164450.d

Injection Date: 22-May-2023 16:29:00

Instrument ID: WC\_IonChrom10

Operator ID: wetchemd

Lims ID: mrl

Worklist Smp#: 3

Client ID:


Injection Vol: 5.0 ul

Dil. Factor: 1.0000

ALS Bottle#: 0

Method: Anions\_IC10

Limit Group: Wet - Anions



Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230522-121704.b\Info 2\_DENPC179\_Anions\_20230522-16  
 Lims ID: lcs  
 Client ID:  
 Sample Type: LCS  
 Inject. Date: 22-May-2023 16:44:00 ALS Bottle#: 0 Worklist Smp#: 4  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121704-004  
 Misc. Info.: 280-0121704-004  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230522-121704.b\Anions\_IC10.m  
 Limit Group: Wet - Anions  
 Last Update: 23-May-2023 10:59:17 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1623

First Level Reviewer: LVW8

Date: 22-May-2023 17:51:17

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response  | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|-----------|------------------|--------------------|-------|
| 1 Fluoride            |              | 3.140            |                  |           | ND               | ND                 |       |
| 2 Chloride            | 4.317        | 4.317            | 0.000            | 922923313 | NC               | NC                 |       |
| 3 Nitrite as N        | 4.998        | 4.998            | 0.000            | 103618457 | 2.50             | 2.44               |       |
| 4 Bromide             | 6.153        | 6.153            | 0.000            | 18372126  | NC               | NC                 |       |
| 5 Nitrate as N        | 6.907        | 6.907            | 0.000            | 107220163 | 2.50             | 2.41               |       |
| 6 Orthophosphate as P | 9.163        | 9.163            | 0.000            | 1342060   | 2.50             | -0.1984            |       |
| 7 Sulfate             | 10.288       | 10.288           | 0.000            | 697393055 | NC               | NC                 |       |

**QC Flag Legend**

Processing Flags

NC - Not Calibrated

ND - Not Detected or Marked ND

**Reagents:**

IC LCS\_01954

Amount Added: 10.00

Units: mL

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230522-121704.b\Info 2\_DENPC179\_Anions\_20230522-17  
 Lims ID: lcsd  
 Client ID:  
 Sample Type: LCSD  
 Inject. Date: 22-May-2023 16:59:00 ALS Bottle#: 0 Worklist Smp#: 5  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121704-005  
 Misc. Info.: 280-0121704-005  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230522-121704.b\Anions\_IC10.m  
 Limit Group: Wet - Anions  
 Last Update: 23-May-2023 10:59:18 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1623

First Level Reviewer: LVW8

Date: 22-May-2023 17:51:22

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response  | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|-----------|------------------|--------------------|-------|
| 1 Fluoride            |              | 3.140            |                  |           | ND               | ND                 |       |
| 2 Chloride            | 4.313        | 4.313            | 0.000            | 926415322 | NC               | NC                 |       |
| 3 Nitrite as N        | 4.997        | 4.997            | 0.000            | 103572545 | 2.50             | 2.44               |       |
| 4 Bromide             | 6.150        | 6.150            | 0.000            | 18444382  | NC               | NC                 |       |
| 5 Nitrate as N        | 6.903        | 6.903            | 0.000            | 107568099 | 2.50             | 2.42               |       |
| 6 Orthophosphate as P | 9.147        | 9.147            | 0.000            | 2150895   | 2.50             | -0.1637            |       |
| 7 Sulfate             | 10.287       | 10.287           | 0.000            | 696494797 | NC               | NC                 |       |

**QC Flag Legend**

Processing Flags

NC - Not Calibrated

ND - Not Detected or Marked ND

**Reagents:**

IC LCS\_01954

Amount Added: 10.00

Units: mL

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230522-121704.b\Info 2\_DENPC179\_Anions\_20230522-171447.d

Injection Date: 22-May-2023 16:59:00

Instrument ID: WC\_IonChrom10

Operator ID: wetchemd

Lims ID: lcsd

Worklist Smp#: 5

Client ID:


Injection Vol: 5.0 ul

Dil. Factor: 1.0000

ALS Bottle#: 0

Method: Anions\_IC10

Limit Group: Wet - Anions



Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230522-121704.b\Info 2\_DENPC179\_Anions\_20230522-172  
 Lims ID: mb  
 Client ID:  
 Sample Type: MB  
 Inject. Date: 22-May-2023 17:14:00 ALS Bottle#: 0 Worklist Smp#: 6  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121704-006  
 Misc. Info.: 280-0121704-006  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230522-121704.b\Anions\_IC10.m  
 Limit Group: Wet - Anions  
 Last Update: 23-May-2023 10:59:18 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-132  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1623

First Level Reviewer: LVW8

Date: 22-May-2023 17:51:27

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 1 Fluoride            |              | 3.140            |                  |          |                  | ND                 |       |
| 2 Chloride            |              | 4.313            |                  |          |                  | ND                 |       |
| 3 Nitrite as N        | 4.755        | 4.755            | 0.000            | 76667    |                  | 0.0244             |       |
| 4 Bromide             | 6.325        | 6.325            | 0.000            | 701320   |                  | NC                 |       |
| 5 Nitrate as N        |              | 6.903            |                  |          |                  | ND                 |       |
| 6 Orthophosphate as P | 9.148        | 9.148            | 0.000            | 2761007  |                  | -0.1375            |       |
| 7 Sulfate             | 10.372       | 10.372           | 0.000            | 3219626  |                  | NC                 |       |

**QC Flag Legend**

Processing Flags

NC - Not Calibrated



Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230522-121704.b\Info 2\_DENPC179\_Anions\_20230522-19  
 Lims ID: 280-176674-A-3  
 Client ID: LL12mw-187-230401-GW  
 Sample Type: Client  
 Inject. Date: 22-May-2023 18:59:00 ALS Bottle#: 0 Worklist Smp#: 13  
 Injection Vol: 5.0 ul Dil. Factor: 2000.0000  
 Sample Info: 280-0121704-013  
 Misc. Info.: 280-0121704-013  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230522-121704.b\Anions\_IC10.m  
 Limit Group: Wet - Anions  
 Last Update: 23-May-2023 10:59:18 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1623

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|----------|--------------------|-------|
| 1 Fluoride            |              | 3.140            |                  |          | ND                 |       |
| 2 Chloride            | 4.300        | 4.313            | -0.013           | 322362   | NC                 |       |
| 3 Nitrite as N        | 4.983        | 4.755            | 0.228            | 160405   | 0.0263             |       |
| 4 Bromide             | 6.148        | 6.325            | -0.177           | 743682   | NC                 |       |
| 5 Nitrate as N        | 6.908        | 6.903            | 0.005            | 32807745 | 0.7550             |       |
| 6 Orthophosphate as P | 9.085        | 9.148            | -0.063           | 10142747 | 0.1793             |       |
| 7 Sulfate             | 10.263       | 10.372           | -0.109           | 7818321  | NC                 |       |

### QC Flag Legend

Processing Flags

NC - Not Calibrated

## Eurofins Denver

|                 |                                                                                                                |                |                      |                |    |
|-----------------|----------------------------------------------------------------------------------------------------------------|----------------|----------------------|----------------|----|
| Data File:      | \\chromfs\Denver\ChromData\WC_IonChrom10\20230522-121704.b\Info 2_DENPC179_Anions_20230522-20230522-121704.b   |                |                      |                |    |
| Lims ID:        | ccv                                                                                                            |                |                      |                |    |
| Client ID:      |                                                                                                                |                |                      |                |    |
| Sample Type:    | CCV                                                                                                            |                |                      |                |    |
| Inject. Date:   | 22-May-2023 19:59:00                                                                                           | ALS Bottle#:   | 0                    | Worklist Smp#: | 17 |
| Injection Vol:  | 5.0 ul                                                                                                         | Dil. Factor:   | 1.0000               |                |    |
| Sample Info:    | 280-0121704-017                                                                                                |                |                      |                |    |
| Misc. Info.:    | 280-0121704-017                                                                                                |                |                      |                |    |
| Operator ID:    | wetchemd                                                                                                       | Instrument ID: | WC_IonChrom10        |                |    |
| Sublist:        | chrom-Anions_IC10*sub5                                                                                         |                |                      |                |    |
| Method:         | \\chromfs\Denver\ChromData\WC_IonChrom10\20230522-121704.b\Anions_IC10.m                                       |                |                      |                |    |
| Limit Group:    | Wet - Anions                                                                                                   |                |                      |                |    |
| Last Update:    | 23-May-2023 10:59:26                                                                                           | Calib Date:    | 18-May-2023 13:10:00 |                |    |
| Integrator:     | Falcon                                                                                                         |                |                      |                |    |
| Quant Method:   | External Standard                                                                                              | Quant By:      | Initial Calibration  |                |    |
| Last ICal File: | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230518-1320230518-121594.b |                |                      |                |    |
| Column 1 :      | Det: Info 2_091554_1                                                                                           |                |                      |                |    |
| Process Host:   | CTX1623                                                                                                        |                |                      |                |    |

Date:

23-May-2023 10:15:37

| Compound              | RT (min.) | Exp RT (min.) | Dlt RT (min.) | Response  | Cal Amt ug/ml | OnCol Amt ug/ml | Flags |
|-----------------------|-----------|---------------|---------------|-----------|---------------|-----------------|-------|
| 1 Fluoride            |           | 3.140         |               |           | ND            | ND              |       |
| 2 Chloride            | 4.307     | 4.313         | -0.006        | 932944574 | NC            | NC              |       |
| 3 Nitrite as N        | 4.988     | 4.755         | 0.233         | 110845562 | 2.50          | 2.61            |       |
| 4 Bromide             | 6.117     | 6.325         | -0.208        | 18494660  | NC            | NC              |       |
| 5 Nitrate as N        | 6.877     | 6.903         | -0.026        | 107747893 | 2.50          | 2.42            |       |
| 6 Orthophosphate as P | 9.103     | 9.148         | -0.045        | 1912491   | 2.50          | -0.1739         |       |
| 7 Sulfate             | 10.242    | 10.372        | -0.130        | 666383474 | NC            | NC              |       |

## QC Flag Legend

## Processing Flags

NC - Not Calibrated

ND - Not Detected or Marked ND

## Reagents:

Amount Added: 10.00

Units: mL

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230522-121704.b\Info 2\_DENPC179\_Anions\_20230522-201457.d

Injection Date: 22-May-2023 19:59:00

Instrument ID: WC\_IonChrom10

Operator ID: wetchemd

Lims ID: ccv

Worklist Smp#: 17

Client ID:


Injection Vol: 5.0 ul

Dil. Factor: 1.0000

ALS Bottle#: 0

Method: Anions\_IC10

Limit Group: Wet - Anions



Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230522-121704.b\Info 2\_DENPC179\_Anions\_20230522-20230522-121704.b  
 Lims ID: ccb  
 Client ID:  
 Sample Type: CCB  
 Inject. Date: 22-May-2023 20:15:00 ALS Bottle#: 0 Worklist Smp#: 18  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121704-018  
 Misc. Info.: 280-0121704-018  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230522-121704.b\Anions\_IC10.m  
 Limit Group: Wet - Anions  
 Last Update: 23-May-2023 10:59:26 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-121594.b  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1623

First Level Reviewer: XAY4

Date: 23-May-2023 10:21:55

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 1 Fluoride            |              | 3.140            |                  |          |                  | ND                 |       |
| 2 Chloride            |              | 4.313            |                  |          |                  | ND                 |       |
| 3 Nitrite as N        |              | 4.755            |                  |          |                  | ND                 | U     |
| 4 Bromide             |              | 6.325            |                  |          |                  | ND                 | U     |
| 5 Nitrate as N        |              | 6.903            |                  |          |                  | ND                 |       |
| 6 Orthophosphate as P | 9.115        | 9.148            | -0.033           | 2602938  |                  | -0.1443            |       |
| 7 Sulfate             | 10.322       | 10.372           | -0.050           | 2963452  |                  | NC                 |       |

**QC Flag Legend**

Processing Flags

NC - Not Calibrated

Review Flags

U - Marked Undetected

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-121594.b  
 Lims ID: STD L1  
 Client ID:  
 Sample Type: IC Calib Level: 1  
 Inject. Date: 18-May-2023 11:56:00 ALS Bottle#: 0 Worklist Smp#: 1  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121594-001  
 Misc. Info.: 280-0121594-001  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Sublist: chrom-Anions\_IC10\*sub5  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
 Limit Group: Wet - Anions 28D  
 Last Update: 19-May-2023 11:53:45 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-121594.b  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1603

First Level Reviewer: LVW8

Date: 18-May-2023 12:57:44

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 1 Fluoride            |              | 3.145            |                  |          |                  | ND                 |       |
| 2 Chloride            | 4.302        | 4.302            | 0.000            | 8390361  | 0.5000           | 0.5100             |       |
| 3 Nitrite as N        |              | 4.998            |                  |          |                  | ND                 |       |
| 4 Bromide             |              | 6.305            |                  |          |                  | ND                 | U     |
| 5 Nitrate as N        |              | 6.913            |                  |          |                  | ND                 |       |
| 6 Orthophosphate as P |              | 9.127            |                  |          |                  | ND                 | U     |
| 7 Sulfate             | 10.243       | 10.243           | 0.000            | 11312645 | 0.5000           | 0.4824             |       |

**QC Flag Legend**

Processing Flags

Review Flags

U - Marked Undetected

**Reagents:**

IC CAL cl/so4\_00480

Amount Added: 0.04

Units: mL

IC Cal low\_00709

Amount Added: 0.00

Units: mL

Report Date: 19-May-2023 11:53:45

Chrom Revision: 2.3 16-May-2023 16:37:30

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-121101.d

Injection Date: 18-May-2023 11:56:00

Instrument ID: WC\_IonChrom10

Operator ID: wetchemd

Lims ID: STD L1

Worklist Smp#: 1

Client ID:


Injection Vol: 5.0 ul

Dil. Factor: 1.0000

ALS Bottle#: 0

Method: Anions\_IC10

Limit Group: Wet - Anions 28D



Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-121594.b  
 Lims ID: STD L2  
 Client ID:  
 Sample Type: IC Calib Level: 2  
 Inject. Date: 18-May-2023 12:11:00 ALS Bottle#: 0 Worklist Smp#: 2  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121594-002  
 Misc. Info.: 280-0121594-002  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Sublist: chrom-Anions\_IC10\*sub5  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
 Limit Group: Wet - Anions 28D  
 Last Update: 19-May-2023 11:53:46 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-121594.b  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1603

First Level Reviewer: LVW8

Date: 18-May-2023 12:58:14

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 1 Fluoride            | 3.133        | 3.133            | 0.000            | 5835197  | 0.2000           | 0.1959             | Ma    |
| 2 Chloride            | 4.300        | 4.300            | 0.000            | 20218551 | 1.25             | 1.14               |       |
| 3 Nitrite as N        | 5.002        | 5.002            | 0.000            | 7062392  | NC               | NC                 |       |
| 4 Bromide             | 6.165        | 6.165            | 0.000            | 2255477  | 0.2000           | 0.1957             |       |
| 5 Nitrate as N        | 6.953        | 6.953            | 0.000            | 7411099  | NC               | NC                 |       |
| 6 Orthophosphate as P | 9.093        | 9.093            | 0.000            | 8660779  | NC               | NC                 |       |
| 7 Sulfate             | 10.245       | 10.245           | 0.000            | 20026015 | 1.25             | 1.11               |       |

**QC Flag Legend**

Processing Flags

NC - Not Calibrated

Review Flags

M - Manually Integrated

a - User Assigned ID

**Reagents:**

IC CAL c/so4\_00480

Amount Added: 0.10

Units: mL

IC Cal low\_00709

Amount Added: 0.08

Units: mL

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-122558.d

Injection Date: 18-May-2023 12:11:00

Instrument ID: WC\_IonChrom10

Operator ID: wetchemd

Lims ID: STD L2

Worklist Smp#: 2

Client ID:


Injection Vol: 5.0 ul

Dil. Factor: 1.0000

ALS Bottle#: 0

Method: Anions\_IC10

Limit Group: Wet - Anions 28D





## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-121594.d  
Injection Date: 18-May-2023 12:11:00 Instrument ID: WC\_IonChrom10  
Lims ID: STD L2  
Client ID:  
Operator ID: wetchemd ALS Bottle#: 0 Worklist Smp#: 2  
Injection Vol: 5.0 ul Dil. Factor: 1.0000  
Method: Anions\_IC10 Limit Group: Wet - Anions 28D  
Column: Detector Info 2\_091554\_1

**1 Fluoride, CAS: 16984-48-8**

Signal: 1


Not Detected  
Expected RT: 3.13

## Processing Integration Results



RT: 3.13  
Area: 5835197  
Amount: 0.195947  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LVW8, 18-May-2023 12:58:04 -06:00:00 (UTC)

Audit Action: Manually Integrated/Assigned Compound ID Audit Reason: Baseline Smoothing

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-121594.b  
 Lims ID: STD L3  
 Client ID:  
 Sample Type: IC Calib Level: 3  
 Inject. Date: 18-May-2023 12:26:00 ALS Bottle#: 0 Worklist Smp#: 3  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121594-003  
 Misc. Info.: 280-0121594-003  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Sublist: chrom-Anions\_IC10\*sub5  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
 Limit Group: Wet - Anions 28D  
 Last Update: 19-May-2023 11:53:47 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-121594.b  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1603

First Level Reviewer: LVW8

Date: 18-May-2023 12:58:34

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 1 Fluoride            | 3.137        | 3.137            | 0.000            | 18179755 | 0.5000           | 0.4967             | Ma    |
| 2 Chloride            | 4.300        | 4.300            | 0.000            | 48348641 | 2.50             | 2.64               |       |
| 3 Nitrite as N        | 5.000        | 5.000            | 0.000            | 21399863 | NC               | NC                 |       |
| 4 Bromide             | 6.162        | 6.162            | 0.000            | 4768597  | 0.5000           | 0.5351             |       |
| 5 Nitrate as N        | 6.938        | 6.938            | 0.000            | 22979043 | NC               | NC                 |       |
| 6 Orthophosphate as P | 9.075        | 9.075            | 0.000            | 20921591 | NC               | NC                 |       |
| 7 Sulfate             | 10.247       | 10.247           | 0.000            | 43450426 | 2.50             | 2.81               |       |

**QC Flag Legend**

Processing Flags

NC - Not Calibrated

Review Flags

M - Manually Integrated

a - User Assigned ID

**Reagents:**

IC CAL c/so4\_00480

Amount Added: 0.20

Units: mL

IC Cal low\_00709

Amount Added: 0.20

Units: mL

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-124056.d

Injection Date: 18-May-2023 12:26:00

Instrument ID: WC\_IonChrom10

Operator ID: wetchemd

Lims ID: STD L3

Worklist Smp#: 3

Client ID:


Injection Vol: 5.0 ul

Dil. Factor: 1.0000

ALS Bottle#: 0

Method: Anions\_IC10

Limit Group: Wet - Anions 28D



## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-121594.d  
Injection Date: 18-May-2023 12:26:00 Instrument ID: WC\_IonChrom10  
Lims ID: STD L3  
Client ID:  
Operator ID: wetchemd ALS Bottle#: 0 Worklist Smp#: 3  
Injection Vol: 5.0 ul Dil. Factor: 1.0000  
Method: Anions\_IC10 Limit Group: Wet - Anions 28D  
Column: Detector Info 2\_091554\_1

**1 Fluoride, CAS: 16984-48-8**

Signal: 1


Not Detected  
Expected RT: 3.14

## Processing Integration Results



RT: 3.14  
Area: 18179755  
Amount: 0.496661  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LVW8, 18-May-2023 12:58:27 -06:00:00 (UTC)

Audit Action: Manually Integrated/Assigned Compound ID Audit Reason: Baseline Smoothing

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-121594.b  
 Lims ID: STD L4  
 Client ID:  
 Sample Type: IC Calib Level: 4  
 Inject. Date: 18-May-2023 12:40:00 ALS Bottle#: 0 Worklist Smp#: 4  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121594-004  
 Misc. Info.: 280-0121594-004  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Sublist: chrom-Anions\_IC10\*sub5  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
 Limit Group: Wet - Anions 28D  
 Last Update: 19-May-2023 11:53:47 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-121594.b  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1603

First Level Reviewer: LVW8

Date: 18-May-2023 13:29:45

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response  | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|-----------|------------------|--------------------|-------|
| 1 Fluoride            | 3.138        | 3.138            | 0.000            | 39638142  | 1.00             | 1.02               | Ma    |
| 2 Chloride            | 4.305        | 4.305            | 0.000            | 570289992 | 30.0             | 30.5               |       |
| 3 Nitrite as N        | 5.000        | 5.000            | 0.000            | 43146865  | NC               | NC                 |       |
| 4 Bromide             | 6.162        | 6.162            | 0.000            | 7900278   | 1.00             | 0.9581             |       |
| 5 Nitrate as N        | 6.928        | 6.928            | 0.000            | 42854846  | NC               | NC                 |       |
| 6 Orthophosphate as P | 9.068        | 9.068            | 0.000            | 33022328  | NC               | NC                 |       |
| 7 Sulfate             | 10.242       | 10.242           | 0.000            | 432139991 | 30.0             | 31.0               |       |

**QC Flag Legend**

Processing Flags

NC - Not Calibrated

Review Flags

M - Manually Integrated

a - User Assigned ID

**Reagents:**

IC CAL c/so4\_00480

Amount Added: 2.40

Units: mL


IC Cal low\_00709

Amount Added: 0.40

Units: mL

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-125554.d  
Injection Date: 18-May-2023 12:40:00 Instrument ID: WC\_IonChrom10 Operator ID: wetchemd  
Lims ID: STD L4 Worklist Smp#: 4  
Client ID:  
Injection Vol: 5.0 ul Dil. Factor: 1.0000 ALS Bottle#: 0  
Method: Anions\_IC10 Limit Group: Wet - Anions 28D



## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-121594.d  
Injection Date: 18-May-2023 12:40:00 Instrument ID: WC\_IonChrom10  
Lims ID: STD L4  
Client ID:  
Operator ID: wetchemd ALS Bottle#: 0 Worklist Smp#: 4  
Injection Vol: 5.0 ul Dil. Factor: 1.0000  
Method: Anions\_IC10 Limit Group: Wet - Anions 28D  
Column: Detector Info 2\_091554\_1

**1 Fluoride, CAS: 16984-48-8**

Signal: 1


Not Detected  
Expected RT: 3.14

## Processing Integration Results



RT: 3.14  
Area: 39638142  
Amount: 1.019388  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LVW8, 18-May-2023 13:29:42 -06:00:00 (UTC)

Audit Action: Manually Integrated/Assigned Compound ID Audit Reason: Baseline Smoothing

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13  
 Lims ID: STD L5  
 Client ID:  
 Sample Type: IC Calib Level: 5  
 Inject. Date: 18-May-2023 12:55:00 ALS Bottle#: 0 Worklist Smp#: 5  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121594-005  
 Misc. Info.: 280-0121594-005  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Sublist: chrom-Anions\_IC10\*sub5  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
 Limit Group: Wet - Anions 28D  
 Last Update: 19-May-2023 11:53:48 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13

Column 1 : Det: Info 2\_091554\_1

Process Host: CTX1603

First Level Reviewer: LVW8

Date: 18-May-2023 13:30:00

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response   | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|------------|------------------|--------------------|-------|
| 1 Fluoride            | 3.140        | 3.140            | 0.000            | 102461401  | 2.50             | 2.55               | Ma    |
| 2 Chloride            | 4.313        | 4.313            | 0.000            | 1123523133 | 60.0             | 60.0               |       |
| 3 Nitrite as N        | 4.995        | 4.995            | 0.000            | 105836223  | NC               | NC                 |       |
| 4 Bromide             | 6.152        | 6.152            | 0.000            | 18980516   | 2.50             | 2.45               |       |
| 5 Nitrate as N        | 6.903        | 6.903            | 0.000            | 110713949  | NC               | NC                 |       |
| 6 Orthophosphate as P | 9.043        | 9.043            | 0.000            | 66560323   | NC               | NC                 |       |
| 7 Sulfate             | 10.227       | 10.227           | 0.000            | 831159688  | 60.0             | 59.9               |       |

### QC Flag Legend

Processing Flags

NC - Not Calibrated

Review Flags

M - Manually Integrated

a - User Assigned ID

### Reagents:

IC CAL c/so4\_00480 Amount Added: 4.80 Units: mL  
 IC Cal low\_00709 Amount Added: 1.00 Units: mL



Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-131053.d

Injection Date: 18-May-2023 12:55:00

Instrument ID: WC\_IonChrom10

Operator ID: wetchemd

Lims ID: STD L5

Worklist Smp#: 5

Client ID:


Injection Vol: 5.0 ul

Dil. Factor: 1.0000

ALS Bottle#: 0

Method: Anions\_IC10

Limit Group: Wet - Anions 28D



## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-131053.d  
Injection Date: 18-May-2023 12:55:00 Instrument ID: WC\_IonChrom10  
Lims ID: STD L5  
Client ID:  
Operator ID: wetchemd ALS Bottle#: 0 Worklist Smp#: 5  
Injection Vol: 5.0 ul Dil. Factor: 1.0000  
Method: Anions\_IC10 Limit Group: Wet - Anions 28D  
Column: Detector Info 2\_091554\_1

**1 Fluoride, CAS: 16984-48-8**

Signal: 1


Not Detected  
Expected RT: 3.14

## Processing Integration Results



RT: 3.14  
Area: 102461401  
Amount: 2.549765  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LVW8, 18-May-2023 13:29:56 -06:00:00 (UTC)

Audit Action: Manually Integrated/Assigned Compound ID Audit Reason: Baseline Smoothing

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-132  
 Lims ID: STD L6  
 Client ID:  
 Sample Type: IC Calib Level: 6  
 Inject. Date: 18-May-2023 13:10:00 ALS Bottle#: 0 Worklist Smp#: 6  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121594-006  
 Misc. Info.: 280-0121594-006  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Sublist: chrom-Anions\_IC10\*sub5  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
 Limit Group: Wet - Anions 28D  
 Last Update: 19-May-2023 11:53:49 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-132  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1603

First Level Reviewer: LVW8

Date: 18-May-2023 13:42:47

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response   | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|------------|------------------|--------------------|-------|
| 1 Fluoride            | 3.145        | 3.145            | 0.000            | 200510278  | 5.00             | 4.94               | Ma    |
| 2 Chloride            | 4.333        | 4.333            | 0.000            | 1864000362 | 100.0            | 99.5               |       |
| 3 Nitrite as N        | 4.997        | 4.997            | 0.000            | 212086804  | NC               | NC                 |       |
| 4 Bromide             | 6.150        | 6.150            | 0.000            | 38244681   | 5.00             | 5.06               |       |
| 5 Nitrate as N        | 6.890        | 6.890            | 0.000            | 224161262  | NC               | NC                 |       |
| 6 Orthophosphate as P | 9.035        | 9.035            | 0.000            | 115018456  | NC               | NC                 |       |
| 7 Sulfate             | 10.220       | 10.220           | 0.000            | 1371227441 | 100.0            | 99.0               |       |

**QC Flag Legend**

Processing Flags

NC - Not Calibrated

Review Flags

M - Manually Integrated

a - User Assigned ID

**Reagents:**

IC CAL c/so4\_00480

Amount Added: 8.00

Units: mL


IC Cal low\_00709

Amount Added: 2.00

Units: mL

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-132551.d  
Injection Date: 18-May-2023 13:10:00 Instrument ID: WC\_IonChrom10 Operator ID: wetchemd  
Lims ID: STD L6 Worklist Smp#: 6  
Client ID:  
Injection Vol: 5.0 ul Dil. Factor: 1.0000 ALS Bottle#: 0  
Method: Anions\_IC10 Limit Group: Wet - Anions 28D



## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-132551.d  
Injection Date: 18-May-2023 13:10:00 Instrument ID: WC\_IonChrom10  
Lims ID: STD L6  
Client ID:  
Operator ID: wetchemd ALS Bottle#: 0 Worklist Smp#: 6  
Injection Vol: 5.0 ul Dil. Factor: 1.0000  
Method: Anions\_IC10 Limit Group: Wet - Anions 28D  
Column: Detector Info 2\_091554\_1

**1 Fluoride, CAS: 16984-48-8**

Signal: 1


Not Detected  
Expected RT: 3.15

## Processing Integration Results



RT: 3.15  
Area: 200510278  
Amount: 4.938239  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LVW8, 18-May-2023 13:30:11 -06:00:00 (UTC)

Audit Action: Manually Integrated/Assigned Compound ID Audit Reason: Baseline Smoothing

## IC Instrument Information

**WL:** 121594 **Inst ID:** 1010 **Analysis Date:** 5/18/23 **Analyst:** MEC

| Rush                     | Job No.       | Samples   | Anions                                                     | QC Req | HT Exp |
|--------------------------|---------------|-----------|------------------------------------------------------------|--------|--------|
| <input type="checkbox"/> | <u>176637</u> | <u>1</u>  | F Cl NO2 Br <u>NO3</u> PO4 SO4                             | MS/D   | _____  |
| <input type="checkbox"/> | <u>176589</u> | <u>1</u>  | F Cl <u>NO2</u> Br NO3 PO4 SO4                             | MS/D   | _____  |
| <input type="checkbox"/> | <u>176670</u> | <u>7</u>  | F <u>Cl</u> NO2 Br <u>NO3</u> PO4 <u>SO4</u>               | MS/D   | _____  |
| <input type="checkbox"/> | <u>176676</u> | <u>15</u> | <u>F</u> <u>Cl</u> <u>NO2</u> Br <u>NO3</u> PO4 <u>SO4</u> | MS/D   | _____  |
| <input type="checkbox"/> | <u>176674</u> | <u>4</u>  | F Cl <u>NO2</u> Br <u>NO3</u> PO4 <u>SO4</u>               | MS/D   | _____  |
| <input type="checkbox"/> | <u>176678</u> | <u>5</u>  | F <u>Cl</u> <u>NO2</u> Br <u>NO3</u> PO4 <u>SO4</u>        | MS/D   | _____  |
| <input type="checkbox"/> | <u>176669</u> | <u>4</u>  | F <u>Cl</u> NO2 Br <u>NO3</u> PO4 <u>SO4</u>               | MS/D   | _____  |
| <input type="checkbox"/> | <u>176683</u> | <u>3</u>  | F <u>Cl</u> NO2 Br <u>NO3</u> PO4 <u>SO4</u>               | MS/D   | _____  |
| <input type="checkbox"/> | <u>176668</u> | <u>1</u>  | F Cl NO2 Br <u>NO3</u> PO4 SO4                             | MS/D   | _____  |
| <input type="checkbox"/> | _____         | _____     | F Cl NO2 Br NO3 PO4 SO4                                    | MS/D   | _____  |
| <input type="checkbox"/> | _____         | _____     | F Cl NO2 Br NO3 PO4 SO4                                    | MS/D   | _____  |
| <input type="checkbox"/> | _____         | _____     | F Cl NO2 Br NO3 PO4 SO4                                    | MS/D   | _____  |
| <input type="checkbox"/> | _____         | _____     | F Cl NO2 Br NO3 PO4 SO4                                    | MS/D   | _____  |
| <input type="checkbox"/> | _____         | _____     | F Cl NO2 Br NO3 PO4 SO4                                    | MS/D   | _____  |

### Dilutions

| Job No. | Samples | Anions                  | Dilution | Reason |
|---------|---------|-------------------------|----------|--------|
| _____   | _____   | F Cl NO2 Br NO3 PO4 SO4 | _____    | _____  |
| _____   | _____   | F Cl NO2 Br NO3 PO4 SO4 | _____    | _____  |
| _____   | _____   | F Cl NO2 Br NO3 PO4 SO4 | _____    | _____  |
| _____   | _____   | F Cl NO2 Br NO3 PO4 SO4 | _____    | _____  |
| _____   | _____   | F Cl NO2 Br NO3 PO4 SO4 | _____    | _____  |
| _____   | _____   | F Cl NO2 Br NO3 PO4 SO4 | _____    | _____  |
| _____   | _____   | F Cl NO2 Br NO3 PO4 SO4 | _____    | _____  |
| _____   | _____   | F Cl NO2 Br NO3 PO4 SO4 | _____    | _____  |
| _____   | _____   | F Cl NO2 Br NO3 PO4 SO4 | _____    | _____  |

Eurofins Environment Testing America  
Initial Calibration Report

Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230519-121638.b\Anions\_IC10.m

Instrument: WC\_IonChrom10

Lims Location: 280

Lock State: Initial Calib Locked

Cpnd Order: Retention Time

Integrator: Falcon

Last Modified: 19-May-2023 12:31:06

No. Compounds: 7

Sublist: chrom-Anions\_IC10\*sub5

Limit Group: Wet - Anions

## Detectors

Detector: 1, Info 2\_091554\_1

Data Type: ic Spec Type: none

Supports Extracted Chromatograms: False

Run Time: 0.000-14.900

No. Points: 7196

## Calibration File Names

|          |                                                                                                  |
|----------|--------------------------------------------------------------------------------------------------|
| Level: 1 | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230518-12110 |
|          | Inj Date: 18-May-2023 11:56:00 Worklist: 121594 Sample#: 1                                       |
| Level: 2 | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230518-12255 |
|          | Inj Date: 18-May-2023 12:11:00 Worklist: 121594 Sample#: 2                                       |
| Level: 3 | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230518-12405 |
|          | Inj Date: 18-May-2023 12:26:00 Worklist: 121594 Sample#: 3                                       |
| Level: 4 | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230518-12555 |
|          | Inj Date: 18-May-2023 12:40:00 Worklist: 121594 Sample#: 4                                       |
| Level: 5 | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230518-13105 |
|          | Inj Date: 18-May-2023 12:55:00 Worklist: 121594 Sample#: 5                                       |
| Level: 6 | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230518-13255 |
|          | Inj Date: 18-May-2023 13:10:00 Worklist: 121594 Sample#: 6                                       |

Start Cal Date: 18-May-2023 11:56:00

End Cal Date: 19-May-2023 12:01:00

## Individual Compound Calibration Parameters

Quant Method: ESTD

RF Calibration: Replace

Rule Name: Linear1

Curve: Linear

Weighting: Conc

Origin: None

Error: raw\_COD

Error Limit: 1.00

RF %Dif: 0.0

SPCC Limit: 0.0

CCC Limit: 0.0

Dependent Variable: Resp

Number of Compounds: 3

| RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | b | m1 | m2 | Error |
|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|---|----|----|-------|
|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|---|----|----|-------|

| RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | b | m1 | m2 | Error |
|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|---|----|----|-------|
|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|---|----|----|-------|

## 3 Nitrite as N

Signal: 1

|             |             |            |            |            |  |         |          |  |       |
|-------------|-------------|------------|------------|------------|--|---------|----------|--|-------|
| 35311960    | 42799726    | 43146865   | 42334489   | 42417361   |  | -967573 |          |  | 1.000 |
| 0.200000(2) | 0.500000(3) | 1.0000 (4) | 2.5000 (5) | 5.0000 (6) |  |         | 42866306 |  |       |
| 7062392     | 21399863    | 43146865   | 105836223  | 212086804  |  |         |          |  |       |
| 121594(2)   | 121594(3)   | 121594(4)  | 121594(5)  | 121594(6)  |  |         |          |  |       |

## 5 Nitrate as N

Signal: 1

|             |             |            |            |            |  |          |          |  |       |
|-------------|-------------|------------|------------|------------|--|----------|----------|--|-------|
| 37055495    | 45958086    | 42854846   | 44285580   | 44832252   |  | -1159043 |          |  | 0.999 |
| 0.200000(2) | 0.500000(3) | 1.0000 (4) | 2.5000 (5) | 5.0000 (6) |  |          | 44990806 |  |       |
| 7411099     | 22979043    | 42854846   | 110713949  | 224161262  |  |          |          |  |       |
| 121594(2)   | 121594(3)   | 121594(4)  | 121594(5)  | 121594(6)  |  |          |          |  |       |

## 6 Orthophosphate as P

Signal: 1

|             |             |            |            |            |  |         |          |  |        |
|-------------|-------------|------------|------------|------------|--|---------|----------|--|--------|
| 43303895    | 41843182    | 33022328   | 26624129   | 23003691   |  | 5964839 | R1       |  | *0.980 |
| 0.200000(2) | 0.500000(3) | 1.0000 (4) | 2.5000 (5) | 5.0000 (6) |  |         | 23299922 |  |        |
| 8660779     | 20921591    | 33022328   | 66560323   | 115018456  |  |         |          |  |        |
| 121594(2)   | 121594(3)   | 121594(4)  | 121594(5)  | 121594(6)  |  |         |          |  |        |

ICalib Error Legend

R1, Curve Fit Fail Error Limit Test



Eurofins Environment Testing America  
Initial Calibration Report

Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230519-121638.b\Anions\_IC10.m

Instrument: WC\_IonChrom10

Lims Location: 280

Lock State: Initial Calib Locked

Cpnd Order: Retention Time

Integrator: Falcon

Last Modified: 19-May-2023 12:31:06

No.Compounds:7

Sublist: chrom-Anions\_IC10\*sub5

Limit Group: Wet - Anions 28D

## Detectors

Detector: 1, Info 2\_091554\_1

Data Type: ic Spec Type: none

Supports Extracted Chromatograms: False

Run Time: 0.000-14.900

No. Points: 7196

## Calibration File Names

|          |                                                                                                  |
|----------|--------------------------------------------------------------------------------------------------|
| Level: 1 | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230518-12110 |
|          | Inj Date: 18-May-2023 11:56:00 Worklist: 121594 Sample#: 1                                       |
| Level: 2 | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230518-12255 |
|          | Inj Date: 18-May-2023 12:11:00 Worklist: 121594 Sample#: 2                                       |
| Level: 3 | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230518-12405 |
|          | Inj Date: 18-May-2023 12:26:00 Worklist: 121594 Sample#: 3                                       |
| Level: 4 | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230518-12555 |
|          | Inj Date: 18-May-2023 12:40:00 Worklist: 121594 Sample#: 4                                       |
| Level: 5 | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230518-13105 |
|          | Inj Date: 18-May-2023 12:55:00 Worklist: 121594 Sample#: 5                                       |
| Level: 6 | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230518-13255 |
|          | Inj Date: 18-May-2023 13:10:00 Worklist: 121594 Sample#: 6                                       |

Start Cal Date: 18-May-2023 11:56:00

End Cal Date: 19-May-2023 12:01:00

## Individual Compound Calibration Parameters

Quant Method: ESTD

RF Calibration: Replace

Rule Name: Linear1

Curve: Linear

Weighting: Conc

Origin: None

Error: raw\_COD

Error Limit: 1.00

RF %Dif: 0.0

SPCC Limit: 0.0

CCC Limit: 0.0

Dependent Variable: Resp

Number of Compounds: 4

| RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | b | m1 | m2 | Error |
|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|---|----|----|-------|
|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|---|----|----|-------|

| RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | RF/Amt(Lvl)<br>Response<br>WL(Smp) | b         | m1       | m2 | Error |
|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|-----------|----------|----|-------|
| 1 Fluoride                         |                                    |                                    |                                    |                                    |                                    | Signal: 1 |          |    |       |
| 29175985                           | 36359510                           | 39638142                           | 40984560                           | 40102056                           |                                    | -2208602  |          |    | 1.000 |
| 0.200000(2)                        | 0.500000(3)                        | 1.0000 (4)                         | 2.5000 (5)                         | 5.0000 (6)                         |                                    |           | 41050846 |    |       |
| M5835197                           | M18179755                          | M39638142                          | M102461401                         | M200510278                         |                                    |           |          |    |       |
| 121594(2)                          | 121594(3)                          | 121594(4)                          | 121594(5)                          | 121594(6)                          |                                    |           |          |    |       |
| 2 Chloride                         |                                    |                                    |                                    |                                    |                                    | Signal: 1 |          |    |       |
| 16780722                           | 16174841                           | 19339456                           | 19009666                           | 18725386                           | 18640004                           | -1171132  |          |    | 1.000 |
| 0.500000(1)                        | 1.2500 (2)                         | 2.5000 (3)                         | 30.0 (4)                           | 60.0 (5)                           | 100.0 (6)                          |           | 18747994 |    |       |
| 8390361                            | 20218551                           | 48348641                           | 570289992                          | 1123523133                         | 1864000362                         |           |          |    |       |
| 121594(1)                          | 121594(2)                          | 121594(3)                          | 121594(4)                          | 121594(5)                          | 121594(6)                          |           |          |    |       |
| 4 Bromide                          |                                    |                                    |                                    |                                    |                                    | Signal: 1 |          |    |       |
| 11277385                           | 9537194                            | 7900278                            | 7592206                            | 7648936                            |                                    | 806494    |          |    | 0.999 |
| 0.200000(2)                        | 0.500000(3)                        | 1.0000 (4)                         | 2.5000 (5)                         | 5.0000 (6)                         |                                    |           | 7404030  |    |       |
| 2255477                            | 4768597                            | 7900278                            | 18980516                           | 38244681                           |                                    |           |          |    |       |
| 121594(2)                          | 121594(3)                          | 121594(4)                          | 121594(5)                          | 121594(6)                          |                                    |           |          |    |       |
| 7 Sulfate                          |                                    |                                    |                                    |                                    |                                    | Signal: 1 |          |    |       |
| 22625290                           | 16020812                           | 17380170                           | 14404666                           | 13852661                           | 13712274                           | 4654296   |          |    | 0.999 |
| 0.500000(1)                        | 1.2500 (2)                         | 2.5000 (3)                         | 30.0 (4)                           | 60.0 (5)                           | 100.0 (6)                          |           | 13803812 |    |       |
| 11312645                           | 20026015                           | 43450426                           | 432139991                          | 831159688                          | 1371227441                         |           |          |    |       |
| 121594(1)                          | 121594(2)                          | 121594(3)                          | 121594(4)                          | 121594(5)                          | 121594(6)                          |           |          |    |       |

## Preliminary Report

Eurofins Denver  
ICV, ICal Verification Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13  
Lims ID: ICV  
Client ID:  
Sample Type: ICV  
Inject. Date: 18-May-2023 13:25:00 ALS Bottle#: 0 Worklist Smp#: 7  
Injection Vol: 5.0 ul Dil. Factor: 1.0000  
Sample Info: 280-0121594-007  
Misc. Info.: 280-0121594-007  
Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
Sublist: chrom-Anions\_IC7\*sub4

Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
Limit Group: Wet - Anions 28D  
Last Update: 19-May-2023 12:32:45 Calib Date: 18-May-2023 13:10:00  
Integrator: Falcon  
Quant Method: External Standard Quant By: Initial Calibration  
Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13

Column 1 : Det: Info 2\_091554\_1

Process Host: CTX1603

First Level Reviewer: LVW8

Date: 19-May-2023 12:38:30

Start Cal Date: 18-May-2023 11:56:00

End Cal Date: 18-May-2023 13:10:00

| Compound   | Amount Added | Amount Detected | %Drift | Max. %Drift | %Rec  | %Rec Limits |
|------------|--------------|-----------------|--------|-------------|-------|-------------|
| 1 Fluoride | 2.00         | 2.34            | * 16.8 | 10.0        | 116.8 |             |
| 2 Chloride | 80.0         | 83.0            | 3.8    | 10.0        | 103.8 |             |
| 4 Bromide  | 2.00         | 1.91            | -4.7   | 10.0        | 95.3  |             |
| 7 Sulfate  | 80.0         | 81.8            | 2.3    | 10.0        | 102.3 |             |

## Preliminary Report

Eurofins Denver  
ICV, ICal Verification Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13  
Lims ID: ICV  
Client ID:  
Sample Type: ICV  
Inject. Date: 18-May-2023 13:25:00 ALS Bottle#: 0 Worklist Smp#: 7  
Injection Vol: 5.0 ul Dil. Factor: 1.0000  
Sample Info: 280-0121594-007  
Misc. Info.: 280-0121594-007  
Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
Sublist: chrom-Anions\_IC7\*sub4  
Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
Limit Group: Wet - Anions  
Last Update: 19-May-2023 12:38:46 Calib Date: 18-May-2023 13:10:00  
Integrator: Falcon  
Quant Method: External Standard Quant By: Initial Calibration  
Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13  
Column 1 : Det: Info 2\_091554\_1  
Process Host: CTX1603  
First Level Reviewer: LVW8 Date: 19-May-2023 12:38:46  
Start Cal Date: 18-May-2023 11:56:00  
End Cal Date: 18-May-2023 13:10:00

| Compound            | Amount Added | Amount Detected | %Drift | Max. %Drift | %Rec  | %Rec Limits |
|---------------------|--------------|-----------------|--------|-------------|-------|-------------|
| 3 Nitrite as N      | 2.00         | 1.99            | -0.3   | 10.0        | 99.7  |             |
| 5 Nitrate as N      | 2.00         | 1.95            | -2.3   | 10.0        | 97.7  |             |
| 6 Orthophosphate as | 2.00         | 2.30            | * 14.8 | 10.0        | 114.8 |             |

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13  
 Lims ID: ICV  
 Client ID:  
 Sample Type: ICV  
 Inject. Date: 18-May-2023 13:25:00 ALS Bottle#: 0 Worklist Smp#: 7  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121594-007  
 Misc. Info.: 280-0121594-007  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Sublist: chrom-Anions\_IC7\*sub4  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
 Limit Group: Wet - Anions 28D  
 Last Update: 19-May-2023 12:38:46 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1603

First Level Reviewer: LVW8

Date: 19-May-2023 12:38:46

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response   | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|------------|------------------|--------------------|-------|
| 1 Fluoride            | 3.145        | 3.145            | 0.000            | 93724806   | 2.00             | 2.34               | Ma    |
| 2 Chloride            | 4.328        | 4.328            | 0.000            | 1555403505 | 80.0             | 83.0               |       |
| 3 Nitrite as N        | 5.000        | 5.000            | 0.000            | 84485309   | NC               | NC                 |       |
| 4 Bromide             | 6.162        | 6.162            | 0.000            | 14920856   | 2.00             | 1.91               |       |
| 5 Nitrate as N        | 6.917        | 6.917            | 0.000            | 86743393   | NC               | NC                 |       |
| 6 Orthophosphate as P | 9.058        | 9.058            | 0.000            | 59448017   | NC               | NC                 |       |
| 7 Sulfate             | 10.230       | 10.230           | 0.000            | 1134012406 | 80.0             | 81.8               |       |

**QC Flag Legend**

Processing Flags

NC - Not Calibrated

Review Flags

M - Manually Integrated


a - User Assigned ID

**Reagents:**

|                  |                    |           |
|------------------|--------------------|-----------|
| IC SO4 ICV_00024 | Amount Added: 0.80 | Units: mL |
| CI ICV Std_00006 | Amount Added: 0.80 | Units: mL |
| IC ICV 5_00405   | Amount Added: 0.80 | Units: mL |

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-134050.d  
Injection Date: 18-May-2023 13:25:00 Instrument ID: WC\_IonChrom10 Operator ID: wetchemd  
Lims ID: ICV Worklist Smp#: 7  
Client ID:  
Injection Vol: 5.0 ul Dil. Factor: 1.0000 ALS Bottle#: 0  
Method: Anions\_IC10 Limit Group: Wet - Anions 28D



Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13  
 Lims ID: ICB  
 Client ID:  
 Sample Type: ICB  
 Inject. Date: 18-May-2023 13:40:00 ALS Bottle#: 0 Worklist Smp#: 8  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121594-008  
 Misc. Info.: 280-0121594-008  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
 Limit Group: Wet - Anions 28D  
 Last Update: 19-May-2023 12:31:30 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1603

First Level Reviewer: LVW8

Date: 19-May-2023 11:57:30

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 1 Fluoride            |              | 3.145            |                  |          |                  | ND                 |       |
| 2 Chloride            |              | 4.328            |                  |          |                  | ND                 |       |
| 3 Nitrite as N        |              | 5.000            |                  |          |                  | ND                 |       |
| 4 Bromide             |              | 6.162            |                  |          |                  | ND                 | U     |
| 5 Nitrate as N        |              | 6.917            |                  |          |                  | ND                 |       |
| 6 Orthophosphate as P | 9.123        | 9.058            | 0.065            | 12987322 |                  | NC                 |       |
| 7 Sulfate             | 10.285       | 10.230           | 0.055            | 1071622  |                  | -0.2595            | M     |

**QC Flag Legend**

Processing Flags

NC - Not Calibrated

Review Flags

M - Manually Integrated

U - Marked Undetected

Report Date: 19-May-2023 12:31:31

Chrom Revision: 2.3 16-May-2023 16:37:30

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-135549.d

Injection Date: 18-May-2023 13:40:00

Instrument ID: WC\_IonChrom10

Operator ID: wetchemd

Lims ID: ICB

Worklist Smp#: 8

Client ID:


Injection Vol: 5.0 ul

Dil. Factor: 1.0000

ALS Bottle#: 0

Method: Anions\_IC10

Limit Group: Wet - Anions 28D





## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-135549.d  
Injection Date: 18-May-2023 13:40:00 Instrument ID: WC\_IonChrom10  
Lims ID: ICB  
Client ID:  
Operator ID: wetchemd ALS Bottle#: 0 Worklist Smp#: 8  
Injection Vol: 5.0 ul Dil. Factor: 1.0000  
Method: Anions\_IC10 Limit Group: Wet - Anions 28D  
Column: Detector Info 2\_091554\_1

**7 Sulfate, CAS: 14808-79-8**

Signal: 1


RT: 10.30  
Area: 7155425  
Amount: 0.181191  
Amount Units: ug/ml

## Processing Integration Results



RT: 10.29  
Area: 1071622  
Amount: -0.259542  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LVW8, 19-May-2023 11:57:25 -06:00:00 (UTC)

Audit Action: Manually Integrated

Audit Reason: Baseline Smoothing

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-14  
 Lims ID: MRL  
 Client ID:  
 Sample Type: MRL  
 Inject. Date: 18-May-2023 14:10:00 ALS Bottle#: 0 Worklist Smp#: 10  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121594-010  
 Misc. Info.: 280-0121594-010  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
 Limit Group: Wet - Anions 28D  
 Last Update: 19-May-2023 12:31:32 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1603

First Level Reviewer: LVW8

Date: 18-May-2023 14:58:41

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 1 Fluoride            |              | 3.145            |                  |          | ND               | ND                 |       |
| 2 Chloride            | 4.302        | 4.313            | -0.011           | 43926233 | 2.50             | 2.41               |       |
| 3 Nitrite as N        | 5.003        | 4.998            | 0.005            | 8516946  | NC               | NC                 |       |
| 4 Bromide             | 6.178        | 6.163            | 0.015            | 2736457  | 0.2500           | 0.2607             |       |
| 5 Nitrate as N        | 6.958        | 6.913            | 0.045            | 9978227  | NC               | NC                 |       |
| 6 Orthophosphate as P | 9.105        | 9.060            | 0.045            | 20685657 | NC               | NC                 |       |
| 7 Sulfate             | 10.265       | 10.247           | 0.018            | 45774744 | 2.50             | 2.98               |       |

**QC Flag Legend**

Processing Flags

NC - Not Calibrated

ND - Not Detected or Marked ND

**Reagents:**

IC CAL cl/so4\_00480

Amount Added: 0.10

Units: mL

IC Cal low\_00709

Amount Added: 0.05

Units: mL

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-142550.d

Injection Date: 18-May-2023 14:10:00

Instrument ID: WC\_IonChrom10

Operator ID: wetchemd

Lims ID: MRL

Worklist Smp#: 10

Client ID:


Injection Vol: 5.0 ul

Dil. Factor: 1.0000

ALS Bottle#: 0

Method: Anions\_IC10

Limit Group: Wet - Anions 28D



Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-14  
 Lims ID: LCS  
 Client ID:  
 Sample Type: LCS  
 Inject. Date: 18-May-2023 14:25:00 ALS Bottle#: 0 Worklist Smp#: 11  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121594-011  
 Misc. Info.: 280-0121594-011  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
 Limit Group: Wet - Anions 28D  
 Last Update: 19-May-2023 12:31:32 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1603

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response  | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|-----------|------------------|--------------------|-------|
| 1 Fluoride            |              | 3.145            |                  |           | ND               | ND                 |       |
| 2 Chloride            | 4.313        | 4.313            | 0.000            | 924355606 | 50.0             | 49.4               |       |
| 3 Nitrite as N        | 4.998        | 4.998            | 0.000            | 102973098 | NC               | NC                 |       |
| 4 Bromide             | 6.163        | 6.163            | 0.000            | 18068153  | 2.50             | 2.33               |       |
| 5 Nitrate as N        | 6.912        | 6.913            | -0.001           | 109333247 | NC               | NC                 |       |
| 6 Orthophosphate as P | 9.062        | 9.060            | 0.002            | 71396675  | NC               | NC                 |       |
| 7 Sulfate             | 10.248       | 10.247           | 0.001            | 701445617 | 50.0             | 50.5               |       |

### QC Flag Legend

Processing Flags

NC - Not Calibrated

ND - Not Detected or Marked ND

### Reagents:

IC LCS\_01954

Amount Added: 5.00

Units: mL

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-144050.d

Injection Date: 18-May-2023 14:25:00

Instrument ID: WC\_IonChrom10

Operator ID: wetchemd

Lims ID: LCS

Worklist Smp#: 11

Client ID:


Injection Vol: 5.0 ul

Dil. Factor: 1.0000

ALS Bottle#: 0

Method: Anions\_IC10

Limit Group: Wet - Anions 28D



Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-14  
 Lims ID: LCSD  
 Client ID:  
 Sample Type: LCSD  
 Inject. Date: 18-May-2023 14:40:00 ALS Bottle#: 0 Worklist Smp#: 12  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121594-012  
 Misc. Info.: 280-0121594-012  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
 Limit Group: Wet - Anions 28D  
 Last Update: 19-May-2023 12:31:32 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1603

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response  | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|-----------|------------------|--------------------|-------|
| 1 Fluoride            |              | 3.145            |                  |           | ND               | ND                 |       |
| 2 Chloride            | 4.313        | 4.313            | 0.000            | 921265063 | 50.0             | 49.2               |       |
| 3 Nitrite as N        | 4.998        | 4.998            | 0.000            | 102379982 | NC               | NC                 |       |
| 4 Bromide             | 6.163        | 6.163            | 0.000            | 17988727  | 2.50             | 2.32               |       |
| 5 Nitrate as N        | 6.912        | 6.913            | -0.001           | 108961704 | NC               | NC                 |       |
| 6 Orthophosphate as P | 9.062        | 9.060            | 0.002            | 71915979  | NC               | NC                 |       |
| 7 Sulfate             | 10.250       | 10.247           | 0.003            | 695470892 | 50.0             | 50.0               |       |

### QC Flag Legend

Processing Flags

NC - Not Calibrated

ND - Not Detected or Marked ND

### Reagents:


IC LCS\_01954

Amount Added: 5.00

Units: mL

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-145552.d  
Injection Date: 18-May-2023 14:40:00 Instrument ID: WC\_IonChrom10 Operator ID: wetchemd  
Lims ID: LCSD Worklist Smp#: 12  
Client ID:  
Injection Vol: 5.0 ul Dil. Factor: 1.0000 ALS Bottle#: 0  
Method: Anions\_IC10 Limit Group: Wet - Anions 28D



Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-15  
 Lims ID: MB  
 Client ID:  
 Sample Type: MB  
 Inject. Date: 18-May-2023 14:55:00 ALS Bottle#: 0 Worklist Smp#: 13  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121594-013  
 Misc. Info.: 280-0121594-013  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
 Limit Group: Wet - Anions 28D  
 Last Update: 19-May-2023 12:31:35 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1603

First Level Reviewer: LVW8

Date: 19-May-2023 11:57:53

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 1 Fluoride            |              | 3.145            |                  |          |                  | ND                 |       |
| 2 Chloride            |              | 4.313            |                  |          |                  | ND                 |       |
| 3 Nitrite as N        |              | 4.998            |                  |          |                  | ND                 |       |
| 4 Bromide             |              | 6.163            |                  |          |                  | ND                 | U     |
| 5 Nitrate as N        |              | 6.913            |                  |          |                  | ND                 |       |
| 6 Orthophosphate as P | 9.127        | 9.127            | 0.000            | 13698318 |                  | NC                 |       |
| 7 Sulfate             | 10.300       | 10.300           | 0.000            | 1160087  |                  | -0.2531            | M     |

**QC Flag Legend**

Processing Flags

NC - Not Calibrated

Review Flags

M - Manually Integrated

U - Marked Undetected



Report Date: 19-May-2023 12:31:35

Chrom Revision: 2.3 16-May-2023 16:37:30

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-151052.d

Injection Date: 18-May-2023 14:55:00

Instrument ID: WC\_IonChrom10

Operator ID: wetchemd

Lims ID: MB

Worklist Smp#: 13

Client ID:


Injection Vol: 5.0 ul

Dil. Factor: 1.0000

ALS Bottle#: 0

Method: Anions\_IC10

Limit Group: Wet - Anions 28D



## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-151052.d  
Injection Date: 18-May-2023 14:55:00 Instrument ID: WC\_IonChrom10  
Lims ID: MB  
Client ID:  
Operator ID: wetchemd ALS Bottle#: 0 Worklist Smp#: 13  
Injection Vol: 5.0 ul Dil. Factor: 1.0000  
Method: Anions\_IC10 Limit Group: Wet - Anions 28D  
Column: Detector Info 2\_091554\_1

**7 Sulfate, CAS: 14808-79-8**

Signal: 1


RT: 10.30  
Area: 7420754  
Amount: 0  
Amount Units: ug/ml

## Processing Integration Results



RT: 10.30  
Area: 1160087  
Amount: -0.253134  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LVW8, 19-May-2023 11:57:49 -06:00:00 (UTC)

Audit Action: Manually Integrated

Audit Reason: Baseline Smoothing

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-179  
 Lims ID: ccv  
 Client ID:  
 Sample Type: CCV  
 Inject. Date: 18-May-2023 17:41:00 ALS Bottle#: 0 Worklist Smp#: 24  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121594-024  
 Misc. Info.: 280-0121594-024  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Sublist: chrom-Anions\_IC10\*sub5  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
 Limit Group: Wet - Anions 28D  
 Last Update: 19-May-2023 12:31:42 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-132  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1603

First Level Reviewer: XAY4

Date: 19-May-2023 09:47:24

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response  | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|-----------|------------------|--------------------|-------|
| 1 Fluoride            |              | 3.145            |                  |           | ND               | ND                 |       |
| 2 Chloride            | 4.312        | 4.313            | -0.001           | 924998506 | 50.0             | 49.4               |       |
| 3 Nitrite as N        | 5.000        | 4.998            | 0.002            | 107797836 | NC               | NC                 |       |
| 4 Bromide             | 6.140        | 6.163            | -0.023           | 17602305  | 2.50             | 2.27               |       |
| 5 Nitrate as N        | 6.903        | 6.913            | -0.010           | 108124722 | NC               | NC                 |       |
| 6 Orthophosphate as P | 9.045        | 9.127            | -0.082           | 53503464  | NC               | NC                 |       |
| 7 Sulfate             | 10.223       | 10.300           | -0.077           | 680174954 | 50.0             | 48.9               |       |

**QC Flag Legend**

Processing Flags

NC - Not Calibrated

ND - Not Detected or Marked ND

**Reagents:**


IC LCS\_01954

Amount Added: 10.00

Units: mL

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-175619.d  
Injection Date: 18-May-2023 17:41:00 Instrument ID: WC\_IonChrom10 Operator ID: wetchemd  
Lims ID: ccv Worklist Smp#: 24  
Client ID:  
Injection Vol: 5.0 ul Dil. Factor: 1.0000 ALS Bottle#: 0  
Method: Anions\_IC10 Limit Group: Wet - Anions 28D



Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-18  
 Lims ID: ccb  
 Client ID:  
 Sample Type: CCB  
 Inject. Date: 18-May-2023 17:56:00 ALS Bottle#: 0 Worklist Smp#: 25  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121594-025  
 Misc. Info.: 280-0121594-025  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
 Limit Group: Wet - Anions 28D  
 Last Update: 19-May-2023 12:31:42 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1603

First Level Reviewer: LVW8

Date: 19-May-2023 12:00:26

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 1 Fluoride            |              | 3.145            |                  |          |                  | ND                 |       |
| 2 Chloride            |              | 4.313            |                  |          |                  | ND                 |       |
| 3 Nitrite as N        |              | 4.998            |                  |          |                  | ND                 | U     |
| 4 Bromide             |              | 6.163            |                  |          |                  | ND                 | U     |
| 5 Nitrate as N        |              | 6.913            |                  |          |                  | ND                 |       |
| 6 Orthophosphate as P | 9.133        | 9.127            | 0.006            | 5344408  |                  | NC                 |       |
| 7 Sulfate             | 10.305       | 10.300           | 0.005            | 810551   |                  | -0.2785            | M     |

**QC Flag Legend**

## Processing Flags

NC - Not Calibrated

## Review Flags

M - Manually Integrated

U - Marked Undetected

Report Date: 19-May-2023 12:31:43

Chrom Revision: 2.3 16-May-2023 16:37:30

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-181124.d

Injection Date: 18-May-2023 17:56:00

Instrument ID: WC\_IonChrom10

Operator ID: wetchemd

Lims ID: ccb

Worklist Smp#: 25

Client ID:


Injection Vol: 5.0 ul

Dil. Factor: 1.0000

ALS Bottle#: 0

Method: Anions\_IC10

Limit Group: Wet - Anions 28D



## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-181124.d  
Injection Date: 18-May-2023 17:56:00 Instrument ID: WC\_IonChrom10  
Lims ID: ccb  
Client ID:  
Operator ID: wetchemd ALS Bottle#: 0 Worklist Smp#: 25  
Injection Vol: 5.0 ul Dil. Factor: 1.0000  
Method: Anions\_IC10 Limit Group: Wet - Anions 28D  
Column: Detector Info 2\_091554\_1

**7 Sulfate, CAS: 14808-79-8**

Signal: 1


RT: 10.31  
Area: 4417823  
Amount: -0.017131  
Amount Units: ug/ml

## Processing Integration Results



RT: 10.31  
Area: 810551  
Amount: -0.278455  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LVW8, 19-May-2023 12:00:19 -06:00:00 (UTC)

Audit Action: Manually Integrated

Audit Reason: Baseline Smoothing

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-19  
 Lims ID: 280-176674-A-7  
 Client ID: FBQmw-173-230401-GW-R  
 Sample Type: Client  
 Inject. Date: 18-May-2023 19:26:00 ALS Bottle#: 0 Worklist Smp#: 31  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121594-031  
 Misc. Info.: 280-0121594-031  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
 Limit Group: Wet - Anions 28D  
 Last Update: 19-May-2023 12:31:42 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1603

First Level Reviewer: LVW8

Date: 19-May-2023 12:02:08

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response  | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|-----------|--------------------|-------|
| 1 Fluoride            |              | 3.145            |                  |           | ND                 |       |
| 2 Chloride            | 4.303        | 4.313            | -0.010           | 19574865  | 1.11               | M     |
| 3 Nitrite as N        |              | 4.998            |                  |           | ND                 |       |
| 4 Bromide             | 6.212        | 6.163            | 0.049            | 1197331   | 0.0528             |       |
| 5 Nitrate as N        | 6.975        | 6.913            | 0.062            | 7707821   | NC                 | M     |
| 6 Orthophosphate as P | 9.138        | 9.127            | 0.011            | 8994927   | NC                 |       |
| 7 Sulfate             | 10.243       | 10.300           | -0.057           | 580169865 | 41.7               | M     |

**QC Flag Legend**

Processing Flags

NC - Not Calibrated


Review Flags

M - Manually Integrated



Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-194147.d  
Injection Date: 18-May-2023 19:26:00 Instrument ID: WC\_IonChrom10 Operator ID: wetchemd  
Lims ID: 280-176674-A-7 Lab Sample ID: 280-176674-7 Worklist Smp#: 31  
Client ID: FBQmw-173-230401-GW-R  
Injection Vol: 5.0 ul Dil. Factor: 1.0000 ALS Bottle#: 0  
Method: Anions\_IC10 Limit Group: Wet - Anions 28D



## Eurofins Denver


|                 |                                                                                                     |                |                  |                |    |
|-----------------|-----------------------------------------------------------------------------------------------------|----------------|------------------|----------------|----|
| Data File:      | \\chromfs\Denver\ChromData\WC_IonChrom10\20230518-121594.b\Info 2_DENPC179_Anions_20230518-194147.d |                |                  |                |    |
| Injection Date: | 18-May-2023 19:26:00                                                                                | Instrument ID: | WC_IonChrom10    |                |    |
| Lims ID:        | 280-176674-A-7                                                                                      | Lab Sample ID: | 280-176674-7     |                |    |
| Client ID:      | FBQmw-173-230401-GW-R                                                                               |                |                  |                |    |
| Operator ID:    | wetchemd                                                                                            | ALS Bottle#:   | 0                | Worklist Smp#: | 31 |
| Injection Vol:  | 5.0 ul                                                                                              | Dil. Factor:   | 1.0000           |                |    |
| Method:         | Anions_IC10                                                                                         | Limit Group:   | Wet - Anions 28D |                |    |
| Column:         |                                                                                                     | Detector       | Info 2_091554_1  |                |    |

**7 Sulfate, CAS: 14808-79-8**

Signal: 1


RT: 10.24  
Area: 593904893  
Amount: 0  
Amount Units: ug/ml

## Processing Integration Results



RT: 10.24  
Area: 580169865  
Amount: 41.692511  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LVW8, 19-May-2023 12:02:04 -06:00:00 (UTC)

Audit Action: Manually Integrated

Audit Reason: Baseline Smoothing

Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-19  
 Lims ID: 280-176674-A-8  
 Client ID: FBQmw-173-230402-GW-R  
 Sample Type: Client  
 Inject. Date: 18-May-2023 19:41:00 ALS Bottle#: 0 Worklist Smp#: 32  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121594-032  
 Misc. Info.: 280-0121594-032  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
 Limit Group: Wet - Anions 28D  
 Last Update: 19-May-2023 12:31:42 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1603

First Level Reviewer: LVW8

Date: 19-May-2023 12:02:17

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response  | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|-----------|--------------------|-------|
| 1 Fluoride            |              | 3.145            |                  |           | ND                 |       |
| 2 Chloride            | 4.302        | 4.313            | -0.011           | 20099453  | 1.13               | M     |
| 3 Nitrite as N        |              | 4.998            |                  |           | ND                 |       |
| 4 Bromide             | 6.210        | 6.163            | 0.047            | 1338653   | 0.0719             |       |
| 5 Nitrate as N        | 6.977        | 6.913            | 0.064            | 8138449   | NC                 | M     |
| 6 Orthophosphate as P | 9.148        | 9.127            | 0.021            | 5690064   | NC                 |       |
| 7 Sulfate             | 10.245       | 10.300           | -0.055           | 627177033 | 45.1               |       |

**QC Flag Legend**

Processing Flags


NC - Not Calibrated

Review Flags

M - Manually Integrated

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-195651.d  
Injection Date: 18-May-2023 19:41:00 Instrument ID: WC\_IonChrom10 Operator ID: wetchemd  
Lims ID: 280-176674-A-8 Lab Sample ID: 280-176674-8 Worklist Smp#: 32  
Client ID: FBQmw-173-230402-GW-R  
Injection Vol: 5.0 ul Dil. Factor: 1.0000 ALS Bottle#: 0  
Method: Anions\_IC10 Limit Group: Wet - Anions 28D



Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-20230518-131594.b  
 Lims ID: ccv  
 Client ID:  
 Sample Type: CCV  
 Inject. Date: 18-May-2023 20:42:00 ALS Bottle#: 0 Worklist Smp#: 36  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121594-036  
 Misc. Info.: 280-0121594-036  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Sublist: chrom-Anions\_IC10\*sub5  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
 Limit Group: Wet - Anions 28D  
 Last Update: 19-May-2023 12:31:51 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-131594.b

Column 1 : Det: Info 2\_091554\_1

Process Host: CTX1603

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response  | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|-----------|------------------|--------------------|-------|
| 1 Fluoride            |              | 3.145            |                  |           | ND               | ND                 |       |
| 2 Chloride            | 4.315        | 4.313            | 0.002            | 926185158 | 50.0             | 49.5               |       |
| 3 Nitrite as N        | 5.003        | 4.998            | 0.005            | 107572491 | NC               | NC                 |       |
| 4 Bromide             | 6.145        | 6.163            | -0.018           | 17667283  | 2.50             | 2.28               |       |
| 5 Nitrate as N        | 6.910        | 6.913            | -0.003           | 108155910 | NC               | NC                 |       |
| 6 Orthophosphate as P | 9.048        | 9.127            | -0.079           | 53833009  | NC               | NC                 |       |
| 7 Sulfate             | 10.232       | 10.300           | -0.068           | 686565829 | 50.0             | 49.4               |       |

### QC Flag Legend

Processing Flags

NC - Not Calibrated

ND - Not Detected or Marked ND

### Reagents:


IC LCS\_01954

Amount Added: 10.00

Units: mL

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-205701.d  
Injection Date: 18-May-2023 20:42:00 Instrument ID: WC\_IonChrom10 Operator ID: wetchemd  
Lims ID: ccv Worklist Smp#: 36  
Client ID:  
Injection Vol: 5.0 ul Dil. Factor: 1.0000 ALS Bottle#: 0  
Method: Anions\_IC10 Limit Group: Wet - Anions 28D



Eurofins Denver  
Target Compound Quantitation Report

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-21.  
 Lims ID: ccb  
 Client ID:  
 Sample Type: CCB  
 Inject. Date: 18-May-2023 20:57:00 ALS Bottle#: 0 Worklist Smp#: 37  
 Injection Vol: 5.0 ul Dil. Factor: 1.0000  
 Sample Info: 280-0121594-037  
 Misc. Info.: 280-0121594-037  
 Operator ID: wetchemd Instrument ID: WC\_IonChrom10  
 Method: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Anions\_IC10.m  
 Limit Group: Wet - Anions 28D  
 Last Update: 19-May-2023 12:31:51 Calib Date: 18-May-2023 13:10:00  
 Integrator: Falcon  
 Quant Method: External Standard Quant By: Initial Calibration  
 Last ICal File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-13.  
 Column 1 : Det: Info 2\_091554\_1  
 Process Host: CTX1603

First Level Reviewer: LVW8

Date: 19-May-2023 12:03:41

| Compound              | RT<br>(min.) | Exp RT<br>(min.) | Dlt RT<br>(min.) | Response | Cal Amt<br>ug/ml | OnCol Amt<br>ug/ml | Flags |
|-----------------------|--------------|------------------|------------------|----------|------------------|--------------------|-------|
| 1 Fluoride            |              | 3.145            |                  |          |                  | ND                 |       |
| 2 Chloride            | 4.305        | 4.313            | -0.008           | 7127336  |                  | 0.4426             | M     |
| 3 Nitrite as N        | 5.003        | 4.998            | 0.005            | 116037   |                  | NC                 | M     |
| 4 Bromide             | 6.167        | 6.163            | 0.004            | 189030   |                  | -0.0834            | M     |
| 5 Nitrate as N        | 6.973        | 6.913            | 0.060            | 76369    |                  | NC                 | M     |
| 6 Orthophosphate as P | 9.052        | 9.127            | -0.075           | 6191477  |                  | NC                 |       |
| 7 Sulfate             | 10.250       | 10.300           | -0.050           | 11085723 |                  | 0.4659             | M     |

**QC Flag Legend**

Processing Flags

NC - Not Calibrated

Review Flags

M - Manually Integrated

Report Date: 19-May-2023 12:31:52

Chrom Revision: 2.3 16-May-2023 16:37:30

Eurofins Denver

Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-211202.d

Injection Date: 18-May-2023 20:57:00

Instrument ID: WC\_IonChrom10

Operator ID: wetchemd

Lims ID: ccb

Worklist Smp#: 37

Client ID:


Injection Vol: 5.0 ul

Dil. Factor: 1.0000

ALS Bottle#: 0

Method: Anions\_IC10

Limit Group: Wet - Anions 28D





## Eurofins Denver


Data File: \\chromfs\Denver\ChromData\WC\_IonChrom10\20230518-121594.b\Info 2\_DENPC179\_Anions\_20230518-211202.d  
Injection Date: 18-May-2023 20:57:00 Instrument ID: WC\_IonChrom10  
Lims ID: ccb  
Client ID:  
Operator ID: wetchemd ALS Bottle#: 0 Worklist Smp#: 37  
Injection Vol: 5.0 ul Dil. Factor: 1.0000  
Method: Anions\_IC10 Limit Group: Wet - Anions 28D  
Column: Detector Info 2\_091554\_1

**7 Sulfate, CAS: 14808-79-8**

Signal: 1


RT: 10.25  
Area: 14926462  
Amount: 0  
Amount Units: ug/ml

## Processing Integration Results



RT: 10.25  
Area: 11085723  
Amount: 0.465917  
Amount Units: ug/ml

## Manual Integration Results



Reviewer: LVW8, 19-May-2023 12:03:21 -06:00:00 (UTC)

Audit Action: Manually Integrated

Audit Reason: Baseline Smoothing

|    | Analysis | Sample Name        | ObjectID                                  |
|----|----------|--------------------|-------------------------------------------|
| 1  | NPOC     | ICV                | 0L-10900100866-10101000-134B17724A86-0000 |
| 2  | NPOC     | ICB                | 0L-10900100866-10101000-134B17724A86-0001 |
| 3  | NPOC     | LCS                | 0L-10900100866-10101000-134B17724A86-0002 |
| 4  | NPOC     | MB                 | 0L-10900100866-10101000-134B17724A86-0003 |
| 5  | NPOC     | TIC                | 0L-10900100866-10101000-134B17724A86-0004 |
| 6  | NPOC     | 280-176674-d-7     | 0L-10900100866-10101000-134B17724A86-0005 |
| 7  | NPOC     | MS 280-176674-d-7  | 0L-10900100866-10101000-134B17724A87-0000 |
| 8  | NPOC     | MSD 280-176674-d-7 | 0L-10900100866-10101000-134B17724A87-0001 |
| 9  | NPOC     | 580-127209-p-1     | 0L-10900100866-10101000-134B17724A87-0002 |
| 10 | NPOC     | 580-127209-p-3     | 0L-10900100866-10101000-134B17724A87-0003 |
| 11 | NPOC     | 580-127209-p-5     | 0L-10900100866-10101000-134B17724A87-0004 |
| 12 | NPOC     | 280-176349-d-8     | 0L-10900100866-10101000-134B17724A87-0005 |
| 13 | NPOC     | 280-176349-c-11    | 0L-10900100866-10101000-134B17724A87-0006 |
| 14 | NPOC     | 280-176123-a-35    | 0L-10900100866-10101000-134B17724A87-0007 |
| 15 | NPOC     | CCV                | 0L-10900100866-10101000-134B17724A87-0008 |
| 16 | NPOC     | CCB                | 0L-10900100866-10101000-134B17724A87-0009 |
| 17 | NPOC     | 280-176363-c-1     | 0L-10900100866-10101000-134B17724A87-000A |
| 18 | NPOC     | 280-176363-a-2     | 0L-10900100866-10101000-134B17724A87-000B |
| 19 | NPOC     | 280-176363-c-3     | 0L-10900100866-10101000-134B17724A87-000C |
| 20 | NPOC     | 280-176363-d-5     | 0L-10900100866-10101000-134B17724A87-000D |
| 21 | NPOC     | MS 280-176363-d-5  | 0L-10900100866-10101000-134B17724A87-000E |
| 22 | NPOC     | MSD 280-176363-d-5 | 0L-10900100866-10101000-134B17724A87-000F |
| 23 | NPOC     | 280-176363-d-6     | 0L-10900100866-10101000-134B17724A87-0010 |
| 24 | NPOC     | 280-176363-c-7     | 0L-10900100866-10101000-134B17724A87-0011 |
| 25 | NPOC     | 280-176363-d-8     | 0L-10900100866-10101000-134B17724A87-0012 |
| 26 | NPOC     | 280-176325-d-1     | 0L-10900100866-10101000-134B17724A87-0013 |
| 27 | NPOC     | CCV                | 0L-10900100866-10101000-134B17724A87-0014 |
| 28 | NPOC     | CCB                | 0L-10900100866-10101000-134B17724A87-0015 |
| 29 | NPOC     | 280-176325-d-2     | 0L-10900100866-10101000-134B17724A87-0016 |
| 30 | NPOC     | 280-176274-d-2     | 0L-10900100866-10101000-134B17724A87-0017 |
| 31 | NPOC     | 280-176296-g-4     | 0L-10900100866-10101000-134B17724A87-0018 |
| 32 | NPOC     | 590-20532-a-9      | 0L-10900100866-10101000-134B17724A87-0019 |
| 33 | NPOC     | 590-20532-a-10     | 0L-10900100866-10101000-134B17724A87-001A |
| 34 | NPOC     | LCS                | 0L-10900100866-10101000-134B17724A87-001B |
| 35 | NPOC     | LCSD               | 0L-10900100866-10101000-134B17724A87-001C |
| 36 | NPOC     | MB                 | 0L-10900100866-10101000-134B17724A87-001D |
| 37 | NPOC     | TIC                | 0L-10900100866-10101000-134B17724A87-001E |
| 38 | NPOC     | 280-176296-g-3     | 0L-10900100866-10101000-134B17724A87-001F |
| 39 | NPOC     | CCV                | 0L-10900100866-10101000-134B17724A87-0020 |
| 40 | NPOC     | CCB                | 0L-10900100866-10101000-134B17724A87-0021 |
| 41 | NPOC     | MS 280-176296-g-3  | 0L-10900100866-10101000-134B17724A87-0022 |
| 42 | NPOC     | MSD 280-176296-g-3 | 0L-10900100866-10101000-134B17724A87-0023 |
| 43 | NPOC     | 280-176296-g-5     | 0L-10900100866-10101000-134B17724A87-0024 |
| 44 | NPOC     | 280-176296-g-7     | 0L-10900100866-10101000-134B17724A87-0025 |
| 45 | NPOC     | 280-176296-g-8     | 0L-10900100866-10101000-134B17724A87-0026 |
| 46 | NPOC     | 280-176296-g-9     | 0L-10900100866-10101000-134B17724A87-0027 |
| 47 | NPOC     | 590-20532-a-1      | 0L-10900100866-10101000-134B17724A87-0028 |
| 48 | NPOC     | 590-20532-a-2      | 0L-10900100866-10101000-134B17724A87-0029 |
| 49 | NPOC     | 590-20532-a-3      | 0L-10900100866-10101000-134B17724A87-002A |
| 50 | NPOC     | 590-20532-b-5      | 0L-10900100866-10101000-134B17724A87-002B |
| 51 | NPOC     | CCV                | 0L-10900100866-10101000-134B17724A87-002C |
| 52 | NPOC     | CCB                | 0L-10900100866-10101000-134B17724A87-002D |
| 53 | NPOC     | 590-20532-a-6      | 0L-10900100866-10101000-134B17724A87-002E |

|    | Origin                                                         | Manual Dilution | Result          | Status    |
|----|----------------------------------------------------------------|-----------------|-----------------|-----------|
| 1  | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:19.72mg/L  | Completed |
| 2  | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:0.2162mg/L | Completed |
| 3  | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:24.75mg/L  | Completed |
| 4  | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:0.2335mg/L | Completed |
| 5  | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:0.2530mg/L | Completed |
| 6  | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:0.8469mg/L | Completed |
| 7  | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:25.05mg/L  | Completed |
| 8  | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:25.05mg/L  | Completed |
| 9  | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:1.076mg/L  | Completed |
| 10 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:0.9633mg/L | Completed |
| 11 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:0.7734mg/L | Completed |
| 12 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 10.00           | NPOC:299.3mg/L  | Completed |
| 13 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 25.00           | NPOC:526.1mg/L  | Completed |
| 14 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 25.00           | NPOC:521.9mg/L  | Completed |
| 15 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:24.55mg/L  | Completed |
| 16 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:0.2058mg/L | Completed |
| 17 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 2.000           | NPOC:45.94mg/L  | Completed |
| 18 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:0.3234mg/L | Completed |
| 19 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 10.00           | NPOC:234.1mg/L  | Completed |
| 20 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 10.00           | NPOC:110.2mg/L  | Completed |
| 21 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 10.00           | NPOC:340.1mg/L  | Completed |
| 22 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 10.00           | NPOC:349.6mg/L  | Completed |
| 23 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:0.3406mg/L | Completed |
| 24 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 2.000           | NPOC:47.96mg/L  | Completed |
| 25 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 10.00           | NPOC:219.6mg/L  | Completed |
| 26 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 500.0           | NPOC:8359mg/L   | Completed |
| 27 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:24.68mg/L  | Completed |
| 28 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:0.2412mg/L | Completed |
| 29 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 500.0           | NPOC:4521mg/L   | Completed |
| 30 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 100.0           | NPOC:226.3mg/L  | Completed |
| 31 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:1.785mg/L  | Completed |
| 32 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:2.131mg/L  | Completed |
| 33 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:1.902mg/L  | Completed |
| 34 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:24.69mg/L  | Completed |
| 35 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:24.82mg/L  | Completed |
| 36 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:0.2179mg/L | Completed |
| 37 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:0.2542mg/L | Completed |
| 38 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:2.625mg/L  | Completed |
| 39 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:24.31mg/L  | Completed |
| 40 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:0.2449mg/L | Completed |
| 41 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:26.70mg/L  | Completed |
| 42 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:27.01mg/L  | Completed |
| 43 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:1.581mg/L  | Completed |
| 44 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:1.497mg/L  | Completed |
| 45 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:0.8961mg/L | Completed |
| 46 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:0.8706mg/L | Completed |
| 47 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:0.4120mg/L | Completed |
| 48 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:0.3433mg/L | Completed |
| 49 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:0.3826mg/L | Completed |
| 50 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:1.602mg/L  | Completed |
| 51 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:24.52mg/L  | Completed |
| 52 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:0.2184mg/L | Completed |
| 53 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:1.944mg/L  | Completed |

|    | Date / Time           | Vial |
|----|-----------------------|------|
| 1  | 5/19/2023 6:09:32 PM  | 1    |
| 2  | 5/19/2023 6:23:44 PM  | 2    |
| 3  | 5/19/2023 6:38:30 PM  | 3    |
| 4  | 5/19/2023 6:52:58 PM  | 4    |
| 5  | 5/19/2023 7:07:18 PM  | 5    |
| 6  | 5/19/2023 7:23:58 PM  | 6    |
| 7  | 5/19/2023 7:38:21 PM  | 7    |
| 8  | 5/19/2023 7:52:58 PM  | 8    |
| 9  | 5/19/2023 8:09:22 PM  | 9    |
| 10 | 5/19/2023 8:25:39 PM  | 10   |
| 11 | 5/19/2023 8:39:43 PM  | 11   |
| 12 | 5/19/2023 8:53:56 PM  | 12   |
| 13 | 5/19/2023 9:08:35 PM  | 13   |
| 14 | 5/19/2023 9:22:44 PM  | 14   |
| 15 | 5/19/2023 9:37:19 PM  | 15   |
| 16 | 5/19/2023 9:52:12 PM  | 16   |
| 17 | 5/19/2023 10:06:44 PM | 17   |
| 18 | 5/19/2023 10:22:59 PM | 18   |
| 19 | 5/19/2023 10:37:37 PM | 19   |
| 20 | 5/19/2023 10:51:46 PM | 20   |
| 21 | 5/19/2023 11:06:24 PM | 21   |
| 22 | 5/19/2023 11:20:45 PM | 22   |
| 23 | 5/19/2023 11:37:03 PM | 23   |
| 24 | 5/19/2023 11:51:23 PM | 24   |
| 25 | 5/20/2023 12:05:34 AM | 25   |
| 26 | 5/20/2023 12:19:50 AM | 26   |
| 27 | 5/20/2023 12:35:02 AM | 27   |
| 28 | 5/20/2023 12:49:11 AM | 28   |
| 29 | 5/20/2023 1:03:40 AM  | 29   |
| 30 | 5/20/2023 1:19:54 AM  | 30   |
| 31 | 5/20/2023 1:36:04 AM  | 31   |
| 32 | 5/20/2023 1:52:23 AM  | 32   |
| 33 | 5/20/2023 2:08:42 AM  | 33   |
| 34 | 5/20/2023 2:22:58 AM  | 34   |
| 35 | 5/20/2023 2:37:44 AM  | 35   |
| 36 | 5/20/2023 2:51:54 AM  | 36   |
| 37 | 5/20/2023 3:06:04 AM  | 37   |
| 38 | 5/20/2023 3:22:19 AM  | 38   |
| 39 | 5/20/2023 3:37:17 AM  | 39   |
| 40 | 5/20/2023 3:51:27 AM  | 40   |
| 41 | 5/20/2023 4:06:01 AM  | 41   |
| 42 | 5/20/2023 4:20:18 AM  | 42   |
| 43 | 5/20/2023 4:36:49 AM  | 43   |
| 44 | 5/20/2023 4:52:59 AM  | 44   |
| 45 | 5/20/2023 5:09:09 AM  | 45   |
| 46 | 5/20/2023 5:25:27 AM  | 46   |
| 47 | 5/20/2023 5:39:55 AM  | 47   |
| 48 | 5/20/2023 5:54:06 AM  | 48   |
| 49 | 5/20/2023 6:08:16 AM  | 49   |
| 50 | 5/20/2023 6:24:31 AM  | 50   |
| 51 | 5/20/2023 6:39:12 AM  | 51   |
| 52 | 5/20/2023 6:53:23 AM  | 52   |
| 53 | 5/20/2023 7:09:37 AM  | 53   |

|    | Analysis | Sample Name        | ObjectID                                  |
|----|----------|--------------------|-------------------------------------------|
| 54 | NPOC     | 280-176326-b-1     | 0L-10900100866-10101000-134B17724A87-002F |
| 55 | NPOC     | MS 280-176326-b-1  | 0L-10900100866-10101000-134B17724A87-0030 |
| 56 | NPOC     | MSD 280-176326-b-1 | 0L-10900100866-10101000-134B17724A87-0031 |
| 57 | NPOC     | 280-176326-b-2     | 0L-10900100866-10101000-134B17724A87-0032 |
| 58 | NPOC     | 280-176326-b-3     | 0L-10900100866-10101000-134B17724A87-0033 |
| 59 | NPOC     | 280-176326-b-4     | 0L-10900100866-10101000-134B17724A87-0034 |
| 60 | NPOC     | 280-176326-a-5     | 0L-10900100866-10101000-134B17724A87-0035 |
| 61 | NPOC     | 590-20532-b-7      | 0L-10900100866-10101000-134B17724A87-0036 |
| 62 | NPOC     | 590-20532-a-8      | 0L-10900100866-10101000-134B17724A87-0037 |
| 63 | NPOC     | CCV                | 0L-10900100866-10101000-134B17724A87-0038 |
| 64 | NPOC     | CCB                | 0L-10900100866-10101000-134B17724A87-0039 |
| 65 | NPOC     | 590-20532-a-4      | 0L-10900100866-10101000-134B17724A87-003A |
| 66 | NPOC     | 590-20532-a-5      | 0L-10900100866-10101000-134B17724A87-003B |
| 67 | NPOC     | 590-20532-b-6      | 0L-10900100866-10101000-134B17724A87-003C |
| 68 | NPOC     | LCS                | 0L-10900100866-10101000-134B17724A87-003D |
| 69 | NPOC     | MB                 | 0L-10900100866-10101000-134B17724A87-003E |
| 70 | NPOC     | 590-20532-a-7      | 0L-10900100866-10101000-134B17724A87-003F |
| 71 | NPOC     | MS 590-20532-a-7   | 0L-10900100866-10101000-134B17724A87-0040 |
| 72 | NPOC     | MSD 590-20532-a-7  | 0L-10900100866-10101000-134B17724A87-0041 |
| 73 | NPOC     | 590-20532-b-8      | 0L-10900100866-10101000-134B17724A87-0042 |
| 74 | NPOC     | 590-20532-b-9      | 0L-10900100866-10101000-134B17724A87-0043 |
| 75 | NPOC     | CCV                | 0L-10900100866-10101000-134B17724A87-0044 |
| 76 | NPOC     | CCB                | 0L-10900100866-10101000-134B17724A87-0045 |
| 77 | NPOC     | 590-20532-b-10     | 0L-10900100866-10101000-134B17724A87-0046 |
| 78 | NPOC     | 280-176081-j-3     | 0L-10900100866-10101000-134B17724A87-0047 |
| 79 | NPOC     | 280-176081-j-4     | 0L-10900100866-10101000-134B17724A87-0048 |
| 80 | NPOC     | 280-176081-j-5     | 0L-10900100866-10101000-134B17724A87-0049 |
| 81 | NPOC     | 280-176081-j-6     | 0L-10900100866-10101000-134B17724A87-004A |
| 82 | NPOC     | 280-176081-j-7     | 0L-10900100866-10101000-134B17724A87-004B |
| 83 | NPOC     | 280-176081-j-8     | 0L-10900100866-10101000-134B17724A87-004C |
| 84 | NPOC     | 280-176081-j-9     | 0L-10900100866-10101000-134B17724A87-004D |
| 85 | NPOC     | MS 280-176081-j-9  | 0L-10900100866-10101000-134B17724A87-004E |
| 86 | NPOC     | MSD 280-176081-j-9 | 0L-10900100866-10101000-134B17724A87-004F |
| 87 | NPOC     | CCV                | 0L-10900100866-10101000-134B17724A87-0050 |
| 88 | NPOC     | CCB                | 0L-10900100866-10101000-134B17724A87-0051 |
| 89 | NPOC     | 280-176081-j-10    | 0L-10900100866-10101000-134B17724A87-0052 |
| 90 | NPOC     | 280-176081-j-11    | 0L-10900100866-10101000-134B17724A87-0053 |
| 91 | NPOC     | 280-176081-j-12    | 0L-10900100866-10101000-134B17724A87-0054 |
| 92 | NPOC     | CCV                | 0L-10900100866-10101000-134B17724A87-0055 |
| 93 | NPOC     | CCB                | 0L-10900100866-10101000-134B17724A87-0056 |

|    | Origin                                                         | Manual Dilution | Result          | Status    |
|----|----------------------------------------------------------------|-----------------|-----------------|-----------|
| 54 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:0.8675mg/L | Completed |
| 55 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:25.16mg/L  | Completed |
| 56 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:24.99mg/L  | Completed |
| 57 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:1.143mg/L  | Completed |
| 58 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:1.043mg/L  | Completed |
| 59 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:0.5668mg/L | Completed |
| 60 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:0.8970mg/L | Completed |
| 61 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:1.996mg/L  | Completed |
| 62 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:1.972mg/L  | Completed |
| 63 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:24.33mg/L  | Completed |
| 64 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:0.2680mg/L | Completed |
| 65 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:0.4453mg/L | Completed |
| 66 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:1.620mg/L  | Completed |
| 67 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:1.911mg/L  | Completed |
| 68 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:24.81mg/L  | Completed |
| 69 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:0.2352mg/L | Completed |
| 70 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:2.123mg/L  | Completed |
| 71 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:26.53mg/L  | Completed |
| 72 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:26.40mg/L  | Completed |
| 73 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:1.986mg/L  | Completed |
| 74 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:2.082mg/L  | Completed |
| 75 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:24.87mg/L  | Completed |
| 76 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:0.2183mg/L | Completed |
| 77 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:1.906mg/L  | Completed |
| 78 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 10.00           | NPOC:11.35mg/L  | Completed |
| 79 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 2.000           | NPOC:8.942mg/L  | Completed |
| 80 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 2.000           | NPOC:6.013mg/L  | Completed |
| 81 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 2.000           | NPOC:6.591mg/L  | Completed |
| 82 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 2.000           | NPOC:0.9532mg/L | Completed |
| 83 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 2.000           | NPOC:14.64mg/L  | Completed |
| 84 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 2.000           | NPOC:5.626mg/L  | Completed |
| 85 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 2.000           | NPOC:54.65mg/L  | Completed |
| 86 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 2.000           | NPOC:54.46mg/L  | Completed |
| 87 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:24.75mg/L  | Completed |
| 88 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:0.2674mg/L | Completed |
| 89 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 2.000           | NPOC:1.446mg/L  | Completed |
| 90 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 2.000           | NPOC:6.315mg/L  | Completed |
| 91 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 2.000           | NPOC:4.506mg/L  | Completed |
| 92 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:24.82mg/L  | Completed |
| 93 | Calibration_10_31_22_Auto Dilution_Cal.2023_05_04_14_32_31.cal | 1.000           | NPOC:0.2411mg/L | Completed |

|    | Date / Time           | Vial |
|----|-----------------------|------|
| 54 | 5/20/2023 7:25:48 AM  | 54   |
| 55 | 5/20/2023 7:40:14 AM  | 55   |
| 56 | 5/20/2023 7:54:37 AM  | 56   |
| 57 | 5/20/2023 8:11:17 AM  | 57   |
| 58 | 5/20/2023 8:27:35 AM  | 58   |
| 59 | 5/20/2023 8:43:47 AM  | 59   |
| 60 | 5/20/2023 8:57:58 AM  | 60   |
| 61 | 5/20/2023 9:14:16 AM  | 61   |
| 62 | 5/20/2023 9:28:23 AM  | 62   |
| 63 | 5/20/2023 9:43:09 AM  | 63   |
| 64 | 5/20/2023 9:57:24 AM  | 64   |
| 65 | 5/20/2023 10:11:35 AM | 65   |
| 66 | 5/20/2023 10:25:46 AM | 66   |
| 67 | 5/20/2023 10:40:05 AM | 67   |
| 68 | 5/20/2023 10:54:21 AM | 68   |
| 69 | 5/20/2023 11:10:37 AM | 69   |
| 70 | 5/20/2023 11:26:47 AM | 70   |
| 71 | 5/20/2023 11:41:23 AM | 71   |
| 72 | 5/20/2023 11:55:33 AM | 72   |
| 73 | 5/20/2023 12:11:49 PM | 73   |
| 74 | 5/20/2023 12:28:10 PM | 74   |
| 75 | 5/20/2023 12:42:57 PM | 75   |
| 76 | 5/20/2023 12:57:08 PM | 76   |
| 77 | 5/20/2023 1:11:19 PM  | 77   |
| 78 | 5/20/2023 1:25:30 PM  | 78   |
| 79 | 5/20/2023 1:42:05 PM  | 79   |
| 80 | 5/20/2023 1:58:53 PM  | 80   |
| 81 | 5/20/2023 2:13:21 PM  | 81   |
| 82 | 5/20/2023 2:27:32 PM  | 82   |
| 83 | 5/20/2023 2:41:43 PM  | 83   |
| 84 | 5/20/2023 2:57:59 PM  | 84   |
| 85 | 5/20/2023 3:12:34 PM  | 85   |
| 86 | 5/20/2023 3:27:21 PM  | 86   |
| 87 | 5/20/2023 3:42:03 PM  | 87   |
| 88 | 5/20/2023 3:56:15 PM  | 88   |
| 89 | 5/20/2023 4:10:27 PM  | 89   |
| 90 | 5/20/2023 4:26:44 PM  | 90   |
| 91 | 5/20/2023 4:43:04 PM  | 91   |
| 92 | 5/20/2023 4:57:54 PM  | 92   |
| 93 | 5/20/2023 5:12:10 PM  | 93   |

Date of Creation 5/4/2023 3:58:47 PM  
User ABW  
System TOC-LCSN

## Cal. Curve

Sample Name: SHI5\_NPOC Auto Dilution Calibration\_05\_04\_2023  
Sample ID: SHI5\_NPOC Auto Dilution Calibration\_05\_04\_2023  
Object ID: 0L-10900100866-10101000-134B16820239-0000  
Cal. Curve: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Comment:

| Type     | Anal. |
|----------|-------|
| Standard | NPOC  |

Conc: 0.000mg/L

| No. | Area  | Inj. Vol. | Aut. Dil. | Rem.  | Ex. | Date / Time         |
|-----|-------|-----------|-----------|-------|-----|---------------------|
| 1   | 0.000 | 50uL      | 1.000     | ***** |     | 5/4/2023 2:40:33 PM |
| 2   | 0.000 | 50uL      | 1.000     | ***** |     | 5/4/2023 2:42:44 PM |
| 3   | 0.000 | 50uL      | 1.000     | ***** |     | 5/4/2023 2:44:53 PM |
| 4   | 0.000 | 50uL      | 1.000     | ***** |     | 5/4/2023 2:47:01 PM |

Acid Add. 1.500%  
Spurge Gas Flow 80mL/min  
Sp. Time 90.00sec  
Mean Area 0.000  
SD Area 0.000  
CV Area 0.00%  
Vial 92  
No Washes 2

Conc: 1.000mg/L

| No. | Area  | Inj. Vol. | Aut. Dil. | Rem.  | Ex. | Date / Time         |
|-----|-------|-----------|-----------|-------|-----|---------------------|
| 1   | 1.771 | 50uL      | 50.00     | ***** |     | 5/4/2023 2:55:42 PM |
| 2   | 1.714 | 50uL      | 50.00     | ***** |     | 5/4/2023 2:57:50 PM |
| 3   | 1.595 | 50uL      | 50.00     | ***** |     | 5/4/2023 2:59:59 PM |
| 4   | 1.786 | 50uL      | 50.00     | ***** |     | 5/4/2023 3:02:12 PM |

Acid Add. 1.500%  
Spurge Gas Flow 80mL/min  
Sp. Time 90.00sec  
Mean Area 1.716  
SD Area 0.08674  
CV Area 5.05%  
Vial 93  
No Washes 2

Conc: 5.000mg/L

| No. | Area  | Inj. Vol. | Aut. Dil. | Rem.  | Ex. | Date / Time         |
|-----|-------|-----------|-----------|-------|-----|---------------------|
| 1   | 10.60 | 50uL      | 10.00     | ***** |     | 5/4/2023 3:09:45 PM |
| 2   | 10.29 | 50uL      | 10.00     | ***** |     | 5/4/2023 3:11:54 PM |
| 3   | 10.25 | 50uL      | 10.00     | ***** |     | 5/4/2023 3:14:04 PM |
| 4   | 10.32 | 50uL      | 10.00     | ***** |     | 5/4/2023 3:16:14 PM |

Acid Add. 1.500%  
Spurge Gas Flow 80mL/min  
Sp. Time 90.00sec  
Mean Area 10.37  
SD Area 0.1593  
CV Area 1.54%  
Vial 93  
No Washes 2

Conc: 10.00mg/L



| No. | Area  | Inj. Vol. | Aut. Dil. | Rem.  | Ex. | Date / Time         |
|-----|-------|-----------|-----------|-------|-----|---------------------|
| 1   | 20.19 | 50uL      | 5.000     | ***** |     | 5/4/2023 3:23:47 PM |
| 2   | 19.43 | 50uL      | 5.000     | ***** |     | 5/4/2023 3:25:57 PM |
| 3   | 19.44 | 50uL      | 5.000     | ***** |     | 5/4/2023 3:28:06 PM |
| 4   | 19.84 | 50uL      | 5.000     | ***** |     | 5/4/2023 3:30:14 PM |

Acid Add. 1.500%  
 Sparge Gas Flow 80mL/min  
 Sp. Time 90.00sec  
 Mean Area 19.73  
 SD Area 0.3641  
 CV Area 1.85%  
 Vial 93  
 No Washes 2

Conc: 25.00mg/L

| No. | Area  | Inj. Vol. | Aut. Dil. | Rem.  | Ex. | Date / Time         |
|-----|-------|-----------|-----------|-------|-----|---------------------|
| 1   | 51.32 | 50uL      | 2.000     | ***** |     | 5/4/2023 3:38:03 PM |
| 2   | 49.72 | 50uL      | 2.000     | ***** |     | 5/4/2023 3:40:19 PM |
| 3   | 49.18 | 50uL      | 2.000     | ***** |     | 5/4/2023 3:42:38 PM |
| 4   | 50.28 | 50uL      | 2.000     | ***** |     | 5/4/2023 3:45:02 PM |


Acid Add. 1.500%  
 Sparge Gas Flow 80mL/min  
 Sp. Time 90.00sec  
 Mean Area 50.13  
 SD Area 0.9145  
 CV Area 1.82%  
 Vial 93  
 No Washes 2

Conc: 50.00mg/L

| No. | Area  | Inj. Vol. | Aut. Dil. | Rem.  | Ex. | Date / Time         |
|-----|-------|-----------|-----------|-------|-----|---------------------|
| 1   | 103.1 | 50uL      | 1.000     | ***** |     | 5/4/2023 3:51:58 PM |
| 2   | 100.5 | 50uL      | 1.000     | ***** |     | 5/4/2023 3:54:13 PM |
| 3   | 99.70 | 50uL      | 1.000     | ***** |     | 5/4/2023 3:56:25 PM |
| 4   | 101.6 | 50uL      | 1.000     | ***** |     | 5/4/2023 3:58:47 PM |

Acid Add. 1.500%  
 Sparge Gas Flow 80mL/min  
 Sp. Time 90.00sec  
 Mean Area 101.2  
 SD Area 1.473  
 CV Area 1.45%  
 Vial 93  
 No Washes 2

Slope: 2.024  
 Intercept -0.1771  
 r<sup>2</sup> 0.9999  
 r 1.0000  
 RSE(%) N/A  
 Zero Shift No



# TOC-Control L Report

ABW  
2023\_05\_19\_001.tx

## Instr.Information

Instrument Options  
Catalyst

TOC/ASI/  
Regular Sensitivity

## Sample

Sample Name: ICV  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 1.000           | NPOC:19.72mg/L |

1. Det

Anal.: NPOC

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 40.95 | 20.32mg/L | 50uL      | 1.000     | E   | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 6:00:51 PM |
| 2   | 39.32 | 19.51mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 6:03:08 PM |
| 3   | 39.05 | 19.38mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 6:05:17 PM |
| 4   | 40.02 | 19.86mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 6:07:27 PM |
| 5   | 40.58 | 20.13mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 6:09:32 PM |

Mean Area 39.74  
Mean Conc. 19.72mg/L



## Sample

Sample Name: ICB  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result

| Type    | Anal. | Manual Dilution | Result          |
|---------|-------|-----------------|-----------------|
| Unknown | NPOC  | 1.000           | NPOC:0.2162mg/L |

1. Det


Anal.: NPOC

# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx

| No. | Area   | Conc.      | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|--------|------------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 0.2875 | 0.2295mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 6:17:21 PM |
| 2   | 0.2897 | 0.2306mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 6:19:30 PM |
| 3   | 0.2334 | 0.2028mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 6:21:39 PM |
| 4   | 0.2315 | 0.2019mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 6:23:44 PM |

Mean Area 0.2605  
Mean Conc. 0.2162mg/L



## Sample

Sample Name: LCS  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 1.000           | NPOC:24.75mg/L |

1. Det

Anal.: NPOC

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 50.52 | 25.04mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 6:31:44 PM |
| 2   | 49.42 | 24.50mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 6:34:00 PM |
| 3   | 49.46 | 24.52mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 6:36:16 PM |
| 4   | 50.34 | 24.95mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 6:38:30 PM |

Mean Area 49.94  
Mean Conc. 24.75mg/L



## Sample

Sample Name: MB  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result

# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx


| Type    | Anal. | Manual Dilution | Result          |
|---------|-------|-----------------|-----------------|
| Unknown | NPOC  | 1.000           | NPOC:0.2335mg/L |

1. Det

Anal.: NPOC

| No. | Area   | Conc.      | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|--------|------------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 0.3467 | 0.2588mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 6:46:24 PM |
| 2   | 0.2552 | 0.2136mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 6:48:39 PM |
| 3   | 0.2324 | 0.2023mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 6:50:47 PM |
| 4   | 0.3480 | 0.2594mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 6:52:58 PM |

Mean Area 0.2956  
Mean Conc. 0.2335mg/L



## Sample

Sample Name: TIC  
Sample ID: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Origin: Completed  
Status:  
Chk. Result:


| Type    | Anal. | Manual Dilution | Result          |
|---------|-------|-----------------|-----------------|
| Unknown | NPOC  | 1.000           | NPOC:0.2530mg/L |

1. Det

Anal.: NPOC

| No. | Area   | Conc.      | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|--------|------------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 0.3570 | 0.2639mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 7:00:52 PM |
| 2   | 0.3151 | 0.2432mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 7:03:02 PM |
| 3   | 0.3526 | 0.2617mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 7:05:10 PM |
| 4   | 0.3152 | 0.2432mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 7:07:18 PM |

Mean Area 0.3350  
Mean Conc. 0.2530mg/L



# TOC-Control L Report

ABW  
2023\_05\_19\_001.tx

## Sample

Sample Name: 280-176674-d-7  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result          |
|---------|-------|-----------------|-----------------|
| Unknown | NPOC  | 1.000           | NPOC:0.8469mg/L |

1. Det

Anal.: NPOC

| No. | Area  | Conc.      | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|------------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 1.732 | 0.9431mg/L | 50uL      | 1.000     | E   | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 7:15:30 PM |
| 2   | 1.583 | 0.8695mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 7:17:39 PM |
| 3   | 1.464 | 0.8107mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 7:19:48 PM |
| 4   | 1.559 | 0.8576mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 7:21:53 PM |
| 5   | 1.543 | 0.8497mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 7:23:58 PM |

Mean Area 1.537  
Mean Conc. 0.8469mg/L



## Sample

Sample Name: MS 280-176674-d-7  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result

| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 1.000           | NPOC:25.05mg/L |

1. Det


Anal.: NPOC

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 50.84 | 25.20mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 7:32:00 PM |
| 2   | 50.55 | 25.06mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 7:34:06 PM |
| 3   | 50.15 | 24.86mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 7:36:16 PM |
| 4   | 50.57 | 25.07mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 7:38:21 PM |

# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx

Mean Area 50.53  
Mean Conc. 25.05mg/L



## Sample

Sample Name: MSD 280-176674-d-7  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 1.000           | NPOC:25.05mg/L |

1. Det

Anal.: NPOC

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 51.41 | 25.48mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 7:46:20 PM |
| 2   | 49.70 | 24.64mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 7:48:29 PM |
| 3   | 50.10 | 24.84mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 7:50:45 PM |
| 4   | 50.94 | 25.25mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 7:52:58 PM |

Mean Area 50.54  
Mean Conc. 25.05mg/L



## Sample

Sample Name: 580-127209-p-1  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result

| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 1.000           | NPOC:1.076mg/L |

1. Det


Anal.: NPOC

# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 2.285 | 1.216mg/L | 50uL      | 1.000     | E   | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 8:00:52 PM |
| 2   | 1.959 | 1.055mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 8:02:57 PM |
| 3   | 2.022 | 1.086mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 8:05:07 PM |
| 4   | 1.965 | 1.058mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 8:07:17 PM |
| 5   | 2.057 | 1.104mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 8:09:22 PM |

Mean Area 2.001  
Mean Conc. 1.076mg/L



## Sample

Sample Name: 580-127209-p-3  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result          |
|---------|-------|-----------------|-----------------|
| Unknown | NPOC  | 1.000           | NPOC:0.9633mg/L |

1. Det

Anal.: NPOC

| No. | Area  | Conc.      | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|------------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 1.868 | 1.010mg/L  | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 8:17:11 PM |
| 2   | 1.607 | 0.8813mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 8:19:20 PM |
| 3   | 1.532 | 0.8443mg/L | 50uL      | 1.000     | E   | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 8:21:29 PM |
| 4   | 1.844 | 0.9984mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 8:23:34 PM |
| 5   | 1.773 | 0.9633mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 8:25:39 PM |

Mean Area 1.773  
Mean Conc. 0.9633mg/L



## Sample

# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx

Sample Name: 580-127209-p-5  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result          |
|---------|-------|-----------------|-----------------|
| Unknown | NPOC  | 1.000           | NPOC:0.7734mg/L |

1. Det

Anal.: NPOC

| No. | Area  | Conc.      | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|------------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 1.459 | 0.8082mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 8:33:28 PM |
| 2   | 1.383 | 0.7707mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 8:35:33 PM |
| 3   | 1.365 | 0.7618mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 8:37:38 PM |
| 4   | 1.347 | 0.7529mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 8:39:43 PM |

Mean Area 1.389  
Mean Conc. 0.7734mg/L



## Sample

Sample Name: 280-176349-d-8  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result

| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 10.00           | NPOC:299.3mg/L |

1. Det

Anal.: NPOC


| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 61.59 | 305.1mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 8:47:37 PM |
| 2   | 59.61 | 295.3mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 8:49:46 PM |
| 3   | 59.74 | 296.0mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 8:51:51 PM |
| 4   | 60.68 | 300.6mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 8:53:56 PM |



# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx

Mean Area 60.41  
Mean Conc. 299.3mg/L



## Sample

Sample Name: 280-176349-c-11  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 25.00           | NPOC:526.1mg/L |

1. Det

Anal.: NPOC

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 42.85 | 531.4mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 9:01:50 PM |
| 2   | 41.96 | 520.4mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 9:03:59 PM |
| 3   | 42.13 | 522.5mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 9:06:04 PM |
| 4   | 42.75 | 530.1mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 9:08:35 PM |

Mean Area 42.42  
Mean Conc. 526.1mg/L



## Sample

Sample Name: 280-176123-a-35  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result

| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 25.00           | NPOC:521.9mg/L |

1. Det


Anal.: NPOC

# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 42.91 | 532.1mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 9:16:29 PM |
| 2   | 41.60 | 515.9mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 9:18:34 PM |
| 3   | 41.49 | 514.6mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 9:20:39 PM |
| 4   | 42.35 | 525.2mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 9:22:44 PM |

Mean Area 42.09  
Mean Conc. 521.9mg/L



## Sample

Sample Name: CCB  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 1.000           | NPOC:24.55mg/L |

1. Det

Anal.: NPOC

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 50.47 | 25.02mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 9:30:48 PM |
| 2   | 48.96 | 24.27mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 9:33:07 PM |
| 3   | 48.99 | 24.29mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 9:35:14 PM |
| 4   | 49.63 | 24.60mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 9:37:19 PM |

Mean Area 49.51  
Mean Conc. 24.55mg/L



## Sample

Sample Name: CCB  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result

# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx


| Type    | Anal. | Manual Dilution | Result          |
|---------|-------|-----------------|-----------------|
| Unknown | NPOC  | 1.000           | NPOC:0.2058mg/L |

1. Det

Anal.: NPOC

| No. | Area   | Conc.      | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|--------|------------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 0.2882 | 0.2299mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 9:45:13 PM |
| 2   | 0.2427 | 0.2074mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 9:47:18 PM |
| 3   | 0.2121 | 0.1923mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 9:49:53 PM |
| 4   | 0.2147 | 0.1936mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 9:52:12 PM |

Mean Area 0.2394  
Mean Conc. 0.2058mg/L



## Sample

Sample Name: 280-176363-c-1  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 2.000           | NPOC:45.94mg/L |

1. Det

Anal.: NPOC

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time           |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|-----------------------|
| 1   | 46.90 | 46.51mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 10:00:17 PM |
| 2   | 45.93 | 45.55mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 10:02:25 PM |
| 3   | 46.01 | 45.63mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 10:04:32 PM |
| 4   | 46.45 | 46.07mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 10:06:44 PM |

Mean Area 46.32  
Mean Conc. 45.94mg/L



# TOC-Control L Report

ABW  
2023\_05\_19\_001.tx

## Sample

Sample Name: 280-176363-a-2  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result          |
|---------|-------|-----------------|-----------------|
| Unknown | NPOC  | 1.000           | NPOC:0.3234mg/L |

1. Det

Anal.: NPOC

| No. | Area   | Conc.      | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time           |
|-----|--------|------------|-----------|-----------|-----|-------------------------------------------------------------------|-----------------------|
| 1   | 1.424  | 0.7909mg/L | 50uL      | 1.000     | E   | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 10:14:38 PM |
| 2   | 0.6734 | 0.4201mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 10:16:43 PM |
| 3   | 0.4585 | 0.3140mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 10:18:49 PM |
| 4   | 0.3988 | 0.2845mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 10:20:54 PM |
| 5   | 0.3799 | 0.2752mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 10:22:59 PM |

Mean Area 0.4777  
Mean Conc. 0.3234mg/L



## Sample

Sample Name: 280-176363-c-3  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result

| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 10.00           | NPOC:234.1mg/L |

1. Det


Anal.: NPOC

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time           |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|-----------------------|
| 1   | 47.82 | 237.1mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 10:31:05 PM |
| 2   | 46.59 | 231.0mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 10:33:18 PM |
| 3   | 47.20 | 234.0mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 10:35:29 PM |
| 4   | 47.26 | 234.3mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 10:37:37 PM |

# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx

Mean Area 47.22  
Mean Conc. 234.1mg/L



## Sample

Sample Name: 280-176363-d-5  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 10.00           | NPOC:110.2mg/L |

1. Det

Anal.: NPOC

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time           |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|-----------------------|
| 1   | 22.58 | 112.4mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 10:45:31 PM |
| 2   | 21.96 | 109.4mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 10:47:36 PM |
| 3   | 21.91 | 109.1mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 10:49:41 PM |
| 4   | 22.11 | 110.1mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 10:51:46 PM |

Mean Area 22.14  
Mean Conc. 110.2mg/L



## Sample

Sample Name: MS 280-176363-d-5  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result

| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 10.00           | NPOC:340.1mg/L |

1. Det


Anal.: NPOC

# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time           |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|-----------------------|
| 1   | 69.67 | 345.0mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 10:59:58 PM |
| 2   | 67.96 | 336.6mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 11:02:10 PM |
| 3   | 67.92 | 336.4mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 11:04:15 PM |
| 4   | 69.11 | 342.3mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 11:06:24 PM |

Mean Area 68.67  
Mean Conc. 340.1mg/L



## Sample

Sample Name: MSD 280-176363-d-5  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 10.00           | NPOC:349.6mg/L |

1. Det

Anal.: NPOC

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time           |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|-----------------------|
| 1   | 71.47 | 353.9mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 11:14:26 PM |
| 2   | 70.24 | 347.8mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 11:16:35 PM |
| 3   | 69.86 | 346.0mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 11:18:40 PM |
| 4   | 70.79 | 350.6mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 11:20:45 PM |

Mean Area 70.59  
Mean Conc. 349.6mg/L



## Sample

Sample Name: 280-176363-d-6  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result

# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx


| Type    | Anal. | Manual Dilution | Result          |
|---------|-------|-----------------|-----------------|
| Unknown | NPOC  | 1.000           | NPOC:0.3406mg/L |

1. Det

Anal.: NPOC

| No. | Area   | Conc.      | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time           |
|-----|--------|------------|-----------|-----------|-----|-------------------------------------------------------------------|-----------------------|
| 1   | 0.6960 | 0.4313mg/L | 50uL      | 1.000     | E   | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 11:28:39 PM |
| 2   | 0.5866 | 0.3773mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 11:30:48 PM |
| 3   | 0.4924 | 0.3307mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 11:32:53 PM |
| 4   | 0.4737 | 0.3215mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 11:34:58 PM |
| 5   | 0.4968 | 0.3329mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 11:37:03 PM |

Mean Area 0.5124  
Mean Conc. 0.3406mg/L



## Sample

Sample Name: 280-176363-c-7  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 2.000           | NPOC:47.96mg/L |

1. Det

Anal.: NPOC

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time           |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|-----------------------|
| 1   | 49.06 | 48.64mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 11:44:59 PM |
| 2   | 48.07 | 47.67mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 11:47:10 PM |
| 3   | 47.98 | 47.58mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 11:49:18 PM |
| 4   | 48.38 | 47.97mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 11:51:23 PM |

Mean Area 48.37  
Mean Conc. 47.96mg/L



# TOC-Control L Report

ABW  
2023\_05\_19\_001.tx

## Sample

Sample Name: 280-176363-d-8  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 10.00           | NPOC:219.6mg/L |

1. Det

Anal.: NPOC

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time           |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|-----------------------|
| 1   | 45.10 | 223.7mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/19/2023 11:59:18 PM |
| 2   | 43.82 | 217.3mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 12:01:23 AM |
| 3   | 43.84 | 217.4mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 12:03:29 AM |
| 4   | 44.34 | 219.9mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 12:05:34 AM |

Mean Area 44.27  
Mean Conc. 219.6mg/L



## Sample

Sample Name: 280-176325-d-1  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result

| Type    | Anal. | Manual Dilution | Result        |
|---------|-------|-----------------|---------------|
| Unknown | NPOC  | 500.0           | NPOC:8359mg/L |

1. Det

Anal.: NPOC


| No. | Area  | Conc.    | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time           |
|-----|-------|----------|-----------|-----------|-----|-------------------------------------------------------------------|-----------------------|
| 1   | 34.43 | 8548mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 12:13:28 AM |
| 2   | 33.36 | 8283mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 12:15:37 AM |
| 3   | 33.05 | 8207mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 12:17:42 AM |
| 4   | 33.82 | 8397mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 12:19:50 AM |



# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx

Mean Area 33.66  
Mean Conc. 8359mg/L



## Sample

Sample Name: CCV  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 1.000           | NPOC:24.68mg/L |

1. Det

Anal.: NPOC

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time           |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|-----------------------|
| 1   | 50.21 | 24.89mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 12:28:01 AM |
| 2   | 49.21 | 24.40mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 12:30:20 AM |
| 3   | 50.02 | 24.80mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 12:32:42 AM |
| 4   | 49.73 | 24.65mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 12:35:02 AM |

Mean Area 49.79  
Mean Conc. 24.68mg/L



## Sample

Sample Name: CCB  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result

| Type    | Anal. | Manual Dilution | Result          |
|---------|-------|-----------------|-----------------|
| Unknown | NPOC  | 1.000           | NPOC:0.2412mg/L |

1. Det


Anal.: NPOC

# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx

| No. | Area   | Conc.      | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time           |
|-----|--------|------------|-----------|-----------|-----|-------------------------------------------------------------------|-----------------------|
| 1   | 0.2843 | 0.2279mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 12:42:56 AM |
| 2   | 0.2469 | 0.2095mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 12:45:01 AM |
| 3   | 0.2716 | 0.2217mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 12:47:06 AM |
| 4   | 0.4420 | 0.3058mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 12:49:11 AM |

Mean Area 0.3112  
Mean Conc. 0.2412mg/L



## Sample

Sample Name: 280-176325-d-2  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result        |
|---------|-------|-----------------|---------------|
| Unknown | NPOC  | 500.0           | NPOC:4521mg/L |

1. Det

Anal.: NPOC

| No. | Area  | Conc.    | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time           |
|-----|-------|----------|-----------|-----------|-----|-------------------------------------------------------------------|-----------------------|
| 1   | 18.33 | 4571mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 12:57:05 AM |
| 2   | 17.99 | 4487mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 12:59:14 AM |
| 3   | 18.00 | 4490mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 1:01:35 AM  |
| 4   | 18.19 | 4536mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 1:03:40 AM  |

Mean Area 18.13  
Mean Conc. 4521mg/L



## Sample

Sample Name: 280-176274-d-2  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result

# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx

| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 100.0           | NPOC:226.3mg/L |


1. Det

Anal.: NPOC

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 4.606 | 236.3mg/L | 50uL      | 1.000     | E   | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 1:11:34 AM |
| 2   | 4.430 | 227.6mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 1:13:39 AM |
| 3   | 4.359 | 224.1mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 1:15:44 AM |
| 4   | 4.363 | 224.3mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 1:17:49 AM |
| 5   | 4.466 | 229.4mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 1:19:54 AM |

Mean Area  
Mean Conc.

4.405  
226.3mg/L



## Sample

Sample Name: 280-176296-g-4  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result

| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 1.000           | NPOC:1.785mg/L |


1. Det

Anal.: NPOC

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 3.446 | 1.790mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 1:27:44 AM |
| 2   | 3.274 | 1.705mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 1:29:49 AM |
| 3   | 3.164 | 1.650mg/L | 50uL      | 1.000     | E   | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 1:31:54 AM |
| 4   | 3.510 | 1.821mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 1:33:59 AM |
| 5   | 3.518 | 1.825mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 1:36:04 AM |

Mean Area  
Mean Conc.

3.437  
1.785mg/L



# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx

## Sample

Sample Name: 590-20532-a-9  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 1.000           | NPOC:2.131mg/L |

1. Det

Anal.: NPOC

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 4.847 | 2.482mg/L | 50uL      | 1.000     | E   | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 1:43:54 AM |
| 2   | 4.324 | 2.223mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 1:45:59 AM |
| 3   | 4.113 | 2.119mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 1:48:04 AM |
| 4   | 4.095 | 2.110mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 1:50:13 AM |
| 5   | 4.015 | 2.071mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 1:52:23 AM |

Mean Area 4.137  
Mean Conc. 2.131mg/L



## Sample

Sample Name: 590-20532-a-10  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result

| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 1.000           | NPOC:1.902mg/L |

1. Det


Anal.: NPOC

# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 3.972 | 2.050mg/L | 50uL      | 1.000     | E   | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 2:00:13 AM |
| 2   | 3.706 | 1.918mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 2:02:18 AM |
| 3   | 3.706 | 1.918mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 2:04:23 AM |
| 4   | 3.712 | 1.921mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 2:06:28 AM |
| 5   | 3.569 | 1.851mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 2:08:42 AM |

Mean Area 3.673  
Mean Conc. 1.902mg/L



## Sample

Sample Name: LCS  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 1.000           | NPOC:24.69mg/L |

1. Det

Anal.: NPOC

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 51.11 | 25.33mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 2:16:34 AM |
| 2   | 49.01 | 24.30mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 2:18:43 AM |
| 3   | 49.44 | 24.51mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 2:20:53 AM |
| 4   | 49.63 | 24.60mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 2:22:58 AM |

Mean Area 49.80  
Mean Conc. 24.69mg/L



## Sample

# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx

Sample Name: LCSD  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 1.000           | NPOC:24.82mg/L |

1. Det

Anal.: NPOC

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 51.30 | 25.43mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 2:31:09 AM |
| 2   | 49.43 | 24.50mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 2:33:26 AM |
| 3   | 49.63 | 24.60mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 2:35:39 AM |
| 4   | 49.92 | 24.75mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 2:37:44 AM |

Mean Area 50.07  
Mean Conc. 24.82mg/L



## Sample

Sample Name: MB  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result

| Type    | Anal. | Manual Dilution | Result          |
|---------|-------|-----------------|-----------------|
| Unknown | NPOC  | 1.000           | NPOC:0.2179mg/L |

1. Det


Anal.: NPOC

| No. | Area   | Conc.      | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|--------|------------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 0.2723 | 0.2220mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 2:45:39 AM |
| 2   | 0.2274 | 0.1998mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 2:47:44 AM |
| 3   | 0.2572 | 0.2146mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 2:49:49 AM |
| 4   | 0.2991 | 0.2353mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 2:51:54 AM |

# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx

Mean Area 0.2640  
Mean Conc. 0.2179mg/L



## Sample

Sample Name: TIC  
Sample ID: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Origin: Completed  
Status:  
Chk. Result


| Type    | Anal. | Manual Dilution | Result          |
|---------|-------|-----------------|-----------------|
| Unknown | NPOC  | 1.000           | NPOC:0.2542mg/L |

1. Det

Anal.: NPOC

| No. | Area   | Conc.      | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|--------|------------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 0.3668 | 0.2687mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 2:59:49 AM |
| 2   | 0.3002 | 0.2358mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 3:01:54 AM |
| 3   | 0.3010 | 0.2362mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 3:03:59 AM |
| 4   | 0.3819 | 0.2762mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 3:06:04 AM |

Mean Area 0.3375  
Mean Conc. 0.2542mg/L



## Sample

Sample Name: 280-176296-g-3  
Sample ID: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Origin: Completed  
Status:  
Chk. Result

| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 1.000           | NPOC:2.625mg/L |

1. Det


Anal.: NPOC

# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 5.127 | 2.620mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 3:13:59 AM |
| 2   | 4.908 | 2.512mg/L | 50uL      | 1.000     | E   | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 3:16:04 AM |
| 3   | 5.065 | 2.589mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 3:18:09 AM |
| 4   | 5.175 | 2.644mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 3:20:14 AM |
| 5   | 5.181 | 2.647mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 3:22:19 AM |

Mean Area 5.137  
Mean Conc. 2.625mg/L



## Sample

Sample Name: CCV  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 1.000           | NPOC:24.31mg/L |

1. Det

Anal.: NPOC

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 49.96 | 24.77mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 3:30:27 AM |
| 2   | 48.59 | 24.09mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 3:32:40 AM |
| 3   | 48.51 | 24.05mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 3:35:03 AM |
| 4   | 49.08 | 24.33mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 3:37:17 AM |

Mean Area 49.04  
Mean Conc. 24.31mg/L



## Sample



# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx

Sample Name: CCB  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result          |
|---------|-------|-----------------|-----------------|
| Unknown | NPOC  | 1.000           | NPOC:0.2449mg/L |

1. Det

Anal.: NPOC

| No. | Area   | Conc.      | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|--------|------------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 0.3573 | 0.2640mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 3:45:11 AM |
| 2   | 0.3310 | 0.2510mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 3:47:17 AM |
| 3   | 0.2215 | 0.1969mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 3:49:22 AM |
| 4   | 0.3651 | 0.2679mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 3:51:27 AM |

Mean Area 0.3187  
Mean Conc. 0.2449mg/L



## Sample

Sample Name: MS 280-176296-g-3  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result

| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 1.000           | NPOC:26.70mg/L |

1. Det


Anal.: NPOC

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 55.05 | 27.28mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 3:59:42 AM |
| 2   | 53.44 | 26.49mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 4:01:50 AM |
| 3   | 53.15 | 26.34mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 4:03:56 AM |
| 4   | 53.85 | 26.69mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 4:06:01 AM |

# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx

Mean Area 53.87  
Mean Conc. 26.70mg/L



## Sample

Sample Name: MSD 280-176296-g-3  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 1.000           | NPOC:27.01mg/L |

1. Det

Anal.: NPOC

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 55.33 | 27.42mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 4:14:03 AM |
| 2   | 54.07 | 26.80mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 4:16:08 AM |
| 3   | 53.75 | 26.64mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 4:18:13 AM |
| 4   | 54.82 | 27.17mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 4:20:18 AM |

Mean Area 54.49  
Mean Conc. 27.01mg/L



## Sample

Sample Name: 280-176296-g-5  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result

| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 1.000           | NPOC:1.581mg/L |

1. Det


Anal.: NPOC

# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 3.279 | 1.707mg/L | 50uL      | 1.000     | E   | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 4:28:13 AM |
| 2   | 2.991 | 1.565mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 4:30:34 AM |
| 3   | 2.938 | 1.539mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 4:32:39 AM |
| 4   | 3.099 | 1.618mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 4:34:44 AM |
| 5   | 3.068 | 1.603mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 4:36:49 AM |

Mean Area 3.024  
Mean Conc. 1.581mg/L



## Sample

Sample Name: 280-176296-g-7  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 1.000           | NPOC:1.497mg/L |

1. Det

Anal.: NPOC

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 3.035 | 1.587mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 4:44:39 AM |
| 2   | 2.743 | 1.442mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 4:46:43 AM |
| 3   | 2.527 | 1.336mg/L | 50uL      | 1.000     | E   | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 4:48:49 AM |
| 4   | 2.766 | 1.454mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 4:50:54 AM |
| 5   | 2.867 | 1.504mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 4:52:59 AM |

Mean Area 2.853  
Mean Conc. 1.497mg/L



## Sample

# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx

Sample Name: 280-176296-g-8  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result          |
|---------|-------|-----------------|-----------------|
| Unknown | NPOC  | 1.000           | NPOC:0.8961mg/L |

1. Det

Anal.: NPOC

| No. | Area  | Conc.      | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|------------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 1.917 | 1.034mg/L  | 50uL      | 1.000     | E   | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 5:00:49 AM |
| 2   | 1.604 | 0.8798mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 5:02:54 AM |
| 3   | 1.592 | 0.8739mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 5:04:59 AM |
| 4   | 1.660 | 0.9075mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 5:07:04 AM |
| 5   | 1.692 | 0.9233mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 5:09:09 AM |

Mean Area 1.637  
Mean Conc. 0.8961mg/L



## Sample

Sample Name: 280-176296-g-9  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result

| Type    | Anal. | Manual Dilution | Result          |
|---------|-------|-----------------|-----------------|
| Unknown | NPOC  | 1.000           | NPOC:0.8706mg/L |

1. Det


Anal.: NPOC

| No. | Area  | Conc.      | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|------------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 1.781 | 0.9673mg/L | 50uL      | 1.000     | E   | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 5:16:59 AM |
| 2   | 1.557 | 0.8566mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 5:19:11 AM |
| 3   | 1.596 | 0.8759mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 5:21:17 AM |
| 4   | 1.558 | 0.8571mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 5:23:22 AM |
| 5   | 1.630 | 0.8927mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 5:25:27 AM |

# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx

Mean Area 1.585  
Mean Conc. 0.8706mg/L



## Sample

Sample Name: 590-20532-a-1  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result          |
|---------|-------|-----------------|-----------------|
| Unknown | NPOC  | 1.000           | NPOC:0.4120mg/L |

1. Det

Anal.: NPOC

| No. | Area   | Conc.      | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|--------|------------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 0.7423 | 0.4542mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 5:33:17 AM |
| 2   | 0.6393 | 0.4033mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 5:35:26 AM |
| 3   | 0.5976 | 0.3827mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 5:37:50 AM |
| 4   | 0.6481 | 0.4076mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 5:39:55 AM |

Mean Area 0.6568  
Mean Conc. 0.4120mg/L



## Sample

Sample Name: 590-20532-a-2  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result

| Type    | Anal. | Manual Dilution | Result          |
|---------|-------|-----------------|-----------------|
| Unknown | NPOC  | 1.000           | NPOC:0.3433mg/L |

1. Det


Anal.: NPOC

# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx

| No. | Area   | Conc.      | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|--------|------------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 0.5870 | 0.3775mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 5:47:50 AM |
| 2   | 0.5106 | 0.3397mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 5:49:56 AM |
| 3   | 0.4813 | 0.3253mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 5:52:00 AM |
| 4   | 0.4928 | 0.3309mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 5:54:06 AM |

Mean Area 0.5179  
Mean Conc. 0.3433mg/L



## Sample

Sample Name: 590-20532-a-3  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result          |
|---------|-------|-----------------|-----------------|
| Unknown | NPOC  | 1.000           | NPOC:0.3826mg/L |

1. Det

Anal.: NPOC

| No. | Area   | Conc.      | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|--------|------------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 0.5870 | 0.3775mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 6:02:01 AM |
| 2   | 0.6129 | 0.3903mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 6:04:05 AM |
| 3   | 0.5721 | 0.3701mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 6:06:10 AM |
| 4   | 0.6175 | 0.3925mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 6:08:16 AM |

Mean Area 0.5974  
Mean Conc. 0.3826mg/L



## Sample

Sample Name: 590-20532-b-5  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result

# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx


| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 1.000           | NPOC:1.602mg/L |

1. Det

Anal.: NPOC

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 3.037 | 1.588mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 6:16:11 AM |
| 2   | 2.995 | 1.567mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 6:18:16 AM |
| 3   | 3.019 | 1.579mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 6:20:21 AM |
| 4   | 3.222 | 1.679mg/L | 50uL      | 1.000     | E   | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 6:22:27 AM |
| 5   | 3.214 | 1.675mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 6:24:31 AM |

Mean Area 3.066  
Mean Conc. 1.602mg/L



## Sample

Sample Name: CCV  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 1.000           | NPOC:24.52mg/L |

1. Det

Anal.: NPOC

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 50.36 | 24.96mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 6:32:41 AM |
| 2   | 49.20 | 24.39mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 6:34:54 AM |
| 3   | 48.63 | 24.11mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 6:37:05 AM |
| 4   | 49.63 | 24.60mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 6:39:12 AM |

Mean Area 49.45  
Mean Conc. 24.52mg/L



# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx

## Sample

Sample Name: CCB  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result          |
|---------|-------|-----------------|-----------------|
| Unknown | NPOC  | 1.000           | NPOC:0.2184mg/L |

1. Det

Anal.: NPOC

| No. | Area   | Conc.      | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|--------|------------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 0.3185 | 0.2448mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 6:47:07 AM |
| 2   | 0.2386 | 0.2054mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 6:49:12 AM |
| 3   | 0.2079 | 0.1902mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 6:51:17 AM |
| 4   | 0.2946 | 0.2330mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 6:53:23 AM |

Mean Area 0.2649  
Mean Conc. 0.2184mg/L



## Sample

Sample Name: 590-20532-a-6  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result

| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 1.000           | NPOC:1.944mg/L |

1. Det

Anal.: NPOC


| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 3.842 | 1.985mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 7:01:18 AM |
| 2   | 3.779 | 1.954mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 7:03:23 AM |
| 3   | 3.581 | 1.856mg/L | 50uL      | 1.000     | E   | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 7:05:27 AM |
| 4   | 3.668 | 1.899mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 7:07:33 AM |
| 5   | 3.744 | 1.937mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 7:09:37 AM |



# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx

Mean Area 3.758  
Mean Conc. 1.944mg/L



## Sample

Sample Name: 280-176326-b-1  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result          |
|---------|-------|-----------------|-----------------|
| Unknown | NPOC  | 1.000           | NPOC:0.8675mg/L |

1. Det

Anal.: NPOC

| No. | Area  | Conc.      | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|------------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 1.686 | 0.9203mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 7:17:28 AM |
| 2   | 1.513 | 0.8349mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 7:19:32 AM |
| 3   | 1.715 | 0.9347mg/L | 50uL      | 1.000     | E   | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 7:21:37 AM |
| 4   | 1.538 | 0.8472mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 7:23:42 AM |
| 5   | 1.579 | 0.8675mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 7:25:48 AM |

Mean Area 1.579  
Mean Conc. 0.8675mg/L



## Sample

Sample Name: MS 280-176326-b-1  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result

| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 1.000           | NPOC:25.16mg/L |

1. Det


Anal.: NPOC

# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 51.88 | 25.71mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 7:33:50 AM |
| 2   | 50.49 | 25.03mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 7:36:04 AM |
| 3   | 50.05 | 24.81mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 7:38:09 AM |
| 4   | 50.59 | 25.08mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 7:40:14 AM |

Mean Area 50.75  
Mean Conc. 25.16mg/L



## Sample

Sample Name: MSD 280-176326-b-1  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 1.000           | NPOC:24.99mg/L |

1. Det

Anal.: NPOC

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 51.38 | 25.47mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 7:48:17 AM |
| 2   | 49.66 | 24.62mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 7:50:22 AM |
| 3   | 49.89 | 24.73mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 7:52:29 AM |
| 4   | 50.68 | 25.12mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 7:54:37 AM |

Mean Area 50.40  
Mean Conc. 24.99mg/L



## Sample

Sample Name: 280-176326-b-2  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result

# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx

| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 1.000           | NPOC:1.143mg/L |


1. Det

Anal.: NPOC

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 2.388 | 1.267mg/L | 50uL      | 1.000     | E   | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 8:02:32 AM |
| 2   | 2.085 | 1.117mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 8:04:49 AM |
| 3   | 2.146 | 1.148mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 8:06:58 AM |
| 4   | 2.204 | 1.176mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 8:09:07 AM |
| 5   | 2.110 | 1.130mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 8:11:17 AM |

Mean Area  
Mean Conc.

2.136  
1.143mg/L



## Sample

Sample Name: 280-176326-b-3  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result

| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 1.000           | NPOC:1.043mg/L |


1. Det

Anal.: NPOC

| No. | Area  | Conc.      | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|------------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 2.319 | 1.233mg/L  | 50uL      | 1.000     | E   | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 8:19:07 AM |
| 2   | 2.016 | 1.083mg/L  | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 8:21:16 AM |
| 3   | 1.832 | 0.9925mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 8:23:22 AM |
| 4   | 1.961 | 1.056mg/L  | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 8:25:30 AM |
| 5   | 1.931 | 1.041mg/L  | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 8:27:35 AM |

Mean Area  
Mean Conc.

1.935  
1.043mg/L



# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx

## Sample

Sample Name: 280-176326-b-4  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result          |
|---------|-------|-----------------|-----------------|
| Unknown | NPOC  | 1.000           | NPOC:0.5668mg/L |

1. Det

Anal.: NPOC

| No. | Area   | Conc.      | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|--------|------------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 1.179  | 0.6699mg/L | 50uL      | 1.000     | E   | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 8:35:26 AM |
| 2   | 0.9629 | 0.5632mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 8:37:31 AM |
| 3   | 1.159  | 0.6600mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 8:39:36 AM |
| 4   | 0.8843 | 0.5243mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 8:41:42 AM |
| 5   | 0.8751 | 0.5198mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 8:43:47 AM |

Mean Area 0.9703  
Mean Conc. 0.5668mg/L



## Sample

Sample Name: 280-176326-a-5  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result

| Type    | Anal. | Manual Dilution | Result          |
|---------|-------|-----------------|-----------------|
| Unknown | NPOC  | 1.000           | NPOC:0.8970mg/L |

1. Det


Anal.: NPOC

| No. | Area  | Conc.      | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|------------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 1.674 | 0.9144mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 8:51:38 AM |
| 2   | 1.646 | 0.9006mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 8:53:43 AM |
| 3   | 1.662 | 0.9085mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 8:55:49 AM |
| 4   | 1.573 | 0.8645mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 8:57:58 AM |

# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx

Mean Area 1.639  
Mean Conc. 0.8970mg/L



## Sample

Sample Name: 590-20532-b-7  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 1.000           | NPOC:1.996mg/L |

1. Det

Anal.: NPOC

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 4.218 | 2.171mg/L | 50uL      | 1.000     | E   | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 9:05:55 AM |
| 2   | 3.755 | 1.942mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 9:07:59 AM |
| 3   | 3.835 | 1.982mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 9:10:04 AM |
| 4   | 3.875 | 2.002mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 9:12:10 AM |
| 5   | 3.987 | 2.057mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 9:14:16 AM |

Mean Area 3.863  
Mean Conc. 1.996mg/L



## Sample

Sample Name: 590-20532-a-8  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result

| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 1.000           | NPOC:1.972mg/L |

1. Det


Anal.: NPOC

# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 3.829 | 1.979mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 9:22:06 AM |
| 2   | 3.880 | 2.004mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 9:24:11 AM |
| 3   | 3.688 | 1.909mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 9:26:17 AM |
| 4   | 3.863 | 1.996mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 9:28:23 AM |

Mean Area 3.815  
Mean Conc. 1.972mg/L



## Sample

Sample Name: CCB  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 1.000           | NPOC:24.33mg/L |

1. Det

Anal.: NPOC

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 49.74 | 24.66mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 9:36:35 AM |
| 2   | 48.60 | 24.09mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 9:38:53 AM |
| 3   | 48.54 | 24.07mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 9:40:58 AM |
| 4   | 49.39 | 24.48mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 9:43:09 AM |

Mean Area 49.07  
Mean Conc. 24.33mg/L



## Sample

Sample Name: CCB  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result

# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx


| Type    | Anal. | Manual Dilution | Result          |
|---------|-------|-----------------|-----------------|
| Unknown | NPOC  | 1.000           | NPOC:0.2680mg/L |

1. Det

Anal.: NPOC

| No. | Area   | Conc.      | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|--------|------------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 0.3777 | 0.2741mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 9:51:04 AM |
| 2   | 0.3175 | 0.2443mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 9:53:10 AM |
| 3   | 0.2903 | 0.2309mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 9:55:19 AM |
| 4   | 0.4760 | 0.3226mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 9:57:24 AM |

Mean Area 0.3654  
Mean Conc. 0.2680mg/L



## Sample

Sample Name: 590-20532-a-4  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result          |
|---------|-------|-----------------|-----------------|
| Unknown | NPOC  | 1.000           | NPOC:0.4453mg/L |

1. Det

Anal.: NPOC

| No. | Area   | Conc.      | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time           |
|-----|--------|------------|-----------|-----------|-----|-------------------------------------------------------------------|-----------------------|
| 1   | 0.7756 | 0.4706mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 10:05:20 AM |
| 2   | 0.6559 | 0.4115mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 10:07:26 AM |
| 3   | 0.7025 | 0.4345mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 10:09:31 AM |
| 4   | 0.7629 | 0.4644mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 10:11:35 AM |

Mean Area 0.7242  
Mean Conc. 0.4453mg/L



# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx

## Sample

Sample Name: 590-20532-a-5  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 1.000           | NPOC:1.620mg/L |

1. Det

Anal.: NPOC

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time           |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|-----------------------|
| 1   | 3.177 | 1.657mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 10:19:31 AM |
| 2   | 3.177 | 1.657mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 10:21:36 AM |
| 3   | 3.063 | 1.601mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 10:23:41 AM |
| 4   | 2.996 | 1.567mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 10:25:46 AM |

Mean Area 3.103  
Mean Conc. 1.620mg/L



## Sample

Sample Name: 590-20532-b-6  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result

| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 1.000           | NPOC:1.911mg/L |

1. Det

Anal.: NPOC


| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time           |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|-----------------------|
| 1   | 3.755 | 1.942mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 10:33:41 AM |
| 2   | 3.734 | 1.932mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 10:35:51 AM |
| 3   | 3.651 | 1.891mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 10:38:00 AM |
| 4   | 3.629 | 1.880mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 10:40:05 AM |



# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx

Mean Area 3.692  
Mean Conc. 1.911mg/L



## Sample

Sample Name: LCS  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 1.000           | NPOC:24.81mg/L |

1. Det

Anal.: NPOC

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time           |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|-----------------------|
| 1   | 50.23 | 24.90mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 10:48:01 AM |
| 2   | 49.93 | 24.75mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 10:50:11 AM |
| 3   | 49.61 | 24.59mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 10:52:16 AM |
| 4   | 50.46 | 25.01mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 10:54:21 AM |

Mean Area 50.06  
Mean Conc. 24.81mg/L



## Sample

Sample Name: MB  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result

| Type    | Anal. | Manual Dilution | Result          |
|---------|-------|-----------------|-----------------|
| Unknown | NPOC  | 1.000           | NPOC:0.2352mg/L |

1. Det


Anal.: NPOC

# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx

| No. | Area   | Conc.      | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time           |
|-----|--------|------------|-----------|-----------|-----|-------------------------------------------------------------------|-----------------------|
| 1   | 0.3045 | 0.2379mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 11:02:16 AM |
| 2   | 0.5415 | 0.3550mg/L | 50uL      | 1.000     | E   | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 11:04:22 AM |
| 3   | 0.2745 | 0.2231mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 11:06:27 AM |
| 4   | 0.2332 | 0.2027mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 11:08:32 AM |
| 5   | 0.3837 | 0.2770mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 11:10:37 AM |

Mean Area 0.2990  
Mean Conc. 0.2352mg/L



## Sample

Sample Name: 590-20532-a-7  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 1.000           | NPOC:2.123mg/L |

1. Det

Anal.: NPOC

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time           |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|-----------------------|
| 1   | 4.218 | 2.171mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 11:18:27 AM |
| 2   | 4.156 | 2.140mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 11:20:32 AM |
| 3   | 3.877 | 2.003mg/L | 50uL      | 1.000     | E   | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 11:22:37 AM |
| 4   | 3.954 | 2.041mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 11:24:42 AM |
| 5   | 4.155 | 2.140mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 11:26:47 AM |

Mean Area 4.121  
Mean Conc. 2.123mg/L



## Sample

# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx

Sample Name: MS 590-20532-a-7  
Sample ID: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Origin: Completed  
Status:  
Chk. Result:


| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 1.000           | NPOC:26.53mg/L |

1. Det

Anal.: NPOC

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time           |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|-----------------------|
| 1   | 54.61 | 27.06mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 11:34:49 AM |
| 2   | 53.18 | 26.36mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 11:37:01 AM |
| 3   | 52.65 | 26.10mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 11:39:15 AM |
| 4   | 53.69 | 26.61mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 11:41:23 AM |

Mean Area 53.53  
Mean Conc. 26.53mg/L



## Sample

Sample Name: MSD 590-20532-a-7  
Sample ID: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Origin: Completed  
Status:  
Chk. Result:

| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 1.000           | NPOC:26.40mg/L |

1. Det


Anal.: NPOC

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time           |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|-----------------------|
| 1   | 54.13 | 26.83mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 11:49:18 AM |
| 2   | 53.00 | 26.27mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 11:51:23 AM |
| 3   | 52.22 | 25.88mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 11:53:28 AM |
| 4   | 53.68 | 26.60mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 11:55:33 AM |

# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx

Mean Area 53.26  
Mean Conc. 26.40mg/L



## Sample

Sample Name: 590-20532-b-8  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 1.000           | NPOC:1.986mg/L |

1. Det

Anal.: NPOC

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time           |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|-----------------------|
| 1   | 3.961 | 2.044mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 12:03:29 PM |
| 2   | 3.774 | 1.952mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 12:05:34 PM |
| 3   | 3.732 | 1.931mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 12:07:39 PM |
| 4   | 3.903 | 2.015mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 12:09:44 PM |
| 5   | 4.036 | 2.081mg/L | 50uL      | 1.000     | E   | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 12:11:49 PM |

Mean Area 3.843  
Mean Conc. 1.986mg/L



## Sample

Sample Name: 590-20532-b-9  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result

| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 1.000           | NPOC:2.082mg/L |

1. Det


Anal.: NPOC

# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time           |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|-----------------------|
| 1   | 4.145 | 2.135mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 12:19:40 PM |
| 2   | 3.791 | 1.960mg/L | 50uL      | 1.000     | E   | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 12:21:45 PM |
| 3   | 3.952 | 2.040mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 12:23:50 PM |
| 4   | 3.940 | 2.034mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 12:25:55 PM |
| 5   | 4.114 | 2.120mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 12:28:10 PM |

Mean Area 4.038  
Mean Conc. 2.082mg/L



## Sample

Sample Name: CCV  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 1.000           | NPOC:24.87mg/L |

1. Det

Anal.: NPOC

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time           |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|-----------------------|
| 1   | 51.23 | 25.39mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 12:36:23 PM |
| 2   | 49.86 | 24.72mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 12:38:41 PM |
| 3   | 49.51 | 24.54mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 12:40:46 PM |
| 4   | 50.08 | 24.83mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 12:42:57 PM |

Mean Area 50.17  
Mean Conc. 24.87mg/L



## Sample

# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx

Sample Name: CCB  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result          |
|---------|-------|-----------------|-----------------|
| Unknown | NPOC  | 1.000           | NPOC:0.2183mg/L |

1. Det

Anal.: NPOC

| No. | Area   | Conc.      | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time           |
|-----|--------|------------|-----------|-----------|-----|-------------------------------------------------------------------|-----------------------|
| 1   | 0.2804 | 0.2260mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 12:50:53 PM |
| 2   | 0.2777 | 0.2247mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 12:52:58 PM |
| 3   | 0.2529 | 0.2124mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 12:55:03 PM |
| 4   | 0.2482 | 0.2101mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 12:57:08 PM |

Mean Area 0.2648  
Mean Conc. 0.2183mg/L



## Sample

Sample Name: 590-20532-b-10  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result

| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 1.000           | NPOC:1.906mg/L |

1. Det


Anal.: NPOC

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 3.773 | 1.951mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 1:05:04 PM |
| 2   | 3.641 | 1.886mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 1:07:09 PM |
| 3   | 3.609 | 1.870mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 1:09:14 PM |
| 4   | 3.702 | 1.916mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 1:11:19 PM |

# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx

Mean Area 3.681  
Mean Conc. 1.906mg/L



## Sample

Sample Name: 280-176081-j-3  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 10.00           | NPOC:11.35mg/L |

1. Det

Anal.: NPOC

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 2.219 | 11.84mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 1:19:15 PM |
| 2   | 2.084 | 11.17mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 1:21:20 PM |
| 3   | 2.065 | 11.08mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 1:23:25 PM |
| 4   | 2.111 | 11.30mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 1:25:30 PM |

Mean Area 2.120  
Mean Conc. 11.35mg/L



## Sample

Sample Name: 280-176081-j-4  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result

| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 2.000           | NPOC:8.942mg/L |

1. Det


Anal.: NPOC

# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 9.472 | 9.533mg/L | 50uL      | 1.000     | E   | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 1:33:26 PM |
| 2   | 8.916 | 8.984mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 1:35:50 PM |
| 3   | 8.838 | 8.907mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 1:37:55 PM |
| 4   | 8.702 | 8.772mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 1:40:00 PM |
| 5   | 9.038 | 9.104mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 1:42:05 PM |

Mean Area 8.873  
Mean Conc. 8.942mg/L



## Sample

Sample Name: 280-176081-j-5  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 2.000           | NPOC:6.013mg/L |

1. Det

Anal.: NPOC

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 6.048 | 6.150mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 1:49:56 PM |
| 2   | 5.829 | 5.934mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 1:52:01 PM |
| 3   | 5.730 | 5.836mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 1:54:06 PM |
| 4   | 6.030 | 6.132mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 1:56:36 PM |
| 5   | 6.271 | 6.370mg/L | 50uL      | 1.000     | E   | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 1:58:53 PM |

Mean Area 5.909  
Mean Conc. 6.013mg/L



## Sample



# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx

Sample Name: 280-176081-j-6  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 2.000           | NPOC:6.591mg/L |

1. Det

Anal.: NPOC

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 6.661 | 6.756mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 2:06:44 PM |
| 2   | 6.372 | 6.470mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 2:08:50 PM |
| 3   | 6.504 | 6.601mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 2:11:16 PM |
| 4   | 6.441 | 6.538mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 2:13:21 PM |

Mean Area 6.495  
Mean Conc. 6.591mg/L



## Sample

Sample Name: 280-176081-j-7  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result

| Type    | Anal. | Manual Dilution | Result          |
|---------|-------|-----------------|-----------------|
| Unknown | NPOC  | 2.000           | NPOC:0.9532mg/L |

1. Det


Anal.: NPOC

| No. | Area   | Conc.      | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|--------|------------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 0.7817 | 0.9473mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 2:21:17 PM |
| 2   | 0.6991 | 0.8657mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 2:23:22 PM |
| 3   | 0.8321 | 0.9971mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 2:25:27 PM |
| 4   | 0.8379 | 1.003mg/L  | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 2:27:32 PM |

# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx

Mean Area 0.7877  
Mean Conc. 0.9532mg/L



## Sample

Sample Name: 280-176081-j-8  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 2.000           | NPOC:14.64mg/L |

1. Det

Anal.: NPOC

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 14.81 | 14.81mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 2:35:28 PM |
| 2   | 14.49 | 14.49mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 2:37:33 PM |
| 3   | 14.53 | 14.53mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 2:39:38 PM |
| 4   | 14.72 | 14.72mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 2:41:43 PM |

Mean Area 14.64  
Mean Conc. 14.64mg/L



## Sample

Sample Name: 280-176081-j-9  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result

| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 2.000           | NPOC:5.626mg/L |

1. Det


Anal.: NPOC

# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 5.688 | 5.794mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 2:49:39 PM |
| 2   | 5.430 | 5.540mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 2:51:44 PM |
| 3   | 5.462 | 5.571mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 2:53:49 PM |
| 4   | 5.488 | 5.597mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 2:55:54 PM |
| 5   | 5.991 | 6.094mg/L | 50uL      | 1.000     | E   | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 2:57:59 PM |

Mean Area 5.517  
Mean Conc. 5.626mg/L



## Sample

Sample Name: MS 280-176081-j-9  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 2.000           | NPOC:54.65mg/L |

1. Det

Anal.: NPOC

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 56.48 | 55.97mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 3:06:13 PM |
| 2   | 54.68 | 54.20mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 3:08:22 PM |
| 3   | 54.38 | 53.90mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 3:10:29 PM |
| 4   | 55.03 | 54.54mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 3:12:34 PM |

Mean Area 55.14  
Mean Conc. 54.65mg/L



## Sample

# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx

Sample Name: MSD 280-176081-j-9  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 2.000           | NPOC:54.46mg/L |

1. Det

Anal.: NPOC

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 56.38 | 55.88mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 3:20:58 PM |
| 2   | 54.48 | 54.00mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 3:23:07 PM |
| 3   | 54.08 | 53.60mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 3:25:16 PM |
| 4   | 54.84 | 54.35mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 3:27:21 PM |

Mean Area 54.95  
Mean Conc. 54.46mg/L



## Sample

Sample Name: CCV  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result

| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 1.000           | NPOC:24.75mg/L |

1. Det


Anal.: NPOC

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 50.96 | 25.26mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 3:35:31 PM |
| 2   | 49.22 | 24.40mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 3:37:49 PM |
| 3   | 49.39 | 24.48mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 3:39:58 PM |
| 4   | 50.14 | 24.86mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 3:42:03 PM |

# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx

Mean Area 49.93  
Mean Conc. 24.75mg/L



## Sample

Sample Name: CCB  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result          |
|---------|-------|-----------------|-----------------|
| Unknown | NPOC  | 1.000           | NPOC:0.2674mg/L |

1. Det

Anal.: NPOC

| No. | Area   | Conc.      | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|--------|------------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 0.3353 | 0.2531mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 3:50:00 PM |
| 2   | 0.4765 | 0.3229mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 3:52:05 PM |
| 3   | 0.3211 | 0.2461mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 3:54:10 PM |
| 4   | 0.3236 | 0.2474mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 3:56:15 PM |

Mean Area 0.3641  
Mean Conc. 0.2674mg/L



## Sample

Sample Name: 280-176081-j-10  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result

| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 2.000           | NPOC:1.446mg/L |

1. Det


Anal.: NPOC

# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 1.318 | 1.477mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 4:04:12 PM |
| 2   | 1.258 | 1.418mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 4:06:17 PM |
| 3   | 1.295 | 1.454mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 4:08:22 PM |
| 4   | 1.275 | 1.435mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 4:10:27 PM |

Mean Area 1.287  
Mean Conc. 1.446mg/L



## Sample

Sample Name: 280-176081-j-11  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 2.000           | NPOC:6.315mg/L |

1. Det

Anal.: NPOC

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 6.667 | 6.762mg/L | 50uL      | 1.000     | E   | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 4:18:24 PM |
| 2   | 6.103 | 6.204mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 4:20:29 PM |
| 3   | 6.068 | 6.170mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 4:22:34 PM |
| 4   | 6.299 | 6.398mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 4:24:39 PM |
| 5   | 6.390 | 6.488mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 4:26:44 PM |

Mean Area 6.215  
Mean Conc. 6.315mg/L



## Sample

# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx

Sample Name: 280-176081-j-12  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 2.000           | NPOC:4.506mg/L |

1. Det

Anal.: NPOC

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 4.520 | 4.641mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 4:34:36 PM |
| 2   | 4.289 | 4.412mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 4:36:41 PM |
| 3   | 4.451 | 4.572mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 4:38:54 PM |
| 4   | 4.274 | 4.398mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 4:40:59 PM |
| 5   | 4.676 | 4.795mg/L | 50uL      | 1.000     | E   | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 4:43:04 PM |

Mean Area 4.384  
Mean Conc. 4.506mg/L



## Sample

Sample Name: CCV  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result

| Type    | Anal. | Manual Dilution | Result         |
|---------|-------|-----------------|----------------|
| Unknown | NPOC  | 1.000           | NPOC:24.82mg/L |

1. Det


Anal.: NPOC

| No. | Area  | Conc.     | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|-------|-----------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 51.03 | 25.30mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 4:51:00 PM |
| 2   | 49.62 | 24.60mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 4:53:11 PM |
| 3   | 49.55 | 24.56mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 4:55:35 PM |
| 4   | 50.08 | 24.83mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 4:57:54 PM |

# TOC-Control L Report

ABW  
2023\_05\_19\_001.tlx

Mean Area 50.07  
Mean Conc. 24.82mg/L



## Sample

Sample Name: CCB  
Sample ID:  
Origin: Calibration\_10\_31\_22\_Auto Dilution\_Cal.2023\_05\_04\_14\_32\_31.cal  
Status: Completed  
Chk. Result


| Type    | Anal. | Manual Dilution | Result          |
|---------|-------|-----------------|-----------------|
| Unknown | NPOC  | 1.000           | NPOC:0.2411mg/L |

1. Det

Anal.: NPOC

| No. | Area   | Conc.      | Inj. Vol. | Aut. Dil. | Ex. | Cal. Curve                                                        | Date / Time          |
|-----|--------|------------|-----------|-----------|-----|-------------------------------------------------------------------|----------------------|
| 1   | 0.4046 | 0.2874mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 5:05:51 PM |
| 2   | 0.2701 | 0.2209mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 5:08:00 PM |
| 3   | 0.2519 | 0.2119mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 5:10:05 PM |
| 4   | 0.3168 | 0.2440mg/L | 50uL      | 1.000     |     | Calibration_10_31_22_Auto<br>Dilution_Cal.2023_05_04_14_32_31.cal | 5/20/2023 5:12:10 PM |

Mean Area 0.3109  
Mean Conc. 0.2411mg/L





Shi 5 TOC

9060

|  | Analysis | Sample Name | Result             | Notes        | Vial |
|--|----------|-------------|--------------------|--------------|------|
|  | 1*       | NPOC        | ICV                |              | 1    |
|  | 2*       | NPOC        | ICB                |              | 2    |
|  | 3*       | NPOC        | LCS                |              | 3    |
|  | 4*       | NPOC        | MB                 |              | 4    |
|  | 5*       | NPOC        | TIC                |              | 5    |
|  | 6*       | NPOC        | 280-176674-d-7     |              | 6    |
|  | 7*       | NPOC        | MS 280-176674-d-7  |              | 7    |
|  | 8*       | NPOC        | MSD 280-176674-d-7 |              | 8    |
|  | 9*       | NPOC        | 580-127209-p-1     |              | 9    |
|  | 10*      | NPOC        | 580-127209-p-3     |              | 10   |
|  | 11*      | NPOC        | 580-127209-p-5     |              | 11   |
|  | 12*      | NPOC        | 280-176349-d-8     | 10x odor RR  | 12   |
|  | 13*      | NPOC        | 280-176349-c-11    | 25x odor RR  | 13   |
|  | 14*      | NPOC        | 280-176123-a-35    | 25x RR       | 14   |
|  | 15*      | NPOC        | CCV                |              | 15   |
|  | 16*      | NPOC        | CCB                |              | 16   |
|  | 17*      | NPOC        | 280-176363-c-1     | 2x color     | 17   |
|  | 18*      | NPOC        | 280-176363-a-2     |              | 18   |
|  | 19*      | NPOC        | 280-176363-c-3     | 10x 2x color | 19   |
|  | 20*      | NPOC        | 280-176363-d-5     | 10x color    | 20   |
|  | 21*      | NPOC        | MS 280-176363-d-5  | 10x          | 21   |
|  | 22*      | NPOC        | MSD 280-176363-d-5 | 10x ↓        | 22   |
|  | 23*      | NPOC        | 280-176363-d-6     |              | 23   |
|  | 24*      | NPOC        | 280-176363-c-7     | 2x color     | 24   |
|  | 25*      | NPOC        | 280-176363-d-8     | 10x color    | 25   |
|  | 26*      | NPOC        | 280-176325-d-1     | 500x color   | 26   |
|  | 27*      | NPOC        | CCV                |              | 27   |
|  | 28*      | NPOC        | CCB                |              | 28   |
|  | 29*      | NPOC        | 280-176325-d-2     | 500x color   | 29   |
|  | 30*      | NPOC        | 280-176274-d-2     | 100x odor    | 30   |
|  | 31*      | NPOC        | 280-176296-g-4     |              | 31   |
|  | 32*      | NPOC        | 590-20532-a-9      |              | 32   |
|  | 33*      | NPOC        | 590-20532-a-10     |              | 33   |
|  | 34*      | NPOC        | LCS                |              | 34   |
|  | 35*      | NPOC        | LCSD               |              | 35   |
|  | 36*      | NPOC        | MB                 |              | 36   |
|  | 37*      | NPOC        | TIC                |              | 37   |
|  | 38*      | NPOC        | 280-176296-g-3     |              | 38   |
|  | 39*      | NPOC        | CCV                |              | 39   |
|  | 40*      | NPOC        | CCB                |              | 40   |
|  | 41*      | NPOC        | MS 280-176296-g-3  |              | 41   |
|  | 42*      | NPOC        | MSD 280-176296-g-3 |              | 42   |
|  | 43*      | NPOC        | 280-176296-g-5     |              | 43   |
|  | 44*      | NPOC        | 280-176296-g-7     |              | 44   |
|  | 45*      | NPOC        | 280-176296-g-8     |              | 45   |
|  | 46*      | NPOC        | 280-176296-g-9     |              | 46   |
|  | 47*      | NPOC        | 590-20532-a-1      |              | 47   |
|  | 48*      | NPOC        | 590-20532-a-2      |              | 48   |
|  | 49*      | NPOC        | 590-20532-a-3      |              | 49   |
|  | 50*      | NPOC        | 590-20532-b-5      |              | 50   |
|  | 51*      | NPOC        | CCV                |              | 51   |
|  | 52*      | NPOC        | CCB                |              | 52   |
|  | 53*      | NPOC        | 590-20532-a-6      |              | 53   |

5310

5310


Shi 5 TOC

|     | Analysis | Sample Name        | Result | Notes | Vial |
|-----|----------|--------------------|--------|-------|------|
| 54* | NPOC     | 280-176326-b-1     | PH 6.2 |       | 54   |
| 55* | NPOC     | MS 280-176326-b-1  | /      |       | 55   |
| 56* | NPOC     | MSD 280-176326-b-1 | /      |       | 56   |
| 57* | NPOC     | 280-176326-b-2     | /      |       | 57   |
| 58* | NPOC     | 280-176326-b-3     | /      |       | 58   |
| 59* | NPOC     | 280-176326-b-4     | /      |       | 59   |
| 60* | NPOC     | 280-176326-a-5     | /      |       | 60   |
| 61* | NPOC     | 590-20532-b-7      | /      |       | 61   |
| 62* | NPOC     | 590-20532-a-8      | /      |       | 62   |
| 63* | NPOC     | CCV                | /      |       | 63   |
| 64* | NPOC     | CCB                | /      |       | 64   |
| 65* | NPOC     | 590-20532-a-4      | /      |       | 65   |
| 66* | NPOC     | 590-20532-a-5      | /      |       | 66   |
| 67* | NPOC     | 590-20532-b-6      | /      |       | 67   |
| 68* | NPOC     | LCS                | /      |       | 68   |
| 69* | NPOC     | MB                 | /      |       | 69   |
| 70* | NPOC     | 590-20532-a-7      | /      |       | 70   |
| 71* | NPOC     | MS 590-20532-a-7   | /      |       | 71   |
| 72* | NPOC     | MSD 590-20532-a-7  | /      |       | 72   |
| 73* | NPOC     | 590-20532-b-8      | /      |       | 73   |
| 74* | NPOC     | 590-20532-b-9      | /      |       | 74   |
| 75* | NPOC     | CCV                | /      |       | 75   |
| 76* | NPOC     | CCB                | /      |       | 76   |
| 77* | NPOC     | 590-20532-b-10     | /      |       | 77   |
| 78* | NPOC     | 280-176081-j-3     | /      |       | 78   |
| 79* | NPOC     | 280-176081-j-4     | /      |       | 79   |
| 80* | NPOC     | 280-176081-j-5     | /      |       | 80   |
| 81* | NPOC     | 280-176081-j-6     | /      |       | 81   |
| 82* | NPOC     | 280-176081-j-7     | /      |       | 82   |
| 83* | NPOC     | 280-176081-j-8     | /      |       | 83   |
| 84* | NPOC     | 280-176081-j-9     | /      |       | 84   |
| 85* | NPOC     | MS 280-176081-j-9  | /      |       | 85   |
| 86* | NPOC     | MSD 280-176081-j-9 | /      |       | 86   |
| 87* | NPOC     | CCV                | /      |       | 87   |
| 88* | NPOC     | CCB                | /      |       | 88   |
| 89* | NPOC     | 280-176081-j-10    | /      |       | 89   |
| 90* | NPOC     | 280-176081-j-11    | /      |       | 90   |
| 91* | NPOC     | 280-176081-j-12    | /      |       | 91   |
| 92* | NPOC     | CCV                | /      |       | 92   |
| 93* | NPOC     | CCB                | /      |       | 93   |

BWH 5000, 1000 MS.

# Calibration Report

**Calibration Record # 664**



## Calibration Settings

|                |              |               |                  |
|----------------|--------------|---------------|------------------|
| Calibration ID | PH           | Date          | 5/19/2023        |
| Channel        | 1            | Time          | 10:48 AM         |
| Probe Type     | pH           | Temperature   | 294.24 K 21.09 C |
| Probe ID       | PH ELECTRODE | Analysis Type | Single Line Fit  |

## Calibration Results

|           |         |           |                             |
|-----------|---------|-----------|-----------------------------|
| Slope     | -59.835 | CorrCoeff | 0.9999                      |
| Intercept | 2.124   | Equation: | $Y = (-59.835) X + (2.124)$ |

**Calibration Validity**      **True**      Operator

|                                | Result  | Minimum | Maximum |
|--------------------------------|---------|---------|---------|
| <b>Slope</b>                   | -59.835 | -65.00  | -53.00  |
| <b>Intercept</b>               | 2.124   | -100.00 | 100.00  |
| <b>Correlation Coefficient</b> | 0.9999  | 0.99    | 1.00    |

Note: "True" means the calibration was within the specified ranges  
"False" means the calibration was NOT within the specified ranges

| Calibration Data | Standard | Reading |
|------------------|----------|---------|
|                  | 2.00     | 300.41  |
|                  | 4.00     | 181.76  |
|                  | 7.00     | 2.07    |
|                  | 10.00    | -173.33 |
|                  | 12.00    | -300.29 |

Test America  
Water Analysis Report

AT4 Analyst:\_\_\_\_\_

| Run Number        |           |          |       | 1902      |       | Order Number |          | 20230519-2 |          |          |               |               |               |          |
|-------------------|-----------|----------|-------|-----------|-------|--------------|----------|------------|----------|----------|---------------|---------------|---------------|----------|
| SampleID          | RunDate   | RunTime  | Temp  | cond (uS) | pH    | palk-ppm     | talk-ppm | bcarb-ppm  | carb-ppm | hydr-ppm | (mL) @<br>8.3 | (mL) @<br>4.5 | (mL) @<br>4.2 | Conc (N) |
| RINSE             | 5/19/2023 | 12:32 PM | 20.77 | -1.00     | 5.27  | -1.00        | -1.00    | -1.00      | -1.00    | -1.00    | -1.00         | -1.00         | -1.00         | -1.00    |
| INITIAL CHECK     | 5/19/2023 | 12:40 PM | 20.49 | -1.00     | 10.51 | 136.58       | 180.71   | .00        | 88.28    | 92.44    | 1.37          | 1.81          | -1.00         | .02      |
| BUFFER 7          | 5/19/2023 | 12:43 PM | 20.11 | -1.00     | 7.03  | -1.00        | -1.00    | -1.00      | -1.00    | -1.00    | -1.00         | -1.00         | -1.00         | -1.00    |
| LCS               | 5/19/2023 | 12:51 PM | 20.15 | -1.00     | 10.61 | 141.01       | 182.57   | .00        | 83.12    | 99.45    | 1.41          | 1.83          | -1.00         | .02      |
| LCSD              | 5/19/2023 | 12:58 PM | 20.16 | -1.00     | 10.63 | 139.13       | 182.96   | .00        | 87.66    | 95.30    | 1.39          | 1.83          | -1.00         | .02      |
| MB                | 5/19/2023 | 1:04 PM  | 20.34 | -1.00     | 8.57  | 1.57         | 5.06     | 1.93       | 3.13     | .00      | .02           | .06           | .07           | .02      |
| 280-176527-a-3    | 5/19/2023 | 1:12 PM  | 20.40 | -1.00     | 7.15  | .00          | 456.58   | 456.58     | .00      | .00      | .00           | 4.57          | -1.00         | .02      |
| DU 280-176527-a-3 | 5/19/2023 | 1:21 PM  | 20.59 | -1.00     | 7.29  | .00          | 464.21   | 464.21     | .00      | .00      | .00           | 4.64          | -1.00         | .02      |
| 280-176527-a-11   | 5/19/2023 | 1:26 PM  | 20.74 | -1.00     | 6.56  | .00          | 6.22     | 6.22       | .00      | .00      | .00           | .07           | .07           | .02      |
| 280-176527-a-19   | 5/19/2023 | 1:33 PM  | 20.96 | -1.00     | 6.53  | .00          | 234.97   | 234.97     | .00      | .00      | .00           | 2.35          | -1.00         | .02      |
| 280-176527-a-28   | 5/19/2023 | 1:39 PM  | 20.96 | -1.00     | 6.70  | .00          | 125.13   | 125.13     | .00      | .00      | .00           | 1.25          | -1.00         | .02      |
| 280-176607-a-3    | 5/19/2023 | 1:47 PM  | 20.81 | -1.00     | 6.98  | .00          | 373.61   | 373.61     | .00      | .00      | .00           | 3.74          | -1.00         | .02      |
| 280-176572-c-4    | 5/19/2023 | 1:53 PM  | 20.76 | -1.00     | 7.03  | .00          | 104.66   | 104.66     | .00      | .00      | .00           | 1.05          | -1.00         | .02      |
| 280-176572-d-6    | 5/19/2023 | 2:00 PM  | 20.70 | -1.00     | 7.00  | .00          | 256.97   | 256.97     | .00      | .00      | .00           | 2.57          | -1.00         | .02      |
| 280-176572-c-7    | 5/19/2023 | 2:04 PM  | 20.63 | -1.00     | 4.79  | .00          | .68      | .68        | .00      | .00      | .00           | .04           | .08           | .02      |
| 280-176572-a-15   | 5/19/2023 | 2:11 PM  | 20.66 | -1.00     | 6.89  | .00          | 254.83   | 254.83     | .00      | .00      | .00           | 2.55          | -1.00         | .02      |
| CCV               | 5/19/2023 | 2:19 PM  | 20.78 | -1.00     | 10.57 | 145.36       | 181.23   | .00        | 71.76    | 109.48   | 1.45          | 1.81          | -1.00         | .02      |
| CCB               | 5/19/2023 | 2:25 PM  | 20.89 | -1.00     | 8.36  | .02          | 3.84     | 3.79       | .05      | .00      | .00           | .05           | .07           | .02      |
| 280-176607-a-12   | 5/19/2023 | 2:31 PM  | 20.98 | -1.00     | 6.20  | .00          | 129.23   | 129.23     | .00      | .00      | .00           | 1.29          | -1.00         | .02      |
| 280-176607-a-19   | 5/19/2023 | 2:37 PM  | 21.21 | -1.00     | 5.91  | .00          | .96      | .96        | .00      | .00      | .00           | .03           | .05           | .02      |
| 280-176628-a-2    | 5/19/2023 | 2:42 PM  | 21.14 | -1.00     | 6.17  | .00          | 61.33    | 61.33      | .00      | .00      | .00           | .61           | -1.00         | .02      |
| 280-176628-a-9    | 5/19/2023 | 2:47 PM  | 20.99 | -1.00     | 5.85  | .00          | 4.46     | 4.46       | .00      | .00      | .00           | .05           | .07           | .02      |
| 280-176628-a-16   | 5/19/2023 | 2:52 PM  | 20.84 | -1.00     | 6.18  | .00          | 66.41    | 66.41      | .00      | .00      | .00           | .66           | -1.00         | .02      |
| 280-176642-a-5    | 5/19/2023 | 2:59 PM  | 20.81 | -1.00     | 6.50  | .00          | 287.41   | 287.41     | .00      | .00      | .00           | 2.87          | -1.00         | .02      |
| 280-176642-a-12   | 5/19/2023 | 3:06 PM  | 20.89 | -1.00     | 6.62  | .00          | 260.58   | 260.58     | .00      | .00      | .00           | 2.61          | -1.00         | .02      |
| 280-176696-a-2    | 5/19/2023 | 3:15 PM  | 20.91 | -1.00     | 7.49  | .00          | 455.40   | 455.40     | .00      | .00      | .00           | 4.55          | -1.00         | .02      |
| 280-176696-a-10   | 5/19/2023 | 3:26 PM  | 21.03 | -1.00     | 7.39  | .00          | 890.86   | 890.86     | .00      | .00      | .00           | 8.91          | -1.00         | .02      |

| Run Number        |           |         | 1902  |           | Order Number |          | 20230519-2 |           |          |          |               |               |               |          |
|-------------------|-----------|---------|-------|-----------|--------------|----------|------------|-----------|----------|----------|---------------|---------------|---------------|----------|
| SampleID          | RunDate   | RunTime | Temp  | cond (uS) | pH           | palk-ppm | talk-ppm   | bcarb-ppm | carb-ppm | hydr-ppm | (mL) @<br>8.3 | (mL) @<br>4.5 | (mL) @<br>4.2 | Conc (N) |
| 280-176696-a-17   | 5/19/2023 | 3:35 PM | 21.21 | -1.00     | 7.62         | .00      | 434.13     | 434.13    | .00      | .00      | .00           | 4.34          | -1.00         | .02      |
| CCV               | 5/19/2023 | 3:42 PM | 21.38 | -1.00     | 10.54        | 150.49   | 182.39     | .00       | 63.80    | 118.58   | 1.50          | 1.82          | -1.00         | .02      |
| CCB               | 5/19/2023 | 3:49 PM | 21.60 | -1.00     | 8.28         | .00      | 3.43       | 3.43      | .00      | .00      | .00           | .05           | .07           | .02      |
| LCS               | 5/19/2023 | 3:57 PM | 21.51 | -1.00     | 10.56        | 141.23   | 181.27     | .00       | 80.09    | 101.19   | 1.41          | 1.81          | -1.00         | .02      |
| LCSD              | 5/19/2023 | 4:04 PM | 21.32 | -1.00     | 10.60        | 140.38   | 183.06     | .00       | 85.36    | 97.70    | 1.40          | 1.83          | -1.00         | .02      |
| MB                | 5/19/2023 | 4:11 PM | 21.21 | -1.00     | 8.33         | .00      | 3.58       | 3.58      | .00      | .00      | .00           | .05           | .07           | .02      |
| 280-176698-a-2    | 5/19/2023 | 4:22 PM | 21.12 | -1.00     | 7.37         | .00      | 933.54     | 933.54    | .00      | .00      | .00           | 9.34          | -1.00         | .02      |
| DU 280-176698-a-2 | 5/19/2023 | 4:34 PM | 21.21 | -1.00     | 7.35         | .00      | 935.40     | 935.40    | .00      | .00      | .00           | 9.35          | -1.00         | .02      |
| 280-176698-a-9    | 5/19/2023 | 4:40 PM | 21.24 | -1.00     | 6.77         | .00      | 2.08       | 2.08      | .00      | .00      | .00           | .04           | .06           | .02      |
| 280-176705-a-2    | 5/19/2023 | 4:48 PM | 21.19 | -1.00     | 6.99         | .00      | 336.19     | 336.19    | .00      | .00      | .00           | 3.36          | -1.00         | .02      |
| 280-176591-c-2    | 5/19/2023 | 4:55 PM | 21.24 | -1.00     | 7.40         | .00      | 271.98     | 271.98    | .00      | .00      | .00           | 2.72          | -1.00         | .02      |
| 280-176591-c-4    | 5/19/2023 | 5:01 PM | 21.38 | -1.00     | 6.76         | .00      | 176.44     | 176.44    | .00      | .00      | .00           | 1.76          | -1.00         | .02      |
| 280-176591-c-6    | 5/19/2023 | 5:06 PM | 21.51 | -1.00     | 6.49         | .00      | 48.52      | 48.52     | .00      | .00      | .00           | .49           | -1.00         | .02      |
| 280-176591-c-7    | 5/19/2023 | 5:11 PM | 21.58 | -1.00     | 5.94         | .00      | 4.71       | 4.71      | .00      | .00      | .00           | .06           | .07           | .02      |
| 280-176674-a-7    | 5/19/2023 | 5:16 PM | 21.38 | -1.00     | 6.01         | .00      | 17.97      | 17.97     | .00      | .00      | .00           | .20           | .21           | .02      |
| 280-176674-a-8    | 5/19/2023 | 5:22 PM | 21.27 | -1.00     | 6.02         | .00      | 17.77      | 17.77     | .00      | .00      | .00           | .19           | .21           | .02      |
| CCV               | 5/19/2023 | 5:30 PM | 21.21 | -1.00     | 10.51        | 140.19   | 182.03     | .00       | 83.67    | 98.36    | 1.40          | 1.82          | -1.00         | .02      |
| CCB               | 5/19/2023 | 5:36 PM | 21.14 | -1.00     | 8.34         | .00      | 3.58       | 3.58      | .00      | .00      | .00           | .05           | .07           | .02      |
| 280-176514-g-1    | 5/19/2023 | 5:43 PM | 21.09 | -1.00     | 7.26         | .00      | 257.23     | 257.23    | .00      | .00      | .00           | 2.57          | -1.00         | .02      |
| 280-176514-g-2    | 5/19/2023 | 5:50 PM | 21.12 | -1.00     | 7.23         | .00      | 181.33     | 181.33    | .00      | .00      | .00           | 1.81          | -1.00         | .02      |
| 280-176699-b-1    | 5/19/2023 | 5:55 PM | 21.15 | -1.00     | 5.93         | .00      | 5.65       | 5.65      | .00      | .00      | .00           | .07           | .09           | .02      |
| 280-176699-b-2    | 5/19/2023 | 6:01 PM | 21.18 | -1.00     | 6.01         | .00      | 79.85      | 79.85     | .00      | .00      | .00           | .80           | -1.00         | .02      |
| 280-176699-b-3    | 5/19/2023 | 6:14 PM | 21.32 | -1.00     | 7.28         | .00      | 1,161.09   | 1,161.09  | .00      | .00      | .00           | 11.61         | -1.00         | .02      |
| 280-176699-b-4    | 5/19/2023 | 6:21 PM | 21.61 | -1.00     | 8.12         | .00      | 192.21     | 192.21    | .00      | .00      | .00           | 1.92          | -1.00         | .02      |
| 280-176699-b-5    | 5/19/2023 | 6:28 PM | 21.69 | -1.00     | 8.63         | 37.92    | 142.54     | 66.70     | 75.84    | .00      | .38           | 1.43          | -1.00         | .02      |
| 280-176363-c-10   | 5/19/2023 | 6:38 PM | 21.38 | -1.00     | 6.87         | .00      | 802.17     | 802.17    | .00      | .00      | .00           | 8.02          | -1.00         | .02      |
| 280-176363-b-19   | 5/19/2023 | 6:50 PM | 21.25 | -1.00     | 6.92         | .00      | 905.30     | 905.30    | .00      | .00      | .00           | 9.05          | -1.00         | .02      |
| 280-176363-b-23   | 5/19/2023 | 7:01 PM | 21.21 | -1.00     | 7.36         | .00      | 801.32     | 801.32    | .00      | .00      | .00           | 8.01          | -1.00         | .02      |
| CCV               | 5/19/2023 | 7:09 PM | 21.15 | -1.00     | 10.44        | 142.53   | 184.29     | .00       | 83.51    | 100.78   | 1.43          | 1.84          | -1.00         | .02      |
| CCB               | 5/19/2023 | 7:14 PM | 21.08 | -1.00     | 8.00         | .00      | 3.97       | 3.97      | .00      | .00      | .00           | .05           | .07           | .02      |

| <u>Run Number</u> |                |                |             | 1902             |           | <u>Order Number</u> |                 |                  |                 |                 | 20230519-2            |                       |                       |                 |
|-------------------|----------------|----------------|-------------|------------------|-----------|---------------------|-----------------|------------------|-----------------|-----------------|-----------------------|-----------------------|-----------------------|-----------------|
| <u>SampleID</u>   | <u>RunDate</u> | <u>RunTime</u> | <u>Temp</u> | <u>cond (uS)</u> | <u>pH</u> | <u>palk-ppm</u>     | <u>talk-ppm</u> | <u>bcarb-ppm</u> | <u>carb-ppm</u> | <u>hydr-ppm</u> | <u>(mL) @<br/>8.3</u> | <u>(mL) @<br/>4.5</u> | <u>(mL) @<br/>4.2</u> | <u>Conc (N)</u> |
| LCS               | 5/19/2023      | 7:22 PM        | 21.02       | -1.00            | 10.48     | 138.48              | 182.98          | .00              | 89.00           | 93.98           | 1.38                  | 1.83                  | -1.00                 | .02             |
| MB                | 5/19/2023      | 7:28 PM        | 21.06       | -1.00            | 8.13      | .00                 | 3.06            | 3.06             | .00             | .00             | .00                   | .05                   | .07                   | .02             |
| 280-176478-c-1    | 5/19/2023      | 7:40 PM        | 21.11       | -1.00            | 7.47      | .00                 | 834.79          | 834.79           | .00             | .00             | .00                   | 8.35                  | -1.00                 | .02             |
| DU 280-176478-c-1 | 5/19/2023      | 7:51 PM        | 21.27       | -1.00            | 7.60      | .00                 | 846.30          | 846.30           | .00             | .00             | .00                   | 8.46                  | -1.00                 | .02             |
| 280-176590-e-1    | 5/19/2023      | 8:05 PM        | 21.43       | -1.00            | 8.53      | 51.83               | 1,041.95        | 938.28           | 103.66          | .00             | .52                   | 10.42                 | -1.00                 | .02             |
| 280-176590-e-2    | 5/19/2023      | 8:21 PM        | 21.76       | -1.00            | 8.43      | 48.98               | 1,389.68        | 1,291.73         | 97.96           | .00             | .49                   | 13.90                 | -1.00                 | .02             |
| 280-176590-f-3    | 5/19/2023      | 8:34 PM        | 21.53       | -1.00            | 8.47      | 44.06               | 970.28          | 882.16           | 88.12           | .00             | .44                   | 9.70                  | -1.00                 | .02             |
| 280-176590-f-4    | 5/19/2023      | 8:49 PM        | 21.34       | -1.00            | 8.83      | 105.39              | 1,252.70        | 1,041.92         | 210.78          | .00             | 1.05                  | 12.53                 | -1.00                 | .02             |
| 280-176593-a-1    | 5/19/2023      | 8:55 PM        | 21.25       | -1.00            | 6.70      | .00                 | 135.22          | 135.22           | .00             | .00             | .00                   | 1.35                  | -1.00                 | .02             |
| 280-176596-a-2    | 5/19/2023      | 9:00 PM        | 21.11       | -1.00            | 6.68      | .00                 | 57.26           | 57.26            | .00             | .00             | .00                   | .57                   | -1.00                 | .02             |
| 280-176601-a-1    | 5/19/2023      | 9:07 PM        | 21.02       | -1.00            | 7.29      | .00                 | 158.20          | 158.20           | .00             | .00             | .00                   | 1.58                  | -1.00                 | .02             |
| 280-176606-a-1    | 5/19/2023      | 9:15 PM        | 20.99       | -1.00            | 7.67      | .00                 | 347.01          | 347.01           | .00             | .00             | .00                   | 3.47                  | -1.00                 | .02             |
| CCV               | 5/19/2023      | 9:23 PM        | 21.02       | -1.00            | 10.41     | 146.52              | 185.93          | .00              | 78.83           | 107.10          | 1.47                  | 1.86                  | -1.00                 | .02             |
| CCB               | 5/19/2023      | 9:29 PM        | 21.03       | -1.00            | 8.17      | .00                 | 3.45            | 3.45             | .00             | .00             | .00                   | .05                   | .07                   | .02             |
| 280-176606-a-2    | 5/19/2023      | 9:36 PM        | 21.08       | -1.00            | 7.68      | .00                 | 353.68          | 353.68           | .00             | .00             | .00                   | 3.54                  | -1.00                 | .02             |
| 280-176606-a-3    | 5/19/2023      | 9:45 PM        | 21.22       | -1.00            | 7.92      | .00                 | 449.21          | 449.21           | .00             | .00             | .00                   | 4.49                  | -1.00                 | .02             |
| 280-176606-a-4    | 5/19/2023      | 9:53 PM        | 21.41       | -1.00            | 7.76      | .00                 | 359.19          | 359.19           | .00             | .00             | .00                   | 3.59                  | -1.00                 | .02             |
| 280-176606-a-5    | 5/19/2023      | 10:02 PM       | 21.58       | -1.00            | 8.05      | .00                 | 530.27          | 530.27           | .00             | .00             | .00                   | 5.30                  | -1.00                 | .02             |
| 280-176606-a-6    | 5/19/2023      | 10:13 PM       | 21.38       | -1.00            | 7.95      | .00                 | 708.56          | 708.56           | .00             | .00             | .00                   | 7.09                  | -1.00                 | .02             |
| 280-176606-a-7    | 5/19/2023      | 10:18 PM       | 21.25       | -1.00            | 6.92      | .00                 | 6.12            | 6.12             | .00             | .00             | .00                   | .07                   | .08                   | .02             |
| 280-176606-a-8    | 5/19/2023      | 10:25 PM       | 21.11       | -1.00            | 7.38      | .00                 | 242.17          | 242.17           | .00             | .00             | .00                   | 2.42                  | -1.00                 | .02             |
| 280-176634-a-1    | 5/19/2023      | 10:32 PM       | 21.02       | -1.00            | 7.60      | .00                 | 97.78           | 97.78            | .00             | .00             | .00                   | .98                   | -1.00                 | .02             |
| 280-176634-a-2    | 5/19/2023      | 10:52 PM       | 20.99       | -1.00            | 7.56      | .00                 | .00             | .00              | .00             | .00             | .00                   | .00                   | .00                   | .02             |
| 280-176634-b-3    | 5/19/2023      | 11:12 PM       | 21.05       | -1.00            | 7.65      | .00                 | .00             | .00              | .00             | .00             | .00                   | .00                   | .00                   | .02             |
| CCV               | 5/19/2023      | 11:21 PM       | 21.02       | -1.00            | 10.33     | 162.94              | 193.27          | .00              | 60.66           | 132.61          | 1.63                  | 1.93                  | -1.00                 | .02             |
| CCB               | 5/19/2023      | 11:27 PM       | 20.96       | -1.00            | 7.74      | .00                 | 3.18            | 3.18             | .00             | .00             | .00                   | .06                   | .08                   | .02             |
| LCS               | 5/19/2023      | 11:35 PM       | 20.98       | -1.00            | 10.37     | 142.27              | 185.64          | .00              | 86.74           | 98.90           | 1.42                  | 1.86                  | -1.00                 | .02             |
| MB                | 5/19/2023      | 11:41 PM       | 21.12       | -1.00            | 8.03      | .00                 | 2.82            | 2.82             | .00             | .00             | .00                   | .05                   | .07                   | .02             |
| 280-176658-c-1    | 5/19/2023      | 11:47 PM       | 21.47       | -1.00            | 6.93      | .00                 | 118.76          | 118.76           | .00             | .00             | .00                   | 1.19                  | -1.00                 | .02             |
| DU 280-176658-c-1 | 5/19/2023      | 11:53 PM       | 21.25       | -1.00            | 7.09      | .00                 | 119.55          | 119.55           | .00             | .00             | .00                   | 1.20                  | -1.00                 | .02             |

Run Number

1902

Order Number

20230519-2

| <u>SampleID</u>   | <u>RunDate</u> | <u>RunTime</u> | <u>Temp</u> | <u>cond (uS)</u> | <u>pH</u> | <u>palk-ppm</u> | <u>talk-ppm</u> | <u>bcarb-ppm</u> | <u>carb-ppm</u> | <u>hydr-ppm</u> | <u>(mL) @<br/>8.3</u> | <u>(mL) @<br/>4.5</u> | <u>(mL) @<br/>4.2</u> | <u>Conc (N)</u> |
|-------------------|----------------|----------------|-------------|------------------|-----------|-----------------|-----------------|------------------|-----------------|-----------------|-----------------------|-----------------------|-----------------------|-----------------|
| 280-176658-b-2    | 5/19/2023      | 11:59 PM       | 21.09       | -1.00            | 7.12      | .00             | 122.09          | 122.09           | .00             | .00             | .00                   | 1.22                  | -1.00                 | .02             |
| 280-176658-b-3    | 5/20/2023      | 12:05 AM       | 20.98       | -1.00            | 7.13      | .00             | 122.93          | 122.93           | .00             | .00             | .00                   | 1.23                  | -1.00                 | .02             |
| 280-176658-b-4    | 5/20/2023      | 12:11 AM       | 20.88       | -1.00            | 7.12      | .00             | 122.57          | 122.57           | .00             | .00             | .00                   | 1.23                  | -1.00                 | .02             |
| 280-176658-b-5    | 5/20/2023      | 12:17 AM       | 20.85       | -1.00            | 7.16      | .00             | 122.91          | 122.91           | .00             | .00             | .00                   | 1.23                  | -1.00                 | .02             |
| 280-176658-b-6    | 5/20/2023      | 12:23 AM       | 20.81       | -1.00            | 7.15      | .00             | 125.25          | 125.25           | .00             | .00             | .00                   | 1.25                  | -1.00                 | .02             |
| 280-176658-b-7    | 5/20/2023      | 12:29 AM       | 20.77       | -1.00            | 7.20      | .00             | 128.46          | 128.46           | .00             | .00             | .00                   | 1.28                  | -1.00                 | .02             |
| 280-176634-b-4    | 5/20/2023      | 12:37 AM       | 20.80       | -1.00            | 7.87      | .00             | 314.73          | 314.73           | .00             | .00             | .00                   | 3.15                  | -1.00                 | .02             |
| 280-176634-a-5    | 5/20/2023      | 12:43 AM       | 20.84       | -1.00            | 7.12      | .00             | 73.82           | 73.82            | .00             | .00             | .00                   | .74                   | -1.00                 | .02             |
| CCV               | 5/20/2023      | 12:51 AM       | 20.92       | -1.00            | 10.33     | 138.77          | 186.94          | .00              | 96.34           | 90.60           | 1.39                  | 1.87                  | -1.00                 | .02             |
| CCB               | 5/20/2023      | 12:57 AM       | 21.02       | -1.00            | 8.08      | .00             | 4.30            | 4.30             | .00             | .00             | .00                   | .06                   | .07                   | .02             |
| 280-176634-a-6    | 5/20/2023      | 1:12 AM        | 21.14       | -1.00            | 9.92      | .00             | .00             | .00              | .00             | .00             | .00                   | .00                   | .00                   | .02             |
| 280-176634-b-7    | 5/20/2023      | 1:20 AM        | 20.94       | -1.00            | 8.64      | 35.16           | 142.29          | 71.98            | 70.31           | .00             | .35                   | 1.42                  | -1.00                 | .02             |
| 280-176634-b-8    | 5/20/2023      | 1:41 AM        | 20.74       | -1.00            | 7.66      | .00             | .00             | .00              | .00             | .00             | .00                   | .00                   | .00                   | .02             |
| 280-176634-b-9    | 5/20/2023      | 2:03 AM        | 20.58       | -1.00            | 7.71      | .00             | .00             | .00              | .00             | .00             | .00                   | .00                   | .00                   | .02             |
| 280-176634-b-10   | 5/20/2023      | 2:10 AM        | 20.44       | -1.00            | 8.04      | .00             | 129.94          | 129.94           | .00             | .00             | .00                   | 1.30                  | -1.00                 | .02             |
| 280-176634-b-11   | 5/20/2023      | 2:16 AM        | 20.33       | -1.00            | 7.71      | .00             | 84.59           | 84.59            | .00             | .00             | .00                   | .85                   | -1.00                 | .02             |
| 280-176634-b-12   | 5/20/2023      | 2:24 AM        | 20.27       | -1.00            | 7.85      | .00             | 316.05          | 316.05           | .00             | .00             | .00                   | 3.16                  | -1.00                 | .02             |
| 280-176669-a-1    | 5/20/2023      | 2:34 AM        | 20.27       | -1.00            | 8.21      | .00             | 631.75          | 631.75           | .00             | .00             | .00                   | 6.32                  | -1.00                 | .02             |
| 280-176669-a-2    | 5/20/2023      | 2:42 AM        | 20.30       | -1.00            | 8.60      | 52.79           | 413.71          | 308.13           | 105.58          | .00             | .53                   | 4.14                  | -1.00                 | .02             |
| 280-176669-a-3    | 5/20/2023      | 2:51 AM        | 20.31       | -1.00            | 8.10      | .00             | 396.08          | 396.08           | .00             | .00             | .00                   | 3.96                  | -1.00                 | .02             |
| CCV               | 5/20/2023      | 2:59 AM        | 20.40       | -1.00            | 10.24     | 153.19          | 190.80          | .00              | 75.22           | 115.58          | 1.53                  | 1.91                  | -1.00                 | .02             |
| CCB               | 5/20/2023      | 3:05 AM        | 20.77       | -1.00            | 8.02      | .00             | 4.67            | 4.67             | .00             | .00             | .00                   | .07                   | .09                   | .02             |
| LCS               | 5/20/2023      | 3:12 AM        | 20.65       | -1.00            | 10.28     | 143.49          | 187.96          | .00              | 88.94           | 99.02           | 1.43                  | 1.88                  | -1.00                 | .02             |
| MB                | 5/20/2023      | 3:18 AM        | 20.52       | -1.00            | 8.06      | .00             | 3.93            | 3.93             | .00             | .00             | .00                   | .06                   | .07                   | .02             |
| 280-176669-a-4    | 5/20/2023      | 3:26 AM        | 20.45       | -1.00            | 8.03      | .00             | 280.02          | 280.02           | .00             | .00             | .00                   | 2.80                  | -1.00                 | .02             |
| DU 280-176669-a-4 | 5/20/2023      | 3:33 AM        | 20.40       | -1.00            | 8.17      | .00             | 284.51          | 284.51           | .00             | .00             | .00                   | 2.85                  | -1.00                 | .02             |
| 280-176670-a-1    | 5/20/2023      | 3:39 AM        | 20.38       | -1.00            | 7.36      | .00             | 122.30          | 122.30           | .00             | .00             | .00                   | 1.22                  | -1.00                 | .02             |
| 280-176670-a-2    | 5/20/2023      | 3:45 AM        | 20.37       | -1.00            | 7.27      | .00             | 127.39          | 127.39           | .00             | .00             | .00                   | 1.27                  | -1.00                 | .02             |
| 280-176670-a-3    | 5/20/2023      | 3:51 AM        | 20.36       | -1.00            | 7.29      | .00             | 127.82          | 127.82           | .00             | .00             | .00                   | 1.28                  | -1.00                 | .02             |
| 280-176670-a-4    | 5/20/2023      | 3:57 AM        | 20.37       | -1.00            | 7.18      | .00             | 112.54          | 112.54           | .00             | .00             | .00                   | 1.13                  | -1.00                 | .02             |

| Run Number     |           |         |       | 1902      |       | Order Number |          |           |          |          | 20230519-2    |               |               |          |
|----------------|-----------|---------|-------|-----------|-------|--------------|----------|-----------|----------|----------|---------------|---------------|---------------|----------|
| SampleID       | RunDate   | RunTime | Temp  | cond (uS) | pH    | palk-ppm     | talk-ppm | bcarb-ppm | carb-ppm | hydr-ppm | (mL) @<br>8.3 | (mL) @<br>4.5 | (mL) @<br>4.2 | Conc (N) |
| 280-176670-a-6 | 5/20/2023 | 4:03 AM | 20.41 | -1.00     | 6.56  | .00          | 80.86    | 80.86     | .00      | .00      | .00           | .81           | -1.00         | .02      |
| 280-176670-a-5 | 5/20/2023 | 4:09 AM | 20.45 | -1.00     | 7.11  | .00          | 129.46   | 129.46    | .00      | .00      | .00           | 1.29          | -1.00         | .02      |
| 280-176670-a-8 | 5/20/2023 | 4:15 AM | 20.54 | -1.00     | 7.53  | .00          | 131.08   | 131.08    | .00      | .00      | .00           | 1.31          | -1.00         | .02      |
| 280-176722-e-1 | 5/20/2023 | 4:22 AM | 20.77 | -1.00     | 7.93  | .00          | 346.50   | 346.50    | .00      | .00      | .00           | 3.46          | -1.00         | .02      |
| CCV            | 5/20/2023 | 4:30 AM | 20.67 | -1.00     | 10.22 | 135.42       | 190.65   | .00       | 110.46   | 80.19    | 1.35          | 1.91          | -1.00         | .02      |
| CCB            | 5/20/2023 | 4:36 AM | 20.59 | -1.00     | 7.97  | .00          | 4.22     | 4.22      | .00      | .00      | .00           | .06           | .08           | .02      |
| 280-176722-e-2 | 5/20/2023 | 4:44 AM | 20.51 | -1.00     | 7.70  | .00          | 337.46   | 337.46    | .00      | .00      | .00           | 3.37          | -1.00         | .02      |
| 280-176683-a-1 | 5/20/2023 | 4:51 AM | 20.48 | -1.00     | 7.71  | .00          | 246.60   | 246.60    | .00      | .00      | .00           | 2.47          | -1.00         | .02      |
| 280-176683-a-2 | 5/20/2023 | 4:58 AM | 20.48 | -1.00     | 7.70  | .00          | 248.05   | 248.05    | .00      | .00      | .00           | 2.48          | -1.00         | .02      |
| 280-176683-a-3 | 5/20/2023 | 5:04 AM | 20.47 | -1.00     | 6.66  | .00          | 4.46     | 4.46      | .00      | .00      | .00           | .07           | .09           | .02      |
| 280-176439-c-2 | 5/20/2023 | 5:15 AM | 20.47 | -1.00     | 9.16  | 177.06       | 798.65   | 444.54    | 354.12   | .00      | 1.77          | 7.99          | -1.00         | .02      |
| 280-176526-j-1 | 5/20/2023 | 5:35 AM | 20.47 | -1.00     | 7.97  | .00          | 2,487.05 | 2,487.05  | .00      | .00      | .00           | 24.87         | -1.00         | .02      |
| 280-176663-m-1 | 5/20/2023 | 5:56 AM | 21.12 | -1.00     | 8.46  | 204.56       | .00      | -409.12   | 409.12   | .00      | 2.05          | .00           | .00           | .02      |
| 280-176663-l-2 | 5/20/2023 | 6:15 AM | 21.51 | -1.00     | 8.29  | .67          | .00      | -1.35     | 1.35     | .00      | .01           | .00           | .00           | .02      |
| 280-176663-k-3 | 5/20/2023 | 6:35 AM | 21.87 | -1.00     | 8.43  | 163.55       | .00      | -327.09   | 327.09   | .00      | 1.64          | .00           | .00           | .02      |
| 280-176663-k-4 | 5/20/2023 | 6:55 AM | 21.90 | -1.00     | 8.61  | 244.19       | .00      | -488.38   | 488.38   | .00      | 2.44          | .00           | .00           | .02      |
| RINSE          | 5/20/2023 | 7:02 AM | 21.90 | -1.00     | 8.31  | .00          | 10.02    | 10.02     | .00      | .00      | .00           | .13           | .16           | .02      |
| CCV            | 5/20/2023 | 7:10 AM | 21.82 | -1.00     | 10.16 | 139.25       | 190.81   | .00       | 103.13   | 87.68    | 1.39          | 1.91          | -1.00         | .02      |
| CCB            | 5/20/2023 | 7:16 AM | 21.77 | -1.00     | 7.96  | .00          | 3.35     | 3.35      | .00      | .00      | .00           | .06           | .08           | .02      |



# Shipping and Receiving Documents

RVAAP-337-TA

COC No.: RVAAP-


Page 1 of 1

[illegible]

0.1, 3.2, 0.3°C IR14CFO.1



280-176674 Chain of Custody



COC No.: RVAAP-558-TA

Page ~~1~~ of 7 Date: 05/17/2023

[illegible]

**Chain of Custody Record - 2 of 7**

COC No.: RVAAP-57-TA

Date: 5/17/23

[illegible]

**Chain of Custody Record** 3067

COC No.: **RVAAP-S90TA**

Date: 5/17/23

Page of 15

[illegible]



Page 1 of 1

5/31/2023 9:08 PM



## Chain of Custody Record

COC No.: **RVAAP-545-TA**

Date: 5/17/2023

Name Leidos  
Address: 8866 Commons Blvd, Suite 201, Twinsburg, OH 44087  
Phone Number: (330) 405-5802  
Project Manager: Jed Thomas  
Project: RVAAP FWGW Sampling Event Spring 2023  
Job/P.O. No.: P010216426  
Sampler (Signature) *Michelle Rowe* (Printed Name) M

Sampler (Signature) *Melissa Taylor* (Printed Name) *Melissa Taylor*

**Laboratory Name:** TA- Denver  
**Address:** 4955 Yarrow Street  
Arvada, CO 80002  
**Phone:** 303-736-0107  
**Contact:** Patrick McEntee

|                                                                                                                                                                                                                                                                                                        |  |  |  |  |  |  |  |  |  |                                                                                                                                   |  |  |  |  |  |  |  |  |  |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|--|--|--|--|--|--|--|--|-----------------------------------------------------------------------------------------------------------------------------------|--|--|--|--|--|--|--|--|--|
| Name Leidos<br>Address: 8866 Commons Blvd, Suite 201, Twinsburg, OH 44087<br>Phone Number: (330) 405-5802<br>Project Manager: Jed Thomas<br>Project: RVAAP FWGW Sampling Event Spring 2023<br>Job/P.O. No.: P010216426<br>Sampler (Signature) <i>Michelle Page</i> (Printed Name) <i>Michelle Page</i> |  |  |  |  |  |  |  |  |  | Laboratory Name: TA- Denver<br>Address: 4955 Yarrow Street<br>Arvada, CO 80002<br>Phone: 303-736-0107<br>Contact: Patrick McEntee |  |  |  |  |  |  |  |  |  |
| Requested Parameters                                                                                                                                                                                                                                                                                   |  |  |  |  |  |  |  |  |  | OBSERVATIONS, COMMENTS<br>SPECIAL INSTRUCTIONS                                                                                    |  |  |  |  |  |  |  |  |  |
| Total Number of Containers                                                                                                                                                                                                                                                                             |  |  |  |  |  |  |  |  |  | Temperature Blank                                                                                                                 |  |  |  |  |  |  |  |  |  |
| AMMONIA 15 (E)                                                                                                                                                                                                                                                                                         |  |  |  |  |  |  |  |  |  | Nitrate (10)(A)                                                                                                                   |  |  |  |  |  |  |  |  |  |
| 1                                                                                                                                                                                                                                                                                                      |  |  |  |  |  |  |  |  |  | 1                                                                                                                                 |  |  |  |  |  |  |  |  |  |
| 2                                                                                                                                                                                                                                                                                                      |  |  |  |  |  |  |  |  |  | 2                                                                                                                                 |  |  |  |  |  |  |  |  |  |
| 3                                                                                                                                                                                                                                                                                                      |  |  |  |  |  |  |  |  |  | 3                                                                                                                                 |  |  |  |  |  |  |  |  |  |
| 4                                                                                                                                                                                                                                                                                                      |  |  |  |  |  |  |  |  |  | 4                                                                                                                                 |  |  |  |  |  |  |  |  |  |
| 5                                                                                                                                                                                                                                                                                                      |  |  |  |  |  |  |  |  |  | 5                                                                                                                                 |  |  |  |  |  |  |  |  |  |
| 6                                                                                                                                                                                                                                                                                                      |  |  |  |  |  |  |  |  |  | 6                                                                                                                                 |  |  |  |  |  |  |  |  |  |
| 7                                                                                                                                                                                                                                                                                                      |  |  |  |  |  |  |  |  |  | 7                                                                                                                                 |  |  |  |  |  |  |  |  |  |
| 8                                                                                                                                                                                                                                                                                                      |  |  |  |  |  |  |  |  |  | 8                                                                                                                                 |  |  |  |  |  |  |  |  |  |
| 9                                                                                                                                                                                                                                                                                                      |  |  |  |  |  |  |  |  |  | 9                                                                                                                                 |  |  |  |  |  |  |  |  |  |
| 10                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 10                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 11                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 11                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 12                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 12                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 13                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 13                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 14                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 14                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 15                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 15                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 16                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 16                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 17                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 17                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 18                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 18                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 19                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 19                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 20                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 20                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 21                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 21                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 22                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 22                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 23                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 23                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 24                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 24                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 25                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 25                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 26                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 26                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 27                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 27                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 28                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 28                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 29                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 29                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 30                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 30                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 31                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 31                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 32                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 32                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 33                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 33                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 34                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 34                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 35                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 35                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 36                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 36                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 37                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 37                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 38                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 38                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 39                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 39                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 40                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 40                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 41                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 41                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 42                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 42                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 43                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 43                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 44                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 44                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 45                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 45                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 46                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 46                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 47                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 47                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 48                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 48                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 49                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 49                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 50                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 50                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 51                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 51                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 52                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 52                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 53                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 53                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 54                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 54                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 55                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 55                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 56                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 56                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 57                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 57                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 58                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 58                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 59                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 59                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 60                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 60                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 61                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 61                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 62                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 62                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 63                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 63                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 64                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 64                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 65                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 65                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 66                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 66                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 67                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 67                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 68                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 68                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 69                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 69                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 70                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 70                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 71                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 71                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 72                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 72                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 73                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 73                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 74                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 74                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 75                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 75                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 76                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 76                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 77                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 77                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 78                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 78                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 79                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 79                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 80                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 80                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 81                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 81                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 82                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 82                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 83                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 83                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 84                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 84                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 85                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 85                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 86                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 86                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 87                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 87                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 88                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 88                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 89                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 89                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 90                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 90                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 91                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 91                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 92                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 92                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 93                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 93                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 94                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 94                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 95                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 95                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 96                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 96                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 97                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 97                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 98                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 98                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 99                                                                                                                                                                                                                                                                                                     |  |  |  |  |  |  |  |  |  | 99                                                                                                                                |  |  |  |  |  |  |  |  |  |
| 100                                                                                                                                                                                                                                                                                                    |  |  |  |  |  |  |  |  |  | 100                                                                                                                               |  |  |  |  |  |  |  |  |  |
| 101                                                                                                                                                                                                                                                                                                    |  |  |  |  |  |  |  |  |  | 101                                                                                                                               |  |  |  |  |  |  |  |  |  |
| 102                                                                                                                                                                                                                                                                                                    |  |  |  |  |  |  |  |  |  | 102                                                                                                                               |  |  |  |  |  |  |  |  |  |
| 103                                                                                                                                                                                                                                                                                                    |  |  |  |  |  |  |  |  |  | 103                                                                                                                               |  |  |  |  |  |  |  |  |  |
| 104                                                                                                                                                                                                                                                                                                    |  |  |  |  |  |  |  |  |  | 104                                                                                                                               |  |  |  |  |  |  |  |  |  |
| 105                                                                                                                                                                                                                                                                                                    |  |  |  |  |  |  |  |  |  | 105                                                                                                                               |  |  |  |  |  |  |  |  |  |
| 106                                                                                                                                                                                                                                                                                                    |  |  |  |  |  |  |  |  |  | 106                                                                                                                               |  |  |  |  |  |  |  |  |  |
| 107                                                                                                                                                                                                                                                                                                    |  |  |  |  |  |  |  |  |  | 107                                                                                                                               |  |  |  |  |  |  |  |  |  |
| 108                                                                                                                                                                                                                                                                                                    |  |  |  |  |  |  |  |  |  | 108                                                                                                                               |  |  |  |  |  |  |  |  |  |
| 109                                                                                                                                                                                                                                                                                                    |  |  |  |  |  |  |  |  |  | 109                                                                                                                               |  |  |  |  |  |  |  |  |  |
| 110                                                                                                                                                                                                                                                                                                    |  |  |  |  |  |  |  |  |  | 110                                                                                                                               |  |  |  |  |  |  |  |  |  |
| 111                                                                                                                                                                                                                                                                                                    |  |  |  |  |  |  |  |  |  | 111                                                                                                                               |  |  |  |  |  |  |  |  |  |
| 112                                                                                                                                                                                                                                                                                                    |  |  |  |  |  |  |  |  |  | 112                                                                                                                               |  |  |  |  |  |  |  |  |  |
| 113                                                                                                                                                                                                                                                                                                    |  |  |  |  |  |  |  |  |  | 113                                                                                                                               |  |  |  |  |  |  |  |  |  |
| 114                                                                                                                                                                                                                                                                                                    |  |  |  |  |  |  |  |  |  | 114                                                                                                                               |  |  |  |  |  |  |  |  |  |
| 115                                                                                                                                                                                                                                                                                                    |  |  |  |  |  |  |  |  |  | 115                                                                                                                               |  |  |  |  |  |  |  |  |  |
| 116                                                                                                                                                                                                                                                                                                    |  |  |  |  |  |  |  |  |  | 116                                                                                                                               |  |  |  |  |  |  |  |  |  |
| 117                                                                                                                                                                                                                                                                                                    |  |  |  |  |  |  |  |  |  | 117                                                                                                                               |  |  |  |  |  |  |  |  |  |
| 118                                                                                                                                                                                                                                                                                                    |  |  |  |  |  |  |  |  |  | 118                                                                                                                               |  |  |  |  |  |  |  |  |  |
| 119                                                                                                                                                                                                                                                                                                    |  |  |  |  |  |  |  |  |  | 119                                                                                                                               |  |  |  |  |  |  |  |  |  |
| 120                                                                                                                                                                                                                                                                                                    |  |  |  |  |  |  |  |  |  | 120                                                                                                                               |  |  |  |  |  |  |  |  |  |
| 121                                                                                                                                                                                                                                                                                                    |  |  |  |  |  |  |  |  |  | 121                                                                                                                               |  |  |  |  |  |  |  |  |  |
| 122                                                                                                                                                                                                                                                                                                    |  |  |  |  |  |  |  |  |  | 122                                                                                                                               |  |  |  |  |  |  |  |  |  |
| 123                                                                                                                                                                                                                                                                                                    |  |  |  |  |  |  |  |  |  | 123                                                                                                                               |  |  |  |  |  |  |  |  |  |
| 124                                                                                                                                                                                                                                                                                                    |  |  |  |  |  |  |  |  |  | 124                                                                                                                               |  |  |  |  |  |  |  |  |  |
| 125                                                                                                                                                                                                                                                                                                    |  |  |  |  |  |  |  |  |  | 125                                                                                                                               |  |  |  |  |  |  |  |  |  |
| 126                                                                                                                                                                                                                                                                                                    |  |  |  |  |  |  |  |  |  | 126                                                                                                                               |  |  |  |  |  |  |  |  |  |
| 127                                                                                                                                                                                                                                                                                                    |  |  |  |  |  |  |  |  |  | 127                                                                                                                               |  |  |  |  |  |  |  |  |  |
| 128                                                                                                                                                                                                                                                                                                    |  |  |  |  |  |  |  |  |  | 128                                                                                                                               |  |  |  |  |  |  |  |  |  |
| 129                                                                                                                                                                                                                                                                                                    |  |  |  |  |  |  |  |  |  | 129                                                                                                                               |  |  |  |  |  |  |  |  |  |
| 130                                                                                                                                                                                                                                                                                                    |  |  |  |  |  |  |  |  |  | 130                                                                                                                               |  |  |  |  |  |  |  |  |  |
| 131                                                                                                                                                                                                                                                                                                    |  |  |  |  |  |  |  |  |  | 131                                                                                                                               |  |  |  |  |  |  |  |  |  |
| 132                                                                                                                                                                                                                                                                                                    |  |  |  |  |  |  |  |  |  | 132                                                                                                                               |  |  |  |  |  |  |  |  |  |
| 133                                                                                                                                                                                                                                                                                                    |  |  |  |  |  |  |  |  |  | 133                                                                                                                               |  |  |  |  |  |  |  |  |  |
| 134                                                                                                                                                                                                                                                                                                    |  |  |  |  |  |  |  |  |  | 134                                                                                                                               |  |  |  |  |  |  |  |  |  |
| 135                                                                                                                                                                                                                                                                                                    |  |  |  |  |  |  |  |  |  | 135                                                                                                                               |  |  |  |  |  |  |  |  |  |
| 136                                                                                                                                                                                                                                                                                                    |  |  |  |  |  |  |  |  |  | 136                                                                                                                               |  |  |  |  |  |  |  |  |  |
| 137                                                                                                                                                                                                                                                                                                    |  |  |  |  |  |  |  |  |  | 137                                                                                                                               |  |  |  |  |  |  |  |  |  |
| 138                                                                                                                                                                                                                                                                                                    |  |  |  |  |  |  |  |  |  | 138                                                                                                                               |  |  |  |  |  |  |  |  |  |
| 139                                                                                                                                                                                                                                                                                                    |  |  |  |  |  |  |  |  |  | 139                                                                                                                               |  |  |  |  |  |  |  |  |  |
| 140                                                                                                                                                                                                                                                                                                    |  |  |  |  |  |  |  |  |  | 140                                                                                                                               |  |  |  |  |  |  |  |  |  |
| 141                                                                                                                                                                                                                                                                                                    |  |  |  |  |  |  |  |  |  | 141                                                                                                                               |  |  |  |  |  |  |  |  |  |
| 142                                                                                                                                                                                                                                                                                                    |  |  |  |  |  |  |  |  |  | 142                                                                                                                               |  |  |  |  |  |  |  |  |  |
| 143                                                                                                                                                                                                                                                                                                    |  |  |  |  |  |  |  |  |  | 143                                                                                                                               |  |  |  |  |  |  |  |  |  |
|                                                                                                                                                                                                                                                                                                        |  |  |  |  |  |  |  |  |  |                                                                                                                                   |  |  |  |  |  |  |  |  |  |

Leidos

White Laboratory

Pink. Project Manager

Yellow Project QAO

*Goldenrod Field Project Manager*








1 From  
Date  
Sender's Name  
Company  
Address  
City  
State  
ZIP  
2 Your Internal Billing Reference

3 To  
Recipient's Name  
Company  
Address  
City  
State  
ZIP  
Phone  
Address  
We cannot deliver to P.O. boxes or P.O. ZIP codes.  
Address  
Use this line for the HOLD location address or for continuation of your shipping address.  
City  
State  
ZIP




8178 1968 1523

1 of 3  
TRK# 8178 1968 1523  
0215  
## MASTER ##  
NX LAAA  
CO-US DEN  
DSR AHS 80002  
PRIORITY OVERNIGHT  
THU - 18 MAY 10:30A



14 A  
DEN



Part # 156297-435 RHD EXP 11/23  
SHIP DATE: 17MAY23  
ACTWT: 51.15 LB  
CAD: 6993680/SSFE2420  
DIMS: 26x14x14 IN  
BILL RECIPIENT  
ORIGIN ID: CAKA (330) 998-4246  
RYAN LAURICH  
LEIDOS  
8866 COMMONS BLVD  
TWINSBURG, OH 44087  
UNITED STATES US  
TO EUROFINS ENVIRONMENT TE  
TESTAMERICA DENVER LAB  
4955 YARROW ST  
ARVADA CO 80002  
(303) 736-0100 REF: DEPT: PO:

4 Express Package Service  
To most locations.  
\* For packages over 150 lbs., use the  
FedEx Express Freight US Airbill.  
2 or 3 Business Days  
Next Business Day  
FedEx 2Day A.M.  
Second business morning.  
Saturday Delivery NOT available.  
FedEx 2Day  
Second business afternoon.  
Thursday shipments  
will be delivered on Monday unless Saturday  
Delivery is selected.  
FedEx Express Saver  
Third business day.  
Saturday Delivery NOT available.  
5 Packaging  
Declared value limit \$500.  
FedEx Envelope\*  
FedEx Pak\*  
Box  
Tube  
Other  
6 Special Handling and Delivery Signature Options  
Fees may apply. See the FedEx Service Guide.  
Saturday Delivery  
NOT available for FedEx Standard Overnight, FedEx 2Day A.M., or FedEx Express Saver.  
No Signature Required  
Package may be left without  
obtaining a signature for delivery.  
Direct Signature  
Someone at recipient's address  
may sign for delivery.  
Indirect Signature  
If no one is available at recipient's  
address, someone at a neighboring  
residential delivers only.  
Does this shipment contain dangerous goods?  
One box must be checked.  
Yes  
As per associated  
Shipper's Declaration  
not required.  
Dry Ice  
Dry Ice, 9 UN 1845  
kg  
7 Payment Bill to:  
Sender  
Acct. No. I will  
be billed.  
Section 1 will  
be billed.  
Recipient  
Third Party  
Enter FedEx Acct. No. below  
Obtain recip.  
FedEx Acct. No.  
Total Packages  
Total Weight  
lbs.  
Your liability is limited to US\$100 unless you declare a higher value. See the current  
FedEx Service Guide for details.  
Rev. Date 4/22 - Part #163134 - ©1994-2022 FedEx - PRINTED IN U.S.A.

For packages up to 150 lbs.  
FedEx Express Freight US Airbill.

**4 Express Package Service** \* To most locations.

**Next Business Day**

☐ FedEx First Overnight  
Earliest next business morning delivery to select locations. Friday shipments will be delivered on Monday unless Saturday Delivery is selected.

☐ FedEx Priority Overnight  
Next business morning. Friday shipments will be delivered on Monday unless Saturday Delivery is selected.

☐ FedEx Standard Overnight  
Next business afternoon.

☐ FedEx Express Saver  
Third business day.

☐ Saturday Delivery NOT available.

**2 or 3 Business Days**

☐ FedEx 2Day  
Second business afternoon. Thursday shipments will be delivered on Monday unless Saturday Delivery is selected.

☐ Saturday Delivery NOT available.

**5 Packaging** \*Declared value limit \$500.

☐ FedEx Envelope\* ☐ FedEx Pak\* ☐ FedEx Box ☐ FedEx Tube ☐ Other

**6 Special Handling and Delivery Signature Options** Fees may apply. See the FedEx Service Guide.

☐ Saturday Delivery  
NOT available for FedEx Standard Overnight, FedEx 2Day A.M., or FedEx Express Saver.

☐ No Signature Required  
Package may be left without obtaining a signature for delivery.

☐ Direct Signature  
Someone at recipient's address may sign for delivery.

☐ Indirect Signature  
If no one is available at recipient's address, someone at a neighboring residential delivery only.

**Does this shipment contain dangerous goods?**

☐ Yes ☐ No

☐ Shipper's Declaration ☐ Shipper's Declaration not required.

☐ Dry Ice ☐ Dry Ice, 9, UN 1845

☐ Cargo Aircraft Only

**7 Payment Bill to:**

☐ Sender ☐ Recipient ☐ Third Party

Enter FedEx Acct. No. below

Obtain recip. FedEx Acct No.

lbs. Total Packages Total Weight

Four liability is limited to \$500 unless you declare a higher value. See the current FedEx Service Guide for details.

Rev. Date 4/22 - Part 163134 - ©1994-2022 FedEx - PRINTED IN U.S.A.

Sender's Name  
Company  
Address  
City  
State  
ZIP

Recipient's Name  
Phone  
Company  
Address  
City  
State  
ZIP

Dept./Floor/Suite/Room

Use this line for the HOLD location address or for continuation of your shipping address.

Address  
City  
State  
ZIP

SHIP DATE: 17 MAY 23  
ACTWT: 49.75 LB  
CAD: 6998880/SSFE2420  
DIMS: 26x14x14 IN  
BILL RECIPIENT

ORIGIN TO: CACA (330) 998-4246  
RYAN LAURICH  
LEIDOS  
8866 COMMONS BLVD  
TWINSBURG, OH 44087  
UNITED STATES US

TO EUROFINS ENVIRONMENT TE  
TESTAMERICA DENVER LAB  
4955 YARROW ST  
ARVADA CO 80002  
(303) 738-0100  
REF: 0215


DEPT: 0215

Part # 156297-435 RRD8 EXP 11/23

3 of 3  
MPS# 3984 4714 4476  
OBR1  
Matr# 8178 1968 1523  
0215  
NX LAAA  
DSR AHS  
80002  
CO - US  
DEN

THU - 18 MAY 10:30A  
PRIORITY OVERNIGHT

FedEx Express



## Login Sample Receipt Checklist

Client: Leidos, Inc.

Job Number: 280-176674-1

**Login Number: 176674**

**List Source: Eurofins Denver**

**List Number: 1**

**Creator: Held, Wesley**

| Question                                                                                 | Answer | Comment |
|------------------------------------------------------------------------------------------|--------|---------|
| Radioactivity wasn't checked or is $\leq$ background as measured by a survey meter.      | N/A    |         |
| The cooler's custody seal, if present, is intact.                                        | True   |         |
| Sample custody seals, if present, are intact.                                            | True   |         |
| The cooler or samples do not appear to have been compromised or tampered with.           | True   |         |
| Samples were received on ice.                                                            | True   |         |
| Cooler Temperature is acceptable.                                                        | True   |         |
| Cooler Temperature is recorded.                                                          | True   |         |
| COC is present.                                                                          | True   |         |
| COC is filled out in ink and legible.                                                    | True   |         |
| COC is filled out with all pertinent information.                                        | True   |         |
| Is the Field Sampler's name present on COC?                                              | True   |         |
| There are no discrepancies between the containers received and the COC.                  | True   |         |
| Samples are received within Holding Time (excluding tests with immediate HTs)            | True   |         |
| Sample containers have legible labels.                                                   | True   |         |
| Containers are not broken or leaking.                                                    | True   |         |
| Sample collection date/times are provided.                                               | True   |         |
| Appropriate sample containers are used.                                                  | True   |         |
| Sample bottles are completely filled.                                                    | True   |         |
| Sample Preservation Verified.                                                            | N/A    |         |
| There is sufficient vol. for all requested analyses, incl. any requested MS/MSDs         | True   |         |
| Containers requiring zero headspace have no headspace or bubble is $<6\text{mm}$ (1/4"). | N/A    |         |
| Multiphasic samples are not present.                                                     | True   |         |
| Samples do not require splitting or compositing.                                         | True   |         |
| Residual Chlorine Checked.                                                               | N/A    |         |